

**REVIEWED**

By Melissa (Nikka) Bradley at 4:08 pm, Aug 09, 2019

Per GG, extraction performed 8/8/19

MB

MB

MB

Quantitative Analysis for Ethanol & Qualitative Analysis for Other Volatiles

Analytical Method(s): 1.0

Device: Hamilton MICROLAB 600A Liquid Processor/Dilutor Serial Number: ML600HC11378

Volatiles Quality Assurance Controls

Run Date(s): 08/08/19

Calibration date: 08/07/19 MB


Control level	Expiration	Lot #	Target Value	Acceptable Range	Overall Results
Level 1	Jan-22	1801036	0.0812	0.0731-0.0893	0.0792 g/100cc 0.0839 g/100cc g/100cc
Level 2	Mar-22	1803028	0.2035	0.1832-0.2238	0.2081 g/100cc g/100cc
Multi-Component mixture:			Lot #	FN06041502	OK
Curve Fit:			Column 1	Column 2	0.99999

Ethanol Calibration Reference Material		
Calibrator level	Target Value	Acceptable Range
50	0.050	0.045 - 0.055
100	0.100	0.090 - 0.110
200	0.200	0.180 - 0.220
300	0.300	0.270 - 0.330
500	0.500	0.450 - 0.550

Aqueous Controls		
Control level	Target Value	Acceptable Range
80	0.080	0.076 - 0.084
		Overall Results
		0.081 g/100cc

**Worklist: 3598**

<u>LAB CASE</u>	<u>ITEM</u>	<u>TASK ID</u>	<u>DESCRIPTION</u>
M2019-3380	1	158184	Alcohol Analysis
M2019-3391	1	158370	Alcohol Analysis
M2019-3431	1	158575	Alcohol Analysis
M2019-3437	1	158758	Alcohol Analysis
M2019-3440	1	158761	Alcohol Analysis
M2019-3467	1	158821	Alcohol Analysis
M2019-3468	1	158822	Alcohol Analysis
M2019-3470	1	158855	Alcohol Analysis
M2019-3472	1	158880	Alcohol Analysis
M2019-3477	1	158888	Alcohol Analysis
M2019-3502	1	159106	Alcohol Analysis
M2019-3539	1	159237	Alcohol Analysis
M2019-3540	1	159238	Alcohol Analysis
M2019-3541	1	159239	Alcohol Analysis
M2019-3556	1	159298	Alcohol Analysis
M2019-3557	1	159299	Alcohol Analysis
M2019-3575	1	159346	Alcohol Analysis
M2019-3587	2	159363	Alcohol Analysis
P2019-2239	3	158367	Alcohol Analysis



=====  
Calibration Table  
=====

-----  
General Calibration Setting  
-----

Calib. Data Modified : Wednesday, August 07, 2019 4:28:30 PM  
Signals calculated separately : No

Rel. Reference Window : 0.000 %  
Abs. Reference Window : 0.100 min  
Rel. Non-ref. Window : 0.000 %  
Abs. Non-ref. Window : 0.100 min  
Uncalibrated Peaks : not reported  
Partial Calibration : Yes, identified peaks are recalibrated  
Correct All Ret. Times: No, only for identified peaks

Curve Type : Linear  
Origin : Ignored  
Weight : Equal

Recalibration Settings:  
Average Response : Average all calibrations  
Average Retention Time: Floating Average New 75%

Calibration Report Options :  
Printout of recalibrations within a sequence:  
    Calibration Table after Recalibration  
    Normal Report after Recalibration  
If the sequence is done with bracketing:  
    Results of first cycle (ending previous bracket)

Default Sample ISTD Information (if not set in sample table):

ISTD #	ISTD Amount [g/100cc]	Name
1	1.00000	n-propanol
2	1.00000	n-propanol

-----  
Signal Details  
-----

Signal 1: FID1 A, Front Signal  
Signal 2: FID2 B, Back Signal  
-----

-----  
Overview Table  
-----

RT	Sig	Lvl	Amount [g/100cc]	Area	Rsp.Factor	Ref	ISTD #	Compound
2.586	1	1	1.00000	3.69669	2.70512e-1	No	No 1	methanol
2.809	1	1	1.00000	4.26100	2.34687e-1	No	No 2	Acetaldehyde
2.977	2	1	1.00000	4.26100	2.34687e-1	No	No 2	Acetaldehyde
3.075	1	1	5.00000e-2	4.48455	1.11494e-2	No	No 1	ethanol
		2	1.00000e-1	9.14104	1.09397e-2			
		3	2.00000e-1	18.45197	1.08390e-2			
		4	3.00000e-1	27.49827	1.09098e-2			
		5	5.00000e-1	46.02060	1.08647e-2			
3.388	2	1	1.00000	4.26062	2.34707e-1	No	No 2	methanol
3.628	1	1	1.00000	9.73055	1.02769e-1	No	No 1	isopropyl alcohol
4.285	2	1	5.00000e-2	4.67916	1.06857e-2	No	No 2	ethanol
		2	1.00000e-1	9.51962	1.05046e-2			
		3	2.00000e-1	19.37570	1.03222e-2			
		4	3.00000e-1	29.07134	1.03194e-2			
		5	5.00000e-1	49.06320	1.01909e-2			
4.308	1	1	1.00000	6.49940	1.53860e-1	No	No 1	acetone
4.620	1	1	1.00000	46.14289	2.16718e-2	No	Yes 1	n-propanol
		2	1.00000	47.04701	2.12553e-2			
		3	1.00000	47.25787	2.11605e-2			
		4	1.00000	46.80367	2.13658e-2			
		5	1.00000	47.08571	2.12379e-2			
4.661	2	1	1.00000	6.89301	1.45075e-1	No	No 2	acetone
4.969	2	1	1.00000	10.70642	9.34019e-2	No	No 2	isopropyl alcohol
7.550	2	1	1.00000	48.53417	2.06040e-2	No	Yes 2	n-propanol
		2	1.00000	49.19382	2.03278e-2			
		3	1.00000	49.32095	2.02754e-2			
		4	1.00000	48.84929	2.04711e-2			
		5	1.00000	49.17508	2.03355e-2			


Peak Sum Table

\*\*\*No Entries in table\*\*\*


1 Warnings or Errors :

Warning : Curve requires more calibration points., (methanol)


Calibration Curves




methanol at exp. RT: 2.586  
 FID1 A, Front Signal  
 Correlation: 1.00000  
 Residual Std. Dev.: 0.00000  
 Formula:  $y = mx + b$   
 m: 8.01141e-2  
 b: 0.00000  
 x: Amount Ratio  
 y: Area Ratio




Acetaldehyde at exp. RT: 2.809  
 FID1 A, Front Signal  
 Correlation: 1.00000  
 Residual Std. Dev.: 0.00000  
 Formula:  $y = mx + b$   
 m:  $8.77938e-2$   
 b: 0.00000  
 x: Amount Ratio  
 y: Area Ratio




Acetaldehyde at exp. RT: 2.977  
 FID2 B, Back Signal  
 Correlation: 1.00000  
 Residual Std. Dev.: 0.00000  
 Formula:  $y = mx + b$   
 m:  $8.77938e-2$   
 b: 0.00000  
 x: Amount Ratio  
 y: Area Ratio




ethanol at exp. RT: 3.075  
 FID1 A, Front Signal  
 Correlation: 1.00000  
 Residual Std. Dev.: 0.00081  
 Formula:  $y = mx + b$   
 m: 1.95745  
 b:  $-8.44568e-4$   
 x: Amount Ratio  
 y: Area Ratio




methanol at exp. RT: 3.388  
 FID2 B, Back Signal  
 Correlation: 1.00000  
 Residual Std. Dev.: 0.00000  
 Formula:  $y = mx + b$   
 m:  $8.77861e-2$   
 b: 0.00000  
 x: Amount Ratio  
 y: Area Ratio




isopropyl alcohol at exp. RT: 3.628  
 FID1 A, Front Signal  
 Correlation: 1.00000  
 Residual Std. Dev.: 0.00000  
 Formula:  $y = mx + b$   
 m: 2.10879e-1  
 b: 0.00000  
 x: Amount Ratio  
 y: Area Ratio




ethanol at exp. RT: 4.285  
 FID2 B, Back Signal  
 Correlation: 0.99999  
 Residual Std. Dev.: 0.00198  
 Formula:  $y = mx + b$   
 m: 2.00585  
 b: -6.22115e-3  
 x: Amount Ratio  
 y: Area Ratio




acetone at exp. RT: 4.308  
 FID1 A, Front Signal  
 Correlation: 1.00000  
 Residual Std. Dev.: 0.00000  
 Formula:  $y = mx + b$   
 m: 1.40854e-1  
 b: 0.00000  
 x: Amount Ratio  
 y: Area Ratio




n-propanol at exp. RT: 4.620  
 FID1 A, Front Signal  
 Correlation: 1.00000  
 Residual Std. Dev.: 0.00000  
 Formula:  $y = mx + b$   
 m: 1.00000  
 b: 0.00000  
 x: Amount Ratio  
 y: Area Ratio



acetone at exp. RT: 4.661  
FID2 B, Back Signal  
Correlation: 1.00000  
Residual Std. Dev.: 0.00000  
Formula:  $y = mx + b$   
m: 1.42024e-1  
b: 0.00000  
x: Amount Ratio  
y: Area Ratio



isopropyl alcohol at exp. RT: 4.969  
FID2 B, Back Signal  
Correlation: 1.00000  
Residual Std. Dev.: 0.00000  
Formula:  $y = mx + b$   
m: 2.20595e-1  
b: 0.00000  
x: Amount Ratio  
y: Area Ratio




n-propanol at exp. RT: 7.550  
FID2 B, Back Signal  
Correlation: 1.00000  
Residual Std. Dev.: 0.00000  
Formula:  $y = mx + b$   
m: 1.00000  
b: 0.00000  
x: Amount Ratio  
y: Area Ratio

=====

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.050 FN04271601  
 Laboratory : Meridian  
 Injection Date : Aug 7, 2019  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	4.48455	0.0501	g/100cc
2.	Ethanol	Column 2:	4.67916	0.0512	g/100cc
3.	n-Propanol	Column 1:	46.14289	1.0000	g/100cc
4.	n-Propanol	Column 2:	48.53417	1.0000	g/100cc

*Handwritten signature*



ISP Forensic Services Blood Alcohol Report

Sample Name : 0.100 FN02271802  
 Laboratory : Meridian  
 Injection Date : Aug 7, 2019  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	9.14104	0.0997	g/100cc
2.	Ethanol	Column 2:	9.51962	0.0996	g/100cc
3.	n-Propanol	Column 1:	47.04701	1.0000	g/100cc
4.	n-Propanol	Column 2:	49.19382	1.0000	g/100cc

*[Handwritten signature]*

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.200 FN03301601  
 Laboratory : Meridian  
 Injection Date : Aug 7, 2019  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	18.45197	0.1999	g/100cc
2.	Ethanol	Column 2:	19.37570	0.1990	g/100cc
3.	n-Propanol	Column 1:	47.25787	1.0000	g/100cc
4.	n-Propanol	Column 2:	49.32095	1.0000	g/100cc

*Handwritten signature*

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.300 FN07311804  
 Laboratory : Meridian  
 Injection Date : Aug 7, 2019  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	27.49827	0.3006	g/100cc
2.	Ethanol	Column 2:	29.07134	0.2998	g/100cc
3.	n-Propanol	Column 1:	46.80367	1.0000	g/100cc
4.	n-Propanol	Column 2:	48.84929	1.0000	g/100cc

68

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.500 FN08031602  
 Laboratory : Meridian  
 Injection Date : Aug 7, 2019  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	46.02060	0.4997	g/100cc
2.	Ethanol	Column 2:	49.06320	0.5005	g/100cc
3.	n-Propanol	Column 1:	47.08571	1.0000	g/100cc
4.	n-Propanol	Column 2:	49.17508	1.0000	g/100cc

*SR*

ISP Forensic Services Blood Alcohol Report

Sample Name : INTERNAL STANDARD BLANK  
 Laboratory : Meridian  
 Injection Date : Aug 7, 2019  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	0.00000	0.0000	g/100cc
2.	Ethanol	Column 2:	0.00000	0.0000	g/100cc
3.	n-Propanol	Column 1:	46.62593	1.0000	g/100cc
4.	n-Propanol	Column 2:	48.73176	1.0000	g/100cc

6r

S a m p l e S u m m a r y

Sequence table: C:\Chem32\1\Data\08-07-19\_CAL\08-07-19\_CAL 2019-08-07 15-23-20\08-07-19\_CAL.S  
 Data directory path: C:\Chem32\1\Data\08-07-19\_CAL\08-07-19\_CAL 2019-08-07 15-23-20\  
 Logbook: C:\Chem32\1\Data\08-07-19\_CAL\08-07-19\_CAL 2019-08-07 15-23-20\08-07-19\_CAL.LOG  
 Sequence start: 8/7/2019 3:37:58 PM  
 Sequence Operator: SYSTEM  
 Operator: SYSTEM

Method file name: C:\Chem32\1\Data\08-07-19\_CAL\08-07-19\_CAL 2019-08-07 15-23-20\ALCOHOL.M


Run #	Location #	Inj #	Sample Name	Sample Amt [g/100cc]	Multip.* Dilution	File name	Cal #	# Cmp
1	1	1	0.050 FN04271601	-	1.0000	001F0101.D	*	4
2	2	1	0.100 FN02271802	-	1.0000	002F0201.D	*	4
3	3	1	0.200 FN03301601	-	1.0000	003F0301.D	*	4
4	4	1	0.300 FN07311804	-	1.0000	004F0401.D	*	4
5	5	1	0.500 FN08031602	-	1.0000	005F0501.D	*	4
6	6	1	INTERNAL STANDAR	-	1.0000	006F0601.D		2

Method file name: C:\Chem32\1\Data\08-07-19\_CAL\08-07-19\_CAL 2019-08-07 15-23-20\SHUTDOWN.M

Run #	Location #	Inj #	Sample Name	Sample Amt [g/100cc]	Multip.* Dilution	File name	Cal #	# Cmp
7	7	1	SHUTDOWN	-	1.0000	007F0701.D		0

ISP Forensic Services Blood Alcohol Report

Sample Name : INTERNAL STD BLK  
 Laboratory : Meridian  
 Injection Date : Aug 8, 2019  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	0.00000	0.0000	g/100cc
2.	Ethanol	Column 2:	0.00000	0.0000	g/100cc
3.	n-Propanol	Column 1:	45.02013	1.0000	g/100cc
4.	n-Propanol	Column 2:	46.89166	1.0000	g/100cc

62

ISP Forensic Services Blood Alcohol Report

Sample Name : MIX VOL FN06041502  
 Laboratory : Meridian  
 Injection Date : Aug 8, 2019  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	6.53519	0.1350	g/100cc
2.	Ethanol	Column 2:	6.74381	0.1359	g/100cc
3.	n-Propanol	Column 1:	24.80456	1.0000	g/100cc
4.	n-Propanol	Column 2:	25.31432	1.0000	g/100cc

*Handwritten signature*



## VOLATILES DETERMINATION CASEFILE WORKSHEET

Laboratory No.: QC1-1

Analysis Date(s): 08 Aug 2019

	Column 1 FID A	Column 2 FID B	Column Precision	Mean Value	Over-all Mean	
Sample Results	0.0790	0.0794	0.0004	0.0792	0.0792	
(g/100cc)	0.0790	0.0796	0.0006	0.0793		

### Analysis Method

Refer to Blood Alcohol Method #1

### Instrument Information

*Instrument method is stored centrally.*

Refer to Instrument Method: Alcohol.m  
Hamilton Auto-Dilutor Serial Number: ML600HC11378

### Reporting of Results

Uncertainty of Measurement (UM%): 5.00%


Overall Mean (g/100cc)	Low	High	5% of Mean
0.079	0.075	0.083	0.004

	Reported Result	
	0.079	

*Calibration and control data are stored centrally.*

ISP Forensic Services Blood Alcohol Report

Sample Name : QC1-1-A  
 Laboratory : Meridian  
 Injection Date : Aug 8, 2019  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	6.86247	0.0790	g/100cc
2.	Ethanol	Column 2:	7.07989	0.0794	g/100cc
3.	n-Propanol	Column 1:	44.61438	1.0000	g/100cc
4.	n-Propanol	Column 2:	46.23821	1.0000	g/100cc

*W*

ISP Forensic Services Blood Alcohol Report

Sample Name : QC1-1-B  
 Laboratory : Meridian  
 Injection Date : Aug 8, 2019  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014 - CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	6.86053	0.0790	g/100cc
2.	Ethanol	Column 2:	7.08202	0.0796	g/100cc
3.	n-Propanol	Column 1:	44.60306	1.0000	g/100cc
4.	n-Propanol	Column 2:	46.18317	1.0000	g/100cc

*Handwritten signature*

## VOLATILES DETERMINATION CASEFILE WORKSHEET

Laboratory No.: 0.08 FN04171701

Analysis Date(s): 08 Aug 2019

	Column 1 FID A	Column 2 FID B	Column Precision	Mean Value	Over-all Mean	
Sample Results	0.0808	0.0813	0.0005	0.0810	0.0815	
(g/100cc)	0.0818	0.0823	0.0005	0.0820		

### Analysis Method

Refer to Blood Alcohol Method #1

### Instrument Information

*Instrument method is stored centrally.*

Refer to Instrument Method: Alcohol.m  
Hamilton Auto-Dilutor Serial Number: ML600HC11378

### Reporting of Results

Uncertainty of Measurement (UM%): 5.00%


Overall Mean (g/100cc)	Low	High	5% of Mean
0.081	0.076	0.086	0.005

	Reported Result	
	0.081	

*Calibration and control data are stored centrally.*

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.08 FN04171701-A  
 Laboratory : Meridian  
 Injection Date : Aug 8, 2019  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	7.15766	0.0808	g/100cc
2.	Ethanol	Column 2:	7.42998	0.0813	g/100cc
3.	n-Propanol	Column 1:	45.49435	1.0000	g/100cc
4.	n-Propanol	Column 2:	47.35122	1.0000	g/100cc

*W*

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.08 FN04171701-B  
 Laboratory : Meridian  
 Injection Date : Aug 8, 2019  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	7.02751	0.0818	g/100cc
2.	Ethanol	Column 2:	7.27327	0.0823	g/100cc
3.	n-Propanol	Column 1:	44.13321	1.0000	g/100cc
4.	n-Propanol	Column 2:	45.78617	1.0000	g/100cc

br

## VOLATILES DETERMINATION CASEFILE WORKSHEET

Laboratory No.: QC2-1

Analysis Date(s): 08 Aug 2019

	Column 1 FID A	Column 2 FID B	Column Precision	Mean Value	Over-all Mean	
Sample Results	0.2073	0.2076	0.0003	0.2074	0.2081	
(g/100cc)	0.2088	0.2088	0.0000	0.2088		

### Analysis Method

Refer to Blood Alcohol Method #1

### Instrument Information

*Instrument method is stored centrally.*

Refer to Instrument Method: Alcohol.m  
Hamilton Auto-Dilutor Serial Number: ML600HC11378

### Reporting of Results

Uncertainty of Measurement (UM%): 5.00%

Overall Mean (g/100cc)	Low	High	5% of Mean
0.208	0.197	0.219	0.011

	Reported Result	
	0.208	

*Calibration and control data are stored centrally.*

Revision: 1


Issue Date: 01/04/2019

Issuing Authority: Quality Manager



ISP Forensic Services Blood Alcohol Report

Sample Name : QC2-1-A  
 Laboratory : Meridian  
 Injection Date : Aug 8, 2019  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	18.40351	0.2073	g/100cc
2.	Ethanol	Column 2:	19.24864	0.2076	g/100cc
3.	n-Propanol	Column 1:	45.45421	1.0000	g/100cc
4.	n-Propanol	Column 2:	46.93238	1.0000	g/100cc

*SW*



ISP Forensic Services Blood Alcohol Report

Sample Name : QC2-1-B  
 Laboratory : Meridian  
 Injection Date : Aug 8, 2019  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	18.75864	0.2088	g/100cc
2.	Ethanol	Column 2:	19.62537	0.2088	g/100cc
3.	n-Propanol	Column 1:	46.00233	1.0000	g/100cc
4.	n-Propanol	Column 2:	47.56442	1.0000	g/100cc

*Handwritten signature*

**VOLATILES DETERMINATION CASEFILE WORKSHEET**

Laboratory No.: QC1-2

Analysis Date(s): 08 Aug 2019

	Column 1 FID A	Column 2 FID B	Column Precision	Mean Value	Over-all Mean	
Sample Results	0.0830	0.0832	0.0002	0.0831	0.0839	
(g/100cc)	0.0847	0.0847	0.0000	0.0847		

**Analysis Method**

Refer to Blood Alcohol Method #1

**Instrument Information**

*Instrument method is stored centrally.*

Refer to Instrument Method: Alcohol.m  
Hamilton Auto-Dilutor Serial Number: ML600HC11378

**Reporting of Results**

Uncertainty of Measurement (UM%): 5.00%

Overall Mean (g/100cc)	Low	High	5% of Mean
0.083	0.078	0.088	0.005

	Reported Result	
	0.083	

*Calibration and control data are stored centrally.*




Revision: 1

Issue Date: 01/04/2019

Issuing Authority: Quality Manager

ISP Forensic Services Blood Alcohol Report

Sample Name : QC1-2-A  
 Laboratory : Meridian  
 Injection Date : Aug 8, 2019  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	7.49067	0.0830	g/100cc
2.	Ethanol	Column 2:	7.71703	0.0832	g/100cc
3.	n-Propanol	Column 1:	46.36748	1.0000	g/100cc
4.	n-Propanol	Column 2:	48.03090	1.0000	g/100cc

BT

ISP Forensic Services Blood Alcohol Report

Sample Name : QC1-2-B  
 Laboratory : Meridian  
 Injection Date : Aug 8, 2019  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	7.75443	0.0847	g/100cc
2.	Ethanol	Column 2:	7.95675	0.0847	g/100cc
3.	n-Propanol	Column 1:	47.01603	1.0000	g/100cc
4.	n-Propanol	Column 2:	48.63589	1.0000	g/100cc

*Handwritten signature*

ISP Forensic Services Blood Alcohol Report

Sample Name : INTERNAL STD BLK  
 Laboratory : Meridian  
 Injection Date : Aug 8, 2019  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	0.00000	0.0000	g/100cc
2.	Ethanol	Column 2:	0.00000	0.0000	g/100cc
3.	n-Propanol	Column 1:	45.71836	1.0000	g/100cc
4.	n-Propanol	Column 2:	47.23723	1.0000	g/100cc

67

S a m p l e S u m m a r y

Sequence table: C:\Chem32\1\Data\08-08-19\_SAMPLES\08-08-19t\_SAMPLES 2019-08-08 10-42-55\0-08-19t\_SAMPLES.S  
 Data directory path: C:\Chem32\1\Data\08-08-19\_SAMPLES\08-08-19t\_SAMPLES 2019-08-08 10-42-55\  
 Logbook: C:\Chem32\1\Data\08-08-19\_SAMPLES\08-08-19t\_SAMPLES 2019-08-08 10-42-55\0-08-19t\_SAMPLES.LOG  
 Sequence start: 8/8/2019 10:57:41 AM  
 Sequence Operator: SYSTEM  
 Operator: SYSTEM  
 Method file name: C:\Chem32\1\Data\08-08-19\_SAMPLES\08-08-19t\_SAMPLES 2019-08-08 10-42-55\ALCOHOL.M

Run #	Location #	Inj #	Sample Name	Sample Amt [g/100cc]	Multip.* Dilution	File name	Cal #	# Cmp
1	1	1	INTERNAL STD BLK	-	1.0000	001F0101.D		2
2	2	1	BOT0235LOT19803-	-	1.0000	002F0201.D		4
3	3	1	BOT0235LOT19803-	-	1.0000	003F0301.D		4
4	4	1	INTERNAL STD BLK	-	1.0000	004F0401.D		2
5	5	1	BOT0070LOT19803-	-	1.0000	005F0501.D		4
6	6	1	BOT0070LOT19803-	-	1.0000	006F0601.D		4
7	7	1	INTERNAL STD BLK	-	1.0000	007F0701.D		2
8	8	1	MIX VOL FN060415	-	1.0000	008F0801.D		10
9	9	1	QC1-1-A	-	1.0000	009F0901.D		4
10	10	1	QC1-1-B	-	1.0000	010F1001.D		4
11	11	1	0.08 FN04171701-	-	1.0000	011F1101.D		4
12	12	1	0.08 FN04171701-	-	1.0000	012F1201.D		4
13	13	1	M2019-3380-1-A	-	1.0000	013F1301.D		2
14	14	1	M2019-3380-1-B	-	1.0000	014F1401.D		2
15	15	1	M2019-3391-1-A	-	1.0000	015F1501.D		4
16	16	1	M2019-3391-1-B	-	1.0000	016F1601.D		4
17	17	1	M2019-3431-1-A	-	1.0000	017F1701.D		4
18	18	1	M2019-3431-1-B	-	1.0000	018F1801.D		4
19	19	1	M2019-3437-1-A	-	1.0000	019F1901.D		4
20	20	1	M2019-3437-1-B	-	1.0000	020F2001.D		4
21	21	1	M2019-3440-1-A	-	1.0000	021F2101.D		4
22	22	1	M2019-3440-1-B	-	1.0000	022F2201.D		4
23	23	1	M2019-3467-1-A	-	1.0000	023F2301.D		4
24	24	1	M2019-3467-1-B	-	1.0000	024F2401.D		4
25	25	1	M2019-3468-1-A	-	1.0000	025F2501.D		4
26	26	1	M2019-3468-1-B	-	1.0000	026F2601.D		4
27	27	1	M2019-3470-1-A	-	1.0000	027F2701.D		4
28	28	1	M2019-3470-1-B	-	1.0000	028F2801.D		4
29	29	1	M2019-3472-1-A	-	1.0000	029F2901.D		4
30	30	1	M2019-3472-1-B	-	1.0000	030F3001.D		4
31	31	1	QC2-1-A	-	1.0000	031F3101.D		4
32	32	1	QC2-1-B	-	1.0000	032F3201.D		4
33	33	1	M2019-3477-1-A	-	1.0000	033F3301.D		4
34	34	1	M2019-3477-1-B	-	1.0000	034F3401.D		4
35	35	1	M2019-3502-1-A	-	1.0000	035F3501.D		4
36	36	1	M2019-3502-1-B	-	1.0000	036F3601.D		4
37	37	1	M2019-3539-1-A	-	1.0000	037F3701.D		4
38	38	1	M2019-3539-1-B	-	1.0000	038F3801.D		4
39	39	1	M2019-3540-1-A	-	1.0000	039F3901.D		4
40	40	1	M2019-3540-1-B	-	1.0000	040F4001.D		4
41	41	1	M2019-3541-1-A	-	1.0000	041F4101.D		4
42	42	1	M2019-3541-1-B	-	1.0000	042F4201.D		4
43	43	1	M2019-3556-1-A	-	1.0000	043F4301.D		4

Run #	Location	Inj #	Sample Name	Sample Amt [g/100cc]	Multip.* Dilution	File name	Cal #
44	44	1	M2019-3556-1-B	-	1.0000	044F4401.D	4
45	45	1	M2019-3557-1-A	-	1.0000	045F4501.D	4
46	46	1	M2019-3557-1-B	-	1.0000	046F4601.D	4
47	47	1	M2019-3575-1-A	-	1.0000	047F4701.D	4
48	48	1	M2019-3575-1-B	-	1.0000	048F4801.D	4
49	49	1	M2019-3587-2-A	-	1.0000	049F4901.D	4
50	50	1	M2019-3587-2-B	-	1.0000	050F5001.D	4
51	51	1	P2019-2239-3-A	-	1.0000	051F5101.D	2
52	52	1	P2019-2239-3-B	-	1.0000	052F5201.D	2
53	53	1	QC1-2-A	-	1.0000	053F5301.D	4
54	54	1	QC1-2-B	-	1.0000	054F5401.D	4
55	55	1	INTERNAL STD BLK	-	1.0000	055F5501.D	2

Method file name: C:\Chem32\1\Data\08-08-19\_SAMPLES\08-08-19t\_SAMPLES 2019-08-08 10-42-55  
 \SHUTDOWN.M

Run #	Location	Inj #	Sample Name	Sample Amt [g/100cc]	Multip.* Dilution	File name	Cal #
56	56	1	EMPTY	-	1.0000	056F5601.D	0