


**Worklist: 1847**

*Bylee*

<u>LAB CASE</u>	<u>ITEM</u>	<u>TASK ID</u>	<u>DESCRIPTION</u>
C2017-1459	2	92399	AM 27 Blood THC Quant by LC
C2017-1467	1	92400	AM 27 Blood THC Quant by LC
C2017-1488	1	92403	AM 27 Blood THC Quant by LC
C2017-1493	1	92401	AM 27 Blood THC Quant by LC
C2017-1495	1	92402	AM 27 Blood THC Quant by LC
C2017-1520	1	92404	AM 27 Blood THC Quant by LC
M2017-3455	1	92405	AM 27 Blood THC Quant by LC
P2017-1759	1	92406	AM 27 Blood THC Quant by LC



*[Handwritten signature]*

# Quantitation of THC and Metabolites in Blood by LC-MS/MS

Extraction Date: 8-8-17  
A

Analyst: Anne Nord

## PRE-ANALYTIC

Plate Lot# Custom - 499102 Plate Exp. 1/29/2018 External QC Lot 91317, exp 9-13-17

✓ 1. Ensure all solutions are within expiration date.

- Mobile Phase A: **0.1% Formic Acid in LCMS Water**
- Mobile Phase B: **0.1% Formic Acid in LCMS Acetonitrile**
- **LCMS Methanol**
- **Blank/Negative Blood: Lot 321632-1**
- **0.1% Formic Acid in water**
- **MTBE**
- **Hexane**

✓ Column: **UCT Selectra DA 100 x 2.1 mm 3um**

2. Check levels of mobile phases and needle wash and refill as necessary. Ensure waste is not full.

✓ 3. Purge Pump and Load appropriate Acq. Method, allow system to equilibrate for approx. 30 min.

✓ 4. Create worklist. Data path name: 8917 cann quant

## ANALYTIC

✓ 1. Remove standards plate, blood, and samples from cold storage. Allow to reach room temperature.

✓ 2. Add **1000 µL blood** to wells of analytical (standards) plate. Place cover on Plate

Blank blood for locations containing standards/QCs and internal standards

Sample blood for locations containing only internal standards

✓ 3. Place on shaking incubator at ambient temp., **900rpm for 15 minutes**. *Shaker ID 66759*

✓ 4. Pipette **500µL 0.1% formic acid** to all wells of standards plate.

✓ 5. Place on shaking incubator at ambient temp., **900rpm for 15 minutes**.

✓ 6. Transfer **800µL of blood+acid mixture** to corresponding wells of SLE+ plate.

✓ 7. Apply positive pressure for approx. 10-15 seconds (or until no liquid remains on top of sorbent). Wait 5 min.  
**(Load blood samples at 85- 100 PSI- Selector to Right)**

✓ 8. Add **2.25mL MTBE** and allow to flow under gravity for 5 minutes. **(add in 3 increments of 750uL)**

✓ 9. Apply positive pressure for approx. 15 seconds **(10-15 PSI- Selector to left -)**.

✓ 10. Add **2.25mL Hexane** and allow to flow under gravity for 5 minutes. **(add in 3 increments of 750uL)**

✓ 11. Apply positive pressure for approx. 15 seconds. **(10-15 PSI Selector to the left)**

✓ 12. Remove collection plate containing eluate.

✓ 13. Place collection plate on SPE Dry and evaporate to dryness at approx. 35°C. *SPE Dry ID 66819*

✓ 14. Reconstitute in **100 µL MeOH** and heat seal plate with foil. Place in autosampler and run worklist.

## POST-ANALYTIC

✓ 1. Open quantitation software and create a new quantitation batch.

Batch name: 8917 cann quant

✓ 2. Make any necessary integration changes. Limit curves based on validated linear ranges (3-50ng/mL).

✓ 3. Were all appropriate standards used in the curve for each analyte? **(Y)**/N

Are r<sup>2</sup> values ≥0.98 for each analyte? **(Y)**/N

✓ 4. Did all QCs pass for each analyte? **(Y)**/N Were QCs entered into QC charting? **(Y)**/N

5. Central File Packet to include: ✓ LIMS Worklist: ✓ Method Checklist ✓ Calibration and Control Reports

COMMENTS

# ISP FORENSICS - Cd'A Instrument # 62340


## Cannabinoids Analysis Report

**Batch Data Path** D:\2017 Data\8917 cann quant\QuantResults\8917 cann quant.batch.bin  
**Analysis Time** 8/10/2017 3:47 PM **Analyst Name** ISP Tox  
**Report Time** 8/10/2017 3:48 PM **Reporter Name** ISP Tox  
**Last Calib Update** 8/10/2017 3:47 PM **Batch State** Processed

### Analysis Info

**Acq Time** 2017-08-10 11:04 **Data File** Negative Control.d  
**Sample Type** Sample **Sample Name** Negative Control  
**Dilution** 1 **Acq Method** AM 27 Quant THC 7-2017.m  
**Position** P1-A2 **Sample Info**  
**Inj Vol** -1 **Comment** AM 27 Cannabinoid Confirmation

### Sample Chromatogram



### Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-COOH	THC-COOH-d9	2.366	6243	74804	0.0835	1.3723

# ISP FORENSICS - Cd'A Instrument # 62340


## Cannabinoids Analysis Report

**Batch Data Path** D:\2017 Data\8917 cann quant\QuantResults\8917 cann quant.batch.bin  
**Analysis Time** 8/10/2017 3:47 PM **Analyst Name** ISP Tox  
**Report Time** 8/10/2017 3:48 PM **Reporter Name** ISP Tox  
**Last Calib Update** 8/10/2017 3:47 PM **Batch State** Processed

### Analysis Info

**Acq Time** 2017-08-10 11:16 **Data File** QC - 10ng.d  
**Sample Type** QC **Sample Name** QC - 10ng/mL  
**Dilution** 1 **Acq Method** AM 27 Quant THC 7-2017.m  
**Position** P1-H1 **Sample Info**  
**Inj Vol** -1 **Comment** AM 27 Cannabinoid Confirmation

### Sample Chromatogram



### Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.336	20548	206921	0.0993	10.5906
THC-COOH	THC-COOH-d9	2.466	20976	70373	0.2981	11.8664
THC	THC-d3	6.493	11413	88303	0.1293	10.8835

# ISP FORENSICS - Cd'A Instrument # 62340


## Cannabinoids Analysis Report

**Batch Data Path** D:\2017 Data\8917 cann quant\QuantResults\8917 cann quant.batch.bin  
**Analysis Time** 8/10/2017 3:47 PM **Analyst Name** ISP Tox  
**Report Time** 8/10/2017 3:48 PM **Reporter Name** ISP Tox  
**Last Calib Update** 8/10/2017 3:47 PM **Batch State** Processed

### Analysis Info

**Acq Time** 2017-08-10 11:28 **Data File** External Control lot 91317.d  
**Sample Type** Sample **Sample Name** External Control lot 91317  
**Dilution** 1 **Acq Method** AM 27 Quant THC 7-2017.m  
**Position** P1-B2 **Sample Info**  
**Inj Vol** -1 **Comment** AM 27 Cannabinoid Confirmation

### Sample Chromatogram



### Results


Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.336	21155	239785	0.0882	9.3779
THC-COOH	THC-COOH-d9	2.446	14889	70815	0.2103	7.5723
THC	THC-d3	6.493	8206	94704	0.0867	7.1800

# ISP Forensics Calibration Curve Report

**Batch Data Path** D:\2017 Data\8917 cann quant\QuantResults\8917 cann quant.batch.bin

**Last Calib Update** 8/10/2017 3:47 PM **Analyst Name** ISP TOX

**Target Compound** *THC-OH*  
**Internal Standard** *THC-OH-d3*




Sample	Level	Enabled	Exp Conc	Final Conc	Accuracy
Cal 1 - 3ng/mL	1	<input checked="" type="checkbox"/>	3	2.9	97.4
Cal 2 - 5ng/mL	2	<input checked="" type="checkbox"/>	5	4.7	94.6
Cal 3 - 10ng/mL	3	<input checked="" type="checkbox"/>	10	10.7	107.4
QC - 10ng/mL	3	<input checked="" type="checkbox"/>	10	10.6	105.9
Cal 4 - 25ng/mL	4	<input checked="" type="checkbox"/>	25	25.6	102.2
Cal 5 - 50ng/mL	5	<input checked="" type="checkbox"/>	50	49.3	98.6
Cal 6 - 100ng/mL	6	<input checked="" type="checkbox"/>	100	99.9	99.9
Cal 7 - 250ng/mL	7	<input checked="" type="checkbox"/>	250	249.9	100.0

# ISP Forensics Calibration Curve Report

**Batch Data Path** D:\2017 Data\8917 cann quant\QuantResults\8917 cann quant.batch.bin

**Last Calib Update** 8/10/2017 3:47 PM **Analyst Name** ISP TOX

**Target Compound** *THC-COOH*  
**Internal Standard** *THC-COOH-d9*




Sample	Level	Enabled	Exp Conc	Final Conc	Accuracy
Cal 1 - 3ng/mL	1	<input checked="" type="checkbox"/>	3	3.0	100.2
Cal 2 - 5ng/mL	2	<input checked="" type="checkbox"/>	5	4.6	92.1
Cal 3 - 10ng/mL	3	<input checked="" type="checkbox"/>	10	11.4	114.1
QC - 10ng/mL	3	<input checked="" type="checkbox"/>	10	11.9	118.7
Cal 4 - 25ng/mL	4	<input checked="" type="checkbox"/>	25	26.4	105.5
Cal 5 - 50ng/mL	5	<input checked="" type="checkbox"/>	50	49.1	98.2
Cal 6 - 100ng/mL	6	<input checked="" type="checkbox"/>	100	95.1	95.1
Cal 7 - 250ng/mL	7	<input checked="" type="checkbox"/>	250	236.9	94.8

# ISP Forensics Calibration Curve Report

**Batch Data Path** D:\2017 Data\8917 cann quant\QuantResults\8917 cann quant.batch.bin

**Last Calib Update** 8/10/2017 3:47 PM **Analyst Name** ISP TOX

**Target Compound** *THC*  
**Internal Standard** *THC-d3*



Sample	Level	Enabled	Exp Conc	Final Conc	Accuracy
Cal 1 - 3ng/mL	1	<input checked="" type="checkbox"/>	3	2.9	96.8
Cal 2 - 5ng/mL	2	<input checked="" type="checkbox"/>	5	5.0	100.2
Cal 3 - 10ng/mL	3	<input checked="" type="checkbox"/>	10	10.1	101.0
QC - 10ng/mL	3	<input checked="" type="checkbox"/>	10	10.9	108.8
Cal 4 - 25ng/mL	4	<input checked="" type="checkbox"/>	25	25.9	103.6
Cal 5 - 50ng/mL	5	<input checked="" type="checkbox"/>	50	50.1	100.1
Cal 6 - 100ng/mL	6	<input checked="" type="checkbox"/>	100	97.9	97.9
Cal 7 - 250ng/mL	7	<input checked="" type="checkbox"/>	250	251.1	100.5



# ISP FORENSICS - Cd'A Instrument # 62340


## Cannabinoids Analysis Report

**Batch Data Path** D:\2017 Data\8917 cann quant\QuantResults\8917 cann quant.batch.bin  
**Analysis Time** 8/10/2017 3:47 PM **Analyst Name** ISP Tox  
**Report Time** 8/10/2017 3:48 PM **Reporter Name** ISP Tox  
**Last Calib Update** 8/10/2017 3:47 PM **Batch State** Processed

### Analysis Info

**Acq Time** 2017-08-10 09:29 **Data File** Cal 1 - 3ng.d  
**Sample Type** Calibration **Sample Name** Cal 1 - 3ng/mL  
**Dilution** 1 **Acq Method** AM 27 Quant THC 7-2017.m  
**Position** P1-A1 **Sample Info**  
**Inj Vol** -1 **Comment** AM 27 Cannabinoid Confirmation

### Sample Chromatogram



### Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.315	5500	188138	0.0292	2.9217
THC-COOH	THC-COOH-d9	2.426	6944	59420	0.1169	3.0054
THC	THC-d3	6.453	2888	77093	0.0375	2.9027

# ISP FORENSICS - Cd'A Instrument # 62340


## Cannabinoids Analysis Report

**Batch Data Path** D:\2017 Data\8917 cann quant\QuantResults\8917 cann quant.batch.bin  
**Analysis Time** 8/10/2017 3:47 PM **Analyst Name** ISP Tox  
**Report Time** 8/10/2017 3:48 PM **Reporter Name** ISP Tox  
**Last Calib Update** 8/10/2017 3:47 PM **Batch State** Processed

### Analysis Info

**Acq Time** 2017-08-10 09:41 **Data File** Cal 2 - 5ng.d  
**Sample Type** Calibration **Sample Name** Cal 2 - 5ng/mL  
**Dilution** 1 **Acq Method** AM 27 Quant THC 7-2017.m  
**Position** P1-B1 **Sample Info**  
**Inj Vol** -1 **Comment** AM 27 Cannabinoid Confirmation

### Sample Chromatogram



### Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.336	10839	236903	0.0458	4.7297
THC-COOH	THC-COOH-d9	2.446	10718	71674	0.1495	4.6031
THC	THC-d3	6.513	5456	88443	0.0617	5.0094

# ISP FORENSICS - Cd'A Instrument # 62340


## Cannabinoids Analysis Report

**Batch Data Path** D:\2017 Data\8917 cann quant\QuantResults\8917 cann quant.batch.bin  
**Analysis Time** 8/10/2017 3:47 PM **Analyst Name** ISP Tox  
**Report Time** 8/10/2017 3:48 PM **Reporter Name** ISP Tox  
**Last Calib Update** 8/10/2017 3:47 PM **Batch State** Processed

### Analysis Info

**Acq Time** 2017-08-10 09:53 **Data File** Cal 3 - 10ng.d  
**Sample Type** Calibration **Sample Name** Cal 3 - 10ng/mL  
**Dilution** 1 **Acq Method** AM 27 Quant THC 7-2017.m  
**Position** P1-C1 **Sample Info**  
**Inj Vol** -1 **Comment** AM 27 Cannabinoid Confirmation

### Sample Chromatogram



### Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.356	21735	215931	0.1007	10.7390
THC-COOH	THC-COOH-d9	2.486	21256	73613	0.2888	11.4113
THC	THC-d3	6.573	10938	90968	0.1202	10.0999

# ISP FORENSICS - Cd'A Instrument # 62340


## Cannabinoids Analysis Report

**Batch Data Path** D:\2017 Data\8917 cann quant\QuantResults\8917 cann quant.batch.bin  
**Analysis Time** 8/10/2017 3:47 PM **Analyst Name** ISP Tox  
**Report Time** 8/10/2017 3:48 PM **Reporter Name** ISP Tox  
**Last Calib Update** 8/10/2017 3:47 PM **Batch State** Processed

### Analysis Info

**Acq Time** 2017-08-10 10:05 **Data File** Cal 4 - 25ng.d  
**Sample Type** Calibration **Sample Name** Cal 4 - 25ng/mL  
**Dilution** 1 **Acq Method** AM 27 Quant THC 7-2017.m  
**Position** P1-D1 **Sample Info**  
**Inj Vol** -1 **Comment** AM 27 Cannabinoid Confirmation

### Sample Chromatogram



### Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.376	56470	239256	0.2360	25.5543
THC-COOH	THC-COOH-d9	2.506	45992	77313	0.5949	26.3813
THC	THC-d3	6.633	28669	94974	0.3019	25.8909

# ISP FORENSICS - Cd'A Instrument # 62340


## Cannabinoids Analysis Report

**Batch Data Path** D:\2017 Data\8917 cann quant\QuantResults\8917 cann quant.batch.bin  
**Analysis Time** 8/10/2017 3:47 PM **Analyst Name** ISP Tox  
**Report Time** 8/10/2017 3:48 PM **Reporter Name** ISP Tox  
**Last Calib Update** 8/10/2017 3:47 PM **Batch State** Processed

### Analysis Info

**Acq Time** 2017-08-10 10:17 **Data File** Cal 5 - 50ng.d  
**Sample Type** Calibration **Sample Name** Cal 5 - 50ng/mL  
**Dilution** 1 **Acq Method** AM 27 Quant THC 7-2017.m  
**Position** P1-E1 **Sample Info**  
**Inj Vol** -1 **Comment** AM 27 Cannabinoid Confirmation

### Sample Chromatogram



### Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.376	111141	245391	0.4529	49.2931
THC-COOH	THC-COOH-d9	2.506	81575	76973	1.0598	49.1158
THC	THC-d3	6.653	56244	96969	0.5800	50.0748

# ISP FORENSICS - Cd'A Instrument # 62340


## Cannabinoids Analysis Report

**Batch Data Path** D:\2017 Data\8917 cann quant\QuantResults\8917 cann quant.batch.bin  
**Analysis Time** 8/10/2017 3:47 PM **Analyst Name** ISP Tox  
**Report Time** 8/10/2017 3:48 PM **Reporter Name** ISP Tox  
**Last Calib Update** 8/10/2017 3:47 PM **Batch State** Processed

### Analysis Info

**Acq Time** 2017-08-10 10:29 **Data File** Cal 6 - 100ng.d  
**Sample Type** Calibration **Sample Name** Cal 6 - 100ng/mL  
**Dilution** 1 **Acq Method** AM 27 Quant THC 7-2017.m  
**Position** P1-F1 **Sample Info**  
**Inj Vol** -1 **Comment** AM 27 Cannabinoid Confirmation

### Sample Chromatogram



### Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.376	219507	239913	0.9149	99.8623
THC-COOH	THC-COOH-d9	2.506	149959	74959	2.0006	95.1208
THC	THC-d3	6.633	106968	94663	1.1300	97.8910

# ISP FORENSICS - Cd'A Instrument # 62340


## Cannabinoids Analysis Report

**Batch Data Path** D:\2017 Data\8917 cann quant\QuantResults\8917 cann quant.batch.bin  
**Analysis Time** 8/10/2017 3:47 PM **Analyst Name** ISP Tox  
**Report Time** 8/10/2017 3:48 PM **Reporter Name** ISP Tox  
**Last Calib Update** 8/10/2017 3:47 PM **Batch State** Processed

### Analysis Info

**Acq Time** 2017-08-10 10:40 **Data File** Cal 7 - 250ng.d  
**Sample Type** Calibration **Sample Name** Cal 7 - 250ng/mL  
**Dilution** 1 **Acq Method** AM 27 Quant THC 7-2017.m  
**Position** P1-G1 **Sample Info**  
**Inj Vol** -1 **Comment** AM 27 Cannabinoid Confirmation

### Sample Chromatogram



### Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.356	590058	258144	2.2858	249.8999
THC-COOH	THC-COOH-d9	2.486	373860	76296	4.9001	236.9152
THC	THC-d3	6.553	295584	102190	2.8925	251.1312