

Worklist: 1915

Boyle

<u>LAB CASE</u>	<u>ITEM</u>	<u>TASK ID</u>	<u>DESCRIPTION</u>	
C2017-1759	1	96035	AM 27 Blood THC Quant by LC	
C2017-1761	1	96036	AM 27 Blood THC Quant by LC	
C2017-1774	1	96037	AM 27 Blood THC Quant by LC	
C2017-1832	1	96352	AM 27 Blood THC Quant by LC	
C2017-1919	1	96042	AM 27 Blood THC Quant by LC	
M2017-2498	2	96043	AM 27 Blood THC Quant by LC	
M2017-3654	1	96038	AM 27 Blood THC Quant by LC	
M2017-3719	1	96039	AM 27 Blood THC Quant by LC	
M2017-3773	1	96044	AM 27 Blood THC Quant by LC	
M2017-3900	2	96040	AM 27 Blood THC Quant by LC	
M2017-3999	1	96041	AM 27 Blood THC Quant by LC	

[Handwritten signature]

Quantitation of THC and Metabolites in Blood by LC-MS/MS

Extraction Date: 9-26-17

Analyst: Anne Nord

Plate lot#: 0499102

Plate Expiration: 1/29/2018

Mobile phase A: 0.1% Formic Acid in LCMS Water
MTBE

Mobile phase B: 0.1% Formic acid in Acetonitrile
LCMS Methanol Hexane

Blank Blood Lot: 321632-1

Column: UCT Selectra DA 100 x 2.1mm 3um

LCMS-QQQ ID: 62340

Pre-Analytic:

- 1. Check levels of mobile phases and needle wash refill as needed. Ensure waste is not full.
- 2. Ensure correct column is installed and begin mobile phase flow allow to equilibrate ~ 30 minutes.
- 3. Create worklist:

Analytic:

- 1. Remove standards, plate, controls, and samples from cold storage. Allow to reach room temperature.
- 2. Pipette **1000µL blood (calibrated pipette) Pipette ID: 2609543** in wells of analytical (standards) plate.
- 3. Place on shaking incubator at ambient temp., 900rpm for 15 minutes. *Shaker ID: 66759*
- 4. Pipette **500µL 0.1% formic acid in water** in wells of analytical plate.
- 5. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 6. Transfer **800µL of blood+acid** mixture to corresponding wells of SLE+ plate.
- 7. Apply positive pressure for approx. 10-15 seconds (or until no liquid remains on top of sorbent).
(Load at 85-100 PSI- Selector to the right) Manifold ID: 66729
- 8. Wait 5 minutes.
- 9. Add **2.25mL MTBE. (Add in 3 increments of 750uL)**
- 10. Wait 5 minutes.
- 11. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 12. Add **2.25mL Hexane. (Add in 3 increments of 750uL)**
- 13. Wait 5 minutes.
- 14. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 15. Remove plate containing eluate. Place on SPE Dry and evaporate to dryness at approx. 35°C.
SPE Dry ID: 66819
- 16. Reconstitute in **100µL 100% MeOH** and heat seal plate with foil. Place in autosampler and run worklist.

Post-Analytic

- 1. Create batch and process data.
Worklist path: 92617 cann quant Batch Name: 92617 cann quant
- 2. Make any necessary integration changes, r^2 values ≥ 0.98 for each analyte.
- 3. Did all QCs pass for each analyte? Y/N Enter QCs into control charting?
- 4. Central File Packet to include: LIMS Worklist, Method Checklist, Calibration and Control Reports

COMMENTS: [Click here to enter text.](#)

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2017 Data\92617 cann quant\QuantResults\92617 cann quant.batch.bin
Analysis Time 9/28/2017 9:52 AM **Analyst Name** ISP Tox
Report Time 9/28/2017 9:54 AM **Reporter Name** ISP Tox
Last Calib Update 9/28/2017 9:52 AM **Batch State** Processed

Analysis Info

Acq Time 2017-09-27 13:07 **Data File** Negative Control.d
Sample Type Sample **Sample Name** Negative Control
Dilution 1 **Acq Method** AM 27 Quant THC 7-2017.m
Position P1-A2 **Sample Info**
Inj Vol -1 **Comment** AM 27 Cannabinoid Confirmation

Sample Chromatogram



ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2017 Data\92617 cann quant\QuantResults\92617 cann quant.batch.bin
Analysis Time 9/28/2017 9:52 AM **Analyst Name** ISP Tox
Report Time 9/28/2017 9:54 AM **Reporter Name** ISP Tox
Last Calib Update 9/28/2017 9:52 AM **Batch State** Processed

Analysis Info

Acq Time 2017-09-27 13:19 **Data File** QC - 10ng.d
Sample Type QC **Sample Name** QC - 10ng
Dilution 1 **Acq Method** AM 27 Quant THC 7-2017.m
Position P1-H1 **Sample Info**
Inj Vol -1 **Comment** AM 27 Cannabinoid Confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.456	19274	203142	0.0949	8.9572
THC-COOH	THC-COOH-d9	2.626	14384	74586	0.1929	8.4252
THC	THC-d3	6.913	10541	84505	0.1247	9.5295

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2017 Data\92617 cann quant\QuantResults\92617 cann quant.batch.bin
Analysis Time 9/28/2017 9:52 AM **Analyst Name** ISP Tox
Report Time 9/28/2017 9:54 AM **Reporter Name** ISP Tox
Last Calib Update 9/28/2017 9:52 AM **Batch State** Processed

Analysis Info

Acq Time 2017-09-27 13:31 **Data File** External Control lot 21718.d
Sample Type Sample **Sample Name** External Control lot 21718
Dilution 1 **Acq Method** AM 27 Quant THC 7-2017.m
Position P1-B2 **Sample Info**
Inj Vol -1 **Comment** AM 27 Cannabinoid Confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.456	14836	133669	0.1110	10.5661
THC-COOH	THC-COOH-d9	2.626	11620	44397	0.2617	11.7636
THC	THC-d3	6.933	11396	71999	0.1583	12.2422

ISP Forensics Calibration Curve Report


Batch Data Path D:\2017 Data\92617 cann quant\QuantResults\92617 cann quant.batch.bin

Last Calib Update 9/28/2017 9:52 AM

Analyst Name ISP TOX

Target Compound *THC-OH*

Internal Standard *THC-OH-d3*




Sample	Level	Enabled	Exp Conc	Final Conc	Accuracy
Cal 1 - 3ng	1	<input checked="" type="checkbox"/>	3	3.0	100.6
Cal 2 - 5ng	2	<input checked="" type="checkbox"/>	5	4.9	98.8
Cal 3 - 10ng	3	<input checked="" type="checkbox"/>	10	9.3	93.1
QC - 10ng	3	<input checked="" type="checkbox"/>	10	9.0	89.6
Cal 4 - 25ng	4	<input checked="" type="checkbox"/>	25	26.2	104.7
Cal 5 - 50ng	5	<input checked="" type="checkbox"/>	50	50.9	101.7
Cal 6 - 100ng	6	<input checked="" type="checkbox"/>	100	102.8	102.8
Cal 7 - 250ng	7	<input checked="" type="checkbox"/>	250	245.9	98.4

ISP Forensics Calibration Curve Report

Batch Data Path D:\2017 Data\92617 cann quant\QuantResults\92617 cann quant.batch.bin

Last Calib Update 9/28/2017 9:52 AM **Analyst Name** ISP TOX

Target Compound *THC-COOH*
Internal Standard *THC-COOH-d9*




Sample	Level	Enabled	Exp Conc	Final Conc	Accuracy
Cal 1 - 3ng	1	<input checked="" type="checkbox"/>	3	3.2	107.2
Cal 2 - 5ng	2	<input checked="" type="checkbox"/>	5	4.5	90.7
Cal 3 - 10ng	3	<input checked="" type="checkbox"/>	10	9.3	92.7
QC - 10ng	3	<input checked="" type="checkbox"/>	10	8.4	84.3
Cal 4 - 25ng	4	<input checked="" type="checkbox"/>	25	25.4	101.5
Cal 5 - 50ng	5	<input checked="" type="checkbox"/>	50	51.4	102.9
Cal 6 - 100ng	6	<input checked="" type="checkbox"/>	100	105.4	105.4
Cal 7 - 250ng	7	<input checked="" type="checkbox"/>	250	249.0	99.6

ISP Forensics Calibration Curve Report

Batch Data Path D:\2017 Data\92617 cann quant\QuantResults\92617 cann quant.batch.bin

Last Calib Update 9/28/2017 9:52 AM **Analyst Name** ISP TOX

Target Compound *THC*
Internal Standard *THC-d3*



Sample	Level	Enabled	Exp Conc	Final Conc	Accuracy
Cal 1 - 3ng	1	<input checked="" type="checkbox"/>	3	2.5	82.0
Cal 2 - 5ng	2	<input checked="" type="checkbox"/>	5	5.5	110.4
Cal 3 - 10ng	3	<input checked="" type="checkbox"/>	10	9.9	98.8
QC - 10ng	3	<input checked="" type="checkbox"/>	10	9.5	95.3
Cal 4 - 25ng	4	<input checked="" type="checkbox"/>	25	25.5	102.0
Cal 5 - 50ng	5	<input checked="" type="checkbox"/>	50	52.6	105.3
Cal 6 - 100ng	6	<input checked="" type="checkbox"/>	100	104.4	104.4
Cal 7 - 250ng	7	<input checked="" type="checkbox"/>	250	242.6	97.0

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2017 Data\92617 cann quant\QuantResults\92617 cann quant.batch.bin
Analysis Time 9/28/2017 9:52 AM **Analyst Name** ISP Tox
Report Time 9/28/2017 9:54 AM **Reporter Name** ISP Tox
Last Calib Update 9/28/2017 9:52 AM **Batch State** Processed

Analysis Info

Acq Time 2017-09-27 11:32 **Data File** Cal 1 - 3ng.d
Sample Type Calibration **Sample Name** Cal 1 - 3ng
Dilution 1 **Acq Method** AM 27 Quant THC 7-2017.m
Position P1-A1 **Sample Info**
Inj Vol -1 **Comment** AM 27 Cannabinoid Confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.396	3452	97526	0.0354	3.0186
THC-COOH	THC-COOH-d9	2.506	3281	38417	0.0854	3.2171
THC	THC-d3	6.733	1097	29405	0.0373	2.4602

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2017 Data\92617 cann quant\QuantResults\92617 cann quant.batch.bin
Analysis Time 9/28/2017 9:52 AM **Analyst Name** ISP Tox
Report Time 9/28/2017 9:54 AM **Reporter Name** ISP Tox
Last Calib Update 9/28/2017 9:52 AM **Batch State** Processed

Analysis Info

Acq Time 2017-09-27 11:44 **Data File** Cal 2 - 5ng.d
Sample Type Calibration **Sample Name** Cal 2 - 5ng
Dilution 1 **Acq Method** AM 27 Quant THC 7-2017.m
Position P1-B1 **Sample Info**
Inj Vol -1 **Comment** AM 27 Cannabinoid Confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.396	4795	87786	0.0546	4.9383
THC-COOH	THC-COOH-d9	2.526	3855	34222	0.1126	4.5372
THC	THC-d3	6.753	2212	29420	0.0752	5.5220

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2017 Data\92617 cann quant\QuantResults\92617 cann quant.batch.bin
Analysis Time 9/28/2017 9:52 AM **Analyst Name** ISP Tox
Report Time 9/28/2017 9:54 AM **Reporter Name** ISP Tox
Last Calib Update 9/28/2017 9:52 AM **Batch State** Processed

Analysis Info

Acq Time 2017-09-27 11:56 **Data File** Cal 3 - 10ngr.d
Sample Type Calibration **Sample Name** Cal 3 - 10ng
Dilution 1 **Acq Method** AM 27 Quant THC 7-2017.m
Position P1-C1 **Sample Info**
Inj Vol -1 **Comment** AM 27 Cannabinoid Confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.376	9133	92839	0.0984	9.3059
THC-COOH	THC-COOH-d9	2.526	7283	34633	0.2103	9.2698
THC	THC-d3	6.733	3972	30761	0.1291	9.8837

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2017 Data\92617 cann quant\QuantResults\92617 cann quant.batch.bin
Analysis Time 9/28/2017 9:52 AM **Analyst Name** ISP Tox
Report Time 9/28/2017 9:54 AM **Reporter Name** ISP Tox
Last Calib Update 9/28/2017 9:52 AM **Batch State** Processed

Analysis Info

Acq Time 2017-09-27 12:08 **Data File** Cal 4 - 25ng.d
Sample Type Calibration **Sample Name** Cal 4 - 25ng
Dilution 1 **Acq Method** AM 27 Quant THC 7-2017.m
Position P1-D1 **Sample Info**
Inj Vol -1 **Comment** AM 27 Cannabinoid Confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.396	27189	101704	0.2673	26.1764
THC-COOH	THC-COOH-d9	2.526	20829	38394	0.5425	25.3729
THC	THC-d3	6.753	11449	35546	0.3221	25.4889

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2017 Data\92617 cann quant\QuantResults\92617 cann quant.batch.bin
Analysis Time 9/28/2017 9:52 AM **Analyst Name** ISP Tox
Report Time 9/28/2017 9:54 AM **Reporter Name** ISP Tox
Last Calib Update 9/28/2017 9:52 AM **Batch State** Processed

Analysis Info

Acq Time 2017-09-27 12:20 **Data File** Cal 5 - 50ng.d
Sample Type Calibration **Sample Name** Cal 5 - 50ng
Dilution 1 **Acq Method** AM 27 Quant THC 7-2017.m
Position P1-E1 **Sample Info**
Inj Vol -1 **Comment** AM 27 Cannabinoid Confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.376	83367	162029	0.5145	50.8564
THC-COOH	THC-COOH-d9	2.526	61669	57083	1.0803	51.4396
THC	THC-d3	6.733	37876	57565	0.6580	52.6492

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2017 Data\92617 cann quant\QuantResults\92617 cann quant.batch.bin
Analysis Time 9/28/2017 9:52 AM **Analyst Name** ISP Tox
Report Time 9/28/2017 9:54 AM **Reporter Name** ISP Tox
Last Calib Update 9/28/2017 9:52 AM **Batch State** Processed

Analysis Info

Acq Time 2017-09-27 12:32 **Data File** Cal 6 - 100ng.d
Sample Type Calibration **Sample Name** Cal 6 - 100ng
Dilution 1 **Acq Method** AM 27 Quant THC 7-2017.m
Position P1-F1 **Sample Info**
Inj Vol -1 **Comment** AM 27 Cannabinoid Confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.396	131632	127256	1.0344	102.7623
THC-COOH	THC-COOH-d9	2.546	100049	45626	2.1928	105.3596
THC	THC-d3	6.773	59294	45666	1.2984	104.4413

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2017 Data\92617 cann quant\QuantResults\92617 cann quant.batch.bin
Analysis Time 9/28/2017 9:52 AM **Analyst Name** ISP Tox
Report Time 9/28/2017 9:54 AM **Reporter Name** ISP Tox
Last Calib Update 9/28/2017 9:52 AM **Batch State** Processed

Analysis Info

Acq Time 2017-09-27 12:43 **Data File** Cal 7 - 250ng.d
Sample Type Calibration **Sample Name** Cal 7 - 250ng
Dilution 1 **Acq Method** AM 27 Quant THC 7-2017.m
Position P1-G1 **Sample Info**
Inj Vol -1 **Comment** AM 27 Cannabinoid Confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.396	298248	120826	2.4684	245.9422
THC-COOH	THC-COOH-d9	2.546	219345	42543	5.1559	248.9758
THC	THC-d3	6.793	134655	44790	3.0064	242.5547

Quantitation of THC and Metabolites in Blood by LC-MS/MS

Extraction Date: 9-13-17

Analyst: Anne Nard

Plate lot#: 0499102

Plate Expiration: 1/29/2018

Mobile phase A: 0.1% Formic Acid in LCMS Water
MTBE

Mobile phase B: 0.1% Formic acid in Acetonitrile
LCMS Methanol Hexane

Blank Blood Lot: 321632-1

Column: UCT Selectra DA 100 x 2.1mm 3um

LCMS-QQQ ID: 62340

Pre-Analytic:

- 1. Check levels of mobile phases and needle wash refill as needed. Ensure waste is not full.
- 2. Ensure correct column is installed and begin mobile phase flow allow to equilibrate ~ 30 minutes.
- 3. Create worklist:

Analytic:

- 1. Remove standards, plate, controls, and samples from cold storage. Allow to reach room temperature.
- 2. Pipette **1000µL blood (calibrated pipette) Pipette ID: 2609543** in wells of analytical (standards) plate.
- 3. Place on shaking incubator at ambient temp., 900rpm for 15 minutes. *Shaker ID: 66759*
- 4. Pipette **500µL 0.1% formic acid in water** in wells of analytical plate.
- 5. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 6. Transfer **800µL of blood+acid** mixture to corresponding wells of SLE+ plate.
- 7. Apply positive pressure for approx. 10-15 seconds (or until no liquid remains on top of sorbent).
(Load at 85-100 PSI- Selector to the right) Manifold ID: 66729
- 8. Wait 5 minutes.
- 9. Add **2.25mL MTBE. (Add in 3 increments of 750uL)**
- 10. Wait 5 minutes.
- 11. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 12. Add **2.25mL Hexane. (Add in 3 increments of 750uL)**
- 13. Wait 5 minutes.
- 14. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 15. Remove plate containing eluate. Place on SPE Dry and evaporate to dryness at approx. 35°C.
SPE Dry ID: 66819
- 16. Reconstitute in **100µL 100% MeOH** and heat seal plate with foil. Place in autosampler and run worklist.

Post-Analytic

- 1. Create batch and process data.
Worklist path: 9/13/17 cann quant Batch Name: 091317 cann quant
- 2. Make any necessary integration changes, r^2 values ≥ 0.98 for each analyte. — fail
- 3. Did all QCs pass for each analyte? Y (N) Enter QCs into control charting?
- 4. Central File Packet to include: LIMS Worklist, Method Checklist, Calibration and Control Reports

COMMENTS: [Click here to enter text.](#)

This run failed. Samples were re-extracted and run with the next batch.

C2017-1759-1 m2017-3654-1 m2017-3999-1
C2017-1761-1 m2017-3719-1
C2017-1774-1 m2017-3900-2

Sample		THC-OH Results										Qualifier (331.2 -> 193.0)...				THC-OH-d3 (I.					
✓	Ⓢ	Name	Type	L e v el	Acq. Date-Time	Pos.	Exp. Conc.	RT	Resp.	S/N	Ml	Calc. Conc.	Accuracy	Transition	Ratio	Ml	Area	S/N	RT	Resp.	
✓	Ⓢ	Cal 1 - 3ng	Cal	1	9/15/2017 6:42 P...	P1-A1	3.00...	2.295	6582	7.91	1.6397	54.7	331.2 -> 193...								
		Cal 2 - 5ng	Cal	2	9/15/2017 6:54 P...	P1-B1	5.00...	2.315	10173	18.40	4.1359	82.7	331.2 -> 193...		12.4		1259	20.28	2.292	2.292	1867.
		Cal 3 - 10ng	Cal	3	9/16/2017 9:08 A...	P1-C1	10.0...						331.2 -> 193...								
✓		Cal 4 - 25ng	Cal	4	9/16/2017 9:20 A...	P1-D1	25.0...	2.315	33450	109.57	16.6161	66.5	331.2 -> 193...		13.1		4378	36.15	2.292	2.292	1871.
✓		Cal 5 - 50ng	Cal	5	9/16/2017 9:31 A...	P1-E1	50.0...	2.315	1962	Infinity	113.6681	227.3	331.2 -> 193...		12.4		24395	287.88	2.312	2.312	1713.
		Cal 6 - 100ng	Cal	6	9/16/2017 9:43 A...	P1-F1	100...	2.315	1631...	501.99	83.4133	83.4	331.2 -> 193...		12.5		20353	235.49	2.292	2.292	1932.
		Cal 7 - 250ng	Cal	7	9/16/2017 9:55 A...	P1-G1	250...	2.315	4172...	Infinity	213.5268	85.4	331.2 -> 193...		12.2		50752	335.27	2.292	2.292	1950.
		Negative Control	Sa...		9/16/2017 10:19...	P1-A2							331.2 -> 193...								
✓		QC - 10ng	QC	3	9/16/2017 10:31...	P1-H1	10.0...	2.336	25385	40.85	16.1562	161.6	331.2 -> 193...		13.0		3302	53.80	2.332	2.332	1457.
		External Control lot 217...	Sa...		9/16/2017 10:43...	P1-B2		2.315	22288	28.07	5.7621		331.2 -> 193...		13.3		2964	23.13	2.312	2.312	3153.
		blank c2017-1759-1	Sa...		9/16/2017 10:54...	Vial 2							331.2 -> 193...						2.292	2.292	1397.
		C2017-1759-1	Sa...		9/16/2017 11:06...	P1-C2							331.2 -> 193...						2.292	2.292	2377.
		blank c2017-1761-1	Sa...		9/16/2017 11:18...	Vial 2							331.2 -> 193...						2.292	2.292	1410.
		c2017-1761-1	Sa...		9/16/2017 11:30...	P1-D2		2.315	30294	143.34	10.9961		331.2 -> 193...		12.4		3752	13.89	2.292	2.292	2467.
		blank c2017-1774-1	Sa...		9/16/2017 11:42...	Vial 2							331.2 -> 193...						2.352	2.352	1355.
		c2017-1774-1	Sa...		9/16/2017 11:54...	P1-E2		2.356	16362	17.91	3.8619		331.2 -> 193...		12.9		2104	15.58	2.352	2.352	3161.
		blank m2017-3654-1	Sa...		9/16/2017 12:06...	Vial 2							331.2 -> 193...						2.352	2.352	1442.
		m2017-3654-1	Sa...		9/16/2017 12:17...	P1-F2		2.356	1466	3.78	0.0000		331.2 -> 193...						2.352	2.352	2058.
		blank m2017-3719-1	Sa...		9/16/2017 12:29...	Vial 2							331.2 -> 193...						2.352	2.352	1373.
		m2017-3719-1	Sa...		9/16/2017 12:41...	P1-G2		2.356	3177	7.09	0.0973		331.2 -> 193...						2.332	2.332	2227.
		blank m2017-3900-2	Sa...		9/16/2017 12:53...	Vial 2							331.2 -> 193...						2.352	2.352	1368.
		m2017-3900-2	Sa...		9/16/2017 1:05 P...	P1-H2		2.336	11940	95.93	3.8448		331.2 -> 193...		13.2		1573	13.60	2.352	2.352	2314.
		blank m2017-3999-1	Sa...		9/16/2017 1:17 P...	Vial 2							331.2 -> 193...						2.352	2.352	1398.
		m2017-3999-1	Sa...		9/16/2017 1:28 P...	P1-A3							331.2 -> 193...								

Sample	THC...				THC Results							Qualifier (315.2 -> 12...				THC-d3 (IST ...		
	Type	Name	Level	Acq. Date-Time	Pos.	Exp. Conc.	RT	Resp.	S/N	MI	Calc. Conc.	Accuracy	Transition	Ratio	MI	Area	S/N	RT
☉	Cal	Cal 1 - 3ng	1	9/15/2017 6:42...	P1-A1	3.0...	6.312	3248	11.73	☐	1.5046	50.2	315.2 -> 123...	30.2	982	4...	6.306	85142
	Cal	Cal 2 - 5ng	2	9/15/2017 6:54...	P1-B1	5.0...	6.312	5087	30.89	☐	4.1841	83.7	315.2 -> 123...	31.6	1609	3...	6.266	71164
☉	Cal	Cal 3 - 10ng	3	9/16/2017 9:08...	P1-C1	10...				☐			315.2 -> 123...					
☉	Cal	Cal 4 - 25ng	4	9/16/2017 9:20...	P1-D1	25...	6.353	19250	152.56	☐	17.1594	68.6	315.2 -> 123...	28.5	5490	1...	6.346	82656
☉	Cal	Cal 5 - 50ng	5	9/16/2017 9:31...	P1-E1	50...	6.373	1031...	343.49	☐	114.9119	229.8	315.2 -> 123...	33.8	34843	2...	6.346	71221
	Cal	Cal 6 - 100ng	6	9/16/2017 9:43...	P1-F1	100...	6.353	87478	Infinity	☐	82.6775	82.7	315.2 -> 123...	34.1	29806	2...	6.346	83476
	Cal	Cal 7 - 250ng	7	9/16/2017 9:55...	P1-G1	250...	6.353	2183...	Infinity	☐	212.5625	85.0	315.2 -> 123...	33.5	73171	In...	6.346	81957
☉	Sa...	Negative Control		9/16/2017 10:1...	P1-A2					☐			315.2 -> 123...					
☉	QC	QC - 10ng	3	9/16/2017 10:3...	P1-H1	10...	6.433	13909	60.79	☐	16.7848	167.8	315.2 -> 123...	31.4	4373	1...	6.426	60944
	Sa...	External Control...		9/16/2017 10:4...	P1-B2		6.353	14140	39.89	☐	6.3900		315.2 -> 123...	29.2	4122	1...	6.346	1429...
☉	Sa...	blank c2017-17...		9/16/2017 10:5...	Vial 2					☐			315.2 -> 123...				6.226	72497
☉	Sa...	C2017-1759-1		9/16/2017 11:0...	P1-C2		6.373	1167	7.94	☐	0.0000		315.2 -> 123...				6.326	96423
☉	Sa...	blank c2017-176...		9/16/2017 11:1...	Vial 2					☐			315.2 -> 123...				6.206	70839
☉	Sa...	c2017-1761-1		9/16/2017 11:3...	P1-D2		6.272	13783	21.89	☐	17.4746		315.2 -> 123...	36.4	5014	3...	6.226	58204
☉	Sa...	blank c2017-177...		9/16/2017 11:4...	Vial 2					☐			315.2 -> 123...				6.406	69694
☉	Sa...	c2017-1774-1		9/16/2017 11:5...	P1-E2		6.473	15575	96.82	☐	12.7043		315.2 -> 123...	32.4	5043	1...	6.406	87764
☉	Sa...	blank m2017-36...		9/16/2017 12:0...	Vial 2					☐			315.2 -> 123...				6.406	69595
☉	Sa...	m2017-3654-1		9/16/2017 12:1...	P1-F2		6.453	1178	5.67	☐	0.0000		315.2 -> 123...				6.426	64390
☉	Sa...	blank m2017-37...		9/16/2017 12:2...	Vial 2					☐			315.2 -> 123...				6.406	69306
☉	Sa...	m2017-3719-1		9/16/2017 12:4...	P1-G2		6.433	1998	9.69	☐	0.9788		315.2 -> 123...	34.0	680	2...	6.426	63209
☉	Sa...	blank m2017-39...		9/16/2017 12:5...	Vial 2					☐			315.2 -> 123...				6.426	68660
☉	Sa...	m2017-3900-2		9/16/2017 1:05...	P1-H2		6.553	5585	15.91	☐	9.5329		315.2 -> 123...	26.9	1502	1...	6.406	40465
☉	Sa...	blank m2017-39...		9/16/2017 1:17...	Vial 2					☐			315.2 -> 123...				6.406	67994
☉	Sa...	m2017-3999-1		9/16/2017 1:28...	P1-A3					☐			315.2 -> 123...					