


Worklist: 2555

<u>LAB_CASE</u>	<u>ITEM</u>	<u>TASK_ID</u>	<u>DESCRIPTION</u>
M2018-2977	1	120819	AM 27 Blood THC Quant by LC-QQQ
M2018-3061	1	120820	AM 27 Blood THC Quant by LC-QQQ
P2018-1849	1	120821	AM 27 Blood THC Quant by LC-QQQ
P2018-1873	1	120822	AM 27 Blood THC Quant by LC-QQQ
P2018-1875	1	120823	AM 27 Blood THC Quant by LC-QQQ
P2018-1890	1	120824	AM 27 Blood THC Quant by LC-QQQ
P2018-1892	1	120825	AM 27 Blood THC Quant by LC-QQQ
P2018-1894	1	120826	AM 27 Blood THC Quant by LC-QQQ



AM# 27: Quantitation of THC and Metabolites in Blood by LC-MS/MS

PS

Extraction Date: 07/10/18
Plate lot#: 0515037

Analyst: Sarah Pickle
Plate Expiration: 09/28/18

Mobile phase A: 0.1% Formic Acid in LCMS Water
MTBE LCMS Methanol

Mobile phase B: 0.1% Formic acid in Acetonitrile
Hexane

Blank Blood Lot: 361331-1

Column: UCT Selectra DA 100 x 2.1mm 3um

Blank Urine: POC062718

LCMS-QQQ ID: 59740

Pre-Analytic:

- 1. Check levels of mobile phases and needle wash refill as needed. Ensure waste is not full.
- 2. Ensure correct column is installed and begin mobile phase flow allow to equilibrate ~ 30 minutes.
- 3. Create worklist: Data Path: _____

Urine Hydrolysis:

- 1. Pipette 1.5 mL urine into empty 48 well plate.
- 2. Add 250ul 1N KOH to urine samples.
- 3. Place on shaking incubator at 40 degrees for 15 mins.

Analytic:

- 1. Remove standards, plate, controls, and samples from cold storage. Allow to reach room temperature.
- 2. Pipette **1000µL blood/urine (calibrated pipette) Pipette ID: 3382167** in wells of analytical (standards) plate.
- 3. Place on shaking incubator at ambient temp., 900rpm for 15 minutes. *Shaker ID: 067105*
- 4. Pipette **500µL 0.1% formic acid in water** in wells of analytical plate for blood samples. Pipette **500ul saturated phosphate buffer** for urine samples.
- 5. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 6. Transfer **800µL of blood+acid** mixture to corresponding wells of SLE+ plate.
- 7. Apply positive pressure for approx. 10-15 seconds (or until no liquid remains on top of sorbent).
(Load at 85-100 PSI- Selector to the right) Manifold ID: 067104
- 8. Wait 5 minutes.
- 9. Add **2.25mL MTBE. (Add in 3 increments of 750uL)**
- 10. Wait 5 minutes.
- 11. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 12. Add **2.25mL Hexane. (Add in 3 increments of 750uL)**
- 13. Wait 5 minutes.
- 14. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 15. Remove plate containing eluate. Place on SPE Dry and evaporate to dryness at approx. 35°C.
SPE Dry ID: 067103
- 16. Reconstitute in **100µL 100% MeOH** and heat seal plate with foil. Place in autosampler and run worklist.

Post-Analytic

- 1. Create batch and process data.
Worklist path: 071018 THC Quant Urine Validation SP worklist 2555
Batch Name: 071018 THCQ SP
- 2. Make any necessary integration changes, Curve weighting of Linear 1/x with r^2 values ≥ 0.98 for each analyte
- 3. RT +/- 3% or 0.100 min, whichever is greater, +/- 20% Accuracy for greater than (+/- 30% for 10ng/ml or less).
Ion ratios must be within +/- 20% of the averaged calibrators
- 4. Case sample response for THC and OH-THC 3ng/mL (quantitative), Carboxy-THC: 10ng/mL (qualitative only) will be reported. Samples with a THC or OH-THC response over 50 ng/mL will be reported out as greater than 50 ng/mL.

15

- 5. Did all QCs pass for each analyte? Y / N
- 6. Enter QCs into control charting.
- 7. Central File Packet to include: LIMS Worklist, Method Checklist, Calibration and Control Reports

COMMENTS: *Curve Range Limited: THC-COOH 5-250*



Idaho State Police Forensic Services

P

AM #27 Quantitative Analysis of THC and Metabolites in Blood by LCMS-QQQ

Analyst: Sarah Pickle
Extraction Date: 07/10/18
Worklist Number: 2555

<i>Reagent</i>	<i>Lot Number</i>	<i>Expiration Date</i>	<i>Date in Service</i>	<i>Date Out of Service</i>	<i>Initials</i>
ToxBox THC/THC Metabolite Plate	0515037	09/28/18			
Negative Blood	361331-1		05/25/18		
Negative Urine	POC062718		06/27/18		
Methanol External Control Solution	WS020718	02/07/19	02/07/18		
Blood External Control Solution	061718	06/17/19	06/17/18		
Urine External Controls	062818	06/28/19	06/28/18		
Methyl Tert-Butyl Ether (MTBE) 99.9%	A0375555		06/26/17		
Hexanes (ACS)	101642		10/26/17		
Methanol (LCMS Grade)	177145		04/11/18		
1 N KOH	091817		09/18/18		
Saturated Phosphate Buffer pH 1.8	020118		02/01/18		
0.1% Formic Acid in Water (Mobile Phase A)	166541		06/26/17		
0.1% Formic Acid in Acetonitrile (Mobile Phase B)	176190		02/08/18		
Needle Rinse--75% LCMS MeOH in LCMS Water	070518		07/05/18		

Methanol External Control Solution (Lot: WS020718)
10 ul of 1mg/mL THC, 100 ul of 100 ug/mL THC-OH, C-THC in 9790 ul MeOH

<i>Component</i>	<i>Source</i>	<i>Source Lot Number</i>	<i>Expiration Date</i>
Methanol (LCMS)	Fisher	172516	
THC	Cerilliant	FE04231406	04/30/2019
C-THC	Cayman	0497429	02/08/2019
THC-OH	Cerilliant	FE01121503	01/31/2020
Prepared:	02/07/18		
Prepared By:	Tamara Salazar		
Expires:	02/07/19		

Blood External Control Solution (Lot: 061718)
100 ul of methanol external control solution was added to 9900 ul of blood.

<i>Component</i>	<i>Source</i>	<i>Source Lot Number</i>
Negative Blood	Hemostat	361331-1
Methanol External Control Solution		WS020718

AM #27: Quantitative Analysis of THC and Metabolites in Blood by LCMS-QQQ

AS

Prepared:	06/17/18
Prepared by:	Tamara Salazar
Expires:	06/17/19

Urine External Control Solutions (Lot: 062818)

Three Controls were made by adding 50 ul, 100 ul, and 200 ul of methanol external control to 4.5 mL of negative urine.

Component	Source	Source Lot Number
Negative Urine		POC062718
Methanol External Control Solution		WS020718
Prepared:	06/28/18	
Prepared by:	Sarah Pickle	
Expires:	06/28/19	

1 N KOH (Lot: 091817)

Component	Source	Source Lot Number
Potassium Hydroxide	Fisher	034727
DI Water	-	-
Prepared:	09/18/17	
Prepared By:	Celena Shrum	

Saturated Phosphate Buffer pH 1.8 (Lot: 020118)

Component	Source	Source Lot Number
Potassium Phosphate monobasic	Fisher	L10-021-61
DI Water	-	-
Prepared:	02/01/18	
Prepared By:	Celena Shrum	

Needle Rinse (75% LCMS MeOH in LCMS Water) (Lot: 070518)

Component	Source	Source Lot Number
MeOH (LCMS Grade)	Fisher	177145
Water (LCMS Grade)	Fisher	177528
Prepared:	07/05/18	
Prepared By:	Tamara Salazar	

ISP FORENSICS - Pocatello Instrument # 59740


Cannabinoids Analysis Report

Batch Data Path C:\MassHunter\Data\2018\THC Quant\071018 THC Quant Urine Validation SP worklist 2555\QuantResults\071018
Analysis Time 7/11/2018 7:39 AM **Analyst Name** ISPUser
Report Time 7/11/2018 7:41 AM **Reporter Name** ISPUser
Last Calib Update 7/11/2018 7:39 AM **Batch State** Processed

Analysis Info

Acq Time 2018-07-10 13:16 **Data File** QC Control.d
Sample Type Sample **Sample Name** QC Control
Dilution 1 **Acq Method** THC Quant 051517 workingmm.m
Position P1-A6 **Sample Info**
Inj Vol -1 **Comment**

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-D3	2.319	140200	1578708	0.0888	8.9710
THC-COOH	THC-COOH-D9	2.432	97725	415636	0.2351	10.1395
THC	THC-D3	6.199	44720	439788	0.1017	9.5897

ISP FORENSICS - Pocatello Instrument # 59740

Cannabinoids Analysis Report


P

Batch Data Path C:\MassHunter\Data\2018\THC Quant\071018 THC Quant Urine Validation SP worklist 2555\QuantResults\071
Analysis Time 7/11/2018 7:39 AM **Analyst Name** ISPUser
Report Time 7/11/2018 7:41 AM **Reporter Name** ISPUser
Last Calib Update 7/11/2018 7:39 AM **Batch State** Processed

Analysis Info

Acq Time 2018-07-10 13:39 **Data File** Negative.d
Sample Type Sample **Sample Name** Negative
Dilution 1 **Acq Method** THC Quant 051517 workingmm.m
Position P1-H5 **Sample Info**
Inj Vol -1 **Comment** Hemostat 361331-1

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-COOH	THC-COOH-D9	2.419	21044	362106	0.0581	2.4494 <10 μ

SP

ISP FORENSICS - Pocatello Instrument # 59740


Cannabinoids Analysis Report

Batch Data Path C:\MassHunter\Data\2018\THC Quant\071018 THC Quant Urine Validation SP worklist 2555\QuantResults\071018
Analysis Time 7/11/2018 7:48 AM **Analyst Name** ISPUser
Report Time 7/11/2018 7:49 AM **Reporter Name** ISPUser
Last Calib Update 7/11/2018 7:48 AM **Batch State** Processed

Analysis Info

Acq Time 2018-07-10 14:03 **Data File** External Control.d
Sample Type Sample **Sample Name** External Control
Dilution 1 **Acq Method** THC Quant 051517 workingmm.m
Position P1-G5 **Sample Info**
Inj Vol -1 **Comment** Hemostat 361331-1 + WS 020718

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-D3	2.332	110760	1330180	0.0833	8.4843
THC-COOH	THC-COOH-D9	2.432	81199	359728	0.2257	9.7312
THC	THC-D3	6.199	34786	359985	0.0966	9.1661


ISP Forensics Calibration Curve Report

S

Batch Data Path C:\MassHunter\Data\2018\THC Quant\071018 THC Quant Urine Validation SP worklist
2555\QuantResults\071018 THCQ SP.batch.bin

Last Calib Update 7/11/2018 7:39 AM **Analyst Name** ISP TOX

Target Compound *THC-OH*
Internal Standard *THC-OH-D3*




Sample	Level	Enabled	Exp Conc	Final Conc	Accuracy
Cal 1	1	<input checked="" type="checkbox"/>	3	3.5	115.7
Cal 2	2	<input checked="" type="checkbox"/>	5	5.0	100.4
Cal 3	3	<input checked="" type="checkbox"/>	10	10.0	100.1
Cal 4	4	<input checked="" type="checkbox"/>	25	21.0	83.9
Cal 5	5	<input checked="" type="checkbox"/>	50	48.7	97.3
Cal 6	6	<input checked="" type="checkbox"/>	100	101.1	101.1
Cal 7	7	<input checked="" type="checkbox"/>	250	253.8	101.5

ISP Forensics Calibration Curve Report

Batch Data Path C:\MassHunter\Data\2018\THC Quant\071018 THC Quant Urine Validation SP worklist
2555\QuantResults\071018 THCQ SP.batch.bin

Last Calib Update 7/11/2018 7:39 AM **Analyst Name** ISP TOX

Target Compound *THC-COOH*
Internal Standard *THC-COOH-D9*



Sample	Level	Enabled	Exp Conc	Final Conc	Accuracy
Cal 1	1	<input type="checkbox"/>	3	4.4	147.7
Cal 2	2	<input checked="" type="checkbox"/>	5	5.4	107.5
Cal 3	3	<input checked="" type="checkbox"/>	10	11.1	111.2
Cal 4	4	<input checked="" type="checkbox"/>	25	20.7	82.8
Cal 5	5	<input checked="" type="checkbox"/>	50	48.3	96.6
Cal 6	6	<input checked="" type="checkbox"/>	100	100.2	100.2
Cal 7	7	<input checked="" type="checkbox"/>	250	254.3	101.7


ISP Forensics Calibration Curve Report

D

Batch Data Path C:\MassHunter\Data\2018\THC Quant\071018 THC Quant Urine Validation SP worklist
2555\QuantResults\071018 THCQ SP.batch.bin

Last Calib Update 7/11/2018 7:39 AM **Analyst Name** ISP TOX

Target Compound *THC*
Internal Standard *THC-D3*



Sample	Level	Enabled	Exp Conc	Final Conc	Accuracy
Cal 1	1	<input checked="" type="checkbox"/>	3	3.6	120.4
Cal 2	2	<input checked="" type="checkbox"/>	5	5.3	105.4
Cal 3	3	<input checked="" type="checkbox"/>	10	9.4	94.3
Cal 4	4	<input checked="" type="checkbox"/>	25	20.7	82.6
Cal 5	5	<input checked="" type="checkbox"/>	50	47.1	94.2
Cal 6	6	<input checked="" type="checkbox"/>	100	100.5	100.5
Cal 7	7	<input checked="" type="checkbox"/>	250	256.5	102.6

ISP FORENSICS - Pocatello Instrument # 59740

Cannabinoids Analysis Report




Batch Data Path C:\MassHunter\Data\2018\THC Quant\071018 THC Quant Urine Validation SP worklist 2555\QuantResults\071
Analysis Time 7/11/2018 7:39 AM **Analyst Name** ISPUser
Report Time 7/11/2018 7:40 AM **Reporter Name** ISPUser
Last Calib Update 7/11/2018 7:39 AM **Batch State** Processed

Analysis Info

Acq Time 2018-07-10 11:41 **Data File** Cal 1.d
Sample Type Calibration **Sample Name** Cal 1
Dilution 1 **Acq Method** THC Quant 051517 workingmm.m
Position P1-H6 **Sample Info**
Inj Vol -1 **Comment**

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-D3	2.332	38047	1451456	0.0262	3.4724
THC-COOH	THC-COOH-D9	2.432	41312	398396	0.1037	4.4297
THC	THC-D3	6.226	12492	411452	0.0304	3.6125

ISP FORENSICS - Pocatello Instrument # 59740

Cannabinoids Analysis Report


P

Batch Data Path C:\MassHunter\Data\2018\THC Quant\071018 THC Quant Urine Validation SP worklist 2555\QuantResults\071018
Analysis Time 7/11/2018 7:39 AM **Analyst Name** ISPUser
Report Time 7/11/2018 7:40 AM **Reporter Name** ISPUser
Last Calib Update 7/11/2018 7:39 AM **Batch State** Processed

Analysis Info

Acq Time 2018-07-10 11:53 **Data File** Cal 2.d
Sample Type Calibration **Sample Name** Cal 2
Dilution 1 **Acq Method** THC Quant 051517 workingmm.m
Position P1-G6 **Sample Info**
Inj Vol -1 **Comment**

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-D3	2.319	71816	1639341	0.0438	5.0180
THC-COOH	THC-COOH-D9	2.419	57248	456278	0.1255	5.3756
THC	THC-D3	6.199	23278	464124	0.0502	5.2714

ISP FORENSICS - Pocatello Instrument # 59740

Cannabinoids Analysis Report


P

Batch Data Path C:\MassHunter\Data\2018\THC Quant\071018 THC Quant Urine Validation SP worklist 2555\QuantResults\071018
Analysis Time 7/11/2018 7:39 AM **Analyst Name** ISPUser
Report Time 7/11/2018 7:40 AM **Reporter Name** ISPUser
Last Calib Update 7/11/2018 7:39 AM **Batch State** Processed

Analysis Info

Acq Time 2018-07-10 12:05 **Data File** Cal 3.d
Sample Type Calibration **Sample Name** Cal 3
Dilution 1 **Acq Method** THC Quant 051517 workingmm.m
Position P1-F6 **Sample Info**
Inj Vol -1 **Comment**

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-D3	2.332	160903	1599698	0.1006	10.0055
THC-COOH	THC-COOH-D9	2.432	107707	417988	0.2577	11.1196
THC	THC-D3	6.212	43957	440528	0.0998	9.4303

D

ISP FORENSICS - Pocatello Instrument # 59740


Cannabinoids Analysis Report

Batch Data Path C:\MassHunter\Data\2018\THC Quant\071018 THC Quant Urine Validation SP worklist 2555\QuantResults\071018
Analysis Time 7/11/2018 7:39 AM **Analyst Name** ISPUser
Report Time 7/11/2018 7:40 AM **Reporter Name** ISPUser
Last Calib Update 7/11/2018 7:39 AM **Batch State** Processed

Analysis Info

Acq Time 2018-07-10 12:16 **Data File** Cal 4.d
Sample Type Calibration **Sample Name** Cal 4
Dilution 1 **Acq Method** THC Quant 051517 workingmm.m
Position P1-E6 **Sample Info**
Inj Vol -1 **Comment**

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-D3	2.332	351070	1557201	0.2254	20.9744
THC-COOH	THC-COOH-D9	2.432	204411	427492	0.4782	20.6985
THC	THC-D3	6.212	99419	425390	0.2337	20.6540

ISP FORENSICS - Pocatello Instrument # 59740

Cannabinoids Analysis Report


D

Batch Data Path C:\MassHunter\Data\2018\THC Quant\071018 THC Quant Urine Validation SP worklist 2555\QuantResults\071018
Analysis Time 7/11/2018 7:39 AM **Analyst Name** ISPUser
Report Time 7/11/2018 7:40 AM **Reporter Name** ISPUser
Last Calib Update 7/11/2018 7:39 AM **Batch State** Processed

Analysis Info

Acq Time 2018-07-10 12:28 **Data File** Cal 5.d
Sample Type Calibration **Sample Name** Cal 5
Dilution 1 **Acq Method** THC Quant 051517 workingmm.m
Position P1-D6 **Sample Info**
Inj Vol -1 **Comment**

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-D3	2.332	857242	1585268	0.5408	48.6725
THC-COOH	THC-COOH-D9	2.432	470175	422265	1.1135	48.2990
THC	THC-D3	6.212	239484	436082	0.5492	47.0906

ISP FORENSICS - Pocatello Instrument # 59740

Cannabinoids Analysis Report


P

Batch Data Path	C:\MassHunter\Data\2018\THC Quant\071018 THC Quant Urine Validation SP worklist 2555\QuantResults\071		
Analysis Time	7/11/2018 7:39 AM	Analyst Name	ISPUser
Report Time	7/11/2018 7:40 AM	Reporter Name	ISPUser
Last Calib Update	7/11/2018 7:39 AM	Batch State	Processed

Analysis Info

Acq Time	2018-07-10 12:40	Data File	Cal 6.d
Sample Type	Calibration	Sample Name	Cal 6
Dilution	1	Acq Method	THC Quant 051517 workingmm.m
Position	P1-C6	Sample Info	
Inj Vol	-1	Comment	

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-D3	2.319	1836042	1614264	1.1374	101.0838
THC-COOH	THC-COOH-D9	2.419	973235	421823	2.3072	100.1617
THC	THC-D3	6.159	544539	459111	1.1861	100.4646

ISP FORENSICS - Pocatello Instrument # 59740

Cannabinoids Analysis Report


D

Batch Data Path C:\MassHunter\Data\2018\THC Quant\071018 THC Quant Urine Validation SP worklist 2555\QuantResults\071
Analysis Time 7/11/2018 7:39 AM **Analyst Name** ISPUser
Report Time 7/11/2018 7:41 AM **Reporter Name** ISPUser
Last Calib Update 7/11/2018 7:39 AM **Batch State** Processed

Analysis Info

Acq Time 2018-07-10 12:52 **Data File** Cal 7.d
Sample Type Calibration **Sample Name** Cal 7
Dilution 1 **Acq Method** THC Quant 051517 workingmm.m
Position P1-B6 **Sample Info**
Inj Vol -1 **Comment**

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-D3	2.319	4596310	1598412	2.8755	253.7733
THC-COOH	THC-COOH-D9	2.419	2343351	400153	5.8561	254.3455
THC	THC-D3	6.186	1399097	459064	3.0477	256.4766