


**Worklist: 2749**

<u>LAB CASE</u>	<u>ITEM</u>	<u>TASK ID</u>	<u>DESCRIPTION</u>
M2018-4837	1	129665	AM 27 Blood THC Quant by LC-QQQ
M2018-4940	1	129666	AM 27 Blood THC Quant by LC-QQQ
P2018-2841	1	129667	AM 27 Blood THC Quant by LC-QQQ
P2018-2950	1	129668	AM 27 Blood THC Quant by LC-QQQ
P2018-2954	1	129669	AM 27 Blood THC Quant by LC-QQQ
P2018-2955	1	129670	AM 27 Blood THC Quant by LC-QQQ
P2018-3004	1	129671	AM 27 Blood THC Quant by LC-QQQ



# AM# 27: Quantitation of THC and Metabolites in Blood by LC-MS/MS

TS  
A

Extraction Date: 10/24/18  
Plate lot#: 0539904

Analyst: Sarah Pickle  
Plate Expiration: 09/10/19

**Mobile phase A:** 0.1% Formic Acid in LCMS Water  
MTBE  
**Blank Blood Lot:** ~~361331-1~~ <sup>18G207D7</sup> <sub>1</sub> <sup>5</sup>  
**LCMS-QQQ ID:** 59740

**Mobile phase B:** 0.1% Formic acid in Acetonitrile  
LCMS Methanol  
**Column:** UCT Selectra DA 100 x 2.1mm 3um  
Hexane

## Pre-Analytic:

- 1. Check levels of mobile phases and needle wash refill as needed. Ensure waste is not full.
- 2. Ensure correct column is installed and begin mobile phase flow allow to equilibrate ~ 30 minutes.
- 3. Create worklist: Data Path: \_\_\_\_\_

## Analytic:

- 1. Remove standards, plate, controls, and samples from cold storage. Allow to reach room temperature.
- 2. Pipette **1000µL blood (calibrated pipette) Pipette ID: 3382167** in wells of analytical (standards) plate.
- 3. Place on shaking incubator at ambient temp., 900rpm for 15 minutes. *Shaker ID: 067105*
- 4. Pipette **500µL 0.1% formic acid in water** in wells of analytical plate.
- 5. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 6. Transfer **800µL of blood+acid** mixture to corresponding wells of SLE+ plate.
- 7. Apply positive pressure for approx. 10-15 seconds (or until no liquid remains on top of sorbent).  
*(Load at 85-100 PSI- Selector to the right) Manifold ID: 067104*
- 8. Wait 5 minutes.
- 9. Add **2.25mL MTBE. (Add in 3 increments of 750uL)**
- 10. Wait 5 minutes.
- 11. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 12. Add **2.25mL Hexane. (Add in 3 increments of 750uL)**
- 13. Wait 5 minutes.
- 14. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 15. Remove plate containing eluate. Place on SPE Dry and evaporate to dryness at approx. 35°C.  
*SPE Dry ID: 067103*
- 16. Reconstitute in **100µL 100% MeOH** and heat seal plate with foil. Place in autosampler and run worklist.

## Post-Analytic

- 1. Create batch and process data.  
Worklist path: 102418 THCQ SP TS Batch Name: THCQ SP TS
- 2. Make any necessary integration changes, Curve weighting of Linear 1/x with r<sup>2</sup> values ≥0.98 for each analyte
- 3. RT +/- 3% or 0.100 min, whichever is greater, +/- 20% Accuracy for greater than (+/- 30% for 10ng/ml or less).  
Ion ratios must be within +/- 20% of the averaged calibrators
- 4. Case sample response for THC and OH-THC 3ng/mL (quantitative), Carboxy-THC: 10ng/mL (qualitative only) will be reported. Samples with a THC or OH-THC response over 50 ng/mL will be reported out as greater than 50 ng/mL.
- 5. Did all QCs pass for each analyte? Y / N
- 6. Enter QCs into control charting.
- 7. Central File Packet to include: LIMS Worklist, Method Checklist, Calibration and Control Reports

COMMENTS: Cases M2018-3955-1, M2018-4522-1, M2018-4572-1, M2018-4670-1, M2018-4887-4, P2018-2754-1, P2018-2800-1, P2018-2818-1 were also ran in this batch due to a possible interferant in carboxy-THC in the original run on 10/5/18. **Curve range limited Carboxy-THC: 5-100**

Sarah Pickle acted as the primary analyst on this run. I witnessed and approved of all methods used. - TS



# Idaho State Police Forensic Services

B

P

## AM #27 Quantitative Analysis of THC and Metabolites in Blood by LCMS-QQQ

**Analyst:** Sarah Pickle and Tamara Salazar

**Extraction Date:** 10/24/18

**Worklist Number:** 2749

<i>Reagent</i>	<i>Lot Number</i>	<i>Expiration Date</i>	<i>Date in Service</i>	<i>Date Out of Service</i>	<i>Initials</i>
ToxBox THC/THC Metabolite Plate	0539904	09/10/19			
Negative Blood	18G207D7		10/24/18		
Methanol External Control Solution	WS102418	02/07/19	10/24/18		
Blood External Control Solution	102418	02/07/19	10/24/18		
Methyl Tert-Butyl Ether (MTBE) 99.9%	A0375555		6/26/17		
Hexanes (ACS)	101642		10/26/17		
Methanol (LCMS Grade)	177145		4/11/18		
0.1% Formic Acid in Water (Mobile Phase A)	100518		10/05/18		
0.1% Formic Acid in Acetonitrile (Mobile Phase B)	176190		2/6/18		
Needle Rinse--75% LCMS MeOH in LCMS Water	100918		10/09/18		

**Methanol External Control Solution (Lot: WS102418)**

*10 ul of 1mg/mL THC, 100 ul of 100 ug/mL THC-OH, C-THC in 9790 ul MeOH*

<i>Component</i>	<i>Source</i>	<i>Source Lot Number</i>	<i>Expiration Date</i>
Methanol (LCMS)	Fisher	177145	
THC	Cerilliant	FE04231406	04/30/2019
C-THC	Cayman	0497429	02/08/2019
THC-OH	Cerilliant	FE01121503	01/31/2020
Prepared:	10/24/18		
Prepared By:	Tamara Salazar		
Expires:	02/07/19		

**Blood External Control Solution (Lot: 102418)**

*100 ul of methanol external control solution was added to 9900 ul of blood.*

<i>Component</i>	<i>Source</i>	<i>Source Lot Number</i>
Negative Blood	Lampire	18G207D7
Methanol External Control Solution		WS102418
Prepared:	10/24/18	
Prepared by:	Sarah Pickle	
Expires:	02/07/19	

TS  
R**0.1% Formic Acid in LCMS Water (Mobile Phase A) (Lot: 100518)**

<i>Component</i>	<i>Source</i>	<i>Source Lot Number</i>
Formic Acid (LCMS Grade)	Fisher	095180B
Water (LCMS Grade)	Fisher	182702
Prepared:	10/05/18	
Prepared By:	Sarah Pickle	

**Needle Rinse (75% LCMS MeOH in LCMS Water) (Lot: 100918)**

<i>Component</i>	<i>Source</i>	<i>Source Lot Number</i>
MeOH (LCMS Grade)	Fisher	177145
Water (LCMS Grade)	Fisher	182702
Prepared:	10/09/18	
Prepared By:	Sarah Pickle	

Handwritten initials: B and S

# ISP FORENSICS - Pocatello Instrument # 59740


## Cannabinoids Analysis Report

**Batch Data Path** C:\MassHunter\Data\2018\THC Quant\102418 THCQ SP TS\QuantResults\THCQ SP TS.batch.bin  
**Analysis Time** 10/26/2018 8:48 AM **Analyst Name** ISPUser  
**Report Time** 10/26/2018 8:50 AM **Reporter Name** ISPUser  
**Last Calib Update** 10/26/2018 8:48 AM **Batch State** Processed

**Analysis Info**

**Acq Time** 2018-10-24 13:18 **Data File** Negative.d  
**Sample Type** Sample **Sample Name** Negative  
**Dilution** 1 **Acq Method** THC Quant 051517 workingmm.m  
**Position** P1-H6 **Sample Info**  
**Inj Vol** -1 **Comment** Hemostat 361331-1

**Sample Chromatogram**



**Results**

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-D3	2.586	1315	467951	0.0028	0.8656
THC-COOH	THC-COOH-D9	2.299	3925	180318	0.0218	1.5203

TS JS

# ISP FORENSICS - Pocatello Instrument # 59740


## Cannabinoids Analysis Report

**Batch Data Path** C:\MassHunter\Data\2018\THC Quant\102418 THCQ SP TS\QuantResults\THCQ SP TS.batch.bin  
**Analysis Time** 10/26/2018 8:48 AM **Analyst Name** ISPUser  
**Report Time** 10/26/2018 8:50 AM **Reporter Name** ISPUser  
**Last Calib Update** 10/26/2018 8:48 AM **Batch State** Processed

### Analysis Info

**Acq Time** 2018-10-24 13:42 **Data File** External Control.d  
**Sample Type** Sample **Sample Name** External Control  
**Dilution** 1 **Acq Method** THC Quant 051517 workingmm.m  
**Position** P1-G6 **Sample Info**  
**Inj Vol** -1 **Comment** Hemostat 361331-1 + WS 020718

### Sample Chromatogram



### Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-D3	2.332	45393	465690	0.0975	8.6796
THC-COOH	THC-COOH-D9	2.432	31201	181566	0.1718	9.3820
THC	THC-D3	6.199	13404	165627	0.0809	9.8142

TS A

# ISP FORENSICS - Pocatello Instrument # 59740


## Cannabinoids Analysis Report

**Batch Data Path** C:\MassHunter\Data\2018\THC Quant\102418 THCQ SP TS\QuantResults\THCQ SP TS.batch.bin  
**Analysis Time** 10/26/2018 8:48 AM **Analyst Name** ISUser  
**Report Time** 10/26/2018 8:50 AM **Reporter Name** ISUser  
**Last Calib Update** 10/26/2018 8:48 AM **Batch State** Processed

**Analysis Info**

**Acq Time** 2018-10-24 12:55 **Data File** QC Control.d  
**Sample Type** Sample **Sample Name** QC Control  
**Dilution** 1 **Acq Method** THC Quant 051517 workingmm.m  
**Position** P1-A7 **Sample Info**  
**Inj Vol** -1 **Comment**

**Sample Chromatogram**



**Results**

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-D3	2.332	29892	565531	0.0529	4.9966
THC-COOH	THC-COOH-D9	2.432	38859	215452	0.1804	9.8280
THC	THC-D3	6.252	7837	203666	0.0385	4.7264

# ISP Forensics Calibration Curve Report


TS

D

**Batch Data Path** C:\MassHunter\Data\2018\THC Quant\102418 THCQ SP TS\QuantResults\THCQ SP  
TS.batch.bin

**Last Calib Update** 10/26/2018 8:48 AM      **Analyst Name** ISP TOX

**Target Compound** *THC-OH*  
**Internal Standard** *THC-OH-D3*



Sample	Level	Enabled	Exp Conc	Final Conc	Accuracy
Cal 1-1ng	1	<input type="checkbox"/>	1	2.2	216.9
Cal 2-3ng	2	<input checked="" type="checkbox"/>	3	3.4	112.9
Cal 3-5ng	3	<input checked="" type="checkbox"/>	5	4.9	98.5
Cal 4-10ng	4	<input checked="" type="checkbox"/>	10	9.3	93.5
Cal 5-25ng	5	<input checked="" type="checkbox"/>	25	23.4	93.5
Cal 6-50ng	6	<input checked="" type="checkbox"/>	50	49.6	99.2
Cal 7-100ng	7	<input checked="" type="checkbox"/>	100	102.4	102.4

# ISP Forensics Calibration Curve Report

TS  
D

**Batch Data Path** C:\MassHunter\Data\2018\THC Quant\102418 THCQ SP TS\QuantResults\THCQ SP


TS.batch.bin

**Last Calib Update** 10/26/2018 8:48 AM

**Analyst Name** ISP TOX

**Target Compound** *THC-COOH*

**Internal Standard** *THC-COOH-D9*



Sample	Level	Enabled	Exp Conc	Final Conc	Accuracy
Cal 1-1ng	1	<input type="checkbox"/>	1	2.1	214.7
Cal 2-3ng	2	<input type="checkbox"/>	3	4.0	133.6
Cal 3-5ng	3	<input checked="" type="checkbox"/>	5	5.3	106.1
Cal 4-10ng	4	<input checked="" type="checkbox"/>	10	10.0	100.5
Cal 5-25ng	5	<input checked="" type="checkbox"/>	25	23.4	93.6
Cal 6-50ng	6	<input checked="" type="checkbox"/>	50	48.6	97.2
Cal 7-100ng	7	<input checked="" type="checkbox"/>	100	102.6	102.6


# ISP Forensics Calibration Curve Report

TS  
D

**Batch Data Path** C:\MassHunter\Data\2018\THC Quant\102418 THCQ SP TS\QuantResults\THCQ SP  
TS.batch.bin

**Last Calib Update** 10/26/2018 8:48 AM      **Analyst Name** ISP TOX

**Target Compound** *THC*  
**Internal Standard** *THC-D3*



Sample	Level	Enabled	Exp Conc	Final Conc	Accuracy
Cal 1-1ng	1	<input type="checkbox"/>	1	1.2	122.0
Cal 2-3ng	2	<input checked="" type="checkbox"/>	3	3.2	106.7
Cal 3-5ng	3	<input checked="" type="checkbox"/>	5	4.9	97.8
Cal 4-10ng	4	<input checked="" type="checkbox"/>	10	9.9	99.5
Cal 5-25ng	5	<input checked="" type="checkbox"/>	25	23.8	95.1
Cal 6-50ng	6	<input checked="" type="checkbox"/>	50	49.8	99.6
Cal 7-100ng	7	<input checked="" type="checkbox"/>	100	101.4	101.4

TS  
D

# ISP FORENSICS - Pocatello Instrument # 59740


## Cannabinoids Analysis Report

**Batch Data Path** C:\MassHunter\Data\2018\THC Quant\102418 THCQ SP TS\QuantResults\THCQ SP TS.batch.bin  
**Analysis Time** 10/26/2018 8:48 AM **Analyst Name** ISPUser  
**Report Time** 10/26/2018 8:49 AM **Reporter Name** ISPUser  
**Last Calib Update** 10/26/2018 8:48 AM **Batch State** Processed

### Analysis Info

**Acq Time** 2018-10-24 11:20 **Data File** Cal 1-1ng.d  
**Sample Type** QC **Sample Name** Cal 1-1ng  
**Dilution** 1 **Acq Method** THC Quant 051517 workingmm.m  
**Position** P1-H7 **Sample Info**  
**Inj Vol** -1 **Comment**

### Sample Chromatogram



### Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-D3	2.345	9337	502107	0.0186	2.1686
THC-COOH	THC-COOH-D9	2.419	6496	192573	0.0337	2.1470
THC	THC-D3	6.239	1653	179163	0.0092	1.2201

15 18

# ISP FORENSICS - Pocatello Instrument # 59740


## Cannabinoids Analysis Report

**Batch Data Path** C:\MassHunter\Data\2018\THC Quant\102418 THCQ SP TS\QuantResults\THCQ SP TS.batch.bin  
**Analysis Time** 10/26/2018 8:48 AM **Analyst Name** ISUser  
**Report Time** 10/26/2018 8:49 AM **Reporter Name** ISUser  
**Last Calib Update** 10/26/2018 8:48 AM **Batch State** Processed

### Analysis Info

**Acq Time** 2018-10-24 11:32 **Data File** Cal 2-3ng.d  
**Sample Type** Calibration **Sample Name** Cal 2-3ng  
**Dilution** 1 **Acq Method** THC Quant 051517 workingmm.m  
**Position** P1-G7 **Sample Info**  
**Inj Vol** -1 **Comment**

### Sample Chromatogram



### Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-D3	2.345	18140	543733	0.0334	3.3874
THC-COOH	THC-COOH-D9	2.446	14281	206191	0.0693	4.0083
THC	THC-D3	6.252	5003	194247	0.0258	3.2016

# ISP FORENSICS - Pocatello Instrument # 59740

## Cannabinoids Analysis Report

TS


D

**Batch Data Path** C:\MassHunter\Data\2018\THC Quant\102418 THCQ SP TS\QuantResults\THCQ SP TS.batch.bin  
**Analysis Time** 10/26/2018 8:48 AM **Analyst Name** ISPUser  
**Report Time** 10/26/2018 8:49 AM **Reporter Name** ISPUser  
**Last Calib Update** 10/26/2018 8:48 AM **Batch State** Processed

**Analysis Info**

**Acq Time** 2018-10-24 11:43 **Data File** Cal 3-5ng.d  
**Sample Type** Calibration **Sample Name** Cal 3-5ng  
**Dilution** 1 **Acq Method** THC Quant 051517 workingmm.m  
**Position** P1-F7 **Sample Info**  
**Inj Vol** -1 **Comment**

**Sample Chromatogram**



**Results**

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-D3	2.345	28369	545466	0.0520	4.9266
THC-COOH	THC-COOH-D9	2.446	19902	211715	0.0940	5.3044
THC	THC-D3	6.239	7892	198171	0.0398	4.8879

# ISP FORENSICS - Pocatello Instrument # 59740

## Cannabinoids Analysis Report

TS


P

<b>Batch Data Path</b>	C:\MassHunter\Data\2018\THC Quant\102418 THCQ SP TS\QuantResults\THCQ SP TS.batch.bin		
<b>Analysis Time</b>	10/26/2018 8:48 AM	<b>Analyst Name</b>	ISPUser
<b>Report Time</b>	10/26/2018 8:49 AM	<b>Reporter Name</b>	ISPUser
<b>Last Calib Update</b>	10/26/2018 8:48 AM	<b>Batch State</b>	Processed

### Analysis Info

<b>Acq Time</b>	2018-10-24 11:55	<b>Data File</b>	Cal 4-10ng.d
<b>Sample Type</b>	Calibration	<b>Sample Name</b>	Cal 4-10ng
<b>Dilution</b>	1	<b>Acq Method</b>	THC Quant 051517 workingmm.m
<b>Position</b>	P1-E7	<b>Sample Info</b>	
<b>Inj Vol</b>	-1	<b>Comment</b>	

### Sample Chromatogram



### Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-D3	2.345	52147	494029	0.1056	9.3465
THC-COOH	THC-COOH-D9	2.446	35332	191493	0.1845	10.0452
THC	THC-D3	6.279	14086	171720	0.0820	9.9464

AS  
D

# ISP FORENSICS - Pocatello Instrument # 59740


## Cannabinoids Analysis Report

<b>Batch Data Path</b>	C:\MassHunter\Data\2018\THC Quant\102418 THCQ SP TS\QuantResults\THCQ SP TS.batch.bin		
<b>Analysis Time</b>	10/26/2018 8:48 AM	<b>Analyst Name</b>	ISPUser
<b>Report Time</b>	10/26/2018 8:50 AM	<b>Reporter Name</b>	ISPUser
<b>Last Calib Update</b>	10/26/2018 8:48 AM	<b>Batch State</b>	Processed

### Analysis Info

<b>Acq Time</b>	2018-10-24 12:07	<b>Data File</b>	Cal 5-25ng.d
<b>Sample Type</b>	Calibration	<b>Sample Name</b>	Cal 5-25ng
<b>Dilution</b>	1	<b>Acq Method</b>	THC Quant 051517 workingmm.m
<b>Position</b>	P1-D7	<b>Sample Info</b>	
<b>Inj Vol</b>	-1	<b>Comment</b>	

### Sample Chromatogram



### Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-D3	2.332	156465	567649	0.2756	23.3855
THC-COOH	THC-COOH-D9	2.446	93432	212649	0.4394	23.3962
THC	THC-D3	6.252	40685	206135	0.1974	23.7713

TS AS

# ISP FORENSICS - Pocatello Instrument # 59740


## Cannabinoids Analysis Report

**Batch Data Path** C:\MassHunter\Data\2018\THC Quant\102418 THCQ SP TS\QuantResults\THCQ SP TS.batch.bin  
**Analysis Time** 10/26/2018 8:48 AM **Analyst Name** ISPUser  
**Report Time** 10/26/2018 8:50 AM **Reporter Name** ISPUser  
**Last Calib Update** 10/26/2018 8:48 AM **Batch State** Processed

**Analysis Info**

**Acq Time** 2018-10-24 12:19 **Data File** Cal 6-50ng.d  
**Sample Type** Calibration **Sample Name** Cal 6-50ng  
**Dilution** 1 **Acq Method** THC Quant 051517 workingmm.m  
**Position** P1-C7 **Sample Info**  
**Inj Vol** -1 **Comment**

**Sample Chromatogram**



**Results**

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-D3	2.345	324380	546910	0.5931	49.5911
THC-COOH	THC-COOH-D9	2.446	188818	205038	0.9209	48.6200
THC	THC-D3	6.266	80519	194325	0.4144	49.7798

TS  
P

# ISP FORENSICS - Pocatello Instrument # 59740


## Cannabinoids Analysis Report

<b>Batch Data Path</b>	C:\MassHunter\Data\2018\THC Quant\102418 THCQ SP TS\QuantResults\THCQ SP TS.batch.bin		
<b>Analysis Time</b>	10/26/2018 8:48 AM	<b>Analyst Name</b>	ISPUser
<b>Report Time</b>	10/26/2018 8:50 AM	<b>Reporter Name</b>	ISPUser
<b>Last Calib Update</b>	10/26/2018 8:48 AM	<b>Batch State</b>	Processed

### Analysis Info

<b>Acq Time</b>	2018-10-24 12:31	<b>Data File</b>	Cal 7-100ng.d
<b>Sample Type</b>	Calibration	<b>Sample Name</b>	Cal 7-100ng
<b>Dilution</b>	1	<b>Acq Method</b>	THC Quant 051517 workingmm.m
<b>Position</b>	P1-B7	<b>Sample Info</b>	
<b>Inj Vol</b>	-1	<b>Comment</b>	

### Sample Chromatogram



### Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-D3	2.345	673072	546132	1.2324	102.3627
THC-COOH	THC-COOH-D9	2.446	388415	198982	1.9520	102.6341
THC	THC-D3	6.279	165126	195387	0.8451	101.4131