








REVIEWED
By Sarah Collins at 12:32 pm, Sep 18, 2023

Worklist: 6497

| <u>LAB_CASE</u> | <u>ITEM</u> | <u>ITEM_TYPE</u> | <u>DESCRIPTION</u> | |
|-----------------|-------------|------------------|---|---|
| C2023-1953 | 2 | BCK | AM 27 Blood THC Quant by LC-QQQ |  |
| C2023-1959 | 1 | BCK | AM 27 Blood THC Quant by LC-QQQ |  |
| C2023-1974 | 1 | UCK | AM 27 Urine Cannabinoids Confirmation by LC-QQQ |  |
| C2023-1979 | 1 | BCK | AM 27 Blood THC Quant by LC-QQQ |  |
| C2023-2003 | 1 | BCK | AM 27 Blood THC Quant by LC-QQQ |  |
| C2023-2077 | 1 | UCK | AM 27 Urine Cannabinoids Confirmation by LC-QQQ |  |
| C2023-2088 | 1 | UCK | AM 27 Urine Cannabinoids Confirmation by LC-QQQ |  |



AM# 27: Quantitation of THC and Metabolites in Blood and Urine by LC-MS/MS

Extraction Date 9/14/23
Plate lot#: 230627

Analyst: Anne Nord
Plate re-test: 12/27/2023

Mobile phase A: 0.1% Formic Acid in LCMS Water
MTBE

Mobile phase B: 0.1% Formic acid in Acetonitrile
Hexane

Blank Blood Lot: 23C57106 **Urine Blank:** 8423 **Column:** UCT Selectra DA 100 x 2.1mm 3um
LCMS-QQQ ID: 69679

Pre-Analytic:

- 1. Check levels of mobile phases and needle wash refill as needed. Ensure waste is not full.
- 2. Ensure correct column is installed and begin mobile phase flow allow to equilibrate ~ 30 minutes.

Analytic:

- 1. Remove standards, plate, controls, and samples from cold storage. Allow to reach room temperature.
- 2. Urine hydrolysis: add 1.5 ml urine to blank plate, add 250 ul 1N KOH mix and incubate at 40 degrees for 15 minutes.
Pipette 1000µL (calibrated pipette) blood or 1000µL hydrolyzed urine Pipette ID: K52558G in wells of analytical (standards) plate.
- 3. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 4. Pipette 500µL 0.1% formic acid in water blood sample, 500 ul saturated phosphate buffer in urine in wells of analytical plate.
- 5. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 6. Transfer 800µL of blood+acid or urine acid mixture to corresponding wells of SLE+ plate.
- 7. Apply positive pressure for approx. 10-15 seconds (or until no liquid remains on top of sorbent).
(Load at 85-100 PSI- Selector to the right) Manifold ID: 66792
- 8. Wait 5 minutes.
- 9. Add 2.25mL MTBE. (Add in 3 increments of 750uL)
- 10. Wait 5 minutes.
- 11. Apply positive pressure for approx. 15 seconds. (10-15 PSI- Selector to the left).
- 12. Add 2.25mL Hexane. (Add in 3 increments of 750uL)
- 13. Wait 5 minutes.
- 14. Apply positive pressure for approx. 15 seconds. (10-15 PSI- Selector to the left).
- 15. Remove plate containing eluate. Place on SPE Dry and evaporate to dryness at approx. 35°C.
SPE Dry ID: 66819
- 16. Reconstitute in 100µL 100% MeOH and heat seal plate with foil. Place in autosampler and run worklist.

Post-Analytic

- 1. Create batch and process data.
- 2. Make any necessary integration changes, Curve weighting of Linear 1/x with r^2 values ≥ 0.98 for each analyte
- 3. RT +/- 3% or 0.100 min, whichever is greater, +/- 20% Accuracy for greater than (+/- 30% for 10ng/ml or less). Ion ratios must be within +/- 20% of the averaged calibrators. SN > 10
- 4. Case sample response for THC 1ng/ml LOD 3ng/ml LOQ, OH-THC 3ng/mL LOD and LOQ, Carboxy-THC: 5 ng/mL (qualitative only). Samples with a THC or OH-THC response over 50 ng/mL will be reported out as greater than 50 ng/mL.
- 5. Did all QCs pass for each analyte? (if not is it describe in comments section)
- 6. Central File Packet to include: LIMS Worklist, Method Checklist, Calibration and Control Reports

COMMENTS:



| | 1 | 2 | 3 | 4 | 5 | 6 |
|---|--------------------------|------------------------|--------|---|---|---|
| a | cal 1 | Internal control urine | 2077-1 | | | |
| b | cal 2 | negative blood | 2088-1 | | | |
| c | cal 3 | 1953-2 | | | | |
| d | cal 4 | 1959-1 | | | | |
| e | cal 5 | 1979-1 | | | | |
| f | cal 6 | 2003-1 | | | | |
| g | cal 7 | negative urine | | | | |
| h | Internal control (blood) | 1974-1 | | | | |

Plate position 3


c2023-____-__

AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2023\am 27-28\091423\QuantResults\cann.batch.bin
Calibration Last Update 9/15/2023 8:26:00 AM

| | | | |
|-------------------------|----------------------|------------------|--|
| Instrument | 69679 | Data File | mj internal control blood.d |
| Type | QC | Sample | mj internal control blood |
| Acq. Method | AM 27 THC quant.m | Operator | Anne Nord |
| Sample Position | P3-H1 | Comment | Only drugs and concentrations listed on the laboratory report itself are appropriate to be used for interpretation purposes. Any drugs or values included in the notes but not included on the report are used by laboratory personnel to make determinations/reach conclusions within the confines of the methods |
| Injection Volume | 10 | | |
| Acq. Date-Time | 9/14/2023 5:05:48 PM | | |
| Sample Info. | | | |

Sample Chromatogram




| Name | RT | Resp. | S/N | Ratio | S/N | ISTD Resp. | Final Conc. |
|----------|-------|--------|--------|--------|-------|------------|--------------|
| THC-OH | 1.334 | 25835 | ∞ | 996.16 | ∞ | 2076723 | 4.197 ng/ml |
| THC-COOH | 1.358 | 41087 | 257.9 | 286.41 | ∞ | 555065 | 14.461 ng/ml |
| THC | 2.926 | 120835 | 1739.4 | 24.93 | 199.1 | 995094 | 4.147 ng/ml |

AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2023\am 27-28\091423\QuantResults\cann.batch.bin
Calibration Last Update 9/15/2023 8:26:00 AM

| | | | |
|-------------------------|----------------------|------------------|--|
| Instrument | 69679 | Data File | mj negative blood.d |
| Type | Sample | Sample | mj negative blood |
| Acq. Method | AM 27 THC quant.m | Operator | Anne Nord |
| Sample Position | P3-B2 | Comment | Only drugs and concentrations listed on the laboratory report itself are appropriate to be used for interpretation purposes. Any drugs or values included in the notes but not included on the report are used by laboratory personnel to make determinations/reach conclusions within the confines of the methods |
| Injection Volume | 10 | | |
| Acq. Date-Time | 9/14/2023 5:12:25 PM | | |
| Sample Info. | | | |

Sample Chromatogram




AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2023\am 27-28\091423\QuantResults\cann.batch.bin
Calibration Last Update 9/15/2023 8:26:00 AM

| | | | |
|-------------------------|----------------------|------------------|--|
| Instrument | 69679 | Data File | mj negative urine.d |
| Type | Sample | Sample | mj negative urine |
| Acq. Method | AM 27 THC quant.m | Operator | Anne Nord |
| Sample Position | P3-G2 | Comment | Only drugs and concentrations listed on the laboratory report itself are appropriate to be used for interpretation purposes. Any drugs or values included in the notes but not included on the report are used by laboratory personnel to make determinations/reach conclusions within the confines of the methods |
| Injection Volume | 10 | | |
| Acq. Date-Time | 9/14/2023 6:11:52 PM | | |
| Sample Info. | | | |

Sample Chromatogram




AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2023\am 27-28\091423\QuantResults\cann.batch.bin
Calibration Last Update 9/15/2023 8:26:00 AM

| | | | |
|-------------------------|----------------------|------------------|--|
| Instrument | 69679 | Data File | Internal control urine.d |
| Type | Sample | Sample | Internal control urine |
| Acq. Method | AM 27 THC quant.m | Operator | Anne Nord |
| Sample Position | P3-A2 | Comment | Only drugs and concentrations listed on the laboratory report itself are appropriate to be used for interpretation purposes. Any drugs or values included in the notes but not included on the report are used by laboratory personnel to make determinations/reach conclusions within the confines of the methods |
| Injection Volume | 10 | | |
| Acq. Date-Time | 9/14/2023 6:58:07 PM | | |
| Sample Info. | | | |

Sample Chromatogram




| Name | RT | Resp. | S/N | Ratio | S/N | ISTD Resp. | Final Conc. |
|----------|-------|-------|-------|--------|-------|------------|--------------|
| THC-OH | 1.334 | 26180 | ∞ | 819.64 | ∞ | 1983234 | 4.440 ng/ml |
| THC-COOH | 1.358 | 33996 | 966.7 | 273.14 | 957.7 | 465956 | 14.275 ng/ml |
| THC | 2.911 | 22803 | 485.0 | 23.39 | 47.4 | 194173 | 4.022 ng/ml |

AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2023\am 27-28\091423\QuantResults\cann.batch.bin
Calibration Last Update 9/15/2023 8:26:00 AM

| | | | |
|-------------------------|----------------------|------------------|--|
| Instrument | 69679 | Data File | end of run blood control blood h1.d |
| Type | Sample | Sample | end of run blood control blood h1 |
| Acq. Method | AM 27 THC quant.m | Operator | Anne Nord |
| Sample Position | P3-H1 | Comment | Only drugs and concentrations listed on the laboratory report itself are appropriate to be used for interpretation purposes. Any drugs or values included in the notes but not included on the report are used by laboratory personnel to make determinations/reach conclusions within the confines of the methods |
| Injection Volume | 10 | | |
| Acq. Date-Time | 9/14/2023 7:04:43 PM | | |
| Sample Info. | | | |

Sample Chromatogram




| Name | RT | Resp. | S/N | Ratio | S/N | ISTD Resp. | Final Conc. |
|----------|-------|--------|-------|--------|-------|------------|--------------|
| THC-OH | 1.334 | 30781 | ∞ | 910.44 | ∞ | 2247192 | 4.598 ng/ml |
| THC-COOH | 1.358 | 45170 | 488.2 | 283.54 | 598.1 | 603086 | 14.614 ng/ml |
| THC | 2.926 | 137190 | ∞ | 23.99 | 600.2 | 1153408 | 4.069 ng/ml |

Compound Calibration Report



Batch results D:\MassHunter\Data\2023\am 27-28\091423\QuantResults\cann.batch.bin
Last Cal. Update 9/15/2023 8:26 AM
Analyst Name ISP\datastor
Analyte THC **Internal Standard** THC-d3

THC - 7 Levels, 7 Levels Used, 7 Points, 7 Points Used, 1 QCs




| Sample | Level | Enabled | Expected Concentration | Final Concentration | Accuracy |
|----------|-------|---------|------------------------|---------------------|----------|
| mj cal 1 | 1 | ✓ | 1.0 | 1.2 | 118.3 |
| mj cal 2 | 2 | ✓ | 3.0 | 2.9 | 97.6 |
| mj cal 3 | 3 | ✓ | 5.0 | 4.7 | 94.2 |
| mj cal 4 | 4 | ✓ | 10.0 | 9.2 | 91.7 |
| mj cal 5 | 5 | ✓ | 25.0 | 24.3 | 97.1 |
| mj cal 6 | 6 | ✓ | 50.0 | 49.4 | 98.8 |
| mj cal 7 | 7 | ✓ | 100.0 | 102.3 | 102.3 |

Compound Calibration Report



Batch results D:\MassHunter\Data\2023\am 27-28\091423\QuantResults\cann.batch.bin
Last Cal. Update 9/15/2023 8:26 AM
Analyst Name ISP\datastor
Analyte THC-COOH **Internal Standard** THC-COOH-d9

THC-COOH - 7 Levels, 7 Levels Used, 7 Points, 7 Points Used, 1 QCs




| Sample | Level | Enabled | Expected Concentration | Final Concentration | Accuracy |
|----------|-------|---------|------------------------|---------------------|----------|
| mj cal 1 | 1 | ✓ | 5.0 | 5.6 | 112.2 |
| mj cal 2 | 2 | ✓ | 10.0 | 9.8 | 97.8 |
| mj cal 3 | 3 | ✓ | 20.0 | 18.9 | 94.5 |
| mj cal 4 | 4 | ✓ | 50.0 | 48.4 | 96.7 |
| mj cal 5 | 5 | ✓ | 75.0 | 73.0 | 97.3 |
| mj cal 6 | 6 | ✓ | 100.0 | 99.6 | 99.6 |
| mj cal 7 | 7 | ✓ | 250.0 | 254.8 | 101.9 |

Compound Calibration Report



Batch results D:\MassHunter\Data\2023\am 27-28\091423\QuantResults\cann.batch.bin
Last Cal. Update 9/15/2023 8:26 AM
Analyst Name ISP\datastor
Analyte THC-OH **Internal Standard** THC-OH-d3

THC-OH - 7 Levels, 7 Levels Used, 7 Points, 7 Points Used, 1 QCs



| Sample | Level | Enabled | Expected Concentration | Final Concentration | Accuracy |
|----------|-------|---------|------------------------|---------------------|----------|
| mj cal 1 | 1 | ✓ | 1.0 | 1.1 | 110.0 |
| mj cal 2 | 2 | ✓ | 3.0 | 3.0 | 99.5 |
| mj cal 3 | 3 | ✓ | 5.0 | 4.9 | 97.1 |
| mj cal 4 | 4 | ✓ | 10.0 | 9.6 | 96.1 |
| mj cal 5 | 5 | ✓ | 25.0 | 23.4 | 93.6 |
| mj cal 6 | 6 | ✓ | 50.0 | 51.6 | 103.3 |
| mj cal 7 | 7 | ✓ | 100.0 | 100.4 | 100.4 |


AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2023\am 27-28\091423\QuantResults\cann.batch.bin
Calibration Last Update 9/15/2023 8:26:00 AM

Instrument 69679
Type Cal
Acq. Method AM 27 THC quant.m
Sample Position P3-A1
Injection Volume 10
Acq. Date-Time 9/14/2023 4:19:27 PM
Sample Info.

Data File mj cal 1.d
Sample mj cal 1
Operator Anne Nord
Comment Only drugs and concentrations listed on the laboratory report itself are appropriate to be used for interpretation purposes. Any drugs or values included in the notes but not included on the report are used by laboratory personnel to make determinations/reach conclusions within the confines of the methods

Sample Chromatogram



| Name | RT | Resp. | S/N | Ratio | S/N | ISTD Resp. | Final Conc. | |
|----------|-------|-------|-------|-------------|-------|------------|-------------|-----|
| THC-OH | 1.334 | 5502 | ∞ | 1060.6 4 | ∞ | 2013799 | 1.100 ng/ml | Low |
| THC-COOH | 1.358 | 10995 | 162.8 | 307.13 | 469.2 | 471875 | 5.608 ng/ml | |
| THC | 2.926 | 20933 | 575.3 | 23.85 | 135.2 | 778372 | 1.183 ng/ml | |


AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2023\am 27-28\091423\QuantResults\cann.batch.bin
Calibration Last Update 9/15/2023 8:26:00 AM

Instrument 69679
Type Cal
Acq. Method AM 27 THC quant.m
Sample Position P3-B1
Injection Volume 10
Acq. Date-Time 9/14/2023 4:26:14 PM
Sample Info.

Data File mj cal 2.d
Sample mj cal 2
Operator Anne Nord
Comment Only drugs and concentrations listed on the laboratory report itself are appropriate to be used for interpretation purposes. Any drugs or values included in the notes but not included on the report are used by laboratory personnel to make determinations/reach conclusions within the confines of the methods

Sample Chromatogram



| Name | RT | Resp. | S/N | Ratio | S/N | ISTD Resp. | Final Conc. | |
|----------|-------|-------|-------|--------|-------|------------|-------------|-----|
| THC-OH | 1.334 | 17670 | ∞ | 866.69 | ∞ | 2044747 | 2.985 ng/ml | Low |
| THC-COOH | 1.358 | 25220 | 270.8 | 314.86 | ∞ | 534203 | 9.782 ng/ml | |
| THC | 2.911 | 70064 | 520.5 | 25.58 | 212.9 | 849201 | 2.927 ng/ml | |


AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2023\am 27-28\091423\QuantResults\cann.batch.bin
Calibration Last Update 9/15/2023 8:26:00 AM

Instrument 69679
Type Cal
Acq. Method AM 27 THC quant.m
Sample Position P3-C1
Injection Volume 10
Acq. Date-Time 9/14/2023 4:32:50 PM
Sample Info.

Data File mj cal 3.d
Sample mj cal 3
Operator Anne Nord
Comment Only drugs and concentrations listed on the laboratory report itself are appropriate to be used for interpretation purposes. Any drugs or values included in the notes but not included on the report are used by laboratory personnel to make determinations/reach conclusions within the confines of the methods

Sample Chromatogram



| Name | RT | Resp. | S/N | Ratio | S/N | ISTD Resp. | Final Conc. |
|----------|-------|--------|--------|--------|-----|------------|--------------|
| THC-OH | 1.334 | 29664 | 1042.3 | 831.17 | ∞ | 2044942 | 4.856 ng/ml |
| THC-COOH | 1.358 | 57397 | ∞ | 288.75 | ∞ | 577061 | 18.901 ng/ml |
| THC | 2.911 | 116953 | 1854.9 | 23.90 | ∞ | 838991 | 4.710 ng/ml |


AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2023\am 27-28\091423\QuantResults\cann.batch.bin
Calibration Last Update 9/15/2023 8:26:00 AM

Instrument 69679
Type Cal
Acq. Method AM 27 THC quant.m
Sample Position P3-D1
Injection Volume 10
Acq. Date-Time 9/14/2023 4:39:26 PM
Sample Info.

Data File mj cal 4.d
Sample mj cal 4
Operator Anne Nord
Comment Only drugs and concentrations listed on the laboratory report itself are appropriate to be used for interpretation purposes. Any drugs or values included in the notes but not included on the report are used by laboratory personnel to make determinations/reach conclusions within the confines of the methods

Sample Chromatogram



| Name | RT | Resp. | S/N | Ratio | S/N | ISTD Resp. | Final Conc. |
|----------|-------|--------|----------|--------|----------|------------|--------------|
| THC-OH | 1.334 | 61098 | ∞ | 878.42 | ∞ | 2078066 | 9.608 ng/ml |
| THC-COOH | 1.358 | 146093 | 2984.7 | 275.38 | ∞ | 544612 | 48.360 ng/ml |
| THC | 2.911 | 266146 | 3125.3 | 25.89 | 544.8 | 945319 | 9.166 ng/ml |


AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2023\am 27-28\091423\QuantResults\cann.batch.bin
Calibration Last Update 9/15/2023 8:26:00 AM

Instrument 69679
Type Cal
Acq. Method AM 27 THC quant.m
Sample Position P3-E1
Injection Volume 10
Acq. Date-Time 9/14/2023 4:46:02 PM
Sample Info.

Data File mj cal 5.d
Sample mj cal 5
Operator Anne Nord
Comment Only drugs and concentrations listed on the laboratory report itself are appropriate to be used for interpretation purposes. Any drugs or values included in the notes but not included on the report are used by laboratory personnel to make determinations/reach conclusions within the confines of the methods

Sample Chromatogram



| Name | RT | Resp. | S/N | Ratio | S/N | ISTD Resp. | Final Conc. |
|----------|-------|--------|---------|--------|-------|------------|--------------|
| THC-OH | 1.334 | 144616 | 4579.0 | 923.56 | ∞ | 1990713 | 23.404 ng/ml |
| THC-COOH | 1.358 | 226130 | ∞ | 289.89 | ∞ | 552597 | 72.962 ng/ml |
| THC | 2.911 | 692261 | 12012.2 | 25.32 | 638.3 | 906442 | 24.282 ng/ml |


AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2023\am 27-28\091423\QuantResults\cann.batch.bin
Calibration Last Update 9/15/2023 8:26:00 AM

Instrument 69679
Type Cal
Acq. Method AM 27 THC quant.m
Sample Position P3-F1
Injection Volume 10
Acq. Date-Time 9/14/2023 4:52:38 PM
Sample Info.

Data File mj cal 6.d
Sample mj cal 6
Operator Anne Nord
Comment Only drugs and concentrations listed on the laboratory report itself are appropriate to be used for interpretation purposes. Any drugs or values included in the notes but not included on the report are used by laboratory personnel to make determinations/reach conclusions within the confines of the methods

Sample Chromatogram



| Name | RT | Resp. | S/N | Ratio | S/N | ISTD Resp. | Final Conc. |
|----------|-------|---------|--------|--------|--------|------------|--------------|
| THC-OH | 1.334 | 339933 | ∞ | 856.11 | ∞ | 2109338 | 51.643 ng/ml |
| THC-COOH | 1.358 | 297636 | 6061.9 | 275.98 | 2968.9 | 529817 | 99.588 ng/ml |
| THC | 2.926 | 1494181 | ∞ | 25.79 | 2401.7 | 954520 | 49.414 ng/ml |


AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2023\am 27-28\091423\QuantResults\cann.batch.bin
Calibration Last Update 9/15/2023 8:26:00 AM

Instrument 69679
Type Cal
Acq. Method AM 27 THC quant.m
Sample Position P3-G1
Injection Volume 10
Acq. Date-Time 9/14/2023 4:59:14 PM
Sample Info.

Data File mj cal 7.d
Sample mj cal 7
Operator Anne Nord
Comment Only drugs and concentrations listed on the laboratory report itself are appropriate to be used for interpretation purposes. Any drugs or values included in the notes but not included on the report are used by laboratory personnel to make determinations/reach conclusions within the confines of the methods

Sample Chromatogram



| Name | RT | Resp. | S/N | Ratio | S/N | ISTD Resp. | Final Conc. |
|----------|-------|---------|---------|--------|-----|------------|---------------|
| THC-OH | 1.334 | 589481 | ∞ | 889.60 | ∞ | 1877354 | 100.404 ng/ml |
| THC-COOH | 1.358 | 689015 | ∞ | 271.34 | ∞ | 474829 | 254.799 ng/ml |
| THC | 2.911 | 2980439 | 39275.6 | 25.88 | ∞ | 916242 | 102.318 ng/ml |