

REVIEWED


By Anne Nord at 10:43 am, Apr 09, 2024

4/8/2024

Bylee

Worklist: 6757

<u>LAB_CASE</u>	<u>ITEM</u>	<u>ITEM_TYPE</u>	<u>DESCRIPTION</u>
C2024-0564	1	BCK	AM 27 Blood THC Quant by LC-QQQ
C2024-0595	1	BCK	AM 27 Blood THC Quant by LC-QQQ
C2024-0598	1	BCK	AM 27 Blood THC Quant by LC-QQQ
C2024-0603	1	BCK	AM 27 Blood THC Quant by LC-QQQ
C2024-0620	1	BCK	AM 27 Blood THC Quant by LC-QQQ
C2024-0633	1	BCK	AM 27 Blood THC Quant by LC-QQQ
C2024-0647	1	BCK	AM 27 Blood THC Quant by LC-QQQ



Bylye

AM# 27: Quantitation of THC and Metabolites in Blood and Urine by LC-MS/MS

Extraction Date: 4/8/24

Plate lot#: 231212

Mobile phase A: 0.1% Formic Acid in LCMS Water

Blank Blood Lot: 24C52042 (Btl 3)

Column: UCT Selectra DA 100 x 2.1mm 3um

Analyst: Britany Wylie

Plate Retest Date: 6/12/24

Mobile phase B: 0.1% Formic acid in Acetonitrile

Blank Urine Lot: 1324

LCMS-QQQ ID: 69679

Pre-Analytic:

- 1. Check levels of mobile phases and needle wash refill as needed. Ensure waste is not full.
- 2. Ensure correct column is installed and begin mobile phase flow allow to equilibrate ~ 30 minutes.

Analytic:

- 1. Remove standards, plate, controls, and samples from cold storage. Allow to reach room temperature.
- 2. Urine hydrolysis: add 1.5mL urine to blank plate, add 250µl 1N KOH. Shake and incubate at 40 degrees for 15 minutes. Using a calibrated pipette, add 1000µl blood and urine (if applicable) (calibrated pipette) into the appropriate wells of analytical (standards) plate. Pipette ID: I41142J
- 3. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 4. Add 500µL of 0.1% formic acid in water to blood samples, and 500µL of saturated phosphate buffer to urine samples in the wells of the analytical plate.
- 5. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 6. Transfer 700-800µL of blood+acid or urine+acid mixture to corresponding wells of SLE+ plate. Amount transferred: 800 µL
- 7. Apply positive pressure for approx. 10-15 seconds (or until no liquid remains on top of sorbent). *(Load at 85-100 PSI- Selector to the right)*
- 8. Wait 5 minutes.
- 9. Add 2.25mL MTBE. *(Add in 3 increments of 750uL)*
- 10. Wait 5 minutes.
- 11. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 12. Add 2.25mL Hexane. *(Add in 3 increments of 750uL)*
- 13. Wait 5 minutes.
- 14. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 15. Remove plate containing eluate. Place on SPE Dry and evaporate to dryness at approx. 35°C.
- 16. Reconstitute in 100µL 100% MeOH and heat seal plate with foil. Place in autosampler and run worklist.

Post-Analytic

- 1. Create batch and process data.
- 2. Make any necessary integration changes, Curve weighting of Linear 1/x with r² values ≥0.98 for each analyte
- 3. RT +/- 3% or 0.100 min, whichever is greater, +/- 20% Accuracy for greater than (+/- 30% for 10ng/ml or less). Ion ratios must be within +/- 20% of the averaged calibrators
- 4. Did all QCs pass for each analyte? (if not, describe in comments section)
- 5. Enter QCs into control charting.
- 6. Central File Packet to include: LIMS Worklist, Method Checklist, Calibration and Control Reports

COMMENTS: 4/5/24: Extraction failed for all samples, samples were re-extracted and run on 4/8/24

BWylee

	1	2	3	4	5	6
a	cal 1	Internal control urine	0598-1			
b	cal 2	negative blood				
c	cal 3	0595-1				
d	cal 4	0603-1				
e	cal 5	0620-1				
f	cal 6	0633-1				
g	cal 7	0647-1				
h	Internal control (blood)	0564-1				

Plate position 3

c2024-____-__

additional samples added to extraction plate upon re-extraction


AM #27 Cannabinoids

Wylie

Batch results D:\MassHunter\Data\2024\am 27-28\040824\QuantResults\am 27 4-8-24.batch.bin
Calibration Last Update 4/8/2024 3:12:16 PM

Instrument	69679	Data File	mj negative blood b2.d
Type	Sample	Sample	mj negative blood b2
Acq. Method	thc quant 50 50.m	Operator	Britany Wylie
Sample Position	P3-B2	Comment	
Injection Volume	10		
Acq. Date-Time	4/8/2024 12:55:04 PM		
Sample Info.			

Sample Chromatogram



AM #27 Cannabinoids


BWylie

Batch results D:\MassHunter\Data\2024\am 27-28\040824\QuantResults\am 27 4-8-24.batch.bin
Calibration Last Update 4/8/2024 3:12:16 PM

Instrument	69679	Data File	mj internal control blood.d
Type	QC	Sample	mj internal control blood
Acq. Method	thc quant 50 50.m	Operator	Britany Wylie
Sample Position	P3-H1	Comment	
Injection Volume	10		
Acq. Date-Time	4/8/2024 12:48:30 PM		

Sample Info.

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	2.051	33527	604.4	888.21	∞	2333362	4.597 ng/ml
THC-COOH	2.122	60105	38206.7	273.43	372.2	859693	14.075 ng/ml
THC	3.573	66858	∞	23.85	∞	439391	5.108 ng/ml

AM #27 Cannabinoids


Wylie

Batch results D:\MassHunter\Data\2024\am 27-28\040824\QuantResults\am 27 4-8-24.batch.bin
Calibration Last Update 4/8/2024 3:12:16 PM

Instrument	69679	Data File	mj internal blood control end of run.r.d
Type	Sample	Sample	mj internal blood control end of run.r
Acq. Method	thc quant 50 50.m	Operator	Britany Wylie
Sample Position	P3-H1	Comment	
Injection Volume	10		
Acq. Date-Time	4/8/2024 2:34:01 PM		

Sample Info.

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	2.051	36550	217.0	885.00	∞	2525369	4.629 ng/ml
THC-COOH	2.122	63934	1129.3	272.05	280.5	904891	14.211 ng/ml
THC	3.573	63283	18303.9	25.24	85.3	455255	4.693 ng/ml


Compound Calibration Report

Batch results D:\MassHunter\Data\2024\am 27-28\040824\QuantResults\am 27 4-8-24.batch.bin
Last Cal. Update 4/8/2024 3:12 PM
Analyst Name ISP\datastor
Analyte THC-COOH

Bylee

Internal Standard THC-COOH-d9

THC-COOH - 7 Levels, 7 Levels Used, 7 Points, 7 Points Used, 1 QCs




Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
mj cal 1	1	✓	5.0	5.3	106.4
mj cal 2	2	✓	10.0	10.0	99.8
mj cal 3	3	✓	20.0	19.4	96.8
mj cal 4	4	✓	50.0	48.7	97.4
mj cal 5	5	✓	75.0	73.8	98.4
mj cal 6	6	✓	100.0	100.1	100.1
mj cal 7	7	✓	250.0	252.7	101.1

Compound Calibration Report

Batch results D:\MassHunter\Data\2024\am 27-28\040824\QuantResults\am 27 4-8-24.batch.bin
Last Cal. Update 4/8/2024 3:12 PM
Analyst Name ISP\datastor
Analyte THC-OH **Internal Standard** THC-OH-d3

Byylee

THC-OH - 7 Levels, 7 Levels Used, 7 Points, 7 Points Used, 1 QCs



Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
mj cal 1	1	✓	1.0	1.2	119.5
mj cal 2	2	✓	3.0	2.9	97.7
mj cal 3	3	✓	5.0	4.6	92.7
mj cal 4	4	✓	10.0	9.2	91.9
mj cal 5	5	✓	25.0	24.0	96.0
mj cal 6	6	✓	50.0	50.2	100.4
mj cal 7	7	✓	100.0	101.9	101.9

AM #27 Cannabinoids


BWylie

Batch results D:\MassHunter\Data\2024\am 27-28\040824\QuantResults\am 27 4-8-24.batch.bin
Calibration Last Update 4/8/2024 3:12:16 PM

Instrument	69679	Data File	mj cal 1.d
Type	Cal	Sample	mj cal 1
Acq. Method	thc quant 50 50.m	Operator	Britany Wylie
Sample Position	P3-A1	Comment	
Injection Volume	10		
Acq. Date-Time	4/8/2024 12:02:26 PM		

Sample Info.

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	2.066	6785	146.6	716.29	150.2	2243332	1.195 ng/ml Low
THC-COOH	2.137	19046	127.8	279.56	44073.8	863529	5.318 ng/ml
THC	3.573	15977	492.2	24.68	∞	590881	1.158 ng/ml

AM #27 Cannabinoids


BWylie

Batch results D:\MassHunter\Data\2024\am 27-28\040824\QuantResults\am 27 4-8-24.batch.bin
Calibration Last Update 4/8/2024 3:12:16 PM

Instrument	69679	Data File	mj cal 2.d
Type	Cal	Sample	mj cal 2
Acq. Method	thc quant 50 50.m	Operator	Britany Wylie
Sample Position	P3-B1	Comment	
Injection Volume	10		
Acq. Date-Time	4/8/2024 12:09:03 PM		

Sample Info.

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	2.051	24306	∞	816.34	∞	2759058	2.930 ng/ml Low
THC-COOH	2.122	45322	70655.3	268.80	3227.4	953506	9.980 ng/ml
THC	3.573	70385	∞	28.23	∞	888325	2.806 ng/ml


AM #27 Cannabinoids

BWylie

Batch results D:\MassHunter\Data\2024\am 27-28\040824\QuantResults\am 27 4-8-24.batch.bin
Calibration Last Update 4/8/2024 3:12:16 PM

Instrument	69679	Data File	mj cal 3.d
Type	Cal	Sample	mj cal 3
Acq. Method	thc quant 50 50.m	Operator	Britany Wylie
Sample Position	P3-C1	Comment	
Injection Volume	10		
Acq. Date-Time	4/8/2024 12:15:37 PM		

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	2.051	34768	∞	873.76	∞	2399562	4.634 ng/ml
THC-COOH	2.122	87616	474.3	268.17	786.4	886563	19.365 ng/ml
THC	3.573	95489	42831.1	28.24	∞	698111	4.623 ng/ml


AM #27 Cannabinoids

BWylie

Batch results D:\MassHunter\Data\2024\am 27-28\040824\QuantResults\am 27 4-8-24.batch.bin
Calibration Last Update 4/8/2024 3:12:16 PM

Instrument	69679	Data File	mj cal 4.d
Type	Cal	Sample	mj cal 4
Acq. Method	thc quant 50 50.m	Operator	Britany Wylie
Sample Position	P3-D1	Comment	
Injection Volume	10		
Acq. Date-Time	4/8/2024 12:22:11 PM		

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	2.051	76728	∞	886.38	1084.9	2584149	9.193 ng/ml
THC-COOH	2.122	238309	1395.5	269.79	560.7	919193	48.721 ng/ml
THC	3.558	252455	∞	24.29	8290.7	839972	9.793 ng/ml


AM #27 Cannabinoids

BWylie

Batch results D:\MassHunter\Data\2024\am 27-28\040824\QuantResults\am 27 4-8-24.batch.bin
Calibration Last Update 4/8/2024 3:12:16 PM

Instrument	69679	Data File	mj cal 5.d
Type	Cal	Sample	mj cal 5
Acq. Method	thc quant 50 50.m	Operator	Britany Wylie
Sample Position	P3-E1	Comment	
Injection Volume	10		
Acq. Date-Time	4/8/2024 12:28:48 PM		

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	2.051	190947	3355.4	868.13	∞	2416197	23.989 ng/ml
THC-COOH	2.122	338650	28469.8	267.20	1735.7	854682	73.783 ng/ml
THC	3.573	521621	25634.2	25.54	2651.0	673587	24.751 ng/ml


AM #27 Cannabinoids

Wylie

Batch results D:\MassHunter\Data\2024\am 27-28\040824\QuantResults\am 27 4-8-24.batch.bin
Calibration Last Update 4/8/2024 3:12:16 PM

Instrument	69679	Data File	mj cal 6.d
Type	Cal	Sample	mj cal 6
Acq. Method	thc quant 50 50.m	Operator	Britany Wylie
Sample Position	P3-F1	Comment	
Injection Volume	10		
Acq. Date-Time	4/8/2024 12:35:22 PM		

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	2.051	446307	∞	834.83	∞	2682449	50.186 ng/ml
THC-COOH	2.122	493108	572702.5	271.37	483205.8	913118	100.095 ng/ml
THC	3.573	1126680	764372.2	26.36	∞	710200	50.385 ng/ml

AM #27 Cannabinoids


Wylie

Batch results D:\MassHunter\Data\2024\am 27-28\040824\QuantResults\am 27 4-8-24.batch.bin
Calibration Last Update 4/8/2024 3:12:16 PM

Instrument	69679	Data File	mj cal 7.d
Type	Cal	Sample	mj cal 7
Acq. Method	thc quant 50 50.m	Operator	Britany Wylie
Sample Position	P3-G1	Comment	
Injection Volume	10		
Acq. Date-Time	4/8/2024 12:41:56 PM		

Sample Info.

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	2.051	864904	36265.7	834.40	∞	2553369	101.874 ng/ml
THC-COOH	2.122	1140854	1287.8	268.09	∞	830168	252.737 ng/ml
THC	3.573	1431052	79687.3	25.61	18591.3	450943	100.485 ng/ml