

**Quantitative Analysis for Ethanol & Qualitative Analysis for Other Volatiles**

*Analytical Method(s): 1.0*

*Device: Hamilton MICROLAB 600 Liquid Processor/Dilutor Serial Number: ML600HC11378*

**Volatiles Quality Assurance Controls**

**Run Date(s): 5/3/2018**

Calibration Date: 5/2/18

Control Level	Expiration	Lot #	Target Value	Acceptable Range	Overall Results
Level 1	Jul-18	1407031	0.0780	0.0702-0.0858	0.0771 g/100cc
					0.0795 g/100cc
					0.2035 g/100cc
Level 2	Jul-18	1407032	0.2020	0.1818-0.2222	g/100cc
					g/100cc
<b>Multi-Component mixture:</b>		<b>Exp date: Sept 2020</b>	<b>Lot #</b>	<b>FN06041503</b>	<b>OK</b>
<b>Curve Fit:</b>		<b>Column 1</b>	<b>0.99995</b>	<b>Column2</b>	<b>0.99991</b>

<b>Ethanol Calibration Reference Material</b>								
Calibrator level	Expiration	Cerilliant Lot #	Target Value	Acceptable Range	Column 1	Column 2	Precision	Mean
0.050	Jul-19	FN06231406	0.050	0.045 - 0.055	0.0510	0.0523	0.0013	0.0516
0.080			0.080	0.072 - 0.088			0	#DIV/0!
0.100	Jun-20	FN06181501	0.100	0.090 - 0.110	0.0989	0.0991	0.0002	0.099
0.200	Dec-19	FN12011401	0.200	0.180 - 0.220	0.1982	0.1963	0.0019	0.1972
0.300	Feb-21	FN02121601	0.300	0.270 - 0.330	0.3026	0.3021	0.0005	0.3023
0.400			0.400	0.360 - 0.440			0	#DIV/0!
0.500	Aug-19	FN07031402	0.500	0.450 - 0.550	0.4992	0.5001	0.0009	0.4996

<b>Aqueous Controls</b>					
Control level	Expiration	Cerilliant Lot #	Target Value	Acceptable Range	Overall Results
0.080	Nov-20	FN10281510	0.08000	0.076 - 0.084	0.080 g/100cc

Issued: 4/22/2015


~Any information on this document can be changed for laboratory use, except for the precision and mean determination formulas.

Volatiles QA/QC data spreadsheet Rev 5

Issuing Authority: Quality Manager


**Worklist: 2376**

<u>LAB CASE</u>	<u>ITEM</u>	<u>TASK ID</u>	<u>DESCRIPTION</u>
M2018-2105	1	113883	Alcohol Analysis
M2018-2106	1	113889	Alcohol Analysis
M2018-2107	1	113893	Alcohol Analysis
M2018-2108	1	113894	Alcohol Analysis
M2018-2134	1	114052	Alcohol Analysis
M2018-2135	1	114054	Alcohol Analysis
M2018-2145	1	114095	Alcohol Analysis
M2018-2145	2	114096	Alcohol Analysis
M2018-2156	1	114146	Alcohol Analysis
M2018-2157	1	114147	Alcohol Analysis
M2018-2170	1	114165	Alcohol Analysis
M2018-2171	1	114166	Alcohol Analysis
M2018-2176	1	114217	Alcohol Analysis
M2018-2177	1	114218	Alcohol Analysis
M2018-2206	1	114353	Alcohol Analysis
M2018-2208	1	114356	Alcohol Analysis
M2018-2213	1	114374	Alcohol Analysis
M2018-2214	1	114375	Alcohol Analysis
M2018-2217	1	114379	Alcohol Analysis



ISP Forensic Services Blood Alcohol Report


Sample Name : 0.050 FN06231406  
 Laboratory : Meridian  
 Injection Date : May 2, 2018  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	4.51944	0.0510	g/100cc
2.	Ethanol	Column 2:	4.65226	0.0523	g/100cc
3.	n-Propanol	Column 1:	46.43829	1.0000	g/100cc
4.	n-Propanol	Column 2:	48.38544	1.0000	g/100cc

ISP Forensic Services Blood Alcohol Report


Sample Name : 0.100 FN06181501  
 Laboratory : Meridian  
 Injection Date : May 2, 2018  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	9.08848	0.0989	g/100cc
2.	Ethanol	Column 2:	9.36702	0.0991	g/100cc
3.	n-Propanol	Column 1:	47.58935	1.0000	g/100cc
4.	n-Propanol	Column 2:	49.03825	1.0000	g/100cc

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.200 FN12011401  
 Laboratory : Meridian  
 Injection Date : May 2, 2018  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	18.23830	0.1982	g/100cc
2.	Ethanol	Column 2:	18.96443	0.1963	g/100cc
3.	n-Propanol	Column 1:	47.32145	1.0000	g/100cc
4.	n-Propanol	Column 2:	48.90527	1.0000	g/100cc

JK

ISP Forensic Services Blood Alcohol Report


Sample Name : 0.300 FN02121601  
 Laboratory : Meridian  
 Injection Date : May 2, 2018  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	27.17083	0.3026	g/100cc
2.	Ethanol	Column 2:	28.50989	0.3021	g/100cc
3.	n-Propanol	Column 1:	46.07795	1.0000	g/100cc
4.	n-Propanol	Column 2:	47.36019	1.0000	g/100cc

ISP Forensic Services Blood Alcohol Report


Sample Name : 0.500 FN07031402  
 Laboratory : Meridian  
 Injection Date : May 2, 2018  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014 -CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	45.59841	0.4992	g/100cc
2.	Ethanol	Column 2:	48.28486	0.5001	g/100cc
3.	n-Propanol	Column 1:	46.80014	1.0000	g/100cc
4.	n-Propanol	Column 2:	48.14449	1.0000	g/100cc

ISP Forensic Services Blood Alcohol Report

Sample Name : INTERNAL STANDARD BLANK  
 Laboratory : Meridian  
 Injection Date : May 2, 2018  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	0.00000	0.0000	g/100cc
2.	Ethanol	Column 2:	0.00000	0.0000	g/100cc
3.	n-Propanol	Column 1:	44.40972	1.0000	g/100cc
4.	n-Propanol	Column 2:	45.62799	1.0000	g/100cc

JG



ISP Forensic Services Blood Alcohol Report

Sample Name : BLANK  
 Laboratory : Meridian  
 Injection Date : May 2, 2018  
 Method : SHUTDOWN.M  
 Acq. Instrument: CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	0.00000	0.0000	g/100cc
2.	Ethanol	Column 2:	0.00000	0.0000	g/100cc
3.	n-Propanol	Column 1:	0.00000	0.0000	g/100cc
4.	n-Propanol	Column 2:	0.00000	0.0000	g/100cc

JK

S a m p l e S u m m a r y

Sequence table: C:\Chem32\1\Data\05-02-18\_CAL\05-02-18\_CAL 2018-05-02 14-58-55\05-02-18\_CAL.S  
 Data directory path: C:\Chem32\1\Data\05-02-18\_CAL\05-02-18\_CAL 2018-05-02 14-58-55\  
 Logbook: C:\Chem32\1\Data\05-02-18\_CAL\05-02-18\_CAL 2018-05-02 14-58-55\05-02-18\_CAL.LOG  
 Sequence start: 5/2/2018 3:13:33 PM  
 Sequence Operator: SYSTEM  
 Operator: SYSTEM

Method file name: C:\Chem32\1\Data\05-02-18\_CAL\05-02-18\_CAL 2018-05-02 14-58-55\ALCOHOL.M

Run #	Location #	Inj #	Sample Name	Sample Amt [g/100cc]	Multip.* Dilution	File name	Cal #	# Cmp
1	1	1	0.050 FN06231406	-	1.0000	001F0101.D	*	4
2	2	1	0.100 FN06181501	-	1.0000	002F0201.D	*	4
3	3	1	0.200 FN12011401	-	1.0000	003F0301.D	*	4
4	4	1	0.300 FN02121601	-	1.0000	004F0401.D	*	4
5	5	1	0.500 FN07031402	-	1.0000	005F0501.D	*	4
6	6	1	INTERNAL STANDAR	-	1.0000	006F0601.D		2

Method file name: C:\Chem32\1\Data\05-02-18\_CAL\05-02-18\_CAL 2018-05-02 14-58-55\SHUTDOWN.M

Run #	Location #	Inj #	Sample Name	Sample Amt [g/100cc]	Multip.* Dilution	File name	Cal #	# Cmp
7	6	1	BLANK	-	1.0000	006F0701.D		0

JG

=====  
Calibration Table  
=====

-----  
General Calibration Setting  
-----

Calib. Data Modified : Wednesday, May 02, 2018 4:04:05 PM  
Signals calculated separately : No

Rel. Reference Window : 0.000 %  
Abs. Reference Window : 0.100 min  
Rel. Non-ref. Window : 0.000 %  
Abs. Non-ref. Window : 0.100 min  
Uncalibrated Peaks : not reported  
Partial Calibration : Yes, identified peaks are recalibrated  
Correct All Ret. Times: No, only for identified peaks

Curve Type : Linear  
Origin : Ignored  
Weight : Equal

Recalibration Settings:  
Average Response : Average all calibrations  
Average Retention Time: Floating Average New 75%

Calibration Report Options :  
Printout of recalibrations within a sequence:  
    Calibration Table after Recalibration  
    Normal Report after Recalibration  
If the sequence is done with bracketing:  
    Results of first cycle (ending previous bracket)

Default Sample ISTD Information (if not set in sample table):

ISTD #	ISTD Amount [g/100cc]	Name
1	1.00000	n-propanol
2	1.00000	n-propanol

-----  
Signal Details  
-----

Signal 1: FID1 A, Front Signal  
Signal 2: FID2 B, Back Signal  
-----

-----  
Overview Table  
-----

RT	Sig	Lvl	Amount [g/100cc]	Area	Rsp.Factor	Ref	ISTD #	Compound
2.586	1	1	1.00000	3.69669	2.70512e-1	No	No	1 methanol
2.809	1	1	1.00000	4.26100	2.34687e-1	No	No	2 Acetaldehyde
2.977	2	1	1.00000	4.26100	2.34687e-1	No	No	2 Acetaldehyde
3.075	1	1	5.00000e-2	4.51944	1.10633e-2	No	No	1 ethanol
		2	1.00000e-1	9.08848	1.10029e-2			
		3	2.00000e-1	18.23830	1.09659e-2			
		4	3.00000e-1	27.17083	1.10413e-2			
		5	5.00000e-1	45.59841	1.09653e-2			
3.388	2	1	1.00000	4.26062	2.34707e-1	No	No	2 methanol
3.628	1	1	1.00000	9.73055	1.02769e-1	No	No	1 isopropyl alcohol
4.285	2	1	5.00000e-2	4.65226	1.07475e-2	No	No	2 ethanol
		2	1.00000e-1	9.36702	1.06758e-2			
		3	2.00000e-1	18.96443	1.05461e-2			
		4	3.00000e-1	28.50989	1.05227e-2			
		5	5.00000e-1	48.28486	1.03552e-2			
4.308	1	1	1.00000	6.49940	1.53860e-1	No	No	1 acetone
4.620	1	1	1.00000	46.43829	2.15340e-2	No	Yes	1 n-propanol
		2	1.00000	47.58935	2.10131e-2			
		3	1.00000	47.32145	2.11321e-2			
		4	1.00000	46.07795	2.17024e-2			
		5	1.00000	46.80014	2.13675e-2			
4.661	2	1	1.00000	6.89301	1.45075e-1	No	No	2 acetone
4.969	2	1	1.00000	10.70642	9.34019e-2	No	No	2 isopropyl alcohol
7.550	2	1	1.00000	48.38544	2.06674e-2	No	Yes	2 n-propanol
		2	1.00000	49.03825	2.03922e-2			
		3	1.00000	48.90527	2.04477e-2			
		4	1.00000	47.36019	2.11148e-2			
		5	1.00000	48.14449	2.07708e-2			

Peak Sum Table


\*\*\*No Entries in table\*\*\*

41 Warnings or Errors (10 first messages follow) :


- Warning : Curve requires more calibration points., (methanol)
- Warning : Curve requires more calibration points. at 2.586 min, signal 1
- Warning : Curve requires more calibration points. at 2.809 min, signal 1
- Warning : Curve requires more calibration points. at 2.977 min, signal 2
- Warning : Curve requires more calibration points. at 3.388 min, signal 2
- Warning : Curve requires more calibration points. at 3.628 min, signal 1
- Warning : Curve requires more calibration points. at 4.308 min, signal 1
- Warning : Curve requires more calibration points. at 4.62 min, signal 1
- Warning : Curve requires more calibration points. at 4.661 min, signal 2
- Warning : Curve requires more calibration points. at 4.969 min, signal 2

JG


=====  
 Calibration Curves  
 =====




methanol at exp. RT: 2.586  
 FID1 A, Front Signal  
 Correlation: 1.00000  
 Residual Std. Dev.: 0.00000  
 Formula:  $y = mx + b$   
 m: 7.96045e-2  
 b: 0.00000  
 x: Amount Ratio  
 y: Area Ratio



Acetaldehyde at exp. RT: 2.809  
 FID1 A, Front Signal  
 Correlation: 1.00000  
 Residual Std. Dev.: 0.00000  
 Formula:  $y = mx + b$   
 m: 8.80637e-2  
 b: 0.00000  
 x: Amount Ratio  
 y: Area Ratio




Acetaldehyde at exp. RT: 2.977  
 FID2 B, Back Signal  
 Correlation: 1.00000  
 Residual Std. Dev.: 0.00000  
 Formula:  $y = mx + b$   
 m: 8.80637e-2  
 b: 0.00000  
 x: Amount Ratio  
 y: Area Ratio




ethanol at exp. RT: 3.075  
 FID1 A, Front Signal  
 Correlation: 0.99995  
 Residual Std. Dev.: 0.00406  
 Formula:  $y = mx + b$   
 m: 1.95654  
 b: -2.46237e-3  
 x: Amount Ratio  
 y: Area Ratio


JG




methanol at exp. RT: 3.388  
 FID2 B, Back Signal  
 Correlation: 1.00000  
 Residual Std. Dev.: 0.00000  
 Formula:  $y = mx + b$   
 m: 8.80559e-2  
 b: 0.00000  
 x: Amount Ratio  
 y: Area Ratio



isopropyl alcohol at exp. RT: 3.628  
 FID1 A, Front Signal  
 Correlation: 1.00000  
 Residual Std. Dev.: 0.00000  
 Formula:  $y = mx + b$   
 m: 2.09537e-1  
 b: 0.00000  
 x: Amount Ratio  
 y: Area Ratio




ethanol at exp. RT: 4.285  
 FID2 B, Back Signal  
 Correlation: 0.99991  
 Residual Std. Dev.: 0.00573  
 Formula:  $y = mx + b$   
 m: 2.02464  
 b: -9.69879e-3  
 x: Amount Ratio  
 y: Area Ratio




acetone at exp. RT: 4.308  
 FID1 A, Front Signal  
 Correlation: 1.00000  
 Residual Std. Dev.: 0.00000  
 Formula:  $y = mx + b$   
 m: 1.39958e-1  
 b: 0.00000  
 x: Amount Ratio  
 y: Area Ratio


JG




n-propanol at exp. RT: 4.620  
FID1 A, Front Signal  
Correlation: 1.00000  
Residual Std. Dev.: 0.00000  
Formula:  $y = mx + b$   
m: 1.00000  
b: 0.00000  
x: Amount Ratio  
y: Area Ratio



acetone at exp. RT: 4.661  
FID2 B, Back Signal  
Correlation: 1.00000  
Residual Std. Dev.: 0.00000  
Formula:  $y = mx + b$   
m: 1.42460e-1  
b: 0.00000  
x: Amount Ratio  
y: Area Ratio



isopropyl alcohol at exp. RT: 4.969  
FID2 B, Back Signal  
Correlation: 1.00000  
Residual Std. Dev.: 0.00000  
Formula:  $y = mx + b$   
m: 2.21273e-1  
b: 0.00000  
x: Amount Ratio  
y: Area Ratio




n-propanol at exp. RT: 7.550  
FID2 B, Back Signal  
Correlation: 1.00000  
Residual Std. Dev.: 0.00000  
Formula:  $y = mx + b$   
m: 1.00000  
b: 0.00000  
x: Amount Ratio  
y: Area Ratio

=====

JK

ISP Forensic Services Blood Alcohol Report

Sample Name : INTERNAL STD BLK 1  
 Laboratory : Meridian  
 Injection Date : May 3, 2018  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	0.00000	0.0000	g/100cc
2.	Ethanol	Column 2:	0.00000	0.0000	g/100cc
3.	n-Propanol	Column 1:	31.69026	1.0000	g/100cc
4.	n-Propanol	Column 2:	32.65249	1.0000	g/100cc



ISP Forensic Services Blood Alcohol Report

Sample Name : MIX VOL FN06041503  
 Laboratory : Meridian  
 Injection Date : May 3, 2018  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	9.58423	0.1412	g/100cc
2.	Ethanol	Column 2:	9.90654	0.1410	g/100cc
3.	n-Propanol	Column 1:	35.00378	1.0000	g/100cc
4.	n-Propanol	Column 2:	35.91213	1.0000	g/100cc

# VOLATILES DETERMINATION CASEFILE WORKSHEET

**Laboratory No.:** QC1-1

**Analysis Date(s):** 03 May 2018

	Column 1 FID A	Column 2 FID B	Column Precision	Mean Value	Over-all Mean	
Sample Results	0.0764	0.0767	0.0003	0.0765	0.0771	
(g/100cc)	0.0775	0.0778	0.0003	0.0776		

**Analysis Method**

Refer to Blood Alcohol Method #1

**Instrument Information**

*Instrument method is stored centrally.*

Refer to Instrument Method: ALCOHOL.M  
Hamilton Auto-Dilutor Serial Number: MDL600HC11378

**Reporting of Results**

**Uncertainty of Measurement (UM%): 5.00%**

Overall Mean (g/100cc)	Low	High	5% of Mean
0.077	0.073	0.081	0.004

	<b>Reported Result</b>  0.077	
--	-------------------------------------	--

*Calibration and control data are stored centrally.*


Issued: 12/30/2016

Volatiles BAC Calculation Spreadsheet Rev 4

Issuing Authority: Quality Manager

ISP Forensic Services Blood Alcohol Report


Sample Name : QC1-1-A  
 Laboratory : Meridian  
 Injection Date : May 3, 2018  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	6.77133	0.0764	g/100cc
2.	Ethanol	Column 2:	6.91145	0.0767	g/100cc
3.	n-Propanol	Column 1:	46.08132	1.0000	g/100cc
4.	n-Propanol	Column 2:	47.45350	1.0000	g/100cc

ISP Forensic Services Blood Alcohol Report

Sample Name : QC1-1-B  
 Laboratory : Meridian  
 Injection Date : May 3, 2018  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	6.76506	0.0775	g/100cc
2.	Ethanol	Column 2:	6.87942	0.0778	g/100cc
3.	n-Propanol	Column 1:	45.38018	1.0000	g/100cc
4.	n-Propanol	Column 2:	46.52525	1.0000	g/100cc

## VOLATILES DETERMINATION CASEFILE WORKSHEET

**Laboratory No.: 0.08 FN10281510**

**Analysis Date(s): 03 May 2018**

	Column 1 FID A	Column 2 FID B	Column Precision	Mean Value	Over-all Mean	
Sample Results	0.0807	0.0807	0.0000	0.0807	0.0807	
(g/100cc)	0.0809	0.0808	0.0001	0.0808		

**Analysis Method**

Refer to Blood Alcohol Method #1

**Instrument Information**

*Instrument method is stored centrally.*

Refer to Instrument Method: ALCOHOL.M  
Hamilton Auto-Dilutor Serial Number: MDL600HC11378

**Reporting of Results**

**Uncertainty of Measurement (UM%): 5.00%**

Overall Mean (g/100cc)	Low	High	5% of Mean
0.080	0.076	0.084	0.004

	<b>Reported Result</b>	
	0.080	

*Calibration and control data are stored centrally.*


Issued: 12/30/2016

Volatiles BAC Calculation Spreadsheet Rev 4

Issuing Authority: Quality Manager

ISP Forensic Services Blood Alcohol Report


Sample Name : 0.08 FN10281510-A  
 Laboratory : Meridian  
 Injection Date : May 3, 2018  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	7.28977	0.0807	g/100cc
2.	Ethanol	Column 2:	7.39190	0.0807	g/100cc
3.	n-Propanol	Column 1:	46.89442	1.0000	g/100cc
4.	n-Propanol	Column 2:	48.06841	1.0000	g/100cc

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.08 FN10281510-B  
 Laboratory : Meridian  
 Injection Date : May 3, 2018  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	7.25302	0.0809	g/100cc
2.	Ethanol	Column 2:	7.34385	0.0808	g/100cc
3.	n-Propanol	Column 1:	46.54592	1.0000	g/100cc
4.	n-Propanol	Column 2:	47.73020	1.0000	g/100cc

JG

## VOLATILES DETERMINATION CASEFILE WORKSHEET

**Laboratory No.: QC2-1**

**Analysis Date(s): 03 May 2018**

	Column 1 FID A	Column 2 FID B	Column Precision	Mean Value	Over-all Mean	
Sample Results	0.2041	0.2046	0.0005	0.2043	0.2035	
(g/100cc)	0.2027	0.2029	0.0002	0.2028		

**Analysis Method**

Refer to Blood Alcohol Method #1

**Instrument Information**

*Instrument method is stored centrally.*

Refer to Instrument Method: ALCOHOL.M  
Hamilton Auto-Dilutor Serial Number: MDL600HC11378

**Reporting of Results**

**Uncertainty of Measurement (UM%): 5.00%**

Overall Mean (g/100cc)	Low	High	5% of Mean
0.203	0.192	0.214	0.011

Reported Result	
0.203	


*Calibration and control data are stored centrally.*

JG



ISP Forensic Services Blood Alcohol Report


Sample Name : QC2-1-A  
 Laboratory : Meridian  
 Injection Date : May 3, 2018  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	18.36788	0.2041	g/100cc
2.	Ethanol	Column 2:	19.13937	0.2046	g/100cc
3.	n-Propanol	Column 1:	46.28157	1.0000	g/100cc
4.	n-Propanol	Column 2:	47.30501	1.0000	g/100cc

ISP Forensic Services Blood Alcohol Report

Sample Name : QC2-1-B  
 Laboratory : Meridian  
 Injection Date : May 3, 2018  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	18.29845	0.2027	g/100cc
2.	Ethanol	Column 2:	19.05187	0.2029	g/100cc
3.	n-Propanol	Column 1:	46.43522	1.0000	g/100cc
4.	n-Propanol	Column 2:	47.48826	1.0000	g/100cc

## VOLATILES DETERMINATION CASEFILE WORKSHEET

**Laboratory No.:** QC1-2

**Analysis Date(s):** 03 May 2018

	Column 1 FID A	Column 2 FID B	Column Precision	Mean Value	Over-all Mean
Sample Results	0.0799	0.0807	0.0008	0.0803	0.0795
(g/100cc)	0.0785	0.0792	0.0007	0.0788	

### Analysis Method

Refer to Blood Alcohol Method #1

### Instrument Information

*Instrument method is stored centrally.*

Refer to Instrument Method: ALCOHOL.M  
Hamilton Auto-Dilutor Serial Number: MDL600HC11378

### Reporting of Results

**Uncertainty of Measurement (UM%): 5.00%**

Overall Mean (g/100cc)	Low	High	5% of Mean
0.079	0.075	0.083	0.004

	<b>Reported Result</b>	
	0.079	

*Calibration and control data are stored centrally.*


Issued: 12/30/2016

Volatiles BAC Calculation Spreadsheet Rev 4

Issuing Authority: Quality Manager

ISP Forensic Services Blood Alcohol Report


Sample Name : QC1-2-A  
 Laboratory : Meridian  
 Injection Date : May 3, 2018  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	7.15126	0.0799	g/100cc
2.	Ethanol	Column 2:	7.29636	0.0807	g/100cc
3.	n-Propanol	Column 1:	46.47025	1.0000	g/100cc
4.	n-Propanol	Column 2:	47.48112	1.0000	g/100cc

ISP Forensic Services Blood Alcohol Report

Sample Name : QC1-2-B  
 Laboratory : Meridian  
 Injection Date : May 3, 2018  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	6.86289	0.0785	g/100cc
2.	Ethanol	Column 2:	6.98576	0.0792	g/100cc
3.	n-Propanol	Column 1:	45.43435	1.0000	g/100cc
4.	n-Propanol	Column 2:	46.38639	1.0000	g/100cc

JK

ISP Forensic Services Blood Alcohol Report

Sample Name : INTERNAL STD BLK 1  
 Laboratory : Meridian  
 Injection Date : May 3, 2018  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	0.00000	0.0000	g/100cc
2.	Ethanol	Column 2:	0.00000	0.0000	g/100cc
3.	n-Propanol	Column 1:	43.70244	1.0000	g/100cc
4.	n-Propanol	Column 2:	44.56030	1.0000	g/100cc

S a m p l e S u m m a r y

Sequence table: C:\Chem32\1\Data\05-03-18\_SAMPLES\05-03-18\_SAMPLES 2018-05-03 11-16-44\05-03-18\_SAMPLES.S  
 Data directory path: C:\Chem32\1\Data\05-03-18\_SAMPLES\05-03-18\_SAMPLES 2018-05-03 11-16-44\  
 Logbook: C:\Chem32\1\Data\05-03-18\_SAMPLES\05-03-18\_SAMPLES 2018-05-03 11-16-44\05-03-18\_SAMPLES.LOG  
 Sequence start: 5/3/2018 11:31:30 AM  
 Sequence Operator: SYSTEM  
 Operator: SYSTEM  
 Method file name: C:\Chem32\1\Data\05-03-18\_SAMPLES\05-03-18\_SAMPLES 2018-05-03 11-16-44\ALCOHOL.M

Run #	Location #	Inj #	Sample Name	Sample Amt [g/100cc]	Multip.* Dilution	File name	Cal #	# Cmp
1	1	1	INTERNAL STD BLK	-	1.0000	001F0101.D		2
2	2	1	MIX VOL FN060415	-	1.0000	002F0201.D		10
3	3	1	QC1-1-A	-	1.0000	003F0301.D		4
4	4	1	QC1-1-B	-	1.0000	004F0401.D		4
5	5	1	0.08 FN10281510-	-	1.0000	005F0501.D		4
6	6	1	0.08 FN10281510-	-	1.0000	006F0601.D		4
7	7	1	M2018-2105-1-A	-	1.0000	007F0701.D		6
8	8	1	M2018-2105-1-B	-	1.0000	008F0801.D		6
9	9	1	M2018-2106-1-A	-	1.0000	009F0901.D		5
10	10	1	M2018-2106-1-B	-	1.0000	010F1001.D		6
11	11	1	M2018-2107-1-A	-	1.0000	011F1101.D		6
12	12	1	M2018-2107-1-B	-	1.0000	012F1201.D		6
13	13	1	M2018-2108-1-A	-	1.0000	013F1301.D		4
14	14	1	M2018-2108-1-B	-	1.0000	014F1401.D		4
15	15	1	M2018-2134-1-A	-	1.0000	015F1501.D		6
16	16	1	M2018-2134-1-B	-	1.0000	016F1601.D		6
17	17	1	M2018-2135-1-A	-	1.0000	017F1701.D		6
18	18	1	M2018-2135-1-B	-	1.0000	018F1801.D		6
19	19	1	M2018-2145-1-A	-	1.0000	019F1901.D		2
20	20	1	M2018-2145-1-B	-	1.0000	020F2001.D		3
21	21	1	M2018-2145-2-A	-	1.0000	021F2101.D		2
22	22	1	M2018-2145-2-B	-	1.0000	022F2201.D		2
23	23	1	M2018-2156-1-A	-	1.0000	023F2301.D		4
24	24	1	M2018-2156-1-B	-	1.0000	024F2401.D		4
25	25	1	QC2-1-A	-	1.0000	025F2501.D		4
26	26	1	QC2-1-B	-	1.0000	026F2601.D		4
27	27	1	M2018-2157-1-A	-	1.0000	027F2701.D		4
28	28	1	M2018-2157-1-B	-	1.0000	028F2801.D		5
29	29	1	M2018-2170-1-A	-	1.0000	029F2901.D		4
30	30	1	M2018-2170-1-B	-	1.0000	030F3001.D		4
31	31	1	M2018-2171-1-A	-	1.0000	031F3101.D		2
32	32	1	M2018-2171-1-B	-	1.0000	032F3201.D		3
33	33	1	M2018-2176-1-A	-	1.0000	033F3301.D		4
34	34	1	M2018-2176-1-B	-	1.0000	034F3401.D		6
35	35	1	M2018-2177-1-A	-	1.0000	035F3501.D		6
36	36	1	M2018-2177-1-B	-	1.0000	036F3601.D		6
37	37	1	M2018-2206-1-A	-	1.0000	037F3701.D		6
38	38	1	M2018-2206-1-B	-	1.0000	038F3801.D		6
39	39	1	M2018-2208-1-A	-	1.0000	039F3901.D		2
40	40	1	M2018-2208-1-B	-	1.0000	040F4001.D		2
41	41	1	M2018-2213-1-A	-	1.0000	041F4101.D		2
42	42	1	M2018-2213-1-B	-	1.0000	042F4201.D		2
43	43	1	M2018-2214-1-A	-	1.0000	043F4301.D		6

Run #	Location #	Inj #	Sample Name	Sample Amt [g/100cc]	Multip.* Dilution	File name	Cal #	# Cmp
44	44	1	M2018-2214-1-B	-	1.0000	044F4401.D	6	6
45	45	1	M2018-2217-1-A	-	1.0000	045F4501.D	6	6
46	46	1	M2018-2217-1-B	-	1.0000	046F4601.D	6	6
47	47	1	QC1-2-A	-	1.0000	047F4701.D	4	4
48	48	1	QC1-2-B	-	1.0000	048F4801.D	4	4
49	49	1	INTERNAL STD BLK	-	1.0000	049F4901.D	2	2

Method file name: C:\Chem32\1\Data\05-03-18\_SAMPLES\05-03-18\_SAMPLES 2018-05-03 11-16-44  
 \SHUTDOWN.M

Run #	Location #	Inj #	Sample Name	Sample Amt [g/100cc]	Multip.* Dilution	File name	Cal #	# Cmp
50	50	1	EMPTY	-	1.0000	050F5001.D	0	0