

Quantitative Analysis for Ethanol & Qualitative Analysis for Other Volatiles

Analytical Method(s): 1.0

Device: Hamilton MICROLAB 600 Liquid Processor/Dilutor Serial Number: ML600HC11378

Volatiles Quality Assurance Controls

Run Date: 12/4/18-12/5/18

Calibration Date: 11/27/18

Control level	Expiration	Lot #	Target Value	Acceptable Range	Overall Results
Level 1	Jan-22	1801036	0.0812	0.0731-0.0893	0.0801 g/100cc
					0.0839 g/100cc
					0.2024 g/100cc
Level 2	Mar-22	1803028	0.2035	0.1832-0.2238	g/100cc
					g/100cc
Multi-Component mixture:			Exp date: Sept. 2020	Lot #	FN06041502
Curve Fit:			Column 1	0.99999	Column2
					OK
					0.99995

Ethanol Calibration Reference Material								
Calibrator level	Expiration	Ceriliant Lot #	Target Value	Acceptable Range	Column 1	Column 2	Precision	Mean
0.050	Jul-19	FN06231406	0.050	0.045 - 0.055	0.0509	0.0516	0.0007	0.0512
0.080			0.080	0.072 - 0.088			0	#DIV/0!
0.100	Aug-21	FN08101601	0.100	0.090 - 0.110	0.0996	0.1005	0.0009	0.1
0.200	Dec-19	FN12011401	0.200	0.180 - 0.220	0.1993	0.1977	0.0016	0.1985
0.300	Feb-21	FN02121601	0.300	0.270 - 0.330	0.2997	0.2988	0.0009	0.2992
0.400			0.400	0.360 - 0.440			0	#DIV/0!
0.500	Sep-21	FN08031602	0.500	0.450 - 0.550	0.5004	0.5014	0.001	0.5009

Aqueous Controls					
Control level	Expiration	Ceriliant Lot #	Target Value	Acceptable Range	Overall Results
0.080	May-22	FN04171701	0.08000	0.076 - 0.084	0.080 g/100cc

Issued: 4/22/2015

~Any information on this document can be changed for laboratory use, except for the precision and mean determination formulas.


Volatiles QA/QC data spreadsheet Rev 5
Issuing Authority: Quality Manager

MB

26


Worklist: 2821

<u>LAB_CASE</u>	<u>ITEM</u>	<u>TASK_ID</u>	<u>DESCRIPTION</u>
M2018-5847	1	132799	Alcohol Analysis
M2018-5848	1	132800	Alcohol Analysis
M2018-5858	1	132825	Alcohol Analysis
M2018-5864	1	132867	Alcohol Analysis
M2018-5868	2	132915	Alcohol Analysis
M2018-5869	1	132916	Alcohol Analysis
M2018-5908	2	133075	Alcohol Analysis
M2018-5910	3	133082	Alcohol Analysis
M2018-5914	1	133119	Alcohol Analysis
M2018-5915	1	133120	Alcohol Analysis
M2018-5916	1	133121	Alcohol Analysis
M2018-5927	1	133204	Alcohol Analysis
M2018-5928	1	133205	Alcohol Analysis
M2018-5936	1	133253	Alcohol Analysis
M2018-5953	1	133406	Alcohol Analysis
M2018-5954	1	133407	Alcohol Analysis
P2018-3284	1	132510	Alcohol Analysis
P2018-3295	1	132533	Alcohol Analysis
P2018-3296	1	132534	Alcohol Analysis



ISP Forensic Services Blood Alcohol Report


Sample Name : INTERNAL STD BLK 1
 Laboratory : Meridian
 Injection Date : Dec 4, 2018
 Method : ALCOHOL.M
 Acq. Instrument: CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	0.00000	0.0000	g/100cc
2.	Ethanol	Column 2:	0.00000	0.0000	g/100cc
3.	n-Propanol	Column 1:	48.57389	1.0000	g/100cc
4.	n-Propanol	Column 2:	50.96790	1.0000	g/100cc

ISP Forensic Services Blood Alcohol Report

Sample Name : MIX VOL FN06041502
 Laboratory : Meridian
 Injection Date : Dec 4, 2018
 Method : ALCOHOL.M
 Acq. Instrument: CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	10.11729	0.1402	g/100cc
2.	Ethanol	Column 2:	10.45784	0.1398	g/100cc
3.	n-Propanol	Column 1:	38.84338	1.0000	g/100cc
4.	n-Propanol	Column 2:	40.05091	1.0000	g/100cc

VOLATILES DETERMINATION CASEFILE WORKSHEET

Laboratory No.: QC1-1

Analysis Date(s): 04 Dec 2018

	Column 1 FID A	Column 2 FID B	Column Precision	Mean Value	Over-all Mean	
Sample Results	0.0801	0.0805	0.0004	0.0803	0.0801	
(g/100cc)	0.0800	0.0799	0.0001	0.0799		

Analysis Method

Refer to Blood Alcohol Method #1

Instrument Information

Instrument method is stored centrally.

Refer to Instrument Method: ALCOHOL.M
Hamilton Auto-Dilutor Serial Number: ML600HC11378

Reporting of Results

Uncertainty of Measurement (UM%): 5.00%

Overall Mean (g/100cc)	Low	High	5% of Mean
0.080	0.076	0.084	0.004

	Reported Result	
	0.080	

Calibration and control data are stored centrally.


Issued: 12/30/2016

Volatiles BAC Calculation Spreadsheet Rev 4

Issuing Authority: Quality Manager

ISP Forensic Services Blood Alcohol Report


Sample Name : QC1-1-A
 Laboratory : Meridian
 Injection Date : Dec 4, 2018
 Method : ALCOHOL.M
 Acq. Instrument: CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	7.34528	0.0801	g/100cc
2.	Ethanol	Column 2:	7.52096	0.0805	g/100cc
3.	n-Propanol	Column 1:	49.47753	1.0000	g/100cc
4.	n-Propanol	Column 2:	51.17125	1.0000	g/100cc

ISP Forensic Services Blood Alcohol Report

Sample Name : QC1-1-B
 Laboratory : Meridian
 Injection Date : Dec 4, 2018
 Method : ALCOHOL.M
 Acq. Instrument: CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	7.20236	0.0800	g/100cc
2.	Ethanol	Column 2:	7.34634	0.0799	g/100cc
3.	n-Propanol	Column 1:	48.61231	1.0000	g/100cc
4.	n-Propanol	Column 2:	50.34078	1.0000	g/100cc

VOLATILES DETERMINATION CASEFILE WORKSHEET

Laboratory No.: 0.08 FN04171701

Analysis Date(s): 04 Dec 2018

	Column 1 FID A	Column 2 FID B	Column Precision	Mean Value	Over-all Mean	
Sample Results	0.0799	0.0800	0.0001	0.0799	0.0800	
(g/100cc)	0.0801	0.0803	0.0002	0.0802		

Analysis Method

Refer to Blood Alcohol Method #1

Instrument Information

Instrument method is stored centrally.

Refer to Instrument Method: ALCOHOL.M
Hamilton Auto-Dilutor Serial Number: ML600HC11378

Reporting of Results

Uncertainty of Measurement (UM%): 5.00%

Overall Mean (g/100cc)	Low	High	5% of Mean
0.080	0.076	0.084	0.004

	Reported Result	
	0.080	

Calibration and control data are stored centrally.


Issued: 12/30/2016

Volatiles BAC Calculation Spreadsheet Rev 4

Issuing Authority: Quality Manager

ISP Forensic Services Blood Alcohol Report


Sample Name : 0.08 FN04171701-A
 Laboratory : Meridian
 Injection Date : Dec 4, 2018
 Method : ALCOHOL.M
 Acq. Instrument: CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	7.24585	0.0799	g/100cc
2.	Ethanol	Column 2:	7.38853	0.0800	g/100cc
3.	n-Propanol	Column 1:	48.94521	1.0000	g/100cc
4.	n-Propanol	Column 2:	50.57655	1.0000	g/100cc

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.08 FN04171701-B
 Laboratory : Meridian
 Injection Date : Dec 4, 2018
 Method : ALCOHOL.M
 Acq. Instrument: CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	7.17747	0.0801	g/100cc
2.	Ethanol	Column 2:	7.34524	0.0803	g/100cc
3.	n-Propanol	Column 1:	48.36038	1.0000	g/100cc
4.	n-Propanol	Column 2:	50.11644	1.0000	g/100cc

VOLATILES DETERMINATION CASEFILE WORKSHEET

Laboratory No.: QC2-1

Analysis Date(s): 04 Dec 2018

	Column 1 FID A	Column 2 FID B	Column Precision	Mean Value	Over-all Mean	
Sample Results	0.2029	0.2027	0.0002	0.2028	0.2024	
(g/100cc)	0.2023	0.2020	0.0003	0.2021		

Analysis Method

Refer to Blood Alcohol Method #1

Instrument Information

Instrument method is stored centrally.

Refer to Instrument Method: ALCOHOL.M
Hamilton Auto-Dilutor Serial Number: ML600HC11378

Reporting of Results

Uncertainty of Measurement (UM%): 5.00%


Overall Mean (g/100cc)	Low	High	5% of Mean
0.202	0.191	0.213	0.011

	Reported Result	
	0.202	

Calibration and control data are stored centrally.

ISP Forensic Services Blood Alcohol Report


Sample Name : QC2-1-A
 Laboratory : Meridian
 Injection Date : Dec 4, 2018
 Method : ALCOHOL.M
 Acq. Instrument: CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	17.99335	0.2029	g/100cc
2.	Ethanol	Column 2:	18.64326	0.2027	g/100cc
3.	n-Propanol	Column 1:	47.67779	1.0000	g/100cc
4.	n-Propanol	Column 2:	48.79485	1.0000	g/100cc

ISP Forensic Services Blood Alcohol Report

Sample Name : QC2-1-B
 Laboratory : Meridian
 Injection Date : Dec 4, 2018
 Method : ALCOHOL.M
 Acq. Instrument: CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	18.14673	0.2023	g/100cc
2.	Ethanol	Column 2:	18.80096	0.2020	g/100cc
3.	n-Propanol	Column 1:	48.24213	1.0000	g/100cc
4.	n-Propanol	Column 2:	49.38381	1.0000	g/100cc

VOLATILES DETERMINATION CASEFILE WORKSHEET

Laboratory No.: QC1-2

Analysis Date(s): 04 Dec 2018

	Column 1 FID A	Column 2 FID B	Column Precision	Mean Value	Over-all Mean	
Sample Results	0.0848	0.0850	0.0002	0.0849	0.0839	
(g/100cc)	0.0830	0.0830	0.0000	0.0830		

Analysis Method

Refer to Blood Alcohol Method #1

Instrument Information

Instrument method is stored centrally.

Refer to Instrument Method: ALCOHOL.M
Hamilton Auto-Dilutor Serial Number: ML600HC11378

Reporting of Results

Uncertainty of Measurement (UM%): 5.00%

Overall Mean (g/100cc)	Low	High	5% of Mean
0.083	0.078	0.088	0.005

	Reported Result	
	0.083	

Calibration and control data are stored centrally.


Issued: 12/30/2016

Volatiles BAC Calculation Spreadsheet Rev 4

Issuing Authority: Quality Manager

ISP Forensic Services Blood Alcohol Report


Sample Name : QC1-2-A
 Laboratory : Meridian
 Injection Date : Dec 4, 2018
 Method : ALCOHOL.M
 Acq. Instrument: CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	7.58705	0.0848	g/100cc
2.	Ethanol	Column 2:	7.69225	0.0850	g/100cc
3.	n-Propanol	Column 1:	48.26966	1.0000	g/100cc
4.	n-Propanol	Column 2:	49.41397	1.0000	g/100cc

ISP Forensic Services Blood Alcohol Report


Sample Name : QC1-2-B
 Laboratory : Meridian
 Injection Date : Dec 4, 2018
 Method : ALCOHOL.M
 Acq. Instrument: CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	7.43099	0.0830	g/100cc
2.	Ethanol	Column 2:	7.51779	0.0830	g/100cc
3.	n-Propanol	Column 1:	48.31082	1.0000	g/100cc
4.	n-Propanol	Column 2:	49.52483	1.0000	g/100cc

ISP Forensic Services Blood Alcohol Report

Sample Name : INTERNAL STD BLK
 Laboratory : Meridian
 Injection Date : Dec 4, 2018
 Method : ALCOHOL.M
 Acq. Instrument: CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	0.00000	0.0000	g/100cc
2.	Ethanol	Column 2:	0.00000	0.0000	g/100cc
3.	n-Propanol	Column 1:	46.61297	1.0000	g/100cc
4.	n-Propanol	Column 2:	47.47593	1.0000	g/100cc

S a m p l e S u m m a r y

Sequence table: C:\Chem32\1\Data\12-04-18_SAMPLES\12-04-18_SAMPLES 2018-12-04 15-10-02\12-04-18_SAMPLES.S
 Data directory path: C:\Chem32\1\Data\12-04-18_SAMPLES\12-04-18_SAMPLES 2018-12-04 15-10-02\
 Logbook: C:\Chem32\1\Data\12-04-18_SAMPLES\12-04-18_SAMPLES 2018-12-04 15-10-02\12-04-18_SAMPLES.LOG
 Sequence start: 12/4/2018 3:24:51 PM
 Sequence Operator: SYSTEM
 Operator: SYSTEM
 Method file name: C:\Chem32\1\Data\12-04-18_SAMPLES\12-04-18_SAMPLES 2018-12-04 15-10-02\ALCOHOL.M

Run #	Location #	Inj #	Sample Name	Sample Amt [g/100cc]	Multip.* Dilution	File name	Cal # Cmp
1	1	1	INTERNAL STD BLK	-	1.0000	001F0101.D	2
2	2	1	MIX VOL FN060415	-	1.0000	002F0201.D	10
3	3	1	QC1-1-A	-	1.0000	003F0301.D	4
4	4	1	QC1-1-B	-	1.0000	004F0401.D	4
5	5	1	0.08 FN04171701-	-	1.0000	005F0501.D	4
6	6	1	0.08 FN04171701-	-	1.0000	006F0601.D	4
7	7	1	M2018-5847-1-A	-	1.0000	007F0701.D	4
8	8	1	M2018-5847-1-B	-	1.0000	008F0801.D	4
9	9	1	M2018-5848-1-A	-	1.0000	009F0901.D	4
10	10	1	M2018-5848-1-B	-	1.0000	010F1001.D	4
11	11	1	M2018-5858-1-A	-	1.0000	011F1101.D	4
12	12	1	M2018-5858-1-B	-	1.0000	012F1201.D	4
13	13	1	M2018-5864-1-A	-	1.0000	013F1301.D	4
14	14	1	M2018-5864-1-B	-	1.0000	014F1401.D	4
15	15	1	M2018-5868-2-A	-	1.0000	015F1501.D	2
16	16	1	M2018-5868-2-B	-	1.0000	016F1601.D	2
17	17	1	M2018-5869-1-A	-	1.0000	017F1701.D	4
18	18	1	M2018-5869-1-B	-	1.0000	018F1801.D	4
19	19	1	M2018-5908-2-A	-	1.0000	019F1901.D	2
20	20	1	M2018-5908-2-B	-	1.0000	020F2001.D	2
21	21	1	M2018-5910-3-A	-	1.0000	021F2101.D	2
22	22	1	M2018-5910-3-B	-	1.0000	022F2201.D	2
23	23	1	M2018-5914-1-A	-	1.0000	023F2301.D	4
24	24	1	M2018-5914-1-B	-	1.0000	024F2401.D	4
25	25	1	QC2-1-A	-	1.0000	025F2501.D	4
26	26	1	QC2-1-B	-	1.0000	026F2601.D	4
27	27	1	M2018-5915-1-A	-	1.0000	027F2701.D	4
28	28	1	M2018-5915-1-B	-	1.0000	028F2801.D	4
29	29	1	M2018-5916-1-A	-	1.0000	029F2901.D	4
30	30	1	M2018-5916-1-B	-	1.0000	030F3001.D	4
31	31	1	M2018-5927-1-A	-	1.0000	031F3101.D	2
32	32	1	M2018-5917-1-B <i>S927-1-B</i>	-	1.0000	032F3201.D	2
33	33	1	M2018-5928-1-A <i>OK</i>	-	1.0000	033F3301.D	4
34	34	1	M2018-5928-1-B	-	1.0000	034F3401.D	4
35	35	1	M2018-5936-1-A	-	1.0000	035F3501.D	4
36	36	1	M2018-5936-1-B	-	1.0000	036F3601.D	4
37	37	1	M2018-5953-1-A	-	1.0000	037F3701.D	4
38	38	1	M2018-5953-1-B	-	1.0000	038F3801.D	4
39	39	1	M2018-5954-1-A	-	1.0000	039F3901.D	4
40	40	1	M2018-5954-1-B	-	1.0000	040F4001.D	4
41	41	1	P2018-3284-1-A	-	1.0000	041F4101.D	3
42	42	1	P2018-3284-1-B	-	1.0000	042F4201.D	4
43	43	1	P2018-3295-1-A	-	1.0000	043F4301.D	2


Run #	Location #	Inj #	Sample Name	Sample Amt [g/100cc]	Multip.* Dilution	File name	Cal #	# Cmp
44	44	1	P2018-3295-1-B	-	1.0000	044F4401.D		2
45	45	1	P2018-3296-1-A	-	1.0000	045F4501.D		2
46	46	1	P2018-3296-1-B	-	1.0000	046F4601.D		2
47	47	1	QC1-2-A	-	1.0000	047F4701.D		4
48	48	1	QC1-2-B	-	1.0000	048F4801.D		4
49	49	1	INTERNAL STD BLK	-	1.0000	049F4901.D		2

Method file name: C:\Chem32\1\Data\12-04-18_SAMPLES\12-04-18_SAMPLES 2018-12-04 15-10-02 \SHUTDOWN.M

Run #	Location #	Inj #	Sample Name	Sample Amt [g/100cc]	Multip.* Dilution	File name	Cal #	# Cmp
50	50	1	EMPTY	-	1.0000	050F5001.D		0

ISP Forensic Services Blood Alcohol Report


Sample Name : INTERNAL STD BLK 1
 Laboratory : Meridian
 Injection Date : Dec 5, 2018
 Method : ALCOHOL.M
 Acq. Instrument: CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	0.00000	0.0000	g/100cc
2.	Ethanol	Column 2:	0.00000	0.0000	g/100cc
3.	n-Propanol	Column 1:	48.50282	1.0000	g/100cc
4.	n-Propanol	Column 2:	50.30000	1.0000	g/100cc

ISP Forensic Services Blood Alcohol Report

Sample Name : TFE 111914
 Laboratory : Meridian
 Injection Date : Dec 5, 2018
 Method : ALCOHOL.M
 Acq. Instrument: CN11180014-CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	0.00000	0.0000	g/100cc
2.	Ethanol	Column 2:	0.00000	0.0000	g/100cc
3.	n-Propanol	Column 1:	45.88667	1.0000	g/100cc
4.	n-Propanol	Column 2:	47.22564	1.0000	g/100cc

Jc

ISP Forensic Services Blood Alcohol Report


Sample Name : INTERNAL STD BLK 1
 Laboratory : Meridian
 Injection Date : Dec 5, 2018
 Method : ALCOHOL.M
 Acq. Instrument: CN11180014 -CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	0.00000	0.0000	g/100cc
2.	Ethanol	Column 2:	0.00000	0.0000	g/100cc
3.	n-Propanol	Column 1:	41.75641	1.0000	g/100cc
4.	n-Propanol	Column 2:	42.69651	1.0000	g/100cc

ISP Forensic Services Blood Alcohol Report


Sample Name : DFE 111914OM
 Laboratory : Meridian
 Injection Date : Dec 5, 2018
 Method : ALCOHOL.M
 Acq. Instrument: CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	0.00000	0.0000	g/100cc
2.	Ethanol	Column 2:	0.00000	0.0000	g/100cc
3.	n-Propanol	Column 1:	46.99356	1.0000	g/100cc
4.	n-Propanol	Column 2:	48.12229	1.0000	g/100cc

ISP Forensic Services Blood Alcohol Report

Sample Name : INTERNAL STD BLK
 Laboratory : Meridian
 Injection Date : Dec 5, 2018
 Method : ALCOHOL.M
 Acq. Instrument: CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	0.00000	0.0000	g/100cc
2.	Ethanol	Column 2:	0.00000	0.0000	g/100cc
3.	n-Propanol	Column 1:	41.68595	1.0000	g/100cc
4.	n-Propanol	Column 2:	42.36742	1.0000	g/100cc

S a m p l e s u m m a r y

Sequence table: C:\Chem32\12-05-18_inhale\12-05-18_inhale 2018-12-05 08-16-27\12-05-18_inhale.S
 Data directory path: C:\Chem32\12-05-18_inhale\12-05-18_inhale 2018-12-05 08-16-27\
 Logbook: C:\Chem32\12-05-18_inhale\12-05-18_inhale 2018-12-05 08-16-27\12-05-18_inhale.LOG
 Sequence start: 12/5/2018 8:31:05 AM
 Sequence Operator: SYSTEM
 Operator: SYSTEM

Method file name: C:\Chem32\12-05-18_inhale\12-05-18_inhale 2018-12-05 08-16-27\ALCOHOL.M

Run Location Inj	Sample Name	Sample Amt	Multip.*	File name	Cal #	Cmp
1 1	INTERNAL STD BLK	1.0000	001F0101.D	2		
2 2	1 TFE 111914	1.0000	002F0201.D	2		
3 3	1 INTERNAL STD BLK	1.0000	003F0301.D	2		
4 4	1 DFE 111914OM	1.0000	004F0401.D	2		
5 5	1 INTERNAL STD BLK	1.0000	005F0501.D	2		

Method file name: C:\Chem32\12-05-18_inhale\12-05-18_inhale 2018-12-05 08-16-27\SHUTDOWN.M

Run Location Inj	Sample Name	Sample Amt	Multip.*	File name	Cal #	Cmp
6 6	1 EMPTY	-	1.0000 006F0601.D	0		

26

26

Overview Table

Signal 1: FID1 A, Front Signal
Signal 2: FID2 B, Back Signal

Signal Details

2	1.00000	n-propanol
1	1.00000	n-propanol

Default Sample ISTD Information (if not set in sample table):
ISTD ISTD Amount Name

Calibration Report Options :
Printout of recalibrations within a sequence:
Calibration Table after Recalibration
Normal Report after Recalibration
If the sequence is done with bracketing:
Results of first cycle (ending previous bracket)

Recalibration Settings:
Average Response :
Average Retention Time:
Average all recalibrations
Floating Average New 75%

Curve Type :
Origin :
Weight :
Linear
Ignored
Equal

Rel. Reference Window : 0.000 %
Abs. Reference Window : 0.100 min
Rel. Non-ref. Window : 0.000 %
Abs. Non-ref. Window : 0.100 min
Uncalibrated Peaks : not reported
Partial Calibration : Yes, identified peaks are recalibrated
Correct All Ret. Times : No, only for identified peaks

Calib. Data Modified : Tuesday, November 27, 2018 4:11:30 PM
Signals calculated separately : NO

General Calibration Setting

=====
Calibration Table
=====

RT	Stg	Lvl	Amount	Area	Rsp. Factor	Ref	ISTD	#	Compound
2.586	1	1	1.00000	3.69669	2.70512e-1	No	No	1	methanol
2.809	1	1	1.00000	4.26100	2.34687e-1	No	No	2	Acetaldehyde
2.977	2	1	1.00000	4.26100	2.34687e-1	No	No	2	Acetaldehyde
3.075	1	1	5.00000e-2	4.56194	1.09602e-2	No	No	1	ethanol
3.388	2	1	1.00000	4.26062	2.34707e-1	No	No	2	methanol
3.628	1	1	1.00000	9.73055	1.02769e-1	No	No	1	isopropyl alcohol
4.285	2	1	5.00000e-2	4.60882	1.08488e-2	No	No	2	ethanol
4.308	1	1	1.00000	6.49940	1.53860e-1	No	No	1	acetone
4.620	1	1	1.00000	48.52174	2.06093e-2	No	No	1	n-propanol
4.661	2	1	1.00000	6.89301	1.45075e-1	No	No	2	acetone
4.969	2	1	1.00000	10.70642	9.34019e-2	No	No	2	isopropyl alcohol
7.550	2	1	1.00000	50.37904	1.98495e-2	No	No	2	n-propanol
	2	2	1.00000	50.89652	1.96477e-2				
	3	3	1.00000	50.80914	1.96815e-2				
	4	4	1.00000	51.05636	1.95862e-2				
	5	5	1.00000	50.53856	1.97869e-2				

No Entries in table


1 Warnings or Errors :

Warning : Curve requires more calibration points., (methanol)





Peak Sum Table

Calibration Curves


methanol at exp. RT: 2.586
 FID1 A, Front Signal
 Correlation: 1.00000
 Residual Std. Dev.: 0.00000
 Formula: $y = mx + b$
 m: 7.61864e-2
 b: 0.00000
 x: Amount Ratio
 y: Area Ratio




26




26




isopropyl alcohol at exp. RT: 3.628
 FID1 A, Front Signal
 Correlation: 1.00000
 Residual Std. Dev.: 0.00000
 Formula: $y = mx + b$
 m: 2.00540e-1
 b: 0.00000
 x: Amount Ratio
 y: Area Ratio



ethanol at exp. RT: 4.285
 FID2 B, Back Signal
 Correlation: 0.99995
 Residual Std. Dev.: 0.00380
 Formula: $y = mx + b$
 m: 1.92377
 b: -7.86251e-3
 x: Amount Ratio
 y: Area Ratio




acetone at exp. RT: 4.308
 FID1 A, Front Signal
 Correlation: 1.00000
 Residual Std. Dev.: 0.00000
 Formula: $y = mx + b$
 m: 1.33948e-1
 b: 0.00000
 x: Amount Ratio
 y: Area Ratio




n-propanol at exp. RT: 4.620
 FID1 A, Front Signal
 Correlation: 1.00000
 Residual Std. Dev.: 0.00000
 Formula: $y = mx + b$
 m: 1.00000
 b: 0.00000
 x: Amount Ratio
 y: Area Ratio


36



acetone at exp. RT: 4.661
FID2 B, Back Signal
Correlation: 1.00000
Residual Std. Dev.: 0.00000
Formula: $y = mx + b$
m: 1.36823e-1
b: 0.00000
x: Amount Ratio
y: Area Ratio



isopropyl alcohol at exp. RT: 4.969
FID2 B, Back Signal
Correlation: 1.00000
Residual Std. Dev.: 0.00000
Formula: $y = mx + b$
m: 2.12517e-1
b: 0.00000
x: Amount Ratio
y: Area Ratio




n-propanol at exp. RT: 7.550
FID2 B, Back Signal
Correlation: 1.00000
Residual Std. Dev.: 0.00000
Formula: $y = mx + b$
m: 1.00000
b: 0.00000
x: Amount Ratio
y: Area Ratio

=====

JG

ISP Forensic Services Blood Alcohol Report


Sample Name : 0.050 FN06231406
 Laboratory : Meridian
 Injection Date : Nov 27, 2018
 Method : ALCOHOL.M
 Acq. Instrument: CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	4.56194	0.0509	g/100cc
2.	Ethanol	Column 2:	4.60882	0.0516	g/100cc
3.	n-Propanol	Column 1:	48.52174	1.0000	g/100cc
4.	n-Propanol	Column 2:	50.37904	1.0000	g/100cc

ISP Forensic Services Blood Alcohol Report


Sample Name : 0.100 FN08101601
 Laboratory : Meridian
 Injection Date : Nov 27, 2018
 Method : ALCOHOL.M
 Acq. Instrument: CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	9.10056	0.0996	g/100cc
2.	Ethanol	Column 2:	9.43937	0.1005	g/100cc
3.	n-Propanol	Column 1:	49.23827	1.0000	g/100cc
4.	n-Propanol	Column 2:	50.89652	1.0000	g/100cc

ISP Forensic Services Blood Alcohol Report


Sample Name : 0.200 FN12011401
 Laboratory : Meridian
 Injection Date : Nov 27, 2018
 Method : ALCOHOL.M
 Acq. Instrument: CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	18.23166	0.1993	g/100cc
2.	Ethanol	Column 2:	18.92577	0.1977	g/100cc
3.	n-Propanol	Column 1:	49.18580	1.0000	g/100cc
4.	n-Propanol	Column 2:	50.80914	1.0000	g/100cc

ISP Forensic Services Blood Alcohol Report


Sample Name : 0.300 FN02121601
 Laboratory : Meridian
 Injection Date : Nov 27, 2018
 Method : ALCOHOL.M
 Acq. Instrument: CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	27.76342	0.2997	g/100cc
2.	Ethanol	Column 2:	28.94269	0.2988	g/100cc
3.	n-Propanol	Column 1:	49.76912	1.0000	g/100cc
4.	n-Propanol	Column 2:	51.05636	1.0000	g/100cc

ISP Forensic Services Blood Alcohol Report


Sample Name : 0.500 FN08031602
 Laboratory : Meridian
 Injection Date : Nov 27, 2018
 Method : ALCOHOL.M
 Acq. Instrument: CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	45.93099	0.5004	g/100cc
2.	Ethanol	Column 2:	48.35098	0.5014	g/100cc
3.	n-Propanol	Column 1:	49.28674	1.0000	g/100cc
4.	n-Propanol	Column 2:	50.53856	1.0000	g/100cc

ISP Forensic Services Blood Alcohol Report

Sample Name : INTERNAL STANDARD BLANK
 Laboratory : Meridian
 Injection Date : Nov 27, 2018
 Method : ALCOHOL.M
 Acq. Instrument: CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	0.00000	0.0000	g/100cc
2.	Ethanol	Column 2:	0.00000	0.0000	g/100cc
3.	n-Propanol	Column 1:	43.35644	1.0000	g/100cc
4.	n-Propanol	Column 2:	44.46200	1.0000	g/100cc

S a m p l e S u m m a r y

Sequence table: C:\Chem32\1\Data\11-27-18_CAL\11-28-18_CAL 2018-11-27 15-06-22\11-28-18_CAL.S
 Data directory path: C:\Chem32\1\Data\11-27-18_CAL\11-28-18_CAL 2018-11-27 15-06-22\
 Logbook: C:\Chem32\1\Data\11-27-18_CAL\11-28-18_CAL 2018-11-27 15-06-22\11-28-18_CAL.LOG
 Sequence start: 11/27/2018 3:20:59 PM
 Sequence Operator: SYSTEM
 Operator: SYSTEM
 Method file name: C:\Chem32\1\Data\11-27-18_CAL\11-28-18_CAL 2018-11-27 15-06-22\ALCOHOL.M

Run #	Location #	Inj #	Sample Name	Sample Amt [g/100cc]	Multip.* Dilution	File name	Cal #	# Cmp
1	1	1	0.050 FN06231406	-	1.0000	001F0101.D	*	4
2	2	1	0.100 FN08101601	-	1.0000	002F0201.D	*	4
3	3	1	0.200 FN12011401	-	1.0000	003F0301.D	*	4
4	4	1	0.300 FN02121601	-	1.0000	004F0401.D	*	4
5	5	1	0.500 FN08031602	-	1.0000	005F0501.D	*	4
6	6	1	INTERNAL STANDAR	-	1.0000	006F0601.D		2

JK