

# REVIEWED


By Melissa (Nikka) Bradley at 3:13 pm, Oct 17, 2019

MB

10/15/2019

## Worklist: 3756

<u>LAB CASE</u>	<u>ITEM</u>	<u>ITEM TYPE</u>	<u>DESCRIPTION</u>
C2019-1752	1	BCK	Alcohol Analysis
C2019-1752	2	BLOOD	Alcohol Analysis
C2019-1752	3	BLOOD	Alcohol Analysis
C2019-1752	4	BLOOD	Alcohol Analysis
C2019-1801	1	BCK	Alcohol Analysis
C2019-1829	1	BCK	Alcohol Analysis
C2019-1840	1	BCK	Alcohol Analysis
C2019-1844	1	BCK	Alcohol Analysis
C2019-1863	1	BCK	Alcohol Analysis
C2019-1889	1	AVK	Alcohol Analysis
C2019-1934	1	BCK	Alcohol Analysis
C2019-1957	1	BCK	Alcohol Analysis
C2019-1968	1	BCK	Alcohol Analysis
C2019-1970	2	BCK	Alcohol Analysis
C2019-1973	1	AVK	Alcohol Analysis
C2019-1976	1	BCK	Alcohol Analysis



JA

Quantitative Analysis for Ethanol & Qualitative Analysis for Other Volatiles

Analytical Method(s): 1.0

Device: Hamilton MICROLAB 600A Liquid Processor/Dilutor Serial Number: ML600HC11379

Volatiles Quality Assurance Controls

Run Date(s): 10/15/19

Control Level	Expiration	Lot #	Target Value	Acceptable Range	Overall Results
Level 1	Jan-22	1801036	0.0812	0.0731-0.0893	0.0789 g/100cc
					g/100cc
Level 2	Mar-22	1803028	0.2035	0.1832-0.2238	0.1960 g/100cc
					0.1988 g/100cc
Multi-Component mixture:			Sep-20	Lot # FN06041502	OK
Curve Fit:			Column 1	0.99999	Column 2 0.99999

Ethanol Calibration Reference Material						
Calibrator level	Target Value	Acceptable Range	Column 1	Column 2	Precision	Mean
50	0.050	0.045 - 0.055	0.0498	0.0495	0.0003	0.0496
100	0.100	0.090 - 0.110	0.0995	0.0984	0.0011	0.0989
200	0.200	0.180 - 0.220	0.1986	0.1983	0.0003	0.1984
300	0.300	0.270 - 0.330	0.2990	0.2995	0.0005	0.2992
500	0.500	0.450 - 0.550	0.5013	0.5013	0.0000	0.5013

Aqueous Controls				
Control level	Target Value	Acceptable Range	Overall Results	
80	0.080	0.076 - 0.084	0.080	g/100cc

Revision: 1

Issue Date: 01/03/2019

Sample Summary

Sequence table: C:\Chem32\1\TEMP\AESEQ\QS\_15.10.2019\_04.05.54\10-15-2019.S  
 Data directory path: C:\Chem32\1\Data\10-15-2019-JJ  
 Logbook: C:\Chem32\1\Data\10-15-2019-JJ\10-15-2019.LOG  
 Sequence start: 10/15/2019 4:19:40 PM  
 Sequence Operator: SYSTEM  
 Operator: SYSTEM

Method file name: C:\CHEM32\1\METHODS\ALCOHOL.M

Run #	Location #	Inj #	Sample Name	Sample Amt [g/100cc]	Multip. * Dilution	File name	Cal # Cmp
1	1	1	water-1	-	1.0000	001F0101.D	0
2	2	1	VOL MIX FN-06041	-	1.0000	002F0201.D	10
3	3	1	ISTD BLANK-1	-	1.0000	003F0301.D	2
4	4	1	QC-2(1)-A	-	1.0000	004F0401.D	4
5	5	1	QC-2(1)-B	-	1.0000	005F0501.D	4
6	6	1	0.08 FN04171701-	-	1.0000	006F0601.D	4
7	7	1	0.08 FN04171701-	-	1.0000	007F0701.D	4
8	8	1	19106-63-1-A	-	1.0000	008F0801.D	4
9	9	1	19106-63-1-B	-	1.0000	009F0901.D	4
10	10	1	19106-116-1-A	-	1.0000	010F1001.D	4
11	11	1	19106-116-1-B	-	1.0000	011F1101.D	4
12	12	1	C2019-1752-1-A	-	1.0000	012F1201.D	4
13	13	1	C2019-1752-1-B	-	1.0000	013F1301.D	4
14	14	1	C2019-1752-2-A	-	1.0000	014F1401.D	4
15	15	1	C2019-1752-2-B	-	1.0000	015F1501.D	4
16	16	1	C2019-1752-3-A	-	1.0000	016F1601.D	4
17	17	1	C2019-1752-3-B	-	1.0000	017F1701.D	4
18	18	1	C2019-1752-4-A	-	1.0000	018F1801.D	4
19	19	1	C2019-1752-4-B	-	1.0000	019F1901.D	4
20	20	1	C2019-1801-21A	-	1.0000	020F2001.D	6
21	21	1	C2019-1801-21B	-	1.0000	021F2101.D	6
22	22	1	C2019-1829-1-A	-	1.0000	022F2201.D	2
23	23	1	C2019-1829-1-B	-	1.0000	023F2301.D	2
24	24	1	C2019-1840-1-A	-	1.0000	024F2401.D	4
25	25	1	C2019-1840-1-B	-	1.0000	025F2501.D	4
26	26	1	QC-1(1)-A	-	1.0000	026F2601.D	4
27	27	1	QC-1(1)-B	-	1.0000	027F2701.D	4
28	28	1	C2019-1844-1-A	-	1.0000	028F2801.D	6
29	29	1	C2019-1844-1-B	-	1.0000	029F2901.D	6
30	30	1	C2019-1863-1-A	-	1.0000	030F3001.D	2
31	31	1	C2019-1863-1-B	-	1.0000	031F3101.D	2
32	32	1	C2019-1889-1-A	-	1.0000	032F3201.D	2
33	33	1	C2019-1889-1-B	-	1.0000	033F3301.D	2
34	34	1	C2019-1934-1-A	-	1.0000	034F3401.D	4
35	35	1	C2019-1934-1-B	-	1.0000	035F3501.D	4
36	36	1	C2019-1957-1-A	-	1.0000	036F3601.D	2
37	37	1	C2019-1957-1-B	-	1.0000	037F3701.D	2
38	38	1	C2019-1968-1-A	-	1.0000	038F3801.D	5
39	39	1	C2019-1968-1-B	-	1.0000	039F3901.D	4
40	40	1	C2019-1970-22A	-	1.0000	040F4001.D	2
41	41	1	C2019-1970-22B	-	1.0000	041F4101.D	2
42	42	1	C2019-1973-1-A	-	1.0000	042F4201.D	4
43	43	1	C2019-1973-1-B	-	1.0000	043F4301.D	4
44	44	1	C2019-1976-1-A	-	1.0000	044F4401.D	2
45	45	1	C2019-1976-1-B	-	1.0000	045F4501.D	2
46	46	1	QC-2(2)-A	-	1.0000	046F4601.D	4

99

Run #	Location #	Inj #	Sample Name	Sample Amt [g/100cc]	Multip.* Dilution	File name	Cal #	# Cmp
47	47	1	QC-2(2)-B	-	1.0000	047F4701.D		4
48	48	1	ISTD BLANK-2	-	1.0000	048F4801.D		2
49	49	1	water-2	-	1.0000	049F4901.D		0
50	50	1	0.05 DIAGNOSTIC	-	1.0000	050F5001.D		4
51	51	1	0.100 DIAGNOSTIC	-	1.0000	051F5101.D		4
52	52	1	0.200 DIAGNOSTIC	-	1.0000	052F5201.D		4
53	53	1	0.300 DIAGNOSTIC	-	1.0000	053F5301.D		4
54	54	1	0.500 DIAGNOSTIC	-	1.0000	054F5401.D		4

99

=====  
Calibration Table  
=====

-----  
General Calibration Setting  
-----

Calib. Data Modified : Tuesday, October 15, 2019 3:13:11 PM  
Signals calculated separately : No

Rel. Reference Window : 0.000 %  
Abs. Reference Window : 0.100 min  
Rel. Non-ref. Window : 0.000 %  
Abs. Non-ref. Window : 0.100 min  
Uncalibrated Peaks : not reported  
Partial Calibration : No recalibration if peaks missing

Curve Type : Linear  
Origin : Forced  
Weight : Equal

Recalibration Settings:  
Average Response : Average all calibrations  
Average Retention Time: Floating Average New 75%

Calibration Report Options :  
Printout of recalibrations within a sequence:  
    Calibration Table after Recalibration  
    Normal Report after Recalibration  
If the sequence is done with bracketing:  
    Results of first cycle (ending previous bracket)

Default Sample ISTD Information (if not set in sample table):

ISTD #	ISTD Amount [g/100cc]	Name
1	1.00000	n-Propanol
2	1.00000	n-Propanol

-----  
Signal Details  
-----

Signal 1: FID1 A, Front Signal  
Signal 2: FID2 B, Back Signal  
-----

-----  
Overview Table  
-----


99

RT	Sig	Lvl	Amount [g/100cc]	Area	Rsp.Factor	Ref	ISTD #	Compound
2.000	2	1	1.00000	5.00000	2.00000e-1	No	No 2	Difluoroethane
2.000	1	1	1.00000	5.00000	2.00000e-1	No	No 1	Difluoroethane
2.494	1	1	1.00000	3.69669	2.70512e-1	No	No 1	Methanol
2.772	1	1	1.00000	3.19311	3.13174e-1	No	No 1	Acetaldehyde
2.797	2	1	1.00000	3.10575	3.21983e-1	No	No 2	Acetaldehyde
3.108	1	1	5.00000e-2	8.85689	5.64532e-3	No	No 1	Ethanol
		2	1.00000e-1	17.76771	5.62819e-3			
		3	2.00000e-1	35.45885	5.64034e-3			
		4	3.00000e-1	54.01428	5.55409e-3			
		5	5.00000e-1	90.81020	5.50599e-3			
3.211	2	1	1.00000	4.26062	2.34707e-1	No	No 2	Methanol
3.715	1	1	1.00000	9.73055	1.02769e-1	No	No 1	Isopropyl alcohol
4.180	2	1	5.00000e-2	8.89237	5.62280e-3	No	No 2	Ethanol
		2	1.00000e-1	17.70471	5.64822e-3			
		3	2.00000e-1	35.50703	5.63269e-3			
		4	3.00000e-1	54.17454	5.53766e-3			
		5	5.00000e-1	90.86906	5.50242e-3			
4.530	1	1	1.00000	6.49940	1.53860e-1	No	No 1	Acetone
4.549	2	1	1.00000	6.89301	1.45075e-1	No	No 2	Acetone
4.870	2	1	1.00000	10.70642	9.34019e-2	No	No 2	Isopropyl alcohol
4.942	1	1	1.00000	90.36797	1.10659e-2	No	Yes 1	n-Propanol
		2	1.00000	90.78079	1.10155e-2			
		3	1.00000	90.82334	1.10104e-2			
		4	1.00000	91.87808	1.08840e-2			
		5	1.00000	92.12633	1.08547e-2			
7.619	2	1	1.00000	88.43904	1.13072e-2	No	Yes 2	n-Propanol
		2	1.00000	88.66056	1.12790e-2			
		3	1.00000	88.20879	1.13367e-2			
		4	1.00000	89.09641	1.12238e-2			
		5	1.00000	89.28278	1.12004e-2			

Peak Sum Table


\*\*\*No Entries in table\*\*\*

Calibration Curves




Difluoroethane at exp. RT: 2.000  
 FID2 B, Back Signal  
 Correlation: 1.00000  
 Residual Std. Dev.: 0.00000  
 Formula:  $y = mx$   
 m: 5.65361e-2  
 x: Amount Ratio  
 y: Area Ratio


99




Difluoroethane at exp. RT: 2.000  
 FID1 A, Front Signal  
 Correlation: 1.00000  
 Residual Std. Dev.: 0.00000  
 Formula:  $y = mx$   
 m: 5.53293e-2  
 x: Amount Ratio  
 y: Area Ratio



Methanol at exp. RT: 2.494  
 FID1 A, Front Signal  
 Correlation: 1.00000  
 Residual Std. Dev.: 0.00000  
 Formula:  $y = mx$   
 m: 4.09071e-2  
 x: Amount Ratio  
 y: Area Ratio




Acetaldehyde at exp. RT: 2.772  
 FID1 A, Front Signal  
 Correlation: 1.00000  
 Residual Std. Dev.: 0.00000  
 Formula:  $y = mx$   
 m: 3.53345e-2  
 x: Amount Ratio  
 y: Area Ratio




Acetaldehyde at exp. RT: 2.797  
 FID2 B, Back Signal  
 Correlation: 1.00000  
 Residual Std. Dev.: 0.00000  
 Formula:  $y = mx$   
 m: 3.51174e-2  
 x: Amount Ratio  
 y: Area Ratio


99




Ethanol at exp. RT: 3.108  
 FID1 A, Front Signal  
 Correlation: 0.99999 ✓  
 Residual Std. Dev.: 0.00221  
 Formula:  $y = mx$   
 m: 1.96632  
 x: Amount Ratio  
 y: Area Ratio



Methanol at exp. RT: 3.211  
 FID2 B, Back Signal  
 Correlation: 1.00000  
 Residual Std. Dev.: 0.00000  
 Formula:  $y = mx$   
 m: 4.81758e-2  
 x: Amount Ratio  
 y: Area Ratio




Isopropyl alcohol at exp. RT: 3.715  
 FID1 A, Front Signal  
 Correlation: 1.00000  
 Residual Std. Dev.: 0.00000  
 Formula:  $y = mx$   
 m: 1.07677e-1  
 x: Amount Ratio  
 y: Area Ratio




Ethanol at exp. RT: 4.180  
 FID2 B, Back Signal  
 Correlation: 0.99999 ✓  
 Residual Std. Dev.: 0.00285  
 Formula:  $y = mx$   
 m: 2.03006  
 x: Amount Ratio  
 y: Area Ratio

99







Acetone at exp. RT: 4.530  
 FID1 A, Front Signal  
 Correlation: 1.00000  
 Residual Std. Dev.: 0.00000  
 Formula:  $y = mx$   
 m: 7.19215e-2  
 x: Amount Ratio  
 y: Area Ratio



Acetone at exp. RT: 4.549  
 FID2 B, Back Signal  
 Correlation: 1.00000  
 Residual Std. Dev.: 0.00000  
 Formula:  $y = mx$   
 m: 7.79408e-2  
 x: Amount Ratio  
 y: Area Ratio




Isopropyl alcohol at exp. RT: 4.870  
 FID2 B, Back Signal  
 Correlation: 1.00000  
 Residual Std. Dev.: 0.00000  
 Formula:  $y = mx$   
 m: 1.21060e-1  
 x: Amount Ratio  
 y: Area Ratio



n-Propanol at exp. RT: 4.942  
 FID1 A, Front Signal  
 Correlation: 1.00000  
 Residual Std. Dev.: 0.00000  
 Formula:  $y = mx$   
 m: 1.00000  
 x: Amount Ratio  
 y: Area Ratio

99



n-Propanol at exp. RT: 7.619  
FID2 B, Back Signal  
Correlation: 1.00000  
Residual Std. Dev.: 0.00000  
Formula:  $y = mx$   
m: 1.00000  
x: Amount Ratio  
y: Area Ratio

=====

S a m p l e S u m m a r y


Sequence table: C:\Chem32\1\TEMP\AESEQ\QS\_15.10.2019\_12.36.01\10-15-19cal.S  
 Data directory path: C:\Chem32\1\Data\10-15-19calJJ  
 Logbook: C:\Chem32\1\Data\10-15-19calJJ\10-15-19cal.LOG  
 Sequence start: 10/15/2019 12:49:44 PM  
 Sequence Operator: SYSTEM  
 Operator: SYSTEM

Method file name: C:\CHEM32\1\METHODS\ALCOHOL.M

Run #	Location #	Inj #	Sample Name	Sample Amt [g/100cc]	Multip.* Dilution	File name	Cal #	# Cmp
1	1	1	0.05	-	1.0000	001F0101.D	*	4
2	2	1	0.100	-	1.0000	002F0201.D	*	4
3	3	1	0.200	-	1.0000	003F0301.D	*	4
4	4	1	0.300	-	1.0000	004F0401.D	*	4
5	5	1	0.500	-	1.0000	005F0501.D	*	4
6	6	1	blank	-	1.0000	006F0601.D		2

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.05  
 Laboratory : Coeur d' Alene  
 Injection Date : Oct 15, 2019  
 Method : ALCOHOL.M  
 Acq. Instrument: CN10742044-IT00725005




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	8.85689	0.0498	g/100cc
2.	Ethanol	Column 2:	8.89237	0.0495	g/100cc
3.	n-Propanol	Column 1:	90.36797	1.0000	g/100cc
4.	n-Propanol	Column 2:	88.43904	1.0000	g/100cc

*Handwritten signature*

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.100  
 Laboratory : Coeur d' Alene  
 Injection Date : Oct 15, 2019  
 Method : ALCOHOL.M  
 Acq. Instrument: CN10742044-IT00725005




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	17.76771	0.0995	g/100cc
2.	Ethanol	Column 2:	17.70471	0.0984	g/100cc
3.	n-Propanol	Column 1:	90.78079	1.0000	g/100cc
4.	n-Propanol	Column 2:	88.66056	1.0000	g/100cc

*[Handwritten signature]*

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.200  
 Laboratory : Coeur d' Alene  
 Injection Date : Oct 15, 2019  
 Method : ALCOHOL.M  
 Acq. Instrument: CN10742044-IT00725005




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	35.45885	0.1986	g/100cc
2.	Ethanol	Column 2:	35.50703	0.1983	g/100cc
3.	n-Propanol	Column 1:	90.82334	1.0000	g/100cc
4.	n-Propanol	Column 2:	88.20879	1.0000	g/100cc

99

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.300  
 Laboratory : Coeur d' Alene  
 Injection Date : Oct 15, 2019  
 Method : ALCOHOL.M  
 Acq. Instrument: CN10742044-IT00725005




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	54.01428	0.2990	g/100cc
2.	Ethanol	Column 2:	54.17454	0.2995	g/100cc
3.	n-Propanol	Column 1:	91.87808	1.0000	g/100cc
4.	n-Propanol	Column 2:	89.09641	1.0000	g/100cc

99

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.500  
 Laboratory : Coeur d' Alene  
 Injection Date : Oct 15, 2019  
 Method : ALCOHOL.M  
 Acq. Instrument: CN10742044-IT00725005




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	90.81020	0.5013	g/100cc
2.	Ethanol	Column 2:	90.86906	0.5013	g/100cc
3.	n-Propanol	Column 1:	92.12633	1.0000	g/100cc
4.	n-Propanol	Column 2:	89.28278	1.0000	g/100cc

99



ISP Forensic Services Blood Alcohol Report

Sample Name : blank  
 Laboratory : Coeur d' Alene  
 Injection Date : Oct 15, 2019  
 Method : ALCOHOL.M  
 Acq. Instrument: CN10742044-IT00725005




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	0.00000	0.0000	g/100cc
2.	Ethanol	Column 2:	0.00000	0.0000	g/100cc
3.	n-Propanol	Column 1:	89.98125	1.0000	g/100cc
4.	n-Propanol	Column 2:	87.86541	1.0000	g/100cc

89

ISP Forensic Services Blood Alcohol Report

Sample Name : water-1  
 Laboratory : Coeur d' Alene  
 Injection Date : Oct 15, 2019  
 Method : ALCOHOL.M  
 Acq. Instrument: CN10742044-IT00725005




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	0.00000	0.0000	g/100cc
2.	Ethanol	Column 2:	0.00000	0.0000	g/100cc
3.	n-Propanol	Column 1:	0.00000	0.0000	g/100cc
4.	n-Propanol	Column 2:	0.00000	0.0000	g/100cc

99

ISP Forensic Services Blood Alcohol Report


Sample Name : VOL MIX FN-06041502  
 Laboratory : Coeur d' Alene  
 Injection Date : Oct 15, 2019  
 Method : ALCOHOL.M  
 Acq. Instrument: CN10742044-IT00725005



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	38.70235	0.2210	g/100cc
2.	Ethanol	Column 2:	38.67793	0.2186	g/100cc
3.	n-Propanol	Column 1:	89.04226	1.0000	g/100cc
4.	n-Propanol	Column 2:	87.14028	1.0000	g/100cc

ISP Forensic Services Blood Alcohol Report

Sample Name : ISTD BLANK-1  
 Laboratory : Coeur d' Alene  
 Injection Date : Oct 15, 2019  
 Method : ALCOHOL.M  
 Acq. Instrument: CN10742044-IT00725005



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	0.00000	0.0000	g/100cc
2.	Ethanol	Column 2:	0.00000	0.0000	g/100cc
3.	n-Propanol	Column 1:	90.94904	1.0000	g/100cc
4.	n-Propanol	Column 2:	88.88799	1.0000	g/100cc

99

**VOLATILES DETERMINATION CASEFILE WORKSHEET**

Laboratory No.: QC-2(1)

Analysis Date(s): 15 Oct 2019

	Column 1 FID A	Column 2 FID B	Column Precision	Mean Value	Over-all Mean
Sample Results	0.1961	0.1957	0.0004	0.1959	0.1960
(g/100cc)	0.1966	0.1957	0.0009	0.1961	

**Analysis Method**

Refer to Blood Alcohol Method #1

**Instrument Information**

*Instrument method is stored centrally.*

Refer to Instrument Method: Alcohol.m  
Hamilton Auto-Dilutor Serial Number: ML600HC11379

**Reporting of Results**

**Uncertainty of Measurement (UM%): 5.00%**

Overall Mean (g/100cc)	Low	High	5% of Mean
0.196	0.186	0.206	0.010

Reported Result	
0.196	

*Calibration and control data are stored centrally.*


Revision: 1

Issue Date: 01/04/2019

Issuing Authority: Quality Manager

ISP Forensic Services Blood Alcohol Report

Sample Name : QC-2(1)-A  
 Laboratory : Coeur d' Alene  
 Injection Date : Oct 15, 2019  
 Method : ALCOHOL.M  
 Acq. Instrument: CN10742044-IT00725005




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	35.51124	0.1961	g/100cc
2.	Ethanol	Column 2:	35.56581	0.1957	g/100cc
3.	n-Propanol	Column 1:	92.10008	1.0000	g/100cc
4.	n-Propanol	Column 2:	89.54060	1.0000	g/100cc

99

ISP Forensic Services Blood Alcohol Report

Sample Name : QC-2(1)-B  
 Laboratory : Coeur d' Alene  
 Injection Date : Oct 15, 2019  
 Method : ALCOHOL.M  
 Acq. Instrument: CN10742044-IT00725005



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	35.12582	0.1966	g/100cc
2.	Ethanol	Column 2:	35.09735	0.1957	g/100cc
3.	n-Propanol	Column 1:	90.84682	1.0000	g/100cc
4.	n-Propanol	Column 2:	88.33175	1.0000	g/100cc

99

**VOLATILES DETERMINATION CASEFILE WORKSHEET**

Laboratory No.: 0.08 FN04171701

Analysis Date(s): 15 Oct 2019

	Column 1 FID A	Column 2 FID B	Column Precision	Mean Value	Over-all Mean	
Sample Results	0.0816	0.0810	0.0006	0.0813	0.0806	
(g/100cc)	0.0806	0.0795	0.0011	0.0800		

**Analysis Method**

Refer to Blood Alcohol Method #1

**Instrument Information**

*Instrument method is stored centrally.*

Refer to Instrument Method: Alcohol.m  
Hamilton Auto-Dilutor Serial Number: ML600HC11379

**Reporting of Results**

Uncertainty of Measurement (UM%): 5.00%

Overall Mean (g/100cc)	Low	High	5% of Mean
0.080	0.076	0.084	0.004

	Reported Result	
	0.080	

*Calibration and control data are stored centrally.*

99

Revision: 1


Issue Date: 01/04/2019

Issuing Authority: Quality Manager



ISP Forensic Services Blood Alcohol Report

Sample Name : 0.08 FN04171701-A  
 Laboratory : Coeur d' Alene  
 Injection Date : Oct 15, 2019  
 Method : ALCOHOL.M  
 Acq. Instrument: CN10742044-IT00725005




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	14.59405	0.0816	g/100cc
2.	Ethanol	Column 2:	14.56347	0.0810	g/100cc
3.	n-Propanol	Column 1:	90.93156	1.0000	g/100cc
4.	n-Propanol	Column 2:	88.61271	1.0000	g/100cc

99

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.08 FN04171701-B  
 Laboratory : Coeur d' Alene  
 Injection Date : Oct 15, 2019  
 Method : ALCOHOL.M  
 Acq. Instrument: CN10742044-IT00725005



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	14.30655	0.0806	g/100cc
2.	Ethanol	Column 2:	14.20097	0.0795	g/100cc
3.	n-Propanol	Column 1:	90.25481	1.0000	g/100cc
4.	n-Propanol	Column 2:	88.03085	1.0000	g/100cc

99

**VOLATILES DETERMINATION CASEFILE WORKSHEET**

Laboratory No.: QC-1(1)

Analysis Date(s): 15 Oct 2019

	Column 1 FID A	Column 2 FID B	Column Precision	Mean Value	Over-all Mean	
Sample Results	0.0790	0.0788	0.0002	0.0789	0.0789	
(g/100cc)	0.0791	0.0789	0.0002	0.0790		

**Analysis Method**

Refer to Blood Alcohol Method #1

**Instrument Information**

*Instrument method is stored centrally.*

Refer to Instrument Method: Alcohol.m  
Hamilton Auto-Dilutor Serial Number: ML600HC11379

**Reporting of Results**

Uncertainty of Measurement (UM%): 5.00%

Overall Mean (g/100cc)	Low	High	5% of Mean
0.078	0.074	0.082	0.004

	Reported Result	
	0.078	

*Calibration and control data are stored centrally.*


Revision: 1

Issue Date: 01/04/2019

Issuing Authority: Quality Manager

ISP Forensic Services Blood Alcohol Report

Sample Name : QC-1(1)-A  
 Laboratory : Coeur d' Alene  
 Injection Date : Oct 15, 2019  
 Method : ALCOHOL.M  
 Acq. Instrument: CN10742044-IT00725005




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	14.19905	0.0790	g/100cc
2.	Ethanol	Column 2:	14.27709	0.0788	g/100cc
3.	n-Propanol	Column 1:	91.40876	1.0000	g/100cc
4.	n-Propanol	Column 2:	89.25763	1.0000	g/100cc

99

ISP Forensic Services Blood Alcohol Report

Sample Name : QC-1(1)-B  
 Laboratory : Coeur d' Alene  
 Injection Date : Oct 15, 2019  
 Method : ALCOHOL.M  
 Acq. Instrument: CN10742044-IT00725005



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	14.33776	0.0791	g/100cc
2.	Ethanol	Column 2:	14.41655	0.0789	g/100cc
3.	n-Propanol	Column 1:	92.14104	1.0000	g/100cc
4.	n-Propanol	Column 2:	90.01804	1.0000	g/100cc

99

## VOLATILES DETERMINATION CASEFILE WORKSHEET

Laboratory No.: QC-2(2)

Analysis Date(s): 16 Oct 2019

	Column 1 FID A	Column 2 FID B	Column Precision	Mean Value	Over-all Mean	
Sample Results	0.1989	0.1980	0.0009	0.1984	0.1988	
(g/100cc)	0.1997	0.1988	0.0009	0.1992		

### Analysis Method

Refer to Blood Alcohol Method #1

### Instrument Information

*Instrument method is stored centrally.*

Refer to Instrument Method: Alcohol.m  
Hamilton Auto-Dilutor Serial Number: ML600HC11379

### Reporting of Results

Uncertainty of Measurement (UM%): 5.00%

Overall Mean (g/100cc)	Low	High	5% of Mean
0.198	0.188	0.208	0.010

	Reported Result	
	0.198	

*Calibration and control data are stored centrally.*

Revision: 1


Issue Date: 01/04/2019

Issuing Authority: Quality Manager

99

ISP Forensic Services Blood Alcohol Report

Sample Name : QC-2(2)-A  
 Laboratory : Coeur d' Alene  
 Injection Date : Oct 16, 2019  
 Method : ALCOHOL.M  
 Acq. Instrument: CN10742044-IT00725005




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1;	36.11684	0.1989	g/100cc
2.	Ethanol	Column 2;	36.11486	0.1980	g/100cc
3.	n-Propanol	Column 1;	92.34664	1.0000	g/100cc
4.	n-Propanol	Column 2;	89.84921	1.0000	g/100cc

99

ISP Forensic Services Blood Alcohol Report

Sample Name : QC-2(2)-B  
 Laboratory : Coeur d' Alene  
 Injection Date : Oct 16, 2019  
 Method : ALCOHOL.M  
 Acq. Instrument: CN10742044-IT00725005




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	36.74314	0.1997	g/100cc
2.	Ethanol	Column 2:	36.79988	0.1988	g/100cc
3.	n-Propanol	Column 1:	93.59361	1.0000	g/100cc
4.	n-Propanol	Column 2:	91.20244	1.0000	g/100cc

99



ISP Forensic Services Blood Alcohol Report

Sample Name : ISTD BLANK-2  
 Laboratory : Coeur d' Alene  
 Injection Date : Oct 16, 2019  
 Method : ALCOHOL.M  
 Acq. Instrument: CN10742044-IT00725005




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	0.00000	0.0000	g/100cc
2.	Ethanol	Column 2:	0.00000	0.0000	g/100cc
3.	n-Propanol	Column 1:	91.07959	1.0000	g/100cc
4.	n-Propanol	Column 2:	88.80874	1.0000	g/100cc

99

ISP Forensic Services Blood Alcohol Report

Sample Name : water-2  
 Laboratory : Coeur d' Alene  
 Injection Date : Oct 16, 2019  
 Method : ALCOHOL.M  
 Acq. Instrument: CN10742044-IT00725005




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	0.00000	0.0000	g/100cc
2.	Ethanol	Column 2:	0.00000	0.0000	g/100cc
3.	n-Propanol	Column 1:	0.00000	0.0000	g/100cc
4.	n-Propanol	Column 2:	0.00000	0.0000	g/100cc

99

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.05 DIAGNOSTIC  
 Laboratory : Coeur d' Alene  
 Injection Date : Oct 16, 2019  
 Method : ALCOHOL.M  
 Acq. Instrument: CN10742044-IT00725005




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	9.20663	0.0522	g/100cc
2.	Ethanol	Column 2:	9.09644	0.0511	g/100cc
3.	n-Propanol	Column 1:	89.62101	1.0000	g/100cc
4.	n-Propanol	Column 2:	87.66021	1.0000	g/100cc

99

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.100 DIAGNOSTIC  
 Laboratory : Coeur d' Alene  
 Injection Date : Oct 16, 2019  
 Method : ALCOHOL.M  
 Acq. Instrument : CN10742044-IT00725005




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	18.45908	0.1030	g/100cc
2.	Ethanol	Column 2:	18.48982	0.1028	g/100cc
3.	n-Propanol	Column 1:	91.11005	1.0000	g/100cc
4.	n-Propanol	Column 2:	88.59604	1.0000	g/100cc

99

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.200 DIAGNOSTIC  
 Laboratory : Coeur d' Alene  
 Injection Date : Oct 16, 2019  
 Method : ALCOHOL.M  
 Acq. Instrument: CN10742044-IT00725005




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	37.54866	0.2067	g/100cc
2.	Ethanol	Column 2:	37.43740	0.2056	g/100cc
3.	n-Propanol	Column 1:	92.39802	1.0000	g/100cc
4.	n-Propanol	Column 2:	89.69575	1.0000	g/100cc

*99*

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.300 DIAGNOSTIC  
 Laboratory : Coeur d' Alene  
 Injection Date : Oct 16, 2019  
 Method : ALCOHOL.M  
 Acq. Instrument: CN10742044-IT00725005




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	55.71714	0.3109	g/100cc
2.	Ethanol	Column 2:	55.84768	0.3109	g/100cc
3.	n-Propanol	Column 1:	91.14788	1.0000	g/100cc
4.	n-Propanol	Column 2:	88.48684	1.0000	g/100cc

*[Handwritten signature]*

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.500 DIAGNOSTIC  
 Laboratory : Coeur d' Alene  
 Injection Date : Oct 16, 2019  
 Method : ALCOHOL.M  
 Acq. Instrument: CN10742044-IT00725005



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	94.48638	0.5206	g/100cc
2.	Ethanol	Column 2:	94.58567	0.5206	g/100cc
3.	n-Propanol	Column 1:	92.29834	1.0000	g/100cc
4.	n-Propanol	Column 2:	89.49181	1.0000	g/100cc

99