

**REVIEWED**

By Melissa (Nikka) Bradley at 9:26 am, Jun 03, 2020

MB

Revision: 2  
Issue Date: 12/23/2019  
Issuing Authority: Quality Manager

**Quantitative Analysis for Ethanol & Qualitative Analysis for Other Volatiles**

*Analytical Method(s): 1.0*

**Device: Hamilton MICROLAB Liquid Processor/Dilutor Serial Number: MD96BC1382**

**Volatiles Quality Assurance Controls**

**Run Date(s): 06/01/2020**

Calibration Date: 05/28/2020

Control level	Expiration	Lot #	Target Value	Acceptable Range	Overall Results
Level 1	Jan-22	1801036	0.0812	0.0731-0.0893	0.0814 g/100cc
Level 2	Mar-22	1803028	0.2035	0.1832-0.2238	0.1992 g/100cc 0.2000 g/100cc
Multi-Component mixture:		Lot #	FN06041502		OK
Curve Fit:		Column 1	0.99998	Column2	0.99990




















**Ethanol Calibration Reference Material**

Calibrator level	Target Value	Acceptable Range	Column 1	Column 2	Precision	Mean
50	0.050	0.045 - 0.055	0.0514	0.0531	0.0017	0.0522
100	0.100	0.090 - 0.110	0.0992	0.0993	0.0001	0.0992
200	0.200	0.180 - 0.220	0.1998	0.1981	0.0017	0.1989
300	0.300	0.270 - 0.330	0.2989	0.2973	0.0016	0.2981
400	0.400	0.360 - 0.440				
500	0.500	0.450 - 0.550	0.5008	0.5022	0.0014	0.5015

**Aqueous Controls**

Control level	Target Value	Acceptable Range	Overall Results
80	0.080	0.076 - 0.084	0.082 g/100cc

**Worklist: 4270**

<u>LAB CASE</u>	<u>ITEM</u>	<u>ITEM TYPE</u>	<u>DESCRIPTION</u>	
M2020-1467	3	UCK	Alcohol Analysis	
M2020-1686	2	BCK	Alcohol Analysis	
M2020-1865	1	UCK	Alcohol Analysis	
M2020-1869	1	BCK	Alcohol Analysis	
M2020-1870	1	BCK	Alcohol Analysis	
M2020-1910	1	BCK	Alcohol Analysis	
M2020-1915	1	BCK	Alcohol Analysis	
M2020-1916	1	BCK	Alcohol Analysis	
M2020-1917	1	BCK	Alcohol Analysis	
M2020-1918	1	BCK	Alcohol Analysis	
M2020-1923	1	BCK	Alcohol Analysis	
M2020-1924	1	BCK	Alcohol Analysis	
M2020-1925	1	BCK	Alcohol Analysis	
M2020-1926	1	BCK	Alcohol Analysis	
M2020-1937	1	BCK	Alcohol Analysis	
M2020-1949	1	BCK	Alcohol Analysis	
M2020-1960	1	BCK	Alcohol Analysis	
M2020-1970	1	BCK	Alcohol Analysis	
M2020-1971	1	BCK	Alcohol Analysis	



=====  
Calibration Table  
=====

-----  
General Calibration Setting  
-----

Calib. Data Modified : Thursday, May 28, 2020 10:13:37 AM  
Signals calculated separately : No

Rel. Reference Window : 0.000 %  
Abs. Reference Window : 0.100 min  
Rel. Non-ref. Window : 0.000 %  
Abs. Non-ref. Window : 0.100 min  
Uncalibrated Peaks : not reported  
Partial Calibration : Yes, identified peaks are recalibrated  
Correct All Ret. Times: No, only for identified peaks

Curve Type : Linear  
Origin : Ignored  
Weight : Equal

Recalibration Settings:  
Average Response : Average all calibrations  
Average Retention Time: Floating Average New 75%

Calibration Report Options :  
Printout of recalibrations within a sequence:  
Calibration Table after Recalibration  
Normal Report after Recalibration  
If the sequence is done with bracketing:  
Results of first cycle (ending previous bracket)

Default Sample ISTD Information (if not set in sample table):

ISTD #	ISTD Amount [g/100cc]	Name
1	1.00000	n-propanol
2	1.00000	n-propanol

-----  
Signal Details  
-----

Signal 1: FID1 A, Front Signal  
Signal 2: FID2 B, Back Signal  
-----

-----  
Overview Table  
-----

RT	Sig	Lvl	Amount [g/100cc]	Area	Rsp.Factor	Ref	ISTD #	Compound
2.586	1	1	1.00000	3.69669	2.70512e-1	No	No 1	methanol
2.809	1	1	1.00000	4.26100	2.34687e-1	No	No 2	Acetaldehyde
2.977	2	1	1.00000	4.26100	2.34687e-1	No	No 2	Acetaldehyde
3.075	1	1	5.00000e-2	4.39316	1.13813e-2	No	No 1	ethanol
		2	1.00000e-1	8.62537	1.15937e-2			
		3	2.00000e-1	17.68529	1.13088e-2			
		4	3.00000e-1	26.30473	1.14048e-2			
		5	5.00000e-1	44.86656	1.11442e-2			
3.388	2	1	1.00000	4.26062	2.34707e-1	No	No 2	methanol
3.628	1	1	1.00000	9.73055	1.02769e-1	No	No 1	isopropyl alcohol
4.285	2	1	5.00000e-2	4.43744	1.12678e-2	No	No 2	ethanol
		2	1.00000e-1	8.79797	1.13663e-2			
		3	2.00000e-1	18.28441	1.09383e-2			
		4	3.00000e-1	27.36551	1.09627e-2			
		5	5.00000e-1	47.25839	1.05801e-2			
4.308	1	1	1.00000	6.49940	1.53860e-1	No	No 1	acetone
4.620	1	1	1.00000	43.39005	2.30468e-2	No	Yes 1	n-propanol
		2	1.00000	42.72400	2.34060e-2			
		3	1.00000	42.74154	2.33964e-2			
		4	1.00000	42.25135	2.36679e-2			
		5	1.00000	42.81190	2.33580e-2			
4.661	2	1	1.00000	6.89301	1.45075e-1	No	No 2	acetone
4.969	2	1	1.00000	10.70642	9.34019e-2	No	No 2	isopropyl alcohol
7.550	2	1	1.00000	44.66753	2.23876e-2	No	Yes 2	n-propanol
		2	1.00000	43.66137	2.29035e-2			
		3	1.00000	43.58196	2.29453e-2			
		4	1.00000	42.84920	2.33377e-2			
		5	1.00000	43.30358	2.30928e-2			


Peak Sum Table

\*\*\*No Entries in table\*\*\*

1 Warnings or Errors :


Warning : Curve requires more calibration points., (methanol)

Calibration Curves




methanol at exp. RT: 2.586  
 FID1 A, Front Signal  
 Correlation: 1.00000  
 Residual Std. Dev.: 0.00000  
 Formula:  $y = mx + b$   
 m: 8.51968e-2  
 b: 0.00000  
 x: Amount Ratio  
 y: Area Ratio


*W*




Acetaldehyde at exp. RT: 2.809  
 FID1 A, Front Signal  
 Correlation: 1.00000  
 Residual Std. Dev.: 0.00000  
 Formula:  $y = mx + b$   
 m:  $9.53937e-2$   
 b: 0.00000  
 x: Amount Ratio  
 y: Area Ratio



Acetaldehyde at exp. RT: 2.977  
 FID2 B, Back Signal  
 Correlation: 1.00000  
 Residual Std. Dev.: 0.00000  
 Formula:  $y = mx + b$   
 m:  $9.53937e-2$   
 b: 0.00000  
 x: Amount Ratio  
 y: Area Ratio




ethanol at exp. RT: 3.075  
 FID1 A, Front Signal  
 Correlation: 0.99998  
 Residual Std. Dev.: 0.00262  
 Formula:  $y = mx + b$   
 m: 2.10668  
 b:  $-7.04178e-3$   
 x: Amount Ratio  
 y: Area Ratio




methanol at exp. RT: 3.388  
 FID2 B, Back Signal  
 Correlation: 1.00000  
 Residual Std. Dev.: 0.00000  
 Formula:  $y = mx + b$   
 m:  $9.53853e-2$   
 b: 0.00000  
 x: Amount Ratio  
 y: Area Ratio


*W*




isopropyl alcohol at exp. RT: 3.628  
 FID1 A, Front Signal  
 Correlation: 1.00000  
 Residual Std. Dev.: 0.00000  
 Formula:  $y = mx + b$   
 m: 2.24258e-1  
 b: 0.00000  
 x: Amount Ratio  
 y: Area Ratio



ethanol at exp. RT: 4.285  
 FID2 B, Back Signal  
 Correlation: 0.99990  
 Residual Std. Dev.: 0.00652  
 Formula:  $y = mx + b$   
 m: 2.20859  
 b: -1.79038e-2  
 x: Amount Ratio  
 y: Area Ratio




acetone at exp. RT: 4.308  
 FID1 A, Front Signal  
 Correlation: 1.00000  
 Residual Std. Dev.: 0.00000  
 Formula:  $y = mx + b$   
 m: 1.49790e-1  
 b: 0.00000  
 x: Amount Ratio  
 y: Area Ratio




n-propanol at exp. RT: 4.620  
 FID1 A, Front Signal  
 Correlation: 1.00000  
 Residual Std. Dev.: 0.00000  
 Formula:  $y = mx + b$   
 m: 1.00000  
 b: 0.00000  
 x: Amount Ratio  
 y: Area Ratio


*W*



acetone at exp. RT: 4.661  
FID2 B, Back Signal  
Correlation: 1.00000  
Residual Std. Dev.: 0.00000  
Formula:  $y = mx + b$   
m: 1.54318e-1  
b: 0.00000  
x: Amount Ratio  
y: Area Ratio



isopropyl alcohol at exp. RT: 4.969  
FID2 B, Back Signal  
Correlation: 1.00000  
Residual Std. Dev.: 0.00000  
Formula:  $y = mx + b$   
m: 2.39691e-1  
b: 0.00000  
x: Amount Ratio  
y: Area Ratio




n-propanol at exp. RT: 7.550  
FID2 B, Back Signal  
Correlation: 1.00000  
Residual Std. Dev.: 0.00000  
Formula:  $y = mx + b$   
m: 1.00000  
b: 0.00000  
x: Amount Ratio  
y: Area Ratio

=====

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.050 FN05211804  
 Laboratory : Meridian  
 Injection Date : May 28, 2020  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	4.39316	0.0514	g/100cc
2.	Ethanol	Column 2:	4.43744	0.0531	g/100cc
3.	n-Propanol	Column 1:	43.39005	1.0000	g/100cc
4.	n-Propanol	Column 2:	44.66753	1.0000	g/100cc

*W*



ISP Forensic Services Blood Alcohol Report

Sample Name : 0.100 FN02271802  
 Laboratory : Meridian  
 Injection Date : May 28, 2020  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	8.62537	0.0992	g/100cc
2.	Ethanol	Column 2:	8.79797	0.0993	g/100cc
3.	n-Propanol	Column 1:	42.72400	1.0000	g/100cc
4.	n-Propanol	Column 2:	43.66137	1.0000	g/100cc

*W*

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.200 FN06231704  
 Laboratory : Meridian  
 Injection Date : May 28, 2020  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	17.68529	0.1998	g/100cc
2.	Ethanol	Column 2:	18.28441	0.1981	g/100cc
3.	n-Propanol	Column 1:	42.74154	1.0000	g/100cc
4.	n-Propanol	Column 2:	43.58196	1.0000	g/100cc

*W*

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.300 FN07311804  
 Laboratory : Meridian  
 Injection Date : May 28, 2020  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	26.30473	0.2989	g/100cc
2.	Ethanol	Column 2:	27.36551	0.2973	g/100cc
3.	n-Propanol	Column 1:	42.25135	1.0000	g/100cc
4.	n-Propanol	Column 2:	42.84920	1.0000	g/100cc

*M*

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.500 FN08031602  
 Laboratory : Meridian  
 Injection Date : May 28, 2020  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	44.86656	0.5008	g/100cc
2.	Ethanol	Column 2:	47.25839	0.5022	g/100cc
3.	n-Propanol	Column 1:	42.81190	1.0000	g/100cc
4.	n-Propanol	Column 2:	43.30358	1.0000	g/100cc

*Handwritten signature*

ISP Forensic Services Blood Alcohol Report

Sample Name : INTERNAL STANDARD BLANK  
 Laboratory : Meridian  
 Injection Date : May 28, 2020  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	0.00000	0.0000	g/100cc
2.	Ethanol	Column 2:	0.00000	0.0000	g/100cc
3.	n-Propanol	Column 1:	26.27010	1.0000	g/100cc
4.	n-Propanol	Column 2:	26.09530	1.0000	g/100cc

*W*

S a m p l e S u m m a r y

Sequence table: C:\Chem32\1\Data\05-28-20\_CAL\5-28-20\_CAL 2020-05-28 09-08-26\5-28-20\_CAL  
S  
 Data directory path: C:\Chem32\1\Data\05-28-20\_CAL\5-28-20\_CAL 2020-05-28 09-08-26\  
 Logbook: C:\Chem32\1\Data\05-28-20\_CAL\5-28-20\_CAL 2020-05-28 09-08-26\5-28-20\_CAL  
 LOG  
 Sequence start: 5/28/2020 9:23:04 AM  
 Sequence Operator: SYSTEM  
 Operator: SYSTEM


Method file name: C:\Chem32\1\Data\05-28-20\_CAL\5-28-20\_CAL 2020-05-28 09-08-26\ALCOHOL.M

Run #	Location #	Inj #	Sample Name	Sample Amt [g/100cc]	Multip.* Dilution	File name	Cal #	# Cmp
1	1	1	0.050 FN05211804	-	1.0000	001F0101.D	*	4
2	2	1	0.100 FN02271802	-	1.0000	002F0201.D	*	4
3	3	1	0.200 FN06231704	-	1.0000	003F0301.D	*	4
4	4	1	0.300 FN07311804	-	1.0000	004F0401.D	*	4
5	5	1	0.500 FN08031602	-	1.0000	005F0501.D	*	4
6	6	1	INTERNAL STANDAR	-	1.0000	006F0601.D		2

*W*

ISP Forensic Services Blood Alcohol Report

Sample Name : INTERNAL STD BLK 1  
 Laboratory : Meridian  
 Injection Date : Jun 1, 2020  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	0.00000	0.0000	g/100cc
2.	Ethanol	Column 2:	0.00000	0.0000	g/100cc
3.	n-Propanol	Column 1:	28.13090	1.0000	g/100cc
4.	n-Propanol	Column 2:	28.68274	1.0000	g/100cc

*W*

ISP Forensic Services Blood Alcohol Report

Sample Name : MIX VOL FN06041502  
 Laboratory : Meridian  
 Injection Date : Jun 1, 2020  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	6.81683	0.1378	g/100cc
2.	Ethanol	Column 2:	6.98596	0.1387	g/100cc
3.	n-Propanol	Column 1:	24.05987	1.0000	g/100cc
4.	n-Propanol	Column 2:	24.22479	1.0000	g/100cc

*W*



**VOLATILES DETERMINATION CASEFILE WORKSHEET**

Laboratory No.: QC1-1

Analysis Date(s): 01 Jun 2020

	Column 1 FID A	Column 2 FID B	Column Precision	Mean Value	Sample A-B Difference	Over-all Mean
Sample Results	0.0805	0.0815	0.0010	0.0810	0.0008	0.0814
(g/100cc)	0.0812	0.0825	0.0013	0.0818		

**Analysis Method**

Refer to Blood Alcohol Method #1

**Instrument Information**

*Instrument information is stored centrally.*

Refer to Instrument Method: Alcohol.m

**Reporting of Results**

Uncertainty of Measurement (UM%): 5.00%

Overall Mean (g/100cc)	Low	High	5% of Mean
0.081	0.076	0.086	0.005

Reported Result	
0.081	

*Calibration and control data are stored centrally.*


Revision: 2

Issue Date: 12/23/2019

Issuing Authority: Quality Manager

ISP Forensic Services Blood Alcohol Report

Sample Name : QC1-1-A  
 Laboratory : Meridian  
 Injection Date : Jun 1, 2020  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	6.90812	0.0805	g/100cc
2.	Ethanol	Column 2:	7.04403	0.0815	g/100cc
3.	n-Propanol	Column 1:	42.48580	1.0000	g/100cc
4.	n-Propanol	Column 2:	43.44244	1.0000	g/100cc

W

ISP Forensic Services Blood Alcohol Report

Sample Name : QC1-1-B  
 Laboratory : Meridian  
 Injection Date : Jun 1, 2020  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	6.90945	0.0812	g/100cc
2.	Ethanol	Column 2:	7.08906	0.0825	g/100cc
3.	n-Propanol	Column 1:	42.10252	1.0000	g/100cc
4.	n-Propanol	Column 2:	43.15152	1.0000	g/100cc

W

**VOLATILES DETERMINATION CASEFILE WORKSHEET**

Laboratory No.: 0.08 FN04171701

Analysis Date(s): 01 Jun 2020

	Column 1 FID A	Column 2 FID B	Column Precision	Mean Value	Sample A-B Difference	Over-all Mean
Sample Results	0.0822	0.0830	0.0008	0.0826	0.0002	0.0825
(g/100cc)	0.0819	0.0829	0.0010	0.0824		

**Analysis Method**

Refer to Blood Alcohol Method #1

**Instrument Information**

*Instrument information is stored centrally.*

Refer to Instrument Method: Alcohol.m

**Reporting of Results**

Uncertainty of Measurement (UM%): 5.00%

Overall Mean (g/100cc)	Low	High	5% of Mean
0.082	0.077	0.087	0.005

Reported Result	
0.082	

*Calibration and control data are stored centrally.*


Revision: 2

Issue Date: 12/23/2019

Issuing Authority: Quality Manager

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.08 FN04171701-A  
 Laboratory : Meridian  
 Injection Date : Jun 1, 2020  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	7.00777	0.0822	g/100cc
2.	Ethanol	Column 2:	7.08737	0.0830	g/100cc
3.	n-Propanol	Column 1:	42.16664	1.0000	g/100cc
4.	n-Propanol	Column 2:	42.87310	1.0000	g/100cc

W

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.08 FN04171701-B  
 Laboratory : Meridian  
 Injection Date : Jun 1, 2020  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	6.98846	0.0819	g/100cc
2.	Ethanol	Column 2:	7.09599	0.0829	g/100cc
3.	n-Propanol	Column 1:	42.24070	1.0000	g/100cc
4.	n-Propanol	Column 2:	42.95771	1.0000	g/100cc

W

## VOLATILES DETERMINATION CASEFILE WORKSHEET

Laboratory No.: QC2-1

Analysis Date(s): 01 Jun 2020

	Column 1 FID A	Column 2 FID B	Column Precision	Mean Value	Sample A-B Difference	Over-all Mean
Sample Results	0.1974	0.1965	0.0009	0.1969	0.0047	0.1992
(g/100cc)	0.2016	0.2016	0.0000	0.2016		

### Analysis Method

Refer to Blood Alcohol Method #1

### Instrument Information

*Instrument information is stored centrally.*

Refer to Instrument Method: Alcohol.m

### Reporting of Results

Uncertainty of Measurement (UM%): 5.00%

Overall Mean (g/100cc)	Low	High	5% of Mean
0.199	0.189	0.209	0.010

	Reported Result	
	0.199	

*Calibration and control data are stored centrally.*

Revision: 2


Issue Date: 12/23/2019

Issuing Authority: Quality Manager



ISP Forensic Services Blood Alcohol Report

Sample Name : QC2-1-A  
 Laboratory : Meridian  
 Injection Date : Jun 1, 2020  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	17.39332	0.1974	g/100cc
2.	Ethanol	Column 2:	18.02763	0.1965	g/100cc
3.	n-Propanol	Column 1:	42.54603	1.0000	g/100cc
4.	n-Propanol	Column 2:	43.31979	1.0000	g/100cc

*W*



ISP Forensic Services Blood Alcohol Report

Sample Name : QC2-1-B  
 Laboratory : Meridian  
 Injection Date : Jun 1, 2020  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	17.63523	0.2016	g/100cc
2.	Ethanol	Column 2:	18.28837	0.2016	g/100cc
3.	n-Propanol	Column 1:	42.23222	1.0000	g/100cc
4.	n-Propanol	Column 2:	42.78906	1.0000	g/100cc

*W*

## VOLATILES DETERMINATION CASEFILE WORKSHEET

Laboratory No.: ~~QCI-2~~ <sup>QC 2-2 06/02/2020 BT</sup>

Analysis Date(s): 01 Jun 2020

	Column 1 FID A	Column 2 FID B	Column Precision	Mean Value	Sample A-B Difference	Over-all Mean
Sample Results	0.2022	0.2014	0.0008	0.2018	0.0035	0.2000
(g/100cc)	0.1986	0.1981	0.0005	0.1983		

### Analysis Method

Refer to Blood Alcohol Method #1

### Instrument Information

*Instrument information is stored centrally.*

Refer to Instrument Method: Alcohol.m

### Reporting of Results

Uncertainty of Measurement (UM%): 5.00%

Overall Mean (g/100cc)	Low	High	5% of Mean
0.200	0.190	0.210	0.010

	Reported Result
	0.200

*Calibration and control data are stored centrally.*



Revision: 2


Issue Date: 12/23/2019

Issuing Authority: Quality Manager

ISP Forensic Services Blood Alcohol Report

QC2-2A 06/02/2020 BS

Sample Name : ~~QC1-2-A~~  
 Laboratory : Meridian  
 Injection Date : Jun 1, 2020  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	18.21429	0.2022	g/100cc
2.	Ethanol	Column 2:	18.86411	0.2014	g/100cc
3.	n-Propanol	Column 1:	43.47986	1.0000	g/100cc
4.	n-Propanol	Column 2:	44.18612	1.0000	g/100cc

W

ISP Forensic Services Blood Alcohol Report

QC2-2B 06/02/2020 60

Sample Name : ~~QC1-2-B~~  
 Laboratory : Meridian  
 Injection Date : Jun 1, 2020  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	17.51109	0.1986	g/100cc
2.	Ethanol	Column 2:	18.13487	0.1981	g/100cc
3.	n-Propanol	Column 1:	42.56499	1.0000	g/100cc
4.	n-Propanol	Column 2:	43.21969	1.0000	g/100cc

*W*

ISP Forensic Services Blood Alcohol Report

Sample Name : INTERNAL STD BLK  
 Laboratory : Meridian  
 Injection Date : Jun 1, 2020  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	0.00000	0.0000	g/100cc
2.	Ethanol	Column 2:	0.00000	0.0000	g/100cc
3.	n-Propanol	Column 1:	38.86176	1.0000	g/100cc
4.	n-Propanol	Column 2:	39.39694	1.0000	g/100cc

*W*

S a m p l e S u m m a r y

Sequence table: C:\Chem32\1\Data\06-01-20\_SAMPLES\6-1-20\_SAMPLES 2020-06-01 10-20-39\6-1-20\_SAMPLES.S  
 Data directory path: C:\Chem32\1\Data\06-01-20\_SAMPLES\6-1-20\_SAMPLES 2020-06-01 10-20-39\  
 Logbook: C:\Chem32\1\Data\06-01-20\_SAMPLES\6-1-20\_SAMPLES 2020-06-01 10-20-39\6-1-20\_SAMPLES.LOG  
 Sequence start: 6/1/2020 10:35:23 AM  
 Sequence Operator: SYSTEM  
 Operator: SYSTEM

Method file name: C:\Chem32\1\Data\06-01-20\_SAMPLES\6-1-20\_SAMPLES 2020-06-01 10-20-39  
 \ALCOHOL.M

Run #	Location #	Inj #	Sample Name	Sample Amt [g/100cc]	Multip. * Dilution	File name	Cal #	# Cmp
1	1	1	INTERNAL STD BLK	-	1.0000	001F0101.D		2
2	2	1	MIX VOL FN060415	-	1.0000	002F0201.D		10
3	3	1	QC1-1-A	-	1.0000	003F0301.D		4
4	4	1	QC1-1-B	-	1.0000	004F0401.D		4
5	5	1	0.08 FN04171701-	-	1.0000	005F0501.D		4
6	6	1	0.08 FN04171701-	-	1.0000	006F0601.D		4
7	7	1	M2020-1467-3-A	-	1.0000	007F0701.D		2
8	8	1	M2020-1467-3-B	-	1.0000	008F0801.D		2
9	9	1	M2020-1686-2-A	-	1.0000	009F0901.D		4
10	10	1	M2020-1686-2-B	-	1.0000	010F1001.D		4
11	11	1	M2020-1865-1-A	-	1.0000	011F1101.D		4
12	12	1	M2020-1865-1-B	-	1.0000	012F1201.D		4
13	13	1	M2020-1869-1-A	-	1.0000	013F1301.D		4
14	14	1	M2020-1869-1-B	-	1.0000	014F1401.D		4
15	15	1	M2020-1870-1-A	-	1.0000	015F1501.D		4
16	16	1	M2020-1870-1-B	-	1.0000	016F1601.D		4
17	17	1	M2020-1910-1-A	-	1.0000	017F1701.D		4
18	18	1	M2020-1910-1-B	-	1.0000	018F1801.D		4
19	19	1	M2020-1915-1-A	-	1.0000	019F1901.D		4
20	20	1	M2020-1915-1-B	-	1.0000	020F2001.D		4
21	21	1	M2020-1916-1-A	-	1.0000	021F2101.D		2
22	22	1	M2020-1916-1-B	-	1.0000	022F2201.D		2
23	23	1	M2020-1917-1-A	-	1.0000	023F2301.D		2
24	24	1	M2020-1917-1-B	-	1.0000	024F2401.D		2
25	25	1	QC2-1-A	-	1.0000	025F2501.D		4
26	26	1	QC2-1-B	-	1.0000	026F2601.D		4
27	27	1	M2020-1918-1-A	-	1.0000	027F2701.D		4
28	28	1	M2020-1918-1-B	-	1.0000	028F2801.D		4
29	29	1	M2020-1923-1-A	-	1.0000	029F2901.D		4
30	30	1	M2020-1923-1-B	-	1.0000	030F3001.D		4
31	31	1	M2020-1924-1-A	-	1.0000	031F3101.D		4
32	32	1	M2020-1924-1-B	-	1.0000	032F3201.D		4
33	33	1	M2020-1925-1-A	-	1.0000	033F3301.D		4
34	34	1	M2020-1925-1-B	-	1.0000	034F3401.D		4
35	35	1	M2020-1926-1-A	-	1.0000	035F3501.D		4
36	36	1	M2020-1926-1-B	-	1.0000	036F3601.D		4
37	37	1	M2020-1937-1-A	-	1.0000	037F3701.D		4
38	38	1	M2020-1937-1-B	-	1.0000	038F3801.D		4
39	39	1	M2020-1949-1-A	-	1.0000	039F3901.D		3
40	40	1	M2020-1949-1-B	-	1.0000	040F4001.D		3
41	41	1	M2020-1960-1-A	-	1.0000	041F4101.D		2
42	42	1	M2020-1960-1-B	-	1.0000	042F4201.D		2
43	43	1	M2020-1970-1-A	-	1.0000	043F4301.D		4

Run #	Location #	Inj #	Sample Name	Sample Amt [g/100cc]	Multip.* Dilution	File name	Cal #	# Cmp
44	44	1	M2020-1970-1-B	-	1.0000	044F4401.D		4
45	45	1	M2020-1971-1-A	-	1.0000	045F4501.D		4
46	46	1	M2020-1971-1-B	-	1.0000	046F4601.D		4
47	47	1	<del>QC1-2-A</del> QC2-2A 06/01/2020 SV	-	1.0000	047F4701.D		4
48	48	1	<del>QC1-2-B</del> QC2-2B 06/01/2020 SV	-	1.0000	048F4801.D		4
49	49	1	INTERNAL STD BLK	-	1.0000	049F4901.D		2

Method file name: C:\Chem32\1\Data\06-01-20\_SAMPLES\6-1-20\_SAMPLES 2020-06-01 10-20-39 \SHUTDOWN.M

Run #	Location #	Inj #	Sample Name	Sample Amt [g/100cc]	Multip.* Dilution	File name	Cal #	# Cmp
50	50	1	EMPTY	-	1.0000	050F5001.D		0