

B. Wylie

[Type here]

Quantitation of THC and Metabolites in Blood by LC-MS/MS

Date: 1-3-17 Analyst: Anne Nord

PRE-ANALYTIC

0490364

Std. Lot# 70103 Std. Exp. 9-21-17 QC Lot# 31317 QC Exp. 3-13-17
 Plate A

- 1. Mobile Phase A 0.1% Formic Acid in water
 Mobile Phase B 0.1% Formic Acid in water
 Methanol Hexane MTBE
- 2. Check levels of mobile phases and needle wash and refill as necessary. Ensure waste is not full.
- 3. Begin mobile phase flow and allow system to equilibrate for approx. 1 hour.
- 4. Create worklist. Data path name: 1-3-17 Cann quant

ANALYTIC

- 1. Remove standards plate, blood, and samples from cold storage. Allow to reach room temperature.
- 2. Add 1000 µL blood to wells of analytical (standards) plate. Mix via aspirate and dispense.
 Blank blood for locations containing standards/QCs and internal standards (Blank Blood Used 321632-1)
 Sample blood for locations containing only internal standards
- 3. Place on shaking incubator at ambient temp., 900rpm for 15 minutes. Shaker ID 66759
- 4. Pipette 500µL 0.1% formic acid to all wells of standards plate.
- 5. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 6. Transfer 800µL of blood+acid mixture to corresponding wells of SLE+ plate.
- 7. Apply positive pressure for approx. 4 seconds (or until no liquid remains on top of sorbent). Wait 5 min.
- 8. Add 2.25mL MTBE and allow to flow under gravity for 5 minutes.
- 9. Apply positive pressure for approx. 15 seconds.
- 10. Add 2.25mL hexane and allow to flow under gravity for 5 minutes.
- 11. Apply positive pressure for approx. 15 seconds.
- 12. Remove plate containing eluate. Place on SPE Dry and evaporate to dryness at approx. 35°C.
- 13. SPE Dry ID 66819
- 14. Reconstitute in 100 µL MeOH and seal plate with foil. Place in autosampler and run worklist.

POST-ANALYTIC

- 1. Open quantitation software and create a new quantitation batch. Batch name: 1317 Cann quant
- 2. Make any necessary integration changes. Limit curves based on validated linear ranges.
- 3. Were all appropriate standards used in the curve for each analyte? Y/N
 Are r² values ≥ 0.99 for each analyte? Y/N 0.98 per method A
- 4. For unknown samples, are the response ratios for each analyte within ± 30% of the average ratios of the standards? ± 20%
- 5. Did all QCs pass for each analyte? Were QCs entered into QC charting? (Y) N
- 6. Print the following items: Method Checklist Worklist Data Report

COMMENTS

Worklist: 1456

<u>LAB CASE</u>	<u>ITEM</u>	<u>TASK ID</u>	<u>DESCRIPTION</u>	
C2016-1405	1	73518	AM 27 Blood THC Quant by LC	
C2016-2151	1	73536	AM 27 Blood THC Quant by LC	
M2016-3278	2	73022	AM 27 Blood THC Quant by LC	
M2016-3640	1	73523	AM 27 Blood THC Quant by LC	
M2016-3677	1	73519	AM 27 Blood THC Quant by LC	
M2016-3677	2	73520	AM 27 Blood THC Quant by LC	
M2016-3792	1	73521	AM 27 Blood THC Quant by LC	
M2016-3948	1	73522	AM 27 Blood THC Quant by LC	
M2016-4114	1	73524	AM 27 Blood THC Quant by LC	
M2016-4158	1	73534	AM 27 Blood THC Quant by LC	
M2016-4477	1	73526	AM 27 Blood THC Quant by LC	
M2016-4536	1	73527	AM 27 Blood THC Quant by LC	
M2016-4557	2	73535	AM 27 Blood THC Quant by LC	
M2016-4741	1	73528	AM 27 Blood THC Quant by LC	
M2016-4866	1	73529	AM 27 Blood THC Quant by LC	
M2016-4969	1	73530	AM 27 Blood THC Quant by LC	
M2016-5099	2	73525	AM 27 Blood THC Quant by LC	
P2016-1858	1	73531	AM 27 Blood THC Quant by LC	
P2016-2556	1	73532	AM 27 Blood THC Quant by LC	
P2016-2624	1	73533	AM 27 Blood THC Quant by LC	



ISP FORENSICS - Cd'A


Cannabinoids Analysis Report

Batch Data Path D:\2016 Data\1-3-17 cann quant\QuantResults\1317 cann quant.batch.bin
Analysis Time 1/4/2017 10:17 AM **Analyst Name** ISP Tox
Report Time 1/4/2017 10:46 AM **Reporter Name** ISP Tox
Last Calib Update 1/4/2017 10:17 AM **Batch State** Processed

Analysis Info

Acq Time 2017-01-03 14:45 **Data File** External control lot 31317.d
Sample Type Sample **Sample Name** External control lot 31317
Dilution 1 **Acq Method** Quant THC 12-15-16.m
Position p1b2 **Sample Info**
Inj Vol -1 **Comment** AM 27 cannabinoid confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.476	28680	263853	0.1087	11.6601
THC-COOH	THC-COOH-d9	2.626	16367	96058	0.1704	11.4684
THC	THC-d3	6.913	8411	85315	0.0986	9.6196

ISP FORENSICS - Cd'A


Cannabinoids Analysis Report

Batch Data Path D:\2016 Data\1-3-17 cann quant\QuantResults\1317 cann quant.batch.bin
Analysis Time 1/4/2017 10:17 AM **Analyst Name** ISP Tox
Report Time 1/4/2017 10:46 AM **Reporter Name** ISP Tox
Last Calib Update 1/4/2017 10:17 AM **Batch State** Processed

Analysis Info

Acq Time 2017-01-03 14:33 **Data File** QC 10.d
Sample Type QC **Sample Name** QC 10
Dilution 1 **Acq Method** Quant THC 12-15-16.m
Position P1-H1 **Sample Info**
Inj Vol -1 **Comment** AM 27 cannabinoid confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.476	31590	398667	0.0792	8.9248
THC-COOH	THC-COOH-d9	2.626	20119	126550	0.1590	10.6331
THC	THC-d3	6.913	11542	128091	0.0901	8.9006

ISP FORENSICS - Cd'A


Cannabinoids Analysis Report

Batch Data Path D:\2016 Data\1-3-17 cann quant\QuantResults\1317 cann quant.batch.bin
Analysis Time 1/4/2017 10:17 AM **Analyst Name** ISP Tox
Report Time 1/4/2017 10:46 AM **Reporter Name** ISP Tox
Last Calib Update 1/4/2017 10:17 AM **Batch State** Processed

Analysis Info

Acq Time 2017-01-03 14:21 **Data File** Neg Control.d
Sample Type Sample **Sample Name** Neg Control
Dilution 1 **Acq Method** Quant THC 12-15-16.m
Position P1-a2 **Sample Info**
Inj Vol -1 **Comment** AM 27 cannabinoid confirmation

Sample Chromatogram



ISP Forensics Calibration Curve Report


Batch Data Path D:\2016 Data\1-3-17 cann quant\QuantResults\1317 cann quant.batch.bin

Last Calib Update 1/4/2017 10:17 AM

Analyst Name ISP TOX

Target Compound *THC-OH*

Internal Standard *THC-OH-d3*




Sample	Level	Enabled	Exp Conc	Final Conc	Accuracy
Cal 1	1	<input checked="" type="checkbox"/>	3	3.4	112.1
cal 2	2	<input checked="" type="checkbox"/>	5	5.3	105.4
cal 3	3	<input checked="" type="checkbox"/>	10	8.9	89.2
QC 10	3	<input checked="" type="checkbox"/>	10	8.9	89.2
cal 4	4	<input checked="" type="checkbox"/>	25	21.2	84.8
cal 5	5	<input checked="" type="checkbox"/>	50	54.2	108.5
cal 6	6	<input type="checkbox"/>	100	132.6	132.6
cal 7	7	<input type="checkbox"/>	250	290.7	116.3

ISP Forensics Calibration Curve Report

Batch Data Path D:\2016 Data\1-3-17 cann quant\QuantResults\1317 cann quant.batch.bin

Last Calib Update 1/4/2017 10:17 AM **Analyst Name** ISP TOX

Target Compound *THC-COOH*
Internal Standard *THC-COOH-d9*



Sample	Level	Enabled	Exp Conc	Final Conc	Accuracy
Cal 1	1	<input checked="" type="checkbox"/>	3	2.9	96.4
cal 2	2	<input checked="" type="checkbox"/>	5	5.0	100.6
cal 3	3	<input checked="" type="checkbox"/>	10	10.4	103.7
QC 10	3	<input checked="" type="checkbox"/>	10	10.6	106.3
cal 4	4	<input checked="" type="checkbox"/>	25	24.7	98.8
cal 5	5	<input checked="" type="checkbox"/>	50	50.1	100.1
cal 6	6	<input type="checkbox"/>	100	88.0	88.0
cal 7	7	<input type="checkbox"/>	250	137.2	54.9


ISP Forensics Calibration Curve Report

Batch Data Path D:\2016 Data\1-3-17 cann quant\QuantResults\1317 cann quant.batch.bin

Last Calib Update 1/4/2017 10:17 AM

Analyst Name ISP TOX

Target Compound *THC*
Internal Standard *THC-d3*



Sample	Level	Enabled	Exp Conc	Final Conc	Accuracy
Cal 1	1	<input checked="" type="checkbox"/>	3	3.5	117.2
cal 2	2	<input checked="" type="checkbox"/>	5	5.1	101.2
cal 3	3	<input checked="" type="checkbox"/>	10	8.8	87.5
QC 10	3	<input checked="" type="checkbox"/>	10	8.9	89.0
cal 4	4	<input checked="" type="checkbox"/>	25	21.4	85.5
cal 5	5	<input checked="" type="checkbox"/>	50	54.3	108.6
cal 6	6	<input type="checkbox"/>	100	135.6	135.6
cal 7	7	<input type="checkbox"/>	250	299.7	119.9

ISP FORENSICS - Cd'A


Cannabinoids Analysis Report

Batch Data Path D:\2016 Data\1-3-17 cann quant\QuantResults\1317 cann quant.batch.bin
Analysis Time 1/4/2017 10:17 AM **Analyst Name** ISP Tox
Report Time 1/4/2017 10:46 AM **Reporter Name** ISP Tox
Last Calib Update 1/4/2017 10:17 AM **Batch State** Processed

Analysis Info

Acq Time 2017-01-03 12:47 **Data File** Cal 1.d
Sample Type Calibration **Sample Name** Cal 1
Dilution 1 **Acq Method** Quant THC 12-15-16.m
Position P1-A1 **Sample Info**
Inj Vol -1 **Comment** AM 27 cannabinoid confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.476	5986	309560	0.0193	3.3629
THC-COOH	THC-COOH-d9	2.606	6715	94979	0.0707	2.8916
THC	THC-d3	6.893	2546	95650	0.0266	3.5147

ISP FORENSICS - Cd'A


Cannabinoids Analysis Report

Batch Data Path D:\2016 Data\1-3-17 cann quant\QuantResults\1317 cann quant.batch.bin
Analysis Time 1/4/2017 10:17 AM **Analyst Name** ISP Tox
Report Time 1/4/2017 10:46 AM **Reporter Name** ISP Tox
Last Calib Update 1/4/2017 10:17 AM **Batch State** Processed

Analysis Info

Acq Time 2017-01-03 12:58 **Data File** cal 2.d
Sample Type Calibration **Sample Name** cal 2
Dilution 1 **Acq Method** Quant THC 12-15-16.m
Position P1-B1 **Sample Info**
Inj Vol -1 **Comment** AM 27 cannabinoid confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.496	13503	338737	0.0399	5.2688
THC-COOH	THC-COOH-d9	2.646	9454	102827	0.0919	5.0304
THC	THC-d3	6.953	4931	109921	0.0449	5.0619

ISP FORENSICS - Cd'A


Cannabinoids Analysis Report

Batch Data Path D:\2016 Data\1-3-17 cann quant\QuantResults\1317 cann quant.batch.bin
Analysis Time 1/4/2017 10:17 AM **Analyst Name** ISP Tox
Report Time 1/4/2017 10:46 AM **Reporter Name** ISP Tox
Last Calib Update 1/4/2017 10:17 AM **Batch State** Processed

Analysis Info

Acq Time 2017-01-03 13:10 **Data File** cal 3.d
Sample Type Calibration **Sample Name** cal 3
Dilution 1 **Acq Method** Quant THC 12-15-16.m
Position P1-C1 **Sample Info**
Inj Vol -1 **Comment** AM 27 cannabinoid confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.476	28551	360385	0.0792	8.9235
THC-COOH	THC-COOH-d9	2.626	17220	110758	0.1555	10.3708
THC	THC-d3	6.933	9934	112390	0.0884	8.7543

ISP FORENSICS - Cd'A


Cannabinoids Analysis Report

Batch Data Path D:\2016 Data\1-3-17 cann quant\QuantResults\1317 cann quant.batch.bin
Analysis Time 1/4/2017 10:17 AM **Analyst Name** ISP Tox
Report Time 1/4/2017 10:46 AM **Reporter Name** ISP Tox
Last Calib Update 1/4/2017 10:17 AM **Batch State** Processed

Analysis Info

Acq Time 2017-01-03 13:22 **Data File** cal 4.d
Sample Type Calibration **Sample Name** cal 4
Dilution 1 **Acq Method** Quant THC 12-15-16.m
Position P1-D1 **Sample Info**
Inj Vol -1 **Comment** AM 27 cannabinoid confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.456	79043	373794	0.2115	21.2015
THC-COOH	THC-COOH-d9	2.626	46360	115954	0.3998	24.6895
THC	THC-d3	6.913	26620	112304	0.2370	21.3628

ISP FORENSICS - Cd'A


Cannabinoids Analysis Report

Batch Data Path D:\2016 Data\1-3-17 cann quant\QuantResults\1317 cann quant.batch.bin
Analysis Time 1/4/2017 10:17 AM **Analyst Name** ISP Tox
Report Time 1/4/2017 10:46 AM **Reporter Name** ISP Tox
Last Calib Update 1/4/2017 10:17 AM **Batch State** Processed

Analysis Info

Acq Time 2017-01-03 13:34 **Data File** cal 5.d
Sample Type Calibration **Sample Name** cal 5
Dilution 1 **Acq Method** Quant THC 12-15-16.m
Position P1-E1 **Sample Info**
Inj Vol -1 **Comment** AM 27 cannabinoid confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.476	218453	385057	0.5673	54.2433
THC-COOH	THC-COOH-d9	2.626	130617	118996	1.0977	50.0672
THC	THC-d3	6.913	71747	114719	0.6254	54.3063

ISP FORENSICS - Cd'A


Cannabinoids Analysis Report

Batch Data Path D:\2016 Data\1-3-17 cann quant\QuantResults\1317 cann quant.batch.bin
Analysis Time 1/4/2017 10:17 AM **Analyst Name** ISP Tox
Report Time 1/4/2017 10:46 AM **Reporter Name** ISP Tox
Last Calib Update 1/4/2017 10:17 AM **Batch State** Processed

Analysis Info

Acq Time 2017-01-03 13:46 **Data File** cal 6.d
Sample Type Calibration **Sample Name** cal 6
Dilution 1 **Acq Method** Quant THC 12-15-16.m
Position P1-F1 **Sample Info**
Inj Vol -1 **Comment** AM 27 cannabinoid confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.476	588718	417236	1.4110	132.5772
THC-COOH	THC-COOH-d9	2.626	341712	123360	2.7700	87.9635
THC	THC-d3	6.933	192490	121579	1.5832	135.5525

ISP FORENSICS - Cd'A


Cannabinoids Analysis Report

Batch Data Path D:\2016 Data\1-3-17 cann quant\QuantResults\1317 cann quant.batch.bin
Analysis Time 1/4/2017 10:17 AM **Analyst Name** ISP Tox
Report Time 1/4/2017 10:46 AM **Reporter Name** ISP Tox
Last Calib Update 1/4/2017 10:17 AM **Batch State** Processed

Analysis Info

Acq Time 2017-01-03 13:58 **Data File** cal 7.d
Sample Type Calibration **Sample Name** cal 7
Dilution 1 **Acq Method** Quant THC 12-15-16.m
Position P1-G1 **Sample Info**
Inj Vol -1 **Comment** AM 27 cannabinoid confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.476	1245714	399995	3.1143	290.7304
THC-COOH	THC-COOH-d9	2.626	702765	115728	6.0726	137.2257
THC	THC-d3	6.913	396071	112559	3.5188	299.7294