


ISP FORENSIC BIOLOGY QUALITY MANUAL




Revision 10
11/29/10
Issuing Authority: Quality Manager

Forensic Biology Quality Manual

Revision #10



APPROVED
November 29, 2010



Cynthia Cunnington
Forensic Scientist IV
Forensic Biology Supervisor/Technical Leader

Property of Idaho State Police Forensic Services
Uncontrolled Interim Copy
OBSOLETE DOCUMENT

Idaho State Police Forensic Services

Approval for Quality System Controlled Documents




Discipline/Name of Document: Biology Quality Manual

Revision Number: 10

Issue Date: 11/29/2010

APPROVED BY:


Quality Manager

11/29/2010
Date Signed

Checklist Submitted and Checked

WJG

Forensic Biology Quality Manual
Table of Contents

- i. Signature/Approval page
 - ii. Revision History
 - iii. Table of Contents
 1. Introduction: Statement of Purpose and Objectives
 2. Organization and Management
 3. Personnel Qualifications and Training
 4. Facilities
 5. Evidence Control
 6. Validation
 7. Chemicals and Reagents
 8. Calibration and Maintenance
 9. Proficiency Testing
 10. Corrective Action
 11. Reports
 12. Review
 13. Safety
 14. Audits
 15. Outsourcing
 16. Practices, Methods and Forms
 17. Controlled Forms
 18. Non-Controlled Forms
- Appendix A (Stutter and Peak Height Ratio Studies)

Appendix B (Mixtures- minor component calculations)

Appendix C (DNA Database Implementation)

Appendix D (Staff Abbreviation Lists)

Appendix E (Pre-2009 Revision Histories)

Appendix F (Annual Biology System Review)

Appendix G (Technical Lead Contingency Plan)

Appendix H (FBI Quality Assurance Documents)

Appendix I (Staff Acknowledgements)

Appendix J Extraction Control Keys

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

INTRODUCTION

The Forensic Biology Quality and Procedure Manuals are not public documents. Copies of the manuals, or portions thereof, will be released only to individuals having official business and upon proper discovery requests relating to a specific case(s).

1.0 STATEMENT OF PURPOSE AND OBJECTIVES

1.1 Statement of Purpose: ISP Forensic Biology exists to provide quality, unbiased and cost-effective analyses in the identification of biological substances and their source(s) relevant to the investigation and prosecution of criminal offenses in Idaho. The ISP Forensic Biology Quality Manual, along with the ISP Forensic Services Quality/Procedure Manual, provide the framework for the evaluation of QC (Quality Control) measures utilized in Forensic Biology to achieve that purpose. A system-wide mission and objectives are enumerated in the ISP Forensic Services Quality/Procedure Manual.

1.2 Objectives:

- 1.2.1 To develop and maintain, through annual review and revision (where necessary), a system of quality procedures, analytical methods, and controls to ensure quality up-to-date personnel training, biological screening and DNA analyses.
- 1.2.2 To evaluate (and revise where appropriate) through proficiency testing, audits, and other means of review, the thoroughness and effectiveness of biology personnel training, procedures and QC measures.
- 1.2.3 To remain scientifically neutral by basing case/evidence acceptance and analysis decisions, case reports and testimony solely on sound scientific rationale.
- 1.2.4 To develop and use practices that respect and protect the right of privacy for the genetic profiles developed in forensic casework or for database entry.
- 1.2.5 To provide high quality training, technical and informational assistance, biological analyses, written reports and testimony.
- 1.2.6 To provide all services in a cost-effective and timely manner.

2.0 ORGANIZATION AND MANAGEMENT


2.1 Organizational Chart and Functional Structure

2.1.1 An organizational chart for ISP Forensic Services appears in the ISP Forensic Services Quality/Procedure Manual. The Forensic Biology organization is delineated below.

2.1.2 An organizational chart for the Idaho State Police appears in the ISP Policy Manual.

2.2 Authority and Accountability in Forensic Biology

2.2.1 The Quality Assurance Standards for Forensic DNA Testing Laboratories and Convicted Offender DNA Databasing Laboratories, developed by the DAB, serve as a model for the ISP Forensic Biology QA Program. These standards delineate specific responsibilities and authority for the DNA Technical Manager and DNA CODIS Manager (see standard 4.1 of the FBI quality audit document). A copy of the document may be found in the ISP Forensic Biology Training Manual. Additionally, the ISP Forensic Services Quality/Procedure Manual designates specific authority for the DNA Technical Manager and DNA CODIS Manager.



Note: Changes (personnel) to this page do not require new revision numbers.

3.0 PERSONNEL QUALIFICATIONS AND TRAINING

3.1 Job Descriptions

General personnel qualifications and responsibilities, as well as personnel record retention policies, are described in the ISP Forensic Services Quality/Procedure Manual. Complete job descriptions are available through the Idaho Division of Human Resources web site:

(<http://dhr.idaho.gov/dhrapp/stateJobs/JobDescriptions.aspx>).

3.2 Training

Refer to ISP Forensic Biology Training Manual and the ISP Forensic Services Quality/Procedure Manual for specific training requirements and retention of training and continuing education records.

3.3 Continuing Education

Forensic Biology personnel must stay abreast of developments relevant to forensic DNA analyses through the attendance (and participation) at DNA related presentations, seminars, courses and/or professional meetings, for a minimum of 8 hours per calendar year. The training will also be supplemented through the routine reading of current scientific literature. The DNA technical Manager, or designee, will distribute a DNA-related article to each member of the biology section on a monthly basis. Each staff member will read the article and date/initial the attached sign-off sheet to indicate the completion of the reading. Additionally, the CODIS manager must stay abreast of developments relevant to CODIS/NDIS database management, computer and data security and computer networks through the attendance (personal or that of the Alternate CODIS Manager) at the bi-annual CODIS State Administrators' meetings and annual CODIS conference.

3.4 Qualifications

Education, training and experience for Forensic Biology personnel is formally established in the following minimum requirement specifications (Minimum requirements for individual positions may be reviewed at the time of job announcement and may exceed those delineated below). The minimum degree and education requirements are verified by review of transcripts as well as course descriptions, as necessary, during the application process. The DNA Technical Manager approves the degree and coursework prior to a job offer being extended to any potential hire. Periodic review of continuing education and overall performance is

accomplished during the annual employee evaluation. Opportunities are provided by an FS training budget.

3.4.1 Forensic Biology/DNA Supervisor/Technical Lead

It is assumed for the purposes of this document (and is currently the case), that in a laboratory system of the size of Idaho's, these functions will be served by a single individual.

3.4.1.1 Education

Must have at minimum, a Master of Science degree in a biological science. Successful completion of a minimum of 12 credit hours, including a combination of graduate and undergraduate coursework in genetics, biochemistry, molecular biology and statistics (or population genetics).

3.4.1.2 Training

Training and experience in molecular biology and DNA-based analyses from academic, governmental, private forensic and/or research laboratory(ies). Must also complete the FBI sponsored DNA auditor training within 1 year of appointment, if not already completed (dependant on FBI scheduling).

3.4.1.3 Experience

Must have a minimum of three years forensic human DNA laboratory experience as an analyst.

3.4.2 CODIS Administrator

This function may or may not be served by the Forensic Biology/DNA Supervisor. It is assumed for the purposes of this document (and is currently the case) that in a laboratory system of the size of Idaho's, the functions of casework and database CODIS Administrators will be served by a single individual. An Alternate CODIS Administrator will also be appointed and must meet the same qualifications as the CODIS Manager. The CODIS Administrator is responsible for administering the laboratory's CODIS network, scheduling and documenting the computer training for analysts, as well as assuring the security and quality of data and match dispositions all in accordance with state and/or federal law and NDIS operational procedures.

3.4.2.1 Education

Must have at minimum, a Bachelor of Science degree in a biological science and successfully completed college coursework in genetics, biochemistry, and molecular biology. Must also have completed coursework and/or training in statistics (or population genetics).

3.4.2.2 Training

A combination of training and experience in the use of computers, and database systems in a laboratory/scientific setting. Must also complete the FBI's CODIS software training and the DNA auditor training within six months of appointment if not already completed (dependant on FBI scheduling).

3.4.2.3 Experience

Must possess a working knowledge of computers, computer networks, computer database management and have an understanding of DNA profile interpretation for database and casework functions, to include mixture interpretation. Must be or have been a qualified DNA analyst.

3.4.3 DNA Analyst

The following delineate requirements for a DNA casework or database analyst whose responsibilities include performing genetic analyses on the capillary electrophoresis instruments and data interpretation. DNA extraction, quantification, and amplification set-up may be performed by appropriately trained laboratory technicians and/or those performing the biological screening of evidence following task-specific training and successful completion of a qualifying examination.

3.4.3.1 Education

Must have at minimum, a Bachelor of Science degree in a biological science and successfully completed college coursework in genetics, biochemistry, and molecular biology. Must also have completed coursework and/or training in statistics (or population genetics).

3.4.3.2 Training

Training in DNA analyses through academic, governmental, private forensic and/or research laboratory(ies). If received elsewhere, documented training must meet or exceed that outlined in the ISP Forensic Biology training manual. Must successfully complete a qualifying examination prior to performing analyses on database or forensic casework samples.

3.4.3.3 Experience

Must have a minimum of six months forensic human DNA laboratory experience.

3.4.4 Forensic Biologist

The following delineate requirements for those individuals responsible for the screening of evidence for the presence of biological substances and reporting and giving testimony regarding their findings.

3.4.4.1 Education

Must have a Bachelor of Science in a biological science.

3.4.4.2 Training

Training specific to this job function in a governmental and/or private forensic laboratory. If received elsewhere, documented training must meet or exceed that outlined in the ISP Forensic Biology training manual. Must successfully complete a qualifying examination prior to performing forensic casework.

3.4.4.3 Experience

Prior to participating in independent forensic casework, must have a minimum of six months forensic laboratory experience in the area of biological screening and/or DNA analysis.

3.4.5 Biology Laboratory Technician

3.4.5.1 Education

Minimum of two years of college to include scientific coursework (lecture and lab); Bachelor

of Science degree in a biological science is preferred.

3.4.5.2 Training

Must receive on the job training specific to assigned duties and successfully complete a qualifying examination before participating in forensic DNA typing or forensic casework responsibilities.

3.4.5.3 Experience

Prior to participating in any forensic DNA typing responsibilities or forensic case processing activities, technician must have a minimum of six months forensic laboratory experience in the area of Biology/DNA; one year is preferred.

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

4.0 FACILITIES

4.1 Laboratory Security

Security of the Forensic Services Laboratory is covered in the ISP Forensic Services Quality/Procedure Manual.

4.1.1 Forensic Biology Security

When not under the direct control of Forensic Biology personnel, evidence and in-progress work product will be secured either by closing and locking the Forensic Biology door or by its return to secure storage (one of the locked evidence refrigerators/freezers/file cabinets or the analyst's personal evidence cabinet). Only Forensic Biology Personnel will have access to the locked storage and laboratory areas. Persons outside the Forensic Biology unit will not be allowed access to the Forensic Biology laboratories. Exceptions will be made in case of emergencies, for maintenance, safety, and/or equipment service needs, and for required annual quality and DNA audits. At these times, access will be limited to only required individuals, the individual(s) will be accompanied by biology program personnel, and all evidence will be placed in secured storage for the duration of the individual(s) being present in the laboratory.

4.1.2 CODIS Security

The CODIS workstation is located in the locked CODIS office and the CODIS Server is located in the secured server room in the CJIS Section. The following security measures have been implemented:

4.1.2.1 Only Forensic Biology personnel will have access to the CODIS office. When a biology staff member is not present, the office will be secured by closing and locking the door.

4.1.2.2 Only the CODIS State Administrator, designated Forensic Biology staff and CJIS personnel will have access to the CODIS Server.

4.1.2.3 A differential backup of the CODIS server will be performed each weekday. A full backup will be performed once weekly with the backup tape being stored off-site. At any given time, one month of data will be stored offsite.

- 4.1.2.4 Only Forensic Biology Personnel that have gone through the NDIS application and approval process will have user-names and passwords for CODIS.
- 4.1.2.5 CODIS users must log in each time they use CODIS and log out prior to leaving the CODIS Workstation.
- 4.1.2.6 DNA Tracker, the convicted offender sample-tracking database, resides on the ISP intranet and is accessible only to personnel designated by the Biology/DNA Supervisor.
- 4.1.2.7 Personal and identifying information on convicted offenders (hard and electronic/DNA Tracker copies) are stored separately from the DNA profile (CODIS) obtained. The DNA profiles are directly associated only with a unique Idaho Convicted Offender ID number, assigned by DNA Tracker upon sample entry.
- 4.1.2.8 CODIS samples and corresponding information is released only in accordance with 19-5514 of the Idaho DNA Database Act of 1996, the Privacy Act Notice in Appendix E of NDIS procedures, and the FBI/CODIS Memorandum of Understanding.

4.2 Forensic Biology Laboratory Set-up

The Forensic Biology and Database Laboratories are designed to minimize contamination potential during the processing and analysis of forensic and convicted offender samples. Separate areas for evidence examination, DNA extraction, PCR Amplification Set-up and Amplified DNA processing and storage are delineated. Some steps of the pre-amplification processes may be conducted in the same area of the main laboratory; however, these steps are separated by time.

4.3 Laboratory Cleaning and Decontamination

In order to minimize the potential for sample contamination, careful cleaning of laboratory work areas and equipment must be conducted on a routine basis. The efficacy of the procedures used is monitored through the use of controls within the analysis process (see the interpretation guidelines section in BI-210 and BI-318). It is also important that each analyst use proper 'clean technique' at all times when in the laboratory, which includes but is not limited to, using only disposable barrier pipette tips and autoclaved microcentrifuge tubes, using a tube de-capping tool, and wearing gloves, a labcoat, and masks as appropriate.

- 4.3.1 All working **benchtop** surfaces will be cleaned with 10% bleach or Dispatch solution before and after use and as part of the monthly QC procedure. Clean white paper and/or a KayDry will be placed on the workbench prior to use and changed as appropriate and necessary.
- 4.3.2 All **small tools/instruments** (i.e. forceps, scissors, etc.) will be cleaned/rinsed with ethanol or germicidal instrument cleaner prior to use and between samples. Kimwipes, used to dry the instrument after cleaning/rinsing, will be single use only.
- 4.3.3 **Pipettes** are to be cleaned thoroughly with Dispatch solution as part of the monthly QC procedure and anytime the barrel comes in contact with DNA or any biological fluid.
- 4.3.4 All **centrifuges** are to be wiped down (interior and exterior) with Dispatch solution as part of the monthly QC procedure and in the event of a spill.
- 4.3.5 The **Biomek 3000** work surface trays and holders are to be removed and cleaned with 10% bleach or Dispatch solution as part of the monthly QC procedure or in the event of a spill. Additionally, each of tools are to be wiped down with ethanol, being careful not touch the electronic end.
- 4.3.6 The exterior surfaces of the **BSD600-Duet Puncher** are to be wiped down with a damp cloth, as part of the monthly QC procedure. In addition, the chute and punch mechanism are to be cleaned by removing and separating the inner and outer chutes. The inner chute is to be cleaned with ethanol, followed by compressed air blown through both chutes, the hole in the underside of the manifold, and between the punch guide and die. Do not use ethanol on the outer chute or around any electrical components.
- 4.3.7 The **thermal cyclers**, to include the heating block and exterior surfaces, are to be wiped down with ethanol or Dispatch solution as part of the monthly QC procedure. Individual wells should be cleaned as needed.
- 4.3.8 All work surfaces in the **amplification/post-amp rooms** are to be cleaned with 10% bleach or Dispatch solution before and after analysis and as part of the monthly QC procedure. Clean white paper and/or a KayDry is to be placed on the benchtop prior to use. Additionally, as part of the

monthly QC procedure, the following are to be conducted: the exterior surfaces of the genetic analyzers and real-time instruments wiped down with ethanol or Dispatch solution, top of the refrigerator/freezers and surface underneath each genetic analyzer wiped down/dusted, and floor mopped.

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

5.0 EVIDENCE CONTROL

Evidence, Individual Characteristic Database (Convicted Offender) samples, and in progress work product, that is collected, received, handled, sampled, analyzed and/or stored by ISP Forensic Services is done so in a manner to preserve its identity, integrity, condition and security.

5.1 Laboratory Evidence Control

Procedures detailing evidence handling are contained in the ISP Forensic Services Quality/Procedure Manual. Portions of individual evidence items that are carried through the analysis process (i.e. substrate cuttings, extracts, amplified product and/or portions thereof) are considered work product while in the process of analysis and do not require sealing. Work product will be identified by labeling the individual sample tube with a unique identifier, or documenting the locations of individual samples within a plate of samples.

5.2 Forensic Biology Evidence Control/Sample Retention

5.2.1 DNA Packet

It has become increasingly important to retain evidence for possible future analyses and to secure samples for non-probative casework analyses that are necessary for the validation of any new technology. Therefore, a DNA packet is created for cases submitted for analysis to Forensic Biology, in which reference sample(s) are present, and/or positive Biological screening results are obtained (See BI-102). Any remaining DNA extracts, upon completion of analysis, will be placed into a sealed container (such as a plastic zip bag) and stored in the DNA packet.

5.2.2 Limited Sample

In every case, care should be taken to save ~1/2 of a sample for independent testing. If testing would consume all or nearly all of a sample **and** there is an identified suspect charged in the case, the accused must receive appropriate notification. Written and/or verbal notification will be given to the prosecuting attorney informing him/her of possible consumption and requesting defense counsel be notified of the situation. Before testing will commence, an allowance will be made for testing by another accredited laboratory agreed upon by both parties. Additionally, a letter from the prosecuting attorney must be received

by the laboratory indicating whether or not the sample may be consumed.

5.2.3 Amplified Product

Amplified DNA product will not be retained after 1) the report has been issued in the case or 2) review of the offender sample data has been completed and certified for CODIS entry. In cases where both the evidence and associated DNA extract have been consumed, the amplified product will be retained in a sealed container within the product room freezer.

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

6.0 VALIDATION

Procedures for the validation and/or performance verification of methods used in ISP Forensic Services are outlined in the ISP Forensic Services Quality/Procedure Manual. Validation/performance verification data, results and summaries for those methods employed in Forensic Biology will be maintained in that section.

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

7.0 CHEMICALS/REAGENTS

General laboratory policies and procedures regarding the purchase of chemicals and preparation of reagents are covered in the ISP Forensic Services Quality/Procedure Manual.

7.1 COMMERCIALLY PURCHASED CHEMICALS

7.1.1 Biology Personnel should consult the electronic Chemical Inventory Log (Form 400-QC) prior to ordering. Chemical grade requirements should be checked and ordered as appropriate. The date ordered should be reflected in the log to avoid duplicate orders. An entry for chemicals not currently on the inventory will be made at this time to reflect the chemical, source, and order date. This inventory will be audited annually, at a minimum, and a printout placed in the Forensic Biology Reagent Binder.

Note: An order form/document must be filled out and approved by the section supervisor (indicated by date and initials) prior to placing the order. Reference the forensic services approved chemical list prior to ordering new chemicals.

7.1.2 Upon receipt of a chemical or reagent, the Chemical Inventory Log will be updated to reflect the new lot number, received date, quantity received, and quantity in stock. The order date will be removed at this time. The chemical(s) will be marked with the date received and the individual's initials. If it is an outer container that the chemical/kit remains in until use, the inner container will be labeled with this information when removed for use. The following commercially purchased reagents do not have manufacturer expiration dates: FTA Reagent, Phenol:Chloroform (PCIAA), HiDi Formamide, and 10X Genetic Analyzer buffer. These will additionally be labeled with a laboratory assigned expiration date of 2 years from the date of receipt, with the exception of the FTA which will be 10 years. Packing slips should be checked to ensure appropriate accounting, including proper reagent grade, where applicable (this will be indicated by dating and initialing the packing slip and making notations as necessary). The packing slip and corresponding order document will be retained in the biology section. If an MSDS sheet came with the chemical, the MSDS binder should be checked for the presence of an MSDS sheet for that chemical. If one exists, no additional copy is kept;

however, if a newer version is received, the old one should be replaced. If one does not exist, place one in the binder. For chemicals without an MSDS, consult the manufacturer or one of the following websites for information:

<http://www.hazard.com/msds>

<http://www.msds.com>

<http://www.ilpi.com/msds/>

Note: Critical Reagents listed in 7.3 will be tracked on the individual QC forms, rather than the chemical inventory log.

7.1.3 Expired chemicals will be disposed of in an appropriate manner.

7.2 REAGENTS PREPARED IN-HOUSE

7.2.1 All biology reagents will be made with great care, following all quality and safety procedures. A mask will be worn by analysts during reagent preparation to help avoid the potential for contamination. See 7.4 and 7.5 below for individual reagent recipes.

7.2.2 Each reagent has a corresponding form to document the making of the reagent and components used. This form must be filled out. A reagent label must be made that has the reagent name, the lab lot number (which consists of the first few letters of the reagent name followed by the date prepared, in the form 'MMDDYY'), and the preparer's initials. The NFPA designation will be completed on all labels. Refillable squirt-bottles of water or ethanol will be labeled but need not bear dates or initials.

7.2.3 An effort should be made to use in-house reagents within one year of preparation; however, they do not expire and may continue to be used beyond the one year timeframe.

7.3 CRITICAL REAGENTS

CRITICAL REAGENTS are those reagents that, if improperly functioning, could result in significant loss or destruction of DNA and are not amenable (or it's not practical) to testing immediately before (e.g., use on forensic samples) each use. The reagents listed below have been identified as critical in Forensic Biology/DNA. These reagents must undergo a **QC ASSAY**

BEFORE use on forensic casework and/or Convicted Offender samples. Reagents received at a later date but having the same lot number as those previously tested and determined acceptable need not have a QC check performed. Critical Reagents (in addition to other DNA-related reagents with manufacturer expiration dates) may be used beyond the listed expiration date for training purposes without any further testing, so long as expected results are obtained for all associated controls. The reagent must be labeled 'for training only' if it is to be retained once the expiration date has been reached.

ABACARD[®] HEMATRACE[®] TEST KIT (Form 410-QC)

OneStep ABACARD[®] p30 TEST KIT (Form 412-QC)

Quantifiler[®] Human DNA Quantification Kit (Form 419-QC)

STR Kit (Taq Polymerase checked with kits; Form 420-QC)

7.4 BIOLOGICAL SCREENING REAGENTS

Phenolphthalein (Kastle-Meyer) Reagent

(NFPA: health 3, flammability 1, reactivity 2)

May be a commercial purchase.

Phenolphthalein	2.0g
KOH	20.0g
Zinc (granular)	20.0g

Phenolphthalein, KOH, and 100ml of dH₂O are refluxed, in a fume hood, with Zinc until solution is colorless (producing phenolphthalin in ~4 hours). Store stock solution refrigerated in dark bottle to which ~5g mossy zinc has been added to keep the solution in its reduced form. Remove for working solution as needed.

Working solution: Mix 2ml stock solution with 8ml Ethanol

Caution: Zinc is flammable. The unreacted portions and used filter paper are to be disposed of properly.

Hydrogen Peroxide 3% (v/v)

(NFPA: health 0, flammability 0, reactivity 1)

Generally a commercial purchase, however, may be made from a 30% Solution (which is a commercial purchase) as follows:

Hydrogen Peroxide (30%) 10ml/90ml nanopure dH₂O

Mix the H₂O₂ with 90ml of nanopure dH₂O and store at ~4°C.

Ortho-Tolidine Reagent

(NFPA: health 3, flammability 1, reactivity 2)

O-Tolidine	0.6g
Glacial Acetic Acid	100ml
Ethanol	100ml

Dissolve O-tolidine in Acetic Acid/Ethanol mixture consistent with ratios above. O-tolidine is light sensitive and should be stored in dark reagent bottle and kept refrigerated when not in use.

Ammonium Hydroxide (~3%)

(NFPA: health 3, flammability 1, reactivity 2)

Ammonium Hydroxide (Concentrated ~30%) 10ml/100ml

Add the NH₄OH to 90ml of nanopure dH₂O, mix well and store at RT.

Ouchterlony Destain

(NFPA: health 3, flammability 3, reactivity 2)

Methanol	45ml
Distilled water	45ml
Glacial Acetic Acid	10ml

Mix well and store refrigerated.

Ouchterlony Stain

(NFPA: health 3, flammability 3, reactivity 2)

Ouchterlony Destain	50ml
Coomassie Blue (Brilliant Blue R-250)	0.1g

Mix well (overnight), filter, and store at RT.

10X Brentamine (Sodium Acetate) Buffer

(NFPA: health 2, flammability 2, reactivity 2)

Sodium Acetate (Anhydrous) 1.2g
Acetic Acid(to adjust to pH 5) ~400µl

Dissolve Sodium Acetate in 10ml of nanopure dH₂O. Add Acetic Acid to pH 5. Store refrigerated.

Brentamine Solution A

(NFPA: health 1, flammability 0, reactivity 0)

O-Dianisidine Tetrazotized (Fast Blue B Salt) 50 mg
10X buffer pH 5 5 ml

Dissolve Fast Blue B Salt in 5 ml of 10X Brentamine Buffer. Store refrigerated in a dark container.

Brentamine Solution B

(NFPA: health 2, flammability 0, reactivity 0)

α-Naphthyl Phosphate (Disodium Salt) 50 mg

Dissolve in 5 ml of nanopure dH₂O. Store Refrigerated.

Saline (0.85% NaCl)

(NFPA: health 1, flammability 0, reactivity 0)

NaCl 4.25g/500ml

Dissolve the NaCl in 500 ml nanopure water. Sterilize by autoclaving. Store refrigerated.

1X Phosphate Buffered Saline (PBS)

(NFPA: health 1, flammability 0, reactivity 1)

PBS 1 commercial pre-made packet

Dissolve one packet of powdered PBS in 1ℓ of nanopure dH₂O. Check that pH≈7.4, autoclave and store at RT.

If pre-made packets are not available, PBS may be prepared by dissolving 0.2g KCl, 8.0g NaCl, 0.2g KH₂PO₄, and 2.2g Na₂HPO₄·7H₂O (or 1.1g Na₂HPO₄ anhydrous) in 800ml nanopure dH₂O. Adjust pH to 7.4 if necessary. Q.S. to 1ℓ with nanopure dH₂O, autoclave and store at RT.

X-mas Tree Stain Solution A (Kernechtrot Solution)

(NFPA: health 1, flammability 0, reactivity 0)

May be a commercial purchase.

Aluminum Sulfate	5g
Nuclear Fast Red	0.1g

For 100ml, Dissolve the Aluminum Sulfate in 100ml **HOT** nanopure dH₂O. Immediately add the Nuclear Fast Red, mix, cool and filter (paper or ≥45µm). May be stored at RT.

X-mas Tree Stain Solution B (Picroindigocarmine Solution)

(NFPA: health 2, flammability 2, reactivity 2)

May be a commercial purchase.

Saturated Picric Acid Solution	100ml
Indigo Carmine	0.33g

For 100ml, dissolve the Indigo Carmine in 100ml of the Picric Acid. Mix and filter (paper or ≥45µm). May be stored at RT.

Amylase Diffusion/Phosphate Buffer (pH 6.9)

(NFPA: health 1, flammability 0, reactivity 1)

NaH ₂ PO ₄ , anhydrous	2.7g
Na ₂ HPO ₄ , anhydrous	3.9g
NaCl	0.2g

Mix the above with 500ml dH₂O, adjust pH to 6.9, and store at RT.

Amylase Iodine Reagent

(NFPA: health 3, flammability 0, reactivity 2)

Potassium Iodide (KI)	1.65g
Iodine (I ₂)	2.54g

Dissolve the above in 30ml nanopure dH₂O heated to ~65°C. Mix well, filter and store at 4°C in an amber bottle. Dilute 1:100 for Amylase Diffusion Test.

Mercuric Chloride 10% (w/v)

(NFPA: health 4, flammability 0, reactivity 1)

Mercuric Chloride	10g/100ml 95% EtOH
-------------------	--------------------

Dissolve the Mercuric Chloride in 100ml of 95% Ethanol, mix well and store at RT.

Zinc Chloride 10% (w/v)

(NFPA: health 2, flammability 0, reactivity 2)

Zinc Chloride	10g/100ml 95% EtOH
---------------	--------------------

Dissolve the Zinc Chloride in 100ml of 95% Ethanol, mix well and store at RT.

7.5 DNA REAGENTS

1M Tris-HCl Buffer pH 7.5

(NFPA: health 2, flammability 1, reactivity 1)

Tris Base(tris[Hydroxymethyl]amino methane)	121.1 g
---	---------

Dissolve Tris in ~800 ml nanopure dH₂O. Adjust to pH7.5 at RT by adding concentrated HCl (approximately 65ml). Q.S. to 1l with nanopure dH₂O, autoclave and store at RT.

1M Tris-HCl Buffer pH 8

(NFPA: health 2, flammability 1, reactivity 1)

Tris Base(tris[Hydroxymethyl]amino methane)	121.1 g
---	---------

Dissolve Tris in ~800 ml nanopure dH₂O. Adjust to pH8 at RT by adding concentrated HCl (approximately 45ml). Q.S. to 1l with nanopure dH₂O, autoclave and store at RT.

0.5M Ethylenediamine Tetraacetic Acid (EDTA)

(NFPA: health 1, flammability 1, reactivity 0)

Na₂EDTA·2H₂O 186.1g/l

Slowly add EDTA to 800ml nanopure H₂O while stirring vigorously. Add ~20g of NaOH pellets to bring the pH to near 8.0. When fully dissolved adjust pH to 8.0 and bring final volume to 1l. Autoclave and store at RT.

Note: EDTA will not go into solution without the pH adjustment.

Stain Extraction Buffer pH8 (10mM EDTA/10mM Tris-HCl/50mM NaCl/2% SDS)

(NFPA: health 2, flammability 1, reactivity 1)

1M Tris-HCl, pH7.5	5ml
0.5M EDTA	10ml
5.0M NaCl	5ml
10% SDS	100ml

Mix the Tris-HCl, EDTA, NaCl and SDS with ~380ml nanopure dH₂O. Store at RT.

Note: Reagent contains SDS, do not autoclave.

Proteinase K (20mg/ml)

(NFPA: health 1, flammability 1, reactivity 0)

May be a commercial purchase of 20mg/ml solution.

Proteinase K 0.2g

Dissolve the ProK in 10ml sterile nanopure dH₂O.

Dispense ~500µl (commercial purchase or in-house prep.) each into sterile microfuge tubes and store at ≈20°C.

1M Sodium Acetate pH 5.2

(NFPA: health 3, flammability 2, reactivity 0)

CH₃COONa·3H₂O 13.6g

Dissolve the CH₃COONa·3H₂O in 80ml nanopure dH₂O. Adjust to pH5.2 by adding glacial acetic acid (approximately 2 ml). Q.S. to 100ml with nanopure dH₂O, autoclave and store at RT.

DTT Solution

(NFPA: health 2, flammability 1, reactivity 0)

Dithiothreitol (DTT) 0.77g

Dissolve the DTT in 5ml nanopure dH₂O. Add 50µl 1M Sodium Acetate, pH5.2. Dispense ~500µl each into sterile microcentrifuge tubes and store at ≈20°C.

Note: Do not autoclave.

PCR-TE (TE⁻⁴) Buffer (10mM Tris-HCl/0.1mM EDTA)

(NFPA: health 2, flammability 1, reactivity 0)

1M Tris-HCl, pH8 10ml
0.5M EDTA, pH8 0.2ml

Mix Tris-HCl and EDTA with 990ml nanopure dH₂O. Autoclave and store at RT.

5N Sodium Hydroxide

(NFPA: health 3, flammability 0, reactivity 2)

NaOH 50g

Slowly dissolve the Sodium Hydroxide in 250ml sterile nanopure dH₂O. Allow to cool and store at RT.

Caution: NaOH is highly caustic. This reaction generates heat.

5M Sodium Chloride

(NFPA: health 1, flammability 0, reactivity 0)

May be a commercial purchase of 5M solution.

NaCl 146.1g/500ml

Dissolve the NaCl in 500 ml nanopure water. Sterilize by autoclaving.

Bovine Serum Albumin 4%

(NFPA: health 0, flammability 1, reactivity 0)

BSA 0.4 g
PCR-TE 10 ml

Dissolve the BSA in PCR-TE. Filter sterilize and dispense ~500µl each into 1.5ml microfuge tubes. Store at -20°C.

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

8.0 EQUIPMENT CALIBRATION AND MAINTENANCE

General laboratory procedures for the calibration and maintenance of equipment are covered in the ISP Forensic Services Quality/Procedure Manual.

8.1 BIOLOGY EQUIPMENT/INSTRUMENTATION

8.1.1 Analytical equipment significant to the results of examination and requiring routine calibration and/or performance verification will be listed on the BIOLOGY CRITICAL EQUIPMENT INVENTORY Spreadsheet (Form 401-QC). Information on the spreadsheet includes (as known or appropriate): equipment identity and its software, manufacturer's name, model, property number, serial number and/or unique identifier, and location. The inventory spreadsheet will be maintained in the Instrument QC binder or Biology QC binder as appropriate.

8.1.2 OPERATING MANUALS for most equipment/instrumentation will be maintained in the product information file (Manuals for the ABI PRISM™ 310 and 3130/3130xl Genetic Analyzers, ABI 7500 Real-Time PCR System, Thermal Cyclers, and Driftcon FFC will be maintained in the Amp/PostAmp Room in close proximity to the instruments). Exceptions may be made for manuals referred to for instructions. In these cases, the manual will be maintained in close proximity to the instrument. The Biomek 3000 manual is built into the Biomek software.

8.1.3 MAINTENANCE/REPAIR/CALIBRATION LOGS will be maintained as follows:

The records for the ABI PRISM™ 310 and 3130/3130xl Genetic Analyzers, ABI 7500 Real-Time PCR System, and Thermal Cyclers will be maintained in the instrument QC binder.

Any equipment/instrumentation function (not documented on weekly, monthly, quarterly, or annual QC Check forms) will be recorded on the Equipment Maintenance/Repair form (Form 402-QC). Equipment Failure will also be reported on this form. This form and the QC check forms will be maintained in the Biology QC Binder, except as listed above.

8.1.4 EQUIPMENT FAILURE will result in that equipment being 'taken out of service'; an 'out of service' sign will be placed on the equipment and it will not be returned to service until it has passed appropriate performance testing. Actions are reported on Form 402-QC.

8.1.5 The SCHEDULE of QC/Performance Checks for both critical and non-critical equipment is as follows:

WEEKLY (Form 404A/B-QC)

(once per week with an interval between dates not less than 3 days and not exceeding 10 days)

- Nanopure System Check
- Refrigerator/Freezer Temperature Check
- Heating Block(s) Temperature Check
- Oven Temperature Check

MONTHLY (Form 406A/B-QC)

(once per month with an interval between dates not less than 15 days and not exceeding 45 days)

- Pipettes Cleaned
- Centrifuges Cleaned
- Biomek 3000 Cleaned
- BSD600 Cleaned
- Lab Cleaned
- Eye Wash Station Check
- Autoclave Clean and Check Sterilization
- ABI 7500 Background Assay/Contamination Test, and Function Test/Bulb Check
- BioRobot EZ1 grease D-rings
- 3130/3130xl Water Wash
- 3130/3130xl Water Trap Flush
- 310, 3130/3130xl, (C and E drives) and 7500 computer defragmentation

QUARTERLY

(once per quarter with an interval between dates not less than 30 days and not exceeding 120 days) Note: * denotes critical equipment

- Thermal Cycler* Temperature Verification
- ABI 7500* Temperature Verification

- Biomek 3000 Robotic System* Framing/Calibration Check (Form 408-QC)
- Chemical Shower Check (Form 408-QC)

ANNUALLY (Form 402-QC)

(once per calendar year with an interval between dates not less than 6 months and not exceeding 18 months) Note: * denotes critical equipment

- Mechanical Pipette* Performance Verification Check (outside vendor)
- NIST Traceable Thermometers* (outside vendor)
- Driftcon FFC Temperature Verification System* (outside vendor)
- Biological and Chemical Hoods Test (outside vendor)
- Digital Temperature Recording Devices Calibration Check (outside vendor)
- ABI PRISM™ 310* Genetic Analyzer Preventative Maintenance (outside vendor)
- ABI PRISM™ 3130/3130xl* Genetic Analyzer Preventative Maintenance (outside vendor)
- ABI 7500* Real-Time PCR System Preventative Maintenance (outside vendor)
- ABI 7500* Pure Dye Calibration, Optical Calibration, and Regions of Interest (ROI's) verification (see 7500 Maintenance Guide for procedures/may be part of PM by request)
- Qiagen BioRobot EZ1* Preventative Maintenance (outside vendor)
- Biomek 3000* preventative Maintenance (outside vendor)
- Microscope Cleaning/Preventative Maintenance (outside vendor)
- Centrifuge Calibration Check (outside vendor)
- Balance* Calibration Check (outside vendor)

In addition to the above schedule, personnel should check appropriate parameter function on all instrumentation with each use (including calibration of the pH meter at the time of use; documented on Form 403-QC), and run a matrix for the ABI PRISM™ 310 Genetic Analyzers and a spatial and spectral calibration for the ABI PRISM™ 3130/3130xl Genetic Analyzers as needed or following CCD camera and/or laser replacement/adjustment. Data for each new matrix will be filed in the instrument QC binder (see BI-210 and BI-318).

Following the annual preventative maintenance, a sensitivity panel (previously characterized DNA) should be run on the 310 and 3130/3130xl and included in the QC binder as a verification of performance. A color plate and framing/calibration check are to be run on the Biomek 3000, documented on Form 428-QC, and included in the QC binder as a performance check following the annual preventative maintenance.

Any problems noted with laboratory equipment, during normal usage or as part of a QC check should be brought to the attention of the necessary supervisory personnel and documented on Form 402-QC and/or the respective QC form.

A certified NIST standard will also be run annually or if substantial procedural changes have been made. The QC run will be documented on Form 426-QC and filed in the QC binder.

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

9.0 PROFICIENCY TESTING

General laboratory guidelines and practices for proficiency testing and retention are outlined in the ISP Forensic Services Quality/Procedure Manual. Additional Biology/DNA requirements are delineated below.

9.1 External DNA Proficiency Test Requirement. DNA analysts will participate in external proficiency tests, twice in every calendar year, in accordance with The FBI Quality Assurance Standards and the results reported to NDIS as necessary.

9.2 Inconclusive/Uninterpretable Proficiency Test Results.

Typically, sample size/quantity in PCR DNA Proficiency Tests is sufficient for multiple analyses to be performed. Therefore, results of DNA proficiency tests are not likely to be either inconclusive, or uninterpretable (e.g., not meeting minimal rfu and/or statistical threshold for inclusion/exclusion). However, in the event data obtained in a proficiency test does not meet the standard guidelines for interpretation/conclusion, it will first be determined, by re-testing and communication with the vendor, that this is not an issue with a given sample(s). Once that determination has been made, the analyst obtaining the inconclusive data will be removed from casework/database sample analysis until satisfactory completion of a competency test and review of the analyst's casework/database analysis performed since the last successful proficiency test.

10.0 CORRECTIVE ACTION

Laboratory corrective-action and retention procedures are detailed in the ISP Forensic Services Quality/Procedure Manual.

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

11.0 FILE DOCUMENTATION AND REPORTS

Meticulous documentation is an important aspect of forensic work. In casework, the scientist's knowledge of case circumstance (and therefore their ability to discern potential significance) may be limited. It is also common to be called upon to testify months, or even years, after processing evidence for a given case. Careful observation and detailed note-taking will not only refresh the scientist's memory and provide support for the conclusion in the laboratory report, but might also provide additional information not thought to have been important at the time of evidence processing. General laboratory policies regarding case record and retention are described in the ISP Forensic Services Quality/Procedure Manual. The note packet is considered complete when the analyst signs the report and submits the packet to be reviewed. Electronic documentation (eg. electropherograms and tables of results) are considered stored at this time. Any changes to the electronic documentation required after this point (typically on or after the review date documented in the note packet) will be made either by hand on the hard copy (initialed and dated by the analyst), or by changing the electronic version, reprinting and making a notation on the new hard copy as to the changes made. The new printed copy will bear the date the changes were made/reprinted.

11.1 CASE NOTES

- 11.1.1 Each page of case notes should have the following: Laboratory Case Number, Date, Scientist's Initials and page number (in a form indicating page/total pages).
- 11.1.2 Case notes are associated with a particular report. Case notes for additional submissions (i.e., for supplemental reports) will be reflected in the page numbering as well (e.g. s1, supp. 1, etc.).
- 11.1.3 All evidence submitted for biological screening should be transferred to the scientist (i.e., documented on the chain of custody) and bear the scientist's initials. This is the case regardless of whether or not they analyze the item of evidence (exception may be made in cases where communication with investigator/attorney identified select items of those submitted). A description of the evidence (e.g., packaging and what it is said to contain) should

also appear in the case notes with a notation about not being examined at the time, if that's the case. Those items should also appear in the "not examined" statement of the report.

- 11.1.4** The description of evidence packaging should include the type and condition of seal(s). Differences in the description on a package versus ETS entry and/or accompanying submission form (or what the evidence is once opened) should be noted.
- 11.1.5** Whenever feasible, every attempt should be made to gain entry into the evidence without breaking the original seals. Any seal altered or created by a scientist will bear their initials and date across the seal.
- 11.1.6** Evidence descriptions should be "unique" inasmuch as possible (i.e., one pair blue jeans is **NOT** adequate). They should include, as appropriate and necessary for identification, colors, sizes (measurements where appropriate- e.g., knife and blade), manufacturer, model, brand, serial numbers or other identifiers and condition (e.g., worn, clean, torn, mud-soaked, blood-soaked, etc.).
- 11.1.7** Photography, digital or otherwise, is often useful in documenting the appearance of evidence items. However, it is not meant to completely replace drawing, but instead as a supplement or in cases when drawing may be too difficult to accurately depict the item. Careful drawing and description result in careful and detailed examinations and, in many instances, may be a better choice than photography. Digital photographs will be transferred to, printed as necessary for case notes, and stored within the Mideo System; refer to BI-119 for Mideo instructions.
- 11.1.8** Evidence numbering must be unique for the purpose of possible later CODIS entry. Items should be numbered as follows (or other similar system):

A single item (e.g., a baseball cap; Item 57) for which:

≤ 1 area tested positive for a biological substance ≡ Item 57

≥ 2 areas tested positive for a biological substance(s) (in this instance 3 areas)

≡ Item 57-1, Item 57-2 and Item 57-3, or 57-A, 57-B and 57-C.

An item with multiple sub-items

(e.g., a SAECK; Item 1)

≡ Item 1A, Item 1B, Item 1C, etc., the scientist should begin with the most relevant item if possible. Multiple areas ≡ Item 1A-1, Item 1A-2 etc.

11.1.9 The Biology Screening Case Summary Form (Form 101-BI) may be used for summarizing analyses if the scientist chooses.

11.1.10 If a form is used for more than one case, a copy of the 'completed' form should be made for any additional case files. Each copy should contain a reference regarding the location (case file) of the original document. For each file, the associated case should be listed and case data highlighted. In general, biology subfolders should be organized from front to back as follows: restitution where applicable, report, chronological case notes/forms, SAECK form where applicable, CODIS entry forms where applicable, case review forms where applicable, copy of evidence submission form or ETS property form, phone/info log ('tangerine' paper may be used for ease of identification), followed by agency materials submitted with evidence. Upon completion of review the analyst should bind (e.g. staple) the documentation together, with the exception of the restitution and report, and submit to the Forensic Evidence Specialists for report/restitution distribution.

11.2 DATABASE PACKETS

11.2.1 Each page of the database packet should have the following: Plate Identifier, Date, Scientist's Initials, and page number (in a form indicating page/total pages).

11.2.2 In general, database packets will be arranged from front to back as follows: chronological worksheets/forms, table of results (it is not necessary to print electropherograms for database packets), followed by applicable review forms. Upon completion of review, the analyst should bind (e.g. staple) the documentation together and file it appropriately.

11.3 CASEWORK REPORTS

In the interest of consistency and clarity of reports between individual scientists the following format should be adhered to:

- 11.3.1 The report will contain the title Forensic Biology Report for biology screening reports, or Forensic DNA Report for DNA reports.
- 11.3.2 For clarity, when a statement(s) is about a particular Item (or multiple items listed individually), the "I" will be capitalized as in a name. When writing in general terms (i.e., the following items:) the "i" will remain lowercase.
- 11.3.3 The case submission information will include, at a minimum: case#, report date, case agency, agency case#, principals (victim, suspect, etc.), and offense date.
- 11.3.4 The body of the report will be separated from the case submission information by the following headings in the format below:

RESULTS AND INTERPRETATIONS

Statements (see below) regarding evidence exam, results and conclusions. The order of statements should be, inasmuch as possible: 1) positive statements (detection of body fluid), 2) inconclusive statements, 3) negative statements and 4) statements regarding (i.e. a list of) items not examined.

Disposition of Evidence

Statements (See below) regarding evidence retention and return.

Evidence Description

The following items were received in the laboratory via Federal Express (UPS, US Mail, etc.) on Month day, year. (or) The following items were received in the laboratory from Agency Representative (Agency) on Month day, year.

Description of items submitted for examination.

In the first report, all items should be listed (any items scientist took possession of, including reference samples). In supplemental reports, only those items relevant to the additional examinations need to be listed.

DNA reports, in which a DNA packet is checked out for analysis, will state: A tape sealed DNA packet envelope, created in the laboratory on Month day, year, and containing the following items:

Description of items contained within the DNA packet.

This report does or may contain opinions and/or interpretations, of the undersigned analyst, based on scientific data. The analyst's signature certifies that all of the above are true and accurate. (Note: the interpretations statement does not need to be included in reports where all items submitted are being returned without analysis, or other instances when no conclusions or interpretations are made.)

Signature

Name of Scientist
Title of Scientist

11.3.5 The following results/conclusions statements are to be used in a **biology screening report**, as dictated by the analysis findings (Where appropriate, descriptions, quantity, and/or locations of individual stains may be included in the corresponding statements. Portions of individual statements may be combined as needed.):

Semen Results/Conclusions Statements:

Chemical and microscopic analyses for the detection of semen were conducted on (items). Semen was confirmed by the presence of spermatozoa on (items). (or) Semen was not detected on (items). (or) No identifiable spermatozoa were detected on (items).

Chemical and microscopic analyses for the detection of semen were conducted on (items). Semen was confirmed on (items) by the presence of a single spermatozoon (or limited number of spermatozoa), which is (or may be) insufficient for further testing at this time.

Chemical, microscopic, and serological analyses for the detection of semen were conducted on (items). Semen was detected on (items) by the presence of the semen specific protein, p30; however, no spermatozoa were observed, which is insufficient for further testing at this time.

Results from presumptive chemical tests for the presence of semen were negative on (items).

Blood Results/Conclusion Statements:

Results from chemical and serological tests performed on (items) indicated the presence of human (or non-human) blood.

Results from presumptive chemical tests performed on (items) indicated the presence of blood; however, serological tests to determine the species of origin were not performed (or were inconclusive).

Results from presumptive chemical tests for the presence of blood were negative on (items).

Saliva Results/Conclusions Statements:

Results from chemical tests performed on (items) indicated the presence of an elevated level of amylase, an enzymatic component of saliva.

Results from chemical tests performed on (items) indicated (or did not indicate, or were inconclusive for) the presence of amylase, an enzymatic component of saliva.

Urine Results/Conclusions Statements:

Results from presumptive chemical tests performed on (items) indicated (or did not indicate, or were inconclusive for) the presence of urine.

Feces Results/Conclusions Statements:

Results from presumptive chemical tests performed on (items) indicated (or did not indicate, or were inconclusive for) the presence of feces.

Further Testing Statements (to be included at the end of the Results of Examination Section):

If additional testing is desired, please contact the laboratory.

DNA testing can be performed (or may be attempted) upon request and submission of a known blood sample(s) from [list name(s)]. Please contact the laboratory regarding the analysis request.

11.3.6 The following results/conclusions statements are to be used in an **STR DNA Report** (**Note:** the epithelial cell fraction of intimate samples, such as vaginal/rectal swabs, etc., are not considered probative if the testing results in a single profile matching the individual from which the sample was collected. In these instances, a statement regarding the DNA source of this fraction is not required):

Deoxyribonucleic Acid (DNA) Analysis, employing the Polymerase Chain Reaction (PCR), was used to generate a Short Tandem Repeat (STR) profile from the following items: "list of items".¹

Note: The following footnote will appear in all reports in which DNA testing was attempted.

¹Loci Examined: D3S1358, TH01, D21S11, D18S51, Penta E, D5S818, D13S317, D7S820, D16S539, CSF1PO, Penta D, vWA, D8S1179, TPOX, and FGA.

Profile Match Statement [meeting the 'source attribution' criterion (estimated frequency in population of ≤ 1 in 1.6×10^{10})] for single source and identifiable major contributors of a mixture:

The DNA profile obtained from the "item description (Item #)" matches

that obtained from the blood stain/sample (or reference oral swab/sample, etc.) of/from "name". Therefore, "name" is the source of the "(DNA, blood, semen, saliva etc.) " on this item².

Note: The following footnote will appear in any report containing the above match statement.

²This conclusion is based upon the following: 1) a genetic match at the gender identity locus, Amelogenin, in addition to the "number" polymorphic STR loci listed above that have an expected population frequency of at least less than 1 in "actual (most conservative of the population groups calculated) frequency estimate", 2) a statistical frequency exceeding the source attribution criterion of 1.6×10^{10} (for $N = 1.6 \times 10^7$, $\alpha = 0.01$; *Forensic Science Communications* 2(3) July 2000), and 3) that "name" does not have a genetically identical twin.

Profile match Statement [not meeting the 'source attribution' criterion (estimated frequency in population of greater than 1 in 1.6×10^{10})] for single source and identifiable major contributors of a mixture:

The DNA profile obtained from the "item description (Item #)" matches that obtained from the blood/oral sample of "name". The probability of selecting an unrelated individual at random from the general population having a DNA profile that would match the DNA profile obtained from "item description (Item #)" is at least less than one in "actual (most conservative of the population groups calculated) frequency estimate".

Partial Profile Statement [profile consistent with item(s) in match statement above]:

The DNA profile obtained from the "item description (Item #)" also matches that obtained from the blood/oral sample of "name", however less genetic information was obtained.

The partial DNA profile obtained from the "item description (Item #)" is consistent with that obtained from the blood sample of "name".

Positive Paternity Statement [profiles consistent with being a biological child]:

Based upon evaluation of the DNA profiles obtained from the above individuals, "name" cannot be excluded as being the biological father of "name". The probability of paternity (assuming a prior probability of 0.5) is "X%" relative to an unrelated man randomly selected from the general population. The combined paternity index for the loci examined is "X". At least "X%" of the male population would be expected to be excluded from the possibility of being the biological father of "name".

Note: The most conservative of the population groups calculated is reported for the statement above.

Mixture Statements:

The DNA profile from "item description (Item#)" indicates a mixture of DNA from at least "X" persons. "Name(s)" is a potential contributor(s) to this mixture. "X%" of unrelated individuals randomly selected from the general population would be expected to be eliminated as potential contributors to this mixture.

The DNA profile from "item description (Item#)" indicates a mixture of DNA from at least two persons. "Name(s)" is a potential contributor(s) to this mixture. The DNA profile obtained from "item description (Item#)" is at least "X" times more likely to be seen if it were the result of a mixture of DNA from "name and name" than if it resulted from "name" and an unrelated individual randomly selected from the general population.

Note: The most conservative of the population groups calculated is reported for the statement above.

The DNA profile from "item description (Item#)" indicates a mixture of DNA with a discernable major contributor/profile. (include match, consistent with, or exclusionary statement regarding major profile). "name" is included/excluded/cannot be excluded as a possible contributor to the minor DNA component of this mixture.

The DNA profile from "item description (Item#)" indicates a mixture of DNA from at least "X" persons. "Name(s)" is a potential contributor(s) to this mixture. At least one in "actual (most conservative of the population groups calculated) frequency estimate" of unrelated individuals randomly selected from the general population would be expected to be included as potential contributors to this mixture.

Exclusionary Statement:

The DNA profile obtained from the "item description (Item #)" does not match that obtained from the blood sample of "name". Therefore, "name" is not the source (or "a contributor" in a mixed profile situation) of the "(DNA, blood, semen, saliva etc.)" on this item.

The DNA profile obtained from the "item description (Item #)" was determined to be from an unknown male/female. "name" is not the source of the "(DNA, blood, semen, saliva etc.)" on this item.

Based upon evaluation of the DNA profiles obtained from the above individuals, "name" is not the biological father of "name".

No DNA Profile Obtained Statement:

Due to insufficient quantity or degradation, no DNA profile was obtained from "item description (Item #)".

CODIS Entry Statement:

The unknown male/female (included if source is not identified) DNA profile obtained from the "item description (Item #)" was entered into the Combined DNA Index System (CODIS) to be routinely searched against the database. The case agency will be notified in the event of a profile match.

Note: This statement is included when an eligible DNA profile has been developed, regardless of whether the profile is from a known or unknown source. Eligibility of forensic profiles for entry into CODIS and upload to NDIS is according to current NDIS procedures and include both solved and unsolved cases in which the profile is associated with a crime and believed to be attributable to the putative perpetrator. Profiles matching the victim(s) and any elimination samples (e.g. consensual partner samples) may not be entered.

11.3.7 The following statements are to be used in both biology screening and DNA STR reports:

Evidence Disposition Section Statements:

The following items have been retained in the laboratory [list all items/portions by description and Item# that have been retained in the DNA Packet (see BI-102)]. All remaining items have been returned to the main laboratory evidence vault for return to the submitting agency.

The following items have been forwarded for DNA analysis: [list all items/portions by description and Item# that have been retained in the DNA Packet (see BI-102)]. Results will follow in a separate report. All remaining items have been returned to the main laboratory evidence vault for return to the submitting agency.

Note: Nonsuspect cases (those with no known/identified suspect) in which biological evidence has been detected, will be forwarded for DNA testing and CODIS entry.

The DNA packet, which contains any remaining DNA extracts, has been retained in the laboratory. All remaining items have been returned to the main laboratory evidence vault for return to the submitting agency.

Evidence Description Section Examples:

A tape-sealed Sexual Assault Evidence Collection Kit (SAECK) containing biological samples, said to have been collected from "name".

A tape-sealed brown paper bag/manila evidence envelope/white cardboard box/etc. containing "description", (include the following if collection information is known) said to have been collected from "name" or "location".

A tape-sealed brown paper bag/manila evidence envelope/white cardboard box/etc. said to contain "label on package", (include the following if collection information is known) collected from "name" or "location".

A tape-sealed DNA packet, created in the laboratory on month day, year, and containing the following items:

- Item #) "description"
- Item #) "description"

11.3.8 It should be noted that the statements (in either the Forensic Biology Screening or DNA Reports) regarding

evidence examination, testing and conclusions are not all-inclusive. There may be situations for which none of these statements is optimum.

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

12.0 REVIEW

Technical/administrative, document, and testimony (to include retention) review; as well as conflict resolution is addressed in the ISP Forensic Services Quality/Procedure Manual. See also, forms 214-BI and 306-BI in this manual.

12.1 BIOLOGY/DNA CASEWORK REVIEW

- 12.1.1 100% of the examinations and reports documented and/or issued from Forensic Biology/DNA will be "peer-reviewed". This review must be completed prior to issuing results (including verbal results) and/or entering eligible profiles into CODIS. Exceptions for release of results may be made on a case-by-case basis and with the Biology Supervisor's approval.
- 12.1.2 "Peer-review" in Forensic Biology will encompass both technical and administrative reviews.
- 12.1.3 The individual performing the "peer-review" will be a second scientist who is "qualified" in the area of the review (i.e., Biological Screening and/or STR Analysis).
- 12.1.4 It is not sufficient to have the scientist performing/reporting the analysis to be the sole person performing the administrative review.
- 12.1.5 The second scientist performing the review will initial each page (and date the first and last page at a minimum).
- 12.1.6 The second scientist will also place their initials below the signature of the scientist issuing the report.
- 12.1.7 Additionally, the second scientist will review the CODIS Entry Form (Form 218-BI) and verify that all eligible profiles have been identified for CODIS entry and the correct specimen categories have been assigned. The reviewer will date and initial the form. Eligible specimens will not be entered into CODIS until review/verification is complete. The specimen details report will be reviewed and initialed by the CODIS Administrator (or alternate) following manual data entry

and prior to searching at SDIS and uploading to NDIS to verify correct allele entry and specimen category.

- 12.1.8 Outsourced casework (when applicable) will undergo the same review as listed above, as well as for compliance with contract technical specifications.

12.2 CONVICTED OFFENDER/DATABASE SAMPLE REVIEW

- 12.2.1 100% of Convicted Offender sample data (including outsourced data when applicable) will be technically reviewed prior to CODIS entry and subsequent NDIS upload.
- 12.2.2 The individual performing the technical review will be a second scientist who is "qualified" in the area of STR Analysis.
- 12.2.3 The second scientist performing the review will initial each page of the data package (and date the first and last page at a minimum).
- 12.2.4 The scientist performing the review of outsourced data (when applicable) will document in an appropriate manner, the review of data for compliance with contract technical specifications and that the .cmf file, if present, contains the correct DNA profiles.
- 12.2.5 Additionally, a documented administrative review will be performed on CODIS hit confirmation letters containing an offender's personally identifiable information, prior to release.

12.3 TESTIMONY REVIEW

Review of courtroom testimony of Forensic Biology personnel shall be accomplished at least once in each calendar year. Preferably, this review will be performed by the Biology/DNA Supervisor or another qualified analyst and documented on the Forensic Services courtroom testimony evaluation form. Alternatively, the evaluation may be completed by criminal justice personnel (e.g., the judge, prosecutor or defense counsel).

13.0 SAFETY

Laboratory safety practices are addressed in the ISP Forensic Services Health and Safety Manual. In Forensic Biology, personnel are introduced to these practices in Module 1 of the ISP Forensic Biology Training Manual. In addition, Section 8 of this manual addresses the monitoring of the chemical eye-wash and shower.

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

14.0 AUDITS

Quality audits and retention schedules are delineated in the ISP Forensic Services Quality/Procedure Manual. Specific Biology/DNA audit requirements are delineated below.

14.1 A DNA audit, using the current FBI DNA Quality Assurance Standards Audit Document(s), will be conducted on an annual basis.

14.2 The interval between annual audits will be in accordance with the current FBI Quality Assurance Standards.

14.3 Every other year, at a minimum, the DNA audit must be an external audit.

14.4 The completed audit document(s) (Quality Assurance Standards Audit for Forensic DNA Testing Laboratories and for DNA Databasing Laboratories) and appropriate accompanying documentation will be submitted to NDIS according to NDIS Operational Procedures.

Property of Idaho State Police Forensic Services
Uncontrolled Internal Copy
OBSOLETE DOCUMENT

15.0 OUTSOURCING

Outsourcing/Subcontracting policies and procedures are described in the ISP Forensic Services Quality/Procedure Manual.

- 15.1 Approved vendor laboratories must provide documentation of accreditation and compliance with the Quality Assurance Standards for Forensic DNA and/or Database Testing Laboratories prior to contract award and for the duration of the contract.
- 15.2 Technical specifications will be outlined in the outsourcing agreement/contract and approved (approval will be documented) by the Biology/DNA Technical Manager prior to award.
- 15.3 An on-site visit of the vendor laboratory will be performed, by the technical leader or a qualified DNA analyst, and documented prior to the submission of any samples to that laboratory.
- 15.4 An annual on-site visit will be performed and documented for any contract extending beyond one year.
- 15.5 When outsourcing convicted offender samples, at least one quality control sample shall be included with each batch. Additionally, at least 5% of the total outsourced samples shall be re-tested and compared for consistency and data integrity.

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

16.0 Practices, Methods and Forms

The following is a list of general practices/administrative procedures, analytical methods and forms utilized in Forensic Biology. Each follows the numbering scheme of: Biology Screening (1XX), DNA Casework Analysis (2XX), CODIS/Database Analysis (3XX) and QC Functions (4XX).

MBI≡Schemes, generally encompassing many procedures.

- MBI-100** EXAMINATION OF BLOODSTAINED EVIDENCE
- MBI-102** EXAMINATION OF EVIDENCE FOR SEMEN
- MBI-104** EXAMINATION OF EVIDENCE FOR BODY FLUIDS
- MBI-200** INDIVIDUALIZATION OF DNA SOURCES BY STR ANALYSIS
- MBI-300** INDIVIDUALIZATION OF DNA SOURCES BY STR ANALYSIS

BI≡Analytical Procedures or Individual Processes

- BI-100** PROCESSING LIQUID BLOOD
- BI-102** DNA PACKETS
- BI-104** PHENOLPHTHALEIN TEST FOR BLOOD
- BI-105** O-TOLIDINE TEST FOR BLOOD
- BI-106** HUMAN BLOOD IDENTIFICATION USING ABACARD® HEMATRACE® TEST
- BI-108** SPECIES IDENTIFICATION: OUCHTERLONY DOUBLE DIFFUSION
- BI-110** BIOLOGICAL SCREENING: USE OF ALTERNATE LIGHT SOURCE
- BI-111** BIOLOGICAL SCREENING: USE OF INFRA RED LIGHT
- BI-114** BRENTAMINE TEST FOR ACID PHOSPHATASE
- BI-116** SAMPLE EXTRACTION FOR SEMEN IDENTIFICATION
- BI-118** SEMEN IDENTIFICATION: MICROSCOPIC EXAMINATION
- BI-119** SPERM DOCUMENTATION: MIDEO SYSTEM
- BI-120** IDENTIFICATION OF SEMEN BY P30 DETECTION (ABAcad®)
- BI-122** AMYLASE TEST: PHADEBAS
- BI-124** AMYLASE TEST: STARCH IODIDE
- BI-126** DETECTION OF URINE (UREASE)
- BI-128** DETECTION OF URINE (CREATININE)
- BI-130** DETECTION OF FECAL MATERIAL (UROBILINOGEN)
- BI-200** EXTRACTION PROTOCOLS FOR PCR DNA TYPING TESTS
- BI-207** DNA QUANTIFICATION: REAL-TIME PCR
- BI-208** STR AMPLIFICATION: PP16
- BI-210** STR TYPING: CAPILLARY ELECTROPHORESIS AND DATA ANALYSIS
- BI-301** CODIS SAMPLE RECEIPT AND DNA TRACKER ENTRY
- BI-302** CODIS SAMPLE DATA ENTRY AND UPLOAD
- BI-303** CODIS DATABASE HIT VERIFICATION
- BI-310** CODIS SAMPLE REMOVAL

- BI-312** EXTRACTION PROTOCOLS FOR PCR DNA TYPING TESTS
- BI-314** DNA QUANTIFICATION : REAL-TIME PCR
- BI-316** STR AMPLIFICATION : PP16HS
- BI-318** STR TYPING : CAPILLARY ELECTROPHORESIS AND DATA ANALYSIS
- BI-400** DRIFTCON FFC : TEMPERATURE VERIFICATION

Form BI≡Various forms used in each discipline

* indicates a controlled form

- 100-BI** PHENOLPHTHALEIN REAGENT (KASTLE-MEYER)
- 102-BI** HYDROGEN PEROXIDE 3% (v/v)
- 103-BI** O-TOLIDINE REAGENT
- 104-BI** AMMONIUM HYDROXIDE (~3%)
- 108-BI** OUCHTERLONY DESTAIN
- 110-BI** OUCHTERLONY STAIN
- 114-BI** 10X BRENTAMINE (SODIUM ACETATE) BUFFER
- 116-BI** BRENTAMINE SOLUTION A
- 118-BI** BRENTAMINE SOLUTION B
- 120-BI** SALINE (0.85% NaCl)
- 124 BI** 1X PHOSPHATE BUFFERED SALINE (PBS)
- 126-BI** XMAS TREE STAIN SOLUTION A (KERNECHTROT SOLUTION)
- 128-BI** XMAS TREE STAIN SOLUTION B (PICROINDIGOCARMINE SOLUTION)
- 132-BI** AMYLASE DIFFUSION BUFFER (pH6.9)
- 134-BI** AMYLASE IODINE REAGENT
- 138-BI** MERCURIC CHLORIDE 10% (w/v)
- 140-BI** ZINC CHLORIDE 10% (w/v)
- 201-BI** 1M TRIS-HCl BUFFER pH7.5
- 203-BI** 1M TRIS-HCl BUFFER pH8
- 205-BI** ETHYLENEDIAMINE TETRAACETIC ACID (EDTA) 0.5M
- 207-BI** STAIN EXTRACTION BUFFER pH8
- 211-BI** PROTEINASE K (20 mg/ml)
- 222-BI** 1M SODIUM ACETATE pH5.2
- 223-BI** DTT (1M)
- 229-BI** PCR-TE (TE⁻⁴) BUFFER (10mM TRIS-HCl, 0.1M EDTA)
- 231-BI** NaOH 5N
- 233-BI** SODIUM CHLORIDE (NaCl) 5M
- 249-BI** BOVINE SERUM ALBUMIN (BSA) 4%
- 101-BI** BIOLOGY SCREENING SUMMARY
- 200-BI** DNA EXTRACTION WORKSHEET
- 202-BI** DIFFERENTIAL DNA EXTRACTION WORKSHEET
- 206-BI*** 7500 LOAD SHEET
- 209-BI*** 7500 RESULTS SHEET
- 210-BI** STR AMPLIFICATION SET-UP
- 212-BI** STR EXTRACTION CONTROL GENOTYPE CHECK
- 214-BI** STR TECHNICAL REVIEW CHECKLIST
- 216-BI*** 3130 LOAD SHEET

218-BI CODIS ENTRY FORM
306-BI STR OFFENDER DATABASE REVIEW CHECKLIST
310-BI CODIS SAMPLE REMOVAL CHECKLIST
312-BI* DATABASE WORKSHEETS (A-E)
314-BI OUTSOURCED OFFENDER DATA REVIEW
400-QC FORENSIC BIOLOGY CHEMICAL INVENTORY
401-QC FORENSIC BIOLOGY CRITICAL EQUIPMENT INVENTORY
402-QC FORENSIC BIOLOGY EQUIPMENT MAINTENANCE/REPAIR RECORD
403-QC* FORENSIC BIOLOGY pH CALIBRATION RECORD
404A-QC* FORENSIC BIOLOGY WEEKLY QC
404B-QC* EVIDENCE VAULT WEEKLY QC
406A-QC* FORENSIC BIOLOGY MONTHLY QC
406B-QC* FORENSIC BIOLOGY MONTHLY QC
408-QC FORENSIC BIOLOGY QUARTERLY QC
410-QC* QC ABACARD® HEMATRACE® KIT
412-QC* QC ONESTEP ABACARD® P30 KIT
419-QC* QC QUANTIFILER® HUMAN DNA QUANTIFICATION KIT
420-QC* QC STR KITS
422A-QC 310 INJECTION LOG
422B-QC 3130/3130x1 INJECTION LOG
426-QC* ANNUAL NIST QC RUN
428-QC BIOMEK 3000 QC

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT



DNA Quantitation

7500 Load Sheet

Form 206-BI

Case Number: _____

Analyst: _____

Plate Name: _____

Date: _____

	1	2	3	4	5	6	7	8	9	10	11	12
A	STD. 1	STD. 1										
B	STD. 2	STD. 2										
C	STD. 3	STD. 3										
D	STD. 4	STD. 4										
E	STD. 5	STD. 5										
F	STD. 6	STD. 6										
G	STD. 7	STD. 7										
H	STD. 8	STD. 8										

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

lot #: _____
 expiry date: _____
 total samples: _____
 Quantifier Kit _____
 Master Mix made for: _____
 reaction mix 0 ul
 primer mix 0 ul
 Std. Prep. Date: _____
 TE lot#: _____



DNA Quantitation

7500 Results Sheet

Form 209-BI

Case Number: _____

Analyst: _____

Plate Name: _____

Date: _____

Well	Sample Name	IPC C _T	Quantity ng/ul	ul Sample for Dilution	ul TE to be added	ng/ul Final	ul to be Amplified
A3		0	0	5	0.0	0.1	10.0
B3		0	0	5	0.0	0.1	10.0
C3		0	0	5	0.0	0.1	10.0
D3		0	0	5	0.0	0.1	10.0
E3		0	0	5	0.0	0.1	10.0
F3		0	0	5	0.0	0.1	10.0
G3		0	0	5	0.0	0.1	10.0
H3		0	0	5	0.0	0.1	10.0
A4		0	0	5	0.0	0.1	10.0
B4		0	0	5	0.0	0.1	10.0
C4		0	0	5	0.0	0.1	10.0
D4		0	0	5	0.0	0.1	10.0
E4		0	0	5	0.0	0.1	10.0
F4		0	0	5	0.0	0.1	10.0
G4		0	0	5	0.0	0.1	10.0
H4		0	0	5	0.0	0.1	10.0
A5		0	0	5	0.0	0.1	10.0
B5		0	0	5	0.0	0.1	10.0
C5		0	0	5	0.0	0.1	10.0
D5		0	0	5	0.0	0.1	10.0
E5		0	0	5	0.0	0.1	10.0
F5		0	0	5	0.0	0.1	10.0
G5		0	0	5	0.0	0.1	10.0
H5		0	0	5	0.0	0.1	10.0
A6		0	0	5	0.0	0.1	10.0
B6		0	0	5	0.0	0.1	10.0
C6		0	0	5	0.0	0.1	10.0
D6		0	0	5	0.0	0.1	10.0
E6		0	0	5	0.0	0.1	10.0
F6		0	0	5	0.0	0.1	10.0
G6		0	0	5	0.0	0.1	10.0
H6		0	0	5	0.0	0.1	10.0
A7		0	0	5	0.0	0.1	10.0
B7		0	0	5	0.0	0.1	10.0
C7		0	0	5	0.0	0.1	10.0
D7		0	0	5	0.0	0.1	10.0
E7		0	0	5	0.0	0.1	10.0
F7		0	0	5	0.0	0.1	10.0
G7		0	0	5	0.0	0.1	10.0
H7		0	0	5	0.0	0.1	10.0
A8		0	0	5	0.0	0.1	10.0
B8		0	0	5	0.0	0.1	10.0
C8		0	0	5	0.0	0.1	10.0
D8		0	0	5	0.0	0.1	10.0
E8		0	0	5	0.0	0.1	10.0
F8		0	0	5	0.0	0.1	10.0
G8		0	0	5	0.0	0.1	10.0

Well	Sample Name	IPC C _T	Quantity ng/ul	ul Sample for Dilution	ul TE to be added	ng/ul Final	ul to be Amplified
H8		0	0	5	0.0	0.1	10.0
A9		0	0	5	0.0	0.1	10.0
B9		0	0	5	0.0	0.1	10.0
C9		0	0	5	0.0	0.1	10.0
D9		0	0	5	0.0	0.1	10.0
E9		0	0	5	0.0	0.1	10.0
F9		0	0	5	0.0	0.1	10.0
G9		0	0	5	0.0	0.1	10.0
H9		0	0	5	0.0	0.1	10.0
A10		0	0	5	0.0	0.1	10.0
B10		0	0	5	0.0	0.1	10.0
C10		0	0	5	0.0	0.1	10.0
D10		0	0	5	0.0	0.1	10.0
E10		0	0	5	0.0	0.1	10.0
F10		0	0	5	0.0	0.1	10.0
G10		0	0	5	0.0	0.1	10.0
H10		0	0	5	0.0	0.1	10.0
A11		0	0	5	0.0	0.1	10.0
B11		0	0	5	0.0	0.1	10.0
C11		0	0	5	0.0	0.1	10.0
D11		0	0	5	0.0	0.1	10.0
E11		0	0	5	0.0	0.1	10.0
F11		0	0	5	0.0	0.1	10.0
G11		0	0	5	0.0	0.1	10.0
H11		0	0	5	0.0	0.1	10.0
A12		0	0	5	0.0	0.1	10.0
B12		0	0	5	0.0	0.1	10.0
C12		0	0	5	0.0	0.1	10.0
D12		0	0	5	0.0	0.1	10.0
E12		0	0	5	0.0	0.1	10.0
F12		0	0	5	0.0	0.1	10.0
G12		0	0	5	0.0	0.1	10.0
H12		0	0	5	0.0	0.1	10.0

Property of Idaho State Police Forensic Services
 Uncontrolled Internet Copy
 OBSOLETE DOCUMENT



3130 Load Sheet

Form 216-BI

Case Number: _____

Analyst: _____

Plate Name: _____

Date: _____


	1	2	3	4	5	6	7	8	9	10	11	12
A												
B												
C												
D												
E												
F												
G												
H												

Property of Idaho State Police Forensic Services
 Uncontrolled Internal Copy
 OBSOLETE DOCUMENT

total samples: _____
 Master Mix made for: _____
 HiDi Formamide 0 ul
 Internal Lane Standard 0 ul

3130 POP4

Lot# _____
 HiDi Formamide Lot# _____
 Expiration Date _____
 Buffer Lot# _____



BSD/Biomek 3000 Load Sheet

Form 312A-BI

Plate Name: _____

Analyst: _____

Test Name: _____

Date: _____

	1	2	3	4	5	6	7	8	9	10	11	12
A												
B												
C												
D												
E												
F												
G												
H												

Property of Idaho State Police Forensic Services
 Uncontrolled Internet Copy
 OBSOLETE DOCUMENT

Punch Size: _____

Elution Volume (µl): _____

DTT Lot# / Date Added: _____

Number of Punches: _____

Biomek Program: _____

Ethanol Lot# / Date Added: _____

2-Propanol Lot# / Date Added: _____

DNA IQ Kit

Lot #: _____

Exp. Date _____



DNA Quantitation

7500 Load Sheet

Form 312B-BI

Plate Name: _____ Analyst: 0

7500: _____ Date: _____

	1	2	3	4	5	6	7	8	9	10	11	12
A	0	0	0	0	0	0	0	0	0	0	STD. 1	STD. 1
B	0	0	0	0	0	0	0	0	0	0	STD. 2	STD. 2
C	0	0	0	0	0	0	0	0	0	0	STD. 3	STD. 3
D	0	0	0	0	0	0	0	0	0	0	STD. 4	STD. 4
E	0	0	0	0	0	0	0	0	0	0	STD. 5	STD. 5
F	0	0	0	0	0	0	0	0	0	0	STD. 6	STD. 6
G	0	0	0	0	0	0	0	0	0	0	STD. 7	STD. 7
H	0	0	0	0	0	0	0	0	0	0	STD. 8	STD. 8

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

Quantifier Kit _____ Total Samples: _____ Master Mix made for: _____

Lot #: _____ Std. Prep. Date: _____ reaction mix 0 ul

Exp. Date: _____ TE lot#: _____ primer mix 0 ul



DNA Quantitation

7500 Results Sheet

Form 312C-BI

Plate Name: 0

Analyst: 0

Date: 1/0/1900

Well	Sample Name	IPC C _T	Quantity ng/ul	ul Sample for Dilution	ul TE to be added	ng/ul Final	ul to be Amplified
A1		0	0.00	5	0.0	0.1	2.5
B1		0	0.00	5	0.0	0.1	2.5
C1		0	0.00	5	0.0	0.1	2.5
D1		0	0.00	5	0.0	0.1	2.5
E1		0	0.00	5	0.0	0.1	2.5
F1		0	0.00	5	0.0	0.1	2.5
G1		0	0.00	5	0.0	0.1	2.5
H1		0	0.00	5	0.0	0.1	2.5
A2		0	0.00	5	0.0	0.1	2.5
B2		0	0.00	5	0.0	0.1	2.5
C2		0	0.00	5	0.0	0.1	2.5
D2		0	0.00	5	0.0	0.1	2.5
E2		0	0.00	5	0.0	0.1	2.5
F2		0	0.00	5	0.0	0.1	2.5
G2		0	0.00	5	0.0	0.1	2.5
H2		0	0.00	5	0.0	0.1	2.5
A3		0	0.00	5	0.0	0.1	2.5
B3		0	0.00	5	0.0	0.1	2.5
C3		0	0.00	5	0.0	0.1	2.5
D3		0	0.00	5	0.0	0.1	2.5
E3		0	0.00	5	0.0	0.1	2.5
F3		0	0.00	5	0.0	0.1	2.5
G3		0	0.00	5	0.0	0.1	2.5
H3		0	0.00	5	0.0	0.1	2.5
A4		0	0.00	5	0.0	0.1	2.5
B4		0	0.00	5	0.0	0.1	2.5
C4		0	0.00	5	0.0	0.1	2.5
D4		0	0.00	5	0.0	0.1	2.5
E4		0	0.00	5	0.0	0.1	2.5
F4		0	0.00	5	0.0	0.1	2.5
G4		0	0.00	5	0.0	0.1	2.5
H4		0	0.00	5	0.0	0.1	2.5
A5		0	0.00	5	0.0	0.1	2.5
B5		0	0.00	5	0.0	0.1	2.5
C5		0	0.00	5	0.0	0.1	2.5
D5		0	0.00	5	0.0	0.1	2.5
E5		0	0.00	5	0.0	0.1	2.5
F5		0	0.00	5	0.0	0.1	2.5
G5		0	0.00	5	0.0	0.1	2.5
H5		0	0.00	5	0.0	0.1	2.5
A6		0	0.00	5	0.0	0.1	2.5
B6		0	0.00	5	0.0	0.1	2.5
C6		0	0.00	5	0.0	0.1	2.5
D6		0	0.00	5	0.0	0.1	2.5
E6		0	0.00	5	0.0	0.1	2.5
F6		0	0.00	5	0.0	0.1	2.5
G6		0	0.00	5	0.0	0.1	2.5
H6		0	0.00	5	0.0	0.1	2.5



DNA Quantitation

7500 Results Sheet

Form 312C-BI

Plate Name: 0

Analyst: 0

Date: 1/0/1900

Well	Sample Name	IPC C _T	Quantity ng/ul	ul Sample for Dilution	ul TE to be added	ng/ul Final	ul to be Amplified
A7		0	0.00	5	0.0	0.1	2.5
B7		0	0.00	5	0.0	0.1	2.5
C7		0	0.00	5	0.0	0.1	2.5
D7		0	0.00	5	0.0	0.1	2.5
E7		0	0.00	5	0.0	0.1	2.5
F7		0	0.00	5	0.0	0.1	2.5
G7		0	0.00	5	0.0	0.1	2.5
H7		0	0.00	5	0.0	0.1	2.5
A8		0	0.00	5	0.0	0.1	2.5
B8		0	0.00	5	0.0	0.1	2.5
C8		0	0.00	5	0.0	0.1	2.5
D8		0	0.00	5	0.0	0.1	2.5
E8		0	0.00	5	0.0	0.1	2.5
F8		0	0.00	5	0.0	0.1	2.5
G8		0	0.00	5	0.0	0.1	2.5
H8		0	0.00	5	0.0	0.1	2.5
A9		0	0.00	5	0.0	0.1	2.5
B9		0	0.00	5	0.0	0.1	2.5
C9		0	0.00	5	0.0	0.1	2.5
D9		0	0.00	5	0.0	0.1	2.5
E9		0	0.00	5	0.0	0.1	2.5
F9		0	0.00	5	0.0	0.1	2.5
G9		0	0.00	5	0.0	0.1	2.5
H9		0	0.00	5	0.0	0.1	2.5
A10		0	0.00	5	0.0	0.1	2.5
B10		0	0.00	5	0.0	0.1	2.5
C10		0	0.00	5	0.0	0.1	2.5
D10		0	0.00	5	0.0	0.1	2.5
E10		0	0.00	5	0.0	0.1	2.5
F10		0	0.00	5	0.0	0.1	2.5
G10		0	0.00	5	0.0	0.1	2.5
H10		0	0.00	5	0.0	0.1	2.5



STR Amplification Set-Up

Form 312D-BI

Analyst: 0

Plate Name: _____

Date: _____

Thermal Cycler: _____

	1	2	3	4	5	6	7	8	9	10	11	12
A	0	0	0	0	0	0	0	0	0	0	0	0
B	0	0	0	0	0	0	0	0	0	0	0	0
C	0	0	0	0	0	0	0	0	0	0	0	0
D	0	0	0	0	0	0	0	0	0	0	0	0
E	0	0	0	0	0	0	0	0	0	0	0	0
F	0	0	0	0	0	0	0	0	0	0	0	0
G	0	0	0	0	0	0	0	0	0	0	0	0
H	0	0	0	0	0	0	0	0	0	0	0	0

STR Kit Type: _____

STR Kit Lot#: _____

PCR TE Lot#: _____

Total Samples: _____

Master Mix Made For: _____

Reagent	μl/sample	μl in Master Mix
5x Master Mix	1.25	0
Primers	0.6	0
H ₂ O		0
Master Mix/Sample	1.85	
DNA Template	4.4	
Total Rxn Volume	6.25	



3130XL Load Sheet

Form 312E-BI

Plate Name: _____

Analyst: _____ 0 _____


Date: _____

	1	2	3	4	5	6	7	8	9	10	11	12
A	0	0	0	0	0	0	0	0	0	0	0	0
B	0	0	0	0	0	0	0	0	0	0	0	0
C	0	0	0	0	0	0	0	0	0	0	0	0
D	0	0	0	0	0	0	0	0	0	0	0	0
E	0	0	0	0	0	0	0	0	0	0	0	0
F	0	0	0	0	0	0	0	0	0	0	0	0
G	0	0	0	0	0	0	0	0	0	0	0	0
H	0	0	0	0	0	0	0	0	0	0	0	0

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

Total Samples: _____
 3130 POP4
 Lot#: _____
 Expiration Date: _____
 HiDi Formamide Lot#: _____
 Buffer Lot#: _____

Master Mix made for: _____
 HiDi Formamide 0 ul
 Internal Lane Standard 0 ul



FORENSIC BIOLOGY pH CALIBRATION RECORD
(Oakton pH meter, serial #135212)

DATE	INITIALS	STANDARD BUFFER pH 4.01 Reading/lot #	STANDARD BUFFER pH 7.00 Reading/lot #	STANDARD BUFFER pH 10.01 Reading/lot #

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

A 3-point calibration of the pH meter will be performed at the time of use (See the Oakton Operating Manual for calibration and pH measurement instructions). The analyst will record the date of calibration, their initials, the measured pH value and lot # for each buffer. The measured reading must fall within ±0.50 pH for the calibration to be confirmed by the meter.

FORENSIC BIOLOGY WEEKLY QC

DATE/INITIALS	°C		°C		°C		°C		°C	
	Min	Max	Min	Max	Min	Max	Min	Max	Min	Max
COMBO F/F A										
FRIDGE A										
FREEZER A										
FREEZER B										
FREEZER C										
COMBO F/F B										
COMBO F/F C										

Observed temperatures should fall between $\pm 5^{\circ}\text{C}$ of the target temperature (4°C for refrigerators and -20°C for freezers). The temperature control should be adjusted to correct for minor variations; however, if the temperature is not corrected or if it falls significantly outside the target range, it should be taken out of service and maintenance/repair performed as needed. Note: frost-free freezers will have a greater temperature range ($\pm 10^{\circ}\text{C}$) due to the heating and cooling cycles. Combination fridge/freezers with a single temperature control may not be able to maintain both units within the given range. In this case, temperature sensitive reagents should be stored in appropriate containers (such as cryo-boxes) to maintain the desired state.

DATE/INITIALS					
NANOPIRE SYSTEM					
°C OVEN	set/observed				
°C HEATING BLOCK A	set/observed				
°C HEATING BLOCK B	set/observed				
°C HEATING BLOCK C (prod. rm.)	set/observed				
°C HEATING BLOCK D	set/observed				

Observed temperatures should fall between $\pm 2^{\circ}\text{C}$ of the temperature set point. The temperature control should be adjusted to correct for minor variations; however, if the temperature is not corrected or if it falls significantly outside the target range, it should be taken out of service and maintenance/repair performed as needed.

The observed water purity for the Nanopure system should be a minimum of 18.0 mega-ohms. If the purity falls below this point, the cartridges should be changed and the system sanitized as necessary.

MERIDIAN EVIDENCE VAULT WEEKLY QC

DATE/INITIALS	°C		°C		°C		°C		°C	
	Min	Max	Min	Max	Min	Max	Min	Max	Min	Max
FRIDGE 1										
FRIDGE 2										
FREEZER 1										
FREEZER 2										
FREEZER 3										
FREEZER 4										
FREEZER 5										
FREEZER 6										
FREEZER 7										

Observed temperatures for refrigerators should fall between $\pm 5^{\circ}\text{C}$ of the 4°C target temperature. Freezers should fall between $\pm 10^{\circ}\text{C}$ of the target -20°C . The temperature control should be adjusted to correct for minor variations; however, if the temperature is not corrected or if it falls significantly outside the target range, it should be taken out of service and maintenance/repair performed as needed. Note: frost-free freezers will have a greater temperature range ($\pm 15^{\circ}\text{C}$) due to the heating and cooling cycles.

FORENSIC BIOLOGY MONTHLY QC

Form 406A-QC

AUTOCLAVE											
DATE/INITIALS											
CLEAN											
STERILIZATION (+)											
STERILIZATION (-)											
LABORATORY AND OTHER EQUIPMENT											
BIROBOT EZ1s											
GREASE D-RINGS CLEAN											
CENTRIFUGES											
CLEAN PIPETS											
CLEAN BSD PUNCHER											
CLEAN BIOMEK 3000											
LAB CLEANED											
EYEWASH CHECK											

*Personnel should initial the duties they perform and date separately if necessary.

Autoclave sterilization is checked by the observation of microbial growth in the (+) control (non-sterilized) and a lack of growth in the (-) control (sterilized) samples. See the BTSure product insert for test instructions and growth indicators. If sterilization is not achieved, the autoclave should be serviced.

FORENSIC BIOLOGY MONTHLY QC

Form 406B-QC

DATE/INITIALS																				
ABI 7500 Instrument Maintenance																				
Background Assay/Contamination Check																				
System Function Test																				
Lamp Status Check																				
7500 & 310 Computer Maintenance																				
Disk Cleanup																				
Defragment Hard Drive																				
3130/3130XL Maintenance																				
Water Seal Trap Flush																				
Water Wash Wizard																				
Defragment Hard Drive (C & E)																				

*Personnel should initial the duties they perform and date separately, if necessary.

Note: See the ABI 7500 Maintenance Guide and/or the April 2007 User Bulletin for additional Instrument Maintenance procedures and pass/fail criteria.

A contamination check may be performed by either the background assay, or block check, or a combination of the two. If outliers are observed during the Background Assay (Intensity Value $\geq 72,000$), or fluorescence (red) observed during the block check, the specific well should be identified and cleaned. Rerun the background calibration after wells have been cleaned sufficiently. Note: a 96-well tray with 50ul TE in each well may be used as a background tray.

A block check is performed by selecting Instrument > Calibrate. Set the exposure time to 2048ms, lamp control to idle, select Filter A and click Snapshot. Holding the cursor over the fluorescence will give pixel intensity.

If a component fails the function test, a service call should be placed.

If the lamp fails the function test and/or status check, it should be replaced, followed by calibration of ROI, background, optical, pure dye, and instrument verification in that order. Reset the Lamp Timer when complete.

The 7500 and 310 Disk Cleanup is performed by selecting Start Menu > Programs > Accessories > System Tools > Disk Cleanup.

7500 and 310 Defragmentation is performed by selecting Start Menu > Programs > Accessories > System Tools > Disk Defragmenter.

3130/3130XL Defragmentation is performed by right-clicking on 'My Computer' and selecting 'Manage'. In the tree tab choose Computer Management (local)>Disk Fragmenter>Drive name>Defragment.

QC ABACARD® HEMATRACE® KIT

HEMATRACE® KIT LOT: _____

DATE RECEIVED: _____

SCIENTIST: _____

QC DATE: _____

Perform test as usual with one 2mm² cutting and one 2mm thread from known bloodstain. Record results (include time it took for positive rxn to be visible). If available, attach photo documentation and place in Forensic Biology QC binder.

SAMPLE	RXN	TIME (min. sec.)
2mm ² cutting		
2mm thread		
Neg		

The 2mm² cutting sample must have a positive reaction within 10 minutes for passing. The 2mm thread should ideally be positive within 10 minutes but is used primarily as a sensitivity indicator of the given test lot. The kit may still be deemed as passing without a positive result for the thread.

QA/QC PASSED: YES NO

Comments:

QC OneStep ABACARD® p30 KIT

ABACARD® p30 KIT LOT: _____ DATE RECEIVED: _____

SCIENTIST: _____ QC DATE: _____

Perform test as usual with a known semen extract, as well as ~10ng/ml (10µl of a 1:500 dilution) and ~50ng/ml (10µl of a 1:100 dilution) of Seri Semen Standard. Record results (include time it took for positive rxn to be visible). If available, attach photo documentation and place in Forensic Biology QC binder.

SAMPLE	RXN	TIME (min. sec.)
Semen Extract		
10ng/ml		
50ng/ml		
Neg		
*250ng/ml or 1:10		

The semen extract must have a positive rxn within 10 minutes for passing. The Seri standards are used to estimate the range of sensitivity of the kit lot.

For the semen standard dilutions, if a positive rxn is not obtained at 10 minutes, continue to monitor and record result at the end of 15 minutes. In addition, *run a 250ng/ml (50µl of the 1:100 dilution to 150µl of extraction buffer) or a 1:10 dilution of the semen stain extract to ensure the kit is operating within reasonable limits for forensic identification. In addition to the neat semen extract, this control sample (250ng/ml or 1:10 extract) must result in a positive rxn within 10 minutes.

QA/QC PASSED: YES NO

Comments:



QC QUANTIFILER HUMAN KITS

Form 419-QC

KIT LOT #: _____

DATE RECEIVED: _____

EXPIRATION DATE: _____

SCIENTIST: _____

QA/QC DATE: _____

KIT COMPONENT	LOT NUMBER
PRIMER MIX	
REACTION MIX	
DNA STANDARD	

To check the new kit lot, perform quantification as usual. For samples, run standards from the new kit to be QC'd and equivalent dilutions of the NIST SRM 2372 Quant Standard, as well as 0.5ng and 10ng of 9947A DNA. Analyze using the SRM as standard and the new kit as unknown. Using an average of the results for the new kit standards, calculate the new volume of TE to be added in the preparation of standard 1, per the equation $C_1V_1=C_2V_2$ (where C=average for std 1, and V=total volume). Record the slope obtained for the standard curve.

As a check of the calculation and resulting TE volume, use the new kit, with corresponding new dilution to perform a 9947A DNA quantification, according to standard procedure. Use 0.5ng and 10ng and compare the results to those obtained from above. A pass will be achieved if the slopes for both standard curves are comparable.

SRM 2372 component used: _____


Standard curve slope: _____

Volume TE to be used for Standard 1: _____

QA/QC PASSED: YES NO

Comments:

Attach the 7500 Load Sheets, Standard Curves, and Results Sheets. Record the calculations in the documentation. Mark the new kit with TE volume for Standard 1 preparation.



QC STR KITS

STR KIT: _____ DATE RECEIVED: _____

KIT MANUFACTURER: _____ KIT LOT #: _____

LAB LOT#: _____ SCIENTIST: _____ QA/QC DATE: _____

PP16 KIT COMPONENT	LOT NUMBER	PP16HS KT COMPONENT	LOT NUMBER
PRIMER MIX		10X PRIMER PAIR MIX	
REACTION MIX		5X MASTER MIX	
CONTROL DNA		CONTROL DNA	
TAQ GOLD*		INTERNAL LANE STANDARD	
INTERNAL LANE STANDARD		ALLELIC LADDER	
ALLELIC LADDER			

An Extraction Control, reagent blank and associated controls are to be processed from extraction. A pass will be achieved by obtaining the expected results for each of the samples run and data of acceptable quality (e.g. sufficient RFUs). Comments regarding quality concerns are to be noted as appropriate.

* Taq Gold is purchased separately for PowerPlex 16, but typically at the same time as a new STR kit. If Taq Gold is received separate from an STR kit, the QC of the Taq will be noted on this form (corresponding to the appropriate STR kit lot#) under comments.

The ILS, once QC'd, may be used interchangeably between the PP16 and PP16HS kits, if necessary, as they are the same product. The ladders are not the same product and cannot be switched between kit types.

Run Date: _____ Run Folder: _____

QA/QC PASSED: YES NO

Comments:

Attach the appropriate extraction/amplification/Extraction Control forms used and the GeneMapper ID Electropherograms; place in the appropriate QC Binder.

ANNUAL NIST QC RUN

SCIENTIST: _____

QC DATE: _____

At a minimum of once a year, an 'in-date', certified NIST-SRM standard will be analyzed with our standard procedures. Control or known reference samples may be analyzed simultaneously to 'certify' them for use as NIST QC samples. These samples will be listed in the comments section of this form with lot # and that they were certified. After completion of the QC, the newly 'certified' samples, or their container, will be marked as "NIST Certified" with the corresponding date.

The GeneMapper® ID Data will be analyzed as usual and quality of results will be reflected in the comments or 'passed' areas as appropriate and necessary. Passing results are obtained by achieving the expected results for the given NIST sample(s) and any associated controls. The GeneMapper® ID Electropherograms and an Allele Table will be printed [for the NIST sample(s)] and stored in the Forensic Biology QC binder.

Run Folder: _____

QC PASSED: YES NO

Comments:

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT