

Idaho State Police  
Forensic Services  
Toxicology Section



---

Section Four  
Analysis of Alcohol and Common Volatile Solvents

---

---

**4.2 Analysis of Solutions Containing Ethanol and Common Volatiles**

---

**4.2.1 BACKGROUND**

The need to establish the ethyl alcohol concentration and/or the presence of other commonly encountered volatiles in a beverage or solution may arise from ABC violations (Idaho Code 23-611, 23-1002, 23-1303, ...), under-age consumption (Idaho Code 23-603, 23-604), open-container violations (Idaho Code 23-505, 23-1333), poisonings and/or an endless variety of situations including questionable samples submitted as blood or other physiological fluid. In addition, ethyl concentration must be verified in simulator solutions used for breath testing instruments (IDAPA 11.03.01).

**4.2.2 SCOPE**


This method describes the analysis of alcoholic beverages and solutions said to contain a specified amount of ethyl alcohol and of unknown solutions via a headspace sampling gas chromatographic method. Unknown solutions may also contain other volatiles such as acetone, methanol, isopropanol and toluene, which can be qualitatively identified with this method. Samples, controls and standards are sealed into vials that contain an aqueous n-propanol internal standard solution and heated by the headspace analyzer. As described in Henry's Law, in a closed container at a given temperature, a direct (proportional) relationship exists between the amount of a volatile substance dissolved in a liquid and the amount of the volatile substance in the headspace vapor above the solution. An aliquot of the vapor is injected into a gas chromatograph (GC) in a dual column configuration. The GC serves to separate out the components of the solution as a function of their chemical properties. The separated components are identified qualitatively on the basis of the retention time determined for each of the columns. Quantitation of ethanol is accomplished through area percent data obtained from a flame ionization detector (FID). The quantitative result is based on a minimum of a three-point calibration curve, which uses the peak area ratio between ethanol and the n-propanol internal standard. These solution samples can be included as part of a toxicology alcohol determination run utilizing SOP 4.1 provided that quality assurance requirements are met. In addition, if this method is applied for specifically for the qualitative identification of volatiles other than ethanol, ethanol calibrators and controls need not be run.

### 4.2.3 EQUIPMENT

- 4.2.3.1 Perkin Elmer Auto System XL Gas Chromatograph (GC)
- 4.2.3.2 Columns
- 4.2.3.2.1 Restek Rtx<sup>®</sup>-BAC1 (#18003: 30 meter X 0.32mm inner diameter (ID), 1.8µm film thickness (FT)) or equivalent column
- 4.2.3.2.2 Restek Rtx<sup>®</sup>-BAC2 (#18002: 30 meter X 0.32mm ID, 1.2 µm film thickness (FT)) or equivalent column
- 4.2.3.3 Perkin Elmer HS-40 or HS-110 Headspace Autosampler (figures 2 and 3)



*Figure 2. HS-40*



*Figure 3. HS-110*

- 4.2.3.4 PE Workstation Software, TotalChrom Version 6.2.0 or more recent version/upgrade.
- 4.2.3.5 Hand Crimper (P-E B003-8134 or equivalent)
- 4.2.3.6 Hamilton MICROLAB 503A or equivalent semi-automatic Dilutor/Pipetter equipped with sample and reagent syringes capable of dispensing 250µL and 2000µL, respectively.
- 4.2.3.7 Glassware
- 4.2.3.7.1 GC-Headspace vials (P-E B010-4236 or equivalent)
- 4.2.3.7.2 Safety Closures {PTBE septa, crimp caps and star springs} (P-E BO10-4240 or equivalent)

### 4.2.4 SAFETY CONCERNS

- 4.2.4.1 Samples should be processed according to safety guidelines in the *Chemical Hygiene and Safety Manual*.

### 4.2.5 REAGENTS

- 4.2.5.1 1-Propanol (Acros/Fisher Scientific # 23207-0010, #A996-1 or equivalent)
- 4.2.5.2 Acetone (Fisher #A929-1 or equivalent)

- 4.2.5.3 Acetaldehyde (Fisher #01004-250 or equivalent)
- 4.2.5.4 Isopropanol (2-Propanol) (Fisher #A416-500 or equivalent)
- 4.2.5.5 Methanol (Fisher #A454-1 or equivalent)
- 4.2.5.6 Toluene (Fisher T324-500 or equivalent)
- 4.2.5.7 Ammonium Sulfate (Fisher #A702-500 or equivalent)
- 4.2.5.8 Sodium Fluoride (Fisher #S299-500 or equivalent)

#### 4.2.6 REFERENCE MATERIAL

Record the preparation of all solutions on reagent log.

##### 4.2.6.1 Ethanol Calibration Standards

Aqueous Ethanol Standards (g/100mL)

As available, 1-mL ampules containing solutions at concentrations of  $\cong 0.025$ , 0.05, 0.08, 0.10, 0.20, 0.30, and 0.40 (Cerilliant or equivalent)

##### 4.2.6.2 Ethanol Control Standard

Aqueous Ethanol Standards (g/100mL)

As available, 1-mL ampules containing concentrations of 0.02 to 0.40 (Restek or equivalent source not used to prepare calibration standards)

##### 4.2.6.3 Aqueous Multicomponent Mixture ( $\cong$ g/100cc)

Multicomponent mixture can be obtained commercially and/or prepared from reagents as described below.

###### 4.2.6.3.1 **Multicomponent Kit**

(Cerilliant #A-054 or comparable).

Cerilliant kit includes an aqueous 0.05, 0.10 or 0.40 ethanol for use as a quantitative ethanol control and acetone, methanol and isopropanol standards which this method utilizes qualitatively.

###### 4.2.6.3.2 **Qualitative Volatile Standard Mix Solution**

Add approximately 200mL of DI water to a 250mL volumetric flask.

Add one or more of the following volatiles, as needed:

100 $\mu$ L acetaldehyde

100 $\mu$ L acetone

500 $\mu$ L methanol

500 $\mu$ L isopropanol

500  $\mu$ L ethanol

100 $\mu$ L toluene.

QS to 250mL.

*Solution is stable for 1 year.*

As the need arises, other volatiles can be included in this mixture or in a single constituent solution. The volatile should be analyzed prior to the analysis run to determine chromatography characteristics including the retention time. Retention time should be programmed into analysis methods.

#### 4.2.6.4 Internal Standard Solution

##### **0.03g/dL 1-propanol in 1.0M (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub>**

Dissolve 132.14g (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub> in approximately 800mL of to a 1L volumetric flask.

Add 1g sodium fluoride {optional}.

Add 375μL 1-propanol.

QS to 1L with distilled water.

*Solution is stable for 2 months.*

## 4.2.7 ANALYSIS PROCEDURE

### 4.2.7.1 General

4.2.7.1.1 Solutions covered under this SOP can be included as part of a routine toxicology alcohol analysis run.

4.2.7.1.2 Bring calibrators, controls, internal standard and samples to room temperature.

4.2.7.1.3 If analysis is for volatiles other than ethanol, ethanol calibrators need not be included in the analysis run.

4.2.7.1.4 Gather necessary vials, closures and ancillary supplies in or near laminar flow hood.

4.2.7.1.5 Sample preparation should take place in a laminar flow hood.

### 4.2.7.2 Pipetter/Dilutor Set-up

4.2.7.2.1 Switch on power.

4.2.7.2.2 Display will inquire as to the sizes of installed syringes. Select the correct size for sample syringe [right] and reagent syringe [left].

4.2.7.2.3 Scroll down to volume option. Select 250μL for sample syringe [right] and 2000μL for reagent syringe [left].

4.2.7.2.4 Scroll down to speed option. Verify that syringe speed is on desired setting.

4.2.7.2.5 Prime the fluid path. Continue priming until no bubbles are observed.

### 4.2.7.3 Preparation of Blanks and Mixed Standard

#### 4.2.7.3.1 Water Blank

- 4.2.7.3.1.1 Label test vial with *water blank*.
- 4.2.7.3.1.2 Add 2mL DI water to labeled headspace vial.
- 4.2.7.3.1.3 Seal **immediately** with crimp cap.
- 4.2.7.3.2 Internal Standard Blank
- 4.2.7.3.2.1 Label test vial with *ISTD blank*.
- 4.2.7.3.2.2 Use Pipetter/Dilutor to dispense 2000 $\mu$ L of internal standard (ISTD) into labeled headspace vial.
- 4.2.7.3.2.3 Seal **immediately** with crimp cap.
- 4.2.7.3.3 Aqueous Ethanol Controls
- 4.2.7.3.3.1 Label appropriate number of headspace vials for *aqueous ethanol controls (1, 2,...)*.
- 4.2.7.3.3.2 Use Pipetter/Dilutor to dispense 250 $\mu$ L of aqueous control and 2000 $\mu$ L of internal standard (ISTD) into each labeled headspace vial.
- 4.2.7.3.3.3 Seal **immediately** with crimp cap.
- 4.2.7.3.4 Mixed/Other Volatiles Standard Solution
- 4.2.7.3.4.1 Label headspace vial with qualitative volatiles solution.
- 4.2.7.3.4.2 Use Pipetter/Dilutor to dispense 250 $\mu$ L of mixed volatile solution and 2000 $\mu$ L of internal standard (ISTD) into labeled headspace vial.
- 4.2.7.3.4.3 Seal **immediately** with crimp cap.
- 4.2.7.4 Ethanol Calibration Standard Preparation
- 4.2.7.4.1 Label vials for standards.
- 4.2.7.4.2 Use Pipetter/Dilutor to dispense 250 $\mu$ L of appropriate ethanol concentration and 2000 $\mu$ L of internal standard (ISTD) into each labeled headspace vial.
- 4.2.7.4.3 Seal **immediately** with crimp cap.
- 4.2.7.4.4 Establish ethanol calibration plot with a minimum of three calibration points.
- 4.2.7.5 Initial Processing of Specimens
- 4.2.7.5.1 Inspecting and document the condition of seals.
- 4.2.7.5.2 Remove sample(s) container from packaging and place laboratory number on each sample.

- 4.2.7.6 Preparation of Samples for Analysis
- 4.2.7.6.1 Label two headspace vials with the laboratory number.
- 4.2.7.6.2 Dilute alcoholic beverages and unknown solutions as necessary. The sample should be diluted for the value to fall on calibration curve. Generally, beer and wine should be diluted 50:1 with DI water and distilled beverages ( $\geq 16\%$  w/v or 20% v/v) diluted 100:1. If available, the dilution of unknown solutions should be based on sample history.
- 4.2.7.6.3 Breath testing simulator solutions and samples, which appear to be serum, do not require dilution.
- 4.2.7.7 Addition of sample to headspace vials.
- 4.2.7.7.1 Use Pipetter/Dilutor dispense 250 $\mu$ L of sample and 2000 $\mu$ L of internal standard (ISTD) to a labeled headspace vial.
- 4.2.7.7.2 Seal headspace vials **immediately** with crimp caps.
- 4.2.7.8 Preparation for Run
- 4.2.7.8.1 Open **Sequence Editor**
- 4.2.7.8.2 Into Sequence log table, enter the sample case numbers, ethanol standards, other volatiles mix, blanks and controls.
- 4.2.7.8.3 Load samples, calibration standards, blank and controls into the carousel of the headspace sampler as noted in the sequence table.
- 4.2.7.9 Instrument Acquisition Parameters
- 4.2.7.9.1 Refer to instrument METHOD printouts for gas chromatograph and headspace analyzer analysis parameters.
- 4.2.7.9.2 Parameters are at the discretion of the analyst and should be optimized for the instrument.
- 4.2.7.9.3 GC oven parameters should provide for the baseline separation of commonly encountered volatiles in the test mixtures described in sections 4.2.6.3.
- 4.2.7.9.4 Each laboratory should maintain a centrally stored current METHOD printout.
- 4.2.7.10 Calibration
- 4.2.7.10.1 A minimum of three ethanol calibrators must be included in each run. The calibrators chosen should characterize the entire range of interest.

4.2.7.10.2 Ethanol calibrators should be analyzed in order of increasing concentration.

4.2.7.11 Acceptance Criteria

4.2.7.11.1 **Acceptance of Analysis Run**

4.2.7.11.1.1 The least squares line resulting from the analysis of the ethanol calibrators must have a coefficient of correlation of  $\geq 0.998$ .

4.2.7.11.1.2 If calibration standards are run in duplicate, it is not required that duplicate calibration points are included as long as linearity requirement is met.

4.1.7.11.2 **Qualitative Accuracy Criteria**

The qualitative presence of ethanol or other volatile substance can be established if the relative retention time(s) for a specimen is within  $\pm 0.10$  minutes of the relative retention time of a standard of the compound in question. This rejection criterion should be designated in the data station analysis method.

4.1.7.11.3 **Quantitative Accuracy Criteria**

The quantitative ethanol results for a batch of samples can be accepted if the values obtained for aqueous control samples fall within  $\pm 10\%$  of target value on Certificate of Analysis.

4.2.7.11.4 **Quantitative Precision Criteria**

The results obtained from duplicate analysis must agree within 0.010g/100mL. For breath testing solutions, the results between different bottles of solution must also agree within 0.010g/100mL. If these precision requirements are not met, the sample(s) must be reanalyzed.

4.2.7.11.5 **High Ethanol Values**

When an elevated ethanol value is obtained, appropriate calibrators must bracket the value. When necessary additional dilutions should be made. The dilution factor is incorporated into final calculations.

4.2.7.12 **Reporting of Results**

4.2.7.12.1 **Breath Testing Solutions**

4.2.8.12.1.1 Provide results to the Breath Testing Program Manager for evaluation.

4.2.7.12.2 **Alcohol Beverages**

4.2.7.12.2.1 To obtain the ethanol concentration value the mean results of analysis should be multiplied by the dilution factor. This will provide the ethanol concentration in g/100cc (weight per volume (w/v) percent).

4.2.7.12.2.2 For volume per volume (v/v) value divide w/v value by 0.79.

4.2.7.12.2.3 Value should be reported as both w/v and v/v percent. The mean value should be truncated and reported as a whole number.

4.2.7.12.3 **Unknown Liquids and "Serum" - Ethanol**

4.2.7.12.3.1 Ethanol results should be reported concentration in g/100cc and/or weight per volume (w/v) percent depending on the sample history.

4.2.7.12.3.2 When dilution is necessary the mean results of analysis should be multiplied by the dilution factor.

4.2.7.12.4 **Unknown Liquids and "Serum" - Other Volatiles**

4.2.8.12.4.1 The qualitative presence of other volatiles such as acetone, isopropyl alcohol, methyl alcohol, toluene and formaldehyde should be noted on the analysis report.

**4.2.8 ANALYSIS DOCUMENTATION**

4.2.8.1 Controls and standards are to be included in individual case files and/or a packet for the batch of samples is prepared.

4.2.8.2 The packet containing original data for controls and standards will be prepared for the analysis run and stored centrally in the file designated



for alcohol/volatiles quality assurance data in the laboratory where the analysis was performed until archiving.

- 4.2.8.3 When necessary, a copy of the control and standard printouts can be prepared from the centrally stored document.

## 4.2.9 QUALITY ASSURANCE

### 4.2.9.1 General

- 4.2.9.1.1 Samples are to be refrigerated while at the laboratory.
- 4.2.9.1.2 Refer to toxicology manual section 5.1 for pipette calibration options.
- 4.2.9.1.3 Refer to toxicology manual section 5.2 for balance calibration requirements.
- 4.2.9.1.4 Refer to toxicology manual section 5.3.2 for GC-HS maintenance requirements.
- 4.2.9.1.5 Calibrator solutions should be ordered prior to the current supply running out. This will allow for the analysis of new lots against existing calibrators.

### 4.2.9.2 Quality Control Requirements Per Analysis Run

- 4.2.9.2.1 An internal standard blank should follow the highest ethanol calibrator.
- 4.2.9.2.2 For the analysis of simulator solutions, a minimum of two bottles of particular lot of simulator solution should be sampled.
- 4.2.9.2.3 For an analysis run which involves the quantitation of ethanol, two aqueous controls of the same ethanol concentration/level must be run per batch of 10 samples (20 vials). For each additional 10 samples, a minimum of one additional aqueous ethanol control should be run.
- 4.2.9.2.4 An aqueous control containing ethanol with or without other volatiles substances meets the "per run" requirement as described in 4.2.9.2.3.
- 4.2.9.2.5 In an analysis run which involves the qualitative identification of volatiles other than ethanol a multicomponent mixture, and/or a single component aqueous standard containing each volatile to be identified, should be run.

- 4.2.9.3 Periodic Quality Control Requirements
- 4.2.9.3.1 The aqueous control concentration(s) should be varied periodically.
- 4.2.9.3.2 Periodically run either the Volatile Standard Mix Solution or the Multicomponent Alcohol Calibration Kit solution to determine and monitor the retention time of “other” volatiles of interest.
- 4.2.9.4 Monitoring of Quality Control Values
- 4.2.9.4.1 On a monthly basis, calculate the mean and standard deviation of quality control samples. The data will serve as a continual check of manufacturer-supplied values.
- 4.2.9.4.2 All control data will be provided monthly to the Discipline Leader for the Toxicology Program.

#### 4.2.10 AUTHENTICATION OF REFERENCE MATERIALS

- 4.2.10.1 Aqueous Ethanol/Volatile Standards
- 4.2.10.1.1 Standards for quantitative purposes must be traceable to NIST standards (or comparable).
- 4.2.10.1.2 *Certificate of Analysis* for all standards will be stored centrally.
- 4.2.10.1.3 New lots of Ethanol/Volatile standards should be included in duplicate in a minimum of one analysis run prior to official use.
- 4.2.10.1.4 Standards authenticated prior to the start date of this SOP revision can be used until consumed. The authentication data must be centrally stored.
- 4.2.10.1.5 **Ethanol**  
The Certificate of Analysis for an aqueous ethanol calibration standard together with a comparison of relative retention time and quantitation data, against existing calibrators, will serve as the authentication of ethanol in the standard. The new lot number can be accepted if the mean relative retention time for the new standard is  $\pm 0.10$  minutes and the mean concentration obtained falls within 6% of their target value.

4.2.10.1.6 **Volatile Standards**  
 For standards (acetone, ethanol, methanol, isopropanol, ...) used to prepare single constituent or mixed standard of volatiles, the qualitative authentication is established with the *Certificate of Analysis* and comparison of relative retention times. The new lot number can be accepted if the mean relative retention time (RRT) for the new standard is  $\pm 0.10$  minutes from the RRT of existing the qualitative standard components.

4.2.10.2 Commercially Obtained Aqueous Volatile Mixtures

4.2.10.2.1 *Certificate of Analysis* will be stored centrally.

4.2.10.2.2 The Certificate of Analysis for an aqueous mixed volatile standard along with a comparison to data from previous runs will serve as the qualitative authentication of the standard. The solution prepared with a new lot number of volatile chemical standard can be accepted if the mean relative retention time for the new standard is  $\pm 0.10$  minutes.

4.2.10.2.3 Refer to Certificate of Analysis for purity information.

4.2.10.2.4 Standards authenticated prior to the start date of this SOP revision can be used until consumed. The authentication data must be centrally stored.

4.2.11 REFERENCES

- 4.2.11.1 Stafford, D.T., Chromatography. in: Principles of Forensic Toxicology, edited by Barry Levin, pp. 93-101, 103-114, AACC Press, 1999.
- 4.2.11.2 Levine, B., Alcohol. in: Principles of Forensic Toxicology, edited by Barry Levin, pp. 170-184, AACC Press, 1999.
- 4.2.11.3 Caplan, Y.H., The Determination of Alcohol in Blood and Breath. in: Forensic Science Handbook, edited by Richard Saferstein, pp. 594-648, Prentice-Hall New Jersey, 1982.
- 4.2.11.4 Julien, R.M., Central Nervous System Depressants: Alcohol and the Inhalants of Abuse, in: Primer of Drug Action, pp. 64-92, Freeman-New York, 1998.

- 4.2.11.5 Saker, E.G., Screening and Quantitation by Head Space Technique of Some of the Vapors Most Commonly Found in Forensic Toxicology, in: Current Approaches in Forensic Toxicology, Chapter 11, SOFT Meeting, 1994.
- 4.2.11.6 Perrine, D.M., Depressants: Alcohol, Benzodiazepines, Barbiturates, in: The Chemistry of Mind-Altering Drugs, pp. 113-129, ACS, Washington, DC, 1996.
- 4.2.11.7 Hobbs, W.R., Rall, T.W. and Verdoorn, T.A., Drugs Acting on the Central Nervous System - Hypnotics and Sedatives; Ethanol, in: Goodman and Gilman's The Pharmacological Basis of Therapeutics, pp. 361, 386-393, McGraw-Hill, 1996.
- 4.2.11.8 Idaho Administration Code, IDAPA 11.03.01, Rules Governing Alcohol Testing.
- 4.2.11.9 Christmore, D.S., Kelly, R.C. and Doshier, L.A. Improved Recovery and Stability of Ethanol in Automated Headspace Analysis, J. Forensic Sci. 29(4): 1038-1044; 1984.
- 4.2.11.10 Restek Applications Note #59598, Dual-Column Confirmational GC Analysis of Blood Alcohols Using the Rtx<sup>®</sup>-BAC1 and Rtx<sup>®</sup>-BAC2 Columns Optimized for the Perkin-Elmer HS-40 Headspace Autosampler, 1999.

**Idaho State Police  
Forensic Services  
Toxicology Section**



Section Four  
**Analysis of Alcohol and Common Volatile Solvents**

**4.2 Quantitative Analysis of Ethanol Containing Solutions**

<b>Revision #</b>	<b>Issue Date</b>	<b>History</b>
1	01-03-03	Original issue in SOP format
2	05-03-04	Clarifications, incorporation of serum and other toxicology unknown solutions, added volatiles other than ethanol. Validation issues covered when SOP 4.1 was validated for "other volatiles" therefore no validation necessary.

**Approval**

**Discipline Leader:** \_\_\_\_\_ **Date:** \_\_\_\_\_  
Susan C. Williamson

**Issuance**

**QC Manager:** \_\_\_\_\_ **Date:** \_\_\_\_\_  
Richard D. Groff