


Idaho State Police

Forensic Services

Approval for Quality System Controlled Documents



Discipline/Name of Document: Toxicology


3.6.1 Liquid-Liquid Extraction Procedure for the Recovery of Basic and Neutral Drugs from Blood

Revision Number: 3

Issue Date: 11/21/2006

APPROVED BY:


Quality Manager


Date Signed

**Idaho State Police
Forensic Services
Toxicology Discipline**

**Section Three
Blood Toxicology**

3.6 Qualitative Liquid-Liquid Extraction Methods for GC/MSD Confirmation

3.6.1 Liquid-Liquid Extraction Procedure for the Recovery of Neutral and Basic Drugs from Blood.

3.6.1.1 BACKGROUND

This method is a general blood extraction procedure for a variety of commonly encountered neutral and basic drugs along with their metabolites. This method prepares an extract which will be subjected to confirmatory analysis with a gas chromatograph with a mass selective detector (GC/MSD). This method does not efficiently extract some basic compounds, such as morphine and hydromorphone, due to pKa/pH considerations.

3.6.1.2 SCOPE

The method is based upon the principle of liquid/liquid extraction. Positive controls are spiked for a resulting concentration of 200ng/mL or 500ng/mL of drugs of interest. The sample pH is adjusted with a pH 9.2 borate buffer and extracted with n-butyl chloride. An optional back extraction procedure removes most frequently encountered interfering substances. Two internal standards are used to monitor extraction efficiency and chromatographic performance. Gas chromatography in conjunction with full scan mass spectrometry is used to confirm the presence of basic and neutral analytes of interest.

3.6.1.3 EQUIPMENT AND SUPPLIES

- 3.6.1.3.1 Tube Rocker
- 3.6.1.3.2 Vortex Mixer
- 3.6.1.3.3 Laboratory Oven or Drybath
- 3.6.1.3.4 Evaporative concentrator equipped with nitrogen tank.
- 3.6.1.3.5 16 x 100mm round bottom glass screw-top tubes
- 3.6.1.3.6 Screw Cap for 16mm O.D. tubes
- 3.6.1.3.7 GC/MS Automated Liquid Sampler (ALS) vials
- 3.6.1.3.8 GC/MS vial microinsert
- 3.3.1.3.9 Gas Chromatograph equipped with a Mass Selective Detector
- 3.3.1.3.10 5%-Diphenyl-95%-Dimethyl-siloxane copolymer capillary GC column, 12.5 to 30M.

3.6.1.4 REAGENTS

Refer to Manual section 5.12 for solution preparation instructions.

- 3.6.1.4.1 Methanol (Certified ACS Grade)
- 3.6.1.4.2 n-Butyl chloride (Certified ACS Grade)
- 3.6.1.4.3 Borate Buffer (pH 9.2)
- 3.6.1.4.4 Deionized/Distilled (DI) Water
- 3.6.1.4.5 1% Hydrochloric Acid in Methanol
- 3.6.1.4.6 100mM Hydrochloric Acid

3.6.1.5 REFERENCE MATERIAL**3.6.1.5.1 Positive Control**

Positive Control can be prepared with the working solution described below and/or obtained commercially.

3.6.1.5.1.1 Positive Control Stock Solution

Obtain 1mg/mL stock drug standard solutions through Cerilliant, Alltech, Sigma or other appropriate vendor.

3.6.1.5.1.2 Positive Control Working Solution

Add the designated volume of stock solution to 10mL methanol. A minimum of four compounds must be used.

| Stock Solution (1.0mg/mL) | Volume (μ L) |
|------------------------------|----------------------|
| Amitriptyline | 20 |
| Caffeine | 20 |
| Codeine | 20 |
| Diphenhydramine | 20 |
| Lidocaine | 20 |
| Meperidine | 20 |
| Methadone | 20 |
| Methamphetamine | 20 |
| Nicotine | 20 |
| PCP | 20 |
| Trazodone | 50 |

Solution is stable for 6-months when stored at room temperature.

3.6.1.5.2 Internal Standard**3.6.1.5.2.1 Stock Solutions**

- 1 mg/mL Benzphetamine
- 1mg/mL Papaverine

3.6.1.5.2.2 **Working Internal Standard Solution**
[10ng/ μ L]

Add 100 μ L Benzphetamine and Papaverine stock solutions to 10mL volumetric ball flask. QS with methanol.

Solution is stable for three months when stored at room temperature.

3.6.1.5.3 Negative Control
Negative Whole Blood

3.6.1.5.4 Reference Standards

3.6.1.5.4.1 **Stock Standard Solution**

Obtain 1mg/mL stock drug standard solutions through Cerilliant, Alltech, Sigma or other appropriate vendor.

3.6.1.5.4.2 **Working Drug Standard Solution**

Add the designated volume of stock solution to 10mL methanol.

3.6.1.6 **PROCEDURE**

3.6.1.6.1 Initial set-up

For each control and case sample label two sets of screw-top extraction tubes and an ALS vial with microinserts.

3.6.1.6.2 Positive Control Samples

3.6.1.6.2.1 Prepare two positive control samples by adding 200 μ L mixed working control solution to 2mL negative whole blood (Utak 44600-WB (F) or equivalent) or pipette two 2mL samples of commercially obtained whole blood positive control.

3.6.1.6.2.2 When the optional back extraction is used, prepare 2 additional positive controls to parallel the back extraction process.

3.6.1.6.3 Casework Samples

Transfer 2mL casework and negative control samples to screw top extraction tube.

3.6.1.6.4 Negative Control Sample

Transfer 2mL negative whole blood to extraction tube.

3.6.1.6.5 Internal Standard

- 3.6.1.6.5.1 Add 20µL of internal standard mixture. Vortex.
- 3.6.1.6.5.2 Allow sample to stand 10 minutes.

3.6.1.6.6 Initial Extraction

- 3.6.1.6.6.1 Add 2mL borate buffer. Vortex.
- 3.6.1.6.6.2 Pipet 4mL n-butyl chloride into each tube, cap.
- 3.6.1.6.6.3 Place tube on rocker for a minimum of 10 minutes.
- 3.6.1.6.6.4 Centrifuge for 10 minutes at 3200 – 3400 rpm.
- 3.6.1.6.6.5 Transfer the n-butyl chloride layer (upper) to second screw-top tube.
- 3.6.1.6.6.6 Add 50µL 1% HCl in methanol.
- 3.6.1.6.6.7 Evaporate to dryness under N₂ at ≤37°C.

3.6.1.6.7 Optional Sample Clean-up

- 3.6.1.6.7.1 Reconstitute with 50µl 100mM HCl.
- 3.6.1.6.7.2 Add 1ml of n-Butyl Chloride. Vortex.
- 3.6.1.6.7.3 Rock for 5 minutes.
- 3.6.1.6.7.4 Centrifuge for 5 minutes at 3200-3400 rpm.
- 3.6.1.6.7.5 Discard upper n-Butyl Chloride layer.
- 3.6.1.6.7.6 Add 2ml of pH 9.2 borate solution. Vortex
- 3.6.1.6.7.7 Add 4 ml of n-Butyl Chloride.
- 3.6.1.6.7.8 Rock for 5 minutes.
- 3.6.1.6.7.9 Centrifuge for 5 minutes at 3200 - 3400 rpm.
- 3.6.1.6.7.10 Transfer upper n-Butyl Chloride layer into screw-top tube.

3.6.1.6.7.11 Evaporate to just dryness under N₂ at ≤45°C.

3.6.1.6.8 Reconstitution

3.6.1.6.8.1 Add 50uL Methanol to the residue, vortex.

3.6.1.6.8.2 Transfer extract to labeled ALS vial with microinsert.

3.6.1.6.9 Preparation for Analysis Run

3.6.1.6.9.1 Into Sequence log table, enter the sample case numbers, blanks and controls.

3.6.1.6.9.2 Load samples, standards, blank and controls into the quadrant rack as noted in the sequence table.

3.6.1.6.10 GC-MSD Analysis Parameters

3.6.1.6.10.1 Refer to instrument METHOD printout for current analysis parameters.

3.6.1.6.10.2 Current analysis method must be stored centrally as a hard or electronic copy.

3.6.1.6.11 GC-MSD Qualitative Detection and Identification Criteria

3.6.1.6.11.1 For the identification of compounds not included in positive control, analyze appropriate non-extracted reference standards.

3.6.1.6.11.2 The presence of a drug compound is indicated if the retention time for the sample versus applicable standard does not differ by more than ±0.2 minutes and there are no significant differences in the mass spectral data.

3.6.1.7 **QUALITY ASSURANCE REQUIREMENTS**

3.6.1.7.1 General

3.6.1.7.1.1 Blood samples are to be stored under refrigeration after aliquots are removed for analysis.

3.6.1.7.1.2 Refer to toxicology analytical method 5.2 for balance calibration requirements.

3.6.1.7.1.3 Refer to toxicology analytical method 5.3.1 for GC-MSD maintenance guidelines.

3.6.1.7.1.4 Refer to toxicology analytical methods 5.8 and 5.10 for reference standard authentication and additional GC-MSD quality assurance requirements.

3.6.1.8 ANALYSIS DOCUMENTATION

3.6.1.8.1 A packet containing original data for controls will be prepared for each analysis run and stored centrally in the laboratory where the analysis was performed until archiving.

3.6.1.8.2 A copy of controls need not be included in individual case files. When necessary, a copy of control printouts can be prepared from the centrally stored document.

3.6.1.9 REFERENCES

3.6.1.9.1 Procedure for Basic Drug Analysis, Courtesy of Jim Hutchison, Montana Department of Justice, Forensic Services Division, 2005.

3.6.1.9.2 Procedure for Back Extraction, Courtesy of Jim Hutchison, Montana Department of Justice, Forensic Services Division, 2006.

3.6.1.9.3 Strong Bases Extractions - Screening SOP, Courtesy of Dr. Graham Jones, Office of the Chief Medical Examiner, Edmonton, Canada, 2003.

3.6.1.9.4 Jones, G., *Postmortem Toxicology*. pp. 98-102, *in*: Clarke's Analysis of Drugs and Poisons, 3rd Edition, Moffat, A.C, Osselton, M.D. and Widdop, B., eds., Pharmaceutical Press, 2004.

3.6.1.9.5 Hearn, W.L. and Walls, H.C., Strategies for Postmortem Toxicology Investigation. pp. 937-939. *in*: Drug Abuse Handbook, S.B. Karch, ed., CRC Press, Boca Raton, FL, 1998.

**Idaho State Police
Forensic Services
Toxicology Discipline**

**Section Three
Blood Toxicology**

3.6 Qualitative Liquid-Liquid Extraction Methods for GC/MSD Confirmation

3.6.1 Liquid-Liquid Extraction Procedure for the Recovery of Neutral and Basic Drugs from Blood.

| Revision # | Issue Date | History |
|------------|------------|---|
| 1 | 04-25-2002 | Original Issue in SOP format |
| 2 | 05-27-2003 | Updated, Clarifications |
| 3 | 11-21-2006 | Addition of internal standard, extraction process restructured. |

Approval

Discipline Leader: _____ Date: _____
Susan C. Williamson

Issuance

QA Manager: _____ Date: _____
Alan C. Spanbauer