

Idaho State Police  
Forensic Services  
Toxicology Section

---

Section Three  
Blood Toxicology

---

**3.6 Liquid-Liquid Extraction Methods for GC/MSD Confirmation**  
**3.6.2 Liquid-Liquid Extraction Procedure for the Recovery of Acidic Drugs from Blood.**

---

**3.6.2.1 BACKGROUND**

This method is a general blood extraction procedure for a variety of commonly encountered acid drugs of abuse. This method prepares an extract that will be subject to confirmatory analysis by gas chromatography/mass spectrometry (GC/MS).

**3.6.2.2 PRINCIPLE**

The method is based upon the principle of liquid/liquid extraction. Acidic compounds can be extracted from blood samples under acidic conditions with an organic solvent. The sample is extracted with n-butyl chloride. Following centrifugation, the organic layer is transferred to a new extraction tube and 0.45N sodium hydroxide is added to back extract acidic analytes. The pH is then adjusted to  $\leq 6$  with concentrated HCl to convert analytes back to a non-ionic form for a final extraction with n-butyl chloride. The final extract is reconstituted with 1:1 hexane/ethanol for confirmation on the GC/MS using SIM and/or full scan monitoring. The sample is extracted with n-butyl chloride. For sample clean-up, the sample is then back extracted. The final extract is reconstituted with 1:1 hexane/ethanol for confirmation on the GC/MS using SIM and/or full scan monitoring.

**3.6.2.3 EQUIPMENT AND SUPPLIES**

- 3.6.2.3.1 Drybath (Fisher or equivalent)
- 3.6.2.3.2 Evaporative concentrator (Zymark TurboVap or equivalent) equipped with nitrogen tank.
- 3.6.2.3.3 Glassware
  - 13x100mm Screw top tubes (Fisher 14-959-35C or equivalent)
  - Screw cap for tubes (Fisher 14-930-15E or equivalent).
  - 16X144mm tapered tip centrifuge tubes (Fisher 05-538-41C or equivalent)
  - Snap Caps (Fisher 05-538-41N or equivalent)

- GC/MS Automated Liquid Sampler (ALS) vials (HP 5182-0865 or equivalent)
- GC/MS vial microinsert (HP 5183-2088 or equivalent)
- 3.6.2.3.4 pH paper (Fisher 09-876-17 or equivalent)
- 3.6.2.3.5 Gas chromatograph equipped with a mass selective detector (HP 6890/5973 or equivalent) and a nonpolar capillary column with a phase composition capable of efficiently separating amines, alkaloids, drugs compounds and other analytes encountered in toxicological specimens (e.g. 100%-dimethylpolysiloxane or 95%-dimethyl-polysiloxane with 5% diphenyl)

**3.6.2.4 REAGENTS**

*Refer to Manual section 3.8 for solution preparation*

- 3.6.2.4.1 Methanol (Fisher A412-4 or equivalent)
- 3.6.2.4.2 Deionized/Distilled (DI) Water
- 3.6.2.4.3 n-Butyl chloride (Fisher B416-1 or equivalent)
- 3.6.2.4.4 Concentrated Hydrochloric Acid (Fisher A144-500)
- 3.6.2.4.5 Hexane (Fisher H292-4 or equivalent)
- 3.6.2.4.6 Ethanol (Fisher A995-4 or equivalent)
- 3.6.2.4.7 Hexane/Ethanol 1:1
- 3.6.2.4.8 0.45N Sodium Hydroxide

**3.6.2.5 STANDARDS**

**3.6.2.5.1 Stock Standard Solution**

- 3.6.2.5.1.1 1.0mg/mL Drug standard (obtain as necessary from Cerilliant, Alltech, Sigma or equivalent vendor).

**3.6.2.5.2 Working Standard Solution (5000ng/mL)**

- 3.4.5.5.2.1 Add 50µL Stock Solution to 10mL Methanol.  
Solution is stable for 12 months when stored at 4°C.

**3.6.2.6 CONTROLS**

- 3.6.2.6.1 Liquid Whole Blood Positive Control (Utak 98818 or equivalent)
- 3.6.2.6.2 Liquid Whole Blood (Utak 44600-WB (F) or equivalent) spiked with working standard solution at 50, 100 and/or 500ng/ml (other levels may be used as needed). To 2mL of negative blood add working standard solution as indicated below.

Desired ng/mL	$\mu$ L Working Standard Solution
50	20
100	40
500	200

3.6.2.6.3 Liquid Whole Blood Negative Control (Utak 44600-WB (F) or equivalent)

### 3.6.2.7 PROCEDURE

#### 3.6.2.7.1 Initial set-up

Label test tubes, and GC/MS vials with microinserts for the negative control (NC), positive control (PC), and appropriate laboratory numbers.

#### 3.6.2.7.2 Sample Preparation

- Transfer 1mL sample, negative control and positive control to screw-top extraction tube.

#### 3.6.2.7.3 Initial Extraction

- Pipet 10mL n-butyl chloride into each tube, cap and extract for  $\geq 3$  minutes.
- Centrifuge for  $\geq 5$  minutes/ Transfer the butyl chloride (top) layer to a second tube.

**The following are clean-up steps. If the sample is clean, proceed to**

#### **3.6.2.7.6**

#### 3.6.2.7.4 Back Extraction

- Pipet 2.0mL of 0.45N sodium hydroxide, cap and extract for  $\geq 3$  minutes.
- Centrifuge for  $\geq 5$  minutes.
- Discard butyl chloride (top) layer.

#### 3.6.2.7.5 Final Extraction

- Add concentrated HCl until the pH is acidic ( $\leq 6$ ).
- Pipet 10mL butyl chloride into extraction tube, cap and extract for  $\geq 5$  minutes.
- Centrifuge for  $\geq 5$  minutes.
- Transfer the butyl chloride (top) layer into centrifuge tube.

#### 3.6.2.7.6 Evaporation and reconstitution

- Evaporate under a gentle stream of nitrogen at  $\leq 37^{\circ}\text{C}$ .
- Add 100uL of 1:1 hexane/ethanol to the residue.

- Vortex.
- Transfer extract to labeled GC/MSD ALS vial with microinsert.

### 3.6.2.8 GAS CHROMATOGRAPHY/MASS SPECTROMETRY (GC/MS) ANALYSIS

#### 3.6.2.8.1 Analysis Parameters

- 3.6.2.8.1.1 Inject 1  $\mu$ L into GC/MS using the ALS.
- 3.6.2.8.1.2 Analyze sample extract(s) in full scan acquisition or SIM monitoring the appropriate ions.
- 3.6.2.8.1.3 Refer to attached GC/MSD method printout for current analysis parameters.

#### 3.6.2.8.2 Detection and Identification Criteria

- 3.6.2.8.2.1 The presence of a drug compound can be established if there are no significant differences in the retention time and mass spectra for the sample versus standards.
  - Acceptable retention time window is +/- 2%.

Property of Idaho State Police Forensic Services  
Uncontrolled Interim Copy  
OBSOLETE DOCUMENT

Idaho State Police  
Forensic Services  
Toxicology Section

---

Section Three  
Blood Toxicology

---

**3.6 Liquid-Liquid Extraction Methods for GC/MSD Confirmation**  
**3.6.2 Liquid-Liquid Extraction Procedure for the Recovery of Acidic Drugs from Blood.**

---

<b>Revision #</b>	<b>Issue Date</b>	<b>History</b>
1	04-25-02	Original Issue in SOP format

---

**Approval**

---

**Technical Leader:** \_\_\_\_\_ **Date:** \_\_\_\_\_  
S. C. Williamson

---

**Issuance**

---

**QC Manager:** \_\_\_\_\_ **Date:** \_\_\_\_\_  
Rick D. Groff

Property of Idaho State Police Forensic Services  
Uncontrolled Internet Copy  
OBSOLETE DOCUMENT