

**Idaho State Police  
Forensic Services  
Toxicology Section**

---

**Section Two**

**2.4 Liquid-Liquid Extraction Methods for GC/MSD Confirmation**

**2.4.1 General Extraction of Urine Samples for Qualitative Confirmation of Basic and Neutral or Acidic Drugs**

---

**2.4.1.1 BACKGROUND**

These extraction procedures are extensions of the TOXI-LAB® TOXI-A and TOXI-B thin layer chromatography (TLC) drug detection systems. The samples are extracted as with the TLC system, however, instead of concentrating the extract onto a disc, the solvent extract is concentrated and placed into an automated liquid sampler (ALS) vial for analysis by a gas chromatograph equipped with a mass selective detector (GC/MSD).

**2.4.1.2 PRINCIPLE**

This procedure describes the extraction of drug compounds from urine. Depending upon the  $pK_a$  of a drug compound, either Toxi-A or Toxi-B tubes are used. Basic compounds are extracted with a Toxi-A tube. Addition of urine to the Toxi-A tubes results in the urine becoming alkaline (pH=9) into 1,2-Dichloroethane, dichloromethane, heptane and isopropanol. Acidic compounds are isolated from an acidified solution (pH=4.5) into methylene chloride and heptane with zinc chloride to facilitate the extraction process. The extraction is achieved with an Ansys Toxi-B extraction tube. Either resulting extract is analyzed by full scan GC/MS in EI mode.

**2.4.1.3 EQUIPMENT AND SUPPLIES**

- 2.4.1.3.1 Tube Rocker (Fisher Scientific or equivalent)
- 2.4.1.3.2 Electric Warmer with Omega-12 extraction solvent concentrator (Ansys 118/153)
- 2.4.1.3.3 Laboratory Centrifuge (Fisher Marathon or equivalent)
- 2.4.1.3.4 Disposable Aluminum Concentration Cups (Ansys 152)
- 2.4.1.3.5 Glassware
  - GC/MS vials (HP 5182-0865 or equivalent)
  - GC/MS vial microinsert (HP 5183-2088 or equivalent)
- 2.4.1.3.6 Gas Chromatograph equipped with a mass selective detector (HP 6890/5973) and a HP-5MS Ultra low bleed (5%-Diphenyl-95%-Dimethylsiloxane copolymer) capillary column (25M).

**2.4.1.4 REAGENTS**

2.4.1.4.1 ANSYS TOXI-TUBES A and B (Ansys 109A-100/ 109B-100)

**2.4.1.5 CONTROLS**

2.4.1.5.1 Toxi-Control No. 19 – Morphine, amphetamine, imipramine, methadone, propoxyphene, phenobarbital, secobarbital and benzoylecgonine (Ansys 191AB).

2.4.1.5.2 Toxi-Control No. 2 – Amphetamine, methamphetamine, nicotine and cotinine. (Ansys 170B).

**2.4.1.6 STANDARDS**

2.4.1.6.1 Run necessary analytical standards as indicated by examination of GC/MSD data.

**2.4.1.7 PROCEDURE**

2.4.1.7.1 Initial set-up

Label TOXI-TUBES A or B, and GC/MS vials with microinserts with negative control, TC-19 and or TC-2 and appropriate laboratory numbers.

2.4.1.7.2 Extraction Procedure Toxi-A Extraction (Basic or Neutral Compounds)

Transfer 5 mL of urine specimen, negative urine or appropriate Toxi-Control to a TOXI-TUBE A (pH=9).

Rock TOXI-TUBE A for 15 minutes.

2.4.1.7.3 Centrifuge tube at 2500 rpm for 15 minutes.

2.4.1.7.4 Transfer solvent from tube into concentration cup in Omega-12 extraction solvent concentrator. Allow cups to warm prior to the addition of extract.

2.4.1.7.5 Evaporate solvent to approximately 200 $\mu$ L on electric warmer.

2.4.1.7.6 Transfer solvent to labeled GC/MS ALS vial with microinsert.

2.4.1.7.7 Extraction Procedure Toxi-B Extraction (Acidic Compounds)

2.4.1.7.8 Transfer 4.5 mL of urine specimen, negative urine or Toxi-Control 19, to a TOXI-TUBE B (pH=4.5).

2.4.1.7.9 Rock TOXI-TUBE B for 15 minutes.

2.4.1.7.10 Centrifuge tube at 2500 rpm for 15 minutes.

2.4.1.7.11 Transfer solvent from tube into concentration cup in Omega-12 extraction solvent concentrator.

- 2.4.1.7.12 Evaporate solvent to approximately 200 $\mu$ L on electric warmer.
- 2.4.1.7.13 Transfer solvent to labeled GC/MS ALS vial with microinsert.

- 2.4.1.7.4 Gas Chromatography/Mass Spectrometry (GC/MS) Analysis
  - 2.4.1.7.4.1 Inject 1  $\mu$ L into GC/MS using the ALS.
  - 2.4.1.7.4.2 Analyze sample extract in full scan acquisition. Refer to attached GC/MSD method printout for current analysis parameters.
- 2.4.1.7.5 Detection and Identification Criteria
  - 2.4.1.7.5.1 The presence of a drug compound can be established if there are no significant differences in the retention time and mass spectra for the sample versus standards.
  - 2.4.1.7.5.2 Acceptable retention time window is  $\pm$ 5%.

Property of Idaho State Police Forensic Services  
Uncontrolled Internet Copy  
OBSOLETE DOCUMENT