


ISP FORENSIC BIOLOGY QUALITY MANUAL




Property

ces

Forensic Biology Quality Manual

Revision #11



APPROVED
August 29, 2011

Cynthia Cunnington
Forensic Scientist IV
Forensic Biology Supervisor/Technical Leader

Property of Idaho State Police Forensic Services
Uncontrolled Copy
OBSOLETE DOCUMENT

Forensic Biology Quality Manual
Table of Contents

- i. Signature/Approval page
 - ii. Revision History
 - iii. Table of Contents
 1. Introduction: Statement of Purpose and Objectives
 2. Organization and Management
 3. Personnel Qualifications and Training
 4. Facilities
 5. Evidence Control
 6. Validation
 7. Chemicals and Reagents
 8. Calibration and Maintenance
 9. Proficiency Testing
 10. Corrective Action
 11. Reports
 12. Review
 13. Safety
 14. Audits
 15. Outsourcing
 16. Practices, Methods and Forms
 17. Controlled Forms
 18. Non-Controlled Forms
- Appendix A (Stutter and Peak Height Ratio Studies)

Appendix B (Mixtures- minor component calculations)

Appendix C (DNA Database Implementation)

Appendix D (Staff Abbreviation Lists)

Appendix E (Pre-2009 Revision Histories)

Appendix F (Annual Biology System Review)

Appendix G (Technical Lead Contingency Plan)

Appendix H (FBI Quality Assurance Documents)

Appendix I (Staff Acknowledgements)

Appendix J Extraction Control Keys

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

INTRODUCTION

The Forensic Biology Quality and Procedure Manuals are not public documents. Copies of the manuals, or portions thereof, will be released only to individuals having official business and upon proper discovery requests relating to a specific case(s).

1.0 STATEMENT OF PURPOSE AND OBJECTIVES

1.1 Statement of Purpose: ISP Forensic Biology exists to provide quality, unbiased and cost-effective analyses in the identification of biological substances and their source(s) relevant to the investigation and prosecution of criminal offenses in Idaho. The ISP Forensic Biology Quality Manual, along with the ISP Forensic Services Quality/Procedure Manual, provide the framework for the evaluation of QC (Quality Control) measures utilized in Forensic Biology to achieve that purpose. A system-wide mission and objectives are enumerated in the ISP Forensic Services Quality/Procedure Manual.

1.2 Objectives:

- 1.2.1 To develop and maintain, through annual review and revision (where necessary), a system of quality procedures, analytical methods, and controls to ensure quality up-to-date personnel training, biological screening and DNA analyses.
- 1.2.2 To evaluate (and revise where appropriate) through proficiency testing, audits, and other means of review, the thoroughness and effectiveness of biology personnel training, procedures and QC measures.
- 1.2.3 To remain scientifically neutral by basing case/evidence acceptance and analysis decisions, case reports and testimony solely on sound scientific rationale.
- 1.2.4 To develop and use practices that respect and protect the right of privacy for the genetic profiles developed in forensic casework or for database entry.
- 1.2.5 To provide high quality training, technical and informational assistance, biological analyses, written reports and testimony.
- 1.2.6 To provide all services in a cost-effective and timely manner.

2.0 ORGANIZATION AND MANAGEMENT


2.1 Organizational Chart and Functional Structure

2.1.1 An organizational chart for ISP Forensic Services appears in the ISP Forensic Services Quality/Procedure Manual. The Forensic Biology organization is delineated below.

2.1.2 An organizational chart for the Idaho State Police appears in the ISP Policy Manual.

2.2 Authority and Accountability in Forensic Biology

2.2.1 The Quality Assurance Standards for Forensic DNA Testing Laboratories and Convicted Offender DNA Databasing Laboratories, developed by the DAB, serve as a model for the ISP Forensic Biology QA Program. These standards delineate specific responsibilities and authority for the DNA Technical Manager and DNA CODIS Manager (see standard 4.1 of the FBI quality audit document). A copy of the document may be found in the ISP Forensic Biology Training Manual. Additionally, the ISP Forensic Services Quality/Procedure Manual designates specific authority for the DNA Technical Manager and DNA CODIS Manager.



Note: Changes (personnel) to this page do not require new revision numbers.

3.0 PERSONNEL QUALIFICATIONS AND TRAINING

3.1 Job Descriptions

General personnel qualifications and responsibilities, as well as personnel record retention policies, are described in the ISP Forensic Services Quality/Procedure Manual. Complete job descriptions are available through the Idaho Division of Human Resources web site:

(<http://dhr.idaho.gov/dhrapp/stateJobs/JobDescriptions.aspx>).

3.2 Training

Refer to ISP Forensic Biology Training Manual and the ISP Forensic Services Quality/Procedure Manual for specific training requirements and retention of training and continuing education records.

3.3 Continuing Education

Forensic Biology personnel must stay abreast of developments relevant to forensic DNA analyses through the attendance (and participation) at DNA related presentations, seminars, courses and/or professional meetings, for a minimum of 8 hours per calendar year. The training will also be supplemented through the routine reading of current scientific literature. The DNA technical Manager, or designee, will distribute a DNA-related article to each member of the biology section on a monthly basis. Each staff member will read the article and date/initial the attached sign-off sheet to indicate the completion of the reading. Additionally, the CODIS manager must stay abreast of developments relevant to CODIS/NDIS database management, computer and data security and computer networks through the attendance (personal or that of the Alternate CODIS Manager) at the bi-annual CODIS State Administrators' meetings and annual CODIS conference.

3.4 Qualifications

Education, training and experience for Forensic Biology personnel is formally established in the following minimum requirement specifications (Minimum requirements for individual positions may be reviewed at the time of job announcement and may exceed those delineated below). The minimum degree and education requirements are verified by review of transcripts as well as course descriptions, as necessary, during the application process. The DNA Technical Manager approves the degree and coursework prior to a job offer being extended to any potential hire. Periodic review of continuing education and overall performance is

accomplished during the annual employee evaluation. Opportunities are provided by an FS training budget.

3.4.1 Forensic Biology/DNA Supervisor/Technical Lead

It is assumed for the purposes of this document (and is currently the case), that in a laboratory system of the size of Idaho's, these functions will be served by a single individual.

3.4.1.1 Education

Must have at minimum, a Master of Science degree in a biological science. Successful completion of a minimum of 12 credit hours, including a combination of graduate and undergraduate coursework in genetics, biochemistry, molecular biology and statistics (or population genetics).

3.4.1.2 Training

Training and experience in molecular biology and DNA-based analyses from academic, governmental, private forensic and/or research laboratory(ies). Must also complete the FBI sponsored DNA auditor training within 1 year of appointment, if not already completed (dependant on FBI scheduling).

3.4.1.3 Experience

Must have a minimum of three years forensic human DNA laboratory experience as an analyst.

3.4.2 CODIS Administrator

This function may or may not be served by the Forensic Biology/DNA Supervisor. It is assumed for the purposes of this document (and is currently the case) that in a laboratory system of the size of Idaho's, the functions of casework and database CODIS Administrators will be served by a single individual. An Alternate CODIS Administrator will also be appointed and must meet the same qualifications as the CODIS Manager. The CODIS Administrator is responsible for administering the laboratory's CODIS network, scheduling and documenting the computer training for analysts, as well as assuring the security and quality of data and match dispositions all in accordance with state and/or federal law and NDIS operational procedures.

3.4.2.1 Education

Must have at minimum, a Bachelor of Science degree in a biological science and successfully completed college coursework in genetics, biochemistry, and molecular biology. Must also have completed coursework and/or training in statistics (or population genetics).

3.4.2.2 Training

A combination of training and experience in the use of computers, and database systems in a laboratory/scientific setting. Must also complete the FBI's CODIS software training and the DNA auditor training within six months of appointment if not already completed (dependant on FBI scheduling).

3.4.2.3 Experience

Must possess a working knowledge of computers, computer networks, computer database management and have an understanding of DNA profile interpretation for database and casework functions, to include mixture interpretation. Must be or have been a qualified DNA analyst.

3.4.3 DNA Analyst

The following delineate requirements for a DNA casework or database analyst whose responsibilities include performing genetic analyses on the capillary electrophoresis instruments and data interpretation. DNA extraction, quantification, and amplification set-up may be performed by appropriately trained laboratory technicians and/or those performing the biological screening of evidence following task-specific training and successful completion of a qualifying examination.

3.4.3.1 Education

Must have at minimum, a Bachelor of Science degree in a biological science and successfully completed college coursework in genetics, biochemistry, and molecular biology. Must also have completed coursework and/or training in statistics (or population genetics).

3.4.3.2 Training

Training in DNA analyses through academic, governmental, private forensic and/or research laboratory(ies). If received elsewhere, documented training must meet or exceed that outlined in the ISP Forensic Biology training manual. Must successfully complete a qualifying examination prior to performing analyses on database or forensic casework samples.

3.4.3.3 Experience

Must have a minimum of six months forensic human DNA laboratory experience.

3.4.4 Forensic Biologist

The following delineate requirements for those individuals responsible for the screening of evidence for the presence of biological substances and reporting and giving testimony regarding their findings.

3.4.4.1 Education

Must have a Bachelor of Science in a biological science.

3.4.4.2 Training

Training specific to this job function in a governmental and/or private forensic laboratory. If received elsewhere, documented training must meet or exceed that outlined in the ISP Forensic Biology training manual. Must successfully complete a qualifying examination prior to performing forensic casework.

3.4.4.3 Experience

Prior to participating in independent forensic casework, must have a minimum of six months forensic laboratory experience in the area of biological screening and/or DNA analysis.

3.4.5 Biology Laboratory Technician

3.4.5.1 Education

Minimum of two years of college to include scientific coursework (lecture and lab); Bachelor

of Science degree in a biological science is preferred.

3.4.5.2 Training

Must receive on the job training specific to assigned duties and successfully complete a qualifying examination before participating in forensic DNA typing or forensic casework responsibilities.

3.4.5.3 Experience

Prior to participating in any forensic DNA typing responsibilities or forensic case processing activities, technician must have a minimum of six months forensic laboratory experience in the area of Biology/DNA; one year is preferred.

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

4.0 FACILITIES

4.1 Laboratory Security

Security of the Forensic Services Laboratory is covered in the ISP Forensic Services Quality/Procedure Manual.

4.1.1 Forensic Biology Security

When not under the direct control of Forensic Biology personnel, evidence and in-progress work product will be secured either by closing and locking the Forensic Biology door or by its return to secure storage (one of the locked evidence refrigerators/freezers/file cabinets or the analyst's personal evidence cabinet). Only Forensic Biology Personnel will have access to the locked storage and laboratory areas. Persons outside the Forensic Biology unit will not be allowed access to the Forensic Biology laboratories. Exceptions will be made in case of emergencies, for maintenance, safety, and/or equipment service needs, and for required annual quality and DNA audits. At these times, access will be limited to only required individuals, the individual(s) will be accompanied by biology program personnel, and all evidence will be placed in secured storage for the duration of the individual(s) being present in the laboratory.

4.1.2 CODIS Security

The CODIS workstation is located in the locked CODIS office and the CODIS Server is located in the secured server room in the CJIS Section. The following security measures have been implemented:

4.1.2.1 Only Forensic Biology personnel will have access to the CODIS office. When a biology staff member is not present, the office will be secured by closing and locking the door.

4.1.2.2 Only the CODIS State Administrator, designated Forensic Biology staff and CJIS personnel will have access to the CODIS Server.

4.1.2.3 A differential backup of the CODIS server will be performed each weekday. A full backup will be performed once weekly with the backup tape being stored off-site. At any given time, one month of data will be stored offsite.

- 4.1.2.4 Only Forensic Biology Personnel that have gone through the NDIS application and approval process will have user-names and passwords for CODIS.
- 4.1.2.5 CODIS users must log in each time they use CODIS and log out prior to leaving the CODIS Workstation.
- 4.1.2.6 DNA Tracker, the convicted offender sample-tracking database, resides on the ISP intranet and is accessible only to personnel designated by the Biology/DNA Supervisor.
- 4.1.2.7 Personal and identifying information on convicted offenders (hard and electronic/DNA Tracker copies) are stored separately from the DNA profile (CODIS) obtained. The DNA profiles are directly associated only with a unique Idaho Convicted Offender ID number, assigned by DNA Tracker upon sample entry.
- 4.1.2.8 CODIS samples and corresponding information are released only in accordance with 19-5514 of the Idaho DNA Database Act of 1996, the Privacy Act Notice in Appendix E of NDIS procedures, and the FBI/CODIS Memorandum of Understanding.

4.2 Forensic Biology Laboratory Set-up

The Forensic Biology and Database Laboratories are designed to minimize contamination potential during the processing and analysis of forensic and convicted offender samples. Separate areas for evidence examination, DNA extraction, PCR Amplification Set-up and Amplified DNA processing and storage are delineated. Some steps of the pre-amplification processes may be conducted in the same area of the main laboratory; however, these steps are separated by time.

4.3 Laboratory Cleaning and Decontamination

In order to minimize the potential for sample contamination, careful cleaning of laboratory work areas and equipment must be conducted on a routine basis. The efficacy of the procedures used is monitored through the use of controls within the analysis process (see the interpretation guidelines section in BI-210 and BI-318). It is also important that each analyst use proper 'clean technique' at all times when in the laboratory, which includes but is not limited to, using only disposable barrier pipette tips and autoclaved microcentrifuge tubes, using a tube de-capping tool, and wearing gloves, a labcoat, and masks as appropriate.

- 4.3.1 All working **benchtop** surfaces will be cleaned with 10% bleach or Dispatch solution before and after use and as part of the monthly QC procedure. Clean white paper and/or a KayDry will be placed on the workbench prior to use and changed as appropriate and necessary.
- 4.3.2 All **small tools/instruments** (i.e. forceps, scissors, etc.) will be cleaned/rinsed with ethanol or germicidal instrument cleaner prior to use and between samples. Kimwipes, used to dry the instrument after cleaning/rinsing, will be single use only.
- 4.3.3 **Pipettes** are to be cleaned thoroughly with Dispatch solution as part of the monthly QC procedure and anytime the barrel comes in contact with DNA or any biological fluid.
- 4.3.4 All **centrifuges** are to be wiped down (interior and exterior) with Dispatch solution as part of the monthly QC procedure and in the event of a spill.
- 4.3.5 The **Biomek 3000** work surface trays and holders are to be removed and cleaned with 10% bleach or Dispatch solution as part of the monthly QC procedure or in the event of a spill. Additionally, each of tools are to be wiped down with ethanol, being careful not touch the electronic end.
- 4.3.6 The exterior surfaces of the **BSD600-Duet Puncher** are to be wiped down with a damp cloth, as part of the monthly QC procedure. In addition, the chute and punch mechanism are to be cleaned by removing and separating the inner and outer chutes. The inner chute is to be cleaned with ethanol, followed by compressed air blown through both chutes, the hole in the underside of the manifold, and between the punch guide and die. Do not use ethanol on the outer chute or around any electrical components.
- 4.3.7 The **thermal cyclers**, to include the heating block and exterior surfaces, are to be wiped down with ethanol or Dispatch solution as part of the monthly QC procedure. Individual wells should be cleaned as needed.
- 4.3.8 All work surfaces in the **amplification/post-amp rooms** are to be cleaned with 10% bleach or Dispatch solution before and after analysis and as part of the monthly QC procedure. Clean white paper and/or a KayDry is to be placed on the benchtop prior to use. Additionally, as part of the

monthly QC procedure, the following are to be conducted: the exterior surfaces of the genetic analyzers and real-time instruments wiped down with ethanol or Dispatch solution, top of the refrigerator/freezers and surface underneath each genetic analyzer wiped down/dusted, and floor mopped.

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

5.0 EVIDENCE CONTROL

Evidence, Individual Characteristic Database (Convicted Offender) samples, and in progress work product, that is collected, received, handled, sampled, analyzed and/or stored by ISP Forensic Services is done so in a manner to preserve its identity, integrity, condition and security.

5.1 Laboratory Evidence Control

Procedures detailing evidence handling are contained in the ISP Forensic Services Quality/Procedure Manual. Portions of individual evidence items that are carried through the analysis process (i.e. substrate cuttings, extracts, amplified product and/or portions thereof) are considered work product while in the process of analysis and do not require sealing. Work product will be identified by labeling the individual sample tube with a unique identifier, or documenting the locations of individual samples within a plate of samples.

5.2 Forensic Biology Evidence Control/Sample Retention

5.2.1 DNA Packet

It has become increasingly important to retain evidence for possible future analyses and to secure samples for non-probative casework analyses that are necessary for the validation of any new technology. Therefore, a DNA packet is created for cases submitted for analysis to Forensic Biology, in which reference sample(s) are present, and/or positive Biological screening results are obtained (See BI-102). Any remaining DNA extracts, upon completion of analysis, will be placed into a sealed container (such as a plastic zip bag) and stored in the DNA packet.

5.2.2 Limited Sample

In every case, care should be taken to save ~1/2 of a sample for independent testing. If testing would consume all or nearly all of a sample and there is an identified suspect charged in the case, the accused must receive appropriate notification. Written and/or verbal notification will be given to the prosecuting attorney informing him/her of possible consumption and requesting defense counsel be notified of the situation. Before testing will commence, an allowance will be made for testing by another accredited laboratory agreed upon by both parties. Additionally, a letter from the prosecuting attorney must be received

by the laboratory indicating whether or not the sample may be consumed.

5.2.3 Amplified Product

Amplified DNA product will not be retained after 1) the report has been issued in the case or 2) review of the offender sample data has been completed and certified for CODIS entry. In cases where both the evidence and associated DNA extract have been consumed, the amplified product will be retained in a sealed container within the product room freezer.

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

6.0 VALIDATION

Procedures for the validation and/or performance verification of methods used in ISP Forensic Services are outlined in the ISP Forensic Services Quality/Procedure Manual. Validation/performance verification data, results and summaries for those methods employed in Forensic Biology will be maintained in that section.

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

7.0 CHEMICALS/REAGENTS

General laboratory policies and procedures regarding the purchase of chemicals and preparation of reagents are covered in the ISP Forensic Services Quality/Procedure Manual.

7.1 COMMERCIALLY PURCHASED CHEMICALS

7.1.1 Biology Personnel should consult the section's electronic Chemical Inventory Log (Form 400-QC) prior to ordering. Chemical grade requirements should be checked and ordered as appropriate. The date ordered should be reflected in the log to avoid duplicate orders. An entry for chemicals not currently on the inventory will be made at this time to reflect the chemical, source, and order date. This inventory will be audited annually, at a minimum, and a printout placed in the Forensic Biology Reagent Binder.

Note: An order form/document must be filled out and approved by the section supervisor (indicated by date and initials) prior to placing the order. Reference the forensic services approved chemical list prior to ordering new chemicals.

7.1.2 Upon receipt of a chemical or reagent, the Chemical Inventory Log will be updated to reflect the new lot number, received date, quantity received, and quantity in stock. The order date will be removed at this time. The chemical(s) will be marked with the date received and the individual's initials. If it is an outer container that the chemical/kit remains in until use, the inner container will be labeled with this information when removed for use. The following commercially purchased reagents do not have manufacturer expiration dates: Phenol:Chloroform (PCIAA), HiDi Formamide, and 10X Genetic Analyzer buffer. These will additionally be labeled with a laboratory assigned expiration date of 2 years from the date of receipt. Packing slips should be checked to ensure appropriate accounting, including proper reagent grade, where applicable (this will be indicated by dating and initialing the packing slip and making notations as necessary). The packing slip and corresponding order document will be retained in the biology section. If an MSDS sheet came with the chemical, the MSDS binder should be checked for the presence of an MSDS for that chemical. If one exists, no additional copy is kept; however, if a newer version is received, the old one should be replaced. If one does not

already exist, place the one received in the binder. For chemicals without a hard copy MSDS, the manufacturer or one of the following websites may be consulted for information as needed:

<http://www.hazard.com/msds>

<http://www.msds.com>

<http://www.ilpi.com/msds/>

Note: Critical Reagents listed in 7.3 will be tracked on the individual QC forms, rather than the chemical inventory log.

7.1.3 Expired chemicals will be disposed of in an appropriate manner.

7.2 REAGENTS PREPARED IN-HOUSE

7.2.1 All biology reagents will be made with great care, following all quality and safety procedures. A mask will be worn by analysts during reagent preparation to help avoid the potential for contamination. See 7.4 and 7.5 below for individual reagent recipes.

7.2.2 Each reagent has a corresponding form to document the making of the reagent and components used. This form must be filled out. A reagent label must be made that has the reagent name, the lab lot number (which consists of the first few letters of the reagent name followed by the date prepared, in the form 'MMDDYY'), and the preparer's initials. The NFPA designation will be completed on all labels. Refillable squirt-bottles of water or ethanol will be labeled but need not bear dates or initials.

7.2.3 An effort should be made to use in-house reagents within one year of preparation; however, they do not expire and may continue to be used beyond the one year timeframe.

7.3 CRITICAL REAGENTS

CRITICAL REAGENTS are those reagents that, if improperly functioning, could result in significant loss or destruction of DNA and are not amenable (or it's not practical) to testing immediately before (e.g., use on forensic samples) each use. The reagents listed below have been identified as critical in Forensic Biology/DNA. These reagents must undergo a **QC ASSAY BEFORE** use on forensic casework and/or Convicted Offender

samples. Reagents received at a later date but having the same lot number as those previously tested and determined acceptable need not have a QC check performed. Critical Reagents (in addition to other DNA-related reagents with manufacturer expiration dates) may be used beyond the listed expiration date for training purposes without any further testing, so long as expected results are obtained for all associated controls. The reagent must be labeled 'for training only' if it is to be retained once the expiration date has been reached.

ABACARD® HEMATRACE® TEST KIT (Form 410-QC)

OneStep ABACARD® p30 TEST KIT (Form 412-QC)

Quantifiler® Human DNA Quantification Kit (Form 419-QC)

STR Kit (Taq Polymerase checked with kits, Form 420-QC)

7.4 BIOLOGICAL SCREENING REAGENTS

Phenolphthalein (Kastle-Meyer) Reagent

(NFPA: health 3, flammability 1, reactivity 2)

May be a commercial purchase.

Phenolphthalein 2.0g

KOH 20.0g

Zinc (granular) 20.0g

Phenolphthalein, KOH, and 100ml of dH₂O are refluxed, in a fume hood, with Zinc until solution is colorless (producing phenolphthalin in ~4 hours). Store stock solution refrigerated in dark bottle to which ~5g mossy zinc has been added to keep the solution in its reduced form. Remove for working solution as needed.

Working solution: Mix 2ml stock solution with 8ml Ethanol

Caution: Zinc is flammable. The unreacted portions and used filter paper are to be disposed of properly.

Hydrogen Peroxide 3% (v/v)

(NFPA: health 0, flammability 0, reactivity 1)

Generally a commercial purchase, however, may be made from a 30% Solution (which is a commercial purchase) as follows:

Hydrogen Peroxide (30%) 10ml/90ml nanopure dH₂O

Mix the H₂O₂ with 90ml of nanopure dH₂O and store at ~4°C.

Ortho-Tolidine Reagent

(NFPA: health 3, flammability 1, reactivity 2)

O-Tolidine	0.6g
Glacial Acetic Acid	100ml
Ethanol	100ml

Dissolve O-tolidine in Acetic Acid/Ethanol mixture consistent with ratios above. O-tolidine is light sensitive and should be stored in dark reagent bottle and kept refrigerated when not in use.

Ammonium Hydroxide (~3%)

(NFPA: health 3, flammability 1, reactivity 2)

Ammonium Hydroxide (Concentrated ~30%) 10ml/100ml

Add the NH₄OH to 90ml of nanopure dH₂O, mix well and store at RT.

Ouchterlony Destain

(NFPA: health 3, flammability 3, reactivity 2)

Methanol	45ml
Distilled water	45ml
Glacial Acetic Acid	10ml

Mix well and store refrigerated.

Ouchterlony Stain

(NFPA: health 3, flammability 3, reactivity 2)

Ouchterlony Destain	50ml
Coomassie Blue (Brilliant Blue R-250)	0.1g

Mix well (overnight), filter, and store at RT.

10X Brentamine (Sodium Acetate) Buffer

(NFPA: health 2, flammability 2, reactivity 2)

Sodium Acetate (Anhydrous) 1.2g
Acetic Acid(to adjust to pH 5) ≈400µl

Dissolve Sodium Acetate in 10ml of nanopure dH₂O. Add Acetic Acid to pH 5. Store refrigerated.

Brentamine Solution A

(NFPA: health 1, flammability 0, reactivity 0)

O-Dianisidine Tetrazotized (Fast Blue B Salt) 50 mg
10X buffer pH 5 5 ml

Dissolve Fast Blue B Salt in 5 ml of 10X Brentamine Buffer. Store refrigerated in a dark container.

Brentamine Solution B

(NFPA: health 2, flammability 0, reactivity 0)

α-Naphthyl Phosphate (Disodium Salt) 50 mg

Dissolve in 5 ml of nanopure dH₂O. Store Refrigerated.

Saline (0.85% NaCl)

(NFPA: health 1, flammability 0, reactivity 0)

NaCl 4.25g/500ml

Dissolve the NaCl in 500 ml nanopure water. Sterilize by autoclaving. Store refrigerated.

1X Phosphate Buffered Saline (PBS)

(NFPA: health 1, flammability 0, reactivity 1)

PBS 1 commercial pre-made packet

Dissolve one packet of powdered PBS in 1l of nanopure dH₂O. Check that pH≈7.4, autoclave and store at RT.

If pre-made packets are not available, PBS may be prepared by dissolving 0.2g KCl, 8.0g NaCl, 0.2g KH₂PO₄, and 2.2g Na₂HPO₄·7H₂O (or 1.1g Na₂HPO₄ anhydrous) in 800ml nanopure dH₂O. Adjust pH to 7.4 if necessary. Q.S. to 1l with nanopure dH₂O, autoclave and store at RT.

X-mas Tree Stain Solution A (Kernechtrot Solution)

(NEPA: health 1, flammability 0, reactivity 0)

May be a commercial purchase.

Aluminum Sulfate	5g
Nuclear Fast Red	0.1g

For 100ml, Dissolve the Aluminum Sulfate in 100ml **HOT** nanopure dH₂O. Immediately add the Nuclear Fast Red, mix, cool and filter (paper or ≥45µm). May be stored at RT.

X-mas Tree Stain Solution B (Picroindigocarmine Solution)

(NEPA: health 2, flammability 2, reactivity 2)

May be a commercial purchase.

Saturated Picric Acid Solution	100ml
Indigo Carmine	0.33g

For 100ml, dissolve the Indigo Carmine in 100ml of the Picric Acid. Mix and filter (paper or ≥45µm). May be stored at RT.

Amylase Diffusion/Phosphate Buffer (pH 6.9)

(NEPA: health 1, flammability 0, reactivity 1)

NaH ₂ PO ₄ , anhydrous	2.7g
Na ₂ HPO ₄ , anhydrous	3.9g
NaCl	0.2g

Mix the above with 500ml dH₂O, adjust pH to 6.9, and store at RT.

Amylase Iodine Reagent

(NFPA: health 3, flammability 0, reactivity 2)

Potassium Iodide (KI)	1.65g
Iodine (I ₂)	2.54g

Dissolve the above in 30ml nanopure dH₂O heated to ~65°C. Mix well, filter and store at 4°C in an amber bottle. Dilute 1:100 for Amylase Diffusion Test.

Mercuric Chloride 10% (w/v)

(NFPA: health 4, flammability 0, reactivity 1)

Mercuric Chloride	10g/100ml 95% EtOH
-------------------	--------------------

Dissolve the Mercuric Chloride in 100ml of 95% Ethanol, mix well and store at RT.

Zinc Chloride 10% (w/v)

(NFPA: health 2, flammability 0, reactivity 2)

Zinc Chloride	10g/100ml 95% EtOH
---------------	--------------------

Dissolve the Zinc Chloride in 100ml of 95% Ethanol, mix well and store at RT.

7.5 DNA REAGENTS

1M Tris-HCl Buffer pH 7.5

(NFPA: health 2, flammability 1, reactivity 1)

Tris Base(tris[Hydroxymethyl]amino methane)	121.1 g
---	---------

Dissolve Tris in ~800 ml nanopure dH₂O. Adjust to pH7.5 at RT by adding concentrated HCl (approximately 65ml). Q.S. to 1l with nanopure dH₂O, autoclave and store at RT.

1M Tris-HCl Buffer pH 8

(NFPA: health 2, flammability 1, reactivity 1)

Tris Base(tris[Hydroxymethyl]amino methane)	121.1 g
---	---------

Dissolve Tris in ~800 ml nanopure dH₂O. Adjust to pH8 at RT by adding concentrated HCl (approximately 45ml). Q.S. to 1l with nanopure dH₂O, autoclave and store at RT.

0.5M Ethylenediamine Tetraacetic Acid (EDTA)

(NFPA: health 1, flammability 1, reactivity 0)

Na₂EDTA·2H₂O 186.1g/l

Slowly add EDTA to 800ml nanopure H₂O while stirring vigorously. Add ~20g of NaOH pellets to bring the pH to near 8.0. When fully dissolved adjust pH to 8.0 and bring final volume to 1l. Autoclave and store at RT.

Note: EDTA will not go into solution without the pH adjustment.

Stain Extraction Buffer pH8 (10mM EDTA/10mM Tris-HCl/50mM NaCl/2% SDS)

(NFPA: health 2, flammability 1, reactivity 1)

1M Tris-HCl, pH7.5	5ml
0.5M EDTA	10ml
5.0M NaCl	5ml
10% SDS	100ml

Mix the Tris-HCl, EDTA, NaCl and SDS with ~380ml nanopure dH₂O. Store at RT.

Note: Reagent contains SDS, do not autoclave.

Proteinase K (20mg/ml)

(NFPA: health 1, flammability 1, reactivity 0)

May be a commercial purchase of 20mg/ml solution.

Proteinase K 0.2g

Dissolve the ProK in 10ml sterile nanopure dH₂O.

Dispense ~500µl (commercial purchase or in-house prep.) each into sterile microfuge tubes and store at ≈20°C.

1M Sodium Acetate pH 5.2

(NFPA: health 3, flammability 2, reactivity 0)

CH₃COONa·3H₂O 13.6g

Dissolve the CH₃COONa·3H₂O in 80ml nanopure dH₂O. Adjust to pH5.2 by adding glacial acetic acid (approximately 2 ml). Q.S. to 100ml with nanopure dH₂O, autoclave and store at RT.

DTT Solution

(NFPA: health 2, flammability 1, reactivity 0)

Dithiothreitol (DTT) 0.77g

Dissolve the DTT in 5ml nanopure dH₂O. Add 50µl 1M Sodium Acetate, pH5.2. Dispense ~500µl each into sterile microcentrifuge tubes and store at ≈20°C.

Note: Do not autoclave.

PCR-TE (TE⁻⁴) Buffer (10mM Tris-HCl/0.1mM EDTA)

(NFPA: health 2, flammability 1, reactivity 0)

1M Tris-HCl, pH8 10ml
0.5M EDTA, pH8 0.2ml

Mix Tris-HCl and EDTA with 990ml nanopure dH₂O. Autoclave and store at RT.

5N Sodium Hydroxide

(NFPA: health 3, flammability 0, reactivity 2)

NaOH 50g

Slowly dissolve the Sodium Hydroxide in 250ml sterile nanopure dH₂O. Allow to cool and store at RT.

Caution: NaOH is highly caustic. This reaction generates heat.

5M Sodium Chloride

(NFPA: health 1, flammability 0, reactivity 0)

May be a commercial purchase of 5M solution.

NaCl 146.1g/500ml

Dissolve the NaCl in 500 ml nanopure water. Sterilize by autoclaving.

Bovine Serum Albumin 4%

(NFPA: health 0, flammability 1, reactivity 0)

BSA 0.4 g
PCR-TE 10 ml

Dissolve the BSA in PCR-TE. Filter-sterilize and dispense ~500µl each into 1.5ml microfuge tubes. Store at ~-20°C.

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

8.0 EQUIPMENT CALIBRATION AND MAINTENANCE

General laboratory procedures for the calibration and maintenance of equipment are covered in the ISP Forensic Services Quality/Procedure Manual.

8.1 BIOLOGY EQUIPMENT/INSTRUMENTATION

- 8.1.1 Analytical equipment significant to the results of examination and requiring routine calibration and/or performance verification will be listed on the BIOLOGY CRITICAL EQUIPMENT INVENTORY Spreadsheet (Form 401-QC). Information on the spreadsheet includes (as known or appropriate): equipment identity and its software, manufacturer's name, model, property number, serial number and/or unique identifier, and location. The inventory spreadsheet will be maintained in the instrument QC binder or section QC binder as appropriate.
- 8.1.2 OPERATING MANUALS for most equipment/instrumentation will be maintained in the product information file (Manuals for the ABI PRISM™ 310 and 3130/3130xl Genetic Analyzers, ABI 7500 Real-Time PCR System, Thermal Cyclers, and Driftcon FFC will be maintained in the Amp/PostAmp Room in close proximity to the instruments). Exceptions may be made for manuals referred to for instructions. In these cases, the manual will be maintained in close proximity to the instrument. The Biomek 3000 manual is built into the Biomek software.
- 8.1.3 MAINTENANCE/REPAIR/CALIBRATION LOGS will be maintained as follows:

The records for the ABI PRISM™ 310 and 3130/3130xl Genetic Analyzers, ABI 7500 Real-Time PCR System, and Thermal Cyclers will be maintained in the instrument QC binder.

Any equipment/instrumentation function (not documented on weekly, monthly, quarterly, or annual QC Check forms) will be recorded on the Equipment Maintenance/Repair form (Form 402-QC). Equipment Failure will also be reported on this form. This form and the QC check forms will be maintained in the section QC Binder, except as listed above.

8.1.4 EQUIPMENT FAILURE will result in that equipment being 'taken out of service'; an 'out of service' sign will be placed on the equipment and it will not be returned to service until it has passed appropriate performance testing. Actions are reported on Form 402-QC.

8.1.5 The SCHEDULE of QC/Performance Checks for both critical and non-critical equipment is as follows:

WEEKLY (Form 404A/B/C-QC)

(once per week with an interval between dates not less than 3 days and not exceeding 10 days)

- Nanopure System Check
- Refrigerator/Freezer Temperature Check
- Heating Block(s) Temperature Check
- Oven Temperature Check
- Water Bath Temperature Check

MONTHLY (Form 406A/B/C-QC)

(once per calendar month with an interval between dates not less than 15 days and not exceeding 45 days)

- Pipettes Cleaned
- Centrifuges Cleaned
- Biomek 3000 Cleaned
- BSD600 Cleaned
- Lab Cleaned
- Eye Wash Station Check
- Autoclave Clean and Check Sterilization
- ABI 7500 Background Assay/Contamination Test, and Function Test/Bulb Check
- BioRobot EZ1 grease D-rings
- 3130/3130x1 Water Wash
- 3130/3130x1 Water Trap Flush
- 310, 3130/3130x1, (C and E drives) and 7500 computer defragmentation

QUARTERLY

(once per quarter with an interval between dates not less than 30 days and not exceeding 120 days) Note: * denotes critical equipment

- Thermal Cycler* Temperature Verification

- ABI 7500* Temperature Verification
- Biomek 3000 Robotic System* Framing/Calibration Check (Form 408-QC)
- Chemical Shower Check (Form 408-QC)

ANNUALLY (Form 402-QC)

(once per calendar year with an interval between dates not less than 6 months and not exceeding 18 months) Note: * denotes critical equipment

- Mechanical Pipette* Performance Verification Check (outside vendor)
- NIST Traceable Thermometers* (outside vendor)
- Driftcon FFC Temperature Verification System* (outside vendor)
- Biological and Chemical Hoods Test (outside vendor)
- Digital Temperature Recording Devices Calibration Check (outside vendor)
- ABI PRISM™ 310* Genetic Analyzer Preventative Maintenance (outside vendor)
- ABI PRISM™ 3130/3130xl* Genetic Analyzer Preventative Maintenance (outside vendor)
- ABI 7500* Real-Time PCR System Preventative Maintenance (outside vendor)
- ABI 7500* Pure Dye Calibration, Optical Calibration, and Regions of Interest (ROI's) verification (see 7500 Maintenance Guide for procedures/may be part of PM by request)
- Qiagen BioRobot EZ1* Preventative Maintenance (outside vendor)
- Biomek 3000* preventative Maintenance (outside vendor)
- Microscope Cleaning/Preventative Maintenance (outside vendor)
- Centrifuge Calibration Check (outside vendor)
- Balance* Calibration Check (outside vendor)

In addition to the above schedule, personnel should check appropriate parameter function on all instrumentation with each use (including calibration of the pH meter at the time of use; documented on Form 403-QC), and run a matrix for the ABI PRISM™ 310 Genetic Analyzers and a spatial and spectral calibration for the ABI PRISM™ 3130/3130xl Genetic Analyzers as needed or following CCD camera and/or laser replacement/adjustment. Data for each new matrix will be filed in the instrument QC binder (see BI-210 and BI-318).

Following the annual preventative maintenance, a sensitivity panel (previously characterized DNA) should be run on the 310 and 3130/3130xl and included in the QC binder as a verification of performance. A color plate and framing/calibration check are to be run on the Biomek 3000, documented on Form 428-QC, and included in the Database QC binder as a

performance check following the annual preventative maintenance. The Driftcon FFC will be run on each thermal cycler (including 7500's) following repair and prior to being placed back in to service as a verification of performance. If no repairs were necessary, the pure dye calibration and ROI's will serve as the performance verification for the 7500's following the annual preventative maintenance. Documentation will be maintained in the section QC binder.

Any problems noted with laboratory equipment, during normal usage or as part of a QC check should be brought to the attention of the necessary supervisory personnel and documented on Form 402-QC and/or the respective QC form.

A certified NIST standard will also be run annually or if substantial procedural changes have been made. The QC run will be documented on Form 426-QC and filed in the section QC binder.

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

9.0 PROFICIENCY TESTING

General laboratory guidelines and practices for proficiency testing and retention are outlined in the ISP Forensic Services Quality/Procedure Manual. Additional Biology/DNA requirements are delineated below.

9.1 External DNA Proficiency Test Requirement. DNA analysts will participate in external proficiency tests, twice in every calendar year, in accordance with The FBI Quality Assurance Standards and the results reported to NDIS as necessary.

9.2 Inconclusive/Uninterpretable Proficiency Test Results.

Typically, sample size/quantity in PCR DNA Proficiency Tests is sufficient for multiple analyses to be performed. Therefore, results of DNA proficiency tests are not likely to be either inconclusive, or uninterpretable (e.g., not meeting minimal rfu and/or statistical threshold for inclusion/exclusion). However, in the event data obtained in a proficiency test does not meet the standard guidelines for interpretation/conclusion, it will first be determined, by re-testing and communication with the vendor, that this is not an issue with a given sample(s). Once that determination has been made, the analyst obtaining the inconclusive data will be removed from casework/database sample analysis until satisfactory completion of a competency test and review of the analyst's casework/database analysis performed since the last successful proficiency test.

10.0 CORRECTIVE ACTION

Laboratory corrective-action and retention procedures are detailed in the ISP Forensic Services Quality/Procedure Manual.

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

11.0 FILE DOCUMENTATION AND REPORTS

Meticulous documentation is an important aspect of forensic work. In casework, the scientist's knowledge of case circumstance (and therefore their ability to discern potential significance) may be limited. It is also common to be called upon to testify months, or even years, after processing evidence for a given case. Careful observation and detailed note-taking will not only refresh the scientist's memory and provide support for the conclusion in the laboratory report, but might also provide additional information not thought to have been important at the time of evidence processing. General laboratory policies regarding case record and retention are described in the ISP Forensic Services Quality/Procedure Manual. The note packet is considered complete when the analyst signs the report and submits the packet to be reviewed. Electronic documentation (eg. electropherograms and tables of results) are considered stored at this time. Any changes to the electronic documentation required after this point (typically on or after the review date documented in the note packet) will be made either by hand on the hard copy (initialed and dated by the analyst), or by changing the electronic version, reprinting and making a notation on the new hard copy as to the changes made. The new printed copy will bear the date the changes were made/reprinted.

11.1 CASE NOTES

- 11.1.1 Each page of case notes should have the following:
Laboratory Case Number, Date, Scientist's Initials and page number (in a form indicating page/total pages).
- 11.1.2 Case notes are associated with a particular report. Case notes for additional submissions (i.e., for supplemental reports) will be reflected in the page numbering as well (e.g. s1, supp. 1, etc.).
- 11.1.3 All evidence submitted for biological screening should be transferred to the scientist (i.e., documented on the chain of custody) and bear the scientist's initials. This is the case regardless of whether or not they analyze the item of evidence (exception may be made in cases where communication with investigator/attorney identified select items of those submitted). A description of the evidence (e.g., packaging and what it is said to contain) should

also appear in the case notes with a notation about not being examined at the time, if that's the case. Those items should also appear in the "not examined" statement of the report.

- 11.1.4** The description of evidence packaging should include the type and condition of seal(s). Differences in the description on a package versus ETS entry and/or accompanying submission form (or what the evidence is once opened) should be noted.
- 11.1.5** Whenever feasible, every attempt should be made to gain entry into the evidence without breaking the original seals. Any seal altered or created by a scientist will bear their initials and date across the seal.
- 11.1.6** Evidence descriptions should be "unique" inasmuch as possible (i.e., one pair blue jeans is **NOT** adequate). They should include, as appropriate and necessary for identification, colors, sizes (measurements where appropriate- e.g., knife and blade), manufacturer, model, brand, serial numbers or other identifiers and condition (e.g., worn, clean, torn, mud-caked, blood-soaked, etc.).
- 11.1.7** Photography, digital or otherwise, is often useful in documenting the appearance of evidence items. However, it is not meant to completely replace drawing, but instead as a supplement or in cases when drawing may be too difficult to accurately depict the item. Careful drawing and description result in careful and detailed examinations and, in many instances, may be a better choice than photography. Digital photographs will be transferred to, printed as necessary for case notes, and stored within the Mideo System; refer to BI-119 for Mideo instructions.
- 11.1.8** Evidence numbering must be unique for the purpose of possible later CODIS entry and chain of custody tracking. Items should be numbered as follows (or other similar system):

A single item (e.g., a baseball cap; Item 57) for which:

≤ 1 area tested positive for a biological substance and the stain is removed for DNA testing ≡ Item 57A (note: if the entire item is to be retained for DNA testing ≡ Item 57)

≥2 areas tested positive for a biological substance(s) (in this instance 3 areas removed for DNA testing)≡ Item 57-1, Item 57-2 and Item 57-3, or 57-A, 57-B and 57-C.

An item with multiple sub-items (e.g., a SAECK; Item 1) ≡ Item 1A, Item 1B, Item 1C, etc., the scientist should begin with the most relevant item if possible. Multiple areas ≡ Item 1A-1, Item 1A-2 etc.

11.1.9 The Biology Screening Case Summary Form (Form 101-BI) may be used for summarizing analyses if the scientist chooses.

11.1.10 If a form is used for more than one case, a copy of the 'completed' form should be made for any additional case files. A reference regarding the location of the original document(s) will be made in the note packet. For each file, the associated case should be listed and case data highlighted. In general, biology subfolders should be organized from front to back as follows: restitution where applicable, report, chronological case notes/forms, SAECK form where applicable, CODIS entry forms where applicable, case review forms where applicable, copy of evidence submission form or ETS property form, phone/info log ('tangerine' paper may be used for ease of identification), followed by agency materials submitted with evidence. Upon completion of review the analyst should bind (e.g. staple) the documentation together, with the exception of the restitution and report, and submit to the Forensic Evidence Specialists for report/restitution distribution.

11.2 DATABASE PACKETS

11.2.1 Each page of the database packet should have the following: Plate Identifier, Date, Scientist's Initials, and page number (in a form indicating page/total pages).

11.2.2 In general, database packets will be arranged from front to back as follows: chronological worksheets, reinjection summary and table of results (it is not necessary to print electropherograms for database packets). Review forms may be placed at the front of the packet for ease of plate

identification. Upon completion of review, the analyst should bind (e.g. staple) the documentation together and file it appropriately.

11.3 CASEWORK REPORTS

In the interest of consistency and clarity of reports between individual scientists the following format should be adhered to:

- 11.3.1** The report will contain the title Forensic Biology Report for biology screening reports, or Forensic DNA Report for DNA reports.
- 11.3.2** For clarity, when a statement(s) is about a particular Item (or multiple items listed individually), the "I" will be capitalized as in a name. When writing in general terms (i.e., the following items:) the "i" will remain lowercase.
- 11.3.3** The case submission information will include, at a minimum: case#, report date, case agency, agency case#, principals (victim, suspect, etc.), and offense date.
- 11.3.4** The body of the report will be separated from the case submission information by the following headings in the format below:

RESULTS AND INTERPRETATIONS

Statements (see below) regarding evidence exam, results and conclusions. The order of statements should be, inasmuch as possible: 1) positive statements (detection of body fluid), 2) inconclusive statements, 3) negative statements and 4) statements regarding (i.e. a list of) items not examined.

Disposition of Evidence

Statements (See below) regarding evidence retention and return.

Evidence Description

The following items were received in the laboratory via Federal Express (UPS, US Mail, etc.) on Month day, year. (or) The following items were received in the laboratory from Agency Representative (Agency) on Month day, year.

Description of items submitted for examination.

In the first report, all items should be listed (any items scientist took possession of, including reference samples). In supplemental reports, only those items relevant to the additional examinations need to be listed.

DNA reports, in which a DNA packet is checked out for analysis, will state: A tape sealed DNA packet envelope, created in the laboratory on Month day, year, and containing the following items:

Description of items contained within the DNA packet.

This report does or may contain opinions and/or interpretations, of the undersigned analyst, based on scientific data. The analyst's signature certifies that all of the above are true and accurate. (Note: the interpretations statement does not need to be included in reports where all items submitted are being returned without analysis, or other instances when no conclusions or interpretations are made.)

Signature

Name of Scientist
Title of Scientist

11.3.5 The following results/conclusions statements are to be used in a **biology screening report**, as dictated by the analysis findings (Where appropriate, descriptions, quantity, and/or locations of individual stains may be included in the corresponding statements. Portions of individual statements may be combined as needed.):

Semen Results/Conclusions Statements:

Chemical and microscopic analyses for the detection of semen were conducted on (items). Semen was confirmed by the presence of spermatozoa on (items). (or) Semen was not detected on (items). (or) No identifiable spermatozoa were detected on (items).

Chemical and microscopic analyses for the detection of semen were conducted on (items). Semen was confirmed on (items) by the presence of a single spermatozoon (or limited number of spermatozoa), which is (or may be) insufficient for further testing at this time.

Chemical, microscopic, and serological analyses for the detection of semen were conducted on (items). Semen was detected on (items) by the presence of the semen specific protein, p30; however, no spermatozoa were observed, which is insufficient for further testing at this time.

Results from presumptive chemical tests for the presence of semen were negative on (items).

Blood Results/Conclusion Statements:

Results from chemical and serological tests performed on (items) indicated the presence of human (or non-human) blood.

Results from presumptive chemical tests performed on (items) indicated the presence of blood; however, serological tests to determine the species of origin were not performed (or were inconclusive).

Results from presumptive chemical tests for the presence of blood were negative on (items).

Saliva Results/Conclusions Statements:

Results from chemical tests performed on (items) indicated the presence of an elevated level of amylase, an enzymatic component of saliva.

Results from chemical tests performed on (items) indicated (or did not indicate, or were inconclusive for) the presence of amylase, an enzymatic component of saliva.

Urine Results/Conclusions Statements:

Results from presumptive chemical tests performed on (items) indicated (or did not indicate, or were inconclusive for) the presence of urine.

Feces Results/Conclusions Statements:

Results from presumptive chemical tests performed on (items) indicated (or did not indicate, or were inconclusive for) the presence of feces.

Further Testing Statements (to be included at the end of the Results of Examination Section):

If additional testing is desired, please contact the laboratory.

DNA testing can be performed (or may be attempted) upon request and submission of a known reference sample(s) from [list name(s)]. Please contact the laboratory regarding the analysis request.

11.3.6 The following results/conclusions statements are to be used in an **STR DNA Report** (Note: the epithelial cell fraction of intimate samples, such as vaginal/rectal swabs, etc., are not considered probative if the testing results in a single profile matching the individual from which the sample was collected. In these instances, a statement regarding the DNA source of this fraction is not required):

Deoxyribonucleic Acid (DNA) Analysis, employing the Polymerase Chain Reaction (PCR), was used to generate a Short Tandem Repeat (STR) profile from the following items: "list of items".¹

Note: The following footnote will appear in all reports in which DNA testing was attempted.

¹Loci Examined: D3S1358, TH01, D21S11, D18S51, Penta E, D5S818, D13S317, D7S820, D16S539, CSF1PO, Penta D, vWA, D8S1179, TPOX, and FGA.

Profile Match Statement [meeting the 'source attribution' criterion (estimated frequency in population of ≤ 1 in 1.6×10^{10})] for single source and identifiable major contributors of a mixture:

The DNA profile obtained from the "item description (Item #)" matches that obtained from the blood stain/sample (or reference oral swab/sample, etc.) of/from "name". Therefore, "name" is the source of the "(DNA, blood, semen, saliva etc.) " on this item².

Note: The following footnote will appear in any report containing the above match statement.

²This conclusion is based upon the following: 1) a genetic match at the gender identity locus, Amelogenin, in addition to the "number" polymorphic STR loci listed above that have an expected population frequency of at least less than 1 in "actual (most conservative of the population groups calculated) frequency estimate", 2) a statistical frequency exceeding the source attribution criterion of 1.6×10^{10} (for $N=1.6 \times 10^7$, $\alpha=0.01$; *Forensic Science Communications* 2(3)July 2000), and 3) that "name" does not have a genetically identical twin.

Profile match Statement [not meeting the 'source attribution' criterion (estimated frequency in population of greater than 1 in 1.6×10^{10})] for single source and identifiable major contributors of a mixture:

The DNA profile obtained from the "item description (Item #)" matches that obtained from the blood/oral sample of "name". The probability of selecting an unrelated individual at random from the general population having a DNA profile that would match the DNA profile obtained from "item description (Item #)" is at least less than one in "actual (most conservative of the population groups calculated) frequency estimate".

Partial Profile Statement [profile consistent with item(s) in match statement above]:

The DNA profile obtained from the "item description (Item #)" also matches that obtained from the blood/oral sample of "name", however less genetic information was obtained.

The partial DNA profile obtained from the "item description (Item #)" is consistent with that obtained from the blood sample of "name".

Positive Paternity Statement [profiles consistent with being a biological child]:

Based upon evaluation of the DNA profiles obtained from the above individuals, "name" cannot be excluded as being the biological father of "name". The probability of paternity (assuming a prior probability of 0.5) is "X%" relative to an unrelated man randomly selected from the general population. The combined paternity index for the loci examined is "X". At least "X%" of the male population would be expected to be excluded from the possibility of being the biological father of "name".

Note: The most conservative of the population groups calculated is reported for the statement above.

Mixture Statements:

The DNA profile from "item description (Item#)" indicates a mixture of DNA from at least "X" persons. "Name(s)" is a potential contributor(s) to this mixture. "X%" of unrelated individuals randomly selected from the general population would be expected to be eliminated as potential contributors to this mixture.

The DNA profile from "item description (Item#)" indicates a mixture of DNA from at least two persons. "Name(s)" is a potential contributor(s) to this mixture. The DNA profile obtained from "item description (Item#)" is at least "X" times more likely to be seen if it were the result of a mixture of DNA from "name and name" than if it resulted from "name" and an unrelated individual randomly selected from the general population.

Note: The most conservative of the population groups calculated is reported for the statement above.

The DNA profile from "item description (Item#)" indicates a mixture of DNA with a discernable major contributor/profile. (include match, consistent with, or exclusionary statement regarding major profile). "name" is included/excluded/cannot be excluded as a possible contributor to the minor DNA component of this mixture.

The DNA profile from "item description (Item#)" indicates a mixture of DNA from at least "X" persons. "Name(s)" is a potential contributor(s) to this mixture. At least one in "actual (most conservative of the population groups calculated) frequency estimate" unrelated individuals randomly selected from the general population would be expected to be included as potential contributors to this mixture.

Exclusionary Statement:

The DNA profile obtained from the "item description (Item #)" does not match that obtained from the blood sample of "name". Therefore, "name" is not the source (or "a contributor" in a mixed profile situation) of the "(DNA, blood, semen, saliva etc.)" on this item.

The DNA profile obtained from the "item description (Item #)" was determined to be from an unknown male/female. "name" is not the source of the "(DNA, blood, semen, saliva etc.)" on this item.

Based upon evaluation of the DNA profiles obtained from the above individuals, "name" is not the biological father of "name".

No DNA Profile Obtained Statement:

Due to insufficient quantity or degradation, no DNA profile was obtained from "item description (Item #)".

CODIS Entry Statement:

The unknown male/female (included if source is not identified) DNA profile obtained from the "item description (Item #)" was entered into the Combined DNA Index System (CODIS) to be routinely searched against the database. The case agency will be notified in the event of a profile match.

Note: This statement is included when an eligible DNA profile has been developed, regardless of whether the profile is from a known or unknown source. Eligibility of forensic profiles for entry into CODIS and upload to NDIS is according to current NDIS procedures and include both solved and unsolved cases in which the profile is associated with a crime and believed to be attributable to the putative perpetrator. Profiles matching the victim(s) and any elimination samples (e.g. consensual partner samples) may not be entered.

11.3.7 The following statements are to be used in both biology screening and DNA STR reports:

Evidence Disposition Section Statements:

The following items have been retained in the laboratory [list all items/portions by description and Item# that have been retained in the DNA Packet (see BI-102)]. All remaining items (include empty evidence containers when the entire item was repackaged and retained) have been returned to the main laboratory evidence vault for return to the submitting agency.

The following items have been forwarded for DNA analysis. [list all items/portions by description and Item# that have been retained in the DNA Packet (see BI-102)]. Results will follow in a separate report. All remaining items (include empty evidence containers when the entire item was repackaged and retained) have been returned to the main laboratory evidence vault for return to the submitting agency.

Note: Nonsuspect cases (those with no known/identified suspect) in which biological evidence has been detected, will be forwarded for DNA testing and CODIS entry.

The DNA packet, which contains any remaining DNA extracts, has been retained in the laboratory. All remaining items have been returned to the main laboratory evidence vault for return to the submitting agency.

Evidence Description Section Examples:

A tape-sealed Sexual Assault Evidence Collection Kit (SAECK) containing biological samples, said to have been collected from "name".

A tape-sealed brown paper bag/manila evidence envelope/white cardboard box/etc. containing "description", (include the following if collection information is known) said to have been collected from "name" or "location".

A tape-sealed brown paper bag/manila evidence envelope/white cardboard box/etc. said to contain "label on package", (include the following if collection information is known) collected from "name" or "location".

A tape-sealed DNA packet, created in the laboratory on month day, year, and containing the following items:

- Item #) "description"
- Item #) "description"

11.3.8 It should be noted that the statements (in either the Forensic Biology Screening or DNA Reports) regarding evidence examination, testing and conclusions are not all-inclusive. There may be situations for which none of these statements is optimum.

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

12.0 REVIEW

Technical/administrative, document, and testimony (to include retention) review; as well as conflict resolution is addressed in the ISP Forensic Services Quality/Procedure Manual. See also, forms 214-BI and 306-BI in this manual.

12.1 BIOLOGY/DNA CASEWORK REVIEW

- 12.1.1 100% of the examinations and reports documented and/or issued from Forensic Biology/DNA will be "peer-reviewed". This review must be completed prior to issuing results (including verbal results) and/or entering eligible profiles into CODIS. Exceptions for release of results may be made on a case-by-case basis and with the Biology Supervisor's approval.
- 12.1.2 "Peer-review" in Forensic Biology will encompass both technical and administrative reviews.
- 12.1.3 The individual performing the "peer-review" will be a second scientist who is "qualified" in the area of the review (i.e., Biological Screening and/or STR Analysis).
- 12.1.4 It is not sufficient to have the scientist performing/reporting the analysis to be the sole person performing the administrative review.
- 12.1.5 The second scientist performing the review will initial each page (and date the first and last page at a minimum).
- 12.1.6 The second scientist will also place their initials below the signature of the scientist issuing the report.
- 12.1.7 Additionally, the second scientist will review the CODIS Entry Form (Form 218-BI) and verify that all eligible profiles have been identified for CODIS entry and the correct specimen categories have been assigned. The reviewer will date and initial the form. Eligible specimens will not be entered into CODIS until review/verification is complete. The specimen details report will be reviewed and initialed by the CODIS Administrator (or alternate) following manual data entry

and prior to searching at SDIS and uploading to NDIS to verify correct allele entry and specimen category.

- 12.1.8 Outsourced casework (when applicable) will undergo the same review as listed above, as well as for compliance with contract technical specifications.

12.2 CONVICTED OFFENDER/DATABASE SAMPLE REVIEW

- 12.2.1 100% of Convicted Offender sample data (including outsourced data when applicable) will be technically reviewed prior to CODIS entry and subsequent NDIS upload.
- 12.2.2 The individual performing the technical review will be a second scientist who is "qualified" in the area of STR Analysis.
- 12.2.3 The second scientist performing the review will initial each page of the data package (and date the first and last page at a minimum)
- 12.2.4 The scientist performing the review of outsourced data (when applicable) will document in an appropriate manner, the review of data for compliance with contract technical specifications and that the .cmf file, if present, contains the correct DNA profiles.
- 12.2.5 Additionally, a documented administrative review will be performed on CODIS hit confirmation letters containing an offender's personally identifiable information, prior to release.

12.3 TESTIMONY REVIEW

Review of courtroom testimony of Forensic Biology personnel shall be accomplished at least once in each calendar year. Preferably, this review will be performed by the Biology/DNA Supervisor or another qualified analyst and documented on the Forensic Services courtroom testimony evaluation form. Alternatively, the evaluation may be completed by criminal justice personnel (e.g., the judge, prosecutor or defense counsel).

13.0 SAFETY

Laboratory safety practices are addressed in the ISP Forensic Services Health and Safety Manual. In Forensic Biology, personnel are introduced to these practices in Module 1 of the ISP Forensic Biology Training Manual. In addition, Section 8 of this manual addresses the monitoring of the chemical eye-wash and shower.

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

14.0 AUDITS

Quality audits and retention schedules are delineated in the ISP Forensic Services Quality/Procedure Manual. Specific Biology/DNA audit requirements are delineated below.

- 14.1 A DNA audit, using the current FBI DNA Quality Assurance Standards Audit Document(s), will be conducted on an annual basis.
- 14.2 The interval between annual audits will be in accordance with the current FBI Quality Assurance Standards.
- 14.3 Every other year, at a minimum, the DNA audit must be an external audit.
- 14.4 The completed audit document(s) (Quality Assurance Standards Audit for Forensic DNA Testing Laboratories and for DNA Databasing Laboratories) and appropriate accompanying documentation will be submitted to NDIS according to NDIS Operational Procedures.

Property of Idaho State Police Forensic Services
Uncontrolled Internal Copy
OBSOLETE DOCUMENT

15.0 OUTSOURCING

Outsourcing/Subcontracting policies and procedures are described in the ISP Forensic Services Quality/Procedure Manual.

- 15.1 Approved vendor laboratories must provide documentation of accreditation and compliance with the Quality Assurance Standards for Forensic DNA and/or Database Testing Laboratories prior to contract award and for the duration of the contract.
- 15.2 Technical specifications will be outlined in the outsourcing agreement/contract and approved (approval will be documented) by the Biology/DNA Technical Manager prior to award.
- 15.3 An on-site visit of the vendor laboratory will be performed, by the technical leader or a qualified DNA analyst, and documented prior to the submission of any samples to that laboratory. Alternatively, the technical leader may review and accept (the review and acceptance will be documented) an on-site visit conducted by designated FBI personnel.
- 15.4 An annual on-site visit will be performed and documented for any contract extending beyond one year.
- 15.5 When outsourcing convicted offender samples, at least one quality control sample shall be included with each batch. Additionally, at least 5% of the total outsourced samples shall be re-tested and compared for consistency and data integrity.

16.0 Practices, Methods and Forms

The following is a list of general practices/administrative procedures, analytical methods and forms utilized in Forensic Biology. Each follows the numbering scheme of: Biology Screening (1XX), DNA Casework Analysis (2XX), CODIS/Database Analysis (3XX) and QC Functions (4XX).

MBI≡Schemes, generally encompassing many procedures.

- MBI-100** EXAMINATION OF BLOODSTAINED EVIDENCE
- MBI-102** EXAMINATION OF EVIDENCE FOR SEMEN
- MBI-104** EXAMINATION OF EVIDENCE FOR BODY FLUIDS
- MBI-200** INDIVIDUALIZATION OF DNA SOURCES BY STR ANALYSIS
- MBI-300** INDIVIDUALIZATION OF DNA SOURCES BY STR ANALYSIS

BI≡Analytical Procedures or Individual Processes

- BI-100** PROCESSING LIQUID BLOOD
- BI-102** DNA PACKETS
- BI-104** PHENOLPHTHALEIN TEST FOR BLOOD
- BI-105** O-TOLIDINE TEST FOR BLOOD
- BI-106** HUMAN BLOOD IDENTIFICATION USING ABACARD[®] HEMATRACE[®] TEST
- BI-108** SPECIES IDENTIFICATION: OUCHTERLONY DOUBLE DIFFUSION
- BI-110** BIOLOGICAL SCREENING: USE OF ALTERNATE LIGHT SOURCE
- BI-111** BIOLOGICAL SCREENING: USE OF INFRA RED LIGHT
- BI-114** BRENTAMINE TEST FOR ACID PHOSPHATASE
- BI-116** SAMPLE EXTRACTION FOR SEMEN IDENTIFICATION
- BI-118** SEMEN IDENTIFICATION: MICROSCOPIC EXAMINATION
- BI-119** SPERM DOCUMENTATION: MIDEO SYSTEM
- BI-120** IDENTIFICATION OF SEMEN BY P30 DETECTION (ABAcad[®])
- BI-122** AMYLASE TEST: PHADEBAS
- BI-124** AMYLASE TEST: STARCH IODIDE
- BI-126** DETECTION OF URINE (UREASE)
- BI-128** DETECTION OF URINE (CREATININE)
- BI-130** DETECTION OF FECAL MATERIAL (UROBILINOGEN)
- BI-200** EXTRACTION PROTOCOLS FOR PCR DNA TYPING TESTS
- BI-207** DNA QUANTIFICATION: REAL-TIME PCR
- BI-208** STR AMPLIFICATION: PP16
- BI-210** STR TYPING: CAPILLARY ELECTROPHORESIS AND DATA ANALYSIS
- BI-301** CODIS SAMPLE RECEIPT AND DNA TRACKER ENTRY
- BI-302** CODIS SAMPLE DATA ENTRY AND UPLOAD
- BI-303** CODIS DATABASE HIT VERIFICATION
- BI-310** CODIS SAMPLE REMOVAL

- BI-312 EXTRACTION PROTOCOLS FOR PCR DNA TYPING TESTS
- BI-314 DNA QUANTIFICATION : REAL-TIME PCR
- BI-316 STR AMPLIFICATION : PP16HS
- BI-318 STR TYPING : CAPILLARY ELECTROPHORESIS AND DATA ANALYSIS
- BI-400 DRIFTCON FFC : TEMPERATURE VERIFICATION

Form BI≡Various forms used in each discipline
 * indicates a controlled form

- 100-BI PHENOLPHTHALEIN REAGENT (KASTLE-MEYER)
- 102-BI HYDROGEN PEROXIDE 3% (v/v)
- 103-BI O-TOLIDINE REAGENT
- 104-BI AMMONIUM HYDROXIDE (~3%)
- 108-BI OUCHTERLONY DESTAIN
- 110-BI OUCHTERLONY STAIN
- 114-BI 10X BRENTAMINE (SODIUM ACETATE) BUFFER
- 116-BI BRENTAMINE SOLUTION A
- 118-BI BRENTAMINE SOLUTION B
- 120-BI SALINE (0.85% NaCl)
- 124 BI 1X PHOSPHATE BUFFERED SALINE (PBS)
- 126-BI XMAS TREE STAIN SOLUTION A (KERNECHTROT SOLUTION)
- 128-BI XMAS TREE STAIN SOLUTION B (PICROINDIGOCARMIN SOLUTION)
- 132-BI AMYLASE DIFFUSION BUFFER (pH6.9)
- 134-BI AMYLASE IODINE REAGENT
- 138-BI MERCURIC CHLORIDE 10% (w/v)
- 140-BI ZINC CHLORIDE 10% (w/v)
- 201-BI 1M TRIS-HCl BUFFER pH7.5
- 203-BI 1M TRIS-HCl BUFFER pH8
- 205-BI ETHYLENEDIAMINE TETRAACETIC ACID (EDTA) 0.5M
- 207-BI STAIN EXTRACTION BUFFER pH8
- 211-BI BROTEINASE K (20 mg/ml)
- 222-BI 1M SODIUM ACETATE pH5.2
- 223-BI DTT (1M)
- 229-BI PCR-TE (TE⁻⁴) BUFFER (10mM TRIS-HCl, 0.1M EDTA)
- 231-BI NaOH 5N
- 233-BI SODIUM CHLORIDE (NaCl) 5M
- 249-BI BOVINE SERUM ALBUMIN (BSA) 4%
- 101-BI BIOLOGY SCREENING SUMMARY
- 200-BI DNA EXTRACTION WORKSHEET
- 202-BI DIFFERENTIAL DNA EXTRACTION WORKSHEET
- 206-BI* 7500 LOAD SHEET
- 209-BI* 7500 RESULTS SHEET
- 210-BI STR AMPLIFICATION SET-UP
- 212-BI STR EXTRACTION CONTROL GENOTYPE CHECK
- 214-BI STR TECHNICAL REVIEW CHECKLIST
- 216-BI* 3130 LOAD SHEET

218-BI CODIS ENTRY FORM
306-BI STR OFFENDER DATABASE REVIEW CHECKLIST
310-BI CODIS SAMPLE REMOVAL CHECKLIST
312-BI* DATABASE WORKSHEETS (A-E)
314-BI OUTSOURCED OFFENDER DATA REVIEW
316-BI DATABASE REINJECTION SUMMARY
400-QC FORENSIC BIOLOGY CHEMICAL INVENTORY
401-QC FORENSIC BIOLOGY CRITICAL EQUIPMENT INVENTORY
402-QC FORENSIC BIOLOGY EQUIPMENT MAINTENANCE/REPAIR RECORD
403-QC* FORENSIC BIOLOGY pH CALIBRATION RECORD
404A-QC* BIOLOGY/DNA CASEWORK WEEKLY QC
404B-QC* EVIDENCE VAULT WEEKLY QC
404C-QC* DNA DATABASE WEEKLY QC
406A-QC* BIOLOGY/DNA CASEWORK MONTHLY QC
406B-QC* FORENSIC BIOLOGY MONTHLY QC
406C-QC* DNA DATABASE MONTHLY QC
408-QC FORENSIC BIOLOGY QUARTERLY QC
410-QC* QC ABACARD® HEMATRACE® KIT
412-QC* QC ONESTEP ABACARD® P30 KIT
419-QC* QC QUANTIFILER® HUMAN DNA QUANTIFICATION KIT
420-QC* QC STR KITS
422A-QC 310 INJECTION LOG
422B-QC 3130/3130x1 INJECTION LOG
426-QC* ANNUAL NIST QC RUN
428-QC BIOMEK 3000 QC

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT



DNA Quantitation

7500 Load Sheet

Form 206-BI

Case Number: _____

Analyst: _____

Plate Name: _____

Date: _____

	1	2	3	4	5	6	7	8	9	10	11	12
A	STD. 1	STD. 1										
B	STD. 2	STD. 2										
C	STD. 3	STD. 3										
D	STD. 4	STD. 4										
E	STD. 5	STD. 5										
F	STD. 6	STD. 6										
G	STD. 7	STD. 7										
H	STD. 8	STD. 8										

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

Quantifier Kit

lot #: _____

Std. Prep. Date: _____

expiry date: _____

TE lot#: _____

Master Mix made for: _____

total samples: _____

reaction mix 0 ul

primer mix 0 ul



DNA Quantitation

7500 Results Sheet

Form 209-BI

Case Number: _____

Analyst: _____

Plate Name: _____

Date: _____

Well	Sample Name	IPC C _T	Quantity ng/ul	ul Sample for Dilution	ul TE to be added	ng/ul Final	ul to be Amplified
A3		0	0	0	5	0.0	0.1
B3		0	0	0	5	0.0	0.1
C3		0	0	0	5	0.0	0.1
D3		0	0	0	5	0.0	0.1
E3		0	0	0	5	0.0	0.1
F3		0	0	0	5	0.0	0.1
G3		0	0	0	5	0.0	0.1
H3		0	0	0	5	0.0	0.1
A4		0	0	0	5	0.0	0.1
B4		0	0	0	5	0.0	0.1
C4		0	0	0	5	0.0	0.1
D4		0	0	0	5	0.0	0.1
E4		0	0	0	5	0.0	0.1
F4		0	0	0	5	0.0	0.1
G4		0	0	0	5	0.0	0.1
H4		0	0	0	5	0.0	0.1
A5		0	0	0	5	0.0	0.1
B5		0	0	0	5	0.0	0.1
C5		0	0	0	5	0.0	0.1
D5		0	0	0	5	0.0	0.1
E5		0	0	0	5	0.0	0.1
F5		0	0	0	5	0.0	0.1
G5		0	0	0	5	0.0	0.1
H5		0	0	0	5	0.0	0.1
A6		0	0	0	5	0.0	0.1
B6		0	0	0	5	0.0	0.1
C6		0	0	0	5	0.0	0.1
D6		0	0	0	5	0.0	0.1
E6		0	0	0	5	0.0	0.1
F6		0	0	0	5	0.0	0.1
G6		0	0	0	5	0.0	0.1
H6		0	0	0	5	0.0	0.1
A7		0	0	0	5	0.0	0.1
B7		0	0	0	5	0.0	0.1
C7		0	0	0	5	0.0	0.1
D7		0	0	0	5	0.0	0.1
E7		0	0	0	5	0.0	0.1
F7		0	0	0	5	0.0	0.1
G7		0	0	0	5	0.0	0.1
H7		0	0	0	5	0.0	0.1
A8		0	0	0	5	0.0	0.1
B8		0	0	0	5	0.0	0.1
C8		0	0	0	5	0.0	0.1
D8		0	0	0	5	0.0	0.1
E8		0	0	0	5	0.0	0.1
F8		0	0	0	5	0.0	0.1
G8		0	0	0	5	0.0	0.1

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

Well	Sample Name	IPC C _T	Quantity ng/ul	ul Sample for Dilution	ul TE to be added	ng/ul Final	ul to be Amplified
H8		0	0	5	0.0	0.1	10.0
A9		0	0	5	0.0	0.1	10.0
B9		0	0	5	0.0	0.1	10.0
C9		0	0	5	0.0	0.1	10.0
D9		0	0	5	0.0	0.1	10.0
E9		0	0	5	0.0	0.1	10.0
F9		0	0	5	0.0	0.1	10.0
G9		0	0	5	0.0	0.1	10.0
H9		0	0	5	0.0	0.1	10.0
A10		0	0	5	0.0	0.1	10.0
B10		0	0	5	0.0	0.1	10.0
C10		0	0	5	0.0	0.1	10.0
D10		0	0	5	0.0	0.1	10.0
E10		0	0	5	0.0	0.1	10.0
F10		0	0	5	0.0	0.1	10.0
G10		0	0	5	0.0	0.1	10.0
H10		0	0	5	0.0	0.1	10.0
A11		0	0	5	0.0	0.1	10.0
B11		0	0	5	0.0	0.1	10.0
C11		0	0	5	0.0	0.1	10.0
D11		0	0	5	0.0	0.1	10.0
E11		0	0	5	0.0	0.1	10.0
F11		0	0	5	0.0	0.1	10.0
G11		0	0	5	0.0	0.1	10.0
H11		0	0	5	0.0	0.1	10.0
A12		0	0	5	0.0	0.1	10.0
B12		0	0	5	0.0	0.1	10.0
C12		0	0	5	0.0	0.1	10.0
D12		0	0	5	0.0	0.1	10.0
E12		0	0	5	0.0	0.1	10.0
F12		0	0	5	0.0	0.1	10.0
G12		0	0	5	0.0	0.1	10.0
H12		0	0	5	0.0	0.1	10.0

Property of Idaho State Police Forensic Services
 Uncontrolled Internet Copy
 OBSOLETE DOCUMENT

3130 Load Sheet

Form 216-BI



Case Number: _____

Analyst: _____

Plate Name: _____

Date: _____

	1	2	3	4	5	6	7	8	9	10	11	12
A												
B												
C												
D												
E												
F												
G												
H												

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

Master Mix made for: _____

total samples: _____

3130 POP4

HiDi Formamide 0 ul

HiDi Formamide

HiDi Formamide Lot# _____

Lot# _____

Internal Lane Standard 0 ul

Internal Lane Standard

Buffer Lot# _____

Expiration Date _____

Revision 11

BSD/Biomex 3000 Load Sheet

Form 312A-BI



Analyst: _____

Plate Name: _____

Date: _____

Test Name: _____

	1	2	3	4	5	6	7	8	9	10	11	12
A												
B												
C												
D												
E												
F												
G												
H												

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

Punch Size: _____
Number of Punches: _____

Elution Volume (µl): _____
Biomek Program: _____

Lot #: _____
Exp. Date: _____

DNA IQ Kit

DTT Lot# / Date Added: _____
Ethanol Lot# / Date Added: _____
2-Propanol Lot# / Date Added: _____

Revision 11
8/29/2011
Issuing Authority: Quality Manager

DNA Quantitation

7500 Load Sheet

Form 312B-BI



Plate Name: _____

Analyst: 0

Date: _____

7500: _____

	1	2	3	4	5	6	7	8	9	10	11	12
A	0	0	0	0	0	0	0	0	0	0	STD. 1	STD. 1
B	0	0	0	0	0	0	0	0	0	0	STD. 2	STD. 2
C	0	0	0	0	0	0	0	0	0	0	STD. 3	STD. 3
D	0	0	0	0	0	0	0	0	0	0	STD. 4	STD. 4
E	0	0	0	0	0	0	0	0	0	0	STD. 5	STD. 5
F	0	0	0	0	0	0	0	0	0	0	STD. 6	STD. 6
G	0	0	0	0	0	0	0	0	0	0	STD. 7	STD. 7
H	0	0	0	0	0	0	0	0	0	0	STD. 8	STD. 8

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

Master Mix made for: _____

Total Samples: _____

Quantifier Kit

reaction mix 0 ul

Std. Prep. Date: _____


Lot #: _____

primer mix 0 ul

TE lot#: _____

Exp. Date: _____

Revision 11
8/29/2011
Issuing Authority: Quality Manager



DNA Quantitation

7500 Results Sheet

Form 312C-BI

Plate Name: 0

Analyst: 0

Date: 1/0/1900

Well	Sample Name	IPC C _T	Quantity ng/ul	ul Sample for Dilution	ul TE to be added	ng/ul Final	ul to be Amplified
A1		0	0.00	5	0.0	0.1	2.5
B1		0	0.00	5	0.0	0.1	2.5
C1		0	0.00	5	0.0	0.1	2.5
D1		0	0.00	5	0.0	0.1	2.5
E1		0	0.00	5	0.0	0.1	2.5
F1		0	0.00	5	0.0	0.1	2.5
G1		0	0.00	5	0.0	0.1	2.5
H1		0	0.00	5	0.0	0.1	2.5
A2		0	0.00	5	0.0	0.1	2.5
B2		0	0.00	5	0.0	0.1	2.5
C2		0	0.00	5	0.0	0.1	2.5
D2		0	0.00	5	0.0	0.1	2.5
E2		0	0.00	5	0.0	0.1	2.5
F2		0	0.00	5	0.0	0.1	2.5
G2		0	0.00	5	0.0	0.1	2.5
H2		0	0.00	5	0.0	0.1	2.5
A3		0	0.00	5	0.0	0.1	2.5
B3		0	0.00	5	0.0	0.1	2.5
C3		0	0.00	5	0.0	0.1	2.5
D3		0	0.00	5	0.0	0.1	2.5
E3		0	0.00	5	0.0	0.1	2.5
F3		0	0.00	5	0.0	0.1	2.5
G3		0	0.00	5	0.0	0.1	2.5
H3		0	0.00	5	0.0	0.1	2.5
A4		0	0.00	5	0.0	0.1	2.5
B4		0	0.00	5	0.0	0.1	2.5
C4		0	0.00	5	0.0	0.1	2.5
D4		0	0.00	5	0.0	0.1	2.5
E4		0	0.00	5	0.0	0.1	2.5
F4		0	0.00	5	0.0	0.1	2.5
G4		0	0.00	5	0.0	0.1	2.5
H4		0	0.00	5	0.0	0.1	2.5
A5		0	0.00	5	0.0	0.1	2.5
B5		0	0.00	5	0.0	0.1	2.5
C5		0	0.00	5	0.0	0.1	2.5
D5		0	0.00	5	0.0	0.1	2.5
E5		0	0.00	5	0.0	0.1	2.5
F5		0	0.00	5	0.0	0.1	2.5
G5		0	0.00	5	0.0	0.1	2.5
H5		0	0.00	5	0.0	0.1	2.5
A6		0	0.00	5	0.0	0.1	2.5
B6		0	0.00	5	0.0	0.1	2.5
C6		0	0.00	5	0.0	0.1	2.5
D6		0	0.00	5	0.0	0.1	2.5
E6		0	0.00	5	0.0	0.1	2.5
F6		0	0.00	5	0.0	0.1	2.5
G6		0	0.00	5	0.0	0.1	2.5
H6		0	0.00	5	0.0	0.1	2.5



DNA Quantitation

7500 Results Sheet

Form 312C-BI

Plate Name: 0

Analyst: 0

Date: 1/0/1900

Well	Sample Name	IPC C _T	Quantity ng/ul	ul Sample for Dilution	ul TE to be added	ng/ul Final	ul to be Amplified
A7		0	0.00	5	0.0	0.1	2.5
B7		0	0.00	5	0.0	0.1	2.5
C7		0	0.00	5	0.0	0.1	2.5
D7		0	0.00	5	0.0	0.1	2.5
E7		0	0.00	5	0.0	0.1	2.5
F7		0	0.00	5	0.0	0.1	2.5
G7		0	0.00	5	0.0	0.1	2.5
H7		0	0.00	5	0.0	0.1	2.5
A8		0	0.00	5	0.0	0.1	2.5
B8		0	0.00	5	0.0	0.1	2.5
C8		0	0.00	5	0.0	0.1	2.5
D8		0	0.00	5	0.0	0.1	2.5
E8		0	0.00	5	0.0	0.1	2.5
F8		0	0.00	5	0.0	0.1	2.5
G8		0	0.00	5	0.0	0.1	2.5
H8		0	0.00	5	0.0	0.1	2.5
A9		0	0.00	5	0.0	0.1	2.5
B9		0	0.00	5	0.0	0.1	2.5
C9		0	0.00	5	0.0	0.1	2.5
D9		0	0.00	5	0.0	0.1	2.5
E9		0	0.00	5	0.0	0.1	2.5
F9		0	0.00	5	0.0	0.1	2.5
G9		0	0.00	5	0.0	0.1	2.5
H9		0	0.00	5	0.0	0.1	2.5
A10		0	0.00	5	0.0	0.1	2.5
B10		0	0.00	5	0.0	0.1	2.5
C10		0	0.00	5	0.0	0.1	2.5
D10		0	0.00	5	0.0	0.1	2.5
E10		0	0.00	5	0.0	0.1	2.5
F10		0	0.00	5	0.0	0.1	2.5
G10		0	0.00	5	0.0	0.1	2.5
H10		0	0.00	5	0.0	0.1	2.5

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT



STR Amplification Set-Up

Form 312D-BI

Plate Name: _____

Analyst: _____ 0

Date: _____

Thermal Cycler: _____

	1	2	3	4	5	6	7	8	9	10	11	12
A	0	0	0	0	0	0	0	0	0	0	0	0
B	0	0	0	0	0	0	0	0	0	0	0	0
C	0	0	0	0	0	0	0	0	0	0	0	0
D	0	0	0	0	0	0	0	0	0	0	0	0
E	0	0	0	0	0	0	0	0	0	0	0	0
F	0	0	0	0	0	0	0	0	0	0	0	0
G	0	0	0	0	0	0	0	0	0	0	0	0
H	0	0	0	0	0	0	0	0	0	0	0	0

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

STR Kit Type: _____

STR Kit Lot#: _____

PCR TE Lot#: _____

Total Samples: _____

Master Mix Made For: _____

Reagent	μl/sample	μl in Master Mix
5x Master Mix	1.25	0
Primers	0.6	0
H ₂ O		0
Master Mix/Sample	1.85	
DNA Template	4.4	
Total Rxn Volume	6.25	

3130XL Load Sheet

Form 312E-BI



Plate Name: _____ Analyst: 0

Date: _____

	1	2	3	4	5	6	7	8	9	10	11	12
A	0	0	0	0	0	0	0	0	0	0	0	0
B	0	0	0	0	0	0	0	0	0	0	0	0
C	0	0	0	0	0	0	0	0	0	0	0	0
D	0	0	0	0	0	0	0	0	0	0	0	0
E	0	0	0	0	0	0	0	0	0	0	0	0
F	0	0	0	0	0	0	0	0	0	0	0	0
G	0	0	0	0	0	0	0	0	0	0	0	0
H	0	0	0	0	0	0	0	0	0	0	0	0

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

Master Mix made for: _____

Total Samples: _____

HiDI Formamide 0 ul
Internal Lane Standard 0 ul

3130 POP4
HiDI Formamide Lot#: _____
Buffer Lot#: _____

Lot#: _____
Expiration Date: _____



FORENSIC BIOLOGY pH CALIBRATION RECORD
(Oakton pH meter, serial #135212)

DATE	INITIALS	STANDARD BUFFER pH 4.01 Reading/lot #	STANDARD BUFFER pH 7.00 Reading/lot #	STANDARD BUFFER pH 10.01 Reading/lot #

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

A 3-point calibration of the pH meter will be performed at the time of use (See the Oakton Operating Manual for calibration and pH measurement instructions). The analyst will record the date of calibration, their initials, the measured pH value and lot # for each buffer. The measured reading must fall within ± 0.50 pH for the calibration to be confirmed by the meter.

BIOLOGY/DNA CASEWORK WEEKLY QC

DATE/INITIALS	°C		°C		°C		°C		°C	
	Min	Max	Min	Max	Min	Max	Min	Max	Min	Max
COMBO F/F CW1										
FRIDGE CW2										
FREEZER CW3										
FREEZER CW4										
FREEZER CW5										
COMBO F/F CW6										
COMBO F/F RP1										

Observed temperatures should fall between $\pm 5^{\circ}\text{C}$ of the target temperature (4°C for refrigerators and -20°C for freezers). The temperature control should be adjusted to correct for minor variations; however, if the temperature is not corrected or if it falls significantly outside the target range, it should be taken out of service and maintenance/repair performed as needed. Note: frost-free freezers will have a greater temperature range ($\pm 10^{\circ}\text{C}$) due to the heating and cooling cycles. Combination fridge/freezers with a single temperature control may not be able to maintain both units within the given range. In this case, temperature sensitive reagents should be stored in appropriate containers (such as cryo-boxes) to maintain the desired state.

DATE/INITIALS					
NANOPIRE SYSTEM					
°C OVEN A	set/observed				
°C OVEN B	set/observed				
°C HEAT BLOCK A	set/observed				
°C HEAT BLOCK B	set/observed				
°C HEAT BLOCK C (prod. rm.)	set/observed				
°C HEAT BLOCK D	set/observed				

Observed temperatures should fall between $\pm 2^{\circ}\text{C}$ of the temperature set point. The temperature control should be adjusted to correct for minor variations; however, if the temperature is not corrected or if it falls significantly outside the target range, it should be taken out of service and maintenance/repair performed as needed.

The observed water purity for the Nanopure system should be a minimum of 18.0 mega-ohms. If the purity falls below this point, the cartridges should be changed and the system sanitized as necessary.

MERIDIAN EVIDENCE VAULT WEEKLY QC

DATE/INITIALS	°C		°C		°C		°C	
	Min	Max	Min	Max	Min	Max	Min	Max
FRIDGE VR1								
FRIDGE VR2								
FREEZER VF1								
FREEZER VF2								
FREEZER VF3								
FREEZER VF4								
FREEZER VF5								
FREEZER VF6								
FREEZER VF7								

Observed temperatures for refrigerators should fall between $\pm 5^{\circ}\text{C}$ of the 4°C target temperature. Freezers should fall between $\pm 10^{\circ}\text{C}$ of the target -20°C . The temperature control should be adjusted to correct for minor variations; however, if the temperature is not corrected or if it falls significantly outside the target range, it should be taken out of service and maintenance/repair performed as needed. Note: frost-free freezers will have a greater temperature range ($\pm 15^{\circ}\text{C}$) due to the heating and cooling cycles.

DNA DATABASE WEEKLY QC

DATE/INITIALS	°C		°C		°C		°C	
	Min	Max	Min	Max	Min	Max	Min	Max
FRIDGE DB1								
COMBO FF DB2								
FREEZER DB3								

Observed temperatures should fall between $\pm 5^{\circ}\text{C}$ of the target temperature (4°C for refrigerators and -20°C for freezers). The temperature control should be adjusted to correct for minor variations; however, if the temperature is not corrected or if it falls significantly outside the target range, it should be taken out of service and maintenance/repair performed as needed. Note: frost-free freezers will have a greater temperature range ($\pm 10^{\circ}\text{C}$) due to the heating and cooling cycles. Combination fridge/freezers with a single temperature control may not be able to maintain both units within the given range. In this case, temperature sensitive reagents should be stored in appropriate containers (such as cryo-boxes) to maintain the desired state.

DATE/INITIALS	
°C WATER BATH	set/observed

Observed temperatures should fall between $\pm 2^{\circ}\text{C}$ of the temperature set point. The temperature control should be adjusted to correct for minor variations; however, if the temperature is not corrected or if it falls significantly outside the target range, it should be taken out of service and maintenance/repair performed as needed.

BIOLOGY/DNA CASEWORK MONTHLY QC

AUTOCLAVE											
DATE/INITIALS											
CLEAN											
STERILIZATION (+)											
STERILIZATION (-)											
LABORATORY AND OTHER EQUIPMENT											
BIROBOT EZ1s											
GREASE D-RINGS											
CLEAN CENTRIFUGES											
CLEAN PIPETS											
LAB CLEANED											
EYEWASH CHECK											

*Personnel should initial the duties they perform and date separately, if necessary.

Autoclave sterilization is checked by the observation of microbial growth in the (+) control (non-sterilized) and a lack of growth in the (-) control (sterilized) samples. See the BTSure product insert for test instructions and growth indicators. If sterilization is not achieved, the autoclave should be serviced.

FORENSIC BIOLOGY MONTHLY QC

DATE/INITIALS																				
ABI 7500 Instrument Maintenance																				
Background Assay/Contamination Check																				
System Function Test																				
Lamp Status Check																				
7500/310 Computer Maintenance																				
Disk Cleanup																				
Defragment Hard Drive																				
3130/3130XL Maintenance																				
Water Seal Trap Flush																				
Water Wash Wizard																				
Defragment Hard Drive (C & E)																				

*Personnel should initial the duties they perform and date separately, if necessary.

Note: See the ABI 7500 Maintenance Guide and/or the April 2007 User Bulletin for additional Instrument Maintenance procedures and pass/fail criteria.

A contamination check will be performed by the background assay. If outliers are observed during the Background Assay (Intensity value $\geq 72,000$), or fluorescence (red) observed during the block check, the specific well should be identified and cleaned. Rerun the background calibration after wells have been cleaned sufficiently. Note: a 96-well tray with 50ul TE in each well may be used as a background tray.

A block check may be performed prior to the background assay by selecting Instrument > Calibrate. Set the exposure time to 2048ms, lamp control to idle, select Filter A and click Snapshot. Holding the cursor over the fluorescence will give pixel intensity.


If a component fails the function test, a service call should be placed.

If the lamp fails the function test and/or status check, it should be replaced, followed by calibration of ROI, background, optical, pure dye, and instrument verification in that order. Reset the Lamp Timer when complete.

The 7500 and 310 Disk Cleanup is performed by selecting Start Menu > Programs > Accessories > System Tools > Disk Cleanup.

7500 and 310 Defragmentation is performed by selecting Start Menu > Programs > Accessories > System Tools > Disk Defragmenter.

3130/3130XL Defragmentation is performed by right-clicking on 'My Computer' and selecting 'Manage'. In the tree tab choose Computer Management (local)>Disk Fragmenter>Drive name>Defragment.




FORM 406C-QC

DNA DATABASE MONTHLY QC

DATE/INITIALS										
CLEAN CENTRIFUGES										
CLEAN PIPETS										
CLEAN BSD PUNCHER										
CLEAN BIOMEK 3000										
LAB CLEANED										
EYEWASH CHECK										

*Personnel should initial the duties they perform and date separately, if necessary.

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT



QC ABACARD® HEMATRACE® KIT

HEMATRACE® KIT LOT: _____

DATE RECEIVED: _____

SCIENTIST: _____

QC DATE: _____


Perform test as usual with one 2mm² cutting and one 2mm thread from known bloodstain. Record results (include time it took for positive rxn to be visible). If available, attach photo documentation and place in Biology Casework QC binder.

SAMPLE	RXN	TIME (min. sec.)
2mm ² cutting		
2mm thread		
Neg		

The 2mm² cutting sample must have a positive reaction within 10 minutes for passing. The 2mm thread should ideally be positive within 10 minutes but is used primarily as a sensitivity indicator of the given test lot. The kit may still be deemed as passing without a positive result for the thread.

QA/QC PASSED: YES NO

Comments:



QC OneStep ABACARD® p30 KIT

ABACARD® p30 KIT LOT: _____ DATE RECEIVED: _____

SCIENTIST: _____ QC DATE: _____

Perform test as usual with a known semen extract, as well as ~10ng/ml (10µl of a 1:500 dilution) and ~50ng/ml (10µl of a 1:100 dilution) of Seri Semen Standard. Record results (include time it took for positive rxn to be visible). If available, attach photo documentation and place in Biology Casework QC binder.


SAMPLE	RXN	TIME (min. sec.)
Semen Extract		
10ng/ml		
50ng/ml		
Neg		
*250ng/ml or 1:10		

The semen extract must have a positive rxn within 10 minutes for passing. The Seri standards are used to estimate the range of sensitivity of the kit lot.

For the semen standard dilutions, if a positive rxn is not obtained at 10 minutes, continue to monitor and record result at the end of 15 minutes. In addition, *run a 250ng/ml (50µl of the 1:100 dilution to 150µl of extraction buffer) or a 1:10 dilution of the semen stain extract to ensure the kit is operating within reasonable limits for forensic identification. In addition to the neat semen extract, this control sample (250ng/ml or 1:10 extract) must result in a positive rxn within 10 minutes.

QA/QC PASSED: YES NO

Comments:



QC QUANTIFILER HUMAN KITS

KIT LOT #: _____

DATE RECEIVED: _____

EXPIRATION DATE: _____

SCIENTIST: _____

QA/QC DATE: _____

KIT COMPONENT	LOT NUMBER
PRIMER MIX	
REACTION MIX	
DNA STANDARD	

To check the new kit lot, perform quantification as usual. For samples, run standards from the new kit to be QC'd and equivalent dilutions of the NIST SRM 2372 Quant Standard, as well as 0.5ng and 2ng of 9947A DNA. Analyze using the SRM as standard and the new kit as unknown. Using an average of the results for the new kit standards, calculate the new volume of TE to be added in the preparation of standard 1, per the equation $C_1V_1=C_2V_2$ (where C=average for std 1, and V=total volume). Record the slope obtained for the standard curve.

As a check of the calculation and resulting TE volume, use the new kit, with corresponding new dilution to perform a 9947A DNA quantification, according to standard procedure. Use 0.5ng and 2ng and compare the results to those obtained from above. A pass will be achieved if the slopes for both standard curves are comparable.

SRM 2372 component used: _____
 Standard curve slope: _____
 Volume TE to be used for Standard 1: _____

QA/QC PASSED: YES NO

Comments:

Attach the 7500 Load Sheets, Standard Curves, and Results Sheets. Record the calculations in the documentation. Mark the new kit with TE volume for Standard 1 preparation.



QC STR KITS

STR KIT: _____ DATE RECEIVED: _____

KIT MANUFACTURER: _____ KIT LOT #: _____

LAB LOT#: _____ SCIENTIST: _____ QA/QC DATE: _____

PP16 KIT COMPONENT	LOT NUMBER	PP16HS KT COMPONENT	LOT NUMBER
PRIMER MIX		10X PRIMER PAIR MIX	
REACTION MIX		5X MASTER MIX	
CONTROL DNA		CONTROL DNA	
TAQ GOLD*		INTERNAL LANE STANDARD	
INTERNAL LANE STANDARD		ALLELIC LADDER	
ALLELIC LADDER			

An Extraction Control, reagent blank and associated controls are to be processed from extraction. A pass will be achieved by obtaining the expected results for each of the samples run and data of acceptable quality (e.g. sufficient RFUs). Comments regarding quality concerns are to be noted as appropriate.

* Taq Gold is purchased separately for PowerPlex 16, but typically at the same time as a new STR kit. If Taq Gold is received separate from an STR kit, the QC of the Taq will be noted on this form (corresponding to the appropriate STR kit lot#) under comments.

The ILS, once QC'd, may be used interchangeably between the PP16 and PP16HS kits, if necessary, as they are the same product. The ladders are not the same product and cannot be switched between kit types.

Run Date: _____ Run Folder: _____

QA/QC PASSED: YES NO

Comments:

Attach the appropriate extraction/amplification/Extraction Control forms used and the GeneMapper ID Electropherograms; place in the appropriate QC Binder.



ANNUAL NIST QC RUN

SCIENTIST: _____

QC DATE: _____

At a minimum of once a year, an 'in-date', certified NIST-SRM standard will be analyzed with our standard procedures. Control or known reference samples may be analyzed simultaneously to 'certify' them for use as NIST QC samples. These samples will be listed in the comments section of this form with lot # and that they were certified. After completion of the QC, the newly 'certified' samples, or their container, will be marked as "NIST Certified" with the corresponding date.

The GeneMapper® ID Data will be analyzed as usual and quality of results will be reflected in the comments or 'passed' areas as appropriate and necessary. Passing results are obtained by achieving the expected results for the given NIST sample(s) and any associated controls. The GeneMapper® ID Electropherograms and an Allele Table will be printed [for the NIST sample(s)] and stored in the section QC binder.

Run Folder: _____

QC PASSED: YES NO

Comments:

310 INJECTION LOG

DATE	ANALYST	FORMAMIDE (LOT#)	H ₂ O CHANGED	BUFFER CHANGED ✓ (LOT#)	310 POP-4 CHANGED ✓ (LOT#/EXCP)	# INJECTIONS RUN/TOTAL	CAPILLARY CHANGED ✓ (LOT#)	CASES RUN

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

3130/3130XL INJECTION LOG

DATE	ANALYST	H ₂ O WASH ✓	H ₂ O TRAP FLUSH ✓	H ₂ O CHANGED (✓)	BUFFER CHANGED ✓ (LOT#)	3130 POP-4 CHANGED ✓ (LOT#/EXP)	CAPILLARY ARRAY CHANGED ✓ (LOT#)	#INJECTIONS RUN/TOTAL	CASES RUN

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

1 Biology Section, Digital Work Place Training Plan

- 1.1 All analysts utilizing imaging technologies shall be trained and tested for competency in the standard operating procedures and the operation of the relevant imaging technologies.
- 1.2 Formal training may be modified at the discretion of the Latent Section Supervisor dependent upon previous training and/or experience.
- 1.3 Recommended formal training consists of:
 - 1.3.1 Reviewing the ISP-FS Latent Print Section Digital Imaging Procedure.
 - 1.3.2 Reviewing the Digital Workplace User manual or equivalent.
 - 1.3.3 Review of relevant chapters of the Adobe Photoshop Users Manual and/or completion of a digital imaging course that utilizes Adobe Photoshop.
 - 1.3.4 Satisfactory creation and digital processing of a mock-case using Digital Workplace and Adobe Photoshop software or equivalent software.
 - 1.3.6 Satisfactory completion of a written test.
- 1.4 Competency testing shall be repeated when significant changes in hardware or software are made (e.g. manufacturer/vendor changes).

2 Training record, Required Reading	Trainee / Completion Date
2.1 Foray ADAMS User Manual Forensic Image Tracking System Pages 1- 61	_____ / _____
2.2 Latent Print Section AM Section 11.	_____ / _____
2.3 Review Current Adobe Photoshop Users Manual.	_____ / _____
2.4 Practical Exercises	Supervisor / Date / P or F
2.4.1 Digital Acquisition form folder	_____ / _____ / _____
2.4.2 Digital Image Enhancement	_____ / _____ / _____
2.5.2 Competency test on Digital Imaging System. The analyst will independently capture, calibrate, enhance, and document _____ digital images.	_____ / _____ / _____

Idaho State Police Forensic Services

Approval for Quality System Controlled Documents




Discipline/Name of Document: Biology CODIS and Database AM Manual

Revision Number: 11

Issue Date: 08/29/2011


APPROVED BY:


Quality Manager

8/29/11
Date Signed


Checklist Submitted and Checked WAG

ISP FORENSIC BIOLOGY
CODIS AND DATABASE ANALYTICAL
METHODS MANUAL




Forensic Biology CODIS and Database Analytical Methods Manual

Revision #11



APPROVED
August 29, 2011



Cynthia Cunningham
Forensic Scientist IV
Forensic Biology Supervisor/Technical Leader

**Forensic Biology CODIS and Database Analytical Methods
Manual**

Table of Contents

- i. Signature/Approval page
- ii. Revision History
- iii. Table of Contents
- 1-2 Individualization by STR Analysis
- 3-9 Offender Sample Receipt and Tracker Entry
- 10-14 CODIS Sample Data Entry and Upload
- 15-18 CODIS Database Hit Verification
- 19-23 CODIS Sample Removal
- 24-31 Extraction Protocols
- 32-37 DNA Quantification: Real-Time PCR
- 38-42 STR Amplification: PP16 HotStart
- 43-64 STR Typing: Capillary Electrophoresis and Data Analysis
- 65-68 Driftcon FFC: Temperature Verification

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT



INDIVIDUALIZATION OF DNA SOURCES BY STR ANALYSIS

1.0 BACKGROUND:

Once a DNA source has been detected, and identified as to 'source type' where applicable and feasible, it is often important to attribute the DNA sample to a particular individual inasmuch as possible. Current DNA technology, in the analysis of STR loci, offers individualization potential. However, the individualization of a particular sample occurs through a comparative process. This process requires a DNA profile from a 'known' sample to which the evidence sample profile can be evaluated. Often times, this 'known' sample is identified through a search of the CODIS (either SDIS, or NDIS) database in cases where the perpetrator of the particular crime is originally unknown. DNA samples are collected from individuals convicted of qualifying crimes under Idaho code and are then analyzed at STR loci for entry into and subsequent searching against the DNA database.

2.0 SCOPE:


To provide uniform processing of DNA samples to achieve high quality data and consistent interpretation.

3.0 EQUIPMENT/REAGENTS:

As listed in individual analytical procedures.

4.0 PROCEDURE:

See Flow Chart on following page.



5.0 COMMENTS:

5.1 Careful scrutiny at each step will ensure insufficiencies are identified, and compensated for where feasible, at the earliest possible point (see BI-318 for specifics).

**OFFENDER SAMPLE RECEIPT AND DNA TRACKER ENTRY****1.0 BACKGROUND:**

The implementation of the Combined DNA Index System (CODIS) in forensic DNA laboratories has provided an additional tool in assisting law enforcement agencies in solving or linking crimes that otherwise may not have resulted in the identification of a suspect. It is important; however, that samples entered into the database be given a unique identifier, which does not include any personal or identifying information, in order to maintain the confidentiality of the individual. Each laboratory must develop a method of identifier assignment so that each sample may be tracked and identified at a later time, if the need arises.

Idaho Statutes: Title 19, Criminal Procedure, Chapter 55
"The Idaho DNA Database Act of 1996"

ISP Forensic Biology Quality/Procedure Manual, Appendix C

2.0 SCOPE:

To provide a method for tracking offender database samples submitted for STR testing and CODIS entry, while ensuring individual confidentiality.

3.0 EQUIPMENT/REAGENTS:

Computer Workstation with ISP Intranet Access
Barcode Equipment
Court Orders, Database Samples, and Report Forms

4.0 PROCEDURE:**4.1 SAMPLE RECEIPT:**

4.1.1 Offender DNA samples and their corresponding DNA Collection Report Forms received by the

laboratory are to be marked with the date of receipt and the initials of the scientist who received them. The sample and report form may be mailed to the laboratory separately; in the event that a sample has not been received, the submitting agency should be notified.

Note: The Idaho Department of Corrections (IDOC) may submit an offender information printout in lieu of the collection report. The printout will be treated the same as the report throughout the accessioning process.

4.1.2 Where possible, compare the DNA sample card information to that of the Collection Report Form to ensure accuracy. Data for a sample may be entered in absence of a DNA Collection Report form. Additionally, an IDOC# is not necessary for data entry into DNA Tracker, but it is preferred. Contact necessary IDOC personnel or search the corrections website (www.corrections.state.id.us) for a number if one has not been recorded. The DNA Collection Report Form will be retained with the sample after DNA Tracker data entry.

4.1.3 The offense listed on the DNA sample card and/or Collection Report must be a qualifying offense under Idaho Code 19-5506 in order for the sample to be entered into DNA Tracker and subsequently CODIS/NDIS. The corrections website listed above, Idaho Court Repository, and/or ILETS will be consulted for more information if the offense listed is non-qualifying or unclear. The scientist relies, in good faith, on the information provided by the submitting agency for entry into Tracker. It is not necessary, nor is it the scientist's responsibility, to verify the offense(s) of every Offender if the information provided by the submitting agency qualifies as listed. The ISP Attorney General's Office will be consulted on a case by case basis regarding the release of information if a non-qualifying offender is involved in a candidate match.

4.1.4 Samples received in the laboratory that do not have an associated qualifying offense will not be entered into DNA Tracker. The collection report (if present) will be marked to indicate that the sample is a non-qualifying offense and returned to the submitting agency. If there is no collection report, a copy of the sample card will be made, marked in the same manner, and returned to the submitting agency. The sample(s) will then be destroyed.

4.1.5 A thumbprint verification stamp will be placed near the thumbprint on the DNA sample card for all samples to be retained and entered into DNA Tracker (duplicate samples may be excluded). The stamp will be used to record the thumbprint verification for hit confirmations and in order to enter a DNA flag in the offender's criminal history. If no thumbprint is present, the submitting agency will be contacted and a new sample with print will be requested. It is preferable to have a state identification number (SID) to facilitate the thumbprint verification. If no SID is listed on the card or report, the scientist should consult ILETS and record the SID, if found. The latent fingerprint section should be consulted, as necessary, to determine if a new sample and print needs to be requested due to insufficient print quality.

4.2 COURT ORDER RECEIPT:

4.2.1 Court order forms received by the laboratory are to be marked with the date of receipt and the initials of the scientist who received them.

4.2.2 The offense listed on the Court order must be a qualifying offense under Idaho Code 19-5506 in order for the court order information to be entered into DNA Tracker.

4.2.3 If the offense on the court order is not a qualifying offense, no further action will be taken. The court order will be marked to indicate it is a non-qualifying offense and returned to the submitting agency.

4.3 DNA TRACKER PRE-ENTRY SEARCH:

4.3.1 Prior to data entry for any new sample or court order, a database search is performed to eliminate duplicate offender entry. Log on to the DNA Tracker database program, located under Forensics on the ISP Intranet.

4.3.2 A duplicate offender search will be performed using the 'Name' field, followed by at least one of the 'ID' fields to maximize the potential for locating an offender. Note: the 'DOB' field may only be searched in combination with a name.

4.3.3 If all of the searches return 'No matches found', the data for the new sample or court order may be entered as a 'New Offender' (see 4.4).

4.3.4 If a record(s) is returned that meets the criteria, the data is examined and compared with the new sample or court order received. If it is determined that the Offender already exists in DNA Tracker, the new sample or court order and any additional Offender information will be entered under the appropriate tabs for the already existing Offender. Each sample received for an offender will be assigned a unique barcode number and will be retained in the laboratory (see 4.4).

4.4 DNA TRACKER ENTRY:

4.4.1 Enter basic Offender information from the court order, sample card, and/or DNA Collection Report form as follows:

- 1) For Offenders not currently in Tracker, click 'New Offender' at the top of the screen and fill in each of the appropriate fields with the Offender's primary information. Additional information for Offenders already in Tracker may be entered by clicking the 'Edit Basic Details' tab for that Offender.
- 2) Verify all of the information is correct and press the 'Save' button in the top right corner of the screen.
- 3) Add any additional alias names, DOB's, SSN's, and State Identification Numbers on the appropriate alias tab(s) that become available after saving the new offender. The SID verification date will be filled in following thumbprint/SID confirmation. Click the 'Save' button after each entry.

4.4.2 Enter each offense and its associated information from the court order, sample card, and/or DNA Collection Report form as follows:

- 1) Click the 'Offenses' tab for the saved offender.
- 2) Choose the appropriate offense from the pull down menu and enter the corresponding information into the remaining fields.
- 3) Verify the information is correct and press the 'Save' button in the top right corner of the window. Multiple offenses for the same Offender must be entered and saved individually.
- 4) Offense information may be updated/edited if additional data is received at a later time (e.g. when a sample arrives, fulfilling a court order or an additional sample for the same offense is received). Under the 'Offenses' tab, click on the appropriate offense code listed in the table of offenses. Enter the appropriate information and click the 'Save' button.

4.4.3 Enter Court Orders for an Offender as follows:

- 1) Click the 'Court Orders' tab for the saved offender.
- 2) Enter the court order issue and received dates.
Note: the order received date will automatically populate in the sample history, once a sample has been received fulfilling that court order.
- 3) Click the 'Add Offense' link and choose the appropriate offense/court case from the pull down menu. Note: the offenses in the pull down list are those previously entered in the 'Offenses' tab. A court order must have an associated offense in order to be entered; therefore, the offense information must be entered prior to the court order.
- 4) Verify the information is correct and click the 'Save' button.
- 5) Once the court order has been entered it is retained in the laboratory.


4.4.4 Enter samples for an Offender and print sample barcodes as follows:

- 1) Click on the 'Samples' tab for the saved Offender.
- 2) Enter the sample information into the appropriate fields, leaving the barcode field blank. (Tracker will automatically generate a unique barcode number if the field is not filled in).
- 3) If the submitting agency is not listed in the agency pull down menu, it may be added by clicking on the 'edit agencies' button to the right of the agency field.
- 4) If the sample received is pursuant to a previously entered court order, click on the 'Show Unfulfilled Court Orders' link. Click the box next to the appropriate court order to mark it as fulfilled by the sample being entered.
- 5) Enter the number of barcode labels to be printed. One label is to be placed on the DNA sample card/folder, one on the collection report, and one inserted behind the FTA card

envelope for placement on the FTA card at the time of analysis.

- 6) Verify the information is correct and click the 'Save' button. Barcode labels will automatically print upon saving the sample information.
- 7) Staple the DNA Collection Report (if present) to the inside of the DNA sample folder and attach the appropriate barcode labels. The labeled DNA sample card, with collection report, will be stored in the DNA Database Laboratory.
- 8) Barcode labels may be printed at a later time for individual samples, by clicking on the appropriate sample from the 'Samples' tab or for multiple samples, by using the 'Batch Samples' function at the top of the screen.

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT



CODIS SAMPLE DATA ENTRY AND UPLOAD

1.0 BACKGROUND:

The implementation of the Combined DNA Index System (CODIS) in forensic DNA laboratories has provided an additional tool in assisting law enforcement agencies in solving or linking crimes that otherwise may not have resulted in the identification of a suspect. Accurate data entry for upload to NDIS (National DNA Index System) is essential.

It is the responsibility of the Forensic Scientist to generate profiles from convicted offender and/or forensic samples, determine CODIS eligibility, review other scientist's profiles, enter their profiles into CODIS and notify the CODIS Administrator of any potential problems associated with their CODIS DNA profiles. The CODIS Manager is responsible, in part, to serve as the laboratory's chief point of contact for all CODIS communications, maintain CODIS software updates and security, update/add/remove CODIS users, and upload profiles to NDIS.

DNA profiles entered into CODIS will be done so according to CODIS and NDIS Operational Procedures. Idaho currently allows for data entry into the following Indices: Offender (convicted offender); Unidentified Human Remains (deduced missing person, unidentified person); Missing Persons (alleged father/mother, biological child/father/mother, sibling, deduced missing person, missing person, maternal/paternal relative); Forensic Mixture; and Forensic (known and unknown). Profiles entered into the Forensic Index include both solved and unsolved cases in which the profile is associated with a crime and believed to be attributable to the putative perpetrator. Suspect reference samples and profiles matching the victim(s) and/or any elimination samples (e.g. consensual partner samples) may not be entered.

2.0 SCOPE:

To provide a mechanism to ensure accurate data entry for all offender and forensic database samples

3.0 EQUIPMENT/REAGENTS:

CODIS Computer Workstation

Allele Tables and/or CMF files generated from sample analyses.

CODIS Training Manual

4.0 PROCEDURE:

4.1 SAMPLE DATA VERIFICATION:

- 4.1.1 Genotypic data (allele calls) are checked for accuracy and verified during the CODIS/Casework Review process and documented on the appropriate form (Form 306-BI and Form 214-BI, respectively).
- 4.1.2 Forensic and Missing Person/Unidentified Human profiles are verified for CODIS eligibility and Index and documented on the CODIS entry form (Form 218-BI).
- 4.1.3 Additionally, when 'STR Data Entry' is used to enter individual sample data (generally forensic samples) verification of allele calls and specimen category is performed prior to an SDIS search and NDIS upload.

4.2 SAMPLE DATA ENTRY IN CODIS:

- 4.2.1 Once sample data accuracy (forensic and offender) has been verified, data may be entered into the CODIS database either by use of the 'STR Data Entry' module to enter specimen ID/info and allele calls for individual samples or by using the 'Import' function for the entry of batches of data contained in a .cmf (or equivalent) file. Refer to the CODIS Training Manual and course documentation for specific steps in accomplishing these tasks.

4.2.2 A second 'read' must also be entered for the individual samples through 'STR Data Entry'. This may be done by the analyst entering the sample data, or by another analyst logging on and entering the allelic data in the second read box. If a single individual enters both reads, the entire profile (i.e. all loci) for the first read must be entered prior to entering any data into the second read. A 'check' indicates agreement between readers at individual loci and discrepancies in entry must be rectified before saving to SDIS.

4.2.3 When using the 'Import' function, the scientist should validate the .cmf file prior to executing. This will ensure that any typos or inconsistencies (i.e. variant allele definitions/equivalencies) will be identified and corrected prior to importing the batch file data into SDIS.

4.3 CODIS DATA UPLOAD:

4.3.1 Prior to uploading eligible samples to NDIS, a search of SDIS is to be performed to identify any potential matches at the state level. This is typically done by running the 'Autosearcher' function in the software. See BI-303 for verification of any candidate matches identified during the autosearch.

4.3.2 There are various reasons that some samples present at SDIS should not be uploaded to NDIS. Prior to NDIS upload, these samples will be selected in Specimen Manager and 'unmarked for upload'. Generally speaking, an incremental upload will be performed. In Specimen Manager, 'incremental upload' is checked on the 'upload' pull-down menu and 'send upload' is selected. The upload is sent to NDIS as a message attachment via DNACOMM. If any 'candidate matches' are identified at NDIS, a match message will appear in DNACOMM and they will also be reflected in Match Manager. For hit verification see BI-303.

5.0 COMMENTS:

- 5.1 Refer to CODIS Training Manual and course documentation for more specifics if necessary.
- 5.2 The CODIS software is redundant and there is generally more than one way to accomplish many tasks. Using a mechanism other than that listed here is acceptable.
- 5.3 The CODIS software is updated periodically and any necessary changes in procedure provided with new updates supercedes those in procedures written prior to update, if appropriate.
- 5.4 The sample history for convicted offender samples in DNA Tracker will be updated to reflect dates of sample analysis (in progress, complete, failed as appropriate), date sample was outsourced for testing, and/or date profile was entered into CODIS. This is accomplished for individual samples under the 'Samples' tab for the Offender by clicking on the appropriate sample and choosing 'Add History Event'. Multiple samples may be updated using the 'Batch Samples' function at the top of the screen. The current disposition box must be checked as appropriate for the history(s) added.
- 5.5 Eligible forensic profiles of unknown origin (i.e. profiles not matching any of the known reference samples) identified in the course of casework analysis must first be compared against the 'staff' index before entering the profile into the database. The 'searcher' function is used by hand entering the profile and choosing the 'staff' index. This profile will not be saved upon completion of the search and will need to be entered as outlined above for storage in the database.
- 5.6 'Searcher' is also used for performing one-time keyboard searches of forensic unknown profiles. A one-time search will only be performed under exigent circumstances and only on single source samples with a minimum of 8 loci present. The DNA profile must have been generated by or on behalf of a criminal justice agency, by an accredited laboratory, in compliance with the FBI Director's Quality Assurance Standards (QAS), and must be maintained by a Federal, State, or Local NDIS participating laboratory. Additionally, requests

for one-time searches of non-ISPFS profiles must be received from the CODIS Administrator of the laboratory that maintains the profile. When performing a one-time search, the requesting agency will be notified either verbally, or in writing in the event of a negative search result. See BI-303 for verification of candidate matches identified during the one-time search.

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT



CODIS DATABASE HIT VERIFICATION

1.0 BACKGROUND:

The implementation of the Combined DNA Index System (CODIS) in forensic DNA laboratories has provided an additional tool in assisting law enforcement agencies in solving or linking crimes that otherwise may not have resulted in the identification of a suspect. This is accomplished by the electronic storage and maintenance of DNA profiles at the local, state and national levels. Hits are obtained when a candidate match(es) is identified through a database search at any level. Hit verification involves evaluating the candidate match to determine if it is a true match and verification of CODIS offender sample data where necessary and possible.

NDIS CODIS Hit Disposition Reporting & Confirm an Interstate Candidate Match - Operational Procedures

2.0 SCOPE:

To provide a method of sample verification to be performed prior to law enforcement agency notification of a database hit to ensure reporting of only true, confirmed matches.

3.0 EQUIPMENT/REAGENTS:

CODIS
NDIS Procedure Manual; CODIS Training Manual/documentation
DNA Tracker
AFIS
Offender Database Sample(s)
Equipment/Reagents for STR Analysis

4.0 PROCEDURE:

4.1 MATCH VERIFICATION (Forensic):

4.1.1 For 'hits'/matches involving an ISP Forensic Biology evidentiary sample (either case-to-case or case-to-offender and including one-time keyboard searches) the primary responsibility for

match verification follow-up and disposition lies with the Idaho CODIS Administrator.

- 4.1.2** The CODIS Administrator or designee (typically the case analyst) will first evaluate the 'Candidate Match' in Match Manager to determine if there is a basis for exclusion and, therefore, disposition as 'No Match'. A 'No Match' disposition doesn't require further verification or, where applicable, notification of the other laboratory involved with the match.
- 4.1.3** If evaluation in Match Manager demonstrates that the candidate match consists of potential high stringency (e.g., exclusions attributable to different typing systems, Promega v. Applied Biosystems), or possibly moderate stringency in the event of a forensic mixture or degraded sample, the disposition is changed from 'Candidate Match' to 'Pending' until the verification process is complete. In general, for case-to-case matches, the verification will consist of communication between scientists regarding the data and case status, while case-to-offender matches typically necessitate sample verification at the 'offender lab'.
- 4.1.4** Once the status of the 'candidate match' has been resolved, the disposition is set accordingly (e.g., 'No Match', 'Offender Hit', 'Forensic Hit', 'Conviction Match', 'Investigative Information', etc.) and 'Investigations Aided' filled in as appropriate and as outlined in the NDIS 'CODIS Hit Disposition Reporting' procedure.
- 4.1.5** If verification results in a 'hit', a hit report is issued by the case analyst. A copy of the hit report, along with the CODIS match report, is filed in the CODIS file. The original hit report and a copy of the CODIS match report are placed in the associated case file. The appropriate law enforcement agency is notified of the 'hit'. If the law enforcement agency submits a sample from the identified offender, appropriate analysis and issuance of a supplemental report will be performed as in **4.2.5**.

4.2 MATCH VERIFICATION (Offender):

- 4.2.1** For 'hits'/matches involving an ISP Forensic Biology convicted offender sample the primary responsibility for match verification follow-up lies with the CODIS Administrator for the laboratory with the forensic (evidentiary) sample. However, the initial evaluation in Match Manager (see 4.1.2-4) and AFIS sample verification (see 4.2.2) will be initiated as soon as feasible once a verification request has been received from the forensic laboratory.
- 4.2.2** Once a potential match has been confirmed and a verification request received, the associated offender sample will be retrieved from the laboratory. If fingerprint verification was not already done, the sample card will be taken to the latent fingerprint section for an AFIS search of the thumbprint to verify identification of the offender.
- 4.2.3** Following AFIS verification of the thumbprint, re-analysis of the offender sample will be performed as appropriate (i.e., if duplicate analysis has already been performed either as a QC function or as the result of a duplicate sample, analysis will not be repeated) prior to agency notification. In situations where a thumbprint was not received with the DNA sample, or is of insufficient quality for verification, a notation will be made and re-analysis for confirmation may proceed. The forensic laboratory (or law enforcement agency for Idaho cases) will be notified that the Offender could not be verified through thumbprint confirmation.
- 4.2.4** Following sample verification (AFIS and re-analysis as appropriate) the forensic case laboratory, in the case of an interstate hit, or submitting law enforcement agency will be notified of the confirmed hit. Laboratory notification may be made verbally and relevant documentation will be provided to the forensic


case laboratory as requested. In Idaho, initial notification as well as the request for a new DNA sample from the identified offender, may be made verbally. However, written notification and a formal request for a new DNA sample, in the form of a hit report, will be sent to the appropriate law enforcement agency.

4.2.5 For intrastate offender hits, where possible, a newly obtained DNA sample from the offender will be analyzed with all deliberate speed. The analysis of the forensic sample may also be repeated, though this is not required. A supplemental report will be issued delineating the match in the usual manner (i.e., same treatment as for matching suspect sample submitted with case evidence).

4.2.6 The sample history in DNA Tracker will be updated to reflect the date the hit was confirmed for the offender sample.

5.0 COMMENTS:

5.1 A good faith effort will be made to resolve matches within 30 business days from the date of the match or from the date of the verification request for interstate matches involving an Idaho offender.

CODIS SAMPLE REMOVAL**1.0 BACKGROUND:**

Participation in the National DNA database, in accordance with the DNA Analysis Backlog Elimination Act of 2000, necessitates provisions for DNA profile expungement in the event that a qualifying offender's conviction is overturned. Additionally, the Idaho DNA Database Act of 1996 addresses court-granted expungement requests (I.C. §19-5513). Removal of DNA profile data and/or destruction of biological samples obtained from Convicted Offenders may be necessary as a result of conviction reversal or sample collection/submission errors.

Expungement is defined as the removal of DNA profile data from local (LDIS), state (SDIS) and national (NDIS) databases in response to a court order overturning the offender's conviction of a qualifying offense. Expungement will include the removal of identifying information from other laboratory documentation and destruction of the biological sample from which the offender database DNA profile was generated.

Administrative removal is defined as the destruction of a DNA sample and removal of any records relating to that sample. Examples for which administrative removal may be warranted include, but are not limited to, the collection of a sample from a non-qualifying offender, or the notification by the collection agency that removal is warranted. Generally, the determination that a DNA sample does not qualify for inclusion in the database occurs prior to entry of the profile into DNA Tracker and subsequently CODIS (see BI-301); however, there may be instances when the collection agency provides notification of an error after the profile has been generated. In these circumstances, the profile will be removed from the local, state, and national databases as part of the administrative removal.

NDIS Expunge a DNA Profile - Operational Procedures

Idaho Statutes: Title 19, Criminal Procedure, Chapter 55 "The Idaho DNA Database Act of 1996"

2.0 SCOPE:

To provide a protocol for CODIS sample expungement and administrative removal that protects an individual's rights to privacy and maintains the integrity of the Idaho DNA database program. All procedures will be documented on the CODIS Sample Removal Checklist (Form 310-BI).

3.0 EQUIPMENT/REAGENTS:

CODIS Workstation
DNA Tracker Database

4.0 PROCEDURE:

4.1 EXPUNGEMENT

Prior to removal of any DNA profile data, source identification, or biological sample destruction, the CODIS Administrator or designee will verify: 1) the request for expungement is accompanied by a certified court order that meets the requirements outlined in I.C. §19-5513, 2) that the offender has no other qualifying offense(s) and 3) the identification of the original DNA database sample where possible.

4.1.1 A search of DNA Tracker is performed to establish whether or not the specified sample has been received in the laboratory and if so, whether or not a DNA profile has been generated and/or entered into CODIS.

4.1.2 If the sample has been received in the laboratory, a criminal history check in ILETS will be performed to ensure that the offense for which the expungement is requested is the only qualifying offense. If the Offender has a separate qualifying offense, the sample will not be expunged and a copy of the criminal history check will be retained for documentation.

4.1.3 The associated offender sample folder will be retrieved from the laboratory and taken to the latent fingerprint section for an AFIS search of the thumbprint to verify identification of the offender, if not already done. The absence of a thumbprint or a poor quality thumbprint does not preclude the sample expungement from proceeding.

- 4.1.4 The offender record will be removed from DNA Tracker and a sample removal report generated.
- 4.1.5 The DNA profile, if applicable, will be deleted from CODIS, followed by an upload to NDIS for removal at that level. A written notification of the expungement will be made to the NDIS Custodian and a request for verification of deletion at the national level.
- 4.1.6 The original DNA Collection Report received with the sample (or a copy of the sample folder when no report is received) will be marked expunged, along with the scientist's initials and date of expungement.
- 4.1.7 The DNA sample will be destroyed and witnessed by a second scientist.
- 4.1.8 Administrative paperwork documenting the event will be retained in the CODIS file. The paperwork may contain some administrative information about the sample/individual it was collected from and will include the following, as applicable: the court order/request for expungement, criminal history check and thumbprint verification, DNA Tracker sample removal report, CODIS deletion report(s), NDIS correspondence, CODIS sample removal checklist, and copies of any correspondence with the requesting party.
- 4.1.9 An expungement notification letter (or letter indicating why expungement did not occur), the DNA collection report or folder copy, and a copy of the CODIS sample removal checklist will be sent to the requesting party.

4.2 ADMINISTRATIVE REMOVAL

Prior to removal of any DNA profile data, source identification, or biological sample destruction, the CODIS Administrator or designee will verify: 1) the request for administrative removal is in writing from the collecting agency and includes a description of the error resulting in the removal request and 2) that the offender has no other qualifying offense(s).

Removal resulting from an entry error on the part of ISP Forensic Biology staff does not require a written request; although, the nature of the error should be documented. Requests for administrative removal from a party other than the collecting agency or due to biology staff error will be referred to ISP legal staff for a determination of the appropriate action to be taken.

- 4.2.1 A search of DNA Tracker is performed to establish whether or not the specified sample has been received in the laboratory and if so, whether or not a DNA profile has been generated and/or entered into CODIS.
- 4.2.2 If the sample has been received in the laboratory, a criminal history check in ILETS will be performed to ensure that the offense for which the administrative removal is requested is the only qualifying offense. If the Offender has a separate qualifying offense, the sample will not be removed and a copy of the criminal history check will be retained for documentation.
- 4.2.3 The offender record will be removed from DNA Tracker and a sample removal report generated.
- 4.2.4 The DNA profile, if applicable, will be deleted from CODIS, followed by an upload to NDIS for removal at that level.
- 4.2.5 The original DNA Collection Report received with the sample (or a copy of the sample folder when no report is received) will be marked expunged, along with the scientist's initials and date of removal.
- 4.2.6 The DNA sample will be destroyed and witnessed by a second scientist.
- 4.2.7 Administrative paperwork documenting the event will be retained in the CODIS file. The paperwork may contain some administrative information about the sample/individual it was collected from and will include the following, as applicable: the request for administrative removal, criminal history check, DNA Tracker sample removal report, CODIS deletion report(s), CODIS sample removal checklist, and copies of any correspondence with the collection agency.

4.2.8 A sample removal notification letter (or letter indicating why removal did not occur), the DNA collection report or folder copy, and a copy of the CODIS sample removal checklist will be sent to the requesting party, or submitting agency in the event of an internal entry error.

5.0 COMMENTS:

- 5.1 An Offender cannot be deleted from DNA Tracker if there are any samples and/or court orders associated with the offender.
- 5.2 A sample cannot be expunged or deleted from DNA Tracker if there are qualifying offenses associated with the Offender.
- 5.3 The 'Delete Sample' function will be used to remove a sample from DNA Tracker for an Administrative Removal.
- 5.4 The 'Expunge Sample' function will be used to remove a sample from DNA Tracker for an Expungement.

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

EXTRACTION PROTOCOLS FOR PCR DNA TYPING TESTS**1.0 BACKGROUND:**

Many methods exist to obtain DNA, suitable for amplification, from a variety of sources. It has become increasingly important to streamline the extraction process through automation. A variety of extraction methods and robotics systems have become available for the successful extraction and downstream processing of known database DNA samples, using both blood and buccal collection methods. The Promega DNA IQ Extraction System and the Beckman Coulter Biomek 3000 robotic workstation are both examples of systems that have proven successful in the forensic community.

Comey, CT et al. "DNA Extraction Strategies for Amplified Fragment Length Polymorphism Analysis." J For Sci, Vol. 39, 1994, pp. 1254-1269.

McLaren, B et al. "Automating the DNA IQ System on the Biomek 3000 Laboratory Automation Workstation." Profiles in DNA, 2006 Vol. 9(1): 11-13

Tereba, A et al. "High-Throughput Processing of Samples on Solid Supports Using the SlicPrep 96 Device." Profiles in DNA, 2005 Vol. 8(2): 3-5.

DNA IQ System-Database Protocol, Technical Bulletin (revision date 4/06), Promega.

Automated DNA IQ System Protocol for the Beckman Coulter Biomek 3000, Technical Bulletin (revision 6/09), Promega.

2.0 SCOPE:

To provide appropriate protocols for the extraction of database sample DNA, both manual and automated, suitable for PCR amplification and subsequent analyses.

3.0 EQUIPMENT:

Pipettes and tips
Harris Punch
BSD600-Duet Semi-Automated Sample Punch Instrument w/ Scanner

96-well Reaction Plate
MicroAmp Tubes
1.5ml microcentrifuge Tubes (1.5ml tubes)
65°C Oven
SlicPrep Device
Adhesive Foil Seal
70°C Water Bath
Centrifuge with Swinging Plate Rotor
Biomek 3000 Robotics System with MagnaBot 96 Magnetic Separation
Device, 1/4 Inch Foam Spacer, Variomag Teleshake, and Heating
Block
Reagent Reservoirs
2.2ml Square-Well, Deep Well Plate
1.2 ml Round-Well, Deep Well Plate

4.0 REAGENTS:

FTA Purification Reagent
PCR TE (TE, 10mM Tris-HCl; 0.1mM EDTA, pH 8.0)
DNA IQ System Reagents
99% Isopropyl Alcohol
95-100% Ethanol
1M Dithiothreitol (DTT)

5.0 DNA EXTRACTION PROCEDURES:

5.1 FTA EXTRACTION FROM OFFENDER DATABASE SAMPLES:

Note: Since the DNA remains bound to the FTA card, regular, non-filter pipette tips may be used throughout and a single tip may be used for each reagent. A multi-channel pipettor and 96-well reaction plate may be used for larger sample batches.

5.1.1 Remove 1-3 "punches" from the FTA card using a 1.2mm Harris punch (this is accomplished by placing punch firmly on card and twisting 1/2 turn clockwise and 1/2 turn counterclockwise). Eject/punch sample(s) into microAmp tube(s)/plate. Alternatively, the BSD puncher may be used to punch 1 2mm punch directly into the plate (see **5.2.1**).

5.1.2 Add 150µl FTA reagent to microAmp tube(s)/plate, gently mix and incubate at RT for ~5 minutes.

5.1.3 Remove and discard FTA reagent from sample(s).

5.1.4 Repeat 5.1.2-5.1.3 twice.

5.1.5 Add 150µl TE to microamp tube(s)/plate, gently mix and incubate at RT for ~5 minutes.

5.1.6 Remove and discard TE from sample(s).

5.1.7 Repeat 5.1.5-5.1.6 twice.

5.1.8 Make sure the punch is at the bottom of the microAmp tube(s)/plate, using a sterile pipette tip if necessary. Place tubes/plate, uncovered in 65°C oven for ≥2 hours.

5.1.9 Proceed to PCR Amplification (see BI-316).

5.2 AUTOMATED DNA IQ EXTRACTION FROM OFFENDER DATABASE SAMPLES

5.2.1 SAMPLE PUNCHING USING THE BSD600-DUET

5.2.1.1 Fill out the BSD/Biomek Load Sheet on the 'BSD Load Sheet' tab of the Database worksheets template (Form 312-BI). Print a copy for the batch record/file. Choose the 'BSD Input File' tab and ensure the information is correct and corresponds to the Load Sheet information entered. Delete any unused wells, making sure to keep 'End File' as the last cell. Perform a 'Save As' of the worksheet to disc (e.g. USB drive) for subsequent transfer to the BSD600 puncher. The document must be saved as a .txt file. Save the template as samples are automatically populated on the other tabs for subsequent processing steps.

Note: The last two columns of the plate should be left empty for control placement in subsequent steps. This allows for 78 samples, an extraction control, and a reagent blank to be tested per plate.

5.2.1.2 Ensure the BSD tray table is adjusted properly by placing the appropriate plates in both positions and adjusting the height, as necessary. The bottom of the chute should be

~2mm from the top of the plates, measured at each corner.

5.2.1.3 Turn on the BSD puncher and computer and plug in the air/pressure pump. Make sure water is present in both filter bottles (~50-100ml). Double click the 'BSD Duet Main Menu' icon and login using the appropriate user name and password.

5.2.1.4 Choose **Configure System** and open the **Files** tab. Ensure that **File Input** is selected from the **Sample Number Mode** dropdown list. Click the **Browse** button and locate your saved input file to fill in the **Input File Name** box. Click **Save and Exit**.

Note: Using the input file option only allows for a single plate to be punched at a time. To punch multiple plates, the process must be repeated for each new set of samples.

5.2.1.5 Choose **Distribute Spots** and click **Continue**. When prompted to select test group to punch, highlight **All Available Tests** and click **Continue**. Check the box next to **Input File test** and check the samples box on the right. Press **Continue**. This test file will punch 3 of the 3.2mm spots (verify before beginning). If punching for FTA extraction, the number and size of punches will need to be changed under the Edit Test Sequences option prior to beginning.

5.2.1.6 Place the appropriate plate into position 1 of the BSD tray table with well A1 located in the upper left corner. Click **Continue** to begin punching. When using the SlicPrep 96 Device, remove the white collar, push the basket down into the plate, and set the collar aside before processing.

5.2.1.7 Scan the sample barcode using the scanner attached to the BSD and place the card under the card clamps on top of the instrument. Use the light guide dots to properly position

the card for punching. Press the foot switch when ready to punch or use the delay switch for automatic punching. Repeat until all samples have been punched.

Note: The scanned barcode must match that on the input file or the instrument will not punch the card.

- 5.2.1.8 Remove the plate and proceed to SlicPrep pre-processing (see 5.2.2) or FTA extraction (see 5.1).

5.2.2 SLICPREP PRE-PROCESSING

- 5.2.2.1 Add 400ul DNA IQ Lysis Buffer with DTT (see 5.2.3.1) to each well of the device plate. Seal the top of the plate with a foil seal and place into a 70°C water bath for 15-60 min.
- 5.2.2.2 Without removing the foil, raise the filter basket and place the white collar back onto the device by sliding it under the basket and aligning the detents on the device.
- 5.2.2.3 Centrifuge the device at 1500 x g (~3000 RPM) for 5 min. in a swinging plate rotor. Remove and discard the collar and filter basket upon completion of centrifugation.
- 5.2.2.4 Proceed to DNA IQ Extraction on the Biomek 3000 (see 5.2.4).

5.2.3 DNA IQ REAGENT PREPARATION

- 5.2.3.1 When opening a new DNA IQ kit, add 1.5ml 1M DTT to the bottle of Lysis Buffer (1ul DTT/100ul lysis buffer) and mix. Mark the bottle with initials and the date the DTT was added. The 'activated' Lysis buffer may be used for up to 2 months following DTT addition. Alternatively, the lysis buffer may be prepared as needed based on the number of samples to be processed.

5.2.3.2 Prepare the new wash buffer (70ml 2X bottle) by diluting with 35ml Isopropyl Alcohol and 35ml Ethanol for a final 1X concentration. Mark the bottle with initials and the date of alcohol addition. The 1X wash buffer may be used for a period of 1 month after 'activation'. Alternatively, the wash buffer may be prepared as needed based on the number of samples to be processed.

5.2.3.3 Prepare/aliquot each of the DNA IQ reagents at the time of use, based on the following calculations and reservoir positions (instructions listed in the 'Reagents for Solid Support' comments window):

- 1: Elution Buffer = 1.2ml + (#samples x 100ul)
- 2: Empty
- 3: 1X Wash Buffer = 1.5ml + (#samples x 0.3ml)
- 4: Lysis Buffer w/ DTT = 1.5ml + (#samples x 0.1ml)
- 5: Lysis Buffer w/ DTT + Resin = [860ul + (#samples x 43ul)lysis] + [140ul + (#samples x 7ul)resin]
- 6: Empty

5.2.4 DNA IQ EXTRACTION USING THE BIOMEK 3000 ROBOTICS SYSTEM

5.2.4.1 Turn on the Biomek 3000, computer, and Watlow heater. Make sure that the heater is set to 85°C and that the shaker electrical box dial is turned to the red at the far left stop.

5.2.4.2 Check the teleshaker connection by choosing **start > programs > teleshaker v.1.2**. In the test program window, choose **options > interface > com 2** and click **Ok**. Select **file > connect**.

5.2.4.3 Open the Biomek Software and choose the 'All Others Method'.

5.2.4.4 Set up the Biomek 3000 deck as follows (outlined in the 'Tools & Labware for Solid Support Method' comments window):

Rack 1 (A1): Tool rack with MP200 (1) and gripper (3)
ML1 (A2): P250 barrier pipette tips
ML2 (A3): P250 and/or Axygen 165ul barrier pipette tips (up to 48 samples)
ML3 (A4): P250 and/or Axygen 165ul barrier pipette tips (up to 88 samples)
ML4 (A5): P250 and/or Axygen 165ul barrier pipette tips (up to 96 samples)
P6 (A6): Heating system
P1 (B1): PCR plate and base
P2 (B2): empty
P3 (B3): Reservoir frame and reservoirs with reagents
P4 (B4): Pre-processed deep square-well plate
P5 (B5): Deep round-well plate stacked on the Magnabot
P7 (B6): Varionag Teleshaker

Note: The number of tip boxes will depend on the number of samples to be extracted. The instrument will only extract full columns on the plate. If using the Axygen tips, do not use them as the first box (ML1).

5.2.4.5 Click the green run instrument arrow button to begin the run. When prompted, enter the number of columns to be extracted, choose 0 for solid support method, and either 50ul, or 100ul for elution volume. The 50ul volume may be preferable for lower concentration samples. Verify the reagent and deck layout and click ok.

5.2.4.6 When the run has finished remove the PCR/elution plate. Discard the deep-well plates, remaining buffer, and pipette tips. Clean the reservoirs between runs.

5.2.4.7 At the completion of all runs for the day, wipe down the deck with 10% bleach or Dispatch solution. Tools may also be cleaned with 70% Ethanol if necessary.

5.2.4.8 Proceed to realtime PCR (BI-314) for quantification of the purified DNA obtained in 5.2.4.6.

6.0 DNA EXTRACTS:

6.1 After a plate of samples has been extracted and during subsequent analyses (e.g. quantification and amplification), the DNA extracts may be stored sealed with foil at 4°C. For longer storage periods, the extracts should be frozen at approximately -20°C. These extracts are in-progress work product during this stage(s).

6.2 Any extract remaining, following the successful completion of analysis and data review, will be discarded.

7.0 Comments:

- 7.1 A reagent blank shall be carried through all extraction steps to check the purity of the reagents being used. There need only be one reagent blank per extraction plate, but the analyst may choose to run additional blanks on a given plate.
- 7.2 These procedures may be used for blood or buccal samples on FTA paper.
- 7.3 The DNA IQ resin must be thoroughly mixed/resuspended prior to aliquoting.
- 7.4 The final elution volume attained on the Biomek 3000 will be ~8-10ul less than the volume chosen, due to evaporation during the processing.
- 7.5 The Biomek 3000 performance verification color plate (Form 428-QC) is run by choosing the 'Color Plate/Performance Check' method in the Biomek Software. Deck set up and dye preparation is described in the 'Instrument Setup' step and the 'Deck Setup Instructions' comment window.
- 7.6 The DNA IQ extraction procedure may be performed manually, if necessary, following the procedure outlined in the DNA IQ System-Database Protocol, Technical Bulletin (revision date 4/06). This procedure follows the same steps as performed on the Biomek 3000.

**DNA QUANTIFICATION: REAL-TIME PCR****1.0 BACKGROUND:**

DNA methodologies that employ the PCR, such as STR analysis, necessitate consistent quantification of human DNA to obtain optimum data.

"Developmental Validation of the Quantifiler™ Real-Time PCR Kits for the Quantification of Human Nuclear DNA Samples," Green, R.L., et al, Journal of Forensic Science, Vol. 50, No. 4, pp. 809-825.

"Improving Efficiency of a Small Forensic DNA Laboratory: Validation of Robotic Assays and Evaluation of Microcapillary Array Device," Crouse, C., et al, Croat Med J 2005 Vol. 46, No. 4, pp. 563-577.

Quantifiler™ Kits (Quantifiler™ Human DNA Quantification Kit and Quantifiler™ Y Human Male DNA Quantification Kit) User's Manual, Applied Biosystems.

7500/7500 Fast Real-Time PCR Systems Maintenance Guide, Applied Biosystems.

2.0 SCOPE:

To provide a reliable method for the consistent quantification of variable amounts of human DNA isolated from database samples.

3.0 EQUIPMENT/REAGENTS:

ABI 7500/Computer	96-well Reaction Plate
ABI 7500 SDS Software	96-well Reaction Plate Base
Pipettes and Tips	Optical Adhesive Covers
Reservoir (optional)	Centrifuge (optional)
Quantifiler™ Human Kit	Microcentrifuge Tubes/Strip Tubes
20 µg/ml Glycogen (optional)	PCR-TE

4.0 PROCEDURE:

4.1 PREPARATION OF DNA STANDARDS:

- 4.1.1 Label 8 sterile microfuge tubes, or strip tubes A through H or 1-8.
- 4.1.2 Dispense 30 μ l (or adjusted amount according to the kit QC results: Form 419-QC) of PCR-TE into tube A (Std. 1) and 20 μ l of PCR-TE into tubes B-H (Std. 2-8).
- 4.1.3 Mix the Quantifiler Human DNA Standard thoroughly by vortexing 3-5 seconds. Transfer 10 μ l to tube A (Std. 1). Mix the dilution thoroughly.
- 4.1.4 Prepare Std. 2-8 via a serial dilution by mixing and subsequent 10 μ l transfers from tubes A through H. The dilution series consists of 50, 16.7, 5.56, 1.85, 0.62, 0.21, 0.068, and 0.023 ng/ μ l, respectively.

4.2 REACTION PREPARATION:

- 4.2.1 Determine the number of samples to be quantified (including, at minimum, 2 sets of DNA standards).
- 4.2.2 Finish filling out the 7500 Load Sheet on the Database worksheets template (Form 312-BI). Print a copy for the batch record/file. Choose the '7500 Plate Document' tab and ensure the information is correct and corresponds to the Load Sheet information. Perform a 'Save As' of the Plate Document Worksheet to disc (e.g. USB drive) for subsequent transfer to the ABI 7500. The document must be saved as a .txt file.
- 4.2.3 Calculate the volume of reaction components needed, based upon the number of samples to be quantified and adding 2 or 3 reactions to compensate for loss and variability due to pipetting. The following are the volumes needed per reaction.

Quantifiler PCR Reaction Mix	12.5 μ l
Quantifiler Human Primer Mix	10.5 μ l

Note: The volume of reaction components necessary to prepare the Master Mix will be automatically calculated upon Load Sheet data entry.

- 4.2.4 Thaw the primer mix and vortex 3-5 seconds. Pulse-spin prior to opening the tube. Mix the PCR Reaction Mix by gently swirling the bottle prior to use.
- 4.2.5 Place a 96-well reaction plate into a base, being careful not to touch the top or individual wells. Do not place the plate directly onto the counter or any surface other than its base or the ABI 7500 thermal block.
- 4.2.6 Prepare the Master Mix by pipetting the required volumes of primer and reaction mixes into an appropriately sized microcentrifuge tube or reservoir. Mix thoroughly and pulse spin if mixing in a tube.
- 4.2.7 Carefully pipet 23 μ l of the PCR Master Mix into the bottom of each reaction well to be used. 'Blowing-out' the pipette is not recommended to avoid splashing and/or bubbles in the well.
- 4.2.8 Add 2 μ l of sample or standard to the appropriate reaction well, being careful to avoid bubbles as much as possible.
- 4.2.9 Seal the reaction plate with an Optical Adhesive Cover. Proceed to 4.3.

4.3 RUNNING THE REACTION:

- 4.3.1 Turn on the 7500 computer and login with the appropriate user name and password. After the computer has completely started up, power on the 7500 instrument, allowing it to warm up at least ~30 seconds. Launch the ABI 7500 SDS Software.
- 4.3.2 Open the instrument tray by pushing on the tray door. Place the plate into the tray holder so that well A1 is in the upper-left corner and the notched corner of the plate is in the upper-right corner.
- 4.3.3 Close the instrument tray by gently pushing the right side of the tray door.

4.3.4 In the SDS software, select **File>New** and choose **Absolute Quantitation** for Assay, **96-Well Clear** for Container, and **Quantifiler Human Database** for Template.

4.3.5 Import the previously saved plate document by selecting **File>Import Sample Set-Up**. Browse to locate the saved .txt file and choose **OK**.

4.3.6 Review the plate document to ensure the appropriate detectors and tasks have been applied to each sample. Change the task for any unused wells to NTC in **View>Well Inspector** or by highlighting and double clicking on a well(s). Make any other changes, as necessary. Select the **Instrument** tab and review the thermal cycler conditions [Stage 1: 1 cycle, 95°C, 10:00 min.; Stage 2: 40 cycles, 95°C, 00:15 min, 60°C, 1:00 min.; 25µl sample volume; 9600 emulation; Data collection: Stage 2, Step 2 (60.00@1:00)]

Note: Detectors are created during the initial instrument set-up and/or kit usage. Refer to the Quantifiler Kits User's Manual (page 2-11) for instructions on creating detectors if needed.

4.3.7 Save the plate document as a .sds file with the appropriate plate name.

4.3.8 Under the **Instrument** tab, click **Start** to begin the run. When the run has completed, proceed to 4.4.

4.4 ANALYSIS AND RESULTS:

4.4.1 Open the plate document to be analyzed.

4.4.2 Select **Analysis>Analysis Settings** and verify the settings are set as follows: All for Detector, Manual Ct, 0.200000 for Threshold, Manual Baseline, 3 for Baseline Start (cycle), and 15 for Baseline End (cycle). Click **OK**.

4.4.3 Select **Analysis>Analyze**.

4.4.4 In the **Results** tab, select the **Standard Curve** tab and choose Quantifiler Human as the detector. Review the data for inconsistencies from the following:

An R^2 value of >0.99 indicates a close fit between the standard curve regression line and the individual C_T data points of quantification standard reactions.

An R^2 value of <0.98 needs further analysis of the standard curve for problems. Refer to the Quantifiler Kits User's Manual (page 5-6) for troubleshooting guidelines.

The slope should fall within the typical slope range of -2.9 to -3.3 . A slope of -3.3 indicates 100% amplification efficiency.

- 4.4.5 Select the **Amplification Plot** tab (in the Results tab) and choose either the Quantifiler Human, or the IPC detector. Ensure the Threshold is set to 0.20 before proceeding (**Note:** the threshold bar will be green if the data has been analyzed and red if analysis is needed). Highlight the sample(s) of interest in the table to view the associated plot(s). Review the plots for both detectors for amplification and/or inconsistencies.
- 4.4.6 Select the **Component** tab within the Results tab. The halogen lamp may need replacement if the dye signal lines contain spikes or appear wavy/unstable and/or if the Rox value begins approaching or has fallen below 500. See Comment 2. Note: it is important to use the same sample well each time.
- 4.4.7 In the Results tab, select the **Report** tab and highlight the sample(s) of interest to view the results. Review the Qty column to determine the amount of DNA present in each sample. Review the Internal Positive Control (IPC) C_T value for each sample. It should fall within a range of 20-30. If the value is >30 for a particular sample, there may be an indication of inhibition.
- 4.4.8 Export the report. Within the report tab, select **Tools>Report Settings** and check the appropriate boxes to be displayed in the report and click **OK**. Print a copy of the Standard Curve for the batch record. Select **File>Export** to export the report (e.g. to USB drive) as a .csv file.
- 4.4.9 Import the .csv file into the **7500 Raw Data** tab of the Database Excel worksheet/template. Choose the **7500 Results** tab and review the imported data. Delete any

unused wells from the sheet. Adjust values in the final concentration and ul Sample for dilution columns. Print a copy of the results sheet for the batch record.

5.0 COMMENTS:

- 5.1 Refer to the Quantifiler Kits User's Manual for specific thermal cycler conditions, additional user information, and troubleshooting guidelines.
- 5.2 If the Component Dye signals appear unstable and/or Rox values approach 500, the Halogen Lamp may be checked manually to determine if replacement is needed. Place the Green Calibration Tray in the block. Select **Instrument > Calibrate** and set the exposure time to **4096ms**, lamp control to **Max**, and select **Filter A**. Click Snapshot and observe results. Expected results should consist of red fluorescence displayed in all wells. Lack of fluorescence indicates the need for lamp replacement. The lamp status should be checked as well by selecting **Instrument > Lamp Status/Replacement** and viewing the condition.
- 5.3 In order to extend the life of the Halogen Lamp, the instrument should be turned off anytime it is not in use. Lamp life is approximately 2,000 hours.

**STR AMPLIFICATION: PP16 HS****1.0 BACKGROUND:**

DNA analyses have revolutionized forensic biology. The advent of PCR allowed scientists to analyze evidentiary material present in minute quantities and degraded states. The identification of forensically significant STR loci has allowed scientists to combine the discrimination attainable with the older RFLP technology with the speed and sampling capabilities of other PCR-based methodologies. The PowerPlex™ 16 HS allows the co-amplification of the core CODIS 13 loci, as well as, Amelogenin, and two pentanucleotide-repeat loci, Penta D and Penta E.

Butler, J. *Forensic DNA Typing: Biology and Technology Behind STR Markers*. (2001) Academic Press.

PowerPlex™ 16 HS System Technical Manual

2.0 SCOPE:

To provide a reliable method for consistent, high quality amplification of DNA from offender database samples ensuring the generation of suitable PCR product for capillary electrophoresis and analyses of these STR loci.

3.0 EQUIPMENT/REAGENTS:

Thermocycler
Microcentrifuge
Microcentrifuge tubes and/or reservoir
96 well Reaction plate and base
Pipettes and tips
Foil Seal or Optical Adhesive Cover
Compression pad
PowerPlex™ 16 HS Kit Contents

4.0 PROCEDURE:

4.1 DNA TEMPLATE:

4.1.1 Based upon the quantity of DNA isolated and its initial concentration, the scientist should have all samples at an optimal concentration for amplification (e.g., 0.1ng/ μ l). It is also convenient to have all samples, that are to be amplified at the same time, at the same concentration if possible for ease in the preparation of PCR Master Mix and reaction additions. The maximum amplification volume for low level or undetected samples is 4.4 μ l for PowerPlex 16 HS, quarter volume reactions. For larger volume samples, it may be necessary to concentrate the sample prior to amplification.

4.1.2 The amount of DNA template added to a quarter volume amplification reaction should be targeted at ~0.25ng, but may be adjusted as necessary.

4.2 AMPLIFICATION SET-UP:

4.2.1 Thaw the PowerPlex 16 HS 5X Master Mix, PowerPlex 16 HS 10X Primer Pair Mix, and Amplification Grade Water (optional) contained in the HotStart kit.

4.2.2 Calculate the volume of reaction components needed based upon the number of samples (including extraction and amplification controls) to be amplified and adding 1 or 2 reactions to compensate for loss and variability due to pipetting. Use the 'Amp. Worksheet' tab of the Database Worksheets Template (Form 312-BI) to automatically calculate necessary volumes. Print a copy of the completed worksheet for the batch record/file. The following is a list of the 'fixed' amounts to be added for a 6.25 μ l reaction.

PowerPlex 16 HS 5X Master Mix	1.25 μ l
PowerPlex 16 HS 10X Primer Pair Mix	0.6 μ l
¹ DNA Template + dH ₂ O	4.4 μ l

¹For FTA extracted samples there is no volume for the DNA template so 4.4 μ l of dH₂O will be added to these tubes.

Note: The amount of Reaction Mix added to each sample is dependent on the volume needed to add the DNA template.

- 4.2.3 Prepare the reaction mix by pipetting the required volumes of 5X master mix, 10X primer mix, and water (if needed) into an appropriately sized microcentrifuge tube or reservoir. Mix thoroughly and pulse-spin if using a tube.
- 4.2.4 Pipet the appropriate amount of reaction mix into each well of a labeled 96-well reaction plate. The negative amplification control should be the last sample processed.
- 4.2.5 If DNA concentrations were not the same, add appropriate volume of dH₂O as necessary.
- 4.2.6 Pipet each DNA sample into the appropriate well. Use 9947A control DNA for the positive amplification control and dH₂O for the negative amplification control.
- 4.2.7 Thoroughly seal the plate with an adhesive foil seal. Use a seal applicator to tightly seal between and over each well.
- 4.2.8 Remove lab coat and, touching only the plate, transport the samples to the thermal cycler in the Amp/PostAmp room, using the other hand on the door knob.
- 4.2.9 Place a compression pad over the plate and put it into the thermal cycler. Remove/discard gloves before exiting the Amp/PostAmp room.

4.3 THERMAL CYCLING PARAMETERS:

4.3.1 After the samples have been placed in the thermal cycler, turn on the power and select the appropriate pre-programmed cycling profile.

4.3.1.1 For DNA IQ-extracted and low level FTA-extracted samples, use 'pp16HS' with the following conditions:

96°C for 2 minutes, then:

ramp 100% to:
94°C for 30 seconds,
ramp 29% to:
60°C for 30 seconds
ramp 23% to:
70°C for 45 seconds
for **10 cycles**, then:

ramp 100% to:
90°C for 30 seconds
ramp 29% to:
60°C for 30 seconds
ramp 23% to:
70°C for 45 seconds
for **20 cycles**, then:

60°C for 30 minutes, then:

4°C soak

4.3.1.2 For FTA-extracted samples (other than low template), 'PP16HS17' with the following conditions may be used:

96°C for 2 minutes, then:

ramp 100% to:
94°C for 30 seconds,
ramp 29% to:
60°C for 30 seconds
ramp 23% to:
70°C for 45 seconds
for **10 cycles**, then:

ramp 100% to:
90°C for 30 seconds
ramp 29% to:
60°C for 30 seconds
ramp 23% to:
70°C for 45 seconds
for **17 cycles**, then

60°C for 30 minutes, then:

4°C soak

5.0 AMPLIFIED DNA PRODUCT:

- 5.1 After cycling has concluded remove samples from thermal cycler. Samples should be run on the Genetic Analyzer as soon as possible after amplification. Prior to capillary electrophoresis and/or before analysis is completed the samples may be stored at 4°C. For longer storage periods, samples should be frozen at -20°C. Amplified product is ONLY stored in the Amp/PostAmp room.
- 5.2 At a point in time after STR analysis is completed (i.e., data has been reviewed and approved for upload), the amplified product will be disposed of in a biohazard container in the amp/post-amp room. As needed, this container will be sealed and transported directly to the dishwashing room. The container will be placed into a second biohazard bag, sealed and disposed of with other biohazardous material.

6.0 COMMENTS:

- 6.1 Clean surfaces with freshly made 10% bleach solution or Dispatch® prior to set-up.
- 6.2 Wear gloves at all times during amplification set-up.
- 6.3 Mix all reagents thoroughly (e.g., vortex) and pulse-spin them in microfuge prior to dispensing.
- 6.4 If excessive evaporation is observed following amplification, an optical adhesive cover may be used to seal the plate rather than the foil seal.
- 6.5 When using a ¼ volume reaction, the amplification reagents may be completely absorbed by the 2mm BSD punch, if used for the FTA extraction method. If this happens, it is recommended that 5ul sterile water be added post-amplification, and thoroughly mixed prior to running on the genetic analyzer.

STR TYPING: CAPILLARY ELECTROPHORESIS AND DATA ANALYSIS**1.0 BACKGROUND:**

Any eukaryotic genome is interspersed with repeated DNA sequences that are typically classified by the length of the core repeat sequence, and the range of contiguous repeats typically seen or the overall length of the repeat region. STR (Short Tandem Repeat) loci are scattered throughout the genome occurring every 10,000 nucleotides or so, and have core repeat units of 2-6bp in length with overall lengths of less than 400 bp.

STR loci examined for human identification purposes were selected for the following characteristics: 1) high discriminating power (generally >0.9) with observed heterozygosity of >70%, 2) loci on separate chromosomes to avoid linkage, 3) ability to obtain robust, quality, reproducible data when multiplex amplification is performed, 4) low stutter, 5) low mutation rate and 6) small allele sizes (<500 bp) for enhancement of analysis of degraded samples.

By 1997, as the result of a community-wide forensic science effort, the following 13 STR loci, all tetranucleotide repeats, were selected as the basis for NDIS, the CODIS (COMbined DNA Index System) National Database: D3S1358, TH01, D21S11, D18S51, D5S818, D13S317, D7S820, D16S539, CSF1PO, vWA, D8S1179, TPOX, FGA. When all 13 CODIS core loci were examined, the average random match probability was found to be <1 in 1×10^{12} among unrelated individuals, offering the promise of individualization.

In addition to the 13 core CODIS loci, the PowerPlex 16 HS multiplex includes Amelogenin, a gender identification locus, and two pentanucleotide repeat STR loci, Penta D and Penta E. STR typing, with amplified products generated from this kit, separated by capillary electrophoresis on the 3130xl Genetic Analyzer with data collection and analysis software employed in developing the genetic profiles, will be used to produce STR profiles from convicted offender samples for entry into CODIS.

Butler, J. *Forensic DNA Typing: Biology and Technology Behind STR Markers*. (2001) Academic Press.

PowerPlex™ 16 HS System Technical Manual

PowerPlex® Matrix Standards, 3100/3130, Technical Bulletin

ABI 3130/3130xl Genetic Analyzer Getting Started Guide

ABI 3130/3130xl Genetic Analyzer Maintenance Troubleshooting and Reference Guide

GeneMapper™ ID Software User Guide

2.0 SCOPE:

To provide a reliable method for generating STR genetic profiles from offender DNA database samples.

3.0 EQUIPMENT/REAGENTS:

3130XL Genetic Analyzer with Data Collection Software
GeneMapper™ ID Software
Computers
9700 Thermal Cycler
Pipettes and Tips
Benchtop Cooler
Capillary Arrays
96 Well Reaction Plates and Septa
Buffer Reservoirs and Septa
POP4 Polymer
Genetic Analyzer Buffer
PowerPlex® 16 Kit HS Contents
PowerPlex® 16 Matrix Standards
Deionized Formamide
Nanopure Water

4.0 PROCEDURE:

4.1 AMPLIFIED FRAGMENT DETECTION USING THE 3130xl

Note: Prior to using the ABI 3130xl Genetic Analyzer for samples, a spectral calibration (matrix standards) must be run to achieve proper color separation of the dyes used for the amplification primers, allelic ladders and size standard. To prepare a matrix, four standards are run under the same capillary electrophoresis conditions that

will be used for samples and allelic ladders. Use the 3130 Matrix Standard set, which includes the Fluorescein Matrix, JOE Matrix, TMR Matrix and CXR Matrix for the blue, green, yellow and red matrix standards, respectively. This is performed when necessary due to performance, or after any instrument maintenance/repair that involves adjustment/replacement of the CCD camera or laser.

Additionally, a Spatial Calibration must be performed prior to running any samples. The instrument uses images collected during the spatial calibration to establish a relationship between the signal emitted for each capillary, as well as the position where that signal falls and is detected by the CCD camera. This is performed any time a capillary is installed or replaced (including temporary removal of a capillary) or if the instrument is moved.

- 4.1.1** Turn on the computer, turn on the instrument, start Data Collection Software and wait for green squares to appear for all applications on the service console. Expand the necessary subfolders on the left tree pane of Data Collection. Refer to the ABI 3130/3130xl Genetic Analyzers Getting Started Guide for detailed instructions on instrument set-up (including creation of instrument protocols, results groups, and spatial calibration). Fill-in appropriate information in the 3130xl Injection Log (Form 422B-QC).

Shut down is performed in the opposite order (Data Collection software, 3130xl, then computer). The Data Collection Software must be closed by choosing 'Stop All' and waiting for all red symbols to appear before closing. Never use the 'X' to close while green or yellow symbols are displayed.

- 4.1.2** Create a new plate record:

- 4.1.2.1** For a spectral calibration plate expand the tree pane of the Data Collection Software and click 'Plate Manager,' under 'ga3130xl'. Choose 'New', and fill in the dialog boxes, with 'Spectral Calibration' as the application. Fill in the applicable dialog boxes on the

Spectral Calibration Editor as follows
(clicking 'OK' when complete to save):

Sample Name: date_Spectral

Priority: May optionally be changed to a number <100 for injection priority.

Instrument Protocol 1: Choose the Spectral instrument protocol from the drop down menu

PowerPlex 16 specific run module and protocol settings for Spectral Calibrations are as follows:

Module Type:	Spectral
Template:	Spect36_POP4
Inj. kV:	3
Data Delay Time:	100
Run Time (seconds):	800
Protocol Type:	Spectral
DyeSet:	F
Array Length:	36
Chemistry:	Matrix Standard
Lower condition bound:	4.0
Upper condition bound:	12.0

Inj.Secs: 5

4.1.2.2 For a sample plate finish filling out the 3130XL Load Sheet on the database worksheets template (Form 312-BI). Print a copy for the batch record/file. Choose the '3130xl Plate Document' tab and ensure the information corresponds to the Load Sheet information. Verify the information on the template is as follows:

Container Type: 96-Well

Application Type: regular

GeneMapper: GeneMapper_Generic_Instance

Sample Name:

Allelic Ladder: LADDER (or PP16_LADDER)

Controls: POS [or (+), etc.], NEG [or (-), etc.], RB (reagent blank)

Database samples: ID#####
(e.g., ID2001001412)

Priority: May optionally be changed to a number <100 for injection priority.

Sample Type: Sample Categories of 'Sample', 'Allelic Ladder', 'Positive Control', or 'Negative Control', may optionally be typed in.

Results Group 1: Enter the appropriate results group. These are typically denoted by the analyst initials and should automatically populate from the Load Sheet.

Instrument Protocol 1: Enter the appropriate instrument protocol (e.g., PP16_5 sec, PP16_3 sec, PP16_10 sec).

PowerPlex 16 specific run module and protocol settings for sample plates are as follows:

Module Type: Regular
Template: HIDFragmentAnalysis36_POP4
Inj. kV: 3
Run Time (seconds):2000

Protocol Type: Regular
DyeSet: F

Inj.Secs: 3-10

3 secs - Generally produces good results for Allelic Ladders and

samples \geq .25ng (injection times may be adjusted [3-10 seconds per analyst's discretion]).

5 secs.- Samples $<$.25ng generally produce good results.

- 4.1.2.3** Delete any unused wells. Perform a 'Save As' of the Plate Template Worksheet to disc (i.e. USB drive) for subsequent transfer to the ABI 3130xl. The document must be saved as a .txt file (tab delimited).
- 4.1.2.4** Import the previously saved plate record by selecting 'Import' on the 'Plate Manager' window. Browse to locate the saved .txt file and choose 'OK'.
- 4.1.2.5** Open the imported plate record by highlighting it and clicking 'Edit'. Review the information in the GeneMapper Plate Editor to ensure that it is correct or make changes as necessary. Click 'OK' when complete to save the plate record.
- 4.1.2.6** To perform more than one run of a sample (e.g. multiple injection times), select Edit/Add/Sample Run in the GeneMapper Plate Editor window. This will add additional Results Group and Instrument Protocol columns to the end of the plate record. These additional runs may be added at any point in the run, prior to the last injection, if the scientist notices that a sample would benefit from re-injection (e.g., repeat because of bad injection or to vary injection times [from 3-10 seconds]). Additional Results Groups and Instrument Protocols may also be filled in on the original Load Sheet template prior to importing.
- 4.1.3** In the manual control window, the scientist may choose to set the oven to 60°C so that it will be ready to run. Choose Oven in the 'Send Defined Command for' drop down menu box. In the 'Command Name' box, choose Turn On/Off oven, with a

'Value' of On, and click 'Send Command'. Next, in the 'Command Name' box, choose Set oven temperature, with a 'Value' of 60.0 and click 'Send command'. **Note:** once the oven has been turned on and the temperature set, the oven will only preheat for 45 minutes before shutting itself off.

4.1.4 Prepare samples for capillary electrophoresis:

4.1.4.1 For amplified products (including controls), typically 1 μ l-1.5 μ l rxn is added to 10 μ l of ILS Master Mix (made by adding 0.5 μ l ILS600 size standard/sample; 9.5 μ l deionized formamide/sample and adding quantities for N+2 in Master) that has been dispensed into the wells of a pre-labeled plate. For Allelic Ladders add 1 μ l Ladder to 10 μ l Master Mix. **Note:** The master mix may be altered by adding 0.25 μ l, 0.5 μ l, or 0.75 μ l ILS600 size standard to 9.75 μ l 9.5 μ l, or 9.25 μ l deionized formamide respectively to account for ILS peak height variability.

4.1.4.2 Matrix samples are diluted 1:10 in Nanopure H₂O (a 1:5 dilution may be necessary). 5 μ l of each matrix dye fragment is then added to 480 μ l of deionized formamide (without size standard). Load 25 μ l of the fragment mix into each of 16 wells on the pre-labeled plate, which will include each of the 16 capillaries (e.g. wells A1 through H1 and A2 through H2).

4.1.5 Following sample addition, place a plate septa on the plate and heat denature for ~3 minutes at 95°C. Immediately chill in benchtop cooler (or on ice) for \geq 3 minutes (perform on all sample types - ladders, matrix, controls and samples). **Note:** the plate septa may be cut to cover only those well columns being used on smaller plate runs.

- 4.1.6 Place the sample plate into the plate base and secure the plate retainer clip on top, making sure that no gray is visible through the holes.
- 4.1.7 Place the plate assembly in the instrument in either position A, or position B and close the doors. The plate map on the 'Plate View' window, under 'Run Scheduler' will turn yellow when the plate is in place and has been detected by the instrument. **Note:** the 3130xl will accommodate two plates per run.
- 4.1.8 Prior to running the plate, confirm that dye set F is selected and the correct active calibration for dye set F is set in spectral viewer.
- 4.1.9 Locate the plate record in the 'Plate View' window and highlight it by clicking on it once. With the plate record highlighted, click the appropriate plate map (position A or position B) to link the plate to that specific record. The plate map will turn from yellow to green when it is successfully linked. Verify the correct scheduling of the run in the 'Run View' window. Select a run and confirm that the corresponding wells highlighted in the plate diagram are correct for that run. Make adjustments to the plate record if necessary.
- 4.1.10 Click the green Run Instrument arrow button in the toolbar to start the run. Monitor electrophoresis by observing the run, view, array, or capillaries viewer window. Each injection (set of sixteen samples) will take ~45 minutes. **Note:** to run a duplicate plate record, the plate may need to be unlinked prior to linking the duplicated record. This is done by highlighting the currently linked plate record and clicking 'unlink'.
- 4.1.11 After completion of the spectral calibration run, open the 'Spectral Viewer' window to evaluate the spectral and set the active calibration. Confirm that Dye Set F is selected. Click on individual wells in the plate diagram to see results for each of the sixteen capillaries. For each capillary, verify that four peaks are present in

the spectral profile (upper pane), that the order of the peaks are, from left to right, blue-green-yellow-red, and that the peaks are regular in appearance. Next verify that four peaks are present in the raw data profile (lower pane), that the order of the peaks are, from left to right, red-yellow-green-blue, and that the peak heights are above 750RFU (1,000-4,000 RFU is ideal). If at least 12 capillaries pass, then the calibration should be saved and set as the active calibration.

Note: A minimum of 12 of the 16 capillaries must pass in order to accept a spectral calibration. A passing capillary will be colored green in the plate diagram. Additionally, capillary status may be viewed in the 'Event Log' under 'Instrument Status'. Rerun the spectral calibration as necessary until at least 12 capillaries pass.

4.1.12 After completion of the run finish filling out the 3130xl Injection Log (Form 422B-QC).

4.2 DATA ANALYSIS: GENEMAPPER® ID (GMID)

4.2.1 Data analysis is NOT performed on the instrument computers. Transfer the run folder (including the 3130xl plate record) to an analysis computer using a portable USB drive. After analysis and review are complete, a copy of the run folder and GMID project(s) will be stored on the ISPF network drive. The Run Folder on the instrument computer may be deleted at this point.

Note: prior to data analysis, the appropriate panels and bins must be imported into GeneMapper® ID. Additionally, previously run Macintosh data must first be converted to PC files using the 'Mac to Win' conversion program.

4.2.2 Set up the analysis methods for GMID analysis as follows (analysis methods are created and stored in the 'Analysis Methods' tab in 'GeneMapper Manager'):

'General' Tab: Name the analysis method so that it reflects what the method is (e.g. 3130PP16-150RFU).

'Allele' Tab: Choose the appropriate bin set. Choose 'Use marker-specific stutter ratio if available', and ensure 'minus stutter distances' are from 3.25 to 4.75 for tetra and from 4.25 to 5.75 for penta. All others should be 0.

'Peak Detector' Tab: Advanced Peak Detection Algorithm, partial sizing (80-550 or 600), light smoothing, Local Southern size calling method with baseline window of 51 pfs, min. peak half width = 2, polynomial degree = 3, peak window size = 15, and slope thresholds = 0.

Analysis range may be set to either full or partial and is empirically determined for each run and/or instrument. When using partial range, the start and stop points are determined by a review of the raw data and choosing points that will not include the primer peaks but will cover the size range of 80 to ≥ 500 bases.

Peak Amplitude Thresholds will depend on sample quality. Generally 150 rfu threshold in all colors. Rfu threshold may be raised in Blue, Green and Yellow for Allelic Ladders if necessary. Rfu threshold may be lowered to 100 at the analyst's discretion. With the exception of Red, peaks below 100 rfu are deemed inconclusive (see 4.3.2 RFU Threshold).

'Peak Quality' Tab: The minimum peak height ratio for Heterozygote Balance should be set at 0.5 for database (CODIS) samples. Set the max peak width to 1.5 bp and pull-up ratio to 0.05. The signal level and allele number may be set according to analyst preference and sample type.

'Quality Flags' Tab: The quality flags are only used as a tool to aid in data analysis and review (i.e. to assist in calling attention to potential artifacts or data quality concerns). These flag

settings may be adjusted according to analyst preference and sample quality.

- 4.2.3 Create and store a size standard for GMID analysis, under the 'Size Standards' tab in 'GeneMapper Manager'. Name the size standard so that it reflects what the standard is (e.g. ILS600 80-600).

Data analysis will be performed using the 'Basic or Advanced' size standard. The size standard consists of the following peaks: 60, 80, 100, 120, 140, 160, 180, 200, 225, 250, 275, 300, 325, 350, 375, 400, 425, 450, 475, 500, 550, and 600 (the 60 and 600 peaks may be optionally defined by the analyst).

- 4.2.4 Create a GeneMapper® ID Project:

- 4.2.4.1 From the GMID main menu, select File/Add Samples to project. Highlight the appropriate run folder in the pop-up window and click 'Add to List'. Once the run folder has been copied to the column on the right, click 'Add' to populate the project with the samples in the run folder.

- 4.2.4.2 In the Samples table, for each sample, select the sample type, analysis method, panel, and size standard from the pull-down lists. Ladders must be assigned the sample type of 'Allelic Ladder' for the analysis to occur. In order to use the control concordance quality flag, all controls must be marked appropriately as either 'Positive Control', or 'Negative Control'. All others may be marked as 'Sample'.

- 4.2.4.3 Save the project as the date MMDDYY, followed by the analyst's initials (and any other descriptors that may be necessary). **Note:** the analyzed project will be exported to the run folder at the completion of analysis/review.

4.2.4.4 Analyze the samples by clicking the green Analyze button. If the project has not already been saved, a prompt will appear to enter a project name before analysis will commence.

4.2.5 Evaluate GeneMapper® ID Data:

4.2.5.1 The Raw Data may be reviewed to determine analysis start/stop points, or to identify baseline problems, off-scale data, excessive 'spikes' or other anomalies that may interfere with data analysis and require re-injection or other corrective measures. Expand the run folder located in the navigation pane on the left. Highlight the sample(s) of interest to view the associated sample information, raw data and EPT data. Minimize or highlight the run folder to return to the main project window.

4.2.5.2 Check the 'SQ' (sizing quality) for all samples. A green square indicates that the sample has passed the sizing criteria and need not be manually examined. Examine the size standard of each sample with yellow and/or red 'SQ' to confirm correct assignment of fragment sizes. Highlight the sample(s) of interest and click the Size Match Editor button. If necessary, adjust the peak assignments by right clicking on a peak and deleting, adding, and/or changing values. If all peaks are correctly labeled but the quality score is below 1.0 (may be checked by choosing Tools/Check Sizing Quality), click the 'Override SQ' button to set the SQ to 1.0. Once all edits have been made, click 'OK' to save the changes and close the Size Match Editor (clicking 'Apply' saves the changes but leaves the Size Match Editor open). These samples are ready for reanalysis in the project window.

Note: Data may still be deemed acceptable without the ILS 60 and/or 600 bp peaks present. If additional peaks are assigned because of bleed-through of TMR peaks (typically Amelogenin peaks), the scientist may choose an Analysis Method, with an increased rfu threshold for the red channel to prevent these peaks from being detected, if desired.

- 4.2.5.3** Examine the blue, green, and yellow allelic ladders. Check that correct allelic assignments were made.

Note: GMID automatically averages all valid ladders in a run for genotyping. Genotypes are assigned by comparing the sizings of unknown alleles from samples with the sizings of known alleles contained within the averaged allelic ladders of each locus. A ladder(s) may be omitted from analysis by deleting it from the main project window prior to analysis.

- 4.2.5.4** Data may be examined in various combinations of colors and/or tables to identify bleed-through, spikes, stutter, -A, off-ladder variants, etc. Sample Plots viewed from the 'Samples' tab/window, allows all loci in a given color(s) to be viewed simultaneously. The Sample Plots view from the 'Genotypes' tab/window; however, allows loci to be viewed individually (more than one locus can be viewed at the same time by adjusting the number of panes displayed).

- 4.2.5.5** GeneMapper® ID includes a series of quality flags (POVs) to alert the analyst of potential sample quality concerns. A green square indicates that sample data has passed all of the quality checks, but yellow or red indicate that the data has a problem with one or more of the quality checks. A yellow or red flag does not

necessarily mean that the data is bad or unusable and the flags are not to be relied on solely. The analyst may choose to use the PQVs, in combination with manual data examination to aid in the identification of bleed-through, spikes, stutter, off-ladder variants, -A, etc. Once the data has been evaluated and deemed acceptable, the analyst may choose to override the yellow or red Genotype Quality (GQ) flag by right clicking on the flag in the Genotypes Sample Plots view. **Note:** overriding the GQ flag will cause all other flags to turn from the original color to gray.

4.2.5.6 All negative controls (including reagent blanks) should be examined to verify that each displays a relatively flat baseline in blue, green and yellow.

4.2.5.7 Review all samples (including positive controls) for the above listed 'artifacts' and evaluate: peak height and shape, spectral quality, and individual sample profiles. Compare each sample with the allelic ladder(s) and examine for off-ladder or microvariants, signals that were too low to be genotyped and assignment of genotypes to stutter peaks (or minor peaks that may have been subtracted as 'stutter', etc.)

4.2.5.8 Reanalyze individual samples with different Analysis Methods, as necessary if the rfu cut-off will need to be changed.

4.2.5.9 Edit peaks as necessary, by right clicking on the peak label and selecting 'add allele', 'delete allele' or 'rename allele'. The allele should be labeled, at minimum with the allele call, however the analyst may select up to four allele labels, including peak height and size, from the 'Plot Settings Editor' window. **Note:** labels added to artifact peaks,

such as spike, pull-up, etc. will appear in the Genotypes table as an additional allele.

4.2.5.10 Samples demonstrating an off-ladder (< or > smallest or largest ladder allele, respectively), tri-allele, or microvariant (alleles with incomplete repeats) allele(s) should be re-analyzed for verification. Microvariants will be labeled and reported as "X.Y" (where X is the number of complete repeats and Y is the number of basepairs of the incomplete repeat). Off-ladder will be reported as > or < the largest or smallest ladder allele, respectively. **Note:** the nomenclature for upload to NDIS may necessitate a change in allele designation.

4.2.5.11 GMID automatically flags off-scale (camera saturation) data. This data may still be acceptable if it is limited to a few or a single peak and the overall data for that sample is of good quality (see 4.3.2.3).

4.2.5.12 Export an allele/genotypes table to Excel and save it in the run folder. The table will be printed for the batch record/file. The table should also be exported as a .cmf file for CODIS import. To create a .cmf file, the specimen category must be assigned and the export fields set in the 'CODIS Export Manager' under tools in the main menu.

4.3 STR INTERPRETATION GUIDELINES

4.3.1 CONTROLS

4.3.1.1 The purpose of a **REAGENT BLANK** (RB) is to determine if the reagents used for DNA extraction/isolation were contaminated with human DNA and as a method for monitoring facility

decontamination. In GeneMapper® ID, peaks above threshold should only appear in the CXR (red dye) lane, corresponding to the ILS600 size standard. Electropherograms for the blue, green and yellow dyes should show a relatively flat baseline throughout the range (discounting primer signal, fluorescent 'spikes' or CXR bleed-through). If detectable signal, with characteristic 'peak' shape is visible in the electropherogram of a reagent blank and does not disappear upon re-injection, results for all associated samples may be deemed inconclusive (close examination at 50 rfu is performed on all samples to examine for presence of any alleles seen in the RB). Data may be deemed acceptable if contamination is 'isolated' to the RB. The reagent blank should be treated the same as the least concentrated DNA sample in terms of volume and amount amplified, injection time/amount, and analysis threshold. Additionally, the reagent blank will be reamplified with samples from the set if any of the amplifications conditions are more sensitive than the original.

4.3.1.2 The purpose of the **POSITIVE AMPLIFICATION CONTROL** (9947A DNA supplied with the PP16 HS kit) is to assess the amplification process, ensuring that adequate sample amplified simultaneously would produce an appropriate signal. All expected alleles (see below) must be detected, using standard parameters or all of the samples associated with amplification may be deemed inconclusive. Data may be deemed acceptable if all alleles are present (though some are below 150-rfu threshold) **AND** the other positive control (Extraction Control) appears as expected (i.e. the problem is confined to the 9947A sample).

LOCUS	GENOTYPE	LOCUS	GENOTYPE
D3S1358	14, 15	TH01	8, 9.3
D21S11	30, 30	D18S51	15, 19
Penta E	12, 13	D5S818	11, 11
D13S317	11, 11	D7S820	10, 11
D16S539	11, 12	CSF1PO	10, 12
Penta D	12, 12	AMELOGENIN	X, X
vWA	17, 18	D8S1179	13, 13
TPOX	8, 8	FGA	23, 24

4.3.1.3 The purpose of the **NEGATIVE AMPLIFICATION CONTROL** is to determine if any human DNA contamination occurred in the process of amplification set-up (or beyond that point) and as another method of monitoring facility decontamination. In the GeneMapper® ID electropherograms, peaks above threshold should only appear in the CXR (red dye) lane, corresponding to the ILS600 size standard. Electropherograms for the blue, green and yellow dyes should show a relatively flat baseline throughout the range (discounting primer signal, fluorescent 'spikes' or CXR bleed-through). If detectable signal, with characteristic 'peak' shape, is visible in the electropherogram of a negative amplification control and does not disappear upon re-injection, results for all of the samples associated with that amplification will be examined for the presence of the same peak(s).

If extraneous peaks appear only in this control, the data for other samples associated with that amplification need not be deemed inconclusive. This occurrence should be documented and the

scientist's determination (and basis for it) documented in the note packet.

4.3.1.4 The purpose of an **EXTRACTION CONTROL** sample is primarily to assess correct genotyping, however, it does take measure of all of the steps in the analytical process from extraction through allele designation. The extraction controls consist of previously typed buccal samples on FTA cards prepared in batches. An extraction control must be run on every sample plate. The reviewing scientist will complete an Extraction Control Check Form (Form 212-BI) for verifying correct genotype(s). A copy of this form will be included in each CODIS Data file. Failure of the extraction control, if isolated to that sample, will not deem other samples inconclusive.

4.3.2 RFU THRESHOLD:

4.3.2.1 For Offender database samples, a minimum of 100 rfu should be achieved for data acceptance (the threshold may be lowered to 50 RFU for ILS peaks as needed). If necessary, go back in the process as follows: repeat injection (changing injection time; 3-10 seconds allowable range), or perform re-analysis (i.e., changing amount of amplified product added for fragment analysis), or re-amplification (increase DNA template), or re-extraction.

4.3.2.2 Peaks below the analysis threshold (based on data obtained and signal/baseline) will not be interpreted but should be noted as being present in the notes (eg. on the table of results).

4.3.2.3 Peaks marked as off-scale in GeneMapper® ID (indicating camera saturation) will not be interpreted if multiple peaks are affected and if it causes excessive artifacts (e.g. split peaks, increased stutter, pull-up, etc.) which interfere with data interpretation (see 4.2.5.11). If the overall quality of the data is not acceptable, the sample must be diluted, re-injected (3-10 seconds), reanalyzed (decrease the amount of amplified product added) or re-amplified (decrease DNA template) as deemed appropriate by the scientist.

4.3.2.4 Multiplex amplification kits are designed so that heterozygous loci in single-source samples generally demonstrate relatively balanced peak heights [typically $\geq 70\%$ peak height ratio (p hr)]. Some samples, although single-source, may at times demonstrate greater imbalance due to degradation, stochastic effects, primer binding site mutations, preferential amplification, etc. Peak height ratios for these loci ($< 50\%$ p hr for database samples) will be flagged in GeneMapper® ID.

4.3.3 EXTRA PEAKS (NON-MIXTURES)

4.3.3.1 PCR amplification of STR loci typically produces a minor product peak one core repeat unit shorter than the main allele peak (n-4 for tetranucleotide loci and n-5 for pentanucleotide loci). This minor peak is referred to as the **stutter** peak. Percent stutter generally increases with allele length and does not change significantly with the quantity of input DNA (peak heights within ~150-4500 RFU). The measurement of percent stutter may be unnaturally high for main peaks that are off-scale or due to problems with spectral performance and can be corrected by diluting (or reamplifying less DNA) the

sample and/or running a new spectral. Loci stutter values are listed in Appendix A of the Biology Quality Manual to assess potential contribution to peaks in stutter positions. N+4 stutter may also be seen on occasion.

4.3.3.2 Electronic or fluorescent **spikes** are random events that produce generally spike-shaped peaks in most or all dye colors at the same location (equivalent bp size) within a single injection. Peak heights usually vary between dye colors for a given spike. These anomalies are generally not reproducible and will typically be eliminated upon reinjection. If the spike is above the analysis threshold and falls within an allelic range that could interfere with either computer analysis or scientist's analysis, the scientist will label the spike in the GeneMapper® ID software.

4.3.3.3 Dye "**blobs**" are anomalies that typically occur in the same approximate location in multiple injections and do not always disappear upon reinjection. Blobs generally look like broad or irregular peaks and may occur in a single color or multiple colors at the same approximate location but can vary in height. The blob should be labeled on the electropherogram (in GMID) if it falls within a diagnostic region and is of significant size to potentially interfere with analysis.

4.3.3.4 **Bleed-through** or pull-up peaks are a result of the spectral not correcting for all of the spectral overlap (most common with the PowerPlex 16 HS kit from yellow into red) and may be increased due to off-scale peaks. These pull-up peaks are in the same location (same bp size) as peaks in another color(s) and are easily

recognized. The presence of bleed-through should be labeled on the corresponding electropherogram (in GMID) if it falls within a diagnostic region and is of significant size to potentially interfere with analysis. If excessive bleed-through occurs in a color other than red, and is not due to off-scale data, a new spectral may need to be performed.

4.3.3.5 Taq Polymerase can catalyze the addition of a single nucleotide (predominantly adenosine) to the 3' ends of double stranded PCR product, resulting in product one base pair longer than the actual target sequence (+A). Amplification parameters include a final extension time, so that the reaction is driven to full A addition (i.e. all product is +A). **Split-peaks** may occur as a result of incomplete A addition and appear as a single allele represented by two peaks one base pair apart (-A and +A). This can occur when the amount of template DNA is too great (overloaded sample). In this instance, Taq is unable to add the A nucleotide to the entire amount of product generated in the time allotted. These samples will typically contain off-scale data as well. Split peaks can be alleviated by diluting the sample or by incubating samples at 60°C for an additional 45 minutes, followed by dilution, prior to reinjection. It may be necessary to re-amplify the sample with less template DNA.

4.3.4 MIXTURES

4.3.4.1 If, after the elimination of possible stutter and/or bleed-through, a profile shows at least 3 peaks at 2 or more loci, this is strong evidence of a mixture.

4.3.4.2 Loci that demonstrate only two alleles but have a heterozygous ratio of <50% may also be indicative of a mixture. However, if data are obtained from multiple loci, a scientist should expect to see this or other mixture indications (> 2 alleles) at additional loci.

4.3.4.3 Samples showing evidence of a mixture are to be deemed inconclusive, as offender database samples should be from a single source. An effort should be made to determine the source of the contamination, if possible and to determine if other samples may be affected.

5.0 Comments:

- 5.1 The 310 POP4 Polymer and the 3130/3130x1 POP4 Polymer are different and are not to be used interchangeably.
- 5.2 The 3130x1 Data Collection Software does not allow the entry of spaces or dashes in titles, sample names, etc. An underscore must be used in place of spaces when entering information.
- 5.3 Low level contamination (few peaks below 50 RFU) may be seen periodically with the Biomek 3000 extraction. The contamination may be disregarded if it is isolated to only a few wells on the plate, all peaks are below 50 RFU, and extraneous peaks do not interfere with allele calling in the samples.
- 5.4 Additionally, low level carry-over may be seen occasionally on the 3130x1. The carryover may also be disregarded in samples and negative controls, so long as the extraneous peaks are low level, do not interfere with allele calling, and the peaks can be attributed to the previously injected sample in the corresponding capillary.

DRIFTCON FFC: TEMPERATURE VERIFICATION**1.0 BACKGROUND:**

Successful DNA amplification is dependent on consistent thermal cycling parameters and achieving proper heating and cooling to facilitate the various steps necessary for PCR. It is necessary to employ a method of monitoring the thermal cycler performance and verify that the correct temperature has been achieved in order to have confidence in the amplification process. A variety of temperature probes have been developed to test the temperature of specific wells within the different thermal cycler instruments. However, most of these have proven to be cumbersome, time consuming, and not amenable to testing the newer real-time PCR instruments. The Driftcon FFC is able to test multiple temperatures within a short period of time and can be used on the thermal cyclers, as well as real-time instruments.

Driftcon Operations Manual, version 1.1

Driftcon Quick Start Guide, version 1.4

2.0 SCOPE:

To provide a reliable method for verifying the performance of laboratory thermal cyclers, to include real-time PCR instruments.

3.0 EQUIPMENT/REAGENTS:

Computer with Driftcon Software
Driftcon Hardware Module and Cables
Driftcon FFC Fixture
Driftcon Smart Card
ABI 7500 FFC Adaptor
Cork Leveler

4.0 PROCEDURE:

4.1 9700 THERMAL CYCLER VERIFICATION

- 4.1.1 Set up the Driftcon by connecting the hardware module to the computer and FFC fixture/probe plate with the provided cables.
- 4.1.2 Insert the smart card into the Driftcon hardware module until it stops.
- 4.1.3 Turn on the Driftcon computer. Open the Driftcon Software and login with the appropriate user name and password.
- 4.1.4 Turn on the 9700, place the fixture into the plate with the cable facing out, and select/start the 'Driftcon' protocol. **Do Not** close the 9700 lid as it will damage the fixture and void the warranty.
- 4.1.5 In the Driftcon software, choose the instrument to be tested. If the list of instruments does not appear automatically, click start and follow the wizard prompts.
- 4.1.6 Choose the Driftcon default protocol and make sure the steps match those in the 9700 protocol. Leave the humidity, pressure, and temperature blank. Enter notes as needed (these will appear on the final report).
- 4.1.7 Start the test by clicking the check mark in the lower right corner. Make sure the 9700 has heated completely and the protocol started before starting the test. The protocol will complete in ~25 min.
- 4.1.8 The analyst need not be present once the test has started; however, additional views/information are available during the run by pressing the F6 key for a surface view (heat dispersment), F7 for measurement values (shown in real time), and F8 for measurement locations (probe wells).

4.1.9 Once the run has completed, print the report and place it in the QC binder.

4.2 7500 REAL-TIME INSTRUMENT VERIFICATION

4.2.1 Set up the Driftcon by connecting the hardware module to the computer and FFC fixture/probe plate with the provided cables.

4.2.2 Insert the smart card into the Driftcon hardware module until it stops.

4.2.3 Turn on the Driftcon computer. Open the Driftcon Software and login with the appropriate user name and password.

4.2.4 Turn on the 7500 computer, login with the appropriate user name and password, and open the 7500 SDS software.

4.2.5 Turn on the 7500. Remove the tray from the plate loader and replace it with the FFC adaptor. Place the fixture into the adaptor with the cable facing out. Place the cork leveler on top of the fixture and close the plate loader.

4.2.6 In the 7500 SDS software, choose **File > New** and select **Absolute Quantitation** for Assay, **96-Well Clear** for Container, and **Driftcon** for Template.

4.2.7 Save the plate document as a .sds file with the appropriate plate name and open the Instrument tab.

4.2.8 In the Driftcon software, choose the instrument to be tested. If the list of instruments does not appear automatically, click start and follow the wizard prompts.

4.2.9 Choose the Driftcon default protocol and make sure the steps match those in the 7500 protocol. Leave the humidity, pressure, and temperature blank. Enter notes as needed (these will appear on the final report).

4.2.10 Start the 7500 run and Driftcon test. Start the test by clicking the check mark in the lower right

corner. Make sure the 7500 has heated completely and the protocol started before starting the test. The protocol will complete in ~25 min.

4.2.11 The analyst need not be present once the test has started; however, additional views/information are available during the run by pressing the F6 key for a surface view (heat dispersment), F7 for measurement values (shown in real time), and F8 for measurement locations (probe wells).

4.2.12 Once the run has completed, print the report and place it in the QC binder.

5.0 Comments

5.1 The second page of the report contains information regarding the number of measurements for each probe. The number should be approximately double the protocol time. If significantly higher or lower, the instrument may be heating too slowly or too quickly. The percentage should be 100% for each probe. If a percentage is lower than 100%, the probe may not have been in the well and measurements missed. Adjust the probes and run the protocol again.

5.2 Pass/better than specifications/fail data may be found beginning on page four, for each temperature. The measured value, along with the target specifications are shown in this section.