

IDAHO DEPARTMENT OF LAW ENFORCEMENT

BUREAU OF FORENSIC SERVICES

**DNA LABORATORY
QUALITY ASSURANCE MANUAL**

Date

4/15/98

Laboratory Manager

Richard D. [Signature]

Technical Leader

[Signature]

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

TABLE OF CONTENTS

TABLE OF CONTENTS 1

INTRODUCTION..... 2

I. MISSION STATEMENT..... 2

II. PERSONNEL..... 4

III. DOCUMENTATION..... 6

IV. METHOD VALIDATION 8

V. QC PROCEDURES – EQUIPMENT, REAGENTS, MATERIALS, AND PHYSICAL FACILITIES.. 16

VI. PROCEDURES FOR EVIDENCE HANDLING AND SECURITY 20

VII. ANALYTICAL PROCEDURES 22

VIII. CASE NOTES..... 26

IX. REPORTS..... 27

X. TECHNICAL REVIEW PROCEDURES..... 28

XI. PROCEDURES FOR THE RELEASE OF CASE REPORT INFORMATION..... 31

XII. TRAINING AND PROFICIENCY TESTING..... 31

XIII. AUDITS..... 34

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

INTRODUCTION

The goal of the DNA Laboratory is to provide DNA analysis and interpretation of the resulting data in case evidence samples. The quality assurance program outlined in this manual will ensure that the work product is of the highest quality, integrity, and reliability.

This manual outlines the quality assurance program of the Idaho Department of Law Enforcement Bureau of Forensic Services DNA Laboratory. Quality assurance is a process that will continually evolve with the Laboratory. Therefore, the QA manual will be periodically modified to reflect current professional standards and our own experience.

Other operational and analytical procedures related to quality assurance are contained in the BUREAU OF FORENSIC SERVICES LABORATORY QUALITY MANUAL, BUREAU OF FORENSIC SERVICES LABORATORY POLICY MANUAL, BUREAU OF FORENSIC SERVICES LABORATORY CHEMICAL HYGIENE PLAN, DNA PROCEDURES MANUAL, DNA REAGENT PREPARATION NOTEBOOK, DNA QC DATA NOTEBOOK, and IDAHO DEPARTMENT OF LAW ENFORCEMENT SAFETY MANUAL.

I. MISSION STATEMENT

The mission of the Bureau of Forensic Services is to provide quality and impartial scientific analysis, testimony, crime scene investigation, education, and research to the criminal justice system. To further that mission, the DNA Laboratory will provide DNA analysis and interpretation of case evidence samples to our client agencies and ensure the quality, integrity and reliability of our work through an ongoing quality assurance program. The quality assurance program will provide for us, the users of our services, and the triers of fact, the assurance that our work product meets recognized standards.

1.1 THE QUALITY ASSURANCE PROGRAM ADDRESSES:

- Sample preservation and chain of custody
- Laboratory QC procedures
- Laboratory protocols and interpretation guidelines
- Casework documentation and report writing
- Analyst training and proficiency

1.2 THE FUNCTION OF THE QUALITY ASSURANCE PROGRAM IS TO ENSURE THAT:

- The entire DNA typing procedure is operating within established performance criteria.
- DNA testing and reporting procedures are monitored by means of a quality control (QC) program.
- Problems are noted and corrective actions are taken.

The primary sources upon which this document is based are the TWGDAM "Guidelines for a Quality Assurance Program for DNA Analysis" and the ASCLD "Guidelines for Forensic Laboratory Management Practices". The Idaho Bureau of Forensic Services DNA Laboratory has demonstrated that it meets recognized standards of quality by becoming accredited by the American Society of Crime Laboratory Directors Laboratory Accreditation Board (ASCLD/LAB).

1.3 AUTHORITY AND ACCOUNTABILITY—SEE ORGANIZATIONAL CHART IN BUREAU POLICY MANUAL (TWGDAM 1.3.1)

1.3.1 Laboratory Manager

The Laboratory Manager is responsible for overall quality assurance, for monitoring QC procedures, and for assigning and recording the results of proficiency tests. With regard to the DNA Laboratory, the Laboratory Manager will proceed on the advice and recommendations of the DNA Technical Leader. In particular, the Laboratory Manager will:

- 1.3.1.1 Review and approve new QA guidelines and QC procedures and modifications to existing procedures and guidelines.**
- 1.3.1.2 Assign proficiency tests, collect the results of the analysis and review, and update files to reflect current status.**
- 1.3.1.3 Review the QC records on a regular basis to ensure that the designated tasks are being performed as outlined in the QA and QC manuals. (TWGDAM 1.3.3)**
- 1.3.1.4 Review and approve corrective action proposed or taken regarding analytical or interpretive problems detected by the review process or DNA Laboratory deficiencies detected by QA/QC process.**

1.3.2 DNA Technical Leader (TWGDAM 1.3.1, 1.3.3)

In the current role of consultant, the DNA Technical Leader will advise the Laboratory Manager on all aspects of the Laboratory's DNA program, including procedures, protocols, QA guidelines, QC procedures, Analyst training and proficiency testing, research and development, and case review. In particular, the Technical Leader will, in collaboration with the Laboratory Manager:

- 1.3.2.1 Review and approve modifications to analytical protocols used for casework.**
- 1.3.2.2 Review and approve new QA guidelines and QC procedures and modifications to existing procedures and guidelines.**
- 1.3.2.3 Represent, or designate a substitute to represent, the DNA Laboratory in meetings and policy discussions where Bureau-wide Quality Assurance matters are discussed. (TWGDAM 1.3.3)**
- 1.3.2.4 Review proficiency test results and verify their successful completion and report on such to the Laboratory Manager.**
- 1.3.2.5 Review each case file for technical accuracy and support of conclusions by the data.**
- 1.3.2.6 Review the QC records periodically to ensure that the designated tasks are being performed as outlined in the QA and QC manuals. (TWGDAM 1.3.3)**
- 1.3.2.7 Review and approve corrective action proposed or taken regarding analytical or interpretive problems detected by the review process or DNA Laboratory deficiencies detected by QC process.**
- 1.3.2.8 Actively seek new knowledge and stay abreast with developments in the field. (TWGDAM 2.2.3.4)**

1.3.3 DNA Analyst (TWGDAM 1.3.2)

Each DNA Analyst has primary responsibility for the quality of his/her work. Analysts are expected to:

- 1.3.3.1** **Appropriately evaluate samples before proceeding to DNA analysis.**
- 1.3.3.2** **Understand and follow established procedures for DNA analysis and interpretation in casework, including the proper use of controls.**
- 1.3.3.3** **Interpret results in an objective manner and provide a complete and accurate report of the findings, including a qualitative and a quantitative (when appropriate) statement summarizing the results.**
- 1.3.3.4** **Discuss any analytical or interpretational problems with a supervisor, including suggestions for resolution of the problem. (TWGDAM 1.3.3)**
- 1.3.3.5** **Actively seek new knowledge and remain current with appropriate literature. (TWGDAM 2.2.4.4)**
- 1.3.3.6** **Participate in peer review, proficiency testing and training as assigned. (TWGDAM 1.3.3)**
- 1.3.3.7** **Perform all designated QC tasks, assisted as needed by support staff. Typically these may include:**
 - (a) Routine calibration checks of equipment (freezers, incubators, thermal cyclers)
 - (b) Routine quality assurance of reagents and materials used for DNA analysis
 - (c) Preparing and issuing QC samples

II. PERSONNEL

2.1 JOB DESCRIPTIONS AND QUALIFICATIONS (TWGDAM 2.1)

The general knowledge, skills and abilities required for the technical classifications within the DNA Laboratory are set forth in the Bureau of Forensic Services job descriptions for DNA Technical Leader, Principal Criminalist, Senior Criminalist, Chemist, and Laboratory Technician.

In addition to the specialized knowledge listed below, employees shall read and abide by the IDAHO DEPARTMENT OF LAW ENFORCEMENT SAFETY MANUAL, the Bureau's policies which include the provisions of the section on Laboratory safety, and the BUREAU OF FORENSIC SERVICES CHEMICAL HYGIENE PLAN. (TWGDAM 11.1, 11.2, and 11.4)

Specific additional qualifications are required for individuals to be assigned to analytical work at various levels in the DNA Laboratory. These requirements are designed to ensure that staff have training, education and proficiency commensurate with their duties and that ASCLD/LAB standards and TWGDAM guidelines are met.

2.1.1 Technical Leader

It is the responsibility of the DNA Technical Leader to review interpretation and technical content of casework reports and other analytical work, and to oversee research and protocol development, the DNA training program, and the development and implementation of the QA/QC program for the DNA Laboratory. The DNA Technical Leader may serve as a second reader of DNA testing results. In the role of consultant, the Technical Leader shall perform the foregoing functions and

advise the Laboratory Manager regarding the status of the above, including any recommendations for change or corrective action.

2.1.1.1

Qualifications

To be qualified to perform technical review, the DNA Technical Leader must have knowledge of forensic DNA testing as indicated by:

- (a) A minimum of a BA/BS degree or its equivalent in a biological, chemical or forensic science (TWGDAM 2.1.3.1) and have credited coursework in biochemistry, genetics and molecular biology (molecular genetics, recombinant DNA technology), which provide a basic understanding of the foundation of forensic DNA analysis (TWGDAM 2.2.3.1), and
- (b) Training in the fundamentals of forensic biology. There must be documented training in DNA analysis that includes the methods, procedures, equipment and materials used in forensic DNA analysis and their limitations and applications (TWGDAM 2.2.3.2), and
- (c) A minimum of 2 years of experience as a forensic biology analyst and meet or exceed the experience requirement listed for DNA Analysts for analysis of samples prior to casework (TWGDAM 2.2.1.2, 2.2.3.3.), and
- (d) Demonstrate competence in review and interpretation of DNA typing results in a variety of samples typical of evidentiary material (TWGDAM 2.2.1.4). This may be established by successful completion of proficiency tests prior to assuming the duties of Technical Leader (TWGDAM 2.2.1.3). For new methods implemented after the Technical Leader is no longer performing active lab analyses, competency in review and interpretation may be demonstrated by successful second reading and review of a demonstration case.

2.1.2

DNA Analyst

It is the responsibility of the DNA Analyst to perform all aspects of criminal casework on evidentiary materials, including evaluation and preparation of the evidence, DNA typing, interpretation of analytical results, preparation of case reports and expert testimony. Casework Analysts may serve as second-readers for another Analyst's DNA casework.

2.1.2.1

Qualifications

An Analyst qualified to analyze and interpret DNA in evidentiary materials must have knowledge of forensic DNA testing as indicated by:

- (a) A minimum of a Bachelor's degree in a biological, chemical, or forensic science-related area, and have credited college coursework in biochemistry, genetics and molecular biology (molecular genetics, recombinant DNA technology) or other related subjects, that provides a basic understanding of the foundation of forensic DNA analysis (TWGDAM 2.2.4.1), and
- (b) Training in the fundamentals of forensic biology and training in DNA analysis that includes hands-on DNA Laboratory work in the procedures performed and their applications and limitations (for example, CCI or Perkin-Elmer courses, and/or formalized in-house training). (TWGDAM 2.2.4.2) , and
- (c) One year or more of forensic biology experience including six months of forensic DNA Laboratory experience. This experience shall include successful analysis and interpretation of DNA typing results of a range of samples typically encountered in forensic casework prior to independent casework analysis using DNA technology (TWGDAM 2.2.4.3) , and
- (d) Demonstration of competence in analysis of DNA and interpretation of DNA typing results as established by proficiency testing. (TWGDAM 2.2.1.3).

2.1.3 Support staff

It is the responsibility of the Support Staff to conduct a variety of duties including preparation of reagents, calibration of equipment, conducting other QC checks, ordering supplies, computer data entry, and receipt and logging of evidence samples. After appropriate training, support staff may be responsible for assisting Analysts in analytical work. Analytical support staff will not perform independent casework, interpret analytical results or provide second readings on casework samples, nor will they testify regarding any of the above.

2.1.3.1

Qualifications:

To be qualified to assist in analytical functions in the DNA Laboratory, Support staff must have:

- (a) A Bachelor's Degree in science (TWGDAM 2.2.5.1) , and
- (b) In-service training by a qualified Analyst and demonstration of competence in the DNA Laboratory procedures performed (TWGDAM 2.2.5.1)

III. DOCUMENTATION

3.1. PURPOSE OF DOCUMENTATION

Documentation of all significant aspects of the DNA analysis procedure and other aspects of the DNA Laboratory operation related to the reliability and interpretation of analytical results is necessary to:

- Support the scientific conclusions in the DNA Laboratory report.
- Permit supervisory/peer review of the work product.
- Allow re-evaluation of the data by outside scientific observers.
- Provide a foundation for introduction of the work product into a court of law.
- Provide an audit trail by which management can demonstrate and verify the continued quality of the DNA Laboratory's work.

3.2 TEST METHODS AND PROCEDURES

Protocols for each analytical method routinely used in casework analysis are maintained in the DNA PROCEDURES MANUAL. These protocols will include the standards and controls required and the date(s) on which the protocol was approved by the Technical Leader and Laboratory Manager. Previous protocols must be archived. (TWGDAM 3.1)

3.3 ANALYTICAL NOTES, PROFICIENCY TESTING AND CASE REPORTS

(See *Section VII* for form and content of notes and reports and *Section VIII* for additional requirements regarding proficiency testing)

- "Note-keeping" is defined as "documentation of work as performed".
- "Documentation" means the recording of details and observations sufficient to both support conclusions and duplicate the experiment at another time.
- "As performed" means both that the record is made contemporaneously with the work as done and that what is recorded is what was actually done.
- Notes will be made as appropriate for all analytical work. Notes are the property of the DNA Laboratory, not the Analyst, and are not to be removed from the DNA Laboratory without express

permission of a supervisor. The type of file in which analytical notes are kept depends on the kind of work performed.

3.3.1 Research notes (TWGDAM 3.10)

Records of research, training, validation, and other non-casework projects are kept in notebooks or files as appropriate. Each major experiment will be assigned an experiment number that will appear on all notes related to that experiment. Research records are archived indefinitely. (TWGDAM 3.16)

3.3.2 Casework analysis (TWGDAM 3.4)

A case file with a unique BFS Laboratory number is established for each case in which evidence has been received in the DNA Laboratory. Upon completion of the analysis and report, all analytical notes, chain of custody information, case contact information, and reports are entered into this file. The file will also contain documentation of technical review. Case files are archived according to the established BFS schedule. (TWGDAM 3.16)

3.3.3 Proficiency Testing (TWGDAM 3.8)

A folder is created for each proficiency test; these are filed in the PROFICIENCY TEST FILE and retained according to BFS policy. The PROFICIENCY TEST FILE contains all analytical data (notes, photos, run sheets, etc.) generated in the analysis and a statement of the Analyst's conclusion. It will also contain the second reader's notes and the technical review checksheet. For external tests, it will also contain the summary report of the test provider regarding that particular test. In any situation where the results of the test are not satisfactory, the file will also include documentation of any corrective action taken. A log of each individual Analyst's proficiency tests will be kept by the Laboratory Manager in the individual's personnel file.

3.4 DOCUMENTATION OF REAGENT CHECKS, EQUIPMENT CALIBRATION AND MAINTENANCE LOGS (TWGDAM 3.3 and 3.7)

Equipment maintenance records are entered in the notebook established for each piece of equipment. Calibration records and the required routine quality control checks, including the record of commercial lot numbers and expiration dates, are filed in the DNA QC DATA NOTEBOOK. Periodic equipment checks, and freezer and refrigerator temperature logs, are also kept in the DNA QC DATA NOTEBOOK. These records will be archived for a minimum of five years. (TWGDAM 3.16) Specific instructions for checks and calibrations are found in the DNA QC DATA NOTEBOOK.

3.5 PERSONNEL TRAINING AND QUALIFICATIONS RECORDS (TWGDAM 3.9)

An assessment of the education, training and experience of each new employee is conducted at the time of hire by review of transcripts and certificates, interview of former employers and, where appropriate, review of prior work or proficiency test results.

3.5.1 A training file will be established for each member of the DNA Laboratory as part of his/her personnel record. This file will include:

(a) A current CV listing relevant education, in-service training and other qualifications.

(b) Records of in-service training courses, seminars, and continuing education provided by the DNA Laboratory.

(c) A list of completed proficiency tests.

(d) A current career development plan, including a plan for in-service training. It is required that the Technical Leader and DNA Analysts stay abreast of developments in the field of DNA typing (TWGDAM 2.2.3.4 and 2.2.4.4).

3.5.2 Performance evaluations

Performance evaluations will be made of each employee by a supervisor according to the established schedule. The performance evaluation will include a review of the training file, proficiency tests, and work plan by the supervisor and employee and updating of the file as needed. The required continuing education will be documented. (TWGDAM 2.2.2)

3.6 EQUIPMENT INVENTORY (TWGDAM 3.13, 5.1.1)

An inventory of all DNA Laboratory capital equipment is kept in the administrative files. This inventory is updated as new equipment is received.

3.7 SAFETY MANUALS (TWGDAM 3.14, 3.15, and 11.3)

Laboratory procedures and policies related to facility safety, emergency response and evacuation, illness and accident prevention, hazard communication, and blood borne pathogens are all contained in the IDAHO DEPARTMENT OF LAW ENFORCEMENT SAFETY MANUAL, the Bureau's POLICY MANUAL in its section on laboratory safety, or the LABORATORY CHEMICAL HYGIENE PLAN. Also available to all staff is a binder or CD ROM containing MSDS's for all materials used in the DNA Laboratory. (TWGDAM 11.3)

IV. METHOD VALIDATION

Validation is the process used by the scientific community to assess the accuracy, reproducibility and reliability of a procedure, to determine the conditions under which the test results accurately reflect the nature of the sample, and to determine the limitations of the method. The validation process identifies the critical aspects of the procedure that must be carefully controlled and monitored.

Testing procedures used in the DNA Laboratory will be validated according to the TWGDAM Guidelines before being implemented in casework. Developmental validation work may be done by an outside DNA Laboratory developing the procedure or internally if the procedure was developed in-house. Each locus used will have been characterized as to inheritance, chromosomal location, and nature of the polymorphism. Validation studies will address reproducibility, mixed specimens, environmental effects, substrate effects, human specificity, detection thresholds, and population frequencies.

4.1 INTERNAL VALIDATION

Prior to implementing a new analysis procedure which has been developed and validated in another DNA Laboratory, the DNA Laboratory will conduct internal validation tests using known samples and non-probative evidentiary samples.

If a substantive modification is made in the procedure by the DNA Laboratory, the modified procedure will be compared with the original procedure on identical samples. Supervision and approval over the method modifications shall be the responsibility of the Technical Leader.

4.2 PCR DQA1 VALIDATION (TWGDAM 4.1.1, 4.1.2, and 4.1.4)

The following section lists validation studies applicable to the PCR-based analysis of the HLA DQA1 locus. References are provided to studies conducted by other laboratories.

4.2.1 General Developmental Validation (TWGDAM 4.1.5)

Walsh, P.S., et al. (1991) "Report of the Blind Trial of the Cetus AmpliType HLA DQ α Forensic Deoxyribonucleic Acid (DNA) Amplification and Typing Kit". Journal Forensic Sciences 36 (5):1551.

Comey, C.T., and Budowle, B. (1991) "Validation Studies on the Analysis of the HLA DQ α Locus Using the Polymerase Chain Reaction". Journal Forensic Sciences 36 (6):1633.

Keel, A., Sims, G., and Buoncristiani, M. (1992) "A Collaborative Study of DQ α Typing by PCR on Sexual Assault Evidence". (Abstract), Journal Forensic Science Society 32 (3):270.

Rudin, N., et al. (1992) "A Systematic Study of the Effect of Various Environmental Abuses on RFLP and PCR Analysis of Forensic Samples". (Abstract), Journal of Forensic Science Society 32 (3):274 and in Proceedings from The Third International Symposium on Human Identification. Promega Corp. 1992:421.

Saiki, R., et al. (1986) "Analysis of Enzymatically Amplified B-globin and HLA-DQ α DNA with Allele-Specific Oligonucleotide Probes". Nature 324:163.

Saiki, R., et al. (1989) "Genetic Analysis of Amplified DNA with Immobilized Sequence-Specific Oligonucleotide Probes". Proceedings of the National Academy of Science 86:6230.

4.2.2 Population Studies (TWGDAM 4.1.5.3)

Comey, C.T., and Budowle, B. (1991) "Validation Studies on the Analysis of the HLA DQ α Locus Using the Polymerase Chain Reaction". Journal of Forensic Sciences 36 (6):1633.

Crouse, C.A., et al. (1994) "Analysis of HLA DQ α Allele and Genotype Frequencies in Populations from Florida". Journal of Forensic Sciences 39 (3):731.

Helmuth, R., et al. (1990) "HLA DQ α Allele and Genotype Frequencies in Various Human Populations, Determined by Using Enzymatic Amplification and Oligonucleotide Probes". American Journal of Human Genetics 47:515.

Perkin Elmer (1992) "AmpliType User Guide". Version 2.

Sullivan, K.M., et al. (1992) "Characterisation of HLA DQ α for Forensic Purposes. Allele and Genotype Frequencies in British Caucasian, Afro-Caribbean and Asian Populations". International Journal of Legal Medicine 105:17.

Tamaki, K., et al. (1991) "Frequency of HLA-DQA1 Alleles in the Japanese Population". Hum. Hered 41:209.

4.2.3 Reproducibility (TWGDAM 4.1.5.4)

See above literature on general validation

4.2.4 Mixed Specimen Studies (TWGDAM 4.1.5.5)

Comey, C.T. and Budowle, B. (1991) "Validation Studies on the Analysis of the HLA DQ α Locus Using the Polymerase Chain Reaction". Journal of Forensic Sciences 36 (6):1633.

4.2.5 Environmental Studies (TWGDAM 4.1.5.6)

Rudin, N., et al. (1992) "A Systematic Study of the Effect of Various Environmental Abuses on RFLP and PCR Analysis of Forensic Samples". (Abstract), Journal of Forensic Science Society 32 (3):274, and in Proceedings from The Third International Symposium on Human Identification. Promega Corp. 1992:421.

Comey, C.T., and Budowle, B. (1991) "Validation Studies on the Analysis of the HLA DQ α Locus Using the Polymerase Chain Reaction". Journal of Forensic Sciences 36 (6):1633.

Sims, G., et al. (1992) "The Recovery, Amplification and DQ α typing of DNA from Partially Cremated Human Bones". (Abstract), Journal of Forensic Science Society 32 (2):263.

4.2.6 Matrix Studies (TWGDAM 4.1.5.7)

Rudin, N., et al. (1992) "A Systematic Study of the Effect of Various Environmental Abuses on RFLP and PCR Analysis of Forensic Samples". (Abstract), Journal of Forensic Sciences Society 32 (3):274, and in Proceedings from The Third International Symposium on Human Identification. Promega Corp. 1992:421.

Comey, C.T., and Budowle, B. (1991) "Validation Studies on the Analysis of the HLA DQ α Locus Using the Polymerase Chain Reaction". Journal of Forensic Sciences 36 (6):1633.

Hochmeister, M.N., et al. (1991) "PCR-Based Typing of DNA Extracted from Cigarette Butts". International Journal of Legal Medicine 104:229.

Jung, M.J., et al. (1991) "Extraction Strategy for Obtaining DNA from Bloodstains for PCR Amplification and Typing of the HLA DQ α Gene". International Journal of Legal Medicine 104:145.

Perkin Elmer (1992) "AmpliType User Guide", Version 2.

4.2.7 Non-probative Evidence (TWGDAM 4.1.5.8)

Keel, A., Sims, G., and Buoncristiani, M. (1992) "A Collaborative Study of DQ α Typing by PCR on Sexual Assault Evidence". (Abstract), Journal of Forensic Science Society 32 (3):270.

Comey, C.T., and Budowle, B. (1991) "Validation Studies on the Analysis of the HLA DQ α Locus Using the Polymerase Chain Reaction". Journal of Forensic Sciences 36 (6):1633.

Blake, E., et al. (1992) "Polymerase Chain Reaction (PCR) Amplification and Human Leukocyte Antigen (HLA-DQ α) Oligonucleotide Typing on Biological Evidence Samples: Casework Experience". Journal of Forensic Sciences 37 (3):700.

4.2.8 Non-human Studies (TWGDAM 4.1.5.9)

Gyllenstein, U. and Erlich, H. (1989) "Ancient Roots for Polymorphism at the HLA DQ α Locus in Primates". Proceedings of the National Academy of Sciences 86:9986.

Comey, C.T., and Budowle, B. (1991) "Validation Studies on the Analysis of the HLA DQ α Locus Using the Polymerase Chain Reaction". Journal of Forensic Sciences 36 (6):1633.

Blake, E., et al. (1992) "Polymerase Chain Reaction (PCR) Amplification and Human Leukocyte Antigen (HLA-DQ α) Oligonucleotide Typing on Biological Evidence Samples: Casework Experience". Journal of Forensic Sciences 37 (3):700.

4.2.9 Minimum Sample (TWGDAM 4.1.5.10)

Comey, C.T., and Budowle, B. (1991) "Validation Studies on the Analysis of the HLA DQ α Locus Using the Polymerase Chain Reaction". Journal of Forensic Sciences 36 (6):1633.

Saiki, R., et al. (1986) "Analysis of Enzymatically Amplified B-globin and HLA-DQ α DNA with allele-specific Oligonucleotide probes". Nature 324:163.

4.3 PCR PM + DQA1 VALIDATION (TWGDAM 4.1.1, 4.1.2, AND 4.1.4)

The following section lists validation studies applicable to the PCR-based analysis of loci included in the AmpliType PM + DQA1 kit manufactured by Roche Molecular Systems. References are provided to external studies conducted by other laboratories, as well as to the validation experiments done internally by the DNA Laboratory.

4.3.1 General Developmental Validation (TWGDAM 4.1.5)

Budowle, B., et al. (1995) "Validation and Population Studies of the Loci LDLR, GYPA, HBG, D7S8 and Gc (PM loci), and HLA-DQ α using a Multiplex Amplification and Typing Procedure". *Journal of Forensic Sciences*, 40 (1):45.

Crouse, C.A., Nippes, D.C., and Ritzline, E.L. (1996) "Confirmation of PM Typing Protocols for Consistent and Reliable Results". *Journal of Forensic Sciences* 41 (3):493.

Roy, R. and Reynolds, R. (1995) "AmpliType PM and HLA DQ α Typing from Pap Smear, Semen Smear, and Postcoital Slides". *Journal of Forensic Sciences* 40 (2):266.

Hochmeister, M.N., et al. (1995) "A Method for the Purification and Recovery of Genomic DNA from an HLA DQA1 Amplification Product and Its Subsequent Amplification and Typing with the AmpliType PM PCR Amplification and Typing Kit". *Journal of Forensic Sciences* 40 (4):649.

Woo, K.M. and Budowle, B. (1995) "Korean Population Data on the PCR-Based Loci LDLR, GYPA, HBG, D7S8, Gc, HLA-DQA1, and DIS80". *Journal of Forensic Sciences* 40 (4):645.

Bing, D.H. and Word, C.J. (1995) "PCR Based Forensic Testing with AmpliType PM PCR Amplification and Typing Kit: The Results of Validation Studies from Forensic Laboratories". Presented at the 1995 meeting of the American Academy of Forensic Sciences.

Fildes, N. and Reynolds, R. (1995) "Consistency and Reproducibility of AmpliType PM Results Between Seven Laboratories: Field Trial Results". *Journal of Forensic Sciences* 40 (2):279.

Bille, T.W., et al. (1995) "Validation Studies on the PM System and Population Database for Indiana". *Crime Laboratory Digest* 22:117.

4.3.2 Population Studies (TWGDAM 4.1.5.3)

Perkin Elmer AmpliType PM + DQA1 Product Insert

Budowle B., et al. (1995) "Validation and Population Studies of the Loci LDLR, GYPA, HBG, D7S8 and Gc (PM loci), and HLA-DQ α using a Multiplex Amplification and Typing Procedure". *Journal of Forensic Sciences* 40 (1):45.

Woo, K.M. and Budowle, B. (1995) "Korean Population Data on the PCR-Based Loci LDLR, GYPA, HBG, D7S8, Gc, HLA-DQA1, and DIS80". *Journal of Forensic Sciences* 40 (4):645.

Hausmann, R., Hantschel, M., and Lotterle, J. (1995) "Frequencies of the 5 PCR-Based Genetic Markers LDLR, GYPA, HBG, D7S8, and GC in a North Bavarian Population". *International Journal of Legal Medicine* 107:227.

Hayes, J.M., Budowle, B., and Freund, M. (1995) "Arab Population Data on the PCR-Based Loci: HLA-DQA1, LDLR, GYPA, HBG, D7S8, Gc, and DIS80". *Journal of Forensic Sciences* 40 (5):888.

Bille, T.W., et al. (1995) "Validation Studies on the PM System and Population Database for Indiana". Crime Laboratory Digest 22:117.

Huang, N.E. and Budowle, B., (1995) "Chinese Population Data on the PCR-Based Loci HLA DQ α , Low Density Lipoprotein Receptor, Glycophorin A, Hemoglobin γ G, D7S8, and Group-Specific Component". Human Heredity 45:34.

Scholl, S., et al. (1996) "Navajo, Pueblo, and Sioux Population Data on the Loci HLA-DQA1, LDLR, GYPA, HBG, D7S8, Gc, and DIS80". Journal of Forensic Sciences 41(2):47.

Rodríguez-Calvo, M.S., et al. (1996) "Population Data on the Loci LDLR, GYPA, HBG, D7S8, and GC in Three Southwest European Populations". Journal of Forensic Sciences 41 (2):291.

4.3.3 Reproducibility (TWGDAM 4.1.5.4)

Idaho DNA Laboratory study on liquid and dried specimens from 7 individuals 3/7/97.

Bing, D.H. and Word, C.J. (1995) "PCR Based Forensic Testing with AmpliType PM PCR Amplification and Typing Kit: The Results of Validation Studies from Forensic Laboratories". Presented at the 1995 meeting of the American Academy of Forensic Sciences.

4.3.4 Mixed Specimen Studies (TWGDAM 4.1.5.5)

See above general validation citations.

Idaho DNA Laboratory study on mixed blood specimens 6/19/97.

Idaho DNA Laboratory study on mixed blood and semen specimens 98-T2.

Idaho DNA Laboratory study on mixed saliva and semen specimens 97-T2.

4.3.5 Environmental Studies (TWGDAM 4.1.5.6)

Budowle B., et al. (1995) "Validation and Population Studies of the Loci LDLR, GYPA, HBG, D7S8 and Gc (PM loci), and HLA-DQ α Using a Multiplex Amplification and Typing Procedure". Journal of Forensic Sciences 40 (1):45.

Bing, D.H. and Word, C.J. (1995) "PCR Based Forensic Testing with AmpliType PM PCR Amplification and Typing Kit: The Results of Validation Studies from Forensic Laboratories." Presented at the 1995 meeting of the American Academy of Forensic Sciences.

Bille, T.W., et al. (1995) "Validation Studies on the PM System and Population Database for Indiana". Crime Laboratory Digest 22:117.

4.3.6 Matrix Studies (TWGDAM 4.1.5.7)

Budowle B., et al. (1995) "Validation and Population Studies of the Loci LDLR, GYPA, HBG, D7S8 AND Gc (PM loci), and HLA-DQ α Using a Multiplex Amplification and Typing Procedure". Journal of Forensic Sciences 40 (1):45.

Bing, D.H. and Word, C.J. (1995) "PCR Based Forensic Testing with AmpliType PM PCR Amplification and Typing Kit: The Results of Validation Studies from Forensic Laboratories". Presented at the 1995 meeting of the American Academy of Forensic Sciences.

Bille, T.W., et al. (1995) "Validation Studies on the PM System and Population Database for Indiana". Crime Laboratory Digest 22:117.

4.3.7 Non-probative Evidence (TWGDAM 4.1.5.8)

Idaho DNA Laboratory non-probative sexual assault cases 97-T3, 98-T3.

Fildes, N. and Reynolds, R. (1995) "Consistency and Reproducibility of AmpliType PM Results Between Seven Laboratories: Field Trial Results". Journal of Forensic Sciences 40 (2):279.

Roy R. and Reynolds R. (1995) "AmpliType PM and HLA DQ α Typing from Pap Smear, Semen Smear, and Postcoital Slides". Journal of Forensic Sciences 40 (2):266.

Herrin, G., Fildes, N., and Reynolds, R. (1994) "Evaluation of the AmpliType PM DNA Test System on Forensic Case Samples". Journal of Forensic Sciences 39:1247.

Bing, D.H. and Word, C.J. (1995) "PCR Based Forensic Testing with AmpliType PM PCR Amplification and Typing Kit: The Results of Validation Studies from Forensic Laboratories". Presented at the 1995 meeting of the American Academy of Forensic Sciences.

4.3.8 Nonhuman studies (TWGDAM 4.1.5.9)

Budowle B., et al. (1995) "Validation and Population Studies of the Loci LDLR, GYPA, HBGG, D7S8 and Gc (PM loci), and HLA-DQ α Using a Multiplex Amplification and Typing Procedure". Journal of Forensic Sciences 40 (1):45.

Bing, D.H. and Word, C.J. (1995) "PCR Based Forensic Testing with AmpliType PM PCR Amplification and Typing Kit: The Results of Validation Studies from Forensic Laboratories". Presented at the 1995 meeting of the American Academy of Forensic Sciences.

Bille, T.W. et al. (1995) "Validation Studies on the PM System and Population Database for Indiana". Crime Laboratory Digest 22:117.

4.3.9 Minimum sample (TWGDAM 4.1.5.10)

Idaho DNA Laboratory sensitivity study (5/7/97)

Budowle B. et al. (1995) "Validation and Population Studies of the Loci LDLR, GYPA, HBGG, D7S8 and Gc (PM loci), and HLA-DQ α Using a Multiplex Amplification and Typing Procedure". Journal of Forensic Sciences 40 (1):45.

Bing, D.H. and Word, C.J. (1995) "PCR Based Forensic Testing with AmpliType PM PCR Amplification and Typing Kit: The Results of Validation Studies from Forensic Laboratories". Presented at the 1995 meeting of the American Academy of Forensic Sciences.

Bille, T.W., et al. (1995) "Validation Studies on the PM System and Population Database for Indiana". Crime Laboratory Digest 22:117.

4.3.10 On-site evaluation (TWGDAM 4.1.5.11)

See above Idaho DNA Laboratory experiments

4.3.11 Meetings/papers (TWGDAM 4.1.5.12)

See above listed publications

4.3.12 Inheritance (TWGDAM 4.2.1.)

Roche Molecular Systems unpublished information

4.3.13 Gene mapping (TWGDAM 4.2.2)

Perkin Elmer (1992) "AmpliType User Guide", Version 2.

Yamamoto, T., et al. (1984) "The Human LDL receptor: A Cysteine-Rich Protein with Multiple Alu Sequences in Its mRNA," *Cell*, 39:27.

Siebert, P.D. and Fukuda, M. (1987) "Molecular cloning of a human glycoporphin B cDNA: Nucleotide sequence and genomic relationship to glycoporphin A," *Proc. Natl. Acad. Sci. USA*, 84:6735.

Slightom, J.L., Blechl, A.E., and Smithies, O. (1980) "Human Fetal Agamma and Ggamma Globin Genes: Complete Nucleotide Sequences Suggest that DNA can be Exchanged Between These Duplicated Genes," *Cell*, 21:627.

Horn, G.T., et al. (1990) "Characterization and Rapid Diagnostic Analysis of DNA Polymorphisms Closely Linked to the Cystic Fibrosis Locus," *Clinical Chemistry*, 36:1614.

Yang, F., et al. (1985) "Human group-specific component (Gc) is a member of the albumin family," *Proc. Natl. Acad. Sci. USA*, 82:7994.

4.3.14 Detection (TWGDAM 4.2.3)

Perkin Elmer (1992) "AmpliType User Guide", Version 2.

Saiki et al. (1986) "Analysis of Enzymatically Amplified B-globin and HLA-DQ DNA with allele-specific Oligonucleotide probes," *Nature* 324:13.

Saiki et al. (1989) "Genetic Analysis of amplified DNA with immobilized sequence-specific oligonucleotide probes," *Proceedings of the National Academy of Science* 86:6230.

Roche Molecular Systems unpublished information

4.3.15 Polymorphism (TWGDAM 4.2.4)

Perkin Elmer (1992) "AmpliType User Guide", Version 2

Yamamoto, T., et al. (1984) "The Human LDL receptor: A Cysteine-Rich Protein with Multiple Alu Sequences in Its mRNA," *Cell*, 39:27.

Siebert, P.D. and Fukuda, M. (1987) "Molecular cloning of a human glycoporphin B cDNA: Nucleotide sequence and genomic relationship to glycoporphin A," *Proc. Natl. Acad. Sci. USA*, 84:6735.

Slightom, J.L., Blechl, A.E., and Smithies, O. (1980) "Human Fetal Agamma- and Ggamma-Globin Genes: complete Nucleotide Sequences Suggest that DNA Can Be Exchanged between These Duplicated Genes," *Cell*, 21:627.

Horn, G.T., et al. (1990) "Characterization and Rapid Diagnostic Analysis of DNA Polymorphisms Closely Linked to the Cystic Fibrosis Locus," *Clinical Chemistry*, 36:1614.

Yang, F., et al. (1985) "Human group-specific component (Gc) is a member of the albumin family," *Proc. Natl. Acad. Sci. USA*, 82:7994.

4.3.16 Amplification (TWGDAM 4.4.1)

Roche Molecular Systems unpublished information

4.3.17 Protection from contamination (TWGDAM 4.4.12)

See *Section 5.5*

Perkin Elmer (1992) "AmpliType User Guide", Version 2

Higuchi and Kwok (1992) "Avoiding false positives with PCR," *Nature* 339:237.

4.3.18 Conditions, Cycles (TWGDAM 4.4.1.3 and 4.4.1.4)

Critical reagent concentrations and thermocycling parameters follow the recommendations published in the AmpliType® PM + DQA1 kit insert. Primer sequences used in the PM+DQA1 kit are unpublished, but are on file at the Idaho DNA lab.

4.3.19 Differential Amplification (TWGDAM 4.4.1.5)

Budowle B. et al. (1995) "Validation and Population Studies of the Loci LDLR, GYPA, HBGG, D7S8 and GC (PM loci), and HLA-DQ α Using a Multiplex Amplification and Typing Procedure," *Journal of Forensic Sciences*, 40:45.

Bing, D.H. and Word, C.J. (1995) "PCR Based Forensic Testing with AmpliType PM PCR Amplification and Typing Kit: The Results of Validation Studies from Forensic Laboratories," Presented at the 1995 meeting of the American Academy of Forensic Sciences.

Bille, T.W. et al. (1995) "Validation Studies on the PM System and Population Data Base for Indiana," *Crime Laboratory Digest*, 22:117.

Walsh et al., (1992) "Preferential PCR amplification of Alleles: Mechanisms and Solutions" in *PCR Methods and Applications*. CSHL Press. PCR. 241.

Roche Molecular Systems unpublished information

4.3.20 Multiplex amplification (TWGDAM 4.4.1.6)

Fildes, N. and Reynolds, R. (1995) "Consistency and reproducibility of AmpliType PM Results between Seven laboratories: Field Trial Results," *Journal of Forensic Sciences*, 40:279.

4.3.21 **Characterization with hybridization** (TWGDAM 4.4.2.2 [a] and [c])

Budowle B. et al. (1995) "Validation and Population Studies of the Loci LDLR, GYPA, HBGG, D7S8 and GC (PM loci), and HLA-DQ α Using a Multiplex Amplification and Typing Procedure," *Journal of Forensic Sciences*, 40:45.

Roche Molecular Systems unpublished information

"C" or "S" dot is present as a threshold control on every strip

Product Gel performed to assess amplification efficacy

4.3.22 **Internal Validation** (TWGDAM 4.5)

Idaho DNA Laboratory validation studies

4.3.23 **Known Samples** (TWGDAM 4.5.1)

Concordance studies with Utah Department of Public Safety DNA Laboratory

Proficiency tests performed by Idaho DNA Laboratory Personnel

Positive control DNA samples extracted and typed with each analysis

4.3.24 **Contamination Control** (TWGDAM 4.5.4)

Ongoing proficiency testing

Proper use of controls and standards as described in DNA Laboratory Quality Assurance Manual

4.3.25 **Proficiency Testing** (TWGDAM 4.5.5)

See *Section XII*

V. **QC PROCEDURES – EQUIPMENT, REAGENTS, MATERIALS, AND PHYSICAL FACILITIES**

The purpose of the procedures in this section is to ensure that the parameters critical to the testing process are routinely monitored in the manner necessary to maintain the success and reliability of the testing procedures.

It is possible to verify "after the fact" that the equipment, materials, and reagents used in an analysis have not significantly affected the reliability of the results. For example, the QC allelic and extraction controls described in *Section VII* are designed to signal potential problems in the analysis. If acceptable results are obtained for these controls, it is reasonable to assume that the results from other samples analyzed with them are also reliable. If the controls indicate a problem with the analysis, it may be possible to determine the source of the problem. Depending on the nature of the problem, re-analysis of the samples may be required.

It is highly desirable to minimize the need for repeat analysis due to failure of equipment, materials or reagents. To that end, QC procedures focus as much as possible on preventing problems before they occur rather than on dealing with them after they happen.

5.1 QC MANUALS AND NOTEBOOKS

Protocols for preparing reagents, reagent logs, QC procedures for commercial products and kits, equipment calibration data and logs, run sheets, worksheets, forms, and other logs are maintained in the DNA QC DATA NOTEBOOK or the DNA REAGENT PREPARATION NOTEBOOK.

5.2 EQUIPMENT CALIBRATION CHECKS AND MAINTENANCE SCHEDULE (TWGDAM 5.1.3)

5.2.1 Refrigerators/freezers

The temperatures of the freezers and refrigerators used for reagent and evidence storage are checked and recorded as part of the weekly maintenance routine. Those in which reagents are stored have thermometers traceable to NIST. The refrigerators and freezers checked are:

- Refrigerator and chest freezer in the evidence vault
- Refrigerator and freezer in the serology DNA Laboratory
- Refrigerator/ freezer in the DNA product room

5.2.2 Water baths and heat blocks

The temperature of devices used to maintain critical parameters (for example, denaturation and hybridization conditions) are checked at the time of use by the Analyst. Maintenance will be performed as needed.

5.2.2.1 Heat blocks

The heat block used for digestion is verified at every use by a NIST-traceable thermometer.

5.2.2.2 Shaking water bath (PCR room)

The calibration of the PCR room shaking water bath is verified at every use by a NIST-traceable thermometer in the bath; this is documented on the typing sheet in the Analyst's notes.

5.2.3 Thermal cycler

The thermal cycler is checked in accordance with the manufacturer's specifications. The diagnostic file checks and temperature calibration and uniformity checks are performed monthly. The probe used to check the temperature is calibrated annually by the manufacturer. Maintenance will be performed as needed.

5.2.4 pH meter

5.2.4.1 At the time of use:

Clean electrode, check fluid level (fill if necessary), standardize with buffers.

5.2.4.2 Maintenance:

Weekly: check fluid level (fill if necessary), check storage

Monthly: standardize at pH 7.0 and record values obtained with pH 4.0 and pH 10.0 solutions.

5.2.5 Balances

Maintenance will be performed as needed.

5.2.5.1 Calibration:

Monthly: check for deviation using internal weights dialed on and off beam; check accuracy using an external class II weight.

Annually: balances will be checked by an outside contractor.

5.2.6

Water purification system

Maintenance will be performed as needed per the manufacturer's instructions

5.2.7

Casework Pipettors

5.2.7.1

Calibration check

Pipettors should be checked at least annually.

Maintenance will be performed as needed.

5.3

MAINTENANCE LOGS

Temperature logs for refrigerators and freezers and calibration check records for the thermal cycler are stored in the DNA QC DATA NOTEBOOK. Instrument operation manuals are kept close to the equipment to which they apply. (TWGDAM 5.1.2)

5.4

MATERIALS AND REAGENT PREPARATION, LABELING

5.4.1

Logs of commercial supplies and kits (TWGDAM 5.2.1)

PCR kits and Quantiblot kits are dated upon receipt. The date, manufacturer's kit lot number and the lot number's and expiration dates of the kit's contents are recorded in DNA QC DATA NOTEBOOK.

5.4.2

Formulation, inventory and labeling of reagents (TWGDAM 5.2.2, 5.2.3)

The following records are kept in the DNA REAGENT PREPARATION NOTEBOOK:

- A list of all solutions, their formulas, and method of preparation.
- Documentation of preparation of each lot of reagents, including lot numbers of chemicals used.

Each container will bear a label with its identity, the date the reagent was prepared, and the identity of the preparer (or a lot number which can be used to track that information).

5.4.3

Glassware and plastic supplies preparation (TWGDAM 5.2.6)

Glassware and plastic supplies will either be purchased sterile or be sterilized in-house by autoclaving.

5.4.4

Quality control of materials and reagents

To assess the quality of each lot of reagents and buffers, samples with known DNA profiles are processed through the current casework protocol.

Routine QC of PCR reagents and kits involves the amplification and detection of at least three extracted DNA samples of known type and one negative amplification control. Testing is performed with supplies and reagents from one kit out of each lot of kits received. After confirmation that the known samples produce the expected types with the expected relative dot intensities, and the negative amplification control produces no type, the kits and reagents are placed into use.

5.5

DEDICATED PCR FACILITIES (TWGDAM 5.3)

[See also *Section 6.3* and *6.4* in the Manufacturer's User Guide and the product information sheet regarding evidence handling and detailed DNA Laboratory set-up guidelines.]

Because of the sensitivity of PCR-based procedures, special precautions must be taken to avoid contamination of samples with foreign human-DNA-containing material. Contamination is defined as

the negligent introduction of foreign material to the sample by responsible personnel after they assume control of the evidence. Precautions must be taken to avoid contamination either during collection and transportation of the evidence or in the DNA Laboratory.

Evidence collection technicians must be trained in proper collection and preservation techniques. DNA Laboratory personnel may assist by offering education and training opportunities. However, DNA Analysts cannot be responsible for any contamination introduced before the evidence is received into custody of the DNA Laboratory.

5.5.1 Sources of contamination

- 5.5.1.1 Sample contamination with human genomic DNA from the environment or analyst**
To minimize the potential of chance contamination, Analysts should wear gloves while handling evidence and performing analyses. Gloves should be changed when necessary. Sample tubes should be closed when not immediately in use. Dispersal of aerosols should be minimized. Careful pipetting technique, spinning of tubes prior to opening and controlled release of tube caps will contribute to this effort.

Reagents used for DNA extraction and PCR amplification will be sterilized by either autoclaving or filtration as directed in the DNA Reagent Preparation manual. Supplies such as tubes, pipettes, tips and toothpicks shall be purchased as sterile or sterilized by autoclave as appropriate for their composition.

- 5.5.1.2 Cross contamination between samples during sample extraction/preparation (TWGDAM 5.3.1, 5.3.2, 5.3.3)**
Precautions should be taken during DNA extraction and PCR set-up to prevent transfer of DNA from one sample to another. Analysts must use a fresh pipette tip for each sample, open tubes carefully, and keep sample tubes closed when not in use. In particular two evidence samples, or an evidence and reference sample, should not be open at the same time in the same vicinity. Evidence and reference samples will be processed separately either in time or space.

DNA extraction and PCR reaction set-up will be performed in separate work areas. Typically DNA extractions are conducted on the Analyst's bench and in the chemical fume hoods. PCR set up is carried out in the biological hood. A dedicated pipettor and rack is used for setting up PCR reactions and for aliquoting PCR pre-mix. Barrier tips should be used for PCR amplification set-up. PCR amplification is performed in a separate PCR room with its own dedicated equipment and supplies.

- 5.5.1.3 Contamination of a sample with amplified DNA from a previous PCR reaction (TWGDAM 5.3.4)**
PCR product carry-over is the most troublesome potential source of contamination. Carry-over is defined as the introduction of amplified DNA into a sample that has not yet been amplified. It is crucial to contain amplified PCR product to prevent it from encroaching on samples before they are amplified. To this end, DNA extraction, PCR reaction set-up, and PCR amplification are performed in a physically different rooms. All activities involving the handling of amplified DNA are confined to the dedicated PCR room. This includes amplification, typing (hybridization and color development), gel electrophoresis of amplified DNA waste disposal of amplified DNA solutions, photography, and storage of amplified DNA.

- 5.5.1.4 Procedures to be followed in the PCR room (TWGDAM 5.1.4)**
Dedicated and distinguishable lab coats will be worn in the PCR product room and in no other area of the DNA Laboratory. Dedicated equipment (e.g. scissors, pipettors, and timers) are marked with red tape and are not normally to be taken from the room (except as described below). (TWGDAM 5.2.5.) Gloves will be removed and disposed of before leaving the PCR

room. Waste from the PCR room will be delivered immediately to the waste storage areas without being stored in other DNA Laboratory rooms. Neither equipment, supplies or data are to leave the PCR room unless they have been thoroughly cleaned with ethanol or bleach or subjected to cross-linking by UV irradiation. (TWGDAM 5.3.5) Any PCR product remaining after typing will be stored in the PCR room refrigerator or freezer.

VI. PROCEDURES FOR EVIDENCE HANDLING AND SECURITY (TWGDAM 6)

6.1 CHAIN OF CUSTODY AND SECURITY

The policies that cover procedures to ensure that the chain of custody of the evidence is maintained and that evidence is protected from loss, deterioration or deleterious change at the DNA laboratory are contained in the BUREAU POLICY MANUAL. Laboratory security is addressed there as well.

6.2 STORAGE

6.2.1 Sample containing, or potentially containing, biological evidence

In general, all biological evidence should be thoroughly dried and stored frozen as soon as possible upon receipt. Evidence which might be compromised by condensation upon thawing may be stored at room temperature. Liquid samples should be refrigerated on arrival and processed as soon as possible thereafter.

6.2.2 Evidence from which all usable biological material has been removed

Evidence from which all usable biological material has been removed may be stored at room temperature.

6.3 EXAMINATION

6.3.1 Responsibility and coordination

Once assigned to a case, the Analyst will sign for all evidence he/she is taking and will coordinate with the referring BFS Regional Laboratory or agency as needed.

6.3.2 Documentation

Before altering the evidence by sampling, the examiner will record its condition by written description, diagram, and/or photography. The Analyst must remain alert to other categories of evidence that would be destroyed by serological/DNA testing and should request a consultation if potentially probative non-biological evidence is found or suspected.

6.3.3 Safeguards against DNA Laboratory contamination

In general, the examiner should examine only one item of evidence at a time, marking the evidence with a unique identifier and returning it to its container before opening another item. In particular, items of evidence collected from suspects should be examined separately from victim samples; reference samples should be handled separately from evidence samples. This practice will guard against sample mix-up and cross transfer.

Exam tables and benches are to be covered with clean paper before the examination begins. Clean paper should be laid down under each item examined. This process will prevent cross transfer from one item to another.

Examiners will confine examination of their evidence as much as possible to their own assigned work areas or evidence exam room. Evidence should remain on the lab bench or exam table only long enough to be examined and have appropriate samples removed for analysis. The evidence will be closed and sealed as soon as analytical samples have been removed.

Once the examination of a case is completed, the lab bench/exam table is to be wiped with disinfectant solution (TWGDAM 5.3.5).

6.4 ANALYSIS

6.4.1 Sample retention guidelines (TWGDAM 7.1.3)

Where possible, the Analyst will ensure that sufficient sample is left for possible re-analysis. If the entire sample must be consumed in order to obtain an interpretable result, the examiner should confer with the submitter to ensure that the legal implications of destruction of the sample have been properly considered. This conversation should be documented in the case file.

6.4.2 Safeguards against DNA Laboratory contamination

Care should be taken to prevent inadvertent transfer between evidence samples, in particular those that might contain disparate DNA concentrations, or between evidence and reference samples, or between victim and suspect reference samples. This may include processing the samples at a different time or space and arranging the order of samples during extraction and analysis such that the impact of any carryover would be minimized. (TWGDAM 7.2)

6.5 DISPOSITION OF EVIDENCE

The examiner is responsible for re-packaging and resealing the evidence prior to return and for making appropriate entries in the chain of custody records of transfer to the evidence custodian. The latter is responsible for maintaining appropriate records of UPS shipment or other record of return to the submitting agency.

The original item of evidence and original packaging containing identifying markings will be returned in sealed condition to the appropriate agency together with the report upon completion of the case. Evidence containers should indicate the appropriate conditions under which the returned evidence should be stored. The evidence may be returned in person or via UPS or other appropriately secure means which provides a return receipt.

6.6 DNA CASE ACCEPTANCE AND PRIORITY CRITERIA

The goal of the DNA Laboratory's case acceptance policy is to ensure that its limited resources are focused and utilized as effectively as possible. In general, the DNA Laboratory will rely on the BFS regional laboratories to evaluate the potential probative value of DNA analysis to their cases and to coordinate the application of DNA testing in the investigation of serial crimes.

6.6.1 Case acceptance guidelines (TWGDAM 7.1.1)

6.6.1.1 The case must be referred to the DNA Laboratory through a BFS Regional Laboratory. The referring Laboratory will make the initial contact with the DNA Laboratory and participate in assessing the priority to be given the case.

6.6.1.2 In general, input from the BFS DNA Laboratory, BFS Regional Laboratory, investigating agency, and prosecuting agency will be sought in making the determination that the case is suitable for DNA analysis.

6.6.1.3 In general, cases will be accepted for DNA analysis only if they concern violent crimes against persons, or felony crimes of a sexual nature. The DNA Laboratory will not accept case evidence which has been previously analyzed by another DNA Laboratory except with express approval of the Bureau Chief.

6.6.1.4 Other cases falling outside these guidelines may be accepted with the express approval of the Bureau Chief.

6.6.1.5 The Bureau Chief's approval will be required to accept cases with a date of offense prior to April 1, 1998.

6.6.2 Case priority

The priority given a particular case will be established by the DNA Laboratory based on:

- a) The investigative or probative value (jointly assessed by the BFS Regional Laboratory, DNA Laboratory, and investigator and/or prosecutor) that DNA is expected to have in the case,
- b) The likelihood (assessed by the DNA Laboratory) that the analysis could provide interpretable and useful results within a reasonable time frame,
- c) Court date or other legal deadline, and the order in which the case is received.

VII. ANALYTICAL PROCEDURES (TWGDAM 7)

7.1 ANALYTICAL PROTOCOLS

Protocols for each analytical or assessment method routinely used in casework are maintained in the DNA PROCEDURES MANUAL. These protocols include sample preparation, sample extraction, and specify the required standards and controls for DNA analysis and data interpretation. Each protocol will exhibit the date(s) on which the protocol was approved by the Manager.

7.2 ASSESSMENT OF QUALITY AND QUANTITY

7.2.1 A microscopic examination of suspected sexual assault samples should be conducted. (TWGDAM 7.1.1)

A microscopic examination of suspected sexual assault samples provides an indication of both the number and relative concentrations of sperm and non-sperm cells. Microscopic observations should be included in the case notes.

7.2.2 The quantity of human DNA in the sample should be estimated. (TWGDAM 7.3.4)

A slot blot using a primate-specific probe should be performed routinely to determine the amount of sample to amplify by PCR. Extenuating circumstances, such as extremely limited sample, might preclude a slot blot evaluation. Dilutions of human DNA are used as semi-quantitative standards on slot blots.

7.2.3 The quality of the DNA sample may be determined. (TWGDAM 7.3.1)

On occasion, it may be desirable to obtain an indication of the average size of the DNA fragments in a sample before analysis. In particular, this assessment may factor into the decision as to which typing system should be chosen. A yield gel will assist in this determination. A size marker is run on the yield gel to help assess the amount of degradation of the extracted samples and to control for the electrophoresis. Dilutions of human DNA are used as semi-quantitative standards on yield gels. Yield gel results may support interpretations which include sample degradation as an explanation for missing dots.

7.3 STANDARDS AND CONTROLS

A variety of standards and controls are required to assess the quantity and quality of the sample and the effectiveness, accuracy, and precision of the analytical procedures in a particular case. Evaluation of the controls is essential for the proper interpretation of the test results. (See DNA PROCEDURES MANUAL for additional details on controls required for each type of analysis.)

7.3.1 Standards and Controls for General DNA Analysis

7.3.1.1 Substrate sample (TWGDAM 7.5.1.3)
Where appropriate, an unstained portion of the substrate adjacent to the questioned stain should be sampled and analyzed. This may identify any possible contribution to positive results by human materials present on the background substrate.

7.3.1.2 Extraction control (QC sample):
The purpose of this control is to ensure that the extraction procedure worked properly. An extraction control is included with each extraction set. The control may be either a bloodstain or other biological material from a previously characterized source which is extracted and typed. Poor yield, degradation, an incorrect type, or other problems observed in the extraction control are useful for diagnosing similar problems which may be observed in the evidence samples and may dictate repeating the analysis.

7.3.1.3 Reagent blank (TWGDAM 7.5.1.1):
The reagent blank is a tube containing no sample that is carried through the entire analytical process and is exposed to all the reagents used for the extraction, amplification and typing with each batch of samples. The purpose of this sample is to detect contamination that might occur from the reagents, between the samples being processed, or from previously amplified PCR product. Generally, the reagent blank should be carried through the amplification and typing process for each of the PCR systems in which the evidence is typed unless the blank has been entirely consumed. At a minimum, the reagent blank must be typed in at least one PCR system.

7.3.2 Assessment of Quality and Quantity

A product gel is run after amplification of an SSO dot blot marker to verify that amplified DNA of the expected product size(s) is present. A size marker standard is used to estimate the product size. The absence of larger bands in, e.g. PM, may indicate degradation extending into that size range.

7.3.3 STANDARDS AND CONTROLS FOR PCR DOT BLOT ANALYSIS (TWGDAM 7.5.1)

7.3.3.1 Positive amplification control (TWGDAM 7.5.1.2)
A DNA sample of known type included to verify that the amplification and typing processes are working properly and to detect conditions that might lead to differential amplification. The positive amplification control should include an allele that is sensitive to amplification conditions in that particular system.

7.3.3.2 Negative amplification control (TWGDAM 7.5.1.1)
The negative amplification control contains only the reagents used to prepare the PCR amplification mixture, including sample buffer, for each batch of samples. The purpose of this sample is to detect contamination which might occur during amplification set-up from the reagents, between the samples being processed, or from previously amplified PCR product. This tube may be left open during the preparation of the PCR cocktails as an added control to monitor the environment of the set up area.

7.3.3.3 Threshold control
The "C" dot on DQAI typing strips and the "S" dot on the PM strips are threshold controls designed to ensure that an adequate amount of amplified product is present for reliable genotyping.

7.3.3.4 Duplicate analysis
Monitoring of critical parameters, including negative controls and other indicators of potential contamination, in the DNA Laboratory has not demonstrated a need for duplicate analysis of all samples. Where duplicate analysis is deemed appropriate, samples are split as early as

possible in the process. Replicate analysis (e.g. reference sample from victim vs. epithelial cell fraction of vaginal swab; analysis of two separate spots within a coherent bloodstain pattern) is considered to serve the same purpose as duplicate analysis. If the results between two positive typings differ, the results for that sample are inconclusive. An additional analysis must be conducted to assess the cause of the difference before any conclusion may be drawn.

7.4 INTERPRETATION GUIDELINES

The interpretation of results in casework is necessarily a matter of professional judgment and expertise. Not every situation can or should be covered by a pre-set rule; nor is it expected that competent Analysts will always be in full agreement in a particular case. However, it is important that the DNA Laboratory develop and adhere to minimum criteria for the interpretation of analytical results. These criteria will be based on validation studies, literature references, and casework experience and will be developed with maximum input from Analysts. It is to be expected that interpretation guidelines will continue to evolve as the collective experience of the DNA Laboratory grows with various analytical procedures. Interpretation of the results within the context of the case is the responsibility of the primary Analyst.

The purpose of these guidelines is to establish a general framework and outline nominal standards to ensure that:

- The conclusions in casework reports are supported by the analytical data, including that from appropriate standards and controls.
- The interpretations are objective, consistent from Analyst to Analyst, adhere to defined guidelines, and fall within established tolerance limits.

7.4.1 PCR Dot Blot Evaluation and Interpretation Guidelines (TWGDAM 8.2.1, 8.2.2)

The goal of the evaluation and interpretation of PCR dot blot strips (DQA1, PM) is to deduce the DNA profile(s) of the donor(s) of the questioned samples. The interpretation of the data is to be based on the education, training and experience of a qualified Analyst operating within the following guidelines.

7.4.1.1

Substrate sample

The alleles, if any, present in the substrate sample will assist in the interpretation the DNA results obtained from the evidentiary sample, especially where the results indicate a mixed sample.

7.4.1.2

Extraction control (QC sample) (TWGDAM 8.2.1.1)

The extraction control must produce the correct type and lack evidence of gross contamination. Where the QC sample produces the expected results and questioned samples are clearly typable, faint signals other than the expected type seen in the QC sample may not be serious enough to affect the final interpretation. They should, however, be noted. If the intensity of extraneous alleles is greater than trace or their presence in questioned samples interferes with interpretation, the test is considered inconclusive and the analysis must be repeated for any conclusions to be drawn.

7.4.1.3

Reagent blank (TWGDAM 8.2.1.1)

Where test samples are clearly typable, faint signals seen in the reagent blank, where the threshold control dot is not visible, may not be serious enough to affect the final interpretation. They should, however, be noted. If signals exceeding a visible control dot in intensity occur in any negative control, the test is considered inconclusive and the analysis must be repeated for any conclusions to be drawn.

- 7.4.1.4 Positive amplification control (TWGDAM 8.2.1.1)**
The positive amplification control(s) must give the expected results. If the positive amplification control fails to give the correct result, the analysis must be repeated.
- 7.4.1.5 Negative amplification control (TWGDAM 8.2.1.1)**
No signal should be present in the negative amplification control. However, trace dots in a negative amplification control will not necessarily negate a useful interpretation of the data in that run. Where test samples are clearly typable, faint signals seen in the negative amplification control, where the threshold control dot is not visible, may not be serious enough to affect the final interpretation. They should, however, be noted. If signals exceeding a visible control dot in intensity occur in any negative control, the test is considered inconclusive and the analysis must be repeated for any conclusions to be drawn.
- 7.4.1.6 Threshold control (TWGDAM 8.2.2.1)**
The results of a mixed sample are interpreted with the knowledge that the threshold control was designed for single source samples and does not preclude the meaningful interpretation of dots less than "C" in this circumstance. If the control dot is not visible on the strip, the amplification/hybridization may not have detected one or more alleles present in the sample. Generally, conclusions about complete profiles will not be drawn from strips without a visible control dot.
- 7.4.1.7 Definition of a "dot" (TWGDAM 8.2.2.1)**
A dot is defined as a circular shaped blue area on a membrane strip. Its edge should be distinct all around the circumference.
- 7.4.1.8 Trace Dots (TWGDAM 8.2.2.1)**
Trace dots are defined as dots that are so weak that alterations in lighting conditions, visual medium (strip or photograph), or Analyst judgment might affect their detection. By their very nature, the presence or absence of trace dots would not affect the major or minor types determined for any particular sample. Therefore, a difference of opinion between Analysts regarding the presence or absence of an allele in the trace range will not be interpreted as a disagreement about the sample profile.
- 7.4.1.9 Dot intensity guidelines (TWGDAM 8.2.2.1)**
Dot intensities are recorded in relationship to the threshold control dot (i.e. "C" or "S") present on the same strip. Dots of greater intensity than the control dot will be recorded as "C+" or "S+". Dots of lesser intensity than the control dot, will be recorded as "C-" or "S-". Dots reacting significantly less than the control dot will be recorded as "trace".
- 7.4.2 Population frequency calculations (TWGDAM 8.2.2.3)**
For a given population(s) and/or hypothesis of relatedness, the statistical interpretation shall, as deemed applicable, be made following the recommendations 4.1, 4.2 or 4.3 of the National Research Council report entitled "The Evaluation of Forensic DNA Evidence" (1996) and/or court directed method. These calculations shall be derived from a documented population database appropriate for the calculation.
- 7.4.2.1 Population database**
Frequency estimates will be calculated for at least three major population groups, generally Caucasian, African American and Hispanic, and for additional racial/ethnic groups known to be relevant to the case for which data is available. For PCR dot blot markers DQA1 and LDLR, GYPA, HBGG, D7S8, GC, frequencies will be calculated using expected values from data in "Validation and population studies of the loci LDLR, GYPA, HBGG, D7S8, and GC (PM loci), and HLA-DQ α using a multiplex amplification and typing procedure", Budowle, et al., Journal of Forensic Sciences, Vol. 40, No. 1, Jan 1995.

- 7.4.2.2 Heterozygote frequencies**
Heterozygote frequencies will be calculated using the standard formula $f_{pq} = 2pq$. Frequency calculations for homozygotes in PCR systems will be done using the NRC II formula: $f_{pp} = p^2 + p(1-p)q$ with $q=0.01$. If a small isolated population such as Native Americans is relevant to a particular case, frequency calculations will use 0.03 for q .
- 7.4.2.3 Rare alleles**
A five event minimum allele frequency will be used for rare alleles. For each individual allele, an observed allele count less than five is raised to five. This modified allele count is converted to a frequency and used for all subsequent frequency calculations.
- 7.4.2.4 Significant figures**
Frequency estimates will be reported using two significant figures. Truncation rather than rounding will be used to arrive at the two significant figures.

VIII. CASE NOTES

- 8.1 CASE NOTES AND DOCUMENTATION (TWGDAM 8.1)**
The functions of note taking are to permit internal review of the work product, to allow re-evaluation of the data by independent reviewers, to support the conclusions in the DNA Laboratory report, and to provide a foundation for the introduction of the work product into court.
- 8.1.1 Evidence evaluation**
Notes on the initial examination of evidence items should describe the item, noting method of sealing or packaging, the item's condition, presence of stains, or other unusual features. Evidentiary material should be diagrammed, described, or photographed before sampling. The locations from which stain and control samples were collected for typing should be documented.
- 8.1.2 Written examination notes**
Notes must be recorded at the time the work is done and must accurately reflect the work as actually performed. They must be sufficiently detailed to support the conclusions in the report and to allow duplication of the work at another time.
- For casework, analytical notes should be recorded on DNA case note forms or the appropriate run sheet or checklist. Each page will be numbered, dated, and initialed by the Analyst. Handwritten notes should be made in ink. Corrections should be made by lining through so the original text is still visible and initialing the correction.
- 8.1.3 Photography**
In general, the results of either assessment or analytical procedures will be documented by direct photography, or other appropriate technique. Because the intensity of the dots on PCR dot blot strips fades over time, the strips should not be retained once photographed. If the photograph does not capture what can be seen at the time, that fact should be documented. Conditions of the assessment or analytical procedure will be recorded on the appropriate run sheet (copies of standard DNA analysis run sheets are found in the DNA QC DATA NOTEBOOK) along with notes on the analytical procedures.
- 8.1.4 Case record**
The DNA Laboratory shall maintain, in a case record, all documentation generated by examiners related to case analyses.
- Photographs, run sheets, and other data generated during the analysis will be retained in the case file.

8.1.5 Double reading of analytical results

Double reading helps guard against clerical/typographical errors and also guards against bias in the interpretation of borderline results. Second readers are to make an independent, objective assessment of the data, without prior knowledge of the results obtained by the first reader. They shall date and initial the record to document their findings.

8.1.5.1 Second reader qualifications

Double reading can be performed by the Technical Leader as a part of the technical review of the case or by any qualified peer. A "qualified second reader" is one who has demonstrated, through proficiency tests, that he/she is competent to read and interpret the type of data being reviewed.

8.1.5.2 Time of second reading

Double reading should be done at the time of the first reading if the primary data may be ephemeral (e.g. PCR dot blots). However, as necessary, this type of data may be double read from a photograph if the photograph accurately captures what could be seen initially. Occasionally extremely faint dots visible on the original strips may not be visible in the photograph. This should be documented.

8.1.5.3 Resolution of disagreement

The primary Analyst is also responsible for verifying that consistent results were obtained by the second reader and for identifying and addressing any differences.

8.1.5.4 Double reading

Double reading is required in casework for PCR dot blots. The second reading confirms identification of allelic dots, relative intensity of allelic dots to control dots and to each other, and indicated types.

IX. REPORTS

The function of the DNA Laboratory report is to communicate the analytical results, conclusions, and interpretation of the Analyst, conveying the essence of what the examiner would say if asked for his/her expert opinion in court. The conclusions should clearly state appropriate qualification or limitations on the evidence interpretation.

9.1 REPORTS ACCORDING TO WRITTEN GUIDELINES SHALL INCLUDE:(TWGDAM 8.3)

9.1.1 Case Identifier (TWGDAM 8.3.1)

The unique number assigned to the case when it entered the DNA Laboratory.

9.1.2 Summary and Conclusion (TWGDAM 8.3.7)

This section will contain a concise summary of the principal findings of the examination including a statement of conclusion. All statements made in the summary must be supported by data in the results section of the report or other referenced reports. Any assumptions relevant to the interpretation should be stated. Where statistics are critical to communicate the limitations of the conclusion, they may be included in the summary.

9.1.3 Description of evidence examined

This section will list all evidence received by the DNA Laboratory from the case, whether examined or not. A unique number, description, and reported source (if any) will be listed for each item. A specific statement should be made indicating which (if any) of the items received were not examined by the DNA Laboratory. Chain of custody documentation should be indicated (e.g., from whom and how the items were received). In general, the numbers used by the submitting agency or referring BFS Regional Laboratory will be used, so long as they provide unique identifiers.

- 9.1.4 The examination conducted, including specific loci and a short description of the methodology (TWGDAM 8.3.4, 8.3.5)**
In the report, analytical methods will be described only in general terms (e.g., PCR typing at the DQA1 locus using the AmpliType™ kit). The loci examined will be specified.
- 9.1.5 Results (TWGDAM 8.3.6)**
In this section the results of any tests conducted will be reported. A specific statement should be made about each item examined. It is preferable to summarize analytical data, e.g., DQA1 types, in a table for ready reference. Items will be referenced by both their unique item number and description. Attributes of the sample which affect the interpretation (e.g. degradation of DNA) will be specified.
- 9.1.6 A quantitative interpretative statement (TWGDAM 8.3.8)**
The approach to calculating profile frequencies should be identified. Assumptions about kinship should be explicated.
- 9.1.6.1 Additional information**
When necessary, additional information on complex issues such as statistical calculations or specific analytical procedures may be provided by way of an appendix.
- 9.1.7 A qualitative interpretative statement**
The interpretation of DNA profile patterns will be described, and indistinguishable patterns noted. An accounting should be provided for all significant (greater than trace) alleles.
- 9.1.8 Date issued (TWGDAM 8.3.3)**
The date of issue of the report will be signified.
- 9.1.9 A signature and title, or equivalent identification, of the person(s) accepting responsibility for the content of the report (TWGDAM 8.3.2, 8.3.9)**
The report will be signed by the Analyst(s) expected to testify about the interpretation of the results. If part of the analysis or interpretation was done by an individual other than the primary Analyst, the second Analyst should be identified in the report. The technical reviewer will be identified in the signature block and will initial the report.
- 9.1.10 Disposition of evidence**
The date and manner of return of each item returned to the submitting agency will be recorded in the evidence log. The recommended conditions for long term storage of the evidence, where appropriate, will be marked on the outer packaging of the evidence item. Evidence retained by the DNA Laboratory (including extracted DNA and PCR product) will be documented in the Analyst's notes.

X. TECHNICAL REVIEW PROCEDURES

- 10.1 TECHNICAL REVIEW (TWGDAM 8.4)**
The DNA Laboratory shall conduct administrative and technical reviews of all case files and reports to ensure conclusions and supporting data are reasonable and within the constraints of scientific knowledge.
- Technical Leaders are expected to verify the primary Analyst's reporting and interpretation of the data as a part of the technical review process.
- 10.1.1 Resolution of disagreement between Analysts and reviewer(s) (TWGDAM 8.4)**
Any significant differences between results obtained by each of the two readers must be resolved prior to issuing a final interpretation and report. The two workers will discuss their findings and

determine if there is a basis for agreement. The fact of, and rationale for, any resulting change in interpretation is to be documented in the notes. If, after discussion and review of the data, disagreement remains regarding a particular result, the result is to be reported as inconclusive. (Inconclusive means that insufficient information exists to support any conclusion.)

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

10.1.2 Technical review worksheet

EVALUATION

Sample evaluation

- Description of sample(s)
- Microscopic assessment of sample(s) as appropriate

ANALYSIS

Extraction protocol-sample assessment and QC

- Slot blot: photo, assessment, calculations
- Yield gel photo, assessment, calculations
- Estimation of quantity/concentration of extracted DNA
- Reagent blank(s)
- Substrate control(s)
- Internal QC sample(s)

Analytical protocol

- Amplification worksheet
- Product gel photo, worksheet
- Hybridization/Typing worksheet

Data Interpretation

- Threshold control (C or S dot)
- Internal QC sample(s)
- Amplification (+), (-) controls
- Types called
- Victim/e.cell agreement for differential extraction
- 2nd reader agreement

DOCUMENTATION

- Sample ID continuity
- Lab notes clear, legible, complete
- All pages dated, numbered, initialed
- 2nd reader worksheets dated, initialed
- Chain of custody documentation
- Population frequency printout

REPORT

- Summary and Conclusions
- Conclusions (as conditioned by stated assumptions) are supported by data
- Account of each item received
- All significant examinations in notes included in report
- Chart of DNA profile
- Frequency calculations
- Evidence disposition if other than routine
- Report signed by Analyst, initialed by technical reviewer

XI. PROCEDURES FOR THE RELEASE OF CASE REPORT INFORMATION

Procedures for the release of case report information are found in the BUREAU POLICY MANUAL.

XII. TRAINING AND PROFICIENCY TESTING

12.1 TRAINING OF NEW STAFF

The purpose of this training outline is to provide a consistent approach to training of new staff within the DNA Laboratory. DNA Analysts are expected to undergo training sufficient to acquire and demonstrate an understanding of the principles, usage, and limitations of the procedures and equipment they will use to generate DNA results. Training will culminate in a proficiency test that will demonstrate the Analyst's readiness to perform independent casework. Because of their role as expert witnesses, Analysts must also be familiar with relevant procedures and practices followed in the broader scientific community.

The specific content and length of the initial training period will naturally be a function of the individual Analyst's prior educational background and experience, as well as of the complexity of the assigned tasks. Projects assigned to Analysts in training may vary depending on the DNA Laboratory's current operational and research needs. Training may proceed in stages.

Work done by Analysts in training will be documented and will be periodically reviewed by the designated trainer(s) and supervisor. The designated trainer(s) or supervisor will evaluate the progress of training by assessing the Analyst's knowledge, understanding, and skill in conducting and describing the procedures, and the Analyst's analytical and interpretive competence as demonstrated through practice analysis, proficiency testing, and report writing exercises.

12.2 DNA TRAINING OUTLINE

12.2.1 Safety training

The DNA Analyst will receive safety training in chemical, biological, and emergency procedures before performing independent labwork.

12.2.2 PCR-DNA analysis

The DNA Analyst will receive the following training in PCR-based techniques before being assigned casework.

12.2.2.1 External training

The DNA Analyst will be sent, at Bureau expense, for introductory course (s) in forensic PCR analysis at the FBI Academy, CCI, NFSTC and/or Perkin-Elmer.

12.2.2.2 Internal training

The DNA Analyst will receive in-lab training in PCR procedures and the specific forensic testing systems in use in the DNA Laboratory under the supervision of a qualified Analyst.

12.2.2.3 Practice

The DNA Analyst will practice the relevant PCR-DNA techniques under minimum supervision using known samples to demonstrate analytical accuracy.

12.2.2.4 Independent analysis

The DNA Analyst will analyze non-compromised, non-sperm-containing samples to evaluate the Analyst's ability to obtain reliable, reproducible, and contamination-free results consistently.

12.2.3 Practice analyses for casework

12.2.3.1 Forensic biology, presumptive testing

Any DNA Analyst entering the program without previous serology experience will receive an introduction to serology evidence handling and presumptive tests. The training, conducted under the supervision of a qualified Analyst, will include presumptive screening tests for blood, semen, and saliva and evaluation of case evidence, specifically microscopic evaluation of stains, in particular identification of sperm. The Analyst will receive training in the interpretation and documentation of these procedures, and their utility in evaluating evidence for DNA testing.

12.2.3.2 Sexual assault samples-training

The DNA Analyst will receive training from a qualified Analyst in the differential lysis procedures used to separate sperm from other components in sexual assault sample.

12.2.3.3 Sexual assault samples-practice

The DNA Analyst will practice differential lysis techniques under minimum supervision using known samples to demonstrate analytical accuracy.

12.2.3.4 Non-probative cases

The DNA Analyst will perform evaluation, analysis, documentation and interpretation on samples from previously adjudicated cases, such as sexual assault kits that no longer retain probative value. The analysis will be followed by a formal technical case review.

12.2.4 Other desirable training

The DNA Analyst should participate in the following activities in an ongoing or periodic basis

12.2.4.1 Meetings

The DNA Analyst should attend scientific meetings, in particular those where information about forensic science and forensic DNA is presented.

12.2.4.2 Internal research and development

When possible, the DNA Analyst should participate in method development and validation studies.

12.2.4.3 Continuing education in forensic science

The DNA Analyst shall obtain continuing education and training in all areas of forensic science, and in particular those with a bearing on DNA analysis (e.g. general forensic science, serology, crime scenes, and courtroom testimony).

12.2.4.4 Continuing education in related academic disciplines

The DNA Analyst shall obtain continuing education in molecular biology, population genetics and related academic areas.

12.3 PROFICIENCY TESTING

The purpose of proficiency testing is to monitor the performance both of individual Analysts and the DNA Laboratory as a whole. Proficiency tests provide a mechanism for critical self review and a means by which others, such as the ASCLD/LAB and DNA Laboratory clients, may evaluate the DNA Laboratory's performance on a continual basis. In general, the proficiency testing program of the DNA Laboratory will follow guidelines set forth by the TWGDAM, and ASCLD/LAB.

Because proficiency tests are intended to monitor work as normally performed in the DNA Laboratory, they are to be conducted using the procedures currently approved for casework samples. Work is to be done independently by the Analyst, supported by notes, photographs, and other documentation, and

summarized in a written report as required by the Proficiency test provider. Where allowed, the report format for casework is desirable. Prior to reporting the proficiency test results, the work is to receive the same level of technical review required for casework. (TWGDAM 9.3.1)

12.3.1 Internal proficiency tests

Internal proficiency tests may be prepared from in-house known samples or may use extra samples from previously completed external tests that are re-issued to other Analysts as unknowns. Internal proficiency tests may be completed by Analysts in training before they attempt external proficiency tests to qualify for casework. The "correct results" are to be unknown to both the Analyst and technical reviewer until after the tests are completed and the results are reported.

12.3.2 External proficiency tests

External proficiency tests are obtained from, and reported back to, an independent outside agency. External proficiency samples should reflect the types of analyses the Analyst will be expected to perform, for instance simulated sexual assault samples or other mixtures in addition to single source blood samples.

12.3.2.1 Casework qualification (TWGDAM 2.2.1.3)

The DNA Analyst will successfully complete an appropriate proficiency test prior to being assigned to work independently or serving as a second reader. The tests should reflect the type of sample typically encountered by the DNA casework Analyst, for instance samples requiring differential extraction and interpretation of mixtures.

12.3.2.2 Assignment, frequency and documentation

The DNA Laboratory subscribes to external tests from ASCLD-approved Proficiency Test Providers. (TWGDAM 9.1.3 and 9.1.4) While actively performing DNA case analyses, the DNA Analyst will complete two external tests per year (TWGDAM 9.1.2). Analysts will complete and return results within the applicable deadlines.

The Laboratory Manager is responsible for receipt of proficiency tests from the provider (TWGDAM 9.1.4), assignment of proficiency test samples, and for maintaining documentation of test results, including a master list of proficiency tests conducted in the DNA Laboratory. The Laboratory Manager records whether test results are reported within applicable deadlines. The Technical Leader is responsible for reviewing each proficiency test for technical accuracy and completeness. Any problems noted in a proficiency test will be reported to the Laboratory Manager. The Laboratory Manager, in consultation with the Technical Leader, will administer corrective action should a problem be noted on a proficiency test.

A file is created for each proficiency test (TWGDAM 9.4); these are filed in chronological order by source and retained in the PROFICIENCY TEST FILE according to Bureau policy. The PROFICIENCY TEST FILE contains all analytical data (notes, photos, run sheets, etc.) generated in the analysis, a report of the Analyst's conclusions, and signature and comments of the technical reviewer. In any situation where the results of the test are not satisfactory, significant discrepancies and appropriate corrective action are documented in this file. For external tests, the file also contains the summary report of the test provider for that particular test. (TWGDAM 9.6)

12.3.2.3 Verification of Performance, Corrective Action (TWGDAM 9.5)

The DNA Technical Leader will report in a timely fashion to the Laboratory Manager any substantive discrepancy in a proficiency test, along with recommendations for corrective action. In general, the Analyst will be expected to generate a written response for the file addressing the problem. Upon the DNA Technical Leader's approval of the recommended corrective action, the Laboratory Manager is responsible for ensuring that the corrective action is carried out and documented.

Discrepancies found to be a result of administrative error (e.g. clerical error, sample confusion, or insufficient documentation) may be handled by counseling, remedial training and other supervisory techniques.

Substantive discrepancies found to be the result of a systematic error (equipment, materials, storage) may require a review of related casework since the Laboratory's last successfully completed proficiency test. Once the cause of the problem has been identified, all Analysts should be made aware of any corrective action taken to minimize the recurrence of the discrepancy.

Any substantive discrepancy (e.g. incompetent labwork, incorrect "match", improper interpretation of mixed sample) determined to be the result of an analytical/interpretive problem will prohibit the individual(s) involved from further examination of case evidence until the cause of the problem is identified and corrected. Depending on the nature of the problem, an audit of prior cases may be required. Before resuming casework, the individual(s) responsible for the discrepancy must satisfactorily complete an additional set of proficiency samples.

12.4 PREPARATION FOR COURT TESTIMONY

Effective testimony is an essential component of professional competence in forensic science. Training will include attending court sessions to observe the testimony of others and participation in moot court sessions. An Analyst's preparation for testimony will be assessed prior to court testimony by his/her performance in moot court situations. Subsequent to this initial qualification, testimony will be monitored periodically in person or by transcript review by a supervisor or qualified peer. The DNA Analyst must demonstrate:

- Familiarity with the literature, court transcripts and legal decisions related to forensic DNA testing
- Experience and training in DNA and forensic biology
- Understanding of the scientific principles underlying the procedures used in the DNA Laboratory
- Knowledge of DNA Laboratory quality assurance policies and procedures
- Ability to explain the DNA analysis procedures to a lay audience
- Ability to advise and assist attorneys in presentation of DNA evidence
- Professional appearance and demeanor

XIII. AUDITS

13.1 INTERNAL AUDIT

The Laboratory Manager will conduct and record the results of an annual audit to verify that the DNA Laboratory continues to meet ASCLD/LAB and TWGDAM standards.

13.2 EXTERNAL AUDIT (TWGDAM 10.1)

At least every 2 years an audit should be conducted by individuals separate from the DNA Laboratory. Ideally, at least one auditor should be from an outside agency.

13.3 AUDIT RECORDS (TWGDAM 10.2)

Records of each inspection, including the date, name of inspector, findings and problems, action taken to resolve existing problems, and the next scheduled inspection, will be maintained by the Laboratory Manager.

Report pg. # R ____ of ____
Case # DNA _____
Date _____
Technical Reviewer _____

**Idaho Department of Law Enforcement
Bureau of Forensic Services**

Technical Review Checksheet

EVALUATION

Sample evaluation

- ___ Description of sample(s)
- ___ Microscopic assessment of sample(s) as appropriate

ANALYSIS

Extraction protocol—sample assessment and QC

- ___ Slot blot: photo, assessment, calculations
- ___ Yield gel photo, assessment, calculations
- ___ Estimation of quantity/concentration of extracted DNA
- ___ Reagent blank(s)
- ___ Substrate control(s) as appropriate
- ___ Internal QC sample(s)

Analytical protocol

- ___ Amplification worksheet
- ___ Product gel photo, worksheet
- ___ Hybridization/Typing worksheet

Data Interpretation

- ___ Threshold control assessment (C or S dot)
- ___ Internal QC sample(s) correct
- ___ Amplification (+), (-) controls correct
- ___ Types called correctly
- ___ Victim/e.cell agreement for differential extraction
- ___ 2nd reader agreement

DOCUMENTATION

- ___ Sample ID continuity
- ___ Lab notes clear, legible, complete
- ___ All pages dated, numbered, initialed
- ___ 2nd reader worksheets dated, initialed
- ___ Population frequency printout
- ___ Evidence disposition if other than routine

REPORT

- ___ Summary and Conclusions
- ___ Conclusions (as conditioned by stated assumptions) are supported by data
- ___ Account of each item received
- ___ All significant examinations in notes included in report
- ___ Chart of DNA profile
- ___ Frequency calculations

DQA1 TYPING SHEET

Date/Analyst _____
DQA1 Strip Lot # _____
Hybe Lot # _____
Wash Lot # _____

Enzyme Conj Lot # _____
Hybe °C _____ Hybe Time _____
Wash °C _____ Wash Time _____

Second Reader _____
Chromogen Lot # _____
H₂O₂ Lot # _____
Citrate Buffer Lot # _____

SAMPLE	1	2	3	4	C	1.1	1.2 1.3 4	1.3	All but 1.3	4.1	4.2 4.3	RESULTS	TYPE
												<i>observed</i>	
												<i>all alleles identified</i>	<i>Mixture etc major minor</i>

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

POLYMARKER TYPING

Date/Analyst _____
 PM Strip Lot # _____
 Hybe Lot # _____
 Wash Lot # _____

Enzyme Conj Lot # _____
 Hybe °C _____ Hybe Time _____
 Wash °C _____ Wash Time _____

Second Reader _____
 Chromogen Lot # _____
 H₂O₂ Lot # _____
 Citrate Buffer Lot # _____

SAMPLES	S Dot	LDLR		GYPA		HBGG			D7S8		GC			COMMENTS
		A	B	A	B	A	B	C	A	B	A	B	C	
Amp date / /														
						5+		5-						
Result						A → C								
Type						A major C minor								
Result														
Type														
Result														
Type														
Result														
Type														
Result														
Type														
Result														
Type														
Result														
Type														

Property of Idaho State Police Forensic Services
 Uncontrolled Internet Copy
 OBSOLETE DOCUMENT

QC samples summary for _____

Blind Internal Quality Control Sample (Bloodstain)

Identity	DQA1	LDLR	GYP A	HBGG	D7S8	GC
Correct profile						

Reviewer

Date

Differential extraction positive control

Identity	DQA1	LDLR	GYP A	HBGG	D7S8	GC
RMT #2	1.2,4.1	BB	AB	BB	AB	AC
RMT #4	2.4.1	AB	AA	BB	AB	AC
RMT #5	1.1,1.3	AB	AB	AB	AA	AA

Do Not
use
used 9702
control

Amplification Positive Control

Identity	DQA1	LDLR	GYP A	HBGG	D7S8	GC
Correct profile	1.1,4.1	BB	AB	AA	AB	BB

DQA1 TYPING SHEET

Date/Analyst _____ Strip color development time _____
 DQA1 Strip Lot # _____ Enzyme Conj Lot # _____ Chromogen Lot # _____
 Reagent Lot # _____ Hybe °C _____ Hybe Time _____ H₂O₂ Lot # _____
 Wash Lot # _____ Wash °C _____ Wash Time _____ Citrate Buffer Lot # _____

SAMPLE	1	2	3	4	C	1.1	1.2 1.3 4	1.3	All but 1.3	4.1	4.2 4.3	RESULTS	TYPE
			/	tr	/	+							Report "allele detected 1.1.2"
													To emphasize it's not a type as the
													common says
													Major (type) follow allele
													Don't list names here
													Be careful about
													1.2 caveat
													do intensity of components for two add up?

ASCLD/LAB DNA Proficiency Test Record

Laboratory Name: _____

Contact Person: _____

Laboratory Accreditation Number: _____

Proficiency Test Provider and Test #: _____

Unique Laboratory Test Code #: _____

Number of Analysts Working on Test: _____

Loci Tested (check the boxes)

RFLP

DIS7	D2S44	D4S139	D5S110	D7S467
D10S28	D17S26	D17S79	Other	Other

PCR

DQA1	PM	D18S0	Other	Other

STR

TH01	TPOX	CSF1PO	VWA	FES/FPS
F13A1	VWF	AMELOGENIN	D5S818	D13S317
D7S80	D3S1358	FGA	D8S1179	D18S51
D21S11	Other	Other	Other	Other

Send to PRC Chair: Beth Carpenter
 1111 S.W. 2nd Room 1201
 Portland, OR 97204

For PRC use only: Correct Typing Obtained _____	Date _____
Discrepancies Noted _____	

ASCLD/LAB DNA Proficiency Test Record

Laboratory Name: _____

Contact Person: _____

Laboratory Accreditation Number: _____

Proficiency Test Provider and Test #: _____

Unique Laboratory Test Code #: _____

Number of Analysts Working on Test: _____

Loci Tested (check the boxes)

RFLP

D1S7	D2S44	D4S139	D5S110	D7S467
D10S28	D17S26	D17S79	Other	Other

PCR

DQA1	PM	D1S80	Other	Other

STR

TH01	TPOX	CSF1PO	VWA	FES/FPS
F13A1	VWF	AMELOGENIN	D5S818	D13S317
D7S80	D3S1358	FGA	D8S1179	D18S51
D21S11	Other	Other	Other	Other

Send to PRC Chair: Beth Carpenter
 1111 S.W. 2nd Room 1201
 Portland, OR 97204

For PRC use only: Correct Typing Obtained _____	Date: _____
Discrepancies Noted _____	

ASCLD/LAB DNA Proficiency Test Record

Laboratory Name: _____

Contact Person: _____

Laboratory Accreditation Number: _____

Proficiency Test Provider and Test #: _____

Unique Laboratory Test Code #: _____

Number of Analysts Working on Test: _____

Loci Tested (check the boxes)

RFLP

DIS7	D2S44	D4S139	D5S110	D7S467
D10S28	D17S26	D17S79	Other	Other

PCR

DQA1	PM	D1S80	Other	Other

STR

TH01	TPOX	CSF1PO	VWA	FES/FPS
F13A1	VWF	AMELOGENIN	D5S818	D13S317
D7S80	D3S1358	FGA	D8S1179	D18S51
D21S11	Other	Other	Other	Other

Send to PRC Chair: Beth Carpenter
 1111 S.W. 2nd Room 1201
 Portland, OR 97204

For PRC use only: Correct Typing Obtained _____	Date: _____
Discrepancies Noted _____	

ASCLD/LAB DNA Proficiency Test Record

Laboratory Name: _____

Contact Person: _____

Laboratory Accreditation Number: _____

Proficiency Test Provider and Test #: _____

Unique Laboratory Test Code #: _____

Number of Analysts Working on Test: _____

Loci Tested (check the boxes)

RFLP

DIS7	D2S44	D4S139	D5S110	D7S467
D10S28	D17S26	D17S79	Other	Other

PCR

DQA1	PM	D1S80	Other	Other

STR

TH01	TPOX	CSF1PO	VWA	FES/FPS
F13A1	VWF	AMELOGENIN	D5S818	D13S317
D7S80	D3S1358	FGA	D8S1179	D18S51
D21S11	Other	Other	Other	Other

Send to PRC Chair: Beth Carpenter
 1111 S.W. 2nd Room 1201
 Portland, OR 97204

For PRC use only: Correct Typing Obtained _____	Date: _____
Discrepancies Noted _____	

ASCLD/LAB DNA Proficiency Test Record

Laboratory Name: _____

Contact Person: _____

Laboratory Accreditation Number: _____

Proficiency Test Provider and Test #: _____

Unique Laboratory Test Code #: _____

Number of Analysts Working on Test: _____

Loci Tested (check the boxes)

RFLP

DIS7	D2S44	D4S139	D5S110	D7S467
D10S28	D17S26	D17S79	Other	Other

PCR

DQA1	PM	D1S80	Other	Other

STR

TH01	TPOX	CSF1PO	VWA	FES/FPS
F13A1	VWF	AMELOGENIN	D5S818	D13S317
D7S80	D3S1358	FGA	D8S1179	D18S51
D21S11	Other	Other	Other	Other

Send to PRC Chair: Beth Carpenter
 1111 S.W. 2nd Room 1201
 Portland, OR 97204

For PRC use only: Correct Typing Obtained _____	Date: _____
Discrepancies Noted _____	

From: nrbiocom@uclink4.berkeley.edu ("N. Rudin")
To: DLE.DLE-PC(abradley)
Date: 9/11/98 3:25pm
Subject: 1.2 disclaimer

found it

"Due to limitations in the DQA1 test system, the presence of a 1.1, 1.3 or a 4 allele in a mixed sample may prevent the conclusive determination of a 1.2 allele."

I like this better - among other things, the "may" as opposed to the "cannot" is important. For instance we call a 1.2 in H2 sub.

N

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

POLYMARKER TYPING

Date/Analyst _____
 PM Strip Lot # _____
 Hybe Lot # _____
 Wash Lot # _____

Enzyme Conj Lot # _____
 Hybe °C _____ Hybe Time _____
 Wash °C _____ Wash Time _____

Chromogen Lot ## _____
 H₂O₂ Lot# _____
 Citrate Buffer Lot # _____
 Color develop. time _____

SAMPLES	S Dot	LDLR A B	GYP A A B	HBGG A B C	D7S8 A B	GC A B C	COMMENTS
Amp date / /							
Result							
Type							
Result							
Type			Major type Minor allele				
Result							
Type							
Result							
Type							
Result							
Type							
Result							
Type							
Result							
Type							
Result							
Type							

Property of Idaho State Police Forensic Services
 Uncontrolled Internet Copy
 OBSOLETE DOCUMENT

From: Norah@forensicdna.com ("N. Rudin")
To: DLE.DLE-PC(abradley)
Date: 1/27/99 7:05pm
Subject: 1957

Mornin'

Here are my comments on 1957.

****For immediate consumption****

1. Analysis and general conclusion look fine.
2. Wording in Summary and for Evidence received and Examination and results is very good. You might want to reserve it somewhere for easy cut and paste.
3. Interpretation wording. I'm a bit uncomfortable with the fractions "contributing" DNA to each other ? feels a bit anthropomorphic. You and I know what you mean, but Also, you never mention the major component of the sperm fraction. Some of this may be cut&paste in a hurry. I don't think you've used this particular wording before.

My suggestion for the first two paragraphs:

In the sperm cell fraction from the sheet cutting, DNA from more than one source was detected at the DQA1 locus. The major type is the same as that of Jerry Palmer as are the types detected at all the other loci tested. The minor component at the DQA1 locus corresponds to Wendi Brandon's type and is consistent with incomplete separation of the sperm an non-sperm cells during differential extraction.

The major profile exhibited in the non-sperm fraction is the same as that of Wendi Brandon. The minor alleles seen at the DQA1, LDLR, and GYPA loci are all contained in Jerry Palmer's profile and are consistent with incomplete separation of the sperm an non-sperm cells during differential extraction.

Third paragraph etc. is fine.

****For later consumption**:**

Noting the yellow sticky about 1.2. This is actually a different situation. It occurs to me that I go through a pretty specific process in my head and it might be worth putting into writing. Tell me if this makes sense to you.

A) Is there any evidence of a mixture in the evidence sample? (w/o any comparison to the reference sample) If not, I don't feel the need to mention possible cryptic 1.2 on typing sheet.

B) Is there any evidence of a mixture in the evidence sample? If yes (could be extra alleles, intensity differences, either in DQA1 or PM),

then, depending on the specific alleles showing up in the evidence sample (still w/o any comparison to the reference sample), the it might be prudent to suggest the possibility of a cryptic 1.2 allele. Note that this is instigated by the suggestion of a possible mixture in the evidence, not any reference type.

C) Now, comparing the referene and evidence profiles, note if minor blatent alleles corresponds to some probitive reference sample and if, in reference sample, it is paired with a 1.2. Only then would you say something on the order of "this person would not excluded if the evidence contained a 1.2 allele and ?the intensity of the trio suggests an extra allele? or ?the presence/absence of the 1.2 is inconclusive? or ?etc.?"

Finally, for rush cases, you might want to consider specifying "first morning delivery" w/ FedEx so I get it in the morning. With regular ON delivery, I don't get it unti about 2PM, making it virtually impossible to get it back out the same day.

Keith says Barry Brown is the one to call, so I will make that a priority for the rest of this week as well as getting all the other cases signed and out.

Sooner rather than later

Norah

--

Norah Rudin, Ph.D.
Forensic DNA Consultant
norah@forensicdna.com
<http://www.forensicdna.com>

Criminalistics Bookstore
<http://www.forensicdna.com/Bookstore/index.html>

DOT BLOT EVALUATION AND INTERPRETATION GUIDELINES

The goal of the evaluation and interpretation of dot blot strips is to determine the DNA profile(s) of the donor(s) of the questioned samples. The interpretation of the data is to be based on the education, training and experience of a qualified analyst operating within the following guidelines.

- In general, a dot is defined as a circular-shaped blue area on a membrane strip. Its edge should be distinct all around the circumference.
- Dot intensities are recorded in relationship to the control dot (i.e., "C" or "S") present on the strip. Dots reacting stronger than the control dot will be recorded as "C+" or "S+". Dots less than the control dot, but not significantly so, will be recorded as "C-" or "S-". Dots reacting significantly less than the control dot will be recorded as "trace".
- Signals significantly weaker than the control dot may not be interpreted reliably. If the control dot is not visible on the strip, the amplification/hybridization may not have detected one or more alleles present in the sample. Generally, conclusions as to genotype will not be drawn from strips without a visible control dot.
- The positive amplification control and QC sample must give the expected results. If the control dot on the positive amplification control strip is not visible but the correct alleles are present, and it is shown that the quantitation of the control DNA is low, the results may be reported as long as the QC sample strip has a visible control dot and gives the correct genotype.
- Occasionally a very weak signal, much weaker than the "C" dot, is seen at the "1.1" dot of a DQA1 strip on which other signals are strong, even when there is no dot present for the "1" probe. This may be due to the relatively low-level amplification of some alleles of a related pseudogene, HLA DQA2 (formerly DX α).
- The "1.2, 1.3, 4" dot can be lighter than the "C" dot when the genotype has a DQA1 4.2 or 4.3 allele paired with a 1.1, 2, 3, 4.2, or 4.3 allele. This is because the 4.2 and 4.3 alleles each have a single partially destabilizing mismatch to the "1.2, 1.3, 4" probe.
- The "All but 1.3" dot can be equal to or lighter than the "C" dot when the genotype has a DQA1 1.3 allele paired with a 4.1, 4.2, or 4.3 allele. This is because the 4.1, 4.2, and 4.3 alleles have a single partially destabilizing mismatch to the "All but 1.3" probe.

The typing probes have been designed to bind specifically to PCR product containing the alleles designated on the DQA1 and PM strips under the hybridization conditions. However, weak signals due to cross-hybridization are occasionally observed. This phenomenon is most pronounced at the dots for the sub-types of the DQA1 1 allele, the HBGG A allele, and the GC A allele. If signals greater than "trace" level are thought to be caused by cross-hybridization, the typing should be repeated.

PM + DQA1 PCR AMPLIFICATION PROTOCOL

This analysis is performed using the Amplitype® PM + DQA1 Amplification and Typing Kit from Perkin-Elmer. Using a dedicated pipetter and working in the biological hood, put 40 µL aliquots of Reaction Mix from the kit into sterile GeneAmp® tubes. Store in the refrigerator.

(Turn on the Thermal Cycler. Work in the biological hood to set up the amplification tubes.)

1. Make appropriate dilutions of samples to total 20 µL using sterile water. (Sample prep: for Chelex extract, vortex and centrifuge 3 min. at max. g.) The samples will typically contain 1-10 ng of DNA.
2. Remove tubes containing 40 µL of reaction mixture from the refrigerator. Pulse spin to ensure all liquid contents are in the bottom of the tube.
3. Label these tubes with sample numbers plus amplification ⊕ and ⊖ controls.
4. Prepare a master mix for all tubes plus one:
 - 40 µL of primer reagent
 - 0.4 µL of 4 % BSAAdd 40 µL of the master-mix to each reaction mixture tube.
At this time, the amplification must be completed within 20 minutes.
5. Add two drops of mineral oil to the tubes, cap each tube loosely.
6. Add 20 µL of sample, 20 µL of kit DNA std. for ⊕ amp control, and 20 µL of sterile H₂O for ⊖ amp control to the appropriately labeled tubes. (BELOW MINERAL OIL LAYER)
7. Place tubes into cardboard tray, and take into product room for amplification.
8. Place samples into wells, checking that lids are tightly on tubes. Use worksheet for documentation. Start the amplification program: Push "FILE 40" then "ENTER" and then push "START". Watch the first cycle to ensure instrument is properly functioning. The amplification conditions for this kit are 32 cycles as follows:
 - 94°C denaturation for 1 minute
 - 60°C annealing for 30 seconds
 - 72 °C extension for 30 secondsAfter the last cycle, extend at 72°C for 7 minutes, then hold at 4 °C
The amplification takes approximately 2.5 hours. When the process is completed, the program will hold the samples at 4°C indefinitely.

Adopted: 11/97

APPENDIX B

EXTRACTION REAGENTS

5% Chelex

5 gram Chelex
100 mL sterile H₂O

***0.5 M EDTA**

186.1 g Na₂EDTA·2H₂O dissolve in 800 mL of DI H₂O
adjust pH to 8.0 with 15-20 g NaOH pellets
adjust volume to 1.0 L

Ethidium Bromide (EtBr)

Purchase pre-prepared (10 mg/mL solution)

***5M NaCl**

292.2 g NaCl
dissolve in DI H₂O and adjust volume to 1.0 L

Nuclear Fast Red Stain

5 g aluminum sulfate
0.1 g Nuclear Fast Red
adjust volume with 100 mL hot dH₂O
Refrigerate overnight then filter

Phenol/chloroform/isoamyl alcohol (25/24/1)

Purchase pre-prepared

***PBS (pH 7.4)**

Purchase pre-prepared; prepare according to manufacturer's instructions if concentrated or pre-weighed, OR make according to the following:

8 g sodium chloride
0.2 g potassium chloride
1.44 g Na₂HPO₄, anhydrous
0.24 g KH₂PO₄, anhydrous
adjust pH to 7.4 with HCL
adjust volume to 1.0 L with DI H₂O

*** Solutions that must be autoclaved**

Appendix B Page 2

Picroindigocarminic Stain

0.33 g indigocarminic
100 mL saturated Picric Acid solution
filter

Proteinase K (20 mg/mL)

500 mg proteinase K
dissolve in 25 mL DI H₂O
aliquot into convenient size volumes and freeze

20% SDS (200 g sodium dodecyl sulfate/ 1.0 L H₂O)

Purchase pre-prepared

Stain Extraction Buffer (10 mM Tris; 10 mM EDTA; 100 mM NaCl; 2% SDS)

10 mL of 1M Tris, pH 8.0
20 mL of 0.5M EDTA, pH 8.0
20 mL 5M NaCl
100 mL 20% SDS
adjust volume to 900 mL with DI H₂O
after autoclaving, add 100 mL 20% SDS

***1M Tris, pH 7.5**

121.1 g Tris base
dissolve in 800 mL DI H₂O
adjust pH to 7.5 with concentrated HCl
adjust volume to 1.0 L with DI H₂O

***1M Tris, pH 8.0**

121.1 g Tris base
dissolve in 800 mL DI H₂O
adjust pH to 8.0 with concentrated HCl
adjust volume to 1.0 with DI H₂O

***TE⁻⁴ (10 mM Tris; pH 7.5; 0.1 mM EDTA)**

10 mL 1 M Tris, pH 7.5
200 µL 0.5 M EDTA, pH 8.0
adjust to 1.0 L with DI H₂O

Appendix B Page 3

QUANTITATION REAGENTS (for yield gel and QuantiBlot™ kit)

NOTE: Some reagents used for quantitation have been described in the previous section.

Citrate Buffer (0.1 M Sodium Citrate, pH 5.0)

Purchase pre-prepared concentrate then dilute according to manufacturer's directions OR make according to the following:

18.4 g trisodium citrate dihydrate ($\text{Na}_3\text{C}_6\text{H}_5\text{O}_7 \cdot 2\text{H}_2\text{O}$)

dissolve in 800 mL DI H_2O

adjust pH to 5.0 (± 0.2) by the addition of 6g citric acid monohydrate ($\text{C}_6\text{H}_8\text{O}_7 \cdot \text{H}_2\text{O}$)

adjust volume to 1.0 L with DI H_2O

mix thoroughly

Hybridization Solution (5X SSPE; 0.5% w/v SDS)

250 mL

20X SSPE

25 mL 20% SDS

adjust to 1.0 L with 725 mL DI H_2O

mix thoroughly

NOTE: Hybridization solution solids must be in solution before use; warming may be required to dissolve solids completely. Preparation in a clear glass container is recommended to facilitate visual inspection for solids during warming.

5N NaOH

20 g NaOH pellets

dissolve in 60 mL of DI H_2O

adjust volume to 100 mL

Pre-wetting Solution (0.4N NaOH; 25 mM EDTA)

40 mL of 5N NaOH

25 mL of 0.5 M EDTA

add to 435 mL of DI H_2O

mix thoroughly

QuantiBlot Wash Solution (1.5X SSPE; 0.5% w/v SDS)

150 mL of 20X SSPE

50 mL of 20% SDS

add to 1,800 mL of DI H_2O

mix thoroughly

NOTE: Wash solution solids must be in solution before use; warming may be required to dissolve solids completely. Preparation in a glass container is recommended to facilitate visual inspection for solids during warming.

Appendix B Page 4

Spotting Solution (0.4 N NaOH; 25 mM EDTA; 0.00008% Bromphenol Blue)

2 mL of 5 N NaOH

1.25 mL of 0.5 M EDTA

50 μ L 0.04% Bromphenol Blue (provided in QuantiBlot™ kit)

add to 22 mL DI H₂O

mix thoroughly

20 X SSPE (3.6 M NaCl; 0.2 M NaH₂PO₄; 20 mM EDTA)

Purchase pre-prepared OR make according to the following:

210.4 g NaCl

27.6 g NaH₂PO₄·H₂O

400 mL 0.5 M EDTA or 7.4 g of Na₂EDTA·2H₂O

adjust pH to 7.4 with NaOH

adjust volume to 1.0 L with DI H₂O

1X TBE (0.1 M Tris; 0.09M Boric Acid; 0.001M EDTA; pH 8.4)

Purchase pre-prepared concentrate then dilute according to the manufacturer's instructions

Yield Gel

add .25 g LE agarose to 25 mL of 1X TBE buffer

melt agarose

add 1.25 μ L of Ethidium Bromide (10 mg/mL)

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

AmpliType® PERKIN ELMER

Part No. N808-0057
 Part No. N808-0094

PM PM+DQA1 PCR Amplification and Typing Kits

**FOR FORENSIC OR RESEARCH USE ONLY
 NOT FOR USE IN DIAGNOSTIC PROCEDURES**
 See notice to purchaser

READ ENTIRE PROCEDURE PRIOR TO USE OF KIT

TABLE OF CONTENTS

Notice to Purchaser: Limited License	1
Storage and Stability	2
Introduction	2
List of Kit Components	6
Materials Required But Not Supplied	9
1.0 Reagent Preparation	12
2.0 Precautions	15
3.0 PCR Amplification Protocols	16
4.0 Preparation of PCR Product for Detection	20
5.0 Verification of PCR Amplification	20
6.0 Preparation and Pouring of Agarose Gel	20
7.0 Gel Loading and Electrophoresis	21
8.0 DNA Hybridization	22
9.0 Color Development	27
10.0 Photography and Storage	28
11.0 Disposal or Reuse of Typing Trays	28
12.0 Interpretation of Results	28
Performance Characteristics	32
Troubleshooting	33
References	40
Ordering Information	44

NOTICE TO PURCHASER: LIMITED LICENSE

A license under U.S. Patents 4,683,202, 4,683,195 and 4,965,188 or their foreign counterparts, owned by Hoffmann-La Roche Inc. and F. Hoffmann-La Roche Ltd ("Roche"), for use in research and forensic testing, has an up-front fee component and a running-royalty component. The purchase price of this product includes limited, nontransferable rights under the running-royalty component to use only this amount of the product to practice the Polymerase Chain Reaction ("PCR") and related processes described in said patents solely for the forensic testing and research activities of the purchaser when this product is used in conjunction with a thermal cycler whose use is covered by the up-front fee component. Rights to the up-front fee component must be obtained by the end user in order to have a complete license. These rights under the up-front fee component may be purchased from Perkin-Elmer or obtained by purchasing an Authorized Thermal Cycler. The right to use this product to perform and to offer commercial services for forensic testing, using PCR, including reporting results of purchaser's activities for a fee or other commercial consideration, is also hereby granted. Further information on purchasing licenses to practice the PCR process may be obtained by contacting the Director of Licensing at The Perkin-Elmer Corporation, 850 Lincoln Centre Drive, Foster City, California 94404 or at Roche Molecular Systems, Inc., 1145 Atlantic Avenue, Alameda, California 94501.

ORDERING INFORMATION CTS ARE NOT FOR USE IN DIAGNOSTIC PROCEDURES

<p>In the United States: For PCR Technical Support call toll-free 1-800-762-4001 FAX: 203-791-2542</p> <p>To order PCR consumables and reagents call toll-free 1-800-327-3002</p> <p>To order PCR instruments or to reach your sales representative call toll-free 1-800-345-7547</p> <p>For PCR product literature call toll-free 1-800-345-5224</p>	<p>Outside of the United States: For the number of your local Perkin-Elmer sales representative call:</p> <table border="0"> <tr><td>Australia</td><td>61 3 212 8595</td></tr> <tr><td>Austria</td><td>431 1 602 31010</td></tr> <tr><td>Belgium</td><td>32 2 725 9099</td></tr> <tr><td>Bulgaria</td><td>359 2 59 1118</td></tr> <tr><td>Canada</td><td>800 669 6913</td></tr> <tr><td>Czechoslovakia</td><td>42 2 772507</td></tr> <tr><td>Denmark</td><td>45 48 141391</td></tr> <tr><td>Finland</td><td>358 0880 144</td></tr> <tr><td>France</td><td>33 1 4990 1800</td></tr> <tr><td>Germany</td><td>49 6150 1010</td></tr> <tr><td>Hungary</td><td>36 1 251 1116</td></tr> <tr><td>Italy</td><td>39 39 23 831</td></tr> <tr><td>Japan</td><td>81 4 7380 8500</td></tr> <tr><td>Latin America-Mexico</td><td>52 5 657 7077</td></tr> <tr><td>Netherlands</td><td>31 1820 754 11</td></tr> <tr><td>Norway</td><td>47 22 68 6565</td></tr> <tr><td>Poland</td><td>48 22 33 0996</td></tr> <tr><td>Portugal</td><td>351 1 796 2172</td></tr> <tr><td>Spain</td><td>34 1 803 4210</td></tr> <tr><td>Sweden</td><td>46 87 3300 10</td></tr> <tr><td>Switzerland</td><td>41 42 657 777</td></tr> <tr><td>United Kingdom</td><td>44 192 582 5650</td></tr> </table>	Australia	61 3 212 8595	Austria	431 1 602 31010	Belgium	32 2 725 9099	Bulgaria	359 2 59 1118	Canada	800 669 6913	Czechoslovakia	42 2 772507	Denmark	45 48 141391	Finland	358 0880 144	France	33 1 4990 1800	Germany	49 6150 1010	Hungary	36 1 251 1116	Italy	39 39 23 831	Japan	81 4 7380 8500	Latin America-Mexico	52 5 657 7077	Netherlands	31 1820 754 11	Norway	47 22 68 6565	Poland	48 22 33 0996	Portugal	351 1 796 2172	Spain	34 1 803 4210	Sweden	46 87 3300 10	Switzerland	41 42 657 777	United Kingdom	44 192 582 5650
Australia	61 3 212 8595																																												
Austria	431 1 602 31010																																												
Belgium	32 2 725 9099																																												
Bulgaria	359 2 59 1118																																												
Canada	800 669 6913																																												
Czechoslovakia	42 2 772507																																												
Denmark	45 48 141391																																												
Finland	358 0880 144																																												
France	33 1 4990 1800																																												
Germany	49 6150 1010																																												
Hungary	36 1 251 1116																																												
Italy	39 39 23 831																																												
Japan	81 4 7380 8500																																												
Latin America-Mexico	52 5 657 7077																																												
Netherlands	31 1820 754 11																																												
Norway	47 22 68 6565																																												
Poland	48 22 33 0996																																												
Portugal	351 1 796 2172																																												
Spain	34 1 803 4210																																												
Sweden	46 87 3300 10																																												
Switzerland	41 42 657 777																																												
United Kingdom	44 192 582 5650																																												

GeneAmp® PCR Instrument Systems from Perkin-Elmer

- 1801-0100 DNA Thermal Cycler 480, 120 V, 60 Hz
- 1801-0101 DNA Thermal Cycler 480, 240 V, 50 Hz
- 1801-0102 DNA Thermal Cycler 480, 100 V, 50/60 Hz
- 1801-0001 GeneAmp® PCR System 9600, 120 V, 60 Hz
- 1801-0002 GeneAmp® PCR System 9600, 240 V, 50 Hz
- 1801-0003 GeneAmp® PCR System 9600, 100 V, 50/60 Hz

PERKIN ELMER



Perkin-Elmer PCR reagents are developed and manufactured by Roche Molecular Systems, Inc. Branchburg, New Jersey, U.S.A.

Applied Biosystems Division
 850 Lincoln Centre Drive
 Foster City, CA 94404

STORAGE AND STABILITY

1. Store all of the components of the AmpliType® PM PCR Amplification and Typing Kit and the AmpliType PM+DQA1 PCR Amplification and Typing Kit at 2° to 8°C. Isolate the kits from any sources of contaminating DNA, especially amplified PCR product.
2. Store the AmpliType DNA Probe Strips, protected from light (e.g., wrapped in foil), at 2° to 8°C. Store the DNA probe strips with the desiccant in the glass tube and ensure that the screw-cap is securely tightened.
3. Under these conditions, components of the kits are stable through the control date printed on the label.
4. The Chromogen Solution is stable for 6 months after its preparation when stored at 2° to 8°C.

INTRODUCTION

One of the most important developments in the field of human identity testing is the use of DNA typing to analyze biological evidence.^{1,2} In particular, the powerful GeneAmp Polymerase Chain Reaction (PCR)³ is used to analyze samples which cannot be typed by other methods, such as samples containing minute amounts of human DNA and very old and/or degraded DNA.⁴⁻⁸ The AmpliType HLA DQ α PCR Amplification and Typing Kit (AmpliType HLA DQ α Kit) was the first PCR-based test applied to forensic casework analysis.⁹⁻¹² The AmpliType PM PCR Amplification and Typing Kit (AmpliType PM Kit, Perkin-Elmer Part No. N808-0057) was the second commercially available human identity testing product based on the reverse dot blot typing technology¹³ for forensic casework analysis.¹⁴⁻¹⁷ The AmpliType PM+DQA1 PCR Amplification and Typing Kit (AmpliType PM+DQA1 Kit, Part No. N808-0094) has been added to the product line to enable amplification and typing of all five of the PM loci and the HLA DQA1 locus from the same PCR amplification reaction.

The AmpliType PM Kit and AmpliType PM+DQA1 Kit are optimized for amplification in fixed-volume reactions using the Perkin-Elmer GeneAmp® PCR Instrument System and reaction tube combinations listed in Table 1. Typing of the DNA amplified by these kits is performed by hybridization of the amplified PCR products to DNA probe strips for the PM loci or the HLA DQA1 locus. Specifically hybridized amplified DNA is visualized upon enzymatic conversion of a soluble, colorless substrate to a blue-colored precipitate. Descriptions of the amplification and typing protocols for each of these kits are provided in the individual kit summaries below. Additional information can be found in the AmpliType User Guide.¹⁸

Table 1. Recommended GeneAmp® PCR Instrument System and Reaction Tube Combinations for Amplification.

Instrument	Reaction Tube	Perkin-Elmer Part Number
DNA Thermal Cycler 480	GeneAmp® Thin-Walled Reaction Tubes*	N801-0537
	GeneAmp Autoclaved Thin-Walled Reaction Tubes	N801-0611
	GeneAmp Thin-Walled Reaction Tube with Fat Cap*	N801-0737
GeneAmp PCR System 9600	MicroAmp® Reaction Tube with Cap*	N801-0540
	MicroAmp Autoclaved Reaction Tube with Cap	N801-0612

*Tubes should be autoclaved by the user prior to use in the amplification process.

AMPLIType® PM PCR AMPLIFICATION AND TYPING KIT DESCRIPTION

The AmpliType® PM PCR Amplification and Typing Kit (AmpliType PM Kit, Part No. N808-0057) includes PCR amplification reagents that direct the simultaneous amplification of specific regions of the following six genetic loci: Human Leukocyte Antigen (HLA) DQA1 (previously referred to as HLA DQ α),¹⁹ Low Density Lipoprotein Receptor (LDLR),²⁰ Glycophorin A (GYPA),²¹ Hemoglobin G Gammaglobin (HBGG),²² D7S8,²³ and Group Specific Component (Gc).^{24,25} The AmpliType PM Kit also contains detection reagents and DNA probe strips for typing LDLR, GYPA, HBGG, D7S8, and Gc but not for typing HLA DQA1 (see AmpliType PM+DQA1 Kit description for information on typing all six loci).

While the AmpliType PM Kit uses the amplification and reverse dot blot typing technologies introduced in the original AmpliType HLA DQ α Kit, minor changes to the amplification and typing protocols were made because of the simultaneous analysis of multiple genetic loci. AmpliType PM PCR amplification is performed in a fixed-volume reaction (100 μ L total volume) instead of in the variable-volume reaction (102-140 μ L total volume) used with the AmpliType HLA DQ α Kit. The typing of DNA amplified with AmpliType PM reagents requires a separate step to bind the horseradish peroxidase-streptavidin (HRP-SA) enzyme conjugate to the biotinylated PCR product. HRP-SA Enzyme Conjugate is bound to PCR products that have hybridized to the probes only after the excess biotinylated primers and unbound PCR products are removed. When using DNA amplified with the original AmpliType HLA DQ α reagents, hybridization of the HLA DQ α PCR product to the probes and binding of the HRP-SA Enzyme Conjugate to the biotinylated PCR products are performed simultaneously. PCR products amplified with the AmpliType PM PCR Reaction Mix and Primer Set should not be typed using the original AmpliType HLA DQ α protocol because all dot intensities on the AmpliType PM DNA Probe Strip will be too light.

minor modifications to the original AmpliType® HLA DQα typing protocol are. The addition of disodium EDTA to the PCR product (see Section 4.0); 2) the reduction of the volume of the denatured PCR product that is added to the Hybridization Step from 35 µL to 20 µL (see Section 8.8); and 3) the reduction of the volume of the Wash Solution, the Citrate Buffer, and the Color Development Solution from 10 mL to 5 mL (see Sections 8 and 9).

The AmpliType HLA DQα DNA Probe Strips, supplied in the original AmpliType HLA DQα Kit, should not be used to obtain the HLA DQA1 type of samples amplified with the AmpliType PM PCR Reaction Mix and Primer Set because the "C" dot on these DNA probe strips may not be visible with weakly amplified samples. To obtain both the PM types and the HLA DQA1 type from samples amplified using the AmpliType PM reagents, the AmpliType PM+DQA1 PCR Amplification and Typing Kit should be used (see below).

AMPLIType PM+DQA1 PCR AMPLIFICATION AND TYPING KIT DESCRIPTION

The AmpliType PM+DQA1 PCR Amplification and Typing Kit (AmpliType PM+DQA1 Kit, Part No. N808-0094) contains all of the amplification and detection reagents for typing the five PM loci and the HLA DQA1 locus from a single aliquot of extracted DNA. PCR amplification is performed using the same amplification reagents that are included in the AmpliType PM Kit. These reagents direct the simultaneous amplification of the HLA DQA1, LDLR, GYPA, HBG, D7S8, and GC loci. The AmpliType PM+DQA1 PCR setup and amplification protocols are identical to those provided in the AmpliType PM Kit. The AmpliType PM+DQA1 Kit also contains both the AmpliType PM and AmpliType HLA DQA1 DNA Probe Strips. The AmpliType PM DNA Probe Strips are identical to those included in the AmpliType PM Kit. The AmpliType HLA DQA1 DNA Probe Strips contain the nine probes found on the original AmpliType HLA DQα DNA Probe Strips that distinguish six common HLA DQA1 alleles (DQA1 1.1, 1.2, 1.3, 2, 3, 4).¹⁹ The AmpliType HLA DQA1 DNA Probe Strips contain two additional probes that distinguish the DQA1 4.1 allele from the DQA1 4.2 and 4.3 alleles.²⁸ The AmpliType PM+DQA1 typing protocol is identical to the one used with the AmpliType PM Kit. PCR products amplified with the AmpliType PM PCR Reaction Mix and Primer Set should not be typed using the original AmpliType HLA DQα protocol because all dot intensities on the AmpliType PM and AmpliType HLA DQA1 DNA Probe Strips will be too light.

AMPLIType PM AND HLA DQA1 GENETIC MARKER CHARACTERISTICS AND POPULATION DATA

The chromosomal location, size of the amplified DNA product, number of alleles distinguished (Table 2), and allele frequencies (Table 3) for the HLA DQA1, LDLR, GYPA, HBG, D7S8, and GC genetic markers are listed below. These allele frequencies were obtained by Roche Molecular Systems' scientists by typing unrelated U.S. Caucasian (n = 200), African American (n = 200), U.S. Hispanic (n = 200), and Japanese (n = 89) individuals.²⁷

The currently available population data¹⁵ suggest there are no significant deviations from Hardy-Weinberg equilibrium expectations for these six markers. Given the statistical independence of alleles at these loci (linkage equilibrium), the frequency of a multi-locus genotype can be estimated by multiplying the genotype frequencies at each locus.

Table 2. AmpliType® Genetic Marker Characteristics

Chromosomal location	HLA DQA1	LDLR	GYPA	HBG	D7S8	GC
6p21.3	19p13.1 - 13.3	4q28 - 31	11p15.5	7q22 - 31.1	4q11 - 13	
PCR Product (bp)	239/242	214	172	151	138	
Number of Alleles	7*	2	2**	3	2	3

*The HLA DQA1 4.2 and 4.3 alleles are detected but not distinguished from each other by the AmpliType HLA DQA1 DNA Probe Strip.

**The GYPA A and A' alleles are detected but not distinguished from each other by the AmpliType PM DNA Probe Strip.

Table 3. Allele Frequencies

Genetic Marker	Allele	U.S. Caucasian (n = 200)	African American (n = 200)	U.S. Hispanic (n = 200)	Japanese (n = 89)
HLA DQA1	1.1	0.158	0.125	0.105	0.084
	1.2	0.190	0.329	0.130	0.118
	1.3	0.073	0.058	0.053	0.236
	2	0.145	0.130	0.115	0.006
LDLR	3	0.192	0.090	0.218	0.444
	4.1	0.214	0.185	0.269	0.073
	4.2/4.3	0.028	0.083	0.110	0.039
GYPA	A	0.448	0.235	0.485	0.202
	B	0.552	0.765	0.515	0.798
HBG	A	0.530	0.527	0.615	0.517
	B	0.470	0.473	0.385	0.483
D7S8	A	0.537	0.439	0.375	0.331
	B	0.450	0.228	0.580	0.669
	C	0.013	0.333	0.045	0.000
GC	A	0.610	0.655	0.622	0.612
	B	0.390	0.345	0.378	0.388
GC	A	0.275	0.090	0.203	0.287
	B	0.178	0.720	0.335	0.471
	C	0.547	0.190	0.462	0.242

The power of discrimination (P_d) values calculated from the combined frequencies for the five AmpliType[®] PM markers and for all six of the markers (including HLA DQA1) are listed in Table 4. The P_d value reflects the probability that two unrelated individuals will have different genotypes at one or more of the six loci.

Table 4. Power of Discrimination for AmpliType[®] Genetic Marker Systems.

Combined power of discrimination (Combined $P_d = 1 - \prod_{i=1}^n (P_{i2} + \dots + P_{in})$; where n is the number of combined markers and P_i is $1 - P_d$ for each number)				
	U.S. Caucasian	African American	U.S. Hispanic	Japanese
AmpliType [®] PM	0.9954	0.9948	0.9961	0.9938
AmpliType [®] PM+DQA1	0.9998	0.9997	0.9998	0.9993

LIST OF KIT COMPONENTS

Caution: Reagents used for amplification must be isolated from any source of contaminating DNA, especially from previously amplified AmpliType PCR products.
Caution: Sodium azide is used as a preservative in several reagents. It can react with lead and copper plumbing to form explosive metal azides. On disposal of the reagent, flush with a large volume of water to prevent azide build-up.

AmpliType PM PCR Amplification and Typing Kit

Reagent	Volume	Description
AmpliType PM PCR Reaction Mix	2.4 mL	1 bottle of PCR Reaction Mix containing AmpliType [®] DNA Polymerase, MgCl ₂ , dATP, dGTP, dCTP, and dTTP, and 0.08% sodium azide in buffer and salt. Store at 2° to 8°C.
AmpliType PM Primer Set	1.2 mL/tube	2 tubes each containing twelve biotinylated primers and 0.05% sodium azide in buffer and salt. Store at 2° to 8°C.
Control DNA 1	0.2 mL	1 tube containing 100 ng/mL human genomic DNA and 0.05% sodium azide in buffer. AmpliType PM type: LDLR BB, GYPA AB, HBGGA, D7S8 AB, GC BB. Store at 2° to 8°C.

AmpliType[®] PM PCR Amplification and Typing Kit (continued)

Reagent (continued)	Volume	Description
Mineral Oil	5 mL	Supplied in dropper bottle. Do NOT expose to strong ultraviolet light. Store at 2° to 30°C.
AmpliType [®] PM DNA Probe Strips	50 strips	The strips are provided in a screw-top tube with a packet of desiccant. Store the strips at 2° to 8°C in the tightly capped tube with the desiccant and protected from light.
Enzyme Conjugate: HRP-SA	2.0 mL	1 bottle containing Horseradish Peroxidase-Streptavidin (HRP-SA) Enzyme Conjugate supplied in buffer with preservative. Store at 2° to 8°C.
Chromogen: TMB	60 mg	1 bottle containing powdered 3,3',5,5'-tetramethylbenzidine (TMB). Dissolve Chromogen as instructed before use (see Preparation of Reagents Supplied, Section 1.1). Store at 2° to 8°C.
Package Insert		Containing the AmpliType PM PCR amplification and typing protocols.
AmpliType PM+DQA1 PCR Amplification and Typing Kit		
Reagent	Volume	Description
AmpliType PM PCR Reaction Mix	2.4 mL	1 bottle of PCR Reaction Mix containing AmpliType [®] DNA Polymerase, MgCl ₂ , dATP, dGTP, dCTP, dTTP, and 0.08% sodium azide in buffer and salt. Store at 2° to 8°C.
AmpliType PM Primer Set	1.2 mL/tube	2 tubes each containing twelve biotinylated primers and 0.05% sodium azide in buffer and salt. Store at 2° to 8°C.

AmpliType® PM+DQA1 PCR Amplification and Typing Kit (continued)

Reagent (continued)	Volume	Description
Control DNA 1	0.2 mL	1 tube containing 100 ng/mL human genomic DNA and 0.05% sodium azide in buffer. AmpliType® PM type: LDLR BB, GYPA AB, HBGGA, D7S8 AB, GC BB. AmpliType® HLA DQA1 type: 1.1, 4.1. Store at 2° to 8°C.
Mineral Oil	5 mL	Supplied in dropper bottle. Do NOT expose to strong ultraviolet light. Store at 2° to 30°C.
AmpliType PM DNA Probe Strips	50 strips	The strips are provided in a screw-top tube with a packet of desiccant. Store the strips at 2° to 8°C in the tightly capped tube with the desiccant and protected from light.
AmpliType HLA DQA1 DNA Probe Strips	50 strips	The strips are provided in a screw-top tube with a packet of desiccant. Store the strips at 2° to 8°C in the tightly capped tube with the desiccant and protected from light.
Enzyme Conjugate: HRP-SA	2.0 mL/bottle	2 bottles each containing Horseradish Peroxidase-Streptavidin (HRP-SA) Enzyme Conjugate supplied in buffer with preservative. Store at 2° to 8°C.
Chromogen: TMB	60 mg/bottle	2 bottles each containing powdered 3,3',5,5'-tetramethylbenzidine (TMB). Dissolve Chromogen as instructed before use (see Preparation of Reagents Supplied, Section 1.1). Store at 2° to 8°C.
Package Insert		Containing the AmpliType PM+DQA1 PCR amplification and typing protocols.

MATERIALS REQUIRED BUT NOT SUPPLIED

The items listed are those required for PCR amplification, PCR product gel analysis, hybridization, and color development. This list does not include reagents or equipment required for DNA extraction. Many of the items listed below are available from major laboratory suppliers (MLS) such as Baxter Scientific Products (McGaw Park, IL), Fisher Scientific (Pittsburgh, PA), or VWR (So. Plainfield, NJ) unless otherwise noted. Equivalent sources may be acceptable where noted.

Reagents	Source
Use Reagent Grade unless otherwise noted.	
123 Base Pair Ladder	GIBCO BRL, Gaithersburg, MD (Cat. No. 15613-011).
Agarose: buy the following and mix according to protocol (Section 6.2): NuSieve® GTG® Agarose	FMC, Rockland, ME (Cat. No. 50082).
Seakem® GTG Agarose	FMC (Cat. No. 50072).
Alcohol	MLS
95% ethanol or 70% Isopropanol	
Boric acid	GIBCO BRL (Cat. No. 15583-016), Sigma, St. Louis, MO (Cat. No. B 7901), or equivalent.
Bromophenol blue	Sigma (Cat. No. B 5525) or equivalent.
Citric acid, monohydrate, granular	JT Baker, Phillipsburg, NJ (Cat. No. 01110-05) or equivalent.
Ethanol, 200 proof (100%) in glass container	Quantum Chemical Company, Cincinnati, OH (order by volume) or Gold Shield, Hayward, CA (Cat. No. EL200P24X1).
Ethidium bromide, 10 mg/mL	MLS
Ethylenediaminetetraacetic acid (EDTA), disodium salt, dihydrate	GIBCO BRL (Cat. No. 15575-012), Sigma (Cat. No. E 4884), or equivalent.
Glycerol	Sigma (Cat. No. G 5516) or equivalent.
Hydrochloric acid, concentrated	MLS
Hydrogen peroxide, 30%	Sigma (Cat. No. H 1009) or equivalent.
Sodium chloride	MLS
Sodium dodecyl sulfate (SDS), ultra pure electrophoresis grade	GIBCO BRL (Cat. No. 15525-025), Bio-Rad, Hercules, CA (Cat. No. 161-0302), or equivalent.
Sodium hydroxide (NaOH), 10 N	VWR, South Plainfield, NJ (JT5674-3) or equivalent.
Sodium hydroxide, pellets	MLS
Sodium phosphate, monobasic, monohydrate	Sigma (Cat. No. S 9638) or equivalent.
Trisodium citrate, dihydrate	Sigma (Cat. No. C 8532) or equivalent.
Tris base	GIBCO BRL (Cat. No. 15504-012), Sigma (Cat. No. T 8404), or equivalent.

Equipment	Source
AmpliT [™] DNA Typing Tray	Perkin-Elmer (Part No. N808-0065).
Aspirator Apparatus (See Reference 18 for further description)	Perkin-Elmer (Part No. N830-2697) or SERI, Richmond, CA (Cat. No. E081).
Balance, accurate to 1.0 mg	MLS
Dispensing Pipette [†] (5 or 10 mL)	Baxter Scientific Products, McGaw Park, IL (Cat. No. P4965-5 or P4965-10) or equivalent.
Electrophoresis Power Supply	Bio-Rad (Cat. No. 165-4761 (110V) or Cat. No. 165-4762 (220V)) or equivalent.
Filters, cellulose nitrate membrane, 0.2 and 0.45 µm Nalgene [®]	Fisher, Pittsburgh, PA (0.2 µm filter Cat. No. 127-0020 or 0.45 µm filter Cat. No. 127-0045) or equivalent.
Filter forceps, non-pointed tips	Fisher (Cat. No. 09-753-50).
Filter, orange, Wtatten 22 (for use with Polaroid [®] camera)	Kodak, New Haven, CT (Cat. No. 149 5571).
Freezer (non frost-free, -20°C)	MLS
Horizontal minigel electrophoresis apparatus with Gel Combs (Teflon, 1.2 mm thick)	Aquebogue Machines, Aquebogue, NY (model 750 Mini Horizontal gel apparatus)
Laboratory glassware	MLS
Laminar flow/biological/biosafety cabinet or designated clean area	Labconco Corp., Kansas City, MO or MLS.
Microcentrifuge	MLS
Microplate, 96-wells	VWR, So. Plainfield, NJ (Cat. No. 62408-220) or equivalent.
Micro Tube Opener	Robbins Scientific Corp., Sunnyvale, CA (Cat. No. 1012-20-0) or equivalent.
Microtube racks, 1.5 mL [†]	MLS
Microwave	MLS
pH meter compatible with Tris solutions (need electrodes without silver chloride) and reference buffers	Corning or equivalent.
Pipettors	Rainin, Woburn, MA (P20, P100, P200, adjustable to deliver 2 - 20 µL [†] , 10 - 100 µL [†] , 50 - 200 µL [†] , and 100 - 1000 µL [†]) or equivalent.
Pipettor, positive displacement repeat pipettor	Fisher (Cat. No. 21-380-8) or equivalent.
Pipette bulbs [†]	MLS
Polaroid camera and film	Fotodyne, Hartland, WI, Sigma, or any leather black and white film (type 55, 57, or 667) or color film (type 59 or 559)
Refrigerator (2° - 8°C) [†]	MLS

[†] One is needed for pre-PCR sample preparation; a second is needed for post-PCR sample handling.

Equipment (continued)	Source
Shaker, variable speed, orbital platform (needed for color development step)	Belco, Vineland, NJ (Model No. 7744-00110 (110V) or 7744-00220 (220V)) or Lab-Line, Melrose Park, IL (Model No. 3520 for RPM X 100).
Stir plate and stir bar	MLS
Thermometer, Total Immersion	Recommend: Cole Parmer, Niles, IL (Cat. No. H-08001-34, Graduated to 0.1°C range 49° to 57°C, specify NIST traceable when ordering). Alternative: Scientific Products (Cat. No. T2099-5, range 50° to 80°C).
Timer, 60 minute (± 1 minute) [†]	MLS
UV transilluminator (and UV eye protection)	MLS
Vacuum pump or Vacuum source	Fisher (Cat. No. 13-875-220) or equivalent.
Vortexer	MLS
Water bath with cover, rotating	Belco (Model No. 7746-22110 (110V) or 7746-22220 (220V)).
Water bath, stationary or incubator (37° to 55°C)	MLS
Weight, approximately 1 kg (e.g., lead ring)	VWR (Cat. No. 29700-048), Cole Parmer (Cat. No. G-06137-06), or equivalent.
PCR Amplification Equipment	Source
Choose one of the following GeneAmp [®] PCR Instrument Systems:	
DNA Thermal Cycler 480 and equipment	Perkin-Elmer (Part No. N801-0100 (120V/60Hz), N801-0101 (230V-240V/50Hz), or N801-0102 (100V/50Hz-60Hz)).
DNA Thermal Cycler 480	Perkin-Elmer (Part No. N801-0537).
GeneAmp Thin-Walled Reaction Tubes	Perkin-Elmer (Part No. N801-0611).
GeneAmp Autoclaved Thin-Walled Reaction Tubes	Perkin-Elmer (Part No. N801-0737), with Flat Cap
GeneAmp Thin-Walled Reaction Tubes	Perkin-Elmer (Part No. N801-0434).
Temperature Verification System (for the DNA Thermal Cycler 480)	Perkin-Elmer (Part No. N801-0001 OR
GeneAmp PCR System 9600 and equipment	Perkin-Elmer (Part No. N801-0002 (230V-240V/50Hz), or N801-0003 (100V/50Hz-60Hz)).
GeneAmp PCR System 9600	Perkin-Elmer (Part No. N801-0001 (120V/60Hz), N801-0002 (230V-240V/50Hz), or N801-0003 (100V/50Hz-60Hz)).

[†] One is needed for pre-PCR sample preparation; a second is needed for post-PCR sample handling.

PCR Amplification Equipment (continued) Source

MicroAmp® Autoclaved Reaction Tube with Cap	Perkin-Elmer (Part No. N801-0612).
MicroAmp Reaction Tube with Cap	Perkin-Elmer (Part No. N801-0540).
MicroAmp Base	Perkin-Elmer (Part No. N801-0531).
MicroAmp Tray/Retainer Set	Perkin-Elmer (Part No. N801-0530).
Temperature Verification System (for the GeneAmp® PCR System 9600)	Perkin-Elmer (Part No. N801-0435).

Supplies

	Source
Aluminum foil	MLS
Calculator†	MLS
Chromatography paper, Whatman® grade 3MM, pure cellulose chromatography and blotting paper deionized, glass-distilled, or ultrafiltered water, referred to as DI H ₂ O in this package insert	WVR (Cat. No. 21427) or equivalent.
Gloves,† disposable	MLS
Lab bench surface protector	MLS
Lab coat†	MLS
Lab wipe†	MLS
Parafilm®	MLS
Pen (for labeling probe strips)	Provided with AmpliType™ DNA Typing Trays.
Permanent Ink Marker† (for marking tubes)	Sanford (Sharpie® ultra fine point) or equivalent.
Pipett, disposable serological (5 - 10 mL)	MLS
Pipet tips† for adjustable pipettors with hydrophobic filter plugged tips (sterile)	GIBCO BRL or equivalent.
Pipet CombiTips® (0.5 mL, sterile)	Fisher (Cat. No. 21-380-8A) or equivalent.

Protective eyewear†	MLS
Protective mask†	MLS
Scissors†	MLS

† One is needed for pre-PCR sample preparation; a second is needed for post-PCR sample handling.

1.0 REAGENT PREPARATION

1.1 PREPARATION OF REAGENTS SUPPLIED

1.1.1 Chromogen: TMB Solution.

Warning: Chromogen: TMB is an irritant. Avoid skin contact and inhalation. Wear lab coat, gloves, mask, and protective eyewear when handling Chromogen: TMB.

Bring the bottle of Chromogen: TMB (TMB) to room temperature (15° to 30°C). Before opening the bottle, tap it on the lab bench to shake the TMB to the bottom

of the bottle. Remove the stopper carefully to prevent loss of the powder. Slowly add 30 mL of room temperature reagent grade 100% ethanol to the bottle. Do NOT use ethanol that has been stored in a metal container. Do NOT use 95% ethanol or other alcohols. Recap the bottle. Seal the stopper with Parafilm®. Secure the bottle to an orbital shaker and shake in an upright position for 2 hours or until completely dissolved. Protect from heavy metal contamination, especially iron oxide (rust). Store bottle at 2° to 8°C. Under these conditions the Chromogen Solution is stable for six months after preparation.

1.1.2 AmpliType® PM PCR Reaction Mix.

This step must be performed either in a biological hood or in a room free from amplified DNA.

All of the AmpliType® PM PCR Reaction Mix should be aliquotted at the same time to decrease the possibility of contamination through repeated opening of the AmpliType PM PCR Reaction Mix bottle.

Upon first use of the AmpliType PM Kit or AmpliType PM+DQA1 Kit, remove the bottle of AmpliType PM PCR Reaction Mix and carefully aliquot 40 µL into autoclaved tubes (GeneAmp® Thin-Walled Reaction Tubes for the DNA Thermal Cycler 480 or MicroAmp® Reaction Tube with Cap if using the GeneAmp PCR System 9600) using a dedicated positive displacement pipettor or a pipettor with hydrophobic filter plugged tips. Ensure that the tubes are capped tightly. Place PCR Reaction Mix tubes in a rack not used for DNA preparation or amplified DNA handling. Store tubes separated from any source of DNA at 2° to 8°C.

1.2 PREPARATION OF REAGENTS NOT SUPPLIED

Use reagent grade chemicals unless otherwise noted. Prepare all solutions using deionized, glass-distilled, or ultrafiltered water (identified below as DI H₂O). Wear gloves and follow safety recommendations provided by the manufacturer for handling chemicals. Comply with any and all laws, regulations, or orders with respect to the disposal of any hazardous or toxic chemical, material, substance, or waste. Store all solutions at room temperature (15° to 30°C) unless otherwise noted. Reagents are prepared as follows:

1.2.1 Citrate buffer (0.1 M sodium citrate, pH 5.0 (1 L)).

Dissolve 18.4 g of trisodium citrate, dihydrate (Na₃C₆H₅O₇•2H₂O) in 800 mL of DI H₂O. Adjust the pH to 5.0 (± 0.2) by addition of approximately 6 g of citric acid, monohydrate (C₆H₈O₇•H₂O). Adjust to a final volume of 1 L using DI H₂O and mix thoroughly. Autoclave the solution.

1.2.2 20X SSPE buffer (3.6 M NaCl, 200 mM NaH₂PO₄•H₂O, 20 mM EDTA, pH 7.4 (1 L)).

Warning: EDTA is an irritant. Wear lab coat, gloves, mask, and protective eyewear when handling EDTA.

Warning: NaOH is corrosive and toxic. Wear lab coat, gloves, and protective eyewear. Use caution when handling NaOH.

Dissolve 7.4 g of disodium ethylenediaminetetraacetic acid, dihydrate (Na₂EDTA•2H₂O) in 800 mL of DI H₂O. Adjust the pH to 6.0 (± 0.2) with 10 N sodium

hydroxide (NaOH). Add 210 g of sodium chloride (NaCl) and 27.6 g sodium phosphate, monobasic, monohydrate ($\text{NaH}_2\text{PO}_4 \cdot \text{H}_2\text{O}$). Adjust the pH to 7.4 (± 0.2) with 10 N NaOH (approximately 10 mL). Adjust to a final volume of 1 L using DI H_2O and mix thoroughly.

1.2.3 20% (w/v) Sodium Dodecyl Sulfate SDS (1 L).

Warning: SDS is an irritant. Avoid skin contact and inhalation. Wear lab coat, gloves, mask, and protective eyewear when handling SDS.

Slowly dissolve 200 g of electrophoresis-grade (ultra pure) SDS in 800 mL of DI H_2O . Warning (e.g., in a 37°C water bath) may be required to dissolve solids completely. Adjust to a final volume of 1 L using DI H_2O and mix thoroughly.

1.2.4 Hybridization Solution [5X SSPE, 0.5% w/v SDS (1 L)].

Add 250 mL of 20X SSPE (Section 1.2.2) and 25 mL of 20% w/v SDS (Section 1.2.3) to 725 mL of DI H_2O and mix thoroughly. Hybridization Solution solids must be in solution and the solution must be well mixed before use. Warning (e.g., to between 37° to 55°C in a stationary water bath or an incubator) may be required to dissolve solids completely.

1.2.5 Wash Solution [2.5X SSPE, 0.1% w/v SDS (2 L)].

Add 250 mL of 20X SSPE (Section 1.2.2) and 10 mL 20% w/v SDS (Section 1.2.3) to 1,740 mL of DI H_2O and mix thoroughly. Wash Solution solids must be in solution and the solution must be well mixed before use; warning (e.g., to between 37° to 55°C in a stationary water bath or an incubator) may be required to dissolve solids completely.

1.2.6 3% Hydrogen peroxide (1 mL).

Add 100 μL of 30% H_2O_2 to 900 μL of DI H_2O and vortex to mix. Protect from light. Store at 2° to 8°C. In practice, 3% hydrogen peroxide has a shelf life of approximately 6 months. Alternatively, one-tenth the volume of 30% H_2O_2 can be used directly in the Color Development Step (see Section 8).

1.2.7 1 M Tris-HCl, pH 8.0 (1 L).

Warning: Hydrochloric acid (HCl) causes severe burns and is irritating to the eyes. When preparing this reagent, use a fume hood and avoid inhalation and contact with the skin. Wear a lab coat, gloves, and protective eyewear when handling.

Dissolve 121.1 g of Tris base in 800 mL of DI H_2O . Adjust to pH 8.0 (± 0.2) at room temperature by adding approximately 45 mL of concentrated HCl. Adjust the final volume to 1 L with DI H_2O and mix thoroughly. The solution should be autoclaved or filtered through a 0.2 μm Nalgene® filter.

1.2.8 0.5 M EDTA (1 L).

Warning: EDTA is an irritant. Wear lab coat, gloves, mask, and protective eyewear when handling EDTA.

Warning: NaOH is corrosive and toxic. Wear lab coat, gloves, and protective eyewear. Use caution when handling. NaOH, when combined with water, results in an exothermic reaction. Dissolve the NaOH pellets in the water

gradually with 3 or 4 additions. Cover and allow the NaOH to dissolve completely and cool between each addition.

Add 186.1 g of disodium ethylenediaminetetraacetic acid dihydrate ($\text{Na}_2\text{EDTA} \cdot 2\text{H}_2\text{O}$) to 800 mL of DI H_2O . Stir vigorously on a magnetic stirrer. To dissolve the EDTA powder, adjust the pH to 8.0 (± 0.2) by adding approximately 20 g of NaOH pellets. Check the pH and add 10 N NaOH solution if further pH adjustment is needed (for small pH adjustments, pellet addition is not advised). Bring volume up to 1 L with DI H_2O and mix thoroughly. The solution should be autoclaved or filtered through a 0.2 μm Nalgene® filter.

1.2.9 200 mM EDTA (10 mL).

Add 4.0 mL of 0.5 M disodium EDTA (Section 1.2.8) to 6.0 mL of DI H_2O and mix thoroughly. The solution should be autoclaved or filtered through a 0.2 μm Nalgene filter. The solution may be aliquotted, 300 μL /tube, into microcentrifuge tubes. Store at 2° to 8°C.

1.2.10 10X TBE [890 mM Tris, 890 mM boric acid, 20 mM EDTA (1 L)].

To 40 mL of 0.5 M EDTA, pH 8.0 (Section 1.2.8), add approximately 900 mL of DI H_2O . Add 108 g of Tris base and 55 g of boric acid to the diluted EDTA solution. Stir vigorously on a magnetic stirrer. Adjust the volume to 1 L with DI H_2O and mix thoroughly. Filter the mixture using a 0.2 or 0.45 μm Nalgene filter unit to remove particulate matter and prevent the formation of a precipitate. Storage in a glass container is recommended to facilitate visual inspection for precipitates. If a white precipitate is noted in the 10X TBE, the buffer should be discarded and remade.

1.2.11 Running Buffer, 0.5X TBE [44.5 mM Tris-borate, 1 mM EDTA (1 L)].

Add 50 mL of 10X TBE (Section 1.2.10) to 950 mL of DI H_2O and mix thoroughly.

1.2.12 Gel Loading Buffer [0.2% bromophenol blue, 50% glycerol, 20 mM Tris (6H₉O₄), 2.5 mM EDTA (50.0 mL)].

Dissolve 0.1 g of bromophenol blue in 23.75 mL of DI H_2O . Allow 2 hours for bromophenol blue to go into solution; stir as necessary. Add 1 mL of 1 M Tris-HCl, pH 6.0 (Section 1.2.7), 0.25 mL of 0.5 M disodium EDTA (Section 1.2.8), and 25 mL of glycerol. Mix thoroughly. Store at 2° to 8°C.

Note: Sections 2.0 and 3.0 should be read before preparing AmpliType® PM PCR amplification reactions.

2.0 PRECAUTIONS

The sensitivity of PCR allows minute quantities of DNA to be typed using the AmpliType PCR Amplification and Typing Kits. Contamination of the samples by handling or by exposure to any other source of human DNA is an important concern. Precautions should be taken to prevent the following three types of contamination: 1) "carryover" of PCR product from one amplification to the next; 2) cross-contamination between extracted DNA samples; 3) contamination with exogenous human genomic DNA.¹⁸²⁸ At a minimum, the pre-PCR amplification area must be separated from the post-PCR amplification area.

3.0 PCR AMPLIFICATION PROTOCOLS

The following protocols detail the PCR amplification and typing procedures specific for the AmpliType® PCR Amplification and Typing Kits. Protocols and precautions for DNA extraction are described elsewhere.¹⁸

Note: Where the use of DI H₂O is indicated below, either deionized, glass-distilled, or ultrafiltered water may be used.

A dedicated area, such as a biological hood or a separate room, should be used for preparing AmpliType PCR amplification reactions whenever possible. All equipment and supplies used to prepare amplification reactions should be kept in this dedicated "clean" area at all times. Do not use these items to handle amplified DNA. To prevent contamination from other potential sources of genomic DNA, pipet tips plugged with hydrophobic filters must be used. Trace amounts of amplified DNA, if carried over into other samples before amplification, can lead to results that are not interpretable. Do NOT bring amplified DNA or equipment and supplies used to handle amplified DNA into the designated clean area.

Note: Wear clean disposable laboratory gloves while preparing samples for PCR amplification. Change gloves frequently or whenever there is a chance they have been contaminated with DNA.

3.1 Prior to preparing the DNA test samples and PCR amplification reactions, turn on the GeneAmp® PCR Instrument System and program the appropriate PCR profile times and temperatures provided in Table 5. Refer to the appropriate GeneAmp PCR Instrument System Users Manual for additional information on programming the instruments. If using a stored file, step through the program to verify the parameters.

If using a GeneAmp PCR System 9600, follow Sections 3.1.1 and 3.1.2 before proceeding to Section 3.2. If using the DNA Thermal Cycler 480, proceed to Section 3.2.

3.1.1 AmpliType PM and AmpliType PM+DQA1 amplification, using the GeneAmp PCR System 9600, requires that the block be preheated to 95°C. To preheat the block of the GeneAmp PCR System 9600 instrument, use the following CYCL program flowchart (Figure 1, page 18) to bring the block to 95°C. The settings in the Default Mode column will appear automatically when stepping through the CYCL program. The default settings that must be changed to preheat the block are indicated in **bold-type**. The new settings are indicated in **bold-type** in the Program Setting column. Link this preheat program to the beginning of a separate program for 32 cycles using the parameters stated in Table 5 for the GeneAmp PCR System 9600.

3.1.2 Run the program that was created to preheat the block. The block must be at 95°C before the tubes are placed in the block. It is not necessary to close the heated cover during this step.

3.2 Prepare the DNA test samples for addition to the PCR amplification reactions. It is recommended that the DNA sample be quantitated with the QuantIBlot Human DNA Quantitation Kit (Part No. N808-0114). The final DNA concentration should be in the range of 0.1 to 0.5 ng/μL so that 2 to 10 ng of DNA will be added to the PCR reaction in a volume of 20 μL. If the sample DNA concentration is greater than 0.5 ng/μL, dilute a portion of it with autoclaved DI H₂O. For most samples, 2 to 5 ng is sufficient. If the DNA sample contains degraded DNA, it may be appropriate to add >10 ng of DNA.

Table 5. Perkin-Elmer GeneAmp® PCR Instrument Systems PCR Profile Times and Temperatures

GeneAmp® PCR Instrument System	Tube	Total Reaction Volume in μL Tube	Tube Type	Time and Temperature			
				AmpliType® Kit	Initial Step	Each of 32 Cycles	
					Denature	Anneal	Extend
DNA Thermal Cycler 480	GeneAmp Thin-Walled Reaction Tubes plus 2 drops Mineral Oil	100 μL	PM	None	60 sec	30 sec	30 sec
					94°C	60°C	72°C
GeneAmp PCR System 9600	MicroAmp® Reaction Tube with Cap	100 μL	PM	Preheat Block* (see 3.1.1)	30 sec	30 sec	10 min
					95°C	63°C	72°C

Note: These programs may be stored as User Files for later use (see Users Manual DNA Thermal Cycler 480 or GeneAmp PCR System 9600 Users Manual).

As with any calibrated piece of laboratory equipment, the GeneAmp PCR Instrument System should have a documented temperature verification test performed on a regular basis (e.g., monthly using the Temperature Verification System (Part No. N808-0434) for the DNA Thermal Cycler 480 or the GeneAmp PCR System 9600 Temperature Verification System (Part No. N808-0435) for the GeneAmp PCR System 9600).

*Preheat Block: The sample block is heated to 95°C prior to placing the tubes into the instrument and the start of amplification cycles (see Section 3.1.1).

Default Mode	Change to	Program Setting
3 Temperature PCR Setpt #1 94.0°C	↓	1 Temperature PCR Setpt #1 95.0°C
Ramp 0:00 Hold 0:30	↓	Ramp 0:00 Hold 0:00
Total cycles Pause during run?*	↓	Total cycles Pause during run?*
= 1 NO	↓	= 1 YES
Pause after setpt? Beep while paused?	↓	Pause after setpt? Beep while paused?
1 YES	↓	1 NO
1st Pause at cyc Pause every 1 cyc	↓	1st Pause at cyc Pause every 1 cyc
1 1	↓	1 1
Pause time ↓	↓	Pause time ↓
10:00	↓	9:59
STORE	↓	STORE

Figure 1. CYCL program for preheating the GeneAmp® PCR System 9600 block.

3.3 Determine the number of samples to be amplified, including positive and negative controls. The Control DNA 1 provided in the kits should be amplified each time the kits are used and run on each PCR product gel (Sections 5.0 and 6.0) as a positive control. A negative control consisting of 20 µL of autoclaved DI H₂O in place of the DNA sample should also be included with each set of amplification reactions.

3.4 Transfer the PCR amplification reagents to the designated clean area. Place the required number of reaction tubes containing 40 µL of aliquotted AmpliType® PM PCR Reaction Mix in a rack or MicroAmp® Base not used for the preparation of DNA or the handling of amplified DNA. Use the MicroAmp Base when the AmpliType PM PCR Reaction Mix has been aliquotted into MicroAmp Reaction Tubes. Label the reaction tubes. Do not label the top of the MicroAmp Reaction Tubes because the heated cover may remove the ink.

3.5 Ensure that the solution is at the bottom of each tube. Open the caps with a clean Micro Tube Opener (this de-capping tool should not have been used on tubes containing amplified or extracted DNA). Avoid touching the inside surface of the tube caps.

3.6 Pipet 40 µL of the AmpliType® PM Primer Set into each tube, including control tubes, using a pipettor with hydrophobic filter plugged tips. Pipet carefully against the side of the tube to minimize mixing and to avoid splashing the solution. It is important to begin the cycling process within 20 minutes after the addition of the AmpliType PM Primer Set to minimize the formation of primer-dimers and other non-specific PCR products.

3.7 If using the GeneAmp® PCR System 9600, cap each tube loosely and proceed to Section 3.8. If using the DNA Thermal Cycler 480, carefully add 2 drops of Mineral Oil from the dropper bottle provided in the kits to all tubes, including the control tubes, before proceeding to Section 3.8. Be careful NOT to touch the reaction tubes with the dropper. Cap each tube loosely. Do NOT vortex, mix, or spin.

Note: Each AmpliType PM and AmpliType PM+DOA1 PCR amplification is performed in a final volume of 100 µL. The volume allocated for sample addition is 20 µL.

3.8 For each of the following additions, complete the processing of each tube before proceeding to the next tube. To prevent contamination, no more than one tube should be open at a time. Use a new, autoclaved pipet tip for each addition. Open the tube and carefully add 20 µL of extracted sample, Control DNA 1, or autoclaved DI H₂O (negative control). If using the DNA Thermal Cycler 480, carefully insert the pipet tip through the mineral oil layer. Discard the pipet tip and re-cap the tube before proceeding to the next sample. Do NOT vortex, mix, or spin. Prepare tubes as follows:

3.8.1 DNA Test Sample tubes. Add 20 µL of sample DNA prepared in Section 3.2 to each labelled DNA test sample tube.

3.8.2 Positive Control tube. Vortex the Control DNA 1 tube and spin briefly in a microcentrifuge before each use to remove any liquid from the cap. Add 20 µL of the 100 ng/mL (2 ng) Control DNA 1 to the labelled positive control tube.

3.8.3 Negative Control tube. Add 20 µL of autoclaved DI H₂O to the labelled negative control tube.

3.9 As soon as all samples have been added, place the tubes into the appropriate GeneAmp PCR Instrument System. If using MicroAmp® Reaction Tubes with Caps, place a MicroAmp Tray (not a MicroAmp Base) in the GeneAmp PCR System 9600 block before inserting the tubes. The position of each tube in the block should be recorded. Use caution when handling tubes since the GeneAmp PCR System 9600 block has been preheated to 95°C.

3.9.1. If using the GeneAmp PCR System 9600, proceed to Section 3.9.2. If using the DNA Thermal Cycler 480, push the tubes down completely into the block and start the 32 cycle amplification file. Verify the cycling parameters by monitoring the first cycle. Check the tubes after the first cycle to ensure they are all still sealed tightly in the block. Proceed to Section 4.0.

3.9.2 Position and tighten the GeneAmp PCR System 9600 heated cover over the tubes. Immediately press the "run" key to end the preheat program and start the 32 cycle amplification program. Verify the cycling parameters by monitoring the first cycle. Proceed to Section 4.0.

4.0 PREPARATION OF PCR PRODUCT FOR DETECTION

After the PCR amplification process, remove the tubes from the GeneAmp® PCR Instrument System. Prior to DNA hybridization (Section 8.0), open the tubes one at a time and add 5 µL of 200 mM disodium EDTA (Section 1.2.9). Use a new pipet tip for each addition. If using the DNA Thermal Cycler 480, carefully insert the pipet tip through the mineral oil layer. Discard the pipet tip and re-cap the tube before proceeding to the next tube. It may be convenient to remove a 5 µL aliquot for gel electrophoresis (see Section 7.0) before recapping the tube.

Note: EDTA must be added prior to heat denaturation of the samples but may be added before or after the aliquot for gel electrophoresis has been removed. Samples are now ready for analysis by gel electrophoresis, DNA hybridization, and color development. Amplified samples containing 9.5 mM EDTA may be stored at 2° to 8°C for 2 months or at -20°C for 6 months. The continued acceptable performance of these samples beyond these periods may vary with the sample. **Store amplified DNA samples separate from all PCR amplification reagents, extracted DNA samples, and casework samples.**

5.0 VERIFICATION OF PCR AMPLIFICATION

It is recommended that the presence and size of the PCR products generated in the AmpliType® PM reactions be verified by agarose gel electrophoresis. This step should be performed prior to denaturing the samples for DNA hybridization to ensure sharp product bands on the gel. Alternatively, a portion of the amplified sample can be set aside before proceeding to the typing step. **Do NOT heat denature PCR products prior to gel electrophoresis; additional bands may appear on the gel and the product bands may smear.**

Optimal gel electrophoresis results have been obtained using a 5.5 x 9 x 0.45 cm gel in a minigel electrophoresis apparatus (e.g., Aquebogue model 750 Mini Horizontal Gel Apparatus) with the recommended comb. Other gel sizes and apparatus can give acceptable results if the recommended protocols in Sections 6.0 and 7.0 are followed.

6.0 PREPARATION AND POURING OF AGAROSE GEL

6.1 Set up a minigel apparatus containing the casting tray on a level surface. If the gel casting tray is separate from the apparatus, tape the ends of the tray to retain the liquefied agarose.

6.2 Prepare a 3% NuSieve®/1% SeaKem® agarose gel solution. Add 3.0 g NuSieve GTG® agarose and 1.0 g SeaKem GTG agarose to 100 mL of 0.5X TBE Gel Running Buffer (Section 1.2.10) in a 200 mL heat-tolerant laboratory screw-cap bottle.

6.3 Prior to melting the agarose, weigh the bottle containing the agarose and buffer and record the weight.

6.4 Cover the bottle loosely with the cap during heating. To melt the agarose, heat it in a microwave oven (set at defrost or low) for 3 to 4 minutes or in a boiling water bath until melted. Do not allow the agarose solution to boil over. During and after heating, carefully swirl the bottle to aid in dissolving the agarose. Additional heating may be required with some microwave ovens.

6.5 After the agarose is completely dissolved, reweigh the bottle. If the bottle weighs less (volume lost through evaporation), then add a sufficient volume of DI H₂O to return the bottle to its original weight. Swirl the solution to ensure that the solution is homogeneous. If necessary, reheat the bottle (10 to 30 seconds) to ensure that the agarose remains completely dissolved.

Warning: Ethidium bromide is a mutagen. Avoid contact with skin. Wear lab coat, gloves, and safety glasses when handling ethidium bromide.

6.6 After the agarose is completely dissolved, add 5 µL of a 10 mg/mL stock of ethidium bromide to the 100 mL of gel solution to yield a 0.5 µg/mL ethidium bromide concentration. Swirl to mix thoroughly.

6.7 Cool the solution to 55° to 65°C on the bench top or in a water bath prior to pouring the gel. If the gel will not be poured immediately after cooling, maintain the temperature of the agarose at 55° to 65°C to prevent the gel from solidifying.

6.8 It is important that the gel is poured on a level, horizontal surface. Pour enough agarose for a 0.45 cm thick gel into either a gel casting tray or the casting area of the minigel apparatus (e.g., to form a gel 5.5 x 9 x 0.45 cm, 22 mL of liquefied agarose are needed). Avoid formation of air bubbles when pouring agarose. Immediately insert a gel comb. The teeth of the comb used should be 1.0 to 1.2 mm thick and 4.0 to 5.0 mm wide. The comb should sit approximately 1.0 mm above the bottom of the gel.

Note: The remaining agarose solution can be reheated and reused. However, it is important to weigh the bottle containing the gel both before and after heating, and to correct the difference in weight by adding DI H₂O. Additional ethidium bromide is not needed. Do not allow the agarose containing ethidium bromide to boil over; ethidium bromide is a mutagen. If a microwave oven is used to reheat the agarose, it must be dedicated for laboratory use only.

6.9 Allow the gel to solidify completely (approximately 30 minutes) at room temperature (15° to 30°C) before proceeding.

6.10 After the gel has solidified, carefully remove the comb. If the gel was poured in a separate casting tray, place the gel in the horizontal minigel apparatus on a level surface. Add a sufficient volume of 0.5X TBE Gel Running Buffer (Section 1.2.10) containing 0.5 µg/mL ethidium bromide (add 50 µL of 10 mg/mL ethidium bromide to 1L of 0.5X TBE Gel Running Buffer) to fill the buffer tanks and cover the gel to a depth of about 2 mm. Check that there are no air bubbles in any of the wells.

7.0 GEL LOADING AND ELECTROPHORESIS

7.1 Add 2 µL of Gel Loading Buffer (Section 1.2.11) either to 0.5 mL microcentrifuge tubes or to the wells of a 96-well microplate corresponding to the number of samples to be analyzed. The GIBCO BRL 123 Base Pair (bp) Ladder will also be run in one well of each gel. If a 5 µL aliquot of the amplified DNA sample was removed in Section 4.0, add the 2 µL of Gel Loading Buffer directly to the aliquot.

7.2 Dilute a portion of the stock solution of GIBCO BRL 123 bp Ladder to 50 ng/µL with autoclaved DI H₂O.

7.3 Add 5 µL of the 50 ng/µL GIBCO BRL 123 bp Ladder (molecular weight marker) to the 2 µL of loading buffer designated for the molecular weight marker and mix using the pipet tip.

7.4 If 5 µL of each amplified DNA sample has already been mixed with 2 µL of loading buffer (see Section 7.1), proceed to Section 7.5. If the aliquots have not been previously removed, withdraw 5 µL of the amplified DNA sample from each reaction tube. Avoid pipetting mineral oil along with the sample if amplification was performed in the DNA Thermal Cycler 480 with GeneAmp® Thin-Walled Reaction Tubes. Use a new pipet tip to add each sample to 2 µL of loading buffer. Mix the sample and loading buffer with the pipet tip.

7.5 After all amplified samples have been mixed with the loading buffer, carefully pipet the 123 bp ladder plus loading buffer solution into the first gel well using an adjustable pipet. Add the amplified DNA sample plus loading buffer solutions into the remaining gel wells using a new pipet tip for each sample. Do NOT produce bubbles or overflow wells. Samples should lie in an even layer at the bottom of each well.

7.6 Connect the minigel electrophoresis apparatus to the power supply such that the DNA migrates toward the positive electrode. For a minigel apparatus that is 5.5 x 15 cm, run the gel at 115 volts or 7.5 volts/cm (maximum voltage) for approximately 1 hour or until the bromophenol blue dye from the loading buffer has run approximately 7.5 cm down the gel to provide adequate resolution of the six amplified product bands. Electrophoresis run times and voltage settings can vary depending on the gel apparatus used, but the bromophenol blue dye should always be run approximately 7.5 cm from the wells.

7.7 After the minigel has been run, disconnect the power supply and remove the gel. Wear UV protective eyewear and handle the gel with gloves.

7.8 To photograph the gel, place it on a UV transilluminator box under a stationary Polaroid® camera with a Kodak™ 22 or 23A Wratten filter. Photograph the gel in the dark under UV illumination with type 55, 57, or 667 black and white Polaroid film.

7.9 The following six bands should be present in the Control DNA 1 lane and in samples in which the DNA was amplified successfully for the AmpliType® PM and AmpliType PM+DOA1 Kits: 242/239 bp (HLA DOA1), 214 bp (LDLFI), 190 bp (GYPA), 172 bp (HBGG), 151 bp (D7S8), and 138 bp (GC). Primer-dimer bands and unincorporated primers may appear as broad bands near the bottom of the gel in the region of lower molecular weight. Occasionally, non-specific bands or smearing can be observed above and/or below the six specific bands, but they do not compromise the typing results on the AmpliType PM and AmpliType HLA DOA1 DNA Probe Strips.

7.10 If all six bands are present on the AmpliType PM and AmpliType PM+DOA1 PCR product gel, proceed to Section 8.0. If less than six bands are present, see Troubleshooting Section before proceeding.

8.0 DNA HYBRIDIZATION

Note: Ensure that 5 µL of 200 mM disodium EDTA has been added to each tube (Section 4.0) before proceeding to the DNA Hybridization step.

The AmpliType DNA Hybridization process involves three steps performed sequentially as follows: 1) hybridization of amplified DNA to DNA probe strips; 2) binding of HRP-SA Enzyme Conjugate to hybridized PCR products; and 3) stringent wash to remove non-specifically bound PCR products. Color Development is performed after the stringent wash step.

Certain steps in the following procedures involve the aspiration of solutions containing amplified DNA. Instructions for the assembly of the aspirator apparatus, its use, and waste disposal are described in the Laboratory Setup section of the AmpliType® User Guide.¹⁸

Before starting the DNA Hybridization and Color Development procedures, assemble the required reagents and equipment as follows:

Reagents and equipment required for both AmpliType kits.

HRP-SA Enzyme Conjugate (Included in Kit)

Chromogen: TMB Solution (Included in Kit; prepared in Section 1.1.1)

Hybridization Solution (prepared in Section 1.2.4)

Wash Solution (prepared in Section 1.2.5)

Citrate Buffer (prepared in Section 1.2.1)

30% Hydrogen Peroxide or 3% Hydrogen Peroxide (prepared in Section 1.2.6)

AmpliType™ DNA Typing Trays with pen for writing on strips (Part No. N808-0065)

Additional reagents and equipment required for the AmpliType PM Kit

AmpliType PM DNA Probe Strips (Included in Kit)

Additional Reagents and equipment required for the AmpliType PM+DOA1 Kit

AmpliType PM DNA Probe Strips (Included in Kit)

AmpliType HLA DOA1 DNA Probe Strips (Included in Kit)

Clean disposable gloves should be worn throughout the DNA Hybridization and Color Development steps to avoid soiling the DNA probe strips and to protect fingers from the 95°C block. Gloves should be discarded when leaving the work area to avoid transport of amplified DNA from the work area.

8.1 A rotating water bath, heated at 55°C (± 1°C), is necessary for the hybridization and wash steps. The rotating water bath must be on a level surface. Do NOT use a hot air shaker.

Note: It is essential to check the temperature with a calibrated total immersion thermometer before the hybridization step is performed. To obtain an accurate thermometer reading, minimize the time that the water bath remains uncovered.

8.2 Heat a covered, rotating water bath to 55°C and maintain the temperature between 54°C and 56°C. To ensure that the water bath will be at the correct temperature for the hybridization and wash steps, maintain the rotation at 50 to 70 rpm while performing Sections 8.3 through Section 8.9. The water level should be between 1/4 and 1/2 inch above the shaker platform. The water level should not be higher than 1/2 inch since higher levels may result in water splashing into the wells of the tray. An empty tray can be used to test the water level prior to use (place a 1 kg weight on the cover of the tray; see Section 8.10).

8.3 Warm the Hybridization Solution and the Wash Solution in a heated stationary water bath or an incubator (37° to 55°C) to dissolve the solids. All solids must be in solution and all solutions must be well mixed before use.

8.4 Remove the tube of AmpliType DNA Probe Strips from the 2° to 8°C storage area. Allow the tube to equilibrate to ambient temperature before opening to prevent condensation inside the tube. Using filter forceps, remove the required number of DNA probe strips from the tube. Using the pen included with the AmpliType™ DNA Typing Trays, label each strip in the space at the right edge of the strip. The use of other pens is not recommended because some inks may affect the quality of the typing results. Place one DNA probe strip in each well of the clean AmpliType DNA Typing Tray. Strips should all be in the same orientation.

Note: AmpliType® PM and AmpliType HLA DQA1 DNA Probe Strips can be used to type PCR products from the same PM amplification reaction at the same time but the DNA probe strips must be placed in separate wells of the tray.

8.5 Prepare the GeneAmp® PCR Instrument System to denature the amplified DNA by setting the parameter to hold the temperature at 95°C (see appropriate GeneAmp PCR Instrument System Users Manual). Start the file or program.

8.6 Place the tubes in the GeneAmp PCR Instrument System after it reaches 95°C. Press the tubes down tightly in the block. When using the GeneAmp PCR System 9600, position the heated cover over the tubes during the denaturation step to hold down the tubes to avoid refluxing into the cap. Denature the amplified DNA by incubation at 95°C for 3 to 10 minutes. Keep each tube at 95°C until use.

8.7 As shown in Figure 2, tilt the DNA Typing Tray towards the labelled end of the strips. The DNA Typing Tray lid turned upside down may be used as the solid support. Add 3 mL of pre-warmed Hybridization Solution to each well at the labelled end of each strip. Do not wet the remainder of the strip.




Figure 2. Handling solutions in DNA Typing Trays.

Note: Hybridization solution solids must be completely dissolved before adding to the tray.

8.8 Perform the following steps for each tube of amplified DNA:
NOTE: For each tube, perform 8.8.1 through 8.8.3 within 20 seconds and use a new pipet tip for each addition.

8.8.1 Remove the tube from the 95°C block of the GeneAmp PCR Instrument System. When using the GeneAmp PCR System 9600, it is not necessary to reposition the heated cover over the tubes between each sample removal.

8.8.2 Carefully open the tube using the Micro Tube Opener designated for use only with amplified DNA.

8.8.3 Withdraw 20 µL of amplified DNA and immediately add it below the surface of the Hybridization Solution in the well of the corresponding DNA probe strip (see Figure 2). If using the DNA Thermal Cycler 480, be sure to withdraw from the aqueous (bottom) layer.

8.8.4 Cap the tube after adding the denatured amplified DNA.

8.8.4.1 If using both of the DNA probe strips provided in the AmpliType® PM+DQA1 Kit, proceed to Section 8.8.4.2. If using the AmpliType PM Kit, or processing only the PM or HLA DQA1 DNA Probe Strips provided in the AmpliType PM+DQA1 Kit, set the capped tube aside. Repeat Sections 8.8.1 through 8.8.4.1 until each amplified DNA sample has been added to the corresponding well.

8.8.4.2 When both AmpliType PM and AmpliType HLA DQA1 DNA Probe Strips are used to type PCR product from the same tube, add 20 µL denatured amplified DNA to one of the two strips designated for this sample (e.g., PM DNA Probe Strip). Cap and return the tube to the 95°C block. Add the remaining denatured samples individually to the corresponding DNA probe strips (e.g., PM DNA Probe Strips) and return the capped tubes to the 95°C block. Repeat Sections 8.8.1 through 8.8.4.2 using the other AmpliType DNA Probe Strips (e.g., HLA DQA1 DNA Probe Strips) until each amplified DNA sample has been added to the corresponding wells. The capped tubes can be set aside instead of returned to the 95°C block after the amplified DNA is added to the second DNA probe strip.

The remaining amplified DNA samples can be stored at 2° to 8°C for 2 months or at -20°C for 6 months. The continued acceptable performance of these samples beyond these periods may vary with the sample. Store amplified DNA samples separate from all PCR amplification reagents, extracted DNA samples, and casework samples.

8.9 Put the lid on the tray and mix by carefully rocking the tray. Ensure that each strip is completely wet. Once hybridization has begun, strips should remain wet through the conclusion of the Color Development and Photography steps.

8.10 Stop the rotation of the rotating water bath (do not turn off the temperature control). Confirm that the temperature of the water is 55°C and transfer the tray to the water bath. Place a 1 kg weight (e.g., a lead ring) on the covered tray to prevent the tray from sliding or floating. Resume rotation of water bath at 50 to 70 rpm. Check the tray position and confirm that water does not splash into the wells of the tray.

8.11 Replace the water bath cover to maintain the temperature of the water in the bath at 55°C (± 1°C). Hybridize the amplified samples to the DNA-probe strips by incubating at 55°C for 15 minutes (± 2 minutes).

8.12 Approximately 5 minutes before the end of the hybridization step, prepare the Enzyme Conjugate Solution in a glass flask using the following equations to determine the volume of each component required:

Number of Strips	x	3.3 mL of Hybridization Solution	=	Volume of Hybridization Solution
Number of Strips	x	27 µL of Enzyme Conjugate: HRP-SA	=	Volume of Enzyme Conjugate: HRP-SA

Mix the solution thoroughly and ensure that the solids remain in solution. Leave at room temperature (15° to 30°C) until ready to use.

8.13 After hybridization, stop the rotation of the water bath and remove the tray. Replace the water bath cover and keep the water bath rotating between incubation steps to maintain the temperature at 55°C (± 1°C).

8.14 Aspirate the contents of each well from the labelled end of the strip while tilting the tray slightly. Remove condensation from the tray lid with a clean lab wipe. Use of paper towels to wipe the tray lid is not advised because some paper towels contain bleach which can affect the color development.

Note: Wash Solution solids must be completely dissolved and well mixed before use.

8.15 Dispense 5 mL of pre-warmed (37° to 55°C) Wash Solution into each well (a dispensing re-pipet is useful for this purpose). Rinse by gently rocking the tray for several seconds. Aspirate the solution from each well.

8.16 Dispense 3 mL of the Enzyme Conjugate Solution prepared in Section 8.12 into each well and cover with the lid. Stop the rotating water bath and transfer the tray to the 55°C water bath. Place a 1 kg weight on the covered tray to prevent the tray from sliding or floating. Adjust the rotating water bath to 50 to 70 rpm. Check the tray position and confirm that the water does not splash into the wells of the tray.

8.17 Replace the water bath cover to maintain the temperature at 55°C (± 1°C). Incubate the Enzyme Conjugate Solution with the DNA probe strips at 55°C (± 1°C) for 5 minutes (± 1 minute).

8.18 After incubation, stop the rotation of the water bath and remove the tray. Replace the water bath cover and keep the water bath rotating between incubation steps to maintain the temperature at 55°C (± 1°C).

8.19 Aspirate the contents of each well from the labelled end of the strip while tilting the tray slightly. Remove condensation from the tray lid with a clean lab wipe. Use of paper towels to wipe the tray lid is not advised because some paper towels contain bleach which can affect the color development. Replace the water bath cover and keep the water bath rotating between incubation steps to maintain the temperature at 55°C (± 1°C).

8.20 Dispense 5 mL of pre-warmed Wash Solution into each well. Rinse by gently rocking the tray for several seconds. Aspirate the solution from each well.

8.21 To perform the stringent wash step, dispense 5 mL of pre-warmed Wash Solution into each well. Cover the tray with the lid and place it in the 55°C water bath. Place a 1 kg weight on the covered tray to prevent the tray from sliding or floating. Adjust the rotation to 50 to 70 rpm and check the tray position to confirm that water does not splash into the wells of the tray.

8.22 Replace the water bath cover to maintain the temperature at 55°C (± 1°C). Incubate the DNA probe strips at 55°C (± 1°C) for 12 minutes (± 1 minute). The temperature and timing of the stringent wash step are critical.

8.23 After incubation, remove the tray from the water bath. Take off the lid and aspirate the contents of each well from the labelled end of the strips.

8.24 Dispense 5 mL of Wash Solution into each well. Gently rock the tray for several seconds. Slowly pour off or aspirate the contents from each well.

9.0 COLOR DEVELOPMENT

9.1 Dispense 5 mL of Citrate Buffer into each well. A dispensing re-pipet is useful for this purpose. Cover the tray with the lid and place it on an orbital shaker set at approximately 50 rpm at room temperature (15° to 30°C) for 5 minutes.

9.2 During this wash step, prepare the Color Development Solution. Do not prepare the Color Development Solution more than 10 minutes before use. Add the following reagents in the order listed to a glass flask and mix thoroughly by swirling. Protect from light. Do NOT vortex. Use the following equations to determine the volumes of each component required:

Number of Strips	x	5 mL Citrate Buffer	=	Volume of Citrate Buffer
Number of Strips	x	5 µL 3% Hydrogen Peroxide	=	Volume of Hydrogen Peroxide

OR
0.5 µL 30% Hydrogen Peroxide

Number of Strips	x	0.25 mL Chromogen: TMB Solution	=	Volume of Chromogen: TMB
------------------	---	---------------------------------	---	--------------------------

9.3 Remove the tray from the orbital shaker. Remove the lid and slowly pour off or aspirate the contents from each well. Add 5 mL of the freshly prepared Color Development Solution (Section 9.2) to each well.

Note: Place the lid on the tray and cover the lid with aluminum foil during Steps 9.4 and 9.5 to protect the DNA probe strips from strong light.

9.4 Develop the strips at room temperature (15° to 30°C) by rotating on an orbital shaker set at approximately 50 rpm for 20 to 30 minutes. Develop until the "S" or "C" dot is visible. Some strips may develop in less than 20 minutes. In this case, the color development may be stopped by proceeding immediately to Section 9.5. Color development of the negative control strip should not be stopped until development of all other test strips is stopped.

9.5 Stop the color development by removing the solution from the well. Immediately dispense approximately 5 mL of DI H₂O into the well. Place the tray on an orbital shaker set at approximately 50 rpm for 5 to 10 minutes. Remove the DI H₂O from the well and repeat the wash steps a minimum of three times. Additional 5 to 10 minute washes will reduce the potential for development of background color.

9.6 Record the pattern of blue dots from each wet DNA probe strip prior to photography (see Interpretation of Results).

Note: Keep strips wet throughout the photography steps.

12.0 INTERPRETATION OF RESULTS

Results are interpreted by observing the pattern and relative intensities of blue dots on the wet AmpliType[®] PM and AmpliType HLA DQA1 DNA Probe Strips to determine which alleles are present in the DNA sample.

12.1 Reading and Interpreting the AmpliType PM DNA Probe Strips

The AmpliType PM DNA Probe Strips have been spotted with a total of fourteen sequence-specific oligonucleotide probes to distinguish the alleles of five genetic loci (a mixture of two probes is spotted at the GYPA "A" allele position). Under the AmpliType hybridization conditions, the typing probes will bind specifically to PCR product containing the alleles designated on the AmpliType PM DNA Probe Strip.

To read the developed AmpliType PM DNA Probe Strip, the "S" dot is examined first and then each locus is examined separately. The standard probe "S" on the AmpliType PM DNA Probe Strip is identical in sequence to the control probe "C" on the AmpliType HLA DQA1 DNA Probe Strip and detects all of the HLA DQA1 alleles. The "S" dot is designed to be the lightest typing dot on the PM DNA Probe Strip and acts as a minimum dot intensity control for the remaining probes. It is recommended that a DNA probe strip with no visible "S" dot not be typed for any locus.

When a "S" dot is visible on the AmpliType PM DNA Probe Strip, the intensities of the dots at the remaining twelve positions are compared to the intensity of the "S" dot. Those dots that appear either darker than or equivalent to the "S" dot are considered positive. Each positive dot indicates the presence of the corresponding allele. Dots that are lighter than the "S" dot should be interpreted with care.^{14,15}

The dots on the AmpliType PM DNA Probe Strip correspond to the following alleles:

The "A" dot for each locus is positive in the presence of the A allele.

Note: The "A" dot for the GYPA locus is positive in the presence of both the A allele and the A' allele. Both the GYPA AB and GYPA A'B heterozygous genotypes have balanced dot intensities, but additional GYPA A and B variant alleles (observed in < 8% of the African American population) may produce a slightly imbalanced heterozygous signal.

The "B" dot for each locus is positive in the presence of the B allele.

The "C" dot for the HBGG and GC loci is positive in the presence of the C allele. For LDLR, GYPA, and D7S8, three genotypes are possible (AA, BB, and AB). For HBGG and GC, six genotypes are possible (AA, BB, CC, AB, AG, and BC).

An example of a developed AmpliType PM DNA Probe Strip using PCR product amplified from 2 ng of Control DNA 1 is shown in Figure 3. A sample from a single individual will produce balanced dot intensities within each locus for which the individual is heterozygous.

	LDLR	GYPA	HBGG	D7S8	GC
S	A	A'	A	B	C
	A	A'	A	B	C

Figure 3. The AmpliType PM types for Control DNA 1 are: LDLR BB, GYPA AB, HBGG AA, D7S8 AB, GC BB. The dot intensities of the GYPA and D7S8 loci are balanced (i.e., the intensities of the "A" and "B" dots within each locus are similar).

10.0 PHOTOGRAPHY AND STORAGE

10.1 Photographs should be taken for a permanent record of the results. Photographs must be taken while the DNA probe strips are still wet.

10.2 Place wet strips on a flat non-absorbent surface (a black background, such as a sheet of exposed X-ray film, is recommended to enhance contrast). Keep the strips wet throughout the photography steps. Minimize exposure to strong light.

10.3 Use a Polaroid[®] camera with Type 55, 57, or 667 (black and white) film or Type 59 or 559 (color) film.

Note: For black and white photography, an orange filter (Wratten 22 or 23A) will enhance contrast.

10.4 Follow Polaroid film exposure and development instructions.

10.5 After photography, the DNA probe strips may be air-dried on any hard, non-absorbent surface or Whatman[®] 3MM chromatography paper (do not use paper towels for this step). Protect from light and oxidizing agents (e.g., acid-treated paper, bleach, and nitric acid). The dot intensities fade upon drying.

11.0 DISPOSAL OR REUSE OF TYPING TRAYS

11.1 The AmpliType[™] DNA Typing Trays are designed to be disposable but may be reused. If the trays are to be reused, immediately wash the trays and lids according to the following procedure:

Note: Trays which have become discolored due to repeated use should be discarded.

11.1.1 Add approximately 5 to 10 mL of 95% ethanol or 70% isopropanol to each well of the used AmpliType DNA Typing Tray.

Note: Do NOT use detergent or bleach.

11.1.2 Cover the tray with the lid and carefully agitate for 15 to 30 seconds to dissolve any residual Chromogen: TMB.

11.1.3 Remove the lid. Visually inspect each well for a blue or yellow color which will indicate the presence of Chromogen: TMB. If necessary, repeat Sections 11.1.1 and 11.1.2 until no color is present and then rinse each well in the tray and the tray lid with DI H₂O. If all the color cannot be removed from the walls of the wells, then discard the tray.

11.1.4 Dry tray and lid before reuse.

12.2 Reading and Interpreting the AmpliType® HLA DQA1 DNA Probe Strips

The AmpliType® HLA DQA1 DNA Probe Strips have been spotted with a total of eleven sequence-specific oligonucleotide probes to detect eight alleles of the HLA DQA1 locus. Under the AmpliType hybridization conditions, the typing probes will bind specifically to PCR product containing the alleles designated on the AmpliType HLA DQA1 DNA Probe Strip.

To read the developed AmpliType HLA DQA1 DNA Probe Strip, the "C" dot is examined first and then the remaining dots are examined. The control probe "C" on the AmpliType HLA DQA1 DNA Probe Strip detects all of the HLA DQA1 alleles. The "C" dot is designed to be the lightest typing dot on the strip and it indicates that adequate amplification and typing of the HLA DQA1 alleles in the sample have occurred. If the "C" dot is absent, an accurate determination of the type cannot be made. Additional information on the "C" dot can be found in the AmpliType User Guide.¹⁸

The accurate interpretation of the HLA DQA1 results depends on the presence and intensity of the "C" dot. The intensities of the dots at the remaining ten positions are compared to the intensity of the "C" dot. Those dots that appear either darker than or equivalent to the "C" dot are considered positive. Each positive dot indicates the presence of the corresponding HLA DQA1 allele. Dots with signals less than the "C" dot should be interpreted with care.^{9,11,28}

The dots on the AmpliType HLA DQA1 DNA Probe Strip correspond to the following alleles:

- The "1" dot is positive in the presence of the HLA DQA1 1.1, 1.2, and 1.3 alleles.
- The "2" dot is positive only in the presence of the HLA DQA1 2 allele.
- The "3" dot is positive only in the presence of the HLA DQA1 3 allele.
- The "4" dot is positive in the presence of the HLA DQA1 4.1, 4.2, and 4.3 alleles.
- Four HLA DQA1 sub-typing probes differentiate the HLA DQA1 1.1, 1.2, and 1.3 alleles.

The "1.1" dot is positive only in the presence of the HLA DQA1 1.1 allele.

Note: A faint "1.1" dot will appear with some HLA DQA2 pseudogene alleles.³⁰

The "1.3" dot is positive only in the presence of HLA DQA1 1.3 allele.

Note: There is no probe that detects only the HLA DQA1 1.2 allele.

The "1.2, 1.3, 4" dot is positive in the presence of HLA DQA1 1.2, 1.3, 4.1, 4.2, and 4.3 alleles.

Note: The "1.2, 1.3, 4" dot can be lighter than the "C" dot when the genotype has a HLA DQA1 4.2 or 4.3 allele paired with a HLA DQA1 1.1, 2, 3, 4.2, or 4.3 allele because the HLA DQA1 4.2 and 4.3 alleles each have a single partially destabilizing mismatch to the "1.2, 1.3, 4" probe.¹⁸ The partially destabilizing mismatch allows these two alleles to bind to this probe weakly relative to the HLA DQA1 1.2, 1.3, and 4.1 alleles.

The "All but 1.3" dot is positive in the presence of all HLA DQA1 alleles EXCEPT 1.3. This probe is necessary to differentiate the 1.2, 1.3 genotype from the 1.3, 1.3 genotype.

Note: The "All but 1.3" dot can be equal to or lighter than the "C" dot when the genotype has a HLA DQA1 1.3 allele paired with a HLA DQA1 4.1, 4.2, or 4.3 allele because the HLA DQA1 4.1, 4.2, and 4.3 alleles have a single partially destabilizing mismatch to the "All but 1.3" probe.¹⁸ The partially destabilizing mismatch allows these three alleles to bind to this probe weakly relative to the HLA DQA1 1.1, 1.2, 2, and 3 alleles.

Two additional HLA DQA1 sub-typing probes differentiate the HLA DQA1 4.1 allele from the HLA DQA1 4.2 and 4.3 alleles.

The "4.1" dot is positive only in the presence of the HLA DQA1 4.1 allele. The "4.2, 4.3" dot is positive in the presence of HLA DQA1 4.2 and 4.3 alleles.

These HLA DQA1 allele designations correspond to the World Health Organization (WHO) nomenclature as shown in Table 6.

Table 6. World Health Organization (WHO) Nomenclature for HLA DQA1 alleles

Allele	WHO Designation	Allele	WHO Designation
1.1	0101	3	0301
1.2	0102	4.1	0501*
1.3	0103	4.2	0401*
2	0201	4.3	0601

*Note that the WHO number designations are out of sequence compared to the original allele designations.

An example of a developed AmpliType® HLA DQA1 DNA Probe Strip using PCR product amplified from 2 µg of Control DNA 1 is shown in Figure 4.



Figure 4. The AmpliType HLA DQA1 type for Control DNA 1 is 1.1, 4.1.

PERFORMANCE CHARACTERISTICS

The AmpliType® PM and AmpliType PM+DQA1 PCR Amplification and Typing Kits are developed and manufactured by Roche Molecular Systems, Inc. (RMS). Each lot of the AmpliType PM and AmpliType PM+DQA1 PCR Amplification and Typing Kits is carefully tested by RMS to ensure that the kits perform according to specifications and are free from interfering contaminants.

The user of the AmpliType PCR Amplification and Typing Kits will be able to amplify and type a minimum of two nanograms of Control DNA 1 when employing the protocols and reagents provided in the kits.

In the laboratories of RMS, the kit components have been used successfully to type samples containing less than one nanogram of human DNA.

TROUBLESHOOTING

Observation	Possible Cause	Recommended Action
1. No signal or faint signal from both the Control DNA 1 and the DNA test samples at all loci.	No PCR amplification or insufficient PCR amplification of all markers. Improper hybridization or assay condition.	Check amplified DNA on agarose gel (see Observation 1.1). Repeat test from Hybridization step (see Observation 1.2).
1.1 No amplified product visible on gel.	No DNA added or insufficient DNA added to PCR Reaction Mix. AmpliType® PM Primer Set not added to AmpliType PM PCR Reaction Mix. GeneAmp® PCR Instrument System failure or wrong program.	Quantitate DNA and add 2 - 10 ng DNA; repeat test. Add AmpliType PM Primer Set; repeat test. See GeneAmp PCR Instrument System Manual and check instrument calibration.
1.2 Amplified product visible on gel but no signal or faint signal on AmpliType DNA Probe Strips.	Tubes not sealed tightly in the DNA Thermal Cycler 480 block during amplification. MicroAmp® Base used with tray and tubes in GeneAmp PCR System 9600. Hybridization and/or Stringent Wash temperature too high.	Push tubes firmly into contact with block after first cycle; repeat test. Remove MicroAmp Base; repeat test. Check that the rotating water bath temperature is at 55°C (± 1°C) with an immersible thermometer; repeat test.
	Hybridization and/or Wash Solution salt concentration too low.	Prepare new solutions; repeat test.
	Stringent Wash time too long.	Repeat test, washing for 12 minutes (± 1 minute) only.

Observation	Possible Cause	Recommended Action
1.2 (continued) Amplified product visible on gel but no signal or faint signal on AmpliType® DNA Probe Strips.	Inadequate agitation of the DNA probe strips during hybridization.	Check speed of rotating water bath (50 to 70 rpm); verify that hybridization solution is washing over the strips; repeat test.
	Amplified DNA was not added to DNA probe strips.	Repeat test, adding amplified DNA to DNA probe strips.
	Amplified DNA was not denatured.	Check GeneAmp® PCR Instrument System block temperature is 95°C; leave sample in block >3 minutes using the heated block cover; repeat test following instructions in Section 8.6.
	Enzyme Conjugate: HRP-SA was not added to the Enzyme Conjugate Solution.	Prepare new diluted Enzyme Conjugate: HRP-SA solution; repeat test.
	Hydrogen peroxide was not added or too much was added to the Color Development Solution.	Make new Color Development Solution with correct amount of hydrogen peroxide; repeat test.
	Hydrogen peroxide inactive.	Make new Color Development Solution using new bottle or dilution of hydrogen peroxide; repeat test.
	Chromogen: TMB was not added to the Color Development Solution.	Make new Color Development Solution adding Chromogen: TMB; repeat test.
	The original AmpliType HLA DQA instead of AmpliType PM and PM+DQA1 typing protocol was followed.	Repeat test following the typing protocol in Section 8.0.

2. Positive signal from Control DNA 1, but no signal from DNA test sample.

Quantity of DNA test sample is below the assay sensitivity.

Quantitate DNA and add 2 - 10 ng DNA; repeat test.

Observation	Possible Cause	Recommended Action
2. (continued) Positive signal from Control DNA 1, but no signal from DNA test sample.	Test sample contains PCR inhibitor (e.g., heme compounds, certain dyes).	Any or all of the following actions may be taken: 1. Quantitate DNA and add minimum necessary volume; repeat test. 2. Wash the sample in Centricon® 100 (see Reference 11); repeat test. 3. Add BSA to reaction and see References 15 and 31; repeat test.
	Test sample DNA is degraded.	If possible, evaluate the quality of the DNA sample by running an agarose gel (see Reference 18). If the DNA is degraded, reamplify with an increased amount of DNA.
	Low or lack of SDS in Hybridization and/or Wash Solution.	Prepare new Hybridization and/or Wash Solution with correct amount of SDS; repeat test.
	Inadequate agitation of the DNA probe strips during hybridization and washing steps.	Check speed of rotating water bath (50 to 70 rpm); verify that solutions are washing over strips; repeat test.
	Tray lid not wiped adequately.	Wipe lid (Section 8.19); repeat test.
	Excess amounts of Enzyme Conjugate: HRP-SA added to Enzyme Conjugate Solution.	Prepare new Enzyme Conjugate Solution with correct amount of Enzyme Conjugate: HRP-SA; repeat test.
	Exposure to light during Color Development.	Cover tray lid with foil during Color Development steps; repeat test.

3. High DNA probe strip background color.

Prepare new diluted Enzyme Conjugate: HRP-SA solution; repeat test.

Make new Color Development Solution with correct amount of hydrogen peroxide; repeat test.

Make new Color Development Solution using new bottle or dilution of hydrogen peroxide; repeat test.

Make new Color Development Solution adding Chromogen: TMB; repeat test.

Repeat test following the typing protocol in Section 8.0.

Quantitate DNA and add 2 - 10 ng DNA; repeat test.

Observation

Possible Cause

Recommended Action

3. (continued) High DNA probe strip background color.

Use of water other than deionized or glass distilled water for water rinses.

Repeat test using deionized or glass distilled water for water rinses.

4. High DNA probe strip background color upon storage.

Exposure to strong light and oxidizing agents.

Store strips in the dark away from oxidizing agents.

5. Presence of unexpected or additional dots in the amplified Control DNA 1 sample.

Insufficient water washes after Color Development.

Increase number of water washes in future assays.

6. Signals weaker than "S" or "C" dot on the same strip (also see Observations 12 and 13).

Hybridization and/or Stringent Wash temperature too low.

Check that the rotating water bath temperature is at 55°C (± 1°C) with a total immersion thermometer; repeat test.

7. More than two alleles present on the AmpliType HLA DQA1 DNA Probe Strip or the HBGG and/or GC marker on the AmpliType PM DNA Probe Strip.

Cross-hybridization caused by Hybridization and/or Wash Solution salt concentration being too high.

Prepare new solutions; repeat test.

8. Some, but not all, loci visible on gel of AmpliType PM PCR products.

Cross-hybridization caused by Stringent Wash time being too short.

Repeat test, washing for 12 minutes (± 1 minute).

9. Amplification of HLA DQA2 pseudogene (faint 1.1 dot).

Contamination by amplified product or samples.

See References 18 and 28.

10. Test sample DNA is degraded.

Hybridization and/or Stringent Wash temperature too high or too low.

Check that the rotating water bath temperature is at 55°C (± 1°C) with an immersion thermometer; repeat test.

11. If possible, evaluate the quality of the DNA sample by running an agarose gel (see Reference 18). If the DNA is degraded, reamplify with an increased amount of DNA.

Hybridization and/or Wash Solution salt concentration too high or too low.

Prepare new solutions; repeat test.

12. Repeat test, washing for 12 minutes (± 1 minute).

Stringent wash time too long or too short.

Repeat test, washing for 12 minutes (± 1 minute).

13. See References 14, 15, 17, 18, and 28.

Mixed sample or contamination.

See References 14, 15, 17, 18, and 28.

Observation

Possible Cause

Recommended Action

6. (continued) Signals weaker than "S" or "C" dot on the same strip (also see Observations 12 and 13).

Amplification of HLA DQA2 pseudogene (faint "1.1" dot).

See References 18 and 30.

7. More than two alleles present on the AmpliType HLA DQA1 DNA Probe Strip or the HBGG and/or GC marker on the AmpliType PM DNA Probe Strip.

EDTA not added to the reaction prior to the heat denaturation step of the DNA hybridization protocol.

Add EDTA to amplified sample (Section 4.0); repeat test.

8. Some, but not all, loci visible on gel of AmpliType PM PCR products.

Cross-hybridization caused by Hybridization and/or Stringent Wash temperature too low.

Check that the rotating water bath temperature is at 55°C (± 1°C) with a total immersion thermometer; repeat test.

9. Amplification of HLA DQA2 pseudogene (faint 1.1 dot).

Cross-hybridization caused by Hybridization and/or Wash Solution salt concentration too high.

Prepare new solutions; repeat test.

10. Repeat test, washing for 12 minutes (± 1 minute).

Stringent Wash time too short.

Repeat test washing for 12 minutes (± 1 minute).

11. See References 14, 15, 17, 18, and 28.

Mixed sample or contamination.

See References 14, 15, 17, 18, and 28.

12. See References 18 and 30.

Amplification of HLA DQA2 pseudogene (faint 1.1 dot).

See References 18 and 30.

13. If possible, evaluate the quality of the DNA sample by running an agarose gel (see Reference 18). If the DNA is degraded, reamplify with an increased amount of DNA.

Test sample DNA is degraded.

If possible, evaluate the quality of the DNA sample by running an agarose gel (see Reference 18). If the DNA is degraded, reamplify with an increased amount of DNA.

Observation	Possible Cause	Recommended Action
8. (continued) Some, but not all, loci visible on gel of AmpliType PM PCR products.	Test sample contains PCR inhibitor (e.g., heme compounds, certain dyes).	Any or all of the following actions may be taken: 1. Quantitate DNA and add minimum necessary volume; repeat test. 2. Wash the sample in Centricon® 100 (see Reference 11); repeat test. 3. Add BSA to reaction and see References 15 and 31; repeat test.
9. Some, but not all, loci visible on AmpliType PM DNA Probe Strip.	Not all loci amplified. Amplified DNA was not denatured.	Verify the presence of amplified loci on an agarose gel (see Observation 8). Check GeneAmp PCR Instrument System block temperature is 95°C; leave sample in block > 3 minutes using the heated block cover; repeat test following instructions in Section 8.6.
10. Imbalanced dot intensities within a locus on the AmpliType PM DNA Probe Strip (this kit is designed to produce balanced dot intensities when heterozygous samples are typed, except as described in Section 12.1).	Hybridization and/or Stringent Wash temperature too high or too low. Hybridization and/or Wash Solution salt concentration too high or too low. Stringent Wash time too long or too short.	Check that the rotating water bath temperature is at 55°C (± 1°C) with a total immersion thermometer; repeat test. Prepare new solutions; repeat test. Repeat test, washing for 12 minutes (± 1 minute).
10. (continued) Imbalanced dot intensities within a locus on the AmpliType PM DNA Probe Strip (this kit is designed to produce balanced dot intensities when heterozygous samples are typed, except as described in Section 12.1).	Mixed sample or contamination. EDTA not added to the reaction prior to the heat denaturation step of the DNA hybridization protocol.	See References 14, 15, 17, 18, and 28. Add EDTA to amplified sample (Section 4.0); repeat test.
11. Weak or absent "4.1" dot on the AmpliType HLA DQA1 DNA Probe Strip in the amplified Control DNA 1 sample.	EDTA not added to the reaction prior to the heat denaturation step of the DNA hybridization protocol.	Add EDTA to amplified sample (Section 4.0); repeat test.
12. "1.2, 1.3, 4" dot weaker than "C" dot on AmpliType HLA DQA1 DNA Probe Strip.	Genotype of sample has a HLA DQA1 4.2 or 4.3 allele paired with a HLA DQA1 1.1, 2, 3, 4.2, or 4.3 allele.	See Reference 18.
13. "1.1" dot weaker than "C" dot but no signal for "1" dot on AmpliType HLA DQA1 DNA Probe Strip.	Amplification of HLA DQA2 pseudogene (faint 1.1 dot).	See References 18 and 30.
14. "All but 1.3" signal weaker than "C" dot on AmpliType HLA DQA1 DNA Probe Strip.	Genotype of sample has a HLA DQA1 1.3 allele paired with a HLA DQA1 4.1, 4.2, or 4.3 allele.	See Reference 18.

REFERENCES

1. von Beroldingen, C.H., Blake, E.T., Higuchi, R., Sensabaugh, G.F., and Erlich, H.A. 1989. Applications of PCR to the analysis of biological evidence. In: *PCR Technology, Principles and Applications for DNA Amplification*, ed. Erlich, H.A., Stockton Press, Inc. New York, NY, p 209-223.
2. Reynolds, R., Sensabaugh, G., and Blake, E. 1991. Analysis of genetic markers in forensic DNA samples using the polymerase chain reaction. *Analytical Chemistry* **63**: 2-15.
3. Saiki, R.K., Gelfand, D.H., Stoffel, S., Scharf, S.J., Higuchi, R., Horn, G.T., Mullis, K.B., and Erlich, H.A. 1988. Primer-directed enzymatic amplification of DNA with a thermostable DNA polymerase. *Science* **239**: 487-491.
4. Higuchi, R., von Beroldingen, C.H., Sensabaugh, G.F., and Erlich, H.A. 1988. DNA typing from single hairs. *Nature* **332**: 543-546.
5. Pääbo, S., Gifford, J.A., and Wilson, A.C. 1988. Mitochondrial DNA sequences from a 7000-year old brain. *Nucleic Acids Research* **16**: 9775-9787.
6. Hochmeister, M.N., Budowle, B., Jung, J., Borer, U.V., Comey, C.T., and Dimhofer, R. 1991. PCR-based typing of DNA extracted from cigarette butts. *International Journal of Legal Medicine* **104**: 229-233.
7. Hochmeister, M.N., Budowle, B., Borer, U.V., Eggmann, U., Comey, C.T., and Dimhofer, R. 1991. Typing of deoxyribonucleic acid (DNA) extracted from compact bone from human remains. *Journal of Forensic Sciences* **36**: 1649-1661.
8. Sajantila, A., Strom, M., Budowle, B., Karhunen, P.J., and Peltonen, L. 1991. The polymerase chain reaction and post-mortem forensic identity testing: application of amplified D1S80 and HLA-DQ alpha loci to the identification of fire victims. *Forensic Science International* **51**: 23-34.
9. Blake, E., Mihalovich, J., Higuchi, R., Walsh, P.S., and Erlich, H. 1992. Polymerase Chain Reaction (PCR) amplification and human leukocyte antigen (HLA)-DQ α oligonucleotide typing on biological evidence samples: casework experience. *Journal of Forensic Sciences* **37**: 700-726.
10. Comey, C.T. and Budowle, B. 1991. Validation studies on the analysis of the HLA DQ α locus using the polymerase chain reaction. *Journal of Forensic Sciences* **36**: 1633-1648.
11. Comey, C.T., Budowle, B., Adams, D.E., Baumstark, A.L., Lindsey, J.A., and Presley, L.A. 1993. PCR amplification and typing of the HLA DQ α gene in forensic samples. *Journal of Forensic Sciences* **38**: 239-249.
12. Sullivan, K., et al. 1992. Characterization of HLA DQA1 for forensic purposes. Allele and genotype frequencies in British Caucasian, Afro-Caribbean and Asian populations. *International Journal of Legal Medicine* **105**: 17-20.
13. Saiki, R.K., Walsh, P.S., Levenson, C.H., and Erlich, H.A. 1989. Genetic analysis of amplified DNA with immobilized sequence-specific oligonucleotide probes. *Proceedings of the National Academy of Sciences, USA* **86**: 6230-6234.
14. Herrin, G., Jr., Fildes, N., and Reynolds, R. 1994. Evaluation of the AmpliType[®] PM DNA test system on forensic case samples. *Journal of Forensic Sciences* **39**: 1247-1253.
15. Budowle, B., Lindsey, J.A., DeCout, J.A., Koons, B.W., Giusti, A.M., and Comey, C.T. 1995. Validation and population studies of the loci LDLR, GYPA, HBBG, D7S8, and GC (PM loci), and HLA-DQ α using a multiplex amplification and typing procedure. *Journal of Forensic Sciences* **40**: 45-50.
16. Roy, R. and Reynolds, R. 1995. AmpliType[®] PM and HLA DQ α typing from pap smear, semen smear and post-coital slides. *Journal of Forensic Sciences* **40**: 266-269.
17. Fildes, N. and Reynolds, R. 1995. Consistency and reproducibility of AmpliType[®] PM results between seven laboratories: Field trial results. *Journal of Forensic Sciences* **40**: 279-286.
18. Perkin-Elmer. 1990. AmpliType[®] User Guide, Version 2.
19. Gyllenstein, U.B. and Erlich, H.A. 1988. Generation of single-stranded DNA by the polymerase chain reaction and its application to direct sequencing of the HLA-DQ α locus. *Proceedings of the National Academy of Sciences, USA* **85**: 7652-7656.
20. Yamamoto, T., Davis, C.G., Brown, M.S., Schneider, W.J., Casey, M.L., Goldstein, J.L., and Russell, D.W. 1984. The Human LDL Receptor. A cysteine-rich protein with multiple Alu sequences in its mRNA. *Cell* **39**: 27-38.
21. Kudo, Shirohji and Fukuda, Minoru. 1989. Structural organization of glycoprotein A and B genes. Glycophorin B gene evolved by homologous recombination at Alu repeat sequences. *Proceedings of the National Academy of Sciences, USA* **86**: 4619-4623.
22. Slightom, J.L., Blechl, A.E. and Smithies, O. 1980. Human fetal γ - and γ -globin genes: complete nucleotide sequences suggest that DNA can be exchanged between these duplicated genes. *Cell* **21**: 627-638.
23. Horn, G.T., Richards, B., Merrill, J.J., and Klinger, K.W. 1980. Characterization and rapid diagnostic analysis of DNA polymorphisms closely linked to the cystic fibrosis locus. *Clinical Chemistry* **36**: 1614-1619.
24. Yang, F., Brune, J.L., Naylor, S.L., Cupples, R.L., Naberhaus, K.H., and Bowman, B.H. 1985. Human group-specific component (Gc) is a member of the albumin family. *Proceedings of the National Academy of Sciences, USA* **82**: 7994-7998.
25. Reynolds, R.L. and Sensabaugh, G.F. 1990. Use of the Polymerase Chain Reaction for typing Gc variants. In: *Advances in Forensic Haemogenetics* **3**, ed. Polesky, H.F. and Mayr, W.R., Springer-Verlag Berlin Heidelberg, p 158-161.
26. Erlich, H.A. and Bugawan, T.L. 1989. HLA Class II Gene Polymorphism: DNA typing, evolution, and relationship to disease susceptibility. In: *PCR Technology, Principles and Applications for DNA Amplification*, ed. Erlich, H.A., Stockton Press, Inc. New York, NY, p 193-208.
27. Fildes, N., Reynolds, R., and Erlich, H.A. In preparation, Roche Molecular Systems, Inc.

28. Higuchi, R. and Kwok, S. 1989. Avoiding false positives with PCR. *Nature* 339: 237-238.
29. Watson, R. 1989. The formation of primer artifacts in polymerase chain reactions. *Amplifications - A Forum for PCR Users* (Perkin-Elmer Newsletter) 2: 5-6.
30. Crouse, C.A., Vincsek, V., and Caraballo, B.K. 1994. Analysis and Interpretation of the HLA DQA "1.1 weak-signal" observed during the PCR-based typing method. *Journal of Forensic Sciences* 39: 41-51.
31. Corney, C.T., Koons, B.W., Presley, K.W., Smerick, J.B., Sobieralski, C.A., Stanley, D.M., and Baechtel, F.S. 1994. DNA Extraction Strategies for Amplified Fragment Length Polymorphism Analysis. *Journal of Forensic Sciences* 39: 1254-1269.

Perkin-Elmer® and MicroAmp® are registered trademarks of The Perkin-Elmer Corporation. Roche® is a registered trademark of Hoffmann-La Roche Inc. AmpliTaq® and AmpliType® are registered trademarks and AmpliType™ is a trademark of Roche Molecular Systems, Inc. GeneAmp® is a registered trademark of Roche Molecular Systems Inc., licensed to The Perkin-Elmer Corporation.

Centricor® is a registered trademark of Amicon, Inc. Combitps® is a registered trademark of Eppendorf. Kodak™ is a trademark of Eastman Kodak Company. Nalgen® is a registered trademark of Nalge Company. NuSieve®, Seakem®, and GIT® are registered trademarks of FMC Corporation.

Parafilm® is a registered trademark of American National Can Company.

Polaroid® is a registered trademark of Polaroid Company.

Repipe® is a registered trademark of Labindustries, Inc.

Sharp® is a registered trademark of Santorid Corporation.

Whatman® is a registered trademark of Whatman Inc.

Property of Idaho State Police Forensic Services
 Uncontrolled Internet Copy
 OBSOLETE DOCUMENT

QuantIBlot™
Human DNA Quantitation Kit

PERKIN ELMER

Part No. N808-0114 FOR FORENSIC AND RESEARCH USE ONLY
NOT FOR USE IN DIAGNOSTIC PROCEDURES

READ ENTIRE PROCEDURE PRIOR TO USE OF KIT

TABLE OF CONTENTS

List of Kit Components.....1
Introduction.....2
Materials Required But Not Supplied.....3
1.0 Reagent Preparation.....5
2.0 Information Regarding Protocols.....7
3.0 Slot Blotting/Immobilization of DNA.....7
4.0 DNA Hybridization.....9
5.0 Detection Steps.....10
6.0 Results Interpretation.....12
Performance Characteristics.....12
Troubleshooting.....13
References.....15
Ordering Information.....16

LIST OF KIT COMPONENTS

IMPORTANT: Do NOT freeze any of the kit components.

Store the QuantIBlot™ Human DNA Quantitation Kit components at 2° to 8°C. Under these conditions, components of the kit are stable through the control date printed on the label.

Note: The QuantIBlot™ Human DNA Quantitation Kit contains reagents for at least 10 hybridization reactions. Each hybridization reaction requires the DNA Standards and DNA Calibrators with space for testing up to 38 samples.

Reagent	Volume	Description
QuantIBlot™ D17Z1 Probe	220 µL	1 vial containing 1 pmole/µL biotinylated DNA oligonucleotide probe of the following sequence: 5'-biotin-TAGAAGCATTCTCAGAA ACTACTTTGTGATGATTC-3'.
Enzyme Conjugate: HRP-SA	2.0 mL	1 vial containing Horseradish Peroxidase-Streptavidin (HRP-SA) conjugate supplied in buffer.
Bromothymol Blue Solution	200 µL	1 vial containing 0.04% Bromothymol Blue in water.
QuantIBlot™ DNA Standard A	250 µL	1 vial containing 2 ng/µL human genomic DNA in TE buffer.
QuantIBlot™ DNA Calibrator 1	100 µL	1 vial containing 0.7 ng/µL human genomic DNA in TE buffer.
QuantIBlot™ DNA Calibrator 2	100 µL	1 vial containing 0.1 ng/µL human genomic DNA in TE buffer.
Package Insert		QuantIBlot™ Human DNA Quantitation Kit protocol.

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

INTRODUCTION

The QuantiBiot™ Human DNA Quantitation Kit provides reagents and protocols necessary for the rapid quantitation of human DNA. The procedure is based on the hybridization of a biotinylated oligonucleotide probe to DNA samples immobilized on a nylon membrane.¹ The probe included in this kit is complementary to a primate-specific alpha satellite DNA sequence at the locus D17Z1.² Subsequent binding of Enzyme Conjugate-HRP-SA (horseradish peroxidase-streptavidin) to the biotin moiety of the probe allows for either colorimetric or chemiluminescent detection (protocols are provided for both detection schemes). In the case of colorimetric detection, the oxidation of 3,3',5,5'-tetramethylbenzidine (Chromogen: TMB) catalyzed by the horseradish peroxidase results in the formation of a blue-colored precipitate directly on the membrane. Alternatively, for chemiluminescent detection the oxidation of a luminol based reagent catalyzed by the horseradish peroxidase results in the emission of photons that are detected on standard autoradiography film.³ This process is called Enhanced Chemiluminescence (ECL™). In both cases, the quantity of sample DNA is determined by comparison of the signal intensity to human DNA standards.

Using the protocols provided, 0.15 to 10.0 nanograms of human DNA can be quantitated. The entire protocol can be performed in less than 2 hours. DNA prepared by a variety of extraction methods, including phenol-chloroform, Chelex®, salting-out, or binding to silica particles, can be quantitated.

The QuantiBiot Human DNA Quantitation Kit provides protocols for two alternative detection schemes: colorimetric and chemiluminescent. The advantages of the colorimetric detection scheme are as follows: 1) A darkroom equipped for X-ray film development is not necessary; and 2) The Chromogen:TMB reagent is available directly from Perkin-Elmer (Part No. N808-0092). The advantages of the chemiluminescent (ECL) detection scheme are as follows: 1) If desired, the sensitivity can be increased below 0.15 ng by performing long exposures to the X-ray film (see Section 5.2.5); and 2) The X-ray film result is a permanent record (no photography is required).

MATERIALS REQUIRED BUT NOT SUPPLIED

The items listed are those required for quantitation procedures. This list does not include reagents or equipment required for DNA extraction (for a list of DNA extraction materials see the AmpliType™ User Guide, Version 2, Section 3). Many of the items listed are available from major laboratory suppliers (MLS) such as Baxter Scientific Products, Fisher Scientific, or VWR unless otherwise noted. Equivalent sources may be acceptable where noted.

Materials	Source
AmpliType™ User Guide, Version 2	Perkin-Elmer, Norwalk, CT (Part No. N808-0111).
Reagents	
Use reagent grade unless otherwise noted.	
Crystalline citric acid, monohydrate	JT Baker, Phillipsburg, NJ (Cat. No. 0110-05).
Deionized or distilled water	
Disodium ethylenediaminetetraacetic acid, dihydrate	GIBCO BRL, Gaithersburg, MD (Cat. No. 5575UAA), or Sigma St. Louis, MO (Cat. No. E5513).
Hydrochloric acid, concentrated	JT Baker (Cat. No. 9535-01).
Hydrogen peroxide, 30%	Sigma Chemical Company, (Cat. No. H1009).
Sodium Chloride	MLS
Sodium dodecyl sulfate, ultra pure electrophoresis grade	GIBCO BRL (Cat. No. 5525UAA).
5N Sodium hydroxide solution	VWR, San Francisco, CA (Cat. No. JT5671-2).
Sodium hydroxide pellets	MLS
Sodium phosphate, monobasic, monohydrate	Sigma (Cat. No. S9638).
The base (Trizma base [FW 121.1])	GIBCO BRL (Cat. No. 5504UAA), or Sigma (Cat. No. 8404).
Trisodium citrate, dihydrate	GIBCO BRL (Cat. No. 5584UAA).
Detection Reagents	
<u>For colorimetric detection, the following reagents are needed:</u>	
Chromogen:TMB	Perkin-Elmer (Part No. N808-0092).
100% ethanol (in a glass container)	Quantum Chemical Company, Cincinnati, OH.
<u>For chemiluminescent (ECL™) detection, the following reagents are needed:</u>	
ECL™ Detection Reagents	Amersham, Arlington Heights, IL (Cat. No. RPN-2109).
Film processing GBX Fixer	Sigma (Cat. No. P-7167).
Film processing GBX Developer	Sigma (Cat. No. P-7042).
General Equipment	
Autoclave	MLS
Hybridization Tray	Perkin-Elmer (Part No. N808-0136).
Hybridization Tray Retainer	Perkin-Elmer (Part No. N808-0137).
Balance, accurate to 1.0 mg	MLS
Bellco Hot Shaker Plus™ rotating water bath, adjustable to 65°C	Bellco, Vineland, NJ [Model No. 7746-22110 (110V), or Model No. 7746-22220 (220V)].

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

General Equipment (continued)

Complete immersion thermometer (N.B.S. or equivalent, graduated to 0.2°C, range bracketing 37° to 55°C)
Lead weights, (e.g., 1 kg lead ring)

Source (continued)

Cole Parmer, Niles, IL
(Cat No. G-08001-34).
VWR Brockport, NJ (Cat. No. 29700-048), or Cole Parmer
Chicago, IL (Cat. No. G-06137-06).

Labware and glassware
Magnetic stirrer and stir bars
Polaroid™ Camera and Type 55 or 667
(black and white) film or Type 59 or
559 (color) film.

pH meter
Pipettors
-adjustable to deliver 1-20 µL, 20 - 200 µL
and 200-1000 µL

Refrigerator (2° to 8°C)
Shaker, variable speed, orbital platform
(capable of 100-125 rpm)
Slot blot apparatus, The Convertible™

Timer, 60 minute (±1 minute)
Vortex
Water bath, stationary, adjustable to 37°C
or incubator

Additional equipment needed for
chemiluminescent detection:

Automatic film processor (optional)
Film cassette
Dark room
MLS

Supplies
Biodyne® B nylon membrane, 0.45 µm
Gloves, disposable
Lab wipes or paper towels
Permanent ink marker (for marking tubes)
Pipet tips for adjustable pipettors with
plugged tips (sterilized by
irradiation or autoclaving)

GeneAmp™ PCR Reaction Tubes
(0.5 mL polypropylene tubes)
Additional supplies needed for
chemiluminescent detection:

Saran Wrap®
Whatman Benchkote® (polythene backed)
X-ray film (Hyperfilm ECL or Kodak XAR5 film)

1.0 REAGENT PREPARATION

1.1 Preparation of reagents not supplied

Use reagent grade chemicals unless otherwise noted. Prepare all solutions using deionized or glass distilled water (identified below as DI H₂O). Wear gloves and follow safety recommendations provided by manufacturer for handling chemicals. Comply with any and all laws, regulations, or orders with respect to the disposal of any hazardous or toxic chemical, material, substance or waste. Store all reagents at room temperature (15° to 30°C) unless otherwise noted. Reagents are prepared as follows:

1.1.1 0.5 M EDTA (1 L). Add 186.1 g of disodium ethylenediaminetetraacetic acid dihydrate (Na₂EDTA•2H₂O) to 800 mL of DI H₂O. Stir vigorously on a magnetic stirrer. To dissolve the EDTA powder, adjust the pH to 8.0 (± 0.2) by adding approximately 20 g of NaOH pellets. Check the pH and add 5N NaOH solution if needed (for small pH adjustments, pellet addition cannot be accurately controlled). Adjust the final volume to 1 liter with DI H₂O and mix thoroughly. The solution should be autoclaved or filtered through a 0.2 µm Nalgene™ filter.

1.1.2 20% (w/v) SDS (1 L). **CAUTION!!! Wear protective mask when weighing SDS.** Slowly dissolve 200 g electrophoresis-grade (ultra pure) sodium dodecyl sulfate (SDS) in 800 mL DI H₂O. Warming (e.g., in a 37°C water bath) may be required to dissolve solids completely. Adjust the final volume to 1 liter with DI H₂O and mix thoroughly.

1.1.3 1 M Tris-HCl, pH 8.0 (1 L). Dissolve 121.1 g Tris base in 800 mL DI H₂O. Adjust to pH 8.0 (± 0.2) at room temperature by adding approximately 45 mL of concentrated HCl. Adjust the final volume to 1 liter with DI H₂O and mix thoroughly. Sterilize by autoclaving.

1.1.4 TE Buffer (10 mM Tris-HCl, 0.1 mM EDTA, pH 8.0) (1 L). Add 10 mL of 1 M Tris-HCl, pH 8.0 (Section 1.1.3) and 0.2 mL of 0.5 M EDTA (Section 1.1.1) to 990 mL DI H₂O and mix thoroughly. Dispense 100 mL aliquots and sterilize by autoclaving.

1.1.5 Citrate Buffer, pH 5.0 (1 L). Dissolve 18.4 g trisodium citrate dihydrate (Na₃C₆H₅O₇•2H₂O) in 800 mL DI H₂O. Adjust the pH to 5.0 (± 0.2) by addition of approximately 6 g of citric acid monohydrate (C₆H₈O₇•H₂O). Adjust the final volume to 1 liter with DI H₂O and mix thoroughly.

1.1.6 20X SSPE Buffer (3.6 M NaCl, 200 mM (NaH₂PO₄•H₂O), 20 mM EDTA, pH 7.4) (1 L). Dissolve 7.4 g of disodium ethylenediaminetetraacetic acid dihydrate (Na₂EDTA•2H₂O) in 800 mL DI H₂O. Adjust the pH to 8.0 (± 0.2) with 10 N NaOH solution. Add 210 g Sodium Chloride (NaCl) and 27.6 g Sodium Phosphate, monobasic, monohydrate (NaH₂PO₄•H₂O). Adjust the pH to 7.4 (± 0.2) with 10 N NaOH (about 10 mL). Adjust the final volume to 1 liter with DI H₂O and mix thoroughly.

1.1.7 Hybridization Solution (5X SSPE, 0.5% w/v SDS (1 L)). Add 250 mL 20X SSPE (Section 1.1.6) and 25 mL 20% w/v SDS (Section 1.1.2) to 725 mL DI H₂O and mix thoroughly. Hybridization Solution solids must be in solution before use; warming (e.g., in a 37°C water bath) may be required to dissolve solids completely. Preparation in a clear glass container is recommended to facilitate visual inspection for solids during warming.

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

1.1.8 Wash Solution [1.5X SSPE, 0.5% w/v SDS (2 L)]. Add 150 mL of 20 X SSPE (Section 1.1.6), and 50 mL of 20% w/v SDS (Section 1.1.2) to 1,800 mL of DI H₂O and mix thoroughly. Wash Solution solids must be in solution before use; warming (e.g., in a 37°C water bath) may be required to dissolve solids completely. Preparation in a glass container is recommended to facilitate visual inspection for solids during warming.

1.1.9 Spotting Solution [0.4N NaOH, 25 mM EDTA, 0.00008% Bromothymol Blue (25 mL)]. Add 6 mL of 5N NaOH, 3.75 mL of 0.5 M EDTA (Section 1.1.1) and 150 µL 0.04 % of Bromothymol Blue (provided in Kit) to 65 mL of DI H₂O and mix thoroughly. Spotting Solution is stable for at least three months at room temperature.

1.1.10 Pre-Wetting Solution [0.4N NaOH, 25 mM EDTA (500 mL)]. Add 40 mL of 5N NaOH, 25 mL of 0.5 M EDTA (Section 1.1.1) to 435 mL of DI H₂O and mix thoroughly.

1.1.11 3% Hydrogen Peroxide (1 mL) [For Colorimetric Detection Only]. Add 100 µL of 30% H₂O₂ to 900 µL of DI H₂O and vortex to mix. Protect from light. Store at 2° to 8°C. 3% Hydrogen Peroxide has a shelf life of approximately 4 weeks when stored at 2° to 8°C.

1.1.12 Chromogen/TMB Solution. Bring the Chromogen/TMB (TMB) to room temperature (15° to 30°C). Before opening the bottle, tap it on the lab bench to shake the TMB to the bottom of the bottle. Remove the stopper carefully to prevent loss of the powder. Slowly add 30 mL of room temperature reagent grade 100% ethanol to the bottle. **Do NOT use ethanol that has been stored in a metal container. Do NOT use 95% ethanol or other alcohols.** Recap the bottle. Seal the stopper with Parafilm. Shake in an upright position on an orbital shaker for 30 minutes or until completely dissolved. Store in bottle at 2° to 8°C and protect from rust. Under these conditions the Chromogen Solution is stable for six months after preparation.

1.2 Preparation of reagents supplied

1.2.1 Human DNA Standards. Prepare a two-fold serial dilution of the DNA Standard A (provided in Kit) in TE Buffer (Section 1.1.4) as follows:

1. Label seven 0.5 mL autoclaved GeneAmp™ PCR Reaction Tubes A through G.
2. Vortex the DNA Standard A to mix it thoroughly.
3. Transfer 120 µL of DNA Standard A into the tube labeled A.
4. Aliquot 60 µL of TE Buffer (Section 1.1.4) into each of the six remaining tubes labeled B through G.
5. Add 60 µL of DNA Standard A (tube A) to the 60 µL of TE Buffer in tube B. **Vortex to mix thoroughly.**
6. Add 60 µL of diluted DNA Standard B (tube B) to the 60 µL of TE Buffer in tube C. **Vortex to mix thoroughly.**
7. Add 60 µL of diluted DNA Standard C (tube C) to the 60 µL of TE Buffer in tube D. **Vortex to mix thoroughly.**
8. Continue the serial dilution through tube G.

9. If the dilution steps are performed as described in Section 1.2.1 above, the seven DNA Standard tubes (tubes A through G) will have the concentrations of human DNA listed in Table 1.

Table 1: DNA Standards and Concentrations

DNA Standard	Concentration (ng/ µL)	Quantity DNA per 5 µL (ng)
A	2	10
B	1	5
C	0.5	2.5
D	0.25	1.25
E	0.125	0.625
F	0.0625	0.3125
G	0.03125	0.15625

Note: Store the diluted DNA Standards at 2° to 8°C. The DNA Standards A through G are stable for at least three months at 2° to 8°C.

2.0 INFORMATION REGARDING PROTOCOLS

This entire section should be read before slot blotting/immobilization of DNA.

The Quantiblot™ Human DNA Quantitation Kit contains reagents for at least 10 hybridization reactions. Each hybridization reaction must include the following ten control samples: seven DNA Standards, the two DNA Calibrators and one blank (Spotting Solution (Section 1.1.9) only). An additional 38 samples can be spotted on the membrane, for a total of up to 48 samples per hybridization reaction. DNA Calibrators are provided as an internal control for DNA Standard performance.

3.0 SLOT BLOTTING/IMMOBILIZATION OF DNA

Note: Wear clean disposable laboratory gloves while preparing samples. Follow safety recommendations provided by manufacturer for handling chemicals. Comply with any and all laws, regulations of orders with respect to the disposal of any hazardous or toxic chemical, material, substance or waste.

Before beginning this section, assemble the required reagents, supplies and equipment as follows:

- DNA Standards (A through G)
- DNA Calibrator 1 (provided in Kit)
- DNA Calibrator 2 (provided in Kit)
- Slot Blot Apparatus
- Spotting Solution (Section 1.1.9)
- Pre-Wetting Solution (Section 1.1.10)
- Biohyne® B nylon membrane

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

3.1 Determine the number of samples to be analyzed including the seven Human DNA Standards (A through G), the DNA Calibrators 1 and 2 (provided in Kit), and the one blank (Spotting Solution only). Aliquot 150 μ L of Spotting Solution into a new 0.5 mL GeneAmp PCR Reaction Tube for each sample.

3.2 Label seven of the tubes containing 150 μ L Spotting Solution as follows: A, B, C, D, E, F, and G and label two of the tubes containing 150 μ L of Spotting Solution as follows: DNA Calibrator 1 and DNA Calibrator 2.

3.3 Vortex the seven DNA standards and the two DNA Calibrators. Add 5 μ L of each solution to the corresponding labeled tube containing 150 μ L of Spotting Solution.

Note: Sample DNA should be MgCl₂-free. See Troubleshooting Section.

3.4 Add 1 to 5 μ L of each test sample DNA to the remaining tubes containing 150 μ L of Spotting Solution.

3.5 While wearing clean gloves, cut a piece of Biohyne[®] B membrane to 11.0 cm x 7.9 cm. Cut a small notch in the upper right corner of the membrane to mark orientation. Place the membrane in the Hybridization Tray (Part No. N808-0136) containing 50 mL of Pre-Wetting Solution. Incubate at room temperature for 1 to 30 minutes.

Note: The following protocol is for use with GIBCO BRL The Convertible[™] slot blot apparatus. Refer to GIBCO BRL instructions for additional details. The vacuum source must have a pressure of at least 8 to 10 inches Hg.

3.6 Using forceps, remove the membrane from the Pre-Wetting Solution. Place the membrane on the gasket of the slot blot apparatus, then place the top plate of the slot blot apparatus on top of the membrane. Turn on the vacuum source (i.e., turn on house vacuum line or vacuum pump). Turn off the sample vacuum and turn on the clamp vacuum on the slot blot apparatus. Push down on the top plate to ensure the formation of a tight seal. Pour off the Pre-Wetting solution and rinse the Hybridization Tray thoroughly with DI H₂O.

3.7 Use a new pipette tip for each sample. Pipet each sample (approximately 155 μ L) into a different well of the slot blot apparatus. Slowly dispense each sample directly into the center of each well of the slot blot apparatus ensuring that the pipet tip is approximately 5 mm above the membrane.

3.8 After all samples have been pipetted into the wells of the slot blot apparatus, slowly turn on the sample vacuum. Leave the sample vacuum on until all of the samples have been drawn through the membrane (approximately 30 seconds). Inspect each slot that contains a sample for a uniform blue band. (If a uniform blue band is not visible, refer to the Troubleshooting Section.) Turn off the sample vacuum.

3.9 Turn off the clamp vacuum. Turn off the vacuum source. Disassemble the slot blot apparatus and remove the membrane. Proceed to Section 4 immediately. Do NOT allow the membrane to dry-out.

Note: After each use, soak the slot blot apparatus in a large volume of 0.1% SDS solution (approximately 5 to 15 minutes). Using a disposable lab towel, clean the gasket and the side of the top plate that contacts the membrane. Then rinse the slot blot apparatus with an excess of water and allow to dry at room temperature. Never use bleach.

4.0 DNA HYBRIDIZATION

The following section involves the hybridization of biotinylated QuantiBlot[™] D17Z1 Probe to DNA samples immobilized on the nylon membrane, the binding of Enzyme Conjugate:HRP-SA to the hybridized probe and a stringent wash to remove non-specifically bound probe.

Before starting the DNA Hybridization procedure, assemble the required reagents and equipment as follows:

- QuantiBlot D17Z1 Probe (provided in Kit)
- Enzyme Conjugate:HRP-SA (provided in Kit)
- Hybridization Solution (Section 1.1.7)
- Wash Solution (Section 1.1.8)
- Citrate Buffer (Section 1.1.5)
- 30% Hydrogen Peroxide
- Hybridization Tray and lid
- Hybridization Tray Retainer

Do NOT allow the membrane to dry at any point in the protocol. Minimize the time the membrane is not submerged in solution. Use the Hybridization Tray with lid for all steps.

Warm the Hybridization Solution and the Wash Solution to between 37° and 50°C in either a water bath or an incubator. All solids must be in solution before use. Mix well.

Note: Clean, disposable gloves should be worn throughout the DNA Hybridization (Section 4.0) and Detection steps (Section 5.0).

4.1 Pre-hybridization: Transfer the membrane to 100 mL of pre-warmed Hybridization Solution in the Hybridization Tray. Add 5 mL of 30% H₂O₂. Place the lid on the tray. Use the Hybridization Tray Retainer or a lead weight to keep tray from floating in the water bath. Rotate in a 50°C (\pm 1°C) water bath (50 to 60 rpm) for 15 minutes (\pm 2 minutes). Pour off the solution.

4.2 Hybridization: Add 30 mL of Hybridization Solution to the Hybridization Tray containing the membrane. Tilt the tray to one side and add 20 μ L of QuantiBlot D17Z1 Probe to the Hybridization Solution. Place the lid on the tray. Rotate in a 50°C (\pm 1°C) water bath (50 to 60 rpm) for 20 minutes (\pm 2 minutes). Pour off the solution.

4.3 Rinse the membrane briefly in 100 mL of pre-warmed Wash Solution by rocking the tray for several seconds. Pour off the solution.

Note: Step 4.4 (below) describes the addition of two alternative volumes of Enzyme Conjugate:HRP-SA. The volume of Enzyme Conjugate:HRP-SA used depends on the detection scheme chosen: colorimetric (TMB) or chemiluminescent (ECL[™]) detection. (See Introduction on page 2 for more general information on these two alternative detection schemes.)

4.4 Stringent Wash/Conjugation: Add 30 mL of the pre-warmed Wash Solution to the Hybridization Tray. Tilt the tray to one side and add the Enzyme Conjugate:HRP-SA to the 30 mL of Wash Solution as follows: for colorimetric (TMB) detection add 180 μ L of Enzyme Conjugate:HRP-SA OR for chemiluminescent (ECL) detection add 90 μ L of Enzyme Conjugate:HRP-SA. Place the lid on the tray. Rotate in a 50 C (\pm 1°C) water bath (50 to 60 rpm) for 10 minutes (\pm 1 minute). Pour off the solution.

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

5.2 Chemiluminescent Detection

Before beginning this section, assemble the required reagents, supplies and equipment as follows:

ECL™ Detection Reagents
Whatman® Benchkote
Saran Wrap®
X-ray Film
Film Cassette

Note: ECL Reagents 1 and 2 should be stored separately (at 2° to 8°C) and not allowed to cross-contaminate each other.

5.2.1 To 5 mL of ECL Reagent 2, add 5 mL of ECL Reagent 1. **Do NOT prepare this mixture more than 5 minutes before use.** Add the 10 mL of ECL reagent mixture to the membrane in the Hybridization Tray and shake for exactly 1 minute at room temperature. Pour off the solution.

Note: For maximum sensitivity, expose the membrane to X-ray film within 10 minutes of incubation in ECL reagents.

5.2.2 Cut a piece of Benchkote to approximately 12 x 16 cm. Place the damp membrane DNA-side-up on the plastic-coated side of the Benchkote. Cover the membrane with a piece of Saran Wrap that is a few centimeters larger than the Benchkote. Use a paper towel to smooth out any wrinkles or air bubbles in the Saran Wrap. Fold the Saran Wrap behind the Benchkote on the top and bottom sides. Again, use a paper towel to flatten the Saran Wrap and remove any air bubbles. Fold the Saran Wrap behind the Benchkote on the two remaining sides.

5.2.3 In a darkroom, place a piece of Hyperfilm ECL or Kodak XAR5 film in the film cassette. Carefully place the covered membrane on top of the film such that the DNA side is in contact with the film. **Do NOT move the membrane once it is placed on top of the film; movement may cause blurring of the resulting image or a "double image".** Close the film cassette. It is very important that the film is airtight, uniform contact with the covered membrane.

5.2.4 Expose the film for 15 minutes at room temperature.

5.2.5 Process the film with an automatic film processor. If an automatic film processor is not available, use Kodak GBX Developer and GBX Fixer as follows:

1. Dilute the Developer and Fixer solutions as instructed on the reagent bottles.
2. In the darkroom, incubate the film in the Developer for 90 seconds before rinsing it in DI H₂O for a few seconds. Then incubate it in the Fixer for 90 seconds.
3. Rinse the film in a continuous flow of fresh DI H₂O for 30 minutes. Hang to dry.

Note: If desired, additional X-ray film exposures can be obtained by repeating steps 5.2.3 through 5.2.5 above. Exposures longer than 15 minutes (up to several hours or overnight) will result in increased sensitivity. Preparation of additional further-diluted DNA Standards may be necessary to quantitate samples containing less than 0.15 nanograms of DNA.

4.5 Rinse the membrane thoroughly for 1 minute in 100 mL of pre-warmed Wash Solution by rocking the tray or rotating it on an orbital shaker (100 to 125 rpm) at room temperature. Pour off the solution. Rinse again for 1 minute. Pour off the solution.

4.6 Wash the membrane by adding 100 mL of pre-warmed Wash Solution to the tray. Place the lid on the tray. Rotate at room temperature on an orbital shaker (100 to 125 rpm) at room temperature for 15 minutes. Pour off the solution.

4.7 Rinse the membrane briefly in 100 mL of Citrate Buffer by rocking the tray. Pour off the solution.

5.0 DETECTION STEPS

If 180 µL of Enzyme Conjugate:HRP-SA was used in Section 4.4, use the colorimetric detection steps in Section 5.1; if 90 µL of Enzyme Conjugate:HRP-SA was used in Section 4.4, proceed to Section 5.2 for chemiluminescent detection.

5.1 Colorimetric Detection

5.1.1 Prepare the Color Development Solution not more than 10 minutes before use. Add the following reagents in the order listed to a glass flask and mix thoroughly by swirling. **Do NOT vortex.** To 30 mL of Citrate Buffer (Section 1.1.5), add 1.5 mL of Chromogen:TMB Solution (Section 1.1.12) and 30 µL of 3% H₂O₂ (Section 1.1.11).

5.1.2 Add Color Development Solution to the membrane in the tray. Cover the tray with the lid to protect the membrane from strong light.

5.1.3 Shake at room temperature on an orbital shaker (50 to 60 rpm) for 20 to 30 minutes.

5.1.4 Remove tray from shaker, pour off liquid.

5.1.5 Stop the color development by washing in deionized H₂O (100 mL). Shake for 5 to 10 minutes (50 to 60 rpm) with the lid on the tray. Repeat for a total of three washes.

5.1.6 Photograph the membrane when it is wet. Saran Wrap® may be placed over the membrane during photography to prevent it from drying out.

1. Place the wet membrane on a flat non-absorbent surface. Keep the membrane wet throughout the photographic procedure. Minimize exposure to strong light.
2. Use a Polaroid camera with Type 55 or 667 (black and white) film or Type 59 or 559 (color) film.
3. Follow Polaroid film exposure and development instructions.

Note: For black and white photography, an orange filter (Wratten 22 or 23A) will enhance contrast.

4. Following photography, the membrane may be air-dried on any hard non-absorbent surface. Protect from light and oxidizing agents (e.g., acid treated paper, bleach and nitric acid). The blue color on the membrane will fade upon drying.


Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

6.0 RESULTS INTERPRETATION

Results are interpreted by comparing the signal intensity of the DNA test sample to the signal intensity obtained for the DNA Standards. The signal intensity for a sample reflects the total amount of DNA spotted on the membrane. The seven DNA Standards represent the following quantities of DNA spotted on the membrane: 10, 5, 2.5, 1.25, 0.625, 0.3125 and 0.15625 ng (See Table 1).

The DNA Calibrators are used to provide DNA of a known concentration to verify that the DNA Standards were correctly diluted and are providing correct results for the test samples. For example, the DNA Calibrator 1 has a stock concentration of 0.7 ng/ μ L. Five μ L of this control was added to 150 μ L of Spotting Solution and the entire 155 μ L was spotted on the membrane. Thus, 3.5 ng of this sample was spotted on the membrane (0.7 ng/ μ L \times 5 μ L = 3.5 ng). The signal obtained for this control sample should have an intensity that is between the 2.5 and 5 ng DNA Standards. Likewise, the DNA Calibrator 2 should have an intensity that is between the 0.3125 and 0.625 ng DNA Standards. If not, see the Troubleshooting Section.

Figure 1: Example of Quanti-Blot™ results (using colorimetric (TMB) detection method). The seven DNA Standards represent the following quantities of DNA spotted on the membrane: A = 10 ng; B = 5 ng; C = 2.5 ng; D = 1.25 ng; E = 0.625 ng; F = 0.3125 ng and G = 0.15625 ng. DNA Calibrator 1 should have an intensity that is between DNA Standards B and C. DNA Calibrator 2 should have an intensity that is between DNA Standards E and F. Quantities for the test samples are determined by comparison of signal intensities to the DNA Standards.



The concentration of a DNA test sample is determined as follows:

1. Determine the quantity of DNA test sample spotted on the membrane by comparing its signal intensity to the intensity of the DNA Standards.
2. Divide this quantity by the volume of DNA test sample added to the Spotting Solution (typically 5 μ L of DNA test sample is added to 150 μ L of Spotting Solution).

This calculation gives DNA concentration in ng/ μ L.

PERFORMANCE CHARACTERISTICS

When the recommended protocols are followed, using the colorimetric detection method with Chromogen:TMB, the user of the QuantiBlot™ Human DNA Quantitation Kit will be able to detect and quantitate 0.15 to 10 ng of human DNA per test sample.

TROUBLESHOOTING

Observation	Possible Cause	Recommended Action
1. No signal or low sensitivity (0.15 ng DNA Standard not visible).	Use of a membrane other than Biodyne® B.	Use Biodyne B nylon membrane. Do not use membranes that have a neutral charge.
	Incorrect NaOH or EDTA concentrations in Spotting Solution.	Prepare Spotting Solution correctly (Section 1.1.9).
	Water bath temperature too high.	Water bath temperature should be 50°C (\pm 1°C).
	DNA Probe was not added at hybridization step.	Add QuantiBlot™ D17Z1 Probe.
	Enzyme conjugate was not added.	Add Enzyme Conjugate: HRP-SA at indicated step in protocol. Use 180 μ L of Enzyme Conjugate:HRP-SA for colorimetric (TMB) detection, and use 90 μ L of HRP-SA for chemiluminescent (ECL™) detection.
	Hydrogen peroxide was inactive (colorimetric detection only).	Prepare a new Color Development Solution using a fresh bottle of hydrogen peroxide.
	Presence of MgCl ₂ in the DNA sample.	Concentrations of MgCl ₂ >0.3 mM can result in reduced sensitivity. Prepare all DNA dilutions in TE Buffer. Any MgCl ₂ can be removed from samples by microdialysis using Centricon™ 100 spin units (follow manufacturer's directions).
	Cross-contamination of ECL Reagents (chemiluminescent detection method only).	ECL Reagents 1 and 2 should be stored separately and should not be mixed until 5 minutes prior to use.
	ECL Reagents stored at room temperature (chemiluminescent detection method only).	Order fresh ECL Reagents. Store ECL Reagents at 2° to 8°C.
	Membrane slipped up onto the side of the Hybridization Tray during Hybridization or Stringent Wash steps.	Reduce the rotation rate of the water bath to 50-60 rpm. Check that the membrane is fully submerged in the bottom of the Hybridization Tray before shaking.
	Membrane dried-out significantly at some point in the protocol.	Do not allow the membrane to dry at any point in the protocol.
2. Areas of low sensitivity across the membrane.		

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

Observation	Possible Cause	Recommended Action
2. Areas of low sensitivity across the membrane.	Poor contact between the film and the covered membrane during film exposure for chemiluminescent detection.	Smooth-out any wrinkles or bubbles in the Saran Wrap® covering the membrane and make sure that the film cassette provides tight, uniform contact between the covered membrane and the film. Slowly pipet the Spotting Solution directly over the center of the wells of the slot blot apparatus, with the pipet tip raised approximately 5 mm above the membrane. Turn on the sample vacuum slowly, not all at once. After being drawn through the membrane, the sample should appear as a uniform blue band on the membrane. If the entire sample is not drawn through the membrane, turn off the sample vacuum. Pipet the sample back into the pipet tip; then pipet the sample back into the well of the slot blot apparatus. Turn on the sample vacuum to draw the sample through the membrane.
3. Non-uniform signal intensity within a slot.	Bubble(s) in slot blot wells when sample was pipetted into well, or when vacuum was applied.	Prepare solutions with proper concentrations of SDS. Pre-wet the membrane in Pre-Wetting Solution prior to slot blotting. Use 180 µL of Enzyme Conjugate:HRP-SA for colorimetric detection or 90 µL of Enzyme Conjugate:HRP-SA for chemiluminescent detection. Thoroughly rinse twice, for 1 minute each, using 100 mL of pre-warmed Wash Solution at this step. These two rinse times can be extended beyond 1 minute if necessary. Immediately after each use, soak the slot blot apparatus in a large volume of 0.1% SDS solution. Never use bleach.
4. Filter background	No or low SDS in the Hybridization Solution or in the Wash Solution. Membrane was not pre-wetted prior to slot blotting. Too much Enzyme Conjugate:HRP-SA was added. Lack of thorough rinsing at Section 4.5 of the DNA Hybridization protocol.	Prepare solutions with proper concentrations of SDS. Pre-wet the membrane in Pre-Wetting Solution prior to slot blotting. Use 180 µL of Enzyme Conjugate:HRP-SA for colorimetric detection or 90 µL of Enzyme Conjugate:HRP-SA for chemiluminescent detection. Thoroughly rinse twice, for 1 minute each, using 100 mL of pre-warmed Wash Solution at this step. These two rinse times can be extended beyond 1 minute if necessary. Immediately after each use, soak the slot blot apparatus in a large volume of 0.1% SDS solution. Never use bleach.

Observation	Possible Cause	Recommended Action
5. Blurry or mis-shapen slots (chemiluminescent detection only).	Bubbles or wrinkles in Saran Wrap covering membrane during film exposure.	Carefully smooth-out wrinkles or air pockets in the Saran Wrap that covers the membrane during film exposure.
6. Blurry slots and/or double image (chemiluminescent detection only).	Membrane was moved after it was placed in contact with the film.	Carefully place the membrane directly on top of the film. Do not move the membrane once it has been placed on top of the film.
7. The DNA Calibrators do not quantitate correctly with respect to the DNA Standards.	DNA Standard serial dilutions prepared incorrectly.	Prepare two-fold serial dilutions of DNA Standard A in TE Buffer as described. Add 5 µL of each dilution to 150 µL of Spotting Solution for slot blotting.
8. Signals obtained for non-human DNA samples.*	Water bath temperature too low. SSPE concentration too high in Wash Solution.	Water bath temperature should be 50°C (±1°C). Check that the 20X SSPE solution and the Wash Solution were prepared correctly.

* DNA from primate species may give signals similar to those obtained from equivalent amounts of human DNA. In Roche Molecular Systems (RMS) laboratories, 30 ng to 300 ng quantities of non-primate DNA samples result in either no signals or signals that are less than or equal to the signal obtained for 0.15 ng of human DNA. The following non-primate DNA samples have been tested in RMS laboratories: *E. coli*, yeast, dog, cat, mouse, rat, pig, cow, chicken, fish and turkey.

REFERENCES

- Walsh, P.S., Varfaro, J., and Reynolds, R., 1982. A rapid chemiluminescent method for quantitation of human DNA. *Nucleic Acids Research* 20: 5061-5065.
- Waye, J.S. and Willard, H.F., 1986. Structure, organization, and sequence of alpha satellite DNA from human chromosome 17: evidence for evolution by unequal Crossing-Over and an ancestral pentamer repeat shared with the human X chromosome. *Molecular and Cellular Biology* 6: 3156-3165.
- Whitehead, T.P., Thorpe, G.H.G., Carter, T.J.N., Groucutt, C., and Kricka, L.J., 1983. Enhanced luminescence procedure for sensitive determination of peroxidase-labeled conjugates in immunoassay. *Nature* 305: 158-159.
- Miller, S.A., Dykes, D.D., and Polesky, H.F., 1988. A simple salting out procedure for extracting DNA from human nucleated cells. *Nucleic Acids Research* 16: 1215.

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

**ORDERING INFORMATION
PRODUCTS ARE NOT FOR USE IN DIAGNOSTIC PROCEDURES**

<p>In the United States: For PCR Technical Support call toll-free 1-800-762-4001 FAX: 203-761-2542</p> <p>To order PCR consumables and reagents call toll-free 1-800-327-3002</p> <p>To order PCR instruments or to reach your sales representative call toll-free 1-800-545-7547</p> <p>For PCR product literature call toll-free 1-800-762-4000</p>	<p>Outside of the United States: For the number of your local Perkin-Elmer sales representative call:</p> <table border="0"> <tr><td>Austria</td><td>431 1 602 31010</td></tr> <tr><td>Belgium</td><td>32 2 725 9099</td></tr> <tr><td>Bulgaria</td><td>359 2 59 1118</td></tr> <tr><td>Canada</td><td>514-737-7575</td></tr> <tr><td>Czechoslovakia</td><td>42 2 772507</td></tr> <tr><td>Denmark</td><td>45 42 81 1555</td></tr> <tr><td>Far East-Australia</td><td>61 3 560 4566</td></tr> <tr><td>Finland</td><td>358 0880 144</td></tr> <tr><td>France</td><td>33 1 3085 6363</td></tr> <tr><td>Germany</td><td>49 6150 1010</td></tr> <tr><td>Hungary</td><td>36 1 176 4011</td></tr> <tr><td>Italy</td><td>39 39 23 831</td></tr> <tr><td>Latin America-Mexico</td><td>52 5 651 7077</td></tr> <tr><td>Netherlands</td><td>31 1820 754 81</td></tr> <tr><td>Norway</td><td>47 22 88 6565</td></tr> <tr><td>Poland</td><td>48 22 33 0936</td></tr> <tr><td>Portugal</td><td>351 1 796 2172</td></tr> <tr><td>Spain</td><td>34 1 803 4210</td></tr> <tr><td>Sweden</td><td>46 87 3300 10</td></tr> <tr><td>Switzerland</td><td>41 61 692 5740</td></tr> <tr><td>United Kingdom</td><td>44 494 67 6161</td></tr> </table>	Austria	431 1 602 31010	Belgium	32 2 725 9099	Bulgaria	359 2 59 1118	Canada	514-737-7575	Czechoslovakia	42 2 772507	Denmark	45 42 81 1555	Far East-Australia	61 3 560 4566	Finland	358 0880 144	France	33 1 3085 6363	Germany	49 6150 1010	Hungary	36 1 176 4011	Italy	39 39 23 831	Latin America-Mexico	52 5 651 7077	Netherlands	31 1820 754 81	Norway	47 22 88 6565	Poland	48 22 33 0936	Portugal	351 1 796 2172	Spain	34 1 803 4210	Sweden	46 87 3300 10	Switzerland	41 61 692 5740	United Kingdom	44 494 67 6161
Austria	431 1 602 31010																																										
Belgium	32 2 725 9099																																										
Bulgaria	359 2 59 1118																																										
Canada	514-737-7575																																										
Czechoslovakia	42 2 772507																																										
Denmark	45 42 81 1555																																										
Far East-Australia	61 3 560 4566																																										
Finland	358 0880 144																																										
France	33 1 3085 6363																																										
Germany	49 6150 1010																																										
Hungary	36 1 176 4011																																										
Italy	39 39 23 831																																										
Latin America-Mexico	52 5 651 7077																																										
Netherlands	31 1820 754 81																																										
Norway	47 22 88 6565																																										
Poland	48 22 33 0936																																										
Portugal	351 1 796 2172																																										
Spain	34 1 803 4210																																										
Sweden	46 87 3300 10																																										
Switzerland	41 61 692 5740																																										
United Kingdom	44 494 67 6161																																										
<p>In Europe: For PCR Technical Support call Perkin-Elmer's European Sales Support Center in Vaterstetten, Germany at 49/8106/381-115 FAX: 49/8106/6697</p>																																											

PERKIN ELMER



Perkin-Elmer PCR reagents are developed and
manufactured by Roche Molecular Systems, Inc.
Branchburg, New Jersey, U.S.A.

Perkin-Elmer® is a registered trademark of
The Perkin-Elmer Corporation.
Roche is a registered trademark of Hoffmann-
La Roche Inc.
Quantiblot™ and AmpfliType™ are trademarks of
Roche Molecular Systems, Inc.
GeneAmp is a trademark of Roche Molecular Systems,
Inc., licensed to The Perkin-Elmer Corporation.
Hyperfilm ECL™ and ECL™ are trademarks of
Amersham International.
Biodyne® is a trademark of Pall Biosupport Corporation
Chelex® is a trademark of Bio-Rad Laboratories.
Benchkote® is a trademark of Whatman Paper Ltd.
Bellco Hot Shaker Plus™ is a trademark of Bellco
Glass, Inc.
Saran Wrap® is a trademark of Dow Chemical.
Nalgene™ is a trademark of Nalge Company.
Polaroid® is a registered trademark of Polaroid.
Centron® is a trademark of Amicon, Inc.
Kodak is a trademark of Eastman Kodak.

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT


DNA LABORATORY

PROTOCOL FOR EXTRACTION AND QUANTITATION OF CASE EVIDENCE

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

Original dated _____

Revision dated _____ approved on _____



Technical Leader, DNA Laboratory

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

DNA LABORATORY

PROTOCOL FOR EXTRACTION AND QUANTITATION OF CASE EVIDENCE


Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

Original dated Dec 1, 1997

Revision dated _____ approved on _____



Laboratory Manager



Technical Leader, DNA Laboratory

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

DNA LABORATORY

PROTOCOL FOR EXTRACTION AND QUANTITATION OF CASE EVIDENCE

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

Original dated _____

Revision dated _____ approved on _____

Laboratory Manager

[Signature]

Technical Leader, DNA Laboratory

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

DNA LABORATORY

PROTOCOL FOR EXTRACTION AND QUANTITATION OF CASE EVIDENCE

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

Original dated Dec 1, 1997

Revision dated _____ approved on _____

Richard D. Gruff
Laboratory Manager

Michael...
Technical Leader, DNA Laboratory

EXTRACTION AND QUANTITATION OF DNA FROM BIOLOGICAL MATERIALS

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

The methods delineated here assume prior training in and familiarity with forensic DNA typing.

Because not all forensic situations may be foreseen, changes may occasionally be made to these methods at the analyst's discretion in order to accommodate a particular sample.

Reagents are listed in Appendix B

A. ORGANIC EXTRACTION OF DNA FROM BIOLOGICAL MATERIALS

There may be cases in which these volumes might be altered proportionately to accommodate the size of the sample.

1. Cut a **stain** into medium-sized pieces and place the pieces into a 1.5 mL tube.
 - Biological material deposited on a **hard or bulky substrate** should be scraped or swabbed as appropriate and placed in a 1.5 mL tube.
 - Place approximately 0.5-1.0 cm of **hair root** into a 1.5 mL tube. A similar size portion of the adjacent shaft should also be extracted separately.
 - Pulp from **teeth** should be scraped out, or if necessary, the tooth may be crushed and placed directly into a 1.5 mL tube. Insoluble material should be spun out prior to phenol extraction.
 - **Bone** should be shaved, crushed or sawed and placed directly into a 1.5 mL tube. Insoluble material should be spun out prior to phenol extraction.

Samples not described here should be treated as deemed appropriate by the experience of the analyst.

2. Add 400 μ L stain extraction buffer (SEB).
Add 200 μ g Proteinase K (Pro K). (10 μ L of 20 mg/mL or 20 μ L of 10 mg/mL).
For hairs add 3 mg of DTT.
Ensure that the sample is in the liquid.
3. Incubate at 56°C for ~2 hours.
4. Add an additional 200 μ g Pro K to the tube.
5. Incubate at 56°C for another 4 to 16 hours. For freshly prepared reference bloodstains, this second incubation may be shortened to ~2 hours.

NOTE: Alternately, perform the first digestion overnight, followed by addition of more Pro K and a short incubation the following day.

6. Recover the maximum amount of clear DNA-containing digest solution by removing any substrate: first agitate, squeeze out as much liquid as possible, then "piggyback" * as necessary. If fine insolubles are present, remember to spin this digest solution and save the supernatant.

*"Piggyback": Place fabric pieces or other substrate material in a spin basket, cap and spin ~5 minutes in microfuge. Remove the cuttings and basket and recap tube. (Alternatively, punch several holes in the lid of the microfuge tube and place the pieces in the lid. Cover the lid and contents tightly with a piece of parafilm. Spin 5 minutes in microfuge. Remove the cuttings and cap and place a new cap on the tube.)

7. To this digest solution add an equal volume (typically 400 μ L) phenol/chloroform/isoamyl alcohol. (Draw this organic reagent from the middle of its container, avoiding the separate layer at the meniscus.) Gently invert the tube by hand to achieve a milky emulsion in the tube. Spin the tube to achieve phase separation and compression of the interface (typically 2-5 minutes).
8. Transfer aqueous phase (top layer) to a new tube without disturbing the interface.
9. Do additional organic extractions as in step 7 until the interface is clean and the aqueous phase clear, then do one more extraction (usually 2-3 extractions total is sufficient).

The results of organic extraction may vary depending on the condition of the sample. The analyst should use his/her judgment in continuing to process the sample. Sometimes it is desirable to perform subsequent extractions with pure chloroform or butanol, e.g. to remove residual phenol or reduce the sample volume.

10. Pipette about \sim 1.7 mL TE⁻⁴ into each labeled Centricon 100[†]. Add the sample (aq. phase from the last phenol extract) and spin 20 minutes at 1000 x g.
11. Do two additional washes with 2 mL TE⁻⁴. The middle wash spin time may be shortened to 10 minutes.
12. Collect retentate into cup by inverting Centricon and spinning 2 minutes at 500 x g. (If a Centricon 100 is not available, the sample may be precipitated with EtOH. See Appendix A).
13. Transfer retentate to microfuge tube, approximating the amount recovered. This sample may be used for both RFLP and PCR analysis.

Proceed to:

D. TEST GEL FOR ASSESSING THE QUALITY AND QUANTITY OF DNA ISOLATED FROM BODY FLUID STAIN or

E. SLOT BLOT FOR ASSESSING THE QUANTITY OF HUMAN DNA

[†]Another microfiltration device such as a Microcon filtration unit may be substituted for the Centricon 100.

B. DIFFERENTIAL EXTRACTION PROTOCOL FOR BODY FLUID/SEMEN MIXTURES

There may be cases in which these volumes might be altered to accommodate the size of the sample. At times, it may be acceptable to skip the PBS soak and begin with step 5, e.g. where hospital slides are available for evaluation or where additional material is being extracted.

1. Cut the substrate into medium-sized pieces and place the pieces into a 1.5 mL tube or spin basket bottom.
2. Incubate in 400 μ L PBS for ~1 hour at 4°C.
3. Vortex 15-30 seconds to loosen cells. Agitate substrate and remove. "Piggyback"* as necessary for maximum fluid and loose cell recovery. Save the substrate for step 5. If no "piggyback" is used, spin the fluid with loosened cells in a microfuge for ~5 minutes. Carefully remove all except ~50 μ L of the supernatant (aqueous extract) and store at -20°C.
4. Resuspend the cellular pellet in the ~50 μ L of liquid remaining, and remove 3 μ L to a microscope slide. Heat dry these slides.
5. Replace the substrate from the sample from step 3 into the tube.
Add: 400 mL stain extraction buffer.
10 μ L Pro K (20 mg/mL), or 20 μ L of (10 mg/mL).

Incubate ~1 hour at 56°C. Gentle agitation is permissible.

Meanwhile, stain slides as follows and examine microscopically:

Stain sample area with nuclear fast red for 10-15 minutes.
Rinse gently with water.
Stain with picroindigocarmine solution for ~10 seconds.
Rinse with MeOH, air dry.

6. Vortex 15-30 seconds in order to loosen sperm cells. Treat as in Step 3, to remove substrate and spin down sperm cells. Prepare separate tubes to receive e. cell fraction and to save the substrate.

* "Piggyback": Place fabric pieces or other substrate material in a spin basket, cap and spin ~5 minutes in microfuge. Remove the cuttings and basket and recap tube. (Alternatively, punch several holes in the lid of the microfuge tube and place the pieces in the lid. Cover the lid and contents tightly with a piece of parafilm. Spin 5 minutes in microfuge. Remove the cuttings and cap and place a new cap on the tube.)

7. Carefully transfer all except 50 μ L of the supernatant (epithelial cell fraction) to a separate labeled tube and reserve. The e. cell fraction may be returned to 56°C incubation while steps 8 and 9 are performed.

8. Add 400-1000 μ L PBS, gently resuspend the sperm pellet and microfuge 3-5 minutes. Carefully remove all except ~50 μ L of the supernatant and discard.

Repeat this wash step, leaving 50 μ L of the supernatant each time.

The number of washes may be varied at the analyst's discretion depending on the condition of the sample. A total of 3 washes is usually adequate.*

*In cases where a high e. cell count is observed upon initial microscopic analysis or the ratio of e. cells to sperm is high, a greater number of washes is generally desirable to remove the free e. cell DNA from the remaining sperm cell pellet. This must be balanced against the potential loss of sperm to lysis during each wash. The number of washes will generally be between 3 and 5 and the volume of each wash can vary from about 400 to 1000 μ L, depending on the circumstances described above.

9. After the last wash, mix the remaining ~50 μ L as well as possible and remove 3 μ L to a microscope slide. Heat dry these slides. Slide staining and examination should be performed as in step 5 above. If e. cells remain, steps may be taken to further digest the residual e. cells in the sperm pellet. See Appendix A.

10. Prepare a master-mix containing the following for each sperm pellet sample plus one:
400 μ L stain extraction buffer
200 μ g Pro K (20 μ L of 10 mg/mL or 10 μ L of 20 mg/mL).
3 mg fresh DTT.

Add 410 μ L or 420 μ L of this master-mix to each sperm pellet tube and incubate at 56°C ~2 hours.

11. Add another 200 μ g Pro K to each e. cell fraction tube and sperm fraction tube and continue incubating at 56°C another 4 to 16 hours.

NOTE: Alternately, perform the first digestion overnight, followed by addition of more Pro K and a short incubation the following day.

12. Continue extraction procedure from Section A, Step 7.

C. CHELEX EXTRACTION OF DNA FROM BLOODSTAINS

1. Cut a stain into medium-sized pieces and combine with 1 mL of sterile distilled water in a 1.5 mL microfuge tube. Mix gently.
2. Incubate at room temperature for 15 to 30 minutes. Mix occasionally by inversion or gentle rotation.
3. Spin in a microfuge for 2-3 minutes at 10,000 to 15,000 x g.
4. Carefully remove supernatant (all but 20 to 30 μ L or enough to cover the substrate), and discard. Leave the substrate in the tube with the pellet.
5. Add 200 μ L of 5% Chelex (w/v in sterile distilled water). Use a P-1000 pipette, stirring with the tip to evenly disperse the Chelex resin before pipetting.
6. Incubate at 56°C for 30 minutes.
7. Vortex at high speed for 5 to 10 seconds.
8. Incubate in a boiling water bath for 8 minutes.
9. Vortex at high speed for 5 to 10 seconds.
10. Spin in a microfuge for 3 minutes at maximum speed (10,000 to 15,000 x g).

Proceed to:

E. SLOT BLOT FOR ASSESSING THE QUANTITY OF HUMAN DNA

In case of extremely minute samples, quantitation may be omitted. The sample may be concentrated using either microfiltration or alcohol precipitation (See Appendix A).

APPENDIX A

ETHANOL PRECIPITATION of DNA

Add 2.5 volumes of cold absolute EtOH to the sample and mix by hand.

Incubate at -20°C for 15-30 minutes.

Spin tube in microfuge 15-30 minutes. Decant the alcohol.

Rinse pellet with an equal volume of room temperature 70% EtOH. Spin for 5 minutes. Remove as much of the EtOH as possible.

Spin tube in Speed-Vac the minimum time necessary to remove remaining alcohol. This should take about 2-5 minutes. (Alternatively, rinse with absolute EtOH and drain over paper towels for 15-30 minutes.)

Add 30 μL TE⁴ to the tube and place at 56°C for ~1 hour to dissolve the DNA. Be sure that the liquid is distributed around the tube and subsequently spun to the bottom.

OPTIONAL EXTRA DIGESTION OF EPITHELIAL CELLS IN DIFFERENTIAL EXTRACTION

In protocol "B", insert at the end of Step 9, after microscopic examination of the washed sperm pellet:

Repeat the digestion by adding an additional 400 μL of SEB and 10 μL of Pro K (20 mg/mL) to the tubes containing sperm pellets.

Incubate ~30 minutes. Microfuge 3-5 minutes and transfer all but ~50 μL of the supernatant (e. cell fraction "2") to a new tube. This e. cell fraction "2" should be reserved, but is likely to contain a mixture of e.cell and sperm DNA, so may not be useful for analysis.

Wash sperm pellet several times as described above, leaving 50 μL after the last wash.

Resuspend the sperm pellet in the 50 μL and repeat the microscopic examination to check for the presence of residual e. cells.

Continue the differential extraction protocol at step 10.

EXTRACTION OF SULFUR DYES FROM DNA ISOLATED FROM BIOLOGICAL MATERIALS ON DENIM

Bead Preparation (sufficient for ~3 samples):

1. Hydrate an ~100 μL equivalent of dry, thiopropyl-activated Sepharose 6B (Sigma T-8387 or equivalent) beads in 1 mL of sterile dH_2O for ~5 minutes at room temperature, mixing occasionally by inversion.
2. Microfuge ~2 minutes to pellet the beads and discard the supernatant. Wash three times with 1 mL of sterile dH_2O discarding the supernatant after each wash.
3. Add an equal volume of sterile dH_2O to the pelleted beads. Resuspend thoroughly and transfer a 200 μL aliquot, using a large bore pipette tip, into a separate microfuge tube for each sample. Microfuge the aliquot for 2 minutes and discard the supernatant. The volume of the hydrated beads should be ~100 μL .

Dye Extraction:

4. Add the ~400 μL aqueous phase (after the last organic extraction) to the pelleted beads (from step 3 above), incubate 15 minutes at room temperature, mixing occasionally by inversion. Note: If dye is being extracted from concentrated sample (i.e. post-Centricon), bring volume up to ~400 μL with TE^{-4} .
5. Microfuge for 5 minutes to pellet the beads. Transfer the supernatant to a clean tube making sure to avoid bead carry-over. Note: If the sample is to be Centricon filtered after dye extraction, the supernatant can be transferred directly to Centricon 100 previously loaded with ~1.5 mL TE^{-4} .
6. Wash the beads with an equal volume (~100 μL) of TE^{-4} for 5 minutes at room temperature, mixing occasionally by inversion. Microfuge for 5 minutes to pellet the beads and combine this bead-free supernatant with the supernatant from step 5.
7. If filtering by Centricon 100, continue from Section A, step 10.

Proceed to:

D. YIELD GEL FOR ASSESSING THE QUALITY AND QUANTITY OF DNA

or

E. SLOT BLOT FOR ASSESSING THE QUANTITY OF HUMAN DNA

CHELEX DIFFERENTIAL LYSIS AND EXTRACTION FROM VAGINAL SWABS AND SEMEN STAINS

Page 1

1. Dissect swab or fabric into thirds (when possible) using a clean or disposable scalpel. Use a clean cutting surface for each different sample.
2. Pipette 1 mL of sterile distilled water into a sterile 1.5 mL microcentrifuge tube. Add swab or fabric cutting.
3. Incubate at room temperature for 30 minutes.
4. Vortex for 10 seconds, or twirl for 2 minutes with toothpick to agitate the cells off the substrate.
5. Punch holes in the microcentrifuge tube lid, then use a toothpick or sterile pipette tip to move the swab or fabric to the lid. Centrifuge the sample in a microcentrifuge for 1 minute at maximum speed (10,000 to 15,000 g).
6. It is advisable not to discard the substrate until microscopic analysis (Step 8) shows that the sample contains sperm. Store swab or fabric in a sterile tube. If sperm are not visible microscopically, re-extract substrate more vigorously (Step 2-4).
7. Without disturbing the pellet, remove and discard all but 50 μL (or twice the volume of the pellet, whichever is greater) of the supernatant using a sterile Pasteur pipette or the tip of a sterile 1 mL disposable pipette, and discard. Resuspend the pellet in the remaining 50 μL by stirring it with a sterile pipette tip. **NOTE: This pellet contains epithelial cells and sperm cells, and is called the cell debris pellet.**
8. Remove about 3 μL of the resuspended sample and spot on a glass microscope slide for examination. Perform Christmas tree stain.

NOTE: If epithelial cells are detected, proceed with differential lysis procedure beginning with Step 9. If no epithelial cells are observed, the differential lysis procedure may be omitted, and the sample may be processed beginning with Step 17. If there are no sperm cells present, an option is to extract substrate more vigorously (Steps 2,3 and 4), or use Proteinase K to aid in the initial extraction. (See Below)

Optional treatment of substrate by Cal DOJ: To substrate, add 400 μL digest buffer and 10 μL of 20 mg/ml Proteinase K. Incubate 1 hour at 56°C. Loosen cells, spin using piggy-back arrangement to separate substrate from liquid and cells. Save the supernatant leaving 50 μL on the pellet. Examine 3 μL as in step 8. This pellet gets combined with the original pellet at step 13.

CHELEX DIFFERENTIAL LYSIS AND EXTRACTION FROM VAGINAL SWABS AND SEMEN STAINS

Page 2

9. Add 150 μL TE⁻⁴ to the cell debris pellet. Add 2 μL Proteinase K(10 mg/mL). Mix gently.
10. Incubate at 37°C for 1 hour to lyse the epithelial cells.
11. Spin in a microcentrifuge for 5 minutes at maximum speed (10,000 - 15,000 g). The resultant pellet is called the sperm pellet.
12. Add 150 μL of the supernatant to 50 μL of 20% Chelex in a fresh 1.5 mL microcentrifuge tube. SAVE FOR EPITHELIAL DNA ANALYSIS, BEGINNING WITH STEP 19.
13. Wash the pellet as follows: Resuspend the pellet in 0.5 mL Digest Buffer. Vortex briefly. Spin in a microcentrifuge for 5 minutes at maximum speed. Remove all but 50 μL of the supernatant and discard the supernatant.
14. Repeat wash Step 13 and additional 1 to 2 times. Additional wash steps are recommended when the ratio of sperm to epithelial cells is low.
15. Wash once with sterile distilled water as follows: Resuspend the pellet in 1 mL water. Vortex briefly. Spin in a microcentrifuge for 5 minutes at maximum speed. Remove all but 50 μL of the supernatant and discard the supernatant.
16. Resuspend the pellet by stirring with a sterile pipette tip. Remove 3 μL of the resuspended sample and spot on a glass microscope slide for examination. Perform Christmas tree stain.
17. Add 150 μL of 5% Chelex to the approximately 50 μL of resuspended sperm pellet. Add 2 μL of Proteinase K(10mg/mL) and 7 μL of 1 M DTT.
18. Vortex epithelial fractions and sperm samples for 10 seconds. Spin in a microcentrifuge for 3 seconds.
19. Incubate all fractions at 56°C for 1 hour.
20. Vortex tubes for 10 seconds.
21. Incubate samples in boiling water for 8 minutes.
22. Vortex tubes for 10 seconds.

**CHELEX DIFFERENTIAL LYSIS AND EXTRACTION FROM VAGINAL SWABS
AND SEMEN STAINS**

Page 3

23. Spin in a microcentrifuge for 3 minutes.
24. Estimate the amount of DNA in samples.
25. The samples are now ready for PCR amplification.
26. Store samples at 4°C or frozen. Prior to reuse of these samples for amplification, repeat Steps 22 and 23.

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

3/97

Adopted:

ORGANIC EXTRACTION OF WHOLE BLOOD/BLOODSTAINS

Page 1

1. Pipette 0.5 mL Digest Buffer into a sterile 1.5 mL microcentrifuge tube.

Add one of the following and mix.

- a. 10 to 50 μL whole blood
- b. 2 to 10 μL buffy coat
- c. 1 cm^2 blood stain

NOTE: *When working with bloodstains, an initial soak in 1 mL of deionized water is recommended.*

2. Add 15 μL of 10mg/mL Proteinase K solution, to a final concentration of about 0.3 mg/ml. Mix.
3. Incubate at 56°C for at least 1 hour. For evidence samples it is recommended that digestion continue for a minimum of 6 hours. Digestion may be performed overnight, but more than 24 hours is not recommended.
4. If a cloth or other substrate is present, punch several holes in the lid of the tube using a clean needle. Using a disposable pipette tip, or a fresh sterile toothpick, place the substrate in the lid and spin for 5 minutes to collect fluid remaining in the cutting. Remove the cutting, and place a new cap on the tube.
5. Add 0.5 mL buffered phenol-chloroform solution to the 0.5 mL lysed and digested cells.
6. Vortex for 15 seconds. Spin in a microcentrifuge for 5 minutes. Repeat steps 5 and 6 as necessary; aqueous phase should be clear and the interface clean. For each repeat, transfer upper layer to fresh tube.
7. (Optional) Use a sterile pipette tip to transfer the **upper aqueous** phase to a fresh sterile 1.5 mL microcentrifuge tube. Add 0.5 mL Chloroform. Vortex 10 to 15 seconds. Spin in a microcentrifuge for 5 minutes.
8. Transfer the **upper aqueous** phase to the upper chamber of Centricon tube.
9. Bring final volume to the 2 ml mark on the tube with TE^{-4} .
10. Centrifuge tubes at 1000 x g for 20 minutes in the Hermle Centrifuge.
11. Discard the liquid in the lower chamber; add 2 ml TE^{-4} to the upper chamber of the Centricon device.

ORGANIC EXTRACTION OF WHOLE BLOOD/BLOODSTAINS

Page 2

12. Centrifuge as in Step 10.
13. Repeat TE⁻⁴ washes (Steps 9-12) two times, for a total of 3 washes. After the last wash cycle, collect the approximately 40 μ L concentrated DNA sample (as per Centricon instructions) by inverting the upper reservoir into the provided retentate cup and centrifuging at 500 x g for 2 minutes to transfer concentrate into cup. Label the retentate cup. Transfer to 1.5 ml Eppendorf tube estimating volume as you transfer. Should be 15-20 μ L. (Use 5 μ L pipette, count as you transfer.) The sample is now ready for amplification. Store the sample at 2 to 8°C or freeze at -20° until ready to perform PCR.

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

3/97

Adopted:

CHELEX EXTRACTION FROM WHOLE BLOOD/BLOODSTAINS

1. Add 3 μL whole blood or a bloodstain approximately 3 mm by 3 mm to a sterile 1.5 mL microcentrifuge tube. Pipette 1 mL sterile distilled water (or PBS) into the tube. Vortex 2 seconds.
2. Incubate at room temperature for 15-30 minutes. Mix occasionally by inversion or Vortex. Vortex 5 seconds.
3. Spin in a microcentrifuge for 3 minutes at maximum speed (10,000-15,000 g).
4. Without disturbing the pellet, carefully remove and discard the supernatant leaving enough behind to cover the pellet without disturbing it. If the sample is a bloodstain, leave the fabric substrate in the tube with the pellet.
5. Add 200 μL 5% Chelex. (Use P-1000 pipette, use tip to stir Chelex before pipetting).
6. Incubate at 56°C for 30 minutes.
7. Vortex at high speed 5 to 10 seconds.
8. Incubate in a boiling water bath for 8 minutes.
9. Vortex at high speed for 5 to 10 seconds.
10. Spin in a microcentrifuge for 3 minutes at maximum speed (10,000-15,000 g).
11. Quantitate the DNA (in the supernatant).
12. The sample is now ready for the PCR amplification process. Take supernatant, avoiding Chelex beads in bottom.
13. Store the remainder of the supernatant at 2-8°C or frozen. To re-use, thaw and repeat Steps 9 and 10.

3/97
Adopted:

PCR SETUP -- PM DQA1

Turn on the thermal cycler

1. Make appropriate dilutions of samples to total 20 μL using sterile water. (Sample prep: Chelex ext: vortex and centrifuge 3 min. at max).
2. Remove tubes containing 40 μL of reaction mixture from main serology refrigerator. SPIN TO REMOVE CONDENSATION.
3. Label these tubes with sample numbers plus amplification \oplus and \ominus controls.
4. Add 40 μL of primer reagent (already containing 4% BSA) to the reaction mixture tubes, pipetting against the side to minimize mixing.
5. At this time, you have 20 minutes to complete the sample setup for amplification.
6. Add two drops of mineral oil to the tubes, cap each tube loosely.
7. Add 20 μL of sample, 20 μL of DNA for \oplus amp control, and 20 μL of sterile H_2O for neg. Amp control to the appropriately labeled tubes. (BELOW MINERAL OIL LAYER)
8. Place tubes into cardboard tray, and take into product room for amplification.
9. Place samples into wells, checking that lids are tightly on tubes. Use worksheet for documentation. Push "FILE 40" then "ENTER" and then push "START". Watch the first cycle to ensure instrument is properly functioning. The amplification takes approximately 2.5 hours. When the process is completed, the samples will be held at 4 degrees until they are removed.

Adopted: 5/97

QUANTIBLOT HUMAN DNA QUANTITATION KIT

This kit is available through Perkin-Elmer (part no. N808-0114). The procedure is based on the hybridization of a biotinylated oligonucleotide probe to DNA samples immobilized on a nylon membrane. The probe included in this kit is complementary to a primate specific alpha satellite DNA sequence at the locus D17Z1.

A. PREPARATION OF DNA STANDARDS:

1. Label seven 0.5 mL autoclaved GeneAmp PCR Reaction Tubes A through G.
2. Vortex the DNA Standard A to mix it thoroughly.
3. Transfer 120 uL of DNA Standard A into the tube labeled A.
4. Aliquot 60 uL of TE⁻⁴ Buffer into each of the six remaining tubes labeled B through G.
5. Add 60 uL of DNA Standard A (tube A) to the 60 uL of TE⁻⁴ Buffer in tube B. Vortex to mix thoroughly.
6. Add 60 uL of diluted DNA Standard B (tube B) to the 60 uL of TE⁻⁴ Buffer in tube C. Vortex to mix thoroughly.
7. Add 60 uL of diluted DNA Standard C (tube C) to the 60 uL of TE⁻⁴ Buffer in tube D. Vortex to mix thoroughly.
8. Continue the serial dilution through tube G.

The concentrations of human DNA is listed below:

DNA STANDARD	CONCENTRATION (ng/uL)	QUANTITY DNA/5uL (ng)
A	2	10
B	1	5
C	0.5	2.5
D	0.25	1.25
E	0.125	0.625
F	0.0625	0.3125
G	0.03125	0.15625

Each hybridization reaction must include the following ten control samples: seven DNA Standards, the two DNA Calibrators and one blank (spotting solution only). DNA Calibrators are provided as an internal control for DNA standard performance.

3/97

Adopted:

B. SLOT BLOTTING/IMMOBILIZATION OF DNA

1. Determine the number of samples to be analyzed including the seven Human DNA Standards (A through G), the DNA Calibrators 1 and 2 (provided in Kit), and the one blank (Spotting solution only). Aliquot 150 uL of Spotting solution into a new 0.5 mL GeneAmp PCR Reaction Tube for each sample.
2. Label seven of the tubes containing 150 uL Spotting Solution as follows: A,B,C,D,E,F and G and label two of the tubes containing 150 uL of Spotting Solution as follows: DNA Calibrator 1 and DNA Calibrator 2. Label one tube "blank."
3. Vortex the seven DNA standards and the two DNA Calibrators. Add 5 uL of each solution to the corresponding labeled tube containing 150 uL of Spotting Solution.
4. Add 1 to 5 uL of each test sample DNA to the remaining tubes containing 150 uL of Spotting Solution.
5. Cut a piece of Biodyne B membrane 11.0 x 7.9 cm. Cut a small notch in the upper right corner of the membrane to mark orientation. Place the membrane in approximately 50 mL of Pre-Wetting Solution. Incubate at room temperature for 1 to 30 minutes.
6. Place the membrane on the gasket of the slot blot apparatus, then place the top plate of the slot blot apparatus on top of the membrane. Turn on the vacuum source. Turn off the sample vacuum and turn on the clamp vacuum. Push down on the top to ensure a tight seal.
7. Pipet each sample (approximately 155 uL) into a different well of the slot blot apparatus. Slowly dispense each sample directly into the center of each well of the slot blot apparatus ensuring that the pipet tip is approximately 5 mm above the membrane.
8. After all samples have been pipetted into the wells of the slot blot apparatus, slowly turn on the sample vacuum. Leave the sample vacuum on until all of the samples have been drawn through the membrane (approximately 30 seconds). Inspect each well that contains a sample for a uniform blue band. Turn off the sample vacuum and clamp vacuum.
9. Disassemble the slot blot apparatus and remove the membrane. Proceed without delay to hybridization.

Adopted:

C. DNA HYBRIDIZATION

Warm the hybridization solution and the wash solution to between 37°C and 50°C in the water bath. All solids must be in solution before use. Mix well. Turn on rotating water bath, set to 50°C.

1. **Pre-hybridization:** Transfer the membrane to 100 mL of pre-warmed Hybridization Solution. Add 5 mL of 30% H₂O₂. Place lid on the container. Use a weight to keep tray from floating in the water bath. Rotate in a 50°C (+/-1°) water bath (50 to 60 rpm) for 15 minutes (+/-2 minutes). Pour off the solution.
2. **Hybridization:** Add 30 mL of Hybridization Solution to the tray containing the membrane. Tilt the tray to one side and add 20 uL of QuantiBlot D17Z1 Probe to the Hybridization Solution. Place lid on the tray. Rotate in a 50°C(+/-1°) water bath (50 to 60 rpm) for 20 minutes (+/- 2 minutes). Pour off solution.
3. Rinse membrane briefly in approximately 100 mL of pre-warmed Wash Solution by rocking the tray for several seconds. Pour off the solution.
4. **Stringent Wash/Conjugation:** Add 30 mL of the pre-warmed Wash Solution to the Hybridization Tray. Tilt the tray to one side and add 180 uL Enzyme Conjugate: HRP-SA to the 30 mL of Wash Solution. Place lid on tray. Rotate in a 50°C(+/- 1°) water bath (50 to 60 rpm) for 10 minutes (+/- 1 minute). Pour off the solution.
5. Rinse the membrane thoroughly for 1 minute in approximately 100 mL of pre-warmed Wash Solution by rocking the tray or rotating it on an orbital shaker (100-125 rpm) at room temperature. Pour off the solution. Rinse again for 1 minute. Pour off the solution.
6. Wash the membrane by adding approximately 100 mL of pre-warmed Wash Solution to the tray. Place the lid on the tray. Rotate at room temperature on an orbital shaker (100 to 125 rpm) for 15 minutes. Pour off the solution.
7. Rinse the membrane briefly in 100 mL of Citrate Buffer by rocking the tray. Pour off the solution.

D. DETECTION STEPS

1. Prepare the Color Development Solution no more than 10 minutes before use. Add the following reagents in the order listed to a glass flask and mix thoroughly by swirling. Do not vortex. To 30 mL of Citrate Buffer add 1.5 mL of Chromogen: TMB Solution and 30 uL of 3% H₂O₂ (or 3 uL of 30% H₂O₂).
2. Add Color Development Solution to the membrane in the tray. Cover the tray with lid to protect the membrane from strong light.
3. Shake at room temperature on an orbital shaker (50 to 60 rpm) for 20 to 30 minutes.
4. Remove tray from shaker, pour off liquid.
5. Stop the color development by washing in deionized water (100 mL). Shake for 5 to 10 minutes (50 to 60 rpm) with the lid on the tray. Repeat for a total of three washes.
6. Photograph the membrane when it is wet.

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

PM + DQA1 PCR AMPLIFICATION PROTOCOL

This analysis is performed using the Amplitype® PM + DQA1 Amplification and Typing Kit from Perkin-Elmer. Using a dedicated pipetter and working in the biological hood, put aliquots of 40µL of Reaction Mix from the kit into GeneAmp® tubes. Store in the refrigerator.

(Turn on the Thermal Cycler. Work in the biological hood to set up the amplification tubes.)

1. Make appropriate dilutions of samples to total 20 µL using sterile water. (Sample prep: for Chelex extract vortex and centrifuge 3 min. at max g). The samples will typically contain 1-10 ng of DNA.
2. Remove tubes containing 40 µL of reaction mixture from the refrigerator. Pulse spin to ensure all liquid contents in bottom of tube.
3. Label these tubes with sample numbers plus amplification \oplus and \ominus controls.
4. Prepare master mix for all tubes plus one :
40 µL of primer reagent
0.4 µL of 4 % BSA
Add 40 µL of this master-mix to each reaction mixture tube, pipetting against the side to minimize mixing.
5. At this time, the amplification set-up must be completed within 20 minutes.
6. Add two drops of mineral oil to the tubes, cap each tube loosely.
7. Add 20 µL of each sample, 20 µL of kit DNA std. for \oplus amp control, and 20 µL of sterile H₂O for \ominus amp control to the appropriately labeled tubes. (BELOW MINERAL OIL)
8. Place tubes into cardboard tray, and take into product room for amplification.
9. Place samples into wells, checking that lids are tightly on tubes. Use worksheet for documentation. Start the amplification program: Push "FILE 40" then "ENTER" and then push "START". Watch the first cycle to ensure instrument is properly functioning. The recommended amplification conditions for this kit are 32 cycles as follows:
94°C denaturation for 1 minute
60°C annealing for 30 seconds
72 °C extension for 30 seconds

After the last cycle, extend at 72°C for 7 minutes, then hold at 4°C

The amplification takes approximately 2.5 hours. When the process is completed, the samples will be held at 4 degrees until they are removed.

Adopted: 10/97

PRODUCT GEL PROTOCOL

- 1) Remove 5 μ L amplified DNA and combine with 2 μ L loading buffer. The loading buffer is placed into the tube first. Mix the sample and loading buffer with the pipet tip.
- 2) Add 5 μ L of 100ng/ μ L GIBCO BRL 123 ladder to 2 μ L of loading buffer. (To prepare a 100ng/ μ L ladder, dilute the 1ng/ μ L GIBCO BRL product 10 μ L + 90 μ L DI H₂O).
- 3) Prepare a 4% gel
 - a) Follow directions in reagent preparation book to prepare 4% Nu Sieve agarose gel in 1X TBE buffer. Add EtBr: 5 μ L of 10 mg/mL per 100 mL of buffer. Agarose is dispensed into tubes containing 30-35ml each (one tube is used for each product gel run).
 - b) Melt agarose in microwave on defrost for 2-3 min., checking frequently.
 - c) Pour a 0.45 cm. thick gel with 8 tooth comb at one end. Allow to solidify about ½ hr.
- 4) Pour enough .5X TBE mixed with ethidium bromide into the tank to cover gel to about 2 mm, typically about 250 mL. (50 μ L of 10mg/ml ethidium bromide in 1L of 0.5X TBE).
- 5) Load the 123 bp ladder and the samples from step 1.
- 6) Set the voltage at 115v and run for 45-60 minutes or until the dye has gone 2 cm. Run towards the anode (red).
- 7) Photograph the gel. (Six bands should be present in the control DNA lane and in samples in which the DNA was amplified successfully. If all six bands are present, proceed to DNA hybridization.)

Adopted: 11/97

PRODUCT GEL

- 1) Remove $5\mu\text{l}$ amplified DNA and combine with $2\mu\text{l}$ loading buffer. The loading buffer is placed into the tube first. Mix the sample and loading buffer with the pipet tip.
- 2) Add $5\mu\text{l}$ of the $100\text{ng}/\mu\text{l}$ GIBCO BRL 123 ladder to the $2\mu\text{l}$ of loading buffer. (To prepare a $100\text{ng}/\mu\text{l}$ ladder, dilute the $1\text{ng}/\mu\text{l}$ GIBCO BRL product $10\mu\text{l} + 90\mu\text{l}$ DI H_2O).
- 3) Prepare a gel
 - a) Follow directions in reagent preparation book to prepare product gel. Agarose dispensed into tubes containing a 30-35ml each (one tube is used for each product gel run).
 - b) Heat on defrost for 2-3 min. in microwave, checking frequently.
 - c) Pour a 0.45 cm. thick gel with 8 tooth comb at one end. Allow to solidify about $\frac{1}{2}$ hr.
- 4) Place enough .5X TBE mixed with ethidium bromide into the tank to cover gel to about 2 mm ($50\mu\text{l}$ of $10\text{mg}/\text{ml}$ ethidium bromide in 1L of 0.5X TBE). Use 250mL.
- 5) Load the 123 bp ladder in first well and then the samples. $7\mu\text{l}$ or $7.5\mu\text{l}$ on P20.
- 6) Set the voltage at 115v and run for 45-60 minutes or until dye has gone 2 cm. (Alternative parameters from Washoe Co.: 135v for 30 min.) Run towards red - origin on right.
- 7) Photograph the gel. (Six bands should be present in the control DNA lane and in samples in which the DNA was amplified successfully. If all six bands are present, proceed to DNA hybridization.) Suggested settings: 100% UV lamp intensity; f 5.6, $1/8$ sec

Check mfg. handbook from kit for specific instructions on preparing reagents, etc.

5/97
Adopted:

PRODUCT GEL PROTOCOL

- 1) Remove 5 μ L amplified DNA and combine with 2 μ L loading buffer. Mix the sample and loading buffer with the pipet tip.
- 2) Add 5 μ L of 100ng/ μ L GIBCO BRL 123 ladder to 2 μ L of loading buffer. (To prepare a 100ng/ μ L ladder, dilute the 1 μ g/ μ L GIBCO BRL ladder 10 μ L + 90 μ L DI H₂O).
- 3) Prepare a 4% gel
 - a) Prepare 4% NuSieve 3:1 agarose gel in 1X TBE buffer.
 - b) Heat to dissolve agarose in microwave on defrost for 2-3 min., checking frequently. Add EtBr: 5 μ L of 10 mg/mL per 100 mL of buffer.
 - c) Pour a 0.45 cm. thick gel with 8 tooth comb at one end. Allow to solidify about ½ hr.
- 4) Remove the comb and pour enough 0.5X TBE into the tank to cover the gel to about 2 mm, typically about 250 mL.
- 5) Load the 123 bp ladder and the samples from step 1.
- 6) Set the voltage at 115v and run for 45-60 minutes or until the dye has gone at least 2 cm. Run towards the anode (red).
- 7) Photograph the gel. (Six bands should be present in the control DNA lane and in samples in which the DNA was amplified successfully. If all six bands are present, proceed to DNA hybridization.)

Adopted: 11/97

E. SLOT BLOT FOR ASSESSING THE QUANTITY OF HUMAN DNA USING QUANTITBLOT KIT

PREPARATION OF DNA STANDARDS:

1. Label seven 0.5 mL autoclaved microcentrifuge tubes A through G.
2. Vortex the DNA Standard A to mix it thoroughly, then spin down all fluid.
3. Transfer 20 μL of DNA Standard A into the tube labeled A.
4. Aliquot 10 μL of TE⁺ Buffer into each of the six remaining tubes labeled B through G.
5. Add 10 μL of DNA Standard A (tube A) to the 10 μL of TE⁺ Buffer in tube B. Mix thoroughly.
6. Add 10 μL of diluted DNA Standard B (tube B) to the 10 μL of TE⁺ Buffer in tube C. Mix thoroughly.
7. Continue the serial dilution through tube G.
8. Dilute the K562 human cell line standard using TE⁺ to 100 ng/ μL for storage as a stock solution. For immediate use, dilute in two steps to 0.2 ng/ μL (e.g. 5:100, then 4:100)

The concentrations of human DNA are listed below:

DNA STANDARD	CONCENTRATION (ng/ μL)	QUANTITY DNA/5 μL (ng)
A	2	10
B	1	5
C	0.5	2.5
D	0.25	1.25
E	0.12	0.6
F	0.06	0.3
G	0.03	0.15
K562 cell line	0.2	1.0

Each hybridization reaction must include the following eleven control samples: seven DNA Standards, the two DNA Calibrators, a dilution of K562 cell line standard and one blank (spotting solution only). DNA Calibrators are provided as an internal control for the dilution series of the DNA standard. The K562 solution is an external control on the Quantiblot kit.

SLOT BLOTTING/IMMOBILIZATION OF DNA

1. Determine the number of samples to be analyzed including the seven human DNA standards (A through G), the DNA calibrators 1 and 2 (provided in kit), the K562 standard, and the one blank (spotting solution only). Label sterile 0.5mL centrifuge tubes for these eight standards, two calibrators and blank. Label additional tubes for the samples to be quantitated.
2. To every tube add 150 μ L of spotting solution.
3. Mix the seven DNA standards and the two DNA Calibrators. Spin. Add 5 μ L of each solution to the corresponding labeled tube containing 150 μ L of spotting solution. Add 5 μ L of K562 solution to its tube.
4. Add 1 to 5 μ L of each test sample DNA to the remaining tubes containing 150 μ L of spotting solution.
5. Cut a piece of Biodyne B membrane 11.0 x 7.9 cm. Cut a small notch in the upper right corner of the membrane to mark orientation. Place the membrane in approximately 50 mL of pre-wetting solution. Incubate at room temperature for 1 to 30 minutes.
6. Place the membrane on the gasket of the slot blot apparatus, then place the top plate of the slot blot apparatus on top of the membrane. Turn on the vacuum source. Turn off the sample vacuum and turn on the clamp vacuum. Push down on the top to ensure a tight seal.
7. Pipet each sample (approximately 155 μ L) into a different well of the slot blot apparatus. Slowly dispense each sample directly into the center of each well of the slot blot apparatus ensuring that the pipet tip is approximately 5 mm above the membrane.
8. After all samples have been pipetted into the wells of the slot blot apparatus, slowly turn on the sample vacuum. Leave the sample vacuum on until all of the samples have been drawn through the membrane (approximately 30 seconds). Inspect each well that contains a sample for a uniform blue band. Turn off the sample vacuum and clamp vacuum.
9. Disassemble the slot blot apparatus and remove the membrane. Proceed without delay to hybridization.

DNA HYBRIDIZATION

Warm the hybridization solution and the wash solution to between 37°C and 50°C in the water bath. All solids must be in solution before use. Mix well. Turn on rotating water bath, set to 50°C.

1. **Pre-hybridization:** Transfer the membrane to 100 mL of pre-warmed hybridization solution. Add 5 mL of 30% H₂O₂. Place lid on the container. Use a weight to keep tray from floating in the water bath. Rotate in a 50°C (+/-1°) water bath (50 to 60 rpm) for 15 minutes (+/-2 minutes). Pour off the solution.
2. **Hybridization:** Add 30 mL of hybridization solution to the tray containing the membrane. Tilt the tray to one side and add 20 µL of QuantiBlot D17Z1 Probe to the hybridization solution. Place lid on the tray. Rotate in a 50°C(+/-1°) water bath (50 to 60 rpm) for 20 minutes (+/- 2 minutes). Pour off solution.
3. Rinse membrane briefly in approximately 100 mL of pre-warmed wash solution by rocking the tray for several seconds. Pour off the solution.
4. **Stringent Wash/Conjugation:** Add 30 mL of the pre-warmed wash solution to the hybridization tray. Tilt the tray to one side and add 180 µL Enzyme Conjugate: HRP-SA to the 30 mL of wash solution. Place lid on tray. Rotate in a 50°C(+/- 1°) water bath (50 to 60 rpm) for 10 minutes (+/- 1 minute). Pour off the solution.
5. Rinse the membrane thoroughly for 1 minute in approximately 100 mL of pre-warmed wash solution by rocking the tray or rotating it on an orbital shaker (100-125 rpm) at room temperature. Pour off the solution. Rinse again for 1 minute. Pour off the solution.
6. Wash the membrane by adding approximately 100 mL of pre-warmed wash solution to the tray. Place the lid on the tray. Rotate at room temperature on an orbital shaker (100 to 125 rpm) for 15 minutes. Pour off the solution.
7. Rinse the membrane briefly in 100 mL of citrate buffer by rocking the tray. Pour off the solution.

D. DETECTION STEPS

1. Prepare the color development solution no more than 10 minutes before use. Add the following reagents in the order listed to a glass flask and mix thoroughly by swirling. Do not vortex. To 30 mL of citrate buffer add 1.5 mL of Chromogen: TMB Solution and 30 μ L of 3% H_2O_2 (or 3 μ L of 30% H_2O_2).
2. Add color development solution to the membrane in the tray. Cover the tray with lid to protect the membrane from strong light.
3. Shake at room temperature on an orbital shaker (50 to 60 rpm) for 20 to 30 minutes.
4. Remove tray from shaker, pour off liquid.
5. Stop the color development by washing in deionized water (100 mL). Shake for 5 to 10 minutes (50 to 60 rpm) with the lid on the tray. Repeat for a total of three washes.
6. Photograph the membrane when it is wet.

Property of Idaho State Police Forensic Services
Uncontrolled Interim Copy
OBSOLETE DOCUMENT

DNA HYBRIDIZATION: PM + DQA1

The AmpliType PM + DQA1 DNA Hybridization process involves 3 steps performed sequentially as follows: 1) hybridization of amplified DNA to DNA probe strips, 2) binding of HRP-SA to hybridized PCR products, and 3) stringent wash to remove non-specifically bound PCR products. Color development follows the stringent wash step.

N.B. Ensure that 5ul of 200mM disodium EDTA is added to each tube before hybridization step.

Perform these steps in the Amplified DNA Work Areas using reagents and pipets dedicated for use in this area. Do not transport any of these items to the Extraction/PCR Set Up Work Area.

1. Heat a rotating water bath to 55° C. The water level should be 0.5 to 1 cm above the platform. The water level should not be higher than 1 cm, as higher levels may result in water splashing into the wells. The temperature should not go below 54° C or above 56° C. It is essential to check the temperature before each hybridization is performed. Set the bath for 50 to 70 rpm.
2. Warm (37° - 55° C) the Hybridization solution and the Wash Solution (if necessary). All solids must be in solution before use.
3. Allow the tube strips to come to room temperature before opening it. Remove the strips using flat forceps. With a waterproof marking pen, label each probe strip with the designation of the corresponding amplified DNA.
4. Place the tubes of amplified DNA in the thermal cycler after it reaches 95° C. Denature the DNA by incubation for 3 to 10 minutes. Keep each tube at 95° C until use.
5. Add 3 mL of pre-warmed Hybridization Solution to each well at the labeled end of each strip while tilted. **DO NOT WET THE REMAINDER OF THE STRIP.**
6. Perform the following steps for each tube of amplified DNA:
 - a. Remove the tube from the 95° C block, open using a tube opener.
 - b. Withdraw 20 uL of amplified DNA and immediately (within 20 seconds) add it below the surface of the hybridization solution in the well of the corresponding probe strip.
 - c. Add only one sample of amplified DNA at a time. Replace in 95° C block.
 - d. If using 2 strips with the same amplified DNA sample, replace the tube in the 95° block after taking the first 20 µL.

The remaining amplified DNA samples can be stored at 2° to 8° C for two months or at -20° C for 6 months.

7. Transfer the tray to the 55° C water bath. Place a weight on the cover. Rotate the water bath 50 to 70 rpm for 15 minutes (+/- 2 minutes). Record temperature and time of hybe step.

DNA HYBRIDIZATION: PM + DQA1 (continued)

PAGE 2

8. Approximately 5 minutes before the end of the hybridization step, prepare the Enzyme Conjugate Solution as follows:

Number of strips X 3.3 mL Hybridization Soln = Volume of Hybridization Soln

Number of strips X 27 uL Enzyme Conjugate = Volume of Enzyme Conjugate

Mix the solution thoroughly and ensure that the solids remain in solution. Leave at room temperature until ready to use. Remove the tray from the water bath, leave the water bath covered and rotating.

9. Aspirate the contents of each well from the labeled end of the strip. Wipe the lid with a clean lab wipe.
10. Dispense 5 mL of pre-warmed Wash Solution into each well. Rock the tray for several seconds with the lid on, then aspirate each well.
11. Dispense 3 mL of the Enzyme Conjugate Solution prepared previously into each well. Transfer to the 55° C water bath and rotate at 50 to 70 rpm for 5 minutes (+\ 1 minute).
12. Aspirate the contents of each well. Wipe the lid with a clean lab wipe. Keep water bath rotating and covered.
13. Dispense 5 mL of pre-warmed Wash Solution into each well. Rinse by rocking the tray for several seconds with the lid on then aspirate the solution from each well.
14. Dispense 5 mL of pre-warmed Wash Solution into each well. Place in the 55° C water bath and rotate at 50 to 70 rpm for 12 minutes (+\ 1 minute). Record the time and temperature of this stringent wash step.
15. After incubation, aspirate the contents out of each well. Wipe the lid with a clean lab wipe.
16. Dispense 5 mL of Wash Solution into each well. Cover and rock the tray gently for several seconds.
17. Aspirate the contents of each well. Wipe the lid with a clean lab wipe.

5/97
Adopted:

COLOR DEVELOPMENT OF STRIPS: PM + DQA1

1. Dispense 5 mL of Citrate buffer into each well. Place on an orbital shaker set at approximately 50 rpm at room temperature for 5 minutes.
2. During this wash step prepare the Color Development Solution as follows:

Number of strips X 5 mL Citrate Buffer = Volume of Citrate Buffer

Number of strips X 0.5 uL 30% H₂O₂ = Volume of Hydrogen Peroxide

Number of strips X 0.25 mL Chromogen = Volume of Chromogen

3. Remove the tray from the orbital shaker and aspirate the contents of each well. Add 5 mL of freshly prepared Color Development Solution to each well.

NOTE: PLACE THE CLEAR PLASTIC LID ON THE TRAY AND COVER THE LID WITH ALUMINUM FOIL DURING STEPS 4, 5 & 6.

4. Develop the strips at room temperature by rotating on an orbital shaker set at approximately 50 rpm for 20 to 30 minutes. Develop until the "C" or "S" dot is visible. Individual strips may develop before 20 minutes and be stopped individually (Step 5).
5. Aspirate the contents of each well. Dispense approximately 5 mL DI water into each well. Place tray on orbital shaker at about 50 rpm for 5 to 10 minutes. Pour off water.
6. Repeat step 5 at least two times for a minimum of three DI water washes. Turn off aspirator, release VAC.
7. Record the pattern of the blue dots from each wet strip prior to photography.

NOTE: KEEP STRIPS WET THROUGHOUT PHOTOGRAPHY.

Adopted: 5/97

STRIP PHOTOS

Use the $\frac{1}{2}$:1 frame lens assembly (the one with a square frame and 3 tripod feet). If necessary, plug the flash unit power cord into the side of the lens assembly, and attach the upper section (the part with the fat handle and film holder). Use the twist button marked "lock/open" to unlock the upper section and to secure it if you have to move it between lens units.

Make sure type 665 film is loaded. The white pull tab is marked.

Confirm settings of f22 and 125th of a second. Filter is already in place.

Turn flash power unit on, setting control to "+1/2".

Arrange strips on black background, blot off any standing water. Click shutter, let develop 30 sec, examine print. If OK*, peel away extra paper from negative and drop it in water. Coat positive print, making sure coating coverage is complete.

*If not OK, adjust for more or less light with flash unit.

Process negative as soon as feasible.

D. YIELD GEL FOR ASSESSING THE QUALITY AND QUANTITY OF DNA

A yield gel may be used by itself or in concert with other methods in order to assess the quality and quantity of DNA. DNA at this stage has been through a microfiltration device such as a Centricon 100 or precipitated by EtOH and is in a volume of ~30 μL .

1. Remove 2-3 μL of the DNA and combine with 1 μL bromophenol blue (B ϕ B) loading buffer solution. Or if using a commercial yield gel set of standards, dilute according to package insert.
2. Preparation of test gel: The size used and number of origins depends on the number of specimens that one needs to test. Refer to Appendix B for the appropriate volume of agarose to achieve the appropriate gel thickness for each size tray.

All gels use 1% LE agarose in 1X TBE buffer supplemented with 0.5 $\mu\text{g}/\text{mL}$ ethidium bromide (5 μL of 10 mg/mL EtBr per 100 mL of buffer). Alternatively, EtBr may be omitted and the gel stained after the run in a staining tray in the same concentration of EtBr.

Prepare the appropriate volume of 1X TBE buffer.

Weigh out the appropriate amount of LE agarose into a flask or bottle.

Add the 1X TBE.

Heat to dissolve agarose.

Add EtBr (if using).

Place an appropriate size comb into the gel tray.

When the temperature is ~60°C, pour agarose into gel tray. Let stand until solid.

3. Pour 1X TBE buffer into electrophoresis tank. (This buffer does not contain EtBr).
4. Place the gel into the tank with the well comb at the cathodic end. Enough buffer should be present to cover the gel to a depth of ~3 mm (refer to Appendix B). Remove the comb.
5. Each DNA sample or standard mixed with loading solution is pipetted into a separate well with the gel submerged. Well #1 is defined as the well at the upper left side of the gel. Standards to be included once on every origin include 150-300 ng $\lambda\text{Hind III}$ (previously heated to 65°C for 5 minutes) and 6 μL each of a range of quantitation standards typically spanning 250-6 ng.
6. Run at constant voltage of 100 volts. When the blue tracking dye has moved ~2 cm from the origin, the run can be stopped.
7. Take a photograph of the gel on a UV transilluminator.
8. From the photograph, visually assess the quality and quantity of DNA in test specimens by comparison with the DNA standards.

PROCEDURE MANUAL

- I. Extraction and Quantitation
 - Section A Organic Extraction
 - Section B Differential Extraction
 - Section C Chelex Extraction
 - Section D Yield Gel
 - Section E Slot Blot

 - Appendix A
 - Ethanol Precipitation
 - Optional Extra Digestion of Epithelial Cells
 - Extraction of Sulfa Dyes
 - Appendix B
 - Extraction Reagents
 - Quantitation Reagents
- II. Amplification
- III. Strip Development
- IV. Product Gel
- V. Former Protocols
- VI. Appendix C
 - Product Inserts PM+DQA1 Kit and Quantiblot Kit