



STATE OF IDAHO

DEPARTMENT OF HEALTH AND WELFARE

BUREAU OF LABORATORIES, 2220 Old Penitentiary Road, Boise, ID 83702

Rick Groff
MAR 31 1978
RECEIVED
DEPARTMENT OF HEALTH AND WELFARE
BUREAU OF LABORATORIES

I.R. methodology for methamphetamine and phentermine used in the central lab.

1. ^{making HCl salt} Acid base extraction method: The methamphetamine or phentermine is extracted from a basic aqueous solution into CH_2Cl_2 and this solution can then be dried through a sodium sulfate column in a pasteur pipette. Add a very small drop of concentrated HCl in .5 to 1 milliliter of methanol, or add HCl in methanol to the CH_2Cl_2 drop by drop, until pH 2, or bubble HCl vapor through the solution until you are satisfied all the sample has been converted to the hydrochloride salt. Reduce the solution to a small volume, add more CH_2Cl_2 and reduce in volume again. The purpose of this step is to get rid of water, HCl, and methanol. If when you add petroleum ether, two phases are formed, go through the CH_2Cl_2 volume reduction phase again. Recrystallize with petroleum ether by adding it to the CH_2Cl_2 solution when reduced to a small volume. If necessary reduce to a small volume and add more petroleum ether. Let the crystals grow before decanting or filtering for I.R. This procedure may give different scans for d-methamphetamine and dl-methamphetamine. Be prepared for scans inbetween. Leaving in beam for a couple of hours can enhance scan occasionally.
2. Use the hexane - H_2O method.
3. Extract methamphetamine HCl directly with CH_2Cl_2 , filter through sodium sulfate, and recrystallize with petroleum ether.
4. A celite column method as outlined in the white DEA manual the Analysis of Drugs or in "Minnibennies Purification for IR" separation IV. This will separate methamphetamine from amphetamine and phentermine.
5. Chloroplatinic derivative of methamphetamine.
(To clean up where other methods don't work or to separate from phentermine).
 - a. Shake out into organic solvent and dry down to decrease sample size.
 - b. Add five drops 0.1N HCl.
 - c. Add 3-4 drops 5% aqueous platinum chloride add enough to form all the precipitate possible.
 - d. After crystals have formed filter (with suction).
 - e. Wash crystals with three or four drops 0.1NHCl then two to three drops methanol (be careful not to dissolve crystals).
 - f. Run on KBr.
 - g. To regenerate methamphetamine extract (if you don't like to run complexes) add several drops of NH_4OH extract methamphetamine

I

I. Introduction: why this conference.

1. each lab tends to go its own way without feedback to other labs - learn techniques aren't always shared
2. source manual of commonly used techniques and procedures
3. since each case varies apply methods where they seem to be useable.
4. eventually we could have a method manual that covers all our routine analysis.

IV Discussion of topics not related to subject at hand (very brief) not often all labs doing things linked up at one time

II Discussion of amphetamine a. Boise
b. Reno
c. Pocatello

Points to include methodology, parameters, solvent systems, TLC, development, etc.

III Discussion of methamphetamine

II

Beise -
Amphetamine

1. Quick report & included so you'll know reasoning behind our analysis. Saves us a lot of time on reports on analysis - less
no type written report
no quant unless needed for a special reason.
- what quant & why
2. For our quick report we tend to do things which are reasonable quick, don't require prior separation and almost always work. Gives strong probability of nature of substances - sufficient for preliminary hearings.
3. For our final analysis we do an I.R. if feasible within context of the situation. If we can't get an I.R. except for NH4LSD - special procedures - we go for a variety of tests.
4. Now you can why our analysis is split into two parts. Following our procedure is ST-DIA procedure from about 1972 so some things have certainly changed.

1. Spot tests: many
a. colorimetric
b. colorimetric - why
c. Froede -

2. T.C very straightforward used Fluor.
Mercuric yellow from amphetamine works well - needs basic pH - Ti works well with it
b. simple
c. doesn't interfere with other sprays. a. word about sources on Fluoroscamine

- (1) article from Fred Cherna - general background doesn't fluor. with many of the things it gives high sensitivity for general screening only
- (2) article from Joan et al Forensic Science ignore 5mg/10ml - works first at 5mg/100ml
(b) lists 325 drugs tells which it reacts with & causes fluorescence such as procaine



3. GLC with Apiezon L + KOH

a. originally wanted column for quant of amphetamine - amaly tested amphetamine column ran at too low temp as over acylation is often incomplete for amphetamine. KOH reduces tailing, can run directly at 150°C.

b. found to be such a good column in terms of sensitivity, reproducibility, & separation. ~~So~~ use it for routine qualitative analysis. Other labs with a smaller number of GLC may not have room for it.

end of quick report 908 to 95%
of all amphetamine analysis

First Comment

1. Most of us would then try for an I.R. - sample size providing. ~~proportion~~ - I Ran mandelic salt

1. extract group into methylene
we use this method 2. may want to filter & then wash with CH₂Cl₂
it because ~~it~~ can't eliminate caffeine - shouldn't be necessary but has
differentiated, & d.l. helped with crystals - washed with caffeine can also
it because it works go through hexane - H₂O wash.

well

3. make base, extract with CH₂Cl₂ - through Mandelic column

separate crystals for IR
reduce volume add mandelic acid in CH₂Cl₂
not too much, it suppresses formation of crystals.

4. tricks

using pH

1. reduce solution to volume of mandelic acid in CH₂Cl₂ if only small amount amphetamine present

2. catch crystals directly to be where only small amount of crystals

3. add H₂O to increase number of crystals

4. Purify with 1.5 NaOH - Hexane - H₂O

5. Some other amides will crystallize with mandelic acid

6. ordered R mandelic acid by accident seems purer than any of mandelic acid I use it.

5. I report as amphetamine not dl amphetamine

IV

2. I.R. by extracting amphetamine from 1.5N NaOH into hexane which is then washed 3x which distilled H₂O eliminates epinephrine & caffeine only works with but does not separate amphetamine, methamphetamine, phentermine, etc.

papers on this are quite self explanatory, note there are two papers the second is a modification of the first I recommend that you follow the second.

3. Minibombs Purification for F.M. - Ann Bredby will discuss

4. Vol. test I copied H₂OAC - do both other labs have Fulton? If you don't should have. Goes into same microcrystalline tests.

How we do volatility tests

1. put sample in vial in spot plate

2. add 5% NaOH

3. put clip vial on microscope slide

4. invert over heat to catch vapors in spot plate

5. turn over to view in microscope.

so many put tests excellent for commercial prep. on cross taps. volatile impurities so good. Ann B & I have had some luck with high interfere chicken tracks produced by Phthalic chloride in H₂O very sensitive not all amphetamine such similar crystals.

6. Finally more GC: short. Sm. and amphetamine. Short sample + 20% acetic anhydride in pyridine, DMF, or ethyl acetate. Pyridine & DMF supposed to enhance acylation but ethyl acetate works fine not such broad solvent front. Pure acetic anhydride can change characteristics of column momentarily.

6. 20% DMF at 180°C. Only if you happen to have this column set up which we do for another purpose

conclusion: for final reports we prefer IR if possible, if not possible, I prefer 3 different techniques in combination, TLC, GLC, & microcrystalline test because each one should test different properties.

Add more columns still may be measuring same properties same for TLC. TLC - GLC combinations much more powerful but things similar enough to give some GLC on two columns may give some TLC also. Probably a different microcrystalline test. Entirely different approach. I have used microcrystalline using GLC & TLC only.

run through DEA method for Cocaine - amph.

Wey - amph.

Methamph.

1. spot test: Wagner's
Fry's & methamph test

2. TLC : Use thamp. reacts best with Dragendorff's may take up to one hour

2. Phentermine - very weakly reactive if at all

3. GLC Apizzen L + KUH - already covered

Diagrams 4. IR crystal tests using gold chloride in H₂O ^{direct} I show diagram in separate page

a. dl meth - clathrospins

b. dl meth - X's on clathrospins stick head to head.

c. my diagram was what I observed - you may see it differently. By looking for these

subtle differences improve specificity of test - most clandestine meth is dl meth

this is much more important than differentiating dl, l, & dl

5 to 10%

IR sheet goes along with this

VI

Final analysis on separate sheet

1. IR: The CH_2Cl_2 + HCl in MeOH recrystallized with pet ether is the way we try to recrystallize any base that we don't have another specific I.R. method for. read so people in Boise etc can make comments and/or criticisms.

2. MeOH + hexane used to separate methamphetamine from epinephrine or caffeine

3. simple, quick procedure, works quite often

4. read self explanatory

6. read then ask for discussion and/or comments: sort of a read critique methodology let Ann discuss

7. Variations in IR of methamphetamine

1. properly prepared all meth gives different scan than d meth may have to warm or set in beam for awhile to get into most stable form

- I d meth HCl
- II d meth HCl I made
- III d meth HCl I made and left in IR beam for 11 hrs.

slight changes at 1465, 1365, 1000 to 1080 peaks at 700 & 750 change relative intensity slightly - this change very minor

d meth is quite different from 1000 to 1110 & also from 1400 to 1500

VII

IV & V can measure all with HCl

IV can scan immediately after recrystallization

V a scan obtained by heating overnight at 70°C

differences from 1400 to 1300 m. but more liked

also at 1570 to 1630

& 1000 to 1120

warmup scan of all with HCl from both sources virtually identical except for concentration differences.

Importance of all of all with HCl differences doesn't legally matter but you should be aware of these differences

1. optical isomers can give different scans
2. heating to convert to most stable crystalline form.

F:

2. more GLC - already discussed

3. Quant if you do quant to many possibilities to trust GC except in commercial preparations.

then cover all alone

then Percy

Order of material on analysis
of ~~methamphetamine~~ amphetamine, methamphetamine,
phenetamine, etc.

1. Copy of quick report used by Boise lab.
2. Procedure for amphetamine in central lab, a brief outline
3. outline of ~~ARDA~~ San Fran. DEA procedures for amphetamine 1972
4. article from Jour. of Chem. ~~104~~ 104 (1975) 201-204 mentions reactions of various drugs to Fluorescamine.
5. article from Jour. of Forensic Science, Jan 1976, vol 21, No. 1, p154; more complete description of reaction of various drugs to Fluorescamine.
6. ~~an~~ article from ^{Microgram VII}, No. 2, Feb 1974 describes use of acetone + 2% KOH as liquid phase for GLC of amphetamine, methamphetamine, etc.
7. ~~an~~ article ~~concerning~~ concerning the forming of mandelic salts of amphetamine for I.R.
8. ~~microgram~~ ~~article~~ ^{publication of} ~~function of~~ ~~IR~~ ~~of~~ amphetamine or methamphetamine for I.R. by Hexane-H₂O wash microgram, April, 1974
9. More advice on using Hexane-H₂O to purify amphetamine etc. for I.R. Microgram, April, 1975
10. ~~Purification~~ ~~with~~ Minibunnies purification for I.R., DEA
11. X-ray crystal test for amph. & methamphetamine
12. Procedure often followed in the central Lab for analysis of methamphetamine, phenetamine, etc. an outline

13. outline of DEA methodology for ~~street~~ microscopy
c 1977.

14. microcrystal test for methamphetamine using gold chloride
direct

15. IR methodology for metham, phentermine, used in the
central lab

16. 5 IR scans of methamphetamine HCl

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

San Fran ~~BMD 11~~ P DEA

ASPIRIN/ASPIRE Approx 1972

Identification (including optical activity)

- A. Infrared spectrum of the demandelate salt(s)
- B. Color and filter reaction examination - Barqais, gold chloride (volatility) platinum chloride (volatility) are used. Optical crystallography (also used for diluents), platinum chloride (direct), Berlinia chloride test for sulfate, gold chloride (direct) are used by a minority of chemists.

Assay

- A. Ultraviolet assay by direct dilution in acid and filtration. Solution sometimes washed with ether to remove interfering peak.

TLC

Most of their chemists now quantify

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT



STATE OF IDAHO

DEPARTMENT OF HEALTH AND WELFARE

BUREAU OF LABORATORIES, 2220 Old Penitentiary Road, Boise, ID 83702

Procedure often followed in the Central Lab for the analysis of amphetamine.

A. For preliminary analysis:

1. Spot tests: Marquis
Feigl's secondary amine test
Mecke
Froede
2. TLC: Either basic extract or a methanolic extraction is spotted on a plate and ran in Clark's T_1 (methanol 100/ NH_4OH , 1.5). The plate is visualized by a series of sprays: Fluorescamine, Iodoplatinate, Dragendorff's. (Dragendorff's I use mainly)
3. GLC: 10% Apiezon 1 + 2.5% potassium hydroxide useful for both quantitation and qualitative identification.

B. Additional analysis for a final report.

1. If there is enough material, we always try for an I.R. First we try to recrystallize as a salt of mandelic acid, perhaps with some preliminary cleanup to remove caffeine and/or ephedrine. Can catch crystals right on KBr in pasteur pipette. Might try the hexane - H_2O cleanup and formation of the hydrochloride salt.

If that fails:

2. Volatility microcrystalline tests (not very successful on the common junky crosstops)
 - a. Gold chloride in phosphoric acid.
 - b. Platonic chloride in phosphoric acid (more sensitive - may work more often).
3. More GLC: 3% OV17 3 feet column with and without acylation to observe peak shift. (Take up sample plus acetic anhydride and inject)
:20% OV17 6 feet column without acylation

Bradley (4)

Journal of Chromatography, 104 (1975) 201-204
© Elsevier Scientific Publishing Company, Amsterdam — Printed in The Netherlands

CHROM. 7864

Note

Thin-layer chromatography of 3,4-methylenedioxyamphetamine, 3,4-methylenedioxymethamphetamine and other phenethylamine derivatives

M. A. SHAW

Crime Detection Laboratory, Royal Canadian Mounted Police, Regina, Saskatchewan (Canada) and

H. W. PEEL*

Scientific Services, "L" Directorate, Royal Canadian Mounted Police, Ottawa, Ontario (Canada)

(Received June 13th, 1974)

Many phenethylamine derivatives (amphetamine drugs) are used medicinally for therapeutic treatment. Illicit manufacture and use of many of these derivatives have increased greatly in recent years. The relative ease with which chemical groups can be manipulated on the basic phenethylamine nucleus (moiety) has resulted in the synthesis and use of many new hallucinogenic substances, including such drugs as 3,4-methylenedioxyamphetamine (MDA) and 3,4-methylenedioxymethamphetamine (MMDA). Substances such as these are a problem to the forensic analyst who must have the capability to screen for them. There are often not sufficient differences of physico-chemical properties to readily distinguish popular phenethylamine drugs by chromatographic means.

The following study was carried out to develop a method that could be used to screen for MDA and MMDA in thin-layer systems already in use in our laboratory. As with other current methods using thin-layer chromatography (TLC) as a screening technique, a sequential spraying pattern is carried out^{1,2}. Indeed, the appearance of a particular color with a specific spray reagent, in conjunction with an accurate R_f value, can often serve to distinguish a particular substance from other chromatographically similar compounds. The use of a gallic acid spray, adapted from a procedure reported by DeMayo *et al.*³, served to distinguish MDA from MMDA. The application of the fluorogenic reagent fluorescamine is described for the general detection of phenethylamine substances.

EXPERIMENTAL

Apparatus

The chromatography tanks used were glass, generally of the size 23 × 12 × 23 cm, with glass tops sealed with starch glycerin paste. The tanks were lined with filter paper. Analtech silica gel G plates (250 μm) were activated for 20 min at 115° before use (available from Mandel Scientific Co., Montreal, Canada).

* To whom reprint requests should be addressed.

Property of the State Police Forensic Services
Unauthorized Photocopying Prohibited
PROPERTY DOCUMENT

Reagents

Developing solvents. (I) Methanol-ammonia (100:1.5); (II) Benzene (distilled)-dioxane (distilled)-ethanol-ammonia (150:120:15:15); (III) Chloroform (distilled)-cyclohexane (distilled)-diethylamine (50:40:10).

(V) *Fluorescamine spray.* Fluorescamine (Fluram[®], Roche Diagnostics, Vaudreuil, Canada), 5 mg/100 ml in acetone.

Ninhydrin spray. Ninhydrin 0.1% (w/v) in acetone is prepared fresh.

Gallic acid spray. Gallic acid 1.0% in sulphuric acid-ethanol (1:1). Store in the dark.

Potassium iodoplatinate (acidified) spray. Platinic chloride (0.25 g) and potassium iodide (5 g) are mixed with water to produce 100 ml. Hydrochloric acid (2 ml) is added to this solution.

Drug standards. Solutions of either the free base or salt were prepared in a known concentration of approximately 5 mg/ml in ethanol or methanol.

Procedure

Approximately 50 µg of each drug were spotted and run. The running distance of each plate was 15 cm. Codeine was run on each plate as an internal standard. A reference mixture of amphetamine, benzphetamine and mescaline was included on most of the plates to monitor the tank condition. The plates obtained from using solvent III were sprayed with a 5% ethanolic solution of hydrochloric acid prior to commencing the sequential spraying pattern to remove the effects of the residual diethylamine. After each run, the plate was air dried, sprayed with fluorescamine spray and viewed under ultraviolet light. All developed spots were noted. The plate was then oversprayed with ninhydrin spray, warmed at 70° for approximately 5 min and observed for any spot development. The plate was again oversprayed with gallic acid spray, warmed again for 10 min at 70° and observed. The plate was finally oversprayed with acidified potassium iodoplatinate spray for general development. After each spray sequence, the running distances of the spots were measured and the R_f and $R_{codeine}$ (R_f relative to codeine) values calculated. The values obtained for the phenethylamine substances and some other hallucinogens are given in Table I.

The sensitivity of detection of some of the drugs with ninhydrin spray, fluorescamine spray and gallic acid spray was determined by diluting the appropriate stock solutions and running the drug on the thin-layer plate in developing solvent I or II. These results are given in Table II.

The data in Table I show that it is often difficult to resolve MDA and MMDA using only R_f determinations. The inclusion of gallic acid spray in a sequential spraying pattern readily distinguishes these two compounds from each other with very good sensitivity. MDA appears as a green spot and MMDA appears as a blue spot. The other phenethylamine derivatives did not develop any interferences with this reagent.

A number of other TLC solvent systems had previously been investigated, but were discarded for reasons of poor resolution, poor sensitivity and difficult handling procedures. The three solvent systems used in this study are useful for many other drugs/poisons and generally give good reproducibility. A second spray reagent for MDA and MMDA was also considered—chromotropic acid, 1.0% in sulphuric acid-ethanol (1:1). A similar chromotropic reagent has been previously investigated by

NOTES

TABLE I
TLC DATA

Substance	R_f
MDA	0.1
MMDA	0.1
Mescaline	0.1
Amphetamine	0.1
Methamphetamine	0.1
Phenethylamine	0.1
Ephedrine	0.1
Phenylephrine	0.1
Benzphetamine	0.1
Phenmetrazine	0.1
Phendimetrazine	0.1
Chlorphentermine	0.1
Dimethyltryptamine	0.1
Diethyltryptamine	0.1
PCP	0.1
LBJ	0.1
Codeine	0.1

* S.D._{codeine} isDeMayo *et al.*,
respectively, and
the gallic acid sp

RESULTS AND DISCUSSION

The R_f values
running distance
usually from theTABLE II
DETECTION LIMITS

Substance	Detection Limit
MDA	
MMDA	
Mescaline	
Amphetamine	
Methamphetamine	
Phenethylamine	
Ephedrine	
Phenylephrine	

TABLE I
TLC DATA

Substance	Developing solvent								
	I			II			III		
	R_F	$R_{codeine}$	$S.D._{codeine}$ *	R_F	$R_{codeine}$	$S.D._{codeine}$ *	R_F	$R_{codeine}$	$S.D._{codeine}$ *
MDA	0.50	1.01	0.04	0.43	1.20	0.05	0.42	1.8	0.1
MMDA	0.55	0.98	0.05	0.41	1.10	0.05	0.40	1.7	0.1
Mescaline	0.46	0.81	0.02	0.19	0.57	0.05	0.26	1.10	0.09
Amphetamine	0.61	1.08	0.01	0.46	1.27	0.04	0.41	1.8	0.2
Methamphetamine	0.56	0.96	0.01	0.47	1.22	0.08	0.46	2.0	0.1
Phenethylamine	0.53	0.96	0.07	0.37	0.96	0.01	0.04	0.18	0.05
Ephedrine	0.48	0.85	0.04	0.26	0.71	0.05	0.13	0.60	0.04
Phenylephrine	0.46	1.81	0.03	0.03	0.09	0.02	0.05	0.20	0.03
Benzphetamine	0.78	1.37	0.03	0.76	2.3	0.3	0.71	3.0	0.2
Phenmetrazine	0.60	1.06	0.07	0.50	1.40	0.07	0.38	1.7	0.2
Phendimetrazine	0.02	0.02	0.01	0.01	0.04	0.01	0.18	0.9	0.2
Chlorphentermine	0.62	1.10	0.05	0.51	1.42	0.09	0.41	1.8	0.2
Dimethyltryptamine	0.59	1.04	0.01	0.48	1.34	0.04	0.21	0.89	0.04
Diethyltryptamine	0.65	1.15	0.03	0.61	1.70	0.08	0.31	1.33	0.03
PCP	0.71	1.27	0.05	0.82	2.30	0.20	0.72	3.1	0.3
LBJ	0.76	1.33	0.07	0.70	2.0	0.20	0.52	2.2	0.2
Codeine	0.57	1.00	—	0.36	1.00	—	0.24	1.00	—

* $S.D._{codeine}$ is the standard deviation of $R_{codeine}$.

DeMayo *et al.*³. The colors obtained with MDA and MMDA were pink and purple, respectively, and somewhat easier to differentiate than the gallic acid spray. However, the gallic acid spray was found to be more sensitive.

RESULTS AND DISCUSSION

The R_F values (and subsequent $R_{codeine}$ values) were obtained by measuring the running distances obtained from three separate tanks of the same solvent system, usually from three separate plates within that tank. Thus, most of the values recorded

TABLE II
DETECTION LIMITS (μ g) OF STRAYS

Substance	Fluorescamine	Ninhydrin	Gallic acid
MDA	0.05	20	1
MMDA	0.05	10	1
Mescaline	0.05	5	
Amphetamine	0.05	20	
Methamphetamine	0.5	20	
Phenethylamine	0.05	5	
Ephedrine	3.0	10	
Phenylephrine	0.1	5	

in Table I are the result of nine values. In addition, similar data have been obtained separately in another laboratory.

The application of fluorescamine as a spray was investigated and found to be very useful. After spraying, the plate is viewed under ultraviolet light. Generally the phenethylamines are visible as green fluorescing spots or as dark-blue absorbing spots. The intensity of the latter may be increased by exposing the plate to a vapor of ammonia prior to spraying. Table II shows the detection limits obtained by the fluorescamine spray on some of the drugs.

ACKNOWLEDGEMENT

The authors wish to express their appreciation to the Health Protection Branch, Health and Welfare, Canada for supplying some of the standard substances.

REFERENCES

- 1 G. F. Phillips and J. Gardiner, *J. Pharm. Pharmacol.*, 21 (1969) 793
- 2 M. L. Bastos, D. Jukofsky and S. J. Mule, *J. Chromatogr.*, 81 (1973) 93.
- 3 M. M. DeMayo, E. J. Briglia, Jr. and L. A. Dal Cortivo, *J. Forensic Sci.*, 17 (1972) 444.

Property of Idaho State Police Forensic Services
Uncontrolled Interim Copy
OBSOLETE DOCUMENT

Note

Analysis of topography

PETER JAMES
Home Office
(Received September 1973)

Metals
variety of
morphine,
ephedrine
parison of
identificati
at least th
(HPLC) of
of thin-lay
trapping
of the pres
field and t
samples.

Seve
alkaloids?
diamorphin
extensively
regarded as
analysis.

In t
ysis of
chromato
graphic se
twelve nu

EXPERIMENTAL

Equipment
Chromatography
Zipax SC
eluent was
a Waters
Instrument
by means

ing the b-

(5)

Studies on Fluorescamine: Part I—Applications of Fluorescamine in Forensic Toxicological Analysis

We use at concentration of 5mg/100ml and its sensitivity is just fine

Weigle and co-workers [1,2] reported the structure and synthesis of the reagent 4-phenylspiro [furan-2(3H), 1'-phthalan]-3,3'-dione (fluorescamine), which reacts with substances containing primary amino groups to yield highly fluorescent products. Recently, Undenfriend and co-workers [3-5] published a paper describing the use of fluorescamine in biochemical analysis. This reagent was of interest because of its possible use to improve the detection and differentiation of amphetamine and its relatives in biological samples and solid dosage forms. The main subject of this paper is to present data which clearly demonstrate the use of fluorescamine in forensic toxicological analysis of amphetamine and differentiation from methamphetamine. This reagent also has potential application for analysis of other drugs containing a primary amine group.

Two compounds of primary interest to the law enforcement effort are amphetamine and methamphetamine, although a few other amphetamine-like drugs are also federally controlled. The physiological effects of these two drugs are similar [6-8], methamphetamine having a slightly greater stimulant effect on the central nervous system. Amphetamine and methamphetamine are both optically active compounds. The dextro isomers are about twice as active as the racemic mixtures, however, all isomers are under Federal control and are listed in Schedule II of Public Law 91-513.

An excellent review of existing methods for analysis of amphetamine analogs has been published [9]. Miles and Schenk [10] made use of the natural fluorescence of phenylethylamines to assay these compounds in pharmaceutical preparations. Few attempts have been made to form fluorescent derivatives of amphetamine and related substances.

In this work, fluorescamine is used to form fluorescent derivatives. It has also been found very simple to further analyze the fluorescent derivative formed in the spot test by thin-layer chromatographic analysis.

Methods and Materials

Fluram™ (fluorescamine) was purchased from Roche Diagnostics, Division of Hoffman-La Roche, Inc., Nutley, N.J. The reagent is supplied in vials containing 100-mg fluorescamine crystals. Fluram is stable at room temperature in both solution and powder form; thus, refrigeration is not recommended by the manufacturer.

The fluorescamine is prepared into a working solution by dissolving 50 mg of fluorescamine in 100 ml of acetone. Fluorescamine has minimal solubility in water

Received for publication 30 Aug. 1974; revised manuscript received 31 March 1975; accepted for publication 3 April 1975.

*Forensic toxicological mass spectroscopist, U.S. Army Criminal Investigation Laboratory.

154

and it will decompose in the presence of water. In the assay, excess reagent is hydrolyzed to water-soluble nonfluorescent products [3].

In Fig. 1, fluorescamine (I) reacts with primary amines (II) to form intensely fluorescent substances (III), providing the basis for a rapid and highly sensitive assay for compounds containing a primary amine group, such as amino acids, primary amines, peptides, and proteins [3]. The reagent (I) does not react with secondary or tertiary amines. Thus, it provides a quick method to distinguish different types of amines; this has been found to be very helpful to distinguish between amphetamine and methamphetamine. Amphetamine is a primary amine and yields an intensely blue-green fluorescent product in the fluorescamine test and methamphetamine does not yield a fluorescent product.




FIG. 1—The reaction between Fluram[®] (I) and a primary amine (II).

The spot test procedure used was to add to a tile-welled spot plate two drops of borate buffer, pH 9.0, and check the fluorescence. All fluorescence examinations were performed in a Chromato-Vue[®] box (manufactured by Ultra-Violet Products, Inc., San Gabriel, Calif.,) using the long wave (366 nanometer) for excitation. In the spot test the fluorescence is so strong that it can be detected in room light with a small hand ultraviolet (UV) source for excitation. Then a small amount of the exhibit sample was added to the borate buffer and the fluorescence was determined. It is important to check the cleanliness of the plate and the natural native fluorescence of the substance for proper interpretation of the results. The test was completed by adding one drop of the fluorescamine acetone solution and checking the long-wave fluorescence. A positive test was the formation of a very intense blue-green fluorescent product. A negative test was one that did not form a fluorescent product or had fluorescence of another color due to the native fluorescence of the substance.

Quite often it was found to be helpful to analyze the intensely blue-green product further; this was done by thin-layer chromatography (TLC). Two to three microlitres of the spot test were applied to a silica gel thin-layer chromatogram which did not contain any fluorescent material. Standard fluorescent products were prepared from known substances and also applied to the chromatogram. The chromatograph was then developed in a chloroform:methanol system (90:10). After development the fluorescent products were located with the aid of long-wave UV light and the R_f values compared to that of standard compounds. Since very small amounts of the aqueous reaction mixture were required for thin-layer chromatography analysis, there were no problems with drying of the chromatogram before chromatography.

Amphetamine present in urine specimens (10 ml) was isolated by extraction at pH 11 with two volumes of chloroform:isopropanol (3:1) [11]. The solvent extracts were treated with 0.10 ml of acidified methanol (0.1 mol sulfuric acid per litre of methanol) and evaporated. The residue was redissolved in 50 µl of methanol for application

to a thin layer plate

Simultaneously, in duplicate, drug free urine and amphetamine-containing urine (drug-free urine that was "spiked" with amphetamine ranging from 0.10 to 30.0 μg per 10 ml urine) specimens were prepared. The duplicates were then spotted on separate TLC plates, coated with a 250- μm layer of silica gel without fluorescent indicator. The TLC plates were developed in ethyl acetate:methanol:concentrated ammonium hydroxide (170:20:10 by volume) [12]. The plates were removed when the solvent migrated between 18 and 19 cm and dried with warm air from a hair drier type of air blower. Then the plates were dried 10 min in an oven at 75°C to remove residual ammonia. Residual ammonia interferes with the ninhydrin spray; it does not interfere with the fluorescamine reagent.

One TLC plate was developed with ninhydrin and the duplicate was developed with fluorescamine. When the fluorescamine test was applied to detection of amphetamine on thin-layer chromatograms, the fluorescamine-acetone solution (50 mg%) was first sprayed onto the dried chromatogram and then oversprayed immediately with the pH 9.0 borate buffer. The fluorescent produce was visualized and located with long-wave UV light (366 nm).

Results and Discussion

Theoretically, the fluorescamine test should yield positive results only with compounds containing a primary amine group. In order to verify this and evaluate possible interfering compounds, many common drugs which are seen as exhibits in the forensic chemistry laboratory were tested. The standard compounds contained in this laboratory consist of a combination of pure drugs, tablets, and capsules. The chemical constituents of the trade name compounds can be found in standard references [13,14]. About 10% of the standard compounds tested yielded a positive fluorescamine test; these results are contained in Table 1. From the molecular structures of the compounds examined, it was concluded that a primary amine was a requirement for obtaining a positive fluorescamine test.

Many compounds have native fluorescence; it is important to note this by determining the presence of fluorescence at the stage of addition of borate buffer, before the addition of fluorescamine. It is important to run blanks with this test because it is extremely sensitive and very small contamination will yield false positive results.

Amphetamine and methamphetamine can be distinguished easily by the fluorescamine test. Both of these compounds give the same colors in the Marquis, Meckes, and Froehdes spot tests, and the UV spectra of these compounds are indistinguishable. However, the fluorescamine test is positive with the primary amine, amphetamine, and negative with the secondary amine, methamphetamine.

Sometimes amphetamine preparations have caffeine or other components which interfere with the UV spectrum. Under this circumstance it is convenient to spot 2 to 3 μl of the positive fluorescamine spot test on a thin-layer chromatogram and compare with standard fluorescent derivatives. The chromatographic mobility of the fluorescent derivative of amphetamine has been found to be unique. No other standard compound (see Table 1) has been found to form a fluorescent derivative with fluorescamine and then migrate the same to that of amphetamine fluorescent derivative. The amino acids tyrosine, phenylalanine, and histidine have been studied, as well as Aldomet® (α -methyldihydroxyphenylalanine). All four of the fluorescent derivatives of these compounds remained at the origin in the chloroform:methanol (90:10) mobile phase, while the amphetamine fluorophore moved at an R_f of 0.30.

Often it is necessary to perform urinalysis to determine amphetamine abuse. The classic method of detection of amphetamine in urine is extraction and TLC analysis with ninhydrin spray. The generally accepted minimal detection limit of ninhydrin reagent spray for this

TABLE I—Examination of standard drugs with the fluoroxamine test.

Drug	Test	Drug	Test
1. Acetylcodeine	-	49. Caffeine (PCP)	-
2. Amesec [®] Caps (aminophylline compound)	-	50. Cafilon [®] (calcium lactate)	-
3. Aminophylline + Phenobarbital Tablets	-	51. Calcium cyclobarbital	-
4. Aminopyrine	-	52. Camoquin [®] HCl (amodiaquin)	-
5. Aminosalicyclic acid	-	53. Camphocodine [®]	-
6. Amobarbital	-	54. Carbarsone [®] (N-carbamoylsalicylic acid)	+
7. dl-Amphetamine	+	55. Carbinoxamine maleate	-
8. d-Amphetamine sulfate	+	56. Carisoprodol	-
9. Sodium Amytal [®]	-	57. Cascara sagrada	-
10. Amphojel [®] (aluminum hydroxide)	-	58. Chloracetophenone	-
11. Ampicillin trihydrate	+	59. Chloramphenicol	-
12. Ampicillin	+	60. Chloromycetin [®] powder (chloramphenicol)	-
13. Anacin [®]	-	61. Chloral hydrate	-
14. Anahist [®]	-	62. Chloropylene	-
15. Analexin-AF [®]	-	63. Chloroquimediophosphate	-
16. Ananase [®] (plant protease concentrate)	-	64. Chloroquine phosphate	-
17. Ansional [®]	-	65. Chlorpheniramine	-
18. Antipyrine	-	66. Chlortetracycline HCl	-
19. Anturane [®]	-	67. Chlor-Trimeton [®] (chlorpheniramine maleate)	-
20. Apresoline [®] HCl (hydralazine HCl)	-	68. Cholan [®] (dehydrocholic acid)	-
21. Ascorbic acid	-	69. Chrysarobin	-
22. Atarax [®] (hydroxyzine HCl)	-	70. Cibalgin [®]	-
23. Atraxin [®] (meprobamate)	-	71. CKD Tablet [®] (ephedrine HCl)	-
24. Atropine sulfate	-	72. Clarnil	-
25. Auramin [®]	-	73. Cocaine	-
26. Aureomycin [®]	-	74. Codeine	-
27. Aventyl [®] HCl (nortriptyline HCl)	-	75. Colehicine	-
28. Barbital	-	76. Combid [®]	-
29. Belladonna	-	77. Compazine [®]	-
30. Bellerгал [®] Tablets	-	78. Co-Pyronil [®]	-
31. Bellerгал [®] Spacetabs [®]	-	79. Cortancyl [®] (prednisone)	-
32. Benadryl [®]	-	80. Coricidin D [®]	-
33. Bendectin [®]	-	81. Cortisone acetate	-
34. Benmid [®] (probenecid)	-	82. Creatine	-
35. Benzadrine [®] (amphetamine sulfate)	+	83. Crystoids [®] (hexylresorcinol)	-
36. Benzocaine (ethyl aminobenzoate)	-	84. Crystoserpine [®] (reserpine)	-
37. Beta-Chlor [®] (chloral betaine)	-	85. Cyclandel	-
38. Bicillin [®] (benzathine penicillin G)	-	86. Cyclobarbitol	-
39. Binoctal [®] (amobarbital, secobarbital)	-	87. Dalmane [®] (flurazepam HCl)	-
40. Bismuth subnitrate	-	88. Daprisal [®]	+
41. Bromanautine [®]	-	89. Dapsone [®] (4', 4'-sulfonyldianiline)	+
42. Brovarin [®]	-	90. Darvon [®] (propoxyphene HCl)	-
43. Brucine	-	91. Darvon [®] Compound-65	-
44. Butabarbital	-	92. Darvo-Tran [®]	-
45. Butalbital	-	93. DBI [®] (phenformin HCl 5)	-
46. Butazolidin [®] (phenylbutazone)	-	94. Declomycin [®] (demethylchlorotetracycline HCl)	-
47. Cafergot [®] (ergotamine tartrate)	-	95. Delta-Cortef [®] (prednisolone)	-
		96. Dexamyil [®]	+
		97. Dexedrine [®] (d-amphetamine sulfate)	+
		98. Diabinese [®] (chlorpropamide)	-
		99. Dianabol [®] (methandrostenolone)	-

Property of Idaho State Police Forensic Services
 UNCONTROLLED INTERFERENCE DOCUMENT
 OBSOLETE

Drug	Test	Drug	Test
101. Digitalis	-	150. High Crelan* (secopyrabital)	-
102. Digitoxin	-	151. Homatropine hydrobromide	-
103. Digoxin	-	152. Hygroton* (chlorthalidone)	-
104. Dihydrostreptomycin polymix	-	153. Hyminal* (methaqualone)	-
105. Diiodohydroxyquin	-	154. Hyoseyamine or stramonium	-
106. Dilantin* (diphenylhydantoin)	-	155. Iloxone* (erythromycin ester late)	+
107. Dimetane* (brompheniramine maleate)	-	156. Immenoctal* (amobarbital)	-
108. Dimetapp*	+	157. Indocin* (indomethacin)	-
109. Dioctyl sodium	-	158. Ismelin* sulfate (guanethidine sulfate)	-
110. Diphenhydramine	-	159. Isopropyl meprobamate	-
111. Diphenylhydantoin with phenobarbital	-	160. Isordil* (isosorbide dinitrate)	-
112. Disophrrol Chronotab*	-	161. Kafe-Sol* (caffeine)	-
113. Diuril*	-	162. Kemicefine* (chloramphenicol)	+
114. Doloran* (allobarbitol)	-	163. Keramine	-
115. Donnatal*	-	164. Leukeran* (chlorambucil)	-
116. Doriden* (glutethimide NF)	-	165. Librium* (chloridiazepoxide HCl)	-
117. Dormopan* (hexobarbital)	-	166. Lidocaine HCl (injection, USP 2%)	-
118. Doxidan*	-	167. Lincocin* (lincomycin HCl monohydrate)	-
119. Dramamine* (dimenhydrinate)	-	168. Domotil* (diphenoxylate HCl with atropine sulfate)	-
120. Drixoral*	-	169. Lysergic acid diethylamide	-
121. Ducolax* (bisacodyl)	-	170. Mandelamine* (methenamine mandelate)	-
122. Edrisal*	+	171. Mannitol (Korean)	-
123. Elavil* HCl (amitriptyline HCl)	-	172. Maple Tablet (maple preparation)	-
124. Ephedrine HCl	-	173. MDA (methylenedioxyamphetamine)	+
125. Ephedrine sulfate	-	174. Marezine* (cyclizine)	-
126. Ergotrate* maleate (ergonovine maleate)	-	175. Meclizine	-
127. Erythrocin* (erythromycin ethylsuccinate)	-	176. Meconic acid	-
128. Erythromycin-base Film Tab*	-	177. Medroxyprogesterone acetate	-
129. Erythromycin stearate (Iloxye*)	-	178. Mellaril* (thioridazine)	-
130. Erythromycin-stan*	-	179. Mephesisin	-
131. Eskatrol* Spansule*	+	180. Mephobarbital	-
132. Estinyl* (ethinyl estradiol)	-	181. Mepivacaine HCl injection 1% (carbocaine)	-
133. Ethylmorphine	-	182. Meprobamate	-
134. Ferrous sulfate	-	183. Meprophen* (meprobamate)	-
135. Fiorinal*	-	184. Meridon	-
136. Flagyl*	-	185. Merthiolate* (thimerosal)	-
137. Floraquin* vaginal tablets (diiodohydroxyquin)	+	186. Methadone	-
138. Flurazepam HCl	-	187. Methamphetamine	-
139. Gantrisin* (sulfisoxazole)	+	188. Methapyrilene	-
140. Glutethimide	-	189. Methaqualone	-
141. Grifulvin V* (griseofulvin microsize)	-	190. dl-Methionine and vitamin B ₁₂	+
142. Griseofulvin	-	191. Methylcellulose	-
143. Gynergen* (ergotamine tartrate)	+	192. Methylhexobarbital	-
144. Halotestin* (fluoxymesterone)	-	193. 3-Monoacetylmorphine	-
145. Hashish	-	194. Mycostatin* (nystatin USP)	+
146. Helozid	-	195. Myleran* (busulfan)	-
147. Heptabarbital	-	196. Mylicon*	-
148. Heroin (diacetylmorphine)	-	197. Mysoline* (primidone)	-
149. Hetrazan* (diethylcarbamazine citrate)	-	198. Naron* (cyclopyrabital)	-

Drug	Test	Drug	Test
295. Tetracycline	-	311. Tylenol ¹ (acetaminophen)	-
296. Tetrex ² (tetracycline phosphate)	-	312. Urecholine ¹ (bethanechol chloride)	-
297. Thebaine	-	313. Valium ¹ (diazepam)	-
298. Theophylline	-	314. Vallestril ¹ (methallenestril)	-
299. Thianphenicol	-	315. Vasodilan ¹ (isoxsuprine HCl)	-
300. Thorazine ² (chlorpromazine)	-	316. Vigosan ¹	-
301. Tigan ² (trimethobenzamide HCl)	+	317. Vistaril ¹ (hydroxyzine pamoate)	-
302. Titalac ² (calcium carbonate)	-	318. Vitamin B ₁₂	+
303. Tofranil ² (imipramine HCl)	-	319. Vitamin K (synthetic)	-
304. Tranquinal ² (meprobamate)	-	320. Warfarin	-
305. Triacetyloleandomycin	+	321. Wyamine sulfate ² (mephentermine)	-
306. Trilafon ² (perphenazine)	-	322. Wyanooids ²	-
307. Tri-Span ²	+	323. Zactrin ²	-
308. Tropacocaine HCl	-	324. Zaronin ² (ethosuximide)	-
309. Tuinal ² (sodium amobarbital and sodium secobarbital)	-	325. Zylorin (allopurinol)	-
310. Tuss-Ornade ²	-		

limit of ninhydrin and fluorescamine. It has been found that fluorescamine is a more sensitive method for the detection of amphetamine on thin-layer chromatograms. These results are tabulated in Table 2. The detection limit of amphetamine in urine samples is increased 100 times when the only parameter varied is the method of detection, ninhydrin or fluorescamine.

Other physiological chemicals containing primary amino groups normally present in drug-free urine extracts will also react with fluorescamine. This was evidenced in the non-drug-containing urine (Table 2). As many as eight well-separated fluorescent areas were routinely observed. These naturally occurring substances containing primary amino groups all have mobilities considerably less than that of amphetamine in the developing solvent system used; thus, the naturally occurring compounds do not interfere with the analysis.

While this paper was in the process of publication another has appeared on the use of fluorescamine in amphetamine detection in urine [16]. These authors have also stressed the increased detection limit of fluorescamine compared to ninhydrin. Klein et al quantitated the amphetamine level by extracting the silica gel area containing the amphetamine fluorophore and analyzed the extract in a microfluorimeter. They also reported that after amphetamine was made visible by fluorescamine spray, the plate could be sprayed with other common identification reagents with no interference. Accordingly, it is possible to substitute fluorescamine for ninhydrin in a routine battery of sprays for drug abuse screening in urine samples.

Summary

This paper describes some applications of the fluorescamine spot test to forensic toxicological analysis. The fluorescamine test only reacts with primary amines; thus, this test makes a clear-cut distinction between amphetamine and methamphetamine. Previous common spot tests used reacted the same with these two amines. Fluorescamine is 100 times more sensitive in detecting amphetamine extracted from urine on thin-layer chromatograms than ninhydrin. Thus, it is a more sensitive method of detecting amphetamine abuse in urinalysis screening programs.

Drug	Test	Drug	Test
199. NegGram [®] (nalidixic acid)	-	247. Promacetin [†] (acetosalzone sodium)	-
200. Nembutal [®] (sodium pentobarbital)	-	248. Pronestyl [†] (procainamide HCl)	-
201. Neomycin	+	249. Propadrine [®] (phenylpropanolamine)	+
202. Nicotinic acid	-	250. Propadrine HCl [®]	+
203. Nitrofurantoin	-	251. Propoxyphene HCl	-
204. Nitroglycerin	-	252. Prostaphilin [®] (sodium oxacillin)	-
205. Norflex [®] (orphenadrine citrate)	-	253. Provest [®] Daypak (birth control)	-
206. Norgescic [®]	-	254. Pyridium [®] (phenazopyridine HCl)	-
207. Norinyl [®] (norethindrone with mestranol)	-	255. Pyridoxine HCl (vitamin B ₆)	-
208. Noludar [®] (methyprylon)	-	256. Quinacrine HCl	-
209. Novahistine [®] LP tablets	-	257. Quinidine sulfate	-
210. Novatophen [®] (neocinchophen)	-	258. Quinine sulfate	-
211. Optalidon [®] (allylisobutylbarbital)	-	259. Rarical	-
212. Opium (raw)	-	260. Rela [®] (carisoprodol)	-
213. Oretic [®] (hydrochlorothiazide)	-	261. Rhubarb	-
214. Ornade [®] Spansule [®]	+	262. Riboflavin (vitamin B ₂)	-
215. Oxsoralen [®] (methoxsalen)	-	263. Ritalin [®] HCl (methylphenidate HCl)	-
216. Paraflex [®] (chlorzoxazone)	-	264. Robaxin [®] (methocarbamol)	-
217. Pavatrine [®] HCl	-	265. Romilar [®] (racemethoxyphane HBr)	-
218. Pavron [®]	-	266. Rotoxamine	-
219. Penicillin G potassium	-	267. Salol [®] (phenyl salicylate)	-
220. Penicillin, phenoxymethyl potassium	-	268. Sand	-
221. Pentobarbital	-	269. Sansert [®] (methysergide maleate)	-
222. Periacin [®] HCl (cyproheptadine HCl)	+	270. Saridon [®] (isopropyl antipyrine)	-
223. Peritrate [®] (pentaerythritol tetranitrate)	+	271. Scopalamine hydrobromide	-
224. Persantine [®] (dipyridamole)	-	272. Secobarbital	-
225. Phenaphen [®]	-	273. Sedes [®] (hexobarbital)	-
226. Phencyclidine	-	274. Sedalin [®] (pyrabital)	-
227. Phenergan [®]	-	275. Sernylan [®] (phencyclidine HCl)	-
228. Phenformin HCl	-	276. Sinequan (10 mg)	-
229. Phenobarbital	-	277. Sintrom [®] (acencoumarol)	-
230. Phenobarbital, ephedrine, and theophylline	-	278. Sodium bicarbonate	-
231. Phenylpropanolamine HCl	+	279. Sodium chloride and sodium bicarbonate	-
232. Physostigmine salicylate	-	280. Sodium salicylate	-
233. Pilocarpine nitrate	-	281. Sparine [®] (pronazine HCl)	-
234. Pival [®] (pindone)	-	282. Stelazine [®] (trifluoperazine)	-
235. Placidyl [®] soft capsules (ethchlorvynol)	-	283. Streptomycin	-
236. Polycillin [®] (ampicillin trihydrate)	-	284. Strychnine	-
237. KMnO [®] tablet	-	285. Sudafed [®] (pseudoephedrine HCl)	-
238. Povan [®] (pryvinium pamoate)	-	286. Sulfadiazine	-
239. Prednisolone	-	287. Surfak [®] [calcium bis-(dioctyl sulfosuccinate)]	-
240. Prednisone	-	288. Syntho Tab	-
241. Preludin [®] (phenmetrazine HCl)	-	289. Surgex [®] (pipradral)	-
242. Premarin [®]	-	290. Talwin [®] (pentazocine)	-
243. Pre-Sate [®] (chlorphentermine HCl)	-	291. Tapazole [®] (methimazole)	-
244. Primaquine phosphate	-	292. Tenuate [®] , Dospan [®] (diethylpropion HCl)	-
245. Pro-Banthine [®] (propantheline bromide)	-	293. Terramycin [†] (oxytetracycline)	-
246. Procaine hydrochloride	+	294. Tessalon [®] (benzonatate)	-

Jell.
tes,
ical
nine
974.

Property of the State Police Forensic Services
PROPERTY CONTROLLED INTERNET DOCUMENTS

TABLE 2—Relative detection of amphetamine isolated from urine with ninhydrin and fluorescamine on thin-layer chromatograms.

Amphetamine*, µg	Ninhydrin	Fluorescamine
30.0	+	+
20.0	±	+
10.0	—	+
5.0	—	+
0.5	—	+
0.25	—	+
0.1	—	±
0	—	—

*Amount of amphetamine added to 10 ml of control, non-drug-containing urine.

Acknowledgments

I gratefully acknowledge the forensic chemists of this laboratory, D. O. Browning, R. E. Dolle, Jr., R. P. Erickson, R. Wojciechowski, E. R. Lieber, and J. H. Stopper, II; and L. L. Pyles of the U.S. Army Medical Laboratory, St. Louis, Mo., for making numerous comments and using the fluorescamine test in everyday forensic toxicological chemical analysis this past year. I also acknowledge the commander of this laboratory, Donald B. Jackson, for supporting this research.

References

- [1] Weigele, M., Blount, J. F., Teng, J. P., Czajkowski, R. C., and Leingruber, W., "The Fluorogenic Ninhydrin Reaction. Structure of the Fluorescent Principle," *Journal of the American Chemical Society*, Vol. 94, No. 8, June 1972, pp. 4052-4054.
- [2] Weigele, M., DeMernardo, S. L., Teng, J. P., and Leingruber, W., "A Novel Reagent for the Fluorometric Assay of Primary Amines," *Journal of the American Chemical Society*, Vol. 94, No. 9, July 1972, pp. 5927-5928.
- [3] Udenfriend, S., Stein, S., Bohlen, P., Dairman, W., Leingruber, W., and Wiegele, M., "Fluorescamine: A Reagent for Assay of Amino Acids, Peptides, Proteins, and Primary Amines in the Picomole Range," *Science*, Vol. 178, No. 24, 1972, pp. 871-872.
- [4] Stein, S., Bohlen, P., Stone, J., Dairman, W., and Udenfriend, S., "Amino Acid Analysis with Fluorescamine at the Picomole Level," *Archives of Biochemistry and Biophysics*, Vol. 155, No. 1, March 1973, pp. 203-212.
- [5] Bohlen, P., Stein, S., Dairman, W., and Udenfriend, S., "Fluorometric Assay of Proteins in the Nanogram Range," *Archives of Biochemistry and Biophysics*, Vol. 155, No. 1, March 1973, pp. 213-220.
- [6] Cutting, W. C., *Handbook of Pharmacology: The Actions and Uses of Drugs*, 3rd ed., Appleton Century Crofts, New York, 1967, pp. 406-420.
- [7] Burger, A., *Medicinal Chemistry*, 2nd ed., Interscience Publishers, Inc., New York, 1966, pp. 594-616.
- [8] Goodman, L. S. and Gilman, A., *The Pharmacological Basis of Therapeutics*, 3rd ed., Macmillan Co., New York, 1965, pp. 484-510.
- [9] Kaistha, K. K., "Methods of Analysis of Amphetamine Analogs," *Journal of Pharmaceutical Sciences*, Vol. 61, No. 3, March 1972, pp. 655-690.
- [10] Miles, C. I. and Schenk, G. H., "Use of Natural Fluorescence of Phenylethylamines to Assay These Compounds in Pharmaceutical Preparations," *Analytical Chemistry*, Vol. 45, No. 1, Jan. 1973, pp. 130-134.
- [11] Mule, S. J., "Routine Identification of Drugs of Abuse in Human Urine, Part I," *Journal of Chromatography*, Vol. 55, No. 2, Feb. 1971, p. 255.
- [12] Davidow, B., Li Petri, N., and Quame, B., "A Thin-Layer Chromatographic Screening Procedure for Detecting Drug Abuse," *American Journal of Clinical Pathology*, Vol. 50, No. 1

- [13] Stecher, P. G., Ed., *The Merck Index*, 8th ed., Merck & Co., Inc., Rahway, N.J., 1968.
- [14] Baker, E. C., Jr., Ed., *Physicians' Desk Reference*, 28th ed., Litton Publication, Oradell, N.J., 1974.
- [15] Kaistha, K. K. and Jaffe, J. J., "TLC Techniques for Identification of Narcotics, Barbiturates, and CNS Stimulants in a Drug Abuse Urine Screening Program," *Journal of Pharmaceutical Sciences*, Vol. 61, No. 3, March 1972, p. 679.
- [16] Klein, B., Sheehan, J. F., and Grunberg, E., "Use of Fluorescamine to Detect Amphetamine in Urine by Thin-Layer Chromatography," *Clinical Chemistry*, Vol. 20, No. 2, March 1974, p. 272.

USACIL-Pacific
APO San Francisco, Calif. 96343

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

LABORATORY NOTES

Microgram, VII, No. 2

Feb., 1974

(6)



DATE

-22-

NO.

DRUG TYPE

METHODOLOGY

The Determination of Amphetamine and Methamphetamine Preparations by Gas Chromatography by Stanley Blasof and Jack Fasanello Northeast Regional Laboratory Drug Enforcement Administration

A gas chromatographic procedure is presented for the quantitative determination of Amphetamine and Methamphetamine. This procedure, which primarily involves a direct methanolic dilution of the sample preparation and the use of naphthalene as an internal standard, can be utilized for the analysis of the above drugs when they appear alone or in combination with each other, or other interfering substances.

Introduction

The quantitative determination of amphetamine mixtures by ultraviolet analysis has presented some difficulty. Current separation techniques often prove difficult, time consuming and incomplete. Interferences from other components within the sample mixture becomes a problem due to the relatively low absorbtivity of amphetamine.

The successful separation of amine drugs on alkaline GC columns has been previously reported (1) (2) (3). To quantitate amphetamine and/or methamphetamine a GC procedure employing a 10% Apiezon L + 2% KOH column is presented.

Experimental

Apparatus

Gas Chromatograph:	Packard Model 804 equipped with a flame ionization detector.
Column:	Borosilicate glass. Length - 6' I.D. - 2mm
Packing:	10% Apiezon L plus 2% KOH (Analabs, North Haven, Connecticut) on Chromasorb WHP 80-100 Mesh.
Operating Conditions:	Column Temperature 145°C Inlet Temperature 300°C Detector Temperature 210°C Carrier Gas Nitrogen Flow Rate 40 ml./min. Range 1 x 10 ⁻⁹ amps

Preparation of Standard Solutions:

a. Internal Standard - prepare a stock solution by accurately weighing and transferring 400-600 mgs. of naphthalene into a 250.0 ml. volumetric flask and dilute to volume with methanol.

b. Reference Standard - accurately weigh portions of the reference standards equivalent to 3-5 mgs. of amphetamine and/or methamphetamine and transfer into a 10.0 ml. volumetric flask. Pipet 1.0 ml. of the internal standard solution into the flask and dilute to volume with methanol.

Procedure

Sample Preparation - Uniformly mix the entire sample powder. Time delay pellets must be thoroughly ground and passed through a 60 mesh sieve. Accurately weigh and transfer a portion of the sample mixture equivalent to 3.0 mgs. of amphetamine to a 10.0 ml. volumetric flask. Add 1.0 ml. of internal standard solution and dilute to volume with methanol. Preparations containing time delay pellets should be shaken on a mechanical shaker for 15 minutes to insure all amphetamine is in solution.

Calculations:

One microliter of the standard and sample solutions are injected into the gas chromatograph (See Fig. 1). The concentration of amphetamine and/or methamphetamine is calculated using the following formula:

$$\% \text{ amphetamine} = \frac{H_{SP}}{H_{ST}} \times \frac{H_{IST}}{H_{ISP}} \times \frac{C_{STD}}{WT.SP} \times 100$$

H_{SP} = height of sample peak

H_{ST} = height of standard peak

H_{IST} = height of internal standard peak in standard solution

H_{ISP} = height of internal standard peak in sample solution

C_{STD} = concentration of standard

WT.SP = weight of sample used

Results and Discussion

The ratio of peak height to concentration was found to be linear over the range of 0.1 - 1 mg. per ml. of methamphetamine and amphetamine. Other concentrations were not determined.

This procedure has been utilized routinely by this laboratory for the past year for both illicit as well as commercial preparations. Illicit

1. Phenolic amines are not eluted from this column.

preparations of amphetamine or methamphetamine containing interfering by products resulting from the synthesis were accurately quantitated by this procedure. An example of an unusual clandestine preparation which this procedure readily resolved consisted of a mixture of amphetamine, methamphetamine and ephedrine. The amphetamine and methamphetamine were quantitated at 145°C; the ephedrine content was estimated² using a column temperature of 165°C (See Fig. 2).

The amphetamine content in "Eskatrol" capsules was also quantitated by this procedure. (Eskatrol TDC's is a commercial preparation consisting of a mixture of amphetamine and prochlorperazine.)

The standard solutions employed are stable indefinitely when refrigerated. Therefore, most sample preparations can be determined rapidly and accurately by direct dilution, addition of internal standard and subsequent injection of the sample and standard solutions into the gas chromatograph. An additional advantage of this procedure is that sample handling is maintained at a minimum and derivatization is eliminated.


It should be noted that this laboratory has concurrently been using a Carbowax 20M + 2% KOH column for the analysis. Amphetamine and methamphetamine exhibit a reversal of elution order when chromatographed on this column, thus providing an excellent qualitative check. We have also used the Carbowax column for quantitative analysis, but prefer the Apiezon column because the separation is better and analysis time is shorter. Further qualitative proof may be obtained using the Schiff base-derivative technique.⁴

Reference

1. Parker, K.D., Fontan, C.R., Kirk, P.L., Analytical Chemistry, 34, 1345 (1962).
2. Beckett, A.H., Rowland, M., Journal of Pharmacy and Pharmacology, 1964, 16.
3. Chromatography Lipids Newsletter, Vol. VI, No. 3.
4. Brochmann-Hanssen, E., and Svendsen, A.B., Journal of Pharmaceutical Science, 51, 938 (1962).

2. Some decomposition was noted for ephedrine.


STANDARD SOLUTION AT 145°C



Compound	Retention Time Relative to Naphthalene
A Amphetamine	0.5
B Methamphetamine	0.7
C Naphthalene	1.0

FIGURE #2

AMPHETAMINE, METHAMPHETAMINE,
NAPHTHALENE AND EPHEDRINE AT 165°C



Compound	Retention Time Relative to Naphthalene
A Amphetamine	0.5
B Methamphetamine	0.7
C Naphthalene	1.0
D Ephedrine	1.8

Infrared Method for Distinguishing Optical Isomers of Amphetamine

James A. Heagy

U. S. Bureau of Narcotics and Dangerous Drugs, Box 36075, 450 Golden Gate Ave., San Francisco, Calif. 94102

THE NECESSITY to identify optically active drugs occurs frequently in forensic science. Quite often such drugs are identified as to the basic drug in question without determining the optical sign and indeed many drugs are found in both licit and illicit markets as only one isomer. Amphetamine and some other drugs, however, occur in the drug trade as *d*-, *dl*-, and *l*-isomers. A simple microcrystalline test (1) will distinguish *dl*-amphetamine from the *d*- or *l*-isomers but cannot distinguish *d*- from *l*-. The *d*- or *l*-isomers are distinguished by mixing the sample with the proper proportion of standard *d*- or *l*-amphetamine and observing a positive test for *dl*-amphetamine. The polarimeter will, of course, distinguish the isomers but this instrument is not available in many labora-

EXPERIMENTAL

A water solution of any amphetamine salt (10–50 mg) is made basic and the amphetamine extracted into methylene chloride. The methylene chloride is passed through anhydrous sodium sulfate into a small beaker and concentrated to ca. 2 ml by heating on a steam bath. A saturated solution of *d*-mandelic acid in methylene chloride is added several drops at a time until the amphetamine is neutralized as determined by a drop of solution on pH paper. The beaker is then covered for several minutes, allowing the *d*-mandelate salt to crystallize and the solution is filtered using suction and the crystals washed with a small portion of methylene chloride. After drying, a KBr disk of the crystals is prepared and the infrared spectra are run.




Figure 1. Infrared spectra of amphetamine *d*-mandelate salts in KBr between 800 and 600 cm^{-1}

tories. A gas-liquid chromatographic method has also been developed (2) using *N*-trifluoroacetyl-(*l*)-prolyl chloride to form diastereoisomeric derivatives with *d*- and *l*-amphetamine. This author has developed a simple method by which three distinct infrared spectra can be produced for *d*-, *dl*-, and *l*-amphetamine as the *d*-mandelate salts.

RESULTS AND CONCLUSION

The resulting spectra of the different isomers show several differences, the greatest being in the 800–600 cm^{-1} region as illustrated in Figure 1. The differences in the three spectra are certainly sufficient to distinguish *d*-, *dl*-, and *l*-amphetamine. This method has been used successfully on several samples of illicit amphetamine tablets. Only impurities precipitated by mandelic acid will interfere and none have been encountered in samples. Other optically active drugs can probably be distinguished using *d*-mandelic or other acids.

(1) Methods of Analysis of the A.O.A.C. 10th ed., 1965, William Horwitz, Ed., p 597.

(2) Clyde E. Wells, *J. Ass. Offic. Anal. Chem.*, 53, 113–115 (1970).

(3) Chouli, N. H., *J. Pharm. Sci.* 54, 1257 (1965)

RECEIVED for review April 30, 1970. Accepted July 15, 1970.

vol. ETP No. 4-51-

From Microgram
April, 1974

8

SEPARATION AND IDENTIFICATION OF AMPHETAMINE OR
METHAMPHETAMINE IN COMBINATION WITH EPHEDRINE
OR CAFFEINE

Samuel B. Stinson and Michael R. Berry
Crime Laboratory Analysts
Minnesota Bureau of Criminal Apprehension Laboratory
St. Paul, Minn.

Introduction

"White cross" tablets containing the following mixtures are being encountered in our laboratory: amphetamine with caffeine, amphetamine with ephedrine, and methamphetamine with ephedrine. The following rapid extraction procedure results in the isolation of the controlled substance and the subsequent identification by infrared spectroscopy.

Procedure

Grind up 1 tablet and place in a test tube. Add 0.5N NaOH and extract with an equal volume of hexane. Transfer hexane to a clean test tube. Wash the hexane 3 times with distilled water to remove the ephedrine or caffeine. Transfer the washed hexane layer to a clean test tube. Form the hydrochloride salt by bubbling HCl gas through the hexane. (The HCl gas can be obtained by withdrawing the vapors over conc. HCl with a disposable pipet equipped with a bulb). Precipitate the hydrochloride salt by centrifuging this mixture. Decant the hexane and dry the residue. Obtain the IR spectrum of the residue using a KBr pellet. (Chloroform may be used to transfer the residue to the mortar, followed by drying). Compare with standard amphetamine HCl (or methamphetamine HCl).

Remarks

This procedure is also suitable for routine analyses of suspected amphetamine tablets, but omitting the 3 water washes. In amphetamine samples containing relatively large amounts of caffeine, this one extraction procedure results in IR spectra clearly recognizable as being of amphetamine HCl (although slightly impure). Usually in the case of caffeine mixtures, one water wash is sufficient.

Property of Ohio State Police Forensic Services
Unconfidential Internal Copy
OBSOLETE DOCUMENT

NOTES ON THE TEST TUBE METHOD FOR
SEPARATION OF AMPHETAMINE OR METHAMPHETAMINE
FROM EPHEDRINE OR CAFFEINE

1/10/75

Microgram

Joseph G. Bibonic
Supervising Criminalist
Illinois Bureau of Identification
Pekin, Illinois

Introduction

Problems were encountered in the separation and crystallization of particularly methamphetamine when mixed with ephedrine or caffeine. The critical part of the procedure is to allow the test tube to set for 5 to 10 minutes after the hydrochloride salt forms. The following procedure with the minor changes is recommended.

Procedure

1. Follow the procedure of Stinson and Berry Microgram, Vol. VII, No. 4, April 1974.
2. For the usual white cross tablets - grind to powder one or two tablets, usually two. Shake NaOH and powder in corked test tube.
3. Add the NaOH, Hexane, distilled water with clean, individual pasteur pipets. Also use a clean pipet to transfer the Hexane layer.
4. Each time a shaking is necessary, use a corked test tube. Spin down each time using centrifuge.
5. After Hydrochloride Salt has been formed, spin down, set tube aside for 5 to 10 minutes to allow crystals to form on those spun to sides of test tube.
6. Pour off Hexane and allow the Hydrochloride Salt to dry in the test tube. A warm but not hot oven should be used.
7. KBR pellet can be formed by scraping the dried powder from the tube with a spatula.

Results

1. Resulting IR's are very sharp.
2. Procedure has been successful for suspected Amphetamine, Methamphetamine, MDA, and Methylphenidate. Other related compounds, which are soluble in Hexane, may also work.

Mini Bennies Purification
for I.R , DEA

1. Amphetamine
2. Desoxyephedrine
3. Diphenhydramine
4. Caffeine
5. Ephedrine
6. Meprobamate
7. Amphetamine, caffeine
- 7A. Acetaminophen, diphenhydramine
8. Aspirin, ephedrine
9. Brucine, ephedrine
10. Chloramphenicol, pemoline ~~1mg level~~
11. Caffeine, diphenhydramine
12. Caffeine, ephedrine
13. Caffeine, desoxyephedrine
14. Desoxyephedrine, phentermine (amphetamine)
15. Ephedrine, strychnine
16. Acetaminophen, caffeine, diphenhydramine
17. Amphetamine, caffeine, pemoline 1mg level
18. Amphetamine, caffeine, 1,3-diphenylisopropylamine
19. Caffeine, ephedrine, phentermine
20. Caffeine, ephedrine, erythromycin
21. Caffeine, ephedrine, thonzylamine
22. Caffeine, diphenhydramine, ephedrine
23. Acetaminophen, caffeine, diphenhydramine, ephedrine
24. Caffeine, diphenhydramine, ephedrine, thonzylamine
25. Caffeine, desoxyephedrine, ephedrine, phentermine

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

diphenhydramine follows...
work state... CH...

7. CHCl_3 wash from acid (caffeine) then CH_2Cl_2 wash from base followed by V (amphetamine).
- 7A. II (hexane diphenhydramine); I start with CH_2Cl_2 followed by CHCl_3 (acetaminophen).
8. Dry extract with CHCl_3 (aspirin) then $\text{CHCl}_3:\text{CH}_3\text{OH}$ 3:2 (ephedrine).
9. I petroleum ether (brucine), ether (brucine), $\text{CHCl}_3:\text{CH}_3\text{OH}$ (ephedrine). It may be necessary to put the brucine through II.
10. CH_2Cl_2 dry extract (chloramphenicol), CHCl_3 extract from acid solution evaporate to dryness, residue may have to be washed with ether to remove impurities (pemoline).
11. II hexane (diphenhydramine) then extract with CHCl_3 (caffeine).
12. I ether, CH_2Cl_2 and CHCl_3 (all caffeine), $\text{CHCl}_3:\text{CH}_3\text{OH}$ (ephedrine).
13. II hexane (desoxyephedrine) then extract with CHCl_3 (caffeine), i.e. desoxyephedrine precipitate out using 0.1N HCl in ether.
14. IV CHCl_3 (desoxyephedrine), column retains (phentermine).
15. Same as #9-brucine.
16. II hexane (diphenhydramine), I CH_2Cl_2 (caffeine) use 2nd CHCl_3 extract (acetaminophen).
17. III neutral (caffeine), basic (amphetamine), acid (pemoline). Pemoline residue may have to be washed with ether.
18. CHCl_3 extract from HCl solution, then CH_2Cl_2 extract from basic solution followed by V (amphetamine). Evaporate CHCl_3 extract then II hexane (1,3-diphenylisopropylamine), then extract with CH_2Cl_2 (caffeine).
19. II hexane (phentermine), I ether (caffeine), $\text{CHCl}_3:\text{CH}_3\text{OH}$ (ephedrine).
20. I erythromycin is usually present as the propionate which comes out with the petroleum ether extract (use larger volume), then proceed as in 19.
21. Same as #19 hexane (thonzylamine).
22. Same as #19 hexane (diphenhydramine).
23. I CH_2Cl_2 (caffeine), use 2nd CHCl_3 extract (acetaminophen), $\text{CHCl}_3:\text{CH}_3\text{OH}$ (ephedrine). Use TLC to locate diphenhydramine-varies from petroleum ether to CH_2Cl_2 extract. Use II - hexane (diphenhydramine) on this extract (evaporated) or on another portion of tablets.

24. Use I as in #23 and II on appropriate extracts (evaporated) for separation of diphenhydramine and thonzylamine.
25. I for ephedrine, caffeine -II for phentermine, desoxyephedrine, then IV on hexane extract.

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

MINIBENNIES
PURIFICATION FOR IR

Most of the separation techniques listed are used only to extract an ingredient for IR in a reasonably pure form and are not quantitative.

SEPARATION I

Finely ground tablet is packed in small column and dry washed with petroleum ether, then ether, CH_2Cl_2 , CHCl_3 , $\text{CHCl}_3:\text{CH}_3\text{OH}$ 3:2, CH_3OH . Use about 2ml each. TLC on each fraction with iodoplatinate spray to locate separated unknowns.

SEPARATION II

Ground tablet material is partitioned between basic water and hexane (3ml each). Hexane layer is washed thoroughly with 5ml water.

SEPARATION III

Ground tablet mixed with water and extracted with CHCl_3 from neutral, then basic, then acid solution.

SEPARATION IV

Ground tablets are mixed with 4N HCl celite, or hexane from separation II is poured onto 4N HCl celite column. Column is eluted with CHCl_3 .

SEPARATION V

A basic water solution of the tablets is extracted with CH_2Cl_2 . The CH_2Cl_2 is dried through Na_2SO_4 and d-mandelic acid in CH_2Cl_2 is added. The volume is reduced by boiling and amphetamine mandelate precipitates out. Too much mandelic acid will prevent the amphetamine from precipitating. Most other amines will not precipitate.

36.528 *Characteristics of Microchemical Tests for Alkaloids and Related Amines—Continued*

Alkaloid	Reagent	Description of Crystals
Narcoine (194)	Iodine potassium iodide, or zinc potassium iodide Platinic chloride	1:400. Blue, radiating needles, sometimes with yellow dichroism. Beautiful feathery rosettes develop in all solns.
Nicotine (204)	Mercuric chloride Mercuric chloride-sodium chloride	Radiating, transparent blades form in presence of slight excess of H ₂ SO ₄ ; feather-like blades form in presence of HCl. Radiating, transparent blades.
Noscapine (194) (Narcotine)	Potassium hydroxide or ammonium hydroxide	1:200. White, amorphous ppt that crystallizes slowly; dense rosettes of needles.
Papaverine (205)	Zinc chloride	Thin, rectangular plates in excess HCl.
Physostigmine (206)	Lead iodide Gold bromide in HCl	1:100. Radiating, serrated plates. 1 mg in 1 drop H ₂ O. Brown, dendritic aggregates.
Pilocarpine (192)	Platinic chloride	Crystals form slowly; layers of thin, yellow, triangular plates of delicate structure.
Procaine (205)	Platinic chloride Gold chloride and HCl	Spherical crystals of radiating branches. Irregular, radiating branches.
Quinidine (197)	Potassium iodide	Small, triangular crystals in great numbers; best in 1:1000 diln; sol. in excess reagent.
Quinine (197)	Disodium phosphate	Silvery, sheaf-like crystals.
Racephedrine (207) (<i>dl</i> -Ephedrine)	Bismuth iodide in dil'd sulfuric acid	1:200. Large orange plates and red prisms and grains.
Scopolamine (203) (Hyoscine)	Gold chloride	Clusters of pale yellow, transparent blades, with coarse, saw-toothed edges form immediately on shaking slide. Crystals grow to large size in 1:200 soln.
Sparteine (204)	Gold chloride	Large numbers of blade-like crystals varying in size according to concn.
Strychnine (208)	Platinic chloride Potassium cadmium iodide	Crystals form immediately in clusters and singly in small, wedge-shaped needles that move about field. Silvery masses, slowly forming rosettes.
Yohimbine (189)	Sodium carbonate	In 1:1000 soln heated to 50°. Fine needles in sheaf-like bundles and rosettes.

For Barbiturates (209)—Official Final Action
(See also 36.536.)

36.529 *Reagent*
Iodine-potassium iodide soln.—Dissolve 5 g I and 80 g KI in enough H₂O (ca 78 ml) to make 100 ml. Dil. with 2 parts by vol. of H₃PO₄. Prep. dil'd reagent every 2-3 weeks.

36.530 *Identification*
Dissolve little barbiturate in drop H₂O on slide. If present as Na salt, it dissolves readily; if present as acid, add little droplet 1% NaOH on stirring rod and mix. Add 1 full drop reagent and let stand until crystn occurs (immediate with some compds, 0.5-1 hr with secobarbital). Free acid may ppt or crystallize. I reaction crystals are easily distinguished by their color, often coupled with strong dichroism. Det. birefringence with polarizing microscope over glass is usually not needed but may be used for observation at high magnification and when slide stands >1 hr; on

standing, KI may crystallize as square, colorless, isotropic crystals.

Note crystals formed and compare characteristics with description, 36.531.

36.531 See pages 712-713.

For Sympathomimetics (207)—Official Final Action

36.532 *Reagents*
(a) Bismuth iodide in diluted sulfuric acid soln.—See 36.525(e).

(b) Gold chloride in diluted phosphoric acid soln.—Dissolve 1 g H₂AuCl₄·3H₂O in 20 ml H₃PO₄ (1+2). H₃PO₄ 1 part concn H₃PO₄ 2 parts H₂O

(c) Platinic chloride in diluted phosphoric acid soln.—Dissolve 1 g H₂PtCl₆·6H₂O in 20 ml H₃PO₄ (1+3). 1 part concn H₃PO₄ 2 parts H₂O

(d) Sodium tetraphenylboron soln.—Aq. soln (1+20).

Continued

of base. Stir, add 1

shaped forms.

apod plates; some

edge, resembling

stir, and add small

black rectangular

ing to large size.

cinchonine.

odles, form imme

more burr-shaped

ing heavier in struc.

rosettes of irregu

slowly in yellow

ng.

develop into dense

plates.

0 add minute frag

ms; also in aggre

Interlacing

ng later.

ends; singly and

ndles form slowly

standing.

d ends and united

standing.

ke slide to start

errated edges. In

brown plates with

s rapidly in neut.

ystals.

lar plates, after

00 soln. Shake

esses of fine net

brown ppt, slowly

as extending "

(Continued)

barbiturate Metharbital (5-Diethyl-1-methyl-5-barbituric acid)	Dark needles, small to large, and splinter blades.	Dichroism black to brown.	Good birefringence with crossed nicols.
5-Methyl-5-phenylbarbituric acid (Rutonal®)	Red-brown irregular platy forms appear after free acid is pptd.	Gradually strongly dichroic rods or blades.	Test fairly sensitive for dil. soln.
Phenobarbital (5-Ethyl-5-phenyl-5-barbituric acid)	Soon crystallizes in little dark grains; also a few larger red blades and dark splinter-rods in clusters.	—	Free acid may also crystallize out.

36.531 Characteristics of Microchemical Tests for Barbiturates—Continued

Barbiturate	Crystal Form	Dichroism or Pleochroism	Remarks
Probarbital (5-Ethyl-5-isopropyl-5-barbituric acid)	Scattered iodine-reaction crystals form in various jagged shapes, color dark brown to black dichroism, or red-black with but little dichroism.	Free acid thrown out, forming long rods with pointed ends.	
Secobarbital (5-Allyl-5-(1-methyl-butyl) barbituric acid)	Crystallizes in plates or elongate and rectangular but mostly distorted into any shape after 1 hr.	Light yellow to orange or red dichroism by polarized light.	Distinctly birefringent.
Sodium Pentobarbital (Sodium 5-ethyl-5-(1-methylbutyl) barbiturate)	Crystallizes quickly in great numbers of small red-brown plates.	Minute light-colored flakes exhibit dichroism; dark brown at black to yellow.	—
Talbutal (5-Allyl-5-sec-butylbarbituric acid)	Amorphous ppt crystallizes in large needles and dichroic blades, lighter to deeper brown in dendrites; then gray-black curled sheaves of threads.	—	Excellent test. Both types of crystals have good birefringence.
Vinbarbital (5-Ethyl-5-(1-methyl-1-butenyl) barbituric acid)	Multitudes of small dark crystals, tiny grains and rods with dichroism brown to black, in quite dil. soln possible to get good small crystals, little dark rods with dichroism red to black, and small plates tending to be square, generally appearing red but with same red to black dichroism, and with square extinction (not diagonal).	—	Very sensitive.

* This drug has barbiturate-type formula (although there is only one N) but is central nervous stimulant instead of depressant.

36.533

Identification

(a) *Direct test*.—Add drop of reagent to little of powd solid or crushed tablet and spread out on slide with little stirring. Do not stir to homogeneity as local concns and dilns will assist crystn. Let stand to evap. to higher acid concn if necessary for crystal formation.

(b) *Volatility test*.—Place small amt of substance or crushed tablet in depression of cavity slide, add drop 5% NaOH soln, and stir briefly. Place very small drop of reagent on thin slide, invert over cavity slide, and let stand. As crystals appear, examine with inverted slide in place. After observing crystals or after 1 hr or more exposure, if only few or no crystals form, reinvert thin slide with hanging drop, and let stand for gradual evapn of H₂O from reagent drop. Examine for crystals. Compare with descriptions, 36.534.

36.534

See page 714.

For Synthetics—Official Final Action

36.535

Reagents

- Acetic acid*.—Dil. 6 ml HOAc to 100 ml with H₂O.
- Ammoniacal nickel acetate soln*.—Mix 1 vol. 5% Ni(OAc)₂·4H₂O soln with 1 vol. NH₄OH (2 + 3). Use clear supernatant.
- Ammoniacal silver nitrate soln*.—See 36.525(a).
- Ammonium thiocyanate soln*.—See 36.525(c).
- Barium hydroxide soln*.—Satd aq. soln.
- Benzaldehyde*.—NF quality.
- Bismuth iodide soln*.—See 36.525(d).
- Bromide-bromate soln*.—Dissolve 0.3 g KBrO₃ and 1.2 g KBr in H₂O, and dil. to 100 ml.
- Glycerol-alcohol mixture*.—(1 + 1).
- Gold bromide in hydrochloric acid soln*.—See 36.525(g).
- Gold chloride soln*.—See 36.525(h).
- Iodine-potassium iodide soln*.—See 36.525(j).
- Lead acetate soln*.—Dissolve 5 g Pb(OAc)₂·3H₂O in H₂O and dil. to 100 ml.
- Lead triethanolamine soln*.—Add 1 ml triethanolamine (tech. 90% is satisfactory) to soln of 1 g Pb(OAc)₂·3H₂O in 20 ml H₂O. Slight turbidity does not interfere.
- Magnesia mixture*.—Dissolve 5.5 g MgCl₂·6H₂O and 14.0 g NH₄Cl in H₂O. Add 13.05 ml NH₄OH and dil. to 100 ml with H₂O.
- Mercuric chloride soln*.—See 36.525(l).
- Mercurous nitrate soln*.—Dissolve 15 g HgNO₃·H₂O in mixt. of 90 ml H₂O and 10 ml HNO₃ (1 + 9). Store in dark, amber bottle contg small globule of Hg.
- Nitric acid*.—(1 + 1).

36.534

Characteristics of Microchemical Tests for Sympathomimetics

Sympathomimetic	Reagent	Test	Description of Crystals
Volatile Substances			
<i>dl</i> -Amphetamine	Gold chloride in dilid phosphoric acid	direct or volatility	Very irregular plates, with irregular blade-arms especially after evapn; square if perfect.
	Platinic chloride in dilid phosphoric acid	volatility	Irregular blades and needles, very low birefringence; after evapn, characteristic plates with narrow irregular arms of blades.
<i>d</i> -Amphetamine	Gold chloride in dilid phosphoric acid	direct or volatility	Long yellow rods and blades; with evapn, some crystals as with <i>dl</i> may form.
	Platinic chloride in dilid phosphoric acid	volatility	Long needles, often bent, very little birefringence, after some evapn, long rectangular blades. (<i>d</i> -Ephedrine in direct test gives similar crystals which are more sol.; <i>l</i> is less volatile and does not normally form crystals in hanging drop.)
Epinephrine	Sodium tetraphenylboron	volatility	MeNH ₂ liberated; birefringent X's or 4-arm crystals; also thick blades with central rib, pointed ends, positive elongation.
Isoproterenol	Sodium tetraphenylboron	volatility	Isopropylamine liberated; plates tending to non-regular hexagons; no birefringence where plates lie flat but there are rods which are birefringent.
<i>d</i> - and <i>dl</i> -Methamphetamine (<i>d</i> - and <i>dl</i> -Desoxyephedrine)	Gold chloride in dilid phosphoric acid	direct or volatility	Long blades and jointed crystals, fairly high birefringence.
	Platinic chloride in dilid phosphoric acid	volatility	Grains with sharp edges which aggregate in chains and short prisms. Birefringent.
<i>d</i> -Methamphetamine	Bismuth iodide in dilid sulfuric acid	volatility	Drops, long orange splinters, blades, needles; also deep red angular grains (red prisms only after evapn).
<i>dl</i> -Methamphetamine	Bismuth iodide in dilid sulfuric acid	volatility	Drops, crystg in orange-red prisms with conspicuously slanting ends; inclined extinction ca 20°; also "mossy" formation of grains and some large deep red grains.
Slightly Volatile Substances			
<i>dl</i> -Ephedrine (racephedrine)	Gold chloride in dilid phosphoric acid	direct or volatility	Irregular plates based on the square growing along diagonals in 4 arms; some birefringent, some not.
	Bismuth iodide in dilid sulfuric acid	volatility	Orange rods or sticks, short and stubby, some plates; more irregular plates on evapn.
<i>l</i> -Ephedrine	Gold chloride in dilid phosphoric acid	direct or volatility	Long needles or splinters and long jointed forms, strong birefringence.
	Bismuth iodide in dilid sulfuric acid	volatility	Long brownish-orange needles, often branching or in sheaves; also, especially with evapn, orange irregular blades.
Pseudoephedrine	Gold chloride in dilid phosphoric acid	direct or volatility (2 hr)	Thin branching sticks, many like combs; some broaden to blades or spear-head plates; very high birefringence.
Phenylpropanolamine	Gold chloride in dilid phosphoric acid	direct	Plates and blades of extremely high birefringence, elongate hexagonal or diamonds, very bright colors. Branch into 4 or 6 irregular arms.
		volatility (2 hr)	After definite drying, pyramidal grains to blades and plates with irregular arms, very birefringent.
Phenmetrazine	Gold chloride in dilid phosphoric acid	direct or volatility	Rectangular plates joined in jagged arms of strongly birefringent crystals, often in X forms, very characteristic.
	Bismuth iodide in dilid sulfuric acid	volatility	Orange-red blades, usually pointed ends, often in rosettes; also with needles in branching aggregates; also red prisms.

(e) Phosphotungstic acid
 $P_2O_5 \cdot 24WO_3 \cdot xH_2O$
 (t) Picric acid
 (u) Picrolonic acid
 1-(*p*-nitrophenyl)-3-naphthol.
 (v) Platinic chloride
 (w) Potassium dichromate
 36.525(o).

36.536

Synthetic

Acetanilid (210)

Acetophenetidin (210)

Acetylsalicylic acid (211)

Aminopyrine (212)

Amobarbital (201)
Amytal®

Antipyrine (213)

Barbital (201)

Benzoic acid (211)

Cinchophen (214)

Callybarbituric acid (215)

p-Nitrophenol (212)



STATE OF IDAHO

DEPARTMENT OF HEALTH AND WELFARE
BUREAU OF LABORATORIES, 2220 Old Penitentiary Road, Boise, ID 83702

Procedure often followed in the Central Lab for the analysis of methamphetamine, phentermine, etc.

A. Fore preliminary analysis:

- 1. Spot tests: Marquis
Feigl's secondary amine test (for methamphetamine)
- 2. TLC: Same extraction and solvent system as amphetamine.
Fluorescamine: No reaction (may appear as dark spots under shortwave UV).
Iodoplatinate: May react with methamphetamine, very weak reaction if any with phentermine.
Dragendorff's: Should react with methamphetamine, (within one hour in small amounts) very weak reaction if any with phentermine.

3. GLC: 10% Apiezon L + 2.5% KOH.

4. Gold chloride in phosphoric acid direct microcrystalline test.
Saturated phosphoric acid direct on dry methamphetamine - gives squares and rectangles

B. For a final report:

- 1. IR: methodology we have used is mentioned separately.
- 2. More GLC: OV 17, 3% 3 ft. column, with and without acylation (take up sample into syringe and also 20% acetic anhydride in ethyl acetate, pyridine, or dimethylformamide).
- 3. If needed, quantitation on home-made methamphetamine can be done on the apiezon column. Use UV only for commercially manufactured preparations.

Property of Idaho State Police Forensic Services
UNCONTROLLED INTERNET COPY
PLEASE DELETE DOCUMENT

SF DEA
41972

DESORXYEPHEDRINE

Identification

- A. Infrared Spectrum - usually as the hydrochloride salt recrystallized from methylene chloride-petroleum ether. Different spectra are obtained for d and dl isomers.
- B. Microscopic and color tests - Platinic chloride (volatility), gold chloride (direct and volatility) are used. Optical Crystallography used by a minority of chemists.

Assay

- A. Blue-violet assay by direct dilution or shakeout.

Optical Activity

- A. Polarimeter - on concentrated solution.
- B. Melting Point for dl-isomer
- C. Crystal test - bismuth iodide

Diluents

- A. Infrared spectrum on chloroform-insoluble portion.
- B. Optical Crystallography

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

Differentiation of d-methamphetamine from dl-methamphetamine using gold chloride in H₂PO₄ directly.

d-methamphetamine

1. should have a number of good-sized crystals

shaped  or 

2. should have some plates shaped



,  etc.


no perfect rectangles

3. Not many X's consisting of larger plates

dl-methamphetamine


1. should have a number of good-sized crystals shaped

 these fill to  which may finally fill to

 the area shown by the dotted lines may be distinguished under crossed-polars.

2. should have some longer crystals which are perfect and also some crystals which are not perfect rectangles.

3. tends to form larger plates in X's:



Infrared Method for Distinguishing Optical Isomers of Amphetamine

James A. Heagy

U. S. Bureau of Narcotics and Dangerous Drugs, Box 36075, 450 Golden Gate Ave., San Francisco, Calif. 94102

The necessity to identify optically active drugs occurs frequently in forensic science. Quite often such drugs are identified as to the basic drug in question without determining the optical sign and indeed many drugs are found in both licit and illicit markets as only one isomer. Amphetamine and some other drugs, however, occur in the drug trade as *d*-, *dl*-, and *l*-isomers. A simple microcrystalline test (1) will distinguish *dl*-amphetamine from the *d*- or *l*-isomers but cannot distinguish *d*- from *l*-. The *d*- or *l*-isomers are distinguished by mixing the sample with the proper proportion of standard *d*- or *l*-amphetamine and observing a positive test for *dl*-amphetamine. The polarimeter will, of course, distinguish the isomers but this instrument is not available in many labora-

EXPERIMENTAL

A water solution of any amphetamine salt (10-50 mg) is made basic and the amphetamine extracted into methylene chloride. The methylene chloride is passed through anhydrous sodium sulfate into a small beaker and concentrated to ca. 2 ml by heating on a steam bath. A saturated solution of *d*-mandelic acid in methylene chloride is added several drops at a time until the amphetamine is neutralized as determined by a drop of solution on pH paper. The beaker is then covered for several minutes, allowing the *d*-mandelate salt to crystallize and the solution is filtered using suction and the crystals washed with a small portion of methylene chloride. After drying a KBr disk of the crystals is prepared and the infrared spectra are run.




Figure 1. Infrared spectra of amphetamine *d*-mandelate salts in KBr between 800 and 600 cm^{-1}

tories. A gas-liquid chromatographic method has also been developed (2) using *N*-trifluoroacetyl-*l*-prolyl chloride to form diastereoisomeric derivatives with *d*- and *l*-amphetamine. This author has developed a simple method by which three distinct infrared spectra can be produced for *d*-, *dl*-, and *l*-amphetamine as the *d*-mandelate salts.

RESULTS AND CONCLUSION

The resulting spectra of the different isomers show several differences, the greatest being in the 800-600 cm^{-1} region as illustrated in Figure 1. The differences in the three spectra are certainly sufficient to distinguish *d*-, *dl*-, and *l*-amphetamine. This method has been used successfully on several samples of illicit amphetamine tablets. Only impurities precipitated by mandelic acid will interfere and none have been encountered in samples. Other optically active drugs can probably be distinguished using *d*-mandelic or other acids.

RECEIVED for review April 30, 1970. Accepted July 15, 1970.

(1) Methods of Analysis of the A.O.A.C., 10th ed., 1965, William Horwitz, Ed., p 597

(2) Clyde E. Wells, *J. Ass. Offic. Anal. Chem.*, **53**, 113-115 (1970).

(3) Chouliis, N. H., *J. Pharm. Sci.* **54**, 1367 (1965)

COCAINE

Background Information-

DEA Booklet for purification process, use, etc.

Screening tests-

Cobalt thiocyanate 2% aq - (Reference: IRS Methods of Analysis)
Modified by using stannous chloride to remove false positives
caused by procaine. But cocaine base does not react, so add HCl to
all negatives.

False positive with phenazocine^{ZONE} (antipyrine) (1)
Microgram Vol 6, No 1 pg 14 (Jan 1973)

Richard Ruybal's modification of cobalt thiocyanate reagent (2)
Microgram Vol 5, No 3 pg 28 (March 1972)

Carolyn Ruybal's Modification of the above reagent (3)
Microgram Vol 6, No 2 pg 28 (Feb 1973)

Scott's Test (4)
Microgram Vol 6, No 11, pg 179 (Nov. 1973)

Lorch advocates FPN, for phenothiazine, in addition to Scott's (5)
Microgram Vol 7, No. 8, pg 100 (Aug 1974)

Lorch: Non-phenothiazine false positives with Scott's test (6)
Microgram Vol 7, No 11 page 129 (Nov. 1974)

Prall: List of false positives with Scott's test, detected with Marquis (7)
Microgram Vol 8, No 9 pg 130 (Sept 1975)

Bleach test (8)
Microgram Vol 7, No 6 pg 68 (June 1974)

Methylbenzoate test (9)
Microgram Vol 8, No 1, page 10 (Jan 1975)
Also Bulletin on Narcotics Vol 27, No 2 (Apr-June 1975)

For Procaine: Sanchez test - positive with procaine or benzocaine (10)
Microgram Vol 5, No 5 pg 51

Pocatello also uses Liebermann's, Mayer's and Wagner's for general screening.

TLC:

Comments

- A. Reference: Microgram Vol 5, No 5, May 1972, Page 51 (10)
- | | | |
|--------------------|-----|----------------------------------|
| (1) Cyclohexane | 5 | procaine, cocaine,
benzocaine |
| Chloroform | 4 | |
| Diethylamine | 1 | |
| (2) Ethylacetate | 6 | procaine, cocaine
benzocaine |
| Benzene | 3.5 | |
| Ammonium hydroxide | 1 | |
- B. Reference: Microgram Vol 7 No. 10 pg 122 (Oct 1974) (11)
- | | | |
|--------------------|----|---|
| Ethyl acetate | 60 | Paper gives Rf's:
cocaine, tetracaine
procaine lidocaine,
piprocaine, benzocaine |
| Dimethyl formamide | 20 | |
| Dioxane | 10 | |
| Acetic Acid | 10 | |
- C. Reference Microgram Vol 7, no 10 page 124 (Oct 1974) (12)
- | | | |
|-------------------|---|---|
| Methylcyclohexane | 3 | cocaine, lidocaine
(Apply as base, not as
HCl salt) |
| CHCl ₃ | 5 | |
- D. Reference Microgram Vol 8, No 1, page 4 (Jan 1975) (13)
- Spray plate with 0.5 N NaOH and dry at 140° C
Develop in methyl acetate
(Can be done without spraying plate but
tailing occurs)
- E. Reference Analytical Manual by Sobol and Moore ("The White Book")
- | | | |
|--------------------|----------|------------------------------------|
| 1. Chloroform | 25 | |
| Dioxane | 60 | |
| Ethyl acetate | 10 | |
| Ammonium hydroxide | 5 | |
| 2. Chloroform | 40 | Separate cocaine,
procaine well |
| Ethyl acetate | 10 | |
| Ammonium hydroxide | 10 drops | |
- F. Chloroform 9
MeOH 1
Separates cocaine
procaine. slight
streaking
- G. Chloroform saturated
with NH₄OH 18
MeOH 1
Separates cocaine,
procaine well, high
Rf.
- H. Cyclohexane 80
Benzene 16
Isopropylame 4
Supposed to be good
for cocaine, procaine,
PCP
- I. Clarke's T₁ See picture (14)
- J. 4% diethylamine in Benzene " "

Microcrystal tests

1. Platinic chloride - 5% aq. (Reference: Methods of Analysis of AOAC Eleventh Edition 1970 pg 709, 710). Gives feathers. Does not work if lidocaine present, or some other diluents.
2. Gold chloride - 5% aq. - dissolve powder in 0.1 NHC1, add drop of 5% gold chloride - floating feathers.
Ref. EGC Clarke - Isolation and Identification of Drugs
Note: Pocatello reminds us l-cocaine and d,l-cocaine give different microcrystal forms. (15)
3. Gold bromide - Reference: Modern Microcrystal Tests for Drugs - Fulton, page 385 #4 (See Pam's write up for recipe in plain English). Almost immediately a large number of irregular blades form which appear green with low birefringence. After standing 1-2 minutes see small x's, usually in a different plane (focus down). This test works even in the presence of lidocaine. (Sometimes called Ruybal's Gold Bromide because Richard Ruybal first described the use of this reagent on cocaine at the American Academy meeting in Dallas). (16)
4. Sodium acetate: 20% aq. (Ref: "Isolation and Estimation of Cocaine in the Presence of Tetracaine, Procaine, Benzocaine, and Amphetamine" by Milos and Porto, Alcohol and Tax Laboratory, IRS, New York). Supposed to form cocaine base crystals, unaffected by procaine, tetracaine, benzocaine, and amphetamine. (17)
5. Ecgonine with phosphomolybdic acid (1 gram phosphomolybdic acid in 20 ml of 1:1 HNO₃). Reference-same as #4 above. See paper for a discussion of the test. (17)

Note: I am unaware of anyone having used #4 and #5 above, and so am not sure if they work. they are listed for sake of completeness and in case you are desperate.

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

GLC -

Microgram Vol IV No. 3 pg 29. (March 1971)
6' 3% OV-1 on Chromosorb WHP at 210°C

(18)

Microgram Vol IV No. 6 pg 60 (July 1971)

Procaine, cocaine, PCP 6' 3% OV-17 Chromosorb WHP at 230°

(19)

Microgram Vol V No 5 pg 51 (May 1972)

Procaine, cocaine, benzocaine 3% OV-1, OV-25, 210° 140° respectively

(10)

Microgram Vol V No 12 pg 140 (Dec 1972)

Procaine, cocaine, benzocaine

Internal std of docosane, dissolve sample 9/1 (CHCl₃:MeOH)
3% OV-17, 6' Chromosorb WHP, 210°

(20)

Microgram Vol 9 No 2 pg 18 (Feb. 1976)

Table of retention values on OV-17, 6' on Gas Chrom Q

(21)

Journal of Chromatographic Science Vol 9, July 1971 page 393

Finkle's Table of Relative Retention times on SE-30, 2' and 6'

Journal of Pharmaceutical Sciences Vol 63, No 12 Dec. 1974 pg 1963.

Derivatization of Cocaine w/trimethylanilinium hydroxide.

(22)

Currently most people shoot on SE-30 (or OV-1 or -101) and/or OV-17, at temperatures appropriate to column length and GC performance.

Example: (a). 3' OV-17 at 205° (Boise prefers to SE-30, OV-1 or -101 because less tailing) OV-17

(b) 2' 3% OV-17 at 225°C and 250°C used by Pocatello

For detection of hydrolysis by-products benzoylecgonine and ecgonine
Journal of Chromatography 101 (1974) pg 215-218.

(23)

Property of Idaho State Police Forensic Services
OCCUPATIONAL DOCUMENT

Microgram Vol 2, No 2 pg 47 June 1969 (27)
Cocaine ion-pairing extraction from CHCl_3

Microgram Vol 5, No 7 July 1972 pg 80
Ion pairing - use phase separating paper

Microgram Vol 3 No 3 May 1970, pg 89, and Vol III No. 4 June 1970, pg 121
Celite column:cocaine and procaine - requires sample containing
5-15 mg cocaine (26 and 27)

Microgram Vol 4, No 3, March 1971, pg 29
Celite column:cocaine, caffeine and procaine (18)

Microgram Vol 5, No 5, May 1972 pg 49
Celite column:cocaine, procaine, benzocaine

Microgram Vol 7, No 6 page 70 (June 1974)
Cocaine extracted from NaHCO_3 with CCl_4 (28)

Microgram Vol 7 No 8, page 96 (Aug. 1974)
Clark's & Allen's paper on cocaine - lidocaine (29)
by Celite column. Rick says this did not work for him.

- I. Quick cleanup if cocaine HCl + sugars; or cocaine-procaine mix
(Does not work if lidocaine present)

Make a column using a pasteur pipette (unconstricted) with glass wool and place a portion of the sample in it. Run CH_2Cl_2 through the samples fast. If you wish, dry the CH_2Cl_2 through a separate column of Na_2SO_4 . Put the CH_2Cl_2 on a steam bath to reduce its volume, and precipitate using pet ether. Dry crystals and run on KBr disc. (This method depends on the fact that procaine HCl is not very soluble in CH_2Cl_2 esp. in a short time).

- II. Acid-base shake out.
As noted in the early microgram papers, cocaine-HCl will ion-pair, that is it will be extracted into CHCl_3 if the HCl concentration is 2.8N to 4N. If you shake it out on the basic side, it's best to use NaHCO_3 or Na_2CO_3 to avoid hydrolysis in very basic solutions. After shaking with CHCl_3 or CH_2Cl_2 , dry the organic phase through Na_2SO_4 , reduce its volume and dry directly on KBr. Do not try to make the HCl salt by adding HCl and MeOH because of risk of hydrolysis.

- III. Cobalt-thiocyanate derivative for cocaine+lidocaine (30)
(Reference: Naylor et al from Bullet - Sn of Midwest Association).

- IV. Platinic chloride derivative for cocaine + lidocaine. (31)
Joe Powers' method.

- V. Spot from TLC plate extracted with dilute acid. Make basic keeping in mind the precautions in II, above.

- VI. Celite column - see U.V. section, reference #1. Follow directions given in that paper. Also see Fred's column chromatography notes. (32)

VII. "Lumpology"

If your cocaine is lumpy or contains 2 kinds of material which you can separate under the stereoscope try running an IR direct on the lumps (or whatever component your screening tests indicate may be cocaine).

VIII. Alumina column

As a last resort you might try alumina. Make a microcolumn of alumina and wet it with CH_2Cl_2 . Apply your sample in CH_2Cl_2 or ether. Wash with ether or CH_2Cl_2 to remove the contaminants (esp dyes). Elute using progressively more polar solvents. The cocaine will come off as the base, not HCl salt. Two cautions: different batches of alumina behave differently -- you may have to try different ones to get one that works properly. Also this procedure is far from quantitative; bad recoveries are to be expected.

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

UV-

1. Paper by Nakamura and Parker "Assay of Cocaine in the Presence of Procaine and Quinine by Column Chromatography"
Journal of Chromatography Vol 52 (1970) pg 107-110. (23)
2. Rick's Micro Adaptation of the above
 - A. Column prep: Mix 1 gram Celite + 0.5 ml of (1M KNO_3 in 0.1 N HCl) Put in column on top of glass wool, tamp down, and put glass wool over top.
 - B. Sample prep: Dissolve sufficient sample in 0.5 ml of (1M KNO_3 in 0.1 N HCl) to give absorbance of approximately 0.5A. (Is about 0.14 mg cocaine HCl/ml when 0.25 ml of solution put on column). Apply 0.25 ml of sample solution to column.
 - C. Elution: Use 9 ml of water-washed CHCl_3 on column. Catch in 10 ml volumetric containing 1 ml of MeOH and 1 drop of HCl; bring to volume with CHCl_3 . Scan on U.V.

Quantitation-

1. By external std on GLC
2. By internal std on GLC - have used methadone as internal std.--need to show none present originally.
3. By U.V. using celite column (if no lidocaine). Or Pocatello uses EtOH if no interfering substances present.
4. "Back seat" method - see separate sheet-- for narcotics officers. (34)

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

Isomer problem

- Microgram Vol 9, No 5 page 66 (May 1976) (35)
Review isomer problem, suggest micro X-tal test
- Microgram Vol 9, No 8, page 107 (Aug 1976) (36)
Offer of small amount of d,l cocaine.
- Microgram Vol 10, No 4, page 52 (April 1977) (37)
Use of salt of tartaric acid to give different IR's.
- Microgram Vol 11, No. 3 page 41 (March 1978) (38)
Preparative TLC and mixed melting points
- Microgram Vol 11, No 3 page 45 (March 1978) (39)
Philosophical discussion of wording of law - "derivative".

Also see latest Board of Pharmacy list. (40)

JOURNAL OF CHROMATOGRAPHY Vol 152 (1978) 589-591 - *TLC of derivatives separate isomers*

Scientific Sleuthing Vol 3 No 3.

Reference to People vs. Harper, 562 P2nd 1112 (Colo. 1977)

Colorado Supreme Court said their statute includes natural and synthetically produced cocaine. The chemists' testimony was admissible to identify cocaine even though his tests did not distinguish the two forms.

Stability problems

1. Discoloration effects of diluents
Journal of Forensic Sciences Vol 19, No. 4 (868-872) (41)
2. Stability in dilute H_2SO_4
Microgram Vol 6, No 9, page 135 (Sept. 1973) (42)

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

ABSTRACT OF A PAPER PRESENTED AT THE SPRING 1975 MEETING OF THE
MIDWESTERN ASSOCIATION OF FORENSIC SCIENTISTS.

"A Simple Procedure for the Separation and Identification of Cocaine".
by Jon D. Naylor, Carl R. Phillips, Robert J. McCurdy & Stephen A. Koers.

The selective extraction of the blue complex of cocaine with cobalt thiocyanate from aqueous acid into chloroform in the presence of common street excipients, e.g. lidocaine, tetracaine, procaine, etc., has suggested a very simple purification and identification procedure. The infrared spectrum of the dried blue extract has been found to be quite characteristic for cocaine. The spectrum consists of a simple addition of the absorption peaks for cocaine and cobalt thiocyanate with very few shifts in frequency or distortions of component peak shapes. The procedure is: a.) addition of about 2 ml of cobalt thiocyanate reagent (2g of cobalt thiocyanate in 100ml of water) to enough street sample to contain one to two mg. of cocaine, b.) addition of up to one-half ml. of concentrated hydrochloric acid (excess HCl results in displacement of cocaine by chloride ion in the complex to form a blue solution), c.) dropwise addition of enough water to dissolve all of the blue precipitate (vigorous shaking is necessary), d.) extraction of the blue cocaine-cobalt thiocyanate complex into chloroform, e.) drying the complex in an evaporating dish, and f.) running the infrared spectrum of the dried blue complex as a KBr mull. The resultant spectrum has been found to be totally free of other materials when mixtures of cocaine with procaine, lidocaine, and tetracaine were treated in the above manner.

If pure cocaine instead of the complex is desired, the chloroform extract in (d) above may be washed with aqueous ammonia to displace the cocaine from the complex. ^{shake!} The tan-colored ammonia complex will be extracted into the aqueous layer, and the clear chloroform layer will contain cocaine free base which can be treated in the usual manner for further identification. 2X

The greater speed, more complete separations, and ability to perform the entire operation in a test tube makes this procedure far superior to the traditional multiple extraction and chromatographic separation procedures.



DATE February 8, 1972

NO. 31

-28-

DRUG TYPE

METHODOLOGY

COLOR TEST TO DIFFERENTIATE BETWEEN
COCAINE AND PROCAINE

Richard Ruybal
Forensic Chemist
Bureau of Narcotics and Dangerous Drugs
Dallas Regional Laboratory

Using the present cobalt thiocyanate reagent, blue colors are obtained with both Cocaine and Procaine as well as with many other drugs. The Sanchez reagent gives a positive reaction with Procaine and other primary amines.

As screening field tests, especially to the agent in the field, the combination of the two reagents affords no real indication that Cocaine is present.

The following reagent is submitted for a color test to differentiate between Cocaine and Procaine in illicit samples. Other drugs were also tested and results are noted.

Stock Solutions

- A. Cobalt Thiocyanate - a 2% solution in H₂O
- B. Phosphoric Acid - Syrupy H₃PO₄
- C. Platinum Chloride - 1 gram H₂PtCl₆·6H₂O in 20 ml H₃PO₄ (1+3)

Reagent

Mix by volume 9 parts of A, 3 parts of B, and 1 part of C for test reagent.

Test

Place small amount of sample on a spot plate, add 3-4 drops of reagent and stir with glass rod.

DRUG


RESULTS

1. Cocaine.....Blue, flaky ppt. remaining undissolved.
2. Procaine.....No ppt. or blue color forms.
3. Benzocaine.....No ppt. or blue color forms.
4. Butacaine.....Green color forms, but fades away.
5. Dibucaine.....Slight green color and ppt. forms, but fades away.
6. Lidocaine.....Blue, flaky ppt. remaining undissolved.
7. Mepivacaine.....Blue ppt. forms, but fades away.
8. Tetracaine.....Blue color forms, but fades away.
9. Quinine.....Green ppt. and color remaining undissolved, but different from Cocaine.
10. Methapyrilene.....Green ppt. and color, but fades away.
11. Heroin.....Blue ppt. and color, but fades away.
12. Methadone.....Blue ppt. and color, but fades away.
13. Demerol.....Blue ppt. and color, but fades away.
14. Phencyclidine HCl (PCP).....Blue-green ppt. remaining undissolved.

The only similarity of any of the above drugs tested to Cocaine was Lidocaine. It was the only other one that formed a blue, flaky ppt. which remained unchanged. All others forming blue ppts. or blue colors lost all blue color within a few minutes leaving either a pink-orange or a yellowish color. Quinine formed a green ppt. and color which is quite easily discernible from Cocaine.

The drugs listed above were also tested with a mixture of Stock Solutions (A) and (B) in a ratio of 9.3. The drugs that gave a blue, flaky ppt. in addition to Cocaine and Lidocaine were Dibucaine, Mepivacaine, Heroin, Methadone, and PCP. Demerol gave a blue color.

BNDD LABORATORY NOTES



DATE

-28-

NO. 52

DRUG TYPE

METHODOLOGY

COLOR TEST TO DIFFERENTIATE BETWEEN COCAINE AND LIDOCAINE

Carolyn N. Ruybal
Forensic Chemist
Bureau of Narcotics and Dangerous Drugs
Dallas Regional Laboratory

The reagent for differentiation between cocaine and procaine described by Mr. Ruybal in BNDD Laboratory Notes in Microgram Vol. 5 No. 3 March 1972 is successful for differentiating cocaine from all the "caines" except lidocaine.

The following preparation of this reagent is submitted as a color test to differentiate between cocaine and lidocaine. It is simply a dilution of the "Ruybal" reagent 3 parts to 2 parts water. Lidocaine will not give a color reaction.

STOCK SOLUTIONS:

- A. Cobalt Thiocyanate - a 2% solution in H₂O.
- B. Phosphoric Acid - syrupy H₃PO₄.
- C. Platinum Chloride - 1 gm H₂PtCl₆-6H₂O in 20 ml H₃PO₄ (1+3)

REAGENT:

By volume mix 9 parts A and 3 parts B. Add 1 part C and mix well. Add 9 parts distilled water and mix. When these stock solutions are first mixed the reagent is clear orange in color. Let stand at least 5 to 7 days. A reddish brown precipitate is formed which settles to the bottom and the supernatant is a clear pink color. The reagent is then ready for use. The precipitate should not be removed or separated from the reagent.

TEST:

Place a small amount of sample on a spot plate, add 2-3 drops of reagent.

RESULTS:

Cocaine gives a blue flaky precipitate.

Procaine, benzocaine, butacaine, dibucaine, lidocaine, methapyrilene, heroin and demerol produce no color with this reagent. Methadone forms a slight blue color on standing. Phencyclidine HCl and antipyrine also give a blue color similar to cocaine.

Antipyrine (phenazone) can be differentiated from cocaine with Mandelin's reagent. Antipyrine gives an intense blue-green color whereas cocaine gives a slight red color which quickly fades.

DISCUSSION:

The first approach to modifying this reagent for cocaine-lidocaine differentiation was to dilute the original Ruybal Cobalt Platinum Reagent with a freshly made 5% solution of stannous chloride. But when distilled water was substituted for the stannous solution the reagent worked equally well.

When preparing this reagent it would appear that a 1% solution of cobalt thiocyanate solution in water could be used in lieu of the final dilution step. This was tried with no success. It is important to mix the reagent in the manner described.

When using the reagent for the test use only the pink portion (supernate) of the reagent with the precipitate present but settled to the bottom. Agitation of the reagent before the test will not affect its performance but is undesirable as the dispersed red particles of precipitate make it harder to see the color reaction of the cocaine.

Positive results were obtained with this reagent on samples containing as low as 5 percent cocaine.


The final dilution with water may be done either at the time of mixing the original solution, as described above, or later. Some of an aged original Ruybal cobalt platinum solution was diluted approximately 2 weeks later in the ratio 3 parts reagent plus 2 parts water. The final diluted mixture performed successfully for differentiating cocaine from lidocaine when used immediately

after dilution and when used again 3 weeks later. It is important when diluting the aged reagent in this manner to agitate well first to disperse the precipitate before measuring the aliquot to be diluted. The precipitate must always be present in the reagent.

An alternate approach to performing this test is to add 3 drops of the original Ruybal cobalt platinum solution to the sample on a spot plate. At this point both cocaine and lidocaine give a blue precipitate. Add 2 drops of water and stir. The lidocaine precipitate will dissolve. The cocaine precipitate will partially dissolve but the larger blue flakes of the cocaine precipitate remain.

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

DEA LABORATORY NOTES



DATE June 11, 1975

NO.

DRUG TYPE Cocaine

METHODOLOGY Color Tests (Field Test)

AN ADDITIONAL SCREENING TEST FOR COCAINE

Josef D. Prall
Forensic Chemist
South Central Regional Laboratory
Drug Enforcement Administration

OBJECTIVE

To develop a series of field test reagents for cocaine, readily available to federal, state, and local agents, which will screen out false positive tests with the Scott Field Test for cocaine (1).

BACKGROUND

A sample purchased as cocaine submitted to this laboratory for analysis gave a positive reaction to all three steps of the Scott Field Test for cocaine. Further analysis by gas chromatography/mass spectrometry and infrared spectroscopy showed the sample to be diphenhydramine. Lorch (2,3) has also encountered similar difficulties with suspected cocaine samples which were found to be phencyclidine and phenothiazine derivatives. Lorch suggested using a second test with the "FPN Reagent" (4). However, since the Scott test is primarily intended as a field test and the FPN Reagent is not readily available to agents in the field, an additional test using a reagent which is readily available was deemed advisable. It was found that the Marquis Reagent gives such a test. In addition to the false positives reported by Lorch, several additional substances were found to give false positives; these too were screened out by the Marquis test.

DRUG ENFORCEMENT ADMINISTRATION / U. S. DEPARTMENT OF JUSTICE

REAGENTS

(1) Scott Field Test (1)

Solution #1 2% cobaltous thiocyanate dissolved in water and then dilute 1:1 with 96% USP Glycerine (use 5 drops).

Solution #2 Concentrated hydrochloric acid (use 1 drop, 2nd drop if necessary -- no more).

Solution #3 Chloroform.

(2) Marquis Reagent -

A 4-5% solution of 40% formaldehyde in concentrated sulfuric acid, also packaged and sold commercially as "Marquis Reagent" used as a test for opium alkaloids, heroin, and amphetamines.

METHOD

The three step Scott Field Test is performed on a portion of the sample as previously described by Scott (1). If a positive result is obtained on all three parts, a second portion is added to Marquis Reagent. In general:

Positive Scott, Positive Marquis (Some Color Reaction)

- Indicates presence of some compound other than cocaine, although it does not rule out the presence of cocaine.

Positive Scott, Negative Marquis (No Color Reaction)

- Indicates presence of cocaine, and absence of the compounds listed below in Table 1 which give a positive Marquis Test.

A positive reaction to the Scott Field Test involves

- 154 -

a blue solution in Step 1, a pink solution in Step 2, and a blue extract in the chloroform layer in Step 3. It must be emphasized that the solution should be blue in Step 1; a pink solution with blue specks should not be taken as a positive for this step.

Based on this criterion, some of the compounds which Lorch reported as giving positive reactions to Step 1 were instead found to be negative. This may be a function of the particular form of the compound used for the test; i.e., free base, hydrochloride salt, or other salt. Note the reactions for the various forms of cocaine, diphenhydramine and doxylamine in Table 1. It might be expected that certain salts, especially the hydrochlorides, of other drugs might give positive Scott tests; for example, phenyltoloxamine, which is similar in structure to diphenhydramine.

Combinations of certain drugs may of course give the same results to both tests as cocaine hydrochloride; for example, a mixture of antipyrine and carbinoxamine maleate.

REFERENCES

1. Scott, L. J., Jr., Microgram, Vol. 6, No. 11, Nov. 1973
2. Lorch, S. K., Microgram, Vol. 7, No. 8, August 1974
3. Lorch, S. K., Microgram, Vol. 7, No. 11, November 1974
4. Clarke, E. G. C., "Isolation and Identification of Drugs", The Pharmaceutical Press, London, U. K. (1969).

TABLE 1

Steps 1 SCOTT 2 3

MARQUIS

Substance	Step 1	SCOTT 2	SCOTT 3	Reaction
Cocaine HCl	+	+	+	NR
Cocaine Base	-	+	+	NR
Chlorpromazine HCl	+	+	+	Pink to Red to Violet
Diphenhydramine HCl	+	+	+	Yellow to Orange to Brown
Diphenhydramine Base	SP	NC	+	Orange to Brown
Diphenhydramine Salicylate	NR	+	+	Yellow to Orange to Brown
Diphenhydramine Theoclate	NR	SP	+	Yellow to Orange to Brown
Doxylamine HCl	+	+	+	Yellow to Dirty Pink to Light Purple
Doxylamine Succinate	+	+	-	Yellow to Dirty Pink to Light Purple
Diphenylpyraline HCl	+	+	+	Yellow to Orange to Brown
Dibucaine	+	+	-	NR
Tetracaine HCl	+	+	-	NR
Carbinoxamine Maleate	+	+	-	NR
Promazine HCl	+	+	-	Pink to Red to Dark Red
Dimethindine	+	+	-	Purple
Pyrilamine Maleate	+	+	-	Purple
Homochlorcyclizine HCl	+	+	-	Pale Yellow to Pale Green
Methapyrilene	+	+	-	Purple
Tripeleennamine HCl	+	+	-	Orange-brown to Brown
Lidocaine HCl	SP	NC	+	Faint Pink
Phencyclidine HCl	SP	NC	+	NR
Phenyltoloxamine Citrate	SP	+	+	Wine Red
Promethazine HCl	SP	NC	+	Pink to Red to Violet
Pyrathiazine	Slow SP	NC	+	Pink to Red to Violet
Pyrathiazine HCl	Slow SP	NC	+	Pink to Red to Violet
Prochlorperazine Edisylate	SP	NC	-	Pink to Red to Violet
Pheniramine Maleate	SP	+	-	NR
Orphenadrine Citrate	SP	NC	+	Orange
Chlorcyclizine	NR	NC	-	Pale Yellow
Cyclizine HCl	NR	NC	NC	NR
Brompheniramine Maleate	NR	NC	NC	NR
Methoxyphenamine HCl	Pale Blue w/Lg.Amt.	NC w/5 drops HCl	+	Wine Red
Perphenazine	NR	NC	NC	Pink to Red to Violet

Property of Idaho State Police Forensic Services
 Uncontrolled Internet Copy
 OBSOLETE DOCUMENT

NR - No Color Reaction
 NC - No Change from Previous Step
 SP - Blue Specks in Pink Solution

TABLE 1 (Continued)

	Steps 1	SCOTT 2	3	MARQUIS
Diethylpropion	SP	NC	+	NR
Benzocaine	NR	NC	NC	NR
Procaine HCl	NR	NC	NC	NR
Antipyrine	NR	Flashes Blue then Pink	+	NR

NR - No Color Reaction
 NC - No Change from Previous Step
 SP - Blue Specks in Pink Solution

Property of Idaho State Police Forensic Services
 Uncontrolled Internet Copy
 OBSOLETE DOCUMENT

SPECIFICITY PROBLEM WITH THE COCAINE-SPECIFIC FIELD TEST
II. NON-PHENOTHIAZINE FALSE POSITIVES AND THE
SEPARATION OF PHENCYCLIDINE - PROMAZINE COMBINATIONS

Steven K. Lorch, Ph.D.
Crime Laboratory Scientist
Division of Crime Detection
Bureau of Laboratories
Michigan Department of Public Health
Lansing, Michigan 48914

OBJECTIVE

We have found that certain non-phenothiazine drug combinations will give the "highly specific" field test (1) for cocaine, and that the FPN test (2) will not screen for these false positives.

BACKGROUND

Since our previous report (3) our laboratory has received three more cases of phencyclidine (PCP) - phenothiazine combinations which gave the cocaine-specific field test color sequence and a positive FPN test for phenothiazines. The phenothiazine was identified as promazine by gas chromatography and gas chromatography-mass spectrometry. As reported earlier promazine gives steps one and two of the cocaine-specific field test (3).

The test was run on standard samples of the following drugs and combinations of drugs:

TABLE 1

	Color after step 1	Color after step 2	Color in CHCL₃ CHCL ₃	FPN
1) Phencyclidine	blue	pink	clear	-
2) Promazine	blue	pink	clear	+
3) Cocaine	blue	pink	blue	-
4) 1 + 2	blue	pink	blue	+
5) Dibucaine	blue	pink	clear	-
6) Methapyrilene	blue	pink	clear	-
7) 1 + 5	blue	pink	blue	-
8) 1 + 6	blue	pink	blue	-
9) 5 + 6	blue	pink	clear	-
10) 2 + 5	blue	pink	blue*	+
11) 2 + 6	blue	pink	blue*	+

*depends on the relative amounts of each component.

It is important to note that combinations #7 and 8 give a false positive for cocaine, but cannot be screened out with the FPN test, as neither compound is a phenothiazine.

SEPARATION OF PCP - PROMAZINE COMBINATIONS

Dissolve mixture in 2.8 N HCl and extract with an equal volume of chloroform. Wash chloroform two times with fresh 2.8 N HCl. Extract PCP from chloroform with an equal volume of 0.1 N sulfuric acid. A U.V. spectrum of the 0.1 N sulfuric acid extract will reveal a promazine peak at 252 m μ with a 269 m μ shoulder representing the PCP (promazine 252 m μ E 1% 1 cm 1122, PCP 269 m μ E 1% 1 cm 9.2). Add three drops of concentrated HCl to the sulfuric acid extract, extract with an equal volume of chloroform and dry. The resulting extract will give an IR spectrum clearly recognizable as that of phencyclidine HCl.

- (1) Scott, L.J., Jr., "Specific Field Test for Cocaine", Microgram, 6:11, Nov. 1973, pages 179-181.
- (2) Clarke, E.B.C., "Isolation and Identification of Drugs", The Pharmaceutical Press, London, U.K. (1969).
- (3) Lorch, S.K., "Specificity Problem With the Cocaine Specific Field Test, and Its Solution", Microgram, 7:8, Aug. 1974, pages 100-101.

The author wishes to thank Dr. Fathi M. Saad, Chief, Warren Regional Laboratory, Bureau of Laboratories, Michigan Department of Public Health, for the G.C.-Mass Spectrographic identification of the promazine.

9/10/74

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

SPECIFICITY PROBLEM WITH THE COCAINE SPECIFIC FIELD TEST,
AND ITS SOLUTION

Steven K. Lorch, Ph.D.
Crime Laboratory Scientist
Division of Crime Detection
Bureau of Laboratories
Michigan Department of Public Health
Lansing, Michigan 48914

OBJECTIVE

We have found that certain drug combinations will give the "highly specific" field test (1) for cocaine. An extra test is suggested for screening out these false positives.

REAGENTS

- Solution #1 (1) 2% cobaltous thiocyanate dissolved in water and then dilute 1:1 with 96% USP Glycerine (use 5 drops).
- Solution #2 (1) Concentrated hydrochloric acid (use 1 drop, 2nd drop if necessary -- no more),
- Solution #3 (1) Chloroform.
- FPN Reagent (2) 5 ml of 5% w/v solution of ferric chloride in water.
45 ml of 20% w/w perchloric acid
50 ml 50% w/w nitric acid

BACKGROUND

A purported cocaine sample submitted to this laboratory gave a positive field test for cocaine. However, further analysis showed it to be a combination of phencyclidine (PCP) and a phenothiazine; our first test, a UV spectrum, led us to believe it was promethazine, but we were unable to finally identify it as such.

The field test run on promethazine alone failed to give a blue color in solution 1, and a pink color in solution 2 but gave a blue color in the chloroform layer. Combinations of PCP, dibucaine or methapyrilene (which as reported earlier by Scott (1), gave blue with solution #1 and pink on addition of concentrated HCl), with promethazine gave the full cocaine specific color sequence (Table 1). Butacaine was not tested as we had no available standard. Two other phenothiazines tested, chlorpromazine and promazine, gave the proper color sequence when in combination with promethazine.

METHOD FOR DISTINGUISHING FALSE POSITIVES

We find we can improve test specificity by performing a fourth test. Add a few drops of FPN solution (2) to another portion of the sample. FPN indicates the presence of phenothiazines, and gives a short term orange, pink or red color, depending on which of the many phenothiazines it is. If one of these colors is observed, there is a strong indication of a combination of the nature described above, and further testing is indicated. FPN tests on cocaine, PCP, dibucaine and methapyrilene alone are negative.

It should be noted that the phenothiazine in the actual case was not promethazine, but probably another phenothiazine as yet unidentified.

TABLE 1

	Color after Step 1	Color after Step 2	Color in CHCl ₃	FPN
1) Cocaine	blue	pink	blue	-
2) Phencyclidine	blue	pink	clear	-
3) Dibucaine	blue	pink	clear	-
4) Methapyrilene	blue	pink	clear	-
5) Chlorpromazine	blue	pink	clear	+
6) Promazine	blue	pink	clear	+
7) Perphenazine	pink	pink	clear	+
8) Prochloroperazine	pink	pink	clear	+
9) Promethazine	pink	pink	blue	+
10) 1 + 9	blue	pink	blue	+
11) 2 + 9	blue	pink	blue	+
12) 3 + 9	blue	pink	blue	+
13) 4 + 9	blue	pink	blue	+
14) 5 + 9	blue	pink	blue	+
15) 6 + 9	blue	pink	blue	+

1. Scott, L.J., Jr., Microgram, Vol. 6, No. 11, Nov. 1973.
2. Clarke, E.B.C., "Isolation and identification of Drugs," The Pharmaceutical Press, London, U.K. (1969).

DEA LABORATORY NOTES

(4)

DATE November 2, 1973

NO. 68

-179-

DRUG TYPE COCAINE

METHODOLOGY Color Test

From

Microgram

Vol 6, No 11, (Nov: 1973)

SPECIFIC FIELD TEST FOR COCAINE

L. J. Scott, Jr.
Forensic Chemist
Drug Enforcement Administration
South Central Regional Laboratory
Dallas, Texas

OBJECTIVE

To develop a rapid, sensitive, highly specific color test for cocaine.

BACKGROUND

The use of cobalt thiocyanate reagent for the presumptive detection of cocaine is of limited reliability since many other compounds give the same reaction as cocaine. Lidocaine, procaine, PCP, tetracaine, benzocaine as well as cocaine give a blue precipitate with this reagent. The following test for cocaine has been empirically developed and successfully tested several hundred times by both chemists and agents. The "Ruybal" reagent¹ and its modification² are very useful in the hands of chemists, but field personnel often have experienced difficulty in their interpretation. The method proposed herein is almost impossible to misinterpret, and is highly sensitive and specific.

REAGENTS

Solution #1: 2% cobaltous thiocyanate dissolved in water and then diluted 1:1 with 96% USP Glycerine(*)

Solution #2: Concentrated Hydrochloric Acid

Solution #3: Chloroform

Other: Small disposable test tubes

(*) Available from K&K Laboratories, Plainview, N.Y.

PROCEDURE

Step 1: place a small amount of suspected cocaine in a test tube, add 5 drops Solution #1 and shake. If cocaine is present a blue color develops at once. If a blue color is not seen, add more sample. If a blue color still does not develop, the sample does not contain cocaine.

Step 2: Add 1 drop of Solution 2 and shake. The blue will disappear and a clear pink solution is seen. If all the blue does not disappear, add a second drop (no more) of HCl and shake.

Step 3: Add several drops of Solution 3 (chloroform) and shake. The CHCl₃ layer will develop an intense blue color if cocaine is present.

DISCUSSION

Only cocaine, from the following list of compounds, will give the results outlined above. The following chart illustrates how the test can select cocaine from other compounds listed.

	<u>Blue Color with Solution #1</u>	<u>Pink After Solution #2 Added</u>	<u>Blue CHCl₃ Layer From Solution #3</u>
Cocaine	yes	yes	yes
Phencyclidine	yes	yes	no
Dibucaine	yes	yes	no
Butacaine	yes	yes	no
Methapyrilene	yes	yes	no

The following compounds do not produce a blue color with Solution #1

- | | |
|-------------|---------------------|
| Antipyrine | Quinine |
| Procaine | Methadone |
| Benzocaine | Tetracaine |
| Mepivacaine | Lidocaine |
| Prilocaine | Heroin |
| Beta Eucain | Demerol |
| Aminopyrine | Sugars and starches |

Instructions to agents who use this procedure specifies that they are to use 5 drops solution 1, 1 or 2 drops (not more) of solution 2 and 5 drops solution 3. The amount of solution 1 is not critical, nor is the amount of solution 3. However, the ratio of solution 1 to solution 2 is critical. If excess concentrated Hydrochloric Acid is added to solution 1 after the blue color has developed with cocaine, a blue rather than pink solution will result; this blue will not extract into the CHCl_3 layer. If excess cocaine is used with solution 1, then it is sometimes necessary to add 2 drops of concentrated Hydrochloric Acid, no more should be used.

Consequently, if one adheres to the 5 drops + 1 drop + 5 drops of the three solutions 1, 2 and 3 respectively, the test will function easily and well. This test has not failed to detect cocaine as low as 1% cocaine in some cases. This system has been used in this laboratory on approximately 150 cocaine samples over a period of about 5 months, and about 50 to 75 times by agents on an experimental basis over a period of about 4 months.

No extensive shelf life tests have been run, but the solution #1 has been stored at room temperature under laboratory conditions for a period of six months with no detectable deterioration or loss of effectiveness.

REFERENCES

- (1) Richard Ruybal, Microgram Vol. 5, No. 3, 1972
- (2) Carolyn N. Ruybal, Microgram Vol. 6, No. 2, 1973

Property of Idaho State Police Forensic Services
Unclassified Inmate Copy Services
OBSOLETE DOCUMENT

LABORATORY NOTES



DATE JUNE 4, 1974

-68-

NO.

DRUG TYPE

METHODOLOGY

A "NEW" FIELD TEST REAGENT
By Ferris H. Van Sickle
NORTH CENTRAL REGIONAL LABORATORY
DRUG ENFORCEMENT ADMINISTRATION

INTRODUCTION

Agent Park Kaestner of the Kansas City, Missouri DEA Task Force recently described the use of laundry bleach as a field test reagent for procaine. It has been used in the Kansas City area in an undercover capacity.

Laundry bleach, commonly sold in grocery stores as Purex or Chlorox, is an aqueous solution of sodium hypochlorite, usually 5.25% by volume. As such, in an alkaline state, it is a strong oxidizing agent. A literature search of available material revealed nothing on sodium hypochlorite as a color test for organic compounds of interest in our field of endeavor.

PROCEDURE

Spot tests were performed in the usual manner using spot plates and comparing the colors developed with a standard. Purex solution purchased from a grocery store was used.

RESULTS AND DISCUSSION

Sodium hypochlorite, nicknamed "doper's reagent", when tried on a number of controlled and non-controlled drugs including some of the "caines" behaved much like the Sanchez Test for certain primary amines. "Doper's reagent" gave a chocolate brown with procaine and benzocaine. It was flaky in appearance. The following drugs were tested with "doper's reagent" and no color developed:

Heroin HCl, Cocaine (HCl and base), MDA HCl, PCP HCl, Methadone HCl, Demerol, common barbiturates, Codeine (phosphate and base), Amphetamine sulfate, Methamphetamine HCl, Mescaline sulfate, LSD, Chlordiazepoxide (Librium), Diazepam (Valium), Phenacetin, Acetanilid, Carbromal, Propoxyphene (Darvon), Phenyl propanolamine, Ephedrine sulfate, Mepivacaine, Dibucaine, Hexylcaine, Quinine, Caffeine, Mannitol, Saccharine, Lactose, Dextrose, Starch, Aspirin, Boric Acid, Stearic Acid, and Talcum.

Listed below are color reactions obtained with other drugs not previously mentioned:

Methapyrilene HCl-Light green-flaky
Lidocaine HCl-orange brown, oily, floating drops
Lidocaine sulfate-orange brown, oily, floating drops

Stovaine HCl-orange, oily drops-on standing
Salicylamide-greenish brown on standing
Salicylic Acid-brown on standing
Oxytetracycline-greenish yellow (chartreuse)
Sulfamerazine-greenish orange, then greenish yellow
Sulfadiazine-greenish orange, then greenish yellow
Chloroprocaine-orange
DMT-instant orange, then brown
DET-instant orange, then brown
Para amino benzoic acid-brown-flaky
Meta amino benzoic acid-brown-flaky
Morphine sulfate-brown-flaky

In conclusion, it would seem that sodium hypochlorite would have some usefulness as an easily obtained field test reagent for suspected procaine in cocaine-procaine mixtures; cobalt thiocyanate being used to check for cocaine. DMT, DET, benzocaine, and morphine sulfate should also be kept in mind. If the Marquis is negative for opium alkaloids, then morphine can be eliminated. The test does not work well with brown heroin-procaine mixtures. Since it gives the chocolate brown color with the hydrolysis product (or metabolite) of procaine, it could possibly be used in TLC drug screening procedures as a color development spray for p-amino benzoic acid.

The use of standards for color comparison with the unknown drug is encouraged whenever possible.

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

Microgram Vol. VIII, No. 1
Jan. 1975

A SIMPLE FIELD TEST FOR COCAINE NOT RELYING ON COBALT THIOCYANATE

Fred W. Grant, William C. Martin and Ralph W. Quackenbush
Marcy Psychiatric Center
Research Division
Marcy, N. Y. 13403

OBJECTIVE


Our purpose was to improve the sensitivity and specificity of present methods of cocaine detection in a simple, easily interpreted field test.

BACKGROUND

Efforts to improve the specificity of cobalt thiocyanate in cocaine field testing (1,2,3,4,5) have tended to compromise field adaptability. Recently Lorch (4) had to add a fourth reagent to the already involved three reagent test of Scott (3,5) to further improve specificity. It is our impression that cobalt thiocyanate has been pushed to the limit in cocaine detection and that perhaps a fresh start should be made in another direction.

METHOD

Our approach relates to the fact that cocaine is quite unique among drugs in being a benzoate ester. We have screened several hundred drugs, including virtually all of the common drugs of abuse, by molecular formula and find only the relatively obscure "piperocaine" sharing this property with cocaine. While we knew of no available color test specific for benzoate esters we were aware of the distinctive odors of the lower alkyl benzoate esters. It is an easy matter to cleave off the benzoate portion of cocaine in the form of methyl benzoate by treatment with sodium methoxide in methanol. Metallic sodium need not be used in the formation of sodium methoxide since sodium hydroxide in methanol will serve this purpose. The reaction is indicated below.



PROCEDURE

A few drops of a 5% solution of sodium hydroxide or potassium hydroxide in methyl alcohol are used to moisten the suspected cocaine specimen. A few minutes are allowed for the alcohol to evaporate and any odor is noted.

RESULTS AND DISCUSSION

The wintergreen-like odor of methyl benzoate is strong and distinctive. In presenting a test based on odor detection one faces a longstanding prejudice that rightfully maintains that subtle qualitative or quantitative distinctions based on odor cannot be made reproducibly by the human nose. However, in this connection it is of interest that a New York State Appeals Court recently ruled that the smell of marijuana smoke provided legal cause for Police to search an automobile and its occupants without a search warrant.

The saving grace of this test resides in the absence of any odor whatsoever when the test is applied to the vast majority of drugs likely to be encountered in the field. Occasionally, a very faint fishy odor arises when a low molecular weight amine, such as amphetamine, is released from its salt. Because of its specificity, the test is unaffected by the presence of excipients or other drugs. Methyl acetate is a product of the test when applied to heroin, aspirin, and other acetate esters but it is removed during the evaporation step. Benzoic acid itself fails to give a positive test. WATER WILL INTERFERE WITH THE TEST SO REAGENT AND SPECIMEN SHOULD BE KEPT DRY.

The test has been field tested by local County and State Police who now consider it the test of choice in the field identification of cocaine.

REFERENCES

- 1) R. Ruybal, Microgram, Vol. 5, No. 3 (1972)
- 2) C. N. Ruybal, Ibid., Vol. 6, No. 2 (1973)
- 3) L. J. Scott, Jr., Ibid., Vol. 6, No. 11 (1973)
- 4) S. K. Lorch, Ibid., Vol. 7, No. 8 (1974)
- 5) Drug Enforcement, Vol. 1, No. 3, pp. 26-27 (1974)

DRUG LABORATORY NOTES



DATE

-49-

NO. 37

DRUG TYPE Mixture

METHODOLOGY

Analysis of Cocaine, Procaine, Benzocaine Mixtures
Jeffrey M. Weber, Forensic Chemist
New York Regional Laboratory
Bureau of Narcotics & Dangerous Drugs

Introduction

This laboratory frequently receives illicit cocaine samples which have been adulterated with synthetic local anesthetics and sugar. Two local anesthetics commonly encountered are benzocaine and procaine.

Since cocaine, benzocaine and procaine have similar UV absorption characteristics, a means of separating all three components has been investigated. Milos and Porto² devised a procedure for the simultaneous spectrophotometric determination of cocaine and procaine mixtures. Canaff³ separates cocaine from procaine or tetracaine by using a 2.0N HCl chromatographic column. Moore⁴ separates the same mixture utilizing a 0.1N HNO₃ column.

The method below is a modification of Moore's procedure which can be utilized to quantify and identify cocaine in the presence of benzocaine and/or procaine.

Method

Apparatus

Chromatographic Column
UV Spectrophotometer

Reagents

Nitric Acid - 0.1N
Sulfuric Acid - 0.1N
Chloroform - Reagent Grade
Ethyl Ether - Reagent Grade
Ammonia - conc.
Celite 545 - acid washed

Procedure

Pack a pledget of fine glass wool into the base of a chromatographic column (25 cm long, 2 cm diameter) with a tamping rod. Place 2 grams of diatomaceous earth (Celite 545, acid washed)* into a 100 ml beaker. Add 1 ml 0.1 Nitric Acid, and mix with a spatula until homogeneous. Transfer the mixture to the column, and tamp moderately to compress the material into a uniform mass.

* Johns - Manville, New York, New York

Into another 100 ml beaker weigh accurately a portion of sample equivalent to about 15-20 mgs cocaine. (Illicit cocaine samples are often damp and lumpy, therefore, sample mixture should be carefully ground before proceeding). Add 2 mls 0.1N Nitric Acid to the sample and swirl to wet powder thoroughly. To this solution add 3 grams of diatomaceous earth and mix until fluffy. Quantitatively transfer this mixture to the column and tamp moderately. Place a pledget of glass wool on top of the packing.

The benzocaine is eluted from the column into a 200 ml volumetric flask with 200 ml water washed ether. Dilute the flask to volume with ether. Evaporate a 20 ml aliquot to dryness on a steam bath. Dissolve the residue in a known volume of 0.1N H₂SO₄ to give a final concentration of approximately 10 mcg/ml. Read the absorbance of the solution at 226 mμ and compare against a standard solution of benzocaine.

The cocaine may then be eluted by passing 100 mls water washed CHCl₃ into a 100 ml volumetric flask. Dilute flask to volume with CHCl₃. Evaporate a 10 ml aliquot to dryness on a steam bath. Dissolve residue in 0.1N H₂SO₄ to give a final concentration of approximately 15 mcg/ml. Read the absorbance of the solution at 233 mμ and compare against standard cocaine.

Procaine is eluted from the column into a 250 ml beaker by passing 100 mls of ammoniacal chloroform (prepared by vigorously mixing 2 ml concentrated ammonium hydroxide with 100 ml CHCl₃ and allowing the phases to separate) followed by 100 mls water saturated chloroform. Remove a 20 ml aliquot and evaporate to dryness on a steam bath. Dissolve the residue in 0.1N H₂SO₄ to give a final concentration of approximately 15 mcg/ml. Read the absorbance of the solution at 228 mμ and compare against standard procaine.

Identification

Methods most frequently used for the qualitative determination of mixtures of this type are infrared spectroscopy and thin layer or gas chromatography.

A - Infrared Spectroscopy

To obtain spectra for benzocaine and procaine, the remainder of the stock solutions eluted from the column are evaporated to dryness. The residue may then be deposited as a film on a sodium chloride window.

The cocaine that is eluted from the column must be re-extracted in order to obtain a satisfactory IR curve.

Evaporate the CHCl₃ stock solution to about 20 ml. Transfer to a 60 ml separatory funnel and shake with 20 mls 4% ammonium hydroxide. Evaporate the CHCl₃ extract on a steam bath (prior to doing so add 1 drop conc. HCl to solution). A KBr disk may now be prepared in the usual manner.

B - Thin Layer Chromatography

Mix a portion of powdered sample with methanol and apply to a silica gel G plate (Analtech). Use the solvent systems mentioned below and observe spots by spraying with iodoplatinate spray (prepared by mixing 0.25 gms Platinic Chloride and 5 grams Potassium Iodide in 100 ml water).

Solvent System A⁵ - Cyclohexane: Chloroform: Diethylamine 5: 4: 1
Solvent System B⁶ - Ethyl Acetate: Benzene: Ammonium Hydride (6:3.5:1)

	(System A)	(System B)
	<u>Rf</u>	<u>Rf</u>
Cocaine	0.65	0.68
Procaine	0.34	0.60
Benzocaine	0.26	0.41

C - Gas Chromatography

Sample is dissolved directly in methanol and injected under the following conditions:

Instrument	Packard #7360	
Column	3% OV1	3% OV25
Temperature	210°C	140°C
Carrier Gas	Nitrogen	Nitrogen
Carrier Flow	40ml/min	25ml/min
Detector	Flame Ionization	Flame Ionization
Air Flow	500ml/min	500ml/min
Hydrogen Flow	30ml/min	25ml/min

The following results were obtained:

Retention Time Cocaine	5.9 min	4.7 min
Retention Time Procaine	4.3 min	2.4 min
Retention Time Benzocaine	0.7 min	0.5 min

Discussion

Cocaine samples should be subject to qualitative tests before quantitative analysis is performed. The presence of procaine or benzocaine may be detected using the following spot tests:

1 - Sanchez reagent (prepared by mixing 2.5 ml freshly distilled furfural, 22.5 ml 95% EtOH and 75 ml glacial acetic acid). A red color is obtained in the presence of benzocaine or procaine.

2 - p-dimethylaminobenzaldehyde (1 gm in 50 ml EtOH and 50 ml conc. hydrochloric acid). An intense yellow color is obtained in the presence of benzocaine or procaine.

Furthermore, chromatographic screening tests should be performed prior to quantitation to detect the presence of benzocaine and/or procaine. In cases where no benzocaine is present, the ether wash will not be necessary during the quantitation.


Quantitative results of better than 95% cocaine have been obtained using this procedure. Results obtained for the benzocaine and procaine ranged between 92 and 94%. The results obtained on simulated mixtures are tabulated as follows:

	mgs added			mgs Found		
	Benz.	Pro.	Coc.	Benz.	Pro.	Coc.
1.	9.8	10.1	10.4	9.3	9.2	10.1
2.	10.0	10.3	20.6	9.2	9.8	19.6
3.	15.1	14.8	15.4	14.0	14.1	14.9


References

1. E.G.C. Clarke, Isolation and Identification of Drugs Pharmaceutical Press, London, 1969
2. C. Milos, P. Porto, Spectrophotometric Determination of Cocaine in the Presence of Procaine and Benzocaine, IRS Publication, Alcohol and Tobacco Tax Laboratory, New York, N.Y.
3. R.F. Canaff, Ion Pairing Chromatographic Separation and Determination of Cocaine, Microgram Vol III, #4, 121 June 1970.
4. J. Moore, The Analysis of Cocaine in Illicit Preparations Microgram Vol III, #3, 89-95, May 1970.
5. Sunshine, Handbook of Analytical Toxicology, Chemical Rubber Company 1969, Reference #1A.
6. IRS Publication #341, June 1967, page 94, system S4

4% diethylamine
in Benzine



Tr (7 drops NH_4OH in
10ml CH_3OH)



Sept 73

OTHER MICROCHEMICAL TESTS

and dl-Amphetamine - 1% HCl (volatility) gives "peacock feather" like crystal forms with low birefringence.

Ethylamphetamine - AuCl in dilute H₃PO₄ (volatility) gives single forms, fans, clusters and sheafs of rods with considerable rectangular-like internal structure. Birefringence colors.

ndimetrazine - AuCl in dilute H₃PO₄ (volatility) gives curved needles and rods in extensive branching systems with moderate birefringence.

phenhydramine - AuCl in dilute H₃PO₄ (volatility) gives very small white rods with slanted ends (= narrow parallelograms).

ocaine - AuCl in H₂O (direct) gives bright birefringent needles and rods of regular shape alone, in clusters and in Xs.

aine - AuCl in H₂O (direct) gives rough, irregular forms with moderate to low birefringence.

haqualone - dissolve sample in drop of MeOH + add drop of 5% Sodium Carbonate which gives diamonds, hexagons, parallelograms and trapezoid forms with bright birefringence.

in - dissolve in 5% HCl + add drop of 10% Sodium Acetate which gives grey & colored birefringent hexagons.

in dil.H₃PO₄ = 1g. HAuCl₄·3H₂O in 20ml. H₃PO₄(1+2)

in H₂O = 1 g. HAuCl₄·3H₂O in 20 ml. H₂O

dry PtBr = 1g. platinum chloride in 10ml of 40% HBr; add 10ml. conc. H₂SO₄ + 10ml. H₂O

dry PtBr = 5g. H₂PtCl₆·6H₂O and 10g. NaBr in 100ml. H₂O. allow to age 48 hours

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

Cocaine in the presence of Lidocaine
~~But~~ Ruybal's "Gold Bromide"

from: Modern Microcrystal Tests for Drugs, Fulton.
page 385 #4

Reagent: $H[AuBr_4]$ in $2H_3PO_4 - 1(2+3)H_2SO_4$
In plain English, you make it
up this way

- $H[AuCl_4 \cdot 3H_2O]$ crystals 1.0 gram
- HBr (40%) 1.5 ml.

Dilute above with

~~2 ml of conc. H_2SO_4 to 3 ml H_2O~~
1 ml. of conc. H_2SO_4 dil. with 18 ml. of H_2O

50 ml.

Property of Andhra State Police Forensic Services
Uncontrolled Internet Copy
SOLETE DOCUMENT

Reprinted from

Journal of Chromatography, 101 (1974) 215-218

© Elsevier Scientific Publishing Company, Amsterdam — Printed in The Netherlands

CHROM. 7751

Noté

Gas chromatographic detection of ecgonine and benzoylecgonine in cocaine

J. M. MOORE

U.S. Department of Justice, Drug Enforcement Administration, Special Testing and Research Laboratory, 7704 Old Springhouse Rd., McLean, Va. 22101 (U.S.A.)

(Received June 24th, 1974)

Cocaine hydrochloride is used in medicine as a topical anesthetic. It also has widespread use as a stimulant in the illicit drug market. During recent analytical investigations into illicitly manufactured cocaine, it was discovered that gas chromatographic (GC) methodology to detect the primary cocaine hydrolysis products was lacking. Cocaine can be hydrolyzed in dilute acid to benzoylecgonine, ecgonine, benzoic acid and methanol¹.

The detection of cocaine hydrolysis products is important in the quality control of pharmaceutical cocaine. A number of papers have described the detection of benzoylecgonine and ecgonine in pharmaceutical cocaine and also in botanical and biological samples²⁻⁵. These procedures utilized primarily paper and thin-layer chromatographic techniques. Fish and Wilson⁶ detected benzoylecgonine in urine by methylation with diazomethane prior to GC analysis. However, in this procedure ecgonine was not determined.

GC was selected as an identification tool for ecgonine and benzoylecgonine because it was believed it would surpass present methodology in sensitivity, speed, and specificity. Due to the highly polar and amphoteric nature of the two degradation products, it was determined that their separation from cocaine by liquid-liquid extraction techniques prior to GC analysis was impractical. Subjection of the sample to direct GC analysis was also difficult due to the very poor chromatographic behavior of ecgonine and benzoylecgonine. Therefore, this study describes a procedure whereby both hydrolysis products are silylated prior to GC analysis using N,O-bis(trimethylsilyl)acetamide (BSA). Both ecgonine silyl and benzoylecgonine silyl possess very good chromatographic properties. The derivatized samples are chromatographed on 10% OV-101 on Chromosorb W-HP and 3% OV-25 on Gas-Chrom Q, using temperature programming. Using this procedure ecgonine and benzoylecgonine can be detected in uncut cocaine samples at levels less than 0.1% and 0.3%, respectively.

EXPERIMENTAL

Reagents and chromatographic materials

The BSA silylating reagent used in this study was supplied by Pierce (Rockford, Ill., U.S.A.). The 10% OV-101 on Chromosorb W-HP (100-120 mesh) and 3% OV-25 on Gas-Chrom Q (100-120 mesh) stationary phases were obtained from Applied Science Labs. (State College, Pa., U.S.A.).

Standards

The cocaine hydrochloride used in this study was supplied by S. B. Penick (New York, N.Y., U.S.A.). Hexadecane, eicosane, tetracosane and triacontane internal standards were obtained from Applied Science Labs. Benzoylecgonine and ecgonine were supplied by the Drug Enforcement Administration.

Apparatus

A Packard Model 7400 gas chromatograph was used for all chromatography.

Sample analysis

A 25-mg sample of cocaine hydrochloride is placed in a 1-ml glass-stoppered test tube and 500 μ l of BSA is added. The tube is stoppered loosely and heated at 75° for 10 min with occasional agitation. After derivatization is complete, 3-4 μ l of the solution are injected into the gas chromatograph using the following parameters:

Column: Coiled glass, 4 ft. \times 4 mm I.D. (OV-101)
 Coiled glass, 6 ft. \times 4 mm I.D. (OV-25)
 Stationary phase: (a) 10% OV-101 on Chromosorb W-HP, 100-120 mesh
 (b) 3% OV-25 on Chromosorb Q, 100-120 mesh
 Carrier gas: Nitrogen, 60 ml/min
 Detector: Flame ionization
 Air: 500 ml/min
 Hydrogen: 50 ml/min
 Injection temperature: 275°
 Detector temperature: 275°
 Column temperature: Temperature programmed
 initial temperature: 180° (OV-101), 170° (OV-25)
 initial hold: 10 min (OV-101), 5 min (OV-25)
 program rate: 3°/min (OV-101 and OV-25)
 final temperature: 260°
 final hold: 5 min
 Sensitivity: 3×10^{-9} a.f.s. (OV-101 and OV-25)
 Chart speed: 0.2 in./min (OV-101 and OV-25)

RESULTS AND DISCUSSION

Fig. 1 illustrates a gas chromatogram on OV-101 of cocaine hydrochloride containing 3% benzoylecgonine and 1% ecgonine following BSA treatment. Since temperature programming was used, two internal standards were desirable. It is also apparent from Fig. 1 that the detection level for ecgonine is lower than that for benzoylecgonine. However, benzoylecgonine can still be detected at levels considerably less than 0.3%.

In order to investigate the effect the moisture content of cocaine samples would have on the silylation process, water was added to the sample shown in Fig. 1 and to a commercial cocaine hydrochloride standard at a 20% w/w level. The presence




Fig. 1. Gas chromatogram of cocaine hydrochloride containing 1% ecgonine and 3% benzoyllecgonine on OV-101 following BSA treatment. See GLC parameters under *Sample analysis*.

of moisture causes no detectable cocaine hydrolysis during the derivatization process. Additionally, the silylation of benzoyllecgonine and ecgonine was not affected noticeably. The silylation of both ecgonine and benzoyllecgonine is rapid and the BSA solution is stable for at least several hours.

In addition to using BSA derivatizing reagent and OV-101 stationary phase, other silylating compounds and chromatographic stationary phases of varying polarity were investigated. The other silylating reagents included N-trifluoroacetyl-imidazole and N,O-bis(trimethylsilyl)trifluoroacetamide (BSTFA), the fluorinated analog of BSA. Of these only BSTFA proved of any value. The only advantage in using BSTFA, instead of BSA, was that a slight increase in resolution of the benzoyllecgonine and cocaine chromatographic peaks was noted. This was probably due to the relative insolubility of cocaine in BSTFA. However, because of this insolubility, and subsequent loss in benzoyllecgonine sensitivity using BSTFA, BSA was selected as the reagent of choice.

In addition to OV-101, 3% OV-1 stationary phase was investigated. However, the resolution between cocaine and benzoyllecgonine was less than desirable on this column. More polar phases were also investigated. These included OV-17, OV-25, OV-210, and OV-225. On OV-17, cocaine and benzoyllecgonine silyl had the same retention times. When OV-210 and OV-225 were used in conjunction with temperature programming, a rather unstable baseline resulted. OV-25 proved to be the most suitable of the polar stationary phases. When using this phase the elution order of cocaine and benzoyllecgonine silyl were reversed when compared to OV-101 (see Table I). The OV-25 was used in this study only as a confirmation for the presence of ecgonine and benzoyllecgonine in cocaine. OV-101 was the column of first choice because the resolution, sensitivity and retention times of cocaine and its hydrolysis products were the most favorable. Table I lists retention times of cocaine, benzoyllecgonine silyl, ecgonine silyl and internal standards on OV-101 and OV-25. All internal standards were chromatographed using separate chloroform solutions.

The procedure given in this paper is a rapid and sensitive method for the de-

TABLE I

RETENTION TIMES OF COCAINE, ECGONINE Silyl, BENZOYLECGONINE Silyl AND INTERNAL STANDARDS ON OV-101 AND OV-25 STATIONARY PHASES
GLC operating parameters for both columns are given in the text under *Sample analysis*.

Compound	Retention time (min)	
	OV-101	OV-25
Hexadecane internal standard	3.8	—
Ecgonine silyl	6.7	2.8
Eicosane internal standard	—	5.2
Cocaine	19.7	25.8
Benzoyllecgonine silyl	22.2	24.5
Tetracosane internal standard	26.7	—
Triacotane internal standard	—	30.2

tection of small amounts of ecgonine and benzoyllecgonine in cocaine. It offers the added advantage of being readily adaptable to quantitative work. The procedure is useful not only in pharmaceutical quality control but in the detection of contaminants in illicitly manufactured cocaine.

REFERENCES

- 1 A. Pictet and H. C. Middle, *Vegetable Alkaloids*, Wiley, New York, 1904.
- 2 C. A. Sanchez, *An. Fac. Farm. Bioquim., Univ. Nac. Mayor San Marcos Lima*, 8 (1957) 82.
- 3 J. Buchi and Schumacher, *Pharm. Acta Helv.*, 32 (1955) 194.
- 4 Maria Theresa Campos Neves and Andre da Silva Campos Neves, *Esc. Farm.*, 27 (1967) 67.
- 5 R. V. Ofiz, *An. Fac. Quim. Farm. Univ. Chile*, 18 (1964) 45.
- 6 P. Majlat and J. Bayer, *J. Chromatogr.*, 20 (1965) 187.
- 7 K. C. Guven and T. Altinkurt, *Eczacilik Bul.*, 7 (1965) 11, 119.
- 8 F. Fish and W. D. C. Wilson, *J. Pharm. Pharmacol.*, 21 (1969) 1355.

IDENTIFICATION OF COCAINE BY INFRA-RED SPECTROSCOPY

L. E. WENER, CHICAGO REGIONAL LABORATORY, BNDD

The spectrogram obtained from the infra-red spectroscopic examination of compounds isolated in a pure or nearly pure state is one of the most substantive methods of examination. It is a very important aid to the forensic chemist now that the courts and defense attorneys are beginning to challenge the identification of alkaloids, such as cocaine, by means of crystal and color tests alone.

Cocaine is ordinarily separated from most of its adulterants and/or diluents by extracting it as the hydrochloride salt with chloroform from a solution or suspension in 1+9 HCl. The filtered chloroform extract is then evaporated on the steam bath to obtain the cocaine salt relatively pure. No difficulty was ordinarily encountered in obtaining the characteristic identifying crystals with the platinum chloride or gold chloride reagent solution. However, the cocaine hydrochloride crystals isolated this way were found not to give infra-red spectrograms which, for forensic purposes, sufficiently matched that given by the authentic.

It was theorized that the above difficulty was due to the acid hydrolysis, however much or slight, of the cocaine at steam bath temperatures due to the hydrochloric acid dissolved in the water solubilized in the chloroform. When the chloroform extract was first dried with excess anhydrous sodium sulfate, decanted, and then evaporated no such deterioration of the cocaine occurred and excellent infra-red spectrograms were obtained.



DATE January 6, 1972

-80-

NO.

DRUG TYPE Narcotic

METHODOLOGY Infrared Identification

RAPID SEPARATION OF COCAINE FROM ADULTERANTS
SUCH AS PROCAINE AND QUININE, AND SUBSEQUENT
INFRARED IDENTIFICATIONRoger G. Fuelster
Forensic Chemist
CHICAGO REGIONAL LABORATORY, BNDD

Previous methods of separation of cocaine from adulterants involved lengthy column chromatography. The following procedure is not only rapid but yields a highly pure product.

PROCEDURE

Place sample in a small separatory funnel. Add approximately 15ml 2.8 normal hydrochloric acid. Extract with an equal volume of chloroform. Pass the extract through Whatman Phase Separating Paper #1PS. Evaporate the chloroform to dryness and dry at 105°C for 10 minutes. Obtain the IR spectra in the usual manner.

SUMMARY

This method effectively removes traces of hydrochloric acid which would otherwise hydrolyze the cocaine during the final evaporating and drying stages.

from Bulletin Vol 3 #2, July 1975

ABSTRACT OF A PAPER PRESENTED AT THE SPRING 1975 MEETING OF THE MIDWESTERN ASSOCIATION OF FORENSIC SCIENTISTS.

"A Simple Procedure for the Separation and Identification of Cocaine".
by Jon D. Naylor, Carl R. Phillips, Robert J. McCurdy & Stephen A. Koers.

The selective extraction of the blue complex of cocaine with cobalt thiocyanate from aqueous acid into chloroform in the presence of common street precipitants, e.g. lidocaine, tetracaine, procaine, etc., has suggested a very simple purification and identification procedure. The infrared spectrum of the dried blue extract has been found to be quite characteristic for cocaine. The spectrum consists of a simple addition of the absorption peaks for cocaine and cobalt thiocyanate with very few shifts in frequency or distortions of component peak shapes. The procedure is: a.) addition of about 2 ml of cobalt thiocyanate agent (2g of cobalt thiocyanate in 100ml of water) to enough street sample to contain one to two mg. of cocaine, b.) addition of up to one-half ml. of concentrated hydrochloric acid (excess HCl results in displacement of cocaine by chloride in the complex to form a blue solution), c.) dropwise addition of enough water to dissolve all of the blue precipitate (vigorous shaking is necessary), d.) extraction of the blue cocaine-cobalt thiocyanate complex into chloroform, e.) ~~dry~~ ^{dry} $CHCl_3$ ^{then} ^{powdered} ^{Na_2SO_4} the complex in an evaporating dish, and f.) running the infrared spectrum of ^{the} ^{dried} ^{blue} ^{complex} ^{as} ^a ^{KBr} ^{mull}. ^{by drying the small volume of $CHCl_3$ right with the KBr.} The resultant spectrum has been found to be totally free of other materials when mixtures of cocaine with procaine, lidocaine, and tetracaine were treated in the above manner.

If pure cocaine instead of the complex is desired, the chloroform extract in the above may be washed with aqueous ammonia to displace the cocaine from the complex. The tan-colored ammonia complex will be extracted into the aqueous layer, the clear chloroform layer will contain cocaine free base which can be treated in the usual manner for further identification.

The greater speed, more complete separations, and ability to perform the entire operation in a test tube makes this procedure far superior to the traditional multiple extraction and chromatographic separation procedures.

Cocaine/Lidocaine Mixture

Rapid separation/Infrared Identification
(11.73386 minutes)

Dissolve sufficient sample in minimum amount 0.1N HCl to provide 5 mg cocaine.

Precipitate cocaine by dropwise addition of 5% Platinic Chloride (H_2PtCl_6) until ppt. ceases.

With disposable pipet transfer ppt and supernatant to a little, tiny buchner funnel and, w/vacuum, remove liquid.

Wash 3X w/ 3ml portions of 0.1NHCl, followed by 3X 3 ml portions of methanol.

Continue to dry the precipitate with suction, for

Three or four minutes - cover buchner with dry filter paper if air is wet.

When sufficiently dry, press with KBr.

Compare with spectrum prepared from known cocaine.

Also works with methamphetamine separation from phentermine but use methanol sparingly as the methamphetamine is slightly methanol soluble.

Joe Power - DEA


Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

Column Chromatograph


Celite 545 aqu. stationary bed

D. Ion Pairing
must use a chlorinated solvent
to elute drug from column

Drug mix
2ml 0.5N NaHCO₃
3g Celite
4ml 0.2N HNO₃
+ 6g Celite



2ml 0.1N
H₂SO₄ + 3g Celite



mix of
 cocaine
 procaine
 Caffeine
 Quinine
 Dextrose
 lactose
 corn starch

ISOLATION AND ESTIMATION OF COCAINE IN THE PRESENCE OF
TETRACAINE, PROCAINE, BENZOCAINE AND AMPHETAMINE

By

Charles Milos and Philip V. Porto
Alcohol and Tobacco Tax Laboratory
Internal Revenue Service
New York, New York

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

Isolation and Estimation of Cocaine in the Presence of
Tetracaine, Procaine, Benzocaine and Amphetamine

By

Charles Milos and Philip V. Porto

Illicit samples of cocaine are generally adulterated with sugar and procaine. These present no problem to the analyst (1). Some samples with low cocaine content may contain all of the components listed in the title of this paper. With this type of sample the problem is twofold -- detecting the cocaine and isolating it from the other components. Sometimes the second objective must be accomplished before the first can be established.

Since existing methods failed to resolve this problem, a new procedure was needed. An extraction method was devised which makes use in part of some of the techniques in existing methods (1).

Experimental

Chloroform completely removes benzocaine from aqueous acid solutions. Separation of tetracaine from procaine, amphetamine and cocaine is effected by first dissolving their free bases in a mixed solvent consisting of equal volumes of ethyl ether, ethyl acetate plus 20% by volume of isobutyl alcohol. Extraction of the mixed solvent with a citric acid-phosphate buffer of pH 2.25 removes the latter three components, leaving the tetracaine

in the solvent. By altering the pH of the buffer to 4.5 the cocaine can be removed by extraction with chloroform. Traces of tetracaine remaining in cocaine thus separated cause no interference with the usual cocaine tests.

Preliminary Examination of Sample

1. Dissolve a small portion of sample in 1 or 2 drops of 1:1 hydrochloric acid on a microscopic slide. Add a drop of platinum chloride reagent (1g. platinum chloride in 20 ml. of water) and examine the slide at once under the microscope, using the low power objective. If appreciable amounts of cocaine are present (5% or more), characteristic feathery crystals form. In a brief period of time they become modified beyond recognition. When relatively large amounts of benzocaine or amphetamine are present these crystals may not be observed. When this occurs place another portion of sample on a microscope slide, add several drops of a 20% aqueous solution of sodium acetate, and after several minutes observe for free base cocaine crystals. These crystals are characteristic and procaine, tetracaine, amphetamine and benzocaine do not interfere.
2. Add concentrated nitric acid to a portion of sample on a spot plate. Formation of a yellow to light brown color indicates tetracaine. Treat another portion with Marquis' reagent. A yellow to orange color indicates amphetamine.

3. Place a few milligrams of dry sample on a microscopic slide and while observing under low power, add a drop of 1:1 hydrochloric acid. Cocaine, procaine, tetracaine and amphetamine dissolve rapidly. Benzocaine dissolves very slowly.
4. To another portion of sample on a spot plate, add several drops of Sanchez reagent (2). A red to purple color indicates benzocaine or procaine. NOTE: Each 100 ml. of Sanchez reagent contains 2.5 ml. of freshly distilled furfural, 22.5 ml. of 95% ethyl alcohol and 75 ml. of glacial acetic acid.
5. To several milligrams of dry sample on a spot plate add a few drops of cobalt thiocyanate, stannous chloride reagent. The formation of a blue precipitate indicates cocaine. NOTE: This reagent contains equal volumes of a 2% aqueous solution of cobalt thiocyanate and a 5% stannous chloride in 10% hydrochloric acid (3).

Method

Reagents:

0.5 Molar citric acid

0.5 Molar disodium phosphate

Buffer solution -- To 40 ml. of 0.5 molar disodium phosphate add 15 ml. of 0.5 molar citric acid and mix.

Mixed solvent ---- To 40 ml. of ethyl ether add 40 ml. of ethyl acetate and 20 ml. of isobutyl alcohol and mix.

Separation Procedure

1. Dissolve a known weight of sample (approximately 0.1 g.) in 10 ml. of approximately 0.1 N hydrochloric acid and quantitatively transfer to a separatory funnel. Extract four times with chloroform using 15, 10, 10, and 5 ml. portions. If a quantitative determination of benzocaine is desired, combine and wash the chloroform with 5 ml. of water in a separator. Transfer the solvent to a tared dish and evaporate to dryness on a steam bath.
2. To the acidic aqueous solution add dilute ammonium hydroxide (10%) until just alkaline to litmus and extract three times with 20 ml. portions of solvent. Combine the solvent in a separator, wash with 5 ml. of water and after complete separation of solvent, discard the lower aqueous layer. Extract the solvent three times with 10 ml. portions of prepared buffer and combine the buffer in a separator. For a quantitative estimation of tetracaine wash the mixed solvent with 10 ml. of water and transfer the solvent to a tared dish. Evaporate to dryness on a steam bath. The residue is free base tetracaine.
3. To the combined buffer add 3 ml. of 0.5 molar citric acid, mix and extract with chloroform as already described for benzocaine. Wash the combined chloroform with 5 ml. of water and transfer the chloroform to a tared dish. Evaporate to dryness on a steam bath. The residue is free base cocaine.

Analytical Data

Table 1

S.	Milligrams Added					Milligrams Found		
	B	T	P	A	C	B	T	C
1.	23	19	50	10	4.7	22.6	19.7	4.4
2.	12	15	20	12	8.0	12	14.8	7.7
3.	13	22.5	30		15.0	13	22.8	14.5
4.		9.5	21	5	2.8		9.3	2.3
5.	10	20	30	12	9.0	9.7	19.6	8.4
6.					8.2			7.6

S--Sample

B--Benzocaine

T--Tetracaine HCl

P--Procaine HCl

A--Amphetamine HCl

C--Cocaine HCl

Table 1 shows results obtained by the proposed method with simulated mixtures. All of the results were obtained gravimetrically.

Discussion

Failure of the preliminary tests to disclose the presence of cocaine should not be taken as conclusive evidence of its absence. Further explorations should be conducted to either verify or negate the preliminary findings. Some samples of illicit cocaine are predominately mixtures of local anesthetics containing small amounts of cocaine. By the proposed method a large number of these diluents can be removed. However, if the cocaine residue is contaminated and the usual tests do not give a clear-cut decision with respect to its presence or absence, part of the residue should be reserved for preparing and identifying the derivative ecgonine. Use of the original sample may introduce complications which can be avoided by using part of the cocaine residue. The reaction of ecgonine with phosphomolybdic acid reagent (1g. phosphomolybdic acid in 20 ml. of 1:1 nitric acid) is so sensitive that the concentration of ecgonine should be kept in the neighborhood of about 15 parts per million. With high concentrations of ecgonine an innumerable number of crystals are formed. These adhere to each other and the characteristics of the individual crystal cannot be discerned. To prepare this derivative heat a dilute slightly acidic solution of cocaine to the boiling point. Add 10% sodium hydroxide solution until alkaline and continue the heating for a minute. Cool and transfer to a separator. Extract several times with 5 ml. portions of chloroform and discard the chloroform.

Acidify with acetic acid and extract once with about 15 ml. of ethyl acetate. Transfer the aqueous layer to a dish and heat to expel the dissolved solvent. Test a few drops and either dilute or concentrate, depending on the results of the test.

In the samples shown in the preceding table, the tetracaine residue was examined for cocaine and none was detected. Small amounts of procaine were present when the amount in the sample was large.

Some cocaine remains in the final buffer. By altering the pH the cocaine can be completely removed, but invariably it is contaminated with procaine.

The proposed method gives results of sufficient accuracy for routine analysis.

References

1. C. H. Riley -- Amer. Jour. Pharm., Vol 127, No. 7, July 1935
2. J. A. Sanchez - Anales Assoc. Quinn, Argentina 23, 431, 1925
3. Methods of Analysis for Alkaloids, Opiates, Marihuana and Misc. Drugs, U.S. Treasury Dept. IRS, Publication No. 341, (Rev. 2-60)

$$w = (f_1 w_1 + f_2 w_2 + \dots + f_n w_n) \quad (\text{Eq. 9})$$

f_1 is the weight fraction of particles of weight, etc. For spherical particles:

$$w = \frac{\pi \rho}{6} \sum f_i d_i^3 \quad (\text{Eq. 10})$$

In the case where particles are the same size, Eq. 10 becomes:

$$w = \frac{\pi \rho d^3}{6} \quad (\text{Eq. 11})$$

For Eqs. 10 and 11 to be equivalent:

$$d^3 = \sum f_i d_i^3 \quad (\text{Eq. 12})$$

Substituting the result from Eq. 8 into Eq. 12 gives $d = 96.2 \mu\text{m}$. Hence, the particle-size distribution given in the example, which was determined by trial and error, is equivalent to the calculated value of $d_{(R)}$ of $96 \mu\text{m}$. This distribution is typical of the kind of result encountered in practice. Use of such a particle-size distribution for all active ingredients would allow a safety margin for ingredients *P* and *Q* but, in the use of *R*, would necessitate achieving a truly random mixture to fulfill the desired tolerance range of 10%.

In practice, a random mix is not always achieved and it may be desirable to introduce an additional safety margin for the lowest concentration drug, *R*. *DeVey et al.* (1) did this in effect by setting the calculated effective mean particle-size limit for *R* as the maximum particle-size limit for the mixture. Alternatively, the coefficient of variation used in Eq. 7 could be set at a lower value than that corresponding to the specified tolerance range of $\pm 10\%$. For example, instead of 3.333%, a C_v value of 2.5% could be used which would give $d_{(R)}$ from Eq. 7 equal to 79 μm . An equivalent particle-size distribution corresponding to this value of $d_{(R)}$ would contain a considerable fraction above the proposed maximum limit of 50 μm (1) while still incorporating a safety margin to allow for the occurrence of nonrandomized mixing.

In conclusion, converting the particle-size limit into an equivalent particle-size distribution increases the utility of the calculations and provides a more convenient guideline in the practical situation. Additionally, a particle-size distribution of a drug obtained on recrystallization or precipitation or after milling can be tested for its suitability with regard to content uniformity by evaluating $\sum f_i d_i^3$ and comparing this value with the value of d^3 derived from Eq. 2 or 7.

- (1) J. A. Hersey, P. Cook, M. Smyth, E. A. Bishop, and E. A. Clark, *J. Pharm. Sci.*, 63, 408(1974).
- (2) M. C. R. Johnson, *Pharm. Acta Helv.*, 47, 546(1972).
- (3) K. Stange, *Chem.-Ing.-Tech. Z.*, 26, 331(1954).

M. C. R. Johnson
 Upjohn Limited
 Crawley, Sussex
 United Kingdom

Received May 15, 1974.
 Accepted for publication September 11, 1974.
 I wish to thank Dr. J. Hersey for his helpful comments.

Definitive GLC Method of Identifying Cocaine

Keyphrases □ Cocaine—definitive GLC identification □ Trimethylanilinium hydroxide—on-column methylation of cocaine, GLC identification □ GLC—identification, cocaine

To the Editor:

The identification of underivatized cocaine by GLC can be misinterpreted and erroneously reported as pentazocine, levorphanol, or methaqualone when using programmed or isothermal temperatures on 7% OV-17¹. TLC can also pose problems and lead to the report of a false positive for methadone instead of cocaine (1). Many laboratories are combining mass spectrometry with GLC to provide a more definitive instrumental method for identifying drugs such as cocaine (2); however, many laboratories cannot afford a mass spectrometer and, therefore, more definitive GLC methods of analysis are desirable.

In view of these problems encountered when employing GLC or TLC as a means of identifying cocaine, we wish to report a novel, definitive GLC method of identifying cocaine via an on-column GLC reaction under methylation reaction conditions that is applicable to confirming the presence of cocaine in various legitimate and illegitimate dosage forms. In our laboratory we have routinely used trimethylanilinium hydroxide in methanol as a methylating reagent for GLC analysis of anticonvulsant drugs in body fluids (3, 4). We anticipated that this methylating reagent would have an interesting on-column




Figure 1—Characteristic chromatogram representing an on-column reaction of cocaine (1.5 μg) and trimethylanilinium hydroxide (no time lapse after adding the methylating reagent to cocaine). See Table I for identification of the peaks.

¹ In our laboratory, these drugs have retention times similar to cocaine under programmed and isothermal conditions and are extracted concurrently with cocaine at basic pH.

Table I—Relative GLC Retention Times of Cocaine and Some Commonly Abused Drugs after On-Column Reaction^a

Drug	R _t , min
Cocaine ^{b,c}	Peak 1 (8.0), <i>N,N</i> -dimethylaniline Peak 2 (10.4), possible Hoffman elimination product Peak 3 (12.5), ecgonidine methyl ester Peak 4 (13.9), under investigation Peak 5 (14.2), ecgonine methyl ester (under investigation) Peak 6 (17.2), under investigation Peak 7 (23.3), cocaine
Amphetamine	8.3
Benzocaine	15.7, 16.8, 17.6, 18.2
Codeine	27.5
Heroin	24.1, 27.2
Levorphanol	21.4
Lidocaine	16.6, 18.5
Methaqualone	23.2
Morphine	27.5
Pentazocine	21.9
Phencyclidine	18.2
Phenobarbital	12.8, 14.3, 14.5, 18.4
Procaine	17.5, 20.8, 21.7, 23.1
Quinine	34.0, 37.0
Secobarbital	15.4

^a Under programmed temperature GLC conditions. ^b The number of chromatographic peaks seen after on-column reaction depends on the concentration of cocaine and the time lapse between addition of the methylating reagent and its on-column injection (see Figs. 1-3). ^c Similar results can be obtained on 3% OV-17 by reducing the nitrogen flow rate from 60 to 30 ml/min and changing the programmed column conditions from 50-250° (10°/min) to 50-250° (8°/min).

Reaction with an unusual bicyclic diester tertiary amine structure such as cocaine, through possible ester cleavage and methylation, as well as serve as a definitive confirmatory method by converting it into one or more identifiable derivatives.

All GLC injections were made on a dual-channel instrument² equipped with four hydrogen flame detectors. The 1.83-m (6-ft) U-shaped glass columns (2




Figure 2—Chromatogram representing an on-column reaction of cocaine (21 µg) and trimethylanilinium hydroxide (no time lapse after adding the methylating reagent to cocaine). See Table I for identification of the peaks.

mm i.d.) were packed with 7% OV-17 on 80-100-mesh Chromosorb W³. Operating temperatures were: injector port, 275°; column (isothermal), 250°; column, 50-250° (programmed at 10°/min); and detector, 275°. Flow rates (milliliters per minute) were: nitrogen, 60; air, 300; and hydrogen, 40. Instrumental attenuation was 8×10^{-10} . Under these conditions, reference standard underivatized cocaine had retention times of 23.3 and 2.9 min under programmed and isothermal column temperatures, respectively.

The reaction of cocaine with trimethylanilinium hydroxide solution was initially studied by adding 50 µl of a 2 M methanolic trimethylanilinium hydroxide solution⁴ (3) to 0.15 mg of cocaine and immediately injecting 0.5 µl (1.5 µg) into the chromatograph under programmed temperature conditions. The resulting chromatogram is illustrated in Fig. 1. Peak 1 is *N,N*-dimethylaniline, a product derived from trimethylanilinium hydroxide during the methylation reaction.

Utilizing a preparative gas chromatograph, we trapped peak 3 and obtained its mass spectrum. Peak 3 was identified as ecgonidine methyl ester (molecular ion, *m/e* 181; base peak, *m/e* 152). Its mass spectrum is identical to the corresponding methylated and trapped reference standard ecgonidine⁵. Ecgonidine methyl ester (peak 3) appeared to be the most characteristic product of on-column reaction between cocaine and trimethylanilinium hydroxide and was observed at all concentrations regardless of the length of time elapsing between addition of the methylating reagent and its on-column injection (peak 3 in Figs. 1-3).

Figure 2 is a representative chromatogram depicting the cocaine-trimethylanilinium hydroxide reaction using a larger amount of cocaine (21 µg), in which case the cocaine-trimethylanilinium hydroxide solution was injected immediately after adding the methylating reagent to cocaine. Mass spectral analysis of peak 2 indicates it to be a potential double Hoffman elimination product, although its structural determination is inconclusive at this time.

Figure 3 represents the chromatogram of the same cocaine-trimethylanilinium hydroxide solution after standing at room temperature for 24 hr. An increase in the concentration of peak 2 (possible Hoffman elimination product) occurred when cocaine was allowed to stand in the alkaline trimethylanilinium hydroxide solution for 1-24 hr (compare Figs. 2 and 3). Peaks 4 and 6 are presently under investigation, and peak 7 is cocaine.

Based on preliminary mass spectral analysis, peak 5 appears to be the ecgonine methyl ester derivative (molecular ion, *m/e* 199; base peak, *m/e* 152). Therefore, to minimize the number of products formed, in-

³ OV-17 and Chromosorb W were obtained from Applied Science Laboratories, Inc., State College, Pa. The 7% OV-17 was prepared in our laboratory using conventional methods (2.1 g OV-17/30 g Chromosorb W).

⁴ Trimethylanilinium hydroxide (0.1 M in methanol) (trimethylphenylammonium hydroxide) was obtained from Eastman Kodak Co., Eastman Organic Chemicals Division. Trimethylanilinium hydroxide (2 M in methanol) was prepared in our laboratory by concentrating 100 ml of 0.1 M reagent to 5 ml.

⁵ K & K Laboratories, Inc., Plainview, N.Y.




Figure 3—Chromatogram representing an on-column reaction of cocaine (21 μ g) and trimethylanilinium hydroxide (24-hr lapse after adding the methylating reagent to cocaine). Table I for identification of the peaks.

jection of the solution should be performed immediately after adding the methylating reagent to the suspected cocaine residue.

We also wish to report the relative GLC retention times of other commonly abused drugs frequently extracted at alkaline pH ranges (pH 8-10) which could possibly interfere with the confirmation of cocaine by this method. These values are reported in Table I and represent the relative retention times of the product(s) produced using the same on-column reac-

tion conditions performed with cocaine. Of the drugs examined, phenobarbital is the only one that could possibly interfere with the characteristic ecgonidine methyl ester peak (peak 3) under programmed temperature conditions. However, there is no interference between cocaine and phenobarbital when underivatized cocaine is analyzed under isothermal temperature conditions.

To test this method further, we analyzed a simulated street sample containing 6% cocaine hydrochloride, 19% quinine, and 75% dextrose. This analysis was carried out by performing a conventional alkaline extraction (pH 10) of 2.7 mg of the sample with chloroform-isopropanol (3:1), evaporation of the solvent, addition of 50 μ l of trimethylanilinium hydroxide to the residue, and immediate injection of 1 μ l. The resultant chromatogram resembled Fig. 1, with no interference from quinine.

In summary, we believe this method will be of value as a definitive confirmatory screening test for cocaine after first tentatively identifying underivatized cocaine using isothermal or programmed GLC temperature conditions.

(1) D. Boyse, N. Radin, D. S. Lewis, and G. Guerrent, "Proficiency Testing. Toxicology-Drug Abuse Survey III," U.S. Department of Health, Education, and Welfare, Public Health Service, Aug. 1973, p. 13.

(2) B. S. Finkle, D. M. Taylor, and E. J. Bonelli, *J. Chromatogr. Sci.*, 10, 312(1972).

(3) R. H. Hammer, B. J. Wilder, R. R. Streiff, and A. Mayersdorf, *J. Pharm. Sci.*, 60, 327(1971).

(4) R. J. Perchalski, K. N. Scott, B. J. Wilder, and R. H. Hammer, *ibid.*, 62, 1735(1973).

Richard H. Hammer *
James L. Templeton
H. L. Panzik

Department of Pharmaceutical Chemistry
College of Pharmacy
University of Florida
Gainesville, FL 32610

Received April 11, 1974.

Accepted for publication September 11, 1974.

* To whom inquiries should be directed.

Reprinted from

Journal of Chromatography
Elsevier Publishing Company, Amsterdam - Printed in The Netherlands
CHROM. 4886

ASSAY OF COCAINE IN THE PRESENCE OF PROCAINE AND QUININE BY COLUMN CHROMATOGRAPHY

GEORGE R. NAKAMURA

*Bureau of Narcotics and Dangerous Drugs, Regional Laboratory, San Francisco, Calif. 94102
(U.S.A.)*

AND

BRIAN P. PARKER

University of California, School of Criminology, Berkeley, Calif. 94720 (U.S.A.)

(Received June 18th, 1970)

SUMMARY

A column chromatography method is presented for the separation and quantitation of cocaine in the presence of procaine, quinine and lactose. Separation procedures are described to afford purification of cocaine for infrared absorption spectrum examination.

INTRODUCTION

Illicit cocaine is examined in the forensic laboratories as a fine, white powder, invariably occurring as a hydrochloride salt, occasionally crude-appearing brown samples have been examined. Adulteration and dilution are effected by illicit dealers with procaine and lactose respectively, and less frequently, quinine is used as an adulterant. The dilution has been mostly in the range of 10 to 60%.

The rapid detection tests for cocaine are described in the U.S. Treasury Manual¹. The preferred methods are in the use of cobalt thiocyanate test and in the formation of the characteristic microcrystals with platonic chloride. These tests afford a rapid and convenient procedure for the presumptive identification of cocaine.

Procaine is detected by the Sanchez color test¹ for cyclic aromatic amines resulting in a vivid red color. Quinine possesses an intense blue fluorescence under UV light, and is readily characterized by placing a sample preferably in an acid solution under fluorescent light.

Thin-layer method for chromatography is recommended if materials other than those mentioned, *i.e.*, procaine, quinine, and lactose are suspected. COMER AND COMER² have reviewed a number of papers describing the use of thin-layer chromatography (TLC) for the separation of different kinds of drugs. He has listed some sixteen different TLC systems for separating cocaine from a number of local anesthetics and

from a number of analgesics-antipyretics. In addition, this paper offers a simple TLC system using Eastman Chromatogram sheets for separating cocaine from procaine and quinine. We have not hitherto encountered adulterants other than procaine and quinine in the illicit samples.

YOUNG³ reported a chemical method for the determination of cocaine in the presence of procaine in which cocaine is hydrolyzed to yield methanol as a product. Methanol is distilled and is measured colorimetrically, using the well-known permanganate oxidation procedure. The method proved time-consuming for forensic work and a number of chloroform shake-out methods¹ were adapted to provide a final residual product suitable for volumetric titration. This procedure proved to be not only laborious but was attended by losses in recovery simply because of its many manipulative steps.

The manual¹ also describes a determination of cocaine, procaine, and tetracaine, employing a direct reading of the three components in alkaline solution and calculating the amounts by simultaneous equation.

The method is beset by the errors introduced by overlapping curves and these errors have been compensated by the introduction of experimental factors. The errors introduced for cocaine determination become more pronounced when proportionally more adulterant is present in the sample than cocaine.

The present method employs the principle of ion pair formation and its extraction by partition chromatography. This principle is discussed by HIGUCHI *et al.*⁴ and it has been adapted by LEVINE and co-workers⁵⁻⁸ for the separation of a large number of pharmaceutical amines.

Among the common anions tested, including chloride, phosphate and sulfate, nitrate proved to be the most efficient for quantitative purposes. Potassium nitrate in HCl solution is used as a stationary phase and chloroform for elution. HCl acts to retain quinine and possibly other basic materials which may be present in the cocaine sample. The eluant from the column is acidified and measured by UV absorption.

PROCEDURE

Mix 4 g of Celite 545 (Johns Mansville) with 2 ml of 1 M KNO₃ in 0.1 N HCl and transfer to a column, such as used by LEVINE⁵. Tamp the mixture on a pad of fine glass wool. Pipette 1 ml of aliquot of solution containing 50 mg of sample in 1 M KNO₃ in 0.1 N HCl onto the surface of the column. Irrigate the column with 45 ml of chloroform saturated with water (spectro-grade solvent preferred); collect the eluant in a 50-ml volumetric flask containing 5 ml of methanol and 5 drops of concentrated HCl. Bring the solution to 50 ml mark with chloroform and read its UV spectrum from 340 to 255 nm, max. at 275 nm, using chloroform as a reference.

Prepare a quantitative standard by dissolving 10 mg of cocaine · HCl in 50 ml chloroform containing 5 ml methanol and 5 drops HCl. Absorbance at 1%/1 cm was 32.5 using a Cary 15 spectrophotometer under these conditions. Calculate % cocaine, as hydrochloride, using the following equation:

$$\frac{A_{\text{sample}}}{A_{\text{standard}}} \times \frac{10}{50} \times 100 = \% \text{ cocaine} \cdot \text{HCl}$$

run for A = .5 (~.15 mg/ml)

Clean-up for IR examination

Shake the chloroform eluant saved from the assay with 10 ml of water. Decant the chloroform phase. Adjust the pH of the aqueous extract with dilute NH_4OH just to alkalinity (observe white precipitate) and extract with 50 ml chloroform twice. Filter the chloroform extract into a beaker and then treat the filtrate with sodium sulfate to remove water. Filter and evaporate the extract to dryness. Press the resulting crystalline material into KBr for IR spectrum examination for cocaine base.

DISCUSSION

While most forensic laboratories do not quantitate fillers or adulterants present in illicit narcotic samples, procaine and quinine can be eluted from the column using triethylamine in chloroform according to a procedure by LEVINE⁸. 1 N HNO_3 or 1 M KNO_3 can be used in the column *per se* if procaine is the only adulterant in the illicit cocaine.

In the clean-up procedure, there is the hazard of hydrolyzing cocaine to benzyl-ecgonine and methanol in the presence of water and heat⁹. Therefore, a desiccant such as sodium sulfate is used to dry the chloroform extract prior to evaporation over a hot steam bath.

For the clean-up, 10 mg of cocaine can be recovered by the method. If sufficient amount of material cannot be recovered, ca. 100 mg of sample can be passed through the column to recover enough material for IR examination.

On occasions, this laboratory has analyzed brown colored cocaine powder preparations containing procaine as an adulterant. The interference due to the presence of chloroform-soluble colored component was removed by washing the sample in the column, under standard conditions, with 50 ml of water-saturated 1,1-dichloroethane, which has been used in the clean-up of brown heroin¹⁰. The cocaine which was eluted subsequently with chloroform yielded a UV curve similar to those obtained with standard cocaine; recovery values of cocaine standards under standard conditions are shown in Table I.

TABLE I

RECOVERIES OF COCAINE STANDARDS

mg added	mg recovered	% recovery
15.00	14.89	99.2
12.00	11.90	99.2
10.00	9.95	99.5
8.00	7.91	99.1
5.00	4.92	98.9

HEAGY¹¹ reported that cocaine can be purified sufficiently for IR determination by its relative solubility in methylene chloride over procaine which is sparingly soluble in this solvent at room temperature. To a sample of illicit cocaine mixed with procaine, just enough methylene chloride is added to the mixture to dissolve cocaine. The extract is filtered immediately and evaporated to dryness. The residue is recrystallized with petroleum ether and pressed into KBr for IR spectrum reading.

TABLE II

TLC OF COCAINE, PROCAINE AND OTHER COMPOUNDS OF FORENSIC INTEREST

Solvent mixture: 40 ml chloroform, 10 ml ethyl acetate, and 10 drops conc. NH_4OH . Chromatography sheet: Eastman Chromatogram.

Compounds	Typical R_F
Cocaine	0.84
Procaine	0.45
Quinine	0.05
Heroin	0.34
Codeine	0.11
Morphine	0.00
LSD	0.34
Mescaline	0.05
STP	0.13
PCP	0.89
Nupercaine	0.45
Benzocaine	0.67
Tetracaine	0.11
dl-Amphetamine	0.00
DMT	0.10

In addition to those TLC solvent systems listed by COMERS AND COMERS² the authors suggest a very simple, rapid procedure using a solvent mixture containing 40 ml chloroform, 10 ml ethyl acetate and 10 drops of ammonium hydroxide. Eastman chromatoshets are cut to approximately 1½ in. and 4 in. and irrigated in the solvent mixture in a small bottle accommodating this sized sheet. The spots are revealed by using iodoplatinate reagent. Table II shows the typical R_F values obtained for cocaine, procaine, quinine and other forensically important compounds.

REFERENCES

- 1 *Methods of Analysis*, U. S. Treasury Department, Internal Revenue Service Publication Number 341, (Rev. 6-67).
- 2 J. P. COMER AND I. COMER, *J. Pharm. Sci.*, 56 (1967) 413.
- 3 J. L. YOUNG, *J. Assoc. Offic. Agr. Chemists*, 31 (1948) 781.
- 4 T. HIGUCHI, T. HIGUCHI, A. MICHAELIS, T. TAN AND A. HOROWITZ, *Anal. Chem.*, 39 (1967) 974.
- 5 J. LEVINE, *J. Assoc. Offic. Agr. Chemists*, 54 (1965) 488.
- 6 T. D. DOYLE AND J. LEVINE, *Anal. Chem.*, 39 (1967) 1282.
- 7 J. LEVINE AND R. T. OTTES, *J. Assoc. Offic. Agr. Chemists*, 44 (1961) 291.
- 8 R. T. OTTES AND J. LEVINE, *J. Assoc. Offic. Agr. Chemists*, 44 (1961) 293.
- 9 L. E. WENNER, *Micrograms*, 17 (1969) 47.
- 10 G. R. NAKAMURA AND H. J. MEURON, *Anal. Chem.*, 41 (1969) 993.
- 11 J. A. HEAGY, Bureau of Narcotics and Dangerous Drugs, San Francisco, private communication.


Southcombe

5-6-77

Cocaine

.200 grams (200mg) 1/5 of a gram

1. Place comparable amount in a plugged pipet - compare with known amount
2. Flick with finger to shake ~~shown~~ down into pipet - tap lightly
3. Mark with a "Sharpie" the top of the powder and the bottom interface with the glass wool.



4. Place 2 droppers full of methylene chloride into pipet - let drain into waste beaker.

5. At this point, you may place pipet ^{draining} over a spot test well - let a few drops fall into the well and let dry. Add Cobalt thiocyanate reagent to well, and observe reaction. - OR SCOTT'S test (3 stage)

6. Compare remainder ^{of powder} with previous marks - make decision on %.



Example

- remaining powder approx 50%



DEA LABORATORY NOTES

DATE October 26, 1976
 NO.
 DRUG TYPE Cocaine
 METHODOLOGY Salt Formation

PURIFICATION OF COCAINE BY SALT FORMATION

GARY J. SORGEN
 FORENSIC CHEMIST
 DRUG ENFORCEMENT ADMINISTRATION
 WESTERN REGIONAL LABORATORY
 SAN FRANCISCO, CALIFORNIA

OBJECTIVE

To provide a simple purification of cocaine mixed with other caines.

BACKGROUND

Cocaine is often mixed with other drugs, commonly with lidocaine and procaine. Cocaine and lidocaine can be easily separated with aqueous platinum chloride.¹ However, there is no simple method for separating cocaine from lidocaine and procaine.

METHOD

The cocaine sample is dissolved in water and made basic with ammonia. The caines are then extracted with dichloromethane. Evaporate off the solvent and dissolve the residue in a minimum amount of acetone. Add an equal volume of Di-p-toluoyl-(-)-tartaric acid solution (about 50mg per ml of acetone). Vigorously stir and scratch the bottom of the beaker with a glass capillary tube until crystals form (2-30 sec.) (salt of cocaine and TTA). If no crystals form, concentrate the solution and try again. A small amount of pet ether may also help to initiate crystallization.

Once crystals of the salt form, they are extremely insoluble and can be washed with acetone. The cocaine can be recovered from the salt by dissolving the salt in dilute ammonia and extracting the cocaine with dichloromethane. The purity of the cocaine residue is indicated by how quickly the oily residue solidifies when stirred with a glass rod. The more quickly it solidifies the purer the cocaine.

DRUG ENFORCEMENT ADMINISTRATION / U. S. DEPARTMENT OF JUSTICE
MICROGRAM, VOL. X, NO. 4, APRIL, 1977

PURIFICATION OF COCAINE BY SALT FORMATION

METHOD - continued


If the cocaine is not purified enough, a simple micro alumina column clean up can be used for further purification. To a disposable pipet plugged with cotton add about one inch of chromatographic aluminum oxide (alumina) powder. Wet the alumina with one milliliter of dichloromethane. Dissolve the cocaine residue in a minimum amount of dichloromethane and add the solution to the column. Elute the cocaine with dichloromethane. Collect the eluate in several beakers. Collect one to two milliliters in the first beaker, about three to five in the second, and about five in the third beaker. Evaporate off the solvent. Most of the cocaine will be in the second beaker.

DISCUSSION

Cocaine has been purified from a mixture of cocaine, procaine, and lidocaine hydrochlorides by the formation of the salt of (-)-cocaine and (+)-Di-p-toluoyl-(-)-tartaric acid. The other optical isomer, (-)-Di-p-toluoyl-(+)-tartaric acid, also reacts with (-)-cocaine but not as quickly. The Di-p-toluoyl-l-tartaric acid from Aldrich is the monohydrate and is soluble in acetone. However, the Di-p-toluoyl-d-tartaric acid is anhydrous and as such is insoluble in acetone. To dissolve add a drop of water to the acetone. Di-p-toluoyl-(-)-tartaric acid is available from Aldrich, Fluka, and Pfaltz & Bauer chemical companies.


REFERENCES

- 1 Power, Joseph A., Jr. DEA, personal communication.




(-) - Cocaine - di - p - toluoyl - (-) - tartrate

(-) - Cocaine - di-p-toluoyl (-) -tartrate




Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OPI-DELETE DOCUMENT

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT



Di-p-toluoyl (-)-tartaric acid

DEA LABORATORY NOTES



DATE
NO.
DRUG TYPE
METHODOLOGY

ISOLATION AND DETERMINATION OF l, d, or dl-COCAINE IN SMALL AMOUNTS OF ILLICIT SAMPLES

by

Paul L. Morgan
Vedoster Ingram
Daniel Francois
Mid-Atlantic Regional Laboratory
Washington, D. C.

OBJECTIVE

To devise a method for the determination of the d, l, or dl isomer of cocaine by using prep-TLC and mixed melting points under the following sample conditions:

1. The total sample quantity is very small
i.e. as little as 50mg.
2. The percentage of cocaine is very low
i.e. as low as 3.3%.
3. The cocaine is found singly or in combination
with other drugs.

BACKGROUND

Recently, the analysis of cocaine has received special attention in order to determine its isomeric form. The isomeric form is generally determined by use of a polarimeter or a mixed melting point determination.

When sufficient quantities of sample are available and the percentage of cocaine high enough, the isomeric determination is readily

DRUG ENFORCEMENT ADMINISTRATION / U. S. DEPARTMENT OF JUSTICE
MICROGRAM, VOL. XI, NO. 3 (MARCH 1978)

accomplished. However, in the absence of such quantities and with low percentages, the task is far more difficult. This is especially true in the presence of other drugs. Under these conditions it is difficult to obtain information concerning the isomeric form of cocaine by use of conventional means.

This paper presents a method which allows such a determination under the stated sample conditions. It employs prep-TLC and mixed melting points of the sample and procedural standard l-cocaine.

REAGENTS AND APPARATUS

Chloroform, Methanol, NH_4OH (concentrated), petroleum ether, anhydrous sodium sulfate, Pasteur type pipets, Thomas Hoover Capillary Melting Point Apparatus, glass wool, and TLC plates (Silica Gel GF, 250 microns thick from Analteck, Inc. Newark, Delaware).

PROCEDURE

Dry extract the sample with approximately 0.5ml of methanol and spot across two or more TLC plates using a Pasteur type pipet. Standard cocaine may be spotted next to the sample for a rate of flow (R_f) comparison. Develop the plates in an appropriate solvent system which separates cocaine from other components in the sample.

Dry the plates thoroughly on a warm hot plate aided by a gentle flow of air. Under short wave u.v. light locate and outline the cocaine area.

Prepare a mini column by packing glass wool (1cm) in a 15cm Pasteur type pipet. Add anhydrous sodium sulfate (1cm) on top of the glass wool.

Scrap off the outlined cocaine area from the dried plates. This may be conveniently done by using a flat head type spatula. Transfer the scraping to the mini column on top of the anhydrous sodium sulfate. The transfer to the mini column can be greatly facilitated by use of a vacuum. This mini column preparation is essentially that described in Microgram, Vol. IX, No. 9, pp. 130-135 (1).

Prepare an eluting solvent by adding 2-3 drops of methanol to 15ml of chloroform saturated with concentrated NH_4OH .

Elute the cocaine from the mini column using 1-2ml of the eluting solvent. Collect the eluate in a mortar dish and take to dryness on a warm hot plate aided by a gentle flow of air.

After drying thoroughly, scrape the residue thus forming a powder. Add 1-2 drops of petroleum ether and take to dryness. The sample is now ready for a melting point determination.

Determine and compare the melting points of the sample, procedural standard l-cocaine and a mixture of the two.

* Standard l-cocaine must be treated in the same manner as the sample and may be done concomitantly with the sample. This is referred to as the procedural standard l-cocaine.

RESULTS AND DISCUSSION

<u>Substance</u>	<u>Melting Points</u>
Sample Cocaine	92 - 95°
Sample plus procedural standard l-cocaine	92 - 95°
Procedural standard l-cocaine	92 - 95°
Non-procedural standard l-cocaine	95 - 97°
Literature values	
l-cocaine (3)	96 - 98°
dl-cocaine (2)	79 - 80°

The usefulness of this method can be readily seen whenever the total sample quantity and percentage of cocaine is small enough to obviate use of the polarimeter.

The presence of other drugs in samples limited by quantity and low percentages of cocaine can, and may often dictate the use of prep-TLC as the only feasible means of separation for isomeric determination.

MICROGRAM, VOL. XI, NO. 3 (MARCH 1978)

REFERENCES

1. P. L. Morgan, D. Francois, Microgram, Vol. IX, No. 9 (Sept. 1976), pp. 130-135.
2. A. Allen, D. Cooper, Microgram, Vol. IX, No. 5 (May, 1976), pp. 66-72.
3. E.G.C. Clarke, Isolation and Identification of Drugs, p. 267.

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

James Conrad Roberson

(39)

Introduction

In the great tradition of the cannabis variety debate, we now see a rise in the question of isomers of cocaine in forensic analysis. The following outlines the legal ramifications of that question and analytical steps taken to avoid the problem.


Legal

On January 14, 1977, Mr. John Partain petitioned (Civil Action No. C77-67A) the United States District Court in Atlanta to overturn his prior conviction of possession of cocaine in violation of the Georgia Controlled Substance Act, Section 79A-807, Schedule II, b4. That Act is patterned after the Federal Schedule and reads as follows:

"(4) Coca leaves, any salt, compound, derivative, or preparation of coca leaves, and any salt, compound, derivative, or preparation thereof which is chemically equivalent or identical with any of these substances, but not including decocainized coca leaves or extractions which do not contain cocaine or ecgonine."

The sample consisted of a trace of cocaine on balance pans whose presence was established by Mr. William Price of the Georgia State Crime Laboratory using UV, TLC, GC, and mass spectra. Dr. Robert Shapiro of the University of Colorado testified that this evidence was insufficient to differentiate naturally occurring l-cocaine from d-cocaine, pseudococaine, allococaine, allospseudococaine, or alpha cocaine. The Georgia State Crime Laboratory contended that mass spectra would differentiate cocaine from its diastereomers and that d-cocaine would be covered by the Controlled Substance Act as a derivative of l-cocaine.

According to the definitions (1, 2, 3) of derivative, it must have the same basic structure as and must be theoretically derivable from the parent compound. The tropane structure is common to both d- and l-cocaine. The following are theoretical pathways from l- to d-cocaine.



Literature ("Introduction to Stereochemistry" by Kurt Mislow) was quoted to indicate that IR or mass spectra would differ between diastereomers. The D.E.A. laboratory in Miami has compared IR and mass spectra of cocaine and pseudococaine and found differences. Pseudococaine is the diastereomer most likely to be spectrally similar to cocaine and is the only other isomer manufactured commercially.

Probably the most important legal point was simply the rarity of isomers other than l-cocaine. Dr. Shapiro admitted that if a person wished to obtain these isomers he would have to synthesize them himself by a difficult chemical procedure. Therefore it is unreasonable to consider their possible presence in a sample in this country.

In his order of August 26, 1977, Judge Newell Edenfield denied the petition of Civil Action No. C77-67A.

Analytical

Mixture melting point⁽⁶⁾ and NMR spectra⁽⁷⁾ are methods of differentiating cocaine and its diastereomers that are fully documented in the literature. The multiplet at 5.27 δ (ppm) in the NMR spectrum is characteristic of cocaine free base. The following clean-up procedure is very useful.

1. Wash sample powder with 5 to 10 portions of 1,4-dioxane (removes lidocaine).
2. Wash sample powder with 2 or 3 small portions of acetone (removes procaine or tetracaine).
3. Dissolve sample in 0.1N HCl, add NaHCO₃ to pH 7, and extract cocaine base with CHCl₃.
4. Evaporate CHCl₃ and dissolve residue in CDCl₃ for NMR spectrum.

Very little cocaine is lost in the dioxane wash, but the acetone wash removes a good deal of cocaine and should be avoided if possible. An alternative procedure is to dissolve the sample in 0.1N boric acid and extract the cocaine with CHCl₃, leaving the procaine in the acid solution.

References

1. Webster's New Collegiate Dictionary, 1973, p. 306.
2. "International Encyclopedia of Chemical Science", D. Van Nostrand Co., Inc., Princeton, N. J., 1964, p. 337.
3. "Hackh's Chemical Dictionary", McGraw-Hill Book Company, New York, N.Y., 1969, p. 203.
4. Communication from A. Allen, D.E.A., Miami, Fla.
5. Fodor, G., Tropane Alkaloids in "Chemistry of the Alkaloids", Van Nostrand Reinhold Co., New York, N. Y., 1970, pp. 440-441.
6. Willstatter, R., et al., Ann. Chem., 434, 111, 1923.
7. Sinnema, A., et al., Recl. Trav. Chim. Pays-Bas, 87, 1027, 1968.

From: The Georgia State Crime Laboratory
959 East Confederate Avenue, S. E.
Atlanta, Georgia 30316

COMPREHENSIVE LIST OF CONTROLLED SUBSTANCES

The Uniform Controlled Substances Act, Title 37, Chapter 27, Section 37-2414, Idaho Code, states that the Board shall revise and republish the schedules semi-annually for two years from the effective date of this act and thereafter annually. Therefore, pursuant to the mandate of that section, the Executive Secretary of the Idaho State Board of Pharmacy hereby orders the annual publication of the schedules of controlled substances. Chemical or trade names follow generic names in parentheses.

SCHEDULE I. (a) The controlled substances listed in this section are included in schedule I.

(b) Any of the following opiates, including their isomers, esters, ethers, salts, and salts of isomers, esters, and ethers, unless specifically excepted, whenever the existence of these isomers, esters, ethers and salts is possible within the specific chemical designation:

- (1) Acetylmethadol; (Acemethadone, Amidolacetate)
- (2) Allylprodine; (Alperidine, NIH-7440, Ro2-7113)
- (3) Alphacetylmethadol; (see Acetylmethadol)
- (4) Alphameprodine; (NU-1932)
- (5) Alphamethadol;
- (6) Benzethidine;
- (7) Betacetylmethadol; (see Acetylmethadol)
- (8) Betameprodine; (NU-1932)
- (9) Betamethadol; (Betametadol)
- (10) Betaprodine; (NU-1779)
- (11) Clonitazene;
- (12) Dextromoramide; (Palfium, Jetricum, Pyrrolamidol, R-B75, SKF-d, 5137)
- (13) Diampromide;
- (14) Diethylthiambutene; (Themalon)
- (15) Difenoxin;
- (16) Dimenoxadol; (NIH-7577, Lokarin)
- (17) Dimpheptanol; (Methadol, Pangerin, Amidol, NIH-2933, Amidalgon)
- (18) Dimethylthiambutene; (Kobaton, Ohton, Skikiton, Takaton)
- (19) Dioxaphetyl butyrate; (Amidalgon, Spasmoxale)
- (20) Dipipanone; (Fenpidon, Pamedon, Pipadone, 378C48, HOECHST 10805)
- (21) Ethylmethylthiambutene; (Emethibutin, 1C50, NIH-5145)
- (22) Etonitazene;
- (23) Etoxadine; (Atenorax, Atenos, Cargetidine)
- (24) Furethidine;
- (25) Hydroxypethidine; (Bemidone)
- (26) Ketobemidone; (Ketogan, Cliradon)
- (27) Levomoramide;
- (28) Levophenacilmorphan; (NIH-7525, Ro4-0288)
- (29) Morpheridine;
- (30) Noracymethadol;
- (31) Norlevorphanol; (NIH-7539)
- (32) Normethadone; (Ticarda, HOECHST 10582, Mepidon, Deatussan, Normedon, Phenylidimazone Veryl)
- (33) Norpipanone; (HOECHST 10, 495, Hexalgon, Orfenso)
- (34) Phenadoxone; (CB-11, Hepagin, Heptalgin, Heptalin, Heptan, Heptone, Heptazone, Supralgin, HOECHST 10, 600)
- (35) Phenampromide;
- (36) Phenomorphan; (NIH-7274)
- (37) Phenoperidine; (Lealgin)

- (38) Piritramide; (A65, Dipidolor)
- (39) Proheptazine; (Proheptazone)
- (40) Properidine; (Gevelina, Ipropethidine, Isopedine, Spasmodolosa)
- (41) Propiram;
- (42) Racemoramide; (R-610)
- (43) Trimeperidine; (Promedol).

(c) Any of the following opium derivatives, their salts, isomers and salts of isomers, unless specifically excepted, whenever the existence of these salts, isomers and salts of isomers is possible within the specific chemical designation:

- (1) Acetorphine;
- (2) Acetyldihydrocodeine; (Acetylcodeine)
- (3) Benzylmorphine; (Peronine)
- (4) Codeine methylbromide; (Eucodin)
- (5) Codeine-N-Oxide; (Genocodeine, Codeigene)
- (6) Cyprenorphine;
- (7) Desomorphine; (Dihydrodesoxymorphine-D, Permonid)
- (8) Dihydromorphine; (Paramorfan)
- (9) Drotebanol;
- (10) Etorphine, except hydrochloride salt; (M99, M183)
- (11) Heroin; (Diacetylmorphine)
- (12) Hydromorphanol;
- (13) Methyldesorphine;
- (14) Methyldihydromorphine; (Metopan)
- (15) Morphine methylbromide; (Morphosan)
- (16) Morphine methylsulfonate;
- (17) Morphine-N-Oxide; (Genomorphine)
- (18) Myrophine; (Leucodine, Myricodine, Peronine myristate)
- (19) Nicocodeine;
- (20) Nicomorphine; (Vilan)
- (21) Normorphine; (Desmethymorphine)
- (22) Pholcodine; (Ethnine, Glycodine, Memine, Codylin, Hibernyl, Pectolin, Prodromine, Weifacodine)
- (23) Thebacon.

(d) Any material, compound, mixture or preparation which contains any quantity of the following hallucinogenic substances, their salts, isomers and salts of isomers, unless specifically excepted, whenever the existence of these salts, isomers, and salts of isomers is possible within the specific chemical designation:

- (1) 3,4-methylenedioxy amphetamine;
- (2) 5-methoxy-3, 4-methylenedioxy amphetamine; (MDA)
- (3) 3,4,5 -trimethoxy amphetamine;
- (4) Bufotenine; (Mappine)
- (5) Diethyltryptamine; (DET)
- (6) Dimethyltryptamine; (DMT)
- (7) 4-methyl-2,5-dimethoxyamphetamine; (STP, DOM, DMA)
- (8) Ibogaine;
- (9) Lysergic acid diethylamide; (LSD)
- (10) Marihuana;
- (11) Mescaline;
- (12) Peyote;
- (13) N-ethyl-3-piperidyl benzilate;
- (14) N-methyl-3-piperidyl benzilate;
- (15) Psilocybin;
- (16) Psilocyn;
- (17) Tetrahydrocannabinols; (THC); (Synthetic equivalents of the substances contained in the plant, or in the resinous extractives of Cannabis, and/or synthetic substances, derivatives, and their isomers with similar chemical structure and pharmacological activity;)
- (18) 2,5-dimethoxyamphetamine (2,5-dimethoxy-a-methylphenethylamine: 2,5-DMA);
- (19) 4-bromo-2,5-dimethoxyamphetamine (4-bromo-2,5-dimethoxy-a-methylphenethylamine: 4-bromo-2,5,DMA);

(20) 4-methoxyamphetamine (4-methoxy-a-methylphenethylamine; paramethoxyamphetamine, PMA).

(21) Thiophene analog of phencyclidine (1-(1-(2-thienyl) cyclohexyl) piperidine).

(e) Any material, compound, mixture or preparation which contains any quantity of the following substances having a depressant effect on the central nervous system, including its salts, isomers, and salts of isomers wherever the existence of such salts, isomers, and salts of isomers is possible within the specific chemical designation:

(1) Mecloqualone.

SCHEDULE II. (a) Schedule II shall consist of the drugs and other substances, by whatever official name, common or usual name, chemical name, or brand name designated, listed in this section.

(b) Substances, vegetable origin or chemical synthesis. Unless specifically excepted or unless listed in another schedule, any of the following substances whether produced directly or indirectly by extraction from substances of vegetable origin, or independently by means of chemical synthesis, or by a combination of extraction and chemical synthesis:

(1) Opium and opiate, and any salt, compound, derivative, or preparation of opium or opiate, excluding naloxone and its salts, and naltrexone and its salts, but including the following:

1. Raw Opium;
2. Opium extracts;
3. Opium fluid extracts;
4. Powdered opium;
5. Granulated opium;
6. Tincture of opium;
7. Codeine; (methymorphine)
8. Ethylmorphine;
9. Etorphine hydrochloride;
10. Hydrocodone;
11. Hydromorphone; (Dilaudid)
12. Metopon;
13. Morphine; (Papine, H-M-C #1 & #2)
14. Oxycodone; (Percodan, Percobarb, Nucodan)
15. Oxymorphone; (Numorphan)
16. Thebaine.

(2) Any salt, compound, derivative, or preparation thereof which is chemically equivalent or identical with any of the substances referred to in paragraph (b) (1) of this section, except that these substances shall not include the isoquinoline alkaloids of opium.

(3) Opium poppy and poppy straw.

(4) Coca leaves and any salt, compound, derivative, or preparation of coca leaves, and any salt, compound, derivative, or preparation thereof which is chemically equivalent or identical with any of these substances, but not including decocainized coca leaves or extractions which do not contain cocaine or ecgonine.

Methylbenzoyllecgonine (Cocaine).

(5) Concentrate of poppy straw (the crude extract of poppy straw in either liquid, solid or powder form which contains the phenanthrine alkaloids of the opium poppy).

(c) Any of the following opiates, including their isomers, esters, ethers, salts, and salts of isomers, whenever the existence of these isomers, esters, ethers and salts is possible within the specific chemical designation:

- (1) Alphaprodine; (Nisentil HCL, Nisintil, Prisilidene)
- (2) Anileridine; (Leritine, Nipecotan, Alidine, Apodol)
- (3) Bezitramide;
- (4) Dihydrocodeine; (Codhydrine, Dehacodin, DF118, DH-Codeine, Didrate, Hydrocodone, Dihydrin, Parazone, Hydrocodin, Paracodin)
- (5) Diphenoxylate;

- (6) Fentanyl; (R4263, Phentanyl, Pentanyl, Sublimaze)
- (7) Isomethadone; (Liden, Isoadanone)
- (8) Levomethorphan;
- (9) Levorphanol; (Aromarine, Dromoran, Levo-Dromoran, Levorphan)
- (10) Metazocine; (NIH 7539)
- (11) Methadone; (Amidone, Dolophine, Methadon, Methajade)
- (12) Methadone--Intermediate, 4-cyano-2-dimethylamino-4, 4-diphenyl butane;
- (13) Moramide--Intermediate, 2-methyl-3-morpholino-1, 1-diphenyl propane-carboxylic acid;
- (14) Pantopon (Hydrochlorides of opium alkaloids);
- (15) Pethidine; (Demerol, Isonipecaine, Meperidine, Mepergan, Demerol APAP)
- (16) Pethidine--Intermediate--A, 4-cyano-1-methyl-4-phenylpiperidine;
- (17) Pethidine--Intermediate--B, ethyl-4-phenyl-piperidine-4-carboxylate;
- (18) Pethidine--Intermediate--C, 1-methyl-4-phenylpiperidine-4-carboxylic acid;
- (19) Phenazocine; (NIH-7519, SKF 6574, Narphen, Prinadol)
- (20) Piminodine; (Alvodine, Anopridine, Cimadon, Pimadine, NIH 7590, WIN-14098)
- (21) Racemethorphan;
- (22) Racemorphan; (Citarin, Methorphanin, Levodromoran)

(d) Stimulants. Unless specifically excepted or unless listed in another schedule, any material, compound, mixture, or preparation which contains any quantity of the following substances having a stimulant effect on the central nervous system:

- (1) Amphetamine, its salts, optical isomers, and salts of its optical isomers;
 1. Dextroamphetamine/amphetamine; (Biphentamine)
 2. Amphetamine Sulfate; (Benzedrine).
 3. Dextroamphetamine sulfate, (Dexedrine, Dextro-Amphetamine Phosphate, Dextro-Amphetamine Sulfate, Amphetamine HCl, D.A.S., Synatan, Desarex, Dex-Sule).

Combinations:

1. Dextroamphetamine sulfate & amobarbital; (Amodex, Dexamyl, Dexobarb, Daprisal, Obozell)
2. Dextroamphetamine sulfate & prochlorperazine; (Eskatrol, Bamadex, Appetrol, Amvical, Amvical X).
- (2) Methamphetamine, its salts, isomers, and salts of its isomers; (Desoxya, Fetamin).
- (3) Phemetrazine and its salts; (Preludin).
- (4) Methyphenidate. (Ritalin HCl).

(e) Depressants. Unless specifically excepted or unless listed in another schedule, any material, compound, mixture, or preparation which contains any quantity of the following substances having a depressant effect on the central nervous system, including its salts, isomers, and salts of isomers whenever the existence of such salts, isomers, and salts of isomers is possible within the specific chemical designation:

- (1) Methaqualone; (Quaalude, Parest, Sopor, Optimil, Somnafac-200, Forte-400)
- (2) Amobarbital; (Amytal)
- (3) Secobarbital; (Tuinal)
- (4) Pentobarbital; (Nembutal)

SCHEDULE III. (a) Schedule III shall consist of the drugs and other substances, by whatever official name, common or usual name, chemical name, or brand name designated, listed in this section.

(b) Stimulants. Unless specifically excepted or unless listed in another schedule, any material, compound, mixture, or preparation which contains any quantity of the following substances having a stimulant effect on the central nervous system, including its salts, isomers, and salts of its isomers, (whether optical or geometric), and salts of such isomers whenever the existence of such salts, isomers, and salts of

isomers is possible within the specific chemical designation:

(1) Those compounds, mixtures, or preparations in dosage unit form containing any stimulant substances listed in schedule II which compounds, mixtures, or preparations were listed on August 25, 1971, as excepted compounds under C.F.R. Sec. 308.32, and any other drug of the quantitative composition shown in that list for those drugs or which is the same except that it contains a lesser quantity of controlled substances.

- (2) Benzphetamine; (Didrex)
- (3) Chlorphentermine; (Pre-Sate)
- (4) Clortermine; (Voramil)
- (5) Mazindol; (Sanorex)
- (6) Phendimetrazine; (Apidex, Edrisal, Genegestic, Mediatric, Thora-Dex #1 & #2, Ropledge, Bacarate, Banobese, Bontril PDM, Melfiat)

(c) Depressants. Unless listed in another schedule, any material, compound, mixture, or preparation which contains any quantity of the following substances having a potential for abuse associated with a depressant effect on the central nervous system:

- (1) Any compound, mixture or preparation containing:
 - i. Amobarbital; (Amytal)
 - ii. Secobarbital; (Tuinal, Efed)
 - iii. Pentobarbital; (Nembutal, Carbrital),or any salt thereof and one or more other active medicinal ingredients which are not listed in any schedule.
- (2) Any suppository dosage form containing:
 - i. Amobarbital; (Alurate)
 - ii. Secobarbital;
 - iii. Pentobarbital; (Emersert)or any salt of any of these drugs and approved by the Food and Drug Administration for marketing only as a suppository.
- (3) Any substance which contains any quantity of a derivative of barbituric acid or any salt thereof
- (4) Chlorhexadol; (Lora, Mecoral, Medodorm)
- (5) Glutethimide; (Doriden, Darmtabs, Rolathimide)
- (6) Lysergic acid;
- (7) Lysergic acid amide; (LSD)
- (8) Methyprylon; (Noludar)
- (9) Phencyclidine; (Sernylan)
- (10) Sulfondiethylmethane;
- (11) Sulfonethylmethane; (Trional, Ethyl Sulfonal)
- (12) Sulfonmethane; (Sulfonal)

(d) Nalorphine; (N-Allyl Nor Morphine, N Alline, Norfin, NANM)

(e) Narcotic Drugs. Unless specifically excepted or unless listed in another schedule, any material, compound, mixture, or preparation containing limited quantities of any of the following narcotic drugs, or any salts thereof:

(1) Not more than 1.8 grams of codeine, or any of its salts, per 100 milliliters or not more than 90 milligrams per dosage unit, with an equal or greater quantity of an isoquinoline alkaloid of opium, and

(2) Not more than 1.8 grams of codeine, or any of its salts, per 100 milliliters or not more than 90 milligrams per dosage unit, with one or more active, non-narcotic ingredients in recognized therapeutic amounts;

1. Codempiral
2. Empirin Compound #1, #2, #3, and #4
3. APC with Codeine
4. Copavin
5. Anexia with codeine
6. Phenaphen with codeine

(3) Not more than 300 milligrams of dihydrocodeinone, or any of its salts, per 100 milliliters or not more than 15 milligrams per dosage unit, with a fourfold or greater quantity of an isoquinoline alkaloid of opium;

(4) Not more than 300 milligrams of dihydrocodeinone, or any of its salts, per 100 milliliters or not more than 15 milligrams per dosage unit, with one or more active, nonnarcotic ingredients in recognized therapeutic amounts;

(5) Not more than 1.8 grams of dihydrocodeine, or any of its salts, per 100 milliliters or not more than 90 milligrams per dosage unit, with one or more active, nonnarcotic ingredients in recognized therapeutic amounts;

(6) Not more than 300 milligrams of ethylmorphine, or any of its salts, per 100 milliliters or not more than 15 milligrams per dosage unit, with one or more ingredients in recognized therapeutic amounts;

(7) Not more than 500 milligrams of opium per 100 milliliters or per 100 grams, or not more than 25 milligrams per dosage unit, with one or more active, nonnarcotic ingredients in recognized therapeutic amounts; (Paregoric).

(8) Not more than 50 milligrams of morphine, or any of its salts, per 100 milliliters or per 100 grams with one or more active, nonnarcotic ingredients in recognized therapeutic amounts.

(f) The board may except by rule any compound, mixture, or preparation containing any stimulant or depressant substance listed in sections (b) and (c) of this section from the application of all or any part of this act if the compound, mixture, or preparation contains one or more active medicinal ingredients not having a stimulant or depressant effect on the central nervous system, and if the admixtures are included therein in combinations, quantity, proportion, or concentration that vitiate the potential for abuse of the substances which have a stimulant or depressant effect on the central nervous system.

SCHEDULE IV. (a) Schedule IV shall consist of the drugs and other substances, by whatever official name, common or usual name, chemical name, or brand name designated, listed in this section.

(b) Depressants. Unless specifically excepted or unless listed in another schedule, any material, compound, mixture, or preparation which contains any quantity of the following substances, including its salts, isomers, and salts of isomers whenever the existence of such salts, isomers, and salts of isomers is possible within the specific chemical designation:

- (1) Barbital; (Barbital Sodium, Veronal)
- (2) Chloral betaine; (Beta-Chlor, Somilan)
- (3) Chloral hydrate; (Noctec, Somnos, Felsules)
- (4) Chlordiazepoxide; (Librium, Libritabs)
- (5) Clonazepam; (R05-4023, Clonopin)
- (6) Clorazepate; (Tranxene)
- (7) Dextropropoxyphene; (Darvocet-N50, N100, Darvon Compound, Dolene AP-65)
- (8) Diazepam; (Valium)
- (9) Ethchlorvynol; (Placidyl)
- (10) Ethinamate; (Valmid)
- (11) Flurazepam; (Dalmane)
- (12) Mebutamate;
- (13) Meprobamate; (Equinal, Miltown, Kesso-bamate, Deprol, Meprospan)
- (14) Methohexital; (Brevital Sodium)
- (15) Methylphenobarbital;
- (16) Oxazepam; (Serax)
- (17) Paraldehyde; (Paral)
- (18) Petrichloral; (Periclor)
- (19) Phenobarbital; (Eskabarb, Luminal, Phenobarbital Sodium)

(c) Fenfluramine - Any material, compound, mixture, or preparation which contains any quantity of the following substances, including its salts, isomers (whether optical, position, or geometric), and salts of such isomers, whenever the existence of

such salts, isomers, and salts of isomers is possible:

- (1) Fenfluramine.

(d) Stimulants. Unless specifically excepted or unless listed in another schedule, any material, compound, mixture, or preparation which contains any quantity of the following substances having a stimulant effect on the central nervous system, including its salts, isomers (whether optical, position, or geometric), and salts of such isomers whenever the existence of such salts, isomers, and salts of isomers is possible within the specific chemical designation:

- (1) Diethylpropion; (Ro-diet)
- (2) Phentermine; (Rolapent, Adipex-P, Fastin, Ionamin)
- (3) Pemoline (including organometallic complexes and chelates thereof. (Cylert)

(e) The board may except by rule any compound, mixture, or preparation containing any depressant substance listed in subsection (b) of this section from the application of all or any part of this act if the compound, mixture or preparation contains one or more active medicinal ingredients not having a depressant effect on the central nervous system, and if the admixtures are included therein in combinations, quantity, proportion, or concentration that vitiate the potential for abuse of the substances which have a depressant effect on the central nervous system.

SCHEDULE V. (a) Schedule V shall consist of the drugs and other substances, by whatever official name, common or usual name, chemical name, or brand name designated, listed in this section.

(b) Narcotic drugs containing nonnarcotic active medicinal ingredients. Any compound, mixture, or preparation containing any of the following limited quantities of narcotic drugs or salts thereof, which shall include one or more nonnarcotic active medicinal ingredients in sufficient proportion to confer upon the compound, mixture, or preparation, valuable medicinal qualities other than those possessed by the narcotic drug alone:

- (1) Not more than 200 milligrams of codeine per 100 milliliters or per 100 grams; (Cosanyl, Cosadaine, Cheracol, TH&C, Robitussin AC, Histadyl Ec)
- (2) Not more than 100 milligrams of dihydrocodeine per 100 milliliters or per 100 grams;
- (3) Not more than 100 milligrams of ethylmorphine per 100 milliliters or per 100 grams;
- (4) Not more than 2.5 milligrams of diphenoxylate and not less than 25 micrograms of atropine sulfate per dosage unit; (Lomotil)
- (5) Not more than 100 milligrams of opium per 100 milliliters or per 100 grams.

NOTE: The Drug Enforcement Administration computer list of controlled substances is in excess of 525 pages in length and contains about 7,000 drug products in various dosage forms. This list is on file in the office of the Idaho Board of Pharmacy.

Forensic Sciences
pp. 868-872
Vol. 19
no. 4;
Oct. 1974

Sumer Dugar,¹ Ph.D.; Thomas Catalano,¹ B.S.; and Richard Cerrato,¹ B.S.

Discoloration Effect of Diluents In Contraband Cocaine

In the course of police investigations of confiscated drug samples, cocaine mixtures of various shades from white to brown have been encountered and analyzed at the Crime Laboratory of the New York City Police Department. Slight discoloration of cocaine mixtures is within reason, in view of the reported discoloration of lactose [1,2] and of lactose-amphetamine [3-5] upon storage.

Recently a very interesting case involving cocaine was resubmitted to the laboratory, wherein a drastic change in the physical appearance of the confiscated evidence was observed. The evidence, originally a tan powder, was analyzed in July 1972. It contained aspirin tablets along with the contraband cocaine powder. The other diluents present were not analyzed at that time. This case was reanalyzed in February 1973 and at that time the powder was brown in color, as compared to the original tan powder. Finally, during a court trial in June 1973, the evidence was opened and found to be a dark tarry substance and the aspirin tablets were no longer present. Due to the discrepancy in the physical appearance of the evidence from the time of the original analysis to its appearance in court, the court had reservations in accepting the evidence unless a reasonable scientific explanation was furnished for this change. Since the trial was in progress, the court permitted a maximum of two days to conduct the related experiments with the evidence in question so that an explanation for this drastic change could be ascertained.

This paper presents the experimental work done with evidence in conjunction with the studies on the discoloration effects of known cocaine and the common diluents associated with it. The interaction of the diluents, which causes discoloration and tarry transformation, is also discussed.

Experimental

Studies were conducted in two parts: The first part deals with the analysis of the evidence and the second part deals with a study on the discoloration of various combinations of diluents commonly associated with contraband cocaine.

Analysis of Evidence

1. Positive cocaine identification was determined from the alkaline chloroform extraction of the water-soluble part of the tarry substance. Identification was made by

Received for publication 17 Sept. 1973; revised manuscript received 7 Jan. 1974; accepted for publication 22 Jan. 1974.

¹Chemist, police officer, and police officer, respectively, New York City Police Department Crime Laboratory, New York, N.Y.

color and mic
chromatograp
2. The tarr
acetate ions w
3. 5-hydrox
lactose, espec
these furfural
in acid medii
identification.
4. The tarr
white crystal
methanol, ar
screening test
and GC/MS
5. All org;
repeated extr
tarry substar
nitrogen in tl
6. Furfura
acidic mediu
media. The
quantitation

Discoloration
Procaine,
most commo
have been fi
derivatives a
agent [1,2].
Studies of
of cocaine, 1
of each of tl
bags on a
browning w.
Experime
which unde
of time, we
effect in a s
Table 1
amines whi
effect of
lactose-ami

Discussion
The resu
long perio
may occur
as 5-hydro
lactose due
cocaine, ar

color and microcrystalline tests [6,7]. Additional confirmation was obtained on a gas chromatograph/mass spectrometer (GC/MS).

2. The tarry substance was acidic and had the odor of acetic acid. The presence of acetate ions was determined by FeCl_3 test and ethyl acetate tests [8,9].

3. 5-hydroxymethyl furfural and related compounds are the degradation products of lactose, especially when the spray-dried process lactose is used [1,2]. Identification of these furfural compounds was made from the benzene extraction of the tarry substance in acid media. The McCance method [10] and Sanchez reaction [7] were used for this identification.

4. The tarry substance dissolved in methanol giving a brown solution and leaving a white crystalline solid. The solid was filtered, purified by further washings with methanol, and identified as lactose. Lactose identification was made by primary screening test for reducing sugars with benedict's solution, followed by X-ray diffraction and GC/MS identification.

5. All organic nitrogeneous compounds were isolated from the tarry substance by repeated extractions with water and chloroform in acid and alkaline media. The residual tarry substance was submitted to the classical sodium fusion test, and presence of nitrogen in that tarry substance was detected.

6. Furfural derivatives were isolated by benzene extraction of the tarry substance in an acidic medium. The cocaine was then extracted with benzene in sodium bicarbonate media. The extracted cocaine was passed through an acid-washed celite column and quantitation was done using a Perkin-Elmer 350 spectrophotometer [11,12].

Discoloration Due to Diluents

Procaine, benzocaine, tetracaine, lidocaine, quinine, and sugars such as lactose are the most common diluents mixed with cocaine. On occasion other chemicals, such as aspirin, have been found with the other diluents. A 5-hydroxymethyl furfural and other furfural derivatives are also present as contaminants when spray-dried lactose is used as a cutting agent [1,2].

Studies on discoloration were conducted using USP-grade lactose and hydrochlorides of cocaine, procaine, benzocaine, tetracaine, lidocaine, and quinine. The thermal stability of each of the compounds was studied by heating them up to four hours in sealed plastic bags on a steam bath. Each was found thermally stable and no color change nor browning was observed.

Experiments on the discoloration effects of various combinations of these compounds, which under normal storage conditions at ambient temperatures would take long periods of time, were conducted. To increase the rate of reaction and to study the discoloration effect in a short period of time the samples were heated on a steam bath.

Table 1 shows the discoloration effect due to the interaction of lactose with those amines which are commonly used as cutting agents with cocaine. Tables 2 to 4 show the effect of other contaminants, such as aspirin, furfural, and moisture, on the lactose-amine reaction resulting in the brown tar formation.

Discussion and Results

The results of this study indicate that lactose, if mixed with cocaine and stored for a long period of time, will not react and give discoloration. However, some discoloration may occur due to the condensation reactions of the degradation products of lactose, such as 5-hydroxymethyl furfural and other furfural derivatives which are associated with lactose due to its spray-drying process [1,2]. If other amines, commonly used for cutting cocaine, are present in addition to lactose, browning on storage is likely to occur due to

TABLE 1—*Browning due to lactose-amine interaction. Heating on steam bath was done for 3 h in heat-sealed plastic bags.*

Reactants	Ratio ^a	Discoloration
Lactose + procaine	7:1	gradual discoloration finally resulting in dark brown color in 3 h
Lactose + benzocaine	7:1	very slow and slight discoloration observed in 3 h
Lactose + tetracaine	7:1	no discoloration observed
Lactose + lidocaine	7:1	very slow and slight discoloration observed in 3 h
Lactose + quinine	7:1	no discoloration observed
Lactose + cocaine	7:1	no discoloration observed

^a Ratio of 7:1 was chosen since it was the most realistic ratio pertaining to the evidence in question.

TABLE 2—*Effect of furfural on browning. Heating on steam bath was done for 3 h in heat-sealed plastic bags.*

Reactants	Ratio	Discoloration
Lactose + furfural	7:1	gradual and slight browning only
Lactose + procaine + furfural	7:1:1	immediate red coloration which turns to a brown tarry substance
Lactose + benzocaine + furfural	7:1:1	immediate red coloration which turns to a brown tarry substance
Lactose + tetracaine + furfural	7:1:1	gradually turns brown
Lactose + lidocaine + furfural	7:1:1	gradually turns brown
Lactose + cocaine + furfural	7:1:1	gradual and slight browning only

TABLE 3—*Effect of acetate ions on browning. Heating on steam bath was done for 3 h in heat-sealed plastic bags.*

Reactants	Ratio	Discoloration
Lactose + aspirin + CH ₃ COOH	7:1:1	no color change
Lactose + procaine + CH ₃ COOH + aspirin	7:1:1:1	starts browning almost immediately, 35 min to a dark tarry substance, no procaine detected
Lactose + procaine + aspirin + 1 drop of water	7:1:1	starts browning almost immediately, 35 min to a dark tarry substance, no procaine detected after the physical change ^a
Lactose + procaine + CH ₃ COOH	7:1:1	starts browning almost immediately, 35 min to a dark tarry substance, no procaine detected after the physical change
Lactose + cocaine + aspirin	7:1:1	no change
Lactose + benzocaine + aspirin	7:1:1	browning starts slow, benzocaine still present
Lactose + tetracaine + aspirin	7:1:1	slight browning

^a An acetic-acid-like odor was detected.

TABLE 4—Effect of moisture on lactose-amine reaction. Heating on steam bath was done for 3 h in heat-sealed plastic bags.

Reactants	Ratio	Discoloration
Cocaine + lactose + aspirin (ovendried)	7:1:1	no discoloration
Cocaine + lactose + aspirin + moisture	7:1:1	no discoloration
Procaine + lactose + aspirin (ovendried)	7:1:1	slow discoloration towards browning
Procaine + lactose + aspirin + moisture	7:1:1	starts turning brown within 5 min and completes tar formation in 30 min, no procaine found present after tar formation ^a

^a An acetic-acid-like odor was detected.

lactose-amine reaction, commonly known as the Milliard reaction [1].² The extent of browning would depend on the reactivity of amines with lactose (Table 1). Procaine, (and other primary amines), is very reactive because it requires a low order of initiation energy and exhibits autocatalytic qualities once the reaction has begun [1]. The furfurals associated with lactose will also react with the amines, as can be seen from Table 2. It has been further reported that certain ions such as acetates, stearates, etc catalyze the Milliard reaction [2]. Table 3 indicates the effect of acetate ions on the Milliard reaction. Aspirin brings about the same catalytic effect as that of acetates, since on hydrolysis it releases acetic acid. Moisture has been found to be an equally important factor in bringing about the Milliard reaction, as can be seen from Table 4.

The evidence in question, after transformation, was found to contain lactose, acetate ions, cocaine, furfural derivatives, and moisture. Quantitation of the evidence indicated that the amount of cocaine found in the tarry substance is in agreement with the amount of cocaine reported in the original analysis and, as previously stated, known mixtures of cocaine and lactose do not result in discoloration (Table 1). Therefore, it can be concluded that cocaine was not involved in the formation of the tarry substance. A sodium fusion test on the isolated tarry substance indicated the presence of nitrogen. In view of the presence of nitrogen in the isolated tarry substance, the experimentation with the common diluents (Table 3), and the similarity in the physical appearance between the evidence and the known lactose-procaine mixture, it was deduced that the most probable reaction which resulted in the tar formation was that of procaine with lactose.

The original evidence, which was analyzed in July 1972, was a tan powder which turned brown by February 1973 during storage. This indicates the initiation of the lactose-procaine browning reaction. The tar formation, which occurred in a relatively short period of time between February and June 1973, could be explained by the catalytic effect of acetates resulting from the decomposition of aspirin during storage.

The hydrolysis of procaine in dilute acetic acid did not occur. However, the effects of moisture and acetate ions on the mechanism of browning due to the lactose-procaine reaction, and the isolation and identification of the end products formed, are presently under study and will be presented in a subsequent communication.

Summary

Discoloration in contraband drugs may occur on storage. The extent of this discoloration depends on the time of storage and the type of cutting agent present with

²The Milliard reaction is a Schiff base reaction.

the drug. In the investigation of the evidence in question the main cause of the brown tarlike transformation of the tan powder was attributed to the procaine-lactose reaction, which was catalyzed by the acetic acid due to the hydrolysis of the aspirin present in the contraband seizure. The self-condensation of the degradation products of lactose, such as 5-hydroxymethyl furfural, can also contribute to the discoloration.

Comments

The explanation for the brown tarry transformation of the confiscated cocaine powder, as a result of the experimental work detailed herein, was accepted by the court. The confiscated cocaine was then admitted as evidence and subsequent convictions were obtained.

Acknowledgments

Our thanks go to Captain Charles V. Rorke and Lieutenant Patrick McCarthy for their encouragement during this investigation. Our special thanks go to Chemist James Fava and members of the Chemistry Section for their cooperation and help in this investigation.

References

- [1] Brownley, C. A. and Lachman, Leon, "Browning of Sprayed Processed Lactose," *Journal of Pharmaceutical Sciences*, Vol. 53, No. 4, April 1964, pp. 452-454.
- [2] Koshy, K. T., Duvall, R. N., Troup, A. E., and Pyles, J. W., "Factors Involved in Browning of Spray Dried Lactose," *Journal of Pharmaceutical Sciences*, Vol. 54, April 1965, pp. 549-554.
- [3] Costello, R. A. and Mattacks, A. M., "Discoloration of Tablets Containing Amines and Lactose," *Journal of Pharmaceutical Sciences*, Vol. 51, Feb. 1962, pp. 106-108.
- [4] Duvall, R. N., Koshy, K. T., and Pyles, J. W., "Comparison of Reactivity of Amphetamine, Methamphetamine and Dimethylamphetamine with Lactose and Related Compounds," *Journal of Pharmaceutical Sciences*, Vol. 54, April 1965, pp. 607-611.
- [5] Blaug, S. M. and Huang, Wen Tzy, "Interaction of Dextroamphetamine Sulphate with Spray Dried Lactose," *Journal of Pharmaceutical Sciences*, Vol. 61, Nov. 1972, pp. 1770-1775.
- [6] Clark, E. C. G., *Screening and Identification of Drugs*, The Pharmaceutical Press, London, 1969.
- [7] "Methods of Analysis," Internal Revenue Service, Washington, D.C., 1967, pp. 78 and 124.
- [8] Feigl, F., *Spot Test in Organic Analysis*, American Elsevier, New York, 1966, p. 456.
- [9] *Remington's Pharmaceutical Sciences*, Mack Publishing Co., Easton, Pa., 1965, p. 1438.
- [10] Brownley, C. A. and Zerban, F. W., *Physical and Chemical Methods of Sugar Analysis*, Wiley, New York, 1948, p. 915.
- [11] Canaff, R. F., "Ion Pairing Chromatographic Separation and Determination of Cocaine," *Microgram*, Vol. 3, No. 4, June 1970, p. 121.
- [12] Moore, J., "The Analysis of Cocaine in Illicit Preparations," *Microgram*, Vol. 3, No. 3, May 1970, pp. 89-96.

New York City Police Department Crime Laboratory
235 E. 20th St.
New York, N.Y. 10003

A
lc
E
—
cap
pre
cou
knc
(DI
Tab
T
clan
Adn
are
auth
ager
cont
these
in g
cann
previ
Ta
inclu
vary
sugar
numl
addit
may l
The
betwe
machi
minut
with l

Rece
For
McLea
Chi
nology,

CHROM. 10,706

Notice: This material may be protected
by copyright law (Title 17 U. S. Code)

Note

Thin-layer chromatographic procedure for the differentiation of the optical isomers of cocaine

DERK ESKES

Forensic Science Laboratory, Municipal Police, Elandsgracht 117, P.O. Box 2287, Amsterdam
(The Netherlands)

(Received November 2nd, 1977)

Medicinally used cocaine is prepared from the leaves of the coca plant (*Erythroxylon coca*, Lamarek) by extracting the total alkaloids and converting them to *l*-ecgonine by acid hydrolysis. The isolated *l*-ecgonine is esterified with methanol and benzoic acid to produce the natural product, *l*-cocaine¹. In 1923 a method for the synthesis of *dl*-cocaine was published². The racemic mixture was separated into the *dextro* and *levo* isomers of cocaine with tartaric acid³. Probably all of the currently available illicit cocaine is prepared from coca leaves and therefore is the natural isomer *l*-cocaine. However, there is the possibility of the occurrence of synthetic *dl*- or *d*-cocaine in illegal drugs and, for legal reasons, it may be necessary for the forensic chemist to be able to distinguish between the optical isomers. The following method was developed for this purpose and makes possible the routine analysis of small samples of illicit cocaine and the determination of the optical isomer or isomers present in the sample. Cocaine was hydrolysed to ecgonine and then esterified with the enantiomeric 2-octanols to give the necessary derivatives. The resulting diastereoisomers may be distinguished by thin-layer chromatography (TLC).

EXPERIMENTAL

Materials and methods

The *dl*-cocaine was prepared by benzylation of *dl*-methylecgonine² and the *d*-cocaine by resolution of *dl*-cocaine². The *l*-cocaine and the benzenesulfonyl chloride were obtained commercially (E. Merck, Darmstadt, G.F.R.) as well as the *d*-2-octanol, *l*-2-octanol and *dl*-2-octanol (Fluka, Buchs, Switzerland).

The diastereomeric derivatives were prepared in 0.3 ml Reacti-Vials[®] provided with magnetic stirrers (Pierce Chem. Co., Rockford, Ill., U.S.A.) according to the procedure of Brewster and Ciotti³.

TLC separation of the isomers was carried out on non-activated pre-coated silica gel plates of thickness 0.25 mm (E. Merck) with methanol as the developing solvent.

Procedure

Approximately 0.5 mg of cocaine base or hydrochloride was placed in a 0.3 ml

screw-cap vial and 30 μ l of 2 *N* hydrochloric acid were added. The vial was capped and heated for 30 min in a heating block maintained at 120°. The cap was removed and the heating continued for about 20 min to evaporate the solvent completely. The residue was allowed to cool to room temperature, then 30 μ l of pyridine were added, followed by 3 μ l of *d*-2-octanol and the vial was placed in a small dish with water and ice. After cooling for some minutes 1 μ l of benzenesulfonyl chloride was added. The vial, provided with a magnetic stirrer, was capped and stirred slowly (60 rpm) for 30 min while the solution was kept cold. Then 30 μ l of 2 *N* ammonia were added, followed by 100 μ l of ether and the mixture was stirred rapidly for a short time.

Then 2 μ l of the ether layer were placed on the starting line of the TLC plate. The plates were developed in cylindrical glass vessels (without prior equilibration) with methanol until the solvent front had reached a height of 15 cm. The plates were air dried at room temperature and sprayed with acidified iodoplatinate⁴. The two reference spots were obtained by using the same procedure with *l*-cocaine and *d*-2-octanol and with *l*-cocaine and *l*-2-octanol or, more simply, with *l*-cocaine and *dl*-2-octanol.

RESULTS AND DISCUSSION

The method for the preparation of esters with benzenesulfonyl chloride in pyridine¹ was found to be the most convenient, giving a high yield in a short time. The optically active alcohols, menthol, borneol, 2-methyl-1-butanol and 2-octanol, were used to prepare esters of *dl*-ecgonine in this way. The TLC separation was best when the diastereoisomers prepared from 2-octanol were used. The R_f values of the diastereoisomers of 2-octanyl-ecgonine are shown in Table I. After spraying with acidified iodoplatinate⁴ the spots were intense blue-green, discrete and well defined. In case of double spots of equal amounts of the diastereoisomers, the spot with the high R_f value was always distinctly smaller than the other spot.

TABLE I
SEPARATION OF DIASTEREOMERS OF 2-OCTANYL-ECGONINE BY THIN-LAYER CHROMATOGRAPHY

Ester prepared from	R_f
<i>l</i> -Ecgonine and <i>d</i> -2-octanol	0.36
<i>l</i> -Ecgonine and <i>l</i> -2-octanol	0.28
<i>l</i> -Ecgonine and <i>dl</i> -2-octanol	0.28, 0.36
<i>d</i> -Ecgonine and <i>d</i> -2-octanol	0.28
<i>d</i> -Ecgonine and <i>l</i> -2-octanol	0.36
<i>d</i> -Ecgonine and <i>dl</i> -2-octanol	0.28, 0.36
<i>dl</i> -Ecgonine and <i>d</i> -2-octanol	0.28, 0.36
<i>dl</i> -Ecgonine and <i>l</i> -2-octanol	0.28, 0.36

ACKNOWLEDGEMENTS

I thank Dr. L. Maat and Mr. T. S. Lie (Laboratory of Organic Chemistry, Technische Hogeschool, Delft, The Netherlands) for the gift of *dl*-methylecgonine.

REFERENCES

- 1 V. E. Tyler, L. R. Br... ed., 1976, pp. 243-2...
- 2 R. Willstätter, O. W...
- 3 J. H. Brewster and...
- 4 E. G. C. Clarke, *Isol...*

NOTES

was capped
 removed
 tely. The
 ere added,
 water and
 added. The
 (m) for 30
 added, fol-
 time.
 .LC plate.
 (vibration)
 plates were
 . The two
 e and *d*-2-
 and *dl*-2-

floride in
 fort time.
 octanol,
 y was best
 y of the
 ng with
 ll defined.
 t with the

UN-LAYER

REFERENCES

- 1 V. E. Tyler, L. R. Brady and J. E. Robbins, *Pharmacognosy*, Lea and Febiger, Philadelphia, 7th ed., 1976, pp. 243-246.
- 2 R. Willstätter, O. Wolfes and H. Mäder, *Ann.*, 434 (1923) 111.
- 3 J. H. Brewster and C. J. Clotti, Jr., *J. Amer. Chem. Soc.*, 77 (1955) 6214.
- 4 E. G. C. Clarke, *Isolation and Identification of Drugs*, Vol. 1, Pharm. Press, London, 1969, p. 801.

Property of Idaho State Police Forensic Services
 Uncontrolled Internet Copy
 OBSOLETE DOCUMENT

Chemistry,
 gonine.

TLC METHOD TO SEPARATE LSD
FROM LYSERGIC ACID METHYL PROPYL AMIDE

Ann R. Bradley
Forensic Laboratory, Idaho Department of Health and Welfare
Boise, Idaho

OBJECTIVE

Because IR identification of LSD in small amounts can be difficult, a method for TLC identification, including separation from lysergic acid methyl propyl amide (LAMPA) is desired.

BACKGROUND

The separation of LSD from LAMPA by thin layer chromatography has been the subject of recent investigation (Microgram VI, No. 6, and J. Association of Offic. Anal. Chem. Volume 56, (1973), 88). This laboratory has tried both methods suggested by these authors: (1) Formation of trimethylsilyl derivatives of LSD and LAMPA, followed by separation on Merck Silica Gel F-254 plates using acetone as the solvent; (2) Impregnating pre-coated silica gel G plates with ethyl acetate: N,N-dimethylformamide: ethanol (13; 1.9; 0.1), letting them dry, then spotting samples and developing for 16 cm in the same solvent. The results were not satisfactory in our hands.

Using the former method, we could not form the TMS derivative using the conditions given and in the latter case, the separation achieved was not adequate. This may have been due to the difference between the Brinkmann plates used by those authors and the Merck plates available to us.

An alternate TLC method is proposed.

APPARATUS

Chromatography tank lined with filter paper for optimum saturation.

Thin layer plates; (5 x 20 cm, precoated 0.25 mm, Merck Silica Gel 60 F-254). These plates are scored and broken to 5 x 10 cm prior to use, then oven dried.

PROCEDURE

Extract the sample and standards into CH_2Cl_2 or similar organic solvent before spotting. Do not use methanol, as the spots must be kept small. Apply spots of LSD standard, mixed LSD standard and LAMPA standard, the sample extract and a mixed spot of sample + LAMPA standard. (LAMPA is available from Applied Science, Cat. No. 01810.)

Prepare the following solvent; chloroform freshly saturated with concentrated NH_4OH (1 part) + acetone (2 parts). Moisten the paper lining of the tank with the solvent and develop the plate to its full height, 10 cms. This takes about 20 minutes. Remove the plate, dry it briefly, examine the spots under UV light. If desired, repeat the development one more time and again dry the plate.

PROCEDURE (continued)

Locate the spots using UV light and spray with p-DMAB spray to produce purple color.

RESULTS

The first pass moves the LSD spot ahead of the LAMPA spot by about 3 mm (measured between the two centers). The second pass improves the separation to about 4 mm, with a total distance traveled of approximately 55-60 mm. If iso-LSD is present, it runs slower than LAMPA and does not interfere.

DISCUSSION

1. The chloroform component of the solvent should be freshly saturated with ammonia and the solvent fresh in the tank. The same solvent may be used for the second development, if one is needed, but then it should be discarded.
2. As with any TLC procedure, care should be taken to not overload the plate. Total amounts should be kept as small as possible, consistent with visualization using p-DMAB spray.
3. Often the separation appears better when the plate is viewed from the back side.
4. Tank geometry seems to have an influence on the quality of separation. Best results were obtained with these short plates in rather short, round jars, in which they fit at an angle.
5. Each chemist has his/her own standard for what constitutes separation between two close spots. Some may find one development sufficient; others may decide two passes are needed.
6. Merck improved their silica gel plates in 1971-72. This separation cannot be achieved using the older type plates. The improved ones carry the designation, Silica Gel 60 F-254.
7. This laboratory does not rely on one TLC system alone for the identification of LSD. The other systems in use are 1,1,1 Trichloroethane/Methanol (9/1) and Benzene/Dimethylformamide (13/2), which were recommended by James A. Heagy, San Francisco Regional DEA Laboratory.

PROC. Manual

LSD (Lysergic acid diethylamide)

Analytical

Sample preparation for TLC:

(Applies to blotter acid, tablets, window panes)

Short Method: One hit either ground up or cut into small pieces; add a small amount of methanol or CH_2Cl_2 and shake, spot the solvent.

Long Method: Cut or grind material; add 0.1N H_2SO_4 - shake for 15 minutes; make slightly alkaline with Na_2CO_3 , extract with ether; concentrate in volume in hot water bath under N_2 . Spot this extract.

TLC Systems: (Use three for confirmatory identification)

1. 1,1,1 Trichloroethane/methanol (9/1)
2. Benzene/Dimethylformamide (13/20)
3. CHCl_3 saturated with NH_4OH /acetone (1/2) See attached paper by Bradley.

It is advisable to spot lysergic acid methyl propyl amide (LAMPA) in addition to LSD.

Visualize with p-DMAB spray.

Attached is an IR preparation method.

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

LSD - Micro-Pelleting for IR

1. Using a centrifuge tube, place sample in 2.0 ml. $0.1\text{N H}_2\text{SO}_4$ (for 3 hits or less)
2. Vortex 1 min.
3. Make basic with concentration ^{of} NH_4OH .
4. Add 1 ml. CHCl_3 - spectro or nano grade.
5. Vortex 1 min.
6. Centrifuge.
7. Remove organic layer to 5.0 ml centrifuge.
8. Repeat steps 4 through 7 3X.
9. Reduce in volume to 1/4-1/2 ml. under N_2 .
10. Spot on TLC plate.
11. Rinse centrifuge tube with 1/4 ml CHCl_3 and vortex.
12. Spot same area.
13. System for TLC - CHCl_3 9/MeOH, 1 (Same grade CHCl_3).
14. Scrape off LSD band from plate into mortar and grind with pestle.
15. Use disposable pipet, plug with glass wool.
16. Test on vacuum - see if plug holds.
17. Vacuum silica gel into pipet.
18. Elute LSD with $0.1\text{N H}_2\text{SO}_4$ in 3-4 ml. aliquots. Pipet plugs up; may be necessary to stir silica gel with long pipet and force $0.1\text{N H}_2\text{SO}_4$ through. Catch in 15 ml centrifuge tube.
19. Add NH_4OH til basic.
20. Extract with 0.5 ml CHCl_3 .
21. Vortex.
22. Centrifuge.
23. Transfer organic to 5.0 ml centrifuge tube.
24. Repeat 20-23 3X.
25. Reduce to dryness under N_2 .
26. Place dry tube in vacuum desiccator for 15 min.

27. Remove and add 1/2 ml. CHCl_3 .
28. Drop on 5 mg KBr - mortar & pestle - let dry.
29. Rinse ^{tube} with 0.5 ml. CHCl_3 .
30. Vortex.
31. Drop on KBr.
32. Place mortar in vac. des. for 15 min.

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

PHENCYCLIDINEBACKGROUND

- late 1950's developed for use as an anesthetic, but too many adverse side effects resulted in it being discontinued for human use in 1967; still available as a veterinary immobilizing agent under the brand name "SERNYL" (Parke, Davis, and Co.)
- first reported illicit use was in mid-1960's on the West Coast; it was known as the "PeaCe Pill" or "angel dust" which was sprinkled on parsley and smoked.
- became a Schedule III drug on 4/6/69, then placed in Schedule II 2/24/78
- TCP placed in Schedule I in 1975; PCE (cyclohexamine) and PHP placed in Schedule I 10/25/78; PCC controlled not as a drug itself, but as an immediate precursor/intermediate in the manufacture of PCP. (C-II)

MANUFACTURE

- ① -PCP most commonly produced through nitrile intermediate, PCC, which is formed in high yield (90-98%) and is frequently found in the final product; the entire process for producing PCP takes only about 4 hours.
- ① -NaCN or KCN is used in this process and is converted to HCN in the presence of acid; ventilation is important if one is involved in any clandestine lab seizures.
- ① -the other process often leaves piperidine as an impurity which can often be recognized by its odor.

PHYSICAL PROPERTIES

- ① -refer to list of structures for PCP and PCP analogs.
- ⑫ PCP
 - a) hydrochloride salt $C_{17}H_{25}N-HCl$ MW 279.86 mp 230-231 C
soluble in MeOH, EtOH, $CHCl_3$, CH_2Cl_2 , H_2O
sparingly soluble in dilute HCl
insoluble in hexane, ether
 - b) free base $C_{17}H_{25}N$ MW 243.41 mp 46-46.5 C
soluble in $CHCl_3$,
insoluble in hexane, ether
- ⑧, ⑭ TCP
 - a) hydrochloride salt $C_{15}H_{23}NS-HCl$ MW 285.85 mp 233-236 C
soluble in MeOH, EtOH, $CHCl_3$, H_2O
insoluble in hexane, ether
 - b) free base $C_{15}H_{23}NS$ MW 249.40
soluble in $CHCl_3$, hexane, ether

SPOT TESTS PCP

- ⑧ Marquis---no rxn., effervesces; rule out carbonates by addition of dilute HCl
- ⑧ Meckes---no rxn., effervesces
Froehdes---no rxn., effervesces
- ⑧ Mandelins---no rxn., effervesces
Wagners---pos. rxn. for alkaloids
Mayers---white ppt.
Cobalt thiocyanate---slow, faint blue ppt.
Ruybals---blue ppt.
- ② Scotts---PCP is negative, but can get false positives when PCP in combination with some drugs.
- ③ Zelonis---clear blue benzene layer; works well on plant material, but best to extract material first with CHCl_3 or CH_2Cl_2 , and take to dryness before adding reagents.
Sodium nitroprusside---some samples give pos. rxns. for secondary amines.

TCP

- ⑧, ⑭, ⑰ Marquis---effervesces, gray-orange, then slowly to gray-green
- Meckes---effervesces, orange or yellow/green, then to yellow-blue/green, then to deep blue
- Mandelins---green, effervesces

Zelonis - ⊕

KETAMINE

⑤

- Marquis---no rxn.
- Cobalt thiocyanate---blue ppt.
- Wagners---pos. rxn. for alkaloids
- Liebermans---slow purple

Zelonis ⊖

Ruybal's ⊖

2° any amine ⊖!!

Cyclohexamine

2° any amine ⊖!!

Zelonis ⊖
Ruybal's ⊕

TLC

- (4) 1) Butyl ether:diethyl ether:diethylamine (45:45:10)
- (4) 2) Benzene:diethylamine (95:5) *May be useful for Cocaine + PCP/Amo.*
- (8, 9) 3) Chloroform:methanol (9:1) - Rf's vary from paper to paper
- (5) 4) Clarke's T₁
- (5) 5) Hexane:benzene:diethylamine (75:25:10)
- (5) 6) Chloroform:dioxane:ethyl acetate:ammonium hydroxide (25:60:10:5)
- (6) 7) Benzene:acetone:pyridine (16:8:1)
- (7) 8) Davidows
- (7) 9) Acetone:chloroform (1:1)
- (7) 10) Acetone:chloroform (65:35)
- 11) Ethyl acetate:methanol:ammonium hydroxide:water (29:1:0.25:0.5)
- 12) Same as (11) with ammonium hydroxide removed

	(1)	(2)		(4)	(5)	(6)
PCP	0.91	0.64	Ketamine	0.72	0.55	0.86
THC	0.61	0.41	PCP	0.59	0.90	0.92
CBD	0.69	0.52				
CBN	0.47	0.30				
DET	0.30	0.15				
DMT	0.18	0.09				

	(7)	(8)	(9)	(10)	Rf's for (3)			
PCP	0.24	0.86	0.23	0.32	0.17	0.43	0.37	0.57
TCP	0.49	0.82	0.45	0.49	0.33	0.67	0.70	
PCC	0.70	--	--	--	0.82	0.91		
PHP	--	0.75	0.05	0.13	0.10	0.43		
PPP	--	0.80	0.29	0.40	0.25			
Piperidine	0.00	--	--	--	--			

11) PCP = 0.90 by eliminating ammonium hydroxide in (12) PCP can be
 12) PCP = 0.20 separated from more neutral drugs with same Rf as (11)

GC

- ⑦ 1) 3% OV-17 on Gas Chrom Q 100/120 DMCS, 6' x 1/4" x 2mm (id)
column = 170 C, inj. port = 250 C, detector = 270 C
Nitrogen flow rate = 30 ml/min
- ⑧ 2) 3% OV-101, all other parameters same as (1) *Better System for analysis*
- ⑩ 3) 3% OV-17 on Gas Chrom Q 100/120, 6' x 1/8" (id)
column = 90 C(3 min hold) to 250 C(2 min hold) @ 10 C/min
injector = 275 C, detector = 300 C, nitrogen = 60 ml/min
- ⑫ 4) 3% OV-17, 6' column @ 170 C, no other parameters given
- ⑫ 5) 3% OV-1, 6' column @ 170 C, no other parameters given

(min.)	(1)	(2)	(3)	(4)	(5)
PCP	6.52	6.70	0.96	--	--
TCP	6.47	6.48	9.62	0.63	0.63
PHP	4.83	4.94	10.49	0.78	0.73
PPP	4.06	4.36	16.48	4.40	3.83

- ⑪ -PCH (phenylcyclohexene) is a degradation product of PCP by the elimination of piperidine; caused by high injection port temps. mainly, but some breakdown does occur on-column; degree of degradation is proportional to the amount of sugars present in the sample and also on the salt form of PCP; the HCl salt breaks down more than the free base and is not consistent in the amount that degrades from one injection to the next; the free base is more stable and has a more consistent response.

Property of Idaho State Police Forensic Services
Unconfidential Internet Copy
OBSOLETE DOCUMENT

UV

	0.1N HCl (nm) -- in order of size			
PCP	262	257	268	252
(9, 13) PHP	262	257	268	252
(13) PPP	261	256.5	268	252
(9, 14) TCP	232			
(11) PCE	262	257	268	236
Ketamine	270	277	263	

- (15) -the free base form of PCP has almost no absorbance in the UV region, so it should be converted to the HCl salt before running the scan.
- (18) -PCP/LSD mixture -- expose solution to longwave UV to decrease the absorbance of LSD so that PCP curve may be seen.

MICROCRYSTAL TESTS

- (19)
- 1) $\text{HAuBr}_4 \cdot \text{HOAc} \cdot (2+3)\text{H}_2\text{SO}_4$ cg
 - 2) $\text{HAuBr}_4 \cdot 2\text{H}_3\text{PO}_4 \cdot (2+3)\text{H}_2\text{SO}_4$ cg
 - 3) HOAc (1 drop) + HAuCl_4 (5% aqueous)

-(1) -- worked best for PCP/marihuana sample; chloroform extract evaporated on the slide with a 50ul pipette.

-(3) -- also worked well, but (2) did not

Other crystal tests

- (20) 10% KI (aqueous) , 10% NH_4SCN (aqueous) -- for PCP & TCP
- (7) 10% HCl or 10% HOAc + 2% KMnO_4 -- for PCP, PHP, PPP
- (21) 5% aqueous HAuCl_4 (volatility) -- for PCP and piperidine

CLEAN-UP

1) Direct extraction using CHCl_3 , CH_2Cl_2 , MeOH, EtOH, or H_2O and recrystallize with pet ether

(22) 2) Ion-pairing with CHCl_3 from dilute HCl

3) 0.1N ^{HCl} celite column - *elute w/ H_2O washed CHCl_3*

(23) 4) Alumina column

5) Prep TLC

6) Acid-base shakeout

7) For PCP & TCP refer to physical properties; free bases of each have different solubilities in hexane & ether

8) Cobalt thiocyanate derivative (ask Herbie)


-(2) useful for separating PCP from many other basic drugs; (6) can be used the same way with a few washes to carry the clean-up one step further.

-(3) & (4) good for cleaning up PCP from plant material; use water-washed solvents on celite; CH_2Cl_2 worked well for alumina, but you take your chances with each new lot no.



-(5) if you've fouled up everything else

Property of Idaho State Police Forensic Services
Uncontrolled Internal Copy
OBSOLETE DOCUMENT

PCP THROUGH NITRILE INTERMEDIATE




NaHSO₃



CYCLOHEXANONE


1-PIPERIDINOCYLO-
HEXANECARBONITRILE
(PCC)



Mg
Et₂O



BROMOBENZENE



PCP

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT


PCP THROUGH ENAMINE INTERMEDIATE




CYCLOHEXANONE

PIPERIDINE


BENZENE



PCP




1-(4-CYCLOHEXYNYL)-
PIPERIDINE




Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

IMPURITIES IN ILLICIT PCP SAMPLES




**1-PIPERIDINOCYLO-
HEXANECARBONITRILE**

PCC




**1-(1-CYCLOHEXENYL)-
PIPERIDINE**

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT



1-PHENYLCYCLOHEXENE


GC degradation product




BIPHENYL

10


PCP ANALOGS




KETAMINE




THIOPHENE, TCP
(schedule I)



CYCLOPENTYL
PPP



PYRROLIDINE, PHP
PCPy









N-ETHYL, PCE
CYCLOHEXAMINE

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

Phencyclidine and Some Analogues and Variants Thereof

1. 1-phenylcyclohexylamine
2. (1-phenylcyclohexyl)methylamine
3. (1-phenylcyclohexyl)dimethylamine
4. (1-phenylcyclohexyl)methylethylamine
5. Cyclohexamine - PCE - C-1
(1-phenylcyclohexyl)ethylamine
6. (1-phenylcyclohexyl)diethylamine
7. (1-phenylcyclohexyl)isopropylamine
8. 1-(1-phenylcyclohexyl)pyrrolidine - PHP - PCPy - C-1
9. Phencyclidine
1-(1-phenylcyclohexyl)piperidine
10. 1-(1-phenylcyclohexyl)morpholine
11. 2-[1-(1-amino)-cyclohexyl]thiophene
12. 2-[1-(1-dimethylamino)cyclohexyl]thiophene
13. 2-[1-(1-ethylmethylamino)cyclohexyl]thiophene
14. 2-[1-(1-ethylamino)cyclohexyl]thiophene
15. 2-[1-(1-diethylamino)cyclohexyl]thiophene
16. 2-[1-(1-isopropylamino)cyclohexyl]thiophene
17. 1-[1-(2-thienyl)cyclohexyl]pyrrolidine
2-[1-(1-pyrrolidino)cyclohexyl]thiophene
18. 2-[1-(1-β-methylpiperidino)cyclohexyl]thiophene
19. Thiophene Analog of Phencyclidine
1-[1-(2-thienyl)cyclohexyl]piperidine
20. 1-[1-(2-thienyl)cyclohexyl]morpholine
21. 2-[1-(1-β-methyl-pyrrolidino)-cyclohexyl]-
4-methyl thiophene

22. 2-[1-(1-piperidino)cyclohexyl]-4-methyl thiophene
23. 2-[1-(1-piperidino)cyclohexyl]-5-methyl thiophene
24. 2-[1-(1-dimethylamino)cyclohexyl]furan
25. 2-[1-(1-piperidino)cyclohexyl]furan
26. 2-[1-(1-piperidino)cyclohexyl]-5-methyl furan
27. 2-methylamino-2-phenylcyclohexanone
28. 2-ethylamino-2-phenylcyclohexanone
29. Ketamine
2-(o-chlorophenyl)-2-(methylamino)-cyclohexanone
30. 2-methylamino-2-(o-hydroxyphenyl)-cyclohexanone
31. 2-methylamino-2-(m-hydroxyphenyl)-cyclohexanone
32. 2-amino-2-(2-thienyl)cyclohexanone
33. Tiletamine
2-ethylamino-2-(2-thienyl)cyclohexanone



1. -NH₂

2. -NHCH₃


3. -N(CH₃)₂

4. -N(CH₃)(C₂H₅)


5. -NHC₂H₅

6. -N(C₂H₅)₂

7. -NHCH(CH₃)₂



17.



18.




21.




22.



9.




19.




22.




23.



25.




26.




10.

24. -N(CH₃)₂


Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT




32. -NH₂




31. -NHCH₃



30. -NHCH₃



29. -NHCH₃



27. -NHCH₃

28. -NHC₂H₅

33. -NHC₂H₅

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

SPECIFICITY PROBLEM WITH THE COCAINE-SPECIFIC FIELD TEST
II. NON-PHENOTHIAZINE FALSE POSITIVES AND THE
SEPARATION OF PHENCYCLIDINE - PROMAZINE COMBINATIONS

Steven K. Lorch, Ph.D.
Crime Laboratory Scientist
Division of Crime Detection
Bureau of Laboratories
Michigan Department of Public Health
Lansing, Michigan 48914

OBJECTIVE

We have found that certain non-phenothiazine drug combinations will give the "highly specific" field test (1) for cocaine, and that the FPN test (2) will not screen for these false positives.

BACKGROUND

Since our previous report (3) our laboratory has received three more cases of phencyclidine (PCP) - phenothiazine combinations which gave the cocaine-specific field test color sequence and a positive FPN test for phenothiazines. The phenothiazine was identified as promazine by gas chromatography and gas chromatography-mass spectrometry. As reported earlier promazine gives steps one and two of the cocaine-specific field test (3).

The test was run on standard samples of the following drugs and combinations of drugs:

TABLE 1

	<u>Color after step 1</u>	<u>Color after step 2</u>	<u>Color in CHCL₃</u>	<u>FPN</u>
1) Phencyclidine	blue	pink	clear	-
2) Promazine	blue	pink	clear	+
3) Cocaine	blue	pink	blue	-
4) 1 + 2	blue	pink	blue	+
5) Dibucaine	blue	pink	clear	-
6) Methapyrilene	blue	pink	clear	-
7) 1 + 5	blue	pink	blue	-
8) 1 + 6	blue	pink	blue	-
9) 5 + 6	blue	pink	clear	-
10) 2 + 5	blue	pink	blue*	+
11) 2 + 6	blue	pink	blue*	+

*depends on the relative amounts of each component.

It is important to note that combinations #7 and 8 give a false positive for cocaine, but cannot be screened out with the FPN test, as neither compound is a phenothiazine.

SEPARATION OF PCP - PROMAZINE COMBINATIONS

Dissolve mixture in 2.8 N HCl and extract with an equal volume of chloroform. Wash chloroform two times with fresh 2.8 N HCl. Extract PCP from chloroform with an equal volume of 0.1 N sulfuric acid. A U.V. spectrum of the 0.1 N sulfuric acid extract will reveal a promazine peak at 252 m μ with a 269 m μ shoulder representing the PCP (promazine 252 m μ E 1% 1 cm 1122, PCP 269 m μ E 1% 1 cm 9.2). Add three drops of concentrated HCl to the sulfuric acid extract, extract with an equal volume of chloroform and dry. The resulting extract will give an IR spectrum clearly recognizable as that of phencyclidine HCl.

- (1) Scott, L.J., Jr., "Specific Field Test for Cocaine", Microgram, 6:11, Nov. 1973, pages 179-181.
- (2) Clarke, E.B.C., "Isolation and Identification of Drugs", The Pharmaceutical Press, London, U.K. (1969).
- (3) Lorch, S.K., "Specificity Problem With the Cocaine Specific Field Test, and Its Solution", Microgram, 7:8, Aug. 1974, pages 100-101.

The author wishes to thank Dr. Fathi M. Saad, Chief, Warren Regional Laboratory, Bureau of Laboratories, Michigan Department of Public Health, for the G.C.-Mass Spectrographic identification of the promazine.

9/10/74

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

COLOR TEST FOR PHENCYCLIDINE AND METHADONE

Paul Zelonis
Forensic Chemist
Drug Enforcement Administration
Southwest Regional Laboratory

OBJECTIVE

To develop a field test possessing sensitivity and high specificity for use in powders, solvents, and plant material.

INTRODUCTION

This test was developed at the request of narcotics agents who desired a reliable field test for phencyclidine and its analogs. The commercial field test previously utilized by agents had failed on occasion to indicate the presence of phencyclidine and its analogs in powdered samples. Laboratory evaluation of the phencyclidine field test previously utilized was found to satisfactorily detect only the pure drug, and not the crude, unrefined material commonly encountered.

The Zelonis Test for phencyclidine and its analogs has enabled agents to successfully test varied materials at clandestine laboratory sites ranging from laboratory glassware residues and solvents, to treated plant material. This test can be utilized for field testing powders and liquids for methadone.

REAGENTS

- Solution A - 1 gram Platinum Chloride (reagent grade)
20 ml water (deionized)
5 ml Glacial Acetic Acid
- Solution B - 2 grams Cobaltous Thiocyanate
50 ml water (deionized)
50 ml Glycerin
- Solution C - Benzene

Solutions A and B should be aged 24 hours, and shaken prior to use.

METHODOLOGY

- 1) Add 3 drops of Solution A to a test tube containing 5-10 mg of powder, 0.5-1.0 ml liquid, or 100-200 mg of plant material.
- 2) Add 5 drops (0.25 ml) of Solution B and shake for 5 seconds.
- 3) Add 10 drops of Solution C and shake vigorously for 5 seconds.
 - (A) The top benzene layer will turn a clear blue color for powders and liquids containing Phencyclidine or its analogs, and Methadone.
 - (B) It is only necessary to wet the plant material and turn the test tube on a horizontal plane to note the blue color.

DISCUSSION

The volumes of Solutions B and C cited are approximations and are not critical for obtaining satisfactory results. If desired, larger volumes of Solutions B and C can be utilized when testing plant material without increasing the volume of Solution A. After utilizing the prescribed methodology with 120 controlled and non-controlled substances, it was found that the clear blue color also developed with methadone and dibucaine. A blue cloudy opaque benzene layer also developed with propoxyphene, bromodiphenhydramine, and PCC, a phencyclidine intermediate. Larger amounts of powder sample (100 mg) will turn the phencyclidine intermediate (PCC) from a cloudy blue to a distinct green within 10 minutes. For screening purposes, it should be noted that phencyclidine and methadone do not develop a Marquis reaction, which maybe used to differentiate them from propoxyphene and bromodiphenhydramine.

Although Solution A enhances the specificity of the test, it can be omitted without seriously jeopardizing the degree of specificity. Although this modification of the test does not have the same specificity of the recommended three reagent procedure, it is superior to other phencyclidine field test presently in use.

One ounce of reagent grade Platinum Chloride is sufficient for the preparation of approximately 5000 test. The sensitivity of the three reagent procedure is 0.5 mg of Phencyclidine powder.

CONCLUSION

No problems have been encountered by the Los Angeles Clandestine Laboratory group utilizing the recommended three reagent test. The test has gained agent acceptance and referral as being suitable for use by law enforcement officers.

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

LABORATORY NOTES

-46-

(4)

DATE March 28, 1974
NO.
DRUG TYPE Hallucinogens
METHODOLOGY Thin-layer Chromatography

SCREENING TEST FOR CANNABINOIDS, PHENCYCLIDINE,
DIMETHYLTRYPTAMINE, DIETHYLTRYPTAMINE.

Paul Zelonis, M.S.
Forensic Chemist
Mid-Atlantic Regional Laboratory
Drug Enforcement Administration

INTRODUCTION

Recent samples of plant material, hashish and smoking paraphernalia have been found to contain mixtures of cannabinoids and phencyclidine. A screening test is needed, which can be run in conjunction with routine marihuana thin-layer chromatograms, for other possible hallucinogens.

REAGENTS

95% Petroleum Ether and 5% Absolute Ethanol solution.
Fast Blue B - 1 gram Fast Blue B salt and 100 ml H₂O.

Iodoplatinate - 10 ml of 10% Chloroplatinic acid, 250 ml
of 4% Potassium Iodide diluted to 500 ml H₂O.

SOLVENT SYSTEMS

System #1 - 45% Butyl Ether, 45% Ethyl Ether, 10%
Diethylamine
System #2 - 5% Diethylamine and 95% Benzene

TLC PLATES

Analtech - Uniplate precoated TLC plates, Silica Gel GF 250 microns

DATA

Rf values (using solvent front as 1.00)

Solvent System #1		Solvent System #2	
PCP	0.91	PCP	0.64
CBD	0.69	CBD	0.52
THC	0.61	THC	0.41
CBN	0.47	CBN	0.30
DET	0.30	DET	0.15
DMT	0.18	DMT	0.09

PROCEDURE

Routine extraction of suspected plant material and smoking paraphernalia using a solution of 95% Petroleum Ether and 5% Absolute Ethanol. Development of thin-layer chromatograms using solvent tank systems 1 and 2. Initial visualization after drying is with Fast Blue B spray for screening of cannabinoids, with subsequent spraying with Iodoplatinate to screen for presence of PCP, DMT or DET.

DISCUSSION

The use of the polar solvent Ethanol aids in elution of PCP, DMT and DET off plant material and smoking paraphernalia, which can be then readily detected using the two sprays. The visualization of PCP, DMT and DET are not effected by the initial Fast Blue B spraying, the use of Iodoplatinate has negligible affect on the prior cannabinoid visualization.


VEGETABLE MATTER COATED WITH A BARBITURATE AND A VETERINARY ANESTHETIC

BY

H. HARRIS, R. BEATUS, M. BIANCHI, T. CATALANO, AND E. MC DONNELL

We recently received a rather unusual sample which we feel should be brought to the attention of other laboratories around the country. The sample was vegetable matter shown to be parsley which also gave a positive Wagners test indicating possible presence of an alkaloid. The material was screened for phencyclidine (PCP) using thin layer chromatography with the T1 developing system (see table 3). One spot migrated, but did not match the known PCP standard. On the theory that this might be a PCP analogue, a small amount was placed in methanol and the solution examined using our gas chromatograph-mass spectrometer (GC-MS). Two distinct peaks were observed in the total ion scan at retention times of 3.4 and 3.8 minutes using our general program of 150-260C at 16/min using a 6'x2 mm I.D. glass column packed with 3% OV-1. The mass spectrum obtained for the second peak clearly indicated either amobarbital or pentobarbital. These two isometric barbiturates have very similar mass spectral fragmentation patterns, but can be differentiated by comparison with spectra of authentic samples of the two barbiturates run under identical conditions. This comparison indicated that our peak at 3.8 minutes was definitely pentobarbital (see table 1). Further confirmation was obtained from infra-red and thin layer chromatographic data.

The first peak (retention time 3.4 minutes) showed a mass spectrum (figure 1) rich in detail with a base peak at m/e 180 and a possible parent ion at m/e 239. Further, doublets in the ratio of about 3:1 at m/e 180-182 and 209-211 suggested the presence of one chlorine in the molecule. Both computer and manual searches of the spectral reference libraries available to us failed to disclose a possible identification. Using the data we had gleaned from the mass spectrum, a manual search of Clarke's "Isolation and Identification of Drugs" disclosed one compound which had the proper molecular weight, contained a chlorine atom, and whose structure appeared to fit reasonably well with the observed fragmentation pattern. This material is used primarily as a veterinary anesthetic with the common name ketamine. The chemical name is 2-o-chlorophenyl -2- methylaminocyclohexane and has the structure shown below:




We were able to obtain an authentic sample of ketamine and confirm that the drug present along with pentobarbital was indeed ketamine. Ketamine does not appear in the subject index of "Microgram" and we have not seen it mentioned before in the forensic literature. The presence of barbiturate on vegetable matter is also very unusual in our experience. The "Physicians Desk Reference" in its section on Ketamine warns in a "Special Note": "Emergence reactions have occurred in approximately 12% of patients. The psychological manifestations vary in severity between pleasant dream-like states, vivid imagery, hallucinations, and emergence delirium".

We would like to alert other laboratories to look for this material. There are several points of similarity between ketamine and phencyclidine, i.e. both are primarily veterinary products and both appear to have mild hallucinogenic properties. In view of the very rapid rise of illicit use of PCP one should be alert for a similar occurrence for ketamine. Further data to aid in identification of Ketamine is given in Tables 2 & 3. It should be noted that Ketamine and methaqualone behave quite similarly with a number of spot tests and with several TLC systems. Therefore it is probably necessary to resort to IR, mass spectrometry or other more sophisticated methods to insure proper identification.

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copying
OBSOLETE DOCUMENT

FIGURE 1

Mass Spectrum of Ketamine



KETAMINE

STD

TABLE 1
MASS SPECTRAL DATA ON BARBITUATES

	RATIO m/e 41:43	m/e 83:85	m/e 155:157
Amobarbital	1.5	2.0	0.2
Pentobarbital	0.8	1.0	0.45
Peak at R.T. 3.4	0.8	1.0	0.5

TABLE 2
PHYSICAL AND CHEMICAL PROPERTIES OF KETAMINE
Major Infra-Red bands of Ketamine Hydrochloride in microns
2.95, 3.4, 3.5, 3.55, 5.8, 5.85, 6.85, 7.0, 7.85, 8.1, 8.6,
8.8, 9.25, 9.45, 9.65, 10.35, 10.5, 11.1, 11.85, 12.3, 13.1
13.3, 13.7, 13.9

Spot Tests and Crystal Tests For Ketamine

Wagners (Alkaloid)	+
Marquis	-
Mandelin	-
Co Thiocyanate	+
+SnCl	-
Liebermann	Slow Purple
Permanganate	Crystals

TABLE 3

Thin Layer Chromatography Data

Compound	Solvent System		
	T1	S2	S16
Ketamine	.72	.86(Gray)	.55(Br.)
Methaqualone	.70	.88(Br.)	.50(Purple)
Methylphenidate	.60	.83(Tan)	.55(Yellow)
Cocaine	.60	.89(Tan)	.67(Br.)
Phencyclidine	.59	.92(Purple)	.90(Purple)

T1 Methanol-Conc. Ammonia (100:1.5)
S2 Chloroform-Dioxane-Ethyl Acetate-Ammonia (25:60:10:5)
S16 Hexane-Benzene-Diethylamine (75:25:10)

Colors in () produced by Iodoplatinate spray

6

DIFFERENTIATION OF PCP, TCP, AND A CONTAMINATING
PRECURSOR PCC, BY THIN-LAYER CHROMATOGRAPHY

Alexander T. Shulgin
1483 Shulgin Road
Lafayette, California

Carmen Helisten
PharmChem Research Foundation
Palo Alto, California

BACKGROUND

Analyses of drugs alleged to be mescaline or synthetic THC have often shown them to be phencyclidine (Sernyl, PCP). This base has also been encountered as a component in drug mixtures involving marijuana and LSD. Recently a pharmacologically related analog thienylcyclohexylpiperidine (TCP) has appeared in the place of, or in admixture with, PCP.

In the illicit synthesis of both PCP and TCP, the intermediate 1-piperidinocyclohexanecarbonitrile (PCC) is usually employed. This material has extraction properties that are extremely similar to those of both PCP and TCP, resulting in the preparation of a contaminated product. As a toxicological problem, this contaminant may contribute to the toxic aspects of these two drugs and its presence or absence might be of value in the treatment of cases of overdose. From the forensic point of view, the ability to detect this precursor should be potentially valuable in two ways: 1. It will help establish the synthetic procedure employed in the manufacture of the drug in question, and 2. As PCC is not found as a contaminant in preparations intended for clinical application, its presence will establish the unlicensed nature of this manufacture.

Attempts to effect this analysis by conventional GLC procedures are thwarted by the thermal instability of PCC. Under ordinary chromatographic conditions, HCN is split out of the molecule resulting in the formation of the enamine cyclohexenylpiperidine. This process occurs during actual passage through the GLC column resulting in an inconsistent spectrum that depends upon column conditions, quantities injected, and the nature of the contamination of the injected sample. A further complication is the inherent instability of this enamine in that, if hydrolytic conditions are encountered, it further degrades to its two components piperidine and cyclohexanone. Both of these latter chemicals are invariably lost under the solvent peak on GLC analysis. Attempts to analyze these mixtures by conventional TLC procedures (activated plates, heated spot application) result in immediate and complete degradation of PCC to these same components.

OBJECTIVE

To provide a chromatographic procedure for the separation and localization of PCP, TCP, and PCC in a single extraction and

MATERIALS

TLC plates: E. Merck, 0.25 mm Silica Gel 60 prepared glass plates which have been deactivated by open-air storage for at least 48 hrs prior to use.

Solvent system: Benzene, acetone and pyridine, in the ratio of 16:8:1.

Visualization reagents: Ninhydrin (0.8% in acetone) and iodoplatinate (4 g chloroplatinic acid and 24 g KI in 1 l. water diluted, after standing, with 1 l. methanol).

PROCEDURE

Approximately one dosage unit of the drug to be analysed is suspended in four drops of anhydrous methanol, and broken up to a fine powder. No external heat is employed. An unactivated TLC plate is appropriately marked (origins located for at least four sample applications, i.e., reference PCP, reference TCP, reference PCC, and the unknown). The site of application of the unknown sample is wetted with a drop of the chromatographic solvent mixture. As soon as the apparent dampness has disappeared, but before the residual discoloration is also lost (the spot is still chalky white against a dull white background), about 5 ul of the methanolic extract is applied. Similar care is needed for the reference PCC solution. The solutions of the two standard drugs PCP and TCP may be applied in the conventional manner except that there must be no heating of the plate.

After air drying, the chromatogram is developed for about 10 cm, air dried until largely free of pyridine smell, then heated in an air oven at 110° for 10 min. While still hot it is sprayed with ninhydrin, and the developed colors noted. The plate is allowed to cool, and oversprayed with iodoplatinate. The following Rf's and colors should be observed:

Compound	Rf	Ninhydrin color:	Iodoplatinate color:
Piperidine (decomposition product of PCC)	0.00	purple	unchanged, fading to a bleached spot on standing
PCP	0.24	none	grey, permanent
TCP	0.49	purple	unchanged, fading to a bleached spot on standing
PCC	0.70	purple	unchanged, fading to a bleached spot on standing

COMMONWEALTH OF VIRGINIA
Division of Consolidated Laboratory Services
BUREAU OF FORENSIC SCIENCE

A. W. TIEDMANN, PH. D.
DIVISION DIRECTOR

CHARLES E. O'REAR, PH. D.
DEPUTY DIRECTOR, FORENSIC SCIENCE

CHARLES L. KILLION
NORTHERN REGIONAL DIRECTOR
2714 DORR AVE.
P. O. BOX 486
MERRIFIELD, VA 22116
(703) 572-8236

Identification of 1-Piperidinocyclohexane
Carbonitrile


Robin D. Porter, M.S.
Analytical Chemist
Northern Virginia Branch
Bureau of Forensic Science

OBJECTIVE

To report occurrence of and provide analytical data for the identification of 1-piperidinocyclohexane carbonitrile, an intermediate in the preparation of phencyclidine (PCP).


BACKGROUND

This laboratory has encountered in recent months a material present, not only in numerous exhibits of phencyclidine (powders, tablets and on plant material), but also in large quantities in pure form. The material was identified as 1-piperidinocyclohexane carbonitrile, an intermediate in the two-step preparation of phencyclidine. Although not readily available commercially, the intermediate is easily synthesized by the reaction:



The reaction is not only fairly rapid but produces the intermediate in good yield.

The subsequent reaction in the production of phencyclidine requires greater technique because of the involvement of the Grignard reagent.



The entire synthesis can be done in approximately four hours with substantial yields. However, when less than a stoichiometric amount of Grignard reagent is used, the final product will be a mixture of phencyclidine and the intermediate, as we have been encountering in the Northern Virginia area.


PROCEDURE


The initial sample of the intermediate encountered (a white powder) was in combination with phencyclidine. Thin layer chromatography (CHCl₃ - MeOH (9:1)) showed the presence of one other material (R_f - 0.82) in addition to phencyclidine (R_f - 0.43). In a benzene-methanol (3:1) solvent system, the second material

an Rf of 0.74 as compared to an Rf of 0.46 for phencyclidine. The second material also developed with acidified iodoplatinate as did phencyclidine. Ultra-violet spectroscopy of the mixture showed no variations from the expected phencyclidine spectrum; however, the infrared spectrum of a basic extract of the sample showed differences from the standard phencyclidine spectrum, the most noteworthy being the addition of a $-C\equiv N$ stretching band at 2220 cm^{-1} . On the gas chromatograph (5% OV-17; column temperature 180° C.), the relative retention time of the second material compared to phencyclidine was 0.11 (Figure 1). Solubility similarities precluded separation by solvent/solvent extraction, thus necessitating preparative thin layer chromatography. Using 1000 microns Silica Gel G plates in a $\text{CHCl}_3\text{-MeOH}$ (9:1) solvent system, the intermediate was isolated and identified on the basis of infrared (Figure 2) and mass spectral characteristics (Figure 3). NMR (Figure 4) data was of little aid in structure determination due to similar methylene protons. The intermediate (M.P. $67.0\text{-}68.0^\circ\text{ C.}$) does not absorb in the UV region. The identification was further confirmed by comparison to an independently synthesized authentic sample.


DISCUSSION

1-Piperidinocyclohexane carbonitrile continues to be seen in both the pure form (90-98 percent) and also in combination with phencyclidine. In some cases, the concentration of the intermediate is three to four times that of phencyclidine. Even though the intermediate is not controlled by federal or state codes, the identification of the intermediate in samples of phencyclidine has provided law enforcement officials with useful information as to possible common source of manufacture.





MICROGRAM, VOL. VIII, NO. 10 (October, 1975)



7.0
6.5
6.0
5.5
5.0
4.5
4.0
3.5
3.0
2.5
2.0
1.5
1.0
0.5
0.0

(7)

R. L. Epstein,¹ Ph.D.; Philip Lorimer,¹ B.A.; and
E. J. Sloma,¹ B.A.

Identification of Phencyclidine-Related Drugs

The appearance of new drugs within the crime laboratory necessitates the development and improvement of analytical schemes for their detection. One class of particular interest is the phencyclidine-related drugs. The drugs studied are 1-(1-phenylcyclohexyl) piperidine (PCP), commonly known as phencyclidine, "Angel's Dust," or DOA [4]; the two homologs 1-(1-phenylcyclohexyl) pyrrolidine (PIHP) and 1-(1-phenylcyclopentyl) piperidine (PPP); and an analog 1-(2[thienyl]cyclohexyl) piperidine (TCP).

The analytical procedures evaluated resulted in a positive identification of these drugs. They are microchemical tests, chemical ionization mass spectrometry (CIMS), thin-layer chromatography (TLC), and gas-liquid chromatography.

A comparison of the four drugs by chemical data is presented in Table 3. The ultra-violet (UV) spectrophotometric data, previously compiled [2,3], are listed in Table 2. As expected [4], only the thiophene analog (TCP) is markedly different. Infrared spectra are also available for each drug but identification by this method can be difficult with a sample that is highly adulterated or in trace quantities.

Experimental Methods

Approximately 3 mg of phencyclidine and the two homologs were placed on a glass slide with 10% aqueous hydrochloric acid or 10% aqueous acetic acid with 2% aqueous potassium permanganate. The crystal formations were observed on a compound microscope at $\times 400$ and are depicted in Figs. 1 and 2.

The chemical ionization mass spectra were taken for all four drugs on a Dupont 21-490 single focusing mass spectrometer. The reagent gas was isobutane (99.9%). All drugs analyzed were admitted through direct probe. The instrument operating conditions were the same as described by Sauerstein et al [5]. The data are tabulated in Table 3.

Thin-layer chromatography was conducted on 250- μ m silica gel plates manufactured by Analtech, Inc., Newark, Del. All drugs were extracted as the free base and spotted with chloroform. After development, the drugs were visualized with 5% aqueous potassium iodoplatinate producing blue-gray spots. All chemicals and solvents were reagent grade and supplied by J. T. Baker, Phillipsburg, N.J. The following solvent systems were used with the resulting $R_f \times 100$ values compiled in Table 4:

- (A) cyclohexane:benzene:diethylamine (75:15:15),
- (B) ethyl acetate:methanol:ammonium hydroxide (85:10:5),
- (C) acetone:dimethylformamide:ammonium hydroxide (85:15:0.5).

Received for publication 7 April 1976; revised manuscript received 1 July 1976; accepted for publication 2 July 1976.

¹Supervising forensic chemist, forensic chemist, and senior forensic chemist, respectively, New Jersey State Police, North Regional Laboratory, U.S. Highway 46, Little Falls, N.J. 07424.

TABLE I—Chemical data.

Drug	Abbreviation	Empirical Formula (Free Base)	Molecular Weight (Free Base)	Molecular Weight (HCl Salt)	Source
1-(1-phenylcyclohexyl) piperidine (phencyclidine)	PCP	C ₁₇ H ₂₇ N	243.4	279.9	clandestine
1-(1-phenylcyclohexyl) pyrrolidine	PHP	C ₁₄ H ₂₁ N	229.4	265.9	Lot P*
1-(1-phenylcyclopentyl) piperidine	PPP	C ₁₆ H ₂₃ N	229.4	265.9	Lot P*
1-(1-[2-thienyl]cyclohexyl) piperidine	TCP	C ₁₈ H ₂₅ NS	249.5	275.8	Lot Q*

*Parke Davis & Co., Midland, Mich.

TABLE 2—Ultraviolet spectra.

Drug	Wavelength Maxima,* 0.1N HCl
PCP	<u>262</u> , 257, 268, and 252
PHP	<u>262</u> , 257, 268, and 252
PPP	<u>261.5</u> , <u>256.5</u> , 268, and 252
TCP	<u>232</u>

*Strongest maxima are underlined.

- (D) methanol:acetic acid (90:10),
 (E) chloroform:methanol (90:10),
 (F) chloroform:methanol (80:20),
 (G) acetone:chloroform (50:50), and
 (H) acetone:chloroform (65:35).

A Varian 2700 dual column gas chromatographic with flame ionization detectors at 270°C was used with two 6-ft (1.8-m) by 1/4-in. (6.35-mm) outside diameter and 2-mm inside diameter glass columns from Analabs, Inc., North Haven, Conn. These columns were treated with dimethyldichlorosilane and packed with 3% OV-17 and 3% OV-101 on 100/120 mesh Gas Chrom Q. The column temperature was 70°C, and the injection port temperature was 250°C. A flow rate of 30 ml/min was used with medical-grade nitrogen as the carrier gas. The detector gases were H₂ (99.999%) and compressed air. All samples were injected as the free base with an 80- μ l Hamilton syringe. A graphical presentation was presented on Hewlett Packard Corp. 3380 A integrator. The tabulated data are presented in Table 5 and Fig. 3.

Discussion

The crystal formations produced by PCP and PPP in 10% aqueous HCl and 2% KMnO₄ were indistinguishable. They can be described as violet, H-shaped plates, whereas PHP produced crystals also resembling violet, H-shaped plates, only thinner (Fig. 1). In 10% acetic acid and 2% KMnO₄ the crystals produced for phenacyclidine and its homologs were somewhat dissimilar and can be used for differentiation. Again, PCP formed violet, H-shaped plates, but PHP formed orange, H-shaped plates. The PPP homolog produced definite, violet, X-shaped crystals by this test, as shown in Fig. 2.

The chemical ionization mass spectra for the drugs studied produced the protonated molecular ion (MH⁺) in greatest abundance. Phenacyclidine (PCP) with an *m/e* of 244 and the thiophene analog (TCP) having an *m/e* of 250 were easily distinguished from the two homologs PPP and PHP having an *m/e* of 230.

The molecular ion M⁺ and the M⁻¹ ion were also present. The homologs PHP and PPP (molecular weight 229.3) produced the same protonated molecular ions and molecular ions. However, the fragmented ions for PHP and PPP, *m/e* 159 and 165, respectively, allowed differentiation (Table 3).

A fragmented ion product of PCP (*m/e* 159) postulated by Hauber [6] as phenylcyclohexene was also present for PHP. Furthermore, the fragmented ion (*m/e* 145) produced by PPP probably corresponds to phenylcyclopentene and the *m/e* of 165, resulting from the fragmentation of the thiophene analog, is most likely thienylcyclohexene.

The alkaline TLC Systems A, B, C and the acidic solvent System D failed to produce adequate resolution. The best results were obtained with the neutral Systems E through H, which successfully resolved PCP, PHP, PPP, and TCP with excellent reproducibility. These systems were composed of various proportions of chloroform with methanol

TABLE 4—Thin-layer chromatography $R_f \times 100$ values.

Drug	Solvent Systems							
	A	B	C	D	E	F	G	H
PCP	69	86	98	54	17	29	23	32
PHP	69	75	95	57	10	24	05	13
PPP	69	80	98	54	25	37	29	40
TCP	69	82	98	55	33	51	45	49

TABLE 5—Gas chromatography data.

Drug	OV-101 Gas Chrom Q		OV-17 Gas Chrom Q	
	Retention Time, min	Relative Time ^a	Retention Time, min	Relative Time ^a
PPP	4.36	0.65	4.06	0.62
PHP	4.94	0.74	4.83	0.75
PCP	6.70	1.00	6.52	1.00
TCP	6.48	0.97	6.47	0.99

^a Based on phencyclidine (PCP).

(Systems E and F) and acetone with chloroform (Systems G and H) (Table 4). There was little difficulty in reproducing the chromatographic study of Shulgin [7] for separating PCP and TCP.

Gas chromatography (GC) completely resolved the homologs in the order of PPP, PHP, and PCP using either 3% OV-17 or 3% OV-101 (Table 5). The liquid phases chosen represent the divergent polarities as exhibited by their McReynold's constants [8]. Better resolution of PCP and TCP could be obtained by lowering the column temperature but these drugs can be easily distinguished by their UV spectra as well as by TLC, as previously discussed. A GC/MS interface could probably serve as an efficient tool for rapidly identifying these structurally related compounds.

Some of these phencyclidine-related drugs were studied by Bailey et al [9]; the study resulted in a separation scheme using absorption spectrometry, electron impact MS, protonated magnetic resonance spectra, and several TLC and GC techniques. The drugs of most interest to forensic laboratories, namely PCP, TCP, and PHP, were only partially resolved in the chromatography techniques described. In Bailey's study PPP was not included.

Conclusion and Summary

The analytical techniques presented allow the forensic chemist to readily differentiate phencyclidine (PCP) from the two homologs 1-(1-phenylcyclohexyl) pyrrolidine (PHP) and 1-(1-phenylcyclopentyl) piperidine (PPP) as well as the analog 1-(2-[thienyl] cyclohexyl) piperidine (TCP).

The UV spectra in dilute mineral acid for PCP, PHP, and PPP are indistinguishable, but the TCP spectrum is markedly different because of the thiophene moiety. The $KMnO_4$ crystal test using HCl can only suggest the presence of a phencyclidine-related drug, but acetic acid does offer more distinguishability.




FIG. 3—Gas chromatographic data.

Chemical ionization mass spectra of all four drugs studied are different, except that the protonated molecular ion (MH^+) for the homologs PHP and PPP are the same. Both TLC and GC contribute to the separation and confirmation of these drugs after preliminary testing indicated their presence.

Acknowledgments

The authors wish to thank Lt. Joseph Barry and the New Jersey State Police for their cooperation as well as Ms. Michele Senko for her assistance in preparing this paper. The mass spectral data was compiled with the cooperation of the Forensic Science Bureau, New Jersey State Police, West Trenton, N.J. 08625.

References

- [1] *Micrograms*, Vol. 2, No. 1, Jan. 1969, p. 3.
- [2] DeZan, P. and Bianchi, P., "Analysis and Identification of Phencyclidine Hydrochloride," *Micrograms*, Vol. 1, No. 12, 1968, pp. 14-18.
- [3] Dobres, H., "Identification of PCP Analogs," *Micrograms*, Vol. 8, No. 2, 1975, pp. 27-30.
- [4] Heagy, J. A., "Thiophene Analog of PCP," *Micrograms*, Vol. 5, No. 11, 1972, pp. 122-124.
- [5] Saferstein, R., Chao, J., and Manura, J., "Identification of Drugs by Chemical Ionization Mass Spectroscopy—Part II," *Journal of Forensic Sciences*, Vol. 19, No. 3, 1974, pp. 463-485.
- [6] Hauber, D., "Identification of a Degradation Peak in Phencyclidine by Gas Chromatography," *Micrograms*, Vol. 8, No. 7, 1975, p. 100.
- [7] Shulgin, A. T., "Differentiation of PCP, TCP and a Contaminating Precursor PCC by Thin-Layer Chromatography," *Micrograms*, Vol. 8, No. 10, 1975, pp. 149-150.
- [8] Supina, W. R., *Handbook of McReynolds Constants*, Supleco, Inc., Bellefonte, Pa., 1971.
- [9] Bailey, K., Gagne, D. R., and Pike, R. K., "Identification of Some Analogs of the Hallucinogen Phencyclidine," *Journal of the Association of Official Analytical Chemists*, Vol. 59, No. 1, 1976, pp. 81-89.

Robert L. Epstein, Ph.D.
New Jersey State Police
North Regional Laboratory
U.S. Highway 46
Little Falls, N.J. 07424


Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

8

ANALYSIS AND IDENTIFICATION OF
1-[1-(2-THIENYL)CYCLOHEXYL]PIPERIDINE (TCP)

by David R. Picard

Wisconsin Department of Justice
Division of Law Enforcement Services
Crime Laboratory Bureau
4706 University Avenue
Madison, Wisconsin 53702



1-[1-(2-thienyl)cyclohexyl]piperidine (TCP), the thiophene analogue of phencyclidine (PCP) has been reported on the illicit drug market.^{1,2,3} It was also encountered on several occasions in 1975 from northwest Wisconsin.

A search of the literature revealed a paucity of spectral data, although Heagy⁴ reported the ultraviolet maximum and minimum and the infrared spectrum of the hydrochloride salt.

Since the occurrence of this drug appears to be on the increase nation wide, a method of analysis including spectroscopic methods of identification is hereby reported.

EXPERIMENTAL

Synthesis

1-[1-(2-thienyl)cyclohexyl]piperidine was synthesized according to Parell⁵, purified by extraction and thin-layer chromatography, and subjected to the following examinations. This was the source of the spectra material in figures 1-3.

Physical Properties:

a. Hydrochloride Salt $C_{15}H_{23}NS \cdot HCl$ Mol. Wt. 285.85

Appearance: fine, white, crystalline powder

Solubility: soluble in chloroform and water,
insoluble in hexane and ether

Melting Point: 233-6° with transition at 182-3°.5

b. Free base $C_{15}H_{23}NS$ Mol. Wt. 249.40

Appearance: viscous, light yellow oil

Solubility: soluble in chloroform, ether, and hexane

Boiling Point: 110-111° at 0.2 mm⁵

Color Tests

Color tests were performed on pure hydrochloride salts of both TCP and PCP for purposes of differentiation.

Color Reagent

TCP

PCP

Marquis gas evolved - faint - grey-brown gas evolved - faint pink salmon

Mecke's gas evolved - orange gas evolved

Mandelin's gas evolved - green gas evolved

All color reactions occur within one minute after reagent addition.

Extraction

Illicit samples were purified by adding a sample of the crude powder to 5 ml. of dilute NaOH. The free base just formed was converted to the chloroform insoluble sulfate by addition of 10 ml 0.5NH₂SO₄. Two 10 ml. chloroform extracts removed neutral and

acidic impurities. The aqueous phase was then rendered alkaline by the addition of NaOH and extracted with 10 ml. of CHCl_3 . The organic phase was dried with Na_2SO_4 , filtered and evaporated to yield a pale yellow oil.

Gas Chromatography

Column Conditions: 200° C; N_2 flow 60 ml/min.

Column: 6 ft. x 2 mm i.d. 3% OV-101 on H.P.

Chromosorb W (80-100 mesh)

Retention time relative to phencyclidine (PCP) - 1.0

Using the above conditions it was not possible to differentiate TCP from PCP.

Thin-Layer Chromatography (TLC)

Analytical TLC was carried out using commercial precoated silica gel GF plates (thickness 200 microns, available from Analtech, Inc., Newark, Delaware). CHCl_3 , MeOH (9:1) was used as the eluting solvent. The spots were visualized with acidified iodoplatinate. Separation of TCP, PCP and 1-piperidinocyclohexane carbonitrile (PCC), the precursor to both of these compounds, was obtained.

<u>Compound</u>	<u>Rf x 100</u>
TCP	67
PCP	37
PCC	91

Preparative TLC was carried out using plates prepared in our laboratory. Silica gel GF-254 (Acc. to Stahl) supplied by Brinkman Instruments, Inc., Des Plaines, Illinois, was spread to a thickness of 0.5 mm on 20 x 20 cm. glass plates and oven dried at 110° C until used.

These plates, when used with CHCl_3 , MeOH (9:1) as solvent,

were suitable for purification of street samples of TCP.

Ultraviolet Spectrophotometry

Purified TCP was run on a Perkin-Elmer 402 ultraviolet-visible spectrophotometer using 0.5N H₂SO₄ as solvent. (figure 1)

Infrared Spectrophotometry

Both the free base and hydrochloride salt of TCP were subjected to infrared spectrophotometry using a Perkin-Elmer 467 Grating Infrared Spectrophotometer (figure 2).

Gas Chromatography/Mass Spectrometry

GC/MS worked well as a means of identifying suspected TCP, (figure 3). Gas chromatography was performed on a 6 ft. x 2 mm i.d. 3% OV-101 on H.P. Chromosorb W (80-100 mesh) column at 210° with a helium flow of 45 ml/min. The transfer line between the gc and the mass spectrometer was kept at 225°.

A Hewlett-Packard 5930A dodecapole mass spectrometer with an electron energy of 70 e.v. was used. The ionization and analyzer areas were maintained at a temperature of 150°.




Comparison of the spectrum of TCP with that of PCP (see table 1) indicates the similarity in the structures of these two compounds.

NOTE: During the preparation of this manuscript Bailey, Gagne and Pike⁶ published spectroscopic and chromatographic data on TCP and other analogs of phencyclidine.

REFERENCES

- ¹Microgram, Vol. VIII, No. 5 (May, 1975), pp. 59-60.
- ²Microgram, Vol. VIII, No. 9 (September, 1975), p. 126.
- ³Microgram, Vol. IX, No. 1 (January, 1976), p. 1.
- ⁴Heagy, J. A., Microgram, Vol. V, No. 11 (Nov., 1972), pp. 122-124.
- ⁵Parell, U.S. Patent 2921076 (1960); Chem. Abstr., Vol. 54, 12159c (1960).
- ⁶Bailey, K., Gagne, D., Pike, R., JAOAC, Vol. 59, No. 1 (Jan., 1976) pp. 81-89.

T A B L E 1

<u>m/e</u>	<u>TCP</u>	<u>Probable Ion Composition</u>	<u>PCP</u>
249		M^+	243
220		$(M-C_2H_5)^+$	214
206		$(M-C_3H_7)^+$	200
192		$(M-C_4H_9)^+$	186
166			166
165		$(M-$  $)^+$	---
164		$(M-C_5H_{11}N)^+$	158
149		$(M-C_6H_{14}N)^+$	143
136		$(M-C_7H_{15}N)^+$	130
123		$(M-C_8H_{16}N)^+$	117
110		$(M-C_9H_{17}N)^+$	104
97		$ArCH_2^+$	91
84			84

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT




figure 1







figure 2

Property of Idaho State Police Forensic Services
UP controlled Intelnet Copy
OBSOLETE DOCUMENT

figure 3



1-[1-(2-thienyl)cyclohexyl]piperidine
GCMS @ 210° 70 e.v.



Mol. Wt.

249

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

DRUGS

Identification of Some Analogs of the Hallucinogen Phencyclidine

KEITH BAILEY, DENISE R. GAGNÉ, and RICHARD K. PIKE
 Drug Research Laboratories, Health Protection Branch, Tunney's Pasture,
 Ottawa, Ontario, Canada K1A 0L2

The drugs 1-[1-(2-thienyl)cyclohexyl]piperidine, 1-[1-(2-thienyl)cyclohexyl]morpholine, 1-[1-(2-thienyl)cyclohexyl]pyrrolidine, 1-(1-phenylcyclohexyl) morpholine, and 1-(1-phenylcyclohexyl) pyrrolidine are identified by spectroscopic techniques. The ultraviolet and proton magnetic resonance spectra of analogs are similar, but mass and infrared spectra are distinctly different, and reference spectra are provided. Gas-liquid and thin layer chromatographic systems for the analysis are discussed.

Phencyclidine (I) hydrochloride, 1-(1-phenylcyclohexyl)piperidine hydrochloride, is marketed as Sernylan,® an anesthetic used in veterinary medicine; it is also subject to abuse (1). It has appeared in preparations that are sold on the illicit market and which are said to contain lysergic acid diethylamide (LSD), cocaine, mescaline, or tetrahydrocannabinol rather than phencyclidine (2). Phencyclidine and its salts and derivatives are scheduled under the Narcotic Control Act in Canada and phencyclidine is designated a depressant in Schedule III of the Controlled Substances Act of the United States. However, analogs of phencyclidine have recently been encountered by narcotic control agencies (3, 4). There is obviously a need for methods which provide an unequivocal identification of these compounds.

This paper describes chromatographic and spectroscopic methods for identifying 5 analogs of phencyclidine: 1-[1-(2-thienyl)cyclohexyl]piperidine (II), 1-[1-(2-thienyl)cyclohexyl]morpholine (III), 1-[1-(2-thienyl)cyclohexyl]pyrrolidine (IV), 1-(1-phenylcyclohexyl)morpholine (V), and 1-(1-phenylcyclohexyl)pyrrolidine (VI).

Experimental

Compound I was obtained by recrystallization of an authentic sample supplied by the Division of Pharmaceutical Chemistry, Health and Welfare Canada. Compound II was made by the reaction of 2-magnesium bromothiophene on 1-(1-cyano-

cyclohexyl)piperidine; the latter intermediate resulted from the addition of the elements of hydrogen cyanide to the piperidine enamine of cyclohexanone. Compounds III-VI were obtained from the corresponding enamine of cyclohexanone by formation of its *p*-toluenesulfonic acid quaternary salt and subsequent attack by the appropriate aryl Grignard reagent. The compounds were purified as the hydrochloride salts by recrystallization from mixtures of isopropanol with ether. Melting points were measured on a Koffler hot stage and are uncorrected, and elemental analysis indicated that the salts were essentially anhydrous (Table 1). The spectra of free bases were recorded on the analog regenerated from the salt with Na₂CO₃ solution and extracted into CHCl₃; the CHCl₃ was removed by warming the solution under a stream of nitrogen. Thin layer chromatograms were developed 15 cm, in saturated tanks under ambient conditions, using precoated plates and sheets as received. They were examined under 254 nm ultraviolet (UV) light and sprayed with a solution of





Table 1. Melting point and analytical data for some analogs of phenacyclidine as their hydrochlorides

Compound	Mp, °C	Found, %			Theoretical, %		
		C	H	N	C	H	N
II	200-203 subl.	62.98	8.43	4.91	63.02	8.46	4.90
III	162-165; 174	58.47	7.68	4.94	58.42	7.71	4.87
IV	187-188	61.87	8.09	5.17	61.85	8.16	5.15
V	180-190 subl.	68.05	8.63	4.08	68.19	8.58	4.97
VI	221-222 subl.	72.29	9.16	5.24	72.29	9.10	5.24

potassium iodoplatinate (5). Mass spectra were determined on a Hitachi Perkin-Elmer Model RMU-61, magnetic deflection instrument, operating at 160-180°C, ionization voltage 70 ev, and acceleration voltage 4-5 v. Samples were introduced via the probe. Infrared (IR) spectra were recorded on a Beckman IR 20A and UV spectra were recorded on a Beckman BD GT spectrophotometer. Proton magnetic resonance spectra were recorded on a Varian A-60A spectrometer. Gas-liquid chromatograms were obtained on a Hydro-flow Series 3000 instrument.

Results and Discussion

Mass Spectra

The compounds as salts give characteristic fragmentation patterns on electron impact; their normalized spectra are shown in Fig. 1. The molecular ions vary widely in relative intensity (12-60% of base peak) but can be assigned. All the compounds gave a brace of neighboring peaks corresponding to the molecular ion minus (a) the mass of the nitrogen heterocycle radical fragment and (b) the mass of the nitrogen heterocycle. This phenomenon is much more pronounced in the thiophene-containing series, one of these peaks usually being the base peak. As would be expected, cleavage of the aromatic moiety leads to a significant peak ($M-83$ in the thiophene series and $M-77$ in the phenyl series). The peak is more noticeable in the phenyl series, a reflection of the stronger aromaticity of the phenyl ring and the consequently greater stability of the phenyl radical. All members of the thiophene series give a strong peak at m/e 97 (70-90% of base peak) while those of the phenyl series give a strong peak at m/e 91 (80-95% of base peak). These are highly characteristic of the mass spectra of thiophene and benzene derivatives (6). All of the analogs and phenacyclidine itself gave a prominent peak at $M-43$. In the thiophene series this varied in intensity

(50-75% of base peak), while in the phenyl series this was the base peak in each case.

The mass spectra are very informative; they can distinguish between the compounds presented in this paper and could be used to elucidate the structure of uncharacterized isomers, analogs, and homologs of phenacyclidine.

Ultraviolet Spectra

The wavelengths and molar absorptivities of the maxima in the spectra of the compounds (Table 2) closely resemble those of the corresponding aromatic moieties, enabling facile distinction between the 2 classes with reference to their aromatic substituent but not among analogs within each of these 2 classes.

Proton Magnetic Resonance Spectra

The spectra were examined as the free bases in $CDCl_3$ solution and as their hydrochloride salts in D_2O . The presence of a multiplet of absorbances (6.85-7.30 δ) due to the thiophyl protons immediately distinguishes compounds containing this ring and their corresponding phenyl relatives whose phenyl protons give much simpler absorption (7.2-7.5 δ) at 60 MHz.

Infrared Spectra

IR spectra of the free bases (films between NaCl plates) and of the hydrochlorides (0.75% in KBr disks) are presented in Figs. 2 and 3. The spectra of phenacyclidine have been published (5) but are included for comparison purposes.

Table 2. Ultraviolet data^a for some analogs of phenacyclidine

Compound	$\lambda_{max}(\epsilon)$	$\lambda_{max}(\epsilon)$	$\lambda_{max}(\epsilon)$	$\lambda_{max}(\epsilon)$
I	251(301)	256(371)	261(413)	268(322)
II	232(7933)			
III	232(7944)			
IV	232(7837)			
V	251(288)	256(363)	261(416)	268(325)
VI	251(247)	256(327)	261(373)	268(293)

^a Solutions of the hydrochloride in ethanol, λ_{max} in nm.

Received August 5, 1975.






FIG. 1—Normalized mass spectra of phenacyclidine (I) and its analogs: III, 1-(1-

Morphides	
%	
N	
4.90	
4.87	
5.15	
4.97	
5.24	

in the phenyl
b case.
mative; they
compounds pre-
used to eluci-
rized isomers,
lidine.

absorptivities of
the compounds
of the corre-
g facie dis-
h reference to
at among ana-

the free bases
hydrochloride
multiplet of ab-
the thienyl pro-
compounds con-
sponding phenyl
ave much sim-
MHz.

(films between
chlorides (0.75%
Figs. 2 and 3.
been published
rison purposes.

or some

$\nu_{max}(c)$	$\lambda_{max}(e)$
1(413)	268(322)

11(416)	268(325)
11(373)	260(293)

ethanol, λ_{max} in nm.




FIG. 1—Normalized mass spectra of I, 1-(1-phenylcyclohexyl)piperidine HCl; II, 1-(1-(2-thienyl)cyclohexyl)piperidine HCl; III, 1-[1-(2-thienyl)cyclohexyl]morpholine HCl; IV, 1-[1-(2-thienyl)cyclohexyl]pyrrolidine HCl; V, 1-(1-phenylcyclohexyl)morpholine HCl; and VI, 1-(1-phenylcyclohexyl)pyrrolidine HCl.







FIG. 2—IR spectra of 1,1-(2-thienyl)ethane

15






FIG. 2.—IR spectra of some analogs of phencyclidine as their bases, NaCl film. I, 1-(1-phenylcyclohexyl)piperidine; II, 1-[1-(2-thienyl)cyclohexyl]piperidine; III, 1-[1-(2-thienyl)cyclohexyl]morpholine; IV, 1-[1-(2-thienyl)cyclohexyl]pyrrolidine; V, 1-(1-phenylcyclohexyl)morpholine; and VI, 1-(1-phenylcyclohexyl)pyrrolidine.

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

m/e 235
111

m/e 245
111

m/e 229
111







FIG. 3—IR sp
hexyl)piperidin
cyclohe

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT







FIG. 3—IR spectra of some analogs of phencyclidine as their hydrochloride salts, KBr disks. I, 1-(1-phenylcyclohexyl)piperidine; II, 1-[1-(2-thienyl)cyclohexyl]piperidine; III, 1-[1-(2-thienyl)cyclohexyl]morpholine; IV, 1-[1-(2-thienyl)cyclohexyl]pyrrolidine; V, 1-(1-phenylcyclohexyl)morpholine; and VI, 1-(1-phenylcyclohexyl)pyrrolidine.



Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

Table 4
Comp

* Columns were
w (HP), exceptin

Table 4

System ^a	Plate
A	Br.
A	Ea.
B	Br.
B	Ea.
C	Br.
C	Ea.
D	Br.
D	Ea.
E	Br.
E	Ea.
F	Br.
F	Ea.
G	Br.
G	Ea.
H	Br.
H	Ea.
I	Br.
I	Ea.

^a A = ethyl
water-ammonium
B = ethyl ac
methanol-water
monia (9+1);
chloroform-me
triethylamine
ethyl ketone-m
ane-methanol
anol (8+10+2)
^b Br. = Brink
Eastman Chrom
Indicator.

The spectra
bases are cle
from those o
logs (4).

Gas-Liquid

The resul
presented in
analog em
ably from t

6.1, 1976)

Table 3. Retention times (min) of some phencyclidine analogs for column packings and oven temperatures indicated*

Compound	3% OV-17			2.5% OV-225		5% OV-7		3% SE-30		
	225°C	200°C	175°C	175°C	150°C	200°C	175°C	150°C	125°C	100°C
I	4.0	7.4	19.5	5.2	14.6	2.8	5.0	13.5	4.3	16.5
II	4.0	8.0	19.3	5.2	14.9	1.1	1.5	2.6	1.0	2.0
III	5.2	10.7	30.3	10.6	32.0	1.1	1.5	2.6	1.1	2.0
IV	3.2	5.9	14.4	4.4	11.1	1.1	1.5	2.6	1.1	2.0
V	5.2	10.7	30.4	10.5	31.2	3.2	6.2	18.4	6.0	23.5
VI	3.2	5.9	14.5	4.3	11.2	2.2	3.9	10.0	3.3	11.8

* Columns were glass, 6' long, injector 275°C, nitrogen flow 30 ml/min. Support material was 80-100 mesh Chromosorb W (HP), excepting 3% SE-30 when 60-80 mesh was used.

Table 4. R_f values ($\times 100$) of some analogs of phencyclidine

System ^a	Plate ^b	I	II	III	IV	V	VI	LSD
A	Br.	82	85	79	72	78	58	21
A	Ea.	78	79	77	77	78	71	44
B	Br.	85	86	81	77	82	70	40
B	Ea.	70	75	75	68	73	67	60
C	Br.	82	85	84	76	84	68	76
C	Ea.	76	74	75	72	71	60	62
D	Br.	84	89	86	73	85	58	52
D	Ea.	79	81	77	76	79	63	62
E	Br.	12	30	79	10	73	12	38
E	Ea.	57	70	82	65	81	43	66
F	Br.	85	87	82	80	82	78	40
F	Ea.	74	71	75	71	75	74	58
G	Br.	24	42	65	26	63	12	21
G	Ea.	62	68	70	58	69	47	45
H	Br.	20	53	82	82	77	20	33
H	Ea.	64	72	74	64	75	55	51
I	Br.	4	8	60	8	42	5	22
I	Ea.	20	34	72	34	69	22	41

^a A = ethyl acetate-cyclohexane-dioxane-methanol-water-ammonium hydroxide (50+50+10+10+1.5+0.5); B = ethyl acetate-cyclohexane-ammonium hydroxide-methanol-water (70+15+2+8+0.5); C = ethanol-5N ammonia (9+1); D = acetone-12N ammonia (99+1); E = chloroform-methanol (9+1); F = chloroform-acetone-triethylamine (3+4+1); G = 1,1,1-trichloroethane-methyl ethyl ketone-methanol (7+2+1); H = n-hexane-nitromethane-methanol (7+13+5); I = toluene-nitromethane-methanol (9+10+2).

^b Br. = Brinkmann silica gel G glass plates and Ea. = Eastman Chromagram 6080 gel sheets with a fluorescent indicator.

The spectra of the hydrochlorides and the free bases are clearly distinct from one another and from those of closely related isomers and homologs (4).

Gas-Liquid Chromatography of Bases

The results obtained with several phases are presented in Table 3. The thiophene-containing analogs emerged more quickly and indistinguishably from the OV-7 and SE-30 columns than did

the phenyl containing analogs, whereas with the OV-17 and OV-225 columns the phenyl and thienyl analogs were indistinguishable and the order of emergence was determined by the nitrogen heterocycle. The chromatograms, particularly of the thienyl analogs on OV-7 and SE-30 columns, are indicative of decomposition, and the nature of the emergent peaks remains to be investigated.

Thin Layer Chromatography

Nine systems were investigated (Table 4). LSD and phencyclidine were included in order to compare data for these compounds. The use of iodoplatinate spray reagent showed the compounds as magenta-colored spots.

Acknowledgments

The assistance of J. C. Ethier, D. Legault, and B. Lodge is appreciated.

REFERENCES

- (1) Reed, A., & Kane, A. W. (1970) *Capsules 2*, 1-2
- (2) National Clearinghouse for Drug Abuse Information: *Phencyclidine (PCP)* (1973) Series 14, No. 1, 1-11
- (3) Marshman, J. A. (Ed.) (1974) *Street Drug Analysis and Its Social and Clinical Implications*, Addiction Research Foundation of Ontario, Canada, p. 26
- (4) Dobres, H. (1975) *Microgram* 8, 27-30
- (5) Beckstead, H. D., & French, W. N. (1971) *Some Analytical Methods for Drugs Subject to Abuse*, Department of National Health and Welfare, Ottawa, Canada
- (6) Budzikiewicz, H., Djerassi, C., & Williams, D. H. (1964) *Interpretation of Mass Spectra of Organic Compounds*, Holden-Day, Inc., San Francisco, pp. 162-165 and 231-235

DEA LABORATORY NOTES

10



DATE

NO.

DRUG TYPE

Phencyclidine

METHODOLOGY

Gas Liquid Chromatography

GAS LIQUID CHROMATOGRAPHIC SCREENING PROCEDURE FOR COMPONENTS FOUND IN CLANDESTINE PHENCYCLIDINE REACTION MIXTURES

BY

JACK K. RANEY
GERALD T. SKOWRONSKI
AND
RONALD J. WAGENHOFER
FORENSIC CHEMISTS,
DRUG ENFORCEMENT ADMINISTRATION
NORTH CENTRAL REGIONAL LABORATORY
CHICAGO, ILLINOIS

INTRODUCTION

Recently there have been numerous clandestine phencyclidine laboratories found operating throughout the country. Many of the samples originating from these laboratories are contaminated with starting materials, reaction intermediates and by-products. It is the scope of this paper to present a GLC temperature program that would enable the chemist to screen for these components. Since most of the PCP clandestine laboratories utilize the "carbonitrile intermediate" method, the data in this paper will reflect this synthetic route.

APPARATUS

Instrumentation: Hewlett-Packard 5840A Gas Chromatograph with a Flame Ionization Detector

Columns: 1) 6 ft. x 1/8 in. I.D. 3% OV-17 on 100/120 gas-chrom Q.
2) 6 ft. x 1/8 in. I.D. 10% OV-101 on 100/120 gas-chrom Q.

DRUG ENFORCEMENT ADMINISTRATION / U. S. DEPARTMENT OF JUSTICE

Parameters:

Carrier Gas: Nitrogen @ approximately 60 ml/min.
Injector Temperature: 275°
Detector Temperature: 300°
Column Temperature: Initial temperature of 90°c for
3 min. with a rate of 10°c/min
for 16 min. Final temperature
of 250°c for 2 min.

PROCEDURAL STANDARD

Starting Reagents: Cyclohexanone; Piperidine and Bromobenzene
Reaction Intermediate: 1-Piperidinocyclohexane Carbonitrile
(1-PCC)

By-Products: 1-Cyclohexylpiperidine; Biphenyl; 1-Phenylcyclo-
hexene;* 1-Phenylcyclohexanol; Phenol and 1-Phenyl-
ethanol.

Final Product: Phencyclidine (PCP) (See table 1 for structures)

The procedural standard is composed of the 11 individual compo-
nents listed above with each being diluted to a concentration of
approximately 1 mg/ml in methanol.

Although listed as a byproduct, 1-phenylcyclohexene may be partially
caused by the thermal degradation of PCP. 1,2

PROCEDURE

A portion of the clandestine sample is dissolved in methanol and the
solution is then injected into the GC using the above stated conditions.
A chromatogram such as figure 1 is obtained and the retention times are
compared to those of a procedural standard run under identical condi-
tions (figure 2).

DISCUSSION

Comparison of the chromatograms in figures 1 and 2 indicate the
possible presence of six or seven of the procedural standard compo-
nents in the clandestine laboratory sample. The variations in reten-
tion times between the two chromatograms for four of these compounds
(1-PCC, 1-phenylcyclohexene, 1-phenylcyclohexanol and PCP) fall with-
in acceptable limits. Although the first two compounds eluting after
the solvent front do not match the retention times of the procedural
standard exactly, they have been identified as piperidine and cyclo-
hexanone by "spiking" the sample with the respective standard (see
figures 3 and 4). The retention times of the piperidine and cyclohexa-
ne in the spiked sample likewise changed; this suggests that the re-
tention times of these components can vary somewhat according to their
concentrations in a sample.

Many of the smaller peaks in the sample have been ignored as being due to minor contaminants. Two of the peaks shown in the procedural standard, phenol and 1-phenylethanol, are rarely encountered in clandestine samples. They appear to be the results of an excessive use of phenyl magnesium bromide.

Since the 3% OV-17 column does not separate cyclohexanone and bromobenzene well, a second column, 10% OV-101, has been used to resolve these peaks (figures 5 and 6). In addition a better separation between 1-phenylcyclohexene and biphenyl is effected on the OV-101 column. Figure 7 is a clandestine laboratory sample "spiked" with piperidine and injected on the 10% OV-101 column. Again, as in figure 3, the "spiked" sample results in one peak which is indicative of piperidine. Although the OV-101 column separates cyclohexanone and bromobenzene, it now has the disadvantage that 1-PCC and biphenyl have the same retention time. We suggest the use of the OV-17 column as the primary column since it is often important to identify the presence of the reaction intermediate, 1-PCC.

The identity of each of the individual components in standard and sample was confirmed by GC-MS, and a subsequent paper will list the normalized mass spectra.

This two-column separation technique has been utilized in the North Central Regional Laboratory for the last year with very favorable results. The authors wish to express their gratitude to all members of the laboratory for suggestions as well as contributions of data obtained from the analysis of clandestine PCP samples.

REFERENCES


1. "Identification of Degradation Peak in Phencyclidine Gas Chromatography", by David Hauber, Microgram, Vol. VIII, No. 7, 1975, p. 100.
2. "Metabolites of Phencyclidine", by L.K. Wong and K. Biemann, Clinical Toxicology, 9 (4), 1976, pp. 583-591.

TABLE 1


COMPOUNDS FOUND IN PCP REACTION MIXTURES




1-CYCLOHEXYLPYPERIDINE



1-PHENYLCYCLOHEXANOL




BIPHENYL




1-PHENYLCYCLOHEXENE




PHENOL



1-PHENYLETHANOL



PCP



PCC




FIGURE 1 PCP CLANDESTINE LABORATORY SAMPLE (3% OV-17 Column)

RETENTION TIME IN MINUTES	COMPONENT PEAK
1.08	Piperidine
2.39	Cyclohexanone and/or Bromobenzene
9.61	1-PCC
10.47	1-Phenylcyclohexene
12.07	1-Phenylcyclohexanol
16.46	PCP




FIGURE 2 PROCEDURAL STANDARD (3% OV-17 Column)

RETENTION TIME IN MINUTES	COMPONENT PEAK
0.96	Piperidine
2.43	Cyclohexanone and Bromobenzene
3.56	Phenol
5.41	1-Phenylethanol
8.23	1-Cyclohexylpiperidine
9.62	1-PCC
10.49	1-Phenylcyclohexene
10.61	Biphenyl
12.07	1-Phenylcyclohexanol
16.48	PCP




FIGURE 3 PCP CLANDESTINE LABORATORY SAMPLE SPIKED WITH PIPERIDINE (3% OV-17 Column)

RETENTION TIME IN MINUTES	COMPONENT PEAK
0.99	Piperidine
2.42	Cyclohexanone and/or Bromobenzene
9.62	1-PCC
10.48	1-Phenylcyclohexene
12.08	1-Phenylcyclohexanol
16.47	PCP




FIGURE 4 PCP CLANDESTINE LABORATORY SAMPLE SPIKED WITH CYCLOHEXANONE (3% OV=17 Column)

RETENTION TIME IN MINUTES	COMPONENT PEAK
1.11	Piperidine
2.36	Cyclohexanone and/or Bromobenzene
9.63	1-PCC
10.49	1-Phenylcyclohexene
12.09	1-Phenylcyclohexanol
16.47	PCP




FIGURE 5 CLANDESTINE LABORATORY SAMPLE (10% OV-101 Column)

RETENTION TIME IN MINUTES	COMPONENT TIME
2.85	Piperidine
4.96	Cyclohexanone
13.84	1-PCC and/or Biphenyl
14.15	1-Phenylcyclohexene
15.27	1-Phenylcyclohexanol
20.48	PCP




FIGURE 6 PROCEDURAL STANDARD (10% OV-101 Column)

RETENTION TIME IN MINUTES	COMPONENT PEAK
2.76	Piperidine
4.93	Cyclohexanone
6.04	Bromobenzene
6.71	Phenol
8.44	1-Phenylethanol
13.04	1-Cyclohexylpiperidine
13.82	1-PCC and Biphenyl
14.16	1-Phenylcyclohexene
15.27	1-Phenylcyclohexanol
20.50	PCP




FIGURE 7 CLANDESTINE LABORATORY SAMPLE SPIKED WITH PIPERIDINE (10% OV-101 Column)

RETENTION TIME IN MINUTES	COMPONENTS PEAK
2.83	Piperidine
4.97	Cyclohexanone
13.88	1-PCC and/or Biphenyl
14.19	1-Phenylcyclohexene
15.30	1-Phenylcyclohexanol
20.52	PCP

(11)

Identification of Degradation Peak in Phencyclidine Gas Chromatography



David Hauber, Kentucky State Police Crime Laboratory
Frankfort, Kentucky

An extra peak in gas chromatograms of phencyclidine has been a constant occurrence in our lab with any PCP sample run at lower column temperatures. On our 5% SE-30 column (6 ft., M80-100), elution of PCP in 6 min. at constant temperature (250°C) will result in an extra peak at approximately 1.6 min. Elution of PCP in 5 min. on our 3% OV-1 column (8 ft., M80-100) at 200°C results in an extra peak at 1.8 min. The mass spectrum of this extra peak eluting from an OV-1 column (see figure) suggests phenylcyclohexene - a logical PCP degradation product resulting from the elimination of piperidine, which would be buried in the solvent elution. There is less degradation in systems with glass columns and glass-lined injection ports. Hot catalytic surfaces of injection ports (approximately 265°C) cause most of the degradation, but an elevated base line between the GC peaks suggests more degradation during the journey through the column. Some evidence exists for more degradation of the hydrochloride than the free base in certain glass systems.

This extra peak gives one more identifying mark in PCP gas chromatography, but the degradation precludes easy quantitation in many GC systems since the relative amounts of this degradation varies.

Computer normalized mass spectrum of the PCP degradation peak obtained from a Finnigan GCMS (Model 3100D).

Suggested structure:



STATE OF MICHIGAN



WILLIAM G. MILLIKEN, GOVERNOR

COL. GEORGE L. HALVERSON, DIRECTOR, DEPARTMENT OF STATE POLICE
MAURICE S. REIZEN, M.D., DIRECTOR, DEPARTMENT OF PUBLIC HEALTH

REGIONAL CRIME DETECTION LABORATORY

XXXXXXXXXXXXXXXXXXXXXXXXXXXX
XXXXXXXXXXXXXXXXXXXXXXXXXXXX
XXXXXXXXXXXXXXXXXXXXXXXXXXXX

30303 Stephenson Hwy., Madison Hgts, Mich.
585-7521

12

1 - PIPERIDINOCYCLOHEXANE CARBONITRILE

A PHENCYCLIDINE PRECURSOR

by

John H. Siefert, Ph.D.
Crime Laboratory Scientist
Michigan Department of Public Health
Madison Heights Crime Detection Laboratory
Madison Heights, Michigan 48071

Our laboratory has recently received numerous submissions of phencyclidine (PCP) containing an impurity. GC/MS data of this impurity suggested that it might be 1 - Piperidinocyclohexane Carbonitrile (PCC), a precursor in the synthesis of PCP^{1,2}.

Concurrently we also received a fairly pure substance which was later proven to be PCC and exhibited the same GC/MS data as the impurity encountered in the above PCP submissions.

This substance was identified as PCC via its IR, UV, GC, GC/MS and conversion to PCP.

Discussion

The IR spectrum of the pure substance was identical to that of PCC³ and showed a nitrile band at 2210 cm⁻¹ as well as major bands at 2975, 1447, 1107, 872 and 787 cm⁻¹.

As expected from its chemical structure, PCC did not show significant absorption in the UV region of 350-200 nm.

A GC/MS and a solid probe MS of PCC differed significantly due to the apparent loss of HCN (27 amu) on the GC column giving N-(1-cyclohexenyl) piperidine.

Gas chromatographic data of PCC using three different liquid phases at 170° are listed in Table 1.



Conversion of PCC to PCP

Phenylmagnesium bromide was prepared in the usual manner⁴ by placing 1.26 g of magnesium turnings, 5 ml of anhydrous diethyl ether and a crystal of iodine in a 100 ml round bottom flask, then adding 8.18 g (5.5 ml) of bromobenzene in 17 ml of anhydrous ether. The mixture was then refluxed until the Mg had completely dissolved. To the ethereal solution of phenylmagnesium bromide was added 8.0 g of PCC in 50 ml of cyclohexane. After refluxing the mixture for 1 hr, the ether was distilled off and the remaining solution was extracted with 25 ml of 2.8 N HCl. The aqueous layer was separated and evaporated to yield a crude crystalline product. Subsequent purification and analysis of this product showed it to be PCP, thus confirming that the starting material was indeed PCC.

Table 1: Gas Chromatographic data using:

a) 3% OV-1, 6', 170°C

	$\frac{R}{t}$
PCC	0.63 min
1-phenylcyclohexene	0.73 min
PCP	3.83 min

b) 3% OV-17, 6', 170°C

PCC	0.63 min
1-Phenylcyclohexene	0.78 min
PCP	4.40 min

c) 3% Ultraphase, 5', 170°C

PCC	0.85 min
1-Phenylcyclohexene	0.92 min
PCP	4.50 min

References

- 1) V.H. Maddox, E.F. Godefroi, R.F. Parcell;
J. Med. Chem., 8, 230 (1965)
- 2) A.T. Shulgin, D.E. MacLean,
Clinical Toxicology, 9, 553 (1976)
- 3) R.D. Porter, Microgram, 8 (10), 151 (1975)
- 4) A.I. Vogel; A Textbook of Practical Organic Chemistry, Longman Group Limited, London, 1972, D756





FIG 1. 1-PIPERIDINOCYCLOHEXANE CARBONITRILE

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

DEA LABORATORY NOTES

- 27 -



DATE

NO.

DRUG TYPE Hallucinogen

METHODOLOGY U.V., I.R., and G.C-M.S Identification

Identification of the PCP analogs,

- (a) 1-(1-phenylcyclohexyl) pyrrolidine, and
- (b) 1-(1-phenylcyclopentyl) piperadine.

Howard Dobres, Forensic Chemist
Mid-Atlantic Regional Laboratory

Compound (a) was encountered in a clandestine PCP laboratory in January, 1974 in Rockville, Maryland. It is believed to have been synthesized by substituting pyrrolidine for piperadine in the PCP synthesis. Compound (b) was alleged to have been synthesized, but not found, in the same laboratory.

EXPERIMENTAL:

Standard samples of 1-(1-phenylcyclohexyl) pyrrolidine monohydrochloride, Lot P, and 1-(1-phenylcyclopentyl) piperadine, Lot P, were supplied by Park-Davis & Company, Research and Development Division.

Samples were recrystallized from chloroform-methanol (plus a drop of conc. HCl) for the I.R. spectra.

Instruments used were the Beckman Acta-5 U.V. spectrophotometer, Perkin-Elmer 457 U.R. spectrophotometer, and Finigan 3000 G.C-M.S with a 4ft. 3% Q.V-1 column.

Operating conditions:

- U.V.: 1-(1-phenylcyclohexyl) pyrrolidine hydrochloride-1.6 mg/3ml 0.1N HCl
- 1-(1-phenylcyclopentyl) piperadine base-1.2mg/3ml.
- I.R.: 2% 1-(1-phenylcyclohexyl) pyrrolidine hydrochloride, and 1-(1-phenylcyclopentyl) piperadine hydrochloride-KBr pellets.
- Scan speed - Fast
- Slit - 7

(See attached spectra)

DRUG ENFORCEMENT ADMINISTRATION / U. S. DEPARTMENT OF JUSTICE


G.C-M.S - 1 mcg. of each compound.
Column temperature = 200°C
Carrier gas flow rate = 28 cc/min.
Sensitivity 10⁻⁶

IDENTIFICATION:


- (a) 1-(1-phenylcyclohexyl) pyrrolidine
U.V.: max. 262nm, 257nm, 268nm, (0.1NHC1)
I.R. principal peaks: 1447, 700, 1353, 1040, 1031,
1305, 920, 769 cm⁻¹
G.C.-M.S - Parent peak - M/e = 229
Base peak - M/e = 186
- (b) 1-(1-phenylcyclopentyl) piperadine
U.V.: max. 261.5nm, 256.5nm, 268 nm
I.R. principal peaks: 707, 1453, 749, 1249, 1023,
765, 975, 1338 cm⁻¹
G.C. - M.S - Parent peak - M/e = 229
Base peak - M/e = 200

Property of Idaho State Police Forensic Services
Uncontrolled Internal Copy
OBSOLETE DOCUMENT

1-(1-phenylcyclopentyl) piperadine
base 1.2 mg/3 ml




1-(1-phenylcyclohexyl) pyrrolidine
hydrochloride 1.6 mg/3 ml




Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

1-(1-phenylcyclohexyl) piperidine HCl
(redissolved CHCl₃-MeOH)



SAMPLE 1-(1-phenylcyclohexyl) piperidine HCl (redissolved from CHCl ₃ -MeOH)	SOLVENT DMF CONCENTRATION 2% CELL PATH REFERENCE	WAVENUMBER (CM ⁻¹) 3500 3000 2500 2000 1500 1000 500	TRANSMITTANCE (%) 0 10 20 30 40 50 60 70 80 90 100
OPERATOR J. J. [unclear] DATE 5-23-74 INSTRUMENT 600 MA PILOT NO. 55-1-507	REMARKS		

1-(1-phenylcyclopentyl) piperidine HCl
(recrystallized from CHCl₃-MeOH-HCl)



SAMPLE 1-(1-phenylcyclopentyl) piperidine HCl (recrystallized from CHCl ₃ -MeOH-HCl)	SOLVENT DMF CONCENTRATION 2% CELL PATH REFERENCE	WAVENUMBER (CM ⁻¹) 3500 3000 2500 2000 1500 1000 500	TRANSMITTANCE (%) 0 10 20 30 40 50 60 70 80 90 100
OPERATOR J. J. [unclear] DATE 5-23-74 INSTRUMENT 600 MA PILOT NO. 55-1-507	REMARKS		

Property of Idaho State Police Forensic Services
UNCONTROLLED INTERNET COPY SERVICE
OBSOLETE DOCUMENT

BND LABORATORY NOTES

DATE November 22, 1972 -122-
NO. 43
DRUG TYPE Hallucinogen
METHODOLOGY I. R. Identification

14

THIOPHENE ANALOG OF PCP

by

James A. Heagy
Forensic Chemist
San Francisco Regional Laboratory
Bureau of Narcotics and Dangerous Drugs
San Francisco, California 94102

A new drug has recently appeared on the illicit market in several parts of the country. It has been identified as 1-[1-(2-thienyl) cyclohexyl] piperidine, is synthesized from 2-bromothiophene and piperidinocyclohexanecarbonitrile (CA 54:12159c) in a manner similar to the synthesis of phencyclidine. All samples encountered so far contain several amine impurities including, in some cases, large amounts of piperidine. A strong cyanide-like odor is also generally present.

Identification characteristics

U. V. max. 232.5 nm, min. 204 nm. (in 0.1 HCl)

I. R. peaks at 708, 1444, 1453, 1251, 1292, 854, 846 cm^{-1} .

(See attachment)

Marquis reagent: gas expelled, turns to gray-orange.

Mecke reagent: gas expelled, color turns to yellow-green, then blue-green, and slowly to deep blue.

Mass spectrum: Molecular ion is not present.

MICRONS

5.0

4.0

3.0

1800

2000

2500

3000


3500


WAVENUMBER (CM⁻¹)

1-[(2-thienyl)cyclohexyl] piperidine HCl
Parke-Davis, Lot No. Q
in KBr

James A. Heagy
October 19, 1972

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT





DEA LABORATORY NOTES

- 96 -

(15)

DATE OCT 23 1974

NO. 88

DRUG TYPE HALLUCINOGEN

METHODOLOGY UV SPECTRA OF PHENCYCLIDINE BASE

BY


ROGER G. FUELSTER
FORENSIC CHEMIST, NORTH CENTRAL REGIONAL LABORATORY
DRUG ENFORCEMENT ADMINISTRATION

A recent sample received in the North Central Regional Laboratory consisted of rock-candy like light yellow chunks of material. Screening with acid CoSCN produced the pale blue of PCP. The material was then dissolved in Spectral Grade MeOH and a UV scan obtained. The scan was featureless (see attached scans A & C). Further analysis of the material by GC - MS and IR showed nearly pure PCP base.

The UV scans were then re-run in acidified spectral grade MeOH (see attached scans B&D), with the expected benzenoid curve of PCP.

It is apparent from these results that care must be exercised when using UV as a screening technique for unknowns. Curves must be run in acid as well as neutral solution. Although PCP is the only example of this behavior encountered to date, other organic bases may exhibit similar aberrations.

DRUG ENFORCEMENT ADMINISTRATION / U. S. DEPARTMENT OF JUSTICE




ⓑ
 SAME AS
 "A" WITH
 1 DROP HCl
 (1+1)

240 260 280 300 320 340 360

Ⓐ

PCP BASE
 ca 5mg/100ml
 SPECTRAL GRADE
 MeOH

240 260 280 300 320 340 360



ⓓ
 SAME AS "C"
 WITH 1 DROP
 HCl (1+1).

ⓐ
 PCP BASE
 ca 50mg/100ml
 SPECTRAL GRADE
 MeOH

240 260 280 300 320 340

(16)

Compound: PHENCYCLIDINE HYDROCHLORIDE
Mol. Formula: $C_{17}H_{25}N \cdot HCl$
Mol. Weight: 279.86



UV DATA:

UV Maxima: 262, 252, 257, 268 m μ

Concentration: 52.4 mgs per 100.0 mls

Abs²⁶² : .628 + .012 = .640

E_{1%}^{1cm} : 12.2



Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

LABORATORY NOTES

(17)

DATE

-127-

NO.

DRUG TYPE Hallucinogen

METHODOLOGY

The Identification of N-ethyl-1-phenylcyclohexylamine hydrochloride
(Cyclohexamine) P.C.C.

by

R. P. Barron
Chemist

Special Testing and Research Laboratory
DRUG ENFORCEMENT ADMINISTRATION

This communication reports data utilized in the structural elucidation and confirmation of a recently encountered drug, N-ethyl-1-phenylcyclohexamine hydrochloride (Figure 1) on the illicit market. The compound is described in U.S. Patent 3,097,136 issued to Parke, Davis and Company¹ and included in "Psychotropic Drugs and Related Compounds"² under the names CI 440 and Cyclohexamine.



FIGURE 1

Identification Characteristics

UV: maxima at 236, 257, 262 and 268 millimicrons

IR: absorption peaks at 2940, 1582, 1463, 780, 701
and 624 cm^{-1} by KBr. See Figure 2 for spectrum.

MS: Molecular ion observed at m/e 203 and base peak at
 m/e 160. See Table 1 for spectrum.

NMR: See Figure 3 for spectrum.


AN OBSERVATION ON THE ULTRAVIOLET ABSORPTION CURVE OF MIXTURES OF LSD AND PCP


Otto Schales
Elyria Memorial Hospital
Elyria, Ohio 44035

Submitted for examination were several tablets, light tan-pink in color, biconvex, diameter 6.4 mm, thickness 3.2 mm at center, average weight 160 mg.

One of the tablets was powdered and extracted three times with 10 ml chloroform each. The combined extracts were filtered and evaporated to dryness at 25° C. by blowing air against the surface of the liquid in a porcelain dish. The residue was taken up in 3 ml 0.1 N sulfuric acid and the UV spectrum was recorded with a Zeiss DMR 21. In addition to the LSD maximum at 310 nm, there were two peaks at 268 and 262 nm. The solution was exposed to long wave UV light for 4 min. and the UV scan was repeated. As expected, there was a lowering of the LSD maximum and a shift toward 292 nm. Unexpectedly, however, there occurred also a marked decrease in the LSD absorbance in the shorter wavelength UV range, resulting in the emergence of a complete PCP absorbance curve. This curve was brought out even better after an additional 4 min. UV exposure.


The spectral behavior of LSD after UV radiation offers a simple means of demonstrating the presence of PCP in preparations containing LSD.






Phencyclidine HCl
HAuCh 52 Agueous soln. Cg.
(HDAI)


white-gray to
light yellow




Cocaine HCl
HAuCh 52 Agueous Cg.



Cocaine HCl
HAuCh 52 Agueous Cg.
(19)




Phencyclidine HCl - Cocaine HCl
HAuCh 52 Agueous soln Cg.
(HDAI)




Phencyclidine HCl
HAuBr - 2H₂PO₄ - (203)H₂SO₄ Cg.

yellow orange brown




Phencyclidine HCl
HAuBr - HCl - (273)H₂SO₄

pyramidal
yellow to



Phencyclidine HCl
HAuCh - (101)H₂SO₄

clear,
to - clear to
multi colored.



Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

DEA LABORATORY NOTES



DATE

NO.

DRUG TYPE, Thiophene Analog of Phencyclidine


METHODOLOGY Microcrystalline, GC, TLC, UV, IR, X-RAY

PART I: ANALYTICAL DATA ON THE THIOPHENE ANALOG OF PHENCYCLIDINE

by Jose J. Alvarez

INTRODUCTION

Infrequently encountered controlled substances, especially those which are sporadically reported in the chemical literature, are difficult to analyze and identify routinely. If the usual color tests are employed and the concentration of the controlled substance is low, the controlled substance may be overlooked. An example of such a controlled substance is the thiophene analog of phencyclidine, 1-[1-(2-thienyl)cyclohexyl] piperidine (TCP):



The following presentation on TCP was initiated to report the information in the literature and to report the results obtained during the laboratory investigations. Due to the length of one comprehensive paper, the study will be presented in two parts. This paper will cover microcrystalline tests, background, and various instrumental and chemical analyses; part two will cover the mass spectrometric analyses.

BACKGROUND

TCP is prepared in a manner similar to phencyclidine, PCP (1-5). It is synthesized from 2-thienylmagnesium bromide (prepared by reaction of 163 g 2-bromothiophene and 24 g Mg in 500 ml EtOH) and 1-piperidinocyclohexanecarbonitrile (4). The latter compound may be prepared as described by Kalir, et al. (5).

MICROGRAM, Vol. X, No. 9 (September, 1977)

It appears that the synthesis was initiated to study its pharmacological activity, as related to PCP, which is used as a veterinary tranquilizer. Various papers have appeared in the literature in regards to its pharmacological activity. For brevity, some Chem. Abstr. references are cited (6). Kalir, et. al. (5) studied a series of N,N-substituted 1-arylcyclohexylamines and found TCP to be the most active in tests conducted for psychopharmacological activity. Also, to be noted is the presence of a positional isomer of TCP as reported by Mousseron, et. al. (7): 1-[1-(3-thienyl)cyclohexyl] piperidine.

TCP was added to the Controlled Substances Act in 1975 - as a Schedule One drug (8). It has been encountered on several occasions; Microgram has reported its appearance throughout the country, as input was received from the various laboratories (9).

EQUIPMENT AND REAGENTS EMPLOYED

Microcrystalline tests were performed with two reagents: 10% solution of KI in distilled water and 10% solution of NH_4SCN in distilled water. Both of these reagents are quite sensitive to PCP and TCP. Both reagents are very stable and can be used for a long period of time. The photomicrographs designated Figure 1 through Figure 4 were taken with the Nikon Microflex Model AFM (automatic photomicrographic attachment) Microscope, using ordinary transmitted light and a magnification of 150x.

Gas chromatographic data was obtained with the Perkin-Elmer 3920 Gas Chromatograph. Parameters are as follows: carrier gas - nitrogen @ 55 ml/min, injector and detector temperatures - 250 degrees C, chart speed - 1/2 in/min, GC interfaced with a Perkin-Elmer M-1 Computing Integrator, attenuation - x32, range - x10. Two columns were employed. Both were 6' x 1/4" glass columns; packing as follows: 3% OV-1 on 100/120 Gas Chrom Q and 10% OV-101 on 100/120 Chromasorb WHP. The oven temperatures employed for the columns were 180 and 270 degrees C, respectively.

Ultraviolet spectra were obtained with the Cary 15 UV Spectrophotometer (consult Figure 6 for absorbance and concentration data).

Infrared spectra (Figures 7-10) were obtained with the Perkin-Elmer 467 Grating Infrared Spectrophotometer.

The X-ray diffraction pattern was obtained with the Siemens X-Ray Diffraction System, which is equipped with a scintillation counter, using nickel-filtered copper radiation ($\text{Cu K}\alpha_1 = 1.540$ Angstroms). Recordings were obtained on a precision line recorder scanning at 1 cm/min.

The TCP standard was obtained from Parke, Davis & Company.

SOLUBILITY AND MELTING POINT

- A. As the HCl salt - soluble in water, chloroform, methanol/ethanol
- B. Melting points, as the HCl salt (as reported by the literature):
- 1) 233 - 236 degrees C with transition at 182-3 degrees C (3,4)
 - 2) 230 - 235 degrees C (5)
 - 3) 200 - 203 degrees C subl. (2)

CHEMICAL COLOR TESTS

Several reagents were tested for chemical color sensitivity, but few produced a discernible color. Other possible tests were evaluated, but few were of total value because of the time and equipment involved (10).

The two color tests, used extensively for opium alkaloids/derivatives, are the preferential tests - the Marquis (10 drops formaldehyde solution in 10 ml conc. sulfuric acid) and the Mecke (0.25 g selenious acid in 25 ml conc. sulfuric acid), as reported by Heagy (1). Of the two, the Mecke's reagent is the most sensitive.

- A. Marquis - effervescence, turns to a gray-orange and slowly to gray-green.
- B. Mecke - effervescence, turns to a yellow-green, slowly to a yellow-blue/green, and then to a deep blue.

MICROCRYSTALLINE TESTS

Grind and mix a portion of the sample well; on a microscope slide place a small drop of distilled/deionized water. Dissolve a minute portion of the sample in the water and place a small drop of the reagent adjacent to the sample solution. Run the solution and the reagent drops together with a spatula. For best results, do not stir the mixture.

The crystal formation will occur within a minute; view the crystals through a microscope, preferably a polarizing microscope at 125x. As with all other microcrystalline tests, it is best to run a test with a standard right before or after the sample run. This should be done until enough familiarization is obtained with the particular crystal formations.

Figures 1 and 2 are the crystal formations of PCP and TCP, respectively, with NH_4SCN . Figures 3 and 4 are those of TCP and PCP, respectively, with KI. Of the two reagents, KI gave the more consistent crystal formation with the two compounds. With NH_4SCN , the PCP crystal formation was also consistent. Successful results have been obtained with PCP and TCP samples which were quantitated at three to five percent.

With NH_4SCN , the crystal formation of TCP was not consistent - the concentration is more critical. Aside from the crystal formation in Figure 2, the formation of long slender needles is prevalent. These needles are similar to Fulton's large, coarse, Class I needles (11). To enhance the proper microcrystal formation in low percentage TCP samples, a preliminary CHCl_3 dissolution was performed. This was done by dissolving a few milligrams of sample in CHCl_3 , filtering the insolubles, and evaporating the CHCl_3 to dryness. The residue was then used for the microcrystalline work. This preliminary cleanup improved the consistency of the crystal formation.

GAS CHROMATOGRAPHY

The results of the gas chromatographic investigations are represented in Figure 5. The quantitation of PCP can be achieved successfully with an OV-1 column. There is noticeable decomposition, but negligible as compared to the OV-101. TCP is a compound difficult to quantitate by GC because of instability to the phases. Major decomposition occurs in both columns. These decomposition products will be discussed in part two of these TCP papers.

Bailey, et. al. investigated the employment of SE-30 and OV-17, 225, and 7 columns (2). They report decomposition in all four, particularly on OV-7 and SE-30.

THIN-LAYER CHROMATOGRAPHY

TLC investigations were not initiated in this laboratory.

The major work by Bailey, et. al. reports the investigations of nine solvent systems (2). Shulgin and Helisten have also investigated aspects of TLC analyses (12). Their work covered a differentiation of PCP, TCP, and the precursor, 1-piperidinocyclohexanecarbonitrile, using a benzene, acetone, and pyridine (16:8:1) solvent system.

ULTRAVIOLET SPECTRA

Figure 6 is the ultraviolet spectrum of TCP in the respective solvents. It is a strong absorber with a maxima at 232.5 nm (EtOH, 0.1N HCl). This is in good agreement with the works of Heagy (1) and Bailey, et. al. (2).

INFRARED SPECTRA

Figure 7 is the infrared spectrum of the Parke-Davis TCP hydrochloride standard. Figures 8-10 are spectra of a sample received at the laboratory.

Figure 8 represents the spectrum of the isolated TCP, converted to the HCl salt. It was isolated by dissolving the sample in CHCl_3 and filtering of the insoluble material. A 0.1N HCl extraction was performed on the CHCl_3 solution. The CHCl_3 layer was again filtered and CHCl_3 evaporated to approximately two ml. A mini-column of alumina was prepared and the two ml. were eluted with an additional ten ml. CHCl_3 . The free base was converted to the HCl salt and pressed into a KBr pellet.

Figure 9 is the spectrum of the same TCP sample with a preliminary cleanup. The sample was dissolved in distilled water, made alkaline with sodium bicarbonate, and extracted with CHCl_3 . The CHCl_3 extract was filtered and approx. five ml. of a 1% HCl solution of MeOH added; the mixture was taken to dryness. The HCl salt in a KBr pellet was scanned. The "impurity"/TCP infrared spectrum was obtained.

Figure 10 is the spectrum of the "impurity" in the sample. The 0.1N HCl solution or layer (obtained in the extractions to isolate the TCP for the spectrum of Figure 8) was made alkaline and extracted with CHCl_3 . The HCl salt was prepared in the same manner as that for the spectrum of Figure 9. The infrared spectrum identifies piperidine as the "impurity".

X-RAY

Figure 11 depicts the diffraction pattern of TCP hydrochloride. It is presented in order of increasing 2θ degrees versus peak intensity, as is usually obtained on the recorder. TCP hydrochloride was scanned from 3 to 40 2θ degrees and that which is depicted in Figure 11 represents the area of interest for this compound. The 2θ degrees of interest are as follows: 10.3, 14.0, 14.8, 15.4, 19.2, 20.7, 21.7, 22.8, 23.5, 24.5, 27.0, and 28.1. This is in good agreement with the X-ray study conducted by Berens (13).

ACKNOWLEDGMENTS

I would like to thank Mr. Carlos L. Castillo and Mr. Wayne Moody for their assistance in the drawings and in the photomicrographic work, respectively. Also, my sincere thanks to Mr. Buddy Goldston and Mr. Charles Teer for their valuable assistance and constructive criticism during the preparation of the paper.

REFERENCES

- (1) Heagy, J. A., Microgram, V(10), 122(1972)
- (2) Bailey, K., Gagne, D. R., and Pike, R. K., JAOAC, 59(1), 81(1976)
- (3) Chem. Abst., 54, 12159c (Parcell, R. F., Parke, Davis & Co., U.S. 2,921,076)
- (4) Chem. Abst., 55, 2692c (Parke, Davis & Co., Brit. 838,748)
- (5) Kalir, A.; Eder, H.; Pelah, Z.; Balderman, D.; and Porath, G., J. Med. Chem., 12(3), 473(1969)
- (6) Chem. Abst.: 63 - 16959c, 79 - 73734h, 80 - 66595f, 81 - 114434x, and 81 - 85956m
- (7) Mousseron, M., et. al., Chim. Ther., 1968, 3(4), 241-7 (Chem. Abst.: 70,57767e)
- (8) Microgram, VIII(7), 101(1975)
- (9) Microgram: V(10), 105(1972); VI(2), 16(1973); VI(6), 80(1973); VI(11), 168(1973); VII(5), 52(1974); VII(10), 114(1974); VIII(5), 59, 60(1975); VIII(9), 126(1975); IX(1), 1(1976); and IX(4), 53(1976)
- (10) Feigl, F., Spot Tests In Organic Analysis, Elsevier Publishing Co., 1966, New York
- (11) Fulton, C.C., Modern Microcrystal Tests for Drugs, Wiley-Interscience, 1969, New York, p. 425 (1-02)
- (12) Shulgin, A. T. and Helisten, C., Microgram, VIII(10), 149(1975) and VIII(11), 171(1975)
- (13) Berens, L. J., personal communications




FIGURE 1: PCP/NH₄SCN 150x




FIGURE 2: TCP/NH₄SCN 150x




FIGURE 3: TCP/KI 150x




FIGURE 4: PCP/KI 150x






FIGURE 5: Chromatograms A and B were obtained with the 3% OV-1 column at 180 degrees C; C and D, with the 10% OV-101 at 270 degrees C. Numbers above the peaks denote the retention time in seconds.

FIGURE 6

THIOPHENE ANALOG OF PCP

Inst: Cary 15 UV
Abs: 0 - 1
Conc: 0.032 mg./ml.






Property of Idaho State Police Forensic Services
 Uncontrolled Internet Copy
 OBSOLETE DOCUMENT

FIGURE 7 (I of II)

THIOPHENE ANALOG OF PCP (TCP)

Inst: Perkin-Elmer 467
 Origin: Parke-Davis
 Salt: Hydrochloride (KBr pellet)
 Scan: Medium Speed
 Slit: Normal


WAVENUMBER (CM-1)



Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

FIGURE 7 (II of II)

TCP






FIGURE 11: X-RAY DIFFRACTION PATTERN OF TCP.HCl

BND LABORATORY NOTES

DATE

NO. 51

-22-

(21)

DRUG TYPE Mixture

METHODOLOGY IR

Identification of Mixtures of Phencyclidine Hydrochloride and Piperidine Hydrochloride

by

Victor A. Folen
Forensic Chemist
Special Testing and Research Laboratory

Several recent exhibits have been found to contain phencyclidine hydrochloride in combination with piperidine hydrochloride. In some of the mixtures phencyclidine was present both as the hydrochloride and as the free base. Piperidine was probably present as an unreacted precursor in the synthesis of phencyclidine. The material was either pale yellow or tan in color, and some of the mixtures were hygroscopic. Such samples gave no color with Marquis solution (conc. H_2SO_4 -formaldehyde 10:1) but there was a rapid effervescence. There was no effervescence, however, with dilute HCl, ruling out the possible presence of carbonate.

Crystal Tests

For the tentative identification of phencyclidine, one of the routine crystal tests may be used. Dissolve a very small amount of sample (ca. 0.1 mg.) in 6-10% acetic acid on a microscope slide. Touch the solution to an adjacent drop of gold chloride solution (5% aq.). There is an immediate dense formation of microscopic droplets, followed by the gradual development of typical crystals¹, rectangular, envelope-shaped, and of low birefringence. If the solution in acetic acid is too concentrated, atypical crystals will result, being in the form of sheaves, or a branched, serrate mosaic, both types of low birefringence. The crystal test for piperidine involves volatilization when treated with alkali.¹ Place a small amount of sample in a depression slide and add one-two drops of 10% NaOH. Place a microscope slide holding a drop of gold chloride solution (5% aq.) over the depression. Piperidine volatilizes quickly, with rapid crystallization in the drop of gold chloride solution of small, highly birefringent rods and clusters of three or four blades of unequal length.

Confirmation of the Presence of Phencyclidine and Piperidine

To separate the hydrochloride salt of phencyclidine from piperidine hydrochloride, advantage may be taken of their partitioning characteristics in water and chloroform.²

BUREAU OF NARCOTICS AND DANGEROUS DRUGS / U.S. DEPARTMENT OF JUSTICE

Dissolve about 200 mg. of sample in water in a separatory funnel. Extract phencyclidine hydrochloride three times using 10 ml. portions of chloroform, collecting the chloroform in a 50 ml. beaker; make the solution in the separatory funnel basic with Na_2CO_3 . Extract piperidine base using three 10 ml. portions of chloroform, again collecting in a 50 ml. beaker. Add concentrated HCl (5-10 drops) to the beaker containing piperidine base and evaporate the contents of both beakers to dryness on a water bath under a stream of air. Obtain infra-red spectra for the isolated components.

Notes


Attached is an infra-red spectrum for piperidine hydrochloride.

After evaporation of the chloroform from the extract containing piperidine hydrochloride, some moisture will persist, which can be removed by drying overnight in a dessicator.

References

1. Koles, Joseph E., Special Testing and Research Laboratory, BNDD - private communication.
2. Kram, Theodore C., Special Testing and Research Laboratory, BNDD - unpublished data.

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT



S-17-E FTIR/DRIFT SPECTROSCOPY		SAMPLE NUMBER (CS-1)	
CONCENTRATION: 2 mg./100 mg. KBr		SCANNING SPEED: 4	
CELL PATH: REFERENCE		SCAN SPEED: 4	
		DATE	
		PERSON ELUMER	
		PART NO. 437-5001	
		REF. NO.	

Property of Idaho State Police Forensic Services
 Unclassified Internet Copy
 OBSOLETE DOCUMENT

VI 511
22

DRUG PURIFICATION BY ION PAIR EXTRACTION

James M. White

Melvyn C. Kong

Orange County Sheriff's Department

Laboratory of Criminalistics

Santa Ana, California

Many drug preparations submitted to Forensic Laboratories require separation of the restricted substance present before identification and quantitation can be made by many of the conventional tests used for this purpose (e.g., Ultra-violet absorption, Infra-red absorption, crystal tests, etc.)

Methods routinely used to effect this separation include solvent-solvent extraction, column chromatography, thin-layer chromatography and gas chromatography.

Gas and thinlayer chromatographic separations are often disadvantageous because of analysis time and their inability to separate effectively the milligram quantities of substance often needed for subsequent analysis.

Column chromatographic methods, although highly selective, do not provide the rapid separation often needed in a high volume forensic drug laboratory.

Certain drugs form an ion pair in an HCl Solution. This ion pair is more soluble in chloroform than in the aqueous acid. Since aqueous acid and chloroform phases do not appreciably mix, this can provide an effective

means of separating certain drugs from adulterants, contaminants or other drugs. (See Analytical Chemistry, Vol. 39, No. 11, September 1967, pgs. 1283-1287).

Abstract

A Beckman DK-2A Spectrophotometer was used to quantitate recovery of drugs extracted under ion-pair and pH (saturated bicarbonate) conditions. Overall, bicarbonate extraction was found more efficient but less selective. The ion pair extraction is less efficient but more selective, lending itself to mixtures of drugs.

Procedure

I Ion Pair

The drugs to be evaluated were prepared in 10% HCl (v/v) solution to obtain an absorbance between .3 and .8. Five ml of this drug solution were pipetted into a 60 ml separatory funnel. Twenty five ml of chloroform were added; the mixture shaken for 1 minute and then allowed to separate for 2 minutes. The chloroform layer was drawn off and filtered. (Whatman No. 1 filter paper was used except for "LSD", "PCP" and "TCP" where glass wool was used). A second 25 ml of CHCl_3 was added and the process repeated. To the resulting 50 ml CHCl_3 was added 5 ml of 5% NaHCO_3 , shaken for 1 minute and allowed to separate for 2 minutes. The CHCl_3 was drawn off to a 50 ml stoppered graduate cylinder. To this CHCl_3 was added 5 ml of 0.2 N

H₂SO₄ and shaken for 1 minute and allowed to separate for 2 minutes. This 0.2 N H₂SO₄ was then measured for absorbance in the appropriate ultra-violet region on the DK-2A.

II Bicarbonate Extraction

Similarly, after obtaining and recording an absorbance between 0.3 and 0.8, 5 ml of drug solution (0.2 N H₂SO₄) were pipetted into a 60 ml separatory funnel. The solution was neutralized to pH 8.4 with excess solid NaHCO₃. To the solution was added 25 ml CHCl₃, the mixture shaken for 30 seconds and allowed to separate for 15 seconds. The CHCl₃ was drawn off. Another 25 ml of CHCl₃ was added and the process repeated. To the approximately 50 ml of CHCl₃ was added 5 ml 0.2 N H₂SO₄. This mixture was shaken for 30 seconds and allowed to separate for 2 minutes. The 0.2 N H₂SO₄ was then measured for absorbance in the appropriate ultraviolet region of the DK-2A.

Calculations

The wavelength of maximum absorbance was chosen on both plots, original and extracted.

For single drugs $\frac{\text{height extracted peak}}{\text{height original peak}} = \% \text{ recovery} = \text{recovery coefficient}$

For mixtures
1 : 1 $\frac{\text{height extracted peak}}{\text{height original solution peak} \times \frac{1}{2}} = \% \text{ recovery} = \text{recovery coefficient}$

Standard Deviation

$$S = \sqrt{\frac{\sum (\bar{x} - x)^2}{N-1}} \quad \bar{x} = \frac{\sum x}{N} \quad N = \text{Number of trials}$$

Error Analysis

There appears to be a small amount of CHCl_3 loss by evaporation, absorption into filter paper and solubility into aqueous phase. Adjustments were not made for the loss. In the direct comparison of absorbances, care was taken so the baseline coincided for original and extracted solutions. The same procedure was used in zeroing the instrument each time and the same sample cells were balanced, so the difference attributed to correcting for absolute zero absorbance is negligible.

Results

Ion pair extraction will separate cocaine at about a 25% efficiency from tetracaine, procaine and benzocaine. Other methods must be used to separate hexylcaine or lidocaine from cocaine. Ion pair will not separate "PCP", "LSD" or "TCP" from each other. Ion pair will separate meperidine, cocaine and "PCP" from nicotine.

Bicarbonate extraction in general provided higher yields but less selectivity. Apparently, the only advantage of HNO_3 over HCl in ion pair extractions would be the increased yield of cocaine (50% vs 27%).

The complete testing of HCl , H_2SO_4 , NaHCO_3 , and CHCl_3 solutions used in the extractions provided valuable information, such as the fact that most of the benzocaine remains in the CHCl_3 phase in this ion pair extraction scheme.

Table 1: Extraction Recoveries: Ion Pair vs pH Extraction

Drug	Ion Pair (10% HCl) Extraction		Bicarb Extraction	
	\bar{x}	s	\bar{x}	s
Cocaine	.272	.011	.906	.028
Hexylcaine	.822	.041	-	-
Lidocaine	.227	.013	-	-
Tetracaine	Negligible *	-	-	-
Procaine	Negligible *	-	-	-
Benzocaine	Negligible **	-	-	-
Phencyclidine ("PCP")	.851	.028	.827	.027
LSD	.791	.059	.897	.006
1-[1-(2-thienyl) cyclohexyl] piperidine ("TCP")	.830	.059	.821	.020
Nicotine	Negligible *	-	-	-
Meperidine	.887	.017	.901	.016
Amphetamine	Negligible *	-	.903	.034
Codeine	-	-	.906	.028
Methamphetamine	-	-	.917	.017
Mescaline	-	-	.915	.016
Heroin	.938	.023	-	-

Table 2: Separation of Mixtures by Ion Pair Extraction

Mixture	Results
Cocaine/Hexylcaine	No separation
Cocaine/Lidocaine	No separation
Cocaine/Tetracaine	Separation 26.9% cocaine
Cocaine/Procaine	Separation 24.0% cocaine
Cocaine/Benzocaine	Separation 23.6% cocaine
PCP/LSD	No separation
PCP/TCP	No separation
TCP/LSD	No separation
Meperidine/Nicotine	Separation 80.6% meperidine
PCP/Nicotine	Separation 87.1% PCP
Cocaine/Nicotine	Separation 23.6% cocaine
Amphetamine/Nicotine	No separation

* Drug remains in 10% HCl

** 18% stays in 10% HCl, remainder in CHCl₃

PHENCYCLIDINE (PCP) ON MINT

J. W. Smith and G. P. Chastean

December 1974

501 D. Wilson St, 3rd Floor
Los Angeles, CA 90012

213-674-4611

23

SOUTHCOMBE

An effective extraction procedure for PCP on mint has been found and is herein reported. This method is more complex than the method previously suggested by J. Smith, however, it is recommended due to the possibility of obtaining an infrared spectrum. The spectrum obtained is comparable to that of the primary standard.

PROCEDURE

Place approximately 1 gram of plant material in a 13 X100 test tube. The method should work well with smaller amounts. Add enough conc. NH_4OH to wet the material, approx. 0.5 ml. Fill the remainder of the test tube with methylene chloride and thoroughly mix by shaking or stirring with a wooden stick. Remove the methylene chloride and extract the plant material again with methylene chloride. The methylene chloride solution will be green in color.

Prepare an alumina column in either a column or pipette (both have been used successfully). Run the methylene chloride through the column and wash with excess methylene chloride until the colored bands are well into the alumina. The methylene chloride should be either amber or colorless at this point.


Evaporate the methylene chloride to approximately 10 ml and extract with 0.1 N HCl. Obtain an ultraviolet spectrum on the acid layer. The UV spectrum should be fairly clean and typical for PCP.

Add Na_2CO_3 to make basic and extract with methylene chloride. Evaporate the methylene chloride adding approx. 2 ml methanol with 1 drop conc. HCl at first and again when near dryness. The previous mentioned infrared spectrum was run directly on this product.

Crystals, GC and TLC can also be run using this product. It will be noted that our present systems 3, 5, and 6 do not separate PCP and its thiophene analog ~~anor~~ with the utilization of multiple development.

WHEN REORDERING SPECIFY CHART NO. 104411

BECKMAN INSTRUMENTS, INC., FULLER



DEA LABORATORY NOTES


IDENTIFICATION OF PROPOXYPHENE IN MIXTURES CONTAINING APC OR ACETAMINOPHEN

Katherine T. Churchill
Southeast Regional Laboratory
Miami, Florida

Dissolve sample in a small amount of 2N HCl (e.g. 1 capsule of Darvon Compound/10 ml 2N HCl) and filter. Make the filtrate strongly alkaline with 5N NaOH (approximately 30 ml). Allow the precipitated propoxyphene base to flocculate and settle out (approximately 45 minutes minimum). Filter and discard the filtrate. Wash the residue from the filter paper with methylene chloride, dry the methylene chloride by passing through anhydrous sodium sulfate. Evaporate to dryness, cool and speed crystallization by scratching the beaker. The crystalline propoxyphene base produced may be used for IR or NMR and polarimetry.


IR and NMR spectra of the resulting crystalline propoxyphene base are presented.

Propoxyphene base is soluble in ether - to remove any last contaminants I added petroleum ether - removed it from



IR Spectrum of Propoxyphene base

the small amt. of ppt. that resulted. (Pet. ether added after CH₂Cl₂ evaporated.)



solv CDCl₃ NMR Spectrum of Propoxyphene base

Poentele



STATE OF IDAHO

DEPARTMENT OF HEALTH AND WELFARE

BUREAU OF LABORATORIES, 2220 Old Penitentiary Road, Boise, ID 83702
FORENSIC SECTION (208) 334-2231

March 5, 1985

SOUTHEASTERN IDAHO FORENSIC
LABORATORY RECEIVED

DATE 3-11-85

TIME 8:00

BY [Signature]

MEMO TO: Forensic Section

FROM: Richard D. Groff, Supervisor

SUBJECT: Sampling Technique to be Used in the Analysis of Solid Dosage Drug Analysis. (Please place in Policy Section of Methods Manual)

- I. Do not use more than half of the sample. (The defendant should have an amount equal to what we used for analysis by an expert of his own choosing.)
 - a. If it is obvious that, with the macro-techniques available in Idaho forensic labs, more than half will be used, consider forwarding the sample to the D.E.A.
 - b. If all the material has to be manipulated during the examination, some way should be found to return at least half of the original specimen. For example, if residue off a spoon is extracted by washing it with methanol, then half the methanol extract should be dried in a glass vial and returned with the rest of the evidence.

- II. Single samples - powder only.
 - a. Under 1 gram: Mix powder thoroughly and take one representative sample, using your own best judgment as to amount.
 - b. A gram or more: Take several portions from different locations in the sample to be mixed thoroughly and used for analysis. This should total .15 gram or more.

- III. Multiple samples (powders, tablets, capsules and plant material.) For this sampling technique to be applicable, each item must have the same appearance. (For instance: In dealing with baggies of marijuana, each baggie should have the same general appearance as to contents, color, size, container, etc.)
 - a. Determine the total net weight of the question material based on the tare weight of the packaging material from one or more items.
 - b. For 1 to 2500 items, form a composite sample by combining portions of the square root of the number of items.
 - c. For 26 to 2500 items, analyze composite sample of the square root of the total number of items. For 2500 items or more, analyze a composite sample of 50 items. If the number of items exceeds 125,000, analyze a composite sample of the cube root.

[Signature]
Richard Groff
Supervising Criminalist

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

PROCEDURES/IDEAS
from
DEA LABORATORY
SAN FRANCISCO
JANUARY 22-24, 1985

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

Susan C. Carre'
Pocatello Crime Lab

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

BENZODIAZEPINES
IR PROCEDURE

BENZODIAZEPINES & DERIVATIVES; Chlordiazepoxide (Librium)
Diazepam (Valium)
Flurazepam (Dalmane)
Prazepam (Centrax)
Oxazepam (Serax)

Diazepam: 1. Place powder on top of a plug of cotton in a pipet.
(base)

2. Extract with ethyl ether

3. Dry down & run IR scan

Others: 1. Place powder on top of plug of cotton in a pipet.
(HCl)

2. Extract with CH_2Cl_2 - reduce volume

3. Recrystallize with pet ether/hexane

4. Run IR scan (beware of polymorphism)

Sometimes it may be necessary to do an acid wash to clean-up a pharmaceutical.

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

COCAINE
TARTARIC ACID

1. Dissolve sample in water in a small test tube.
2. Make basic with a drop of NH_4OH
3. Add one drop of sat'd NaCl solution
4. Extract with pet ether
5. Run pet ether through cotton plugged pipet
6. Reduce volume to a couple of drops
7. Estimate amount of cocaine present and add corresponding amount of Di-p-toluoyl-l-tartaric acid monohydrate (end of spatula)
8. Add approx. $\frac{1}{2}$ ml acetone. Drip down side to rinse in tartaric acid.
9. Let stand until crystals form. Pipet off excess liquid & use blotter paper to remove all liquid.
10. Rinse crystals with acetone - dry
11. Run IR scan


Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

p-DMAB
FLYING SAUCER

When using p-DMAB to visualize LSD and psilocybin/cybin the following is safer than spraying the reagent.


1. Prepare a sat'd solution of p-DMAB in Pet Ether. (Evaporates off)
2. apply p-DMAB to plate with cotton or Q-tip
3. place plate in a watch plate that has a squirt of conc. HCl in it and place an additional plate over the top,
4. Allow to sit until spots visualize.



Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

Do-it-yourself
GAS CELL



RUN IR ON SLOW
CUT DOWN SLIT TO 2


This is especially good for use in analyzing materials from lab raids.
e.g. methylamine & other organics

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT



Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

HEROIN
IR PROCEDURES

Especially for small amounts



Or:



*The alumina will absorb part of the heroin that passes through it which is why this technique is less effective for small amounts.

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

LSD
IR Procedure

1. Extract with CH_2Cl_2 from a basic soln. (NH_4OH)
2. Run through an alumina column
3. Run through $\text{CH}_2\text{Cl}_2/\text{MeOH}$ - watch for LSD band
4. Reduce volume to 1 ml
5. Run an IR using a very small amount of KBr and a pin hole in a card.
6. If the above IR needs clean-up, do so on a TLC plate

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

METHAMPHETAMINE
PIT DERIVATIVE

1. Dissolve sample in water. If sample was obtained from a lab raid, use dilute H_2SO_4 and wash with CH_2Cl_2 . (To remove phenylacetic acid and P-2-P.)
2. Make aqueous basic (drop of NH_4OH)
3. Extract with pet ether
4. Filter through cotton
5. Reduce volume (preferably over a steam bath) Heating drives off Methylamine.
6. Place melting point tube into PIT and place it in pet ether.
7. Let solution stand. White ppt. should form within 15 minutes.
8. Suck off liquid.
Wash with more pet ether.
Suck off pet ether and use blotter paper to remove all of the liquid.

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

MUSHROOMS
GC CLEAN-UP

1. Wash material with ether.
2. Grind material in MeOH
3. Filter with Buchner & glass fiber filter paper.
4. Evaporate MeOH over steam with air.

n-BuOH sat'd water:HOAc (9:1) is an excellent system. TLC $R_f = .82$

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

d-PROPOXYPHENE (Darvon)
IR PROCEDURE


1. Add water to ground tablet in tt.
2. Filter through cotton plugged pipet
3. Make basic with drop of NH_4OH
4. Extract with Pet Ether - Suck-up with pipet to mix.
5. Pipet off Pet Ether into another tt
6. Add water (equal volume) to wash
7. Filter Pet Ether through cotton
8. Evaporate down - tease until crystals form or
- run on salt plates in CH_2Cl_2
9. Run IR scan

****THIS PROCEDURE CAN BE USED FOR MANY PHARMACEUTICAL DRUGS****

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

MIXTURES
SEPARATION OF AMPHETAMINE
& METHAMPHETAMINE



Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

MIXTURES
SEPARATION SCHEME

DRY EXTRACTION WITH:

1. Pet ether (will dissolve the least)
2. Ethyl ether
3. CH_2Cl_2
4. CH_3Cl
5. CH_3/MeOH 3:1 (best @ dissolving)

This procedure will indicate: What the substance is soluble in & suggests clean-up scheme.

What TLC system to use to separate

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT

Property of Idaho State Police Forensic Services
Uncontrolled Internet Copy
OBSOLETE DOCUMENT