

## Section Two

### Urine Toxicology

---

#### 2.2 ANSYS<sup>®</sup> Thin Layer Chromatography (TLC) Methods

##### 2.2.4 Toxi-Lab<sup>®</sup> THC II-PLUS 11-nor- $\Delta$ 9-tetrahydrocannabinol-9-carboxylic acid ( $\Delta$ 9-THC-COOH) Detection System

---

###### 2.2.4.1 BACKGROUND

The TOXI-LAB<sup>®</sup> THC thin layer chromatography (TLC) drug detection system provides for extraction, concentration, inoculation, elution, and visualization steps for the detection of  $\Delta$ 9-THC-COOH in urine specimens.<sup>1</sup> The preliminary identification is based on matching the position of a drug (*R<sub>f</sub>*) and visualization color characteristics with that of corresponding reference material.

###### 2.2.4.2 SCOPE

This method is an option to screen for the presence of  $\Delta$ 9-THC-COOH in urine. The results serve to support the results of the enzyme immunoassay (EIA) screen or used in lieu of EIA screen. The TOXI-LAB THC II system provides a preliminary result that must be confirmed by GC-MSD.

###### 2.2.4.3 EQUIPMENT AND SUPPLIES

- 2.2.4.3.1 Tube rocker
- 2.2.4.3.2 Laboratory centrifuge
- 2.2.4.3.3 Solvent concentrator with appropriate concentration cups
- 2.2.4.3.4 Electric (plate) warmer
- 2.2.4.3.5 Fixed and adjustable volume single channel air displacement pipettors, and appropriate tips, capable of accurate and precise dispensing of volumes indicated.
- 2.2.4.3.6 Forceps
- 2.2.4.3.7 Disc handling pins
- 2.2.4.3.8 Index cards for use as disc press cards
- 2.2.4.3.9 TOXI-GRAMS Blank THC-II-PLUS
- 2.2.4.3.10 TOXI-GRAMS Blank THC-II
- 2.2.4.3.11 TOXI-DISCS THC
- 2.2.4.3.12 SPEC-C18-1 Extraction Cartridges
- 2.2.4.3.13 THC II Wash Reagent 1 Bottle
- 2.2.4.3.14 THC II Wash Reagent 2 Bottle
- 2.2.4.3.15 TOXI-DIP THC-1 Reagent Tank with Lid
- 2.2.4.3.16 TOXI-DIP THC-2 Reagent Tank with Lid
- 2.2.4.3.17 HCl Reagent Tank with Lid
- 2.2.4.3.18 Chromatography Tank with Lid (THC-II-PLUS)
- 2.2.4.3.19 Chromatography Jar with Lid (THC-II)
- 2.2.4.3.20 TOXI-LAB THC Elution Solvent Bottle

## 2.2.4.4 REAGENTS

*Refer to manual section 5.12 for solution preparation instructions not listed below.*

- 2.2.4.4.1 11.8N KOH
- 2.2.4.4.2 Methanol (ACS Grade)
- 2.2.4.4.3 Wash Reagent 1  
**20% Acetic Acid**
- 2.2.4.4.4 Wash Reagent 2  
20% **Methylene Chloride (ACS Grade)** in **n-Heptane (ACS Grade)**.
- 2.2.4.4.5 TOXI-DIP THC-1 Fast Blue BB  
Prepared with approximately 1g **Fast Blue BB Salt (Purified Grade)**. Add Fast Blue BB to reagent tank for TOXI-DIP THC-1. Add approximately 700mL **Methylene Chloride (ACS Grade)**. Solution should be pale yellow in color. Mix well.  
Store at room temperature. Solution stable for 2 - 3 months.
- 2.2.4.4.6 TOXI-DIP THC-2 Diethylamine Fuming  
Pipet 40mL **Diethylamine (DEA) (ACS Grade)** through an opening in the standoff to the bottom of the tank. Remove any DEA on standoff surface.  
Store at room temperature. Replace DEA weekly.
- 2.2.4.4.7 Hydrochloric Acid Fuming  
Pipet 40mL **Concentrated Hydrochloric Acid (ACS Grade)** through an opening in the standoff to the bottom of the tank. Remove any HCl on standoff surface.  
Store at room temperature. Replace HCl weekly.
- 2.2.4.4.8 THC II Stock Elution Solvent  
In THC II Elution Solvent Bottle, mix 50mL **n-Heptane (ACS Grade)**, 50mL **Acetone (ACS Grade)** and 1mL **Glacial Acetic Acid (ACS Grade)**. Cap tightly and mix.  
Store at room temperature.

## 2.2.4.5 REFERENCE MATERIAL

- 2.2.4.5.1 Positive Control  
Positive Control can be prepared by adding specified amount of Working Control Solution to negative urine and/or obtained commercially. Use the same lot of negative urine to prepare positive control as used to prepare negative control.
- 2.2.4.5.1.1 **Stock Reference Solution**  
100µg/mL (+) 11-nor-9-carboxy- $\Delta^9$ -THC
- 2.2.4.5.1.2 **Working Reference Solution (1800ng/mL)**

Add 900 $\mu$ L Stock Solution to 49.1mL Methanol. Solution is stable for six months when stored at 4°C.

- 2.2.4.5.2 Negative Control  
Negative Urine  
May be either in-house verified or commercially obtained.

## 2.2.4.6 PROCEDURE

- 2.2.4.6.1 Initial set-up  
Label extraction tubes and extraction cartridges for the negative control, positive control, and appropriate laboratory numbers.
- 2.2.4.6.2 60ng/mL Carboxy-THC Positive Control
- 2.2.4.6.2.1 Transfer 6mL of negative urine to extraction tube.
- Use the same lot of urine used for negative control.
- 2.2.4.6.2.2 200 $\mu$ L of working reference solution. Vortex.
- 2.2.4.6.3 Negative Control  
Transfer 6mL of negative urine to extraction tube.
- 2.2.4.6.4 Casework Sample Preparation  
Transfer 6 mL of casework urine specimen to extraction tube.
- 2.2.4.6.5 Sample Hydrolysis
- 2.2.4.6.5.1 To 6mL of urine, add 12 drops 11.8N KOH. Vortex.
- 2.2.4.6.5.2 Allow to hydrolyze for 10 minutes.
- 2.2.4.6.5.3 Add 1.5mL glacial acetic acid. Vortex.
- 2.2.4.6.6 Extraction
- 2.2.4.6.6.1 Condition cartridge with 1mL methanol. Aspirate at approximately 5 in. Hg. **Do not allow the disc to dry.**
- 2.2.4.6.6.2 Add acidified samples to cartridge reservoirs. Aspirate such that the sample passes through the column no faster than 2mL/min.
- 2.2.4.6.6.3 Once the sample is completely through the reservoir, remove filter.
- 2.2.4.6.6.4 Add 1mL 20% acetic acid. Aspirate  $\geq$ 2 minutes at 10-12 in. Hg.

- 2.2.4.6.6.5 Add 500 $\mu$ L wash reagent. Aspirate at 10-12 in. Hg.
- 2.2.4.6.6.6 After solvent has past through, allow to aspirate  $\geq 2$  minutes.
- 2.2.4.6.6.7 Remove disc from cartridge and place into a pre-heated concentrated cup to remove all residual moisture.
- 2.2.4.6.7 TLC
- 2.2.4.6.7.1 Place disc into labeled three or 10-channel TOXI-GRAM for THC-II.
- 2.2.4.6.7.2 Add THC-DISC THC disc.
- 2.2.4.6.7.3 Note: If not all channels are used, remove excess with razor blade or scissors.
- 2.2.4.6.7.4 Heat the GRAM, with the disc end slightly off the warmer edge, for 30-60 seconds.
- 2.2.4.6.7.5 Add 12.5mL of developing solution to chromatography tank (10-channel) or 3mL solution to chromatography jar (3-channel).
- 2.2.4.6.7.6 Place gram into chromatography tank or jar. Allow dye marker to migrate to  $\cong 4$  cm. [***This only takes 2-3 minutes***]
- 2.2.4.6.7.7 Remove GRAM from tank/jar and place face down on warmer for 1-2 minutes.
- 2.2.4.6.7.8 Dip GRAM into TOXI-DIP THC-1, hold to dry until GRAM becomes speckled.
- 2.2.4.6.7.9 Place GRAM into TOXI-DIP THC-2 until scarlet spots develop.
- 2.2.4.6.7.10 Place GRAM in hood so that the Diethylamine (DEA) can evaporate. If any DEA is present when the HCl is added, fuming will occur.
- 2.2.4.6.7.11 For HCl fuming, either of the following options may be pursued. Option 2 will produce a more intense color which photocopies better.
- Option 1: Place GRAM into HCl Fuming Tank until a purple spot develops. Place GRAM into page protector, label and photocopy.
- Option 2: Place GRAM on to a page protector. With bulb pipet, add concentrated HCl to just cover GRAM. Note desired color change to deep purple. Label and photocopy GRAM.
- 2.2.4.6.7.12 Place a copy of GRAM into each associated casefile.

**2.2.4.7 DETECTION AND IDENTIFICATION CRITERIA**

The position (*R<sub>f</sub>*) and color characteristics at each stage of visualization of a spot noted for a specimen must correspond to that of reference material.

**2.2.4.8 REFERENCES AND RECOMMENDED READING**

2.2.4.8.1 Toxi-Lab<sup>®</sup> THC II 11-nor- $\Delta$ 9-Tetrahydrocannabinol-9-COOH Detection System Instruction Manual, © 1998.

2.2.4.8.2 Toxi-Lab<sup>®</sup> THC II-PLUS 11-nor- $\Delta$ 9-Tetrahydrocannabinol-9-COOH Detection System Instruction Manual, © 1998.

Property of Idaho State Police Forensic Services  
Uncontrolled Internet Copy  
OBSOLETE DOCUMENT 01/07/2013

## *Revision History*

---



---

### Section Two

#### Urine Toxicology

---

#### **2.2 ANSYS® Thin Layer Chromatography (TLC) Methods**

##### 2.2.4 Toxi-Lab® THC II-PLUS 11-nor- $\Delta$ 9-tetrahydrocannabinol-9-carboxylic acid ( $\Delta$ 9-THC-COOH) Detection System

---



---

<b>Revision No.</b>	<b>Issue Date</b>	<b>Comments</b>
0	10/1991	Original Issue
1	11-27-2001	Introduction into Reformatted SOP Manual
2	04-25-2002	THC-II Method Summary Added
3	10-18-2002	Refinements, other TLC methods added to alternative method binder
4	05-07-2007	Reformat, c-THC TLC only
5	07-28-2008	Clarified that negative urine used to prepare positive control must be the same lot as used for negative control.