

# *Idaho State Police*

## *Forensic Services*

### *Approval for Quality System Controlled Documents*



Discipline/Name of Document: Toxicology – BAC Calculation Excel  
Spreadsheet (casefile, calibrators and controls tabs)

Revision Number: 0

Issue Date: 6/5/2007

APPROVED BY: *Coinna C. Orsley*  
Quality Manager

6/5/07  
Date Signed

**VOLATILES DETERMINATION CASEFILE WORKSHEET**

Laboratory No.:

Analysis Date(s):

Sample Type:  Blood  Vitreous Humor  Urine

	Column 1 Channel A	Column 2 Channel B	Column Precision	Mean Value	Over-all Mean	Reported Result
Sample Results	0.0000	0.0000	0.0000	0.0000	0.0000	0.00
(g/100cc)	0.0000	0.0000	0.0000	0.0000		

**Analysis Method**

Refer to Toxicology Analytical Method 4.1.

**Instrument Information**

Perkin Elmer HS-40 Automatic Headspace Sampler  
AutoSystem XL Gas Chromatograph equipped with dual FID detectors.

**Columns**

Vendor	Column	Channel	Film thickness	Length	ID	
Restek	Rtx®-BAC1	A	1.80mm	30M	0.32mm	
Restek	Rtx®-BAC2	B	1.20mm	30M	0.32mm	

**Gas Chromatograph Instrumental Temperature Conditions**

Initial Oven	Ramp Rate	Final Oven	Hold	Injection Port	Detector A	Detector B
35°C	8°C/min	100°C	None	200°C	150°C	150°C

**Headspace Analyzer Instrumental Conditions**

Needle Temperature	90°C		Sample Oven	60°C
Transfer Line	90°C		Cycle Time	10 minutes
Thermostating Time	12 minutes		Pressurization Time	0.5 minute
Injection Time	0.08 minutes		Withdrawal Time	0.2 minutes

Calibration and control data are stored centrally.

Analyst: \_\_\_\_\_

Page \_\_\_\_ of \_\_\_\_

Quantitative Analysis for Ethanol

Ethanol Calibration Standards

Run Date(s):

	0.05g/100cc	0.08g/100cc	0.10g/100cc
Source of Standard	Cerilliant	Cerilliant	Cerilliant
Lot Number	FN072105-01	FN082905-01	FN041204-01
Results for Column 1 Channel (A)	0.0494	0.0796	0.0979
Results for Column 2 Channel (B)	0.0494	0.0804	0.0982
Column Precision	0.0000	-0.0008	-0.0003
Mean	0.0494	0.0800	0.0980
	0.20g/100cc	0.30g/100cc	0.40g/100cc
Source of Standard	Cerilliant	Cerilliant	Cerilliant
Lot Number	FN012204-01	FN031204-01	FN033006-02
Results for Column 1 Channel (A)	0.1999	0.2994	0.4012
Results for Column 2 Channel (B)	0.1996	0.2994	0.4011
Column Precision	0.0003	0.0000	0.0001
Mean	0.1997	0.2994	0.4011
Analyst: _____ Page _____			

**Quantitative Analysis for Ethanol & Qualitative Analysis for Other Volatiles**

**Volatiles Quality Assurance Controls**

**Run Date(s):**

Property of Idaho State Police Forensic Services Uncontrolled Internet Copy OBSOLETE DOCUMENT								
Control Type	Whole blood				Whole blood			
Mean Value	Level 1: 0.0565g/100cc				Level 2: 0.1396 g/100cc			
Acceptable Range	0.0508 - 0.0621				0.1256 - 0.1535			
Source of Control	BIORAD				BIORAD			
Lot Number	44591				44592			
Results for Column 1 Channel (A)	0.0551	0.0551			0.1415	0.1408		
Results for Column 2 Channel (B)	0.0545	0.0544			0.1401	0.1400		
Column Precision	0.0006	0.0007	0.0000	0.0000	0.0014	0.0008	0.0000	0.0000
Sample Mean	0.0548	0.0547	#DIV/0!	#DIV/0!	0.1408	0.1404	#DIV/0!	#DIV/0!
Control Type	Multi-Component Aqueous						Water Blank	ISTD Blank
Mean Value	0.40g/100cc							
Acceptable Range								
Source of Control	Cerilliant							
Lot Number	FN082905-01							
Results for Column 1 Channel (A)	0.4415	0.4233					No peaks	0.0000
Results for Column 2 Channel (B)	0.4162	0.4245					No peaks	0.0000
Column Precision	0.0253	-0.0012	0.0000	0.0000	0.0000	0.0000		
Sample Mean	0.4288	0.4239	#DIV/0!	#DIV/0!	#DIV/0!	#DIV/0!		

Analyst: \_\_\_\_\_

Page \_\_\_\_\_