

**Quantitative Analysis for Ethanol & Qualitative Analysis for Other Volatiles**

*Analytical Method(s): 1.0*

*Device: Hamilton MICROLAB 503A Liquid Processor/Dilutor Serial Number: MD-96BC1382/MD94AM10010*

**Volatiles Quality Assurance Controls**

**Run Date(s): 03/22/2017-03/23/2017**

Control level	Expiration	Lot #	Target Value	Acceptable Range	Overall Results	
Level 1	Jul-18	1407031	0.0780	0.0702 - 0.0858	0.0793 g/100cc	
					0.0811 g/100cc	
					g/100cc	
Level 2	Jul-18	1407032	0.2020	0.1818 - 0.2222	0.2098 g/100cc	
					0.2113 g/100cc	
<b>Multi-Component Mixture</b>		<b>Exp: Oct 2019</b>	<b>Lot #</b>	FN09231404	OK	
<b>Curve Fit:</b>			<b>Column 1</b>	1.00000	<b>Column 2</b>	0.99982

<b>Ethanol Calibration Reference Material</b>								
Calibrator level	Expiration	Cerilliant Lot #	Target Value	Acceptable Range	Column 1	Column 2	Precision	Mean
0.050	Jul-19	FN06231406	0.050	0.045 - 0.055	0.0506	0.0538	0.0032	0.0522
0.080			0.080	0.072 - 0.088			0	#DIV/0!
0.100	Jun-20	FN06181501	0.100	0.090 - 0.110	0.0995	0.0998	0.0003	0.0996
0.200	Oct-20	FN07201502	0.200	0.180 - 0.220	0.2000	0.1972	0.0028	0.1986
0.300	Jun-20	FN06051501	0.300	0.270 - 0.330	0.2995	0.2962	0.0033	0.2978
0.400			0.400	0.360 - 0.440			0	#DIV/0!
0.500	Aug-19	FN07031402	0.500	0.450 - 0.550	0.5003	0.5031	0.0028	0.5017

<b>Aqueous Controls</b>					
Control level	Expiration	Cerilliant Lot #	Target Value	Acceptable Range	Overall Results
0.080	Nov-20	FN10281510	0.08000	0.076 - 0.084	0.082 g/100cc

Issued: 4/22/2015






-Any information on this document can be changed for laboratory use, except for the precision and mean determination formulas.

Volatiles QA/QC data spreadsheet Rev 5

Issuing Authority: Quality Manager


JG

**Worklist: 1620**

<u>LAB CASE</u>	<u>ITEM</u>	<u>TASK ID</u>	<u>DESCRIPTION</u>	
M2017-1038	1	78739	Alcohol Analysis	
M2017-1053	1	78764	Alcohol Analysis	
M2017-1054	1	78765	Alcohol Analysis	
M2017-1055	1	78766	Alcohol Analysis	
M2017-1056	1	78770	Alcohol Analysis	
M2017-1057	1	78771	Alcohol Analysis	
M2017-1058	1	78772	Alcohol Analysis	
M2017-1077	1	78879	Alcohol Analysis	
M2017-1095	1	78973	Alcohol Analysis	
M2017-1096	1	78977	Alcohol Analysis	
M2017-1110	1	79060	Alcohol Analysis	
M2017-1115	1	79068	Alcohol Analysis	
M2017-1116	1	79069	Alcohol Analysis	
M2017-1117	1	79073	Alcohol Analysis	
M2017-1118	1	79074	Alcohol Analysis	
M2017-1119	1	79075	Alcohol Analysis	
M2017-1122	1	79154	Alcohol Analysis	
M2017-1148	2	79297	Alcohol Analysis	
M2017-1169	1	79405	Alcohol Analysis	
M2017-1170	1	79406	Alcohol Analysis	
M2017-1171	1	79410	Alcohol Analysis	
M2017-1173	1	79420	Alcohol Analysis	
M2017-1174	1	79421	Alcohol Analysis	

**Worklist: 1620**

<u>LAB CASE</u>	<u>ITEM</u>	<u>TASK ID</u>	<u>DESCRIPTION</u>
M2017-1188	1	79476	Alcohol Analysis
M2017-1190	1	79497	Alcohol Analysis
P2017-0449	1	77963	Alcohol Analysis
P2017-0501	1	78439	Alcohol Analysis
P2017-0548	1	79040	Alcohol Analysis



=====  
Calibration Table  
=====

-----  
General Calibration Setting  
-----

Calib. Data Modified : Wednesday, March 22, 2017 11:27:59 AM  
Signals calculated separately : No

Rel. Reference Window : 0.000 %  
Abs. Reference Window : 0.100 min  
Rel. Non-ref. Window : 0.000 %  
Abs. Non-ref. Window : 0.100 min  
Uncalibrated Peaks : not reported  
Partial Calibration : Yes, identified peaks are recalibrated  
Correct All Ret. Times: No, only for identified peaks

Curve Type : Linear  
Origin : Ignored  
Weight : Equal

Recalibration Settings:  
Average Response : Average all calibrations  
Average Retention Time: Floating Average New 75%

Calibration Report Options :  
Printout of recalibrations within a sequence:  
Calibration Table after Recalibration  
Normal Report after Recalibration  
If the sequence is done with bracketing:  
Results of first cycle (ending previous bracket)

Default Sample ISTD Information (if not set in sample table):

ISTD #	ISTD Amount [g/100cc]	Name
1	1.00000	n-propanol
2	1.00000	n-propanol

-----  
Signal Details  
-----

Signal 1: FID1 A, Front Signal  
Signal 2: FID2 B, Back Signal  
-----

-----  
Overview Table  
-----

JG

RT	Sig	Lvl	Amount [g/100cc]	Area	Rsp.Factor	Ref	ISTD #	Compound
2.586	1	1	1.00000	3.69669	2.70512e-1	No	No 1	methanol
3.072	1	1	5.00000e-2	4.44027	1.12606e-2	No	No 1	ethanol
		2	1.00000e-1	9.04249	1.10589e-2			
		3	2.00000e-1	17.89574	1.11758e-2			
		4	3.00000e-1	27.30493	1.09870e-2			
		5	5.00000e-1	44.96452	1.11199e-2			
3.388	2	1	1.00000	4.26062	2.34707e-1	No	No 2	methanol
3.628	1	1	1.00000	9.73055	1.02769e-1	No	No 1	isopropyl alcohol
4.282	2	1	5.00000e-2	4.47271	1.11789e-2	No	No 2	ethanol
		2	1.00000e-1	9.08370	1.10087e-2			
		3	2.00000e-1	18.09025	1.10557e-2			
		4	3.00000e-1	27.94878	1.07339e-2			
		5	5.00000e-1	46.74831	1.06956e-2			
4.308	1	1	1.00000	6.49940	1.53860e-1	No	No 1	acetone
4.619	1	1	1.00000	42.17959	2.37081e-2	No	Yes 1	n-propanol
		2	1.00000	43.06885	2.32186e-2			
		3	1.00000	42.11362	2.37453e-2			
		4	1.00000	42.79330	2.33681e-2			
		5	1.00000	42.10326	2.37511e-2			
4.661	2	1	1.00000	6.89301	1.45075e-1	No	No 2	acetone
4.969	2	1	1.00000	10.70642	9.34019e-2	No	No 2	isopropyl alcohol
7.553	2	1	1.00000	42.95251	2.32815e-2	No	Yes 2	n-propanol
		2	1.00000	43.51038	2.29830e-2			
		3	1.00000	42.03962	2.37871e-2			
		4	1.00000	42.63439	2.34552e-2			
		5	1.00000	41.51231	2.40892e-2			

Peak Sum Table


\*\*\*No Entries in table\*\*\*

41 Warnings or Errors (10 first messages follow) :


- Warning : Curve requires more calibration points., (methanol)
- Warning : Curve requires more calibration points. at 2.586 min, signal 1
- Warning : Curve requires more calibration points. at 3.388 min, signal 2
- Warning : Curve requires more calibration points. at 3.628 min, signal 1
- Warning : Curve requires more calibration points. at 4.308 min, signal 1
- Warning : Curve requires more calibration points. at 4.619 min, signal 1
- Warning : Curve requires more calibration points. at 4.661 min, signal 2
- Warning : Curve requires more calibration points. at 4.969 min, signal 2
- Warning : Curve requires more calibration points. at 7.553 min, signal 2
- Warning : Curve requires more calibration points. at 2.586 min, signal 1

JG


=====  
 Calibration Curves  
 =====




methanol at exp. RT: 2.586  
 FID1 A, Front Signal  
 Correlation: 1.00000  
 Residual Std. Dev.: 0.00000  
 Formula:  $y = mx + b$   
 m: 8.76418e-2  
 b: 0.00000  
 x: Amount Ratio  
 y: Area Ratio



ethanol at exp. RT: 3.072  
 FID1 A, Front Signal  
 Correlation: 1.00000  
 Residual Std. Dev.: 0.00121  
 Formula:  $y = mx + b$   
 m: 2.14080  
 b: -3.14635e-3  
 x: Amount Ratio  
 y: Area Ratio




methanol at exp. RT: 3.388  
 FID2 B, Back Signal  
 Correlation: 1.00000  
 Residual Std. Dev.: 0.00000  
 Formula:  $y = mx + b$   
 m: 9.91939e-2  
 b: 0.00000  
 x: Amount Ratio  
 y: Area Ratio




isopropyl alcohol at exp. RT: 3.628  
 FID1 A, Front Signal  
 Correlation: 1.00000  
 Residual Std. Dev.: 0.00000  
 Formula:  $y = mx + b$   
 m: 2.30693e-1  
 b: 0.00000  
 x: Amount Ratio  
 y: Area Ratio


JG




ethanol at exp. RT: 4.282  
 FID2 B, Back Signal  
 Correlation: 0.99982  
 Residual Std. Dev.: 0.00893  
 Formula:  $y = mx + b$   
 m: 2.27461  
 b: -1.81820e-2  
 x: Amount Ratio  
 y: Area Ratio



acetone at exp. RT: 4.308  
 FID1 A, Front Signal  
 Correlation: 1.00000  
 Residual Std. Dev.: 0.00000  
 Formula:  $y = mx + b$   
 m: 1.54089e-1  
 b: 0.00000  
 x: Amount Ratio  
 y: Area Ratio




n-propanol at exp. RT: 4.619  
 FID1 A, Front Signal  
 Correlation: 1.00000  
 Residual Std. Dev.: 0.00000  
 Formula:  $y = mx + b$   
 m: 1.00000  
 b: 0.00000  
 x: Amount Ratio  
 y: Area Ratio




acetone at exp. RT: 4.661  
 FID2 B, Back Signal  
 Correlation: 1.00000  
 Residual Std. Dev.: 0.00000  
 Formula:  $y = mx + b$   
 m: 1.60480e-1  
 b: 0.00000  
 x: Amount Ratio  
 y: Area Ratio

JG



isopropyl alcohol at exp. RT: 4.969  
FID2 B, Back Signal  
Correlation: 1.00000  
Residual Std. Dev.: 0.00000  
Formula:  $y = mx + b$   
m: 2.49262e-1  
b: 0.00000  
x: Amount Ratio  
y: Area Ratio



n-propanol at exp. RT: 7.553  
FID2 B, Back Signal  
Correlation: 1.00000  
Residual Std. Dev.: 0.00000  
Formula:  $y = mx + b$   
m: 1.00000  
b: 0.00000  
x: Amount Ratio  
y: Area Ratio


=====

JG



ISP Forensic Services Blood Alcohol Report

Sample Name : 0.050 FN06231406  
 Laboratory : Meridian  
 Injection Date : Mar 22, 2017  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	4.44027	0.0506	g/100cc
2.	Ethanol	Column 2:	4.47271	0.0538	g/100cc
3.	n-Propanol	Column 1:	42.17959	1.0000	g/100cc
4.	n-Propanol	Column 2:	42.95251	1.0000	g/100cc

JG

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.100 FN06181501  
 Laboratory : Meridian  
 Injection Date : Mar 22, 2017  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	9.04249	0.0995	g/100cc
2.	Ethanol	Column 2:	9.08370	0.0998	g/100cc
3.	n-Propanol	Column 1:	43.06885	1.0000	g/100cc
4.	n-Propanol	Column 2:	43.51038	1.0000	g/100cc

JG

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.200 FN07201502  
 Laboratory : Meridian  
 Injection Date : Mar 22, 2017  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	17.89574	0.2000	g/100cc
2.	Ethanol	Column 2:	18.09025	0.1972	g/100cc
3.	n-Propanol	Column 1:	42.11362	1.0000	g/100cc
4.	n-Propanol	Column 2:	42.03962	1.0000	g/100cc

UG

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.300 FN06051501  
 Laboratory : Meridian  
 Injection Date : Mar 22, 2017  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	27.30493	0.2995	g/100cc
2.	Ethanol	Column 2:	27.94878	0.2962	g/100cc
3.	n-Propanol	Column 1:	42.79330	1.0000	g/100cc
4.	n-Propanol	Column 2:	42.63439	1.0000	g/100cc

JG

ISP Forensic Services Blood Alcohol Report


Sample Name : 0.500 FN07031402  
 Laboratory : Meridian  
 Injection Date : Mar 22, 2017  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	44.96452	0.5003	g/100cc
2.	Ethanol	Column 2:	46.74831	0.5031	g/100cc
3.	n-Propanol	Column 1:	42.10326	1.0000	g/100cc
4.	n-Propanol	Column 2:	41.51231	1.0000	g/100cc

ISP Forensic Services Blood Alcohol Report

Sample Name : INTERNAL STANDARD BLANK  
 Laboratory : Meridian  
 Injection Date : Mar 22, 2017  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	0.00000	0.0000	g/100cc
2.	Ethanol	Column 2:	0.00000	0.0000	g/100cc
3.	n-Propanol	Column 1:	42.90772	1.0000	g/100cc
4.	n-Propanol	Column 2:	42.60557	1.0000	g/100cc

UG

S a m p l e S u m m a r y

Sequence table: C:\Chem32\1\Data\03-22-17\_CAL\03-22-17\_CAL 2017-03-22 10-22-49\03-22-17\_CAL.S  
 Data directory path: C:\Chem32\1\Data\03-22-17\_CAL\03-22-17\_CAL 2017-03-22 10-22-49\  
 Logbook: C:\Chem32\1\Data\03-22-17\_CAL\03-22-17\_CAL 2017-03-22 10-22-49\03-22-17\_CAL.LOG  
 Sequence start: 3/22/2017 10:37:25 AM  
 Sequence Operator: SYSTEM  
 Operator: SYSTEM

Method file name: C:\Chem32\1\Data\03-22-17\_CAL\03-22-17\_CAL 2017-03-22 10-22-49\ALCOHOL.M

Run #	Location #	Inj #	Sample Name	Sample Amt [g/100cc]	Multip.* Dilution	File name	Cal #	Cmp
1	1	1	0.050 FN06231406	-	1.0000	001F0101.D	*	4
2	2	1	0.100 FN06181501	-	1.0000	002F0201.D	*	4
3	3	1	0.200 FN07201502	-	1.0000	003F0301.D	*	4
4	4	1	0.300 FN06051501	-	1.0000	004F0401.D	*	4
5	5	1	0.500 FN07031402	-	1.0000	005F0501.D	*	4
6	6	1	INTERNAL STANDAR	-	1.0000	006F0601.D		2

JG

## VOLATILES DETERMINATION CASEFILE WORKSHEET

Laboratory No.: QC1-1

Analysis Date(s): 22 Mar 2017

	Column 1 FID A	Column 2 FID B	Column Precision	Mean Value	Over-all Mean	
Sample Results	0.0783	0.0802	0.0019	0.0792	0.0793	
(g/100cc)	0.0784	0.0806	0.0022	0.0795		

### Analysis Method

Refer to Blood Alcohol Method #1

### Instrument Information

*Instrument method is stored centrally.*

Refer to Instrument Method: ALCOHOL.M  
Hamilton Auto-Dilutor Serial Number:  
MD96BC1382/MD94AM10010

### Reporting of Results

Uncertainty of Measurement (UM%): 5.00%

Overall Mean (g/100cc)	Low	High	5% of Mean
0.079	0.075	0.083	0.004

	<b>Reported Result</b>	
	0.079	

*Calibration and control data are stored centrally.*

Issued: 12/30/2016

Volatiles BAC Calculation Spreadsheet Rev 4


Issuing Authority: Quality Manager

JG



ISP Forensic Services Blood Alcohol Report

Sample Name : QC1-1-A  
 Laboratory : Meridian  
 Injection Date : Mar 22, 2017  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	6.75568	0.0783	g/100cc
2.	Ethanol	Column 2:	6.69796	0.0802	g/100cc
3.	n-Propanol	Column 1:	41.09568	1.0000	g/100cc
4.	n-Propanol	Column 2:	40.78588	1.0000	g/100cc

JG

ISP Forensic Services Blood Alcohol Report

Sample Name : QC1-1-B  
 Laboratory : Meridian  
 Injection Date : Mar 22, 2017  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	6.81249	0.0784	g/100cc
2.	Ethanol	Column 2:	6.73674	0.0806	g/100cc
3.	n-Propanol	Column 1:	41.34389	1.0000	g/100cc
4.	n-Propanol	Column 2:	40.80205	1.0000	g/100cc

## VOLATILES DETERMINATION CASEFILE WORKSHEET

Laboratory No.: QC1-2

Analysis Date(s): 23 Mar 2017

	Column 1 FID A	Column 2 FID B	Column Precision	Mean Value	Over-all Mean	
Sample Results	0.0813	0.0839	0.0026	0.0826	0.0811	
(g/100cc)	0.0786	0.0806	0.0020	0.0796		

### Analysis Method

Refer to Blood Alcohol Method #1

### Instrument Information

*Instrument method is stored centrally.*

Refer to Instrument Method: ALCOHOL.M  
Hamilton Auto-Dilutor Serial Number:  
MD96BC1382/MD94AM10010

### Reporting of Results

Uncertainty of Measurement (UM%): 5.00%

Overall Mean (g/100cc)	Low	High	5% of Mean
0.081	0.076	0.086	0.005

	<b>Reported Result</b>	
	0.081	

*Calibration and control data are stored centrally.*

Issued: 12/30/2016


Volatiles BAC Calculation Spreadsheet Rev 4

Issuing Authority: Quality Manager

SG

ISP Forensic Services Blood Alcohol Report

Sample Name : QC1-2-A  
 Laboratory : Meridian  
 Injection Date : Mar 23, 2017  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	7.18279	0.0813	g/100cc
2.	Ethanol	Column 2:	7.06120	0.0839	g/100cc
3.	n-Propanol	Column 1:	42.00691	1.0000	g/100cc
4.	n-Propanol	Column 2:	40.90062	1.0000	g/100cc

JG

ISP Forensic Services Blood Alcohol Report

Sample Name : QC1-2-B  
 Laboratory : Meridian  
 Injection Date : Mar 23, 2017  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	6.70230	0.0786	g/100cc
2.	Ethanol	Column 2:	6.54951	0.0806	g/100cc
3.	n-Propanol	Column 1:	40.59548	1.0000	g/100cc
4.	n-Propanol	Column 2:	39.68166	1.0000	g/100cc

JG

## VOLATILES DETERMINATION CASEFILE WORKSHEET

Laboratory No.: QC2-1

Analysis Date(s): 22 Mar 2017

	Column 1 FID A	Column 2 FID B	Column Precision	Mean Value	Over-all Mean	
Sample Results	0.2081	0.2093	0.0012	0.2087	0.2098	
(g/100cc)	0.2102	0.2118	0.0016	0.2110		

### Analysis Method

Refer to Blood Alcohol Method #1

### Instrument Information

*Instrument method is stored centrally.*

Refer to Instrument Method: ALCOHOL.M  
Hamilton Auto-Dilutor Serial Number:  
MD96BC1382/MD94AM10010

### Reporting of Results

Uncertainty of Measurement (UM%): 5.00%

Overall Mean (g/100cc)	Low	High	5% of Mean
0.209	0.198	0.220	0.011

	<b>Reported Result</b>	
	0.209	

*Calibration and control data are stored centrally.*

Issued: 12/30/2016


Volatiles BAC Calculation Spreadsheet Rev 4

Issuing Authority: Quality Manager

JG

ISP Forensic Services Blood Alcohol Report

Sample Name : QC2-1-A  
 Laboratory : Meridian  
 Injection Date : Mar 22, 2017  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	18.05231	0.2081	g/100cc
2.	Ethanol	Column 2:	18.29752	0.2093	g/100cc
3.	n-Propanol	Column 1:	40.81908	1.0000	g/100cc
4.	n-Propanol	Column 2:	39.95520	1.0000	g/100cc

SG

ISP Forensic Services Blood Alcohol Report

Sample Name : QC2-1-B  
 Laboratory : Meridian  
 Injection Date : Mar 22, 2017  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	18.64679	0.2102	g/100cc
2.	Ethanol	Column 2:	18.90981	0.2118	g/100cc
3.	n-Propanol	Column 1:	41.72538	1.0000	g/100cc
4.	n-Propanol	Column 2:	40.79625	1.0000	g/100cc

JG



## VOLATILES DETERMINATION CASEFILE WORKSHEET

Laboratory No.: QC2-2

Analysis Date(s): 23 Mar 2017

	Column 1 FID A	Column 2 FID B	Column Precision	Mean Value	Over-all Mean	
Sample Results	0.2134	0.2161	0.0027	0.2147	0.2113	
(g/100cc)	0.2076	0.2084	0.0008	0.2080		

### Analysis Method

Refer to Blood Alcohol Method #1

### Instrument Information

*Instrument method is stored centrally.*

Refer to Instrument Method: ALCOHOL.M  
Hamilton Auto-Dilutor Serial Number:  
MD96BC1382/MD94AM10010

### Reporting of Results

Uncertainty of Measurement (UM%): 5.00%

Overall Mean (g/100cc)	Low	High	5% of Mean
0.211	0.200	0.222	0.011

	<b>Reported Result</b>	
	0.211	

*Calibration and control data are stored centrally.*

Issued: 12/30/2016


Volatiles BAC Calculation Spreadsheet Rev 4

Issuing Authority: Quality Manager

JF

ISP Forensic Services Blood Alcohol Report

Sample Name : QC2-2-A  
 Laboratory : Meridian  
 Injection Date : Mar 23, 2017  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	18.87840	0.2134	g/100cc
2.	Ethanol	Column 2:	19.11911	0.2161	g/100cc
3.	n-Propanol	Column 1:	41.61736	1.0000	g/100cc
4.	n-Propanol	Column 2:	40.39738	1.0000	g/100cc

JG

ISP Forensic Services Blood Alcohol Report

Sample Name : QC2-2-B  
 Laboratory : Meridian  
 Injection Date : Mar 23, 2017  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	18.50614	0.2076	g/100cc
2.	Ethanol	Column 2:	18.68274	0.2084	g/100cc
3.	n-Propanol	Column 1:	41.93214	1.0000	g/100cc
4.	n-Propanol	Column 2:	40.97992	1.0000	g/100cc

SG

## VOLATILES DETERMINATION CASEFILE WORKSHEET

Laboratory No.: 0.08 FN10281510

Analysis Date(s): 22 Mar 2017

	Column 1 FID A	Column 2 FID B	Column Precision	Mean Value	Over-all Mean	
Sample Results	0.0817	0.0833	0.0016	0.0825	0.0824	
(g/100cc)	0.0816	0.0832	0.0016	0.0824		

### Analysis Method

Refer to Blood Alcohol Method #1

### Instrument Information

*Instrument method is stored centrally.*

Refer to Instrument Method: ALCOHOL.M  
Hamilton Auto-Dilutor Serial Number:  
MD96BC1382/MD94AM10010

### Reporting of Results

Uncertainty of Measurement (UM%): 5.00%

Overall Mean (g/100cc)	Low	High	5% of Mean
0.082	0.077	0.087	0.005

	<b>Reported Result</b>  0.082	
--	-------------------------------------	--

*Calibration and control data are stored centrally.*

Issued: 12/30/2016


Volatiles BAC Calculation Spreadsheet Rev 4

Issuing Authority: Quality Manager

JG

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.08 FN10281510-A  
 Laboratory : Meridian  
 Injection Date : Mar 22, 2017  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	6.99458	0.0817	g/100cc
2.	Ethanol	Column 2:	6.89366	0.0833	g/100cc
3.	n-Propanol	Column 1:	40.73009	1.0000	g/100cc
4.	n-Propanol	Column 2:	40.26515	1.0000	g/100cc

JG

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.08 FN10281510-B  
 Laboratory : Meridian  
 Injection Date : Mar 22, 2017  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	7.07359	0.0816	g/100cc
2.	Ethanol	Column 2:	7.01517	0.0832	g/100cc
3.	n-Propanol	Column 1:	41.21565	1.0000	g/100cc
4.	n-Propanol	Column 2:	41.00610	1.0000	g/100cc

JG

ISP Forensic Services Blood Alcohol Report

Sample Name : INTERNAL STD BLK  
 Laboratory : Meridian  
 Injection Date : Mar 23, 2017  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	0.00000	0.0000	g/100cc
2.	Ethanol	Column 2:	0.00000	0.0000	g/100cc
3.	n-Propanol	Column 1:	41.75019	1.0000	g/100cc
4.	n-Propanol	Column 2:	40.67232	1.0000	g/100cc

UG

ISP Forensic Services Blood Alcohol Report

Sample Name : MIX VOL FN09231404  
 Laboratory : Meridian  
 Injection Date : Mar 22, 2017  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	6.87512	0.1618	g/100cc
2.	Ethanol	Column 2:	6.84842	0.1643	g/100cc
3.	n-Propanol	Column 1:	20.02605	1.0000	g/100cc
4.	n-Propanol	Column 2:	19.26696	1.0000	g/100cc

JG



ISP Forensic Services Blood Alcohol Report

Sample Name : INTERNAL STD BLK 1  
 Laboratory : Meridian  
 Injection Date : Mar 22, 2017  
 Method : ALCOHOL.M  
 Acq. Instrument: CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	0.00000	0.0000	g/100cc
2.	Ethanol	Column 2:	0.00000	0.0000	g/100cc
3.	n-Propanol	Column 1:	41.25090	1.0000	g/100cc
4.	n-Propanol	Column 2:	41.51320	1.0000	g/100cc

S a m p l e S u m m a r y

Sequence table: C:\Chem32\1\Data\03-01-17\_SAMPLES\03-22-17\_SAMPLES 2017-03-22 16-05-36\03-22-17\_SAMPLES.S  
 Data directory path: C:\Chem32\1\Data\03-01-17\_SAMPLES\03-22-17\_SAMPLES 2017-03-22 16-05-36\  
 Logbook: C:\Chem32\1\Data\03-01-17\_SAMPLES\03-22-17\_SAMPLES 2017-03-22 16-05-36\03-22-17\_SAMPLES.LOG  
 Sequence start: 3/22/2017 4:20:28 PM  
 Sequence Operator: SYSTEM  
 Operator: SYSTEM  
 Method file name: C:\Chem32\1\Data\03-01-17\_SAMPLES\03-22-17\_SAMPLES 2017-03-22 16-05-36\ALCOHOL.M

Run #	Location #	Inj #	Sample Name	Sample Amt [g/100cc]	Multip.* Dilution	File name	Cal #	Cmp
1	1	1	INTERNAL STD BLK	-	1.0000	001F0101.D		2
2	2	1	MIX VOL FN092314	-	1.0000	002F0201.D		10
3	3	1	QC1-1-A	-	1.0000	003F0301.D		4
4	4	1	QC1-1-B	-	1.0000	004F0401.D		4
5	5	1	0.08 FN10281510-	-	1.0000	005F0501.D		4
6	6	1	0.08 FN10281510-	-	1.0000	006F0601.D		4
7	7	1	M2017-1038-1-A	-	1.0000	007F0701.D		4
8	8	1	M2017-1038-1-B	-	1.0000	008F0801.D		4
9	9	1	M2017-1053-1-A	-	1.0000	009F0901.D		4
10	10	1	M2017-1053-1-B	-	1.0000	010F1001.D		4
11	11	1	M2017-1054-1-A	-	1.0000	011F1101.D		4
12	12	1	M2017-1054-1-B	-	1.0000	012F1201.D		4
13	13	1	M2017-1055-1-A	-	1.0000	013F1301.D		4
14	14	1	M2017-1055-1-B	-	1.0000	014F1401.D		4
15	15	1	M2017-1056-1-A	-	1.0000	015F1501.D		4
16	16	1	M2017-1056-1-B	-	1.0000	016F1601.D		4
17	17	1	M2017-1057-1-A	-	1.0000	017F1701.D		4
18	18	1	M2017-1057-1-B	-	1.0000	018F1801.D		4
19	19	1	M2017-1058-1-A	-	1.0000	019F1901.D		2
20	20	1	M2017-1058-1-B	-	1.0000	020F2001.D		2
21	21	1	M2017-1077-1-A	-	1.0000	021F2101.D		4
22	22	1	M2017-1077-1-B	-	1.0000	022F2201.D		4
23	23	1	M2017-1095-1-A	-	1.0000	023F2301.D		2
24	24	1	M2017-1095-1-B	-	1.0000	024F2401.D		2
25	25	1	QC2-1-A	-	1.0000	025F2501.D		4
26	26	1	QC2-1-B	-	1.0000	026F2601.D		4
27	27	1	M2017-1096-1-A	-	1.0000	027F2701.D		2
28	28	1	M2017-1096-1-B	-	1.0000	028F2801.D		2
29	29	1	M2017-1110-1-A	-	1.0000	029F2901.D		6
30	30	1	M2017-1110-1-B	-	1.0000	030F3001.D		6
31	31	1	M2017-1115-1-A	-	1.0000	031F3101.D		2
32	32	1	M2017-1115-1-B	-	1.0000	032F3201.D		2
33	33	1	M2017-1116-1-A	-	1.0000	033F3301.D		2
34	34	1	M2017-1116-1-B	-	1.0000	034F3401.D		2
35	35	1	M2017-1117-1-A	-	1.0000	035F3501.D		4
36	36	1	M2017-1117-1-B	-	1.0000	036F3601.D		4
37	37	1	M2017-1118-1-A	-	1.0000	037F3701.D		4
38	38	1	M2017-1118-1-B	-	1.0000	038F3801.D		4
39	39	1	M2017-1119-1-A	-	1.0000	039F3901.D		4
40	40	1	M2017-1119-1-B	-	1.0000	040F4001.D		4
41	41	1	M2017-1122-1-A	-	1.0000	041F4101.D		2
42	42	1	M2017-1122-1-B	-	1.0000	042F4201.D		2
43	43	1	M2017-1148-1-A	-	1.0000	043F4301.D		2

*JG*

Run #	Location #	Inj #	Sample Name	Sample Amt [g/100cc]	Multip.* Dilution	File name	Cal #	# Cmp
44	44	1	M2017-1148-1-B	-	1.0000	044F4401.D		2
45	45	1	M2017-1169-1-A	-	1.0000	045F4501.D		4
46	46	1	M2017-1169-1-B	-	1.0000	046F4601.D		4
47	47	1	QC1-2-A	-	1.0000	047F4701.D		4
48	48	1	QC1-2-B	-	1.0000	048F4801.D		4
49	49	1	M2017-1170-1-A	-	1.0000	049F4901.D		4
50	50	1	M2017-1170-1-B	-	1.0000	050F5001.D		4
51	51	1	M2017-1171-1-A	-	1.0000	051F5101.D		4
52	52	1	M2017-1171-1-B	-	1.0000	052F5201.D		4
53	53	1	M2017-1173-1-A	-	1.0000	053F5301.D		4
54	54	1	M2017-1173-1-B	-	1.0000	054F5401.D		4
55	55	1	M2017-1174-1-A	-	1.0000	055F5501.D		4
56	56	1	M2017-1174-1-B	-	1.0000	056F5601.D		4
57	57	1	M2017-1188-1-A	-	1.0000	057F5701.D		4
58	58	1	M2017-1188-1-B	-	1.0000	058F5801.D		4
59	59	1	M2017-1190-1-A	-	1.0000	059F5901.D		4
60	60	1	M2017-1190-1-B	-	1.0000	060F6001.D		4
61	61	1	P2017-0449-1-A	-	1.0000	061F6101.D		4
62	62	1	P2017-0449-1-B	-	1.0000	062F6201.D		4
63	63	1	P2017-0501-1-A	-	1.0000	063F6301.D		4
64	64	1	P2017-0501-1-B	-	1.0000	064F6401.D		4
65	65	1	P2017-0548-1-A	-	1.0000	065F6501.D		4
66	66	1	P2017-0548-1-B	-	1.0000	066F6601.D		4
67	67	1	QC2-2-A	-	1.0000	067F6701.D		4
68	68	1	QC2-2-B	-	1.0000	068F6801.D		4
69	69	1	INTERNAL STD BLK	-	1.0000	069F6901.D		2

Method file name: C:\Chem32\1\Data\03-01-17\_SAMPLES\03-22-17\_SAMPLES 2017-03-22 16-05-36  
 \SHUTDOWN.M

Run #	Location #	Inj #	Sample Name	Sample Amt [g/100cc]	Multip.* Dilution	File name	Cal #	# Cmp
70	70	1	EMPTY	-	1.0000	070F7001.D		0