

**Quantitative Analysis for Ethanol & Qualitative Analysis for Other Volatiles**

*Analytical Method(s): 1.0*

**Device: Hamilton MICROLAB 600 Liquid Processor/Dilutor Serial Number: ML600HC11378**

**Volatiles Quality Assurance Controls**

**Run Date(s): 8/3/18-8/4/18**  
Calibration Date: 08/01/18

Control level	Expiration	Lot #	Target Value	Acceptable Range	Overall Results
Level 1	Jan-22	1801036	0.0731	0.0731-0.0893	0.0758 g/100cc 0.0798 g/100cc g/100cc
Level 2	Mar-22	1803028	0.2035	0.1832-0.2238	0.2019 g/100cc 0.2120 g/100cc g/100cc
<b>Multi-Component mixture:</b>		<b>Exp date: Sept 2020</b>	<b>Lot #</b>	<b>FN06041502</b>	<b>OK</b>
<b>Curve Fit:</b>		<b>Column 1</b>	<b>Column 2</b>	<b>0.99991</b>	<b>0.99981</b>

**Ethanol Calibration Reference Material**

Calibrator level	Expiration	Cerilliant Lot #	Target Value	Acceptable Range	Column 1	Column 2	Precision	Mean
0.050	Jul-19	FN06231406	0.050	0.045 - 0.055	0.0536	0.0550	0.0014	0.0543
0.080			0.080	0.072 - 0.088			0	#DIV/0!
0.100	Aug-21	FN08101601	0.100	0.090 - 0.110	0.0975	0.0974	0.0001	0.0974
0.200	Apr-21	FN03301601	0.200	0.180 - 0.220	0.1985	0.1970	0.0015	0.1977
0.300	Feb-21	FN02121601	0.300	0.270 - 0.330	0.2991	0.2985	0.0006	0.2988
0.400			0.400	0.360 - 0.440			0	#DIV/0!
0.500	Sep-21	FN08031602	0.500	0.450 - 0.550	0.5013	0.5021	0.0008	0.5017

**Aqueous Controls**

Control level	Expiration	Cerilliant Lot #	Target Value	Acceptable Range	Overall Results
0.080	<del>Nov-20</del> May-22	FN04171701	0.08000	0.076 - 0.084	0.078 g/100cc

~Any information on this document can be changed for laboratory use, except for the precision and mean determination formulas.

Issued: 4/22/2015

Volatiles QA/QC data spreadsheet Rev 5


Issuing Authority: Quality Manager



**Worksheet: 2616**


LAB CASE	ITEM	TASK ID	DESCRIPTION
M2018-3615	4	122369	Alcohol Analysis
M2018-3639	1	121932	Alcohol Analysis
M2018-3640	1	121936	Alcohol Analysis
M2018-3641	1	121937	Alcohol Analysis
M2018-3641	2	121941	Alcohol Analysis
M2018-3641	3	121945	Alcohol Analysis
M2018-3642	1	122021	Alcohol Analysis
M2018-3643	1	122027	Alcohol Analysis
M2018-3644	1	122032	Alcohol Analysis
M2018-3648	1	122062	Alcohol Analysis
M2018-3649	1	122067	Alcohol Analysis
M2018-3650	1	122071	Alcohol Analysis
M2018-3651	1	122081	Alcohol Analysis
M2018-3652	1	122085	Alcohol Analysis
M2018-3653	1	122095	Alcohol Analysis
M2018-3658	1	122171	Alcohol Analysis
M2018-3659	1	122172	Alcohol Analysis
M2018-3660	1	122176	Alcohol Analysis
M2018-3673	1	122244	Alcohol Analysis
M2018-3678	1	122349	Alcohol Analysis
M2018-3679	1	122350	Alcohol Analysis
M2018-3694	1	122368	Alcohol Analysis
M2018-3702	1	122403	Alcohol Analysis

8/6/2018



**Worklist: 2616**

<u>LAB CASE</u>	<u>ITEM</u>	<u>TASK ID</u>	<u>DESCRIPTION</u>
M2018-3708	1	122416	Alcohol Analysis
M2018-3709	1	122420	Alcohol Analysis
M2018-3759	1	122517	Alcohol Analysis
M2018-3762	1	122561	Alcohol Analysis
M2018-3773	1	122602	Alcohol Analysis
P2018-1239	1	121180	Alcohol Analysis



96

=====  
Calibration Table  
=====

-----  
General Calibration Setting  
-----

Calib. Data Modified : Wednesday, August 01, 2018 11:12:59 AM  
Signals calculated separately : No

Rel. Reference Window : 0.000 %  
Abs. Reference Window : 0.100 min  
Rel. Non-ref. Window : 0.000 %  
Abs. Non-ref. Window : 0.100 min  
Uncalibrated Peaks : not reported  
Partial Calibration : Yes, identified peaks are recalibrated  
Correct All Ret. Times : No, only for identified peaks

Curve Type : Linear  
Origin : Ignored  
Weight : Equal

Recalibration Settings : Average Response : Average Retention Time :  
Average all calibrations  
Floating Average New 75%

Calibration Report Options :

Printout of recalibrations within a sequence:  
Calibration Table after Recalibration  
Normal Report after Recalibration  
If the sequence is done with bracketing:  
Results of first cycle (ending previous bracket)

Default Sample ISTD Information (if not set in sample table):  
ISTD ISTD Amount Name # [g/100cc]  
-----|-----  
1 1.00000 n-propanol  
2 1.00000 n-propanol





-----  
Signal Details  
-----





Signal 1: FID1 A, Front Signal  
Signal 2: FID2 B, Back Signal

-----  
Overview Table  
-----












26




=====

50



ISP Forensic Services Blood Alcohol Report

Sample Name : 0.050 FN06231406  
 Laboratory : Meridian  
 Injection Date : Aug 1, 2018  
 Method : ALCOHOL.M  
 Acq. Instrument : CN11180014-CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	4.50304	0.0536	g/100cc
2.	Ethanol	Column 2:	4.60733	0.0550	g/100cc
3.	n-Propanol	Column 1:	42.45125	1.0000	g/100cc
4.	n-Propanol	Column 2:	43.90464	1.0000	g/100cc

26

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.100 FN08101601  
 Laboratory : Meridian  
 Injection Date : Aug 1, 2018  
 Method : ALCOHOL.M  
 Acq. Instrument : CN1180014-CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	9.55978	0.0975	g/100cc
2.	Ethanol	Column 2:	9.83303	0.0974	g/100cc
3.	n-Propanol	Column 1:	50.31847	1.0000	g/100cc
4.	n-Propanol	Column 2:	52.12196	1.0000	g/100cc

JK

ISP Forensic Services Blood Alcohol Report


Sample Name : 0.200 FN03301601  
 Laboratory : Meridian  
 Injection Date : Aug 1, 2018  
 Method : ALCOHOL.M  
 Acq. Instrument : CN1180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	19.01844	0.1985	g/100cc
2.	Ethanol	Column 2:	19.78436	0.1970	g/100cc
3.	n-Propanol	Column 1:	49.67036	1.0000	g/100cc
4.	n-Propanol	Column 2:	51.30992	1.0000	g/100cc

ISP Forensic Services Blood Alcohol Report


Sample Name : 0.300 FN02121601  
 Laboratory : Meridian  
 Injection Date : Aug 1, 2018  
 Method : ALCOHOL.M  
 Acq. Instrument : CN1180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	28.44850	0.2991	g/100cc
2.	Ethanol	Column 2:	29.87221	0.2985	g/100cc
3.	n-Propanol	Column 1:	49.47174	1.0000	g/100cc
4.	n-Propanol	Column 2:	50.96388	1.0000	g/100cc

ISP Forensic Services Blood Alcohol Report


Sample Name : 0.500 FN08031602  
 Laboratory : Meridian  
 Injection Date : Aug 1, 2018  
 Method : ALCOHOL.M  
 Acq. Instrument : CN1180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	47.58224	0.5013	g/100cc
2.	Ethanol	Column 2:	50.32014	0.5021	g/100cc
3.	n-Propanol	Column 1:	49.49506	1.0000	g/100cc
4.	n-Propanol	Column 2:	50.89685	1.0000	g/100cc

ISP Forensic Services Blood Alcohol Report

Sample Name : INTERNAL STANDARD BLANK  
 Laboratory : Meridian  
 Injection Date : Aug 1, 2018  
 Method : ALCOHOL.M  
 Acq. Instrument : CN1180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	0.00000	0.0000	g/100cc
2.	Ethanol	Column 2:	0.00000	0.0000	g/100cc
3.	n-Propanol	Column 1:	48.96428	1.0000	g/100cc
4.	n-Propanol	Column 2:	50.36459	1.0000	g/100cc

S a m p l e    s u m m a r y

Sequence table: C:\Chem32\1\Data\08-01-18\_CAL\08-01-18\_CAL 2018-08-01 10-07-47\08-01-18\_CAL.S  
Data directory path: C:\Chem32\1\Data\08-01-18\_CAL\08-01-18\_CAL 2018-08-01 10-07-47\  
Logbook: C:\Chem32\1\Data\08-01-18\_CAL\08-01-18\_CAL 2018-08-01 10-07-47\08-01-18\_CAL.LOG  
Sequence start: 8/1/2018 10:22:27 AM  
Sequence Operator: SYSTEM  
Operator: SYSTEM

Method file name: C:\Chem32\1>Data\08-01-18\_CAL\08-01-18\_CAL 2018-08-01 10-07-47\ALCOHOL.M

Run Location Inj	#	Sample Name	Sample Amt	Multlp.*	File name	Cal #	Cmp #
1 1	1	0.050 FN06231406	1.0000	001F0101.D		4	
2 2	1	0.100 FN08101601	1.0000	002F0201.D		4	
3 3	1	0.200 FN03301601	1.0000	003F0301.D		4	
4 4	1	0.300 FN02121601	1.0000	004F0401.D		4	
5 5	1	0.500 FN08031602	1.0000	005F0501.D		4	*
6 6	1	INTERNAL STANDARD	1.0000	006F0601.D		2	

✓

ISP Forensic Services Blood Alcohol Report


Sample Name : INTERNAL STD BLK 1

Laboratory : Meridian

Injection Date : Aug 3, 2018

Method : ALCOHOL.M

Acq. Instrument : CN1180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	0.00000	0.0000	g/100cc
2.	Ethanol	Column 2:	0.00000	0.0000	g/100cc
3.	n-Propanol	Column 1:	48.08860	1.0000	g/100cc
4.	n-Propanol	Column 2:	50.06171	1.0000	g/100cc



ISP Forensic Services Blood Alcohol Report


Sample Name : MIX VOL FN06041502

Laboratory : Meridian

Injection Date : Aug 3, 2018

Method : ALCOHOL.M

Acq. Instrument : CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	7.99926	0.1618	g/100cc
2.	Ethanol	Column 2:	8.37904	0.1666	g/100cc
3.	n-Propanol	Column 1:	25.57914	1.0000	g/100cc
4.	n-Propanol	Column 2:	25.74659	1.0000	g/100cc

2c

*Calibration and control data are stored centrally.*


Reporting of Results						
Overall Mean (g/100cc)						
0.075						
Low						
High						
5% of Mean						
0.004						
Reported Result						
0.075						
Uncertainty of Measurement (UM%): 5.00%						
Refer to Instrument Method: ALCOHOL.M						
Hamilton Auto-Dilutor Serial Number: ML600HC11378						
Instrument Information						
<i>Instrument method is stored centrally.</i>						
Refer to Blood Alcohol Method #1						
Analysis Method						
Sample Results						
0.0752						
0.0750						
0.0760						
0.0770						
0.0008						
0.0020						
0.0756						
0.0760						
0.0758						
Column 1	FID A	Column 2	FID B	Column	Precision	Mean Value
Over-all Mean						

Laboratory No.: QC1-1      Analysis Date(s): 03 Aug 2018

**VOLATILES DETERMINATION CASEFILE WORKSHEET**

ISP Forensic Services Blood Alcohol Report


Sample Name : Q1-1-A  
 Laboratory : Meridian  
 Injection Date : Aug 3, 2018  
 Method : ALCOHOL.M  
 Acq. Instrument : CN1180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	7.22514	0.0752	g/100cc
2.	Ethanol	Column 2:	7.40212	0.0760	g/100cc
3.	n-Propanol	Column 1:	49.05032	1.0000	g/100cc
4.	n-Propanol	Column 2:	50.56157	1.0000	g/100cc

ISP Forensic Services Blood Alcohol Report

Sample Name : Q1-1-B  
 Laboratory : Meridian  
 Injection Date : Aug 3, 2018  
 Method : ALCOHOL.M  
 Acq. Instrument : CN1180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	7.21546	0.0750	g/100cc
2.	Ethanol	Column 2:	7.50752	0.0770	g/100cc
3.	n-Propanol	Column 1:	49.10493	1.0000	g/100cc
4.	n-Propanol	Column 2:	50.60512	1.0000	g/100cc

**VOLATILES DETERMINATION CASEFILE WORKSHEET**

Laboratory No.: 08 FN04171701 Analysis Date(s): 03 Aug 2018

Sample Results	Column 1 FID A	Column 2 FID B	Column Precision	Mean Value	Over-all Mean
(g/100cc)	0.0785	0.0798	0.0013	0.0791	0.0789
	0.0782	0.0793	0.0011	0.0787	

**Analysis Method**

Refer to Blood Alcohol Method #1

**Instrument Information**

Refer to Instrument Method: ALCOHOL.M  
Hamilton Auto-Dilutor Serial Number: ML600HC1378

Reporting of Results				Uncertainty of Measurement (UM%): 5.00%	
Overall Mean (g/100cc)	Low	High	5% of Mean		
0.078	0.074	0.082	0.004		
Reported Result			0.078		

*Calibration and control data are stored centrally.*

ISP Forensic Services Blood Alcohol Report


Sample Name : 0.08 FN04171701-A

Laboratory : Meridian

Injection Date : Aug 3, 2018

Method : ALCOHOL.M


Acq. Instrument : CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	7.41553	0.0782	g/100cc
2.	Ethanol	Column 2:	7.60016	0.0793	g/100cc
3.	n-Propanol	Column 1:	48.45267	1.0000	g/100cc
4.	n-Propanol	Column 2:	49.71123	1.0000	g/100cc

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.08 FN04171701-B  
 Laboratory : Meridian  
 Injection Date : Aug 3, 2018  
 Method : ALCOHOL.M  
 Acq. Instrument : CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	7.50238	0.0785	g/100cc
2.	Ethanol	Column 2:	7.72650	0.0798	g/100cc
3.	n-Propanol	Column 1:	48.81833	1.0000	g/100cc
4.	n-Propanol	Column 2:	50.21921	1.0000	g/100cc

**VOLATILES DETERMINATION CASEFILE WORKSHEET**

Laboratory No.: QC2-1

Analysis Date(s): 03 Aug 2018

Sample Results	Column 1 FID A	Column 2 FID B	Column Precision	Mean Value	Over-all Mean	
					0.2019	0.2019
(g/100cc)	0.2003	0.2024	0.0021	0.2013		
Sample Results	0.2022	0.2029	0.0007	0.2025		

**Analysis Method**

Refer to Blood Alcohol Method #1

**Instrument Information**

Refer to Instrument Method: ALCOHOL.M  
Hamilton Auto-Dilutor Serial Number: ML600HC11378


Reporting of Results				Uncertainty of Measurement (UM%): 5.00%	
Overall Mean (g/100cc)	Low	High	5% of Mean		
0.201	0.190	0.212	0.011		
Reported Result			0.201		

Calibration and control data are stored centrally.



ISP Forensic Services Blood Alcohol Report


Sample Name : QC2-1-A  
 Laboratory : Meridian  
 Injection Date : Aug 3, 2018  
 Method : ALCOHOL.M  
 Acq. Instrument : CN1180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	19.19377	0.2022	g/100cc
2.	Ethanol	Column 2:	19.95613	0.2029	g/100cc
3.	n-Propanol	Column 1:	49.21776	1.0000	g/100cc
4.	n-Propanol	Column 2:	50.23718	1.0000	g/100cc

ISP Forensic Services Blood Alcohol Report

Sample Name : QC2-1-B  
 Laboratory : Meridian  
 Injection Date : Aug 3, 2018  
 Method : ALCOHOL.M  
 Acq. Instrument : CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	18.85115	0.2003	g/100cc
2.	Ethanol	Column 2:	19.71850	0.2024	g/100cc
3.	n-Propanol	Column 1:	48.78908	1.0000	g/100cc
4.	n-Propanol	Column 2:	49.75737	1.0000	g/100cc

**VOLATILES DETERMINATION CASEFILE WORKSHEET**

Laboratory No.: QC1-2

Analysis Date(s): 03 Aug 2018

	Column 1 FID A	Column 2 FID B	Column Precision	Mean Value	Over-all Mean
Sample Results	0.0790	0.0803	0.0013	0.0796	0.0798
(g/100cc)	0.0795	0.0804	0.0009	0.0799	

**Analysis Method**

Refer to Blood Alcohol Method #1

**Instrument Information**

*Instrument method is stored centrally.*

Refer to Instrument Method: ALCOHOL.M  
Hamilton Auto-Dilutor Serial Number: ML600HC11378


Reporting of Results				Uncertainty of Measurement (UM%): 5.00%	
Overall Mean (g/100cc)	Low	High	5% of Mean	0.079	
0.079	0.075	0.083	0.004	Reported Result	
0.079			0.079		

*Calibration and control data are stored centrally.*

26

ISP Forensic Services Blood Alcohol Report


Sample Name : Q1-2-A  
 Laboratory : Meridian  
 Injection Date : Aug 3, 2018  
 Method : ALCOHOL.M  
 Instrument : CN1180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	7.65681	0.0790	g/100cc
2.	Ethanol	Column 2:	7.83465	0.0803	g/100cc
3.	n-Propanol	Column 1:	49.54123	1.0000	g/100cc
4.	n-Propanol	Column 2:	50.58554	1.0000	g/100cc

ISP Forensic Services Blood Alcohol Report

Sample Name : Q11-2-B  
 Laboratory : Meridian  
 Injection Date : Aug 3, 2018  
 Method : ALCOHOL.M  
 Acq. Instrument : CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	7.67616	0.0795	g/100cc
2.	Ethanol	Column 2:	7.82422	0.0804	g/100cc
3.	n-Propanol	Column 1:	49.34435	1.0000	g/100cc
4.	n-Propanol	Column 2:	50.47264	1.0000	g/100cc

26

*Calibration and control data are stored centrally.*


Reporting of Results						
Overall Mean (g/100cc)			Low	High	5% of Mean	
0.212			0.201	0.223	0.011	
Reported Result			0.212			
Refer to Blood Alcohol Method #1 Refer to Instrument Method: ALCOHOL.M Hamilton Auto-Dilutor Serial Number: ML600HC11378						
Instrument Information						
<i>Instrument method is stored centrally.</i>						
Analysis Method						
Refer to Blood Alcohol Method #1						
(g/100cc)						
Column 1	Column 2	Column Precision	Mean Value	Over-all Mean	Sample Results	
FID A	FID B	0.0001	0.2117	0.2120	0.2117	
0.2116	0.2130	0.0014	0.2123	0.2120	0.2130	

Laboratory No.: QC2-2      Analysis Date(s): 04 Aug 2018

**VOLATILES DETERMINATION CASEFILE WORKSHEET**

ISP Forensic Services Blood Alcohol Report

Sample Name : QC2-2-A  
 Laboratory : Meridian  
 Injection Date : Aug 4, 2018  
 Method : ALCOHOL.M  
 Acq. Instrument : CN1180014-CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	19.88034	0.2117	g/100cc
2.	Ethanol	Column 2:	20.57656	0.2118	g/100cc
3.	n-Propanol	Column 1:	48.71247	1.0000	g/100cc
4.	n-Propanol	Column 2:	49.60549	1.0000	g/100cc

2

ISP Forensic Services Blood Alcohol Report

Sample Name : Q2-2-B  
 Laboratory : Meridian  
 Injection Date : Aug 4, 2018  
 Method : ALCOHOL.M  
 Acq. Instrument : CN1180014-CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	20.14013	0.2116	g/100cc
2.	Ethanol	Column 2:	21.00805	0.2130	g/100cc
3.	n-Propanol	Column 1:	49.37965	1.0000	g/100cc
4.	n-Propanol	Column 2:	50.35028	1.0000	g/100cc



ISP Forensic Services Blood Alcohol Report

Sample Name : INTERNAL STD BLK  
 Laboratory : Meridian  
 Injection Date : Aug 4, 2018  
 Method : ALCOHOL.M  
 Acq. Instrument : CN1180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	0.00000	0.0000	g/100cc
2.	Ethanol	Column 2:	0.00000	0.0000	g/100cc
3.	n-Propanol	Column 1:	47.77999	1.0000	g/100cc
4.	n-Propanol	Column 2:	48.53879	1.0000	g/100cc

Sample Summary

Sequence table: C:\chem32\1\data\08-03-18\_SAMPLES\08-03-18\_SAMPLES 2018-08-03 15-17-27\08-03-18\_SAMPLES.S  
 Data directory path: C:\chem32\1\data\08-03-18\_SAMPLES\08-03-18\_SAMPLES 2018-08-03 15-17-27\08-03-18\_SAMPLES.LOG  
 Logbook: C:\chem32\1\data\08-03-18\_SAMPLES\08-03-18\_SAMPLES 2018-08-03 15-17-27\08-03-18\_SAMPLES.S  
 Sequence start: 8/3/2018 3:32:12 PM  
 Sequence Operator: SYSTEM  
 Operator: SYSTEM  
 Method file name: C:\chem32\1\data\08-03-18\_SAMPLES\08-03-18\_SAMPLES 2018-08-03 15-17-27\ALCOHOL.M

Run Location Inj #	Sample Name	Sample Amt [g/100cc]	Multip. Dilution	File name	Cal #	Cmp
--------------------	-------------	----------------------	------------------	-----------	-------	-----

1	INTERNAL STD BLK	1.0000	001F0101.D		2	
2	MIX VOL FN060415	1.0000	002F0201.D		10	
3	QC1-1-A	1.0000	003F0301.D		4	
4	QC1-1-B	1.0000	004F0401.D		4	
5	0.08 FN04171701	1.0000	005F0501.D		4	
6	0.08 FN04171701	1.0000	006F0601.D		4	
7	M2018-3639-1-A	1.0000	007F0701.D		2	
8	M2018-3639-1-B	1.0000	008F0801.D		2	
9	M2018-3640-1-A	1.0000	009F0901.D		6	
10	M2018-3640-1-B	1.0000	010F1001.D		6	
11	M2018-3641-1-A	1.0000	011F1101.D		4	
12	M2018-3641-1-B	1.0000	012F1201.D		4	
13	M2018-3641-2-A	1.0000	013F1301.D		2	
14	M2018-3641-2-B	1.0000	014F1401.D		2	
15	M2018-3641-3-A	1.0000	015F1501.D		2	
16	M2018-3641-3-B	1.0000	016F1601.D		2	
17	M2018-3644-1-A	1.0000	017F1701.D		6	
18	M2018-3644-1-B	1.0000	018F1801.D		6	
19	M2018-3648-1-A	1.0000	019F1901.D		2	
20	M2018-3648-1-B	1.0000	020F2001.D		2	
21	M2018-3649-1-A	1.0000	021F2101.D		4	
22	M2018-3649-1-B	1.0000	022F2201.D		4	
23	M2018-3650-1-A	1.0000	023F2301.D		2	
24	M2018-3650-1-B	1.0000	024F2401.D		2	
25	QC2-1-A	1.0000	025F2501.D		4	
26	QC2-1-B	1.0000	026F2601.D		4	
27	M2018-3651-1-A	1.0000	027F2701.D		6	
28	M2018-3651-1-B	1.0000	028F2801.D		5	
29	M2018-3652-1-A	1.0000	029F2901.D		6	
30	M2018-3652-1-B	1.0000	030F3001.D		4	
31	M2018-3653-1-A	1.0000	031F3101.D		4	
32	M2018-3653-1-B	1.0000	032F3201.D		4	
33	M2018-3658-1-A	1.0000	033F3301.D		4	
34	M2018-3658-1-B	1.0000	034F3401.D		4	
35	M2018-3659-1-A	1.0000	035F3501.D		2	
36	M2018-3659-1-B	1.0000	036F3601.D		2	
37	M2018-3660-1-A	1.0000	037F3701.D		4	
38	M2018-3660-1-B	1.0000	038F3801.D		4	
39	M2018-3673-1-A	1.0000	039F3901.D		6	
40	M2018-3673-1-B	1.0000	040F4001.D		6	
41	M2018-3615-1-A	1.0000	041F4101.D		2	
42	M2018-3615-1-B	1.0000	042F4201.D		2	
43	M2018-3642-1-A	1.0000	043F4301.D		4	



Run Location Inj	#	Sample Name	Sample Amt	Multipl.*	File name	Cal #	Cmp #
44	44	1 M2018-3642-1-B	1.0000	-	044F4401.D	4	4
45	45	1 M2018-3643-1-A	1.0000	-	045F4501.D	4	4
46	46	1 M2018-3643-1-B	1.0000	-	046F4601.D	4	4
47	47	1 QCI-2-A	1.0000	-	047F4701.D	4	4
48	48	1 QCI-2-B	1.0000	-	048F4801.D	4	4
49	49	1 M2018-3678-1-A	1.0000	-	049F4901.D	4	4
50	50	1 M2018-3678-1-B	1.0000	-	050F5001.D	4	4
51	51	1 M2018-3679-1-A	1.0000	-	051F5101.D	4	4
52	52	1 M2018-3679-1-B	1.0000	-	052F5201.D	4	4
53	53	1 M2018-3694-1-A	1.0000	-	053F5301.D	5	4
54	54	1 M2018-3694-1-B	1.0000	-	054F5401.D	4	4
55	55	1 M2018-3702-1-A	1.0000	-	055F5501.D	4	4
56	56	1 M2018-3702-1-B	1.0000	-	056F5601.D	4	4
57	57	1 M2018-3708-1-A	1.0000	-	057F5701.D	4	4
58	58	1 M2018-3708-1-B	1.0000	-	058F5801.D	4	4
59	59	1 M2018-3709-1-A	1.0000	-	059F5901.D	6	4
60	60	1 M2018-3709-1-B	1.0000	-	060F6001.D	4	4
61	61	1 M2018-3759-1-A	1.0000	-	061F6101.D	4	4
62	62	1 M2018-3759-1-B	1.0000	-	062F6201.D	4	4
63	63	1 M2018-3762-1-A	1.0000	-	063F6301.D	6	6
64	64	1 M2018-3762-1-B	1.0000	-	064F6401.D	6	6
65	65	1 M2018-3773-1-A	1.0000	-	065F6501.D	4	4
66	66	1 M2018-3773-1-B	1.0000	-	066F6601.D	5	5
67	67	1 P2018-1239-1-A	1.0000	-	067F6701.D	2	2
68	68	1 P2018-1239-1-B	1.0000	-	068F6801.D	2	2
69	69	1 QCI-2-A	1.0000	-	069F6901.D	4	4
70	70	1 QCI-2-B	1.0000	-	070F7001.D	4	4
71	71	1 INTERNAL STD BLK	-	-	071F7101.D	2	2

Method file name:

C:\Chem32\1\Data\08-03-18\_SAMPLES\08-03-18\_SAMPLES 2018-08-03 15-17-27\SHUTDOWN.M

Run Location Inj	#	Sample Name	Sample Amt	Multipl.*	File name	Cal #	Cmp #
72	72	1 EMPTY	-	-	072F7201.D	0	0

Due to a gas line failure, the following samples were previously extracted on 8/1/18 and re-extracted on 8/3/18.

- M2018-3639-1
- M2018-3640-1
- M2018-3641-1
- M2018-3641-2
- M2018-3641-3
- M2018-3644-1