

Quantitative Analysis for Ethanol & Qualitative Analysis for Other Volatiles

Analytical Method(s): 1.0

Device: Hamilton MICROLAB 600A Liquid Processor/Dilutor Serial Number: ML600HC11378

Volatiles Quality Assurance Controls Run Date(s): 3/28/19

Calibration Date: 3/28/19

Control level	Expiration	Lot #	Target Value	Acceptable Range	Overall Results
Level 1	Jan-22	1801036	0.0812	0.0731-0.0893	0.0788 g/100cc
					0.0829 g/100cc
Level 2	Mar-22	1803028	0.2035	0.1832-0.2238	0.2068 g/100cc
					g/100cc
Multi-Component mixture:		Sep-20	Lot #	FN06041502	ok
Curve Fit:		Column 1	1.00000	Column 2	0.99999

Ethanol Calibration Reference Material

Calibrator level	Target Value	Acceptable Range	Column 1	Column 2	Precision	Mean
50	0.050	0.045 - 0.055	0.0501	0.0512	0.0011	0.0506
100	0.100	0.090 - 0.110	0.0994	0.0989	0.0005	0.0991
200	0.200	0.180 - 0.220	0.2001	0.1992	0.0009	0.1996
300	0.300	0.270 - 0.330	0.3007	0.3006	0.0001	0.3006
500	0.500	0.450 - 0.550	0.4996	0.5000	0.0004	0.4998

Aqueous Controls

Control level	Target Value	Acceptable Range	Overall Results
80	0.080	0.076 - 0.084	0.081 g/100cc


REVIEWED
By Rachel Cutler at 3:55 pm, Mar 29, 2019

jc

12

Worklist: 3177

LAB CASE	ITEM	TASK ID	DESCRIPTION
M2019-1361	1	145581	Alcohol Analysis
M2019-1372	1	145652	Alcohol Analysis
M2019-1373	1	145656	Alcohol Analysis
M2019-1383	1	145697	Alcohol Analysis
M2019-1399	1	145963	Alcohol Analysis
M2019-1402	1	146049	Alcohol Analysis
M2019-1403	1	146050	Alcohol Analysis
M2019-1404	1	146051	Alcohol Analysis
M2019-1406	1	146059	Alcohol Analysis
M2019-1414	1	146097	Alcohol Analysis
M2019-1415	2	146099	Alcohol Analysis



JK

Overview Table

Signal 1: FID1 A, Front Signal
Signal 2: FID2 B, Back Signal

Signal Details

ISTD Amount	ISTD Name	# [g/100cc]
1.00000	n-propanol	1
1.00000	n-propanol	2

Default Sample ISTD Information (if not set in sample table):

Printout of recalibrations within a sequence:
 Calibration Table after Recalibration
 Normal Report after Recalibration
 If the sequence is done with bracketing:
 Results of first cycle (ending previous bracket)

Calibration Report Options :

Recalibration Settings:
 Average Response : Average all recalibrations
 Average Retention Time: Floating Average New 75%

Curve Type : Linear
 Origin : Ignored
 Weight : Equal

Rel. Reference Window : 0.000 %
 Abs. Reference Window : 0.100 min
 Rel. Non-ref. Window : 0.000 %
 Abs. Non-ref. Window : 0.100 min
 Uncalibrated Peaks : not reported
 Partial Calibration : Yes, identified peaks are recalibrated
 Correct All Ret. Times: No, only for identified peaks

Calib. Data Modified : Thursday, March 28, 2019 3:48:31 PM
 Signals calculated separately : NO

General Calibration Setting

=====
Calibration Table
=====

RT	Sig	Lvl	Amount	Area	Rsp.Factor	Ref	ISTD	#	Compound
2.586	1	1	1.00000	3.69669	2.70512e-1	No	No	1	methanol
2.809	1	1	1.00000	4.26100	2.34687e-1	No	No	2	Acetaldehyde
2.977	2	1	1.00000	4.26100	2.34687e-1	No	No	2	Acetaldehyde
3.075	1	1	5.00000e-2	4.60040	1.08686e-2	No	No	1	ethanol
3.388	2	1	1.00000	4.26062	2.34707e-1	No	No	2	methanol
3.628	1	1	1.00000	9.73055	1.02769e-1	No	No	1	isopropyl alcohol
4.285	2	1	5.00000e-2	4.75733	1.05101e-2	No	No	2	ethanol
4.308	1	1	1.00000	6.49940	1.53860e-1	No	No	1	acetone
4.620	1	1	1.00000	44.73034	2.23562e-2	No	Yes	1	n-propanol
4.661	2	1	1.00000	6.89301	1.45075e-1	No	No	2	acetone
4.969	2	1	1.00000	10.70642	9.34019e-2	No	No	2	isopropyl alcohol
7.550	2	1	1.00000	46.76783	2.13822e-2	No	Yes	2	n-propanol
				46.62830	2.14462e-2				
				46.86003	2.13401e-2				
				46.73500	2.13972e-2				
				46.98234	2.12846e-2				

Peak Sum Table


No Entries in table

1 Warnings or Errors :





Warning : Curve requires more calibration points., (methanol)

Calibration Curves





methanol at exp. RT: 2.586
 FID1 A, Front Signal
 Correlation: 1.00000
 Residual Std. Dev.: 0.00000
 Formula: $Y = mx + b$
 m: 8.26440e-2
 b: 0.00000
 x: Amount Ratio
 y: Area Ratio






Handwritten mark resembling a stylized 'K' or '6'.



26



96




=====

26

ISP Forensic Services Blood Alcohol Report


Sample Name : 0.050 FN04271601
 Laboratory : Meridian
 Injection Date : Mar 28, 2019
 Method : ALCOHOL.M
 Acq. Instrument : CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	4.60040	0.0501	g/100cc
2.	Ethanol	Column 2:	4.75733	0.0512	g/100cc
3.	n-Propanol	Column 1:	44.73034	1.0000	g/100cc
4.	n-Propanol	Column 2:	46.76783	1.0000	g/100cc

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.100 FN08101601
 Laboratory : Meridian
 Injection Date : Mar 28, 2019
 Method : ALCOHOL.M
 Acq. Instrument : CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	9.14620	0.0994	g/100cc
2.	Ethanol	Column 2:	9.50124	0.0989	g/100cc
3.	n-Propanol	Column 1:	44.66777	1.0000	g/100cc
4.	n-Propanol	Column 2:	46.62830	1.0000	g/100cc

86

ISP Forensic Services Blood Alcohol Report


Sample Name : 0.200 FN03301601

Laboratory : Meridian

Injection Date : Mar 28, 2019

Method : ALCOHOL.M

Acq. Instrument : CN11180014-CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	18.68408	0.2001	g/100cc
2.	Ethanol	Column 2:	19.59364	0.1992	g/100cc
3.	n-Propanol	Column 1:	45.21098	1.0000	g/100cc
4.	n-Propanol	Column 2:	46.86003	1.0000	g/100cc

JK

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.300 FN02121601
 Laboratory : Meridian
 Injection Date : Mar 28, 2019
 Method : ALCOHOL.M
 Acq. Instrument: CN11180014-CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	28.06101	0.3007	g/100cc
2.	Ethanol	Column 2:	29.66684	0.3006	g/100cc
3.	n-Propanol	Column 1:	45.15730	1.0000	g/100cc
4.	n-Propanol	Column 2:	46.73500	1.0000	g/100cc

26

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.500 FN08031602
 Laboratory : Meridian
 Injection Date : Mar 28, 2019
 Method : ALCOHOL.M
 Acq. Instrument : CN11180014-CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	46.89567	0.4996	g/100cc
2.	Ethanol	Column 2:	49.84876	0.5000	g/100cc
3.	n-Propanol	Column 1:	45.39736	1.0000	g/100cc
4.	n-Propanol	Column 2:	46.98234	1.0000	g/100cc

26

ISP Forensic Services Blood Alcohol Report

Sample Name : INTERNAL STANDARD BLANK
 Laboratory : Meridian
 Injection Date : Mar 28, 2019
 Method : ALCOHOL.M
 Instrument: CN1180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	0.00000	0.0000	g/100cc
2.	Ethanol	Column 2:	0.00000	0.0000	g/100cc
3.	n-Propanol	Column 1:	43.61692	1.0000	g/100cc
4.	n-Propanol	Column 2:	45.03482	1.0000	g/100cc

S a m p l e S u m m a r y

Sequence table: C:\Chem32\1\data\03-28-19\CAL\03-28-19\CAL 2019-03-28 14-43-22\03-28-19\CAL.S
 Data directory path: C:\Chem32\1\data\03-28-19\CAL\03-28-19\CAL 2019-03-28 14-43-22\
 Logbook: C:\Chem32\1\data\03-28-19\CAL\03-28-19\CAL 2019-03-28 14-43-22\03-28-19\CAL.LOG
 Sequence start: 3/28/2019 2:57:59 PM
 Sequence Operator: SYSTEM
 Operator: SYSTEM
 Method file name: C:\Chem32\1\data\03-28-19\CAL\03-28-19\CAL 2019-03-28 14-43-22\ALCOHOL.M


Run Location Inj	#	Sample Name	Sample Amt	Multip.*	File name	Cal #	Cmp #
1 1	1	FN04271601	1.0000	001F0101.D	1.0000	4	*
2 2	1	FN08101601	1.0000	002F0201.D	1.0000	4	*
3 3	1	FN03301601	1.0000	003F0301.D	1.0000	4	*
4 4	1	FN02121601	1.0000	004F0401.D	1.0000	4	*
5 5	1	FN08031602	1.0000	005F0501.D	1.0000	4	*
6 6	1	INTERNAL STANDARD	1.0000	006F0601.D	1.0000	2	

JK

26

ISP Forensic Services Blood Alcohol Report


Sample Name : INTERNAL STD BLK 1
 Laboratory : Meridian
 Injection Date : Mar 28, 2019
 Method : ALCOHOL.M
 Acq. Instrument : CN1180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	0.00000	0.0000	g/100cc
2.	Ethanol	Column 2:	0.00000	0.0000	g/100cc
3.	n-Propanol	Column 1:	43.17646	1.0000	g/100cc
4.	n-Propanol	Column 2:	44.75474	1.0000	g/100cc

ISP Forensic Services Blood Alcohol Report

Sample Name : MIX VOL FN06041502
 Laboratory : Meridian
 Injection Date : Mar 28, 2019
 Method : ALCOHOL.M
 Acq. Instrument : CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	8.90980	0.1126	g/100cc
2.	Ethanol	Column 2:	9.24349	0.1133	g/100cc
3.	n-Propanol	Column 1:	38.36989	1.0000	g/100cc
4.	n-Propanol	Column 2:	39.42739	1.0000	g/100cc

DL

VOLATILES DETERMINATION CASEFILE WORKSHEET

Laboratory No.: QCI-1

Analysis Date(s): 28 Mar 2019

Sample Results	Column 1 FID A	Column 2 FID B	Column Precision	Mean Value	Over-all Mean
0.0783	0.0793	0.0010	0.0788	0.0788	0.0788
0.0787	0.0789	0.0002	0.0788	0.0788	

Analysis Method

Refer to Blood Alcohol Method #1

Instrument Information
Instrument method is stored centrally.

Refer to Instrument Method: Alcohol.m
Hamilton Auto-Dilutor Serial Number: ML600HC11378


Reporting of Results			
Overall Mean (g/100cc)	Low	0.074	0.078
High	0.082	0.004	
Uncertainty of Measurement (UM%): 5.00%			
Reported Result		0.078	

Calibration and control data are stored centrally.

16

ISP Forensic Services Blood Alcohol Report

Sample Name : Q1-1-A
 Laboratory : Meridian
 Injection Date : Mar 28, 2019
 Method : ALCOHOL.M
 Acq. Instrument : CN1180014-CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	7.05676	0.0783	g/100cc
2.	Ethanol	Column 2:	7.28550	0.0793	g/100cc
3.	n-Propanol	Column 1:	43.79214	1.0000	g/100cc
4.	n-Propanol	Column 2:	45.02147	1.0000	g/100cc

2

DL

ISP Forensic Services Blood Alcohol Report

Sample Name : Q1-1-B
 Laboratory : Meridian
 Injection Date : Mar 28, 2019
 Method : ALCOHOL.M
 Acq. Instrument : CN1180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	7.10325	0.0787	g/100cc
2.	Ethanol	Column 2:	7.29320	0.0789	g/100cc
3.	n-Propanol	Column 1:	43.85428	1.0000	g/100cc
4.	n-Propanol	Column 2:	45.28424	1.0000	g/100cc

VOLATILES DETERMINATION CASEFILE WORKSHEET

Laboratory No.: 08 FN04171701 Analysis Date(s): 28 Mar 2019

	Column 1	Column 2	Column Precision	Mean Value	Over-all Mean
Sample Results	0.0812	0.0813	0.0001	0.0812	0.0810
(g/100cc)	0.0803	0.0813	0.0010	0.0808	

Analysis Method

Refer to Blood Alcohol Method #1

Instrument Information

Refer to Instrument Method: Alcohol.m
Hamilton Auto-Dilutor Serial Number: ML600HC11378

Reporting of Results

Uncertainty of Measurement (UM%): 5.00%

Overall Mean (g/100cc)	Low	High	5% of Mean
0.081	0.076	0.086	0.005

Reported Result


0.081

Calibration and control data are stored centrally.

JK

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.08 FN04171701-A
 Laboratory : Meridian
 Injection Date : Mar 28, 2019
 Method : ALCOHOL.M
 Acq. Instrument : CN11180014-CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	7.25831	0.0812	g/100cc
2.	Ethanol	Column 2:	7.43800	0.0813	g/100cc
3.	n-Propanol	Column 1:	43.45200	1.0000	g/100cc
4.	n-Propanol	Column 2:	44.80337	1.0000	g/100cc

JK

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.08 FN04171701-B
 Laboratory : Meridian
 Injection Date : Mar 28, 2019
 Method : ALCOHOL.M
 Acq. Instrument: CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	7.28260	0.0803	g/100cc
2.	Ethanol	Column 2:	7.51392	0.0813	g/100cc
3.	n-Propanol	Column 1:	44.09018	1.0000	g/100cc
4.	n-Propanol	Column 2:	45.24621	1.0000	g/100cc

26

VOLATILES DETERMINATION CASEFILE WORKSHEET

Laboratory No.: QC2-1 Analysis Date(s): 28 Mar 2019

	Column 1	Column 2	Column Precision	Mean Value	Over-all Mean
Sample Results	0.2059	0.2053	0.0006	0.2056	0.2068
(g/100cc)	0.2075	0.2087	0.0012	0.2081	

Analysis Method

Refer to Blood Alcohol Method #1

Instrument Information

Refer to Instrument Method: Alcohol.m
Hamilton Auto-Dilutor Serial Number: ML600HC11378

Reporting of Results

Overall Mean (g/100cc)	Low	High	5% of Mean
0.206	0.195	0.217	0.011

Uncertainty of Measurement (UM%): 5.00%

Reported Result	0.206

Calibration and control data are stored centrally.

26


Revision: 1

Issue Date: 01/04/2019

Issuing Authority: Quality Manager

ISP Forensic Services Blood Alcohol Report

Sample Name : QC2-1-A
 Laboratory : Meridian
 Injection Date : Mar 28, 2019
 Method : ALCOHOL.M
 Acq. Instrument : CN11180014-CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	18.48496	0.2059	g/100cc
2.	Ethanol	Column 2:	19.22343	0.2053	g/100cc
3.	n-Propanol	Column 1:	43.47303	1.0000	g/100cc
4.	n-Propanol	Column 2:	44.59336	1.0000	g/100cc

26

ISP Forensic Services Blood Alcohol Report

Sample Name : QC2-1-B
 Laboratory : Meridian
 Injection Date : Mar 28, 2019
 Method : ALCOHOL.M
 Instrument : CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	18.89463	0.2075	g/100cc
2.	Ethanol	Column 2:	19.73756	0.2087	g/100cc
3.	n-Propanol	Column 1:	44.10130	1.0000	g/100cc
4.	n-Propanol	Column 2:	45.03465	1.0000	g/100cc

JK

VOLATILES DETERMINATION CASEFILE WORKSHEET

Laboratory No.: QC1-2 Analysis Date(s): 28 Mar 2019

	Column 1 FID A	Column 2 FID B	Column Precision	Mean Value	Over-all Mean
Sample Results	0.0828	0.0843	0.0015	0.0835	0.0829
(g/100cc)	0.0821	0.0827	0.0006	0.0824	

Analysis Method
 Refer to Blood Alcohol Method #1

Instrument Information
 Refer to Instrument Method: Alcohol.m
 Hamilton Auto-Dilutor Serial Number: ML600HC11378
Instrument method is stored centrally.


Reporting of Results			
Overall Mean (g/100cc)	Low	High	5% of Mean
0.082	0.077	0.087	0.005
Reported Result		0.082	

Calibration and control data are stored centrally.

7P

ISP Forensic Services Blood Alcohol Report


Sample Name : Q1-2-A
 Laboratory : Meridian
 Injection Date : Mar 28, 2019
 Method : ALCOHOL.M
 Acq. Instrument : CN1180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	7.66010	0.0828	g/100cc
2.	Ethanol	Column 2:	7.91551	0.0843	g/100cc
3.	n-Propanol	Column 1:	44.96146	1.0000	g/100cc
4.	n-Propanol	Column 2:	45.88955	1.0000	g/100cc

ISP Forensic Services Blood Alcohol Report

Sample Name : Q1-2-B
 Laboratory : Meridian
 Injection Date : Mar 28, 2019
 Method : ALCOHOL.M
 Acq. Instrument: CN1180014-CN11041167




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	7.61887	0.0821	g/100cc
2.	Ethanol	Column 2:	7.81172	0.0827	g/100cc
3.	n-Propanol	Column 1:	45.07674	1.0000	g/100cc
4.	n-Propanol	Column 2:	46.18452	1.0000	g/100cc

2

ISP Forensic Services Blood Alcohol Report

Sample Name : INTERNAL STD BLK
 Laboratory : Meridian
 Injection Date : Mar 28, 2019
 Method : ALCOHOL.M
 Acq. Instrument : CN11180014-CN11041167



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	0.00000	0.0000	g/100cc
2.	Ethanol	Column 2:	0.00000	0.0000	g/100cc
3.	n-Propanol	Column 1:	43.01557	1.0000	g/100cc
4.	n-Propanol	Column 2:	43.86008	1.0000	g/100cc

16

S a m p l e S u m m a r y

Sequence table: C:\Chem32\1\Data\03-28-19_SAMPLES\03-28-19_SAMPLES\2019-03-28 16-13-56\03-28-19_SAMPLES.S
 Data directory path: C:\Chem32\1\Data\03-28-19_SAMPLES\03-28-19_SAMPLES\2019-03-28 16-13-56\03-28-19_SAMPLES.LOG
 Logbook: C:\Chem32\1\Data\03-28-19_SAMPLES\03-28-19_SAMPLES\2019-03-28 16-13-56\03-28-19_SAMPLES.S
 Sequence start: 3/28/2019 4:28:40 PM
 Sequence Operator: SYSTEM
 Operator: SYSTEM
 Method file name: C:\Chem32\1\Data\03-28-19_SAMPLES\03-28-19_SAMPLES\2019-03-28 16-13-56\ALCOHOL.M

Run Location	Inj #	Sample Name	Sample Amt	Multip.*	File name	Cal #	Cmp #
--------------	-------	-------------	------------	----------	-----------	-------	-------

1	1	INTERNAL STD BLK	1.0000		001F0101.D	2	2
2	2	MIX VOL FN060415	1.0000		002F0201.D	10	10
3	3	QC1-1-A	1.0000		003F0301.D	4	4
4	4	QC1-1-B	1.0000		004F0401.D	4	4
5	5	0.08 FN04171701-	1.0000		005F0501.D	4	4
6	6	0.08 FN04171701-	1.0000		006F0601.D	4	4
7	7	M2019-1361-1-A	1.0000		007F0701.D	4	4
8	8	M2019-1361-1-B	1.0000		008F0801.D	4	4
9	9	M2019-1372-1-A	1.0000		009F0901.D	4	4
10	10	M2019-1372-1-B	1.0000		010F1001.D	4	4
11	11	M2019-1373-1-A	1.0000		011F1101.D	2	2
12	12	M2019-1373-1-B	1.0000		012F1201.D	2	2
13	13	M2019-1383-1-A	1.0000		013F1301.D	4	4
14	14	M2019-1383-1-B	1.0000		014F1401.D	4	4
15	15	M2019-1399-1-A	1.0000		015F1501.D	4	4
16	16	M2019-1399-1-B	1.0000		016F1601.D	4	4
17	17	M2019-1402-1-A	1.0000		017F1701.D	4	4
18	18	M2019-1402-1-B	1.0000		018F1801.D	4	4
19	19	M2019-1403-1-A	1.0000		019F1901.D	4	4
20	20	M2019-1403-1-B	1.0000		020F2001.D	4	4
21	21	M2019-1404-1-A	1.0000		021F2101.D	2	2
22	22	M2019-1404-1-B	1.0000		022F2201.D	2	2
23	23	M2019-1405-1-A	1.0000		023F2301.D	4	4
24	24	M2019-1405-1-B	1.0000		024F2401.D	4	4
25	25	QC2-1-A	1.0000		025F2501.D	4	4
26	26	QC2-1-B	1.0000		026F2601.D	4	4
27	27	M2019-1406-1-A	1.0000		027F2701.D	4	4
28	28	M2019-1406-1-B	1.0000		028F2801.D	4	4
29	29	M2019-1414-1-A	1.0000		029F2901.D	4	4
30	30	M2019-1414-1-B	1.0000		030F3001.D	4	4
31	31	M2019-1415-1-A	1.0000		031F3101.D	2	2
32	32	M2019-1415-1-B	1.0000		032F3201.D	2	2
33	33	QC1-2-A	1.0000		033F3301.D	4	4
34	34	QC1-2-B	1.0000		034F3401.D	4	4
35	35	INTERNAL STD BLK	1.0000		035F3501.D	2	2

Handwritten notes in blue ink:
 -2-A (next to row 31)
 -2-B (next to row 32)
 10 (next to row 26)
 10X+ (next to row 24)
 10X+ (next to row 23)

Handwritten blue mark resembling the number '16'.

Sequence File C:\Chem32\...9_SAMPLES\03-28-19_SAMPLES 2019-03-28 16-13-56
 Method file name: C:\Chem32\1\Data\03-28-19_SAMPLES\03-28-19_SAMPLES 2019-03-28 16-13-56
 \SHUTDOWN.M

Run Location Inj	Sample Name	Sample Amt	Multip.*	File name	Cal #	Cmp #
36 36	1 EMPTY	-	1.0000	036F3601.D	0	0

JK