


REVIEWED

By Rachel Cutler at 4:28 pm, Aug 14, 2019

8/12/2019

Worklist: 3608

<u>LAB CASE</u>	<u>ITEM</u>	<u>TASK ID</u>	<u>DESCRIPTION</u>
C2019-1463	1	158176	Alcohol Analysis
C2019-1485	1	158836	Alcohol Analysis
C2019-1489	1	158865	Alcohol Analysis
C2019-1521	1	159337	Alcohol Analysis
C2019-1531	1	159377	Alcohol Analysis
C2019-1532	1	159382	Alcohol Analysis
C2019-1556	1	159703	Alcohol Analysis
C2019-1557	1	159850	Alcohol Analysis



99

Worklist: 3609

<u>LAB CASE</u>	<u>ITEM</u>	<u>TASK ID</u>	<u>DESCRIPTION</u>
C2019-1512	1	159209	Alcohol Analysis



99

Quantitative Analysis for Ethanol & Qualitative Analysis for Other Volatiles

Analytical Method(s): 1.0

Device: Hamilton MICROLAB 600A Liquid Processor/Dilutor Serial Number: ML600HC11379

Volatiles Quality Assurance Controls

Run Date(s): 8/12/19

Control level	Expiration	Lot #	Target Value	Acceptable Range	Overall Results
Level 1	Jan-22	1801036	0.0812	0.0731-0.0893	0.0793 g/100cc
					0.0799 g/100cc
					g/100cc
Level 2	Mar-22	1803028	0.2035	0.1832-0.2238	0.1999 g/100cc
					g/100cc
					g/100cc
Multi-Component mixture:		Sep-20	Lot #	FN06041502	OK
Curve Fit:		Column 1	1.00000	Column 2	0.99999

Ethanol Calibration Reference Material						
Calibrator level	Target Value	Acceptable Range	Column 1	Column 2	Precision	Mean
50	0.050	0.045 - 0.055	0.0501	0.0505	0.0004	0.0503
100	0.100	0.090 - 0.110	0.0993	0.0991	0.0002	0.0992
200	0.200	0.180 - 0.220	0.2000	0.1995	0.0005	0.1997
300	0.300	0.270 - 0.330	0.3002	0.3006	0.0004	0.3004
500	0.500	0.450 - 0.550			0.0000	#DIV/0!

Aqueous Controls			
Control level	Target Value	Acceptable Range	Overall Results
80	0.080	0.076 - 0.084	0.081 g/100cc

Revision: 1

Issue Date: 01/03/2019

Issuing Authority: Quality Manager

Sample Summary

Sequence table: C:\Chem32\1\TEMP\AESEQ\QS_12.08.2019_02.34.44\8-12-2019.S
 Data directory path: C:\Chem32\1\Data\8-12-2019-JJ
 Logbook: C:\Chem32\1\Data\8-12-2019-JJ\8-12-2019.LOG
 Sequence start: 8/12/2019 2:48:28 PM
 Sequence Operator: SYSTEM
 Operator: SYSTEM

Method file name: C:\CHEM32\1\METHODS\ALCOHOL.M

Run #	Location #	Inj #	Sample Name	Sample Amt [g/100cc]	Multip.* Dilution	File name	Cal #	# Cmp
1	1	1	water	-	1.0000	001F0101.D		0
2	2	1	VOL MIX FN-06041	-	1.0000	002F0201.D		10
3	3	1	ISTD BLANK	-	1.0000	003F0301.D		2
4	4	1	QC-1-A	-	1.0000	004F0401.D		4
5	5	1	QC-1-B	-	1.0000	005F0501.D		4
6	6	1	0.08 FN04171701-	-	1.0000	006F0601.D		4
7	7	1	0.08 FN04171701-	-	1.0000	007F0701.D		4
8	8	1	19803-180-A	-	1.0000	008F0801.D		4
9	9	1	19803-180-B	-	1.0000	009F0901.D		4
10	10	1	19803-022-A	-	1.0000	010F1001.D		4
11	11	1	19803-022-B	-	1.0000	011F1101.D		4
12	12	1	C2019-1463-1-A	-	1.0000	012F1201.D		4
13	13	1	C2019-1463-1-B	-	1.0000	013F1301.D		4
14	14	1	C2019-1485-1-A	-	1.0000	014F1401.D		2
15	15	1	C2019-1485-1-B	-	1.0000	015F1501.D		2
16	16	1	C2019-1489-1-A	-	1.0000	016F1601.D		4
17	17	1	C2019-1489-1-B	-	1.0000	017F1701.D		4
18	18	1	C2019-1512-1-A	-	1.0000	018F1801.D		2
19	19	1	C2019-1512-1-B	-	1.0000	019F1901.D		2
20	20	1	C2019-1521-1-A	-	1.0000	020F2001.D		4
21	21	1	C2019-1521-1-B	-	1.0000	021F2101.D		4
22	22	1	C2019-1531-1-A	-	1.0000	022F2201.D		2
23	23	1	C2019-1531-1-B	-	1.0000	023F2301.D		2
24	24	1	C2019-1532-1-A	-	1.0000	024F2401.D		4
25	25	1	C2019-1532-1-B	-	1.0000	025F2501.D		4
26	26	1	QC-2-A	-	1.0000	026F2601.D		4
27	27	1	QC-2-B	-	1.0000	027F2701.D		4
28	28	1	C2019-1556-1-A	-	1.0000	028F2801.D		2
29	29	1	C2019-1556-1-B	-	1.0000	029F2901.D		2
30	30	1	C2019-1557-1-A	-	1.0000	030F3001.D		4
31	31	1	C2019-1557-1-B	-	1.0000	031F3101.D		4
32	32	1	QC-1-A	-	1.0000	032F3201.D		4
33	33	1	QC-1-B	-	1.0000	033F3301.D		4
34	34	1	ISTD BLANK	-	1.0000	034F3401.D		2
35	35	1	water	-	1.0000	035F3501.D		0

0.500 calibrator was found to be expired after running the sequence and processing the data. Calibration curve was reprocessed after the removal of the expired calibrator and all of the data was reprocessed. None of the reported results changed in value. The data was reprocessed on 8/13/2019.

Jeremy Johnson

99

=====
 Calibration Table
 =====

 General Calibration Setting

Calib. Data Modified : Tuesday, August 13, 2019 1:31:41 PM
 Signals calculated separately : No

Rel. Reference Window : 0.000 %
 Abs. Reference Window : 0.100 min
 Rel. Non-ref. Window : 0.000 %
 Abs. Non-ref. Window : 0.100 min
 Uncalibrated Peaks : not reported
 Partial Calibration : No recalibration if peaks missing

Curve Type : Linear
 Origin : Forced
 Weight : Equal

Recalibration Settings:
 Average Response : Average all calibrations
 Average Retention Time: Floating Average New 75%

Calibration Report Options :
 Printout of recalibrations within a sequence:
 Calibration Table after Recalibration
 Normal Report after Recalibration
 If the sequence is done with bracketing:
 Results of first cycle (ending previous bracket)

Default Sample ISTD Information (if not set in sample table):

ISTD #	ISTD Amount [g/100cc]	Name
1	1.00000	n-Propanol
2	1.00000	n-Propanol

 Signal Details

Signal 1: FID1 A, Front Signal
 Signal 2: FID2 B, Back Signal


 Overview Table

RT	Sig	Lvl	Amount [g/100cc]	Area	Rsp.Factor	Ref	ISTD #	Compound
2.000	2	1	1.00000	5.00000	2.00000e-1	No	No 2	Difluoroethane
2.000	1	1	1.00000	5.00000	2.00000e-1	No	No 1	Difluoroethane
2.494	1	1	1.00000	3.69669	2.70512e-1	No	No 1	Methanol
2.772	1	1	1.00000	3.19311	3.13174e-1	No	No 1	Acetaldehyde
2.797	2	1	1.00000	3.10575	3.21983e-1	No	No 2	Acetaldehyde
3.109	1	1	5.00000e-2	8.90005	5.61794e-3	No	No 1	Ethanol
		2	1.00000e-1	17.97227	5.56413e-3			
		3	2.00000e-1	35.37431	5.65382e-3			
		4	3.00000e-1	54.59329	5.49518e-3			
3.211	2	1	1.00000	4.26062	2.34707e-1	No	No 2	Methanol
3.715	1	1	1.00000	9.73055	1.02769e-1	No	No 1	Isopropyl alcohol
4.181	2	1	5.00000e-2	9.07075	5.51223e-3	No	No 2	Ethanol
		2	1.00000e-1	18.25399	5.47825e-3			
		3	2.00000e-1	35.69363	5.60324e-3			
		4	3.00000e-1	54.75084	5.47937e-3			
4.530	1	1	1.00000	6.49940	1.53860e-1	No	No 1	Acetone
4.549	2	1	1.00000	6.89301	1.45075e-1	No	No 2	Acetone
4.870	2	1	1.00000	10.70642	9.34019e-2	No	No 2	Isopropyl alcohol
4.944	1	1	1.00000	87.61784	1.14132e-2	No	Yes 1	n-Propanol
		2	1.00000	89.28069	1.12006e-2			
		3	1.00000	87.26619	1.14592e-2			
		4	1.00000	89.72559	1.11451e-2			
		5	1.00000	90.23061	1.10827e-2			
7.624	2	1	1.00000	87.80907	1.13883e-2	No	Yes 2	n-Propanol
		2	1.00000	90.06332	1.11033e-2			
		3	1.00000	87.50938	1.14273e-2			
		4	1.00000	89.09099	1.12245e-2			
		5	1.00000	89.18039	1.12132e-2			

Peak Sum Table


No Entries in table

Calibration Curves




Difluoroethane at exp. RT: 2.000
 FID2 B, Back Signal
 Correlation: 1.00000
 Residual Std. Dev.: 0.00000
 Formula: $y = mx$
 m: 5.69417e-2
 x: Amount Ratio
 y: Area Ratio


99




Difluoroethane at exp. RT: 2.000
FID1 A, Front Signal
Correlation: 1.00000
Residual Std. Dev.: 0.00000
Formula: $y = mx$
m: 5.70660e-2
x: Amount Ratio
y: Area Ratio



Methanol at exp. RT: 2.494
FID1 A, Front Signal
Correlation: 1.00000
Residual Std. Dev.: 0.00000
Formula: $y = mx$
m: 4.21911e-2
x: Amount Ratio
y: Area Ratio




Acetaldehyde at exp. RT: 2.772
FID1 A, Front Signal
Correlation: 1.00000
Residual Std. Dev.: 0.00000
Formula: $y = mx$
m: 3.64436e-2
x: Amount Ratio
y: Area Ratio




Acetaldehyde at exp. RT: 2.797
FID2 B, Back Signal
Correlation: 1.00000
Residual Std. Dev.: 0.00000
Formula: $y = mx$
m: 3.53694e-2
x: Amount Ratio
y: Area Ratio


99




Ethanol at exp. RT: 3.109
 FID1 A, Front Signal
 Correlation: 1.00000 ✓
 Residual Std. Dev.: 0.00084
 Formula: $y = mx$
 m: 2.02677
 x: Amount Ratio
 y: Area Ratio



Methanol at exp. RT: 3.211
 FID2 B, Back Signal
 Correlation: 1.00000
 Residual Std. Dev.: 0.00000
 Formula: $y = mx$
 m: 4.85215e-2
 x: Amount Ratio
 y: Area Ratio




Isopropyl alcohol at exp. RT: 3.715
 FID1 A, Front Signal
 Correlation: 1.00000
 Residual Std. Dev.: 0.00000
 Formula: $y = mx$
 m: 1.11057e-1
 x: Amount Ratio
 y: Area Ratio




Ethanol at exp. RT: 4.181
 FID2 B, Back Signal
 Correlation: 0.99999 ✓
 Residual Std. Dev.: 0.00150
 Formula: $y = mx$
 m: 2.04473
 x: Amount Ratio
 y: Area Ratio


99




Acetone at exp. RT: 4.530
 FID1 A, Front Signal
 Correlation: 1.00000
 Residual Std. Dev.: 0.00000
 Formula: $y = mx$
 m: $7.41790e-2$
 x: Amount Ratio
 y: Area Ratio



Acetone at exp. RT: 4.549
 FID2 B, Back Signal
 Correlation: 1.00000
 Residual Std. Dev.: 0.00000
 Formula: $y = mx$
 m: $7.85000e-2$
 x: Amount Ratio
 y: Area Ratio




Isopropyl alcohol at exp. RT: 4.870
 FID2 B, Back Signal
 Correlation: 1.00000
 Residual Std. Dev.: 0.00000
 Formula: $y = mx$
 m: $1.21928e-1$
 x: Amount Ratio
 y: Area Ratio



n-Propanol at exp. RT: 4.944
 FID1 A, Front Signal
 Correlation: 1.00000
 Residual Std. Dev.: 0.00000
 Formula: $y = mx$
 m: 1.00000
 x: Amount Ratio
 y: Area Ratio

99



n-Propanol at exp. RT: 7.624
FID2 B, Back Signal
Correlation: 1.00000
Residual Std. Dev.: 0.00000
Formula: $y = mx$
m: 1.00000
x: Amount Ratio
y: Area Ratio

=====

99

S a m p l e S u m m a r y

Sequence table: C:\Chem32\1\TEMP\AESEQ\QS_12.08.2019_12.46.31\8-12-19cal.S
 Data directory path: C:\Chem32\1\Data\8-12-19calJJ
 Logbook: C:\Chem32\1\Data\8-12-19calJJ\8-12-19cal.LOG
 Sequence start: 8/12/2019 1:00:14 PM
 Sequence Operator: SYSTEM
 Operator: SYSTEM


Method file name: C:\CHEM32\1\METHODS\ALCOHOL.M

Run #	Location #	Inj #	Sample Name	Sample Amt [g/100cc]	Multip.* Dilution	File name	Cal #	# Cmp
1	1	1	0.05	-	1.0000	001F0101.D	*	4
2	2	1	0.100	-	1.0000	002F0201.D	*	4
3	3	1	0.200	-	1.0000	003F0301.D	*	4
4	4	1	0.300	-	1.0000	004F0401.D	*	4
5	5	1	0.500	-	1.0000	005F0501.D	*	4
6	6	1	blank	-	1.0000	006F0601.D		2

99

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.05
 Laboratory : Coeur d' Alene
 Injection Date : Aug 12, 2019
 Method : ALCOHOL.M
 Acq. Instrument: CN10742044-IT00725005




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	8.90005	0.0501	g/100cc
2.	Ethanol	Column 2:	9.07075	0.0505	g/100cc
3.	n-Propanol	Column 1:	87.61784	1.0000	g/100cc
4.	n-Propanol	Column 2:	87.80907	1.0000	g/100cc

99

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.100
 Laboratory : Coeur d' Alene
 Injection Date : Aug 12, 2019
 Method : ALCOHOL.M
 Acq. Instrument: CN10742044-IT00725005




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	17.97227	0.0993	g/100cc
2.	Ethanol	Column 2:	18.25399	0.0991	g/100cc
3.	n-Propanol	Column 1:	89.28069	1.0000	g/100cc
4.	n-Propanol	Column 2:	90.06332	1.0000	g/100cc

99

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.200
 Laboratory : Coeur d' Alene
 Injection Date : Aug 12, 2019
 Method : ALCOHOL.M
 Acq. Instrument: CN10742044-IT00725005




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	35.37431	0.2000	g/100cc
2.	Ethanol	Column 2:	35.69363	0.1995	g/100cc
3.	n-Propanol	Column 1:	87.26619	1.0000	g/100cc
4.	n-Propanol	Column 2:	87.50938	1.0000	g/100cc

99

ISP Forensic Services Blood Alcohol Report


Sample Name : 0.300
 Laboratory : Coeur d' Alene
 Injection Date : Aug 12, 2019
 Method : ALCOHOL.M
 Acq. Instrument: CN10742044-IT00725005



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	54.59329	0.3002	g/100cc
2.	Ethanol	Column 2:	54.75084	0.3006	g/100cc
3.	n-Propanol	Column 1:	89.72559	1.0000	g/100cc
4.	n-Propanol	Column 2:	89.09099	1.0000	g/100cc

ISP Forensic Services Blood Alcohol Report

Sample Name : blank
 Laboratory : Coeur d' Alene
 Injection Date : Aug 12, 2019
 Method : ALCOHOL.M
 Acq. Instrument: CN10742044-IT00725005




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	0.00000	0.0000	g/100cc
2.	Ethanol	Column 2:	0.00000	0.0000	g/100cc
3.	n-Propanol	Column 1:	90.61393	1.0000	g/100cc
4.	n-Propanol	Column 2:	89.14573	1.0000	g/100cc

99

ISP Forensic Services Blood Alcohol Report

Sample Name : water
 Laboratory : Coeur d' Alene
 Injection Date : Aug 12, 2019
 Method : ALCOHOL.M
 Acq. Instrument: CN10742044-IT00725005




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	0.00000	0.0000	g/100cc
2.	Ethanol	Column 2:	0.00000	0.0000	g/100cc
3.	n-Propanol	Column 1:	0.00000	0.0000	g/100cc
4.	n-Propanol	Column 2:	0.00000	0.0000	g/100cc

99

ISP Forensic Services Blood Alcohol Report

Sample Name : VOL MIX FN-06041502
 Laboratory : Coeur d' Alene
 Injection Date : Aug 12, 2019
 Method : ALCOHOL.M
 Acq. Instrument: CN10742044-IT00725005




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	37.41499	0.2113	g/100cc
2.	Ethanol	Column 2:	37.00815	0.2110	g/100cc
3.	n-Propanol	Column 1:	87.36483	1.0000	g/100cc
4.	n-Propanol	Column 2:	85.77782	1.0000	g/100cc

99

ISP Forensic Services Blood Alcohol Report

Sample Name : ISTD BLANK
 Laboratory : Coeur d' Alene
 Injection Date : Aug 12, 2019
 Method : ALCOHOL.M
 Acq. Instrument: CN10742044-IT00725005



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	0.00000	0.0000	g/100cc
2.	Ethanol	Column 2:	0.00000	0.0000	g/100cc
3.	n-Propanol	Column 1:	90.23061	1.0000	g/100cc
4.	n-Propanol	Column 2:	89.18039	1.0000	g/100cc

99

VOLATILES DETERMINATION CASEFILE WORKSHEET

Laboratory No.: QC-1-A

Analysis Date(s): 12 Aug 2019

	Column 1 FID A	Column 2 FID B	Column Precision	Mean Value	Over-all Mean
Sample Results	0.0794	0.0791	0.0003	0.0792	0.0793
(g/100cc)	0.0794	0.0795	0.0001	0.0794	

Analysis Method

Refer to Blood Alcohol Method #1

Instrument Information

Instrument method is stored centrally.

Refer to Instrument Method: Alcohol.m
Hamilton Auto-Dilutor Serial Number: ML600HC11379

Reporting of Results

Uncertainty of Measurement (UM%): 5.00%


Overall Mean (g/100cc)	Low	High	5% of Mean
0.079	0.075	0.083	0.004

	Reported Result	
	0.079	

Calibration and control data are stored centrally.

ISP Forensic Services Blood Alcohol Report


Sample Name : QC-1-A
 Laboratory : Coeur d' Alene
 Injection Date : Aug 12, 2019
 Method : ALCOHOL.M
 Acq. Instrument: CN10742044-IT00725005



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	14.55525	0.0794	g/100cc
2.	Ethanol	Column 2:	14.49028	0.0791	g/100cc
3.	n-Propanol	Column 1:	90.46012	1.0000	g/100cc
4.	n-Propanol	Column 2:	89.61170	1.0000	g/100cc

ISP Forensic Services Blood Alcohol Report

Sample Name : QC-1-B
 Laboratory : Coeur d' Alene
 Injection Date : Aug 12, 2019
 Method : ALCOHOL.M
 Acq. Instrument: CN10742044-IT00725005



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	14.58807	0.0794	g/100cc
2.	Ethanol	Column 2:	14.60417	0.0795	g/100cc
3.	n-Propanol	Column 1:	90.62993	1.0000	g/100cc
4.	n-Propanol	Column 2:	89.82311	1.0000	g/100cc

VOLATILES DETERMINATION CASEFILE WORKSHEET

Laboratory No.: 0.08 FN04171701-B

Analysis Date(s): 12 Aug 2019

	Column 1 FID A	Column 2 FID B	Column Precision	Mean Value	Over-all Mean	
Sample Results	0.0810	0.0810	0.0000	0.0810	0.0810	
(g/100cc)	0.0810	0.0810	0.0000	0.0810		

Analysis Method

Refer to Blood Alcohol Method #1

Instrument Information

Instrument method is stored centrally.

Refer to Instrument Method: Alcohol.m
Hamilton Auto-Dilutor Serial Number: ML600HC11379

Reporting of Results

Uncertainty of Measurement (UM%): 5.00%


Overall Mean (g/100cc)	Low	High	5% of Mean
0.081	0.076	0.086	0.005

	Reported Result	
	0.081	

Calibration and control data are stored centrally.

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.08 FN04171701-A
 Laboratory : Coeur d' Alene
 Injection Date : Aug 12, 2019
 Method : ALCOHOL.M
 Acq. Instrument: CN10742044-IT00725005




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	14.39749	0.0802	g/100cc
2.	Ethanol	Column 2:	14.43176	0.0806	g/100cc
3.	n-Propanol	Column 1:	88.58889	1.0000	g/100cc
4.	n-Propanol	Column 2:	87.53783	1.0000	g/100cc

99

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.08 FN04171701-B
 Laboratory : Coeur d' Alene
 Injection Date : Aug 12, 2019
 Method : ALCOHOL.M
 Acq. Instrument: CN10742044-IT00725005



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	14.42969	0.0810	g/100cc
2.	Ethanol	Column 2:	14.42987	0.0810	g/100cc
3.	n-Propanol	Column 1:	87.91169	1.0000	g/100cc
4.	n-Propanol	Column 2:	87.12988	1.0000	g/100cc

VOLATILES DETERMINATION CASEFILE WORKSHEET

Laboratory No.: QC-2-A

Analysis Date(s): 12 Aug 2019

	Column 1 FID A	Column 2 FID B	Column Precision	Mean Value	Over-all Mean	
Sample Results	0.2000	0.2014	0.0014	0.2007	0.1999	
(g/100cc)	0.1984	0.1999	0.0015	0.1991		

Analysis Method

Refer to Blood Alcohol Method #1

Instrument Information

Instrument method is stored centrally.

Refer to Instrument Method: Alcohol.m
Hamilton Auto-Dilutor Serial Number: ML600HC11379

Reporting of Results

Uncertainty of Measurement (UM%): 5.00%


Overall Mean (g/100cc)	Low	High	5% of Mean
0.199	0.189	0.209	0.010

Reported Result	
0.199	

Calibration and control data are stored centrally.

ISP Forensic Services Blood Alcohol Report

Sample Name : QC-2-A
 Laboratory : Coeur d' Alene
 Injection Date : Aug 12, 2019
 Method : ALCOHOL.M
 Acq. Instrument: CN10742044-IT00725005




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	36.84611	0.2000	g/100cc
2.	Ethanol	Column 2:	37.06873	0.2014	g/100cc
3.	n-Propanol	Column 1:	90.88922	1.0000	g/100cc
4.	n-Propanol	Column 2:	90.03351	1.0000	g/100cc

99

ISP Forensic Services Blood Alcohol Report

Sample Name : QC-2-B
 Laboratory : Coeur d' Alene
 Injection Date : Aug 12, 2019
 Method : ALCOHOL.M
 Acq. Instrument: CN10742044-IT00725005



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	36.87164	0.1984	g/100cc
2.	Ethanol	Column 2:	37.03089	0.1999	g/100cc
3.	n-Propanol	Column 1:	91.68819	1.0000	g/100cc
4.	n-Propanol	Column 2:	90.59061	1.0000	g/100cc

99

VOLATILES DETERMINATION CASEFILE WORKSHEET

Laboratory No.: QC-1-A

Analysis Date(s): 12 Aug 2019

	Column 1 FID A	Column 2 FID B	Column Precision	Mean Value	Over-all Mean
Sample Results	0.0798	0.0804	0.0006	0.0801	0.0799
(g/100cc)	0.0795	0.0802	0.0007	0.0798	

Analysis Method

Refer to Blood Alcohol Method #1

Instrument Information

Instrument method is stored centrally.

Refer to Instrument Method: Alcohol.m
Hamilton Auto-Dilutor Serial Number: ML600HC11379

Reporting of Results

Uncertainty of Measurement (UM%): 5.00%


Overall Mean (g/100cc)	Low	High	5% of Mean
0.079	0.075	0.083	0.004

Reported Result	
0.079	

Calibration and control data are stored centrally.

ISP Forensic Services Blood Alcohol Report

Sample Name : QC-1-A
 Laboratory : Coeur d' Alene
 Injection Date : Aug 12, 2019
 Method : ALCOHOL.M
 Acq. Instrument: CN10742044-IT00725005




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	14.89281	0.0798	g/100cc
2.	Ethanol	Column 2:	14.97219	0.0804	g/100cc
3.	n-Propanol	Column 1:	92.07324	1.0000	g/100cc
4.	n-Propanol	Column 2:	91.03889	1.0000	g/100cc

99

ISP Forensic Services Blood Alcohol Report

Sample Name : QC-1-B
 Laboratory : Coeur d' Alene
 Injection Date : Aug 12, 2019
 Method : ALCOHOL.M
 Acq. Instrument: CN10742044-IT00725005




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	14.72321	0.0795	g/100cc
2.	Ethanol	Column 2:	14.81103	0.0802	g/100cc
3.	n-Propanol	Column 1:	91.33966	1.0000	g/100cc
4.	n-Propanol	Column 2:	90.36684	1.0000	g/100cc

99

ISP Forensic Services Blood Alcohol Report

Sample Name : ISTD BLANK
 Laboratory : Coeur d' Alene
 Injection Date : Aug 12, 2019
 Method : ALCOHOL.M
 Acq. Instrument: CN10742044-IT00725005




#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	0.00000	0.0000	g/100cc
2.	Ethanol	Column 2:	0.00000	0.0000	g/100cc
3.	n-Propanol	Column 1:	89.71379	1.0000	g/100cc
4.	n-Propanol	Column 2:	88.96717	1.0000	g/100cc

99

ISP Forensic Services Blood Alcohol Report

Sample Name : water
 Laboratory : Coeur d' Alene
 Injection Date : Aug 12, 2019
 Method : ALCOHOL.M
 Acq. Instrument: CN10742044-IT00725005



#	Compound	Column	Area	Amount	Units
1.	Ethanol	Column 1:	0.00000	0.0000	g/100cc
2.	Ethanol	Column 2:	0.00000	0.0000	g/100cc
3.	n-Propanol	Column 1:	0.00000	0.0000	g/100cc
4.	n-Propanol	Column 2:	0.00000	0.0000	g/100cc

99