









REVIEWED

By Galina Giso at 12:59 pm, Jun 01, 2020

5/31/2020

Worklist: 4268

| <u>LAB CASE</u> | <u>ITEM</u> | <u>ITEM TYPE</u> | <u>DESCRIPTION</u> | |
|-----------------|-------------|------------------|--------------------|---|
| C2020-0946 | 1 | BCK | Alcohol Analysis |  |
| C2020-0968 | 1 | BCK | Alcohol Analysis |  |
| C2020-0991 | 1 | BCK | Alcohol Analysis |  |
| C2020-0992 | 1 | BCK | Alcohol Analysis |  |
| C2020-0999 | 1 | BCK | Alcohol Analysis |  |
| C2020-1002 | 1 | BCK | Alcohol Analysis |  |
| C2020-1030 | 1 | BCK | Alcohol Analysis |  |
| C2020-1031 | 1 | BCK | Alcohol Analysis |  |



Quantitative Analysis for Ethanol & Qualitative Analysis for Other Volatiles

Analytical Method(s): 1.0

Device: Hamilton MICROLAB Liquid Processor/Dilutor Serial Number: ML600HC11379

Volatiles Quality Assurance Controls Run Date(s): 5-31-20

Worksheet #4268

| Control level | Expiration | Lot # | Target Value | Acceptable Range | Overall Results |
|--------------------------|------------|----------|--------------|------------------|-----------------|
| Level 1 | Jan-22 | 1801036 | 0.0812 | 0.0731-0.0893 | 0.0767 g/100cc |
| | | | | | g/100cc |
| | | | | | 0.1981 g/100cc |
| Level 2 | Mar-22 | 1803028 | 0.2035 | 0.1832-0.2238 | g/100cc |
| | | | | | g/100cc |
| | | | | | g/100cc |
| Multi-Component mixture: | | Sep-20 | Lot # | FN06041502 | OK |
| Curve Fit: | | Column 1 | 1.00000 | Column2 | 1.00000 |

| Ethanol Calibration Reference Material | | | | | | |
|--|--------------|------------------|----------|----------|-----------|---------|
| Calibrator level | Target Value | Acceptable Range | Column 1 | Column 2 | Precision | Mean |
| 50 | 0.050 | 0.045 - 0.055 | 0.0497 | 0.0491 | 0.0006 | 0.0494 |
| 100 | 0.100 | 0.090 - 0.110 | 0.1003 | 0.0995 | 0.0008 | 0.0999 |
| 200 | 0.200 | 0.180 - 0.220 | 0.2010 | 0.2006 | 0.0004 | 0.2008 |
| 300 | 0.300 | 0.270 - 0.330 | 0.3003 | 0.3000 | 0.0003 | 0.3001 |
| 400 | 0.400 | 0.360 - 0.440 | | | 0 | #DIV/0! |
| 500 | 0.500 | 0.450 - 0.550 | 0.4994 | 0.4999 | 0.0005 | 0.4996 |

| Aqueous Controls | | | |
|------------------|--------------|------------------|-----------------|
| Control level | Target Value | Acceptable Range | Overall Results |
| 80 | 0.080 | 0.076 - 0.084 | 0.080 g/100cc |

99

Sample Summary

Sequence table: C:\Chem32\1\TEMP\AESEQ\QS_31.05.2020_01.18.35\5-31-2020.S
 Data directory path: C:\Chem32\1\Data\5-31-20jj
 Logbook: C:\Chem32\1\Data\5-31-20jj\5-31-2020.LOG
 Sequence start: 5/31/2020 1:32:18 PM
 Sequence Operator: SYSTEM
 Operator: SYSTEM

Method file name: C:\CHEM32\1\METHODS\ALCOHOL.M

| Run # | Location # | Inj # | Sample Name | Sample Amt [g/100cc] | Multip.* Dilution | File name | Cal # | # Cmp |
|-------|------------|-------|------------------|----------------------|-------------------|------------|-------|-------|
| 1 | 1 | 1 | water-1 | - | 1.0000 | 001F0101.D | | 0 |
| 2 | 2 | 1 | VOL MIX FN-06041 | - | 1.0000 | 002F0201.D | | 10 |
| 3 | 3 | 1 | ISTD BLANK-1 | - | 1.0000 | 003F0301.D | | 2 |
| 4 | 4 | 1 | QC-1(1)-A | - | 1.0000 | 004F0401.D | | 4 |
| 5 | 5 | 1 | QC-1(1)-B | - | 1.0000 | 005F0501.D | | 4 |
| 6 | 6 | 1 | 0.08 FN09181807- | - | 1.0000 | 006F0601.D | | 4 |
| 7 | 7 | 1 | 0.08 FN09181807- | - | 1.0000 | 007F0701.D | | 4 |
| 8 | 8 | 1 | C2020-0946-1-A | - | 1.0000 | 008F0801.D | | 4 |
| 9 | 9 | 1 | C2020-0946-1-B | - | 1.0000 | 009F0901.D | | 4 |
| 10 | 10 | 1 | C2020-0968-1-A | - | 1.0000 | 010F1001.D | | 4 |
| 11 | 11 | 1 | C2020-0968-1-B | - | 1.0000 | 011F1101.D | | 4 |
| 12 | 12 | 1 | C2020-0991-1-A | - | 1.0000 | 012F1201.D | | 4 |
| 13 | 13 | 1 | C2020-0991-1-B | - | 1.0000 | 013F1301.D | | 4 |
| 14 | 14 | 1 | C2020-0992-1-A | - | 1.0000 | 014F1401.D | | 4 |
| 15 | 15 | 1 | C2020-0992-1-B | - | 1.0000 | 015F1501.D | | 4 |
| 16 | 16 | 1 | C2020-0999-1-A | - | 1.0000 | 016F1601.D | | 4 |
| 17 | 17 | 1 | C2020-0999-1-B | - | 1.0000 | 017F1701.D | | 4 |
| 18 | 18 | 1 | C2020-1002-1-A | - | 1.0000 | 018F1801.D | | 4 |
| 19 | 19 | 1 | C2020-1002-1-B | - | 1.0000 | 019F1901.D | | 4 |
| 20 | 20 | 1 | C2020-1030-1-A | - | 1.0000 | 020F2001.D | | 4 |
| 21 | 21 | 1 | C2020-1030-1-B | - | 1.0000 | 021F2101.D | | 4 |
| 22 | 22 | 1 | C2020-1031-1-A | - | 1.0000 | 022F2201.D | | 4 |
| 23 | 23 | 1 | C2020-1031-1-B | - | 1.0000 | 023F2301.D | | 4 |
| 24 | 24 | 1 | QC-2(1)-A | - | 1.0000 | 024F2401.D | | 4 |
| 25 | 25 | 1 | QC-2(1)-B | - | 1.0000 | 025F2501.D | | 4 |
| 26 | 26 | 1 | ISTD BLANK-2 | - | 1.0000 | 026F2601.D | | 2 |
| 27 | 27 | 1 | water-2 | - | 1.0000 | 027F2701.D | | 0 |
| 28 | 28 | 1 | 0.05 CHECK | - | 1.0000 | 028F2801.D | | 4 |
| 29 | 29 | 1 | 0.100 CHECK | - | 1.0000 | 029F2901.D | | 4 |
| 30 | 30 | 1 | 0.200 CHECK | - | 1.0000 | 030F3001.D | | 4 |
| 31 | 31 | 1 | 0.300 CHECK | - | 1.0000 | 031F3101.D | | 4 |
| 32 | 32 | 1 | 0.500 CHECK | - | 1.0000 | 032F3201.D | | 4 |

99

=====
Calibration Table
=====

General Calibration Setting

Calib. Data Modified : Sunday, May 31, 2020 1:16:42 PM
Signals calculated separately : No

Rel. Reference Window : 0.000 %
Abs. Reference Window : 0.100 min
Rel. Non-ref. Window : 0.000 %
Abs. Non-ref. Window : 0.100 min
Uncalibrated Peaks : not reported
Partial Calibration : No recalibration if peaks missing

Curve Type : Linear
Origin : Forced
Weight : Equal

Recalibration Settings:
Average Response : Average all calibrations
Average Retention Time: Floating Average New 75%

Calibration Report Options :
Printout of recalibrations within a sequence:
 Calibration Table after Recalibration
 Normal Report after Recalibration
If the sequence is done with bracketing:
 Results of first cycle (ending previous bracket)

Default Sample ISTD Information (if not set in sample table):

| ISTD # | ISTD Amount [g/100cc] | Name |
|--------|-----------------------|------------|
| 1 | 1.00000 | n-Propanol |
| 2 | 1.00000 | n-Propanol |

Signal Details

Signal 1: FID1 A, Front Signal
Signal 2: FID2 B, Back Signal

Overview Table


99

| RT | Sig | Lvl | Amount [g/100cc] | Area | Rsp.Factor | Ref | ISTD # | Compound |
|-------|-----|-----|---------------------|----------|------------|-----|--------|---------------------|
| 1.977 | 2 | 1 | 1.00000 | 1.06794 | 9.36380e-1 | No | No | 2 Difluoroethane |
| 2.000 | 1 | 1 | 1.00000 | 5.00000 | 2.00000e-1 | No | No | 1 Difluoroethane |
| 2.494 | 1 | 1 | 1.00000 | 3.69669 | 2.70512e-1 | No | No | 1 Methanol |
| 2.772 | 1 | 1 | 1.00000 | 3.19311 | 3.13174e-1 | No | No | 1 Acetaldehyde |
| 2.797 | 2 | 1 | 1.00000 | 3.10575 | 3.21983e-1 | No | No | 2 Acetaldehyde |
| 3.109 | 1 | 1 | 5.00000e-2 | 8.95644 | 5.58257e-3 | No | No | 1 Ethanol |
| | | | 1.00000e-1 | 17.97128 | 5.56443e-3 | | | |
| | | | 2.00000e-1 | 35.82807 | 5.58222e-3 | | | |
| | | | 3.00000e-1 | 54.45925 | 5.50871e-3 | | | |
| | | | 5.00000e-1 | 89.91124 | 5.56104e-3 | | | |
| 3.211 | 2 | 1 | 1.00000 | 4.26062 | 2.34707e-1 | No | No | 2 Methanol |
| 3.715 | 1 | 1 | 1.00000 | 9.73055 | 1.02769e-1 | No | No | 1 Isopropyl alcohol |
| 4.181 | 2 | 1 | 5.00000e-2 | 8.94641 | 5.58883e-3 | No | No | 2 Ethanol |
| | | | 1.00000e-1 | 17.97812 | 5.56232e-3 | | | |
| | | | 2.00000e-1 | 35.96552 | 5.56088e-3 | | | |
| | | | 3.00000e-1 | 54.67428 | 5.48704e-3 | | | |
| | | | 5.00000e-1 | 90.19679 | 5.54343e-3 | | | |
| 4.530 | 1 | 1 | 1.00000 | 6.49940 | 1.53860e-1 | No | No | 1 Acetone |
| 4.549 | 2 | 1 | 1.00000 | 6.89301 | 1.45075e-1 | No | No | 2 Acetone |
| 4.870 | 2 | 1 | 1.00000 | 10.70642 | 9.34019e-2 | No | No | 2 Isopropyl alcohol |
| 4.943 | 1 | 1 | 1.00000 | 93.16283 | 1.07339e-2 | No | Yes | 1 n-Propanol |
| | | | 1.00000 | 92.58790 | 1.08005e-2 | | | |
| | | | 1.00000 | 92.07500 | 1.08607e-2 | | | |
| | | | 1.00000 | 93.68510 | 1.06741e-2 | | | |
| | | | 1.00000 | 93.00353 | 1.07523e-2 | | | |
| 7.622 | 2 | 1 | 1.00000 | 90.25261 | 1.10800e-2 | No | Yes | 2 n-Propanol |
| | | | 1.00000 | 89.55220 | 1.11667e-2 | | | |
| | | | 1.00000 | 88.86015 | 1.12536e-2 | | | |
| | | | 1.00000 | 90.32391 | 1.10713e-2 | | | |
| | | | 1.00000 | 89.43030 | 1.11819e-2 | | | |

Peak Sum Table


No Entries in table

Calibration Curves




Difluoroethane at exp. RT: 1.977
 FID2 B, Back Signal
 Correlation: 1.00000
 Residual Std. Dev.: 0.00000
 Formula: $y = mx$
 m: 1.18328e-2
 x: Amount Ratio
 y: Area Ratio


99




Difluoroethane at exp. RT: 2.000
 FID1 A, Front Signal
 Correlation: 1.00000
 Residual Std. Dev.: 0.00000
 Formula: $y = mx$
 m: $5.36695e-2$
 x: Amount Ratio
 y: Area Ratio



Methanol at exp. RT: 2.494
 FID1 A, Front Signal
 Correlation: 1.00000
 Residual Std. Dev.: 0.00000
 Formula: $y = mx$
 m: $3.96799e-2$
 x: Amount Ratio
 y: Area Ratio




Acetaldehyde at exp. RT: 2.772
 FID1 A, Front Signal
 Correlation: 1.00000
 Residual Std. Dev.: 0.00000
 Formula: $y = mx$
 m: $3.42745e-2$
 x: Amount Ratio
 y: Area Ratio




Acetaldehyde at exp. RT: 2.797
 FID2 B, Back Signal
 Correlation: 1.00000
 Residual Std. Dev.: 0.00000
 Formula: $y = mx$
 m: $3.44117e-2$
 x: Amount Ratio
 y: Area Ratio


99




Ethanol at exp. RT: 3.109
 FID1 A, Front Signal
 Correlation: 1.00000 ✓
 Residual Std. Dev.: 0.00124
 Formula: $y = mx$
 m: 1.93581
 x: Amount Ratio
 y: Area Ratio



Methanol at exp. RT: 3.211
 FID2 B, Back Signal
 Correlation: 1.00000
 Residual Std. Dev.: 0.00000
 Formula: $y = mx$
 m: 4.72078e-2
 x: Amount Ratio
 y: Area Ratio




Isopropyl alcohol at exp. RT: 3.715
 FID1 A, Front Signal
 Correlation: 1.00000
 Residual Std. Dev.: 0.00000
 Formula: $y = mx$
 m: 1.04447e-1
 x: Amount Ratio
 y: Area Ratio




Ethanol at exp. RT: 4.181
 FID2 B, Back Signal
 Correlation: 1.00000 ✓
 Residual Std. Dev.: 0.00119
 Formula: $y = mx$
 m: 2.01748
 x: Amount Ratio
 y: Area Ratio


99




Acetone at exp. RT: 4.530
FID1 A, Front Signal
Correlation: 1.00000
Residual Std. Dev.: 0.00000
Formula: $y = mx$
m: 6.97639e-2
x: Amount Ratio
y: Area Ratio



Acetone at exp. RT: 4.549
FID2 B, Back Signal
Correlation: 1.00000
Residual Std. Dev.: 0.00000
Formula: $y = mx$
m: 7.63746e-2
x: Amount Ratio
y: Area Ratio




Isopropyl alcohol at exp. RT: 4.870
FID2 B, Back Signal
Correlation: 1.00000
Residual Std. Dev.: 0.00000
Formula: $y = mx$
m: 1.18627e-1
x: Amount Ratio
y: Area Ratio



n-Propanol at exp. RT: 4.943
FID1 A, Front Signal
Correlation: 1.00000
Residual Std. Dev.: 0.00000
Formula: $y = mx$
m: 1.00000
x: Amount Ratio
y: Area Ratio

99




n-Propanol at exp. RT: 7.622
FID2 B, Back Signal
Correlation: 1.00000
Residual Std. Dev.: 0.00000
Formula: $y = mx$
m: 1.00000
x: Amount Ratio
y: Area Ratio

=====

99

ISP Forensic Services Blood Alcohol Report


Sample Name : WATER
 Laboratory : Coeur d' Alene
 Injection Date : May 31, 2020
 Method : ALCOHOL.M
 Acq. Instrument: CN10742044-IT00725005



| # | Compound | Column | Area | Amount | Units |
|----|------------|-----------|---------|--------|---------|
| 1. | Ethanol | Column 1: | 0.00000 | 0.0000 | g/100cc |
| 2. | Ethanol | Column 2: | 0.00000 | 0.0000 | g/100cc |
| 3. | n-Propanol | Column 1: | 0.00000 | 0.0000 | g/100cc |
| 4. | n-Propanol | Column 2: | 0.00000 | 0.0000 | g/100cc |

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.05
 Laboratory : Coeur d' Alene
 Injection Date : May 31, 2020
 Method : ALCOHOL.M
 Acq. Instrument: CN10742044-IT00725005




| # | Compound | Column | Area | Amount | Units |
|----|------------|-----------|----------|--------|---------|
| 1. | Ethanol | Column 1: | 8.95644 | 0.0497 | g/100cc |
| 2. | Ethanol | Column 2: | 8.94641 | 0.0491 | g/100cc |
| 3. | n-Propanol | Column 1: | 93.16283 | 1.0000 | g/100cc |
| 4. | n-Propanol | Column 2: | 90.25261 | 1.0000 | g/100cc |

99

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.100
 Laboratory : Coeur d' Alene
 Injection Date : May 31, 2020
 Method : ALCOHOL.M
 Acq. Instrument: CN10742044-IT00725005




| # | Compound | Column | Area | Amount | Units |
|----|------------|-----------|----------|--------|---------|
| 1. | Ethanol | Column 1: | 17.97128 | 0.1003 | g/100cc |
| 2. | Ethanol | Column 2: | 17.97812 | 0.0995 | g/100cc |
| 3. | n-Propanol | Column 1: | 92.58790 | 1.0000 | g/100cc |
| 4. | n-Propanol | Column 2: | 89.55220 | 1.0000 | g/100cc |

99

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.200
 Laboratory : Coeur d' Alene
 Injection Date : May 31, 2020
 Method : ALCOHOL.M
 Acq. Instrument: CN10742044-IT00725005




| # | Compound | Column | Area | Amount | Units |
|----|------------|-----------|----------|--------|---------|
| 1. | Ethanol | Column 1: | 35.82807 | 0.2010 | g/100cc |
| 2. | Ethanol | Column 2: | 35.96552 | 0.2006 | g/100cc |
| 3. | n-Propanol | Column 1: | 92.07500 | 1.0000 | g/100cc |
| 4. | n-Propanol | Column 2: | 88.86015 | 1.0000 | g/100cc |

99

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.300
 Laboratory : Coeur d' Alene
 Injection Date : May 31, 2020
 Method : ALCOHOL.M
 Acq. Instrument: CN10742044-IT00725005




| # | Compound | Column | Area | Amount | Units |
|----|------------|-----------|----------|--------|---------|
| 1. | Ethanol | Column 1: | 54.45925 | 0.3003 | g/100cc |
| 2. | Ethanol | Column 2: | 54.67428 | 0.3000 | g/100cc |
| 3. | n-Propanol | Column 1: | 93.68510 | 1.0000 | g/100cc |
| 4. | n-Propanol | Column 2: | 90.32391 | 1.0000 | g/100cc |

99

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.500
 Laboratory : Coeur d' Alene
 Injection Date : May 31, 2020
 Method : ALCOHOL.M
 Acq. Instrument: CN10742044-IT00725005




| # | Compound | Column | Area | Amount | Units |
|----|------------|-----------|----------|--------|---------|
| 1. | Ethanol | Column 1: | 89.91124 | 0.4994 | g/100cc |
| 2. | Ethanol | Column 2: | 90.19679 | 0.4999 | g/100cc |
| 3. | n-Propanol | Column 1: | 93.00353 | 1.0000 | g/100cc |
| 4. | n-Propanol | Column 2: | 89.43030 | 1.0000 | g/100cc |

99

ISP Forensic Services Blood Alcohol Report

Sample Name : blank
 Laboratory : Coeur d' Alene
 Injection Date : May 31, 2020
 Method : ALCOHOL.M
 Acq. Instrument: CN10742044-IT00725005



| # | Compound | Column | Area | Amount | Units |
|----|------------|-----------|----------|--------|---------|
| 1. | Ethanol | Column 1: | 0.00000 | 0.0000 | g/100cc |
| 2. | Ethanol | Column 2: | 0.00000 | 0.0000 | g/100cc |
| 3. | n-Propanol | Column 1: | 93.31607 | 1.0000 | g/100cc |
| 4. | n-Propanol | Column 2: | 90.19881 | 1.0000 | g/100cc |

99

S a m p l e S u m m a r y

Sequence table: C:\Chem32\1\TEMP\AESEQ\QS_31.05.2020_11.40.33\5-31-20cal.S
Data directory path: C:\Chem32\1\Data\5-31-20calJJ
Logbook: C:\Chem32\1\Data\5-31-20calJJ\5-31-20cal.LOG
Sequence start: 5/31/2020 11:54:18 AM
Sequence Operator: SYSTEM
Operator: SYSTEM


Method file name: C:\CHEM32\1\METHODS\ALCOHOL.M

| Run # | Location # | Inj # | Sample Name | Sample Amt [g/100cc] | Multip.* Dilution | File name | Cal # | # Cmp |
|-------|------------|-------|-------------|----------------------|-------------------|------------|-------|-------|
| 1 | 1 | 1 | WATER | - | 1.0000 | 001F0101.D | | 0 |
| 2 | 2 | 1 | 0.05 | - | 1.0000 | 002F0201.D | * | 4 |
| 3 | 3 | 1 | 0.100 | - | 1.0000 | 003F0301.D | * | 4 |
| 4 | 4 | 1 | 0.200 | - | 1.0000 | 004F0401.D | * | 4 |
| 5 | 5 | 1 | 0.300 | - | 1.0000 | 005F0501.D | * | 4 |
| 6 | 6 | 1 | 0.500 | - | 1.0000 | 006F0601.D | * | 4 |
| 7 | 7 | 1 | blank | - | 1.0000 | 007F0701.D | | 2 |

99

ISP Forensic Services Blood Alcohol Report

Sample Name : water-1
 Laboratory : Coeur d' Alene
 Injection Date : May 31, 2020
 Method : ALCOHOL.M
 Acq. Instrument: CN10742044-IT00725005




| # | Compound | Column | Area | Amount | Units |
|----|------------|-----------|---------|--------|---------|
| 1. | Ethanol | Column 1: | 0.00000 | 0.0000 | g/100cc |
| 2. | Ethanol | Column 2: | 0.00000 | 0.0000 | g/100cc |
| 3. | n-Propanol | Column 1: | 0.00000 | 0.0000 | g/100cc |
| 4. | n-Propanol | Column 2: | 0.00000 | 0.0000 | g/100cc |

99

ISP Forensic Services Blood Alcohol Report

Sample Name : VOL MIX FN-06041502
 Laboratory : Coeur d' Alene
 Injection Date : May 31, 2020
 Method : ALCOHOL.M
 Acq. Instrument: CN10742044-IT00725005




| # | Compound | Column | Area | Amount | Units |
|----|------------|-----------|-----------|--------|---------|
| 1. | Ethanol | Column 1: | 19.51800 | 0.0847 | g/100cc |
| 2. | Ethanol | Column 2: | 19.52291 | 0.0836 | g/100cc |
| 3. | n-Propanol | Column 1: | 119.08755 | 1.0000 | g/100cc |
| 4. | n-Propanol | Column 2: | 115.77052 | 1.0000 | g/100cc |

99

ISP Forensic Services Blood Alcohol Report

Sample Name : ISTD BLANK-1
 Laboratory : Coeur d' Alene
 Injection Date : May 31, 2020
 Method : ALCOHOL.M
 Acq. Instrument: CN10742044-IT00725005



| # | Compound | Column | Area | Amount | Units |
|----|------------|-----------|----------|--------|---------|
| 1. | Ethanol | Column 1: | 0.00000 | 0.0000 | g/100cc |
| 2. | Ethanol | Column 2: | 0.00000 | 0.0000 | g/100cc |
| 3. | n-Propanol | Column 1: | 90.88826 | 1.0000 | g/100cc |
| 4. | n-Propanol | Column 2: | 87.93625 | 1.0000 | g/100cc |

99

VOLATILES DETERMINATION CASEFILE WORKSHEET

Laboratory No.: QC-1(1)

Analysis Date(s): 31 May 2020

| | Column 1 FID A | Column 2 FID B | Column Precision | Mean Value | Sample A-B Difference | Over-all Mean |
|----------------|-------------------|-------------------|------------------|------------|--------------------------|---------------|
| Sample Results | 0.0772 | 0.0766 | 0.0006 | 0.0769 | 0.0003 | 0.0767 |
| (g/100cc) | 0.0772 | 0.0761 | 0.0011 | 0.0766 | | |

Analysis Method

Refer to Blood Alcohol Method #1

Instrument Information

Instrument information is stored centrally.

Refer to Instrument Method: Alcohol.m

Reporting of Results

Uncertainty of Measurement (UM%): 5.00%

| Overall Mean (g/100cc) | Low | High | 5% of Mean |
|------------------------|-------|-------|------------|
| 0.076 | 0.072 | 0.080 | 0.004 |

| Reported Result | |
|-----------------|--|
| 0.076 | |

Calibration and control data are stored centrally.

99


Revision: 2

Issue Date: 12/23/2019

Issuing Authority: Quality Manager

ISP Forensic Services Blood Alcohol Report

Sample Name : QC-1(1)-A
 Laboratory : Coeur d' Alene
 Injection Date : May 31, 2020
 Method : ALCOHOL.M
 Acq. Instrument: CN10742044-IT00725005




| # | Compound | Column | Area | Amount | Units |
|----|------------|-----------|----------|--------|---------|
| 1. | Ethanol | Column 1: | 13.88048 | 0.0772 | g/100cc |
| 2. | Ethanol | Column 2: | 13.91205 | 0.0766 | g/100cc |
| 3. | n-Propanol | Column 1: | 92.88947 | 1.0000 | g/100cc |
| 4. | n-Propanol | Column 2: | 89.97144 | 1.0000 | g/100cc |

99

ISP Forensic Services Blood Alcohol Report

Sample Name : QC-1(1)-B
 Laboratory : Coeur d' Alene
 Injection Date : May 31, 2020
 Method : ALCOHOL.M
 Acq. Instrument: CN10742044-IT00725005



| # | Compound | Column | Area | Amount | Units |
|----|------------|-----------|----------|--------|---------|
| 1. | Ethanol | Column 1: | 13.98739 | 0.0772 | g/100cc |
| 2. | Ethanol | Column 2: | 13.93441 | 0.0761 | g/100cc |
| 3. | n-Propanol | Column 1: | 93.61187 | 1.0000 | g/100cc |
| 4. | n-Propanol | Column 2: | 90.76152 | 1.0000 | g/100cc |

99

VOLATILES DETERMINATION CASEFILE WORKSHEET

Laboratory No.: 0.08 FN09181807

Analysis Date(s): 31 May 2020

| | Column 1 FID A | Column 2 FID B | Column Precision | Mean Value | Sample A-B Difference | Over-all Mean |
|----------------|-------------------|-------------------|------------------|------------|--------------------------|---------------|
| Sample Results | 0.0811 | 0.0803 | 0.0008 | 0.0807 | 0.0010 | 0.0802 |
| (g/100cc) | 0.0797 | 0.0797 | 0.0000 | 0.0797 | | |

Analysis Method

Refer to Blood Alcohol Method #1

Instrument Information

Instrument information is stored centrally.

Refer to Instrument Method: Alcohol.m

Reporting of Results

Uncertainty of Measurement (UM%): 5.00%

| Overall Mean (g/100cc) | Low | High | 5% of Mean |
|------------------------|-------|-------|------------|
| 0.080 | 0.076 | 0.084 | 0.004 |


| Reported Result | |
|-----------------|--|
| 0.080 | |

Calibration and control data are stored centrally.

99

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.08 FN09181807-A
 Laboratory : Coeur d' Alene
 Injection Date : May 31, 2020
 Method : ALCOHOL.M
 Acq. Instrument: CN10742044-IT00725005




| # | Compound | Column | Area | Amount | Units |
|----|------------|-----------|----------|--------|---------|
| 1. | Ethanol | Column 1: | 14.59977 | 0.0811 | g/100cc |
| 2. | Ethanol | Column 2: | 14.58305 | 0.0803 | g/100cc |
| 3. | n-Propanol | Column 1: | 92.99659 | 1.0000 | g/100cc |
| 4. | n-Propanol | Column 2: | 89.96109 | 1.0000 | g/100cc |

99

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.08 FN09181807-B
 Laboratory : Coeur d' Alene
 Injection Date : May 31, 2020
 Method : ALCOHOL.M
 Acq. Instrument: CN10742044-IT00725005



| # | Compound | Column | Area | Amount | Units |
|----|------------|-----------|----------|--------|---------|
| 1. | Ethanol | Column 1: | 14.29927 | 0.0797 | g/100cc |
| 2. | Ethanol | Column 2: | 14.41291 | 0.0797 | g/100cc |
| 3. | n-Propanol | Column 1: | 92.63902 | 1.0000 | g/100cc |
| 4. | n-Propanol | Column 2: | 89.68399 | 1.0000 | g/100cc |

99

VOLATILES DETERMINATION CASEFILE WORKSHEET

Laboratory No.: QC-2(1)

Analysis Date(s): 31 May 2020

| | Column 1 FID A | Column 2 FID B | Column Precision | Mean Value | Sample A-B Difference | Over-all Mean |
|----------------|-------------------|-------------------|------------------|------------|--------------------------|---------------|
| Sample Results | 0.1978 | 0.1969 | 0.0009 | 0.1973 | 0.0015 | 0.1981 |
| (g/100cc) | 0.1993 | 0.1984 | 0.0009 | 0.1988 | | |

Analysis Method

Refer to Blood Alcohol Method #1

Instrument Information

Instrument information is stored centrally.

Refer to Instrument Method: Alcohol.m

Reporting of Results

Uncertainty of Measurement (UM%): 5.00%

| Overall Mean (g/100cc) | Low | High | 5% of Mean |
|------------------------|-------|-------|------------|
| 0.198 | 0.188 | 0.208 | 0.010 |

| Reported Result | |
|-----------------|--|
| 0.198 | |

Calibration and control data are stored centrally.


Revision: 2

Issue Date: 12/23/2019

Issuing Authority: Quality Manager

ISP Forensic Services Blood Alcohol Report

Sample Name : QC-2(1)-A
 Laboratory : Coeur d' Alene
 Injection Date : May 31, 2020
 Method : ALCOHOL.M
 Acq. Instrument: CN10742044-IT00725005




| # | Compound | Column | Area | Amount | Units |
|----|------------|-----------|----------|--------|---------|
| 1. | Ethanol | Column 1: | 36.02095 | 0.1978 | g/100cc |
| 2. | Ethanol | Column 2: | 36.06114 | 0.1969 | g/100cc |
| 3. | n-Propanol | Column 1: | 94.08472 | 1.0000 | g/100cc |
| 4. | n-Propanol | Column 2: | 90.77829 | 1.0000 | g/100cc |

29

ISP Forensic Services Blood Alcohol Report

Sample Name : QC-2(1)-B
 Laboratory : Coeur d' Alene
 Injection Date : May 31, 2020
 Method : ALCOHOL.M
 Acq. Instrument: CN10742044-IT00725005




| # | Compound | Column | Area | Amount | Units |
|----|------------|-----------|----------|--------|---------|
| 1. | Ethanol | Column 1: | 36.37951 | 0.1993 | g/100cc |
| 2. | Ethanol | Column 2: | 36.39520 | 0.1984 | g/100cc |
| 3. | n-Propanol | Column 1: | 94.28717 | 1.0000 | g/100cc |
| 4. | n-Propanol | Column 2: | 90.92138 | 1.0000 | g/100cc |

99

ISP Forensic Services Blood Alcohol Report

Sample Name : ISTD BLANK-2
 Laboratory : Coeur d' Alene
 Injection Date : May 31, 2020
 Method : ALCOHOL.M
 Acq. Instrument: CN10742044-IT00725005




| # | Compound | Column | Area | Amount | Units |
|----|------------|-----------|----------|--------|---------|
| 1. | Ethanol | Column 1: | 0.00000 | 0.0000 | g/100cc |
| 2. | Ethanol | Column 2: | 0.00000 | 0.0000 | g/100cc |
| 3. | n-Propanol | Column 1: | 91.54604 | 1.0000 | g/100cc |
| 4. | n-Propanol | Column 2: | 88.66466 | 1.0000 | g/100cc |

99

ISP Forensic Services Blood Alcohol Report

Sample Name : water-2
 Laboratory : Coeur d' Alene
 Injection Date : May 31, 2020
 Method : ALCOHOL.M
 Acq. Instrument: CN10742044-IT00725005




| # | Compound | Column | Area | Amount | Units |
|----|------------|-----------|---------|--------|---------|
| 1. | Ethanol | Column 1: | 0.00000 | 0.0000 | g/100cc |
| 2. | Ethanol | Column 2: | 0.00000 | 0.0000 | g/100cc |
| 3. | n-Propanol | Column 1: | 0.00000 | 0.0000 | g/100cc |
| 4. | n-Propanol | Column 2: | 0.00000 | 0.0000 | g/100cc |

99

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.05 CHECK
 Laboratory : Coeur d' Alene
 Injection Date : May 31, 2020
 Method : ALCOHOL.M
 Acq. Instrument: CN10742044-IT00725005




| # | Compound | Column | Area | Amount | Units |
|----|------------|-----------|----------|--------|---------|
| 1. | Ethanol | Column 1: | 9.03751 | 0.0512 | g/100cc |
| 2. | Ethanol | Column 2: | 8.98211 | 0.0505 | g/100cc |
| 3. | n-Propanol | Column 1: | 91.20735 | 1.0000 | g/100cc |
| 4. | n-Propanol | Column 2: | 88.08714 | 1.0000 | g/100cc |

29

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.100 CHECK
 Laboratory : Coeur d' Alene
 Injection Date : May 31, 2020
 Method : ALCOHOL.M
 Acq. Instrument: CN10742044-IT00725005




| # | Compound | Column | Area | Amount | Units |
|----|------------|-----------|----------|--------|---------|
| 1. | Ethanol | Column 1: | 18.42128 | 0.1019 | g/100cc |
| 2. | Ethanol | Column 2: | 18.33600 | 0.1008 | g/100cc |
| 3. | n-Propanol | Column 1: | 93.38719 | 1.0000 | g/100cc |
| 4. | n-Propanol | Column 2: | 90.12579 | 1.0000 | g/100cc |

99

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.200 CHECK
 Laboratory : Coeur d' Alene
 Injection Date : May 31, 2020
 Method : ALCOHOL.M
 Acq. Instrument: CN10742044-IT00725005




| # | Compound | Column | Area | Amount | Units |
|----|------------|-----------|----------|--------|---------|
| 1. | Ethanol | Column 1: | 36.49109 | 0.2044 | g/100cc |
| 2. | Ethanol | Column 2: | 36.39232 | 0.2030 | g/100cc |
| 3. | n-Propanol | Column 1: | 92.20997 | 1.0000 | g/100cc |
| 4. | n-Propanol | Column 2: | 88.87759 | 1.0000 | g/100cc |

99

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.300 CHECK
 Laboratory : Coeur d' Alene
 Injection Date : May 31, 2020
 Method : ALCOHOL.M
 Acq. Instrument: CN10742044-IT00725005




| # | Compound | Column | Area | Amount | Units |
|----|------------|-----------|----------|--------|---------|
| 1. | Ethanol | Column 1: | 53.82040 | 0.3070 | g/100cc |
| 2. | Ethanol | Column 2: | 53.70264 | 0.3052 | g/100cc |
| 3. | n-Propanol | Column 1: | 90.56096 | 1.0000 | g/100cc |
| 4. | n-Propanol | Column 2: | 87.21037 | 1.0000 | g/100cc |

99

ISP Forensic Services Blood Alcohol Report

Sample Name : 0.500 CHECK
 Laboratory : Coeur d' Alene
 Injection Date : May 31, 2020
 Method : ALCOHOL.M
 Acq. Instrument: CN10742044-IT00725005



| # | Compound | Column | Area | Amount | Units |
|----|------------|-----------|----------|--------|---------|
| 1. | Ethanol | Column 1: | 91.30092 | 0.5106 | g/100cc |
| 2. | Ethanol | Column 2: | 91.59035 | 0.5127 | g/100cc |
| 3. | n-Propanol | Column 1: | 92.36597 | 1.0000 | g/100cc |
| 4. | n-Propanol | Column 2: | 88.54821 | 1.0000 | g/100cc |

99