








B. Wylee

Worklist: 1943

<u>LAB CASE</u>	<u>ITEM</u>	<u>TASK ID</u>	<u>DESCRIPTION</u>	
C2017-1925	1	97412	AM 27 Blood THC Quant by LC	
C2017-1938	1	97413	AM 27 Blood THC Quant by LC	
C2017-1969	1	97415	AM 27 Blood THC Quant by LC	
C2017-1974	1	97414	AM 27 Blood THC Quant by LC	
C2017-2016	1	97416	AM 27 Blood THC Quant by LC	
C2017-2030	1	97417	AM 27 Blood THC Quant by LC	
M2017-4292	1	97418	AM 27 Blood THC Quant by LC	



Quantitation of THC and Metabolites in Blood by LC-MS/MS

Extraction Date: 10-11-2017

Analyst: Anne Nord

PRE-ANALYTIC

Plate Lot# Custom - 499102 Plate Exp. 1/29/2018 External QC Lot 21718 exp 2-17-18


- ✓ 1. Ensure all solutions are within expiration date.
 - Mobile Phase A: *0.1% Formic Acid in LCMS Water* • *0.1% Formic Acid in water*
 - Mobile Phase B: *0.1% Formic Acid in LCMS Acetonitrile* • *MTBE*
 - *LCMS Methanol* • *Hexane*
 - **Blank/Negative Blood: Lot 321632-1**
- Column: *UCT Selectra DA 100 x 2.1 mm 3um*
- ✓ 2. Check levels of mobile phases and needle wash and refill as necessary. Ensure waste is not full.
- ✓ 3. Purge Pump and Load appropriate Acq. Method, allow system to equilibrate for approx. 30 min.
- ✓ 4. Create worklist. Data path name: 101117 cann quant r

ANALYTIC

- ✓ 1. Remove standards plate, blood, and samples from cold storage. Allow to reach room temperature.
- ✓ 2. Add **1000 µL blood** to wells of analytical (standards) plate. Place cover on Plate
 - Blank blood for locations containing standards/QCs and internal standards
 - Sample blood for locations containing only internal standards
- ✓ 3. Place on shaking incubator at ambient temp., **900rpm for 15 minutes**. *Shaker ID 66759*
- ✓ 4. Pipette **500µL 0.1% formic acid** to all wells of standards plate.
- ✓ 5. Place on shaking incubator at ambient temp., **900rpm for 15 minutes**.
- ✓ 6. Transfer **800µL of blood+acid mixture** to corresponding wells of SLE+ plate.
- ✓ 7. Apply positive pressure for approx. 10-15 seconds (or until no liquid remains on top of sorbent). Wait 5 min. (**Load blood samples at 85- 100 PSI- Selector to Right**)
- ✓ 8. Add **2.25mL MTBE** and allow to flow under gravity for 5 minutes. (**add in 3 increments of 750uL**)
- ✓ 9. Apply positive pressure for approx. 15 seconds (**10-15 PSI- Selector to left -**).
- ✓ 10. Add **2.25mL Hexane** and allow to flow under gravity for 5 minutes. (**add in 3 increments of 750uL**)
- ✓ 11. Apply positive pressure for approx. 15 seconds. (**10-15 PSI Selector to the left**)
- ✓ 12. Remove collection plate containing eluate.
- ✓ 13. Place collection plate on SPE Dry and evaporate to dryness at approx. 35°C. *SPE Dry ID 66819*
- ✓ 14. Reconstitute in **100 µL MeOH** and heat seal plate with foil. Place in autosampler and run worklist.

POST-ANALYTIC

- ✓ 1. Open quantitation software and create a new quantitation batch.
Batch name: Cann quant
- ✓ 2. Make any necessary integration changes. Limit curves based on validated linear ranges (3-50ng/mL).
- ✓ 3. Were all appropriate standards used in the curve for each analyte? Y / N
Are r^2 values ≥ 0.98 for each analyte? Y / N
- ✓ 4. Did all QCs pass for each analyte? Y / N Were QCs entered into QC charting? Y / N
- ✓ 5. Central File Packet to include: ✓ LIMS Worklist: ✓ Method Checklist ✓ Calibration and Control Reports

COMMENTS Samples run in data path 101117 cann quant. The blank run did not have internal standard in it. Samples were reconstituted and run 10-13-17. 




ISP Forensics Calibration Curve Report

Batch Data Path D:\2017 Data\101117 cann quant r\QuantResults\cann quant.batch.bin

Last Calib Update 10/15/2017 11:35 AM **Analyst Name** ISP TOX

Target Compound *THC-OH*
Internal Standard *THC-OH-d3*




Sample	Level	Enabled	Exp Conc	Final Conc	Accuracy
Cal 1 - 3ng	1	<input checked="" type="checkbox"/>	3	3.3	109.5
Cal 2 - 5ng	2	<input checked="" type="checkbox"/>	5	4.4	88.3
Cal 3 - 10ng	3	<input checked="" type="checkbox"/>	10	10.4	104.0
QC - 10ng	3	<input checked="" type="checkbox"/>	10	9.3	92.9
Cal 4 - 25ng	4	<input checked="" type="checkbox"/>	25	23.9	95.5
Cal 5 - 50ng	5	<input checked="" type="checkbox"/>	50	51.3	102.6
Cal 6 - 100ng	6	<input checked="" type="checkbox"/>	100	100.1	100.1
Cal 7 - 250ng	7	<input checked="" type="checkbox"/>	250	249.6	99.8

ISP Forensics Calibration Curve Report

Batch Data Path D:\2017 Data\101117 cann quant r\QuantResults\cann quant.batch.bin

Last Calib Update 10/15/2017 11:35 AM **Analyst Name** ISP TOX

Target Compound *THC-COOH*
Internal Standard *THC-COOH-d9*




Sample	Level	Enabled	Exp Conc	Final Conc	Accuracy
Cal 1 - 3ng	1	<input checked="" type="checkbox"/>	3	2.9	97.6
Cal 2 - 5ng	2	<input checked="" type="checkbox"/>	5	5.2	103.5
Cal 3 - 10ng	3	<input checked="" type="checkbox"/>	10	10.2	102.0
QC - 10ng	3	<input checked="" type="checkbox"/>	10	10.9	109.1
Cal 4 - 25ng	4	<input checked="" type="checkbox"/>	25	24.1	96.3
Cal 5 - 50ng	5	<input checked="" type="checkbox"/>	50	51.5	103.0
Cal 6 - 100ng	6	<input checked="" type="checkbox"/>	100	97.0	97.0
Cal 7 - 250ng	7	<input checked="" type="checkbox"/>	250	251.6	100.6

ISP Forensics Calibration Curve Report

Batch Data Path D:\2017 Data\101117 cann quant r\QuantResults\cann quant.batch.bin

Last Calib Update 10/15/2017 11:35 AM **Analyst Name** ISP TOX

Target Compound *THC*
Internal Standard *THC-d3*



Sample	Level	Enabled	Exp Conc	Final Conc	Accuracy
Cal 1 - 3ng	1	<input checked="" type="checkbox"/>	3	2.9	97.3
Cal 2 - 5ng	2	<input checked="" type="checkbox"/>	5	5.6	111.2
Cal 3 - 10ng	3	<input checked="" type="checkbox"/>	10	8.9	88.6
QC - 10ng	3	<input checked="" type="checkbox"/>	10	10.6	106.1
Cal 4 - 25ng	4	<input checked="" type="checkbox"/>	25	24.9	99.8
Cal 5 - 50ng	5	<input checked="" type="checkbox"/>	50	52.8	105.7
Cal 6 - 100ng	6	<input checked="" type="checkbox"/>	100	97.2	97.2
Cal 7 - 250ng	7	<input checked="" type="checkbox"/>	250	250.7	100.3

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2017 Data\101117 cann quant r\QuantResults\cann quant.batch.bin
Analysis Time 10/15/2017 11:35 AM **Analyst Name** ISP Tox
Report Time 10/15/2017 11:37 AM **Reporter Name** ISP Tox
Last Calib Update 10/15/2017 11:35 AM **Batch State** Processed

Analysis Info

Acq Time 2017-10-13 12:05 **Data File** Negative Control.d
Sample Type Sample **Sample Name** Negative Control
Dilution 1 **Acq Method** AM 27 Quant THC 7-2017.m
Position P1-A2 **Sample Info**
Inj Vol -1 **Comment** AM 27 Cannabinoid Confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-COOH	THC-COOH-d9	2.285	4371	29136	0.1500	0.2790

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2017 Data\101117 cann quant r\QuantResults\cann quant.batch.bin
Analysis Time 10/15/2017 11:35 AM **Analyst Name** ISP Tox
Report Time 10/15/2017 11:37 AM **Reporter Name** ISP Tox
Last Calib Update 10/15/2017 11:35 AM **Batch State** Processed

Analysis Info

Acq Time 2017-10-13 12:17 **Data File** QC - 10ng.d
Sample Type QC **Sample Name** QC - 10ng
Dilution 1 **Acq Method** AM 27 Quant THC 7-2017.m
Position P1-H1 **Sample Info**
Inj Vol -1 **Comment** AM 27 Cannabinoid Confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.115	10537	107151	0.0983	9.2875
THC-COOH	THC-COOH-d9	2.245	14259	39525	0.3608	10.9128
THC	THC-d3	5.572	5237	39074	0.1340	10.6066

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2017 Data\101117 cann quant r\QuantResults\cann quant.batch.bin
Analysis Time 10/15/2017 11:35 AM **Analyst Name** ISP Tox
Report Time 10/15/2017 11:37 AM **Reporter Name** ISP Tox
Last Calib Update 10/15/2017 11:35 AM **Batch State** Processed

Analysis Info

Acq Time 2017-10-13 12:29 **Data File** External Control lot 21718.d
Sample Type Sample **Sample Name** External Control lot 21718
Dilution 1 **Acq Method** AM 27 Quant THC 7-2017.m
Position P1-B2 **Sample Info**
Inj Vol -1 **Comment** AM 27 Cannabinoid Confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.135	2076	24928	0.0833	7.8283
THC-COOH	THC-COOH-d9	2.205	3351	9980	0.3358	9.6542
THC	THC-d3	5.412	909	8329	0.1091	8.6099

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2017 Data\101117 cann quant r\QuantResults\cann quant.batch.bin
Analysis Time 10/15/2017 11:35 AM **Analyst Name** ISP Tox
Report Time 10/15/2017 11:37 AM **Reporter Name** ISP Tox
Last Calib Update 10/15/2017 11:35 AM **Batch State** Processed

Analysis Info

Acq Time 2017-10-13 11:18 **Data File** Cal 1 - 3ngr.d
Sample Type Calibration **Sample Name** Cal 1 - 3ng
Dilution 1 **Acq Method** AM 27 Quant THC 7-2017.m
Position P1-A1 **Sample Info**
Inj Vol -1 **Comment** AM 27 Cannabinoid Confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.155	1178	32280	0.0365	3.2864
THC-COOH	THC-COOH-d9	2.265	2659	13130	0.2025	2.9290
THC	THC-d3	5.532	417	10922	0.0382	2.9186

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2017 Data\101117 cann quant r\QuantResults\cann quant.batch.bin
Analysis Time 10/15/2017 11:35 AM **Analyst Name** ISP Tox
Report Time 10/15/2017 11:37 AM **Reporter Name** ISP Tox
Last Calib Update 10/15/2017 11:35 AM **Batch State** Processed

Analysis Info

Acq Time 2017-10-13 11:29 **Data File** Cal 2 - 5ngr.d
Sample Type Calibration **Sample Name** Cal 2 - 5ng
Dilution 1 **Acq Method** AM 27 Quant THC 7-2017.m
Position P1-B1 **Sample Info**
Inj Vol -1 **Comment** AM 27 Cannabinoid Confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.155	2912	60495	0.0481	4.4164
THC-COOH	THC-COOH-d9	2.265	5915	23939	0.2471	5.1763
THC	THC-d3	5.492	1603	22527	0.0711	5.5616

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2017 Data\101117 cann quant r\QuantResults\cann quant.batch.bin
Analysis Time 10/15/2017 11:35 AM **Analyst Name** ISP Tox
Report Time 10/15/2017 11:37 AM **Reporter Name** ISP Tox
Last Calib Update 10/15/2017 11:35 AM **Batch State** Processed

Analysis Info

Acq Time 2017-10-13 11:41 **Data File** Cal 3 - 10ngr.d
Sample Type Calibration **Sample Name** Cal 3 - 10ng
Dilution 1 **Acq Method** AM 27 Quant THC 7-2017.m
Position P1-C1 **Sample Info**
Inj Vol -i **Comment** AM 27 Cannabinoid Confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.135	10340	94146	0.1098	10.4031
THC-COOH	THC-COOH-d9	2.245	12365	35669	0.3467	10.2019
THC	THC-d3	5.592	3937	35095	0.1122	8.8552

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2017 Data\101117 cann quant r\QuantResults\cann quant.batch.bin
Analysis Time 10/15/2017 11:35 AM **Analyst Name** ISP Tox
Report Time 10/15/2017 11:37 AM **Reporter Name** ISP Tox
Last Calib Update 10/15/2017 11:35 AM **Batch State** Processed

Analysis Info

Acq Time 2017-10-13 11:53 **Data File** Cal 4 - 25ngr.d
Sample Type Calibration **Sample Name** Cal 4 - 25ng
Dilution 1 **Acq Method** AM 27 Quant THC 7-2017.m
Position P1-D1 **Sample Info**
Inj Vol -1 **Comment** AM 27 Cannabinoid Confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.135	27135	109116	0.2487	23.8773
THC-COOH	THC-COOH-d9	2.225	26413	42510	0.6213	24.0627
THC	THC-d3	5.552	12218	39061	0.3128	24.9473

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2017 Data\101117 cann quant r\QuantResults\cann quant.batch.bin
Analysis Time 10/15/2017 11:35 AM **Analyst Name** ISP Tox
Report Time 10/15/2017 11:37 AM **Reporter Name** ISP Tox
Last Calib Update 10/15/2017 11:35 AM **Batch State** Processed

Analysis Info

Acq Time 2017-10-13 10:42 **Data File** Cal 5 - 50ng.d
Sample Type Calibration **Sample Name** Cal 5 - 50ng
Dilution 1 **Acq Method** AM 27 Quant THC 7-2017.m
Position P1-E1 **Sample Info**
Inj Vol -1 **Comment** AM 27 Cannabinoid Confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.135	52333	98461	0.5315	51.3230
THC-COOH	THC-COOH-d9	2.225	39503	33912	1.1649	51.4891
THC	THC-d3	5.532	22269	33718	0.6604	52.8346

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2017 Data\101117 cann quant r\QuantResults\cann quant.batch.bin
Analysis Time 10/15/2017 11:35 AM **Analyst Name** ISP Tox
Report Time 10/15/2017 11:37 AM **Reporter Name** ISP Tox
Last Calib Update 10/15/2017 11:35 AM **Batch State** Processed

Analysis Info

Acq Time 2017-10-13 10:54 **Data File** Cal 6 - 100ng.d
Sample Type Calibration **Sample Name** Cal 6 - 100ng
Dilution 1 **Acq Method** AM 27 Quant THC 7-2017.m
Position P1-F1 **Sample Info**
Inj Vol -1 **Comment** AM 27 Cannabinoid Confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.135	78527	75934	1.0342	100.0997
THC-COOH	THC-COOH-d9	2.225	58210	28174	2.0661	96.9645
THC	THC-d3	5.512	28557	23533	1.2135	97.1997

ISP FORENSICS - Cd'A Instrument # 62340


Cannabinoids Analysis Report

Batch Data Path D:\2017 Data\101117 cann quant r\QuantResults\cann quant.batch.bin
Analysis Time 10/15/2017 11:35 AM **Analyst Name** ISP Tox
Report Time 10/15/2017 11:37 AM **Reporter Name** ISP Tox
Last Calib Update 10/15/2017 11:35 AM **Batch State** Processed

Analysis Info

Acq Time 2017-10-13 11:06 **Data File** Cal 7 - 250ng.d
Sample Type Calibration **Sample Name** Cal 7 - 250ng
Dilution 1 **Acq Method** AM 27 Quant THC 7-2017.m
Position P1-G1 **Sample Info**
Inj Vol -1 **Comment** AM 27 Cannabinoid Confirmation

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	2.135	256117	99475	2.5747	249.5940
THC-COOH	THC-COOH-d9	2.205	176125	34332	5.1300	251.5714
THC	THC-d3	5.552	106791	34153	3.1269	250.6831