


*BWylee*

**Worklist: 2524**

<u>LAB_CASE</u>	<u>ITEM</u>	<u>TASK_ID</u>	<u>DESCRIPTION</u>
C2018-1077	1	119750	AM 27 Blood THC Quant by LC-QQQ
C2018-1081	1	119751	AM 27 Blood THC Quant by LC-QQQ
C2018-1136	1	119752	AM 27 Blood THC Quant by LC-QQQ
C2018-1186	1	119753	AM 27 Blood THC Quant by LC-QQQ
C2018-1206	1	119754	AM 27 Blood THC Quant by LC-QQQ
M2018-2862	3	119755	AM 27 Blood THC Quant by LC-QQQ



*A*

# AM# 27: Quantitation of THC and Metabolites in Blood by LC-MS/MS

Extraction Date: 6-26-18

Analyst: Anne Nord

Plate lot#: 0515037

Plate Expiration: 9/28/2018

**Mobile phase A:** 0.1% Formic Acid in LCMS Water  
MTBE

**Mobile phase B:** 0.1% Formic acid in Acetonitrile  
LCMS Methanol Hexane

**Blank Blood Lot:** 17J0718

**Column:** UCT Selectra DA 100 x 2.1mm 3um

**LCMS-QQQ ID:** 62340

## Pre-Analytic:

1. Check levels of mobile phases and needle wash refill as needed. Ensure waste is not full.
2. Ensure correct column is installed and begin mobile phase flow allow to equilibrate ~ 30 minutes.
3. Create worklist:

## Analytic:

1. Remove standards, plate, controls, and samples from cold storage. Allow to reach room temperature.
2. Pipette **1000µL blood (calibrated pipette) Pipette ID: k52558g** in wells of analytical (standards) plate.
3. Place on shaking incubator at ambient temp., 900rpm for 15 minutes. *Shaker ID: 66759*
4. Pipette **500µL 0.1% formic acid in water** in wells of analytical plate.
5. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
6. Transfer **800µL of blood+acid** mixture to corresponding wells of SLE+ plate.
7. Apply positive pressure for approx. 10-15 seconds (or until no liquid remains on top of sorbent).  
*(Load at 85-100 PSI- Selector to the right) Manifold ID: 66792*
8. Wait 5 minutes.
9. Add **2.25mL MTBE. (Add in 3 increments of 750uL)**
10. Wait 5 minutes.
11. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
12. Add **2.25mL Hexane. (Add in 3 increments of 750uL)**
13. Wait 5 minutes.
14. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
15. Remove plate containing eluate. Place on SPE Dry and evaporate to dryness at approx. 35°C.  
*SPE Dry ID: 66819*
16. Reconstitute in **100µL 100% MeOH** and heat seal plate with foil. Place in autosampler and run worklist.

## Post-Analytic

1. Create batch and process data.  
Worklist path: 06272018 cann quant Batch Name: cann quant
2. Make any necessary integration changes, Curve weighting of Linear 1/x with  $r^2$  values  $\geq 0.98$  for each analyte
3. RT +/- 3% or 0.100 min, whichever is greater, +/- 20% Accuracy for greater than (+/- 30% for 10ng/ml or less).  
Ion ratios must be within +/- 20% of the averaged calibrators
4. Case sample response for THC and OH-THC 3ng/mL (quantitative), Carboxy-THC: 10ng/mL (qualitative only) will be reported. Samples with a THC or OH-THC response over 50 ng/mL will be reported out as greater than 50 ng/mL.
5. Did all QCs pass for each analyte? Y/N- *QC internal control THC-OH did not pass. THC-OH will not be evaluated in this run.*
6. Enter QCs into control charting.
7. Central File Packet to include: LIMS Worklist, Method Checklist, Calibration and Control Reports

COMMENTS: *curves limited 3-100 ng THC and Carboxy-THC*


*THC-OH not evaluated internal control failed.*

I started the run and evaluated the initial data. The calibrators, negative control, internal control, and external control were re-injected, THC was not coming out in the acquisition window.

Toxicology AM method 27 external prep information  
working solution 1 ug/ml in meoh C-THC, THC-OH, THC  
Stock solution 1mg/ml 10 ul each THC, THC-OH 100 ug/ml 100 ul C-THC in 9890 ul meOH  
Ppd 6/5/18 Exp: 4/1/19 lot 6518 by AMN

Drug	lot (cerilliant)	expiration
C-THC	FE03121501	3/1/2020
THC-OH	FE01141502	1/1/2020
THC	FE04231406	4/1/2019

AM 27 control 50 ul working solution lot (6518) in 4950 ul blood lot (17J20718)  
ppd 6/5/18 Exp 4/1/19 lot 6518 Concentration 10 ng/ml each by AMN



# ISP FORENSICS - Cd'A Instrument # 62340


## Cannabinoids Analysis Report

Batch Data Path D:\2018 Data\06272018 cann quant\QuantResults\cann quant.batch.bin  
Analysis Time 6/28/2018 10:38 AM Analyst Name ISP Tox  
Report Time 6/28/2018 10:39 AM Reporter Name ISP Tox  
Last Calib Update 6/28/2018 10:38 AM Batch State Processed

### Analysis Info

Acq Time 2018-06-27 16:53 Data File Negative Control.r.d  
Sample Type Sample Sample Name Negative Control  
Dilution 1 Acq Method AM 27 Quant THC 7-2017.m  
Position P1-A2 Sample Info  
Inj Vol -1 Comment AM 27 Cannabinoid Confirmation

### Sample Chromatogram



### Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-COOH	THC-COOH-d9	1.585	3536	118086	0.0299	1.1155 <i>L10</i>

# ISP FORENSICS - Cd'A Instrument # 62340


## Cannabinoids Analysis Report

**Batch Data Path** D:\2018 Data\06272018 cann quant\QuantResults\cann quant.batch.bin  
**Analysis Time** 6/28/2018 10:38 AM **Analyst Name** ISP Tox  
**Report Time** 6/28/2018 10:39 AM **Reporter Name** ISP Tox  
**Last Calib Update** 6/28/2018 10:38 AM **Batch State** Processed

### Analysis Info

**Acq Time** 2018-06-27 17:05 **Data File** QC - 10ngr.d  
**Sample Type** QC **Sample Name** QC - 10ng  
**Dilution** 1 **Acq Method** AM 27 Quant THC 7-2017.m  
**Position** P1-H1 **Sample Info**  
**Inj Vol** -1 **Comment** AM 27 Cannabinoid Confirmation

### Sample Chromatogram



### Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	1.755	45466	464020	0.0980	8.7963 - ratio out
THC-COOH	THC-COOH-d9	1.805	98534	141670	0.6955	29.7347
THC	THC-d3	4.091	16251	125310	0.1297	9.3139

THC-OH will not be evaluated in this run the qualifier ratio is out of range in this control. ✖

# ISP FORENSICS - Cd'A Instrument # 62340


## Cannabinoids Analysis Report

**Batch Data Path** D:\2018 Data\06272018 cann quant\QuantResults\cann quant.batch.bin  
**Analysis Time** 6/28/2018 10:38 AM **Analyst Name** ISP Tox  
**Report Time** 6/28/2018 10:39 AM **Reporter Name** ISP Tox  
**Last Calib Update** 6/28/2018 10:38 AM **Batch State** Processed

### Analysis Info

**Acq Time** 2018-06-27 17:17 **Data File** External Control lot 6518r.d  
**Sample Type** QC **Sample Name** External Control lot 6518  
**Dilution** 1 **Acq Method** AM 27 Quant THC 7-2017.m  
**Position** P1-B2 **Sample Info**  
**Inj Vol** -1 **Comment** AM 27 Cannabinoid Confirmation

### Sample Chromatogram



### Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	1.755	25443	309997	0.0821	7.4257 - Not evaluated
THC-COOH	THC-COOH-d9	1.805	17246	102990	0.1675	7.0283
THC	THC-d3	4.051	12121	95202	0.1273	9.1482


# ISP Forensics Calibration Curve Report

\*Not Evaluated in this batch

**Batch Data Path** D:\2018 Data\06272018 cann quant\QuantResults\cann quant.batch.bin

**Last Calib Update** 6/28/2018 7:47 AM **Analyst Name** ISP TOX

**Target Compound** *THC-OH*  
**Internal Standard** *THC-OH-d3*




Sample	Level	Enabled	Exp Conc	Final Conc	Accuracy
Cal 1 - 3ng	1	<input checked="" type="checkbox"/>	3	3.3	109.6
Cal 2 - 5ng	2	<input checked="" type="checkbox"/>	5	5.0	100.9
Cal 3 - 10ng	3	<input checked="" type="checkbox"/>	10	9.5	95.4
QC - 10ng	3	<input checked="" type="checkbox"/>	10	8.8	88.0
External Control lot 6518	3	<input checked="" type="checkbox"/>	10	7.4	74.3
Cal 4 - 25ng	4	<input checked="" type="checkbox"/>	25	24.0	95.8
Cal 5 - 50ng	5	<input checked="" type="checkbox"/>	50	47.0	94.1
Cal 6 - 100ng	6	<input checked="" type="checkbox"/>	100	104.1	104.1
Cal 7 - 250ng	7	<input type="checkbox"/>	250	189.9	76.0

# ISP Forensics Calibration Curve Report

**Batch Data Path** D:\2018 Data\06272018 cann quant\QuantResults\cann quant.batch.bin

**Last Calib Update** 6/28/2018 7:47 AM **Analyst Name** ISP TOX

**Target Compound** *THC-COOH*  
**Internal Standard** *THC-COOH-d9*



Sample	Level	Enabled	Exp Conc	Final Conc	Accuracy
Cal 1 - 3ng	1	<input checked="" type="checkbox"/>	3	3.0	99.6
Cal 2 - 5ng	2	<input checked="" type="checkbox"/>	5	4.9	98.8
Cal 3 - 10ng	3	<input checked="" type="checkbox"/>	10	11.3	113.0
QC - 10ng	3	<input checked="" type="checkbox"/>	10	30.2	302.1
External Control lot 6518	3	<input checked="" type="checkbox"/>	10	7.0	70.3
Cal 4 - 25ng	4	<input checked="" type="checkbox"/>	25	23.2	92.6
Cal 5 - 50ng	5	<input checked="" type="checkbox"/>	50	45.3	90.6
Cal 6 - 100ng	6	<input checked="" type="checkbox"/>	100	105.3	105.3
Cal 7 - 250ng	7	<input type="checkbox"/>	250	212.6	85.0




# ISP Forensics Calibration Curve Report

**Batch Data Path** D:\2018 Data\06272018 cann quant\QuantResults\cann quant.batch.bin

**Last Calib Update** 6/28/2018 7:47 AM **Analyst Name** ISP TOX

**Target Compound** *THC*  
**Internal Standard** *THC-d3*



Sample	Level	Enabled	Exp Conc	Final Conc	Accuracy
Cal 1 - 3ng	1	<input checked="" type="checkbox"/>	3	3.2	106.4
Cal 2 - 5ng	2	<input checked="" type="checkbox"/>	5	5.1	102.9
Cal 3 - 10ng	3	<input checked="" type="checkbox"/>	10	9.7	96.7
QC - 10ng	3	<input checked="" type="checkbox"/>	10	9.3	93.1
External Control lot 6518	3	<input checked="" type="checkbox"/>	10	9.1	91.5
Cal 4 - 25ng	4	<input checked="" type="checkbox"/>	25	23.9	95.6
Cal 5 - 50ng	5	<input checked="" type="checkbox"/>	50	47.4	94.8
Cal 6 - 100ng	6	<input checked="" type="checkbox"/>	100	103.7	103.7
Cal 7 - 250ng	7	<input type="checkbox"/>	250	201.8	80.7

# ISP FORENSICS - Cd'A Instrument # 62340


## Cannabinoids Analysis Report

**Batch Data Path** D:\2018 Data\06272018 cann quant\QuantResults\cann quant.batch.bin  
**Analysis Time** 6/28/2018 10:38 AM **Analyst Name** ISP Tox  
**Report Time** 6/28/2018 10:39 AM **Reporter Name** ISP Tox  
**Last Calib Update** 6/28/2018 10:38 AM **Batch State** Processed

### Analysis Info

**Acq Time** 2018-06-27 15:18 **Data File** Cal 1 - 3ngr.d  
**Sample Type** Calibration **Sample Name** Cal 1 - 3ng  
**Dilution** 1 **Acq Method** AM 27 Quant THC 7-2017.m  
**Position** P1-A1 **Sample Info**  
**Inj Vol** -1 **Comment** AM 27 Cannabinoid Confirmation

### Sample Chromatogram



### Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	1.755	11347	333268	0.0340	3.2882
THC-COOH	THC-COOH-d9	1.805	8105	110609	0.0733	2.9785
THC	THC-d3	4.071	4099	96776	0.0424	3.1909

# ISP FORENSICS - Cd'A Instrument # 62340


## Cannabinoids Analysis Report

**Batch Data Path** D:\2018 Data\06272018 cann quant\QuantResults\cann quant.batch.bin  
**Analysis Time** 6/28/2018 10:38 AM **Analyst Name** ISP Tox  
**Report Time** 6/28/2018 10:39 AM **Reporter Name** ISP Tox  
**Last Calib Update** 6/28/2018 10:38 AM **Batch State** Processed

### Analysis Info

**Acq Time** 2018-06-27 15:30 **Data File** Cal 2 - 5ngr.d  
**Sample Type** Calibration **Sample Name** Cal 2 - 5ng  
**Dilution** 1 **Acq Method** AM 27 Quant THC 7-2017.m  
**Position** P1-B1 **Sample Info**  
**Inj Vol** -1 **Comment** AM 27 Cannabinoid Confirmation

### Sample Chromatogram



### Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	1.735	18211	334388	0.0545	5.0467
THC-COOH	THC-COOH-d9	1.805	13255	111528	0.1188	4.9383
THC	THC-d3	4.051	6382	90913	0.0702	5.1428

# ISP FORENSICS - Cd'A Instrument # 62340


## Cannabinoids Analysis Report

Batch Data Path D:\2018 Data\06272018 cann quant\QuantResults\cann quant.batch.bin  
Analysis Time 6/28/2018 10:38 AM Analyst Name ISP Tox  
Report Time 6/28/2018 10:39 AM Reporter Name ISP Tox  
Last Calib Update 6/28/2018 10:38 AM Batch State Processed

### Analysis Info

Acq Time 2018-06-27 15:42 Data File Cal 3 - 10ngr.d  
Sample Type Calibration Sample Name Cal 3 - 10ng  
Dilution 1 Acq Method AM 27 Quant THC 7-2017.m  
Position P1-C1 Sample Info  
Inj Vol -1 Comment AM 27 Cannabinoid Confirmation

### Sample Chromatogram



### Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	1.755	29826	279685	0.1066	9.5422
THC-COOH	THC-COOH-d9	1.805	25234	94437	0.2672	11.3174
THC	THC-d3	4.071	10302	76460	0.1347	9.6676

# ISP FORENSICS - Cd'A Instrument # 62340


## Cannabinoids Analysis Report

**Batch Data Path** D:\2018 Data\06272018 cann quant\QuantResults\cann quant.batch.bin  
**Analysis Time** 6/28/2018 10:38 AM **Analyst Name** ISP Tox  
**Report Time** 6/28/2018 10:39 AM **Reporter Name** ISP Tox  
**Last Calib Update** 6/28/2018 10:38 AM **Batch State** Processed

### Analysis Info

**Acq Time** 2018-06-27 15:54 **Data File** Cal 4 - 25ngr.d  
**Sample Type** Calibration **Sample Name** Cal 4 - 25ng  
**Dilution** 1 **Acq Method** AM 27 Quant THC 7-2017.m  
**Position** P1-D1 **Sample Info**  
**Inj Vol** -1 **Comment** AM 27 Cannabinoid Confirmation

### Sample Chromatogram



### Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	1.755	72982	266390	0.2740	23.9577
THC-COOH	THC-COOH-d9	1.805	48918	89976	0.5437	23.2058
THC	THC-d3	4.051	25268	74801	0.3378	23.9045

# ISP FORENSICS - Cd'A Instrument # 62340


## Cannabinoids Analysis Report

**Batch Data Path** D:\2018 Data\06272018 cann quant\QuantResults\cann quant.batch.bin  
**Analysis Time** 6/28/2018 10:38 AM **Analyst Name** ISP Tox  
**Report Time** 6/28/2018 10:39 AM **Reporter Name** ISP Tox  
**Last Calib Update** 6/28/2018 10:38 AM **Batch State** Processed

### Analysis Info

**Acq Time** 2018-06-27 16:06 **Data File** Cal 5 - 50ngr.d  
**Sample Type** Calibration **Sample Name** Cal 5 - 50ng  
**Dilution** 1 **Acq Method** AM 27 Quant THC 7-2017.m  
**Position** P1-E1 **Sample Info**  
**Inj Vol** -1 **Comment** AM 27 Cannabinoid Confirmation

### Sample Chromatogram



### Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	1.755	140632	259519	0.5419	47.0407
THC-COOH	THC-COOH-d9	1.805	93101	87860	1.0597	45.3922
THC	THC-d3	4.071	48173	71600	0.6728	47.3921

# ISP FORENSICS - Cd'A Instrument # 62340


## Cannabinoids Analysis Report

**Batch Data Path** D:\2018 Data\06272018 cann quant\QuantResults\cann quant.batch.bin  
**Analysis Time** 6/28/2018 10:38 AM **Analyst Name** ISP Tox  
**Report Time** 6/28/2018 10:39 AM **Reporter Name** ISP Tox  
**Last Calib Update** 6/28/2018 10:38 AM **Batch State** Processed

### Analysis Info

**Acq Time** 2018-06-27 16:17 **Data File** Cal 6 - 100ngr.d  
**Sample Type** Calibration **Sample Name** Cal 6 - 100ng  
**Dilution** 1 **Acq Method** AM 27 Quant THC 7-2017.m  
**Position** P1-F1 **Sample Info**  
**Inj Vol** -1 **Comment** AM 27 Cannabinoid Confirmation

### Sample Chromatogram



### Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	1.755	310323	257642	1.2045	104.1244
THC-COOH	THC-COOH-d9	1.805	203691	83146	2.4498	105.1678
THC	THC-d3	4.031	107059	72534	1.4760	103.7021

# ISP FORENSICS - Cd'A Instrument # 62340


## Cannabinoids Analysis Report

**Batch Data Path** D:\2018 Data\06272018 cann quant\QuantResults\cann quant.batch.bin  
**Analysis Time** 6/28/2018 10:38 AM **Analyst Name** ISP Tox  
**Report Time** 6/28/2018 10:39 AM **Reporter Name** ISP Tox  
**Last Calib Update** 6/28/2018 10:38 AM **Batch State** Processed

### Analysis Info

<b>Acq Time</b>	2018-06-27 16:29	<b>Data File</b>	Cal 7 - 250ngr.d
<b>Sample Type</b>	Calibration	<b>Sample Name</b>	Cal 7 - 250ng
<b>Dilution</b>	1	<b>Acq Method</b>	AM 27 Quant THC 7-2017.m
<b>Position</b>	P1-G1	<b>Sample Info</b>	
<b>Inj Vol</b>	-1	<b>Comment</b>	AM 27 Cannabinoid Confirmation

### Sample Chromatogram



### Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-d3	1.755	746813	339420	2.2003	189.9149
THC-COOH	THC-COOH-d9	1.805	487631	98268	4.9622	213.2009
THC	THC-d3	4.071	264949	92146	2.8753	201.8096