

**REVIEWED**

By Sarah Pickle at 8:57 am, Dec 19, 2019

12/18/2019





TS

**Worklist: 3889**

<u>LAB CASE</u>	<u>ITEM</u>	<u>ITEM TYPE</u>	<u>DESCRIPTION</u>	
M2019-3507	3	BCK	AM 25/AM 26 Blood MultiDrug/THC Screen by LC-QQQ	
M2019-4718	2	BCK	AM 25/AM 26 Blood MultiDrug/THC Screen by LC-QQQ	
M2019-4835	2	BCK	AM 25/AM 26 Blood MultiDrug/THC Screen by LC-QQQ	
M2019-4927	1	BCK	AM 25/AM 26 Blood MultiDrug/THC Screen by LC-QQQ	
M2019-4938	1	BCK	AM 25/AM 26 Blood MultiDrug/THC Screen by LC-QQQ	
M2019-4938	2	BCK	AM 25/AM 26 Blood MultiDrug/THC Screen by LC-QQQ	
M2019-5117	1	BCK	AM 25/AM 26 Blood MultiDrug/THC Screen by LC-QQQ	
M2019-5141	2	BCK	AM 25/AM 26 Blood MultiDrug/THC Screen by LC-QQQ	
M2019-5163	2	BCK	AM 25/AM 26 Blood MultiDrug/THC Screen by LC-QQQ	
M2019-5191	1	BCK	AM 25/AM 26 Blood MultiDrug/THC Screen by LC-QQQ	
M2019-5282	2	BCK	AM 25/AM 26 Blood MultiDrug/THC Screen by LC-QQQ	
P2019-3510	1	BCK	AM 25/AM 26 Blood MultiDrug/THC Screen by LC-QQQ	
P2019-3521	1	BCK	AM 25/AM 26 Blood MultiDrug/THC Screen by LC-QQQ	
P2019-3530	1	BCK	AM 25/AM 26 Blood MultiDrug/THC Screen by LC-QQQ	
P2019-3532	1	BCK	AM 25/AM 26 Blood MultiDrug/THC Screen by LC-QQQ	
P2019-3534	1	BCK	AM 25/AM 26 Blood MultiDrug/THC Screen by LC-QQQ	
P2019-3549	1	BCK	AM 25/AM 26 Blood MultiDrug/THC Screen by LC-QQQ	
P2019-3592	1	BCK	AM 25/AM 26 Blood MultiDrug/THC Screen by LC-QQQ	
P2019-3593	1	BCK	AM 25/AM 26 Blood MultiDrug/THC Screen by LC-QQQ	
P2019-3604	1	BCK	AM 25/AM 26 Blood MultiDrug/THC Screen by LC-QQQ	
P2019-3631	1	BCK	AM 25/AM 26 Blood MultiDrug/THC Screen by LC-QQQ	

Worklist: 3889

15

<u>LAB CASE</u>	<u>ITEM</u>	<u>ITEM TYPE</u>	<u>DESCRIPTION</u>	
P2019-3644	1	BCK	AM 25/AM 26 Blood MultiDrug/THC Screen by LC-QQQ	
P2019-3645	1	BCK	AM 25/AM 26 Blood MultiDrug/THC Screen by LC-QQQ	
P2019-3654	1	BCK	AM 25/AM 26 Blood MultiDrug/THC Screen by LC-QQQ	
P2019-3677	1	BCK	AM 25/AM 26 Blood MultiDrug/THC Screen by LC-QQQ	

# AM# 25: Multi-Drug Screen in Blood by LC-MS/MS

Extraction Date: 12/17/19  
Plate Item #: IDP-107 Plate Lot#: 190725

Analyst: Tamara Salazar  
Plate Expiration: 01/25/20

**Mobile phase A:** 10mM Amm Form  
0.5M Ammonium Hydroxide

**Mobile phase B:** 0.1% Formic Acid in MeOH  
Ethyl Acetate LC Methanol

**Blank Blood Lot:** Hemostat 445283-3  
**LCMS-QQQ ID:** 069901

**Column:** Phenomenex Phenyl Hexyl (4.6x50mm, 2.6um)

## Pre-Analytic:

- 1. Check levels of mobile phases and needle wash refill as needed. Ensure waste is not full.
- 2. Ensure correct column is installed and begin mobile phase flow allow to equilibrate ~ 30 minutes.

## Analytic:

- 1. Remove standards, plate, controls, and samples from cold storage. Allow to reach room temperature.
- 2. Pipette **250µL blood (calibrated pipette) Pipette ID: 16** in wells of analytical (standards) plate.
- 3. Place on shaking incubator at ambient temp., 900rpm for 15 minutes. *Shaker ID: 067105*
- 4. Pipette **250µL 0.5M ammonium hydroxide** in wells of analytical plate.
- 5. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 6. Transfer **300µL of blood+base** mixture to corresponding wells of SLE+ plate.
- 7. Apply positive pressure for approx. 10-15 seconds (or until no liquid remains on top of sorbent).  
*(Load at 85-100 PSI- Selector to the right) Manifold ID: 067104*
- 8. Wait 5 minutes.
- 9. Add **900uL ethyl acetate**.
- 10. Wait 5 minutes.
- 11. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left)*.
- 12. Add **900uL ethyl acetate**.
- 13. Wait 5 minutes.
- 14. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left)*.
- 15. Remove plate containing eluate. Place on SPE Dry and evaporate to dryness at approx. 35°C.  
*SPE Dry ID: 067103*
- 16. Reconstitute in **100µL 100% LC MeOH** and heat seal plate with foil. Place in autosampler and run worklist.

## Post-Analytic

- 1. Open quantitation software and create a new quantitation batch.  
Batch Name: AM 25 MDS  
Worklist path: *D:\MassHunter\Data\2019\AM 25-26\121719 wklst 3889 TS*
- 2. Make necessary changes to integration limits
- 3. Evaluate samples, S/N of primary transition >5 and S/N of secondary transition >3 or evaluation of peak symmetry and resolution. Within +/- 2% or 0.1 min RT of administrative control. Calculated concentration 5 or greater, or 2-5 for discretionary range.
- 4. Did all QCs pass for each analyte? Y / N \_\_\_\_\_
- 5. Central File Packet to include: LIMS Worklist, Method Checklist, Calibration and Control Reports

COMMENTS:



# Idaho State Police Forensic Services

## AM #25 Blood Multi-Drug Screen by LCMS-QQQ

### Methanol External Control Solution (Lot: 042719)

100 ul of 1mg/mL stock was added to each drug to 9600 ul of LC MeOH.

<i>Component</i>	<i>Source</i>	<i>Source Lot Number</i>	<i>Expiration Date</i>
Methanol (LCMS)	Fisher	184782	
Morphine	Cerilliant	FE08141515	November 2020
Metoprolol	Cerilliant	FN06091510	July 2020
Flunitrazepam	Cerilliant	FE08051602	August 2021
Trazodone	Cerilliant	FN12151403	January 2020
Prepared:	04/27/19		
Prepared By:	Tamara Salazar		
Expires:	01/31/2020		

### Blood External Control Solution (Lot: WS042719)

100 ul of methanol external control solution was added to 9900 ul of blood.

Approximately 50ng/mL of each compound.

<i>Component</i>	<i>Source</i>	<i>Source Lot Number</i>
Negative Blood	Hemostat	445283-1
Methanol External Control Solution		042719
Prepared:	04/27/19	
Prepared by:	Tamara Salazar	
Expires:	01/31/2020	




# AM #25 Multi-Drug Screen Results



**Batch results** D:\MassHunter\Data\2019\AM 25-26\121719 wk1st 3889 TS\QuantResults\AM 25 MDS.batch.bin  
**Calibration Last Update** 12/17/2019 3:08:20 PM

<b>Instrument</b>	Falco	<b>Data File</b>	External Control- Blood.d
<b>Type</b>	Sample	<b>Sample</b>	External Control- Blood
<b>Acq. Method</b>	am 25 all.m	<b>Operator</b>	Tamara Salazar
<b>Sample Position</b>	P1-E12	<b>Comment</b>	
<b>Injection Volume</b>	5		
<b>Acq. Date-Time</b>	12/17/2019 11:26:17 AM		
<b>Sample Info.</b>			

## Sample Chromatogram



Name	RT	Resp.	S/N	S/N	ISTD Resp.	Calc. Conc.
Flunitrazepam	4.562	14234263	7434.20	4371.13	501187	61.4331
Metoprolol	3.353	4837523	1962.33	20814.22	17368364	55.0137
Morphine	2.174	1195967	76166.13	20759.32	104777	67.4978
Trazodone	4.021	28595131	2053186.95	37042.61	20964037	55.2117

# AM #25 Multi-Drug Screen Results




**Batch results** D:\MassHunter\Data\2019\AM 25-26\121719 wk1st 3889 TS\QuantResults\AM 25 MDS.batch.bin  
**Calibration Last Update** 12/17/2019 3:08:20 PM

<b>Instrument</b>	Falco	<b>Data File</b>	Negative Blood.d
<b>Type</b>	Sample	<b>Sample</b>	Negative Blood
<b>Acq. Method</b>	am 25 all.m	<b>Operator</b>	Tamara Salazar
<b>Sample Position</b>	P1-F12	<b>Comment</b>	
<b>Injection Volume</b>	5		
<b>Acq. Date-Time</b>	12/17/2019 11:17:59 AM		
<b>Sample Info.</b>			

## Sample Chromatogram

+ TIC MRM (\*\* -> \*\*) Negative Blood.d (Negative Blood)




# AM #25 Multi-Drug Screen Results



**Batch results** D:\MassHunter\Data\2019\AM 25-26\121719 wkst 3889 TS\QuantResults\AM 25 MDS.batch.bin  
**Calibration Last Update** 12/17/2019 3:08:20 PM

<b>Instrument</b>	Falco	<b>Data File</b>	Cal.d
<b>Type</b>	Cal	<b>Sample</b>	Cal
<b>Acq. Method</b>	am 25 all.m	<b>Operator</b>	Tamara Salazar
<b>Sample Position</b>	P1-H12	<b>Comment</b>	
<b>Injection Volume</b>	5		
<b>Acq. Date-Time</b>	12/17/2019 11:09:30 AM		

## Sample Chromatogram



Name	RT	Resp.	S/N	S/N	ISTD Resp.	Calc. Conc.
6-MAM	2.783	43619	74.63	139.41	1096285	10.0000
7-aminoclonazepam	3.551	1522053	45737.03	222.50	6472960	10.0000
7-aminoflunitrazepam	3.765	2872755	1440.15	105.45	19423504	10.0000
Acetyl Fentanyl	3.670	454156	612.37	126.67	32681787	10.0000
Acetyl Norfentanyl	2.808	285213	1118.38	87.00	13503180	10.0000
a-hydroxyalprazolam	4.499	217562	503.41	∞	1141010	10.0000
alpha-hydroxymidazolam	4.513	1656731	294.69	1747.92	10394497	10.0000
alpha-PVP	3.420	5454058	963.24	856.30	24429965	10.0000
Alprazolam	4.609	1914617	∞	187.00	5551707	10.0000
Amitriptyline	4.322	4685510	74.78	345.95	10613029	10.0000
Amphetamine	2.782	2761625	∞	510.16	8011027	10.0000
Benzoylcegonine	3.351	1090610	266.10	339.97	4932660	10.0000
Buprenorphine	4.004	506544	152.74	65220.92	2161465	10.0000
Bupropion	3.619	5168121	20987.85	∞	16166243	10.0000
Carbamazepine	4.218	6211225	∞	∞	36983477	10.0000
Carisoprodol	4.201	1241060	133930.38	103.28	6806442	10.0000
Chlordiazepoxide	4.611	306675	∞	∞	11710003	10.0000
Chlorpheniramine	3.843	27299	39.86	∞	52781037	10.0000
Citalopram	3.977	2549178	1067.23	2529.79	12185708	10.0000
Clonazepam	4.439	1199253	16913.75	375.78	2335475	10.0000
Cocaine	3.458	6387451	859114.42	1882.68	30285944	10.0000
Codeine	2.681	450928	90.92	423.13	1889168	10.0000
Cyclobenzaprine	4.245	3534109	1343.14	101.75	12561905	10.0000
Desipramine	4.277	5266389	∞	4429.14	29935915	10.0000
Dextromethorphan	3.967	1837456	∞	232.87	8960430	10.0000
Dextrorphan	3.293	2873382	783.48	717.43	18684420	10.0000
Diazepam	4.841	1238843	347.35	∞	6416647	10.0000
Dihydrocodeine	2.650	889398	78.09	71.98	4119908	10.0000
Diphenhydramine	3.923	10840285	1693899.77	3197.61	52781037	10.0000
Doxepin	4.043	2545926	∞	∞	16546558	10.0000
Doxylamine	3.522	11989154	∞	231420.25	43586444	10.0000
EDDP	3.981	3943542	2604.31	867.20	27097346	10.0000
Estazolam	4.519	5379694	841.57	556.54	15383606	10.0000
Etizolam	4.619	338234	240.64	∞	15383606	10.0000

Cal



# AM #25 Multi-Drug Screen Results



Name	RT	Resp.	S/N	S/N	ISTD Resp.	Calc. Conc.
Fentanyl	3.899	374889	178.63	257.02	23911519	10.0000
Flunitrazepam	4.562	2429872	2732.39	264.77	525594	10.0000
Fluoxetine	4.240	3481303	3123.65	181.29	15636167	10.0000
Flurazepam	4.020	2356680	1106.00	221.94	525594	10.0000
Hydrocodone	2.864	1452360	39.44	14.10	9418169	10.0000
Hydromorphone	2.353	1191577	172.02	∞	4610723	10.0000
Imipramine	4.274	6174203	23840.74	165.34	23377225	10.0000
Ketamine	3.250	4584052	∞	∞	21945558	10.0000
Lamotrigine	3.417	371649	586.23	983.37	14908852	10.0000
Levamisole	2.824	3508491	416.94	450.74	30285944	10.0000
Lorazepam	4.423	581236	109.90	∞	2335475	10.0000
Maprotiline	4.322	4705003	72.69	311.15	10613029	10.0000
MDA	2.917	2122525	∞	∞	9740654	10.0000
MDEA	3.145	5502003	4176.41	4578.12	25091729	10.0000
MDMA	2.992	5967368	∞	2096.95	4276920	10.0000
Meperidine	3.479	3037170	447.13	2035.47	14908852	10.0000
Meprobamate	3.636	571118	2074.67	110.09	2371331	10.0000
Methadone	4.286	6608098	34519.84	342.53	29069586	10.0000
Methamphetamine	2.887	4078809	784.80	1301.46	21038329	10.0000
Methocarbamol	3.541	533996	330.62	278.71	14908852	10.0000
Methylphenidate	3.404	10817408	645.22	∞	44053979	10.0000
Metoprolol	3.353	754810	560.14	1395.85	14908852	10.0000
Midazolam	4.437	720337	336.49	351806.42	8354661	10.0000
Mirtazapine	3.584	3154389	463.92	242885.26	14908852	10.0000
Mitragynine	4.035	367020	20414.21	57874.40	16546558	10.0000
Morphine	2.174	200358	3465.98	212519.18	118479	10.0000
Norbuprenorphine	3.729	88811	60788.21	28.77	407047	10.0000
Nordiazepam	4.691	1831766	795.65	670.79	5739621	10.0000
Norfentanyl	3.235	7354866	∞	1025.43	31579752	10.0000
Norhydrocodone	2.850	25972	∞	35.04	1019371	10.0000
Normeperidine	3.512	2545841	1605.27	∞	9051289	10.0000
Noroxycodone	2.818	908348	82.59	74.80	2845514	10.0000
Nortriptyline	4.325	2424187	302759.91	210.92	5572577	10.0000
O-desmethyl-tramadol	2.822	8545384	7779.58	37.84	40274443	10.0000
Olanzapine	3.318	554039	124.27	180.63	803524	10.0000
Oxazepam	4.504	2574984	420.18	128.23	16550888	10.0000
Oxycodone	2.815	2265936	∞	204.83	9721213	10.0000
Oxymorphone	2.258	1175764	124.51	576.02	4131675	10.0000
Paroxetine	4.253	456030	165.91	26.66	12643951	10.0000
Phenazepam	4.634	1325005	42.46	248.93	6503260	10.0000
Phencyclidine	3.816	6064147	7146.16	1951.26	27537115	10.0000
Phentermine	3.055	1963109	252.85	20.36	22544043	10.0000
Phenytoin	4.109	123574	8327.13	42.82	803524	10.0000
Promethazine	4.197	10509983	10654.91	364.57	37858222	10.0000
Pseudoephedrine	2.613	50468480	5735.00	12846.40	133712227	10.0000
Quetiapine	4.144	3714172	801180.17	806.31	5662239	10.0000
Sertraline	4.456	2499844	570.59	438.20	12643951	10.0000
Sufentanil	4.159	302007	11186.03	195.69	21059251	10.0000
Tapentadol	3.343	4824322	1867.98	1425.86	24950236	10.0000
Temazepam	4.656	3812457	∞	176.59	18400950	10.0000
Tramadol	3.324	10899406	6886.20	280.82	43442289	10.0000
Trazodone	4.021	4661026	∞	∞	18866680	10.0000
Venlafaxine	3.704	7563925	1165.93	244.77	38738839	10.0000
Zaleplon	4.349	3171443	242546.76	3123.77	7087027	10.0000
Zolpidem	3.732	8768798	534039.14	∞	40816046	10.0000
Zopiclone	3.653	384494	381220.56	274.06	2025612	10.0000

# AM# 26: THC and Metabolites Screen in Blood by LC-MS/MS

5  
B

Extraction Date: 12/17/19

Analyst: Tamara Salazar

Plate lot# IDP-108, 190716

Plate Expiration: 01/16/20

**Mobile phase A:** 10mM Ammonium Formate  
0.1% Formic Acid in Water

**Mobile phase B:** 0.1% Formic acid in MeOH  
Hexane

**Blank Blood Lot:** 445283-3

**Column:** Phenomenex Phenyl Hexyl (4.6x50mm: 2.6 um)

**LCMS-QQQ ID:** 069901

## Pre-Analytic:

- 1. Check levels of mobile phases and needle wash refill as needed. Ensure waste is not full.
- 2. Ensure correct column is installed and begin mobile phase flow allow to equilibrate ~ 30 minutes.
- 3. Create worklist:

## Analytic:

- 1. Remove standards, plate, controls, and samples from cold storage. Allow to reach room temperature.
- 2. Pipette **1000 µL blood (calibrated pipette) Pipette ID: 16** in wells of analytical (standards) plate.
- 3. Place on shaking incubator at ambient temp., 900rpm for 15 minutes. *Shaker ID: 067105*
- 4. Pipette **500 µL 0.1% formic acid** for blood in wells of analytical plate.
- 5. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 6. Transfer **800 µL of blood+base** mixture to corresponding wells of SLE+ plate.
- 7. Apply positive pressure for approx. 4 seconds (or until no liquid remains on top of sorbent).  
*(Load at 85-100 PSI- Selector to the right) Manifold ID: 067104*
- 8. Wait 5 minutes.
- 9. Add **2.25 mL MTBE** (add in 3 increments of 750uL).
- 10. Wait 5 minutes.
- 11. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 12. Add **2.25 mL hexane** (add in 3 increments of 750uL).
- 13. Wait 5 minutes.
- 14. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 15. Remove plate containing eluate. Place on SPE Dry and evaporate to dryness at approx. 35°C.  
*SPE Dry ID: 067103*
- 16. Reconstitute in **100µL 100% MeOH** and heat seal plate with foil. Place in autosampler and run worklist.

## Post-Analytic

- 1. Create batch and process data.  
Worklist path: *D:\MassHunter\Data\2019\AM 25-26\121719 wklst 3889 TS*  
Batch Name: *AM 26 THCS*
- 2. Calculated sample concentration of 3 ng/mL or greater for THC and THC-OH, a calculated sample concentration of 10 ng/mL or greater for Carboxy-THC.
- 3. Retention time within +/- 2% or +/- .100 min whichever is greater of the average retention time of the calibrators.
- 4. Did all QCs pass for each analyte? Y / N
- 5. Central File Packet to include: LIMS Worklist, Method Checklist, Calibration and Control Reports

COMMENTS:


# AM #26 Cannabinoids Screen Results



**Batch results** D:\MassHunter\Data\2019\AM 25-26\121719 wk1st 3889 TS\QuantResults\AM 26 THCS.batch.bin  
**Calibration Last Update** 12/18/2019 1:20:17 PM

<b>Instrument</b>	Falco	<b>Data File</b>	MJ QC Control.d
<b>Type</b>	Sample	<b>Sample</b>	MJ QC Control
<b>Acq. Method</b>	am 26 test.m	<b>Operator</b>	Tamara Salazar
<b>Sample Position</b>	P3-H1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	12/18/2019 9:56:12 AM		
<b>Sample Info.</b>			

## Sample Chromatogram



Name	RT	Resp.	ISTD Resp.	Final Conc.
THC	2.839	1396	50051	4.1220 ng/ml
THC-COOH	2.565	17679	86477	15.7137 ng/ml
THC-OH	2.471	18700	410745	4.9990 ng/ml

TS




# AM #26 Cannabinoids Screen Results

**Batch results** D:\MassHunter\Data\2019\AM 25-26\121719 wkst 3889 TS\QuantResults\AM 26 THCS.batch.bin  
**Calibration Last Update** 12/18/2019 1:20:17 PM

<b>Instrument</b>	Falco	<b>Data File</b>	MJ Negative Blood.d
<b>Type</b>	Sample	<b>Sample</b>	MJ Negative Blood
<b>Acq. Method</b>	am 26 test.m	<b>Operator</b>	Tamara Salazar
<b>Sample Position</b>	P3-A2	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	12/18/2019 10:09:15 AM		
<b>Sample Info.</b>			

## Sample Chromatogram






TS



# AM #26 Cannabinoids Screen Calibration Curve Report

**Batch results** D:\MassHunter\Data\2019\AM 25-26\121719 wklst 3889 TS\QuantResults\AM 26 THCS.batch.bin  
**Last Cal. Update** 12/18/2019 1:20 PM  
**Analyst Name** ISP\Datastor  
**Analyte** THC **Internal Standard** THC-d3




Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
MJ Cal 1	1	✓	1.0	1.2	116.3
MJ Cal 2	2	✓	3.0	2.9	95.9
MJ Cal 3	3	✓	5.0	4.9	99.0
MJ Cal 4	4	✓	10.0	10.5	105.2
MJ Cal 5	5	✓	25.0	21.6	86.5
MJ Cal 6	6	✓	50.0	44.2	88.3
MJ Cal 7	7	✓	100.0	108.7	108.7



# AM #26 Cannabinoids Screen Calibration Curve Report

**Batch results** D:\MassHunter\Data\2019\AM 25-26\121719 wklst 3889 TS\QuantResults\AM 26 THCS.batch.bin  
**Last Cal. Update** 12/18/2019 1:20 PM  
**Analyst Name** ISP\Datastor  
**Analyte** THC-COOH **Internal Standard** THC-COOH-d9




Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
MJ Cal 1	1	✓	5.0	5.2	103.6
MJ Cal 2	2	✓	10.0	9.3	92.8
MJ Cal 3	3	✓	20.0	19.9	99.5
MJ Cal 4	4	✓	50.0	53.9	107.8
MJ Cal 5	5	✓	75.0	65.4	87.2
MJ Cal 6	6	✓	100.0	99.8	99.8
MJ Cal 7	7	✓	250.0	273.1	109.2

TS



# AM #26 Cannabinoids Screen Calibration Curve Report

**Batch results** D:\MassHunter\Data\2019\AM 25-26\121719 wk1st 3889 TS\QuantResults\AM 26 THCS.batch.bin  
**Last Cal. Update** 12/18/2019 1:20 PM  
**Analyst Name** ISP\Datastor  
**Analyte** THC-OH **Internal Standard** THC-OH-d3



Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
MJ Cal 1	1	✓	1.0	0.9	93.0
MJ Cal 2	2	✓	3.0	3.2	107.4
MJ Cal 3	3	✓	5.0	5.5	109.9
MJ Cal 4	4	✓	10.0	9.9	98.8
MJ Cal 5	5	✓	25.0	21.0	83.8
MJ Cal 6	6	✓	50.0	53.5	107.0
MJ Cal 7	7	×	100.0	143.7	143.7


# AM #26 Cannabinoids Screen Results



**Batch results** D:\MassHunter\Data\2019\AM 25-26\121719 wkst 3889 TS\QuantResults\AM 26 THCS.batch.bin  
**Calibration Last Update** 12/18/2019 1:20:17 PM

<b>Instrument</b>	Falco	<b>Data File</b>	MJ Cal 1.d
<b>Type</b>	Cal	<b>Sample</b>	MJ Cal 1
<b>Acq. Method</b>	am 26 test.m	<b>Operator</b>	Tamara Salazar
<b>Sample Position</b>	P3-A1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	12/18/2019 9:10:30 AM		
<b>Sample Info.</b>			

## Sample Chromatogram



Name	RT	Resp.	ISTD Resp.	Final Conc.	
THC	2.839	914	131488	1.1634 ng/ml	Low
THC-COOH	2.565	11918	194652	5.1821 ng/ml	
THC-OH	2.471	11054	961551	0.9298 ng/ml	Low

# AM #26 Cannabinoids Screen Results




15

**Batch results** D:\MassHunter\Data\2019\AM 25-26\121719 wkst 3889 TS\QuantResults\AM 26 THCS.batch.bin  
**Calibration Last Update** 12/18/2019 1:20:17 PM

<b>Instrument</b>	Falco	<b>Data File</b>	MJ Cal 2.d
<b>Type</b>	Cal	<b>Sample</b>	MJ Cal 2
<b>Acq. Method</b>	am 26 test.m	<b>Operator</b>	Tamara Salazar
<b>Sample Position</b>	P3-B1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	12/18/2019 9:17:10 AM		
<b>Sample Info.</b>			

## Sample Chromatogram



Name	RT	Resp.	ISTD Resp.	Final Conc.	
THC	2.819	3202	167795	2.8782 ng/ml	<b>Low</b>
THC-COOH	2.565	26925	230356	9.2750 ng/ml	
THC-OH	2.471	29680	967688	3.2225 ng/ml	




# AM #26 Cannabinoids Screen Results



**Batch results** D:\MassHunter\Data\2019\AM 25-26\121719 wkst 3889 TS\QuantResults\AM 26 THCS.batch.bin  
**Calibration Last Update** 12/18/2019 1:20:17 PM

<b>Instrument</b>	Falco	<b>Data File</b>	MJ Cal 3.d
<b>Type</b>	Cal	<b>Sample</b>	MJ Cal 3
<b>Acq. Method</b>	am 26 test.m	<b>Operator</b>	Tamara Salazar
<b>Sample Position</b>	P3-C1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	12/18/2019 9:23:41 AM		
<b>Sample Info.</b>			

## Sample Chromatogram



Name	RT	Resp.	ISTD Resp.	Final Conc.
THC	2.819	4725	140069	4.9483 ng/ml
THC-COOH	2.565	42929	164235	19.9013 ng/ml
THC-OH	2.471	33738	679085	5.4958 ng/ml

15




# AM #26 Cannabinoids Screen Results

Batch results D:\MassHunter\Data\2019\AM 25-26\121719 wkst 3889 TS\QuantResults\AM 26 THCS.batch.bin  
Calibration Last Update 12/18/2019 1:20:17 PM

<b>Instrument</b>	Falco	<b>Data File</b>	MJ Cal 4.d
<b>Type</b>	Cal	<b>Sample</b>	MJ Cal 4
<b>Acq. Method</b>	am 26 test.m	<b>Operator</b>	Tamara Salazar
<b>Sample Position</b>	P3-D1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	12/18/2019 9:30:11 AM		
<b>Sample Info.</b>			

## Sample Chromatogram



Name	RT	Resp.	ISTD Resp.	Final Conc.
THC	2.819	8829	120683	10.5204 ng/ml
THC-COOH	2.565	115278	159255	53.9108 ng/ml
THC-OH	2.471	58783	680495	9.8844 ng/ml




# AM #26 Cannabinoids Screen Results



**Batch results** D:\MassHunter\Data\2019\AM 25-26\121719 wkst 3889 TS\QuantResults\AM 26 THCS.batch.bin  
**Calibration Last Update** 12/18/2019 1:20:17 PM

<b>Instrument</b>	Falco	<b>Data File</b>	MJ Cal 5.d
<b>Type</b>	Cal	<b>Sample</b>	MJ Cal 5
<b>Acq. Method</b>	am 26 test.m	<b>Operator</b>	Tamara Salazar
<b>Sample Position</b>	P3-E1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	12/18/2019 9:36:41 AM		
<b>Sample Info.</b>			

## Sample Chromatogram



Name	RT	Resp.	ISTD Resp.	Final Conc.
THC	2.839	8015	52799	21.6324 ng/ml
THC-COOH	2.565	85929	97599	65.4246 ng/ml
THC-OH	2.471	74436	415903	20.9563 ng/ml


# AM #26 Cannabinoids Screen Results



**Batch results** D:\MassHunter\Data\2019\AM 25-26\121719 wkst 3889 TS\QuantResults\AM 26 THCS.batch.bin  
**Calibration Last Update** 12/18/2019 1:20:17 PM

<b>Instrument</b>	Falco	<b>Data File</b>	MJ Cal 6.d
<b>Type</b>	Cal	<b>Sample</b>	MJ Cal 6
<b>Acq. Method</b>	am 26 test.m	<b>Operator</b>	Tamara Salazar
<b>Sample Position</b>	P3-F1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	12/18/2019 9:43:12 AM		
<b>Sample Info.</b>			

## Sample Chromatogram



Name	RT	Resp.	ISTD Resp.	Final Conc.
THC	2.819	49736	159820	44.1595 ng/ml
THC-COOH	2.565	330155	244952	99.7972 ng/ml
THC-OH	2.471	488273	1082099	53.5112 ng/ml


# AM #26 Cannabinoids Screen Results



**Batch results** D:\MassHunter\Data\2019\AM 25-26\121719 wkst 3889 TS\QuantResults\AM 26 THCS.batch.bin  
**Calibration Last Update** 12/18/2019 1:20:17 PM

<b>Instrument</b>	Falco	<b>Data File</b>	MJ Cal 7.d
<b>Type</b>	Cal	<b>Sample</b>	MJ Cal 7
<b>Acq. Method</b>	am 26 test.m	<b>Operator</b>	Tamara Salazar
<b>Sample Position</b>	P3-G1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	12/18/2019 9:49:42 AM		
<b>Sample Info.</b>			

## Sample Chromatogram



Name	RT	Resp.	ISTD Resp.	Final Conc.
THC	2.819	73555	95789	108.6978 ng/ml
THC-COOH	2.565	496272	133954	273.1235 ng/ml
THC-OH	2.471	704999	585026	143.6533 ng/ml