


**Worklist: 2990**

<u>LAB_CASE</u>	<u>ITEM</u>	<u>TASK_ID</u>	<u>DESCRIPTION</u>
M2019-0371	1	142963	AM 27 Blood THC Quant by LC-QQQ
M2019-0479	1	142964	AM 27 Blood THC Quant by LC-QQQ
P2019-0474	1	142965	AM 27 Blood THC Quant by LC-QQQ
P2019-0501	1	142966	AM 27 Blood THC Quant by LC-QQQ
P2019-0512	1	142967	AM 27 Blood THC Quant by LC-QQQ
P2019-0513	1	142968	AM 27 Blood THC Quant by LC-QQQ
P2019-0602	1	142969	AM 27 Blood THC Quant by LC-QQQ
P2019-0603	1	142970	AM 27 Blood THC Quant by LC-QQQ



# AM# 27: Quantitation of THC and Metabolites in Blood by LC-MS/MS

D

Extraction Date: 2/27/19  
Plate lot#: 0539904

Analyst: Sarah Pickle  
Plate Expiration: 09/10/19

**Mobile phase A:** 0.1% Formic Acid in LCMS Water  
MTBE  
**Blank Blood Lot:** 445283-1  
**LCMS-QQQ ID:** 59740

**Mobile phase B:** 0.1% Formic acid in Acetonitrile  
LCMS Methanol  
Hexane  
**Column:** UCT Selectra DA 100 x 2.1mm 3um

## Pre-Analytic:

- 1. Check levels of mobile phases and needle wash refill as needed. Ensure waste is not full.
- 2. Ensure correct column is installed and begin mobile phase flow allow to equilibrate ~ 30 minutes.
- 3. Create worklist: Data Path: \_\_\_\_\_

## Analytic:

- 1. Remove standards, plate, controls, and samples from cold storage. Allow to reach room temperature.
- 2. Pipette **1000µL blood (calibrated pipette) Pipette ID: #27** in wells of analytical (standards) plate.
- 3. Place on shaking incubator at ambient temp., 900rpm for 15 minutes. *Shaker ID: 067105*
- 4. Pipette **500µL 0.1% formic acid in water** in wells of analytical plate.
- 5. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 6. Transfer **800µL of blood+acid** mixture to corresponding wells of SLE+ plate.
- 7. Apply positive pressure for approx. 10-15 seconds (or until no liquid remains on top of sorbent).  
*(Load at 85-100 PSI- Selector to the right) Manifold ID: 067104*
- 8. Wait 5 minutes.
- 9. Add **2.25mL MTBE. (Add in 3 increments of 750uL)**
- 10. Wait 5 minutes.
- 11. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 12. Add **2.25mL Hexane. (Add in 3 increments of 750uL)**
- 13. Wait 5 minutes.
- 14. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 15. Remove plate containing eluate. Place on SPE Dry and evaporate to dryness at approx. 35°C.  
*SPE Dry ID: 067103*
- 16. Reconstitute in **100µL 100% MeOH** and heat seal plate with foil. Place in autosampler and run worklist.

## Post-Analytic

- 1. Create batch and process data.  
Worklist path: 022719 THCQ TS SP Batch Name: THCQ SP
- 2. Make any necessary integration changes, Curve weighting of Linear 1/x with  $r^2$  values  $\geq 0.98$  for each analyte
- 3. RT +/- 3% or 0.100 min, whichever is greater, +/- 20% Accuracy for greater than (+/- 30% for 10ng/ml or less).  
Ion ratios must be within +/- 20% of the averaged calibrators
- 4. Case sample response for THC and OH-THC 3ng/mL (quantitative), Carboxy-THC: 10ng/mL (qualitative only) will be reported. Samples with a THC or OH-THC response over 50 ng/mL will be reported out as greater than 50 ng/mL.
- 5. Did all QCs pass for each analyte? Y / N
- 6. Enter QCs into control charting.
- 7. Central File Packet to include: LIMS Worklist, Method Checklist, Calibration and Control Reports

COMMENTS: *Curve range limited: Carboxy-THC 10-100*



# Idaho State Police Forensic Services

D

## AM #27 Quantitative Analysis of THC and Metabolites in Blood by LCMS-QQQ

### Methanol External Control Solution (Lot: WS020419)

10 ul of 1mg/mL THC, 100 ul of 100 ug/mL THC-OH, C-THC in 9790 ul MeOH  
Approximate concentration 1ug/mL.

<i>Component</i>	<i>Source</i>	<i>Source Lot Number</i>	<i>Expiration Date</i>
Methanol (LCMS)	Fisher	184782	
THC	Cerilliant	FE04231406	04/30/2019
C-THC	Cerilliant	FE07171501	09/31/2020
THC-OH	Cerilliant	FE01121503	01/31/2020
Prepared:	02/04/19		
Prepared By:	Tamara Salazar		
Expires:	04/30/2019		

### Blood External Control Solution (Lot: 020419)

100 ul of methanol external control solution was added to 9900 ul of blood.  
Approximately 10ng/mL of each compound.

<i>Component</i>	<i>Source</i>	<i>Source Lot Number</i>
Negative Blood	Hemostat	445283-1
Methanol External Control Solution	-	WS020419
Prepared:	02/04/19	
Prepared by:	Tamara Salazar	
Expires:	04/30/2019	

# ISP FORENSICS - Pocatello Instrument # 59740

## Cannabinoids Analysis Report


P

**Batch Data Path** C:\MassHunter\Data\2019\AM 27\022719 THCQ TS SP\QuantResults\THCQ SP.batch.bin  
**Analysis Time** 2/28/2019 12:45 PM **Analyst Name** datastor  
**Report Time** 3/1/2019 2:15 PM **Reporter Name** datastor  
**Last Calib Update** 2/28/2019 12:45 PM **Batch State** Processed

**Analysis Info**

**Acq Time** 2019-02-27 13:11 **Data File** QC Control.d  
**Sample Type** QC **Sample Name** QC Control  
**Dilution** 1 **Acq Method** THC Quant 051517 workingmm.m  
**Position** P2-H1 **Sample Info**  
**Inj Vol** -1 **Comment**

**Sample Chromatogram**



**Results**

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-D3	2.265	12087	252283	0.0479	4.9329
THC-COOH	THC-COOH-D9	2.379	17713	84323	0.2101	9.5760
THC	THC-D3	6.039	3945	87102	0.0453	4.9408

P

# ISP FORENSICS - Pocatello Instrument # 59740

## Cannabinoids Analysis Report


---

<b>Batch Data Path</b>	C:\MassHunter\Data\2019\AM 27\022719 THCQ TS SP\QuantResults\THCQ SP.batch.bin		
<b>Analysis Time</b>	2/28/2019 12:45 PM	<b>Analyst Name</b>	datastor
<b>Report Time</b>	3/1/2019 2:15 PM	<b>Reporter Name</b>	datastor
<b>Last Calib Update</b>	2/28/2019 12:45 PM	<b>Batch State</b>	Processed

### Analysis Info

<b>Acq Time</b>	2019-02-27 13:34	<b>Data File</b>	Negative.d
<b>Sample Type</b>	Sample	<b>Sample Name</b>	Negative
<b>Dilution</b>	1	<b>Acq Method</b>	THC Quant 051517 workingmm.m
<b>Position</b>	P2-A2	<b>Sample Info</b>	
<b>Inj Vol</b>	-1	<b>Comment</b>	Hemostat 445283-1

### Sample Chromatogram



# ISP FORENSICS - Pocatello Instrument # 59740

## Cannabinoids Analysis Report


DS

**Batch Data Path** C:\MassHunter\Data\2019\AM 27\022719 THCQ TS SP\QuantResults\THCQ SP.batch.bin  
**Analysis Time** 2/28/2019 12:45 PM **Analyst Name** datastor  
**Report Time** 3/1/2019 2:15 PM **Reporter Name** datastor  
**Last Calib Update** 2/28/2019 12:45 PM **Batch State** Processed

**Analysis Info**

**Acq Time** 2019-02-27 13:58 **Data File** External Control.d  
**Sample Type** Sample **Sample Name** External Control  
**Dilution** 1 **Acq Method** THC Quant 051517 workingmm.m  
**Position** P2-B2 **Sample Info**  
**Inj Vol** -1 **Comment** Hemostat 445283-1 + WS020419

**Sample Chromatogram**



**Results**

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-D3	2.265	23851	315259	0.0757	7.8231
THC-COOH	THC-COOH-D9	2.379	20372	101129	0.2014	9.1015
THC	THC-D3	5.999	7229	105231	0.0687	7.6092


# ISP Forensics Calibration Curve Report

P

**Batch Data Path** C:\MassHunter\Data\2019\AM 27\022719 THCQ TS SP\QuantResults\THCQ SP.batch.bin

**Last Calib Update** 2/28/2019 12:45 PM      **Analyst Name** ISP TOX

**Target Compound** *THC-OH*  
**Internal Standard** *THC-OH-D3*



Sample	Level	Enabled	Exp Conc	Final Conc	Accuracy
Cal 1-3ng	1	<input checked="" type="checkbox"/>	3	3.5	115.8
Cal 2-5ng	2	<input checked="" type="checkbox"/>	5	4.6	92.6
QC Control	2	<input checked="" type="checkbox"/>	5	4.9	98.7
Cal 3-10ng	3	<input checked="" type="checkbox"/>	10	9.4	93.6
Cal 4-25ng	4	<input checked="" type="checkbox"/>	25	23.3	93.1
Cal 5-50ng	5	<input checked="" type="checkbox"/>	50	52.7	105.4
Cal 6-100ng	6	<input checked="" type="checkbox"/>	100	99.6	99.6


# ISP Forensics Calibration Curve Report

P

**Batch Data Path** C:\MassHunter\Data\2019\AM 27\022719 THCQ TS SP\QuantResults\THCQ SP.batch.bin

**Last Calib Update** 2/28/2019 12:45 PM      **Analyst Name** ISP TOX

**Target Compound** *THC-COOH*  
**Internal Standard** *THC-COOH-D9*



Sample	Level	Enabled	Exp Conc	Final Conc	Accuracy
Cal 1-3ng	1	<input type="checkbox"/>	3	2.7	89.9
Cal 2-5ng	2	<input type="checkbox"/>	5	5.0	100.7
QC Control	2	<input type="checkbox"/>	5	9.6	191.5
Cal 3-10ng	3	<input checked="" type="checkbox"/>	10	10.3	102.7
Cal 4-25ng	4	<input checked="" type="checkbox"/>	25	23.9	95.5
Cal 5-50ng	5	<input checked="" type="checkbox"/>	50	50.9	101.8
Cal 6-100ng	6	<input checked="" type="checkbox"/>	100	99.9	99.9




# ISP Forensics Calibration Curve Report

B

**Batch Data Path** C:\MassHunter\Data\2019\AM 27\022719 THCQ TS SP\QuantResults\THCQ SP.batch.bin

**Last Calib Update** 2/28/2019 12:45 PM      **Analyst Name** ISP TOX

**Target Compound** THC  
**Internal Standard** THC-D3



Sample	Level	Enabled	Exp Conc	Final Conc	Accuracy
Cal 1-3ng	1	<input checked="" type="checkbox"/>	3	3.4	112.4
Cal 2-5ng	2	<input checked="" type="checkbox"/>	5	5.3	105.7
QC Control	2	<input checked="" type="checkbox"/>	5	4.9	98.8
Cal 3-10ng	3	<input checked="" type="checkbox"/>	10	8.6	86.0
Cal 4-25ng	4	<input checked="" type="checkbox"/>	25	23.3	93.2
Cal 5-50ng	5	<input checked="" type="checkbox"/>	50	50.3	100.6
Cal 6-100ng	6	<input checked="" type="checkbox"/>	100	102.1	102.1

15

# ISP FORENSICS - Pocatello Instrument # 59740


## Cannabinoids Analysis Report

**Batch Data Path** C:\MassHunter\Data\2019\AM 27\022719 THCQ TS SP\QuantResults\THCQ SP.batch.bin  
**Analysis Time** 2/28/2019 12:45 PM **Analyst Name** datastor  
**Report Time** 3/1/2019 2:14 PM **Reporter Name** datastor  
**Last Calib Update** 2/28/2019 12:45 PM **Batch State** Processed

### Analysis Info

**Acq Time** 2019-02-27 11:48 **Data File** Cal 1-3ng.d  
**Sample Type** Calibration **Sample Name** Cal 1-3ng  
**Dilution** 1 **Acq Method** THC Quant 051517 workingmm.m  
**Position** P2-B1 **Sample Info**  
**Inj Vol** -1 **Comment**

### Sample Chromatogram



### Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-D3	2.279	12882	379958	0.0339	3.4737
THC-COOH	THC-COOH-D9	2.379	11233	131844	0.0852	2.6983
THC	THC-D3	6.012	4433	140585	0.0315	3.3714

# ISP FORENSICS - Pocatello Instrument # 59740

## Cannabinoids Analysis Report


P

**Batch Data Path** C:\MassHunter\Data\2019\AM 27\022719 THCQ TS SP\QuantResults\THCQ SP.batch.bin  
**Analysis Time** 2/28/2019 12:45 PM **Analyst Name** datastor  
**Report Time** 3/1/2019 2:14 PM **Reporter Name** datastor  
**Last Calib Update** 2/28/2019 12:45 PM **Batch State** Processed

### Analysis Info

**Acq Time** 2019-02-27 11:59 **Data File** Cal 2-5ng.d  
**Sample Type** Calibration **Sample Name** Cal 2-5ng  
**Dilution** 1 **Acq Method** THC Quant 051517 workingmm.m  
**Position** P2-C1 **Sample Info**  
**Inj Vol** -1 **Comment**

### Sample Chromatogram



### Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-D3	2.279	16953	376788	0.0450	4.6288
THC-COOH	THC-COOH-D9	2.379	16281	127552	0.1276	5.0363
THC	THC-D3	5.985	6589	136376	0.0483	5.2849

# ISP FORENSICS - Pocatello Instrument # 59740

## Cannabinoids Analysis Report


R

**Batch Data Path** C:\MassHunter\Data\2019\AM 27\022719 THCQ TS SP\QuantResults\THCQ SP.batch.bin  
**Analysis Time** 2/28/2019 12:45 PM **Analyst Name** datastor  
**Report Time** 3/1/2019 2:14 PM **Reporter Name** datastor  
**Last Calib Update** 2/28/2019 12:45 PM **Batch State** Processed

**Analysis Info**

**Acq Time** 2019-02-27 12:11 **Data File** Cal 3-10ng.d  
**Sample Type** Calibration **Sample Name** Cal 3-10ng  
**Dilution** 1 **Acq Method** THC Quant 051517 workingmm.m  
**Position** P2-D1 **Sample Info**  
**Inj Vol** -1 **Comment**

**Sample Chromatogram**



**Results**

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-D3	2.279	33615	371937	0.0904	9.3566
THC-COOH	THC-COOH-D9	2.365	27935	125419	0.2227	10.2743
THC	THC-D3	6.025	11816	152743	0.0774	8.5963

# ISP FORENSICS - Pocatello Instrument # 59740

## Cannabinoids Analysis Report


P

<b>Batch Data Path</b>	C:\MassHunter\Data\2019\AM 27\022719 THCQ TS SP\QuantResults\THCQ SP.batch.bin		
<b>Analysis Time</b>	2/28/2019 12:45 PM	<b>Analyst Name</b>	datastor
<b>Report Time</b>	3/1/2019 2:14 PM	<b>Reporter Name</b>	datastor
<b>Last Calib Update</b>	2/28/2019 12:45 PM	<b>Batch State</b>	Processed

**Analysis Info**

<b>Acq Time</b>	2019-02-27 12:23	<b>Data File</b>	Cal 4-25ng.d
<b>Sample Type</b>	Calibration	<b>Sample Name</b>	Cal 4-25ng
<b>Dilution</b>	1	<b>Acq Method</b>	THC Quant 051517 workingmm.m
<b>Position</b>	P2-E1	<b>Sample Info</b>	
<b>Inj Vol</b>	-1	<b>Comment</b>	

**Sample Chromatogram**



**Results**

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-D3	2.279	75286	336054	0.2240	23.2794
THC-COOH	THC-COOH-D9	2.365	54597	116274	0.4696	23.8699
THC	THC-D3	6.025	26572	128763	0.2064	23.3052

P

# ISP FORENSICS - Pocatello Instrument # 59740


## Cannabinoids Analysis Report

**Batch Data Path** C:\MassHunter\Data\2019\AM 27\022719 THCQ TS SP\QuantResults\THCQ SP.batch.bin  
**Analysis Time** 2/28/2019 12:45 PM **Analyst Name** datastor  
**Report Time** 3/1/2019 2:14 PM **Reporter Name** datastor  
**Last Calib Update** 2/28/2019 12:45 PM **Batch State** Processed

### Analysis Info

**Acq Time** 2019-02-27 12:35 **Data File** Cal 5-50ng.d  
**Sample Type** Calibration **Sample Name** Cal 5-50ng  
**Dilution** 1 **Acq Method** THC Quant 051517 workingmm.m  
**Position** P2-F1 **Sample Info**  
**Inj Vol** -1 **Comment**

### Sample Chromatogram



### Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-D3	2.279	170351	336428	0.5064	52.6888
THC-COOH	THC-COOH-D9	2.379	108077	112504	0.9607	50.9214
THC	THC-D3	6.025	55380	124981	0.4431	50.2975

# ISP FORENSICS - Pocatello Instrument # 59740

## Cannabinoids Analysis Report


P

**Batch Data Path** C:\MassHunter\Data\2019\AM 27\022719 THCQ TS SP\QuantResults\THCQ SP.batch.bin  
**Analysis Time** 2/28/2019 12:45 PM **Analyst Name** datastor  
**Report Time** 3/1/2019 2:15 PM **Reporter Name** datastor  
**Last Calib Update** 2/28/2019 12:45 PM **Batch State** Processed

**Analysis Info**

<b>Acq Time</b>	2019-02-27 12:47	<b>Data File</b>	Cal 6-100ng.d
<b>Sample Type</b>	Calibration	<b>Sample Name</b>	Cal 6-100ng
<b>Dilution</b>	1	<b>Acq Method</b>	THC Quant 051517 workingmm.m
<b>Position</b>	P2-G1	<b>Sample Info</b>	
<b>Inj Vol</b>	-1	<b>Comment</b>	

**Sample Chromatogram**



**Results**

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
THC-OH	THC-OH-D3	2.279	327800	342736	0.9564	99.5727
THC-COOH	THC-COOH-D9	2.379	204746	110647	1.8504	99.9344
THC	THC-D3	6.039	109056	121464	0.8978	102.1447