


Worklist: 3445

<u>LAB_CASE</u>	<u>ITEM</u>	<u>TASK_ID</u>	<u>DESCRIPTION</u>
M2019-2094	1	153733	AM 27 Blood THC Quant by LC-QQQ
M2019-2231	3	153734	AM 27 Blood THC Quant by LC-QQQ
P2019-1302	1	153735	AM 27 Blood THC Quant by LC-QQQ
P2019-1424	1	153736	AM 27 Blood THC Quant by LC-QQQ
P2019-1425	1	153737	AM 27 Blood THC Quant by LC-QQQ
P2019-1426	2	153738	AM 27 Blood THC Quant by LC-QQQ
P2019-1463	1	153739	AM 27 Blood THC Quant by LC-QQQ
P2019-1467	1	153740	AM 27 Blood THC Quant by LC-QQQ
P2019-1491	1	153741	AM 27 Blood THC Quant by LC-QQQ
P2019-1498	1	153742	AM 27 Blood THC Quant by LC-QQQ
P2019-1511	1	153743	AM 27 Blood THC Quant by LC-QQQ
P2019-1515	1	153744	AM 27 Blood THC Quant by LC-QQQ



# AM# 27: Quantitation of THC and Metabolites in Blood by LC-MS/MS

Extraction Date: 6/3/19  
Plate lot#: 0539904

Analyst: Sarah Pickle  
Plate Expiration: 9/10/19

**Mobile phase A:** 0.1% Formic Acid in LCMS Water  
MTBE  
**Mobile phase B:** 0.1% Formic acid in Acetonitrile  
LCMS Methanol Hexane  
**Blank Blood Lot:** 445283-2  
**Column:** UCT Selectra DA 100 x 2.1mm 3um  
**LCMS-QQQ ID:** 069901

## Pre-Analytic:

- 1. Check levels of mobile phases and needle wash refill as needed. Ensure waste is not full.
- 2. Ensure correct column is installed and begin mobile phase flow allow to equilibrate ~ 30 minutes.

## Analytic:

- 1. Remove standards, plate, controls, and samples from cold storage. Allow to reach room temperature.
- 2. Pipette **1000 µL blood (calibrated pipette)** in wells of analytical (standards) plate. **Pipette ID: #3**
- 3. Place on shaking incubator at ambient temp., 900rpm for 15 minutes. *Shaker ID: 067105*
- 4. Pipette **500 µL 0.1% formic acid in LCMS water** in wells of analytical plate.
- 5. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 6. Transfer **800 µL of blood+acid** mixture to corresponding wells of SLE+ plate.
- 7. Apply positive pressure for approx. 10-15 seconds (or until no liquid remains on top of sorbent).  
*(Load at 85-95 PSI- Selector to the right)* Manifold ID: 067104
- 8. Wait 5 minutes.
- 9. Add **2.25 mL MTBE. (Add in 3 increments of 750 µL)**
- 10. Wait 5 minutes.
- 11. Apply positive pressure for approx. 15 seconds. *(12-15 PSI- Selector to the left).*
- 12. Add **2.25 mL Hexane. (Add in 3 increments of 750 µL)**
- 13. Wait 5 minutes.
- 14. Apply positive pressure for approx. 15 seconds. *(12-15 PSI- Selector to the left).*
- 15. Remove plate containing eluate. Place on SPE Dry and evaporate to dryness at approx. 35°C.  
*SPE Dry ID: 067103*
- 16. Reconstitute in **100 µL 100% MeOH** and heat seal plate with foil. Place in autosampler and run worklist.

## Post-Analytic

- 1. Create batch and process data.  
Worklist path: D:\MassHunter\Data\2019\AM 27\060319 THCQ SP TS Batch Name: THCQ SP
- 2. Make any necessary integration changes, Curve weighting of Linear 1/x with r<sup>2</sup> values ≥0.98 for each analyte
- 3. RT +/- 3% or 0.100 min, whichever is greater, +/- 20% Accuracy for greater than (+/- 30% for 10ng/ml or less).  
Ion ratios must be within +/- 20% of the averaged calibrators
- 4. Case sample response for THC and OH-THC 3ng/mL (quantitative), Carboxy-THC: 10ng/mL (qualitative only) will be reported. Samples with a THC or OH-THC response over 50 ng/mL will be reported out as greater than 50 ng/mL.
- 5. Did all QCs pass for each analyte? Y / N
- 6. Enter QCs into control charting.
- 7. Central File Packet to include: LIMS Worklist, Method Checklist, Calibration and Control Reports

COMMENTS: *Curve Range limited: THC-COOH 5-100*



# Idaho State Police Forensic Services

## AM #27 Quantitative Analysis of THC and Metabolites in Blood by LCMS-QQQ

### Methanol External Control Solution (Lot: WS041619)

10 ul of 1mg/mL THC, 100 ul of 100 ug/mL THC-OH, C-THC in 9790 ul MeOH  
Approximate concentration 1ug/mL.

Component	Source	Source Lot Number	Expiration Date
Methanol (LCMS)	Fisher	184782	
THC	Cerilliant	FE09101501	11/30/2020
C-THC	Cerilliant	FE07171501	09/30/2019
THC-OH	Cerilliant	FE01121503	01/31/2020
Prepared:	04/16/2019		
Prepared By:	Tamara Salazar		
Expires:	01/31/2020		

### Blood External Control Solution (Lot: 041619)

100 ul of methanol external control solution was added to 9900 ul of blood.  
Approximately 10ng/mL of each compound.

Component	Source	Source Lot Number
Negative Blood	Hemostat	445283-1
Methanol External Control Solution	-	WS041619
Prepared:	04/16/19	
Prepared by:	Tamara Salazar	
Expires:	01/31/2020	

# AM #27 Cannabinoids Quant. Results

S



## Batch results

D:\MassHunter\Data\2019\AM 27\060319 THCQ SP TS\QuantResults\THCQ SP.batch.bin

## Calibration Last Update

6/7/2019 10:20:25 AM

## Instrument

FALCO-LCMS (Property ID 069901)

## Type

Sample

## Data File

Negative.d

## Acq. Method

AM 27 THC quant.m

## Sample

Negative

## Sample Position

P3-A5

## Comment

## Injection Volume

10


## Acq. Date-Time

6/3/2019 12:09:35 PM

## Sample Info.

## Sample Chromatogram

+ TIC MRM (\*\* -> \*\*) Negative.d (Negative)



# AM #27 Cannabinoids Quant. Results

**Batch results**

D:\MassHunter\Data\2019\AM 27\060319 THCQ SP TS\QuantResults\THCQ SP.batch.bin

**Calibration Last Update**

6/7/2019 10:20:25 AM

**Instrument**

FALCO-LCMS (Property ID 069901)

**Data File**

External Control.d

**Type**

Sample

**Sample**

External Control

**Acq. Method**

AM 27 THC quant.m

**Sample Position**

P3-B5

**Comment****Injection Volume**


10

**Acq. Date-Time**

6/3/2019 12:24:46 PM

**Sample Info.****Sample Chromatogram**

+ TIC MRM (\*\* -&gt; \*\*) External Control.d (External Control)



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC	3.506	158066	619.60	28.6	143.57	2569463	7.7360 ng/ml
THC-COOH	1.580	34176	104.49	50.0	87.21	186079	8.4337 ng/ml
THC-OH	1.558	80184	∞	12.4	254.49	602267	8.5866 ng/ml

# AM #27 Cannabinoids Quant. Results

**Batch results**

D:\MassHunter\Data\2019\AM 27\060319 THCQ SP TS\QuantResults\THCQ SP.batch.bin

**Calibration Last Update**

6/10/2019 1:49:13 PM

**Instrument**

FALCO-LCMS (Property ID 069901)

**Data File**

QC Control.d

**Type**

Sample

**Sample**

QC Control

**Acq. Method**

AM 27 THC quant.m

**Sample Position**

P3-H4

**Comment****Injection Volume**


10

**Acq. Date-Time**

6/3/2019 11:54:24 AM

**Sample Info.****Sample Chromatogram**

+ TIC MRM (\*\* -&gt; \*\*) QC Control.d (QC Control)




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC	3.506	95272	661.46	28.4	$\infty$	2432603	4.9282 ng/ml
THC-COOH	1.580	42221	175.16	52.8	425.16	209701	9.3542 ng/ml
THC-OH	1.558	59001	152.02	11.8	71.82	732008	4.7568 ng/ml

SP



# AM #27 Cannabinoids Quant. Calibration Curve Report

**Batch results** D:\MassHunter\Data\2019\AM 27\060319 THCQ SP TS\QuantResults\THCQ SP.batch.bin  
**Last Cal. Update** 6/7/2019 10:20 AM  
**Analyst Name** ISP\datastor  
**Analyte** THC **Internal Standard** THC-D3




Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
Cal 1-3ng	1	✓	3.0	3.1	103.6
Cal 2- 5ng	2	✓	5.0	5.0	99.5
Cal 3 -10ng	3	✓	10.0	9.8	98.0
Cal 4-25ng	4	✓	25.0	24.5	97.8
Cal 5-50ng	5	✓	50.0	50.4	100.8
Cal 6-100ng	6	✓	100.0	100.3	100.3



# AM #27 Cannabinoids Quant. Calibration Curve Report

**Batch results** D:\MassHunter\Data\2019\AM 27\060319 THCQ SP TS\QuantResults\THCQ SP.batch.bin  
**Last Cal. Update** 6/7/2019 10:20 AM  
**Analyst Name** ISP\datastor  
**Analyte** THC-COOH **Internal Standard** THC-COOH-D9




Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
Cal 1-3ng	1	x	3.0	3.8	127.8
Cal 2- 5ng	2	✓	5.0	5.0	100.0
Cal 3 -10ng	3	✓	10.0	10.3	102.9
Cal 4-25ng	4	✓	25.0	24.0	96.0
Cal 5-50ng	5	✓	50.0	50.4	100.7
Cal 6-100ng	6	✓	100.0	100.3	100.3





# AM #27 Cannabinoids Quant. Calibration Curve Report

**Batch results** D:\MassHunter\Data\2019\AM 27\060319 THCQ SP TS\QuantResults\THCQ SP.batch.bin  
**Last Cal. Update** 6/7/2019 10:20 AM  
**Analyst Name** ISP\datastor  
**Analyte** THC-OH **Internal Standard** THC-OH-D3



Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
Cal 1-3ng	1	✓	3.0	3.1	102.3
Cal 2- 5ng	2	✓	5.0	4.8	96.8
Cal 3 -10ng	3	✓	10.0	9.8	97.8
Cal 4-25ng	4	✓	25.0	25.4	101.6
Cal 5-50ng	5	✓	50.0	51.7	103.3
Cal 6-100ng	6	✓	100.0	98.3	98.3


# AM #27 Cannabinoids Quant. Results



**Batch results** D:\MassHunter\Data\2019\AM 27\060319 THCQ SP TS\QuantResults\THCQ SP.batch.bin  
**Calibration Last Update** 6/7/2019 10:20:25 AM

**Instrument** FALCO-LCMS (Property ID 069901) **Data File** Cal 1-3ng.d  
**Type** Cal **Sample** Cal 1-3ng  
**Acq. Method** AM 27 THC quant.m  
**Sample Position** P3-B4 **Comment**  
**Injection Volume** 10  
**Acq. Date-Time** 6/3/2019 11:08:53 AM  
**Sample Info.**

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC	3.506	53914	323.26	29.2	38.88	2184878	3.1083 ng/ml
THC-COOH	1.595	18866	∞	42.7 <b>Low</b>	83.38	197909	3.8332 ng/ml
THC-OH	1.573	39232	∞	11.2	20.53	683007	3.0682 ng/ml

# AM #27 Cannabinoids Quant. Results

**Batch results**

D:\MassHunter\Data\2019\AM 27\060319 THCQ SP TS\QuantResults\THCQ SP.batch.bin

**Calibration Last Update**

6/7/2019 10:20:25 AM

**Instrument**

FALCO-LCMS (Property ID 069901)

**Type**

Cal

**Acq. Method**

AM 27 THC quant.m

**Sample Position**

P3-C4

**Injection Volume**

10

**Acq. Date-Time**

6/3/2019 11:16:28 AM

**Sample Info.****Data File**


Cal 2- 5ng.d

**Sample**

Cal 2- 5ng

**Comment****Sample Chromatogram**

+ TIC MRM (\*\* -&gt; \*\*) Cal 2- 5ng.d (Cal 2- 5ng)



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC	3.506	86688	2671.23	27.5	45.46	2191515	4.9774 ng/ml
THC-COOH	1.595	23215	∞	47.5	120.27	197202	4.9996 ng/ml
THC-OH	1.573	55200	71.38	12.0	50.23	675180	4.8409 ng/ml


# AM #27 Cannabinoids Quant. Results



**Batch results** D:\MassHunter\Data\2019\AM 27\060319 THCQ SP TS\QuantResults\THCQ SP.batch.bin  
**Calibration Last Update** 6/7/2019 10:20:25 AM

**Instrument** FALCO-LCMS (Property ID 069901) **Data File** Cal 3 -10ng.d  
**Type** Cal **Sample** Cal 3 -10ng  
**Acq. Method** AM 27 THC quant.m  
**Sample Position** P3-D4 **Comment**  
**Injection Volume** 10  
**Acq. Date-Time** 6/3/2019 11:24:03 AM  
**Sample Info.**

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC	3.506	183089	695.71	28.1	186.48	2349587	9.7969 ng/ml
THC-COOH	1.595	43308	91.02	51.4	239.91	197415	10.2934 ng/ml
THC-OH	1.558	103508	245.67	12.1	181.01	692552	9.7766 ng/ml


# AM #27 Cannabinoids Quant. Results



**Batch results** D:\MassHunter\Data\2019\AM 27\060319 THCQ SP TS\QuantResults\THCQ SP.batch.bin  
**Calibration Last Update** 6/7/2019 10:20:25 AM

**Instrument** FALCO-LCMS (Property ID 069901) **Data File** Cal 4-25ng.d  
**Type** Cal **Sample** Cal 4-25ng  
**Acq. Method** AM 27 THC quant.m  
**Sample Position** P3-E4 **Comment**  
**Injection Volume** 10  
**Acq. Date-Time** 6/3/2019 11:31:38 AM  
**Sample Info.**

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC	3.506	433507	∞	27.6	∞	2227133	24.4589 ng/ml
THC-COOH	1.595	91994	249.18	56.2	686.75	190633	23.9997 ng/ml
THC-OH	1.558	239479	600.80	13.4	520.79	658631	25.3879 ng/ml

# AM #27 Cannabinoids Quant. Results

**Batch results**

D:\MassHunter\Data\2019\AM 27\060319 THCQ SP TS\QuantResults\THCQ SP.batch.bin

**Calibration Last Update**

6/7/2019 10:20:25 AM

**Instrument**

FALCO-LCMS (Property ID 069901)

**Type**

Cal

**Acq. Method**

AM 27 THC quant.m

**Sample Position**

P3-F4

**Injection Volume**

10

**Acq. Date-Time**

6/3/2019 11:39:13 AM

**Sample Info.****Data File**


Cal 5-50ng.d

**Sample**

Cal 5-50ng

**Comment****Sample Chromatogram**

+ TIC MRM (\*\* -&gt; \*\*) Cal 5-50ng.d (Cal 5-50ng)



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC	3.506	890690	4015.70	27.7	480.34	2221145	50.3799 ng/ml
THC-COOH	1.595	189681	∞	56.8	986.19	191808	50.3681 ng/ml
THC-OH	1.558	502997	∞	12.8	1846.46	694722	51.6633 ng/ml


# AM #27 Cannabinoids Quant. Results



**Batch results** D:\MassHunter\Data\2019\AM 27\060319 THCQ SP TS\QuantResults\THCQ SP.batch.bin  
**Calibration Last Update** 6/7/2019 10:20:25 AM

**Instrument** FALCO-LCMS (Property ID 069901) **Data File** Cal 6-100ng.d  
**Type** Cal **Sample** Cal 6-100ng  
**Acq. Method** AM 27 THC quant.m  
**Sample Position** P3-G4 **Comment**  
**Injection Volume** 10  
**Acq. Date-Time** 6/3/2019 11:46:49 AM  
**Sample Info.**

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC	3.506	1555625	4936.02	27.4	∞	1948802	100.2786 ng/ml
THC-COOH	1.595	323443	661.18	58.5	2074.23	165997	100.3393 ng/ml
THC-OH	1.558	819913	1778.30	14.4	∞	601444	98.2630 ng/ml