









REVIEWED

By Britany Wylie at 2:08 pm, Sep 24, 2019

TS

9/18/2019

Worklist: 3692

<u>LAB CASE</u>	<u>ITEM</u>	<u>TASK ID</u>	<u>DESCRIPTION</u>	
P2019-2490	1	163857	AM 27 Blood THC Quant by LC-QQQ	
P2019-2517	1	163858	AM 27 Blood THC Quant by LC-QQQ	
P2019-2518	1	163864	AM 27 Blood THC Quant by LC-QQQ	
P2019-2533	1	163859	AM 27 Blood THC Quant by LC-QQQ	
P2019-2562	1	163860	AM 27 Blood THC Quant by LC-QQQ	
P2019-2572	1	163861	AM 27 Blood THC Quant by LC-QQQ	
P2019-2594	1	163862	AM 27 Blood THC Quant by LC-QQQ	
P2019-2777	1	163863	AM 27 Blood THC Quant by LC-QQQ	

TS

AM# 27: Quantitation of THC and Metabolites in Blood by LC-MS/MS

Extraction Date: 09/16/19
Plate lot#: IDP-108-190716

Analyst: Tamara Salazar
Plate Expiration: 01/16/2020

Mobile phase A: 0.1% Formic Acid in LCMS Water
MTBE LCMS Methanol

Mobile phase B: 0.1% Formic acid in Acetonitrile
Hexane

Blank Blood Lot: Hemostat 445283-2

Column: UCT Selectra DA 100 x 2.1mm 3um

LCMS-QQQ ID: 069901

Pre-Analytic:

- 1. Check levels of mobile phases and needle wash refill as needed. Ensure waste is not full.
- 2. Ensure correct column is installed and begin mobile phase flow allow to equilibrate ~ 30 minutes.
- 3. Create worklist:

Analytic:

- 1. Remove standards, plate, controls, and samples from cold storage. Allow to reach room temperature.
- 2. Pipette **1000µL blood/urine (calibrated pipette) Pipette ID: 3** in wells of analytical (standards) plate.
- 3. Place on shaking incubator at ambient temp., 900rpm for 15 minutes. *Shaker ID: 067105*
- 4. Pipette **500µL 0.1% formic acid in water** in wells of analytical plate for blood samples.
- 5. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 6. Transfer **800µL of blood+acid** mixture to corresponding wells of SLE+ plate.
- 7. Apply positive pressure for approx. 10-15 seconds (or until no liquid remains on top of sorbent).
(Load at 85-100 PSI- Selector to the right) Manifold ID: 067104
- 8. Wait 5 minutes.
- 9. Add **2.25mL MTBE. (Add in 3 increments of 750uL)**
- 10. Wait 5 minutes.
- 11. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 12. Add **2.25mL Hexane. (Add in 3 increments of 750uL)**
- 13. Wait 5 minutes.
- 14. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 15. Remove plate containing eluate. Place on SPE Dry and evaporate to dryness at approx. 35°C.
SPE Dry ID: 067103
- 16. Reconstitute in **100µL 100% MeOH** and heat seal plate with foil. Place in autosampler and run worklist.

Post-Analytic

- 1. Create batch and process data.
Worklist path: D:\MassHunter\Data\2019\AM 27\091619 MDS SP TS
Batch Name THCQ TS
- 2. Make any necessary integration changes, Curve weighting of Linear 1/x with r² values ≥0.98 for each analyte
- 3. RT +/- 3% or 0.100 min, whichever is greater, +/- 20% Accuracy for greater than (+/- 30% for 10ng/ml or less).
Ion ratios must be within +/- 20% of the averaged calibrators
- 4. Case sample response for THC and OH-THC 3ng/mL (quantitative), Carboxy-THC: 10ng/mL (qualitative only) will be reported. Samples with a THC or OH-THC response over 50 ng/mL will be reported out as greater than 50 ng/mL.
- 5. Did all QCs pass for each analyte? Y / N
- 6. Enter QCs into control charting.
- 7. Central File Packet to include: LIMS Worklist, Method Checklist, Calibration and Control Reports

COMMENTS: *Curves limited: THC 3-100, THC-COOH 10-250, THC-OH 3-100*




AM #27 Cannabinoids Quant. Results

Batch results D:\MassHunter\Data\2019\AM 25\091619 MDS SP TS\QuantResults\THCQ_TS.batch.bin
Calibration Last Update 9/17/2019 8:17:47 AM

Instrument Type	Falco Sample	Data File	TS_MJ_Negative.d
Acq. Method	AM 27 THC quant.m	Sample	TS_MJ_Negative
Sample Position	P4-A2	Comment	
Injection Volume	10		
Acq. Date-Time	9/17/2019 4:11:49 AM		

Sample Info.

Sample Chromatogram



TS




AM #27 Cannabinoids Quant. Results

Batch results D:\MassHunter\Data\2019\AM 27\091619 MDS SP TS\QuantResults\THCQ_TS.batch.bin
Calibration Last Update 9/21/2019 7:29:40 AM

Instrument Type	Falco Sample	Data File Sample	TS_MJ_QC Control.d
Acq. Method	AM 27 THC quant.m	Comment	TS_MJ_QC Control
Sample Position	P4-H1		
Injection Volume	10		
Acq. Date-Time	9/17/2019 3:56:36 AM		

Sample Chromatogram




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC	3.134	92669	620.25	30.1	∞	2628842	4.2905 ng/ml
THC-COOH	1.429	47061	767.00	55.0	844.82	141588	14.6005 ng/ml
THC-OH	1.408	46919	∞	10.1	41.72	525642	4.1229 ng/ml

5



AM #27 Cannabinoids Quant. Calibration Curve Report

Batch results D:\MassHunter\Data\2019\AM 27\091619 MDS SP TS\QuantResults\THCQ_TS.batch.bin
Last Cal. Update 9/21/2019 7:29 AM
Analyst Name ISP\datastor
Analyte THC **Internal Standard** THC-D3




Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
TS MJ Cal 1	1	×	1.0	1.2	116.0
TS MJ Cal 2	2	✓	3.0	3.0	98.7
TS MJ Cal 3	3	✓	5.0	4.8	95.1
TS MJ Cal 4	4	✓	10.0	9.5	95.1
TS MJ Cal 5	5	✓	25.0	28.8	115.2
TS MJ Cal 6	6	✓	50.0	48.9	97.7
TS MJ Cal 7	7	✓	100.0	98.1	98.1



AM #27 Cannabinoids Quant. Calibration Curve Report

Batch results D:\MassHunter\Data\2019\AM 27\091619 MDS SP TS\QuantResults\THCQ_TS.batch.bin
Last Cal. Update 9/21/2019 7:29 AM
Analyst Name ISP\datastor
Analyte THC-COOH **Internal Standard** THC-COOH-D9




Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
TS MJ Cal 1	1	×	5.0	6.3	125.4
TS MJ Cal 2	2	✓	10.0	10.4	103.6
TS MJ Cal 3	3	✓	20.0	18.4	92.2
TS MJ Cal 4	4	✓	50.0	47.6	95.2
TS MJ Cal 5	5	✓	75.0	85.1	113.5
TS MJ Cal 6	6	✓	100.0	96.9	96.9
TS MJ Cal 7	7	✓	250.0	246.6	98.6



AM #27 Cannabinoids Quant. Calibration Curve Report

Batch results D:\MassHunter\Data\2019\AM 27\091619 MDS SP TS\QuantResults\THCQ_TS.batch.bin
Last Cal. Update 9/21/2019 7:29 AM
Analyst Name ISP\datastor
Analyte THC-OH **Internal Standard** THC-OH-D3



Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
TS MJ Cal 1	1	×	1.0	1.0	97.3
TS MJ Cal 2	2	✓	3.0	2.6	88.0
TS MJ Cal 3	3	✓	5.0	4.9	98.1
TS MJ Cal 4	4	✓	10.0	9.9	98.7
TS MJ Cal 5	5	✓	25.0	29.8	119.1
TS MJ Cal 6	6	✓	50.0	50.2	100.3
TS MJ Cal 7	7	✓	100.0	95.6	95.6




AM #27 Cannabinoids Quant. Results

Batch results D:\MassHunter\Data\2019\AM 27\091619 MDS SP TS\QuantResults\THCQ_TS.batch.bin
Calibration Last Update 9/21/2019 7:29:40 AM

Instrument	Falco	Data File	TS_MJ_Cal 1.d
Type	Cal	Sample	TS_MJ_Cal 1
Acq. Method	AM 27 THC quant.m	Comment	
Sample Position	P4-A1		
Injection Volume	10		
Acq. Date-Time	9/17/2019 3:03:08 AM		

Sample Chromatogram

+ TIC: MRM (** -> **) TS_MJ_Cal 1.d (TS_MJ_Cal 1)



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC	3.134	30355	122.26	33.3	46.12	3795239	1.1599 ng/ml
THC-COOH	1.399	28129	∞	43.1 Low	108.70	183669	6.2692 ng/ml
THC-OH	1.393	32441	∞	5.1 Low	6.43 Low	685853	0.9733 ng/ml

TS




AM #27 Cannabinoids Quant. Results

Batch results D:\MassHunter\Data\2019\AM 27\091619 MDS SP TS\QuantResults\THCQ_TS_batch.bin
Calibration Last Update 9/21/2019 7:29:40 AM

Instrument	Falco	Data File	TS_MJ_Cal 2.d
Type	Cal	Sample	TS_MJ_Cal 2
Acq. Method	AM 27 THC quant.m	Comment	
Sample Position	P4-B1		
Injection Volume	10		
Acq. Date-Time	9/17/2019 3:11:00 AM		
Sample Info.			

Sample Chromatogram

+ TIC: MRM (** -> **) TS_MJ_Cal 2.d (TS_MJ_Cal 2)



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC	3.134	94127	203.67	30.5	∞	3973457	2.9624 ng/ml
THC-COOH	1.429	42965	84.85	52.0	141.26	178182	10.3588 ng/ml
THC-OH	1.408	45839	∞	9.7	35.01	659391	2.6410 ng/ml

TS




AM #27 Cannabinoids Quant. Results

Batch results D:\MassHunter\Data\2019\AM 27\091619 MDS SP TS\QuantResults\THCQ_TS_batch.bin
Calibration Last Update 9/21/2019 7:29:40 AM

Instrument Type	Falco Cal	Data File Sample	TS_MJ_Cal 3.d TS_MJ_Cal 3
Acq. Method	AM 27 THC quant.m	Comment	
Sample Position	P4-C1		
Injection Volume	10		
Acq. Date-Time	9/17/2019 3:18:37 AM		

Sample Chromatogram

+ TIC:MRM (**->**) TS_MJ_Cal 3.d (TS_MJ_Cal 3)



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC	3.194	139526	276.25	28.7	45.85	3552354	4.7530 ng/ml
THC-COOH	1.444	71095	149.52	60.5	485.19	171329	18.4391 ng/ml
THC-OH	1.423	63395	∞	10.0	67.56	635783	4.9073 ng/ml

TS




AM #27 Cannabinoids Quant. Results

Batch results D:\MassHunter\Data\2019\AM 27\091619 MDS SP TS\QuantResults\THCQ_TS_batch.bin
Calibration Last Update 9/21/2019 7:29:40 AM

Instrument	Falco	Data File	TS_MJ_Cal 4.d
Type	Cal	Sample	TS_MJ_Cal 4
Acq. Method	AM 27 THC quant.m	Comment	
Sample Position	P4-D1		
Injection Volume	10		
Acq. Date-Time	9/17/2019 3:26:11 AM		
Sample Info.			

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC	3.179	271734	406.55	28.5	401.77	3366741	9.5126 ng/ml
THC-COOH	1.444	177457	366.18	61.9	947.28	170231	47.6070 ng/ml
THC-OH	1.438	104423	∞	12.0	158.13	629645	9.8710 ng/ml

TS




AM #27 Cannabinoids Quant. Results

Batch results D:\MassHunter\Data\2019\AM 27\091619 MDS SP TS\QuantResults\THCQ_TS_batch.bin
Calibration Last Update 9/21/2019 7:29:40 AM

Instrument	Falco	Data File	TS_MJ_Cal 5.d
Type	Cal	Sample	TS_MJ_Cal 5
Acq. Method	AM 27 THC quant.m	Comment	
Sample Position	P4-E1		
Injection Volume	10		
Acq. Date-Time	9/17/2019 3:33:47 AM		

Sample Chromatogram

+ TIC MIRM (** -> **) TS_MJ_Cal 5.d (TS_MJ_Cal 5)



Name	RT	Resp.	S/N	Ratio	ISTD Resp.	Final Conc.
THC	3.149	789159	∞	27.5	3173279	28.8085 ng/ml
THC-COOH	1.429	294342	1514.72	62.5	159207	85.0896 ng/ml
THC-OH	1.408	259511	∞	13.4	601862	29.7866 ng/ml

5




AM #27 Cannabinoids Quant. Results

Batch results D:\MassHunter\Data\2019\AM 27\091619 MDS SP TS\QuantResults\THCQ_TS_batch.bin
Calibration Last Update 9/21/2019 7:29:40 AM

Instrument	Falco	Data File	TS_MJ_Cal 6.d
Type	Cal	Sample	TS_MJ_Cal 6
Acq. Method	AM 27 THC quant.m	Comment	
Sample Position	P4-F1		
Injection Volume	10		
Acq. Date-Time	9/17/2019 3:41:24 AM		
Sample Info.			

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC	3.164	1209231	1416.50	27.6	∞	2856439	48.8705 ng/ml
THC-COOH	1.444	316856	252.75	62.6	4373.37	150660	96.9117 ng/ml
THC-OH	1.423	405611	∞	13.8	1164.80	577189	50.1690 ng/ml

TS




AM #27 Cannabinoids Quant. Results

Batch results D:\MassHunter\Data\2019\AM 27\091619 MDS SP TS\QuantResults\THCQ_TS.batch.bin
Calibration Last Update 9/21/2019 7:29:40 AM

Instrument	Falco	Data File	TS_MJ_Cal 7.d
Type	Cal	Sample	TS_MJ_Cal 7
Acq. Method	AM 27 THC quant.m	Comment	
Sample Position	P4-G1		
Injection Volume	10		
Acq. Date-Time	9/17/2019 3:49:00 AM		

Sample Chromatogram

+ TIC:MRM (** -> **) TS_MJ_Cal 7.d (TS_MJ_Cal 7)



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC	3.119	2365150	8325.08	27.3	∞	2776538	98.0930 ng/ml
THC-COOH	1.414	731354	2182.66	62.8	3436.21	137389	246.5938 ng/ml
THC-OH	1.393	712130	∞	14.2	31122.40	544297	95.6251 ng/ml

5