


REVIEWED
By Sarah Pickle at 8:36 am, Nov 19, 2019

11/13/2019 *TS*
CS

Worklist: 3818

<u>LAB CASE</u>	<u>ITEM</u>	<u>ITEM TYPE</u>	<u>DESCRIPTION</u>
M2019-4750	3	BCK	AM 27 Blood THC Quant by LC-QQQ
M2019-4870	2	BCK	AM 27 Blood THC Quant by LC-QQQ
P2019-3261	1	BCK	AM 27 Blood THC Quant by LC-QQQ
P2019-3277	1	BCK	AM 27 Blood THC Quant by LC-QQQ
P2019-3290	1	BCK	AM 27 Blood THC Quant by LC-QQQ
P2019-3325	1	BCK	AM 27 Blood THC Quant by LC-QQQ



TS
CS

AM# 27: Quantitation of THC and Metabolites in Blood by LC-MS/MS

Extraction Date: 11/08/19
Plate lot#: IDP-108-190716

Analyst: Tamara Salazar
Plate Expiration: 01/16/2020

Mobile phase A: 0.1% Formic Acid in LCMS Water
MTBE LCMS Methanol

Mobile phase B: 0.1% Formic acid in Acetonitrile
Hexane

Blank Blood Lot: Hemostat 445283-3
LCMS-QQQ ID: 069901

Column: UCT Selectra DA 100 x 2.1mm 3um

Pre-Analytic:

- 1. Check levels of mobile phases and needle wash refill as needed. Ensure waste is not full.
- 2. Ensure correct column is installed and begin mobile phase flow allow to equilibrate ~ 30 minutes.
- 3. Create worklist:

Analytic:

- 1. Remove standards, plate, controls, and samples from cold storage. Allow to reach room temperature.
- 2. Pipette **1000µL blood/urine (calibrated pipette) Pipette ID: 3** in wells of analytical (standards) plate.
- 3. Place on shaking incubator at ambient temp., 900rpm for 15 minutes. *Shaker ID: 067105*
- 4. Pipette **500µL 0.1% formic acid in water for blood samples, 500µl saturated phosphate buffer for urine samples** in wells of analytical plate.
- 5. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 6. Transfer **800µL of blood+acid** mixture to corresponding wells of SLE+ plate.
- 7. Apply positive pressure for approx. 10-15 seconds (or until no liquid remains on top of sorbent). *(Load at 85-100 PSI- Selector to the right) Manifold ID: 067104*
- 8. Wait 5 minutes.
- 9. Add **2.25mL MTBE. (Add in 3 increments of 750uL)**
- 10. Wait 5 minutes.
- 11. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 12. Add **2.25mL Hexane. (Add in 3 increments of 750uL)**
- 13. Wait 5 minutes.
- 14. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 15. Remove plate containing eluate. Place on SPE Dry and evaporate to dryness at approx. 35°C. *SPE Dry ID: 067103*
- 16. Reconstitute in **100µL 100% MeOH** and heat seal plate with foil. Place in autosampler and run worklist.

Post-Analytic

- 1. Create batch and process data.
Worklist path: D:\MassHunter\Data\2019\AM 27\110819 THCQ TS Meth Test SP
Batch Name: THCQ TS wklst 3818
- 2. Make any necessary integration changes, Curve weighting of Linear 1/x with r² values ≥0.98 for each analyte
- 3. RT +/- 3% or 0.100 min, whichever is greater, +/- 20% Accuracy for greater than (+/- 30% for 10ng/ml or less).
Ion ratios must be within +/- 20% of the averaged calibrators
- 4. Case sample response for THC and OH-THC 3ng/mL (quantitative), Carboxy-THC: 10ng/mL (qualitative only) will be reported. Samples with a THC or OH-THC response over 50 ng/mL will be reported out as greater than 50 ng/mL.
- 5. Did all QCs pass for each analyte? Y / N
- 6. Enter QCs into control charting.
- 7. Central File Packet to include: LIMS Worklist, Method Checklist, Calibration and Control Reports

COMMENTS: *Curves limited: THC-OH 3-100* THC-OH was not evaluated for the urine cases due to the external control concentration being lower than anticipated. CS
Urine controls were included for the re-extraction of urine samples from another worklist.



Idaho State Police Forensic Services

TS
CS

AM #27 Quantitative Analysis of THC and Metabolites in Blood and Urine by LCMS-QQQ

Methanol External Control Solution (Lot: WS041619)

10 ul of 1mg/mL THC, 100 ul of 100 ug/mL THC-OH, C-THC in 9790 ul MeOH
Approximate concentration 1ug/mL.

Component	Source	Source Lot Number	Expiration Date
Methanol (LCMS)	Fisher	184782	
THC	Cerilliant	FE09101501	11/30/2020
C-THC	Cerilliant	FE07171501	09/30/2020
THC-OH	Cerilliant	FE01121503	01/31/2020
Prepared:	04/16/2019		
Prepared By:	Tamara Salazar		
Expires:	01/31/2020		

Urine External Control Solution (Lot: 100719)

100 ul of methanol external control solution was added to 9900 ul of urine.
Approximately 10ng/mL of each compound.

Component	Source	Source Lot Number
Negative Urine	Pocatello Lab	POC031319
Methanol External Control Solution	-	WS041619
Prepared:	10/07/19	
Prepared by:	Celena Shrum	
Expires:	01/31/2020	




AM #27 Cannabinoids Quant. Results

Batch results D:\MassHunter\Data\2019\AM 27\110819 THCQ TS Meth Test SP\QuantResults\THCQ TS wklist 3818.batch.bin
Calibration Last Update 11/18/2019 8:13:38 AM

Instrument Type	Sample	Data File	
Falco	MJ_Negative.d		
Sample	MJ_Negative		
Acq. Method	AM 27 THC quant.m		
Sample Position	P3-H5		
Injection Volume	10		
Acq. Date-Time	11/8/2019 9:00:05 PM		

Sample Info.

Sample Chromatogram



TS
S




AM #27 Cannabinoids Quant. Results

Batch results
Calibration Last Update

Instrument Type: Falco
Sample: AM 27 THC quant.m
Acq. Method: P3-A6
Sample Position: 10
Injection Volume: 11/8/2019 8:44:54 PM
Acq. Date-Time: MJ_QC Control.d
Sample Info: MJ_QC Control

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	ISTD Resp.	Final Conc.
THC	3.405	158546	191.69	26.4	4650249	4.3496 ng/ml
THC-COOH	1.504	100174	406.41	57.5	297027	14.1836 ng/ml
THC-OH	1.468	90765	89.84	10.3	1184715	4.2538 ng/ml

ST




AM #27 Cannabinoids Quant. Results

Batch results Calibration Last Update D:\MassHunter\Data\2019\AM 27\110819 THCQ TS Meth Test SP\QuantResults\THCQ TS wk1st 3818.batch.bin
 11/18/2019 8:13:38 AM

Instrument	Falco	Data File	MJ_urine ext. ctrl.d
Type	Sample	Sample	MJ_urine ext. ctrl
Acq. Method	AM 27 THC quant.m	Comment	
Sample Position	P3-H4		
Injection Volume	10		
Acq. Date-Time	11/8/2019 11:01:40 PM		
Sample Info.			

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC	3.405	20836	92.14	26.0	25.36	658257	4.0458 ng/ml
THC-COOH	1.504	27995	186.54	61.3	159.75	211927	5.0908 ng/ml
THC-OH	1.468	55272	263.19	14.2	315.03	974283	2.5839 ng/ml


25



AM #27 Cannabinoids Quant. Calibration Curve Report

Batch results D:\MassHunter\Data\2019\AM 27\110819 THCQ TS Meth Test SP\QuantResults\THCQ TS wk1st 3818.batch.bin
Last Cal. Update 11/18/2019 8:13 AM
Analyst Name ISP\datastor
Analyte THC **Internal Standard** THC-D3

THC - 7 Levels, 7 Levels Used, 7 Points, 7 Points Used, 0 QCs




Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
MJ Cal 1	1	✓	1.0	1.1	105.9
MJ Cal 2	2	✓	3.0	3.0	98.8
MJ Cal 3	3	✓	5.0	4.6	92.9
MJ Cal 4	4	✓	10.0	10.4	104.2
MJ Cal 5	5	✓	25.0	24.5	97.9
MJ Cal 6	6	✓	50.0	49.8	99.5
MJ Cal 7	7	✓	100.0	100.7	100.7



AM #27 Cannabinoids Quant. Calibration Curve Report

Batch results D:\MassHunter\Data\2019\AM 27\110819 THCQ TS Meth Test SP\QuantResults\THCQ TS wklst 3818.batch.bin
Last Cal. Update 11/18/2019 8:13 AM
Analyst Name ISP\datastor
Analyte THC-COOH **Internal Standard** THC-COOH-D9

THC-COOH - 7 Levels, 7 Levels Used, 7 Points, 7 Points Used, 0 QCs




Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
MJ Cal 1	1	✓	5.0	5.3	106.7
MJ Cal 2	2	✓	10.0	9.3	93.2
MJ Cal 3	3	✓	20.0	19.2	95.9
MJ Cal 4	4	✓	50.0	52.3	104.6
MJ Cal 5	5	✓	75.0	75.9	101.2
MJ Cal 6	6	✓	100.0	98.6	98.6
MJ Cal 7	7	✓	250.0	249.4	99.7



AM #27 Cannabinoids Quant. Calibration Curve Report

Batch results D:\MassHunter\Data\2019\AM 27\110819 THCQ TS Meth Test SP\QuantResults\THCQ TS wk1st
 3818.batch.bin
Last Cal. Update 11/18/2019 8:13 AM
Analyst Name ISP\datastor
Analyte THC-OH **Internal Standard** THC-OH-D3



Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
MJ Cal 1	1	×	1.0	0.7	69.5
MJ Cal 2	2	✓	3.0	3.3	110.3
MJ Cal 3	3	✓	5.0	4.4	88.6
MJ Cal 4	4	✓	10.0	10.2	102.2
MJ Cal 5	5	✓	25.0	24.7	98.8
MJ Cal 6	6	✓	50.0	49.8	99.6
MJ Cal 7	7	✓	100.0	100.5	100.5




AM #27 Cannabinoids Quant. Results

Batch results
Calibration Last Update

D:\MassHunter\Data\2019\AM 27\110819 THCQ TS Meth Test SP\QuantResults\THCQ TS wklist 3818.batch.bin
11/18/2019 8:13:38 AM

Instrument	Falco	Data File	MJ_Cal 1.d
Type	Cal	Sample	MJ_Cal 1
Acq. Method	AM 27 THC quant.m	Comment	
Sample Position	P3-B6		
Injection Volume	10		
Acq. Date-Time	11/8/2019 7:44:03 PM		

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC	3.435	32264	128.04	30.5	15.27	4210422	1.0595 ng/ml
THC-COOH	1.504	39255	∞	49.3	222.84	285274	5.3350 ng/ml
THC-OH	1.513 High	37994	∞	6.2 Low	10.28	1109545	0.6952 ng/ml

TS
S




AM #27 Cannabinoids Quant. Results

Batch results
 Calibration Last Update
 D:\MassHunter\Data\2019\AM 27\110819 THCQ TS Meth Test SP\QuantResults\THCQ TS wk1st 3818.batch.bin
 11/18/2019 8:13:38 AM

Instrument	Falco	Data File	MJ_Cal 2.d
Type	Cal	Sample	MJ_Cal 2
Acq. Method	AM 27 THC quant.m	Comment	
Sample Position	P3-C6		
Injection Volume	10		
Acq. Date-Time	11/8/2019 7:51:48 PM		

Sample Chromatogram

+ TIC MIRM (** -> **) MJ_Cal 2.d (MJ_Cal 2)



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC	3.450	92710	290.08	28.8	127.01	4035611	2.9653 ng/ml
THC-COOH	1.519	63088	265.51	58.3	221.09	277226	9.3223 ng/ml
THC-OH	1.483	70773	∞	10.2	200.35	1082932	3.3081 ng/ml

TS
 &




AM #27 Cannabinoids Quant. Results

Batch results Calibration Last Update D:\MassHunter\Data\2019\AM 27\110819 THCQ TS Meth Test SP\QuantResults\THCQ TS wk1st 3818.batch.bin
11/18/2019 8:13:38 AM

Instrument	Falco	Data File	MJ_Cal 3.d
Type	Cal	Sample	MJ_Cal 3
Acq. Method	AM 27 THC quant.m	Comment	
Sample Position	P3-D6		
Injection Volume	10		
Acq. Date-Time	11/8/2019 7:59:23 PM		

Sample Chromatogram

+ TIC MRM (** -> **) MJ_Cal 3.d (MJ_Cal 3)



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC	3.435	164127	∞	28.9	218.28	4500799	4.6449 ng/ml
THC-COOH	1.519	129582	∞	57.0	662.01	287964	19.1803 ng/ml
THC-OH	1.483	89521	∞	11.0	73.10	1137274	4.4303 ng/ml

TS
S



AM #27 Cannabinoids Quant. Results


Batch results
Calibration Last Update

Instrument: Falco
Type: Cal
Acq. Method: AM 27 THC quant.m
Sample Position: P3-E6
Injection Volume: 10
Acq. Date-Time: 11/8/2019 8:06:57 PM
Sample Info:

Data File: MJ_Cal 4.d
Sample: MJ_Cal 4
Comment:

Sample Chromatogram

+ TIC MRM (** -> **) MJ_Cal 4.d (MJ_Cal 4)



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC	3.420	346623	711.78	27.6	∞	4183144	10.4202 ng/ml
THC-COOH	1.519	337675	631.84	59.8	2206.33	281936	52.3194 ng/ml
THC-OH	1.468	163966	231.35	12.1	85.78	1110618	10.2189 ng/ml

Handwritten initials 'S' and 'B' in blue ink.



AM #27 Cannabinoids Quant. Results


Batch results
Calibration Last Update
D:\MassHunter\Data\2019\AM 27\110819 THCQ TS Meth Test.SP\QuantResults\THCQ TS wk1st 3818.batch.bin
11/18/2019 8:13:38 AM

Instrument
Type
Acq. Method
Sample Position
Injection Volume
Acq. Date-Time
Sample Info.

Falco
Cal
AM 27 THC quant.m
P3-F6
10
11/8/2019 8:14:32 PM

Data File
Sample
Comment

MJ_Cal 5.d
MJ_Cal 5



Sample Chromatogram
+ TIC-MIRM (** -> **) MJ_Cal 5.d (MJ_Cal 5)

Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC	3.420	747926	∞	26.8	702.49	3820372	24.4752 ng/ml
THC-COOH	1.504	442361	1142.87	59.7	4493.41	255689	75.9147 ng/ml
THC-OH	1.468	324817	521.37	13.2	557.82	1014985	24.6976 ng/ml

Handwritten initials 'S' and 'B' in blue ink.




AM #27 Cannabinoids Quant. Results

Batch results D:\MassHunter\Data\2019\AM 27\110819 THCQ TS Meth Test SP\QuantResults\THCQ TS wk1st 3818.batch.bin
 Calibration Last Update 11/18/2019 8:13:38 AM

Instrument Type	Falco Cal	Data File Sample	MJ_Cal 6.d MJ_Cal 6
Acq. Method	AM 27 THC quant.m	Comment	
Sample Position	P3-G6		
Injection Volume	10		
Acq. Date-Time	11/8/2019 8:22:07 PM		

Sample Chromatogram

+ TIC MRM (** -> **) MJ_Cal 6.d (MJ_Cal 6)



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC	3.420	1459241	13925.71	26.6	2270.90	3657155	49.7738 ng/ml
THC-COOH	1.504	542895	∞	60.9	4213.54	242211	98.5779 ng/ml
THC-OH	1.468	602451	∞	13.9	∞	973449	49.7989 ng/ml

SD




AM #27 Cannabinoids Quant. Results

Batch results Calibration Last Update D:\MassHunter\Data\2019\AM 27\110819 THCQ TS Meth Test SP\QuantResults\THCQ TS wkst 3818.batch.bin
 11/18/2019 8:13:38 AM

Instrument	Falco	Data File	MJ_Cal 7.d
Type	Cal	Sample	MJ_Cal 7
Acq. Method	AM 27 THC quant.m	Comment	
Sample Position	P3-H6		
Injection Volume	10		
Acq. Date-Time	11/8/2019 8:29:42 PM		
Sample Info.			

Sample Chromatogram

+ TIC.MRM (**->**) MJ_Cal 7.d (MJ_Cal 7)



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC	3.420	2946889	6930.58	27.0	∞	3647982	100.6612 ng/ml
THC-COOH	1.504	1329480	∞	60.5	5380.55	235588	249.3504 ng/ml
THC-OH	1.468	1208435	∞	13.7	6452.95	988020	100.5461 ng/ml

Handwritten initials/signature in blue ink.