

**REVIEWED**


By Sarah Pickle at 9:23 am, Dec 12, 2019

12/10/2019

BW

**Worklist: 3874**

<u>LAB CASE</u>	<u>ITEM</u>	<u>ITEM TYPE</u>	<u>DESCRIPTION</u>
C2019-2197	1	BCK	AM 27 Blood THC Quant by LC-QQQ
C2019-2254	1	BCK	AM 27 Blood THC Quant by LC-QQQ
C2019-2256	1	BCK	AM 27 Blood THC Quant by LC-QQQ
C2019-2289	1	BCK	AM 27 Blood THC Quant by LC-QQQ
C2019-2302	1	BCK	AM 27 Blood THC Quant by LC-QQQ
C2019-2334	1	BCK	AM 27 Blood THC Quant by LC-QQQ



BW

# AM# 27: Quantitation of THC and Metabolites in Blood and Urine by LC-MS/MS

Extraction Date: 12/9/19  
Plate lot#: 190716

Analyst: Britany Wylie  
Plate Expiration: 1-16-2020

**Mobile phase A:** 0.1% Formic Acid in LCMS Water MTBE      **Mobile phase B:** 0.1% Formic acid in Acetonitrile Hexane  
LCMS Methanol

**Blank Blood Lot:** 19H52275      **Column:** UCT Selectra DA 100 x 2.1mm 3um  
**LCMS-QQQ ID:** 69679

## Pre-Analytic:

- 1. Check levels of mobile phases and needle wash refill as needed. Ensure waste is not full.
- 2. Ensure correct column is installed and begin mobile phase flow allow to equilibrate ~ 30 minutes.

## Analytic:

- 1. Remove standards, plate, controls, and samples from cold storage. Allow to reach room temperature.
- 2. Pipette **1000µL blood (calibrated pipette) Pipette ID: k52558g** in wells of analytical (standards) plate.
- 3. Place on shaking incubator at ambient temp., 900rpm for 15 minutes. *Shaker ID: 66759*
- 4. Pipette **500µL 0.1% formic acid in water blood sample** in wells of analytical plate.
- 5. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 6. Transfer **800µL of blood+acid** mixture to corresponding wells of SLE+ plate.
- 7. Apply positive pressure for approx. 10-15 seconds (or until no liquid remains on top of sorbent). *(Load at 85-100 PSI- Selector to the right) Manifold ID: 66792*
- 8. Wait 5 minutes.
- 9. Add **2.25mL MTBE. (Add in 3 increments of 750uL)**
- 10. Wait 5 minutes.
- 11. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 12. Add **2.25mL Hexane. (Add in 3 increments of 750uL)**
- 13. Wait 5 minutes.
- 14. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 15. Remove plate containing eluate. Place on SPE Dry and evaporate to dryness at approx. 35°C. *SPE Dry ID: 66819*
- 16. Reconstitute in **100µL 100% MeOH** and heat seal plate with foil. Place in autosampler and run worklist.

## Post-Analytic

- 1. Create batch and process data.
- 2. Make any necessary integration changes, Curve weighting of Linear 1/x with r<sup>2</sup> values ≥0.98 for each analyte
- 3. RT +/- 3% or 0.100 min, whichever is greater, +/- 20% Accuracy for greater than (+/- 30% for 10ng/ml or less). Ion ratios must be within +/- 20% of the averaged calibrators
- 4. Case sample response for THC and OH-THC 3ng/mL (quantitative blood), Carboxy-THC: 10ng/mL (qualitative only) will be reported. Samples with a THC or OH-THC response over 50 ng/mL will be reported out as greater than 50 ng/mL.
- 5. Did all QCs pass for each analyte? (if not is it describe in comments section)
- 6. Enter QCs into control charting.
- 7. Central File Packet to include: LIMS Worklist, Method Checklist, Calibration and Control Reports

COMMENTS: *Blood samples only in this run. Curve ranges: carboxy-THC: 10-250, THC 3-100; Not Evaluated: THC-OH*


# AM #27 Cannabinoids

BW

**Batch results** D:\MassHunter\Data\2019\am27-28 120919\QuantResults\cann quant.batch.bin  
**Calibration Last Update** 12/11/2019 11:49:39 AM

<b>Instrument</b>	69679	<b>Data File</b>	mj negative blood.d
<b>Type</b>	Sample	<b>Sample</b>	mj negative blood
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Britany Wylie
<b>Sample Position</b>	P3-A2	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	12/9/2019 12:47:07 PM		
<b>Sample Info.</b>			

## Sample Chromatogram




# AM #27 Cannabinoids

BW

**Batch results** D:\MassHunter\Data\2019\am27-28 120919\QuantResults\cann quant.batch.bin  
**Calibration Last Update** 12/11/2019 11:49:39 AM

<b>Instrument</b>	69679	<b>Data File</b>	mj internal control.d
<b>Type</b>	QC	<b>Sample</b>	mj internal control
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Britany Wylie
<b>Sample Position</b>	P3-H1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	12/9/2019 12:39:25 PM		
<b>Sample Info.</b>			

## Sample Chromatogram




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-COOH	1.040	103027	455.3	36.6	1581.1	621884	14.439 ng/ml
THC-OH	<del>1.044</del>	<del>173010</del>	∞	<del>7.4</del> <b>Low</b>	∞	<del>1741032</del>	<del>4.420</del> ng/ml <b>NE</b>
THC	1.935	19825	467.1	25.8	∞	691200	4.318 ng/ml

# Compound Calibration Report

**Batch results**      D:\MassHunter\Data\2019\am27-28 120919\QuantResults\cann quant.batch.bin  
**Last Cal. Update**    12/11/2019 11:49 AM  
**Analyst Name**        ISP\datastor  
**Analyte**                THC                                **Internal Standard**      THC-d3

BW

THC - 6 Levels, 6 Levels Used, 6 Points, 6 Points Used, 2 QCs



Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
mj cal 6 <span style="color: blue;">cal 2</span>	2	✓	3.0	3.2	106.7
mj cal 5 <span style="color: blue;">cal 3</span>	3	✓	5.0	5.0	100.5
mj cal 4	4	✓	10.0	9.6	96.3
mj cal 3 <span style="color: blue;">cal 5</span>	5	✓	25.0	23.8	95.0
mj cal 2 <span style="color: blue;">cal 6</span>	6	✓	50.0	50.1	100.2
mj qc1 <span style="color: blue;">cal 7</span>	7	✓	100.0	101.3	101.3

BW 12/11/19


# Compound Calibration Report

**Batch results** D:\MassHunter\Data\2019\am27-28 120919\QuantResults\cann quant.batch.bin  
**Last Cal. Update** 12/11/2019 11:49 AM  
**Analyst Name** ISP\datastor  
**Analyte** THC-COOH

**Internal Standard** THC-COOH-d9

BW

THC-COOH - 6 Levels, 6 Levels Used, 6 Points, 6 Points Used, 2 QCs




Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
mj <del>cal 6</del> cal 2	2	✓	10.0	10.3	103.4
mj <del>cal 5</del> cal 3	3	✓	20.0	19.4	97.2
mj cal 4	4	✓	50.0	49.8	99.6
mj <del>cal 3</del> cal 5	5	✓	75.0	75.5	100.7
mj <del>cal 2</del> cal 6	6	✓	100.0	98.6	98.6
mj <del>cal 1</del> cal 7	7	✓	250.0	251.3	100.5

BW 12/11/19

# Compound Calibration Report

**Batch results** D:\MassHunter\Data\2019\am27-28 120919\QuantResults\cann quant.batch.bin  
**Last Cal. Update** 12/11/2019 11:49 AM  
**Analyst Name** ISP\datastor  
**Analyte** THC-OH **Internal Standard** THC-OH-d3 *BW*

THC-OH - 6 Levels, 6 Levels Used, 6 Points, 6 Points Used, 2 QCs



Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
mj cal 6	2	✓	3.0	3.1	104.1
mj cal 5	3	✓	5.0	4.9	98.5
mj cal 4	4	✓	10.0	9.9	98.7
mj cal 3	5	✓	25.0	24.4	97.8
mj cal 2	6	✓	50.0	50.2	100.4
mj qc1	7	✓	100.0	100.4	100.4

\* compound not evaluated in this run, ratios were out of acceptable range and poor peak shape

*BW*  
12/11/19

# AM #27 Cannabinoids


BW

**Batch results** D:\MassHunter\Data\2019\am27-28 120919\QuantResults\cann quant.batch.bin  
**Calibration Last Update** 12/11/2019 11:49:39 AM

**Instrument** 69679  
**Type** QC  
**Acq. Method** AM 27 THC quant.m  
**Sample Position** P3-A1  
**Injection Volume** 10  
**Acq. Date-Time** 12/9/2019 11:45:27 AM  
**Sample Info.**

**Data File** mj cal 7.d  
**Sample** mj cal 7 qc1  
**Operator** Brittany Wylie  
**Comment** BW 12/11/19

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-COOH	1.040	39250	169.3	31.1	947.7	590868	6.053 ng/ml <b>Low</b>
THC	1.935	4179	191.5	45.7 <b>High</b>	∞	673895	1.508 ng/ml <b>Low</b>




# AM #27 Cannabinoids

BW

**Batch results** D:\MassHunter\Data\2019\am27-28 120919\QuantResults\cann quant.batch.bin  
**Calibration Last Update** 12/11/2019 11:49:39 AM

<b>Instrument</b>	69679	<b>Data File</b>	mj cal 6.d
<b>Type</b>	Cal	<b>Sample</b>	mj cal-6 cal 2 BW
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Britany Wylie 12/11/19
<b>Sample Position</b>	P3-B1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	12/9/2019 11:53:11 AM		
<b>Sample Info.</b>			

**Sample Chromatogram**



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-COOH	1.040	69947	223.2	34.8	591.9	597212	10.337 ng/ml
THC-OH	<del>1.044</del>	<del>134797</del>	∞	<del>6.3</del> <b>Low</b>	∞	<del>1662982</del>	<del>3.124</del> ng/ml NE
THC	1.935	13465	649.1	28.1	240.8	681864	3.201 ng/ml


# AM #27 Cannabinoids

BW

**Batch results** D:\MassHunter\Data\2019\am27-28 120919\QuantResults\cann quant.batch.bin  
**Calibration Last Update** 12/11/2019 11:49:39 AM

<b>Instrument</b>	69679	<b>Data File</b>	mj cal 5.d
<b>Type</b>	Cal	<b>Sample</b>	mj cal 5 cal 3 BW
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Britany Wylie 12/11/19
<b>Sample Position</b>	P3-C1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	12/9/2019 12:00:54 PM		

**Sample Chromatogram**



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-COOH	1.040	133176	701.5	37.7	1872.9	592577	19.431 ng/ml
THC-OH	<del>1.044</del>	<del>173162</del>	$\infty$	<del>8.1</del>	$\infty$	<del>1625445</del>	<del>4.926</del> ng/ml
THC	1.935	22200	1015.8	25.7	169.0	646396	5.025 ng/ml

NE

# AM #27 Cannabinoids


BW

**Batch results** D:\MassHunter\Data\2019\am27-28 120919\QuantResults\cann quant.batch.bin  
**Calibration Last Update** 12/11/2019 11:49:39 AM

<b>Instrument</b>	69679	<b>Data File</b>	mj cal 4.d
<b>Type</b>	Cal	<b>Sample</b>	mj cal 4
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Britany Wylie
<b>Sample Position</b>	P3-D1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	12/9/2019 12:08:36 PM		

**Sample Info.**

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-COOH	1.040	335928	1470.8	37.9	3739.0	574812	49.824 ng/ml
<del>THC-OH</del>	<del>1.044</del>	<del>286608</del>	<del>∞</del>	<del>8.4</del>	<del>∞</del>	<del>1623713</del>	<del>9.874 ng/ml</del> NE
THC	1.950	47466	1308.9	26.9	∞	666798	9.629 ng/ml

# AM #27 Cannabinoids


BW

**Batch results** D:\MassHunter\Data\2019\am27-28 120919\QuantResults\cann quant.batch.bin  
**Calibration Last Update** 12/11/2019 11:49:39 AM

**Instrument** 69679  
**Type** Cal  
**Acq. Method** AM 27 THC quant.m  
**Sample Position** P3-E1  
**Injection Volume** 10  
**Acq. Date-Time** 12/9/2019 12:16:19 PM  
**Sample Info.**

**Data File** mj cal 3.d  
**Sample** mj cal 3 cal 5 BW  
**Operator** Britany Wylie  
**Comment** 12/11/19

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-COOH	1.040	491192	4682.3	37.9	1888.8	552728	75.534 ng/ml
THC-OH	1.044	601987	∞	11.3	∞	1573842	24.440 ng/ml
THC	1.935	122265	887.3	26.7	6125.7	663620	23.757 ng/ml

NE


# AM #27 Cannabinoids

BW

**Batch results** D:\MassHunter\Data\2019\am27-28 120919\QuantResults\cann quant.batch.bin  
**Calibration Last Update** 12/11/2019 11:49:39 AM

<b>Instrument</b>	69679	<b>Data File</b>	mj cal 2.d
<b>Type</b>	Cal	<b>Sample</b>	mj cal 2 cal 6 BW
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Britany Wylie 12/11/19
<b>Sample Position</b>	P3-F1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	12/9/2019 12:24:01 PM		
<b>Sample Info.</b>			

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-COOH	1.040	649850	3375.5	39.0	9488.1	559324	98.618 ng/ml
<del>THC-OH</del>	<del>1.044</del>	<del>1172252</del>	<del>∞</del>	<del>11.3</del>	<del>∞</del>	<del>1569900</del>	<del>50.194 ng/ml</del> NE
THC	1.935	254150	4402.9	24.6	∞	643564	50.083 ng/ml


# AM #27 Cannabinoids

BW

**Batch results** D:\MassHunter\Data\2019\am27-28 120919\QuantResults\cann quant.batch.bin  
**Calibration Last Update** 12/11/2019 11:49:39 AM

<b>Instrument</b>	69679	<b>Data File</b>	mj qc1.d
<b>Type</b>	Cal	<b>Sample</b>	mj qc1 cal 7 BW
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Britany Wylie 12/11/19
<b>Sample Position</b>	P3-G1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	12/9/2019 12:31:43 PM		
<b>Sample Info.</b>			

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-COOH	1.040	1660714	10779.9	39.5	7663.6	559506	251.257 ng/ml
THC-OH	<del>1.044</del>	<del>2434297</del>	$\infty$	<del>11.9</del>	$\infty$	<del>1670437</del>	<del>100.441</del> ng/ml
THC	1.935	530475	6947.3	25.3	$\infty$	659143	101.304 ng/ml