

REVIEWED

By Sarah Pickle at 10:57 am, Jan 08, 2019

TS

12/28/2018

Worklist: 2863

<u>LAB CASE</u>	<u>ITEM</u>	<u>TASK ID</u>	<u>DESCRIPTION</u>
P2018-3029	1	134329	AM 28 Blood Multi-Drug Quant Panel 2 by LC-Q
P2018-3284	1	135620	AM 28 Blood Multi-Drug Quant Panel 2 by LC-Q



AM# 28: Multi-Drug Quantitation in Blood by LC-MS/MS

TS

Extraction Date: 01/02/18
Plate lot#: 0530563

Analyst: Tamara Salazar
Plate Expiration: 05/21/2019

Mobile phase A: 5mM Amm Form + 0.01% FA
0.5M Ammonium Hydroxide

Mobile phase B: 0.01% Formic Acid in MeOH
Ethyl Acetate 20% Methanol in Water

Blank Blood Lot: 445283-1
LCMS-QQQ ID: 59740

Column: Agilent 120 EC-C18 (2.1x 100-4um)

Pre-Analytic:

- 1. Check levels of mobile phases and needle wash refill as needed. Ensure waste is not full.
- 2. Ensure correct column is installed and begin mobile phase flow allow to equilibrate ~ 30 minutes.
- 3. Create worklist: _____

Analytic:

- 1. Remove standards, plate, controls, and samples from cold storage. Allow to reach room temperature.
- 2. Pipette **250µL blood (calibrated pipette) Pipette ID: 3** in wells of analytical (standards) plate.
- 3. Place on shaking incubator at ambient temp., 900rpm for 15 minutes. *Shaker ID: 067105*
- 4. Pipette **250µL 0.5M ammonium hydroxide** in wells of analytical plate.
- 5. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 6. Transfer **300µL of blood+base** mixture to corresponding wells of SLE+ plate.
- 7. Apply positive pressure for approx. 10-15 seconds (or until no liquid remains on top of sorbent).
(Load at 85-100 PSI- Selector to the right) Manifold ID: 067104
- 8. Wait 5 minutes.
- 9. Add **900uL ethyl acetate.**
- 10. Wait 5 minutes.
- 11. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 12. Add **900uL ethyl acetate.**
- 13. Wait 5 minutes.
- 14. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 15. Remove plate containing eluate. Place on SPE Dry and evaporate to dryness at approx. 35°C.
SPE Dry ID: 067103
- 16. Reconstitute in **100µL 20% MeOH** and heat seal plate with foil. Place in autosampler and run worklist.

Post-Analytic

- 1. Create batch and process data.
Worklist path: *C:\MassHunter\Data\2019\AM 28\010219 MDQ P1 P2 wklst 2862 2863 TS*
Batch Name: *P2a MDQ TS*
- 2. Make necessary changes to integration limits
- 3. Integration linear and R² values ≥0.98 for each analyte.
- 4. For unknown samples and controls: response ratio within 20% of average of controls and standards, RT within +/- 5% (tramadol RT +/-2%), S/N for primary transition >10 and secondary transitions >5.
- 5. Did all QCs pass for each analyte? Y / N _____ Add Control data to QC tracking spreadsheet.
- 6. Central File Packet to include: LIMS Worklist, Method Checklist, Calibration and Control Reports.

COMMENTS:

Curves limited: Nortriptyline 10-1000

Compounds evaluated: Chlordiazepoxide, Zopiclone, Nortriptyline, Pseudophedrine, Amitriptyline



Idaho State Police Forensic Services

TS

AM #28 Blood Multi-Drug Confirmatory Analysis by LCMS-QQQ

Analyst: Tamara Salazar
Extraction Date: 01/02/18
Worklist Number: 2863

<i>Reagent</i>	<i>Lot Number</i>	<i>Expiration Date</i>	<i>Date in Service</i>	<i>Date Out of Service</i>	<i>Initials</i>
ToxBox Benzo/Opiate Quant Plate	0530563	05/21/19	06/05/18		
Negative Blood	445283-1		10/25/18		
Methanol External Control Solution	022018	02/20/19	02/20/18		
Blood External Control Solution	WS010219	02/20/19	01/02/19		
0.5M Ammonium Hydroxide	082818	02/28/19	08/28/18		
Formic Acid (LCMS Grade)	095180B		11/15/17		
Water (LCMS Grade)	183905		12/14/18		
Methanol (LCMS Grade)	184782		12/14/18		
2M Ammonium Formate	121418		12/14/18		
Ammonium Formate (LCMS Grade)	A0287238		01/30/17		
Ammonium Hydroxide (ACS grade or better)	053200		08/03/07		
Ethyl Acetate (ACS grade or better)	100362		06/26/17		
20% LCMS Methanol in LCMS Water	082818		08/28/18		
5mM Ammonium Formate + 0.01% Formic Acid in Water (Mobile Phase A)	010318		01/03/18		*Due to an instrument error, the panel 1 analysis took more mobile phase than expected. A fresh lot of mobile phase A was prepared for the panel 2 analysis.
0.01% Formic Acid in Methanol (Mobile Phase B)	121418		12/14/18		
Needle Rinse--75% LCMS MeOH in LCMS Water	122618		12/26/18		

Methanol External Control Solution (Lot: 022018)
100 ul each 1 mg/mL stock solution in 9600 ul MeOH

<i>Component</i>	<i>Source</i>	<i>Source Lot Number</i>	<i>Expiration Date</i>
Methanol (LCMS)	Fisher	172516	
Zopiclone	Cerilliant	FE08291601	12/31/2021
Chlordiazepoxide	Cerilliant	FE07241502	08/31/2020
Flunitrazepam	Cerilliant	FE08051602	08/31/2021
Prepared:	02/20/18		
Prepared By:	Tamara Salazar		
Expires:	02/20/19		

Blood External Control Solution (Lot: WS010219)

50 ul of methanol external control solution was added to 9950 ul of blood.

<i>Component</i>	<i>Source</i>	<i>Source Lot Number</i>
Negative Blood	Hemostat	445283-1
Methanol External Control Solution		022018
Prepared:	01/02/19	
Prepared by:	Tamara Salazar	
Expires:	02/20/19	

0.5M Ammonium Hydroxide (Lot: 082818)

<i>Component</i>	<i>Source</i>	<i>Source Lot Number</i>
Ammonium Hydroxide (ACS grade or better)	Fisher	053200
DI Water	-	-
Prepared:	08/28/18	
Prepared By:	Tamara Salazar	
Expires:	02/28/19	

2M Ammonium Formate Solution (Lot: 121418)

<i>Component</i>	<i>Source</i>	<i>Source Lot Number</i>
Ammonium Formate	Arcos Organic	A0287238
LCMS Water	Fisher	183905
Prepared:	12/14/18	
Prepared By:	Tamara Salazar	
Expires:	Remake as needed	

20% LCMS Methanol in LCMS Water (Lot: 082818)

<i>Component</i>	<i>Source</i>	<i>Source Lot Number</i>
MeOH (LCMS Grade)	Fisher	177145
Water (LCMS Grade)	Fisher	182702
Prepared:	08/28/18	
Prepared By:	Tamara Salazar	

5mM Ammonium Formate + 0.01% Formic Acid in Water (Lot: 010319)

<i>Component</i>	<i>Source</i>	<i>Source Lot Number</i>
2M Ammonium Formate		121418
LCMS Water	Fisher	183905
Formic Acid	Fisher	095180B
Prepared:	01/03/19	
Prepared By:	Tamara Salazar	

0.01% Formic Acid in LCMS Methanol (Mobile Phase B) (Lot: 121418)

<i>Component</i>	<i>Source</i>	<i>Source Lot Number</i>
Formic Acid	Fisher	095180B
MeOH (LCMS Grade)	Fisher	184782
Prepared:	12/14/18	
Prepared By:	Tamara Salazar	

5

AM #28: Blood Multi-Drug Confirmatory Analysis by LCMS-QQQ

Needle Rinse (75% LCMS MeOH in LCMS Water) (Lot: 122618)

<i>Component</i>	<i>Source</i>	<i>Source Lot Number</i>
MeOH (LCMS Grade)	Fisher	177145
Water (LCMS Grade)	Fisher	183905
Prepared:	12/26/18	
Prepared By:	Tamara Salazar	

TS

ISP FORENSICS - Pocatello Instrument # 59740


Multi-Drug Confirmatory Analysis Report

Batch Data Path C:\MassHunter\Data\2019\AM 28\010219 MDQ P1 P2 wklst 2862 2863 TS\QuantResults\P2a MDQ
Analysis Time 1/7/2019 2:14 PM **Analyst Name** ISPUser
Report Time 1/7/2019 2:16 PM **Reporter Name** ISPUser
Last Calib Update 1/7/2019 2:14 PM **Batch State** Processed

Analysis Info

Acq Time 2019-01-04 09:19 **Data File** P2 a Negative control.d
Sample Type Sample **Sample Name** P2 a Negative control.d
Dilution 1 **Acq Method** MDQ Panel 2a 06-22-18.m
Position P1-E6 **Sample Info**
Inj Vol -1 **Comment** AM 28-Panel 2

Sample Chromatogram



TS

ISP FORENSICS - Pocatello Instrument # 59740


Multi-Drug Confirmatory Analysis Report

Batch Data Path C:\MassHunter\Data\2019\AM 28\010219 MDQ P1 P2 wk1st 2862 2863 TS\QuantResults\P2a MDQ TS.batch.bii
Analysis Time 1/7/2019 2:14 PM **Analyst Name** ISPUser
Report Time 1/7/2019 2:16 PM **Reporter Name** ISPUser
Last Calib Update 1/7/2019 2:14 PM **Batch State** Processed

Analysis Info

Acq Time 2019-01-04 09:32 **Data File** P2 a External control-50ng mL.d
Sample Type Sample **Sample Name** P2 a External control-50ng mL.d
Dilution 1 **Acq Method** MDQ Panel 2a 06-22-18.m
Position P1-F6 **Sample Info**
Inj Vol -1 **Comment** AM 28; 445283-1 + WS022018 50 ng flunitrazepam, zopiclone, chlorodia:

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
Zopiclone	Zopiclone-D4	6.130	166343	136617	1.2176	40.1025
Chlordiazepoxide	Chlordiazepoxide-D5	7.667	651734	3237600	0.2013	47.4615

TS

ISP FORENSICS - Pocatello Instrument # 59740


Multi-Drug Confirmatory Analysis Report

Batch Data Path C:\MassHunter\Data\2019\AM 28\010219 MDQ P1 P2 wklst 2862 2863 TS\QuantResults\P2a MDQ TS.batch.bi
Analysis Time 1/7/2019 2:14 PM **Analyst Name** ISPUser
Report Time 1/7/2019 2:16 PM **Reporter Name** ISPUser
Last Calib Update 1/7/2019 2:14 PM **Batch State** Processed

Analysis Info

Acq Time 2019-01-04 09:57 **Data File** P2 a QC 10 (QC1).d
Sample Type QC **Sample Name** P2 a QC 10 (QC1).d
Dilution 1 **Acq Method** MDQ Panel 2a 06-22-18.m
Position P1-A6 **Sample Info**
Inj Vol -1 **Comment** AM 28-Panel 2

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
Pseudoephedrine	Pseudoephedrine-D3	3.715	649898	2887723	0.2251	10.2464
Zopiclone	Zopiclone-D4	6.130	36086	143110	0.2522	9.6763
Amitriptyline	Amitriptyline-D3	7.509	23088	54375	0.4246	8.8444
Nortriptyline	Nortriptyline-D3	7.533	6452	14778	0.4366	10.5965
Chlordiazepoxide	Chlordiazepoxide-D5	7.667	126882	3019086	0.0420	10.0724

15

ISP FORENSICS - Pocatello Instrument # 59740


Multi-Drug Confirmatory Analysis Report

Batch Data Path C:\MassHunter\Data\2019\AM 28\010219 MDQ P1 P2 wk1st 2862 2863 TS\QuantResults\P2a MDQ TS.batch.bii
Analysis Time 1/7/2019 2:14 PM **Analyst Name** ISPUser
Report Time 1/7/2019 2:16 PM **Reporter Name** ISPUser
Last Calib Update 1/7/2019 2:14 PM **Batch State** Processed

Analysis Info

Acq Time 2019-01-04 10:48 **Data File** P2 a QC 100 (QC2).d
Sample Type QC **Sample Name** P2 a QC 100 (QC2).d
Dilution 1 **Acq Method** MDQ Panel 2a 06-22-18.m
Position P1-B6 **Sample Info**
Inj Vol -1 **Comment** AM 28-Panel 2

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
Pseudoephedrine	Pseudoephedrine-D3	3.715	7934092	3182613	2.4929	96.7170
Zopiclone	Zopiclone-D4	6.130	310380	95600	3.2466	104.0497
Amitriptyline	Amitriptyline-D3	7.509	186946	41018	4.5577	102.3223
Nortriptyline	Nortriptyline-D3	7.533	57276	13453	4.2574	90.2096
Chlordiazepoxide	Chlordiazepoxide-D5	7.667	1148091	2606091	0.4405	103.6222

TS

ISP FORENSICS - Pocatello Instrument # 59740


Multi-Drug Confirmatory Analysis Report

Batch Data Path C:\MassHunter\Data\2019\AM 28\010219 MDQ P1 P2 wk1st 2862 2863 TS\QuantResults\P2a MDQ TS.batch.bii
Analysis Time 1/7/2019 2:14 PM **Analyst Name** ISPUser
Report Time 1/7/2019 2:17 PM **Reporter Name** ISPUser
Last Calib Update 1/7/2019 2:14 PM **Batch State** Processed

Analysis Info

Acq Time 2019-01-04 11:38 **Data File** P2 a QC 250 (QC 3).d
Sample Type QC **Sample Name** P2 a QC 250 (QC 3).d
Dilution 1 **Acq Method** MDQ Panel 2a 06-22-18.m
Position P1-C6 **Sample Info**
Inj Vol -1 **Comment** AM 28-Panel 2

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
Pseudoephedrine	Pseudoephedrine-D3	3.695	21068661	3212612	6.5581	251.7143
Zopiclone	Zopiclone-D4	6.130	713742	86209	8.2792	262.6541
Amitriptyline	Amitriptyline-D3	7.509	510843	43874	11.6435	262.5818
Nortriptyline	Nortriptyline-D3	7.533	144794	12425	11.6532	244.3187
Chlordiazepoxide	Chlordiazepoxide-D5	7.667	2034568	1989331	1.0227	240.2910

ISP FORENSICS - Pocatello Instrument # 59740

Multi-Drug Confirmatory Analysis Report


TS

Batch Data Path C:\MassHunter\Data\2019\AM 28\010219 MDQ P1 P2 wklst 2862 2863 TS\QuantResults\P2a MDQ TS.batch.bii
Analysis Time 1/7/2019 2:14 PM **Analyst Name** ISPUser
Report Time 1/7/2019 2:17 PM **Reporter Name** ISPUser
Last Calib Update 1/7/2019 2:14 PM **Batch State** Processed

Analysis Info

Acq Time 2019-01-04 12:04 **Data File** P2 a QC 1000 (QC 4).d
Sample Type QC **Sample Name** P2 a QC 1000 (QC 4).d
Dilution 1 **Acq Method** MDQ Panel 2a 06-22-18.m
Position P1-D6 **Sample Info**
Inj Vol -1 **Comment** AM 28-Panel 2

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
Pseudoephedrine	Pseudoephedrine-D3	3.695	70610957	2729206	25.8723	988.1311
Zopiclone	Zopiclone-D4	6.110	1257964	40556	31.0176	979.2744
Amitriptyline	Amitriptyline-D3	7.509	2681426	64788	41.3876	935.3074
Nortriptyline	Nortriptyline-D3	7.533	796635	14530	54.8278	1143.9521
Chlordiazepoxide	Chlordiazepoxide-D5	7.687	3484276	869892	4.0054	940.4617


ISP Forensics Calibration Curve Report

TS

Batch Data Path C:\MassHunter\Data\2019\AM 28\010219 MDQ P1 P2 wklst 2862 2863 TS\QuantResults\P2a MDQ TS.batch.bin

Last Calib Update 1/7/2019 2:14 PM **Analyst Name** ISP TOX

Target Compound *Pseudoephedrine*
Internal Standard *Pseudoephedrine-D3*




Sample	Level	Enabled	Exp Conc	Final Conc	Accuracy
P2 a cal 1-5ng.d	1	<input checked="" type="checkbox"/>	5	5.8	115.6
P2 a cal 2-10ng.d	2	<input checked="" type="checkbox"/>	10	10.1	101.2
P2 a QC 10 (QC1).d	2	<input checked="" type="checkbox"/>	10	10.2	102.5
P2 a cal 3-25ng.d	3	<input checked="" type="checkbox"/>	25	23.5	93.8
P2 a cal 4-50ng.d	4	<input checked="" type="checkbox"/>	50	46.5	93.1
P2 a cal 5-100ng.d	5	<input checked="" type="checkbox"/>	100	94.4	94.4
P2 a QC 100 (QC2).d	5	<input checked="" type="checkbox"/>	100	96.7	96.7
P2 a cal 6-250ng.d	6	<input checked="" type="checkbox"/>	250	246.0	98.4
P2 a QC 250 (QC 3).d	6	<input checked="" type="checkbox"/>	250	251.7	100.7
P2 a cal 7-500ng.d	7	<input checked="" type="checkbox"/>	500	521.0	104.2
P2 a cal 8-1000ng.d	8	<input checked="" type="checkbox"/>	1000	992.6	99.3
P2 a QC 1000 (QC 4).d	8	<input checked="" type="checkbox"/>	1000	988.1	98.8

ISP Forensics Calibration Curve Report

TS

Batch Data Path C:\MassHunter\Data\2019\AM 28\010219 MDQ P1 P2 wk1st 2862 2863 TS\QuantResults\P2a MDQ TS.batch.bin
Last Calib Update 1/7/2019 2:14 PM **Analyst Name** ISP TOX

Target Compound Zopiclone
Internal Standard Zopiclone-D4




Sample	Level	Enabled	Exp Conc	Final Conc	Accuracy
P2 a cal 1-5ng.d	1	<input checked="" type="checkbox"/>	5	5.5	110.0
P2 a cal 2-10ng.d	2	<input checked="" type="checkbox"/>	10	10.1	100.9
P2 a QC 10 (QC1).d	2	<input checked="" type="checkbox"/>	10	9.7	96.8
P2 a cal 3-25ng.d	3	<input checked="" type="checkbox"/>	25	23.5	94.1
P2 a cal 4-50ng.d	4	<input checked="" type="checkbox"/>	50	46.4	92.9
P2 a cal 5-100ng.d	5	<input checked="" type="checkbox"/>	100	100.2	100.2
P2 a QC 100 (QC2).d	5	<input checked="" type="checkbox"/>	100	104.0	104.0
P2 a cal 6-250ng.d	6	<input checked="" type="checkbox"/>	250	245.5	98.2
P2 a QC 250 (QC 3).d	6	<input checked="" type="checkbox"/>	250	262.7	105.1
P2 a cal 7-500ng.d	7	<input checked="" type="checkbox"/>	500	528.9	105.8
P2 a cal 8-1000ng.d	8	<input checked="" type="checkbox"/>	1000	979.8	98.0
P2 a QC 1000 (QC 4).d	8	<input checked="" type="checkbox"/>	1000	979.3	97.9

ISP Forensics Calibration Curve Report

B

Batch Data Path C:\MassHunter\Data\2019\AM 28\010219 MDQ P1 P2 wklst 2862 2863 TS\QuantResults\P2a MDQ TS.batch.bin
Last Calib Update 1/7/2019 2:14 PM **Analyst Name** ISP TOX

Target Compound Amitriptyline
Internal Standard Amitriptyline-D3




Sample	Level	Enabled	Exp Conc	Final Conc	Accuracy
P2 a cal 1-5ng.d	1	<input checked="" type="checkbox"/>	5	5.1	101.9
P2 a cal 2-10ng.d	2	<input checked="" type="checkbox"/>	10	10.8	107.6
P2 a QC 10 (QC1).d	2	<input checked="" type="checkbox"/>	10	8.8	88.4
P2 a cal 3-25ng.d	3	<input checked="" type="checkbox"/>	25	23.6	94.4
P2 a cal 4-50ng.d	4	<input checked="" type="checkbox"/>	50	44.4	88.8
P2 a cal 5-100ng.d	5	<input checked="" type="checkbox"/>	100	101.9	101.9
P2 a QC 100 (QC2).d	5	<input checked="" type="checkbox"/>	100	102.3	102.3
P2 a cal 6-250ng.d	6	<input checked="" type="checkbox"/>	250	265.5	106.2
P2 a QC 250 (QC 3).d	6	<input checked="" type="checkbox"/>	250	262.6	105.0
P2 a cal 7-500ng.d	7	<input checked="" type="checkbox"/>	500	502.7	100.5
P2 a cal 8-1000ng.d	8	<input checked="" type="checkbox"/>	1000	986.1	98.6
P2 a QC 1000 (QC 4).d	8	<input checked="" type="checkbox"/>	1000	935.3	93.5

ISP Forensics Calibration Curve Report

15

Batch Data Path C:\MassHunter\Data\2019\AM 28\010219 MDQ P1 P2 wk1st 2862 2863 TS\QuantResults\P2a MDQ TS.batch.bin
Last Calib Update 1/7/2019 2:14 PM **Analyst Name** ISP TOX

Target Compound Nortriptyline
Internal Standard Nortriptyline-D3




Sample	Level	Enabled	Exp Conc	Final Conc	Accuracy
P2 a cal 1-5ng.d	1	<input type="checkbox"/>	5	7.3	145.7
P2 a cal 2-10ng.d	2	<input checked="" type="checkbox"/>	10	12.0	120.4
P2 a QC 10 (QC1).d	2	<input checked="" type="checkbox"/>	10	10.6	106.0
P2 a cal 3-25ng.d	3	<input checked="" type="checkbox"/>	25	24.4	97.7
P2 a cal 4-50ng.d	4	<input checked="" type="checkbox"/>	50	47.1	94.1
P2 a cal 5-100ng.d	5	<input checked="" type="checkbox"/>	100	88.8	88.8
P2 a QC 100 (QC2).d	5	<input checked="" type="checkbox"/>	100	90.2	90.2
P2 a cal 6-250ng.d	6	<input checked="" type="checkbox"/>	250	248.5	99.4
P2 a QC 250 (QC 3).d	6	<input checked="" type="checkbox"/>	250	244.3	97.7
P2 a cal 7-500ng.d	7	<input checked="" type="checkbox"/>	500	481.7	96.3
P2 a cal 8-1000ng.d	8	<input checked="" type="checkbox"/>	1000	1032.5	103.2
P2 a QC 1000 (QC 4).d	8	<input checked="" type="checkbox"/>	1000	1144.0	114.4

ISP Forensics Calibration Curve Report

15

Batch Data Path C:\MassHunter\Data\2019\AM 28\010219 MDQ P1 P2 wklst 2862 2863 TS\QuantResults\P2a MDQ TS.batch.bin
Last Calib Update 1/7/2019 2:14 PM **Analyst Name** ISP TOX

Target Compound *Chlordiazepoxide*
Internal Standard *Chlordiazepoxide-D5*



Sample	Level	Enabled	Exp Conc	Final Conc	Accuracy
P2 a cal 1-5ng.d	1	<input checked="" type="checkbox"/>	5	5.2	104.5
P2 a cal 2-10ng.d	2	<input checked="" type="checkbox"/>	10	9.9	98.8
P2 a QC 10 (QC1).d	2	<input checked="" type="checkbox"/>	10	10.1	100.7
P2 a cal 3-25ng.d	3	<input checked="" type="checkbox"/>	25	25.7	102.7
P2 a cal 4-50ng.d	4	<input checked="" type="checkbox"/>	50	51.8	103.5
P2 a cal 5-100ng.d	5	<input checked="" type="checkbox"/>	100	99.2	99.2
P2 a QC 100 (QC2).d	5	<input checked="" type="checkbox"/>	100	103.6	103.6
P2 a cal 6-250ng.d	6	<input checked="" type="checkbox"/>	250	221.2	88.5
P2 a QC 250 (QC 3).d	6	<input checked="" type="checkbox"/>	250	240.3	96.1
P2 a cal 7-500ng.d	7	<input checked="" type="checkbox"/>	500	501.8	100.4
P2 a cal 8-1000ng.d	8	<input checked="" type="checkbox"/>	1000	1025.3	102.5
P2 a QC 1000 (QC 4).d	8	<input checked="" type="checkbox"/>	1000	940.5	94.0

15

ISP FORENSICS - Pocatello Instrument # 59740


Multi-Drug Confirmatory Analysis Report

Batch Data Path C:\MassHunter\Data\2019\AM 28\010219 MDQ P1 P2 wk1st 2862 2863 TS\QuantResults\P2a MDQ TS.batch.bii
Analysis Time 1/7/2019 2:14 PM **Analyst Name** ISPUser
Report Time 1/7/2019 2:15 PM **Reporter Name** ISPUser
Last Calib Update 1/7/2019 2:14 PM **Batch State** Processed

Analysis Info

Acq Time 2019-01-03 18:41 **Data File** P2 a cal 1-5ng.d
Sample Type Calibration **Sample Name** P2 a cal 1-5ng.d
Dilution 1 **Acq Method** MDQ Panel 2a 06-22-18.m
Position P1-A5 **Sample Info**
Inj Vol -1 **Comment** AM 28-Panel 2

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
Pseudoephedrine	Pseudoephedrine-D3	3.715	287918	2668628	0.1079	5.7790
Zopiclone	Zopiclone-D4	6.130	21300	178094	0.1196	5.4986
Amitriptyline	Amitriptyline-D3	7.509	12312	47549	0.2589	5.0972
Nortriptyline	Nortriptyline-D3	7.533	4454	16033	0.2778	7.2871
Chlordiazepoxide	Chlordiazepoxide-D5	7.667	60630	2837272	0.0214	5.2230

15

ISP FORENSICS - Pocatello Instrument # 59740


Multi-Drug Confirmatory Analysis Report

Batch Data Path C:\MassHunter\Data\2019\AM 28\010219 MDQ P1 P2 wk1st 2862 2863 TS\QuantResults\P2a MDQ TS.batch.bi
Analysis Time 1/7/2019 2:14 PM **Analyst Name** ISPUser
Report Time 1/7/2019 2:15 PM **Reporter Name** ISPUser
Last Calib Update 1/7/2019 2:14 PM **Batch State** Processed

Analysis Info

Acq Time 2019-01-03 18:54 **Data File** P2 a cal 2-10ng.d
Sample Type Calibration **Sample Name** P2 a cal 2-10ng.d
Dilution 1 **Acq Method** MDQ Panel 2a 06-22-18.m
Position P1-B5 **Sample Info**
Inj Vol -1 **Comment** AM 28-Panel 2

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
Pseudoephedrine	Pseudoephedrine-D3	3.715	688660	3105400	0.2218	10.1208
Zopiclone	Zopiclone-D4	6.130	51031	192418	0.2652	10.0876
Amitriptyline	Amitriptyline-D3	7.509	22002	43182	0.5095	10.7647
Nortriptyline	Nortriptyline-D3	7.533	5860	11588	0.5057	12.0358
Chlordiazepoxide	Chlordiazepoxide-D5	7.687	134847	3272482	0.0412	9.8798

B

ISP FORENSICS - Pocatello Instrument # 59740


Multi-Drug Confirmatory Analysis Report

Batch Data Path C:\MassHunter\Data\2019\AM 28\010219 MDQ P1 P2 wklst 2862 2863 TS\QuantResults\P2a MDQ TS.batch.bii
Analysis Time 1/7/2019 2:14 PM **Analyst Name** ISPUser
Report Time 1/7/2019 2:15 PM **Reporter Name** ISPUser
Last Calib Update 1/7/2019 2:14 PM **Batch State** Processed

Analysis Info

Acq Time 2019-01-03 19:07 **Data File** P2 a cal 3-25ng.d
Sample Type Calibration **Sample Name** P2 a cal 3-25ng.d
Dilution 1 **Acq Method** MDQ Panel 2a 06-22-18.m
Position P1-C5 **Sample Info**
Inj Vol -1 **Comment** AM 28-Panel 2

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
Pseudoephedrine	Pseudoephedrine-D3	3.715	1673569	2928694	0.5714	23.4533
Zopiclone	Zopiclone-D4	6.130	125671	181737	0.6915	23.5225
Amitriptyline	Amitriptyline-D3	7.509	44881	41670	1.0771	23.6010
Nortriptyline	Nortriptyline-D3	7.533	13549	12310	1.1006	24.4329
Chlordiazepoxide	Chlordiazepoxide-D5	7.667	335577	3093607	0.1085	25.6707

15

ISP FORENSICS - Pocatello Instrument # 59740


Multi-Drug Confirmatory Analysis Report

Batch Data Path C:\MassHunter\Data\2019\AM 28\010219 MDQ P1 P2 wk1st 2862 2863 TS\QuantResults\P2a MDQ TS.batch.bi
Analysis Time 1/7/2019 2:14 PM **Analyst Name** ISPUser
Report Time 1/7/2019 2:15 PM **Reporter Name** ISPUser
Last Calib Update 1/7/2019 2:14 PM **Batch State** Processed

Analysis Info

Acq Time 2019-01-03 19:19 **Data File** P2 a cal 4-50ng.d
Sample Type Calibration **Sample Name** P2 a cal 4-50ng.d
Dilution 1 **Acq Method** MDQ Panel 2a 06-22-18.m
Position P1-D5 **Sample Info**
Inj Vol -1 **Comment** AM 28-Panel 2

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
Pseudoephedrine	Pseudoephedrine-D3	3.715	3428337	2912487	1.1771	46.5467
Zopiclone	Zopiclone-D4	6.130	256281	180638	1.4188	46.4424
Amitriptyline	Amitriptyline-D3	7.509	77612	38866	1.9969	44.4052
Nortriptyline	Nortriptyline-D3	7.533	21445	9808	2.1864	47.0579
Chlordiazepoxide	Chlordiazepoxide-D5	7.667	632730	2880862	0.2196	51.7646

15

ISP FORENSICS - Pocatello Instrument # 59740


Multi-Drug Confirmatory Analysis Report

Batch Data Path C:\MassHunter\Data\2019\AM 28\010219 MDQ P1 P2 wk1st 2862 2863 TS\QuantResults\P2a MDQ TS.batch.bi
Analysis Time 1/7/2019 2:14 PM **Analyst Name** ISPUser
Report Time 1/7/2019 2:15 PM **Reporter Name** ISPUser
Last Calib Update 1/7/2019 2:14 PM **Batch State** Processed

Analysis Info

Acq Time 2019-01-03 19:32 **Data File** P2 a cal 5-100ng.d
Sample Type Calibration **Sample Name** P2 a cal 5-100ng.d
Dilution 1 **Acq Method** MDQ Panel 2a 06-22-18.m
Position P1-E5 **Sample Info**
Inj Vol -1 **Comment** AM 28-Panel 2

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
Pseudoephedrine	Pseudoephedrine-D3	3.715	7293803	2998065	2.4328	94.4251
Zopiclone	Zopiclone-D4	6.130	497615	159265	3.1244	100.1985
Amitriptyline	Amitriptyline-D3	7.509	201485	44410	4.5369	101.8534
Nortriptyline	Nortriptyline-D3	7.533	54188	12930	4.1907	88.8213
Chlordiazepoxide	Chlordiazepoxide-D5	7.667	1098812	2606686	0.4215	99.1607

ISP FORENSICS - Pocatello Instrument # 59740


Multi-Drug Confirmatory Analysis Report

Batch Data Path C:\MassHunter\Data\2019\AM 28\010219 MDQ P1 P2 wk1st 2862 2863 TS\QuantResults\P2a MDQ TS.batch.bii
Analysis Time 1/7/2019 2:14 PM **Analyst Name** ISPUser
Report Time 1/7/2019 2:15 PM **Reporter Name** ISPUser
Last Calib Update 1/7/2019 2:14 PM **Batch State** Processed

Analysis Info

Acq Time 2019-01-03 19:45 **Data File** P2 a cal 6-250ng.d
Sample Type Calibration **Sample Name** P2 a cal 6-250ng.d
Dilution 1 **Acq Method** MDQ Panel 2a 06-22-18.m
Position P1-F5 **Sample Info**
Inj Vol -1 **Comment** AM 28-Panel 2

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
Pseudoephedrine	Pseudoephedrine-D3	3.695	18486773	2884307	6.4094	246.0455
Zopiclone	Zopiclone-D4	6.130	1055535	136440	7.7362	245.5428
Amitriptyline	Amitriptyline-D3	7.509	531010	45110	11.7714	265.4762
Nortriptyline	Nortriptyline-D3	7.533	165288	13946	11.8518	248.4552
Chlordiazepoxide	Chlordiazepoxide-D5	7.667	1799147	1911064	0.9414	221.2055

B

ISP FORENSICS - Pocatello Instrument # 59740


Multi-Drug Confirmatory Analysis Report

Batch Data Path C:\MassHunter\Data\2019\AM 28\010219 MDQ P1 P2 wk1st 2862 2863 TS\QuantResults\P2a MDQ TS.batch.bii
Analysis Time 1/7/2019 2:14 PM **Analyst Name** ISPUser
Report Time 1/7/2019 2:15 PM **Reporter Name** ISPUser
Last Calib Update 1/7/2019 2:14 PM **Batch State** Processed

Analysis Info

Acq Time 2019-01-04 08:29 **Data File** P2 a cal 7-500ng.d
Sample Type Calibration **Sample Name** P2 a cal 7-500ng.d
Dilution 1 **Acq Method** MDQ Panel 2a 06-22-18.m
Position P1-G5 **Sample Info**
Inj Vol -1 **Comment** AM 28-Panel 2

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
Pseudoephedrine	Pseudoephedrine-D3	3.695	40216564	2952707	13.6202	520.9804
Zopiclone	Zopiclone-D4	6.130	1261471	75413	16.7276	528.9118
Amitriptyline	Amitriptyline-D3	7.509	1505229	67619	22.2603	502.7043
Nortriptyline	Nortriptyline-D3	7.533	461431	20022	23.0459	481.7093
Chlordiazepoxide	Chlordiazepoxide-D5	7.687	2671558	1250377	2.1366	501.7657

B

ISP FORENSICS - Pocatello Instrument # 59740


Multi-Drug Confirmatory Analysis Report

Batch Data Path C:\MassHunter\Data\2019\AM 28\010219 MDQ P1 P2 wk1st 2862 2863 TS\QuantResults\P2a MDQ TS.batch.bii
Analysis Time 1/7/2019 2:14 PM **Analyst Name** ISPUser
Report Time 1/7/2019 2:15 PM **Reporter Name** ISPUser
Last Calib Update 1/7/2019 2:14 PM **Batch State** Processed

Analysis Info

Acq Time 2019-01-04 08:42 **Data File** P2 a cal 8-1000ng.d
Sample Type Calibration **Sample Name** P2 a cal 8-1000ng.d
Dilution 1 **Acq Method** MDQ Panel 2a 06-22-18.m
Position P1-H5 **Sample Info**
Inj Vol -1 **Comment** AM 28-Panel 2

Sample Chromatogram



Results

Compound	ISTD Compound	RT	Response	ISTD Resp	Resp Ratio	Final Conc
Pseudoephedrine	Pseudoephedrine-D3	3.695	67920916	2613264	25.9908	992.6490
Zopiclone	Zopiclone-D4	6.130	1585977	51104	31.0342	979.7958
Amitriptyline	Amitriptyline-D3	7.509	3522796	80736	43.6333	986.0980
Nortriptyline	Nortriptyline-D3	7.533	1081965	21867	49.4785	1032.4875
Chlordiazepoxide	Chlordiazepoxide-D5	7.687	3528659	808039	4.3669	1025.3299