


**REVIEWED**

By Tamara Salazar at 10:58 am, Jan 22, 2020

§ 1/21/2020

**Worklist: 3954**

<u>LAB CASE</u>	<u>ITEM</u>	<u>ITEM TYPE</u>	<u>DESCRIPTION</u>
M2019-5642	2	BCK	AM 27 Blood THC Quant by LC-QQQ
M2020-0040	2	BCK	AM 27 Blood THC Quant by LC-QQQ
P2019-3736	1	BCK	AM 27 Blood THC Quant by LC-QQQ
P2020-0003	1	BCK	AM 27 Blood THC Quant by LC-QQQ
P2020-0046	1	BCK	AM 27 Blood THC Quant by LC-QQQ
P2020-0068	1	BCK	AM 27 Blood THC Quant by LC-QQQ
P2020-0086	1	BCK	AM 27 Blood THC Quant by LC-QQQ
P2020-0102	1	BCK	AM 27 Blood THC Quant by LC-QQQ



**Worklist: 3956**

<u>LAB CASE</u>	<u>ITEM</u>	<u>ITEM TYPE</u>	<u>DESCRIPTION</u>
P2019-3618	2	UCK	AM 27 Urine Cannabinoids Confirmation by LC-QQQ



Sample was ran with  
worklist 3954

\$

# AM# 27: Quantitation of THC and Metabolites in Blood and Urine by LC-MS/MS

01/21/20 SP

Extraction Date: 01/16/20

Plate lot#: 190716

**Mobile phase A:** 0.1% Formic Acid in LCMS Water**Blank Blood Lot:** 445283-3**LCMS-QQQ ID:** 069901Analyst: ~~Celena Shrum~~ Sarah Pickle

Plate Expiration: 01/16/2020

**Mobile phase B:** 0.1% Formic acid in Acetonitrile**Column:** UCT Selectra DA 100 x 2.1mm 3um**Blank Urine Lot:** POC031319

## Pre-Analytic:

- 1. Check levels of mobile phases and needle wash refill as needed. Ensure waste is not full.
- 2. Ensure correct column is installed and begin mobile phase flow allow to equilibrate ~ 30 minutes.

## Analytic:

- 1. Remove standards, plate, controls, and samples from cold storage. Allow to reach room temperature.
- 2. Urine hydrolysis: add 1.5mL urine to blank plate, add 250µl 1N KOH. Shake and incubate at 40 degrees for 15 minutes. Using a calibrated pipette, add **1000µl blood and urine (if applicable) (calibrated pipette)** into the appropriate wells of analytical (standards) plate. **Pipette ID: #3382167**
- 3. Place on shaking incubator at ambient temp., 900rpm for 15 minutes. *Shaker ID: 067105*
- 4. Pipette **500µL 0.1% formic acid in water blood sample, 500 µL saturated phosphate buffer in urine** in wells of analytical plate.
- 5. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 6. Transfer **800µL of blood+acid or urine+acid** mixture to corresponding wells of SLE+ plate.
- 7. Apply positive pressure for approx. 10-15 seconds (or until no liquid remains on top of sorbent). **(Load at 85-100 PSI- Selector to the right)** Manifold ID: 067104
- 8. Wait 5 minutes.
- 9. Add **2.25mL MTBE. (Add in 3 increments of 750uL)**
- 10. Wait 5 minutes.
- 11. Apply positive pressure for approx. 15 seconds. **(10-15 PSI- Selector to the left).**
- 12. Add **2.25mL Hexane. (Add in 3 increments of 750uL)**
- 13. Wait 5 minutes.
- 14. Apply positive pressure for approx. 15 seconds. **(10-15 PSI- Selector to the left).**
- 15. Remove plate containing eluate. Place on SPE Dry and evaporate to dryness at approx. 35°C. *SPE Dry ID: 067103*
- 16. Reconstitute in **100µL 100% MeOH** and heat seal plate with foil. Place in autosampler and run worklist.

## Post-Analytic

- 1. Create batch and process data.
- 2. Make any necessary integration changes, Curve weighting of Linear 1/x with R<sup>2</sup> values ≥0.98 for each analyte
- 3. RT +/- 3% or 0.100 min, whichever is greater, +/- 20% Accuracy for greater than (+/- 30% for 10ng/ml or less). Ion ratios must be within +/- 20% of the averaged calibrators
- 4. Case sample response for THC and OH-THC 3ng/mL (quantitative), Carboxy-THC: 10ng/mL (qualitative only) will be reported. Samples with a THC or OH-THC response over 50 ng/mL will be reported out as greater than 50 ng/mL.
- 5. Did all QCs pass for each analyte? (if not, describe in comments section)
- 6. Enter QCs into control charting.
- 7. Central File Packet to include: LIMS Worklist, Method Checklist, Calibration and Control Reports

COMMENTS: Curve ranges are THC: 1-100, THC-COOH: 5-250, THC-OH: 3-100



# Idaho State Police Forensic Services

## AM #26 Urine THC and Metabolites Screen by LCMS-QQQ and AM #27 Quantitative Analysis of THC and Metabolites in Urine by LCMS-QQQ

**Urine External Control Solution (Lot: 011620)**

*8 µL of 100 µg/mL C-THC in 9.992 mL urine*

*Approximate concentration ~80 ng/mL C-THC*

<i>Component</i>	<i>Source</i>	<i>Source Lot Number</i>	<i>Expiration Date</i>
Methanol (LCMS)	Fisher	193941	-
Negative Urine	Pocatello Lab	POC031319	-
C-THC	Cerilliant	FE07171501	09/30/2020
Prepared:	01/16/2020		
Prepared By:	Sarah Pickle		
Expires:	09/30/2020		



# Idaho State Police Forensic Services

## AM #26 Blood THC and Metabolites Screen by LCMS-QQQ and AM #27 Quantitative Analysis of THC and Metabolites in Blood by LCMS-QQQ

**Methanol External Control Solution (Lot: WS011620)**

*10 µL of 1mg/mL THC, 100 µL of 100 µg/mL THC-OH, C-THC in 9790 µL MeOH  
Approximate concentration 1µg/mL.*

<i>Component</i>	<i>Source</i>	<i>Source Lot Number</i>	<i>Expiration Date</i>
Methanol (LCMS)	Fisher	193941	
THC	Cerilliant	FE09101501	11/30/2020
C-THC	Cerilliant	FE07171501	09/30/2020
THC-OH	Cerilliant	FE07221601	07/31/2021
Prepared:	01/16/2020		
Prepared By:	Tamara Salazar		
Expires:	09/30/2020		

**Blood External Control Solution (Lot: 011620)**

*100 µL of methanol external control solution was added to 9900 µL of blood.  
Approximately 10 ng/mL of each compound.*

<i>Component</i>	<i>Source</i>	<i>Source Lot Number</i>
Negative Blood	Hemostat	445283-3
Methanol External Control Solution	-	WS011620
Prepared:	01/16/2020	
Prepared by:	Tamara Salazar	
Expires:	09/30/2020	

**Idaho State Police  
Forensic Services  
Toxicology Discipline**

---

---

**Request for Departure from an Analytical Method**

---

Date of Request

01/13/2020

Forensic Scientist

Celena Shrum

Analytical Methods

Toxicology AM #25, Toxicology AM #26/27, and AM #28

Deviation

The expiration dates listed for the current batch of PinPoint ToxBox extraction plates are as follows:

- \*MDS (batch IDP-107-190725)- Expiration is 1/25/2020
- \*THC (batch IDP-108-190716)- Expiration is 1/16/2020
- \*MDQ P1 (batch IDP-111-190729)- Expiration is 1/29/2020
- \*MDQ P2 (batch IDP-112-190730)- Expiration is 1/30/2020

I am issuing a deviation to allow for the use of the remaining plates of these batches. The controls will be used to evaluate if the plate is working as intended. In addition, at least one external control must be included for each run.

*Celena Shrum*

Date: 01/13/2020

Celena Shrum

Toxicology Discipline Lead

\$




# AM #27 Cannabinoid Quant. Results

**Batch results** D:\MassHunter\Data\2020\AM 27-28\011620 AM 27 28 SP\QuantResults\AM 27.batch.bin  
**Calibration Last Update** 1/17/2020 7:52:16 AM

<b>Instrument</b>	Falco	<b>Data File</b>	MJ_QC Control.d
<b>Type</b>	Sample	<b>Sample</b>	MJ_QC Control
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Sarah Pickle
<b>Sample Position</b>	P3-A6	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	1/16/2020 3:51:34 PM		

**Sample Chromatogram**



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.468	113340	∞	8.9	253.53	1045811	4.2230 ng/ml
THC-COOH	1.489	97038	273.86	59.7	515.85	278760	13.6263 ng/ml
THC	3.360	230414	1713.31	28.4	154.69	6693184	4.2268 ng/ml

# AM #27 Cannabinoid Quant. Results




**Batch results** D:\MassHunter\Data\2020\AM 27-28\011620 AM 27 28 SP\QuantResults\AM 27.batch.bin  
**Calibration Last Update** 1/17/2020 7:52:16 AM

<b>Instrument</b>	Falco	<b>Data File</b>	MJ_Blood External Control.d
<b>Type</b>	Sample	<b>Sample</b>	MJ_Blood External Control
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Sarah Pickle
<b>Sample Position</b>	P3-G5	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	1/16/2020 4:21:56 PM		

**Sample Info.**

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.453	201833	∞	9.9	∞	1090042	10.3907 ng/ml
THC-COOH	1.489	67413	∞	59.4	212.03	286102	9.1696 ng/ml <b>Low</b>
THC	3.375	589520	3621.07	26.8	∞	7182306	9.8698 ng/ml




# AM #27 Cannabinoid Quant. Results



**Batch results** D:\MassHunter\Data\2020\AM 27-28\011620 AM 27 28 SP\QuantResults\AM 27.batch.bin  
**Calibration Last Update** 1/17/2020 7:52:16 AM

<b>Instrument</b>	Falco	<b>Data File</b>	MJ_Negative Blood.d
<b>Type</b>	Sample	<b>Sample</b>	MJ_Negative Blood
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Sarah Pickle
<b>Sample Position</b>	P3-H5	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	1/16/2020 4:06:45 PM		
<b>Sample Info.</b>			

## Sample Chromatogram



\$




# AM #27 Cannabinoid Quant. Results

**Batch results** D:\MassHunter\Data\2020\AM 27-28\011620 AM 27 28 SP\QuantResults\AM 27.batch.bin  
**Calibration Last Update** 1/17/2020 7:52:16 AM

<b>Instrument</b>	Falco	<b>Data File</b>	MJ_Urine External Control.d
<b>Type</b>	Sample	<b>Sample</b>	MJ_Urine External Control
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Sarah Pickle
<b>Sample Position</b>	P3-E4	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	1/16/2020 4:52:19 PM		

**Sample Chromatogram**



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-COOH	1.489	480898	2248.15	60.3	∞	323575	58.7191 ng/ml


# AM #27 Cannabinoid Quant. Results



**Batch results** D:\MassHunter\Data\2020\AM 27-28\011620 AM 27 28 SP\QuantResults\AM 27.batch.bin  
**Calibration Last Update** 1/17/2020 7:52:16 AM

<b>Instrument</b>	Falco	<b>Data File</b>	MJ_Negative Urine.d
<b>Type</b>	Sample	<b>Sample</b>	MJ_Negative Urine
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Sarah Pickle
<b>Sample Position</b>	P3-F4	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	1/16/2020 4:37:08 PM		
<b>Sample Info.</b>			

## Sample Chromatogram






# AM #27 Cannabinoids Quant. Calibration Curve Report

**Batch results** D:\MassHunter\Data\2020\AM 27-28\011620 AM 27 28 SP\QuantResults\AM 27.batch.bin  
**Last Cal. Update** 1/17/2020 7:52 AM  
**Analyst Name** ISP\Datastor  
**Analyte** THC **Internal Standard** THC-D3

THC - 7 Levels, 7 Levels Used, 7 Points, 7 Points Used, 0 QCs



Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
MJ_Cal 1	1	✓	1.0	1.1	108.2
MJ_Cal 2	2	✓	3.0	2.9	97.8
MJ_Cal 3	3	✓	5.0	4.8	95.3
MJ_Cal 4	4	✓	10.0	10.0	99.6
MJ_Cal 5	5	✓	25.0	24.6	98.3
MJ_Cal 6	6	✓	50.0	50.1	100.3
MJ_Cal 7	7	✓	100.0	100.6	100.6


S



# AM #27 Cannabinoids Quant. Calibration Curve Report

**Batch results** D:\MassHunter\Data\2020\AM 27-28\011620 AM 27 28 SP\QuantResults\AM 27.batch.bin  
**Last Cal. Update** 1/17/2020 7:52 AM  
**Analyst Name** ISP\Datastor  
**Analyte** THC-COOH **Internal Standard** THC-COOH-D9

THC-COOH - 7 Levels, 7 Levels Used, 7 Points, 7 Points Used, 0 QCs




Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
MJ_Cal 1	1	✓	5.0	4.9	97.1
MJ_Cal 2	2	✓	10.0	10.7	107.0
MJ_Cal 3	3	✓	20.0	19.2	95.9
MJ_Cal 4	4	✓	50.0	50.9	101.9
MJ_Cal 5	5	✓	75.0	73.8	98.4
MJ_Cal 6	6	✓	100.0	99.2	99.2
MJ_Cal 7	7	✓	250.0	251.3	100.5

\$



# AM #27 Cannabinoids Quant. Calibration Curve Report

**Batch results** D:\MassHunter\Data\2020\AM 27-28\011620 AM 27 28 SP\QuantResults\AM 27.batch.bin  
**Last Cal. Update** 1/17/2020 7:52 AM  
**Analyst Name** ISP\Datastor  
**Analyte** THC-OH **Internal Standard** THC-OH-D3



Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
MJ_Cal 1	1	x	1.0	0.0	0.0
MJ_Cal 2	2	✓	3.0	2.9	97.6
MJ_Cal 3	3	✓	5.0	5.0	99.1
MJ_Cal 4	4	✓	10.0	10.9	109.0
MJ_Cal 5	5	✓	25.0	23.4	93.5
MJ_Cal 6	6	✓	50.0	50.0	100.0
MJ_Cal 7	7	✓	100.0	100.8	100.8

# AM #27 Cannabinoid Quant. Results




**Batch results** D:\MassHunter\Data\2020\AM 27-28\011620 AM 27 28 SP\QuantResults\AM 27.batch.bin  
**Calibration Last Update** 1/17/2020 7:52:16 AM

<b>Instrument</b>	Falco	<b>Data File</b>	MJ_Cal 1.d
<b>Type</b>	Cal	<b>Sample</b>	MJ_Cal 1
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Sarah Pickle
<b>Sample Position</b>	P3-B6	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	1/16/2020 2:50:46 PM		

**Sample Info.**

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-COOH	1.474	35031	181.43	55.4	103.85	276532	4.8531 ng/ml <b>Low</b>
THC	3.345	50808	446.88	31.8	45.82	6455047	1.0824 ng/ml <b>Low</b>

# AM #27 Cannabinoid Quant. Results




**Batch results** D:\MassHunter\Data\2020\AM 27-28\011620 AM 27 28 SP\QuantResults\AM 27.batch.bin  
**Calibration Last Update** 1/17/2020 7:52:16 AM

<b>Instrument</b>	Falco	<b>Data File</b>	MJ_Cal 2.d
<b>Type</b>	Cal	<b>Sample</b>	MJ_Cal 2
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Sarah Pickle
<b>Sample Position</b>	P3-C6	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	1/16/2020 2:58:31 PM		

**Sample Info.**

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.483	96192	∞	8.9	88.73	1042685	2.9280 ng/ml <b>Low</b>
THC-COOH	1.504	75444	137.10	53.2	683.30	275024	10.7026 ng/ml
THC	3.375	153570	788.22	27.5	226.09	6531156	2.9347 ng/ml <b>Low</b>






# AM #27 Cannabinoid Quant. Results

**Batch results** D:\MassHunter\Data\2020\AM 27-28\011620 AM 27 28 SP\QuantResults\AM 27.batch.bin  
**Calibration Last Update** 1/17/2020 7:52:16 AM

<b>Instrument</b>	Falco	<b>Data File</b>	MJ_Cal 3.d
<b>Type</b>	Cal	<b>Sample</b>	MJ_Cal 3
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Sarah Pickle
<b>Sample Position</b>	P3-D6	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	1/16/2020 3:06:05 PM		

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.468	121925	∞	8.7	∞	1037848	4.9542 ng/ml
THC-COOH	1.504	135095	239.28	59.5	651.96	276643	19.1824 ng/ml
THC	3.375	253883	2798.12	27.0	133.18	6516409	4.7639 ng/ml

\$




# AM #27 Cannabinoid Quant. Results

**Batch results** D:\MassHunter\Data\2020\AM 27-28\011620 AM 27 28 SP\QuantResults\AM 27.batch.bin  
**Calibration Last Update** 1/17/2020 7:52:16 AM

<b>Instrument</b>	Falco	<b>Data File</b>	MJ_Cal 4.d
<b>Type</b>	Cal	<b>Sample</b>	MJ_Cal 4
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Sarah Pickle
<b>Sample Position</b>	P3-E6	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	1/16/2020 3:13:40 PM		

**Sample Chromatogram**



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.468	189550	∞	9.7	119.43	989909	10.8985 ng/ml
THC-COOH	1.504	333237	∞	59.7	2123.42	258342	50.9417 ng/ml
THC	3.375	512182	2510.44	27.4	358.76	6185497	9.9555 ng/ml




# AM #27 Cannabinoid Quant. Results

**Batch results** D:\MassHunter\Data\2020\AM 27-28\011620 AM 27 28 SP\QuantResults\AM 27.batch.bin  
**Calibration Last Update** 1/17/2020 7:52:16 AM

<b>Instrument</b>	Falco	<b>Data File</b>	MJ_Cal 5.d
<b>Type</b>	Cal	<b>Sample</b>	MJ_Cal 5
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Sarah Pickle
<b>Sample Position</b>	P3-F6	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	1/16/2020 3:21:14 PM		

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.468	356203	∞	12.7	1098.78	1027484	23.3645 ng/ml
THC-COOH	1.489	494876	2909.54	60.5	3930.50	265083	73.8016 ng/ml
THC	3.375	1350880	7206.95	26.1	888.16	6552198	24.5641 ng/ml


# AM #27 Cannabinoid Quant. Results



**Batch results** D:\MassHunter\Data\2020\AM 27-28\011620 AM 27 28 SP\QuantResults\AM 27.batch.bin  
**Calibration Last Update** 1/17/2020 7:52:16 AM

<b>Instrument</b>	Falco	<b>Data File</b>	MJ_Cal 6.d
<b>Type</b>	Cal	<b>Sample</b>	MJ_Cal 6
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Sarah Pickle
<b>Sample Position</b>	P3-G6	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	1/16/2020 3:28:48 PM		

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.468	667898	∞	12.3	2292.38	984406	50.0165 ng/ml
THC-COOH	1.489	635819	∞	60.1	8092.94	253582	99.1781 ng/ml
THC	3.375	2627958	14125.50	26.5	3628.38	6224605	50.1436 ng/ml




# AM #27 Cannabinoid Quant. Results

**Batch results** D:\MassHunter\Data\2020\AM 27-28\011620 AM 27 28 SP\QuantResults\AM 27.batch.bin  
**Calibration Last Update** 1/17/2020 7:52:16 AM

<b>Instrument</b>	Falco	<b>Data File</b>	MJ_Cal 7.d
<b>Type</b>	Cal	<b>Sample</b>	MJ_Cal 7
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Sarah Pickle
<b>Sample Position</b>	P3-H6	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	1/16/2020 3:36:23 PM		

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.468	1281837	∞	12.3	3249.69	977623	100.8382 ng/ml
THC-COOH	1.504	1538510	∞	60.1	7571.08	242370	251.3405 ng/ml
THC	3.375	5285605	33005.20	26.8	∞	6233643	100.5558 ng/ml