









**Worklist: 4114**

<u>LAB CASE</u>	<u>ITEM</u>	<u>ITEM TYPE</u>	<u>DESCRIPTION</u>	
C2020-0455	1	BCK	AM 27 Blood THC Quant by LC-QQQ	
C2020-0459	1	UCK	AM 27 Urine Cannabinoids Confirmation by LC-QQQ	
C2020-0489	2	BCK	AM 27 Blood THC Quant by LC-QQQ	
C2020-0506	1	BCK	AM 27 Blood THC Quant by LC-QQQ	
C2020-0525	1	BCK	AM 27 Blood THC Quant by LC-QQQ	
C2020-0537	1	UCK	AM 27 Urine Cannabinoids Confirmation by LC-QQQ	
C2020-0554	3	BCK	AM 27 Blood THC Quant by LC-QQQ	
C2020-0562	1	BCK	AM 27 Blood THC Quant by LC-QQQ	

## AM# 27: Quantitation of THC and Metabolites in Blood and Urine by LC-MS/MS

Extraction Date: 3/25/20Analyst: Britany Wylie

Plate lot#: 200303

Plate Expiration: 9-3-2020

**Mobile phase A:** 0.1% Formic Acid in LCMS Water  
MTBE

LCMS Methanol

**Mobile phase B:** 0.1% Formic acid in Acetonitrile  
Hexane

**Blank Blood Lot:** 20A52255**Urine Blank:** 32420**Column:** UCT Selectra DA 100 x 2.1mm 3um**LCMS-QQQ ID:** 69679

### Pre-Analytic:

- 1. Check levels of mobile phases and needle wash refill as needed. Ensure waste is not full.
- 2. Ensure correct column is installed and begin mobile phase flow allow to equilibrate ~ 30 minutes.

### Analytic:

- 1. Remove standards, plate, controls, and samples from cold storage. Allow to reach room temperature.
- 2. Urine hydrolysis: add 1.5 ml urine to blank plate, add 250 ul 1N KOH mix and incubate at 40 degrees for 15 minutes.  
Pipette **1000µL blood (calibrated pipette) Pipette ID: k52558g** in wells of analytical (standards) plate.
- 3. Place on shaking incubator at ambient temp., 900rpm for 15 minutes. *Shaker ID: 66759*
- 4. Pipette **500µL 0.1% formic acid in water blood sample** in wells of analytical plate.
- 5. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 6. Transfer **800µL of blood+acid or urine acid** mixture to corresponding wells of SLE+ plate.
- 7. Apply positive pressure for approx. 10-15 seconds (or until no liquid remains on top of sorbent).  
*(Load at 85-100 PSI- Selector to the right)* Manifold ID: 66792
- 8. Wait 5 minutes.
- 9. Add **2.25mL MTBE. (Add in 3 increments of 750uL)**
- 10. Wait 5 minutes.
- 11. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 12. Add **2.25mL Hexane. (Add in 3 increments of 750uL)**
- 13. Wait 5 minutes.
- 14. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 15. Remove plate containing eluate. Place on SPE Dry and evaporate to dryness at approx. 35°C.  
*SPE Dry ID: 66819*
- 16. Reconstitute in **100µL 100% MeOH** and heat seal plate with foil. Place in autosampler and run worklist.

### Post-Analytic

- 1. Create batch and process data.
- 2. Make any necessary integration changes, Curve weighting of Linear 1/x with  $r^2$  values  $\geq 0.98$  for each analyte
- 3. RT +/- 3% or 0.100 min, whichever is greater, +/- 20% Accuracy for greater than (+/- 30% for 10ng/ml or less).  
Ion ratios must be within +/- 20% of the averaged calibrators
- 4. Case sample response for THC and OH-THC 3ng/mL (quantitative blood), Carboxy-THC: 10ng/mL (qualitative only) will be reported. Samples with a THC or OH-THC response over 50 ng/mL will be reported out as greater than 50 ng/mL.
- 5. Did all QCs pass for each analyte? (if not is it describe in comments section)
- 6. Enter QCs into control charting.
- 7. Central File Packet to include: LIMS Worklist, Method Checklist, Calibration and Control Reports

COMMENTS: Urine cases only evaluated for Carboxy-THC


*Curves limited: 3-100- THC and Hydroxy-THC; 5-250 Carboxy-THC*

# AM #27 Cannabinoids

**Batch results** D:\MassHunter\Data\2020 Data\am 27-28 3-25-2020\QuantResults\thcq.batch.bin  
**Calibration Last Update** 3/25/2020 10:12:04 PM

<b>Instrument</b>	69679	<b>Data File</b>	mj negative blood.d
<b>Type</b>	Sample	<b>Sample</b>	mj negative blood
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Britany Wylie
<b>Sample Position</b>	P3-A2	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	3/25/2020 1:42:28 PM		
<b>Sample Info.</b>			

## Sample Chromatogram




# AM #27 Cannabinoids

**Batch results** D:\MassHunter\Data\2020 Data\am 27-28 3-25-2020\QuantResults\thcq.batch.bin  
**Calibration Last Update** 3/25/2020 10:12:04 PM

<b>Instrument</b>	69679	<b>Data File</b>	mj internal control.d
<b>Type</b>	QC	<b>Sample</b>	mj internal control
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Britany Wylie
<b>Sample Position</b>	P3-H1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	3/25/2020 1:34:44 PM		

**Sample Chromatogram**




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.375	79863	∞	11.2	∞	553546	4.999 ng/ml
THC-COOH	1.400	64671	347.2	256.3	1384.6	413384	15.664 ng/ml
THC	3.138	25934	∞	25.9	∞	655646	4.663 ng/ml

# AM #27 Cannabinoids

**Batch results** D:\MassHunter\Data\2020 Data\am 27-28 3-25-2020\QuantResults\thcq.batch.bin  
**Calibration Last Update** 3/25/2020 10:12:04 PM

<b>Instrument</b>	69679	<b>Data File</b>	mj negative urine.d
<b>Type</b>	Sample	<b>Sample</b>	mj negative urine
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Britany Wylie
<b>Sample Position</b>	P3-H2	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	3/25/2020 3:30:03 PM		

## Sample Chromatogram




# AM #27 Cannabinoids

**Batch results** D:\MassHunter\Data\2020 Data\am 27-28 3-25-2020\QuantResults\thcq.batch.bin  
**Calibration Last Update** 3/25/2020 10:12:04 PM

**Instrument** 69679 **Data File** MJ external control urine.d  
**Type** Sample **Sample** MJ external control urine  
**Acq. Method** AM 27 THC quant.m **Operator** Britany Wylie  
**Sample Position** P3-A3 **Comment**  
**Injection Volume** 10  
**Acq. Date-Time** 3/25/2020 3:37:48 PM  
**Sample Info.**

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.375	184463	∞	12.6	∞	551257	11.314 ng/ml
THC-COOH	1.415	145815	1688.9	233.6	595.0	356689	38.034 ng/ml
THC	3.138	8719	∞	29.3	4256607 436253. 1	345868	3.229 ng/ml

**Toxicology AM method 27/26 external prep information**

working solution 1 ug/ml in meoh C-THC, THC-OH, THC

Stock solution 1mg/ml 7.5 ul each THC, 100 ug/ml 150 ul C-THC, 75 ul THC-OH in 9767.5 ul meOH

Ppd 2/13/20 Exp: 8/13/20 lot 21320 by AMN

Drug	lot	expiration
C-THC	FE07171501	9/1/2020
THC-OH	FE07721601	7/1/2021
THC	FE001041701	3/1/2022

AM 27/26 blood control 100 ul working solution lot (91319) in 9900 ul blood lot (20A52255)

ppd 02/13/20 Exp 08/13/20 lot b81320 Concentration 7.5 ng/ml THC, THC-OH and 15 ng/ml C-THC by AMN


AM 27/26 urine control 400 ul working solution lot (21320) in 9600 ul urine lot (11420)

ppd 02/13/20 Exp 08/13/20 lot u81320 Concentration 30 ng/ml THC, THC-OH and 60 ng/ml C-THC by AMN

# Compound Calibration Report

**Batch results** D:\MassHunter\Data\2020 Data\am 27-28 3-25-2020\QuantResults\thcq.batch.bin  
**Last Cal. Update** 3/25/2020 10:12 PM  
**Analyst Name** ISP\datastor  
**Analyte** THC **Internal Standard** THC-d3

THC - 7 Levels, 6 Levels Used, 7 Points, 6 Points Used, 1 QCs




Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
mj qc1	1	x	1.0	1.6	162.1
mj cal2	2	✓	3.0	3.3	110.3
mj cal 3	3	✓	5.0	5.0	100.2
mj cal 4	4	✓	10.0	9.3	93.0
mj cal 5	5	✓	25.0	24.1	96.3
mj cal 6	6	✓	50.0	48.9	97.8
mj cal 7	7	✓	100.0	102.4	102.4



# Compound Calibration Report

**Batch results** D:\MassHunter\Data\2020 Data\am 27-28 3-25-2020\QuantResults\thcq.batch.bin  
**Last Cal. Update** 3/25/2020 10:12 PM  
**Analyst Name** ISP\datastor  
**Analyte** THC-COOH **Internal Standard** THC-COOH-d9


THC-COOH - 7 Levels, 6 Levels Used, 7 Points, 6 Points Used, 1 QCs



Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
mj qc1	1	x	5.0	5.8	115.3
mj cal2	2	✓	10.0	10.8	107.8
mj cal 3	3	✓	20.0	19.6	98.0
mj cal 4	4	✓	50.0	48.4	96.7
mj cal 5	5	✓	75.0	73.1	97.5
mj cal 6	6	✓	100.0	97.9	97.9
mj cal 7	7	✓	250.0	255.2	102.1

# Compound Calibration Report

**Batch results** D:\MassHunter\Data\2020 Data\am 27-28 3-25-2020\QuantResults\thcq.batch.bin  
**Last Cal. Update** 3/25/2020 10:12 PM  
**Analyst Name** ISP\datastor  
**Analyte** THC-OH **Internal Standard** THC-OH-d3



Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
mj qc1	1	x	1.0	1.6	155.8
mj cal2	2	✓	3.0	3.3	109.5
mj cal 3	3	✓	5.0	5.3	105.2
mj cal 4	4	✓	10.0	8.5	85.4
mj cal 5	5	✓	25.0	24.8	99.4
mj cal 6	6	✓	50.0	49.4	98.8
mj cal 7	7	✓	100.0	101.6	101.6


# AM #27 Cannabinoids

**Batch results** D:\MassHunter\Data\2020 Data\am 27-28 3-25-2020\QuantResults\thcq.batch.bin  
**Calibration Last Update** 3/25/2020 10:12:04 PM

<b>Instrument</b>	69679	<b>Data File</b>	mj qc1.d
<b>Type</b>	Cal	<b>Sample</b>	mj qc1
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Britany Wylie
<b>Sample Position</b>	P3-A1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	3/25/2020 12:40:41 PM		

**Sample Info.**

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.375	23349	∞	8.1 <b>Low</b>	∞	575645	1.558 ng/ml <b>Low</b>
THC-COOH	1.415	17962	37.2	343.3 <b>High</b>	27.2	401087	5.766 ng/ml <b>Low</b>
THC	3.138	5525	∞	25.1	∞	605310	1.621 ng/ml <b>Low</b>


# AM #27 Cannabinoids

**Batch results** D:\MassHunter\Data\2020 Data\am 27-28 3-25-2020\QuantResults\thcq.batch.bin  
**Calibration Last Update** 3/25/2020 10:12:04 PM

<b>Instrument</b>	69679	<b>Data File</b>	mj cal2.d
<b>Type</b>	Cal	<b>Sample</b>	mj cal2
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Britany Wylie
<b>Sample Position</b>	P3-B1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	3/25/2020 12:48:25 PM		

**Sample Info.**

**Sample Chromatogram**




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.375	53560	∞	10.8	∞	578399	3.284 ng/ml
THC-COOH	1.400	40846	106.9	262.7	64.5	403066	10.780 ng/ml
THC	3.123	17054	∞	24.4	∞	655839	3.308 ng/ml

# AM #27 Cannabinoids

**Batch results** D:\MassHunter\Data\2020 Data\am 27-28 3-25-2020\QuantResults\thcq.batch.bin  
**Calibration Last Update** 3/25/2020 10:12:04 PM

**Instrument** 69679 **Data File** mj cal 3.d  
**Type** Cal **Sample** mj cal 3  
**Acq. Method** AM 27 THC quant.m **Operator** Britany Wylie  
**Sample Position** P3-C1 **Comment**  
**Injection Volume** 10  
**Acq. Date-Time** 3/25/2020 12:56:09 PM  
**Sample Info.**

### Sample Chromatogram




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.390	83412	∞	10.6	∞	548080	5.261 ng/ml
THC-COOH	1.400	79093	145.9	249.3	240.3	394004	19.591 ng/ml
THC	3.138	26633	∞	28.3	∞	618873	5.011 ng/ml

# AM #27 Cannabinoids

**Batch results** D:\MassHunter\Data\2020 Data\am 27-28 3-25-2020\QuantResults\thcq.batch.bin  
**Calibration Last Update** 3/25/2020 10:12:04 PM

<b>Instrument</b>	69679	<b>Data File</b>	mj cal 4.d
<b>Type</b>	Cal	<b>Sample</b>	mj cal 4
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Britany Wylie
<b>Sample Position</b>	P3-D1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	3/25/2020 1:03:54 PM		

## Sample Chromatogram




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.375	144063	∞	14.2	∞	573581	8.545 ng/ml
THC-COOH	1.400	217301	3180.4	227.7	1892.8	413640	48.364 ng/ml
THC	3.138	56245	∞	27.1	∞	654893	9.296 ng/ml

# AM #27 Cannabinoids

**Batch results** D:\MassHunter\Data\2020 Data\am 27-28 3-25-2020\QuantResults\thcq.batch.bin  
**Calibration Last Update** 3/25/2020 10:12:04 PM

**Instrument** 69679 **Data File** mj cal 5.d  
**Type** Cal **Sample** mj cal 5  
**Acq. Method** AM 27 THC quant.m **Operator** Britany Wylie  
**Sample Position** P3-E1 **Comment**  
**Injection Volume** 10  
**Acq. Date-Time** 3/25/2020 1:11:36 PM  
**Sample Info.**

## Sample Chromatogram




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.375	428948	∞	12.4	∞	577772	24.843 ng/ml
THC-COOH	1.400	336994	1127.3	229.1	973518.4	418663	73.148 ng/ml
THC	3.138	158437	∞	26.7	∞	677953	24.077 ng/ml

# AM #27 Cannabinoids

**Batch results** D:\MassHunter\Data\2020 Data\am 27-28 3-25-2020\QuantResults\thcq.batch.bin  
**Calibration Last Update** 3/25/2020 10:12:04 PM

<b>Instrument</b>	69679	<b>Data File</b>	mj cal 6.d
<b>Type</b>	Cal	<b>Sample</b>	mj cal 6
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Britany Wylie
<b>Sample Position</b>	P3-F1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	3/25/2020 1:19:20 PM		

**Sample Chromatogram**



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.390	809233	∞	13.0	∞	545580	49.421 ng/ml
THC-COOH	1.415	420614	3753.7	228.4	2272.1	387939	97.906 ng/ml
THC	3.138	332332	∞	25.6	∞	689253	48.923 ng/ml




# AM #27 Cannabinoids

**Batch results** D:\MassHunter\Data\2020 Data\am 27-28 3-25-2020\QuantResults\thcq.batch.bin  
**Calibration Last Update** 3/25/2020 10:12:04 PM

<b>Instrument</b>	69679	<b>Data File</b>	mj cal 7.d
<b>Type</b>	Cal	<b>Sample</b>	mj cal 7
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Britany Wylie
<b>Sample Position</b>	P3-G1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	3/25/2020 1:27:02 PM		

**Sample Chromatogram**



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.375	1666115	∞	13.4	∞	544937	101.647 ng/ml
THC-COOH	1.400	1049964	19892.0	224.4	7134.4	367274	255.211 ng/ml
THC	3.138	637549	∞	26.8	∞	627017	102.384 ng/ml