









**APPROVED**

By Anne Nord at 8:56 am, Apr 21, 2020

4/20/2020 BW

**Worklist: 4176**

<u>LAB CASE</u>	<u>ITEM</u>	<u>ITEM TYPE</u>	<u>DESCRIPTION</u>	
C2020-0606	1	BCK	AM 27 Blood THC Quant by LC-QQQ	
C2020-0607	1	BCK	AM 27 Blood THC Quant by LC-QQQ	
C2020-0622	1	UCK	AM 27 Urine Cannabinoids Confirmation by LC-QQQ	
C2020-0637	1	BCK	AM 27 Blood THC Quant by LC-QQQ	
C2020-0639	2	UCK	AM 27 Urine Cannabinoids Confirmation by LC-QQQ	
C2020-0661	1	UCK	AM 27 Urine Cannabinoids Confirmation by LC-QQQ	
C2020-0665	1	BCK	AM 27 Blood THC Quant by LC-QQQ	
C2020-0676	1	UCK	AM 27 Urine Cannabinoids Confirmation by LC-QQQ	

## AM# 27: Quantitation of THC and Metabolites in Blood and Urine by LC-MS/MS

Extraction Date: 4/17/20Analyst: Britany Wylie

Plate lot#: 200303

Plate Expiration: 9-3-2020

**Mobile phase A:** 0.1% Formic Acid in LCMS Water  
MTBE

LCMS Methanol

**Mobile phase B:** 0.1% Formic acid in Acetonitrile  
Hexane

**Blank Blood Lot:** 20A52255**Urine Blank:** 32420**Column:** UCT Selectra DA 100 x 2.1mm 3um**LCMS-QQQ ID:** 69679

### Pre-Analytic:

- 1. Check levels of mobile phases and needle wash refill as needed. Ensure waste is not full.
- 2. Ensure correct column is installed and begin mobile phase flow allow to equilibrate ~ 30 minutes.

### Analytic:

- 1. Remove standards, plate, controls, and samples from cold storage. Allow to reach room temperature.
- 2. Urine hydrolysis: add 1.5 ml urine to blank plate, add 250 ul 1N KOH mix and incubate at 40 degrees for 15 minutes.  
**Pipette 1000µL blood (calibrated pipette) Pipette ID: k52558g** in wells of analytical (standards) plate.
- 3. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 4. Pipette **500µL 0.1% formic acid in water blood sample** in wells of analytical plate.
- 5. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 6. Transfer **800µL of blood+acid or urine acid** mixture to corresponding wells of SLE+ plate.
- 7. Apply positive pressure for approx. 10-15 seconds (or until no liquid remains on top of sorbent).  
**(Load at 85-100 PSI- Selector to the right)** Manifold ID: 66792
- 8. Wait 5 minutes.
- 9. Add **2.25mL MTBE. (Add in 3 increments of 750uL)**
- 10. Wait 5 minutes.
- 11. Apply positive pressure for approx. 15 seconds. **(10-15 PSI- Selector to the left).**
- 12. Add **2.25mL Hexane. (Add in 3 increments of 750uL)**
- 13. Wait 5 minutes.
- 14. Apply positive pressure for approx. 15 seconds. **(10-15 PSI- Selector to the left).**
- 15. Remove plate containing eluate. Place on SPE Dry and evaporate to dryness at approx. 35°C.  
*SPE Dry ID: 66819*
- 16. Reconstitute in **100µL 100% MeOH** and heat seal plate with foil. Place in autosampler and run worklist.

### Post-Analytic

- 1. Create batch and process data.
- 2. Make any necessary integration changes, Curve weighting of Linear 1/x with  $r^2$  values  $\geq 0.98$  for each analyte
- 3. RT +/- 3% or 0.100 min, whichever is greater, +/- 20% Accuracy for greater than (+/- 30% for 10ng/ml or less).  
Ion ratios must be within +/- 20% of the averaged calibrators
- 4. Case sample response for THC and OH-THC 3ng/mL (quantitative blood), Carboxy-THC: 10ng/mL (qualitative only) will be reported. Samples with a THC or OH-THC response over 50 ng/mL will be reported out as greater than 50 ng/mL.
- 5. Did all QCs pass for each analyte? (if not is it describe in comments section)
- 6. Enter QCs into control charting.
- 7. Central File Packet to include: LIMS Worklist, Method Checklist, Calibration and Control Reports


COMMENTS: Urine cases only evaluated for Carboxy-THC  
*Curves limited: 3-100- Hydroxy-THC; 10-250 Carboxy-THC*

# AM #27 Cannabinoids

**Batch results** D:\MassHunter\Data\2020 Data\am 27-28 4-17-2020\QuantResults\thcq.batch.bin  
**Calibration Last Update** 4/18/2020 10:24:06 AM

<b>Instrument</b>	69679	<b>Data File</b>	mj negative blood.d
<b>Type</b>	Sample	<b>Sample</b>	mj negative blood
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Britany Wylie
<b>Sample Position</b>	P3-A2	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	4/17/2020 2:50:06 PM		
<b>Sample Info.</b>			

**Sample Chromatogram**




# AM #27 Cannabinoids

**Batch results** D:\MassHunter\Data\2020 Data\am 27-28 4-17-2020\QuantResults\thcq.batch.bin  
**Calibration Last Update** 4/18/2020 10:24:06 AM

<b>Instrument</b>	69679	<b>Data File</b>	mj internal control.d
<b>Type</b>	QC	<b>Sample</b>	mj internal control
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Britany Wylie
<b>Sample Position</b>	P3-H1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	4/17/2020 2:42:23 PM		

**Sample Chromatogram**



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.390	93608	∞	9.7	∞	1013531	4.900 ng/ml
THC-COOH	1.415	115769	148.1	201.8	150.9	634011	15.619 ng/ml
THC	3.138	34203	∞	26.7	∞	923211	4.216 ng/ml


BW

# AM #27 Cannabinoids

**Batch results** D:\MassHunter\Data\2020 Data\am 27-28 4-17-2020\QuantResults\thcq.batch.bin  
**Calibration Last Update** 4/18/2020 10:24:06 AM

<b>Instrument</b>	69679	<b>Data File</b>	mj urine neg.d
<b>Type</b>	Sample	<b>Sample</b>	mj urine neg
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Britany Wylie
<b>Sample Position</b>	P3-F2	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	4/17/2020 4:06:57 PM		
<b>Sample Info.</b>			

## Sample Chromatogram




# AM #27 Cannabinoids

**Batch results** D:\MassHunter\Data\2020 Data\am 27-28 4-17-2020\QuantResults\thcq.batch.bin  
**Calibration Last Update** 4/18/2020 10:24:06 AM

<b>Instrument</b>	69679	<b>Data File</b>	mj urine ext ctrl.d
<b>Type</b>	Sample	<b>Sample</b>	mj urine ext ctrl
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Britany Wylie
<b>Sample Position</b>	P3-G2	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	4/17/2020 4:14:41 PM		
<b>Sample Info.</b>			

**Sample Chromatogram**



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.390	316651	∞	10.5	∞	1006449	16.354 ng/ml
THC-COOH	1.415	258449	4819.4	178.2	5471.9	512776	40.441 ng/ml
THC	3.138	123798	∞	24.5	∞	697126	18.986 ng/ml

**Toxicology AM method 27/26 external prep information**

working solution 1 ug/ml in meoh C-THC, THC-OH, THC

Stock solution 1mg/ml 7.5 ul each THC, 100 ug/ml 150 ul C-THC, 75 ul THC-OH in 9767.5 ul meOH

Ppd 2/13/20 Exp: 8/13/20 lot 21320 by AMN

Drug	lot	expiration
C-THC	FE07171501	9/1/2020
THC-OH	FE07721601	7/1/2021
THC	FE001041701	3/1/2022

AM 27/26 blood control 100 ul working solution lot (91319) in 9900 ul blood lot (20A52255)

ppd 02/13/20 Exp 08/13/20 lot b81320 Concentration 7.5 ng/ml THC, THC-OH and 15 ng/ml C-THC by AMN

AM 27/26 urine control 400 ul working solution lot (21320) in 9600 ul urine lot (32420)

ppd 4/17/20 Exp 9120 lot u101720 Concentration 30 ng/ml THC, THC-OH and 60 ng/ml C-THC by BAW


# Compound Calibration Report

**Batch results**      D:\MassHunter\Data\2020 Data\am 27-28 4-17-2020\QuantResults\thcq.batch.bin  
**Last Cal. Update**    4/18/2020 10:24 AM  
**Analyst Name**        ISP\datastor  
**Analyte**                THC

*BW*

**Internal Standard**      THC-d3

THC - 7 Levels, 7 Levels Used, 7 Points, 7 Points Used, 1 QCs



Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
mj qc1	1	✓	1.0	1.2	118.3
mj cal2	2	✓	3.0	2.9	97.9
mj cal 3	3	✓	5.0	4.7	93.9
mj cal 4	4	✓	10.0	9.2	92.0
mj cal 5	5	✓	25.0	24.1	96.3
mj cal 6	6	✓	50.0	49.7	99.3
mj cal 7	7	✓	100.0	102.2	102.2




# Compound Calibration Report

**Batch results** D:\MassHunter\Data\2020 Data\am 27-28 4-17-2020\QuantResults\thcq.batch.bin  
**Last Cal. Update** 4/18/2020 10:24 AM  
**Analyst Name** ISP\datastor  
**Analyte** THC-COOH

*BW*

**Internal Standard**      **THC-COOH-d9**

THC-COOH - 7 Levels, 6 Levels Used, 7 Points, 6 Points Used, 1 QCs




Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
mj qc1	1	x	5.0	6.1	122.3
mj cal2	2	✓	10.0	10.6	106.1
mj cal 3	3	✓	20.0	19.6	97.9
mj cal 4	4	✓	50.0	48.7	97.5
mj cal 5	5	✓	75.0	73.2	97.6
mj cal 6	6	✓	100.0	99.8	99.8
mj cal 7	7	✓	250.0	253.1	101.3

# Compound Calibration Report

**Batch results** D:\MassHunter\Data\2020 Data\am 27-28 4-17-2020\QuantResults\thcq.batch.bin  
**Last Cal. Update** 4/18/2020 10:24 AM  
**Analyst Name** ISP\datastor  
**Analyte** THC-OH **Internal Standard** THC-OH-d3

*BW*

THC-OH - 7 Levels, 6 Levels Used, 7 Points, 6 Points Used, 1 QCs



Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
mj qc1	1	x	1.0	0.9	88.3
mj cal2	2	✓	3.0	3.3	109.3
mj cal 3	3	✓	5.0	4.4	88.0
mj cal 4	4	✓	10.0	10.3	102.7
mj cal 5	5	✓	25.0	24.9	99.5
mj cal 6	6	✓	50.0	50.4	100.7
mj cal 7	7	✓	100.0	99.8	99.8

# AM #27 Cannabinoids


BW

**Batch results** D:\MassHunter\Data\2020 Data\am 27-28 4-17-2020\QuantResults\thcq.batch.bin  
**Calibration Last Update** 4/18/2020 10:24:06 AM

<b>Instrument</b>	69679	<b>Data File</b>	mj qc1.d
<b>Type</b>	Cal	<b>Sample</b>	mj qc1
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Britany Wylie
<b>Sample Position</b>	P3-A1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	4/17/2020 1:48:19 PM		

**Sample Info.**

## Sample Chromatogram




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.390	16442	∞	15.8 <b>High</b>	∞	1141617	0.883 ng/ml <b>Low</b>
THC-COOH	1.415	42673	579.6	236.0 <b>High</b>	377.9	716686	6.116 ng/ml <b>Low</b>
THC	3.138	8267	∞	23.9	∞	1009604	1.183 ng/ml <b>Low</b>

# AM #27 Cannabinoids

**Batch results** D:\MassHunter\Data\2020 Data\am 27-28 4-17-2020\QuantResults\thcq.batch.bin  
**Calibration Last Update** 4/18/2020 10:24:06 AM

<b>Instrument</b>	69679	<b>Data File</b>	mj cal2.d
<b>Type</b>	Cal	<b>Sample</b>	mj cal2
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Britany Wylie
<b>Sample Position</b>	P3-B1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	4/17/2020 1:56:04 PM		

**Sample Chromatogram**




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.390	69404	∞	9.8	∞	1140018	3.278 ng/ml
THC-COOH	1.415	83005	287.8	209.9	1088.7	705203	10.608 ng/ml
THC	3.138	26724	∞	24.8	∞	1074261	2.937 ng/ml <b>Low</b>

# AM #27 Cannabinoids

**Batch results** D:\MassHunter\Data\2020 Data\am 27-28 4-17-2020\QuantResults\thcq.batch.bin  
**Calibration Last Update** 4/18/2020 10:24:06 AM

<b>Instrument</b>	69679	<b>Data File</b>	mj cal 3.d
<b>Type</b>	Cal	<b>Sample</b>	mj cal 3
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Britany Wylie
<b>Sample Position</b>	P3-C1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	4/17/2020 2:03:47 PM		

**Sample Chromatogram**




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.390	91444	∞	13.7	∞	1106161	4.401 ng/ml
THC-COOH	1.415	162997	488.6	194.5	678040.7	697197	19.572 ng/ml
THC	3.138	39465	∞	25.1	∞	948191	4.697 ng/ml

# AM #27 Cannabinoids

**Batch results** D:\MassHunter\Data\2020 Data\am 27-28 4-17-2020\QuantResults\thcq.batch.bin  
**Calibration Last Update** 4/18/2020 10:24:06 AM

<b>Instrument</b>	69679	<b>Data File</b>	mj cal 4.d
<b>Type</b>	Cal	<b>Sample</b>	mj cal 4
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Britany Wylie
<b>Sample Position</b>	P3-D1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	4/17/2020 2:11:30 PM		

**Sample Chromatogram**




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.390	212689	∞	10.6	∞	1082474	10.266 ng/ml
THC-COOH	1.415	411566	178.1	183.6	563.9	673253	48.726 ng/ml
THC	3.138	77378	∞	25.9	∞	915761	9.203 ng/ml

# AM #27 Cannabinoids

**Batch results** D:\MassHunter\Data\2020 Data\am 27-28 4-17-2020\QuantResults\thcq.batch.bin  
**Calibration Last Update** 4/18/2020 10:24:06 AM

<b>Instrument</b>	69679	<b>Data File</b>	mj cal 5.d
<b>Type</b>	Cal	<b>Sample</b>	mj cal 5
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Britany Wylie
<b>Sample Position</b>	P3-E1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	4/17/2020 2:19:14 PM		

**Sample Chromatogram**



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.390	554657	∞	10.9	∞	1155819	24.870 ng/ml
THC-COOH	1.400	616475	410915.3	177.3	23427.7	664240	73.189 ng/ml
THC	3.153	189145	∞	25.8	∞	837030	24.072 ng/ml


# AM #27 Cannabinoids

**Batch results** D:\MassHunter\Data\2020 Data\am 27-28 4-17-2020\QuantResults\thcq.batch.bin  
**Calibration Last Update** 4/18/2020 10:24:06 AM

<b>Instrument</b>	69679	<b>Data File</b>	mj cal 6.d
<b>Type</b>	Cal	<b>Sample</b>	mj cal 6
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Britany Wylie
<b>Sample Position</b>	P3-F1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	4/17/2020 2:26:57 PM		

**Sample Info.**

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.390	1039718	∞	11.5	∞	1066627	50.373 ng/ml
THC-COOH	1.415	816620	9563.1	179.1	8485.2	641865	99.767 ng/ml
THC	3.138	449551	∞	24.9	∞	957498	49.667 ng/ml




# AM #27 Cannabinoids

**Batch results** D:\MassHunter\Data\2020 Data\am 27-28 4-17-2020\QuantResults\thcq.batch.bin  
**Calibration Last Update** 4/18/2020 10:24:06 AM

<b>Instrument</b>	69679	<b>Data File</b>	mj cal 7.d
<b>Type</b>	Cal	<b>Sample</b>	mj cal 7
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Britany Wylie
<b>Sample Position</b>	P3-G1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	4/17/2020 2:34:39 PM		

**Sample Chromatogram**



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.390	2198214	∞	12.0	∞	1136538	99.811 ng/ml
THC-COOH	1.415	2033334	6141.2	174.8	270.2	624043	253.138 ng/ml
THC	3.138	923891	∞	25.2	∞	952709	102.242 ng/ml