






**REVIEWED**  
By Sarah Pickle at 9:05 am, Apr 27, 2020

**Worklist: 4191**

<u>LAB CASE</u>	<u>ITEM</u>	<u>ITEM TYPE</u>	<u>DESCRIPTION</u>	
M2020-0358	3	UCK	AM 27 Urine Cannabinoids Confirmation by LC-QQQ	
M2020-0981	4	UCK	AM 27 Urine Cannabinoids Confirmation by LC-QQQ	
P2020-1042	1	UCK	AM 27 Urine Cannabinoids Confirmation by LC-QQQ	
P2020-1085	1	UCK	AM 27 Urine Cannabinoids Confirmation by LC-QQQ	
P2020-1135	1	UCK	AM 27 Urine Cannabinoids Confirmation by LC-QQQ	

# AM# 27: Quantitation of THC and Metabolites in Blood and Urine by LC-MS/MS

Extraction Date: 04/22/2020

Analyst: Celena Shrum

Plate lot#: IDP-108-2-200303

Plate Expiration: 09/30/2020

**Mobile phase A:** 0.1% Formic Acid in LCMS Water

**Mobile phase B:** 0.1% Formic acid in Acetonitrile

**Blank Blood Lot:** 445283-3

**Column:** UCT Selectra DA 100 x 2.1mm 3um

**LCMS-QQQ ID:** 069901

**Blank Urine Lot:** POC031319

## Pre-Analytic:

- 1. Check levels of mobile phases and needle wash refill as needed. Ensure waste is not full.
- 2. Ensure correct column is installed and begin mobile phase flow allow to equilibrate ~ 30 minutes.

## Analytic:

- 1. Remove standards, plate, controls, and samples from cold storage. Allow to reach room temperature.
- 2. Urine hydrolysis: add 1.5mL urine to blank plate, add 250µl 1N KOH. Shake and incubate at 40 degrees for 15 minutes. Using a calibrated pipette, add **1000µl blood and urine (if applicable) (calibrated pipette)** into the appropriate wells of analytical (standards) plate. **Pipette ID: #42**
- 3. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 4. Pipette **500µL 0.1% formic acid in water blood sample, 500 µL saturated phosphate buffer in urine** in wells of analytical plate.
- 5. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 6. Transfer **800µL of blood+acid or urine+acid** mixture to corresponding wells of SLE+ plate.
- 7. Apply positive pressure for approx. 10-15 seconds (or until no liquid remains on top of sorbent). **(Load at 85-100 PSI- Selector to the right)** Manifold ID: 067104
- 8. Wait 5 minutes.
- 9. Add **2.25mL MTBE. (Add in 3 increments of 750uL)**
- 10. Wait 5 minutes.
- 11. Apply positive pressure for approx. 15 seconds. **(10-15 PSI- Selector to the left).**
- 12. Add **2.25mL Hexane. (Add in 3 increments of 750uL)**
- 13. Wait 5 minutes.
- 14. Apply positive pressure for approx. 15 seconds. **(10-15 PSI- Selector to the left).**
- 15. Remove plate containing eluate. Place on SPE Dry and evaporate to dryness at approx. 35°C. **SPE Dry ID: 067103**
- 16. Reconstitute in **100µL 100% MeOH** and heat seal plate with foil. Place in autosampler and run worklist.

## Post-Analytic

- 1. Create batch and process data.
- 2. Make any necessary integration changes, Curve weighting of Linear 1/x with  $r^2$  values  $\geq 0.98$  for each analyte
- 3. RT +/- 3% or 0.100 min, whichever is greater, +/- 20% Accuracy for greater than (+/- 30% for 10ng/ml or less). Ion ratios must be within +/- 20% of the averaged calibrators
- 4. Case sample response for THC and OH-THC 3ng/mL (quantitative), Carboxy-THC: 10ng/mL (qualitative only) will be reported. Samples with a THC or OH-THC response over 50 ng/mL will be reported out as greater than 50 ng/mL.
- 5. Did all QCs pass for each analyte? (if not, describe in comments section)
- 6. Enter QCs into control charting.
- 7. Central File Packet to include: LIMS Worklist, Method Checklist, Calibration and Control Reports

COMMENTS: Curve ranges are THC: 3-100, THC-COOH: 5-250. THC-OH not evaluated for this run.



# Idaho State Police Forensic Services

## AM #26 Screening of THC and Metabolites and AM #27 Confirmation of THC and Metabolites Urine External Control Prep Sheet

### Methanol External Control Solution (Lot: WS011620)

10  $\mu$ L of 1mg/mL THC, 100  $\mu$ L of 100  $\mu$ g/mL THC-OH, C-THC in 9790  $\mu$ L MeOH  
*Approximate concentration 1ug/mL.*

<i>Component</i>	<i>Source</i>	<i>Source Lot Number</i>	<i>Expiration Date</i>
Methanol (LCMS)	Fisher	193941	
THC	Cerilliant	FE09101501	11/30/2020
C-THC	Cerilliant	FE07171501	09/30/2020
THC-OH	Cerilliant	FE07221601	07/31/2021
Prepared:	01/16/2020		
Prepared By:	Tamara Salazar		
Expires:	09/30/2020		

### Urine External Control Solution (Lot: 042220)

200 ul of methanol external control solution was added to 9800 ul of urine.  
*Approximately 20ng/mL each*

<i>Component</i>	<i>Source</i>	<i>Source Lot Number</i>
Negative Urine	Pocatello Lab	POC031319
Methanol External Control Solution	-	WS011620
Prepared:	04/22/2020	
Prepared by:	Celena Shrum	
Expires:	09/30/2020	


# AM #27 Cannabinoid Quant. Results



**Batch results** D:\MassHunter\Data\2020\AM 27-28\AM 27-28 042220 CS\QuantResults\AM 27 THC-OH removed.batch.bin  
**Calibration Last Update** 4/27/2020 8:43:20 AM

<b>Instrument</b>	Falco	<b>Data File</b>	MJ_QC Control.d
<b>Type</b>	Sample	<b>Sample</b>	MJ_QC Control
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Celena Shrum
<b>Sample Position</b>	P3-A6	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	4/22/2020 3:22:45 PM		
<b>Sample Info.</b>			

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-COOH	1.489	68093	∞	48.4	∞	176193	15.0093 ng/ml
THC	3.375	47125	425.36	32.5	111.00	1226632	4.4008 ng/ml


# AM #27 Cannabinoid Quant. Results



**Batch results** D:\MassHunter\Data\2020\AM 27-28\AM 27-28 042220 CS\QuantResults\AM 27 THC-OH removed.batch.bin  
**Calibration Last Update** 4/27/2020 8:43:20 AM

<b>Instrument</b>	Falco	<b>Data File</b>	MJ_Negative Blood.d
<b>Type</b>	Sample	<b>Sample</b>	MJ_Negative Blood
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Celena Shrum
<b>Sample Position</b>	P3-H5	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	4/22/2020 3:37:57 PM		
<b>Sample Info.</b>			

## Sample Chromatogram




# AM #27 Cannabinoid Quant. Results



**Batch results** D:\MassHunter\Data\2020\AM 27-28\AM 27-28 042220 CS\QuantResults\AM 27 THC-OH removed.batch.bin  
**Calibration Last Update** 4/27/2020 8:43:20 AM

<b>Instrument</b>	Falco	<b>Data File</b>	MJ_Negative Urine.d
<b>Type</b>	Sample	<b>Sample</b>	MJ_Negative Urine
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Celena Shrum
<b>Sample Position</b>	P3-G5	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	4/22/2020 3:53:07 PM		
<b>Sample Info.</b>			

## Sample Chromatogram




# AM #27 Cannabinoid Quant. Results



**Batch results** D:\MassHunter\Data\2020\AM 27-28\AM 27-28 042220 CS\QuantResults\AM 27 THC-OH removed.batch.bin  
**Calibration Last Update** 4/27/2020 8:43:20 AM

<b>Instrument</b>	Falco	<b>Data File</b>	MJ_Urine External Control.d
<b>Type</b>	Sample	<b>Sample</b>	MJ_Urine External Control
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Celena Shrum
<b>Sample Position</b>	P3-F5	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	4/22/2020 4:08:18 PM		
<b>Sample Info.</b>			

## Sample Chromatogram




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-COOH	1.504	58473	∞	55.5	∞	149741	15.1933 ng/ml
THC	3.375	52984	139.78	28.8	35.23	403642	14.3991 ng/ml



# AM #27 Cannabinoids Quant. Calibration Curve Report

**Batch results** D:\MassHunter\Data\2020\AM 27-28\AM 27-28 042220 CS\QuantResults\AM 27 THC-OH removed.batch.bin  
**Last Cal. Update** 4/27/2020 8:43 AM  
**Analyst Name** ISP\Datastor  
**Analyte** THC **Internal Standard** THC-D3




Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
MJ_Cal 1	1	x	1.0	1.2	122.6
MJ_Cal 2	2	✓	3.0	3.0	99.5
MJ_Cal 3	3	✓	5.0	5.2	104.5
MJ_Cal 4	4	✓	10.0	9.9	98.8
MJ_Cal 5	5	✓	25.0	24.5	98.0
MJ_Cal 6	6	✓	50.0	48.8	97.6
MJ_Cal 7	7	✓	100.0	101.6	101.6





# AM #27 Cannabinoids Quant. Calibration Curve Report

**Batch results** D:\MassHunter\Data\2020\AM 27-28\AM 27-28 042220 CS\QuantResults\AM 27 THC-OH removed.batch.bin  
**Last Cal. Update** 4/27/2020 8:43 AM  
**Analyst Name** ISP\Datastor  
**Analyte** THC-COOH **Internal Standard** THC-COOH-D9



Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
MJ_Cal 1	1	✓	5.0	5.4	108.1
MJ_Cal 2	2	✓	10.0	10.4	104.0
MJ_Cal 3	3	✓	20.0	18.8	93.8
MJ_Cal 4	4	✓	50.0	48.4	96.8
MJ_Cal 5	5	✓	75.0	72.2	96.2
MJ_Cal 6	6	✓	100.0	98.5	98.5
MJ_Cal 7	7	✓	250.0	256.4	102.5


# AM #27 Cannabinoid Quant. Results



**Batch results** D:\MassHunter\Data\2020\AM 27-28\AM 27-28 042220 CS\QuantResults\AM 27 THC-OH removed.batch.bin  
**Calibration Last Update** 4/27/2020 8:43:20 AM

<b>Instrument</b>	Falco	<b>Data File</b>	MJ_Cal 1.d
<b>Type</b>	Cal	<b>Sample</b>	MJ_Cal 1
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Celena Shrum
<b>Sample Position</b>	P3-H6	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	4/22/2020 2:21:54 PM		

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-COOH	1.474	32526	∞	45.5	∞	184211	5.4045 ng/ml <b>Low</b>
THC	3.375	11825	∞	30.4	8.13 <b>Low</b>	1323031	1.2262 ng/ml <b>Low</b>


# AM #27 Cannabinoid Quant. Results



**Batch results** D:\MassHunter\Data\2020\AM 27-28\AM 27-28 042220 CS\QuantResults\AM 27 THC-OH removed.batch.bin  
**Calibration Last Update** 4/27/2020 8:43:20 AM

<b>Instrument</b>	Falco	<b>Data File</b>	MJ_Cal 2.d
<b>Type</b>	Cal	<b>Sample</b>	MJ_Cal 2
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Celena Shrum
<b>Sample Position</b>	P3-G6	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	4/22/2020 2:29:39 PM		
<b>Sample Info.</b>			

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-COOH	1.504	50663	∞	46.2	88.26	177314	10.3993 ng/ml
THC	3.390	31059	52.67	35.2	37.46	1229333	2.9844 ng/ml <b>Low</b>


# AM #27 Cannabinoid Quant. Results



**Batch results** D:\MassHunter\Data\2020\AM 27-28\AM 27-28 042220 CS\QuantResults\AM 27 THC-OH removed.batch.bin  
**Calibration Last Update** 4/27/2020 8:43:20 AM

<b>Instrument</b>	Falco	<b>Data File</b>	MJ_Cal 3.d
<b>Type</b>	Cal	<b>Sample</b>	MJ_Cal 3
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Celena Shrum
<b>Sample Position</b>	P3-F6	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	4/22/2020 2:37:13 PM		
<b>Sample Info.</b>			

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-COOH	1.504	80074	∞	52.6	∞	170867	18.7691 ng/ml
THC	3.405	52001	181.79	26.5	30.71	1128983	5.2237 ng/ml


# AM #27 Cannabinoid Quant. Results



**Batch results** D:\MassHunter\Data\2020\AM 27-28\AM 27-28 042220 CS\QuantResults\AM 27 THC-OH removed.batch.bin  
**Calibration Last Update** 4/27/2020 8:43:20 AM

<b>Instrument</b>	Falco	<b>Data File</b>	MJ_Cal 4.d
<b>Type</b>	Cal	<b>Sample</b>	MJ_Cal 4
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Celena Shrum
<b>Sample Position</b>	P3-E6	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	4/22/2020 2:44:47 PM		

**Sample Chromatogram**



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-COOH	1.504	203273	∞	57.3	1200.51	182083	48.4097 ng/ml
THC	3.390	108100	150.12	29.7	197.10	1210602	9.8794 ng/ml


# AM #27 Cannabinoid Quant. Results



**Batch results** D:\MassHunter\Data\2020\AM 27-28\AM 27-28 042220 CS\QuantResults\AM 27 THC-OH removed.batch.bin  
**Calibration Last Update** 4/27/2020 8:43:20 AM

<b>Instrument</b>	Falco	<b>Data File</b>	MJ_Cal 5.d
<b>Type</b>	Cal	<b>Sample</b>	MJ_Cal 5
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Celena Shrum
<b>Sample Position</b>	P3-D6	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	4/22/2020 2:52:22 PM		

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-COOH	1.504	354723	∞	59.2	∞	216888	72.1648 ng/ml
THC	3.390	330163	702.84	26.1	∞	1466622	24.5056 ng/ml


# AM #27 Cannabinoid Quant. Results



**Batch results** D:\MassHunter\Data\2020\AM 27-28\AM 27-28 042220 CS\QuantResults\AM 27 THC-OH removed.batch.bin  
**Calibration Last Update** 4/27/2020 8:43:20 AM

<b>Instrument</b>	Falco	<b>Data File</b>	MJ_Cal 6.d
<b>Type</b>	Cal	<b>Sample</b>	MJ_Cal 6
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Celena Shrum
<b>Sample Position</b>	P3-C6	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	4/22/2020 2:59:58 PM		
<b>Sample Info.</b>			

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-COOH	1.489	405968	∞	58.6	6460.55	183635	98.4868 ng/ml
THC	3.375	583365	3276.78	26.1	164.22	1293675	48.8229 ng/ml


# AM #27 Cannabinoid Quant. Results



**Batch results** D:\MassHunter\Data\2020\AM 27-28\AM 27-28 042220 CS\QuantResults\AM 27 THC-OH removed.batch.bin  
**Calibration Last Update** 4/27/2020 8:43:20 AM

<b>Instrument</b>	Falco	<b>Data File</b>	MJ_Cal 7.d
<b>Type</b>	Cal	<b>Sample</b>	MJ_Cal 7
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Celena Shrum
<b>Sample Position</b>	P3-B6	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	4/22/2020 3:07:32 PM		

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-COOH	1.504	950033	∞	58.7	∞	167823	256.3659 ng/ml
THC	3.390	1168830	4301.33	25.0	124.65	1242254	101.5841 ng/ml