










REVIEWED  
By Briana Wylie at 10:10 am, Jun 15, 2020



6/12/2020

**Worklist: 4301**

<u>LAB CASE</u>	<u>ITEM</u>	<u>ITEM TYPE</u>	<u>DESCRIPTION</u>	
C2020-0959	1	BCK	AM 27 Blood THC Quant by LC-QQQ	
C2020-0982	1	BCK	AM 27 Blood THC Quant by LC-QQQ	
C2020-0994	1	UCK	AM 27 Urine Cannabinoids Confirmation by LC-QQQ	
C2020-0995	1	UCK	AM 27 Urine Cannabinoids Confirmation by LC-QQQ	
C2020-1038	1	BCK	AM 27 Blood THC Quant by LC-QQQ	
C2020-1039	4	BCK	AM 27 Blood THC Quant by LC-QQQ	
C2020-1047	1	BCK	AM 27 Blood THC Quant by LC-QQQ	
C2020-1076	1	BCK	AM 27 Blood THC Quant by LC-QQQ	
C2020-1084	1	BCK	AM 27 Blood THC Quant by LC-QQQ	



# AM# 27: Quantitation of THC and Metabolites in Blood and Urine by LC-MS/MS

Extraction Date: 6/10/2020  
Plate lot#: 200303

Analyst: Anne Nord  
Plate Expiration: 09-03-2020

**Mobile phase A:** 0.1% Formic Acid in LCMS Water MTBE  
**Mobile phase B:** 0.1% Formic acid in Acetonitrile Hexane  
LCMS Methanol

**Blank Blood Lot:** 20A52255 **Urine Blank:** 6920 **Column:** UCT Selectra DA 100 x 2.1mm 3um  
**LCMS-QQQ ID:** 69679

## Pre-Analytic:

- 1. Check levels of mobile phases and needle wash refill as needed. Ensure waste is not full.
- 2. Ensure correct column is installed and begin mobile phase flow allow to equilibrate ~ 30 minutes.

## Analytic:

- 1. Remove standards, plate, controls, and samples from cold storage. Allow to reach room temperature.
- 2. Urine hydrolysis: add 1.5 ml urine to blank plate, add 250 ul 1N KOH mix and incubate at 40 degrees for 15 minutes.  
Pipette 1000µL blood (calibrated pipette) Pipette ID: k52558g in wells of analytical (standards) plate.
- 3. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 4. Pipette 500µL 0.1% formic acid in water blood sample, 500 ul saturated phosphate buffer in urine in wells of analytical plate.
- 5. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 6. Transfer 800µL of blood+acid or urine acid mixture to corresponding wells of SLE+ plate.
- 7. Apply positive pressure for approx. 10-15 seconds (or until no liquid remains on top of sorbent).  
*(Load at 85-100 PSI- Selector to the right)* Manifold ID: 66792
- 8. Wait 5 minutes.
- 9. Add 2.25mL MTBE. *(Add in 3 increments of 750uL)*
- 10. Wait 5 minutes.
- 11. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 12. Add 2.25mL Hexane. *(Add in 3 increments of 750uL)*
- 13. Wait 5 minutes.
- 14. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 15. Remove plate containing eluate. Place on SPE Dry and evaporate to dryness at approx. 35°C.  
*SPE Dry ID: 66819*
- 16. Reconstitute in 100µL 100% MeOH and heat seal plate with foil. Place in autosampler and run worklist.

## Post-Analytic

- 1. Create batch and process data.
- 2. Make any necessary integration changes, Curve weighting of Linear 1/x with  $r^2$  values  $\geq 0.98$  for each analyte
- 3. RT +/- 3% or 0.100 min, whichever is greater, +/- 20% Accuracy for greater than (+/- 30% for 10ng/ml or less).  
Ion ratios must be within +/- 20% of the averaged calibrators
- 4. Case sample response for THC and OH-THC 3ng/mL (quantitative blood), Carboxy-THC: 10ng/mL (qualitative only) will be reported. Samples with a THC or OH-THC response over 50 ng/mL will be reported out as greater than 50 ng/mL.
- 5. Did all QCs pass for each analyte? (if not is it describe in comments section)
- 6. Enter QCs into control charting.
- 7. Central File Packet to include: LIMS Worklist, Method Checklist, Calibration and Control Reports

COMMENTS: *THC curve 3-100*  
*THC-OH not evaluated for urine samples.*



**Toxicology AM method 27/26 external prep information**

working solution 1 ug/ml in meoh C-THC, THC-OH, THC

Stock solution 1mg/ml 7.5 ul each THC, 100 ug/ml 150 ul C-THC, 75 ul THC-OH in 9767.5 ul meOH

Ppd 2/13/20 Exp: 8/13/20 lot 21320 by AMN

Drug	lot	expiration
C-THC	FE07171501	9/1/2020
THC-OH	FE07721601	7/1/2021
THC	FE001041701	3/1/2022

AM 27/26 blood control 100 ul working solution lot (91319) in 9900 ul blood lot (20A52255)

ppd 02/13/20 Exp 08/13/20 lot b81320 Concentration 7.5 ng/ml THC, THC-OH and 15 ng/ml C-THC by AMN

AM 27/26 urine control 400 ul working solution lot (21320) in 9600 ul urine lot (6920)

ppd 4/17/20 Exp 9120	lot u101720	Concentration 30 ng/ml THC, THC-OH and 60 ng/ml C-THC	by BAW	out of use
ppd 6/9/20 exp 8/13/20	lot 6920	Concentration 30 ng/ml THC, THC-OH and 60 ng/ml C-THC	by amn	6/8/2020

# AM #27 Cannabinoids




**Batch results**  
**Calibration Last Update**

D:\MassHunter\Data\2020 Data\am 27-28 6-10-20\QuantResults\mj.batch.bin  
6/11/2020 8:34:00 AM

**Instrument** 69679  
**Type** QC  
**Acq. Method** AM 27 THC quant.m  
**Sample Position** P3-H1  
**Injection Volume** 10  
**Acq. Date-Time** 6/10/2020 12:38:09 PM  
**Sample Info.**

**Data File** mj internal control.d  
**Sample** mj internal control  
**Operator** Anne Nord  
**Comment**

**Sample Chromatogram**




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.390	113407	∞	9.7	∞	1186669	4.985 ng/ml
THC-COOH	1.415	133206	665.1	37.9	1684.1	773839	15.311 ng/ml
THC	3.153	30160	∞	23.4	∞	809359	4.518 ng/ml

# AM #27 Cannabinoids

<b>Batch results</b>	D:\MassHunter\Data\2020 Data\am 27-28 6-10-20\QuantResults\mj.batch.bin		
<b>Calibration Last Update</b>	6/11/2020 8:34:00 AM		
<b>Instrument</b>	69679	<b>Data File</b>	mj negative blood.d
<b>Type</b>	Sample	<b>Sample</b>	mj negative blood
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Anne Nord
<b>Sample Position</b>	P3-A2	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	6/10/2020 12:45:51 PM		
<b>Sample Info.</b>			

## Sample Chromatogram



# AM #27 Cannabinoids

**Batch results**

D:\MassHunter\Data\2020 Data\am 27-28 6-10-20\QuantResults\mj.batch.bin

**Calibration Last Update**

6/11/2020 8:34:00 AM

**Instrument**

69679

**Type**

Sample

**Acq. Method**

AM 27 THC quant.m

**Sample Position**

Vial 4

**Injection Volume**

10

**Acq. Date-Time**

6/10/2020 2:41:15 PM

**Sample Info.****Data File**


mj blank neg urine.d

**Sample**

mj blank neg urine

**Operator**


Anne Nord

**Comment****Sample Chromatogram**

# AM #27 Cannabinoids

<b>Batch results</b>	D:\MassHunter\Data\2020 Data\am 27-28 6-10-20\QuantResults\mj.batch.bin		
<b>Calibration Last Update</b>	6/11/2020 8:34:00 AM		
<b>Instrument</b>	69679	<b>Data File</b>	mj neg urine.d
<b>Type</b>	Sample	<b>Sample</b>	mj neg urine
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Anne Nord
<b>Sample Position</b>	P3-A3	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	6/10/2020 2:48:53 PM		
<b>Sample Info.</b>			

## Sample Chromatogram



# AM #27 Cannabinoids




**Batch results**  
**Calibration Last Update**

D:\MassHunter\Data\2020 Data\am 27-28 6-10-20\QuantResults\mj.batch.bin  
6/11/2020 8:34:00 AM

**Instrument** 69679  
**Type** Sample  
**Acq. Method** AM 27 THC quant.m  
**Sample Position** P3-B3  
**Injection Volume** 10  
**Acq. Date-Time** 6/10/2020 2:56:38 PM  
**Sample Info.**

**Data File** mj Urine external control.d  
**Sample** mj Urine external control  
**Operator** Anne Nord  
**Comment**

**Sample Chromatogram**



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.390	398815	∞	56.8 <b>High</b>	∞	1265314	16.874 ng/ml
THC-COOH	1.415	272270	320.5	37.2	654.1	563156	42.026 ng/ml
THC	3.153	57209	∞	24.5	∞	430973	14.220 ng/ml

Urine evaluated for THC and THC-COOH only.




# Compound Calibration Report




**Batch results** D:\MassHunter\Data\2020 Data\am 27-28 6-10-20\QuantResults\mj.batch.bin  
**Last Cal. Update** 6/11/2020 8:34 AM  
**Analyst Name** ISP\datastor  
**Analyte** THC **Internal Standard** THC-d3

THC - 7 Levels, 6 Levels Used, 7 Points, 6 Points Used, 1 QCs




Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
mj qc1	1	×	1.0	1.6	160.0
mj cal2	2	✓	3.0	3.2	105.1
mj cal3	3	✓	5.0	5.1	102.4
mj cal4	4	✓	10.0	9.3	93.5
mj cal5	5	✓	25.0	24.8	99.4
mj cal6	6	✓	50.0	49.1	98.3
mj cal7	7	✓	100.0	101.4	101.4

# Compound Calibration Report



**Batch results** D:\MassHunter\Data\2020 Data\am 27-28 6-10-20\QuantResults\mj.batch.bin  
**Last Cal. Update** 6/11/2020 8:34 AM  
**Analyst Name** ISP\datastor  
**Analyte** THC-COOH **Internal Standard** THC-COOH-d9




Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
mj qc1	1	✓	5.0	5.2	104.8
mj cal2	2	✓	10.0	10.0	99.7
mj cal 3	3	✓	20.0	20.2	100.9
mj cal 4	4	✓	50.0	48.6	97.2
mj cal 5	5	✓	75.0	72.9	97.2
mj cal 6	6	✓	100.0	98.4	98.4
mj cal 7	7	✓	250.0	254.8	101.9

# Compound Calibration Report



**Batch results** D:\MassHunter\Data\2020 Data\am 27-28 6-10-20\QuantResults\mj.batch.bin  
**Last Cal. Update** 6/11/2020 8:34 AM  
**Analyst Name** ISP\datastor  
**Analyte** THC-OH **Internal Standard** THC-OH-d3

THC-OH - 7 Levels, 7 Levels Used, 7 Points, 7 Points Used, 1 QCs




Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
mj qc1	1	✓	1.0	0.8	81.5
mj cal2	2	✓	3.0	3.4	113.7
mj cal 3	3	✓	5.0	5.4	107.2
mj cal 4	4	✓	10.0	9.9	99.3
mj cal 5	5	✓	25.0	25.0	99.8
mj cal 6	6	✓	50.0	48.9	97.8
mj cal 7	7	✓	100.0	100.6	100.6

# AM #27 Cannabinoids

**Batch results** D:\MassHunter\Data\2020 Data\am 27-28 6-10-20\QuantResults\mj.batch.bin  
**Calibration Last Update** 6/11/2020 8:34:00 AM

<b>Instrument</b>	69679	<b>Data File</b>	mj qc1.d
<b>Type</b>	Cal	<b>Sample</b>	mj qc1
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Anne Nord
<b>Sample Position</b>	P3-A1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	6/10/2020 11:44:01 AM		
<b>Sample Info.</b>			

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.390	24030	∞	11.1	∞	1295639	0.815 ng/ml <b>Low</b>
THC-COOH	1.415	42605	159.8	35.1	236.3	778026	5.240 ng/ml <b>Low</b>
THC	3.153	6358	256.0	28.9	110.1	744521	1.600 ng/ml <b>Low</b>

# AM #27 Cannabinoids




**Batch results**  
**Calibration Last Update**

D:\MassHunter\Data\2020 Data\am 27-28 6-10-20\QuantResults\mj.batch.bin  
6/11/2020 8:34:00 AM

**Instrument** 69679  
**Type** Cal  
**Acq. Method** AM 27 THC quant.m  
**Sample Position** P3-B1  
**Injection Volume** 10  
**Acq. Date-Time** 6/10/2020 11:51:46 AM  
**Sample Info.**

**Data File** mj cal2.d  
**Sample** mj cal2  
**Operator** Anne Nord  
**Comment**

**Sample Chromatogram**



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.390	84162	∞	8.4	∞	1265104	3.412 ng/ml
THC-COOH	1.415	85309	501.1	36.9	375.3	776675	9.966 ng/ml <b>Low</b>
THC	3.153	16789	∞	25.4	69.7	704413	3.154 ng/ml


# AM #27 Cannabinoids



**Batch results** D:\MassHunter\Data\2020 Data\am 27-28 6-10-20\QuantResults\mj.batch.bin  
**Calibration Last Update** 6/11/2020 8:34:00 AM

<b>Instrument</b>	69679	<b>Data File</b>	mj cal 3.d
<b>Type</b>	Cal	<b>Sample</b>	mj cal 3
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Anne Nord
<b>Sample Position</b>	P3-C1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	6/10/2020 11:59:30 AM		
<b>Sample Info.</b>			

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.390	129249	∞	9.9	∞	1261272	5.359 ng/ml
THC-COOH	1.415	180802	226466.2	37.0	565.0	789667	20.187 ng/ml
THC	3.153	31206	∞	27.4	22082445 295712.0	722937	5.118 ng/ml

# AM #27 Cannabinoids

**Batch results**

D:\MassHunter\Data\2020 Data\am 27-28 6-10-20\QuantResults\mj.batch.bin

**Calibration Last Update**

6/11/2020 8:34:00 AM

**Instrument**

69679

**Type**

Cal

**Acq. Method**

AM 27 THC quant.m

**Sample Position**

P3-D1

**Injection Volume**

10

**Acq. Date-Time**

6/10/2020 12:07:14 PM

**Sample Info.****Data File**


mj cal 4.d

**Sample**

mj cal 4

**Operator**

Anne Nord

**Comment****Sample Chromatogram**

Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.390	235140	∞	10.1	∞	1257780	9.932 ng/ml
THC-COOH	1.415	445147	9956.5	38.6	1765.8	795089	48.581 ng/ml
THC	3.153	65649	∞	26.4	∞	774319	9.347 ng/ml

# AM #27 Cannabinoids




**Batch results**  
**Calibration Last Update**

D:\MassHunter\Data\2020 Data\am 27-28 6-10-20\QuantResults\mj.batch.bin  
6/11/2020 8:34:00 AM

**Instrument** 69679  
**Type** Cal  
**Acq. Method** AM 27 THC quant.m  
**Sample Position** P3-E1  
**Injection Volume** 10  
**Acq. Date-Time** 6/10/2020 12:14:58 PM  
**Sample Info.**

**Data File** mj cal 5.d  
**Sample** mj cal 5  
**Operator** Anne Nord  
**Comment**

**Sample Chromatogram**



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.390	537889	∞	11.0	∞	1158262	24.952 ng/ml
THC-COOH	1.400	584958	1215.4	38.9	2940.7	693969	72.868 ng/ml
THC	3.153	150718	∞	25.1	20470718 283739.9	635216	24.841 ng/ml




# AM #27 Cannabinoids



**Batch results** D:\MassHunter\Data\2020 Data\am 27-28 6-10-20\QuantResults\mj.batch.bin  
**Calibration Last Update** 6/11/2020 8:34:00 AM

<b>Instrument</b>	69679	<b>Data File</b>	mj cal 6.d
<b>Type</b>	Cal	<b>Sample</b>	mj cal 6
<b>Acq. Method</b>	AM 27 THC quant.m	<b>Operator</b>	Anne Nord
<b>Sample Position</b>	P3-F1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	6/10/2020 12:22:42 PM		
<b>Sample Info.</b>			

## Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.390	1123282	∞	11.6	∞	1238414	48.915 ng/ml
THC-COOH	1.415	880879	1325.8	39.3	4631.2	772478	98.388 ng/ml
THC	3.153	355600	∞	25.1	∞	746273	49.148 ng/ml

# AM #27 Cannabinoids




**Batch results**  
**Calibration Last Update**

D:\MassHunter\Data\2020 Data\am 27-28 6-10-20\QuantResults\mj.batch.bin  
6/11/2020 8:34:00 AM

**Instrument** 69679  
**Type** Cal  
**Acq. Method** AM 27 THC quant.m  
**Sample Position** P3-G1  
**Injection Volume** 10  
**Acq. Date-Time** 6/10/2020 12:30:27 PM  
**Sample Info.**

**Data File** mj cal 7.d  
**Sample** mj cal 7  
**Operator** Anne Nord  
**Comment**

**Sample Chromatogram**



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.390	2278680	∞	11.7	∞	1223760	100.615 ng/ml
THC-COOH	1.415	2128815	14260.5	39.7	4084.8	718505	254.770 ng/ml
THC	3.153	735859	∞	26.1	∞	742785	101.393 ng/ml