

REVIEWED


By Sarah Pickle at 10:58 am, Aug 28, 2020

8/28/2020

TS

Worklist: 4463

<u>LAB CASE</u>	<u>ITEM</u>	<u>ITEM TYPE</u>	<u>DESCRIPTION</u>
M2020-3036	1	BCK	AM 27 Blood THC Quant by LC-QQQ
M2020-3037	1	BCK	AM 27 Blood THC Quant by LC-QQQ
P2020-1980	1	BCK	AM 27 Blood THC Quant by LC-QQQ
P2020-2273	1	BCK	AM 27 Blood THC Quant by LC-QQQ
P2020-2495	1	BCK	AM 27 Blood THC Quant by LC-QQQ
P2020-2504	1	BCK	AM 27 Blood THC Quant by LC-QQQ
P2020-2506	1	BCK	AM 27 Blood THC Quant by LC-QQQ
P2020-2510	1	BCK	AM 27 Blood THC Quant by LC-QQQ



AM# 27: Quantitation of THC and Metabolites in Blood by LC-MS/MS

Extraction Date: 08/27/2020

Analyst: Tamara Salazar

Plate lot#: IDP-108-200303

Plate Expiration: 09/03/2020

Mobile phase A: 0.1% Formic Acid in LCMS Water
MTBE LCMS Methanol**Mobile phase B:** 0.1% Formic acid in Acetonitrile
Hexane**Blank Blood Lot:** Hemostat 445283-4**Column:** UCT Selectra DA 100 x 2.1mm 3um**LCMS-QQQ ID:** 069901**Pre-Analytic:**

- 1. Check levels of mobile phases and needle wash refill as needed. Ensure waste is not full.
- 2. Ensure correct column is installed and begin mobile phase flow allow to equilibrate ~ 30 minutes.

Analytic:

- 1. Remove standards, plate, controls, and samples from cold storage. Allow to reach room temperature.
- 2. Pipette **1000µL blood (calibrated pipette) Pipette ID: 42** in wells of analytical (standards) plate.
- 3. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 4. Pipette **500µL 0.1% formic acid in water for blood samples** in wells of analytical plate.
- 5. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 6. Transfer **800µL of blood+acid** mixture to corresponding wells of SLE+ plate.
- 7. Apply positive pressure for approx. 10-15 seconds (or until no liquid remains on top of sorbent).
(Load at 85-100 PSI- Selector to the right)
- 8. Wait 5 minutes.
- 9. Add **2.25mL MTBE. (Add in 3 increments of 750uL)**
- 10. Wait 5 minutes.
- 11. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 12. Add **2.25mL Hexane. (Add in 3 increments of 750uL)**
- 13. Wait 5 minutes.
- 14. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 15. Remove plate containing eluate. Place on SPE Dry and evaporate to dryness at approx. 35°C.
- 16. Reconstitute in **100µL 100% MeOH** and heat seal plate with foil. Place in autosampler and run worklist.

Post-Analytic

- 1. Create batch and process data.
- 2. Make any necessary integration changes, Curve weighting of Linear 1/x with r^2 values ≥ 0.98 for each analyte
- 3. RT +/- 3% or 0.100 min, whichever is greater, +/- 20% Accuracy for greater than (+/- 30% for 10ng/ml or less).
Ion ratios must be within +/- 20% of the averaged calibrators
- 4. Case sample response for THC and OH-THC 3ng/mL (quantitative), Carboxy-THC: 5ng/mL (qualitative only) will be reported. Samples with a THC or OH-THC response over 50 ng/mL will be reported out as greater than 50 ng/mL.
- 5. Did all QCs pass for each analyte? Y / N
- 6. Enter QCs into control charting.
- 7. Central File Packet to include: LIMS Worklist, Method Checklist, Calibration and Control Reports

COMMENTS: *Curves limited: THC-OH 3-100*

TS




AM #27 Cannabinoid Quant. Results

Batch results D:\MassHunter\Data\2020\AM 27-28\082720 AM 27 AM 28 P1 and P2 TS\QuantResults\AM 27_THCQ.batch.bin
Calibration Last Update 8/28/2020 10:27:05 AM

Instrument	Falco	Data File	MJ_Negative Blood.d
Type	Sample	Sample	MJ_Negative Blood
Acq. Method	AM 27 THC quant.m	Operator	Tamara Salazar
Sample Position	P3-A2	Comment	
Injection Volume	10		
Acq. Date-Time	8/27/2020 3:02:35 PM		
Sample Info.			

Sample Chromatogram



TS




AM #27 Cannabinoid Quant. Results

Batch results D:\MassHunter\Data\2020\AM 27-28\082720 AM 27 AM 28 P1 and P2 TS\QuantResults\AM 27_THCQ.batch.bin
Calibration Last Update 8/28/2020 10:27:05 AM

Instrument	Falco	Data File	MJ_QC Control.d
Type	Sample	Sample	MJ_QC Control
Acq. Method	AM 27 THC quant.m	Operator	Tamara Salazar
Sample Position	P3-H1	Comment	
Injection Volume	10		
Acq. Date-Time	8/27/2020 2:47:24 PM		

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.468	154233	∞	8.5	217.97	1100631	4.2483 ng/ml
THC-COOH	1.474	119917	∞	58.4	∞	342261	15.1403 ng/ml
THC	3.315	181493	303.18	29.7	100.76	4605121	4.3151 ng/ml



TS


AM #27 Cannabinoids Quant. Calibration Curve Report

Batch results D:\MassHunter\Data\2020\AM 27-28\082720 AM 27 AM 28 P1 and P2 TS\QuantResults\AM 27_THCQ.batch.bin

Last Cal. Update 8/28/2020 10:27 AM

Analyst Name ISP\datastor

Analyte THC **Internal Standard** THC-D3




Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
MJ_Cal 1	1	✓	1.0	1.2	117.1
MJ_Cal 2	2	✓	3.0	2.9	96.2
MJ_Cal 3	3	✓	5.0	4.8	95.6
MJ_Cal 4	4	✓	10.0	9.3	93.3
MJ_Cal 5	5	✓	25.0	24.2	96.7
MJ_Cal 6	6	✓	50.0	49.4	98.7
MJ_Cal 7	7	✓	100.0	102.3	102.3

TS



AM #27 Cannabinoids Quant. Calibration Curve Report

Batch results D:\MassHunter\Data\2020\AM 27-28\082720 AM 27 AM 28 P1 and P2 TS\QuantResults\AM 27_THCQ.batch.bin
Last Cal. Update 8/28/2020 10:27 AM
Analyst Name ISP\datastor
Analyte THC-COOH **Internal Standard** THC-COOH-D9




Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
MJ_Cal 1	1	✓	5.0	5.1	101.5
MJ_Cal 2	2	✓	10.0	9.9	99.4
MJ_Cal 3	3	✓	20.0	20.1	100.5
MJ_Cal 4	4	✓	50.0	49.0	97.9
MJ_Cal 5	5	✓	75.0	75.3	100.4
MJ_Cal 6	6	✓	100.0	100.0	100.0
MJ_Cal 7	7	✓	250.0	250.6	100.2



TS

AM #27 Cannabinoids Quant. Calibration Curve Report

Batch results D:\MassHunter\Data\2020\AM 27-28\082720 AM 27 AM 28 P1 and P2 TS\QuantResults\AM 27_THCQ.batch.bin
Last Cal. Update 8/28/2020 10:27 AM
Analyst Name ISP\datastor
Analyte THC-OH **Internal Standard** THC-OH-D3



Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
MJ_Cal 1	1	x	1.0	1.3	128.4
MJ_Cal 2	2	✓	3.0	2.7	91.3
MJ_Cal 3	3	✓	5.0	5.0	99.6
MJ_Cal 4	4	✓	10.0	11.0	109.6
MJ_Cal 5	5	✓	25.0	25.2	100.7
MJ_Cal 6	6	✓	50.0	49.7	99.4
MJ_Cal 7	7	✓	100.0	99.5	99.5

TS




AM #27 Cannabinoid Quant. Results

Batch results D:\MassHunter\Data\2020\AM 27-28\082720 AM 27 AM 28 P1 and P2 TS\QuantResults\AM 27_THCQ.batch.bin
Calibration Last Update 8/28/2020 10:27:05 AM

Instrument	Falco	Data File	MJ_Cal 1.d
Type	Cal	Sample	MJ_Cal 1
Acq. Method	AM 27 THC quant.m	Operator	Tamara Salazar
Sample Position	P3-A1	Comment	
Injection Volume	10		
Acq. Date-Time	8/27/2020 1:46:37 PM		

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.483	106016	∞	5.5 Low	∞	1228953	1.2844 ng/ml Low
THC-COOH	1.459	39754	∞	57.3	∞	377954	5.0769 ng/ml
THC	3.285	45264	∞	29.3	23.26	4938277	1.1710 ng/ml Low

TS




AM #27 Cannabinoid Quant. Results

Batch results D:\MassHunter\Data\2020\AM 27-28\082720 AM 27 AM 28 P1 and P2 TS\QuantResults\AM 27_THCQ.batch.bin
Calibration Last Update 8/28/2020 10:27:05 AM

Instrument	Falco	Data File	MJ_Cal 2.d
Type	Cal	Sample	MJ_Cal 2
Acq. Method	AM 27 THC quant.m	Operator	Tamara Salazar
Sample Position	P3-B1	Comment	
Injection Volume	10		
Acq. Date-Time	8/27/2020 1:54:21 PM		

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.483	137831	∞	8.5	220.03	1222972	2.7390 ng/ml Low
THC-COOH	1.489	82399	∞	58.8	∞	368529	9.9368 ng/ml
THC	3.330	112792	417.67	26.7	121.81	4396474	2.8851 ng/ml Low

TS




AM #27 Cannabinoid Quant. Results

Batch results D:\MassHunter\Data\2020\AM 27-28\082720 AM 27 AM 28 P1 and P2 TS\QuantResults\AM 27_THCQ.batch.bin
Calibration Last Update 8/28/2020 10:27:05 AM

Instrument	Falco	Data File	MJ_Cal 3.d
Type	Cal	Sample	MJ_Cal 3
Acq. Method	AM 27 THC quant.m	Operator	Tamara Salazar
Sample Position	P3-C1	Comment	
Injection Volume	10		
Acq. Date-Time	8/27/2020 2:01:55 PM		
Sample Info.			

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.468	182998	∞	9.0	∞	1192591	4.9810 ng/ml
THC-COOH	1.489	168655	∞	59.3	∞	357943	20.0990 ng/ml
THC	3.315	171121	714.50	26.7	112.40	3898395	4.7812 ng/ml

TS




AM #27 Cannabinoid Quant. Results

Batch results D:\MassHunter\Data\2020\AM 27-28\082720 AM 27 AM 28 P1 and P2 TS\QuantResults\AM 27_THCQ.batch.bin
Calibration Last Update 8/28/2020 10:27:05 AM

Instrument	Falco	Data File	MJ_Cal 4.d
Type	Cal	Sample	MJ_Cal 4
Acq. Method	AM 27 THC quant.m	Operator	Tamara Salazar
Sample Position	P3-D1	Comment	
Injection Volume	10		
Acq. Date-Time	8/27/2020 2:09:30 PM		

Sample Info.

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.453	309042	∞	9.3	∞	1179431	10.9555 ng/ml
THC-COOH	1.474	414808	∞	59.5	∞	353285	48.9520 ng/ml
THC	3.315	352274	4289.02	26.0	246.00	4016963	9.3345 ng/ml

TS




AM #27 Cannabinoid Quant. Results

Batch results D:\MassHunter\Data\2020\AM 27-28\082720 AM 27 AM 28 P1 and P2 TS\QuantResults\AM 27_THCQ.batch.bin
Calibration Last Update 8/28/2020 10:27:05 AM

Instrument	Falco	Data File	MJ_Cal 5.d
Type	Cal	Sample	MJ_Cal 5
Acq. Method	AM 27 THC quant.m	Operator	Tamara Salazar
Sample Position	P3-E1	Comment	
Injection Volume	10		
Acq. Date-Time	8/27/2020 2:17:04 PM		

Sample Info.

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.453	607278	∞	11.7	480.21	1167030	25.1703 ng/ml
THC-COOH	1.474	626585	∞	60.2	∞	344862	75.3341 ng/ml
THC	3.315	948665	1754.83	25.2	1488.96	4114299	24.1874 ng/ml

TS




AM #27 Cannabinoid Quant. Results

Batch results D:\MassHunter\Data\2020\AM 27-28\082720 AM 27 AM 28 P1 and P2 TS\QuantResults\AM 27_THCQ.batch.bin
Calibration Last Update 8/28/2020 10:27:05 AM

Instrument	Falco	Data File	MJ_Cal 6.d
Type	Cal	Sample	MJ_Cal 6
Acq. Method	AM 27 THC quant.m	Operator	Tamara Salazar
Sample Position	P3-F1	Comment	
Injection Volume	10		
Acq. Date-Time	8/27/2020 2:24:38 PM		

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.453	1124676	∞	12.4	∞	1164300	49.6895 ng/ml
THC-COOH	1.489	825990	∞	60.1	∞	341570	100.0143 ng/ml
THC	3.330	2593448	∞	25.1	2463.28	5484863	49.3711 ng/ml

TS




AM #27 Cannabinoid Quant. Results

Batch results D:\MassHunter\Data\2020\AM 27-28\082720 AM 27 AM 28 P1 and P2 TS\QuantResults\AM 27_THCQ.batch.bin
Calibration Last Update 8/28/2020 10:27:05 AM

Instrument Falco **Data File** MJ_Cal 7.d
Type Cal **Sample** MJ_Cal 7
Acq. Method AM 27 THC quant.m **Operator** Tamara Salazar
Sample Position P3-G1 **Comment**
Injection Volume 10
Acq. Date-Time 8/27/2020 2:32:13 PM
Sample Info.

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.453	2044380	∞	12.6	∞	1092918	99.4647 ng/ml
THC-COOH	1.489	1910789	∞	62.3	∞	313927	250.5868 ng/ml
THC	3.315	3460094	∞	25.9	∞	3524573	102.2696 ng/ml