

REVIEWED





By Celena Shrum at 11:52 am, Feb 08, 2021

2/3/2021

4/2/21- reviewed corrections to control/calibration positions
on am 27 plate map

BW

Worklist: 4771

<u>LAB CASE</u>	<u>ITEM</u>	<u>ITEM TYPE</u>	<u>DESCRIPTION</u>	
C2021-0106	1	BCK	AM 27 Blood THC Quant by LC-QQQ	
C2021-0117	1	UCK	AM 27 Urine Cannabinoids Confirmation by LC-QQQ	
C2021-0143	1	BCK	AM 27 Blood THC Quant by LC-QQQ	
C2021-0216	1	BCK	AM 27 Blood THC Quant by LC-QQQ	

AM# 27: Quantitation of THC and Metabolites in Blood and Urine by LC-MS/MS

Extraction Date 02/02/21
Plate lot#: 201206

Analyst: Anne Nord
Plate Expiration: 06/06/2021

Mobile phase A: 0.1% Formic Acid in LCMS Water
MTBE LCMS Methanol

Mobile phase B: 0.1% Formic acid in Acetonitrile
Hexane

Blank Blood Lot: 20J20793 **Urine Blank:** 10120 **Column:** UCT Selectra DA 100 x 2.1mm 3um
LCMS-QQQ ID: 69679

Pre-Analytic:

- 1. Check levels of mobile phases and needle wash refill as needed. Ensure waste is not full.
- 2. Ensure correct column is installed and begin mobile phase flow allow to equilibrate ~ 30 minutes.

Analytic:

- 1. Remove standards, plate, controls, and samples from cold storage. Allow to reach room temperature.
- 2. Urine hydrolysis: add 1.5 ml urine to blank plate, add 250 ul 1N KOH mix and incubate at 40 degrees for 15 minutes.
Pipette 1000µL blood (calibrated pipette) Pipette ID: k52558g in wells of analytical (standards) plate.
- 3. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 4. Pipette 500µL 0.1% formic acid in water blood sample, 500 ul saturated phosphate buffer in urine in wells of analytical plate.
- 5. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 6. Transfer 800µL of blood+acid or urine acid mixture to corresponding wells of SLE+ plate.
- 7. Apply positive pressure for approx. 10-15 seconds (or until no liquid remains on top of sorbent).
(Load at 85-100 PSI- Selector to the right) Manifold ID: 66792
- 8. Wait 5 minutes.
- 9. Add 2.25mL MTBE. *(Add in 3 increments of 750uL)*
- 10. Wait 5 minutes.
- 11. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 12. Add 2.25mL Hexane. *(Add in 3 increments of 750uL)*
- 13. Wait 5 minutes.
- 14. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 15. Remove plate containing eluate. Place on SPE Dry and evaporate to dryness at approx. 35°C.
SPE Dry ID: 66819
- 16. Reconstitute in 100µL 100% MeOH and heat seal plate with foil. Place in autosampler and run worklist.

Post-Analytic

- 1. Create batch and process data.
- 2. Make any necessary integration changes, Curve weighting of Linear 1/x with r^2 values ≥ 0.98 for each analyte
- 3. RT +/- 3% or 0.100 min, whichever is greater, +/- 20% Accuracy for greater than (+/- 30% for 10ng/ml or less).
Ion ratios must be within +/- 20% of the averaged calibrators
- 4. Case sample response for THC 1ng/ml, OH-THC 3ng/mL (quantitative blood), Carboxy-THC: 5 ng/mL (qualitative only) will be reported. Samples with a THC or OH-THC response over 50 ng/mL will be reported out as greater than 50 ng/mL.
- 5. Did all QCs pass for each analyte? (if not is it describe in comments section)
- 6. Enter QCs into control charting.
- 7. Central File Packet to include: LIMS Worklist, Method Checklist, Calibration and Control Reports

COMMENTS: *THC-OH not evaluated in urine, the ratio in the control is off.*

A

	1	2	3	4	5	6
a	cal 1 cal 100 ng	neg blood				QC 1
b	cal 2 cal 50 ng	106-1				cal 100 ng
c	cal 3 cal 25 ng	143-1				cal 50 ng
d	cal 4 cal 10ng	216-1				cal 25 ng
e	cal 5 cal 5 ng	neg urine				cal 10ng
f	cal 6 cal 3 ng	urine control				cal 5 ng
g	cal 7 cal 1ng	117-1				cal 3 ng
h	QC 1					cal 1ng

C2021-0__-__

Toxicology AM method 27/26 external prep information



working solution 15 ug/ml in meoh C-THC, THC-OH, 7.5 ug/ml THC

Stock solution 1mg/ml 7.5 ul each THC, 100 ug/ml 150 ul C-THC, 150 ul THC-OH in 9692.5 ul meOH

Ppd 8/26/20 Exp: 7/1/21 lot 82620 by AMN

Drug	lot	expiration
C-THC	FE01061702	3/1/2022
THC-OH	FE07221601	7/1/2021
THC	FE01041701	3/1/2022

AM 27/26 blood control 100 ul working solution lot () in 9900 ul blood lot ()

		Concentration 7.5 ng/ml THC, 15 ng/ml C-THC, THC-OH	
--	--	--	--

AM 27/26 urine control 400 ul working solution lot (82620) in 9600 ul urine

out of use

ppd 8/26/20 Exp 7/1/21 neg urine lot 73020	lot u82620	Concentration 30 ng/ml THC, and 60 ng/ml C-THC, THC-OH	by amn	10/4/2020
ppd 10/5/20 Exp 7/1/21 neg urine lot 10120	lot 10520	Concentration 30 ng/ml THC, and 60 ng/ml C-THC, THC-OH	by amn	1/12/2021
ppd 1/13/21 Exp 7/1/21 neg urine lot 10120	lot 11321	Concentration 30 ng/ml THC, and 60 ng/ml C-THC, THC-OH	by amn	


GA

AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2021\am 27-28\020221\QuantResults\mj.batch.bin
Calibration Last Update 2/4/2021 12:06:29 PM

Instrument	69679	Data File	mj internal control blood.d
Type	QC	Sample	mj internal control blood
Acq. Method	AM 27 THC quant.m	Operator	Anne Nord
Sample Position	P3-H1	Comment	
Injection Volume	10		
Acq. Date-Time	2/2/2021 12:01:41 PM		

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.391	111277	∞	9.3	∞	833153	4.793 ng/ml
THC-COOH	1.401	68057	192106.3	34.1	82.8	384701	14.112 ng/ml
THC	3.077	36718	674.0	24.6	74.5	383066	4.052 ng/ml


GA

AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2021\am 27-28\020221\QuantResults\mj.batch.bin
Calibration Last Update 2/4/2021 12:06:29 PM

Instrument	69679	Data File	mj blood neg.d
Type	Sample	Sample	mj blood neg
Acq. Method	AM 27 THC quant.m	Operator	Anne Nord
Sample Position	P3-A2	Comment	
Injection Volume	10		
Acq. Date-Time	2/2/2021 12:08:25 PM		

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.436	27464	∞			889669	0.884 ng/ml Low

GA


AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2021\am 27-28\020221\QuantResults\mj.batch.bin
Calibration Last Update 2/4/2021 12:06:29 PM

Instrument	69679	Data File	mj negative urine.d
Type	Sample	Sample	mj negative urine
Acq. Method	AM 27 THC quant.m	Operator	Anne Nord
Sample Position	P3-E2	Comment	
Injection Volume	10		
Acq. Date-Time	2/2/2021 1:15:19 PM		

Sample Info.

Sample Chromatogram




AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2021\am 27-28\020221\QuantResults\mj.batch.bin
Calibration Last Update 2/4/2021 12:06:29 PM

Instrument	69679	Data File	mj urine control.d
Type	Sample	Sample	mj urine control
Acq. Method	AM 27 THC quant.m	Operator	Anne Nord
Sample Position	P3-F2	Comment	
Injection Volume	10		
Acq. Date-Time	2/2/2021 1:22:03 PM		

Sample Info.

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.391	1494588	∞	14.8 High	∞	1224060	46.192 ng/ml
THC-COOH	1.416	197114	1349.1	34.9	1454.1	408120	36.472 ng/ml
THC	3.137	141398	∞	23.7	1292.3	431653	13.086 ng/ml


Compound Calibration Report



Batch results D:\MassHunter\Data\2021\am 27-28\020221\QuantResults\mj.batch.bin
Last Cal. Update 2/4/2021 12:06 PM
Analyst Name ISP\datastor
Analyte THC

Internal Standard THC-d3

THC - 7 Levels, 7 Levels Used, 7 Points, 7 Points Used, 1 QCs



Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
mj cal1	1	✓	1.0	1.2	116.5
mj cal2	2	✓	3.0	2.9	96.3
mj cal 3	3	✓	5.0	4.8	95.5
mj cal 4	4	✓	10.0	9.6	95.6
mj cal 5	5	✓	25.0	23.6	94.3
mj cal 6	6	✓	50.0	49.7	99.5
mj cal 7	7	✓	100.0	102.3	102.3


Compound Calibration Report



Batch results D:\MassHunter\Data\2021\am 27-28\020221\QuantResults\mj.batch.bin
Last Cal. Update 2/4/2021 12:06 PM
Analyst Name ISP\datastor
Analyte THC-COOH

Internal Standard THC-COOH-d9

THC-COOH - 7 Levels, 7 Levels Used, 7 Points, 7 Points Used, 1 QCs




Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
mj cal1	1	✓	5.0	5.4	108.1
mj cal2	2	✓	10.0	9.4	94.2
mj cal 3	3	✓	20.0	20.1	100.7
mj cal 4	4	✓	50.0	49.4	98.7
mj cal 5	5	✓	75.0	74.0	98.7
mj cal 6	6	✓	100.0	98.1	98.1
mj cal 7	7	✓	250.0	253.5	101.4

AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2021\am 27-28\020221\QuantResults\mj.batch.bin
Calibration Last Update 2/4/2021 12:06:29 PM

Instrument	69679	Data File	mj cal1.d
Type	Cal	Sample	mj cal1
Acq. Method	AM 27 THC quant.m	Operator	Anne Nord
Sample Position	P3-A1	Comment	
Injection Volume	10		
Acq. Date-Time	2/2/2021 11:14:32 AM		

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.391	31089	∞	10.2	∞	974959	0.923 ng/ml Low
THC-COOH	1.416	23132	30.8	37.7	42.8	400604	5.406 ng/ml
THC	3.107	11049	133.4	24.6	29.5	507266	1.165 ng/ml

GA


AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2021\am 27-28\020221\QuantResults\mj.batch.bin
Calibration Last Update 2/4/2021 12:06:29 PM

Instrument	69679	Data File	mj cal2.d
Type	Cal	Sample	mj cal2
Acq. Method	AM 27 THC quant.m	Operator	Anne Nord
Sample Position	P3-B1	Comment	
Injection Volume	10		
Acq. Date-Time	2/2/2021 11:21:18 AM		

Sample Info.

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.391	75814	∞	11.3	∞	913043	2.870 ng/ml
THC-COOH	1.416	44053	56.6	38.0	18087.6	390845	9.422 ng/ml
THC	3.107	29406	175.0	27.2	78.3	445342	2.890 ng/ml


AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2021\am 27-28\020221\QuantResults\mj.batch.bin
Calibration Last Update 2/4/2021 12:06:29 PM

Instrument	69679	Data File	mj cal 3.d
Type	Cal	Sample	mj cal 3
Acq. Method	AM 27 THC quant.m	Operator	Anne Nord
Sample Position	P3-C1	Comment	
Injection Volume	10		
Acq. Date-Time	2/2/2021 11:28:02 AM		

Sample Info.

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.391	142206	∞	9.7	∞	928228	5.541 ng/ml
THC-COOH	1.401	104668	192.2	35.6	56383.7	403508	20.138 ng/ml
THC	3.092	54943	543.2	26.4	402.9	480126	4.777 ng/ml


GA

AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2021\am 27-28\020221\QuantResults\mj.batch.bin
Calibration Last Update 2/4/2021 12:06:29 PM

Instrument	69679	Data File	mj cal 4.d
Type	Cal	Sample	mj cal 4
Acq. Method	AM 27 THC quant.m	Operator	Anne Nord
Sample Position	P3-D1	Comment	
Injection Volume	10		
Acq. Date-Time	2/2/2021 11:34:46 AM		

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.376	267297	∞	10.5	∞	953622	10.380 ng/ml
THC-COOH	1.401	274367	547.2	35.9	1674.3	415966	49.375 ng/ml
THC	3.092	121803	1236.8	25.7	513.7	513671	9.560 ng/ml


GA

AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2021\am 27-28\020221\QuantResults\mj.batch.bin
Calibration Last Update 2/4/2021 12:06:29 PM

Instrument	69679	Data File	mj cal 5.d
Type	Cal	Sample	mj cal 5
Acq. Method	AM 27 THC quant.m	Operator	Anne Nord
Sample Position	P3-E1	Comment	
Injection Volume	10		
Acq. Date-Time	2/2/2021 11:41:29 AM		

Sample Chromatogram




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.376	617471	∞	10.7	∞	943973	24.611 ng/ml
THC-COOH	1.401	407725	525.8	36.7	31127	409010	74.015 ng/ml
THC	3.092	282745	1489.6	24.4	596.8	473958	23.572 ng/ml

AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2021\am 27-28\020221\QuantResults\mj.batch.bin
Calibration Last Update 2/4/2021 12:06:29 PM

Instrument	69679	Data File	mj cal 6.d
Type	Cal	Sample	mj cal 6
Acq. Method	AM 27 THC quant.m	Operator	Anne Nord
Sample Position	P3-F1	Comment	
Injection Volume	10		
Acq. Date-Time	2/2/2021 11:48:13 AM		

Sample Chromatogram




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.376	1227620	∞	11.6	∞	942036	49.319 ng/ml
THC-COOH	1.401	543034	5072.3	37.2	771.2	409188	98.141 ng/ml
THC	3.077	576069	1669.5	24.0	3587.4	454441	49.733 ng/ml

AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2021\am 27-28\020221\QuantResults\mj.batch.bin
Calibration Last Update 2/4/2021 12:06:29 PM

Instrument	69679	Data File	mj cal 7.d
Type	Cal	Sample	mj cal 7
Acq. Method	AM 27 THC quant.m	Operator	Anne Nord
Sample Position	P3-G1	Comment	
Injection Volume	10		
Acq. Date-Time	2/2/2021 11:54:57 AM		

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.376	2476538	∞	11.8	∞	936734	100.357 ng/ml
THC-COOH	1.401	1339998	2591.8	37.8	1801.3	387991	253.503 ng/ml
THC	3.077	1226881	6523.8	23.8	3760.2	468955	102.304 ng/ml