


BW

Worklist: 4789

REVIEWED
By Britany Wylie at 2:29 pm, Feb 17, 2021

<u>LAB CASE</u>	<u>ITEM</u>	<u>ITEM TYPE</u>	<u>DESCRIPTION</u>
C2021-0225	1	BCK	AM 27 Blood THC Quant by LC-QQQ
C2021-0262	1	BCK	AM 27 Blood THC Quant by LC-QQQ
C2021-0281	1	BCK	AM 27 Blood THC Quant by LC-QQQ
C2021-0289	1	BCK	AM 27 Blood THC Quant by LC-QQQ



A

AM# 27: Quantitation of THC and Metabolites in Blood and Urine by LC-MS/MS

Extraction Date 02/11/21
Plate lot#: 201206

Analyst: Anne Nord
Plate Expiration: 06/06/2021

Mobile phase A: 0.1% Formic Acid in LCMS Water MTBE
Mobile phase B: 0.1% Formic acid in Acetonitrile Hexane
LCMS Methanol

Blank Blood Lot: 20J20793 **Urine Blank:** **Column:** UCT Selectra DA 100 x 2.1mm 3um
LCMS-QQ ID: 69679

Pre-Analytic:

- 1. Check levels of mobile phases and needle wash refill as needed. Ensure waste is not full.
- 2. Ensure correct column is installed and begin mobile phase flow allow to equilibrate ~ 30 minutes.

Analytic:

- 1. Remove standards, plate, controls, and samples from cold storage. Allow to reach room temperature.
- 2. Urine hydrolysis: add 1.5 ml urine to blank plate, add 250 ul 1N KOH mix and incubate at 40 degrees for 15 minutes.
Pipette 1000µL blood (calibrated pipette) Pipette ID: k52558g in wells of analytical (standards) plate.
- 3. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 4. Pipette 500µL 0.1% formic acid in water blood sample, 500 ul saturated phosphate buffer in urine in wells of analytical plate.
- 5. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 6. Transfer 800µL of blood+acid or urine acid mixture to corresponding wells of SLE+ plate.
- 7. Apply positive pressure for approx. 10-15 seconds (or until no liquid remains on top of sorbent).
(Load at 85-100 PSI- Selector to the right) Manifold ID: 66792
- 8. Wait 5 minutes.
- 9. Add 2.25mL MTBE. *(Add in 3 increments of 750uL)*
- 10. Wait 5 minutes.
- 11. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 12. Add 2.25mL Hexane. *(Add in 3 increments of 750uL)*
- 13. Wait 5 minutes.
- 14. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 15. Remove plate containing eluate. Place on SPE Dry and evaporate to dryness at approx. 35°C.
SPE Dry ID: 66819
- 16. Reconstitute in 100µL 100% MeOH and heat seal plate with foil. Place in autosampler and run worklist.

Post-Analytic

- 1. Create batch and process data.
- 2. Make any necessary integration changes, Curve weighting of Linear 1/x with r^2 values ≥ 0.98 for each analyte
- 3. RT +/- 3% or 0.100 min, whichever is greater, +/- 20% Accuracy for greater than (+/- 30% for 10ng/ml or less).
Ion ratios must be within +/- 20% of the averaged calibrators
- 4. Case sample response for THC 1ng/ml, OH-THC 3ng/mL (quantitative blood), Carboxy-THC: 5 ng/mL (qualitative only) will be reported. Samples with a THC or OH-THC response over 50 ng/mL will be reported out as greater than 50 ng/mL.
- 5. Did all QCs pass for each analyte? (if not is it describe in comments section)
- 6. Enter QCs into control charting.
- 7. Central File Packet to include: LIMS Worklist, Method Checklist, Calibration and Control Reports

COMMENTS: *Blood only in this run*

Samples originally injected on 2/11/21, response for lower level calcs was poor, samples reconstituted on 2/12/21 and reinjected 2/16/21. The injections on 2/16/21 were evaluated.

GA

	1	2	3	4	5	6
a	cal 1 cal 100 ng	neg blood				QC 1
b	cal 2 cal 50 ng	225-1				cal 100 ng
c	cal 3 cal 25 ng	262-1				cal 50 ng
d	cal 4 cal 10ng	281-1				cal 25 ng
e	cal 5 cal 5 ng	289-1				cal 10ng
f	cal 6 cal 3 ng					cal 5 ng
g	cal 7 cal 1ng					cal 3 ng
h	QC 1					cal 1ng

C2021-0__-__


AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2021\am 27-28\021621\QuantResults\cann quant.batch.bin
Calibration Last Update 2/17/2021 7:46:55 AM

Instrument	69679	Data File	WorklistData-0162.d
Type	QC	Sample	mj internal control blood
Acq. Method	AM 27 THC quant.m	Operator	Anne Nord
Sample Position	P3-H1	Comment	
Injection Volume	10		
Acq. Date-Time	2/16/2021 4:15:37 PM		

Sample Info.

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.394	6015	∞	776.5	∞	78537	4.381 ng/ml
THC-COOH	1.399	37434	∞	23.3	∞	66827	14.614 ng/ml
THC	3.107	58876	4570.4	35.1	195.6	575348	3.817 ng/ml


GA

AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2021\am 27-28\021621\QuantResults\cann quant.batch.bin
Calibration Last Update 2/17/2021 7:46:55 AM

Instrument	69679	Data File	WorklistData-0163.d
Type	Sample	Sample	mj blood neg
Acq. Method	AM 27 THC quant.m	Operator	Anne Nord
Sample Position	P3-A2	Comment	
Injection Volume	10		
Acq. Date-Time	2/16/2021 4:22:21 PM		
Sample Info.			

Sample Chromatogram




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.409	823	∞	1678.5 High	∞	130744	0.618 ng/ml Low

AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2021\am 27-28\021621\QuantResults\cann quant.batch.bin
Calibration Last Update 2/17/2021 7:46:55 AM

Instrument	69679	Data File	mj internal control blood end1.d
Type	QC	Sample	mj internal control blood end
Acq. Method	AM 27 THC quant.m	Operator	Anne Nord
Sample Position	P3-H1	Comment	
Injection Volume	10		
Acq. Date-Time	2/16/2021 5:22:33 PM		
Sample Info.			

Sample Chromatogram




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.379	5123	∞	802.9	∞	66818	4.386 ng/ml
THC-COOH	1.414	33337	∞	22.9	1854.2	59142	14.693 ng/ml
THC	3.107	61964	∞	37.6	∞	578940	3.977 ng/ml

Compound Calibration Report



Batch results D:\MassHunter\Data\2021\am 27-28\021621\QuantResults\cann quant.batch.bin
Last Cal. Update 2/17/2021 7:46 AM
Analyst Name ISP\datastor
Analyte THC **Internal Standard** THC-d3

THC - 7 Levels, 7 Levels Used, 7 Points, 7 Points Used, 2 QCs




Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
mj cal1	1	✓	1.0	1.2	115.9
mj cal2	2	✓	3.0	3.1	102.1
mj cal 3	3	✓	5.0	4.7	94.0
mj cal 4	4	✓	10.0	9.1	91.3
mj cal 5	5	✓	25.0	23.8	95.2
mj cal 6	6	✓	50.0	49.4	98.8
mj cal 7	7	✓	100.0	102.7	102.7

Compound Calibration Report



Batch results D:\MassHunter\Data\2021\am 27-28\021621\QuantResults\cann quant.batch.bin
Last Cal. Update 2/17/2021 7:46 AM
Analyst Name ISP\datastor
Analyte THC-COOH **Internal Standard** THC-COOH-d9

THC-COOH - 7 Levels, 7 Levels Used, 7 Points, 7 Points Used, 2 QCs




Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
mj cal1	1	✓	5.0	6.2	124.7
mj cal2	2	✓	10.0	9.4	94.4
mj cal 3	3	✓	20.0	18.4	91.9
mj cal 4	4	✓	50.0	45.6	91.2
mj cal 5	5	✓	75.0	70.1	93.5
mj cal 6	6	✓	100.0	100.4	100.4
mj cal 7	7	✓	250.0	259.9	103.9

Compound Calibration Report

Batch results D:\MassHunter\Data\2021\am 27-28\021621\QuantResults\cann quant.batch.bin
Last Cal. Update 2/17/2021 7:46 AM
Analyst Name ISP\datastor
Analyte THC-OH **Internal Standard** THC-OH-d3

THC-OH - 7 Levels, 7 Levels Used, 7 Points, 7 Points Used, 2 QCs



Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
mj cal1	1	✓	1.0	1.2	121.8
mj cal2	2	✓	3.0	2.7	91.1
mj cal 3	3	✓	5.0	4.7	93.7
mj cal 4	4	✓	10.0	9.7	96.7
mj cal 5	5	✓	25.0	24.1	96.5
mj cal 6	6	✓	50.0	48.7	97.3
mj cal 7	7	✓	100.0	102.9	102.9


AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2021\am 27-28\021621\QuantResults\cann quant.batch.bin
Calibration Last Update 2/17/2021 7:46:55 AM

Instrument	69679	Data File	WorklistData-0155.d
Type	Cal	Sample	mj cal1
Acq. Method	AM 27 THC quant.m	Operator	Anne Nord
Sample Position	P3-A1	Comment	
Injection Volume	10		
Acq. Date-Time	2/16/2021 3:28:28 PM		

Sample Info.

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.379	1348	∞	739.6	∞	77001	1.218 ng/ml Low
THC-COOH	1.399	12363	∞	23.9	∞	65364	6.233 ng/ml
THC	3.092	9034	345.7	36.7	∞	372407	1.159 ng/ml


GA

AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2021\am 27-28\021621\QuantResults\cann quant.batch.bin
Calibration Last Update 2/17/2021 7:46:55 AM

Instrument	69679	Data File	WorklistData-0156.d
Type	Cal	Sample	mj cal2
Acq. Method	AM 27 THC quant.m	Operator	Anne Nord
Sample Position	P3-B1	Comment	
Injection Volume	10		
Acq. Date-Time	2/16/2021 3:35:13 PM		
Sample Info.			

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.379	3733	∞	733.3	∞	81498	2.733 ng/ml Low
THC-COOH	1.399	24621	∞	26.2	99.9	74316	9.445 ng/ml
THC	3.092	29046	∞	32.8	∞	362200	3.063 ng/ml

GA


AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2021\am 27-28\021621\QuantResults\cann quant.batch.bin
Calibration Last Update 2/17/2021 7:46:55 AM

Instrument	69679	Data File	WorklistData-0157.d
Type	Cal	Sample	mj cal 3
Acq. Method	AM 27 THC quant.m	Operator	Anne Nord
Sample Position	P3-C1	Comment	
Injection Volume	10		
Acq. Date-Time	2/16/2021 3:41:56 PM		

Sample Info.

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.379	9194	∞	782.3	∞	111805	4.683 ng/ml
THC-COOH	1.399	71543	∞	24.5	∞	98455	18.375 ng/ml
THC	3.092	62104	∞	35.3	∞	484328	4.699 ng/ml

GA


AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2021\am 27-28\021621\QuantResults\cann quant.batch.bin
Calibration Last Update 2/17/2021 7:46:55 AM

Instrument	69679	Data File	WorklistData-0158.d
Type	Cal	Sample	mj cal 4
Acq. Method	AM 27 THC quant.m	Operator	Anne Nord
Sample Position	P3-D1	Comment	
Injection Volume	10		
Acq. Date-Time	2/16/2021 3:48:40 PM		

Sample Info.

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.379	10813	∞	616.1	∞	61647	9.672 ng/ml
THC-COOH	1.399	107628	∞	23.6	∞	55685	45.621 ng/ml
THC	3.092	97852	∞	33.1	∞	378861	9.126 ng/ml

GA


AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2021\am 27-28\021621\QuantResults\cann quant.batch.bin
Calibration Last Update 2/17/2021 7:46:55 AM

Instrument	69679	Data File	WorklistData-0159.d
Type	Cal	Sample	mj cal 5
Acq. Method	AM 27 THC quant.m	Operator	Anne Nord
Sample Position	P3-E1	Comment	
Injection Volume	10		
Acq. Date-Time	2/16/2021 3:55:24 PM		

Sample Info.

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.379	30188	∞	612.0	∞	67816	24.114 ng/ml
THC-COOH	1.399	182718	∞	23.2	9680.4	60578	70.095 ng/ml
THC	3.092	265243	∞	34.8	∞	384777	23.802 ng/ml


AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2021\am 27-28\021621\QuantResults\cann quant.batch.bin
Calibration Last Update 2/17/2021 7:46:55 AM

Instrument	69679	Data File	WorklistData-0160.d
Type	Cal	Sample	mj cal 6
Acq. Method	AM 27 THC quant.m	Operator	Anne Nord
Sample Position	P3-F1	Comment	
Injection Volume	10		
Acq. Date-Time	2/16/2021 4:02:08 PM		

Sample Info.

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.379	56558	∞	623.4	∞	62578	48.671 ng/ml
THC-COOH	1.399	229845	∞	22.1	∞	52757	100.375 ng/ml
THC	3.092	603614	∞	31.1	∞	418724	49.413 ng/ml


AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2021\am 27-28\021621\QuantResults\cann quant.batch.bin
Calibration Last Update 2/17/2021 7:46:55 AM

Instrument	69679	Data File	WorklistData-0161.d
Type	Cal	Sample	mj cal 7
Acq. Method	AM 27 THC quant.m	Operator	Anne Nord
Sample Position	P3-G1	Comment	
Injection Volume	10		
Acq. Date-Time	2/16/2021 4:08:52 PM		

Sample Info.

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.379	117829	∞	601.1	∞	61470	102.909 ng/ml
THC-COOH	1.399	565983	∞	22.5	14131.0	49575	259.856 ng/ml
THC	3.092	1443698	∞	31.7	∞	479980	102.738 ng/ml