


Worklist: 4834

REVIEWED

By Brittany Wylie at 12:33 pm, Mar 12, 2021

<u>LAB CASE</u>	<u>ITEM</u>	<u>ITEM TYPE</u>	<u>DESCRIPTION</u>
M2020-5027	1	BCK	AM 27 Blood THC Quant by LC-QQQ
M2021-0533	1	BCK	AM 27 Blood THC Quant by LC-QQQ
M2021-0671	2	BCK	AM 27 Blood THC Quant by LC-QQQ
P2021-0324	1	BCK	AM 27 Blood THC Quant by LC-QQQ
P2021-0390	1	BCK	AM 27 Blood THC Quant by LC-QQQ
P2021-0437	1	BCK	AM 27 Blood THC Quant by LC-QQQ
P2021-0438	1	BCK	AM 27 Blood THC Quant by LC-QQQ
P2021-0471	1	BCK	AM 27 Blood THC Quant by LC-QQQ
P2021-0481	1	BCK	AM 27 Blood THC Quant by LC-QQQ
P2021-0482	1	BCK	AM 27 Blood THC Quant by LC-QQQ
P2021-0513	1	BCK	AM 27 Blood THC Quant by LC-QQQ
P2021-0527	1	BCK	AM 27 Blood THC Quant by LC-QQQ
P2021-0528	1	BCK	AM 27 Blood THC Quant by LC-QQQ



8C

AM# 27: Quantitation of THC and Metabolites in Blood by LC-MS/MS

Extraction Date: 03/09/21
 Plate lot#: IDP-108-2-201206

Analyst: Sarah Collins
 Plate Expiration: 06/06/21

Mobile phase A: 0.1% Formic Acid in LCMS Water
 MTBE
Mobile phase B: 0.1% Formic acid in Acetonitrile
 LCMS Methanol
 Hexane
Blank Blood Lot: Lampire 20L20724
Column: UCT Selectra DA 100 x 2.1mm 3um
LCMS-QQQ ID: 069901

Pre-Analytic:

- 1. Check levels of mobile phases and needle wash refill as needed. Ensure waste is not full.
- 2. Ensure correct column is installed and begin mobile phase flow allow to equilibrate ~ 30 minutes.

Analytic:

- 1. Remove standards, plate, controls, and samples from cold storage. Allow to reach room temperature.
- 2. Pipette **1000 µL blood (calibrated pipette)** in wells of analytical (standards) plate. **Pipette ID: #3382167**
- 3. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 4. Pipette **500 µL 0.1% formic acid in LCMS water** in wells of analytical plate.
- 5. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 6. Transfer **800 µL of blood+acid** mixture to corresponding wells of SLE+ plate.
- 7. Apply positive pressure for approx. 10-15 seconds (or until no liquid remains on top of sorbent).
(Load at 85-95 PSI- Selector to the right)
- 8. Wait 5 minutes.
- 9. Add **2.25 mL MTBE. (Add in 3 increments of 750 µL)**
- 10. Wait 5 minutes.
- 11. Apply positive pressure for approx. 15 seconds. *(12-15 PSI- Selector to the left).*
- 12. Add **2.25 mL Hexane. (Add in 3 increments of 750 µL)**
- 13. Wait 5 minutes.
- 14. Apply positive pressure for approx. 15 seconds. *(12-15 PSI- Selector to the left).*
- 15. Remove plate containing eluate. Place on SPE Dry and evaporate to dryness at approx. 35°C.
- 16. Reconstitute in **100 µL 100% MeOH** and heat seal plate with foil. Place in autosampler and run worklist.

Post-Analytic

- 1. Create batch and process data.
- 2. Make any necessary integration changes, Curve weighting of Linear 1/x with r² values ≥0.98 for each analyte
- 3. RT +/- 3% or 0.100 min, whichever is greater, +/- 20% Accuracy for greater than (+/- 30% for 10ng/ml or less).
 Ion ratios must be within +/- 20% of the averaged calibrators
- 4. Case sample response for THC 1ng/mL and OH-THC 3ng/mL (quantitative), Carboxy-THC: 5ng/mL (qualitative only) will be reported. Samples with a THC or OH-THC response over 50 ng/mL will be reported out as greater than 50 ng/mL.
- 5. Did all QCs pass for each analyte? Y / N
- 6. Enter QCs into control charting.
- 7. Central File Packet to include: LIMS Worklist, Method Checklist, Calibration and Control Reports

COMMENTS: *Curve range limited: THC 3-100 Reconstituted and reinjected Calibrator 1 due retention time shift in initial injection. Reinject data used.*

8C

	1	2	3	4	5	6
A	IS + Cal. 1				p2021-0438-1	IS + QC_1
B	IS + Cal. 2			m2021-0533-1	p2021-0437-1	IS + Cal. 7
C	IS + Cal. 3			p2021-0528-1	p2021-0390-1	IS + Cal. 6
D	IS + Cal. 4			p2021-0527-1	p2021-0324-1	IS + Cal. 5
E	IS + Cal. 5			p2021-0513-1	m2021-0671-2	IS + Cal. 4
F	IS + Cal. 6			p2021-0482-1	m2021-0533-1*	IS + Cal. 3
G	IS + Cal. 7			p2021-0481-1	m2020-5027-1	IS + Cal. 2
H	IS + QC_1			p2021-0471-1	negative	IS + Cal. 1

All wells to contain 100 µl of residual DMSO

*Sample moved during analytical step 6 due to blood clot

**Idaho State Police
Forensic Services
Toxicology Discipline**

Request for Departure from an Analytical Method

Date of Request

02/24/2021

Forensic Scientist

Anne Nord

Analytical Method

Toxicology AM #27: Quantitative Analysis of THC and Metabolites in Blood and Urine by LCMS-QQQ

Request

The method currently reads:

4.3.2.5 If any points are dropped from the approved quantitative range of the curve, the compound will be reported qualitatively. For calibrators and controls 10 ng and below, the accuracy must be within 30%, for calibrators and controls greater than 10 ng/mL the accuracy must be within 20%. If a control falls outside the accuracy range, at the analyst's discretion, the compound may be reported qualitatively.

I would like to add in the following exception:

If the 1ng/ml point is dropped for THC. If the 1 ng/ml point is dropped the quantitative range will be 3-50 ng/ml.

Discipline Leader Review

Departure approved

Comments: This deviation is approved and will remain in effect until it is changed in the actual method.

Departure Not Approved

Comments:



Celena Shrum
Toxicology Discipline Lead
Date: 02/24/2021

8C




AM #27 Cannabinoid Quant. Results

Batch results D:\MassHunter\Data\2021\AM 25-26\030821 AM 25 26 SAEC\QuantResults\AM 27.batch.bin
Calibration Last Update 3/10/2021 8:49:17 AM

Instrument	Instrument 1	Data File	MJQ_Negative Blood.d
Type	Sample	Sample	MJQ_Negative Blood
Acq. Method	AM 27 THCQ.m	Operator	Sarah Collins
Sample Position	P1-H5	Comment	
Injection Volume	10		
Acq. Date-Time	3/9/2021 4:02:24 PM		
Sample Info.			

Sample Chromatogram



8C




AM #27 Cannabinoid Quant. Results

Batch results D:\MassHunter\Data\2021\AM 25-26\030821 AM 25 26 SAEC\QuantResults\AM 27.batch.bin
Calibration Last Update 3/10/2021 8:49:17 AM

Instrument	Instrument 1	Data File	MJQ_QC1.d
Type	Sample	Sample	MJQ_QC1
Acq. Method	AM 27 THCQ.m	Operator	Sarah Collins
Sample Position	P1-A6	Comment	
Injection Volume	10		
Acq. Date-Time	3/9/2021 4:17:38 PM		
Sample Info.			

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.438	122672	∞	10.4	∞	2006622	3.9432 ng/ml
THC-COOH	1.459	184769	720.90	54.4	360.51	537398	13.7469 ng/ml
THC	3.224	108881	∞	32.5	∞	2969039	4.4373 ng/ml

8C




AM #27 Cannabinoid Quant. Results

Batch results D:\MassHunter\Data\2021\AM 25-26\030821 AM 25 26 SAEC\QuantResults\AM 27.batch.bin
Calibration Last Update 3/10/2021 8:49:17 AM

Instrument	Instrument 1	Data File	MJQ_QC-End.d
Type	Sample	Sample	MJQ_QC-End
Acq. Method	AM 27 THCQ.m	Operator	Sarah Collins
Sample Position	P1-A6	Comment	
Injection Volume	10		
Acq. Date-Time	3/9/2021 7:50:47 PM		

Sample Chromatogram




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.438	131261	∞	10.6	∞	2187842	3.8536 ng/ml
THC-COOH	1.459	200398	334.04	55.2	259.56	583757	13.7244 ng/ml
THC	3.209	117438	764.14	34.1	∞	3467077	4.1291 ng/ml

8C



AM #27 Cannabinoids Quant. Calibration Curve Report

Batch results D:\MassHunter\Data\2021\AM 25-26\030821 AM 25 26 SAEC\QuantResults\AM 27.batch.bin
Last Cal. Update 3/12/2021 9:16 AM
Analyst Name ISP\Datastor
Analyte THC **Internal Standard** THC-D3



Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
MJQ_Cal 1R	1	x	1.0	1.4	141.5
MJQ_Cal 2	2	✓	3.0	3.0	100.8
MJQ_Cal 3	3	✓	5.0	5.1	102.5
MJQ_Cal 4	4	✓	10.0	9.9	99.2
MJQ_Cal 5	5	✓	25.0	24.3	97.1
MJQ_Cal 6	6	✓	50.0	49.8	99.5
MJQ_Cal 7	7	✓	100.0	100.9	100.9

Cal 1 dropped due to not meeting retention time requirement


8C



AM #27 Cannabinoids Quant. Calibration Curve Report

Batch results D:\MassHunter\Data\2021\AM 25-26\030821 AM 25 26 SAEC\QuantResults\AM 27.batch.bin
Last Cal. Update 3/12/2021 9:16 AM
Analyst Name ISP\Datastor
Analyte THC-COOH **Internal Standard** THC-COOH-D9

THC-COOH - 7 Levels, 7 Levels Used, 7 Points, 7 Points Used, 0 QCs




Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
MJQ_Cal 1R	1	✓	5.0	5.0	100.3
MJQ_Cal 2	2	✓	10.0	10.0	100.2
MJQ_Cal 3	3	✓	20.0	20.1	100.7
MJQ_Cal 4	4	✓	50.0	49.5	98.9
MJQ_Cal 5	5	✓	75.0	74.8	99.7
MJQ_Cal 6	6	✓	100.0	99.9	99.9
MJQ_Cal 7	7	✓	250.0	250.7	100.3

8C



AM #27 Cannabinoids Quant. Calibration Curve Report

Batch results D:\MassHunter\Data\2021\AM 25-26\030821 AM 25 26 SAEC\QuantResults\AM 27.batch.bin
Last Cal. Update 3/12/2021 9:16 AM
Analyst Name ISP\Datastor
Analyte THC-OH **Internal Standard** THC-OH-D3



Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
MJQ_Cal 1R	1	✓	1.0	1.1	107.3
MJQ_Cal 2	2	✓	3.0	2.8	94.4
MJQ_Cal 3	3	✓	5.0	5.0	99.2
MJQ_Cal 4	4	✓	10.0	10.1	100.9
MJQ_Cal 5	5	✓	25.0	24.4	97.7
MJQ_Cal 6	6	✓	50.0	49.8	99.6
MJQ_Cal 7	7	✓	100.0	100.8	100.8

8C




AM #27 Cannabinoid Quant. Results

Batch results D:\MassHunter\Data\2021\AM 25-26\030821 AM 25 26 SAEC\QuantResults\AM 27.batch.bin
Calibration Last Update 3/10/2021 8:49:17 AM

Instrument Instrument 1
Type Cal
Acq. Method AM 27 THCQ.m
Sample Position P1-H6
Injection Volume 10
Acq. Date-Time 3/9/2021 3:01:25 PM
Sample Info. Due to retention time shifts, sample was reconstituted and reinjected. See reinject data.

Data File MJQ_Cal 1.d
Sample MJQ_Cal 1
Operator Sarah Collins
Comment

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.468 High	55799	∞	7.4 Low	19.87	2301293	1.0376 ng/ml Low
THC-COOH	1.474	94390	∞	45.5	156.61	596652	5.9013 ng/ml
THC	3.315	35844	25.43	69.7 High	117.49	4052319	1.3742 ng/ml

8C




AM #27 Cannabinoid Quant. Results

Batch results D:\MassHunter\Data\2021\AM 25-26\030821 AM 25 26 SAEC\QuantResults\AM 27.batch.bin
Calibration Last Update 3/10/2021 8:49:17 AM

Instrument	Instrument 1	Data File	MJQ_Cal 1R.d
Type	Cal	Sample	MJQ_Cal 1R
Acq. Method	AM 27 THCQ.m	Operator	Sarah Collins
Sample Position	P1-H6	Comment	
Injection Volume	10		
Acq. Date-Time	3/10/2021 8:39:17 AM		
Sample Info.			

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.438	45284	∞	8.3	28.17	1856312	1.0492 ng/ml Low
THC-COOH	1.429	61783	185.97	49.2	337.37	495120	4.4890 ng/ml Low
THC	3.179 Low	28220	10.83	29.1	∞	3061419	1.4152 ng/ml

8C




AM #27 Cannabinoid Quant. Results

Batch results D:\MassHunter\Data\2021\AM 25-26\030821 AM 25 26 SAEC\QuantResults\AM 27.batch.bin
Calibration Last Update 3/10/2021 8:49:17 AM

Instrument	Instrument 1	Data File	MJQ_Cal 2.d
Type	Cal	Sample	MJQ_Cal 2
Acq. Method	AM 27 THCQ.m	Operator	Sarah Collins
Sample Position	P1-G6	Comment	
Injection Volume	10		
Acq. Date-Time	3/9/2021 3:09:11 PM		

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.453	108920	∞	9.6	78.77	2330087	2.8098 ng/ml Low
THC-COOH	1.474	146946	∞	52.2	655.52	602390	9.5252 ng/ml
THC	3.285	90707	138.22	33.1	46.84	3804802	3.0248 ng/ml

8C




AM #27 Cannabinoid Quant. Results

Batch results D:\MassHunter\Data\2021\AM 25-26\030821 AM 25 26 SAEC\QuantResults\AM 27.batch.bin
Calibration Last Update 3/10/2021 8:49:17 AM

Instrument	Instrument 1	Data File	MJQ_Cal 3.d
Type	Cal	Sample	MJQ_Cal 3
Acq. Method	AM 27 THCQ.m	Operator	Sarah Collins
Sample Position	P1-F6	Comment	
Injection Volume	10		
Acq. Date-Time	3/9/2021 3:16:47 PM		

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.438	173221	69.88	10.2	∞	2348201	4.9384 ng/ml
THC-COOH	1.459	291460	∞	56.8	2939.30	601615	19.6913 ng/ml
THC	3.254	161162	260.54	30.1	48.36	3756944	5.1226 ng/ml

8C




AM #27 Cannabinoid Quant. Results

Batch results D:\MassHunter\Data\2021\AM 25-26\030821 AM 25 26 SAEC\QuantResults\AM 27.batch.bin
Calibration Last Update 3/10/2021 8:49:17 AM

Instrument	Instrument 1	Data File	MJQ_Cal 4.d
Type	Cal	Sample	MJQ_Cal 4
Acq. Method	AM 27 THCQ.m	Operator	Sarah Collins
Sample Position	P1-E6	Comment	
Injection Volume	10		
Acq. Date-Time	3/9/2021 3:24:24 PM		

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.423	330458	∞	10.3	374.66	2379222	10.0685 ng/ml
THC-COOH	1.459	702817	∞	57.8	683.94	594220	49.2056 ng/ml
THC	3.224	324379	∞	32.0	∞	3752124	9.9171 ng/ml

8C




AM #27 Cannabinoid Quant. Results

Batch results D:\MassHunter\Data\2021\AM 25-26\030821 AM 25 26 SAEC\QuantResults\AM 27.batch.bin
Calibration Last Update 3/10/2021 8:49:17 AM

Instrument	Instrument 1	Data File	MJQ_Cal 5.d
Type	Cal	Sample	MJQ_Cal 5
Acq. Method	AM 27 THCQ.m	Operator	Sarah Collins
Sample Position	P1-D6	Comment	
Injection Volume	10		
Acq. Date-Time	3/9/2021 3:32:00 PM		
Sample Info.			

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.423	798670	∞	11.1	∞	2486918	24.4249 ng/ml
THC-COOH	1.459	1085447	2194.08	58.5	4260.65	608271	74.6383 ng/ml
THC	3.224	814294	∞	28.1	∞	3753490	24.2814 ng/ml

8C




AM #27 Cannabinoid Quant. Results

Batch results D:\MassHunter\Data\2021\AM 25-26\030821 AM 25 26 SAEC\QuantResults\AM 27.batch.bin
Calibration Last Update 3/10/2021 8:49:17 AM

Instrument	Instrument 1	Data File	MJQ_Cal 6.d
Type	Cal	Sample	MJQ_Cal 6
Acq. Method	AM 27 THCQ.m	Operator	Sarah Collins
Sample Position	P1-C6	Comment	
Injection Volume	10		
Acq. Date-Time	3/9/2021 3:39:35 PM		
Sample Info.			

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.423	1584917	∞	11.4	∞	2462232	49.8319 ng/ml
THC-COOH	1.459	1438075	2288.72	59.0	∞	603707	99.8965 ng/ml
THC	3.224	1701278	∞	24.6	1584.53	3793197	49.7717 ng/ml

8C




AM #27 Cannabinoid Quant. Results

Batch results D:\MassHunter\Data\2021\AM 25-26\030821 AM 25 26 SAEC\QuantResults\AM 27.batch.bin
Calibration Last Update 3/10/2021 8:49:17 AM

Instrument	Instrument 1	Data File	MJQ_Cal 7.d
Type	Cal	Sample	MJQ_Cal 7
Acq. Method	AM 27 THCQ.m	Operator	Sarah Collins
Sample Position	P1-B6	Comment	
Injection Volume	10		
Acq. Date-Time	3/9/2021 3:47:11 PM		
Sample Info.			

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.423	3084816	∞	11.6	∞	2389046	100.8396 ng/ml
THC-COOH	1.459	3395415	26112.38	58.8	5367.90	568503	251.6529 ng/ml
THC	3.224	3335490	∞	24.5	∞	3654060	100.8824 ng/ml