


Worklist: 4890

<u>LAB CASE</u>	<u>ITEM</u>	<u>ITEM TYPE</u>	<u>DESCRIPTION</u>
M2021-1048	2	BCK	AM 27 Blood THC Quant by LC-QQQ
M2021-1387	3	BCK	AM 27 Blood THC Quant by LC-QQQ
P2021-0895	1	BCK	AM 27 Blood THC Quant by LC-QQQ
P2021-0940	1	BCK	AM 27 Blood THC Quant by LC-QQQ
P2021-0956	1	BCK	AM 27 Blood THC Quant by LC-QQQ
P2021-0957	1	BCK	AM 27 Blood THC Quant by LC-QQQ
P2021-0958	1	BCK	AM 27 Blood THC Quant by LC-QQQ
P2021-0960	1	BCK	AM 27 Blood THC Quant by LC-QQQ
P2021-0981	1	BCK	AM 27 Blood THC Quant by LC-QQQ
P2021-0990	1	BCK	AM 27 Blood THC Quant by LC-QQQ
P2021-0991	1	BCK	AM 27 Blood THC Quant by LC-QQQ
P2021-1015	1	BCK	AM 27 Blood THC Quant by LC-QQQ



80

AM# 27: Quantitation of THC and Metabolites in Blood by LC-MS/MS

Extraction Date: 04/07/21
 Plate lot#: IDP-108-2-201206

Analyst: Sarah Collins
 Plate Expiration: 06/06/21

Mobile phase A: 0.1% Formic Acid in LCMS Water
 MTBE
Mobile phase B: 0.1% Formic acid in Acetonitrile
 LCMS Methanol Hexane
Blank Blood Lot: Lampire 20L20724
LCMS-QQQ ID: 069901
Column: UCT Selectra DA 100 x 2.1mm 3um

Pre-Analytic:

- 1. Check levels of mobile phases and needle wash refill as needed. Ensure waste is not full.
- 2. Ensure correct column is installed and begin mobile phase flow allow to equilibrate ~ 30 minutes.

Analytic:

- 1. Remove standards, plate, controls, and samples from cold storage. Allow to reach room temperature.
- 2. Pipette **1000 µL blood (calibrated pipette)** in wells of analytical (standards) plate. **Pipette ID: #3382167**
- 3. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 4. Pipette **500 µL 0.1% formic acid in LCMS water** in wells of analytical plate.
- 5. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 6. Transfer **800 µL of blood+acid** mixture to corresponding wells of SLE+ plate.
- 7. Apply positive pressure for approx. 10-15 seconds (or until no liquid remains on top of sorbent).
(Load at 85-95 PSI- Selector to the right)
- 8. Wait 5 minutes.
- 9. Add **2.25 mL MTBE. (Add in 3 increments of 750 µL)**
- 10. Wait 5 minutes.
- 11. Apply positive pressure for approx. 15 seconds. *(12-15 PSI- Selector to the left).*
- 12. Add **2.25 mL Hexane. (Add in 3 increments of 750 µL)**
- 13. Wait 5 minutes.
- 14. Apply positive pressure for approx. 15 seconds. *(12-15 PSI- Selector to the left).*
- 15. Remove plate containing eluate. Place on SPE Dry and evaporate to dryness at approx. 35°C.
- 16. Reconstitute in **100 µL 100% MeOH** and heat seal plate with foil. Place in autosampler and run worklist.

Post-Analytic

- 1. Create batch and process data.
- 2. Make any necessary integration changes, Curve weighting of Linear 1/x with r² values ≥0.98 for each analyte
- 3. RT +/- 3% or 0.100 min, whichever is greater, +/- 20% Accuracy for greater than (+/- 30% for 10ng/ml or less).
 Ion ratios must be within +/- 20% of the averaged calibrators
- 4. Case sample response for THC 1ng/mL and OH-THC 3ng/mL (quantitative), Carboxy-THC: 5ng/mL (qualitative only) will be reported. Samples with a THC or OH-THC response over 50 ng/mL will be reported out as greater than 50 ng/mL.
- 5. Did all QCs pass for each analyte? Y / N
- 6. Enter QCs into control charting.
- 7. Central File Packet to include: LIMS Worklist, Method Checklist, Calibration and Control Reports

COMMENTS: *Curve range limited: THC 3-100 THC-OH 3-100
 Reinject Calibrator 6 due low internal standard response in initial injection. Reinject data used.*

**Idaho State Police
Forensic Services
Toxicology Discipline**

Request for Departure from an Analytical Method

Date of Request

02/24/2021

Forensic Scientist

Anne Nord

Analytical Method

Toxicology AM #27: Quantitative Analysis of THC and Metabolites in Blood and Urine by LCMS-QQQ

Request

The method currently reads:

4.3.2.5 If any points are dropped from the approved quantitative range of the curve, the compound will be reported qualitatively. For calibrators and controls 10 ng and below, the accuracy must be within 30%, for calibrators and controls greater than 10 ng/mL the accuracy must be within 20%. If a control falls outside the accuracy range, at the analyst's discretion, the compound may be reported qualitatively.

I would like to add in the following exception:

If the 1ng/ml point is dropped for THC. If the 1 ng/ml point is dropped the quantitative range will be 3-50 ng/ml.

Discipline Leader Review

Departure approved

Comments: This deviation is approved and will remain in effect until it is changed in the actual method.

Departure Not Approved

Comments:



Celena Shrum

Toxicology Discipline Lead

Date: 02/24/2021

80

	1	2	3	4	5	6
A	IS + Cal. 1	negative	p2021-0960-1			IS + QC_1
B	IS + Cal. 2	m2021-1048-2	p2021-0981-1			IS + Cal. 7
C	IS + Cal. 3	m2021-1387-3	p2021-0990-1			IS + Cal. 6
D	IS + Cal. 4	p2021-0895-1	p2021-0991-1			IS + Cal. 5
E	IS + Cal. 5	p2021-0940-1	p2021-1015-1*			IS + Cal. 4
F	IS + Cal. 6	p2021-0956-1	p2021-1015-1			IS + Cal. 3
G	IS + Cal. 7	p2021-0957-1				IS + Cal. 2
H	IS + QC_1	p2021-0958-1				IS + Cal. 1

All wells to contain 100 µl of residual DMSO

*Sample moved during analytical step 6 due to blood clot

SC




AM #27 Cannabinoid Quant. Results

Batch results D:\MassHunter\Data\2021\AM 27-28\040721 AM 27 28 SC\QuantResults\AM 27.batch.bin
Calibration Last Update 4/8/2021 9:07:24 AM

Instrument	Instrument 1	Data File	MJQ_Negative Blood.d
Type	Sample	Sample	MJQ_Negative Blood
Acq. Method	AM 27 THCQ.m	Operator	Sarah Collins
Sample Position	P1-A2	Comment	
Injection Volume	10		
Acq. Date-Time	4/7/2021 12:09:58 PM		
Sample Info.			

Sample Chromatogram



SC




AM #27 Cannabinoid Quant. Results

Batch results D:\MassHunter\Data\2021\AM 27-28\040721 AM 27 28 SC\QuantResults\AM 27.batch.bin
Calibration Last Update 4/8/2021 9:07:24 AM

Instrument	Instrument 1	Data File	MJQ_QC Start.d
Type	Sample	Sample	MJQ_QC Start
Acq. Method	AM 27 THCQ.m	Operator	Sarah Collins
Sample Position	P1-H1	Comment	
Injection Volume	10		
Acq. Date-Time	4/7/2021 12:25:10 PM		
Sample Info.			

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.468	181094	∞	8.2	59.57	2360605	4.0092 ng/ml
THC-COOH	1.489	196710	∞	54.0	606.05	523262	14.6748 ng/ml
THC	3.345	117510	493.75	30.3	44.75	3436313	4.0773 ng/ml

SC




AM #27 Cannabinoid Quant. Results

Batch results D:\MassHunter\Data\2021\AM 27-28\040721 AM 27 28 SC\QuantResults\AM 27.batch.bin
Calibration Last Update 4/8/2021 9:07:24 AM

Instrument	Instrument 1	Data File	MJQ_QC End.d
Type	Sample	Sample	MJQ_QC End
Acq. Method	AM 27 THCQ.m	Operator	Sarah Collins
Sample Position	P1-H1	Comment	
Injection Volume	10		
Acq. Date-Time	4/7/2021 3:50:32 PM		

Sample Chromatogram




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.483	207121	∞	8.3	78.71	2837023	3.7189 ng/ml
THC-COOH	1.504	234805	215.22	55.2	699.63	659112	13.8907 ng/ml
THC	3.345	146658	702.07	32.4	∞	4232621	4.1258 ng/ml

SC



AM #27 Cannabinoids Quant. Calibration Curve Report

Batch results D:\MassHunter\Data\2021\AM 27-28\040721 AM 27 28 SC\QuantResults\AM 27.batch.bin
Last Cal. Update 4/8/2021 9:07 AM
Analyst Name ISP\Datastor
Analyte THC **Internal Standard** THC-D3



Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
MJQ_Cal 1	1	x	1.0	1.2	121.2
MJQ_Cal 2	2	✓	3.0	3.0	101.2
MJQ_Cal 3	3	✓	5.0	5.1	101.5
MJQ_Cal 4	4	✓	10.0	10.1	100.5
MJQ_Cal 5	5	✓	25.0	24.0	96.0
MJQ_Cal 6R	6	✓	50.0	49.9	99.7
MJQ_Cal 7	7	✓	100.0	101.0	101.0


Calibrator 1 dropped due to not meeting ratio requirement

SC



AM #27 Cannabinoids Quant. Calibration Curve Report

Batch results D:\MassHunter\Data\2021\AM 27-28\040721 AM 27 28 SC\QuantResults\AM 27.batch.bin
Last Cal. Update 4/8/2021 9:07 AM
Analyst Name ISP\Datastor
Analyte THC-COOH **Internal Standard** THC-COOH-D9




Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
MJQ_Cal 1	1	✓	5.0	5.2	103.1
MJQ_Cal 2	2	✓	10.0	9.9	98.7
MJQ_Cal 3	3	✓	20.0	20.9	104.3
MJQ_Cal 4	4	✓	50.0	48.3	96.7
MJQ_Cal 5	5	✓	75.0	73.1	97.5
MJQ_Cal 6R	6	✓	100.0	97.9	97.9
MJQ_Cal 7	7	✓	250.0	254.8	101.9

SC



AM #27 Cannabinoids Quant. Calibration Curve Report

Batch results D:\MassHunter\Data\2021\AM 27-28\040721 AM 27 28 SC\QuantResults\AM 27.batch.bin
Last Cal. Update 4/8/2021 9:07 AM
Analyst Name ISP\Datastor
Analyte THC-OH **Internal Standard** THC-OH-D3



Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
MJQ_Cal 1	1	x	1.0	0.7	70.0
MJQ_Cal 2	2	✓	3.0	2.9	96.4
MJQ_Cal 3	3	✓	5.0	5.0	100.9
MJQ_Cal 4	4	✓	10.0	10.5	104.7
MJQ_Cal 5	5	✓	25.0	24.5	97.8
MJQ_Cal 6R	6	✓	50.0	50.0	100.0
MJQ_Cal 7	7	✓	100.0	100.1	100.1

Calibrator 1 dropped due to not meeting ratio requirement

SC




AM #27 Cannabinoid Quant. Results

Batch results D:\MassHunter\Data\2021\AM 27-28\040721 AM 27 28 SC\QuantResults\AM 27.batch.bin
Calibration Last Update 4/8/2021 9:07:24 AM

Instrument	Instrument 1	Data File	MJQ_Cal 1.d
Type	Cal	Sample	MJQ_Cal 1
Acq. Method	AM 27 THCQ.m	Operator	Sarah Collins
Sample Position	P1-A1	Comment	
Injection Volume	10		
Acq. Date-Time	4/7/2021 11:09:04 AM		

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.513 High	85761	25.31	5.3 Low	27.80	2490498	0.6998 ng/ml Low
THC-COOH	1.474	74871	∞	49.3	∞	546813	5.1550 ng/ml
THC	3.315	26172	65.33	53.5 High	34.45	3500706	1.2119 ng/ml

SC




AM #27 Cannabinoid Quant. Results

Batch results D:\MassHunter\Data\2021\AM 27-28\040721 AM 27 28 SC\QuantResults\AM 27.batch.bin
Calibration Last Update 4/8/2021 9:07:24 AM

Instrument	Instrument 1	Data File	MJQ_Cal 2.d
Type	Cal	Sample	MJQ_Cal 2
Acq. Method	AM 27 THCQ.m	Operator	Sarah Collins
Sample Position	P1-B1	Comment	
Injection Volume	10		
Acq. Date-Time	4/7/2021 11:16:49 AM		

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.468	144933	∞	8.0	52.01	2320574	2.8930 ng/ml Low
THC-COOH	1.489	138880	∞	51.8	711.18	544043	9.8690 ng/ml
THC	3.345	81136	305.03	33.2	61.53	3311924	3.0372 ng/ml

SC




AM #27 Cannabinoid Quant. Results

Batch results D:\MassHunter\Data\2021\AM 27-28\040721 AM 27 28 SC\QuantResults\AM 27.batch.bin
Calibration Last Update 4/8/2021 9:07:24 AM

Instrument	Instrument 1	Data File	MJQ_Cal 3.d
Type	Cal	Sample	MJQ_Cal 3
Acq. Method	AM 27 THCQ.m	Operator	Sarah Collins
Sample Position	P1-C1	Comment	
Injection Volume	10		
Acq. Date-Time	4/7/2021 11:24:24 AM		
Sample Info.			

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.468	219201	50.65	8.3	304.62	2437274	5.0441 ng/ml
THC-COOH	1.489	280564	∞	54.2	∞	528317	20.8535 ng/ml
THC	3.330	145616	626.98	31.3	∞	3345786	5.0773 ng/ml

SC




AM #27 Cannabinoid Quant. Results

Batch results D:\MassHunter\Data\2021\AM 27-28\040721 AM 27 28 SC\QuantResults\AM 27.batch.bin
Calibration Last Update 4/8/2021 9:07:24 AM

Instrument	Instrument 1	Data File	MJQ_Cal 4.d
Type	Cal	Sample	MJQ_Cal 4
Acq. Method	AM 27 THCQ.m	Operator	Sarah Collins
Sample Position	P1-D1	Comment	
Injection Volume	10		
Acq. Date-Time	4/7/2021 11:32:00 AM		

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.468	372639	∞	9.4	303.60	2339864	10.4700 ng/ml
THC-COOH	1.489	653320	1464.58	56.2	∞	535138	48.3283 ng/ml
THC	3.345	289509	∞	27.3	∞	3219465	10.0533 ng/ml

SC




AM #27 Cannabinoid Quant. Results

Batch results D:\MassHunter\Data\2021\AM 27-28\040721 AM 27 28 SC\QuantResults\AM 27.batch.bin
Calibration Last Update 4/8/2021 9:07:24 AM

Instrument	Instrument 1	Data File	MJQ_Cal 5.d
Type	Cal	Sample	MJQ_Cal 5
Acq. Method	AM 27 THCQ.m	Operator	Sarah Collins
Sample Position	P1-E1	Comment	
Injection Volume	10		
Acq. Date-Time	4/7/2021 11:39:35 AM		
Sample Info.			

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.468	797666	∞	10.7	303.85	2360364	24.4562 ng/ml
THC-COOH	1.489	981000	∞	56.8	∞	532288	73.1086 ng/ml
THC	3.330	752835	∞	27.2	∞	3423240	23.9934 ng/ml

SC




AM #27 Cannabinoid Quant. Results

Batch results D:\MassHunter\Data\2021\AM 27-28\040721 AM 27 28 SC\QuantResults\AM 27.batch.bin
Calibration Last Update 4/8/2021 9:07:24 AM

Instrument	Instrument 1	Data File	MJQ_Cal 6R.d
Type	Cal	Sample	MJQ_Cal 6R
Acq. Method	AM 27 THCQ.m	Operator	Sarah Collins
Sample Position	P1-F1	Comment	
Injection Volume	10		
Acq. Date-Time	4/7/2021 1:48:50 PM		

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.468	1661117	∞	11.1	5074.42	2499535	50.0225 ng/ml
THC-COOH	1.504	1395246	∞	57.3	4708.20	565997	97.8882 ng/ml
THC	3.345	1678588	455.29	25.1	242.79	3639645	49.8670 ng/ml

SC




AM #27 Cannabinoid Quant. Results

Batch results D:\MassHunter\Data\2021\AM 27-28\040721 AM 27 28 SC\QuantResults\AM 27.batch.bin
Calibration Last Update 4/8/2021 9:07:24 AM

Instrument	Instrument 1	Data File	MJQ_Cal 7.d
Type	Cal	Sample	MJQ_Cal 7
Acq. Method	AM 27 THCQ.m	Operator	Sarah Collins
Sample Position	P1-G1	Comment	
Injection Volume	10		
Acq. Date-Time	4/7/2021 11:54:46 AM		
Sample Info.			

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.468	3095918	∞	11.7	1966.21	2373207	100.1142 ng/ml
THC-COOH	1.489	3189628	8122.00	58.3	∞	498028	254.7974 ng/ml
THC	3.345	3231347	∞	25.1	∞	3445818	100.9718 ng/ml