


**REVIEWED**

By Tamara Salazar at 11:42 am, May 11, 2021

5/7/2021

**Worklist: 4956**

<u>LAB CASE</u>	<u>ITEM</u>	<u>ITEM TYPE</u>	<u>DESCRIPTION</u>
P2021-1379	10	URINE	AM 27 Urine Cannabinoids Confirmation by LC-QQQ
P2021-1379	4	BLOOD	AM 27 Blood THC Quant by LC-QQQ
P2021-1379	9	URINE	AM 27 Urine Cannabinoids Confirmation by LC-QQQ



## AM# 27: Quantitation of THC and Metabolites in Blood and Urine by LC-MS/MS

Extraction Date: 5/03/2021

Analyst: Amber Gerheart

Plate lot#: 210412

Plate Expiration: 10/12/2021

**Mobile phase A:** 0.1% Formic Acid in LCMS Water

**Mobile phase B:** 0.1% Formic acid in Acetonitrile

**Blank Blood Lot:** 20L20723

**Column:** UCT Selectra DA 100 x 2.1mm 3um

**LCMS-QQQ ID:** 069901

**Blank Urine Lot:** POC031319

### Pre-Analytic:

- 1. Check levels of mobile phases and needle wash refill as needed. Ensure waste is not full.
- 2. Ensure correct column is installed and begin mobile phase flow allow to equilibrate ~ 30 minutes.

### Analytic:

- 1. Remove standards, plate, controls, and samples from cold storage. Allow to reach room temperature.
- 2. Urine hydrolysis: add 1.5mL urine to blank plate, add 250µl 1N KOH. Shake and incubate at 40 degrees for 15 minutes. Using a calibrated pipette, add **1000µl blood and urine (if applicable) (calibrated pipette)** into the appropriate wells of analytical (standards) plate. **Pipette ID: 42**
- 3. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 4. Pipette **500µL 0.1% formic acid in water blood sample, 500 µL saturated phosphate buffer in urine** in wells of analytical plate.
- 5. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 6. Transfer **700-800µL of blood+acid or urine+acid** mixture to corresponding wells of SLE+ plate. Amount transferred: 800 µL
- 7. Apply positive pressure for approx. 10-15 seconds (or until no liquid remains on top of sorbent). **(Load at 85-100 PSI- Selector to the right)**
- 8. Wait 5 minutes.
- 9. Add **2.25mL MTBE. (Add in 3 increments of 750uL)**
- 10. Wait 5 minutes.
- 11. Apply positive pressure for approx. 15 seconds. **(10-15 PSI- Selector to the left).**
- 12. Add **2.25mL Hexane. (Add in 3 increments of 750uL)**
- 13. Wait 5 minutes.
- 14. Apply positive pressure for approx. 15 seconds. **(10-15 PSI- Selector to the left).**
- 15. Remove plate containing eluate. Place on SPE Dry and evaporate to dryness at approx. 35°C.
- 16. Reconstitute in **100µL 100% MeOH** and heat seal plate with foil. Place in autosampler and run worklist.

### Post-Analytic

- 1. Create batch and process data.
- 2. Make any necessary integration changes, Curve weighting of Linear 1/x with  $r^2$  values  $\geq 0.98$  for each analyte
- 3. RT +/- 3% or 0.100 min, whichever is greater, +/- 20% Accuracy for greater than (+/- 30% for 10ng/ml or less). Ion ratios must be within +/- 20% of the averaged calibrators
- 4. Case samples with calculated concentrations for THC at 1ng/mL or greater and OH-THC at 3ng/mL or greater may be reported quantitatively (blood only). Calculated concentrations for carboxy-THC of 5ng/mL may be reported qualitatively. Samples with a THC or OH-THC response over 50 ng/mL will be reported out as greater than 50 ng/mL.
- 5. Did all QCs pass for each analyte? (if not, describe in comments section)
- 6. Enter QCs into control charting.
- 7. Central File Packet to include: LIMS Worklist, Method Checklist, Calibration and Control Reports

COMMENTS: *Calibration Curve Range: THC 3-100 ng/mL, THC-OH 3-100 ng/mL*

	1	2	3	4	5	6
A	IS + Cal. 1	Negative Blood	IS + Sample	IS + Sample	IS + Sample	IS + QC_1
B	IS + Cal. 2	Negative Urine	IS + Sample	IS + Sample	IS + Sample	IS + Cal. 7
C	IS + Cal. 3	Urine External	IS + Sample	IS + Sample	IS + Sample	IS + Cal. 6
D	IS + Cal. 4	P2021-1379-4	IS + Sample	IS + Sample	IS + Sample	IS + Cal. 5
E	IS + Cal. 5	P2021-1379-9	IS + Sample	IS + Sample	IS + Sample	IS + Cal. 4
F	IS + Cal. 6	P2021-1379-10	IS + Sample	IS + Sample	IS + Sample	IS + Cal. 3
G	IS + Cal. 7	IS + Sample	IS + Sample	IS + Sample	IS + Sample	IS + Cal. 2
H	IS + QC_1	IS + Sample	IS + Sample	IS + Sample	IS + Sample	IS + Cal. 1

All wells to contain 100 µl of residual DMSO

AK



**Idaho State Police  
Forensic Services**

**AM #26 Screening of THC and Metabolites and AM #27  
Confirmation of THC and Metabolites Blood External  
Control Prep Sheet**

**Methanol External Control Solution (Lot: WS03052021)**

10 µL of 1mg/mL THC, 100 µL of 100 µg/mL THC-OH, C-THC in 9790 µL MeOH

*Approximate concentration 1ug/mL.*

<i>Component</i>	<i>Source</i>	<i>Source Lot Number</i>	<i>Expiration Date</i>
Methanol (LCMS)	Fisher	200921	
THC	Cerilliant	FE01041701	03/31/2022
C-THC	Cerilliant	FE08011801	08/31/2023
THC-OH	Cerilliant	FE07221601	07/31/2021
Prepared:	03/05/2021		
Prepared By:	Tamara Salazar/Amber Gerheart		

**Urine External Control Solution (Lot: 04232021)**

200 ul of methanol external control solution was added to 9800 ul of blood.

*Approximately 20ng/mL each*

<i>Component</i>	<i>Source</i>	<i>Source Lot Number</i>
Negative Blood	Lampire	20L20724
Methanol External Control Solution	-	WS03052021
Prepared:	04/23/2021	
Prepared by:	Sarah Collins	

AA




# AM #27 Cannabinoid Quant. Results

**Batch results** D:\MassHunter\Data\2021\AM 27-28\050321 AM 27 28 AG\QuantResults\AM 27.batch.bin  
**Calibration Last Update** 5/10/2021 10:24:13 AM

<b>Instrument</b>	Falco (069901)	<b>Data File</b>	MJQ_Negative Urine.d
<b>Type</b>	Sample	<b>Sample</b>	MJQ_Negative Urine
<b>Acq. Method</b>	AM 27 THCQ.m	<b>Operator</b>	Amber Gerheart
<b>Sample Position</b>	P1-B2	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	5/3/2021 3:01:11 PM		
<b>Sample Info.</b>			

## Sample Chromatogram



AG




# AM #27 Cannabinoid Quant. Results

**Batch results** D:\MassHunter\Data\2021\AM 27-28\050321 AM 27 28 AG\QuantResults\AM 27.batch.bin  
**Calibration Last Update** 5/10/2021 10:24:13 AM

<b>Instrument</b>	Falco (069901)	<b>Data File</b>	MJQ_Negative Blood.d
<b>Type</b>	Sample	<b>Sample</b>	MJQ_Negative Blood
<b>Acq. Method</b>	AM 27 THCQ.m	<b>Operator</b>	Amber Gerheart
<b>Sample Position</b>	P1-A2	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	5/3/2021 2:45:58 PM		
<b>Sample Info.</b>			

## Sample Chromatogram



AA




# AM #27 Cannabinoid Quant. Results

**Batch results** D:\MassHunter\Data\2021\AM 27-28\050321 AM 27 28 AG\QuantResults\AM 27.batch.bin  
**Calibration Last Update** 5/10/2021 10:24:13 AM

<b>Instrument</b>	Falco (069901)	<b>Data File</b>	MJQ_Urine Control.d
<b>Type</b>	Sample	<b>Sample</b>	MJQ_Urine Control
<b>Acq. Method</b>	AM 27 THCQ.m	<b>Operator</b>	Amber Gerheart
<b>Sample Position</b>	P1-C2	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	5/3/2021 3:16:24 PM		

**Sample Chromatogram**



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.468	316894	∞	12.2	∞	1978672	9.5234 ng/ml
THC-COOH	1.489	154450	∞	48.9	∞	459920	12.6682 ng/ml
THC	3.330	215931	243.44	31.7	424.55	1968432	10.7891 ng/ml

AG




# AM #27 Cannabinoid Quant. Results

**Batch results** D:\MassHunter\Data\2021\AM 27-28\050321 AM 27 28 AG\QuantResults\AM 27.batch.bin  
**Calibration Last Update** 5/10/2021 10:24:13 AM

<b>Instrument</b>	Falco (069901)	<b>Data File</b>	MJQ_QC Start.d
<b>Type</b>	Sample	<b>Sample</b>	MJQ_QC Start
<b>Acq. Method</b>	AM 27 THCQ.m	<b>Operator</b>	Amber Gerheart
<b>Sample Position</b>	P1-H1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	5/3/2021 2:30:45 PM		

**Sample Chromatogram**



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.468	189162	∞	10.8	90.92	2185879	4.9867 ng/ml
THC-COOH	1.504	248317	591.55	51.6	1239.43	615242	15.1972 ng/ml
THC	3.330	228211	795.49	31.2	216.57	5136118	4.5755 ng/ml



AG




# AM #27 Cannabinoid Quant. Results

**Batch results** D:\MassHunter\Data\2021\AM 27-28\050321 AM 27 28 AG\QuantResults\AM 27.batch.bin  
**Calibration Last Update** 5/10/2021 10:24:13 AM

<b>Instrument</b>	Falco (069901)	<b>Data File</b>	MJQ_QC End.d
<b>Type</b>	Sample	<b>Sample</b>	MJQ_QC End
<b>Acq. Method</b>	AM 27 THCQ.m	<b>Operator</b>	Amber Gerheart
<b>Sample Position</b>	P1-H1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	5/3/2021 4:09:40 PM		

**Sample Chromatogram**




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.468	174503	∞	11.7	∞	2197104	4.5483 ng/ml
THC-COOH	1.489	248585	∞	52.2	1537.41	622254	15.0436 ng/ml
THC	3.330	243496	505.04	30.3	47.64	5487703	4.5696 ng/ml

AK



# AM #27 Cannabinoids Quant. Calibration Curve Report

**Batch results** D:\MassHunter\Data\2021\AM 27-28\050321 AM 27 28 AG\QuantResults\AM 27.batch.bin  
**Last Cal. Update** 5/10/2021 10:24 AM  
**Analyst Name** ISP\datastor  
**Analyte** THC **Internal Standard** THC-D3



Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
MJQ_Cal 1	1	x	1.0	1.3	131.8
MJQ_Cal 2	2	✓	3.0	3.2	106.0
MJQ_Cal 3	3	✓	5.0	5.1	101.3
MJQ_Cal 4	4	✓	10.0	9.5	94.8
MJQ_Cal 5	5	✓	25.0	24.3	97.3
MJQ_Cal 6	6	✓	50.0	49.6	99.2
MJQ_Cal 7	7	✓	100.0	101.3	101.3


\*Cal 1 dropped due to S/N

AS



# AM #27 Cannabinoids Quant. Calibration Curve Report

**Batch results** D:\MassHunter\Data\2021\AM 27-28\050321 AM 27 28 AG\QuantResults\AM 27.batch.bin  
**Last Cal. Update** 5/10/2021 10:24 AM  
**Analyst Name** ISP\datastor  
**Analyte** THC-COOH **Internal Standard** THC-COOH-D9



Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
MJQ_Cal 1	1	✓	5.0	5.6	112.4
MJQ_Cal 2	2	✓	10.0	9.9	99.1
MJQ_Cal 3	3	✓	20.0	19.2	96.1
MJQ_Cal 4	4	✓	50.0	47.4	94.8
MJQ_Cal 5	5	✓	75.0	71.9	95.9
MJQ_Cal 6	6	✓	100.0	98.9	98.9
MJQ_Cal 7	7	✓	250.0	257.0	102.8


AS



# AM #27 Cannabinoids Quant. Calibration Curve Report

**Batch results** D:\MassHunter\Data\2021\AM 27-28\050321 AM 27 28 AG\QuantResults\AM 27.batch.bin  
**Last Cal. Update** 5/10/2021 10:24 AM  
**Analyst Name** ISP\datastor  
**Analyte** THC-OH **Internal Standard** THC-OH-D3

THC-OH - 7 Levels, 6 Levels Used, 7 Points, 6 Points Used, 0 QCs



Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
MJQ_Cal 1	1	x	1.0	3.4	341.4
MJQ_Cal 2	2	✓	3.0	3.1	102.5
MJQ_Cal 3	3	✓	5.0	5.0	100.2
MJQ_Cal 4	4	✓	10.0	9.7	96.5
MJQ_Cal 5	5	✓	25.0	25.0	100.1
MJQ_Cal 6	6	✓	50.0	50.5	101.0
MJQ_Cal 7	7	✓	100.0	99.7	99.7

\*Cal 1 dropped due to retention time

AA




# AM #27 Cannabinoid Quant. Results

**Batch results** D:\MassHunter\Data\2021\AM 27-28\050321 AM 27 28 AG\QuantResults\AM 27.batch.bin  
**Calibration Last Update** 5/10/2021 10:24:13 AM

<b>Instrument</b>	Falco (069901)	<b>Data File</b>	MJQ_Cal 1.d
<b>Type</b>	Cal	<b>Sample</b>	MJQ_Cal 1
<b>Acq. Method</b>	AM 27 THCQ.m	<b>Operator</b>	Amber Gerheart
<b>Sample Position</b>	P1-A1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	5/3/2021 1:29:49 PM		
<b>Sample Info.</b>			

**Sample Chromatogram**



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.513 <b>High</b>	129894	∞	5.4 <b>Low</b>	32.95	2128914	3.4138 ng/ml
THC-COOH	1.474	88404	∞	46.6	∞	601819	5.6201 ng/ml
THC	3.315	57205	165.18	28.2	5.13 <b>Low</b>	5596434	1.3184 ng/ml

AA




# AM #27 Cannabinoid Quant. Results

**Batch results** D:\MassHunter\Data\2021\AM 27-28\050321 AM 27 28 AG\QuantResults\AM 27.batch.bin  
**Calibration Last Update** 5/10/2021 10:24:13 AM

<b>Instrument</b>	Falco (069901)	<b>Data File</b>	MJQ_Cal 2.d
<b>Type</b>	Cal	<b>Sample</b>	MJQ_Cal 2
<b>Acq. Method</b>	AM 27 THCQ.m	<b>Operator</b>	Amber Gerheart
<b>Sample Position</b>	P1-B1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	5/3/2021 1:37:35 PM		

**Sample Chromatogram**



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.468	96401	∞	9.4	∞	1736547	3.0747 ng/ml
THC-COOH	1.459	119745	∞	46.7	1031.64	457410	9.9064 ng/ml
THC	3.435	96939	305.41	31.2	∞	3254027	3.1815 ng/ml

AG




# AM #27 Cannabinoid Quant. Results

**Batch results** D:\MassHunter\Data\2021\AM 27-28\050321 AM 27 28 AG\QuantResults\AM 27.batch.bin  
**Calibration Last Update** 5/10/2021 10:24:13 AM

<b>Instrument</b>	Falco (069901)	<b>Data File</b>	MJQ_Cal 3.d
<b>Type</b>	Cal	<b>Sample</b>	MJQ_Cal 3
<b>Acq. Method</b>	AM 27 THCQ.m	<b>Operator</b>	Amber Gerheart
<b>Sample Position</b>	P1-C1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	5/3/2021 1:45:11 PM		

**Sample Chromatogram**



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.453	158359	∞	9.7	26.89	1822739	5.0078 ng/ml
THC-COOH	1.459	237515	∞	49.4	970.89	464435	19.2188 ng/ml
THC	3.465	146414	∞	25.2	12.42	2954471	5.0634 ng/ml

AA




# AM #27 Cannabinoid Quant. Results

**Batch results** D:\MassHunter\Data\2021\AM 27-28\050321 AM 27 28 AG\QuantResults\AM 27.batch.bin  
**Calibration Last Update** 5/10/2021 10:24:13 AM

<b>Instrument</b>	Falco (069901)	<b>Data File</b>	MJQ_Cal 4.d
<b>Type</b>	Cal	<b>Sample</b>	MJQ_Cal 4
<b>Acq. Method</b>	AM 27 THCQ.m	<b>Operator</b>	Amber Gerheart
<b>Sample Position</b>	P1-D1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	5/3/2021 1:52:47 PM		

**Sample Chromatogram**



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.468	361264	∞	10.7	584.55	2227109	9.6502 ng/ml
THC-COOH	1.474	698685	∞	52.9	∞	551371	47.4139 ng/ml
THC	3.405	388552	499.83	26.7	133.32	4049130	9.4812 ng/ml



AA




# AM #27 Cannabinoid Quant. Results

**Batch results** D:\MassHunter\Data\2021\AM 27-28\050321 AM 27 28 AG\QuantResults\AM 27.batch.bin  
**Calibration Last Update** 5/10/2021 10:24:13 AM

<b>Instrument</b>	Falco (069901)	<b>Data File</b>	MJQ_Cal 5.d
<b>Type</b>	Cal	<b>Sample</b>	MJQ_Cal 5
<b>Acq. Method</b>	AM 27 THCQ.m	<b>Operator</b>	Amber Gerheart
<b>Sample Position</b>	P1-E1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	5/3/2021 2:00:23 PM		

**Sample Chromatogram**



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.468	839192	∞	10.4	695.02	2038092	25.0287 ng/ml
THC-COOH	1.474	986137	∞	52.5	∞	512576	71.9132 ng/ml
THC	3.405	950072	554.02	26.1	260.90	3772094	24.3248 ng/ml

AG




# AM #27 Cannabinoid Quant. Results

**Batch results** D:\MassHunter\Data\2021\AM 27-28\050321 AM 27 28 AG\QuantResults\AM 27.batch.bin  
**Calibration Last Update** 5/10/2021 10:24:13 AM

<b>Instrument</b>	Falco (069901)	<b>Data File</b>	MJQ_Cal 6.d
<b>Type</b>	Cal	<b>Sample</b>	MJQ_Cal 6
<b>Acq. Method</b>	AM 27 THCQ.m	<b>Operator</b>	Amber Gerheart
<b>Sample Position</b>	P1-F1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	5/3/2021 2:07:58 PM		

**Sample Chromatogram**



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.468	1790164	∞	11.1	∞	2169788	50.4981 ng/ml
THC-COOH	1.489	1534946	∞	54.4	∞	579660	98.9279 ng/ml
THC	3.330	2716277	∞	26.7	∞	5247856	49.6240 ng/ml

AA




# AM #27 Cannabinoid Quant. Results

**Batch results** D:\MassHunter\Data\2021\AM 27-28\050321 AM 27 28 AG\QuantResults\AM 27.batch.bin  
**Calibration Last Update** 5/10/2021 10:24:13 AM

<b>Instrument</b>	Falco (069901)	<b>Data File</b>	MJQ_Cal 7.d
<b>Type</b>	Cal	<b>Sample</b>	MJQ_Cal 7
<b>Acq. Method</b>	AM 27 THCQ.m	<b>Operator</b>	Amber Gerheart
<b>Sample Position</b>	P1-G1	<b>Comment</b>	
<b>Injection Volume</b>	10		
<b>Acq. Date-Time</b>	5/3/2021 2:15:34 PM		

**Sample Chromatogram**



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.468	3328062	∞	11.3	∞	2049192	99.7405 ng/ml
THC-COOH	1.489	3738878	∞	53.3	3076.40	543037	256.9997 ng/ml
THC	3.330	5775716	∞	26.6	3602.78	5445506	101.3252 ng/ml