








REVIEWED

By Sarah Collins at 1:40 pm, Nov 18, 2021



11/17/2021

Worklist: 5391

<u>LAB CASE</u>	<u>ITEM</u>	<u>ITEM TYPE</u>	<u>DESCRIPTION</u>	
C2021-2399		BCK	AM 27 Blood THC Quant by LC-QQQ	
C2021-2434		UCK	AM 27 Urine Cannabinoids Confirmation by LC-QQQ	
C2021-2435	4	UCK	AM 27 Urine Cannabinoids Confirmation by LC-QQQ	
C2021-2443		UCK	AM 27 Urine Cannabinoids Confirmation by LC-QQQ	
C2021-2459		BCK	AM 27 Blood THC Quant by LC-QQQ	
C2021-2461		BCK	AM 27 Blood THC Quant by LC-QQQ	



AM# 27: Quantitation of THC and Metabolites in Blood and Urine by LC-MS/MS

Extraction Date 11/17/21
Plate lot#: 210609

Analyst: Anne Nord
Plate Expiration: 12-9-21

Mobile phase A: 0.1% Formic Acid in LCMS Water MTBE
Mobile phase B: 0.1% Formic acid in Acetonitrile Hexane
LCMS Methanol

Blank Blood Lot: 21D52496 **Urine Blank:** 83121 **Column:** UCT Selectra DA 100 x 2.1mm 3um
LCMS-QQQ ID: 69679

Pre-Analytic:

- 1. Check levels of mobile phases and needle wash refill as needed. Ensure waste is not full.
- 2. Ensure correct column is installed and begin mobile phase flow allow to equilibrate ~ 30 minutes.


Analytic:

- 1. Remove standards, plate, controls, and samples from cold storage. Allow to reach room temperature.
- 2. Urine hydrolysis: add 1.5 ml urine to blank plate, add 250 ul 1N KOH mix and incubate at 40 degrees for 15 minutes.
Pipette 1000µL blood (calibrated pipette) Pipette ID: k52558g in wells of analytical (standards) plate.
- 3. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 4. Pipette 500µL 0.1% formic acid in water blood sample, 500 ul saturated phosphate buffer in urine in wells of analytical plate.
- 5. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 6. Transfer 800µL of blood+acid or urine acid mixture to corresponding wells of SLE+ plate.
- 7. Apply positive pressure for approx. 10-15 seconds (or until no liquid remains on top of sorbent).
(Load at 85-100 PSI- Selector to the right) Manifold ID: 66792
- 8. Wait 5 minutes.
- 9. Add 2.25mL MTBE. *(Add in 3 increments of 750uL)*
- 10. Wait 5 minutes.
- 11. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 12. Add 2.25mL Hexane. *(Add in 3 increments of 750uL)*
- 13. Wait 5 minutes.
- 14. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left).*
- 15. Remove plate containing eluate. Place on SPE Dry and evaporate to dryness at approx. 35°C.
SPE Dry ID: 66819
- 16. Reconstitute in 100µL 100% MeOH and heat seal plate with foil. Place in autosampler and run worklist.

Post-Analytic

- 1. Create batch and process data.
- 2. Make any necessary integration changes, Curve weighting of Linear 1/x with r^2 values ≥ 0.98 for each analyte
- 3. RT +/- 3% or 0.100 min, whichever is greater, +/- 20% Accuracy for greater than (+/- 30% for 10ng/ml or less).
Ion ratios must be within +/- 20% of the averaged calibrators
- 4. Case sample response for THC 1ng/ml, OH-THC 3ng/mL (quantitative blood), Carboxy-THC: 5 ng/mL (qualitative only) will be reported. Samples with a THC or OH-THC response over 50 ng/mL will be reported out as greater than 50 ng/mL.
- 5. Did all QCs pass for each analyte? (if not is it describe in comments section)
- 6. Enter QCs into control charting.
- 7. Central File Packet to include: LIMS Worklist, Method Checklist, Calibration and Control Reports


COMMENTS: *THC-OH 3-100*



	1	2	3	4	5	6
a	cal 1	neg blood	2443-1			
b	cal 2	2399-1				
c	cal 3	2459-1				
d	cal 4	2461-1				
e	Cal 5	neg urine				
f	cal 6	Urine positive				
g	cal 7	2434-1				
h	Internal control	2435-4				

C2021-____-__

Toxicology AM method 27/26 external prep informati



working solution 15 ug/ml in meoh C-THC, THC-OH, 7.5 ug/ml THC

Stock solution 1mg/ml 7.5 ul each THC, 100 ug/ml 150 ul C-THC, 150 ul THC-OH in 9692.5 ul meOH

Ppd 8/26/21 Exp: 8/26/22 lot 82621 by AMN

Drug	lot	expiration
C-THC	FE04151901	6/1/2024
THC-OH	FE06152002	6/1/2025
THC	FE04222001	5/1/2025

AM 27/26 blood control 100 ul working solution lot () in 9900 ul blood lot ()

		Concentration 7.5 ng/ml THC, 15 ng/ml C-THC, THC-OH	
--	--	--------------------------------------------------------	--

AM 27/26 urine control 400 ul working solution lot (82620) in 9600 ul urine

out of use


ppd 8/26/21 Exp 8/26/22 neg urine lot 5621	lot u82621	Concentration 30 ng/ml THC, and 60 ng/ml C-THC, THC-OH	by amn	last used 11/1/21
ppd 11/2/21 Exp 8/26/22 neg urine lot 83121	lot u11221	Concentration 30 ng/ml THC, and 60 ng/ml C-THC, THC-OH	by amn	

AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2021\am 27-28\111721\QuantResults\cann.batch.bin
Calibration Last Update 11/18/2021 8:04:14 AM

Instrument	69679	Data File	mj internal control blood.d
Type	QC	Sample	mj internal control blood
Acq. Method	AM 27 THC quant.m	Operator	Anne Nord
Sample Position	P3-H1	Comment	
Injection Volume	10		
Acq. Date-Time	11/17/2021 2:53:25 PM		
Sample Info.			

Sample Chromatogram




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.394	36920	∞	1072.3	2638.1	2259917	5.247 ng/ml
THC-COOH	1.416	140872	166.3	34.9	307.8	737766	14.970 ng/ml
THC	3.137	326369	∞	23.5	1376.7	2655260	4.666 ng/ml

AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2021\am 27-28\111721\QuantResults\cann.batch.bin
Calibration Last Update 11/18/2021 8:04:14 AM

Instrument	69679	Data File	mj negative blood.d
Type	Sample	Sample	mj negative blood
Acq. Method	AM 27 THC quant.m	Operator	Anne Nord
Sample Position	P3-A2	Comment	
Injection Volume	10		
Acq. Date-Time	11/17/2021 3:00:08 PM		
Sample Info.			

Sample Chromatogram




AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2021\am 27-28\111721\QuantResults\cann.batch.bin
Calibration Last Update 11/18/2021 8:04:14 AM

Instrument	69679	Data File	mj negative urine.d
Type	Sample	Sample	mj negative urine
Acq. Method	AM 27 THC quant.m	Operator	Anne Nord
Sample Position	P3-E2	Comment	
Injection Volume	10		
Acq. Date-Time	11/17/2021 3:53:26 PM		
Sample Info.			

Sample Chromatogram




AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2021\am 27-28\111721\QuantResults\cann.batch.bin
Calibration Last Update 11/18/2021 8:04:14 AM

Instrument	69679	Data File	mj urine positive control.d
Type	Sample	Sample	mj urine positive control
Acq. Method	AM 27 THC quant.m	Operator	Anne Nord
Sample Position	P3-F2	Comment	
Injection Volume	10		
Acq. Date-Time	11/17/2021 4:00:10 PM		
Sample Info.			

Sample Chromatogram




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.394	408807	∞	832.2	∞	2584307	47.008 ng/ml
THC-COOH	1.416	251500	1488.3	36.3	655.2	517931	37.063 ng/ml
THC	3.122	153371	∞	25.0	484.2	509880	10.861 ng/ml

AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2021\am 27-28\111721\QuantResults\cann.batch.bin
Calibration Last Update 11/18/2021 8:04:14 AM

Instrument	69679	Data File	mj internal control blood end of run.d
Type	QC	Sample	mj internal control blood end of run
Acq. Method	AM 27 THC quant.m	Operator	Anne Nord
Sample Position	P3-H1	Comment	
Injection Volume	10		
Acq. Date-Time	11/17/2021 4:53:30 PM		
Sample Info.			

Sample Chromatogram




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.394	38315	∞	1021.1	∞	2386258	5.164 ng/ml
THC-COOH	1.416	154475	2067.3	32.8	244.2	779000	15.521 ng/ml
THC	3.137	380357	99855.1	23.0	1122.5	3029184	4.759 ng/ml

Compound Calibration Report


Batch results D:\MassHunter\Data\2021\am 27-28\111721\QuantResults\cann.batch.bin
Last Cal. Update 11/18/2021 8:04 AM
Analyst Name ISP\datastor
Analyte THC **Internal Standard** THC-d3

THC - 7 Levels, 7 Levels Used, 7 Points, 7 Points Used, 2 QCs




Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
mj cal 1	1	✓	1.0	1.2	116.6
mj cal2	2	✓	3.0	2.9	95.4
mj cal 3	3	✓	5.0	4.7	93.9
mj cal 4	4	✓	10.0	9.4	94.5
mj cal 5	5	✓	25.0	24.6	98.5
mj cal 6	6	✓	50.0	49.9	99.8
mj cal 7	7	✓	100.0	101.3	101.3

Compound Calibration Report



Batch results D:\MassHunter\Data\2021\am 27-28\111721\QuantResults\cann.batch.bin
Last Cal. Update 11/18/2021 8:04 AM
Analyst Name ISP\datastor
Analyte THC-COOH **Internal Standard** THC-COOH-d9

THC-COOH - 7 Levels, 7 Levels Used, 7 Points, 7 Points Used, 2 QCs




Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
mj cal 1	1	✓	5.0	5.2	103.6
mj cal2	2	✓	10.0	10.1	100.6
mj cal 3	3	✓	20.0	19.5	97.7
mj cal 4	4	✓	50.0	48.4	96.8
mj cal 5	5	✓	75.0	76.0	101.4
mj cal 6	6	✓	100.0	99.5	99.5
mj cal 7	7	✓	250.0	251.3	100.5

Compound Calibration Report



Batch results D:\MassHunter\Data\2021\am 27-28\111721\QuantResults\cann.batch.bin
Last Cal. Update 11/18/2021 8:04 AM
Analyst Name ISP\datastor
Analyte THC-OH **Internal Standard** THC-OH-d3

THC-OH - 7 Levels, 6 Levels Used, 7 Points, 6 Points Used, 2 QCs



Sample	Level	Enabled	Expected Concentration	Final Concentration	Accuracy
mj cal 1	1	x	1.0	1.5	153.4
mj cal2	2	✓	3.0	3.3	109.2
mj cal 3	3	✓	5.0	5.0	100.0
mj cal 4	4	✓	10.0	9.4	94.0
mj cal 5	5	✓	25.0	24.0	96.1
mj cal 6	6	✓	50.0	49.3	98.7
mj cal 7	7	✓	100.0	102.0	102.0

Cal 1 dropped ratio out of range


AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2021\am 27-28\111721\QuantResults\cann.batch.bin
Calibration Last Update 11/18/2021 8:04:14 AM

Instrument	69679	Data File	mj cal 1.d
Type	Cal	Sample	mj cal 1
Acq. Method	AM 27 THC quant.m	Operator	Anne Nord
Sample Position	P3-A1	Comment	
Injection Volume	10		
Acq. Date-Time	11/17/2021 1:59:53 PM		

Sample Info.

Sample Chromatogram




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.	
THC-OH	1.394	9830	∞	1273.6 High	∞	2637463	1.534 ng/ml	Low
THC-COOH	1.416	45363	24.7	35.7	37230.9	751427	5.179 ng/ml	
THC	3.137	70055	591.0	25.7	∞	3125405	1.166 ng/ml	

AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2021\am 27-28\111721\QuantResults\cann.batch.bin
Calibration Last Update 11/18/2021 8:04:14 AM

Instrument	69679	Data File	mj cal2.d
Type	Cal	Sample	mj cal2
Acq. Method	AM 27 THC quant.m	Operator	Anne Nord
Sample Position	P3-B1	Comment	
Injection Volume	10		
Acq. Date-Time	11/17/2021 2:06:37 PM		
Sample Info.			

Sample Chromatogram




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.394	24432	∞	1100.3	∞	2532368	3.277 ng/ml
THC-COOH	1.416	90491	138.8	34.9	871.7	721407	10.058 ng/ml
THC	3.137	213236	4459.2	25.2	363.2	2999944	2.861 ng/ml

AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2021\am 27-28\111721\QuantResults\cann.batch.bin
Calibration Last Update 11/18/2021 8:04:14 AM

Instrument	69679	Data File	mj cal 3.d
Type	Cal	Sample	mj cal 3
Acq. Method	AM 27 THC quant.m	Operator	Anne Nord
Sample Position	P3-C1	Comment	
Injection Volume	10		
Acq. Date-Time	11/17/2021 2:13:19 PM		
Sample Info.			

Sample Chromatogram




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.394	40055	7753.4	1043.6	17558 57127 31558. 0	2583393	5.002 ng/ml
THC-COOH	1.416	183743	513.7	34.8	366.0	729610	19.536 ng/ml
THC	3.137	370559	∞	23.4	4855.7	2994172	4.696 ng/ml

AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2021\am 27-28\111721\QuantResults\cann.batch.bin
Calibration Last Update 11/18/2021 8:04:14 AM

Instrument	69679	Data File	mj cal 4.d
Type	Cal	Sample	mj cal 4
Acq. Method	AM 27 THC quant.m	Operator	Anne Nord
Sample Position	P3-D1	Comment	
Injection Volume	10		
Acq. Date-Time	11/17/2021 2:20:01 PM		
Sample Info.			

Sample Chromatogram




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.394	80081	∞	992.2	1907.7	2631117	9.397 ng/ml
THC-COOH	1.416	466366	690.0	36.5	756.2	732749	48.376 ng/ml
THC	3.137	819769	4586237653 743590.0	23.3	6610.7	3149873	9.449 ng/ml

AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2021\am 27-28\111721\QuantResults\cann.batch.bin
Calibration Last Update 11/18/2021 8:04:14 AM

Instrument	69679	Data File	mj cal 5.d
Type	Cal	Sample	mj cal 5
Acq. Method	AM 27 THC quant.m	Operator	Anne Nord
Sample Position	P3-E1	Comment	
Injection Volume	10		
Acq. Date-Time	11/17/2021 2:26:43 PM		
Sample Info.			

Sample Chromatogram




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.394	213763	∞	954.9	∞	2667381	24.030 ng/ml
THC-COOH	1.416	723889	3573.7	36.2	6804.0	720195	76.020 ng/ml
THC	3.137	2157370	∞	23.7	16160. 8	3100376	24.619 ng/ml

AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2021\am 27-28\111721\QuantResults\cann.batch.bin
Calibration Last Update 11/18/2021 8:04:14 AM

Instrument	69679	Data File	mj cal 6.d
Type	Cal	Sample	mj cal 6
Acq. Method	AM 27 THC quant.m	Operator	Anne Nord
Sample Position	P3-F1	Comment	
Injection Volume	10		
Acq. Date-Time	11/17/2021 2:33:25 PM		
Sample Info.			

Sample Chromatogram




Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.394	382017	∞	905.8	∞	2299665	49.343 ng/ml
THC-COOH	1.416	945433	106812.4	36.5	3408.6	716947	99.531 ng/ml
THC	3.137	3953840	∞	24.3	46877 73445 41163. 0	2780092	49.915 ng/ml

AM #27 Cannabinoids

Batch results D:\MassHunter\Data\2021\am 27-28\111721\QuantResults\cann.batch.bin
Calibration Last Update 11/18/2021 8:04:14 AM

Instrument	69679	Data File	mj cal 7.d
Type	Cal	Sample	mj cal 7
Acq. Method	AM 27 THC quant.m	Operator	Anne Nord
Sample Position	P3-G1	Comment	
Injection Volume	10		
Acq. Date-Time	11/17/2021 2:40:07 PM		
Sample Info.			

Sample Chromatogram



Name	RT	Resp.	S/N	Ratio	S/N	ISTD Resp.	Final Conc.
THC-OH	1.394	845131	∞	860.3	26301 8.3	2451005	101.951 ng/ml
THC-COOH	1.416	2148348	2901.3	37.6	2112.0	642689	251.301 ng/ml
THC	3.137	7995943	∞	24.7	11501 33690 41732 0.0	2759583	101.294 ng/ml