


REVIEWED
By Sarah Collins at 2:33 pm, Apr 27, 2021

TS

4/15/2021

Worklist: 4906

| <u>LAB CASE</u> | <u>ITEM</u> | <u>ITEM TYPE</u> | <u>DESCRIPTION</u> |
|-----------------|-------------|------------------|--|
| M2021-1319 | 2 | BCK | AM 28 Blood Multi-Drug Quant Panel 2 by LC-QQQ |
| P2021-0238 | 1 | BCK | AM 28 Blood Multi-Drug Quant Panel 2 by LC-QQQ |
| P2021-0968 | 1 | BCK | AM 28 Blood Multi-Drug Quant Panel 2 by LC-QQQ |
| P2021-1048 | 1 | BCK | AM 28 Blood Multi-Drug Quant Panel 2 by LC-QQQ |



AM# 28: Multi-Drug Quantitation by LC-MS/MS

TS

Extraction Date: 4/15/2021

Analyst: Tamara Salazar – HOA Amber Gerheart

Plate lot#: 201207

Plate Expiration: 6/7/2021

Mobile phase A: 5mM Amm Form + 0.01% FA **Mobile phase B:** 0.01% Formic Acid in MeOH

Blank Blood Lot: 20L20724

Blank Urine Lot: N/A

Column: Agilent 120 EC-C18 (2.1x 100-2.7um) **LCMS-QQQ ID:** 069901

Pre-Analytic:

- 1. Check levels of mobile phases and needle wash refill as needed. Ensure waste is not full.
- 2. Ensure correct column is installed and begin mobile phase flow allow to equilibrate ~ 30 minutes.
- 3. Create worklist

Analytic:

- 1. Remove standards, plate, controls, and samples from cold storage. Allow to reach room temperature.
- 2. Urine hydrolysis pipette 250 ul urine in blank well, add 40 ul BG Turbo, add 100 ul 500 mM sodium phosphate buffer, mix for at least 5 minutes at ambient temperature. Pipette **250µL blood (calibrated pipette) or 250µL hydrolyzed urine** in wells of analytical (standards) plate. **Pipette ID: 42**
- 3. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 4. Pipette **250µL 0.5M ammonium hydroxide** in wells of analytical plate.
- 5. Place on shaking incubator at ambient temp., 900rpm for 15 minutes.
- 6. Transfer **300µL of blood+base/urine+base** mixture to corresponding wells of SLE+ plate.
- 7. Apply positive pressure for approx. 10-15 seconds (or until no liquid remains on top of sorbent). *(Load at 85-100 PSI- Selector to the right)*
- 8. Wait 5 minutes.
- 9. Add **900uL ethyl acetate**.
- 10. Wait 5 minutes.
- 11. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left)*.
- 12. Add **900uL ethyl acetate**.
- 13. Wait 5 minutes.
- 14. Apply positive pressure for approx. 15 seconds. *(10-15 PSI- Selector to the left)*.
- 15. Remove plate containing eluate. **If run contains urine, add 50µL 1% HCl in MeOH to wells and place plate cover on plate before drying (optional).** Place on SPE Dry and evaporate to dryness at approx. 35°C.
- 16. Reconstitute in **100µL 20% MeOH** and heat seal plate with foil.

Post-Analytic

- 1. Create batch and process data.
- 2. Make necessary changes to integration limits
- 3. Integration linear and R² values ≥0.98 for each analyte.
- 4. For unknown samples and controls: response ratio within 20% of average of controls and standards, RT within +/- 5% (tramadol RT +/-2%), S/N for primary transition >10 and secondary transitions >5.
- 5. Did all QCs pass for each analyte? Yes, see comments Add Control data to QC tracking spreadsheet.
- 6. Central File Packet to include: LIMS Worklist, Method Checklist, Calibration and Control Reports.

COMMENTS: Compounds evaluated: Cocaethylene, Doxepin, Methylphenidate, Midazolam, Norketamine
Curve ranges: Norketamine 5-500 ng/mL

| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 |
|---|-------------|----------------|--------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|
| A | IS + Cal. 1 | IS + QC_1 | P2021-0968-1 | IS + Sample | IS + Cal. 1 | IS + QC_1 | IS + Sample | IS + Sample | IS + Sample | IS + Sample | IS + Sample | IS + Cal. 8 |
| B | IS + Cal. 2 | IS + QC_2 | P2021-1048-1 | IS + Sample | IS + Cal. 2 | IS + QC_2 | IS + Sample | IS + Sample | IS + Sample | IS + Sample | IS + Sample | IS + Cal. 7 |
| C | IS + Cal. 3 | IS + QC_3 | IS + Sample | IS + Sample | IS + Cal. 3 | IS + QC_3 | IS + Sample | IS + Sample | IS + Sample | IS + Sample | IS + Sample | IS + Cal. 6 |
| D | IS + Cal. 4 | IS + QC_4 | IS + Sample | IS + Sample | IS + Cal. 4 | IS + QC_4 | IS + Sample | IS + Sample | IS + Sample | IS + Sample | IS + QC_2 | IS + Cal. 5 |
| E | IS + Cal. 5 | IS + QC_2 | IS + Sample | IS + Sample | IS + Cal. 5 | IS + QC_2 | IS + Sample | IS + Sample | IS + Sample | IS + Sample | IS + QC_4 | IS + Cal. 4 |
| F | IS + Cal. 6 | Negative Blood | IS + Sample | IS + Sample | IS + Cal. 6 | IS + Sample | IS + Sample | IS + Sample | IS + Sample | IS + Sample | IS + QC_3 | IS + Cal. 3 |
| G | IS + Cal. 7 | M2021-1319-2 | IS + Sample | IS + Sample | IS + Cal. 7 | IS + Sample | IS + Sample | IS + Sample | IS + Sample | IS + Sample | IS + QC_2 | IS + Cal. 2 |
| H | IS + Cal. 8 | P2021-0238-1 | IS + Sample | IS + Sample | IS + Cal. 8 | IS + Sample | IS + Sample | IS + Sample | IS + Sample | IS + Sample | IS + QC_1 | IS + Cal. 1 |

All wells to contain 60 µl of Trapping Solution

| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 |
|---|-------------|-------------|-------------|-------------|-------------|----------------|--------------|-------------|-------------|-------------|-------------|-------------|
| A | IS + Cal. 1 | IS + QC_1 | IS + Sample | IS + Sample | IS + Cal. 1 | IS + QC_1 | P2021-0968-1 | IS + Sample | IS + Sample | IS + Sample | IS + Sample | IS + Cal. 8 |
| B | IS + Cal. 2 | IS + QC_2 | IS + Sample | IS + Sample | IS + Cal. 2 | IS + QC_2 | P2021-1048-1 | IS + Sample | IS + Sample | IS + Sample | IS + Sample | IS + Cal. 7 |
| C | IS + Cal. 3 | IS + QC_3 | IS + Sample | IS + Sample | IS + Cal. 3 | IS + QC_3 | IS + Sample | IS + Sample | IS + Sample | IS + Sample | IS + Sample | IS + Cal. 6 |
| D | IS + Cal. 4 | IS + QC_4 | IS + Sample | IS + Sample | IS + Cal. 4 | IS + QC_4 | IS + Sample | IS + Sample | IS + Sample | IS + Sample | IS + QC_2 | IS + Cal. 5 |
| E | IS + Cal. 5 | IS + QC_2 | IS + Sample | IS + Sample | IS + Cal. 5 | IS + QC_2 | IS + Sample | IS + Sample | IS + Sample | IS + Sample | IS + QC_4 | IS + Cal. 4 |
| F | IS + Cal. 6 | IS + Sample | IS + Sample | IS + Sample | IS + Cal. 6 | Negative Blood | IS + Sample | IS + Sample | IS + Sample | IS + Sample | IS + QC_3 | IS + Cal. 3 |
| G | IS + Cal. 7 | IS + Sample | IS + Sample | IS + Sample | IS + Cal. 7 | M2021-1319-2 | IS + Sample | IS + Sample | IS + Sample | IS + Sample | IS + QC_2 | IS + Cal. 2 |
| H | IS + Cal. 8 | IS + Sample | IS + Sample | IS + Sample | IS + Cal. 8 | P2021-0238-1 | IS + Sample | IS + Sample | IS + Sample | IS + Sample | IS + QC_1 | IS + Cal. 1 |

All wells to contain 60 µl of Trapping Solution

TS




AM #28 Multi-Drug Quant. Results

Batch results D:\MassHunter\Data\2021\AM 27-28\04152021 AM 28 AG\QuantResults\Evaluated AM 28 P2 AG.batch.bin
Calibration Last Update 4/20/2021 3:29:03 PM

| | | | |
|-------------------------|------------------------|------------------|---------------------|
| Instrument Type | Instrument 1 Sample | Data File | p2 Negative Blood.d |
| Acq. Method | AM 28 MDQ P2.m | Sample | p2 Negative Blood |
| Sample Position | P6-F6 | Operator | Tamara Salazar |
| Injection Volume | 2 | Comment | |
| Acq. Date-Time | 4/15/2021 4:27:33 PM | | |
| Sample Info. | | | |

Sample Chromatogram



TS




AM #28 Multi-Drug Quant. Results

Batch results D:\MassHunter\Data\2021\AM 27-28\04152021 AM 28 AG\QuantResults\Evaluated AM 28 P2 AG.batch.bin
Calibration Last Update 4/20/2021 3:29:03 PM

| | | | |
|-------------------------|----------------------|------------------|----------------|
| Instrument | Instrument 1 | Data File | p2 QC 10.d |
| Type | QC | Sample | p2 QC 10 |
| Acq. Method | AM 28 MDQ P2.m | Operator | Tamara Salazar |
| Sample Position | P6-A6 | Comment | |
| Injection Volume | 2 | | |
| Acq. Date-Time | 4/15/2021 3:24:02 PM | | |
| Sample Info. | | | |

Sample Chromatogram



| Name | RT | Resp. | S/N | Ratio | S/N | ISTD Resp. | Final Conc. |
|-----------------|-------|--------|---------|-------|----------|------------|---------------|
| Cocaethylene | 4.714 | 125080 | 5361.78 | 54.8 | 2229.22 | 241000 | 9.9061 ng/ml |
| Doxepin | 5.356 | 6774 | 106.59 | 49.8 | 12.08 | 29652 | 10.5833 ng/ml |
| Methylphenidate | 4.310 | 216134 | 7561.16 | 21.8 | 46245.88 | 817738 | 9.8642 ng/ml |
| Midazolam | 5.764 | 16641 | 2696.48 | 100.8 | 933.28 | 241000 | 9.9429 ng/ml |
| Norketamine | 4.075 | 30972 | 162.21 | 441.6 | 3151.85 | 669558 | 12.0861 ng/ml |

TS




AM #28 Multi-Drug Quant. Results

Batch results D:\MassHunter\Data\2021\AM 27-28\04152021 AM 28 AG\QuantResults\Evaluated AM 28 P2 AG.batch.bin
Calibration Last Update 4/20/2021 3:29:03 PM

| | | | |
|-------------------------|----------------------|------------------|-----------------------------|
| Instrument | Instrument 1 | Data File | p2 QC 250.d |
| Type | QC | Sample | p2 QC 250 QC 100 |
| Acq. Method | AM 28 MDQ P2.m | Operator | Tamara Salazar |
| Sample Position | P6-B6 | Comment | 4/22/2021 TS |
| Injection Volume | 2 | | |
| Acq. Date-Time | 4/15/2021 3:45:13 PM | | |
| Sample Info. | | | |

Sample Chromatogram



| Name | RT | Resp. | S/N | Ratio | S/N | ISTD Resp. | Final Conc. |
|-----------------|-------|---------|----------|-------|----------|------------|----------------|
| Cocaethylene | 4.714 | 1293761 | 595.50 | 54.6 | 32185.77 | 251812 | 92.1237 ng/ml |
| Doxepin | 5.356 | 68142 | 169.51 | 44.1 | 226.15 | 29523 | 102.8033 ng/ml |
| Methylphenidate | 4.310 | 2255236 | 48511.46 | 21.8 | 22521.43 | 853385 | 97.8896 ng/ml |
| Midazolam | 5.764 | 175502 | 1072.75 | 89.4 | 2477.41 | 251812 | 98.6245 ng/ml |
| Norketamine | 4.075 | 306499 | 2109.01 | 429.1 | 14148.41 | 700492 | 122.2698 ng/ml |

AM #28 Multi-Drug Quant. Results




TS

Batch results D:\MassHunter\Data\2021\AM 27-28\04152021 AM 28 AG\QuantResults\Evaluated AM 28 P2 AG.batch.bin
Calibration Last Update 4/20/2021 3:29:03 PM

| | | | |
|-------------------------|----------------------|------------------|--------------------------------|
| Instrument | Instrument 1 | Data File | p2 QC 100.d |
| Type | QC | Sample | p2 QC 100 QC 250 TS |
| Acq. Method | AM 28 MDQ P2.m | Operator | Tamara Salazar 4/22/2021 |
| Sample Position | P6-C6 | Comment | |
| Injection Volume | 2 | | |
| Acq. Date-Time | 4/15/2021 6:13:27 PM | | |
| Sample Info. | | | |

Sample Chromatogram



| Name | RT | Resp. | S/N | Ratio | S/N | ISTD Resp. | Final Conc. |
|-----------------|-------|---------|------------|-------|------------|------------|----------------|
| Cocaethylene | 4.714 | 3470772 | 1494068.05 | 54.0 | 3013967.64 | 252597 | 245.2547 ng/ml |
| Doxepin | 5.356 | 176029 | 1286.95 | 46.2 | 287.60 | 31611 | 247.3842 ng/ml |
| Methylphenidate | 4.310 | 6249130 | 196872.67 | 21.3 | 15119.73 | 918545 | 251.8755 ng/ml |
| Midazolam | 5.771 | 446759 | 1132.14 | 94.0 | 4680.57 | 252597 | 249.9860 ng/ml |
| Norketamine | 4.081 | 658110 | 4596.77 | 428.8 | 7965.65 | 694009 | 266.0843 ng/ml |

TS




AM #28 Multi-Drug Quant. Results

Batch results D:\MassHunter\Data\2021\AM 27-28\04152021 AM 28 AG\QuantResults\Evaluated AM 28 P2 AG.batch.bin
Calibration Last Update 4/20/2021 3:29:03 PM

| | | | |
|-------------------------|----------------------|------------------|----------------|
| Instrument | Instrument 1 | Data File | p2 QC 1000.d |
| Type | QC | Sample | p2 QC 1000 |
| Acq. Method | AM 28 MDQ P2.m | Operator | Tamara Salazar |
| Sample Position | P6-D6 | Comment | |
| Injection Volume | 2 | | |
| Acq. Date-Time | 4/15/2021 4:06:23 PM | | |

Sample Chromatogram



| Name | RT | Resp. | S/N | Ratio | S/N | ISTD Resp. | Final Conc. |
|-----------------|------------------|--------------------|--------------------|------------------|-----------------------|-------------------|---------------------------|
| Cocaethylene | 4.714 | 13401859 | 71683.99 | 54.0 | 1499829.41 | 241773 | 987.3863 ng/ml |
| Doxepin | 5.356 | 747284 | 8103.93 | 48.8 | 547.79 | 33136 | 1000.4871 ng/ml |
| Methylphenidate | 4.318 | 20208096 | 87027.76 | 24.0 | 65637.61 | 819499 | 912.7286 ng/ml |
| Midazolam | 5.764 | 1678761 | 1712.81 | 88.7 | 914.68 | 241773 | 980.8564 ng/ml |
| * Norketamine | 4.081 | 1960595 | 3846.42 | 442.4 | 2529492.49 | 767339 | 718.5381 ng/ml |

*Outside curve range



TS


AM #28 Multi-Drug Quant. Calibration Curve Report

Batch results D:\MassHunter\Data\2021\AM 27-28\04152021 AM 28 AG\QuantResults\Evaluated AM 28 P2
AG.batch.bin

Last Cal. Update 4/20/2021 3:29 PM

Analyst Name ISP\Datastor

Analyte Doxepin **Internal Standard** Doxepin-D3



| Sample | Level | Enabled | Expected Concentration | Final Concentration | Accuracy |
|-----------------|-------|---------|------------------------|---------------------|----------|
| p2 Cal 1-5ng | 1 | ✓ | 5.0 | 5.0 | 100.5 |
| p2 Cal 2- 10ng | 2 | ✓ | 10.0 | 10.2 | 101.7 |
| p2 Cal 3 -25ng | 3 | ✓ | 25.0 | 23.0 | 92.0 |
| p2 Cal 4-50ng | 4 | ✓ | 50.0 | 51.0 | 101.9 |
| p2 Cal 5-100ng | 5 | ✓ | 100.0 | 101.5 | 101.5 |
| p2 Cal 6-250ng | 6 | ✓ | 250.0 | 253.7 | 101.5 |
| p2 Cal 7-500ng | 7 | ✓ | 500.0 | 505.1 | 101.0 |
| p2 Cal 8-1000ng | 8 | ✓ | 1000.0 | 999.2 | 99.9 |



TS


AM #28 Multi-Drug Quant. Calibration Curve Report

Batch results D:\MassHunter\Data\2021\AM 27-28\04152021 AM 28 AG\QuantResults\Evaluated AM 28 P2
AG.batch.bin

Last Cal. Update 4/20/2021 3:29 PM

Analyst Name ISP\Datastor

Analyte Norketamine **Internal Standard** acetyl-norfentanyl-d5



| Sample | Level | Enabled | Expected Concentration | Final Concentration | Accuracy |
|-----------------|-------|---------|------------------------|---------------------|----------|
| p2 Cal 1-5ng | 1 | ✓ | 5.0 | 4.7 | 94.6 |
| p2 Cal 2- 10ng | 2 | ✓ | 10.0 | 10.6 | 105.8 |
| p2 Cal 3 -25ng | 3 | ✓ | 25.0 | 27.6 | 110.4 |
| p2 Cal 4-50ng | 4 | ✓ | 50.0 | 53.2 | 106.3 |
| p2 Cal 5-100ng | 5 | ✓ | 100.0 | 100.5 | 100.5 |
| p2 Cal 6-250ng | 6 | ✓ | 250.0 | 238.8 | 95.5 |
| p2 Cal 7-500ng | 7 | ✓ | 500.0 | 434.5 | 86.9 |
| p2 Cal 8-1000ng | 8 | x | 1000.0 | 693.0 | 69.3 |


*Cal 8 removed due to accuracy



TS

AM #28 Multi-Drug Quant. Calibration Curve Report

Batch results D:\MassHunter\Data\2021\AM 27-28\04152021 AM 28 AG\QuantResults\Evaluated AM 28 P2
 AG.batch.bin
Last Cal. Update 4/20/2021 3:29 PM
Analyst Name ISP\Datastor
Analyte Cocaethylene **Internal Standard** Midazolam-D4




| Sample | Level | Enabled | Expected Concentration | Final Concentration | Accuracy |
|-----------------|-------|---------|------------------------|---------------------|----------|
| p2 Cal 1-5ng | 1 | ✓ | 5.0 | 5.3 | 105.8 |
| p2 Cal 2- 10ng | 2 | ✓ | 10.0 | 9.3 | 93.5 |
| p2 Cal 3 -25ng | 3 | ✓ | 25.0 | 21.7 | 86.9 |
| p2 Cal 4-50ng | 4 | ✓ | 50.0 | 51.2 | 102.4 |
| p2 Cal 5-100ng | 5 | ✓ | 100.0 | 91.7 | 91.7 |
| p2 Cal 6-250ng | 6 | ✓ | 250.0 | 277.0 | 110.8 |
| p2 Cal 7-500ng | 7 | ✓ | 500.0 | 526.7 | 105.3 |
| p2 Cal 8-1000ng | 8 | ✓ | 1000.0 | 1036.9 | 103.7 |

TS



AM #28 Multi-Drug Quant. Calibration Curve Report

Batch results D:\MassHunter\Data\2021\AM 27-28\04152021 AM 28 AG\QuantResults\Evaluated AM 28 P2
 AG.batch.bin
Last Cal. Update 4/20/2021 3:29 PM
Analyst Name ISP\Datastor
Analyte Methylphenidate **Internal Standard** Methylphenidate-D4



| Sample | Level | Enabled | Expected Concentration | Final Concentration | Accuracy |
|-----------------|-------|---------|------------------------|---------------------|----------|
| p2 Cal 1-5ng | 1 | ✓ | 5.0 | 5.0 | 99.3 |
| p2 Cal 2- 10ng | 2 | ✓ | 10.0 | 10.2 | 101.7 |
| p2 Cal 3 -25ng | 3 | ✓ | 25.0 | 24.7 | 98.7 |
| p2 Cal 4-50ng | 4 | ✓ | 50.0 | 50.0 | 100.0 |
| p2 Cal 5-100ng | 5 | ✓ | 100.0 | 101.3 | 101.3 |
| p2 Cal 6-250ng | 6 | ✓ | 250.0 | 260.6 | 104.2 |
| p2 Cal 7-500ng | 7 | ✓ | 500.0 | 514.2 | 102.8 |
| p2 Cal 8-1000ng | 8 | ✓ | 1000.0 | 919.5 | 91.9 |



TS


AM #28 Multi-Drug Quant. Calibration Curve Report

Batch results D:\MassHunter\Data\2021\AM 27-28\04152021 AM 28 AG\QuantResults\Evaluated AM 28 P2
AG.batch.bin

Last Cal. Update 4/20/2021 3:29 PM

Analyst Name ISP\Datastor

Analyte Midazolam **Internal Standard** Midazolam-D4



| Sample | Level | Enabled | Expected Concentration | Final Concentration | Accuracy |
|-----------------|-------|---------|------------------------|---------------------|----------|
| p2 Cal 1-5ng | 1 | ✓ | 5.0 | 5.0 | 99.2 |
| p2 Cal 2- 10ng | 2 | ✓ | 10.0 | 9.9 | 98.7 |
| p2 Cal 3 -25ng | 3 | ✓ | 25.0 | 27.0 | 107.8 |
| p2 Cal 4-50ng | 4 | ✓ | 50.0 | 50.5 | 101.0 |
| p2 Cal 5-100ng | 5 | ✓ | 100.0 | 96.3 | 96.3 |
| p2 Cal 6-250ng | 6 | ✓ | 250.0 | 245.9 | 98.4 |
| p2 Cal 7-500ng | 7 | ✓ | 500.0 | 504.5 | 100.9 |
| p2 Cal 8-1000ng | 8 | ✓ | 1000.0 | 977.1 | 97.7 |

AM #28 Multi-Drug Quant. Results


TS



Batch results D:\MassHunter\Data\2021\AM 27-28\04152021 AM 28 AG\QuantResults\Evaluated AM 28 P2 AG.batch.bin
Calibration Last Update 4/20/2021 3:29:03 PM

| | | | |
|-------------------------|----------------------|------------------|----------------|
| Instrument | Instrument 1 | Data File | p2 Cal 1-5ng.d |
| Type | Cal | Sample | p2 Cal 1-5ng |
| Acq. Method | AM 28 MDQ P2.m | Operator | Tamara Salazar |
| Sample Position | P6-A5 | Comment | |
| Injection Volume | 2 | | |
| Acq. Date-Time | 4/15/2021 1:37:54 PM | | |
| Sample Info. | | | |

Sample Chromatogram



| Name | RT | Resp. | S/N | Ratio | S/N | ISTD Resp. | Final Conc. |
|-----------------|-------|--------|-----------|-------|----------|------------|--------------|
| Cocaethylene | 4.714 | 65639 | 163175.34 | 55.1 | 403.12 | 252850 | 5.2885 ng/ml |
| Doxepin | 5.363 | 3359 | 41.70 | 52.2 | 10.73 | 32589 | 5.0230 ng/ml |
| Methylphenidate | 4.310 | 119948 | 3315.35 | 22.5 | 48514.81 | 909327 | 4.9641 ng/ml |
| Midazolam | 5.764 | 8541 | 532.16 | 90.3 | 11582.17 | 252850 | 4.9613 ng/ml |
| Norketamine | 4.075 | 13072 | 102.83 | 447.7 | 5916.98 | 649433 | 4.7286 ng/ml |

TS




AM #28 Multi-Drug Quant. Results

Batch results D:\MassHunter\Data\2021\AM 27-28\04152021 AM 28 AG\QuantResults\Evaluated AM 28 P2 AG.batch.bin
Calibration Last Update 4/20/2021 3:29:03 PM

| | | | |
|-------------------------|----------------------|------------------|------------------|
| Instrument | Instrument 1 | Data File | p2 Cal 2- 10ng.d |
| Type | Cal | Sample | p2 Cal 2- 10ng |
| Acq. Method | AM 28 MDQ P2.m | Operator | Tamara Salazar |
| Sample Position | P6-B5 | Comment | |
| Injection Volume | 2 | | |
| Acq. Date-Time | 4/15/2021 1:48:41 PM | | |
| Sample Info. | | | |

Sample Chromatogram



| Name | RT | Resp. | S/N | Ratio | S/N | ISTD Resp. | Final Conc. |
|-----------------|-------|--------|-----------|-------|-----------------|------------|---------------|
| Cocaethylene | 4.714 | 128338 | 281631.94 | 54.5 | 91319.06 | 263115 | 9.3500 ng/ml |
| Doxepin | 5.356 | 6540 | 71.93 | 51.1 | 8.78 Low | 29858 | 10.1659 ng/ml |
| Methylphenidate | 4.317 | 232950 | 14581.23 | 21.5 | 5062.61 | 854757 | 10.1687 ng/ml |
| Midazolam | 5.764 | 18026 | 7680.97 | 96.3 | 7343.81 | 263115 | 9.8664 ng/ml |
| Norketamine | 4.075 | 27222 | 202.25 | 444.4 | 7755.84 | 665192 | 10.5843 ng/ml |

TS




AM #28 Multi-Drug Quant. Results

Batch results D:\MassHunter\Data\2021\AM 27-28\04152021 AM 28 AG\QuantResults\Evaluated AM 28 P2 AG.batch.bin
Calibration Last Update 4/20/2021 3:29:03 PM

| | | | |
|-------------------------|----------------------|------------------|------------------|
| Instrument | Instrument 1 | Data File | p2 Cal 3 -25ng.d |
| Type | Cal | Sample | p2 Cal 3 -25ng |
| Acq. Method | AM 28 MDQ P2.m | Operator | Tamara Salazar |
| Sample Position | P6-C5 | Comment | |
| Injection Volume | 2 | | |
| Acq. Date-Time | 4/15/2021 1:59:16 PM | | |
| Sample Info. | | | |

Sample Chromatogram



| Name | RT | Resp. | S/N | Ratio | S/N | ISTD Resp. | Final Conc. |
|-----------------|-------|--------|-----------|-------|-----------|------------|---------------|
| Cocaethylene | 4.714 | 303905 | 295166.35 | 56.0 | 350450.92 | 256998 | 21.7170 ng/ml |
| Doxepin | 5.356 | 15355 | 43.24 | 47.5 | 40.69 | 30195 | 23.0029 ng/ml |
| Methylphenidate | 4.310 | 525474 | 228.06 | 21.8 | 5593.54 | 791109 | 24.6654 ng/ml |
| Midazolam | 5.764 | 48694 | 63564.42 | 88.9 | 1390.14 | 256998 | 26.9502 ng/ml |
| Norketamine | 4.075 | 68592 | 356.05 | 432.7 | 25680.88 | 677038 | 27.5889 ng/ml |

TS




AM #28 Multi-Drug Quant. Results

Batch results D:\MassHunter\Data\2021\AM 27-28\04152021 AM 28 AG\QuantResults\Evaluated AM 28 P2 AG.batch.bin
Calibration Last Update 4/20/2021 3:29:03 PM

| | | | |
|-------------------------|----------------------|------------------|-----------------|
| Instrument | Instrument 1 | Data File | p2 Cal 4-50ng.d |
| Type | Cal | Sample | p2 Cal 4-50ng |
| Acq. Method | AM 28 MDQ P2.m | Operator | Tamara Salazar |
| Sample Position | P6-D5 | Comment | |
| Injection Volume | 2 | | |
| Acq. Date-Time | 4/15/2021 2:09:52 PM | | |
| Sample Info. | | | |

Sample Chromatogram



| Name | RT | Resp. | S/N | Ratio | S/N | ISTD Resp. | Final Conc. |
|-----------------|-------|---------|------------|-------|-----------|------------|---------------|
| Cocaethylene | 4.714 | 848477 | 2544755.87 | 56.1 | 352355.42 | 299019 | 51.1774 ng/ml |
| Doxepin | 5.356 | 38892 | 640.00 | 46.9 | 391.96 | 34138 | 50.9717 ng/ml |
| Methylphenidate | 4.310 | 1656764 | 1163.62 | 21.3 | 3850.13 | 1227899 | 50.0192 ng/ml |
| Midazolam | 5.764 | 106523 | 770.08 | 94.6 | 722.14 | 299019 | 50.5038 ng/ml |
| Norketamine | 4.075 | 146658 | 1072.16 | 430.9 | 3247.39 | 763399 | 53.1571 ng/ml |

TS




AM #28 Multi-Drug Quant. Results

Batch results D:\MassHunter\Data\2021\AM 27-28\04152021 AM 28 AG\QuantResults\Evaluated AM 28 P2 AG.batch.bin
Calibration Last Update 4/20/2021 3:29:03 PM

| | | | |
|-------------------------|----------------------|------------------|------------------|
| Instrument | Instrument 1 | Data File | p2 Cal 5-100ng.d |
| Type | Cal | Sample | p2 Cal 5-100ng |
| Acq. Method | AM 28 MDQ P2.m | Operator | Tamara Salazar |
| Sample Position | P6-E5 | Comment | |
| Injection Volume | 2 | | |
| Acq. Date-Time | 4/15/2021 2:20:26 PM | | |
| Sample Info. | | | |

Sample Chromatogram



| Name | RT | Resp. | S/N | Ratio | S/N | ISTD Resp. | Final Conc. |
|-----------------|-------|---------|-----------|-------|------------|------------|----------------|
| Cocaethylene | 4.714 | 1501847 | 21839.46 | 54.9 | 1722425.37 | 293765 | 91.6715 ng/ml |
| Doxepin | 5.356 | 78058 | 126.96 | 44.6 | 245.75 | 34252 | 101.5098 ng/ml |
| Methylphenidate | 4.310 | 2818039 | 753704.04 | 21.5 | 14410.01 | 1030169 | 101.3248 ng/ml |
| Midazolam | 5.764 | 199981 | 3272.95 | 98.0 | 1108.72 | 293765 | 96.3356 ng/ml |
| Norketamine | 4.075 | 269396 | 884.23 | 432.3 | 7867.16 | 747943 | 100.4843 ng/ml |

TS




AM #28 Multi-Drug Quant. Results

Batch results D:\MassHunter\Data\2021\AM 27-28\04152021 AM 28 AG\QuantResults\Evaluated AM 28 P2 AG.batch.bin
Calibration Last Update 4/20/2021 3:29:03 PM

| | | | |
|-------------------------|----------------------|------------------|------------------|
| Instrument | Instrument 1 | Data File | p2 Cal 6-250ng.d |
| Type | Cal | Sample | p2 Cal 6-250ng |
| Acq. Method | AM 28 MDQ P2.m | Operator | Tamara Salazar |
| Sample Position | P6-F5 | Comment | |
| Injection Volume | 2 | | |
| Acq. Date-Time | 4/15/2021 2:31:04 PM | | |

Sample Chromatogram



| Name | RT | Resp. | S/N | Ratio | S/N | ISTD Resp. | Final Conc. |
|-----------------|-------|---------|----------|-------|----------|------------|----------------|
| Cocaethylene | 4.714 | 4551541 | 88446.08 | 54.8 | 4970.46 | 293206 | 276.9929 ng/ml |
| Doxepin | 5.356 | 197809 | 712.53 | 46.3 | 640.79 | 34642 | 253.6572 ng/ml |
| Methylphenidate | 4.310 | 7802050 | 34685.54 | 21.3 | 381.21 | 1108556 | 260.5633 ng/ml |
| Midazolam | 5.764 | 510075 | 2170.48 | 94.1 | 6921.07 | 293206 | 245.8879 ng/ml |
| Norketamine | 4.075 | 715189 | 1802.89 | 433.5 | 22105.25 | 840018 | 238.8052 ng/ml |

TS




AM #28 Multi-Drug Quant. Results

Batch results D:\MassHunter\Data\2021\AM 27-28\04152021 AM 28 AG\QuantResults\Evaluated AM 28 P2 AG.batch.bin
Calibration Last Update 4/20/2021 3:29:03 PM

| | | | |
|-------------------------|----------------------|------------------|------------------|
| Instrument | Instrument 1 | Data File | p2 Cal 7-500ng.d |
| Type | Cal | Sample | p2 Cal 7-500ng |
| Acq. Method | AM 28 MDQ P2.m | Operator | Tamara Salazar |
| Sample Position | P6-G5 | Comment | |
| Injection Volume | 2 | | |
| Acq. Date-Time | 4/15/2021 2:41:39 PM | | |
| Sample Info. | | | |

Sample Chromatogram



| Name | RT | Resp. | S/N | Ratio | S/N | ISTD Resp. | Final Conc. |
|-----------------|-------|----------|-----------|-------|-----------|------------|----------------|
| Cocaethylene | 4.714 | 7945780 | 290653.64 | 54.4 | 223575.30 | 268870 | 526.7205 ng/ml |
| Doxepin | 5.356 | 383130 | 949.16 | 47.0 | 615.30 | 33663 | 505.1448 ng/ml |
| Methylphenidate | 4.310 | 14195405 | 234990.21 | 21.5 | 94725.41 | 1021940 | 514.1821 ng/ml |
| Midazolam | 5.764 | 960035 | 1118.88 | 90.6 | 2560.00 | 268870 | 504.4848 ng/ml |
| Norketamine | 4.075 | 1178635 | 13051.77 | 433.6 | 7964.29 | 762124 | 434.5428 ng/ml |

TS




AM #28 Multi-Drug Quant. Results

Batch results D:\MassHunter\Data\2021\AM 27-28\04152021 AM 28 AG\QuantResults\Evaluated AM 28 P2 AG.batch.bin
Calibration Last Update 4/20/2021 3:29:03 PM

Instrument Instrument 1
Type Cal
Acq. Method AM 28 MDQ P2.m
Sample Position P6-H5
Injection Volume 2
Acq. Date-Time 4/15/2021 2:52:15 PM
Sample Info.

Data File p2 Cal 8-1000ng.d
Sample p2 Cal 8-1000ng
Operator Tamara Salazar
Comment

Sample Chromatogram



| Name | RT | Resp. | S/N | Ratio | S/N | ISTD Resp. | Final Conc. |
|--------------------------|------------------|--------------------|--------------------|------------------|----------------------|-------------------|---------------------------|
| Cocaethylene | 4.714 | 13696990 | 260566.77 | 54.3 | 10282.69 | 235279 | 1036.9481 ng/ml |
| Doxepin | 5.356 | 901704 | 1381.71 | 49.3 | 259.23 | 40033 | 999.2500 ng/ml |
| Methylphenidate | 4.317 | 21153577 | 190027.97 | 24.3 | 142961.66 | 851567 | 919.4522 ng/ml |
| Midazolam | 5.764 | 1627485 | 2726.57 | 93.0 | 3264.15 | 235279 | 977.1421 ng/ml |
| * Norketamine | 4.081 | 1915102 | 7918.79 | 443.2 | 213788.10 | 777069 | 693.0441 ng/ml |

*Outside curve range