

~~GA~~

5/6/2021

Worklist: 4955

<u>LAB CASE</u>	<u>ITEM</u>	<u>ITEM TYPE</u>	<u>DESCRIPTION</u>
C2021-0896	1	UCK	AM 3 Urine Carboxy-THC
C2021-0939	1	UCK	AM 3 Urine Carboxy-THC
C2021-0987	1	UCK	AM 3 Urine Carboxy-THC
C2021-1019	1	UCK	AM 3 Urine Carboxy-THC

REVIEWED

By Britany Wylie at 12:47 pm, May 10, 2021



AM 3: Carboxy-THC Urine Extraction

Extraction Date: 5/6/21

Analyst: Anne Nord

Negative Urine Lot: 2121

GC/MS ID: 65198

Hexane Lot: fisher 42712

Ethyl Acetate Lot: fisher 020419

Silylating Agent Lot: Cerilliant FN08181601

Pre-Analytic:

- 1. *Positive Control Working Solution Preparation Instructions:*
Add 180uL of 100ug/mL 11-nor-9-carboxy- Δ 9-THC Stock Solution to 9.82mL Methanol. Other volumes may be prepared. Solution is stable for 1-year or the expiration of the stock reference material (whichever is sooner). Store under refrigeration
- 2. Verify Tune and Tune evaluation completed within the previous 7 days. Tune and Tune evaluation reports initialed and filed.
- 3. Create GCMS sequence to include controls, case blanks and case samples.

Analytic:

- 1. Remove working solutions, external control, negative urine and case samples from cold storage.
- 2. To each labeled round bottomed tube add 3mL sample, using negative urine sample for both negative and positive control. Positive control: added working solution to equal approximately 60 ng/ml Carboxy-THC
- 3. Add 500uL 1N KOH to all tubes. Check pH. (*If pH <12, add additional 500ul 1N KOH*). (*Note: put a mark on the tube or separate the tubes that have a pH <12 as you will need to know this in step 5*).
- 4. Place tubes in 40C water bath for ~15 min. Remove and allow to cool.
- 5. If original pH was >12, add 1.5mL pH 1.8 Saturated Phosphate Buffer and 3mL Hexane/Ethyl Acetate (87:13)
If original pH was <12, add 3mL pH 1.8 Saturated Phosphate Buffer and 4mL Hexane/Ethyl Acetate (87:13)
- 6. Rock at ambient temp for ~ 10 minutes.
- 7. Centrifuge for ~ 10 min at ~3500rpm.
- 8. Transfer solvent to tapered bottom tube and evaporate to dryness under nitrogen @ 37C.
- 9. Add 50uL Ethyl Acetate and 50uL silylating reagent, cap and vortex. Heat @ 95C for 15min, then allow to cool.
- 10. Transfer sample to labeled ALS vial with insert.
- 11. Place ALS Vials in appropriate location on GCMS rack and run using appropriate GCMS method.

Post-Analytic

- 1. Complete Data analysis on all samples and corresponding sample blanks
- 2. Did positive and negative control samples provide intended response? yes
- 3. **Criteria for ID:** RT +/1 0.1 min., Ion Ratio of 347:473 & 371:488 within +/- 20%.
Sample response greater than Min Corrected Area, Diluted samples
- 4. Central File Packet to include: LIMS Worklist, Method Checklist, Working solution prep sheet, and Control sample GCMS data printouts

Comments:



Toxicology AM 3 Method 3 spiked 60 ng/ml

Stock solution 3.6 ul (100 ug/ml) C-THC in 5.996 mls neg urine lot 2121

Ppd 5/6/21 expires 11/6/21 lot 2621

Drug	Lot	expiration
C-THC	FE01061702	3/1/2022

~~A~~

Toxicology AM 3 Commercial control

Pinpoint THC quantitation plate (IDP-108-2) lot 201206 exp 06/06/21

Added 3 mls negative urine (2121) to Cal 5 (approximate concentration 25 ng/ml THC-COOH)



Run Date: 5/6/2021

Worklist #: 4955

Laboratory Case #: Positive control PinPoint

Case Calculations:

EIC Ion 371 Corr. Area: 1027094

EIC Ion 473 Corr. Area: 486085

EIC Ion 488 Corr. Area: 289348

473:371 Ratio: **0.473**

488:371 Ratio: **0.282**

Quality Control Calculations:

1st Control

2nd Control

EIC Ion 371 Corr. Area: 1313508

1524670

EIC Ion 473 Corr. Area: 612754

716147

EIC Ion 488 Corr. Area: 388564

452457

Ratio of 473:371: 0.467

0.470

Average: 0.468

Ratio of 488:371: 0.296

0.297

Average: 0.296

Acceptable Retention Time Range: 4.404

to 4.594

Approximate Minimum Corrected Area of 371 Ion: 262702

Acceptable 473:371 Ratio Range: 0.374

to 0.562

Acceptable 488:371 Ratio Range: 0.237

to 0.356



Data Path : D:\DATA\2021\am 2\050621\
 Data File : 01401012.D
 Acq On : 06 May 2021 18:03
 Operator : Instrument 65198
 Sample : C-THC control pinpoint
 Misc : am 3
 ALS Vial : 14 Sample Multiplier: 1

Integration Parameters: autoint1.e
 Integrator: ChemStation

Method : D:\MassHunter\GCMS\1\methods\TOXI-A 10115.M
 Title :

Signal : EIC Ion 371.00 (370.70 to 371.70): 01401012.D\data.ms

peak #	R.T. min	first scan	max scan	last scan	PK TY	peak height	corr. area	corr. % max.	% of total
1	4.241	3	11	35	BB 2	3469	60280	5.87%	5.544%
2	4.493	51	58	81	VB	60883	1027094	100.00%	94.456%

Sum of corrected areas: 1087374
 Signal : EIC Ion 473.00 (472.70 to 473.70): 01401012.D\data.ms


peak #	R.T. min	first scan	max scan	last scan	PK TY	peak height	corr. area	corr. % max.	% of total
1	4.495	36	58	81	BB	28422	486085	100.00%	84.480%
2	4.931	128	139	154	BB	3832	89302	18.37%	15.520%

Sum of corrected areas: 575387
 Signal : EIC Ion 488.00 (487.70 to 488.70): 01401012.D\data.ms

peak #	R.T. min	first scan	max scan	last scan	PK TY	peak height	corr. area	corr. % max.	% of total
1	4.441	36	48	53	BV	8943	143464	49.58%	33.147%
2	4.497	53	59	69	VV	17683	289348	100.00%	66.853%

Sum of corrected areas: 432812

TOXI-A 10115.M Mon May 10 10:09:55 2021





Data Path : D:\DATA\2021\am 2\050621\
Data File : 00801002.D
Acq On : 06 May 2021 16:29
Operator : Instrument 65198
Sample : 60ng C-THC
Misc : am3
ALS Vial : 8 Sample Multiplier: 1

Integration Parameters: autoint1.e
Integrator: ChemStation

Method : D:\MassHunter\GCMS\1\methods\TOXI-A 10115.M
Title :

Signal : EIC Ion 371.00 (370.70 to 371.70): 00801002.D\data.ms

peak #	R.T. min	first scan	max scan	last scan	PK TY	peak height	corr. area	corr. % max.	% of total
1	4.504	41	60	83	BB	73156	1313508	100.00%	100.000%

Sum of corrected areas: 1313508
Signal : EIC Ion 473.00 (472.70 to 473.70): 00801002.D\data.ms


peak #	R.T. min	first scan	max scan	last scan	PK TY	peak height	corr. area	corr. % max.	% of total
1	4.506	50	60	84	BB	33965	612754	100.00%	100.000%

Sum of corrected areas: 612754
Signal : EIC Ion 488.00 (487.70 to 488.70): 00801002.D\data.ms

peak #	R.T. min	first scan	max scan	last scan	PK TY	peak height	corr. area	corr. % max.	% of total
1	4.452	40	50	55	BV	4911	82258	21.17%	17.471%
2	4.508	55	61	83	VB	21360	388564	100.00%	82.529%

Sum of corrected areas: 470823

TOXI-A 10115.M Mon May 10 10:02:39 2021





Data Path : D:\DATA\2021\am 2\050621\
 Data File : 00801015.D
 Acq On : 06 May 2021 18:31
 Operator : Instrument 65198
 Sample : 60 ng CTHC control
 Misc : am 3 end of run control
 ALS Vial : 8 Sample Multiplier: 1

Integration Parameters: autoint1.e
 Integrator: ChemStation

Method : D:\MassHunter\GCMS\1\methods\TOXI-A 10115.M
 Title :

Signal : EIC Ion 371.00 (370.70 to 371.70): 00801015.D\data.ms

peak #	R.T. min	first scan	max scan	last scan	PK TY	peak height	corr. area	corr. % max.	% of total
1	4.242	3	11	35	BB	2953	55156	3.62%	3.491%
2	4.494	51	58	92	VB	90164	1524670	100.00%	96.509%

Sum of corrected areas: 1579826
 Signal : EIC Ion 473.00 (472.70 to 473.70): 00801015.D\data.ms


peak #	R.T. min	first scan	max scan	last scan	PK TY	peak height	corr. area	corr. % max.	% of total
1	4.495	49	58	80	BB	42372	716147	100.00%	92.398%
2	4.932	126	139	153	BB	2532	58923	8.23%	7.602%

Sum of corrected areas: 775070
 Signal : EIC Ion 488.00 (487.70 to 488.70): 00801015.D\data.ms

peak #	R.T. min	first scan	max scan	last scan	PK TY	peak height	corr. area	corr. % max.	% of total
1	4.441	39	48	53	BV	5939	96909	21.42%	17.640%
2	4.497	53	59	79	VB	26671	452457	100.00%	82.360%

Sum of corrected areas: 549366

TOXI-A 10115.M Mon May 10 10:03:13 2021





Data Path : D:\DATA\2021\am 2\050621\
Data File : 00901003.D
Acq On : 06 May 2021 16:39
Operator : Instrument 65198
Sample : negative control am 3
Misc : am 3
ALS Vial : 9 Sample Multiplier: 1

Integration Parameters: autoint1.e
Integrator: ChemStation

Method : D:\MassHunter\GCMS\1\methods\TOXI-A 10115.M
Title :

Signal : EIC Ion 371.00 (370.70 to 371.70): 00901003.D\data.ms

peak #	R.T. min	first scan	max scan	last scan	PK TY	peak height	corr. area	corr. % max.	% of total
--------	----------	------------	----------	-----------	-------	-------------	------------	--------------	------------

No peaks were detected using the method integration parameters!

Signal : EIC Ion 473.00 (472.70 to 473.70): 00901003.D\data.ms

peak #	R.T. min	first scan	max scan	last scan	PK TY	peak height	corr. area	corr. % max.	% of total
--------	----------	------------	----------	-----------	-------	-------------	------------	--------------	------------

No peaks were detected using the method integration parameters!

Signal : EIC Ion 488.00 (487.70 to 488.70): 00901003.D\data.ms

peak #	R.T. min	first scan	max scan	last scan	PK TY	peak height	corr. area	corr. % max.	% of total
1	4.445	40	49	64	BB	5314	91574	100.00%	100.000%

Sum of corrected areas: 91574

TOXI-A 10115.M Mon May 10 10:04:01 2021

